

Bacampicillin


Molecular formula: $C_{21}H_{27}N_3O_7S$

Molecular weight: 465.53

CAS Registry No.: 50972-17-3, 37661-08-8 (HCl)

Merck Index: 961

Lednicer No.: 3 204



SAMPLE

Matrix: blood

Sample preparation: Add 100000 fold excess of ampicillin to eliminate adsorption of bacampicillin, extract into butyl acetate, re-extract into pH 2 buffer. Wash the aqueous phase with n-hexane, inject an aliquot.

HPLC VARIABLES

Column: $100 \times 2.9 \mu\text{m}$ Nucleosil C18

Mobile phase: MeCN:pH 7.4 phosphate buffer 41:49

Injection volume: 20

Detector: F with post-column reaction. The column effluent mixed with two volumes of sodium borate buffer then with one volume of 150 $\mu\text{g/mL}$ fluorescamine in acetone (Science 1972, 178, 871).

CHROMATOGRAM

Limit of detection: 800 pg/mL

REFERENCE

Sjövall,J.; Westerlund,D.; Alván,G.; Magni,L.; Nord,C.E.; Sörstad,J. Rectal bioavailability of bacampicillin hydrochloride in man as determined by reversed-phase liquid chromatography, *Chemotherapy*, 1984, 30, 137–147.

SAMPLE

Matrix: blood

Sample preparation: 0.5 mL Plasma + 1 mL MeOH, stir for 5 min, centrifuge at 2400 g for 10 min. Remove 1 mL supernatant, add 2 μg cefazolin, inject an aliquot.

HPLC VARIABLES

Column: $300 \times 3.9 \mu\text{m}$ $\mu\text{Bondapak C18}$

Mobile phase: MeOH:67 mM KH_2PO_4 20:80

Flow rate: 1.5

Injection volume: 50

Detector: UV 225

CHROMATOGRAM

Retention time: 9 (measured as ampicillin peak)

Internal standard: cefazolin (14)

Limit of detection: 500 ng/mL

KEY WORDS

plasma

REFERENCE

Marzo,A.; Monti,N.; Ripamonti,M.; Arrigoni Martelli,E.; Picari,M. High-performance liquid chromatographic assay of ampicillin and its prodrug lenampicillin, *J.Chromatogr.*, 1990, 507, 235–239.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A: B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 14.687

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, 1997, 763, 149-163.


Bacitracin

Molecular formula: C₆₆H₁₀₃N₁₇O₁₆S (bacitracin A)

Molecular weight: 1422.61 (bacitracin A)

CAS Registry No.: 1405-87-4, 1405-89-6 (zinc salt), 1405-88-5 (methylenedisalicylate)

Merck Index: 965



SAMPLE

Matrix: bulk

Sample preparation: Dissolve 1 g in 100 mL 20 mM HCl in MeOH:water 80:20, mix, vortex, centrifuge, filter (0.45 µm), inject a 100 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm YMC basic C8 200 Å (YMC)

Mobile phase: Gradient. MeOH:50 mM pH 6.5 KH₂PO₄ from 57:43 to 63:37 over 1 h

Column temperature: 25

Flow rate: 1

Injection volume: 100

Detector: UV 215

CHROMATOGRAM

Retention time: 14, 16, 20, 21, 25, 27, 30 (bioactive fractions)

REFERENCE

Bell, R.G. Preparative high-performance liquid chromatographic separation and isolation of bacitracin components and their relationship to microbiological activity, *J. Chromatogr.*, **1992**, 590, 163–168.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in saline at a concentration of 1-10 mg/mL, inject an aliquot.

HPLC VARIABLES

Column: 300 × 3.9 C18 (Vydac)

Mobile phase: Gradient. MeCN:water (both containing 0.05% trifluoroacetic acid) from 5:95 to 65:35 over 20 min

Flow rate: 2

Detector: UV 210

KEY WORDS

saline

REFERENCE

Drapeau,G.; Petittclerc,E.; Toulouse,A.; Marceau,F. Dissociation of the antimicrobial activity of bacitracin USP from its renovascular effects, *Antimicrob.Agents Chemother.*, **1992**, 36, 955-961.

Baclofen


Molecular formula: $C_{10}H_{12}ClNO_2$

Molecular weight: 213.66

CAS Registry No.: 1134-47-0

Merck Index: 967

Lednicer No.: 2 121



SAMPLE

Matrix: CSF

Sample preparation: Prepare a 100×7 column of Dowex 50X4-400, treat with an excess of aqueous ammonia, wash with water until the washings are neutral, wash with an excess of 4 M HCl, wash with water until the washings are neutral, rinse with 8 mL water. Add 200 μ L CSF to the column, wash with 8 mL water, elute with 2 mL 10% ammonia. Lyophilize the eluate, reconstitute the residue with 1 mL pH 7.4 phosphate buffer, extract twice with 1 mL portions of 1-butanol. Combine the organic layers and evaporate them to dryness, add 10-50 μ L EtOH:water:triethylamine 50:25:25, evaporate to dryness, add 20 μ L reagent, let stand at room temperature for 20 min, evaporate to dryness, reconstitute with pH 7.4 sodium phosphate buffer, inject a 5-25 μ L aliquot. (Reagent was EtOH:triethylamine:water:phenylisothiocyanate 70:10:10:10, store at -20° (Anal.Biochem. 1989, 176, 269).)

HPLC VARIABLES

Column: 150×3.9 Pico-Tag

Mobile phase: Gradient. A was sodium acetate adjusted to pH 6.4 with glacial acetic acid. B was MeCN:water 60:40. A:B from 100:0 to 60:40 over 10 min, to 0:100 over 1 min, maintain at 0:100 over 3 min.

Column temperature: 38

Injection volume: 5-25

Detector: UV 254

CHROMATOGRAM

Retention time: 13.2

Limit of detection: 5-10 ng/mL

KEY WORDS

derivatization; SPE; pharmacokinetics

REFERENCE

Sallerin-Caute,B.; Monsarrat,B.; Lazorthes,Y.; Cros,J.; Bastide,R. A sensitive method for the determination of baclofen in human CSF by high performance liquid chromatography, *J.Liq.Chromatogr.*, 1988, 11, 1753-1761.

SAMPLE

Matrix: blood

Sample preparation: Condition a 1 mL Bond Elut SCX strong cation exchange SPE cartridge with 2 mL hexane, 2 mL MeOH, 2 mL water, and 3 mL saturated NaCl solution. 1 mL Plasma + 100 μ L water + 1 mL citrate buffer, mix, add to the SPE cartridge, wash with 4 mL water, wash with 1 mL saturated NaCl solution, dry the SPE cartridge, elute with 1.5 mL pH 10.4 borate buffer. Mix a 200 μ L aliquot of the eluate and add it to 50 μ L reagent, mix, inject a 20 μ L aliquot. (Prepare citrate buffer by diluting 89 mL 100 mM citric acid to 100 mL with 200 mM Na_2HPO_4 , pH 2.6. Prepare pH 10.4 borate buffer by mixing 54 mL 200 mM boric acid in 100 mM NaOH with 46 mL 100 mM NaOH. Prepare pH 9.3 borate buffer by mixing 87 mL 200 mM boric acid in 100 mM NaOH with

13 mL 100 mM NaOH. Prepare reagent daily by mixing 75 mg o-phthalaldehyde, 5 mL MeOH, 50 μ L tert-butanethiol, and 5 mL pH 9.3 borate buffer.)

HPLC VARIABLES

Column: 150 \times 3.9 4 μ m Novapak

Mobile phase: MeOH:buffer 74:36 (Prepare buffer by adjusting the pH of 60 mM Na₂HPO₄ to 7 with 60 mM KH₂PO₄.)

Flow rate: 0.8

Injection volume: 20

Detector: E, ESA Coulochem II, Model 5020 guard cell +1.2 V, Model 5011 glassy carbon working cell, screen electrode +0.2 V, quantifying electrode +0.7 V

CHROMATOGRAM

Retention time: 26

Limit of detection: 2.5 ng/mL

Limit of quantitation: 10 ng/mL

KEY WORDS

derivatization; plasma; SPE; pharmacokinetics

REFERENCE

Millerioux,L.; Brault,M.; Gualano,V.; Mignot,A. High-performance liquid chromatographic determination of baclofen in human plasma, *J.Chromatogr.A*, **1996**, 729, 309–314.

SAMPLE

Matrix: blood, CSF, dialysate, tissue

Sample preparation: Dilute dialysate and CSF. Homogenize brain tissue with 4 volumes of MeOH:water in an ice bath, let stand at -20° for 8 h. 30 μ L Plasma + 60 μ L MeOH, mix, let stand at -20° for 1 h, centrifuge at 4° at 10000 rpm for 5 min, dilute the supernatant with water. Mix an aliquot of tissue homogenate, deproteinized plasma, diluted dialysate, or diluted CSF with an equal volume of the reagent, let stand for 1.5 min, inject a 30 μ L aliquot. (Prepare reagent by mixing 50 mg o-phthalaldehyde, 900 μ L MeOH, 100 μ L 400 mM pH 9.2 borate buffer, and 50 μ L 2-mercaptoethanol.)

HPLC VARIABLES

Guard column: μ Bondapak C18 Guard-Pak

Column: 250 \times 4.6 5 μ m Finepak SIL C18S ODS (Jasco)

Mobile phase: MeOH:THF:100 mM pH 6.95 acetate buffer 45.5:2:52.5 (dialysate, CSF, plasma) or 43:2:55 (tissue)

Flow rate: 1

Injection volume: 30

Detector: F ex 368 em 434

CHROMATOGRAM

Limit of detection: 100 nM

KEY WORDS

plasma; derivatization; rat; brain; pharmacokinetics

REFERENCE

Deguchi,Y.; Inabe,K.; Tomiyasu,K.; Nozawa,K.; Yamada,S.; Kimura,R. Study on brain interstitial fluid distribution and blood-brain barrier transport of baclofen in rats by microdialysis, *Pharm.Res.*, **1995**, 12, 1838–1844.

SAMPLE

Matrix: blood, urine

Sample preparation: 500 μ L Plasma + 2 mL MeOH, centrifuge. 1 mL Supernatant or 100 μ L urine + 500 μ L pH 9 sodium tetraborate buffer + 250 μ L 0.2% 4-chloro-7-nitrobenzofurazan in MeOH, heat at 60° for 45 min, acidify with 100 mM HCl, extract with 5 mL ethyl acetate. Evaporate 3 mL of the extract to about 500 μ L, dry with anhydrous sodium sulfate, pass through SPE column, elute with ethyl acetate to give a final volume of 2 mL. Evaporate it to dryness under a stream of nitrogen at 50°, reconstitute the residue in 1 mL mobile phase, inject a 10 μ L aliquot. (Prepare the 25 \times 6 SPE column with 0.063-0.200 mm silica gel (Merck), wet it with ethyl acetate.)

HPLC VARIABLES

Column: 300 \times 3.9 10 μ m Bondapak C18

Mobile phase: MeOH:water 45:55

Flow rate: 0.4

Injection volume: 10

Detector: F ex 463 em 524

CHROMATOGRAM

Retention time: 4.0

Limit of detection: 100 ng/mL (urine), 20 ng/mL (plasma)

KEY WORDS

plasma; pharmacokinetics; SPE; derivatization

REFERENCE

Tosunoglu,S.; Ersoy,L. Determination of baclofen in human plasma and urine by high-performance liquid chromatography with fluorescence detection, *Analyst*, **1995**, 120, 373–375.

SAMPLE

Matrix: formulations

Sample preparation: Dilute an aliquot with water to a concentration of 5 μ g/mL, filter (0.22 μ m), inject a 15 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Speri-5 ODS (Applied Biosystems)

Mobile phase: MeCN:50 mM NaH₂PO₄ 20:80, adjusted to pH 3.5 with 85% phosphoric acid

Flow rate: 1

Injection volume: 15

Detector: UV 220

CHROMATOGRAM

Retention time: 5.5

KEY WORDS

syrup; stability-indicating

REFERENCE

Johnson,C.E.; Hart,S.M. Stability of an extemporaneously compounded baclofen oral liquid, *Am.J.Hosp.Pharm.*, **1993**, 50, 2353–2355.

SAMPLE

Matrix: formulations

Sample preparation: Dilute formulation 1:10 with water, inject a 50 μ L aliquot.

HPLC VARIABLES

Guard column: 5 \times 4 35-60 μ m Perisorb RP18

Column: 250 \times 4 10 μ m LiChrosorb RP18

Mobile phase: MeOH:MeCN:2.72 g/L KH₂PO₄ 2:12:86

Injection volume: 50

Detector: UV 220

CHROMATOGRAM

Retention time: 10.3

KEY WORDS

injections; water

REFERENCE

Sadjak,A.; Wintersteiger,R. Compatibility of morphine, baclofen, flouxuridine and fluorouracil in an implantable medication pump, *Arzneimittelforschung*, **1995**, *45*, 93–98.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 5 μ m Adsorbosphere C18

Mobile phase: MeOH:buffer 40:60 (Buffer was 50 mM sodium acetate adjusted to pH 6.4 with glacial acetic acid.)

Flow rate: 1.0

Injection volume: 100

Detector: UV 266

CHROMATOGRAM

Limit of quantitation: 10 μ g/mL

REFERENCE

A simple, rapid and reliable HPLC method for the analysis of baclofen in tablets (Abstract 3350), *Pharm.Res.*, **1997**, *14*, S581–S581.

SAMPLE

Matrix: solutions

Sample preparation: Mix sample:50 (?) mM NaCN in 50 mM pH 9.3 borate buffer: 25 (?) mM naphthalene-2,3-dicarboxaldehyde in MeOH 3:1:1, let stand for 15 min, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 200 \times 3 5 μ m Chromspher ODS-2 C18 (Chrompack)

Mobile phase: Gradient. A was THF:50 mM pH 6.8 potassium phosphate buffer 5:95. B was MeCN:MeOH:50 mM pH 6.8 potassium phosphate buffer 55:10:35. A:B from 70:30 to 0:100 over 1 h, maintain at 0:100 for 20 min.

Flow rate: 0.5

Injection volume: 50

Detector: F ex 420

CHROMATOGRAM

Retention time: 32

OTHER SUBSTANCES

Simultaneous: amphetamine, tranlycypromine

KEY WORDS

derivatization

REFERENCE

Koning,H.; Wolf,H.; Venema,K.; Korf,J. Automated precolumn derivatization of amino acids, small peptides, brain amines and drugs with primary amino groups for reversed-phase high-performance liquid chromatography using naphthalenedialdehyde as the fluorogenic label, *J.Chromatogr.*, **1990**, 533, 171-178.

SAMPLE

Matrix: solutions

Sample preparation: 200 μ L 10 mM Baclofen in 100 mM sodium bicarbonate + 200 μ L 10 mM 1-fluoro-2,4-dinitrophenyl-5-L-alanine amide in acetone (freshly prepared), stir at 40° for 1 h, cool, add 100 μ L 200 mM HCl, inject an aliquot. (Derivatization in Anal.Biochem. 1992, 202, 210.)

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m YMC.GEL C8 (YMC)

Mobile phase: MeOH:5% KH₂PO₄ (pH 4.0) 11:8

Flow rate: 1

Detector: UV 340

CHROMATOGRAM

Retention time: 17 (R(-)), 20 (S(+))

KEY WORDS

chiral; derivatization

REFERENCE

Shimada,K.; Mitamura,K.; Morita,M.; Hirakata,K. Separation of the diastereomers of baclofen by high performance liquid chromatography using cyclodextrin as a mobile phase additive, *J.Liq. Chromatogr.*, **1993**, 16, 3311-3320.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Supelcosil LC-DP (A) or 250 \times 4 5 μ m LiChrospher 100 RP-8 (B)

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 5.71 (A), 3.47 (B)

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordiazepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, fu-

roseamide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazindol, mefenamic acid, meperidine, mephénytoin, mepivacaine, mesoridazine, metaproterenol, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl-dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfipyrazone, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocinide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, trifluoromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

KEY WORDS

also details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J.Chromatogr.A*, **1995**, 692, 103–119.

SAMPLE

Matrix: solutions

Sample preparation: Inject an aliquot of a 100 µg/mL solution in mobile phase.

HPLC VARIABLES

Column: 150 × 4.5 µm Crownpak CR(+) immobilized crown ether

Mobile phase: MeOH:0.1% pH 1.9 perchloric acid 15:85

Column temperature: 40

Flow rate: 1

Detector: UV 210

CHROMATOGRAM

Retention time: 9.32, 14.23

OTHER SUBSTANCES

Simultaneous: levodopa, norephedrine, primaquine

KEY WORDS

chiral; comparison with capillary electrophoresis

REFERENCE

Nishi, H.; Nakamura, K.; Nakai, H.; Sato, T. Separation of enantiomers and isomers of amino compounds by capillary electrophoresis and high-performance liquid chromatography utilizing crown ethers, *J.Chromatogr.A*, **1997**, 757, 225–235.


Bambuterol

Molecular formula: $C_{18}H_{29}N_3O_5$

Molecular weight: 367.45

CAS Registry No.: 81732-65-2, 81732-46-9 (HCl)

Merck Index: 980



SAMPLE

Matrix: blood

Sample preparation: Dilute plasma with buffer containing IS, inject an aliquot corresponding to 0.5 μ L plasma.

HPLC VARIABLES

Guard column: 4 \times 4 LiChrosorb RP-select B

Column: 50 \times 4.5 μ m LiChrospher RP-select B

Mobile phase: Gradient. MeOH:100 mM pH 5 ammonium acetate from 3:97 to 10:90 over 5 min, to 25:75 over 1 min, to 38:62 over 7 min

Flow rate: 1.4

Detector: MS, thermospray, Finnigan 4500 quadrupole, ion source repeller equipped with polyimide sleeve, trap at -90° , analyzer pressure 0.000035 Torr, manifold fore pressure 0.20 Torr, exhaust line pressure 1.1 Torr, repeller 45 V, vaporizer 115° , jet 200° , aerosol 270

CHROMATOGRAM

Retention time: 11

Internal standard: hedxadeuterobambuterol

KEY WORDS

plasma; LC-MS; dog

REFERENCE

Lindberg,C.; Paulson,J.; Blomqvist,A. Evaluation of an automated thermospray liquid chromatography-mass spectrometry system for quantitative use in bioanalytical chemistry, *J.Chromatogr.*, **1991**, *554*, 215-226.

SAMPLE

Matrix: microsomal incubations

Sample preparation: Condition a C18 Sep-Pak SPE cartridge with three 3 mL portions of EtOH, two 3 mL portions of water, and with 3 mL 10 mM pH 7.5 phosphate buffer. Add the incubation mixture to the SPE cartridge, wash with two 3 mL portions of water, elute with two 1 mL portions of EtOH:50 mM pH 8.5 ammonium chloride buffer 95:5. Evaporate the eluate to dryness under a stream of nitrogen at 60° , reconstitute the residue in 500 μ L of the initial mobile phase, inject a 200 μ L aliquot.; SPE

HPLC VARIABLES

Column: 150 \times 5 Nucleosil 10SA

Mobile phase: Gradient. A was 250 mM pH 4.6 ammonium acetate buffer. B was MeCN: 500 mM pH 4.6 ammonium acetate buffer 50:50. From A:B 90:10 to 10:90 over 20 min.

Flow rate: 1

Injection volume: 200

Detector: UV 254

CHROMATOGRAM

Retention time: 16.5

OTHER SUBSTANCES

Extracted: terbutaline, metabolites

KEY WORDS

rat; SPE

REFERENCE

Lindberg,C.; Roos,C.; Tunek,A.; Svensson,L.-Å. Metabolism of bambuterol in rat liver microsomes: identification of hydroxylated and demethylated products by liquid chromatography mass spectrometry, *Drug Metab.Dispos.*, **1989**, *17*, 311-322.

SAMPLE

Matrix: perfusate, tissue

Sample preparation: Perfusate. 400 µL Lung perfusate + 400 µL 5% perchloric acid, mix, centrifuge, inject a 200 µL aliquot of the supernatant. Tissue. Homogenize lung in 2 volumes of water (Polytron Homogenizer), mix 400 µL homogenate with 400 µL 5% perchloric acid, mix, centrifuge, inject a 200 µL aliquot.

HPLC VARIABLES

Column: 150 × 5 Nucleosil 10SA

Mobile phase: Gradient. A was 19.3 g ammonium acetate and 14.4 mL acetic acid in 1 L water. B was 19.3 g ammonium acetate and 14.4 mL acetic acid in 1 L MeCN:water 50:50. A:B from 90:10 to 10:90 over 20 min, stay at 10:90 for 3 min, return to initial conditions over 3 min, re-equilibrate for 7 min.

Flow rate: 1

Injection volume: 200

Detector: UV 254

CHROMATOGRAM

Retention time: 18

OTHER SUBSTANCES

Extracted: metabolites, terbutaline


KEY WORDS

lung; guinea pig

REFERENCE

Ryrfeldt,Å.; Nilsson,E.; Tunek,A.; Svensson,L.-Å. Bambuterol: uptake and metabolism in guinea pig isolated lungs, *Pharm.Res.*, **1988**, *5*, 151-155.

Bamethan



Molecular formula: $C_{12}H_{19}NO_2$

Molecular weight: 209.29

CAS Registry No.: 3703-79-5, 5716-20-1 (sulfate)

Merck Index: 981

Lednicer No.: 2 39

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A: B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 5.91

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 μ g/mL solution in MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 125 \times 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAMRetention time: 1.7

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, benactyzine, benperidol, benzethidine, benzocaine, benzoctamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipanone, diprenorphine, dipyridamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J.Chromatogr.*, **1985**, 323, 191-225.

Bamifylline


Molecular formula: C₂₀H₂₇N₅O₃

Molecular weight: 385.47

CAS Registry No.: 2016-63-9, 20684-06-4 (HCl)

Merck Index: 982

Lednicer No.: 1 426



SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A: B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 207.5

CHROMATOGRAM

Retention time: 10.288


KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

Barbital



Molecular formula: $C_8H_{12}N_2O_3$

Molecular weight: 184.19

CAS Registry No.: 57-44-3, 144-02-5 (sodium salt)

Merck Index: 989

Lednicer No.: 1 267

SAMPLE

Matrix: blood

Sample preparation: Inject a 5-20 μ L aliquot of serum directly onto the column with mobile phase A or B.

HPLC VARIABLES

Column: 100 \times 4.6 5-10 μ m Silicalite (by sieving Silicalite, 3M Co.(?))

Mobile phase: MeCN:20 mM pH 6.9 phosphate buffer 8:92 (A) or Gradient. MeCN:20 mM pH 6.9 phosphate buffer from 5:95 to 20:80 over 2 min, to 25:75 over 2 min, to 30:70 over 4 min, to 50:50 over 2 min, maintain at 50:50 for 10 min (B)

Flow rate: 1

Injection volume: 5 (A), 20 (B)

Detector: UV 254

CHROMATOGRAM

Retention time: 7.60 (A), 8 (B)

OTHER SUBSTANCES

Simultaneous: acetaminophen (B), carbamazepine (B), phenobarbital (B), phenytoin (B), primidone (B), sulfapyridine (B)

KEY WORDS

serum

REFERENCE

Ambrose,D.L.; Fntz,J.S. High-performance liquid chromatographic determination of drugs and metabolites in human serum and urine using direct injection and a unique molecular sieve, *J.Chromatogr.B*, **1998**, 709, 89-96.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A: B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 10.445

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenpropfen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaacal, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephénytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methypylon, metoprolol, mibolerone, morphine, nadolol, na-

lorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nyldrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyriethyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, rescinamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sufadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleennamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.


Barnidipine

Molecular formula: $C_{27}H_{29}N_3O_6$

Molecular weight: 491.54

CAS Registry No.: 104713-75-9

Merck Index: 1031



SAMPLE

Matrix: bile, urine

Sample preparation: Inject urine and bile directly. Hydrolyse urine or bile by heating a 1 mL aliquot with 200 μ L 200 mM pH 5.5 sodium acetate buffer and 20 μ L 100 U/mL β -glucuronidase at 37° for 24 h, add to a Sep-Pak C18 SPE cartridge, elute with MeOH. Evaporate the eluate, reconstitute the residue with MeOH:50 mM pH 5.0 sodium acetate buffer 5:95, inject an aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 Nucleosil 5C18

Mobile phase: Gradient. MeOH:50 mM pH 5.0 sodium acetate buffer 5:95 for 5 min, to 50:50 over 95 min

Flow rate: 0.8

Detector: UV

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

rat; dog; SPE

REFERENCE

Teramura,T.; Tokunaga,T.; Matsumoto,H.; Watanabe,T.; Higuchi,S. Metabolism of barnidipine hydrochloride, a potent calcium antagonist, in rat and dog, *Xenobiotica*, **1996**, 26, 177-187.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Serum or plasma + 500 μ L 500 mM pH 11 phosphate buffer + 3 mL benzene (Caution! Benzene is a carcinogen!), shake mechanically for 5 min, centrifuge at 3000 rpm for 10 min. Remove 2.7 mL of the organic layer and add it to 50 ng IS, evaporate to dryness under a stream of nitrogen, reconstitute the residue in 75 μ L mobile phase, inject a 5 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.5 μ m Aluspher RP-select B (Merck)

Mobile phase: MeCN:100 mM pH 11.8 Britton-Robinson buffer 40:60

Flow rate: 0.9

Injection volume: 5

Detector: E, JASCO Model 840-EC, glassy carbon working electrode 0.6 V, Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 8.5

Internal standard: 2-(N-benzyl-N-methylamino) ethylmethyl 1,4-dihydro-2,6-dimethyl-4-phenyl-3,5-pyridinedicarboxylate hydrochloride (YC-204, Yamanouchi) (7.3)

Limit of detection: 1 ng/mL

OTHER SUBSTANCES

Noninterfering: acetazolamide, allopurinol, amitriptyline, chlorpromazine, clonidine, cyproheptadine, hydrochlorothiazide, imipramine, methyldopa, phenoxybenzamine, propranolol, ascorbic acid

KEY WORDS

human; dog; serum; plasma; protect from light; pharmacokinetics

REFERENCE

Takamura,K.; Abdel-Wadood,H.M.; Kusu,F.; Rafaat,I.H.; Saleh,G.A.; El-Rabbat,N.A.; Otagiri,M. Determination of barnidipine in human serum and dog plasma by HPLC with electrochemical detection, *Biol.Pharm.Bull.*, **1995**, 18, 1311–1314.

SAMPLE

Matrix: blood, enzyme incubations

Sample preparation: 1 mL Enzyme incubation or plasma + 5 mL diethyl ether + 500 μ L saturated sodium bicarbonate, extract, centrifuge at 800 g for 5 min. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 100 μ L MeOH, inject an aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 Nucleosil 5C18

Mobile phase: MeOH:20 mM pH 3.0 bromo-tetra-n-propylammonium phosphate 50:50

Flow rate: 1

Detector: UV

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

rat; dog; plasma

REFERENCE

Teramura,T.; Tokunaga,T.; Matsumoto,H.; Watanabe,T.; Higuchi,S. Metabolism of barnidipine hydrochloride, a potent calcium antagonist, in rat and dog, *Xenobiotica*, **1996**, 26, 177–187.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 500 μ L Microsomal incubation + 500 μ L ice-cold MeOH, vortex, centrifuge at 1000 g for 15 min, filter (0.22 μ m), inject an aliquot of the filtrate.

HPLC VARIABLES

Column: 250 \times 4.6 Nucleosil 5C18

Mobile phase: Gradient. MeCN:20 mM pH 7.0 ammonium acetate buffer 30:70 for 5 min, to 70:30 over 55 min, maintain at 70:30 for 30 min.

Flow rate: 1

Detector: radioactivity

CHROMATOGRAM

Retention time: 66

OTHER SUBSTANCES

Extracted: metabolites


KEY WORDS

rat; dog; liver; 14 C labeled

REFERENCE

Teramura,T.; Tokunaga,T.; Matsumoto,H.; Watanabe,T.; Higuchi,S. Metabolism of barnidipine hydrochloride, a potent calcium antagonist, in rat and dog, *Xenobiotica*, **1996**, 26, 177-187.

Beclobrate



Molecular formula: $C_{20}H_{23}ClO_3$

Molecular weight: 346.85

CAS Registry No.: 55937-99-0

Merck Index: 1046

SAMPLE

Matrix: blood

Sample preparation: 250 μ L Plasma + 250 μ L acetic acid + 1.75 mL 10 μ g/mL IS in hexane, shake for 5 min, centrifuge at 2000 rpm for 20 min, inject a 50 μ L aliquot of the supernatant.

HPLC VARIABLES

Guard column: 30 \times 4 30-40 μ m RP18 (E. Merck)

Column: 125 \times 4 5 μ m LiChrosorb RP18

Mobile phase: MeCN:MeOH:water:acetic acid 40:30:30:0.1

Flow rate: 2

Injection volume: 50

Detector: UV 27

CHROMATOGRAM

Retention time: 4.1 (beclobric acid)

Internal standard: SGD 2774 (5.8)

KEY WORDS

plasma; pharmacokinetics; bioequivalence

REFERENCE

Gikalov,I.; Ifflaender,U. Pharmacokinetik und Bioäquivalenz von zwei peroralen Beclobrat-Zubereitungen [Pharmacokinetics and bioequivalence of two peroral beclobrate preparations], *Arzneimittelforschung*, **1987**, *37*, 1065-1068.

SAMPLE

Matrix: blood, urine

Sample preparation: 200 μ L Plasma or 100 μ L urine + 500 (plasma) or 300 (urine) mg NaCl + 50 μ L 1 μ g/mL clobfibric acid in MeOH + 1 (plasma) or 0.3 (urine) mL pH 4 buffer + 5 mL n-hexane:EtOH 90:10, shake horizontally for 10 min, centrifuge. Remove 4 mL of the organic layer and evaporate it to dryness under a stream of nitrogen at 55°, add 50 μ L toluene and evaporate it to remove traces of water. Reconstitute the residue in 500 μ L dichloromethane, add 50 μ L 1 mg/mL 1-hydroxybenzotriazole in dichloromethane:pyridine 99:1, add 50 μ L 1 mg/mL 1-(3-dimethylaminopropyl)-3-ethylcarbodiimide hydrochloride in dichloromethane, add 50 μ L 1 mg/mL FLOPA, vortex, let stand at room temperature for 2 h, evaporate to dryness, reconstitute in 500 μ L mobile phase, inject a 50 μ L aliquot. (FLOPA is the corresponding amine hydrochloride from (+)-(S)-flunoxapropfen. Synthesis is as follows (protect from light). 500 mg (+)-(S)-Flunoxapropfen in 50 mL dry toluene, slowly add 5 mL freshly distilled thionyl chloride, reflux for 1 h, evaporate to dryness under vacuum, dry the acyl chloride under vacuum over KOH for 2 days. Dissolve 0.5 mmoles acyl chloride in 5 mL acetone, add 600 mg sodium azide dissolved in ice water with stirring, stir at 0° for 30 min, add 10 mL ice-cold water, filter, dry solid in a desiccator under vacuum. Dissolve the solid in 1 mL toluene or dichloromethane (dried over 3 Å molecular sieve), reflux for 10 min, evaporate, store resulting isocyanate under vacuum over a desiccant. Dissolve 0.5 mmole isocyanate in 5 mL acetone, add 20 mL 8.5% phosphoric acid, heat to 80° for 1.5 h, adjust to pH 13, extract with diethyl ether:dichloromethane 4:1. Wash the organic layer twice with water, dry over anhydrous sodium sulfate,

evaporate to dryness, dissolve in 1 mL toluene, evaporate to give crystals (mp 91°). Dissolve in ether, add 0.5 M HCl in ether, filter, dissolve solid in a small volume of MeOH, precipitate with ether, dry FLOPA over phosphorus pentoxide under vacuum (Pharm.Res. 1990, 7, 1262).)

HPLC VARIABLES

Column: 250 × 4.6 5 µm Zorbax Sil

Mobile phase: Gradient. A was n-hexane:chloroform:EtOH 100:10:0.75. B was n-hexane:chloroform:EtOH 100:10:20. A:B 100:0 for 10 min, 50:50 for 5 min, 100:0 for 5 min (stepwise).

Flow rate: 2

Injection volume: 50

Detector: F ex 305 em 355

CHROMATOGRAM

Retention time: 6 (-), 6.5 (+)

Internal standard: clofibril acid (8)

Limit of detection: 25 ng/mL

KEY WORDS

plasma; pharmacokinetics; chiral; derivatization; normal phase

REFERENCE

Mayer,S.; Mutschler,E.; Spahn-Langguth,H. Pharmacokinetic studies with the lipid-regulating agent beclobrate: enantiospecific assay for beclobril acid using a new fluorescent chiral coupling component (S-FLOPA), *Chirality*, **1991**, 3, 35–42.


Beclomethasone

Molecular formula: $C_{22}H_{29}ClO_5$

Molecular weight: 408.92

CAS Registry No.: 4419-39-0, 5534-09-8 (beclomethasone dipropionate)

Merck Index: 1047



SAMPLE

Matrix: blood

Sample preparation: Add 30 mL EtOH to 50 mL plasma, extract with dichloromethane for 15 min, centrifuge at 2500 rpm for 5 min. Evaporate dichloromethane layer to dryness under vacuum at 30°, reconstitute residue in mobile phase and inject a 100 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Alltima C18 (Alltech Associates, Australia) (A) or 250 \times 4.6 5 μ m Econosphere C18 (Alltech Associates, Australia) (B) or 150 \times 4.6 5 μ m Ultrasphere C8 (C)

Mobile phase: MeOH:MeCN:acetic acid:water 30.9:14.6:4.4:0.1 (A) or 25:20.5:4.4:0.1 (B) or 27.5:18:4.4:0.1 (C)

Flow rate: 1.3 (A, C) or 1.2 (B)

Injection volume: 100

Detector: UV 242

CHROMATOGRAM

Retention time: 4.9 (A), 5.5 (B), 3.3 (C) (beclomethasone); 23.0 (A), 25.9 (B), 17.6 (C) (beclomethasone dipropionate)

OTHER SUBSTANCES

Extracted: degradation products

KEY WORDS

plasma

REFERENCE

Foe, K.; Cheung, H.T.A.; Tattam, B.N.; Brown, K.F.; Seale, J.P. Degradation products of beclomethasone dipropionate in human plasma, *Drug Metab. Dispos.*, **1998**, *26*, 132–137.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Plasma + 300 μ L water + 6 mL diethyl ether, shake for 5 min, centrifuge at 1000 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen. Dissolve in 100 μ L MeOH. Add 20 μ L copper acetate solution and let stand at room temperature for 1 h. Add 80 μ L diaminophthalhydrazide solution, heat at 80° for 110 min, cool, inject a 20 μ L aliquot. (Prepare copper acetate solution by dissolving 700 mg copper(II) acetate in 10 mL water, dilute to 100 mL with MeOH, discard after 1 month. Diaminophthalhydrazide solution was 7.5 mM 4,5-diaminophthalhydrazide in 3.5 M hydrochloric acid containing 625 mM β -mercaptoethanol, discard after 5 h. Prepare 4,5-diaminophthalhydrazide dihydrochloride as follows. Reflux 316 g 4-nitrophthalic acid and 50 mL concentrated sulfuric acid in 500 mL MeOH for 10 h, recrystallize the product (dimethyl 4-nitrophthalate) from MeOH (mp 64–65E). Hydrogenate 47.8 g dimethyl 4-nitrophthalate in 300 mL MeOH over 13 g 5% platinum on carbon at an initial hydrogen pressure of 50 psi. When the calculated amount of hydrogen has been absorbed remove the catalyst and evaporate to dryness under reduced pressure, recrystallize the residue from aqueous MeOH to give dimethyl 4-aminophthalate (mp 83–84E). Stir 146.3 g dimethyl 4-aminophthalate in 1.4 L acetic anhydride at 60–70E for 2 h then leave over-

night, precipitate product with MeOH. Dry the product and rinse it with sodium carbonate solution, re-dry, recrystallize from benzene/MeOH (Caution! Benzene is a carcinogen!) to give dimethyl 4-acetamidophthalate (mp 138-140E). Add 100.4 g to 600 mL fuming (90%) nitric acid at 0-5E over 30 min, stir at 5-10E for 2.5 h, mix the reaction mixture with 800 mL cold dichloromethane, shake with crushed ice. Remove the organic layer and extract the aqueous layer with 200 mL cold dichloromethane. Combine the organic layers and wash them with ice water, cold sodium bicarbonate solution, and cold water. Dry over anhydrous magnesium sulfate, evaporate to dryness under reduced pressure and, recrystallize repeatedly from MeOH to give dimethyl 4-acetamido-5-nitrophthalate (mp 123-124.5E). Hydrolyze dimethyl 4-acetamido-5-nitrophthalate to dimethyl 4-amino-5-nitrophthalate. Hydrogenate 20.3 g dimethyl 4-amino-5-nitrophthalate in 250 mL MeOH over 1 g 5% platinum on carbon at an initial hydrogen pressure of 50 psi, remove the catalyst, evaporate to dryness under reduced pressure at 25E, recrystallize from chloroform/dichloromethane to give dimethyl 4,5-diaminophthalate (mp 111.5-113E). Add 1.1 g dimethyl 4,5-diaminophthalate to 3 mL hydrazine hydrate (Caution! Hydrazine hydrate is a carcinogen!) and 3 mL triethylamine in 20 mL MeOH, concentrate the resulting solution, triturate with benzene/MeOH, recrystallize from N,N'-dimethylacetamide/acetic acid to give 4,5-diaminophthalhydrazide (6,7-diamino-2,3-dihydrophthalazine-1,4-dione) (mp 407E) (J. Heterocycl. Chem. 1973, 10, 891), mix 4,5-diaminophthalhydrazide with a small amount of concentrated HCl, recrystallize from EtOH to give 4,5-diaminophthalhydrazide dihydrochloride.)

HPLC VARIABLES

Column: 250 × 4.6 5 µm TSKgel ODS-120T (Tosoh, Japan)

Mobile phase: MeCN:tetrahydrofuran:100 mM pH 7.0 phosphate buffer 24:3:73

Flow rate: 1.0

Injection volume: 20

Detector: Chemiluminiscence following post-column reaction. The column effluent mixed with 20 mM hydrogen peroxide in water pumped at 1.0 mL/min and then with 30 mM potassium hexacyanoferrate(III) in 3.0 M NaOH pumped at 2.0 mL/min and flowed to the detector.

CHROMATOGRAM

Retention time: 52

Internal standard: beclomethasone

OTHER SUBSTANCES

Extracted: dexamethasone

Simultaneous: aldosterone, corticosterone, cortisone, 11-deoxycortisol, hydrocortisone, 18-hydroxycorticosterone, 18-hydroxydeoxycorticosterone, prednisolone, prednisone

Noninterfering: androstendione, cholesterol, estrone, estradiol, estriol, pregnenolone, progesterone

Interfering: betamethasone.

KEY WORDS

plasma; derivatization; beclomethasone is IS

REFERENCE

Ishida,J.; Sonezaki,S.; Yamaguchi,M.; Yoshitake,T., Determination of dexamethasone in plasma by high-performance liquid chromatography with chemiluminiscence detection, *Anal.Sci.*, **1993**, 9, 319-322.

SAMPLE

Matrix: blood

Sample preparation: Mix 1 mL plasma with 750 µL 40 µg/mL IS in EtOH, add 8 mL dichloromethane, extract on a roller mixer for 30 min. Centrifuge at 2500 rpm for 10 min at 25°, collect the organic layer, evaporate to dryness at 30° under a stream of nitrogen. Reconstitute the residue in 1 mL mobile phase, centrifuge at 15000 rpm for 2 min, inject a 100 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 Alltima C18 (Alltech)

Mobile phase: MeCN:MeOH:water:glacial acetic acid 8.8:65:26.2:0.175

Flow rate: 1.3

Injection volume: 100

Detector: UV 242

CHROMATOGRAM

Retention time: 15.4 (dipropionate)

Internal standard: dexamethasone-21-acetate (5.2)

OTHER SUBSTANCES

Extracted: degradation products

KEY WORDS

plasma; serum

REFERENCE

Foe,K.; Brown,K.F.; Seale,J.P. Decomposition of beclomethasone propionate esters in human plasma, *Biopharm.Drug Dispos.*, **1998**, *19*, 1–8.

SAMPLE

Matrix: blood

Sample preparation: 50 µL Plasma + 100 µL 20 µg/mL cloprednol + 3 mL ether, shake 10 min, centrifuge at 3000 g, remove the organic phase and evaporate it to dryness under nitrogen. Take up the residue in 200 µL mobile phase, inject a 50 µL aliquot.

HPLC VARIABLES

Column: Nucleosil R 10 C 18

Mobile phase: MeOH:MeCN:water:acetic acid 400:100:200:1

Injection volume: 50

Detector: UV 254

CHROMATOGRAM

Internal standard: cloprednol

Limit of detection: 500 ng/mL

KEY WORDS

for beclomethasone dipropionate; plasma

REFERENCE

Würthwein,G.; Rohdewald,P. Activation of beclomethasone dipropionate by hydrolysis to beclomethasone-17-monopropionate, *Biopharm.Drug Dispos.*, **1990**, *11*, 381–394.

SAMPLE

Matrix: blood, tissue

Sample preparation: Acidify plasma or lung tissue homogenate to pH 2 with 500 mM HCl, add 100 µL 20 µg/mL IS, extract with 8 mL dichloromethane. Evaporate the organic layer to dryness under vacuum, reconstitute in 120 µL MeOH:5% acetic acid 50:50, inject an 80 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Zorbax ODS C18

Mobile phase: MeCN:MeOH:water 44:11:45

Flow rate: 1

Injection volume: 80

Detector: UV 242 or radioactivity

CHROMATOGRAM

Internal standard: hydrocortisone 21-S-propionate (JO 498)

OTHER SUBSTANCES

Extracted: metabolites, budesonide

KEY WORDS

for beclomethasone dipropionate; plasma; rat; lung; radiolabeled; pharmacokinetics

REFERENCE

Chanoine,F.; Grenot,C.; Heidmann,P.; Junien,J.L. Pharmacokinetics of butixocort 21-propionate, budesonide, and beclomethasone dipropionate in the rat after intratracheal, intravenous, and oral treatments, *Drug Metab.Dispos.*, **1991**, 19, 546–553.

SAMPLE

Matrix: ileostomy effluent

Sample preparation: Dilute ileostomy effluent 1:2 by weight with water and mix with 100 μL 11 $\mu\text{g/mL}$ 17-hydroxyprogesterone. Extract 3 g aliquot three times with 10 mL dichloromethane by shaking for 1 min and centrifuging at 2000 rpm for 2 min. Wash combined extracts successively with 2 mL 0.1 M NaOH and 4 mL water by shaking for 30 s and centrifuging for 1 min then dry the organic layer under air at 40°. Take up the extract in 1 mL MeOH, add 1.1 mL water and apply to C18 Bond Elut SPE cartridge. Wash with 10 mL water, wash with 5 mL MeOH:water 45:55, elute with 2 mL MeOH. Add 50 μL 20 $\mu\text{g/mL}$ progesterone to the eluate, dry at 40°, take up in 100 μL MeOH, inject 10 μL aliquot.

HPLC VARIABLES

Guard column: Bondapak C18/Corasil

Column: 300 \times 3.9 $\mu\text{Bondapak C18}$

Mobile phase: MeOH:50 mM pH 3.0 sodium phosphate buffer 55:45

Flow rate: 3

Injection volume: 10

Detector: UV 254 and 238

CHROMATOGRAM

Retention time: 21.3 (beclomethasone dipropionate)

Internal standard: 17-Hydroxyprogesterone (6.0) and progesterone (11.6)

OTHER SUBSTANCES

Extracted: beclomethasone alcohol, beclomethasone 17-monopropionate

KEY WORDS

SPE

REFERENCE

Levine,D.S.; Raisys,V.A.; Ainardi,V. Coating of oral beclomethasone dipropionate capsules with cellulose acetate phthalate enhances delivery of topically active antiinflammatory drug to the terminal ileum, *Gastroenterology*, **1987**, 92, 1037–1044.

SAMPLE

Matrix: tissue

Sample preparation: 100 mg Tissue + 2 mL Ringer's pH 6.8 phosphate buffer + 2 mL EtOH, centrifuge, wash residue twice. Pool supernatant and washings and evaporate to dryness, take up in 400 μL EtOH, inject an aliquot.

HPLC VARIABLES

Guard column: Used but not specified

Column: μ Bondapak C18

Mobile phase: Gradient. MeOH:water 40:60 to 80:20, time not specified

Flow rate: 1.5

Detector: UV 254

OTHER SUBSTANCES

Extracted: beclomethasone monopropionate, beclomethasone, cyclomethasone

KEY WORDS

for beclomethasone dipropionate; lung

REFERENCE

Ronca-Testoni, S. Hydrolysis of cyclomethasone by the human lung, *Int.J.Clin.Pharmacol.Res.*, **1983**, 3, 17-20.


Befunolol

Molecular formula: C₁₆H₂₁NO₄

Molecular weight: 291.35

CAS Registry No.: 39552-01-7

Merck Index: 1050



SAMPLE

Matrix: perfusate

Sample preparation: 50 μ L Perfusate + 50 μ L pH 7.4 phosphate-buffered saline or 100 mM HCl + 100 μ L 50 μ g/mL salicyl methionine in MeOH, centrifuge at 12000 g for 10 min, inject a 50 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 150 \times 4.6 Cosmosil 5C18-P (Nacalai Tesque)

Mobile phase: MeOH:50 mM NaH₂PO₄ 45:55

Flow rate: 1

Injection volume: 50

Detector: F ex 300 em 500

CHROMATOGRAM

Internal standard: salicyl methionine

KEY WORDS

rabbit

REFERENCE

Sasaki,H.; Igarashi,Y.; Nagano,T.; Nishida,K.; Nakamura,J. Different effects of absorption promoters on corneal and conjunctival penetration of ophthalmic β -blockers, *Pharm.Res.*, **1995**, *12*, 1146–1150.

SAMPLE

Matrix: solutions

Sample preparation: 50 μ L Solution + 50 μ L pH 7.4 PBS + 100 μ L 50 μ g/mL salicylmethionine in MeOH, centrifuge at 12000 g for 10 min, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 Cosmosil 5C18-P (Nacalai Tesque)

Mobile phase: MeOH:50 mM NaH₂PO₄ 45:55

Flow rate: 1

Injection volume: 50

Detector: F ex 300 em 500

CHROMATOGRAM

Internal standard: salicylmethionine

KEY WORDS

buffer; earle's balanced salt solution

REFERENCE

Sasaki,H.; Igarishi,Y.; Nishida,K.; Nakamura,J. Intestinal permeability of ophthalmic β -blockers for predicting ocular permeability, *J.Pharm.Sci.*, **1994**, *83*, 1335–1338.

Benactyzine


Molecular formula: $C_{20}H_{25}NO_3$

Molecular weight: 327.42

CAS Registry No.: 302-40-9, 57-37-4 (HCl)

Merck Index: 1055

Lednicer No.: 1 93



SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 $\mu\text{g/mL}$ solution in MeOH, inject a 20 μL aliquot.

HPLC VARIABLES

Column: 125 \times 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 2.3

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benperidol, benzethidine, benzocaine, benzoctamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipanone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscaphine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl,

protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranylcypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleppamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J.Chromatogr.*, **1985**, *323*, 191–225.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 1 mg/mL solution in MeOH, inject a 5 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Lichrosphere cyanopropyl

Mobile phase: Carbon dioxide:MeOH:isopropylamine 90:10:0.05

Column temperature: 50

Flow rate: 3

Injection volume: 5

Detector: UV 220

CHROMATOGRAM

Retention time: 1.99

OTHER SUBSTANCES

Simultaneous: buclizine, hydroxyzine, perphenazine, thioridazine, amitriptyline, desipramine, imipramine, nortriptyline, protriptyline

KEY WORDS

SFC; pressure 200 bar

REFERENCE

Berger, T.A.; Wilson, W.H. Separation of drugs by packed column supercritical fluid chromatography. 2. Antidepressants, *J.Pharm.Sci.*, **1994**, *83*, 287–290.

Benazepril


Molecular formula: C₂₄H₂₈N₂O₅

Molecular weight: 424.50

CAS Registry No.: 86541-75-5, 86541-74-4 (HCl)

Merck Index: 1058

Lednicer No.: 5 135



SAMPLE

Matrix: formulations

Sample preparation: Add MeOH:water 50:50 to powdered capsules or tablets so as to give a benazepril concentration of ca. 20 µg/mL, stir for 15 min, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.5 5 µm Hypersil ODS

Mobile phase: MeCN:THF:20 mM pH 2.5 sodium heptanesulfonate 45.6:2.4:52

Flow rate: 1

Injection volume: 20

Detector: UV 215

CHROMATOGRAM

Retention time: 15.5

OTHER SUBSTANCES

Simultaneous: ramipril

KEY WORDS

capsules; tablets

REFERENCE

Bonazzi,D.; Gotti,R.; Andrisano,V.; Cavrini,V. Analysis of ACE inhibitors in pharmaceutical dosage forms by derivative UV spectroscopy and liquid chromatography (HPLC), *J.Pharm.Biomed.Anal.*, **1997**, 16, 431-438.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A: B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 205.2

CHROMATOGRAM

Retention time: 17.003

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

SAMPLE

Matrix: formulations

Sample preparation: Grind tablets to a fine powder. Weigh out an amount equivalent to 25 mg benazepril, extract with MeOH, filter. Mix 100-500 μ L filtrate with 200 μ L 4 mg/mL IS in MeOH, make up to 10 mL with MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4 10 μ m LiChrosorb RP-18

Mobile phase: MeCN:buffer 30:70 (Buffer was 67 mM KH_2PO_4 adjusted to pH 2.4 with phosphoric acid.)

Flow rate: 1

Injection volume: 20

Detector: UV 211

CHROMATOGRAM

Retention time: 18.53

Internal standard: enalapril (8.58)

Limit of detection: 5 μ g/mL

Limit of quantitation: 10 μ g/mL

OTHER SUBSTANCES

Simultaneous: cilazapril


KEY WORDS

tablets

REFERENCE

Gumieniczek,A.; Przyborowski,L. Determination of benazepril and cilazapril in pharmaceuticals by high performance liquid chromatography, *J.Liq.Chromatogr.Rel.Technol.*, **1997**, 20, 2135-2142.

Bendroflumethiazide



Molecular formula: C₁₅H₁₄F₃N₃O₄S₂

Molecular weight: 421.42

CAS Registry No.: 73-48-3

Merck Index: 1064

Lednicer No.: 2 358

SAMPLE

Matrix: blood

Sample preparation: 3 mL Plasma + 2 mL 10 mM NaOH + 2 mL 10 mM HCl, mix, add 10 mL diethyl ether, shake gently on a platform shaker for 15 min, centrifuge at -10° at 2200 g for 15 min, freeze in dry ice/acetone. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 70° for 45 min (these conditions are required to remove residual benzyl alcohol that is present as a preservative in the heparin), reconstitute the residue in 50 µL 10 mM NaOH, vortex for 25 s, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 300 × 3.9 10 µm µBondapak C18

Mobile phase: Isopropanol:water:acetic acid 17:82:1

Flow rate: 2

Injection volume: 20

Detector: UV 269

CHROMATOGRAM

Retention time: 14.4

Internal standard: bendroflumethiazide

OTHER SUBSTANCES

Extracted: trichlormethiazide

KEY WORDS

plasma; silanize glassware; bendroflumethiazide is IS

REFERENCE

Meyer, M.C.; Hwang, P.T.R. Determination of trichlormethiazide in human plasma and urine by high-performance liquid chromatography, *J. Chromatogr.*, **1981**, 223, 466–472.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A: B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 208.7

CHROMATOGRAM

Retention time: 18.632

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

SAMPLE

Matrix: bulk

Sample preparation: Dissolve in solvent, inject an aliquot. (Solvent was 750 mg KCl in 10 mL 1 M HCl, add 400 mL water, add 400 mL MeOH, make up to 1 L with water.)

HPLC VARIABLES

Guard column: 5 × 4 7 µm Nucleosil-100 phenyl

Column: 300 × 4 7 µm Nucleosil-100 phenyl

Mobile phase: MeOH:water 40:60

Column temperature: 35

Flow rate: 1.5

Injection volume: 50

Detector: UV 270

CHROMATOGRAM

Retention time: 10.5

OTHER SUBSTANCES

Simultaneous: hydroflumethiazide, degradation products

REFERENCE

Frontini, R.; Mielck, J. B. Determination and quantitation of bendroflumethiazide and its degradation products using HPLC, *J. Liq. Chromatogr.*, **1992**, 15, 2519-2528.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Chirex 3001 (Phenomenex)

Mobile phase: Hexane:1,2-dichloroethane:EtOH/trifluoroacetic acid 55:35:10 (EtOH/trifluoroacetic acid was premixed 20:1.)

Flow rate: 1

Injection volume: 20

Detector: UV 272

CHROMATOGRAM

Retention time: 19, 21.5 (enantiomers)

KEY WORDS

chiral

REFERENCE

Cleveland,T. Pirkle-concept chiral stationary phases for the HPLC separation of pharmaceutical racemates, *J.Liq.Chromatogr.*, **1995**, *18*, 649–671.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a solution in MeOH:water 80:20, inject a 6 μ L aliquot.

HPLC VARIABLES

Guard column: 5 \times 4 10 μ m LiChrosorb RP-8

Column: 100 \times 4.6 5 μ m Spheri RP-18 (Brownlee)

Mobile phase: MeOH:water 80:20 containing 2 g/L lithium perchlorate

Flow rate: 0.5

Injection volume: 6

Detector: E, ESA Model 5100A Coulochem, model 5020 guard cell +950 mV, Model 5010 analytical cell + 400 mV, palladium reference electrode, following post-column photolysis. The effluent from the column flowed through a 20 m \times 0.3 mm coil of PTFE tubing and was irradiated at 254 nm with a Sylvania GTE 8 W low-pressure lamp.

CHROMATOGRAM

Retention time: 5

Limit of detection: 267 ng/mL

OTHER SUBSTANCES

Also analyzed: butizide, chlorthalidone, ethacrynic acid, furosemide, hydrochlorothiazide

KEY WORDS

post-column reaction

REFERENCE

Macher,M.; Wintersteiger,R. Improved electrochemical detection of diuretics in high-performance liquid chromatographic analysis by postcolumn on-line photolysis, *J.Chromatogr.A*, **1995**, *709*, 257–264.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 150 \times 4.6 cellulose 3,5-dimethylphenylcarbamate/10-undecenoate bonded to allylsilica

Mobile phase: Heptane:isopropanol:diethylamine 80:20:0.1

Flow rate: 1

Injection volume: 1000

Detector: UV 254

CHROMATOGRAM

Retention time: k' 8.34

KEY WORDS

chiral; α 1.17

REFERENCE

Oliveros,L.; Lopez,P.; Minguillon,C.; Franco,P. Chiral chromatographic discrimination ability of a cellulose 3,5-dimethylphenylcarbamate/10-undecenoate mixed derivative fixed on several chromatographic matrices, *J.Liq.Chromatogr.*, **1995**, *18*, 1521–1532.

SAMPLE**Matrix:** urine**Sample preparation:** 2 mL Urine + 500 mg solid sodium bicarbonate, mix, add 2 mL 10 mM NaOH, mix, add 10 mL diethyl ether, shake gently on a platform shaker for 15 min, centrifuge at -10° at 2200 g for 15 min, freeze in dry ice/acetone. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 70° for 45 min (these conditions are required to remove residual benzyl alcohol that is present as a preservative in the heparin), reconstitute the residue in 100 μ L MeOH, vortex for 25 s, inject a 20 μ L aliquot.

HPLC VARIABLES**Column:** 300 \times 3.9 10 μ m μ Bondapak C18**Mobile phase:** MeCN:MeOH:water:acetic acid 5:35:59:1**Flow rate:** 2**Injection volume:** 20**Detector:** UV 280

CHROMATOGRAM**Retention time:** 10.7**Internal standard:** bendroflumethiazide

OTHER SUBSTANCES**Extracted:** trichlormethiazide

KEY WORDSsilanize glassware; bendroflumethiazide is IS

REFERENCEMeyer, M.C.; Hwang, P.T.R. Determination of trichlormethiazide in human plasma and urine by high-performance liquid chromatography, *J. Chromatogr.*, **1981**, 223, 466–472.

SAMPLE**Matrix:** urine**Sample preparation:** 2 mL Urine + 2 mL 1 M pH 4.1 NaH_2PO_4 + 4 mL ethyl acetate, vortex for 2 min, centrifuge at 1500 g for 5 min. Remove the organic phase and add it to 5 mL 100 mM pH 7.5 Na_2HPO_4 , vortex for 2 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 60° , reconstitute the residue in 100 μ L MeCN:10 mM pH 3.0 phosphate buffer, inject a 5 μ L aliquot.

HPLC VARIABLES**Column:** 125 \times 4 5 μ m LiChrosorb RP-18**Mobile phase:** Gradient. MeCN:10 mM pH 3.0 phosphate buffer 10:90 for 1.5 min then to 35:65 over 2 min**Column temperature:** 50**Flow rate:** 1.5**Injection volume:** 5**Detector:** UV 271

CHROMATOGRAM**Retention time:** 6.8**Limit of quantitation:** 500 ng/mL

OTHER SUBSTANCES**Extracted:** chlorothiazide, hydrochlorothiazide, quinethazone, chlorthalidone, clopamide, methyclothiazide, furosemide, metolazone, mefruside, cyclopenthiazide, bumetanide**Simultaneous:** indapamide, clorexolone, ethacrynic acid**Noninterfering:** aspirin, albuterol, allopurinol, alprenolol, atenolol, captopril, carbimazole, clonidine, coloxyl, danthron, diazepam, digoxin, doxepin, glibenclamide, hydralazine, in-

domethacin, labetalol, metformin, methyl dopa, metoprolol, mianserin, minoxidil, nifedipine, nitrazepam, oxazepam, oxprenolol, pindolol, prazosin, propranolol, senokot, theophylline, trifluoperazine

REFERENCE

Fullinlaw, R.O.; Bury, R.W.; Moulds, R.F.W. Liquid chromatographic screening of diuretics in urine, *J. Chromatogr.*, **1987**, *415*, 347–356.

SAMPLE

Matrix: urine

Sample preparation: 2 mL Urine + 0.5 g solid buffer I (pH 5-5.5), vortex 15 s, add 4 mL ethyl acetate, agitate for 10 min, centrifuge at 600 g for 5 min. Remove organic layer and vortex it with 2 mL 5% aqueous lead acetate for 10 s, centrifuge at 600 g for 5 min, remove and keep organic phase. 2 mL Urine + 0.5 g solid buffer II (pH 9-9.5), vortex 15 s, add 4 mL ethyl acetate, agitate for 10 min, centrifuge at 600 g for 5 min. Remove organic layer and combine it with previous organic layer. Evaporate to dryness at 50° under a stream of nitrogen, reconstitute in 300 μ L 50 μ g/mL β -hydroxyethyltheophylline in MeOH, inject 5 μ L aliquot. (Solid buffer I was $\text{KH}_2\text{PO}_4\text{:Na}_2\text{HPO}_4$ 99:1, solid buffer II was $\text{NaHCO}_3\text{:K}_2\text{CO}_3$ 3:2.)

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m HP Hypersil ODS (A) or HP LiChrosorb RP-18 (B)

Mobile phase: Gradient. MeCN:buffer from 15:85 at 2 min to 80:20 at 20 min (Buffer was 50 mM NaH_2PO_4 containing 16 mM propylamine hydrochloride, adjusted to pH 3 with concentrated phosphoric acid.)

Flow rate: 1

Injection volume: 5

Detector: UV 230, UV 275

CHROMATOGRAM

Retention time: 15 (A), 15.4 (B)

Internal standard: β -hydroxyethyltheophylline (3.7 (A), 4.4 (B))

Limit of detection: 1000 ng/mL

OTHER SUBSTANCES

Extracted: furosemide, metolazone, amiloride, acetazolamide, chlorothiazide, hydrochlorothiazide, quinethazone, triamterene, hydroflumethiazide, chlorthalidone, dichlorophenamide, trichloromethiazide, methyclothiazide, benzthiazide, cyclothiazide, ethacrynic acid, bumetanide, probenecid, spironolactone, canrenone, flumethiazide

Noninterfering: acetaminophen, aspirin, caffeine, diflunisal, fenoprofen, ibuprofen, indomethacin, methocarbamol, naproxen, phenylbutazone, sulindac, tetracycline, theobromine, theophylline, tolmetin, trimethoprim, verapamil

Interfering: polythiazide

REFERENCE

Cooper, S.F.; Massé, R.; Dugal, R. Comprehensive screening procedure for diuretics in urine by high-performance liquid chromatography, *J. Chromatogr.*, **1989**, *489*, 65–88.

SAMPLE

Matrix: urine

Sample preparation: 2 mL Urine + 1 mL 10 mM HCl + 2000 ng bendroflumethiazide, extract with 5 mL ethyl acetate, centrifuge at 3000 rpm for 5 min. Remove the organic layer and dry it under a stream of nitrogen at 40°. Reconstitute with 100 μ L MeOH, inject a 2 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 2.1 5 μ m Hypersil ODS

Mobile phase: Gradient. MeOH: 50 mM ammonium acetate from 10:90 to 60:40 over 10 min, maintain at 60:40 for 10 min.

Column temperature: 40

Flow rate: 0.3

Injection volume: 2

Detector: F ex 223 em 415 or UV 230

CHROMATOGRAM

Retention time: 8.6

Internal standard: bendroflumethiazide

OTHER SUBSTANCES

Extracted: bumetanide (F ex 231 em 426 or UV), furosemide (UV), piretanide (UV), cyclopenthiazide (UV), etozolin (UV), canrenone (UV)

KEY WORDS

bendroflumethiazide is IS

REFERENCE

Gradeen,C.Y.; Billay,D.M.; Chan,S.C. Analysis of bumetanide in human urine by high-performance liquid chromatography with fluorescence detection and gas chromatography/mass spectrometry, *J.Anal.Toxicol.*, **1990**, *14*, 123–126.

SAMPLE

Matrix: urine

Sample preparation: Make 5 mL urine alkaline (pH 9-10), add 2 g NaCl, extract twice with 6 mL ethyl acetate. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 200 μ L MeCN/water, inject a 10-20 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 4.5 μ m SGE 100 GL-4 C18P (Scientific Glass Engineering)

Mobile phase: MeCN:MeOH:water:trifluoroacetic acid 4.5:10.5:85:0.5

Flow rate: 0.8 or 1

Injection volume: 10-20

Detector: MS, ZAB2-SEQ (VG), PSP source coupled to LC, source 250°, probe 240-260°, scan m/z 200-550 or UV 270

CHROMATOGRAM

Retention time: 9.2

OTHER SUBSTANCES

Extracted: amiloride, chlorthalidone, triamterene, furosemide, benzthiazide

REFERENCE

Ventura,R.; Fraisse,D.; Becchi,M.; Paisse,O.; Segura,J. Approach to the analysis of diuretics and masking agents by high-performance liquid chromatography-mass spectrometry in doping control, *J.Chromatogr.*, **1991**, *562*, 723–736.

SAMPLE

Matrix: urine

Sample preparation: Direct injection.

HPLC VARIABLES

Guard column: 35 \times 4.5 μ m Spherisorb ODS-2

Column: 120 \times 4.5 μ m Spherisorb ODS-2

Mobile phase: MeOH:50 mM sodium dodecyl sulfate 5:95

Column temperature: 50

Flow rate: 1

Injection volume: 20

Detector: UV 224

CHROMATOGRAM

Retention time: 15.2

Limit of detection: 500 ng/mL

OTHER SUBSTANCES

Simultaneous: chlorthalidone

REFERENCE

Bonet Domingo, E.; Medina Hernández, M.J.; Ramis Ramos, G.; Garcia Alvarez-Coque, M.C. High-performance liquid chromatographic determination of diuretics in urine by micellar liquid chromatography, *J. Chromatogr.*, **1992**, 582, 189–194.

SAMPLE

Matrix: urine

Sample preparation: Buffer urine to 4.9 by mixing with an equal volume of pH 4.9 200 mM sodium phosphate buffer. Inject a 40 μ L aliquot onto column A with mobile phase A, after 3 min backflush the contents of column A onto column B with mobile phase B and start the gradient. At the end of the run re-equilibrate for 10 min.

HPLC VARIABLES

Column: A 20 \times 4.5 μ m Hypersil octadecylsilica ODS; B 200 \times 4.6 5 μ m Shiseido SG-120 polymer-based C18

Mobile phase: A water; B Gradient. MeCN:buffer from 7:93 to 15:85 over 3.5 min, to 50:50 over 8.5 min, maintain at 50:50 for 11 min (Buffer was 6.9 g $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$ in 1 L water, pH adjusted to 3.1 with phosphoric acid.)

Flow rate: 1

Injection volume: 40

Detector: UV 270

CHROMATOGRAM

Retention time: 18.8

Limit of detection: 500 ng/mL

OTHER SUBSTANCES

Extracted: acetazolamide, amiloride, benzthiazide, bumetanide, caffeine, carbamazepine, chlorothiazide, chlorthalidone, clopamide, dichlorfenamide, ethacrynic acid, furosemide, hydrochlorothiazide, metyrapone, probenecid, spironolactone, triamterene, trichlormethiazide

KEY WORDS

column-switching; optimum detection wavelengths vary for each drug

REFERENCE

Saarenin, M.; Sirén, H.; Riekkola, M.-L. A column switching technique for the screening of diuretics in urine by high performance liquid chromatography, *J. Liq. Chromatogr.*, **1993**, 16, 4063–4078.

SAMPLE

Matrix: urine

Sample preparation: 5 mL Urine + 50 μ L 100 μ g/mL 7-propyltheophylline in MeOH + 200 μ L ammonium chloride buffer + 2 g NaCl, extract with 6 mL ethyl acetate by rocking at 40 movements/min for 20 min and centrifuging at 800 g for 5 min, repeat extraction, combine organic layers, evaporate to dryness at 40° under a stream of nitrogen. Recon-

stitute in 200 μL MeCN:water 15:85 and inject 20 μL aliquots. (Ammonium chloride buffer was 28 g ammonium chloride in 100 mL water with the pH adjusted to 9.5 with concentrated ammonia solution.)

HPLC VARIABLES

Column: 75 \times 4.6 3 μm Ultrasphere ODS

Mobile phase: Gradient. MeCN:100 mM ammonium acetate adjusted to pH 3 with concentrated phosphoric acid. From 10:90 to 15:85 over 2 min to 55:45 over 3 min to 60:40 over 3 min. Kept at 60:40 for 1 min, decreased to 10:90 over 1 min and equilibrated at 10:90 for 2 min.

Flow rate: 1

Injection volume: 20

Detector: UV 270

CHROMATOGRAM

Retention time: 7.0

Internal standard: 7-propyltheophylline (4.5)

Limit of detection: 50 ng/mL

OTHER SUBSTANCES

Simultaneous: xipamide, bumetanide, acetazolamide, amiloride, benzthiazide, buthiazide, caffeine, canrenone, chlorthalidone, clopamide, cyclothiazide, diclofenamide, furosemide, hydrochlorothiazide, mesocarb, morazone, piretanide, probenecid, spironolactone, torsemide, triamterene

Interfering: polythiazide, ethacrynic acid

REFERENCE

Ventura,R.; Nadal,T.; Alcalde,P.; Pascual,J.A.; Segura,J. Fast screening method for diuretics, probenecid and other compounds of doping interest, *J.Chromatogr.A*, **1993**, 655, 233–242.

SAMPLE

Matrix: urine

Sample preparation: Direct injection into column A with mobile phase A for 1 min then back flush onto column B with mobile phase B.

HPLC VARIABLES

Column: A 20 \times 2.1 30 μm Hypersil ODS-C18; B 250 \times 4 5 μm Hypersil ODS-C18

Mobile phase: A Water; B Gradient. MeCN:buffer 15:85 for 1.5 min then to 80:20 over 8 min. Keep at 80:20 for 2.5 min then re-equilibrate with 15:85. (Buffer was 50 mM NaH_2PO_4 + 1.4 mL propylamine hydrochloride per liter adjusted to pH 3 with concentrated phosphoric acid.)

Flow rate: 1

Injection volume: 50

Detector: UV 230

CHROMATOGRAM

Retention time: 9.8

Limit of detection: 20 ng/mL.

OTHER SUBSTANCES

Simultaneous: bumetanide, ethacrynic acid, acetazolamide, amiloride, chlorthalidone, cyclothiazide, furosemide, hydrochlorothiazide, probenecid, spironolactone, triamterene

REFERENCE

Campíns-Falco,P.; Herráez-Hernández,R.; Sevillano-Cabeza,A. Column-switching techniques for screening of diuretics and probenecid in urine samples, *Anal.Chem.*, **1994**, 66, 244–248.


Benfluorex

Molecular formula: C₁₉H₂₀F₃NO₂

Molecular weight: 351.37

CAS Registry No.: 23602-78-0, 23642-66-2 (HCl)

Merck Index: 1066



SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A: B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 16.378

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.


Benidipine

Molecular formula: C₂₈H₃₁N₃O₆

Molecular weight: 505.57

CAS Registry No.: 105979-17-7

Merck Index: 1071



SAMPLE

Matrix: bile, blood, urine

Sample preparation: Urine, bile. Inject urine and bile directly. Plasma. Add an equal volume of MeCN, centrifuge, remove supernatant and evaporate it to dryness. Reconstitute residue in MeCN:10 mM pH 5 ammonium acetate 20:80.

HPLC VARIABLES

Column: 300 × 8 Unisil Q5C18 (Gaschro Kogyo)

Mobile phase: Gradient. A was MeCN:10 mM pH 5 ammonium acetate 20:80. B was MeCN:10 mM pH 5 ammonium acetate 80:20. A:B from 100:0 to 92:8 over 4 min, to 83:17 over 4 min, to 75:25 over 4 min, to 63:37 over 4 min, to 55:45 over 8 min, to 50:50 over 8 min, to 42:58 over 4 min, to 0:100 over 9 min, hold at 0:100 for 25 min.

Flow rate: 2 for 36 min, then 3

Detector: Radioactivity

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; semi-preparative; rat; dog; radiolabelled

REFERENCE

Kobayashi,H.; Okumura,S.; Kosaka,Y.; Kobayashi,S.; Inoue,A.; Oka,T.; Nakamizo,N. Identification of benidipine hydrochloride metabolites in rats and dogs, *Arzneimittelforschung*, **1988**, *38*, 1753–1756.

SAMPLE

Matrix: bulk

Sample preparation: Dissolve 50 mg in 50 mL MeOH:water 70:30, remove a 10 mL aliquot and add it to 10 mL 1 mg/mL diphenylamine in MeOH, make this mixture up to 100 mL with MeOH:water 70:30, inject a 10 µL aliquot.

HPLC VARIABLES

Column: 300 × 3.9 µBondapak C18

Mobile phase: MeOH:100 mM pH 5 ammonium acetate 70:30 (MeOH:50 mM pH 3.5 phosphate buffer:diisopropylamine 60:40:1 at 0.7 mL/min to detect other diastereomer)

Flow rate: 1

Injection volume: 10

Detector: UV 236

CHROMATOGRAM

Retention time: 10

Internal standard: diphenylamine (6)

OTHER SUBSTANCES

Simultaneous: impurities

REFERENCE

Suzuki,H.; Ono,E.; Ueno,H.; Takemoto,Y.; Nakamizo,N. Physico-chemical properties and stabilities of the highly potent calcium antagonist benidipine hydrochloride, *Arzneimittelforschung*, **1988**, 38, 1671-1676.

SAMPLE

Matrix: formulations

Sample preparation: Weigh out solid dispersion or solid dispersion granules equivalent to 20 mg benidipine hydrochloride, add 200 mL water, 100 mM hydrochloric solution, or McIlvaine buffer (pH 4.0, 6.0), filter (0.2 μ m), dilute with MeOH, inject an aliquot.

HPLC VARIABLES

Column: C18 (YMC-Pack ODS-A)

Mobile phase: MeOH:THF:50 mM pH 3.0 phosphate buffer 27:8:65

Flow rate: 1

Detector: UV 237

CHROMATOGRAM

Internal standard: benzoïn

KEY WORDS

solid dispersions; solid dispersion granules

REFERENCE

Suzuki,H.; Miyamoto,N.; Masada,T.; Hayakawa,E.; Ito,K. Solid dispersions of benidipine hydrochloride. I. Preparations using different solvent systems and dissolution properties, *Chem.Pharm.Bull.*, **1996**, 44, 364-371.

SAMPLE

Matrix: solutions

Sample preparation: Direct injection of a MeOH solution containing 200-1000 ng.

HPLC VARIABLES

Column: 250 \times 4.6 Sumchiral OA-4500 (Sumika Chemical Analysis Service)

Mobile phase: n-Hexane:1,2-dichloroethane:MeOH:trifluoroacetic acid 250:140:10:1

Flow rate: 1

Detector: UV 254

CHROMATOGRAM

Retention time: 33 (+), 42(-) (α = 1.30)


KEY WORDS

chiral

REFERENCE

Ohkubo,T.; Uno,T.; Sugawara,K. Enantiomer separation of dihydropyridine derivative calcium antagonists by high-performance liquid chromatography with chiral stationary phases, *J.Chromatogr.A*, **1994**, 659, 467-471.

Benoxaprofen



Molecular formula: C₁₆H₁₂ClNO₃

Molecular weight: 301.73

CAS Registry No.: 67434-14-4

Merck Index: 1075

SAMPLE

Matrix: bile, blood

Sample preparation: 10 µL Plasma, or bile + 180 µL MeCN + 30 µL 100 µL/mL IS in DMSO + 30 µL water, vigorously mix. Centrifuge at 15000 g for 10 min at 4°, inject a 10 µL aliquot of the supernatant.

HPLC VARIABLES

Guard column: 10 × 4.6 SUMICHIRAL OA 2500 (Sumika Chemical Analysis Service, Japan)

Column: 250 × 4.6 SUMICHIRAL OA 2500 (Sumika Chemical Analysis Service, Japan)

Mobile phase: MeOH containing 40 mM ammonium acetate (plasma) or MeCN:MeOH: water 15:85:5 containing 10 mM ammonium acetate (bile)

Flow rate: 1.0

Injection volume: 10

Detector: F ex 315, em 365

CHROMATOGRAM

Retention time: 10.3 (R, plasma) 11.8 (S, plasma), 22.5 (R, bile) 26.0 (S, bile)

Internal standard: naproxen methyl ester (4.4) (plasma), (4.3) (bile)

Limit of detection: 10 pg/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; rat; chiral; pharmacokinetics

REFERENCE

Mohri,K.; Okada,K.; Benet,L.Z. Stereoselective metabolism of benoxaprofen in rats. Biliary excretion of benoxaprofen taurine conjugate and glucuronide, *Drug Metab.Dispos.*, **1998**, 26, 332–337.

SAMPLE

Matrix: bile, blood, urine

Sample preparation: 10 µL Plasma, bile, or urine + 180 µL MeCN + 30 µL 100 µL/mL IS in DMSO + 30 µL water, mix vigorously. Centrifuge at 15000 g for 10 min at 4°, inject a 10 µL aliquot of the supernatant.

HPLC VARIABLES

Guard column: 10 × 4.6 5 µm Capcell Pak C18 (Shiseido, Japan)

Column: 250 × 4.6 5 µm Capcell Pak C18 (Shiseido, Japan)

Mobile phase: MeCN:THF:10 mM tetrabutylammonium hydrogen sulfate buffer 35:35:100

Flow rate: 1.3

Injection volume: 10

Detector: F ex 315 em 365

CHROMATOGRAM

Retention time: 17.5

Internal standard: naproxen methyl ester (13.2)

Limit of detection: 10 pg/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; rat; pharmacokinetics

REFERENCE

Mohri,K.; Okada,K.; Benet,L.Z. Stereoselective metabolism of benoxaprofen in rats. Biliary excretion of benoxaprofen taurine conjugate and glucuronide, *Drug Metab.Dispos.*, **1998**, *26*, 332–337.