

Benzonatate

Molecular formula: C₃₀H₅₃NO₁₁

Molecular weight: 603.75

CAS Registry No.: 104-31-4

Merck Index: 1127



SAMPLE

Matrix: tissue

Sample preparation: Homogenize 10 g tissue in 40 mL water, make alkaline with NaOH, extract twice with ether. Combine extracts and evaporate them to dryness, reconstitute the residue in 1 mL EtOH, inject a 50 µL aliquot.

HPLC VARIABLES

Column: Reverse phase ODS

Mobile phase: MeCN:100 mM phosphoric acid 50:50

Column temperature: 40

Flow rate: 1

Injection volume: 50

Detector: UV 205

CHROMATOGRAM

Limit of detection: 200 ng/mL


KEY WORDS

liver; brain; blood; kidney

REFERENCE

Cohan, J.A.; Manning, T.J.; Lukash, L.; Long, C.; Ziminski, K.R.; Conradi, S.E. Two fatalities resulting from Tessalon (benzonatate), *Vet. Hum. Toxicol.*, **1986**, 28, 543-544.

Benzoyl peroxide



Molecular formula: C₁₄H₁₀O₄

Molecular weight: 242.23

CAS Registry No.: 94-36-0

Merck Index: 1149

SAMPLE

Matrix: blood

Sample preparation: Condition a 100 mg Bond Elut C18 SPE cartridge with 2 mL MeCN and 2 mL 50 mM pH 7.5 phosphate buffer. Add 1 mL serum to the SPE cartridge, wash with 500 µL 50 mM pH 7.5 phosphate buffer, elute with 1 mL MeCN, inject a 20 µL aliquot of the eluate.

HPLC VARIABLES

Column: Capcell Pak C-18 SG-120 (Shiseido)

Mobile phase: MeCN:water 50:50

Flow rate: 1.2

Injection volume: 20

Detector: UV 235

CHROMATOGRAM

Retention time: 2.4

OTHER SUBSTANCES

Extracted: methyl methacrylate, N,N-dimethyl-p-toluidine

KEY WORDS

serum; dental material; methyl methacrylate polymer; plexiglass; horse; SPE

REFERENCE

Shintani,H.; Tsuchiya,T.; Hata,Y.; Nakamura,A. Solid phase extraction and HPLC analysis of toxic components eluted from methyl methacrylate dental materials, *J.Anal.Toxicol.*, **1993**, 17, 73-78.

SAMPLE

Matrix: food

Sample preparation: Mix 1 g flour with 8.0 mL EtOH, shake vigorously for 1 min and extract by sonication for 15 min at room temperature. Centrifuge at 4000 rpm for 20 min, remove 4 mL supernatant, add 2 mL 100 mM KOH, mix in an ultrasonic water bath for 2 min at room temperature, dilute to 25 mL with water and set aside overnight. Filter (0.45 µm) and inject a 50 µL aliquot.

HPLC VARIABLES

Column: Dionex OmniPac PAX-100 anion exchange column

Mobile phase: MeOH:3 mM sodium carbonate 2:98

Flow rate: 1.0

Injection volume: 50

Detector: UV 222

CHROMATOGRAM

Retention time: 2.5

Limit of detection: 19 ng/mL

Limit of quantitation: 120 ng/mL

OTHER SUBSTANCES

Interfering: benzoic acid

KEY WORDS

wheat flour; derivatization

REFERENCE

Chen, Q.-C.; Mou, S.-F.; Hou, X.-P.; Ni, Z.-M. Determination of benzoyl peroxide in wheat flour by ion chromatography with precolumn derivatization, *J. Liq. Chromatogr. Rel. Technol.*, **1998**, *21*, 705–716.

SAMPLE

Matrix: plastic

Sample preparation: Leach with 20 mL water, MeOH, acetone, or THF at room temperature for 24 h, repeat. Combine solutions and add an equal volume of MeCN:water 50:50, filter (0.45 μm) the supernatant, inject a 20 μl aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 Capcell Pak SG-120 C18 (Shiseido)

Mobile phase: MeCN:water 48:52

Flow rate: 1

Injection volume: 20

Detector: UV 235

CHROMATOGRAM

Retention time: 29.6


OTHER SUBSTANCES

Simultaneous: methyl methacrylate, N,N-dimethyl-p-toluidine

REFERENCE

Shintani, H. HPLC analysis of toxic additives and residual monomer from dental plate, *J. Liq. Chromatogr.*, **1995**, *18*, 613–626.

Benzphetamine



Molecular formula: C₁₇H₂₁N

Molecular weight: 239.36

CAS Registry No.: 156-08-1, 5411-22-3 (HCl)

Merck Index: 1151

Lednicher No.: 1 70

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in MeOH at a concentration of 1 mg/mL, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 5 Spherisorb S5W

Mobile phase: MeOH:buffer 90:10 (Buffer was 94 mL 35% ammonia and 21.5 mL 70% nitric acid in 884 mL water, adjust the pH to 10.1 with ammonia.)

Flow rate: 2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 1.52

OTHER SUBSTANCES

Simultaneous: phendimetrazine, methylphenidate, phenelzine, epinephrine, pipradol, phenylpropanolamine, fencamfamin, chlorphentermine, norpseudoephedrine, phentermine, fenfluramine, methylenedioxymphetamine, amphetamine, normetanephine, 4-hydroxyamphetamine, bromo-STP, STP, prolintane, 2-phenethylamine, tyramine, trimethoxyamphetamine, phenylephrine, pseudoephedrine, ephedrine, methylephedrine, dimethylamphetamine, methamphetamine, mescaline, mephentermine, norpipanone, levallorphan, hydroxypethidine, normethadone, meperidine, dipipanone, diamorphine, pentazocine, acetylcodeine, monoacetylmorphine, thebacon, oxycodone, thebaine, norlevorphanol, methadone, benzylmorphine, ethylmorphine, morphine-N-oxide, codeine, codeine-N-oxide, morphine, ethoheptazine, morphine-3-glucuronide, pholcodeine, norpethidine, hydrocodone, dihydrocodeine, dihydromorphine, levorphanol, norcodeine, normorphine

Noninterfering: dopamine, levodopa, methyl dopa, methyl dopate, norepinephrine

Interfering: pemoline, diethylpropion, mazindol, tranlycypromine, caffeine, fenethyline, buprenorphine, dextromoramide, phenoperidine, fentanyl, etorphine, piritramide, noscapine, papaverine, naloxone, dextropropoxyphene, nalorphine, phenazocine

REFERENCE

Law,B.; Gill,R.; Moffat,A.C. High-performance liquid chromatography retention data for 84 basic drugs of forensic interest on a silica column using an aqueous methanol eluent, *J.Chromatogr.*, **1984**, 301, 165-172.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 μ g/mL solution in MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 125 \times 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2**Injection volume:** 20**Detector:** E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM**Retention time:** 1.9

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzoctamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazepine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipanone, diprenorphine, dipyridamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserin, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypropazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenylglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J.Chromatogr.*, **1985**, 323, 191-225.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlordiazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenoprofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephentyoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methypylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nyldrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidone, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sufadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxizole, sulfanilamide, sulfapyridine, sulfasoxizole, sulindac, tamoxifen, temazepam, testosterone, thiacaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, 1994, 18, 233-242.

Benzquinamide


Molecular formula: C₂₂H₃₂N₂O₅

Molecular weight: 404.51

CAS Registry No.: 63-12-7

Merck Index: 1154

Lednicer No.: 1 350



SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 125 × 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 1.2

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipanone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl,

protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranylcypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

- Jane, I.; McKinnon, A.; Flanagan, R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J.Chromatogr.*, **1985**, 323, 191-225.


Benzthiazide

Molecular formula: C₁₅H₁₄ClN₃O₄S₃

Molecular weight: 431.94

CAS Registry No.: 91-33-8

Merck Index: 1155



SAMPLE

Matrix: blood, feces, urine

Sample preparation: Plasma. 3 mL Plasma + 500 μ L 1 μ g/mL polythiazide in 10 mM NaOH + 1 mL 10 mM NaOH + 800 μ L 100 mM HCl + 10 mL dichloromethane, shake on a platform shaker for 20 min, centrifuge at -10° at 3000 rpm for 15 min, repeat extraction twice more. Combine all organic layers and evaporate them to dryness under a stream of nitrogen at 50°, reconstitute the residue in 50 μ L, vortex, inject whole amount. Urine. 5 mL Urine + 1 mL 2 μ g/mL polythiazide in 10 mM NaOH + 1 mL 10 mM NaOH + 2 mL 0.68% KH₂PO₄ adjusted to pH 6.1 + 10 mL dichloromethane, shake on a platform shaker for 20 min, centrifuge at -10° at 3000 rpm for 15 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 50°, reconstitute the residue in 100 μ L, vortex, inject a 20 μ L aliquot. Feces. Extract 20 g feces with 350 mL acetone, filter (Whatman No. 1 paper), evaporate to ca. 50 mL, make up to 100 mL with acetone. Remove a 10 mL aliquot and add it to 1 mL 300 μ g/mL polythiazide in acetone, evaporate to dryness under nitrogen, dissolve in 10 mL MeOH, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: μ Bondapak C18

Mobile phase: MeCN:glacial acetic acid:water 35:2:63 (plasma, feces) or MeCN:water 40:60 (urine)

Flow rate: 2

Injection volume: 10-50

Detector: UV 280

CHROMATOGRAM

Retention time: 7 (plasma), 5 (urine)

Internal standard: polythiazide (9 (plasma), 8 (urine))

Limit of detection: 10 ng/mL (plasma)

Limit of quantitation: 50 ng/mL (urine), 20 ng/mL (plasma)

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Meyer,M.C.; Hwang,P.; Straughn,A.B.; Rotenberg,K. HPLC determination of benzthiazide in biologic material, *Biopharm.Drug Dispos.*, **1982**, 3, 1-9.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: not specified

Mobile phase: MeOH:water containing 200 mM ammonia and 200 mM glycolic acid 25:75

Flow rate: 0.8

Detector: MS, Hewlett-Packard 9000-300 quadrupole, thermospray, stem 108-115°, tip 180-205°, ion source 276°, filament on, negative ion mode

KEY WORDS

LC-MS

REFERENCE

Kim,Y.; Park,S.; Park,J.; Lee,W. Detection of benzthiazide by high-performance liquid chromatography-thermospray mass spectrometry, *J.Chromatogr.A*, **1995**, 689, 170-174.

SAMPLE

Matrix: urine

Sample preparation: 2 mL Urine + 0.5 g solid buffer I (pH 5-5.5), vortex 15 s, add 4 mL ethyl acetate, agitate for 10 min, centrifuge at 600 g for 5 min. Remove organic layer and vortex it with 2 mL 5% aqueous lead acetate for 10 s, centrifuge at 600 g for 5 min, remove and keep organic phase. 2 mL Urine + 0.5 g solid buffer II (pH 9-9.5), vortex 15 s, add 4 mL ethyl acetate, agitate for 10 min, centrifuge at 600 g for 5 min. Remove organic layer and combine it with previous organic layer. Evaporate to dryness at 50° under a stream of nitrogen, reconstitute in 300 μ L 50 μ g/mL β -hydroxyethyltheophylline in MeOH, inject 5 μ L aliquot. (Solid buffer I was $\text{KH}_2\text{PO}_4\text{:Na}_2\text{HPO}_4$ 99:1, solid buffer II was $\text{NaHCO}_3\text{:K}_2\text{CO}_3$ 3:2.)

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m HP Hypersil ODS (A) or HP LiChrosorb RP-18 (B)

Mobile phase: Gradient. MeCN:buffer from 15:85 at 2 min to 80:20 at 20 min (Buffer was 50 mM NaH_2PO_4 containing 16 mM propylamine hydrochloride, adjusted to pH 3 with concentrated phosphoric acid.)

Flow rate: 1

Injection volume: 5

Detector: UV 230, UV 275

CHROMATOGRAM

Retention time: 13.5 (A), 14.3 (B)

Internal standard: β -hydroxyethyltheophylline (3.7 (A), 4.4 (B))

Limit of detection: 1000 ng/mL

OTHER SUBSTANCES

Extracted: furosemide, metolazone, amiloride, acetazolamide, chlorothiazide, hydrochlorothiazide, quinethazone, triamterene, hydroflumethiazide, chlorthalidone, dichlorphenamide, trichloromethiazide, methyclothiazide, cyclothiazide, polythiazide, bendroflumethiazide, ethacrynic acid, bumetanide, probenecid, spironolactone, canrenone, flumethiazide

Noninterfering: acetaminophen, aspirin, caffeine, diflunisal, fenoprofen, ibuprofen, indomethacin, methocarbamol, naproxen, phenylbutazone, sulindac, tetracycline, theobromine, theophylline, tolmetin, trimethoprim, verapamil

REFERENCE

Cooper,S.F.; Massé,R.; Dugal,R. Comprehensive screening procedure for diuretics in urine by high-performance liquid chromatography, *J.Chromatogr.*, **1989**, 489, 65-88.

SAMPLE

Matrix: urine

Sample preparation: Make 5 mL urine alkaline (pH 9-10), add 2 g NaCl, extract twice with 6 mL ethyl acetate. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 200 μ L MeCN/water, inject a 10-20 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 4 5 μ m SGE 100 GL-4 C18P (Scientific Glass Engineering)

Mobile phase: MeCN:MeOH:water:trifluoroacetic acid 4.5:10.5:85:0.5

Flow rate: 0.8 or 1

Injection volume: 10-20

Detector: MS, ZAB2-SEQ (VG), PSP source coupled to LC, source 250°, probe 240-260°, scan m/z 200-550 or UV 270

CHROMATOGRAM**Retention time:** 7.2**Limit of detection:** 50 ng (by MS)

OTHER SUBSTANCES**Extracted:** amiloride, chlorthalidone, triamterene, furosemide, bendroflumethiazide

REFERENCE

Ventura,R.; Fraisse,D.; Becchi,M.; Paisse,O.; Segura,J. Approach to the analysis of diuretics and masking agents by high-performance liquid chromatography-mass spectrometry in doping control, *J.Chromatogr.*, **1991**, 562, 723–736.

SAMPLE**Matrix:** urine

Sample preparation: Buffer urine to 4.9 by mixing with an equal volume of pH 4.9 200 mM sodium phosphate buffer. Inject a 40 μ L aliquot onto column A with mobile phase A, after 3 min backflush the contents of column A onto column B with mobile phase B and start the gradient. At the end of the run re-equilibrate for 10 min.

HPLC VARIABLES

Column: A 20 \times 4 5 μ m Hypersil octadecylsilica ODS; B 200 \times 4.6 5 μ m Shiseido SG-120 polymer-based C18

Mobile phase: A water; B Gradient. MeCN:buffer from 7:93 to 15:85 over 3.5 min, to 50:50 over 8.5 min, maintain at 50:50 for 11 min (Buffer was 6.9 g $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$ in 1 L water, pH adjusted to 3.1 with phosphoric acid.)

Flow rate: 1**Injection volume:** 40**Detector:** UV 230

CHROMATOGRAM**Retention time:** 17.6**Limit of detection:** 500 ng/mL

OTHER SUBSTANCES

Extracted: acetazolamide, amiloride, bendroflumethiazide, bumetanide, caffeine, carbamazepine, chlorothiazide, chlorthalidone, clopamide, dichlorfenamide, ethacrynic acid, furosemide, hydrochlorothiazide, metyrapone, probenecid, spironolactone, triamterene, trichlormethiazide

KEY WORDS

column-switching; optimum detection wavelengths vary for each drug

REFERENCE

Saarinén,M.; Sirén,H.; Riekkola,M.-L. A column switching technique for the screening of diuretics in urine by high performance liquid chromatography, *J.Liq.Chromatogr.*, **1993**, 16, 4063–4078.

SAMPLE**Matrix:** urine

Sample preparation: 5 mL Urine + 50 μ L 100 μ g/mL 7-propyltheophylline in MeOH + 200 μ L ammonium chloride buffer + 2 g NaCl, extract with 6 mL ethyl acetate by rocking at 40 movements/min for 20 min and centrifuging at 800 g for 5 min, repeat extraction, combine organic layers, evaporate to dryness at 40° under a stream of nitrogen. Reconstitute in 200 μ L MeCN:water 15:85 and inject 20 μ L aliquots. (Ammonium chloride buffer was 28 g ammonium chloride in 100 mL water with the pH adjusted to 9.5 with concentrated ammonia solution.)

HPLC VARIABLES

Column: 75 \times 4.6 3 μ m Ultrasphere ODS

Mobile phase: Gradient. MeCN:100 mM ammonium acetate adjusted to pH 3 with concentrated phosphoric acid. From 10:90 to 15:85 over 2 min to 55:45 over 3 min to 60:40 over 3 min. Kept at 60:40 for 1 min, decreased to 10:90 over 1 min and equilibrated at 10:90 for 2 min.

Flow rate: 1

Injection volume: 20

Detector: UV 270

CHROMATOGRAM

Retention time: 6.5

Internal standard: 7-propyltheophylline (4.5)

Limit of detection: 100 ng/mL

OTHER SUBSTANCES

Simultaneous: xipamide, bumetanide, acetazolamide, amiloride, bendroflumethiazide, buthiazide, caffeine, canrenone, chlorthalidone, clopamide, cyclothiazide, diclofenamide, ethacrynic acid, furosemide, hydrochlorothiazide, mesocarb, morazone, piretanide, polythiazide, probenecid, spironolactone, torsemide, triamterene

REFERENCE

Ventura,R.; Nadal,T.; Alcalde,P.; Pascual,J.A.; Segura,J. Fast screening method for diuretics, probenecid and other compounds of doping interest, *J.Chromatogr.A*, **1993**, 655, 233–242.


Bzntropine mesylate

Molecular formula: $C_{22}H_{29}NO_4S$

Molecular weight: 403.54

CAS Registry No.: 132-17-2, 86-13-5 (free base)

Merck Index: 1156



SAMPLE

Matrix: blood

Sample preparation: 10 mL Plasma or whole blood + 1 mL 1 M NaOH, extract twice with 10 mL hexane for 30 min. Remove the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 1 mL 100 mM HCl, add 5 mL chloroform, vortex for 1 min, centrifuge. Remove a 4.5 mL aliquot of the organic layer and evaporate it to dryness, reconstitute the residue in 100 μ L mobile phase, inject a 50 μ L aliquot. (It is implied, but not explicitly stated in the paper, that this extraction procedure works for this compound.)

HPLC VARIABLES

Column: 10 μ m Micropak CN (Varian)

Mobile phase: MeCN:20 mM ammonium acetate 90:10

Flow rate: 2.5

Injection volume: 50

Detector: UV 254

CHROMATOGRAM

Retention time: 13.2

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Simultaneous: acetophenazine, amitriptyline, butaperazine, carphenazine, chlorpromazine, fluphenazine, haloperidol, imipramine, mesoridazine, nortriptyline, orphenadrine, piperacetazine, promazine, promethazine, thioridazine, thiothixene, trifluoperazine, tri-flupromazine, trihexyphenidyl, trimeprazine

KEY WORDS

plasma; whole blood

REFERENCE

Curry, S.H.; Brown, E.A.; Hu, O.Y.-P.; Perrin, J.H. Liquid chromatographic assay of phenothiazine, thioxanthene and butyrophenone neuroleptics and antihistamines in blood and plasma with conventional and radial compression columns and UV and electrochemical detection, *J. Chromatogr.*, **1982**, 231, 361-376.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Plasma + 150 μ L 200 ng/mL desipramine hydrochloride + 150 μ L 5 M NaOH, vortex, add 1 mL ethylene glycol, vortex, add 10 mL hexane, shake on a rotary shaker at 30 rpm for 30 min, centrifuge at 1000 g at 4°. Remove the organic layer and add it to 300 μ L 100 mM HCl, shake at high speed for 20 min, centrifuge, inject a 200 μ L aliquot of the aqueous layer.

HPLC VARIABLES

Column: 150 \times 3.9 5 μ m Spherisorb C8

Mobile phase: MeCN:buffer 60:40 (Buffer was 1.5 mL triethylamine in 1 L water adjusted to pH 3.0 with 85% phosphoric acid.)

Flow rate: 1.5

Injection volume: 200

Detector: UV 199

CHROMATOGRAM

Retention time: 7.0

Internal standard: desipramine hydrochloride (5.3)

Limit of quantitation: 0.25 ng/mL

OTHER SUBSTANCES

Simultaneous: hyoscyamine, orphenadrine, bromocriptine, biperiden

Noninterfering: amantadine, carbidopa, levodopa

KEY WORDS

plasma

REFERENCE

Selinger,K.; Lebel,G.; Hill,H.M.; Discenza,C. High-performance liquid chromatographic method for the analysis of benztropine in human plasma, *J.Chromatogr.*, **1989**, 491, 248–252.

SAMPLE

Matrix: blood, tissue

Sample preparation: Blood or serum. 1 mL Blood or serum + 1 µg cianopramine + 1 mL water, vortex, add 1 mL 200 mM sodium carbonate, vortex, add 6 mL hexane:1-butanol 95:5, gently agitate for 30 min, centrifuge at 2500 g for 5 min. Remove the organic layer and add it to 100 µL 0.2% phosphoric acid, agitate gently for 30 min, centrifuge for 5 min. Remove the organic layer and inject a 30 µL aliquot of the aqueous layer. Liver homogenate. 0.5 mL Liver homogenate + 10 µg cianopramine + 500 µL 2% sodium tetraborate + 8 mL hexane:1-butanol 95:5, gently agitate for 30 min, centrifuge at 2500 g for 5 min. Remove the organic layer and add it to 400 µL 0.2% phosphoric acid, agitate gently for 30 min, centrifuge for 5 min. Remove the organic layer and inject a 30 µL aliquot of the aqueous layer.

HPLC VARIABLES

Guard column: 15 × 3.2 7 µm RP-18 Newguard (Applied Biosystems)

Column: 100 × 4.6 5 µm Brownlee Spheri-5 RP-18

Mobile phase: MeCN:100 mM NaH₂PO₄:diethylamine 40:57.5:2.5

Flow rate: 2

Injection volume: 30

Detector: UV 220

CHROMATOGRAM

Retention time: 22.75

Internal standard: cianopramine (8.93)

OTHER SUBSTANCES

Simultaneous: amoxapine, brompheniramine, chlorpheniramine, chlorpromazine, clomipramine, cyproheptadine, desipramine, diphenhydramine, dothiepin, doxepin, fluoxetine, haloperidol, imipramine, loxapine, maprotiline, meperidine, mesoridazine, methadone, metoclopramide, mianserin, moclobemide, nomifensine, nordoxepin, norfluoxetine, norpropoxyphene, northiaden, nortriptyline, pentobarbital, pheniramine, propoxyphene, propranolol, protriptyline, quinidine, quinine, sulfonidazine, thioridazine, thiothixene, translycypromine, trazodone, trihexiphenidyl, triprolidine

Noninterfering: dextromethorphan, norphethidine, phenoxybenzamine, prochlorperazine, trifluoperazine

Interfering: promethazine, trimipramine, amitriptyline

KEY WORDS

serum; whole blood; liver

REFERENCE

McIntyre, I.M.; King, C.V.; Skafidis, S.; Drummer, O.H. Dual ultraviolet wavelength high-performance liquid chromatographic method for the forensic or clinical analysis of seventeen antidepressants and some selected metabolites, *J. Chromatogr.*, **1993**, 621, 215–223.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in MeOH:water 1:1 at a concentration of 50 µg/mL, inject a 10 µL aliquot.

HPLC VARIABLES

Column: 300 × 3.9 10 µm µBondapak C18

Mobile phase: MeOH:acetic acid:triethylamine:water 60:1.5:0.5:38

Flow rate: 1.5

Injection volume: 10

Detector: UV

CHROMATOGRAM

Retention time: k' 2.10

REFERENCE

Roos, R.W.; Lau-Cam, C.A. General reversed-phase high-performance liquid chromatographic method for the separation of drugs using triethylamine as a competing base, *J. Chromatogr.*, **1986**, 370, 403–418.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlordiazeoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenoprofen, fenpro-

porex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiacol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephenetermine, mephentyoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methyprylon, metoprolol, mi-boleron, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sufadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripelennamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.

Benzydamine


Molecular formula: $C_{19}H_{23}N_3O$

Molecular weight: 309.41

CAS Registry No.: 642-72-8, 132-69-4 (HCl)

Merck Index: 1157

Lednicer No.: 1 323



SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A: B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 215.8

CHROMATOGRAM

Retention time: 14.955


KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

Benzyl alcohol



Molecular formula: C₇H₈O

Molecular weight: 108.14

CAS Registry No.: 100-51-6

Merck Index: 1159

SAMPLE

Matrix: bulk

Sample preparation: Dissolve 600 µg 7-[(imidazolemethanoyl)methoxy]-4-methylcoumarin in 2.5 mL toluene and 750 µL MeCN, add 20 µL benzyl alcohol, add 150 µL 93.3 µg/mL 4-dimethylaminopyridine in MeCN, mix vigorously, heat at 60° for 1.5 h, cool, evaporate to dryness under a stream of nitrogen at room temperature, reconstitute the residue in 300 µL mobile phase, inject a 5 µL aliquot. (Preparation of 7-[(imidazolemethanoyl)methoxy]-4-methylcoumarin is as follows. Stir 102.3 mg 7-(carboxymethoxy)-4-methylcoumarin in 7 mL THF, add 70.9 mg 1,1'-carbonyldiimidazole in one portion, reflux for 30 min, stir at room temperature for 5 h. Filter and dry the solid under reduced pressure to obtain 7-[(imidazolemethanoyl)methoxy]-4-methylcoumarin as a white solid (mp 161-162°). Fluorescence detection can also be used.)

HPLC VARIABLES

Column: 150 × 3.9 4 µm Nova-Pak C18

Mobile phase: MeCN:MeOH:100 mM pH 5.5 ammonium acetate buffer 50:1.5:48.5

Flow rate: 0.7

Injection volume: 5

Detector: MS, Hewlett-Packard 5989A, thermospray interface, filament-assisted ionization mode, ion source 280°, probe stem 112°, probe tip 235-245°

CHROMATOGRAM

Retention time: 7

Limit of detection: 0.8 ng


KEY WORDS

derivatization

REFERENCE

Phillips, L.R.; Supko, J.G.; Wolfe, T.L.; Malspeis, L. Precolumn derivatization of hydroxy compounds with 7-[(imidazolemethanoyl)methoxy]-4-methylcoumarin (IMMC) and LC/TSP-MS of the resulting esters, *Proc. Am. Soc. Mass Spectrom.*, **1995**, *43*, 163-164.

Benzyl benzoate



Molecular formula: $C_{14}H_{12}O_2$

Molecular weight: 212.25

CAS Registry No.: 120-51-4

Merck Index: 1162

SAMPLE

Matrix: formulations

Sample preparation: Dilute 0.5 mL of nanocapsules suspension 1:200 with MeCN, filter, inject an aliquot. Alternatively, evaporate 5 mL of a nanocapsule suspension to dryness and dissolve the residue in 150 mL dichloromethane or ethyl acetate, dry over anhydrous sodium sulfate. Evaporate to dryness under reduced pressure, take up the residue in 50 mL MeOH, dilute 1:20 with MeOH, inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.5 µm Nucleosil C18

Mobile phase: MeCN: water 75:25

Flow rate: 1

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 4.71

OTHER SUBSTANCES

Simultaneous: progesterone

Noninterfering: poly-ε-caprolactone

KEY WORDS

nanocapsules

REFERENCE

Benali,S.; Tharasse-Bloch,C.; André; Vérité,P.; Duclos,R.; Lafont,O. Determination of progesterone in nanocapsules by high performance liquid chromatography, *J.Liq.Chromatogr.Rel.Technol.*, **1997**, 20, 3233–3243.

SAMPLE

Matrix: formulations

Sample preparation: Injections. Extract 2 mL with EtOH:water 85:15, make up extracts to 100 mL with EtOH:water 85:15, remove a 2 mL aliquot and add it to 1 mL 1 mg/mL hydrocortisone in EtOH. Dilute this mixture to 50 mL with EtOH:water 50:50, inject an aliquot. Suspensions. Dilute 2 mL suspension to 100 mL with EtOH, filter (if necessary), remove a 2 mL aliquot and add it to 1 mL 1 mg/mL hydrocortisone in EtOH. Dilute this mixture to 50 mL with EtOH, inject an aliquot.

HPLC VARIABLES

Column: 300 × 4 µm Bondapak CN

Mobile phase: MeOH:20 mM KH_2PO_4 30:70

Flow rate: 2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 6.5

Internal standard: hydrocortisone (2.5)

OTHER SUBSTANCES

Simultaneous: medroxyprogesterone acetate

Noninterfering: polyethylene glycol 4000, myristyl-gamma-picolinium chloride, methyl-cellulose, thimerosal

Interfering: progesterone

REFERENCE

Das Gupta, V. Quantitation of hydroxyprogesterone caproate, medroxyprogesterone acetate, and progesterone by reversed-phase high-pressure liquid chromatography, *J.Pharm.Sci.*, **1982**, *71*, 294–297.

SAMPLE

Matrix: formulations

Sample preparation: Oils. 1 mL Sample + 25 mL MeOH:water 90:10, shake vigorously for 5 min, centrifuge, inject a 10 μ L aliquot of the supernatant. Tablets. Grind a tablet to a fine powder, add 25 mL MeOH, sonicate for 5–10 min, filter (0.45 μ m), discard first 5 mL of filtrate, inject a 10 μ L aliquot of the remaining filtrate. Suspensions (aqueous). Make up 5 mL to 50 mL with MeOH, filter (0.45 μ m), discard first 5 mL of filtrate, inject a 10 μ L aliquot of the remaining filtrate.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Zorbax ODS

Mobile phase: MeOH:water 75:25

Flow rate: 1.5

Injection volume: 10

Detector: UV 240

CHROMATOGRAM

Retention time: 11.8

Limit of detection: 5 μ g/mL

OTHER SUBSTANCES

Simultaneous: aspirin, caffeine, formebolone, benzyl alcohol, testolactone, cortisone, fluoxymesterone, norethindrone, oxandrolone (UV 210), boldenone, ethisterone, methandrostenedione, nandrolone, norgestrel, testosterone, dehydroepiandrosterone (UV 210), mibolone, methyltestosterone, methandriol (UV 210), norethindrone acetate, calusterone, mesterolone (UV 210), norethandrolone, trenbolone acetate, nandrolone acetate, testosterone acetate, stanozolol, oxymetholone, nandrolone propionate, methenolone acetate, testosterone propionate

KEY WORDS

oils; tablets; suspensions

REFERENCE

Walters, M.J.; Ayers, R.J.; Brown, D.J. Analysis of illegally distributed anabolic steroid products by liquid chromatography with identity confirmation by mass spectrometry or infrared spectrophotometry, *J.Assoc.Off.Anal.Chem.*, **1990**, *73*, 904–926.

SAMPLE

Matrix: formulations

Sample preparation: Weigh out 50 mg formulation, add 5 mL 1.5 mg/mL benzophenone in MeOH, make up to 50 mL with MeOH. Dilute 1 mL of this solution to 10 mL with MeOH, filter (0.45 μ m PTFE membrane), inject a 10 μ L aliquot.

HPLC VARIABLES

Guard column: 30 mm long Brownlee guard column

Column: 220 × 4.6 5 µm C18 (Brownlee)

Mobile phase: MeCN:water 60:40

Flow rate: 2

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 4.5

Internal standard: benzophenone (3.3)

Limit of quantitation: 7000 ng/mL

OTHER SUBSTANCES

Simultaneous: benzocaine

KEY WORDS

dermatological preparations

REFERENCE

Gigante,B.; Barros,A.M.V.; Teixeira,A.; Marcelo-Curto,M.J. Separation and simultaneous high-performance liquid chromatographic determination of benzocaine and benzyl benzoate in a pharmaceutical preparation, *J.Chromatogr.*, **1991**, 549, 217–220.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in MeOH at a concentration of 100 µg/mL, inject a 5 µL aliquot.

HPLC VARIABLES

Guard column: 70 × 2.1 CO:Pell ODS

Column: 300 × 3.9 Bondex C18 (Phenomenex)

Mobile phase: MeOH:water 85:15

Flow rate: 1

Injection volume: 5

Detector: UV 254

CHROMATOGRAM

Retention time: 5

OTHER SUBSTANCES

Also analyzed: boldenone, testosterone, nandrolone, and their esters

REFERENCE

Noggle,F.T.,Jr.; Clark,C.R.; DeRuiter,J. Liquid chromatographic and mass spectral analysis of the anabolic 17-hydroxy steroid esters, *J.Chromatogr.Sci.*, **1990**, 28, 263–268.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 150 × 3.9 5 µm Zorbax BP-ODS

Mobile phase: MeCN:50 mM pH 7.2 sodium phosphate buffer 40:60

Flow rate: 1

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 8

Internal standard: benzyl benzoate

OTHER SUBSTANCES

Extracted: HO-221

KEY WORDS

benzyl benzoate is IS

REFERENCE

Kondo,N.; Iwao,T.; Hirai,K.-I.; Fukuda,M.; Yamanouchi,K.; Yokoyama,K.; Miyaji,M.; Ishihara,Y.; Kon,K.; Ogawa,Y.; Mayumi,T. Improved oral absorption of enteric coprecipitates of a poorly soluble drug, *J.Pharm.Sci.*, **1994**, *83*, 566–570.

Bepridil


Molecular formula: C₂₄H₃₄N₂O

Molecular weight: 366.55

CAS Registry No.: 64706-54-3, 74764-40-2 (HCl monohydrate)

Merck Index: 1188

Lednicer No.: 3 46



SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform: isopropanol:n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 249

CHROMATOGRAM

Retention time: 18.30

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds(all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfipyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; mocllobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vindesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phencyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol;

aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrodine; phenylbutazone; demoxiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, 1995, 40, 254–262.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A: B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10–30

Detector: UV 202.8

CHROMATOGRAM

Retention time: 19.503

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, 1997, 763, 149–163.

SAMPLE

Matrix: formulations

Sample preparation: Weigh out ground tablets or capsules equivalent to 150 mg bepridil.HCl, add 150 mL MeCN, shake for 30 min, dilute to 200 mL with MeCN, filter (Schleicher & Schüll paper, grade 588), inject a 20 µL aliquot.

HPLC VARIABLES

Column: 300 × 4.6 10 μm μBondapak C18

Mobile phase: MeCN:buffer 580:405 (Buffer was 1.1 g sodium 1-heptanesulfonate in 405 mL water, adjust to pH 2.37 with glacial acetic acid (ca. 15 mL).)

Column temperature: 35

Flow rate: 1.3

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 6.8

OTHER SUBSTANCES

Simultaneous: impurities, benzaldehyde, N-benzylaniline, benzoic acid

KEY WORDS

stability-indicating; rugged; capsules; tablets

REFERENCE

Renzi, N.L.; Fronheiser, M.E.; Duong, H.T.; Fulton, D.J.; Rabinowitz, M. Stability-indicating high-performance liquid chromatography assay for bepridil hydrochloride drug substance and drug products, *J. Chromatogr.*, **1989**, 462, 398–405.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Chirex 3014 (Phenomenex)

Mobile phase: Hexane:1,2-dichloroethane:EtOH/trifluoroacetic acid 62:35:3 (EtOH/trifluoroacetic acid was premixed 20:1.)

Flow rate: 0.7–1

Injection volume: 20

Detector: UV 250

KEY WORDS

chiral; $\alpha = 1.22$ for enantiomers

REFERENCE

Cleveland, T. Pirkle-concept chiral stationary phases for the HPLC separation of pharmaceutical racemates, *J. Liq. Chromatogr.*, **1995**, 18, 649–671.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 100 μM solution in buffer, inject a 20 μL aliquot.

HPLC VARIABLES

Column: 100 × 4.6 column containing riboflavin binding proteins (Prepare as follows. Add riboflavin to saturate protein of egg yolk, homogenize with 3 volumes buffer, centrifuge, add the supernatant to a 500 × 30 column of DEAE-cellulose (Whatman) equilibrated with buffer, wash extensively with buffer to remove bound protein, elute riboflavin binding proteins (RFBP) with buffer containing 200 mM NaCl (RFBP has intense yellow color, absorption at 455 nm). Purify RFBP on a Sephadex G-100 column with 50 mM pH 7.5 Tris-HCl buffer as eluent, remove the bound riboflavin by extensive dialysis at pH 3.0. Add 4.5 g N,N-disuccinylimidyl carbonate to 3 g Nucleosil 5NH₂ slurried in MeCN, filter, wash with MeCN, wash with 50 mM pH 7.5 phosphate buffer. Suspend 300 mg RFBP in 50 mM phosphate buffer, add the activated silica, mix gently for 2 h using a rotary evaporator, filter, wash with sterile water, wash with isopropanol:water 1:2, pack in a 100 × 4.6 column.) (Buffer was 100 mM pH 5.3 sodium acetate.)

Mobile phase: EtOH:50 mM pH 5.5 KH₂PO₄ 5:95

Flow rate: 0.8

Injection volume: 20

Detector: UV

CHROMATOGRAM

Retention time: k' 11.96

OTHER SUBSTANCES

Simultaneous: lorazepam, manidipine, nicardipine, oxazepam

KEY WORDS

chiral; $\alpha = 1.21$

REFERENCE

Massolini,G.; De Lorenzi,E.; Ponci,M.C.; Gandini,C.; Caccialanza,G.; Monaco,H.L. Egg yolk riboflavin binding protein as a new chiral stationary phase in high-performance liquid chromatography, *J.Chromatogr.A*, **1995**, 704, 55-65.


Berberine

Molecular formula: $C_{20}H_{18}NO_4$

Molecular weight: 336.37

CAS Registry No.: 2086-83-1

Merck Index: 1192



SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbomal, chloramphenicol, chlordiazepoxide, chloroquine, chlorothiazide, chloroxylonol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenpropofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiacal, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarboxazid, isonicarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methylidopa, methylidopamine, methylphenidate, methylprednisolone, methyltestosterone, methypylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nyldrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, py-

rilamine, pyriethyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleennamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

- Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.


Besipirdine

Molecular formula: C₁₆H₁₇N₃

Molecular weight: 251.33

CAS Registry No.: 119257-34-0, 130953-69-4 (HCl), 119257-40-8 (maleate)

Merck Index: 1223



SAMPLE

Matrix: cell cultures

Sample preparation: Mix 1 mL cell culture with 50 μ L 1 mg/mL IS in MeOH and 1 mL saturated aqueous sodium bicarbonate, add 3 mL ethyl acetate:cyclohexane 50:50, invert at 18 cycles/min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 2 mL mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 μ m Hypersil-Phenyl

Mobile phase: MeCN:buffer 85:15 (Buffer was 5.8 mM triethylammonium formate adjusted to pH 2.5 with 90% formic acid.)

Flow rate: 1

Injection volume: 20

Detector: UV 270

CHROMATOGRAM

Retention time: 10

Internal standard: 3-ethyl-N-methyl-N-(4-pyridinyl)-1H-indol-1-amine hydrochloride

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

details of semipreparative HPLC

REFERENCE

Rao,G.P.; Davis,P.J. Microbial models of mammalian metabolism. Biotransformations of HP 749 (besipirdine) using *Cunninghamella elegans*, *Drug Metab.Dispos.*, **1997**, *25*, 709–715.


Beta-carotene

Molecular formula: C₄₀H₅₆

Molecular weight: 536.88

CAS Registry No.: 7235-40-7

Merck Index: 1902



SAMPLE

Matrix: blood

Sample preparation: Centrifuge 200 μ L serum, add 200 μ L IS and 200 μ L EtOH, mix on orbital shaker for 5 min, add 200 μ L water and 500 μ L hexane, mix for 10 min, centrifuge at 2000 g for 10 min at 17°, remove 300 μ L upper organic layer. Re-extract with 300 μ L hexane, mix for 10 min, centrifuge at 4000 g for 10 min at 17°, remove 300 μ L upper organic layer. Combine the organic layers, and evaporate them to dryness under vacuum in 15 min. Reconstitute the residue with 300 μ L MeOH:EtOH:hexane 88:10:2, vortex for 10 min, inject a 40 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 4.6 3 μ m Adsorbosphere HS C18 + 150 \times 4.6 3 μ m Adsorbosphere HS C18 in series

Mobile phase: Gradient. A was MeCN:MeOH 60:40 containing 0.05% acetic acid. B was MeCN:MeCN:dichloromethane 45.6:30.4:24 containing 0.04% acetic acid. A:B 100:0 for 7 min then 0:100 for 10.4 min (step gradient), re-equilibrate at initial conditions for 5.6 min.

Column temperature: 37

Flow rate: 0.9

Injection volume: 40

Detector: UV 450

CHROMATOGRAM

Retention time: 16.5

Internal standard: tocol (UV 292) (10.1), echinenone (UV 450) (12.8)

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Extracted: canthaxanthine (UV 473), α -carotene (UV 450), β -cryptoxanthine (UV 450), lutein (UV 450), lycopene (UV 473), vitamin A (UV 325), vitamin E (UV 292), zeaxanthin (UV 450), nonidentified carotenoids

KEY WORDS

serum

REFERENCE

Steghens,J.-P.; van Kappel,A.L.; Riboli,E.; Collombel,C. Simultaneous measurement of seven carotenoids, retinol and α -tocopherol in serum by high-performance liquid chromatography, *J.Chromatogr.B*, **1997**, 694, 71–81.

SAMPLE

Matrix: blood

Sample preparation: Add 1 mL 1 μ g/mL retinyl palmitate, 1 μ g/mL retinyl palmitate, and 25 μ g/mL α -tocopheryl acetate in EtOH to 1 mL serum or plasma while continuously vortexing, add 3 mL hexane, vortex for 2 min, centrifuge at 2500 g for 2 min, remove the upper phase, add 2 mL hexane to the lower layer, repeat extraction. Combine the upper layers and evaporate them to dryness under a stream of nitrogen at 37°, reconstitute the residue in 200 μ L mobile phase, inject a 40 μ L aliquot.

HPLC VARIABLES**Guard column:** C18 (Waters)**Column:** 5 μ m Biophase ODS C18 (Bioanalytical Systems)**Mobile phase:** MeCN:chloroform:isopropanol:water 78:16:3.5:2.5**Flow rate:** 2**Injection volume:** 40**Detector:** UV 460 (UV 292 for tocopherol)

CHROMATOGRAM**Retention time:** 16.46**Internal standard:** retinyl acetate (3.07), retinyl palmitate (18.66), α -tocopheryl acetate (8.33)

OTHER SUBSTANCES**Extracted:** vitamin A (retinol), vitamin E (α -tocopherol), gamma-tocopherol, α -carotene, lycopene, cryptoxanthin

KEY WORDS

serum; plasma

REFERENCEKaplan, L.A.; Miller, J.A.; Stein, E.A.; Stampfer, M.J. Simultaneous, high-performance liquid chromatographic analysis of retinol, tocopherols, lycopene, and α - and β -carotene in serum and plasma, *Methods Enzymol.*, **1990**, 189, 155–167.

SAMPLE**Matrix:** blood**Sample preparation:** 250 μ L Serum + 25 μ L 80 μ g/mL tocol in EtOH + 250 μ L 20 μ g/mL BHT (butylated hydroxytoluene) in EtOH + 1.5 mL hexane, vortex for 1 min, remove 1 mL of upper layer, add 500 μ L hexane, vortex for 1 min, remove 300 μ L of upper layer. Combine the hexane extracts, evaporate to dryness under a stream of inert gas. Reconstitute in 250 μ L 20 μ g/mL BHT in EtOH, sonicate, centrifuge if necessary, inject a 25 μ L aliquot.

HPLC VARIABLES**Column:** 250 \times 4.6 5 μ m Vydac 201TP54 (wide pore, polymerically bonded C18)**Mobile phase:** Gradient. A was MeOH:n-butanol:water 75:10:15 containing 50 mM ammonium acetate, pH 5.5. B was MeOH:n-butanol:water 88:10:2 containing 50 mM ammonium acetate, pH 5.5. A:B 100:0 for 3 min, to 0:100 over 15 min, maintain at 0:100 for 17 min**Injection volume:** 25**Detector:** UV 325 for 7 min, UV 295 for 13 min, UV 450 for 14 min or E, glassy carbon electrode, Ag/AgCl reference electrode +1050 mV for retinol, +900 mV for tocol, +750 mV for α -tocopherol, +700 mV for β -carotene

CHROMATOGRAM**Retention time:** 31**Internal standard:** tocol (13)**Limit of detection:** 2.1 μ g/mL (E), 29 μ g/mL (UV)

OTHER SUBSTANCES**Extracted:** vitamin A (retinol), vitamin E (α -tocopherol), gamma-tocopherol, lutein, zeaxanthin, cryptoxanthin, α -carotene, 9-cis- β -carotene

KEY WORDS

serum

REFERENCE

MacCrehan, W.A. Determination of retinol, α -tocopherol, and β -carotene in serum by liquid chromatography, *Methods Enzymol.*, **1990**, 189, 172–181.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Serum + 100 μ L EtOH + 100 μ L α -tocopheryl acetate in EtOH, vortex for 5 s, add 500 μ L hexane, vortex for 2 min, centrifuge at 700 g for 5 min. Remove 250 μ L of the hexane layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 200 μ L mobile phase, mix for 2 min, inject a 50 μ L aliquot.

HPLC VARIABLES

Guard column: 30 \times 4.6 10 μ m Spheri-10 RP18

Column: 150 \times 4.6 5 μ m Ultrasphere ODS

Mobile phase: MeCN:dichloromethane:MeOH 70:20:10

Flow rate: 1.2

Injection volume: 50

Detector: UV 325 for 3.5 min, UV 291 for 4.5 min, UV 450 for 6 min

CHROMATOGRAM

Retention time: 11.79

Internal standard: α -tocopheryl acetate (6.30)

Limit of detection: 50 nM

OTHER SUBSTANCES

Extracted: vitamin A, vitamin E

KEY WORDS

protect from light; serum

REFERENCE

Arnaud, J.; Fortis, I.; Blachier, S.; Kia, D.; Favier, A. Simultaneous determination of retinol, α -tocopherol and β -carotene in serum by isocratic high-performance liquid chromatography, *J. Chromatogr.*, **1991**, 572, 103–116.

SAMPLE

Matrix: blood

Sample preparation: 2.5 mL Plasma + 2.5 mL 18 ng/mL IS1 and 10 ng/mL IS2 in EtOH, shake vigorously for 20 s, centrifuge at 1200 g for 5 min, add 5 mL diethyl ether, shake vigorously, centrifuge for 5 min, extract twice more with 5 mL ether. Combine ether layers, wash with 15 mL 5% NaCl, dry over sodium sulfate, evaporate to dryness under vacuum at 35°. Dissolve residue in 1–2 mL dichloromethane, filter (0.45 μ m). Evaporate to dryness under a stream of nitrogen, make up to 100 μ L with MeCN:MeOH:dichloromethane:hexane 45:10:22.5:22.5, inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: 30 \times 4.6 5 μ m Spheri-5-C18 (Brownlee)

Column: 250 \times 4.6 5 μ m Microsorb C18 (Rainin)

Mobile phase: Gradient. MeCN:MeOH:dichloromethane:hexane 85:10:2.5:2.5 for 10 min then to 45:10:22.5:22.5 over 30 min, re-equilibrate for 15 min

Flow rate: 0.7

Injection volume: 20

Detector: UV 470

CHROMATOGRAM

Retention time: 34

Internal standard: IS1 ethyl β -apo-8'-carotenate (18), IS2 (3R)-8'-apo- β -carotene-3,8'-diol (5)

OTHER SUBSTANCES

Extracted: carotenoids, vitamin A (retinol), vitamin E (α -tocopherol)

KEY WORDS

plasma; handle under yellow lights

REFERENCE

Khachik, F.; Beecher, G.R.; Goli, M.B.; Lusby, W.R.; Smith, J.C., Jr. Separation and identification of carotenoids and their oxidation products in the extracts of human plasma, *Anal.Chem.*, **1992**, *64*, 2111–2122.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Serum or plasma + 200 μ L 25 μ g/mL tocopheryl acetate in EtOH, vortex, add 400 μ L butanol:ethyl acetate 50:50, mix for 1 min, add 20 mg sodium sulfate, vortex for 1 min, let stand at -20° for 20 min, centrifuge at 15000 g for 2 min, inject a 10 μ L aliquot of the upper organic layer.

HPLC VARIABLES

Guard column: 5 μ m C18

Column: 110 \times 4.7 5 μ m Partisphere 5 C18 (Whatman)

Mobile phase: MeOH:butanol:water 89.5:5:5.5

Column temperature: 45

Flow rate: 1.5

Injection volume: 10

Detector: UV 340 for 3 min, UV 290 for 1.5 min, UV 280 for 10.5 min, UV 450 for 7 min

CHROMATOGRAM

Retention time: 20.1

Internal standard: tocopheryl acetate (5.3)

Limit of detection: 100 ng/mL

OTHER SUBSTANCES

Extracted: α -carotene, lycopene, δ -tocopherol, gamma-tocopherol, vitamin A, vitamin E, xanthophyll

KEY WORDS

serum; plasma; protect from light

REFERENCE

Lee, B.L.; Chua, S.C.; Ong, H.Y.; Ong, C.N. High-performance liquid chromatographic method for routine determination of vitamins A and E and β -carotene in plasma, *J.Chromatogr.*, **1992**, *581*, 41–47.

SAMPLE

Matrix: blood

Sample preparation: Dilute 1 mL serum 0.5-5 times with saline. Add 1 mL EtOH to 1 mL diluted serum dropwise while vortexing, add 1.5 mL n-heptane, vortex for 1 min, centrifuge at 3000 rpm (Labofuge) for 15 min. Remove 1.3 mL of the organic layer and evaporate it to dryness under a stream of nitrogen at 40° , reconstitute the residue in 40 μ L MeCN:THF 50:50, inject a 5 μ L aliquot.

HPLC VARIABLES

Column: 200 \times 2.1 5 μ m ODS Hypersil

Mobile phase: MeCN:water:THF 81.3:5.7:13

Column temperature: 40
Flow rate: 0.4
Injection volume: 5
Detector: UV 450

CHROMATOGRAM

Retention time: 14.230
Limit of detection: 1 ng

OTHER SUBSTANCES

Extracted: Vitamin A (retinol), Vitamin E (α -tocopherol), probucol, gamma-tocopherol, lycopene, α -carotene, metabolites

KEY WORDS

serum

REFERENCE

Schäfer Elinder,L.; Walldius,G. Simultaneous measurement of serum probucol and lipid-soluble antioxidants, *J.Lipid Res.*, **1992**, 33, 131–137.

SAMPLE

Matrix: blood

Sample preparation: 20-500 μ L Serum + 2 volumes EtOH + 1 mL ethyl acetate + 4-7 μ L of a solution containing 16 mg/mL tocopheryl acetate, 2-3 μ g/mL canthaxanthin, and 10 μ g/mL retinoic acid, vortex for 30 s, centrifuge for 30 s, extract the pellet twice with 0.5-1 mL portions of ethyl acetate, extract the pellet with 0.5-1 mL hexane. Combine the supernatants, add 500 μ L water, vortex, centrifuge. Remove the upper organic layer and evaporate it to dryness under a stream of argon, reconstitute the residue in 100 μ L MeOH:dichloromethane 2:1, inject a 10-90 μ L aliquot.

HPLC VARIABLES

Guard column: C18 (Upchurch)
Column: 300 \times 3.9 5 μ m Resolve C18 (Waters)
Mobile phase: MeCN:dichloromethane:MeOH:1-octanol 90:15:10:0.1
Flow rate: 1
Injection volume: 10-90
Detector: UV 450

CHROMATOGRAM

Retention time: 21
Internal standard: tocopheryl acetate, canthaxanthin, retinoic acid

OTHER SUBSTANCES

Extracted: carotenoids, vitamin A (UV 325), vitamin E (UV 290)

KEY WORDS

protect from light; serum

REFERENCE

Barua,A.B.; Kostic,D.; Olson,J.A. New simplified procedures for the extraction and simultaneous high-performance liquid chromatographic analysis of retinol, tocopherols and carotenoids in human serum, *J.Chromatogr.*, **1993**, 617, 257–264.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Serum or plasma + 500 μ L EtOH containing 4.27 μ M retinyl acetate and 0.31 μ M echinenone, rotamix for 30 s, add 2 mL n-hexane, rotamix for 30 s,

centrifuge at 2000 g for 2 min, repeat extraction with 2 mL n-hexane. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 50 μ L THF, make up to 200 μ L with EtOH, inject a 50 μ L aliquot.

HPLC VARIABLES

Guard column: 50 \times 4.6 5 μ m Spherisorb ODS1

Column: 250 \times 4.6 5 μ m Spherisorb ODS1

Mobile phase: Gradient. A was MeCN:MeOH 20:80 containing 100 mM ammonium acetate. B was 100 mM ammonium acetate in water. A:B from 90:10 to 100:0 over 12 min, maintain at 100:0 for 10 min, re-equilibrate at initial conditions for 5 min

Flow rate: 2

Injection volume: 50

Detector: UV 325 for 7.5 min, UV 292 for 5.5 min, then UV 450

CHROMATOGRAM

Retention time: 19.50

Internal standard: retinyl acetate (5.96), echinenone (15.15)

Limit of detection: 0.13 μ M

OTHER SUBSTANCES

Extracted: cryptoxanthin, lutein, lycopene, vitamin A, vitamin E

KEY WORDS

plasma; protect from light; serum

REFERENCE

Zaman,Z.; Fielden,P.; Frost,P.G. Simultaneous determination of vitamins A and E and carotenoids in plasma by reversed-phase HPLC in elderly and younger subjects, *Clin.Chem.*, **1993**, 39, 2229–2234.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Plasma or serum + 200 μ L 850 ng/mL retinyl acetate in EtOH, mix for 1 min, add 1 mL 0.4 g/L BHT (2,6-di-tert-butyl-4-methylphenol) in n-hexane, shake on a mechanical shaker for 10 min, centrifuge at 2000 g for 5 min, remove 800 μ L of the supernatant, evaporate to dryness at 40° under a stream of nitrogen, re-constitute in 100 μ L MeCN:THF:MeOH 68:22:7, inject a 15 μ L aliquot.

HPLC VARIABLES

Guard column: 15 \times 3.2 7 μ m Lichrosorb RP18

Column: 250 \times 4.6 5 μ m Nucleosil 100-5 C18

Mobile phase: MeCN:THF:MeOH 68:22:7 made up to 100 with 1% ammonium acetate

Flow rate: 1.5

Injection volume: 15

Detector: UV 325 for 3 min, UV 450 for 1.9 min, UV 290 for 2.5 min, UV 470 for 4.6 min, UV 450 for 3 min, then UV 325 for rest of run

CHROMATOGRAM

Retention time: 13

Internal standard: retinyl acetate (2.7)

Limit of detection: 20 ng/mL

OTHER SUBSTANCES

Extracted: vitamin A (retinol), vitamin E (α -tocopherol), lutein, lycopene, α -carotene, zeaxanthin, trans β -carotene, δ -tocopherol

KEY WORDS

plasma; serum; protect from sunlight

REFERENCE

Bui, M.H. Simple determination of retinol, α -tocopherol and carotenoids (lutein, all-*trans*-lycopene, α - and β -carotenes) in human plasma by isocratic liquid chromatography, *J. Chromatogr. B*, **1994**, 654, 129–133.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Serum + 200 μ L nonapreno- β -carotene and retinyl butyrate in EtOH, vortex for 10 s, add 1 mL hexane, vortex for 30 s, centrifuge at 1500 g for 5 min. Remove 900 μ L of the hexane layer and evaporate it to a waxy or glassy consistency (not dryness) under vacuum, dissolve in 100 μ L EtOH, add 100 μ L MeCN, vortex, filter (0.45 μ m), inject a 30 μ L aliquot.

HPLC VARIABLES

Column: 1540 \times 4.6 5 μ m Ultramex C18 (Phenomenex)

Mobile phase: MeCN:EtOH 50:50 containing 0.1 mL/L diethylamine

Column temperature: 29

Flow rate: 0.9

Injection volume: 30

Detector: UV 450

CHROMATOGRAM

Retention time: 8.10

Internal standard: nonapreno- β -carotene (9.5, UV 450), retinyl butyrate (3.5, UV 300)

Limit of detection: 13 nM

OTHER SUBSTANCES

Extracted: vitamin A (retinol), vitamin E (α -tocopherol), lutein, zeaxanthin, β -cryptoxanthin, lycopene, α -carotene, retinyl linoleate, retinyl oleate, retinyl palmitate, retinyl stearate

KEY WORDS

serum; use gold fluorescent lamps; hold sample at 4° before injection

REFERENCE

Sowell, A.L.; Huff, D.L.; Yeager, P.R.; Caudill, S.P.; Gunter, E.W. Retinol, α -tocopherol, lutein/zeaxanthin, β -cryptoxanthin, lycopene, α -carotene, trans- β -carotene, and four retinyl esters in serum determined simultaneously by reversed-phase HPLC with multiwavelength detection, *Clin. Chem.*, **1994**, 40, 411–416.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Serum + 2.5 mL EtOH, mix for 5 min, add 5 mL n-hexane, mix vigorously, centrifuge at 2000 g for 5 min, repeat extraction with 3 mL n-hexane. Combine the n-hexane layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in dichloromethane, inject an aliquot.

HPLC VARIABLES

Guard column: μ Bondapak C18 Guard-Pak + 50 mm long C18 ODS(4) (Shimadzu)

Column: 250 \times 4.6 5 μ m Vydac 201 TP 54 C18

Mobile phase: MeOH:MeCN 90:10 (Every 100 injections wash column with MeOH:MeCN: dichloromethane 8:1:1.)

Flow rate: 1

Detector: UV 451

CHROMATOGRAM

Retention time: 18.5 (all-*trans*), 20.7 (9-*cis*), 21.8 (15-*cis*), 19.7 (9,15-*dici*s)

OTHER SUBSTANCES

Extracted: vitamin A (UV 324), vitamin E (UV 291)

KEY WORDS

serum

REFERENCE

Ben-Amotz, A. Simultaneous profiling and identification of carotenoids, retinols, and tocopherols by high performance liquid chromatography equipped with three-dimensional photodiode array detection, *J.Liq.Chromatogr.*, **1995**, 18, 2813–2825.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Serum + 200 μ L 650 ng/mL tocopheryl acetate in MeOH, vortex for 30 s, add 200 μ L n-hexane, shake for 15 min, centrifuge at 3000 rpm for 10 min. Remove 120 μ L of the hexane layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 20 μ L dichloromethane, add 100 μ L MeCN:MeOH 50:50, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.5 μ m Nucleosil C18

Mobile phase: MeCN:MeOH:dichloromethane 50:45:5

Flow rate: 0.7

Injection volume: 50

Detector: UV 450

CHROMATOGRAM

Retention time: 18.4 (all trans), 19.5 (15,15'-cis)

Internal standard: tocopheryl acetate (F ex 295 em 330) (9.3)

OTHER SUBSTANCES

Extracted: vitamin A (F ex 325 em 480), vitamin E (F ex 295 em 330), γ -tocopherol (F ex 295 em 330), retinyl palmitate (F ex 325 em 480), α -carotene

KEY WORDS

serum

REFERENCE

Yakushina, L.; Taranova, A. Rapid HPLC simultaneous determination of fat-soluble vitamins, including carotenoids, in human serum, *J.Pharm.Biomed.Anal.*, **1995**, 13, 715–718.

SAMPLE

Matrix: blood, tissue

Sample preparation: Plasma. 100 μ L Plasma + 100 μ L 0.9% NaCl + 200 μ L MeOH, vortex for 30 s, let stand for 10 min, add 400 μ L chloroform, vortex for 4 min. Remove the chloroform layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in MeOH, inject an aliquot. Tissue. Homogenize (liver with Mikro-dismembrator II in liquid nitrogen; other tissue with Ultra-Turrax) tissue with 3 mL 1% acetic acid containing 1 mg/mL ascorbic acid and 10 mM EDTA, add 2 mL MeOH, vortex for 30 s, let stand for 10 min, add 4 mL chloroform, vortex for 4 min. Remove the chloroform layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in MeOH, inject an aliquot.

HPLC VARIABLES

Column: 125 \times 4.6 μ m Hypersil ODS

Mobile phase: MeCN:dichloromethane:MeOH:water 70:10:15:5

Flow rate: 0.5 for 13 min, to 1 over 1 min, maintain at 1 for 10 min, to 1.5 over 1 min, maintain at 1.5 for 21 min, to 2 over 1 min, maintain at 2 for 10 min, return to 0.5 over 1 min, maintain at 0.5 for 2 min.

Injection volume: 50

Detector: UV 445

CHROMATOGRAM

Retention time: 39

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Extracted: vitamin A (UV 350), vitamin E (UV 292)

KEY WORDS

rat; protect from light; liver; plasma; lung

REFERENCE

Van Vliet,T.; Van Schaik,F.; Van Schoonhoven,J.; Schrijver,J. Determination of several retinoids, carotenoids and E vitamers by high-performance liquid chromatography. Application to plasma and tissues of rats fed a diet rich in either β -carotene or canthaxanthin, *J.Chromatogr.*, **1991**, 553, 179–186.

SAMPLE

Matrix: cheese

Sample preparation: 500 mg Cheese + 2 mL 60% KOH + 2 mL 95% EtOH + 1 mL 1% NaCl + 5 mL 6% pyrogallol in EtOH, flush tube with nitrogen, seal, heat at 70° for 30 min, cool in ice water, add 15 mL 1% NaCl, extract twice with 15 mL portions of n-hexane:ethyl acetate 90:10. Combine the organic layers and evaporate them to dryness, dissolve the residue in 2 mL mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Ultrasphere Si

Mobile phase: Gradient. A was n-hexane:isopropanol 99:1. B was n-hexane. A:B 50:50 for 7 min; to 90:10 over 4 min, maintain at 90:10 for 7 min, to 50:50 over 1 min, maintain at 50:50 for 4 min. (About every 100 injections recondition column with 50 mL dichloromethane, 50 mL isopropanol, and 50 mL dichloromethane.)

Flow rate: 1.5

Injection volume: 20

Detector: UV 450 (β -carotene) and F ex 325 em 475 for 3.5 min, ex 280 em 475 for 10.5 min, ex 325 em 475 for 9 min (others)

CHROMATOGRAM

Retention time: 2

Limit of detection: 0.16 ng

OTHER SUBSTANCES

Extracted: vitamin E (α -tocopherol), vitamin A (all-trans-retinol), β -tocopherol, gamma-tocopherol, δ -tocopherol, 13-cis-retinol

KEY WORDS

normal phase; cheese

REFERENCE

Panfili,G.; Manzi,P.; Pizzoferrato,L. High-performance liquid chromatographic method for the simultaneous determination of tocopherols, carotenes, and retinol and its geometric isomers in Italian cheeses, *Analyst*, **1994**, 119, 1161–1165.

SAMPLE

Matrix: food

Sample preparation: Dissolve 10 g margarine in 50 mL dichloromethane, add 3 g anhydrous magnesium sulfate, let stand for 2 h with frequent agitation, filter (fritted glass), make up filtrate to 100 mL with dichloromethane, inject a 250 μ L aliquot on to four μ Styragel 100 Å GPC columns (Waters) in series, elute with dichloromethane at 1 mL/min, monitor at 313 nm, collect the fraction corresponding to β -carotene (at 23.5-26.5 min) and evaporate it to dryness, reconstitute with mobile phase, inject a 100 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Zorbax C18

Mobile phase: MeCN:dichloromethane:MeOH 70:30:0.2

Flow rate: 1

Injection volume: 100

Detector: UV 436

CHROMATOGRAM

Retention time: 15

OTHER SUBSTANCES

Simultaneous: vitamin A palmitate

KEY WORDS

margarine

REFERENCE

Chase, G.W., Jr.; Akoh, C.C.; Eitenmiller, R.R.; Landen, W.O. Liquid chromatographic method for the concurrent analysis of sucrose polyester, vitamin A palmitate, and β -carotene in margarine, *J.Liq.Chromatogr.*, **1995**, 18, 3129-3138.

SAMPLE

Matrix: fruit, juice, oil, vegetables

Sample preparation: Squash, peach. Homogenize (Waring blender) 10-30 g finely chopped squash or peach with 100 mL acetone and 10 g Hyflo supercel for 3 min, filter. repeat extraction of the solids until all the pigment was removed, extract with petroleum ether, wash the organic layer with water, pass it through a short column containing anhydrous sodium sulfate, concentrate under reduced pressure, inject an aliquot. Orange juice, palm tree oil. Homogenize (Waring blender) 10-30 g orange juice or palm tree oil with 100 mL acetone and 10 g Hyflo supercel for 3 min, filter. repeat extraction of the solids until all the pigment was removed, extract with petroleum ether, add 10% KOH in MeOH, let stand overnight, wash the organic layer with water, pass it through a short column containing anhydrous sodium sulfate, concentrate under reduced pressure, inject an aliquot.

HPLC VARIABLES

Column: 250 \times 2.2 Ca(OH)₂ laboratory packed (Ca(OH)₂ from Nakarri Chemicals.)

Mobile phase: Isooctane (squash, palm tree oil) or Gradient. Isooctane:acetone from 100:0 to 80:20 over 1 h. (orange juice, peach)

Flow rate: 0.5

Injection volume: 10

Detector: UV 450

CHROMATOGRAM

Retention time: 10 (13-cis), 15 (trans), 22 (9-cis) (isocratic mobile phase)

KEY WORDS

squash; peach; orange juice; palm tree oil; normal phase

REFERENCE

Carvalho, C.R.L.; Carvalho, P.R.N.; Collins, C.H. High-performance liquid chromatographic determination of the geometrical isomers of β -carotene in several foodstuffs, *J.Chromatogr.A*, **1995**, 697, 289-294.

SAMPLE**Matrix:** solutions**Sample preparation:** Inject a 2-25 μL aliquot of a solution on mobile phase.

HPLC VARIABLES**Column:** 250 \times 4.6 5 μm C30 silica (Anal.Chem. 1994, 66, 1667)**Mobile phase:** MTBE:MeOH 89:11**Flow rate:** 1**Injection volume:** 2-25**Detector:** UV 453

CHROMATOGRAM**Retention time:** 50 (all-trans)

OTHER SUBSTANCES**Simultaneous:** isomers

REFERENCE

Emenhiser,C.; Sander,L.C.; Schwartz,S.J. Capability of a polymeric C₃₀ stationary phase to resolve *cis-trans* carotenoid isomers in reversed-phase liquid chromatography, *J.Chromatogr.A*, **1995**, 707, 205–216.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 150 \times 4.6 3 μm Adsorbosphere-HS C18**Mobile phase:** MeCN:isopropanol:MeOH 60:30:10 containing 0.1% ammonium acetate**Flow rate:** 1**Detector:** UV 234, UV 295, UV 450

OTHER SUBSTANCES**Simultaneous:** vitamin E

REFERENCE

Maitra,I.; Marcocci,L.; Droy-Lefaix,M.T.; Packer,L. Peroxyl radical scavenging activity of *Ginkgo biloba* extract EGb 761, *Biochem.Pharmacol.*, **1995**, 49, 1649–1655.

SAMPLE**Matrix:** vegetables

Sample preparation: Macerate vegetables (Waring blender), remove 5 g of this material and homogenize it with 10 mL EtOH (Biohomogenizer) for 3 min, add 2 mL pentane, homogenize for 2 min (purge homogenizer motor with nitrogen), centrifuge at 7000 g for 3 min, remove pentane layer, add 2 g NaCl and 5 mL water to lower layer, shake gently, add 8 mL pentane, shake vigorously for 2 min, centrifuge, remove pentane layer. Combine pentane layers and evaporate (if necessary) under a stream of nitrogen to reduce the volume, make up to 5 or 10 mL with pentane, inject a 20 μL aliquot.

HPLC VARIABLES**Column:** 250 \times 4.6 5 μm Vydac 201 TP54 C18**Mobile phase:** MeOH:MeCN:dichloromethane:hexane 65:27:4:4**Flow rate:** 1.5**Injection volume:** 20**Detector:** UV 450

CHROMATOGRAM**Retention time:** 14

KEY WORDS

details for SFE also given; squash; broccoli; carrots; collard greens; turnip; kale; mustard greens; zucchini

REFERENCE

Marsili,R.; Callahan,D. Comparison of a liquid solvent extraction technique and supercritical fluid extraction for the determination of α - and β -carotene in vegetables, *J.Chromatogr.Sci.*, **1993**, 31, 422–428.

SAMPLE

Matrix: vegetables

Sample preparation: Blanch vegetables in water for 90 s. Stir 20 g homogenized and mashed vegetables with 50 mL ethyl ether for about 4 h, filter, extract the residue again. Combine the filtrates, add 5 mL 20% KOH in MeOH, mix, place in a refrigerator for about 12 h, add 100 mL water, shake vigorously, let stand for about 5 h. Remove the ether layer and concentrate it under vacuum, evaporate traces of ether with a stream of nitrogen, reconstitute with EtOH (?), inject an aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Wakosil-II 5C18 AR polymeric ODS (Wako Chemicals)

Mobile phase: THF:MeOH 10:90

Column temperature: 20

Flow rate: 1

Injection volume: 100

Detector: UV 450

CHROMATOGRAM

Retention time: 15

OTHER SUBSTANCES

Extracted: α -carotene, β -cryptoxanthin, lutein, lycopene, zeaxanthin


KEY WORDS

carrots; tomatoes; pumpkins

REFERENCE

Jinno,K.; Lin,Y. Separation of carotenoids by high-performance liquid chromatography with polymeric and monomeric octadecylsilica stationary phases, *Chromatographia*, **1995**, 41, 311–317.

Betahistine



Molecular formula: C₈H₁₂N₂

Molecular weight: 136.20

CAS Registry No.: 5638-76-6, 5579-84-0 (2HCL)

Merck Index: 1224

Lednicer No.: 2 279

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A: B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 3.155


KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

Betaine



Molecular formula: $C_5H_{11}NO_2$

Molecular weight: 117.15

CAS Registry No.: 107-43-7, 590-46-5 (HCl), 590-47-6 (monohydrate), 93227-64-6 (sodium aspartate)

Merck Index: 1225

SAMPLE

Matrix: solutions

Sample preparation: Mix 600 μ L solution with 60 μ L IS solution, vortex, add a 550 μ L aliquot to a 70×5 column of Dowex 1-X8 (OH⁻) form in a Pasteur pipette, elute with 2 mL water, evaporate the eluate to dryness under a stream of air, reconstitute with 100 μ L 2 mM N,N-diisopropylethylamine in MeCN, add 100 μ L 5 mM 4'-bromophenacyl trifluoromethanesulfonate in MeCN, vortex for 10 min, add 100 μ L 10 mM hydroxyacetic acid N,N-diisopropylethylammonium salt in MeCN, vortex for 2 min, inject a 15 min aliquot. (Prepare 4'-bromophenacyl trifluoromethanesulfonate as follows. Add 8.8 g p-bromobenzoyl chloride in 40 mL dry ether over 20-30 min to 100 mmoles diazomethane stirred in an ice bath, stir in an ice bath for 8-9 h, let stand at room temperature for 3 h, evaporate the solvent under reduced pressure, recrystallize 4'-bromo-2-diazoacetophenone from ether/hexane (mp 123.5-124° d) (J. Am. Chem. Soc. 1951, 73, 5301). Condense 50 mL anhydrous sulfur dioxide in a flask fitted with a calcium sulfate drying tube, cool in a dry ice/acetone bath, add 2.25 g 4'-bromo-2-diazoacetophenone, stir for 5 min, add 900 μ L anhydrous trifluoromethanesulfonic acid from a freshly opened bottle in one portion, stir for 15 min, remove the cooling bath, after 30 min use an ice/water bath to evaporate the solvent. Dissolve the residue in 100 mL boiling dichloromethane, treat twice with 5 g portions of decolorizing carbon, filter, evaporate the filtrate, recrystallize the residue from pentane:dichloromethane 80:20 to give 4'-bromophenacyl trifluoromethanesulfonate as colorless plates (mp 137-8°) (J. Chromatogr. 1984, 299, 365).)

HPLC VARIABLES

Guard column: 50 \times 4 Co:Pell ODS

Column: 100 \times 8 5 μ m Radial-PAK C18

Mobile phase: MeCN:buffer 70:30 (Prepare by dissolving 70 mg sodium dodecyl sulfate, 140 mg $NaH_2PO_4 \cdot H_2O$, and 300 μ L 3-dimethylamino-1-propanol in 150 mL water, adjusting pH to 6.5 with 85% phosphoric acid, and adding 350 mL MeCN.)

Flow rate: 5

Injection volume: 15

Detector: UV 254

CHROMATOGRAM

Retention time: 7

Internal standard: (4-bromophenyl)carboxymethyl (6-trimethylammonium)hexanoate (Add an aqueous solution of 6-(trimethylammonium)hexanoic acid (Cl⁻ salt, prepared by methylation of 6-aminohexanoic acid) to a column of Dowex 1-X8 (OH⁻) form, elute with 4 column volumes of water. Evaporate the eluate to dryness, dissolve the residue in DMF, add 1.1 equivalents of triethylamine, add 1.1 equivalents of 2,4'-dibromoacetophenone, stir at 40° for 3 h, add ethyl acetate. Collect the precipitate by filtration and recrystallize it from ethanol/acetone to obtain (4-bromophenyl)carboxymethyl (6-trimethylammonium)hexanoate.) (8)

Limit of quantitation: 10 μ M

KEY WORDS

derivatization

REFERENCE

Minkler,P.E.; Ingalls,S.T.; Kormos,L.S.; Weir,D.E.; Hoppel,C.L. Determination of carnitine, butyrobetaine, and betaine as 4'-bromophenacyl ester derivatives by high-performance liquid chromatography, *J.Chromatogr.*, **1984**, 336, 271-283.