

Budesonide


Molecular formula: C₂₅H₃₄O₆

Molecular weight: 430.54

CAS Registry No.: 51333-22-3 (β,16α), 51372-29-3 (11β,16α[R]), 51372-28-2 (11β, 16α[S])

Merck Index: 1490

Lednicer No.: 3 95



SAMPLE

Matrix: blood

Sample preparation: Add 5 mL diethyl ether/dichloromethane (ratio not given) to 1 mL serum, mix for 30 min. Centrifuge and freeze at -80° for 15 min. Evaporate the supernatant under a stream of nitrogen, reconstitute the residue in 200 µL mobile phase. Inject a 50 µL aliquot.

HPLC VARIABLES

Column: Lichrospher RP Select B

Mobile phase: MeCN:20 mM ammonium acetate buffer 80:20

Flow rate: 1

Injection volume: 50

Detector: MS, PE-Sciex API 300, negative ion mode

CHROMATOGRAM

Internal standard: budesonide

OTHER SUBSTANCES

Extracted: flunisolide

KEY WORDS

serum; budesonide is IS

REFERENCE

Möllmann,H.; Derendorf,H.; Barth,J.; Meibohm,B.; Wagner,M.; Krieg,M.; Weisser,H.; Knöller,J.; Möllmann,A.; Hochhaus,G. Pharmacokinetic/pharmacodynamic evaluation of systemic effects of flunisolide after inhalation, *J.Clin.Pharmacol.*, **1997**, 37, 893–903.

SAMPLE

Matrix: blood

Sample preparation: Condition a C18 SPE cartridge twice with 3 mL portions of EtOH and twice with 3 mL portions of water. Add 1 mL EtOH:water 30:70 to 1 mL plasma, vortex, let stand for 15 min. Centrifuge at 1800 g for 15 min. Add the supernatant dropwise to the SPE cartridge, wash with 3 mL EtOH:water 25:75, 3 mL water, and twice with 2 mL ethyl acetate:n-heptane 2:98. Dry the SPE cartridge under vacuum, elute with 2 mL ethyl acetate:heptane 35:65. Evaporate the eluate to dryness under a stream of nitrogen at 35°. Reconstitute the residue with 100 µL mobile phase, let stand 15 min. Inject a 20 µL aliquot.

HPLC VARIABLES

Guard column: 20 × 2.1 5 µm Hypersil ODS C18

Column: 100 × 2.1 5 µm Hypersil ODS

Mobile phase: EtOH:water 50:50

Flow rate: 0.45

Injection volume: 20

Detector: MS, Finnigan Mat TSQ 7000, ESI mode, positive ion APCI mode, m/z 473.2

CHROMATOGRAM**Retention time:** 3.35 (budesonide acetate 22R epimer)**Internal standard:** budesonide acetate 22R epimer

OTHER SUBSTANCES**Extracted:** fluticasone

KEY WORDSbudesonide is IS; plasma; SPE

REFERENCE

Li,Y.N.; Tattam,B.N.; Brown,K.F.; Seale,J.P. A sensitive method for the quantification of fluticasone propionate in human plasma by high-performance liquid chromatography/atmospheric pressure chemical ionisation mass spectrometry, *J.Pharm.Biomed.Anal.*, **1997**, 16, 447–452.

SAMPLE**Matrix:** blood

Sample preparation: Condition a 6 mL 500 mg Extract Clean C18 SPE cartridge (Alltech) with two 3 mL portions of EtOH and two 3 mL portions of water. Add 50 µL 200 ng/mL IS solution and 1 mL EtOH:water 30:70 to 1 mL plasma, vortex carefully, let stand for 15 min, centrifuge at 1200 g for 20 min. Add the supernatant to the SPE cartridge attached to a vacuum manifold operating at 5.1 kPa, wash with 3 mL EtOH:water 25:75, 3 mL water, and 2 mL heptane:ethyl acetate 98:2. Elute with 2 mL heptane:ethyl acetate 65:35 and with 2 mL EtOH, evaporate the eluate to dryness under a stream of nitrogen at 35°, add 100 µL MeCN:triethylamine:acetic anhydride 75:12.5:12.5 to the residue, let react for 15 min, evaporate to dryness under a stream of nitrogen, reconstitute the residue with 100 µL mobile phase, let stand for at least 15 min, inject a 20 µL aliquot.

HPLC VARIABLES**Guard column:** 20 × 2.1 5 µm ODS Hypersil C18**Column:** 100 × 2.1 5 µm ODS Hypersil C18**Mobile phase:** EtOH:water 43:57**Flow rate:** 0.5**Injection volume:** 20**Detector:** MS, Finnigan MAT TSQ 7000, APCI source, m/z 431, 476

CHROMATOGRAM**Retention time:** 8.0 (22R), 9.0 (22S)**Internal standard:** ²H₃ budesonide-21-acetate (preparation described in paper) (9-10.5)**Limit of detection:** 700fg**Limit of quantitation:** 250 pg/mL

KEY WORDSderivatization; plasma; pharmacokinetics; SPE

REFERENCE

Li,Y.N.; Tattam,B.; Brown,K.F.; Seale,J.P. Determination of epimers 22R and 22S of budesonide in human plasma by high-performance liquid chromatography-atmospheric pressure chemical ionization mass spectrometry, *J.Chromatogr.B*, **1996**, 683, 259–268.

SAMPLE**Matrix:** blood

Sample preparation: 500 µL Plasma + 1.5 mL water + 4 mL dichloromethane, shake gently for 30 min, centrifuge at 2000 rpm for 15 min. Remove 3 mL of the organic layer and evaporate it to dryness under a stream of nitrogen with gentle heating, reconstitute the residue in 200 µL mobile phase, inject a 100 µL aliquot.

HPLC VARIABLES

Column: 100 × 3.9 10 µm Spherisorb ODS

Mobile phase: EtOH:water 40:60

Flow rate: 0.5

Injection volume: 100

Detector: UV 254

CHROMATOGRAM

Retention time: 7

KEY WORDS

plasma; dog; radiolabeled; pharmacokinetics

REFERENCE

Ryrfeldt, Å.; Tönnesson, M.; Nilsson, E.; Wikby, A. Pharmacokinetic studies of a potent glucocorticoid (budesonide) in dogs by high-performance liquid chromatography, *J. Steroid Biochem.*, **1979**, *10*, 317–324.

SAMPLE

Matrix: blood

Sample preparation: 250–500 µL blood + 10 mL EtOH:10 mM acetic acid 10:90, mix, add 9 mL to a conditioned Sep-Pak C18 SPE cartridge, wash with 10 mL EtOH:10 mM acetic acid 8:92, elute with 4 mL EtOH:10 mM acetic acid 60:40. Dilute the eluate to an EtOH content of 8% and add it to another conditioned C18 SPE cartridge, wash with 10 mL EtOH:10 mM acetic acid 8:92, elute with 4 mL EtOH. Evaporate the eluate to dryness under a stream of nitrogen at 40°, reconstitute the residue in 250 µL mobile phase, inject a 100 µL aliquot.

HPLC VARIABLES

Column: 250 × 5 5 µm Nucleosil NO2 + 200 × 5 5 µm Nucleosil NO2 in series

Mobile phase: dichloromethane:isopropanol:water 99:1:0.2 for 35 min then 75:25:0.2 for 15 min, re-equilibrate at initial conditions for 45 min

Flow rate: 3

Injection volume: 100

Detector: UV 254

CHROMATOGRAM

Retention time: 22–28

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

mouse; pharmacokinetics; SPE

REFERENCE

Andersson, P.; Appelgren, L.-E.; Ryrfeldt, Å. Tissue distribution and fate of budesonide in the mouse, *Acta Pharmacol. Toxicol. (Copenh)*, **1986**, *59*, 392–402.

SAMPLE

Matrix: blood

Sample preparation: Condition a Bond Elut C18 SPE cartridge. Add plasma + IS to SPE cartridge, wash with aqueous EtOH, wash with water, wash with heptane, elute with ethyl acetate in heptane, esterify with acetic anhydride and triethylamine in MeCN, evaporate, reconstitute in mobile phase, inject an aliquot.

HPLC VARIABLES

Guard column: 10 × 3 Chromguard

Column: 33 × 4.6 3 μm Supelcosil LC-8-DB

Mobile phase: MeOH:100 mM pH 5 ammonium acetate 64:36

Flow rate: 1.4

Injection volume: 100

Detector: MS, Finnigan 4500 quadrupole, thermospray, scan time 40 ms, source block 220°, repeller 45 V, vaporizer 105°, jet block 180°, aerosol 220

CHROMATOGRAM

Retention time: 3.6 (as budesonide 21-acetate)

Internal standard: octadeutero budesonide

Limit of quantitation: 0.1 nM

KEY WORDS

plasma; LC-MS; SPE; derivatization

REFERENCE

Lindberg,C.; Paulson,J.; Blomqvist,A. Evaluation of an automated thermospray liquid chromatography-mass spectrometry system for quantitative use in bioanalytical chemistry, *J.Chromatogr.*, **1991**, *554*, 215–226.

SAMPLE

Matrix: blood

Sample preparation: Centrifuge plasma at 2500 g for 10 min, mix the supernatant with an equal volume of 1 M pH 2.5 glycine buffer containing 0.2% Tween 20, centrifuge at 2500 g for 10 min, inject a 50 μL aliquot of the supernatant on to column A and elute to waste with mobile phase A, after 1.4 min backflush the contents of column A on to column B with mobile phase B, after 3.6 min remove column A from the circuit, elute column B with mobile phase B and monitor the effluent from column B. Flush column A with mobile phase A for 9 min and mobile phase B for 16 min.

HPLC VARIABLES

Column: A 30 × 2.1 Apex II aminopropyl (Jones Chromatography); B 50 × 2.1 Spherisorb C1 pH stable

Mobile phase: A 10 mM pH 2.5 glycine buffer containing 0.1% Tween 20; B Isopropanol: 10 mM pH 2.5 glycine buffer containing 0.1% Tween 20 3:97

Column temperature: 40

Flow rate: 0.5

Injection volume: 50

Detector: UV

CHROMATOGRAM

Retention time: 22, 25 (epimers)

KEY WORDS

plasma; column-switching

REFERENCE

Lövgren,U.; Johansson,M.; Kronkvist,K.; Edholm,L.-E. Biocompatible sample pretreatment for immunochemical techniques using micellar liquid chromatography for separation of corticosteroids, *J.Chromatogr.B*, **1995**, *672*, 33–44.

SAMPLE

Matrix: blood, tissue

Sample preparation: Acidify plasma or lung tissue homogenate to pH 2 with 500 mM HCl, add 100 μL 20 μg/mL IS, extract with 8 mL dichloromethane. Evaporate the organic

layer to dryness under vacuum, reconstitute in 120 μ L MeOH:5% acetic acid 50:50, inject an 80 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Zorbax ODS C18

Mobile phase: MeCN:MeOH:water 44:11:45

Flow rate: 1

Injection volume: 80

Detector: UV 242 or radioactivity

CHROMATOGRAM

Internal standard: hydrocortisone 21-S-propionate (JO 498)

OTHER SUBSTANCES

Extracted: metabolites, beclomethasone dipropionate

KEY WORDS

plasma; rat; lung; radiolabeled; pharmacokinetics; epimers are separated

REFERENCE

Chanoine,F.; Grenot,C.; Heidmann,P.; Junien,J.L. Pharmacokinetics of butixocort 21-propionate, budesonide, and beclomethasone dipropionate in the rat after intratracheal, intravenous, and oral treatments, *Drug Metab.Dispos.*, **1991**, 19, 546–553.

SAMPLE

Matrix: broncho-alveolar lavage fluid

Sample preparation: Centrifuge bronchoalveolar lavage fluid at 1150 g for 10 min. 1 mL Supernatant + 5 mL dichloromethane, shake on an alternating agitator for 10 min, centrifuge at 1150 g for 10 min, repeat the extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 50°, reconstitute the residue in 100 μ L mobile phase, centrifuge, inject an 80 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Spherisorb ODS C18

Mobile phase: MeOH:buffer 69:31 (Buffer was 0.1% acetic acid, pH 3.)

Flow rate: 1

Injection volume: 80

Detector: UV 250

CHROMATOGRAM

Retention time: 10.51

Limit of quantitation: 5 ng/mL

REFERENCE

Faouzi,M.A.; Dine,T.; Luyckx,M.; Brunet,C.; Gressier,B.; Cazin,M.; Wallaert,B. High-performance liquid chromatographic method for the determination of budesonide in bronchoalveolar lavage of asthmatic patients, *J.Chromatogr.B*, **1995**, 664, 463–467.

SAMPLE

Matrix: bulk

Sample preparation: Inject a 5 μ L aliquot of a solution in EtOH.

HPLC VARIABLES

Column: 300 \times 3.9 μ Bondapak C18

Mobile phase: EtOH:water 48:52 or 43:57

Flow rate: 1

Injection volume: 5

Detector: UV 254

CHROMATOGRAM

Retention time: 10, 11 (48:52), 18, 25 (43:57) (two epimers)

OTHER SUBSTANCES

Simultaneous: impurities, 16 α -hydroxyprednisolone

REFERENCE

Roth,G.; Wikby,A.; Nilsson,L.; Thalén,A. High-performance liquid chromatographic determination of epimers, impurities, and content of the glucocorticoid budesonide and preparation of primary standard, *J.Pharm.Sci.*, **1980**, 69, 766–770.

SAMPLE

Matrix: solutions

Sample preparation: Mix 250 μ L of a solution of budesonide in MeOH with 625 μ L 9% trifluoromethanesulfonic acid in MeOH, immediately add 625 μ L 160 μ m dansyl hydrazine in MeOH, after 12 h inject a 100 μ L aliquot onto column A and elute to waste with mobile phase A, after 100 s backflush the contents of column A onto column B with mobile phase B, elute with mobile phase B, monitor the effluent from column B. (Caution! Trifluoromethanesulfonic acid is highly toxic! Decontaminate surplus solutions with an equal volume of 2 M aqueous ammonia for 1 h!)

HPLC VARIABLES

Column: A 10 \times 4.6 5 μ m Vydac C18; B 150 \times 4.6 5 μ m Nucleosil C18

Mobile phase: A MeCN:7.7 mM pH 7 phosphate buffer 6:94; B MeCN:7.7 mM pH 7 phosphate buffer 65:35

Flow rate: 1

Injection volume: 100

Detector: F ex 350 em 520

CHROMATOGRAM

Retention time: 8, 10 (syn and anti isomers of 3-keto derivative)

Limit of detection: 1.5 pmole

KEY WORDS

derivatization; column-switching

REFERENCE

Hyytiäinen,M.; Appelblad,P.; Pontén,E.; Stigbrand,M.; Irgum,K.; Jaegfeldt,H. Trifluoromethanesulfonic acid as a catalyst for the formation of dansylhydrazone derivatives, *J.Chromatogr.A*, **1996**, 740, 279–283.

SAMPLE

Matrix: tissue

Sample preparation: Homogenize tissue in ice cold saline for 30 s. 100 μ L Homogenate + 4 mL EtOH:10 mM acetic acid 70:30, shake for 30 min, centrifuge at 1000 g for 15 min. Remove the supernatant and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 250 μ L EtOH:water 30:70, add 5 mL 10 mM acetic acid, add to a conditioned Sep-Pak C18 SPE cartridge, wash with 5 mL 10 mM acetic acid, elute with 5 mL EtOH:water 70:30. Evaporate the eluate to dryness under a stream of nitrogen at 40°, reconstitute in 250 μ L EtOH:water 30:70, inject a 200 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 5 5 μ m Nucleosil C18

Mobile phase: EtOH:water 50:50 for 12 min then EtOH for 10 min

Flow rate: 1

Injection volume: 200

Detector: UV 254

CHROMATOGRAM

Retention time: 11-13

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

mouse; pharmacokinetics; liver; lung; spleen; brain; SPE

REFERENCE

Andersson,P.; Appelgren,L.-E.; Ryrfeldt,Å. Tissue distribution and fate of budesonide in the mouse, *Acta Pharmacol.Toxicol.(Copenh)*, **1986**, 59, 392-402.


Buflomedil

Molecular formula: C₁₇H₂₅NO₄

Molecular weight: 307.39

CAS Registry No.: 55837-25-7, 35543-24-9 (HCl)

Merck Index: 1498



SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A: B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 201.7

CHROMATOGRAM

Retention time: 11.3

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.


Bufotenine

Molecular formula: C₁₂H₁₆N₂O

Molecular weight: 204.27

CAS Registry No.: 487-93-4

Merck Index: 1502



SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 125 × 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 3.5

OTHER SUBSTANCES


Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazepine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipanone, diprenorphine, dipyrizidamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypropazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phen-tolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, pir-tramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolin-tane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, pro-thipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine,

quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocinide, tolpropamine, tolycaine, tran-ylcypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleennamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J.Chromatogr.*, **1985**, *323*, 191-225.

Bufuralol



Molecular formula: C₁₆H₂₃NO₂, C₁₆H₂₃NO₂.HCl (hydrochloride)

Molecular weight: 261.36

CAS Registry No.: 54340-62-4, 59652-29-8 (HCl), 57704-10-6 ((-), HCl), 57704-11-7 ((+), HCl)

Merck Index: 1504

SAMPLE

Matrix: blood

Sample preparation: Mix 1 mL plasma with 500 µL immobilized antibody (details of preparation in paper), wash with 5 mL 1 M NaCl, wash with 10 mL water, elute with 5 mL MeOH:10 mM pH 5.0 ammonium acetate buffer 95:5, inject an aliquot of the eluate.

HPLC VARIABLES

Column: 150 × 4.6 5 µm Ultron ES-OVM (Shinwa, Osaka, Japan)

Mobile phase: MeCN:0.3% pH 6.7 ammonium acetate buffer from 1:14 to 1:2 over 40 (?) min

Flow rate: 1

Detector: UV 248

CHROMATOGRAM

Retention time: 32

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; chiral separation of the metabolites but not necessarily of bufuralol

REFERENCE

Ikegawa,S.; Matsuura,K.; Sato,T.; Made,N.; Isriyanthi,R.; Niwa,T.; Miyairi,S.; Takashima,H.; Kawashima,Y.; Goto,J. Enantioselective immunoaffinity extraction for simultaneous determination of optically active bufuralol and its metabolites in human plasma by HPLC, *J.Pharm.Biomed.Anal.*, **1998**, 17, 1-9.


Bumadizon

Molecular formula: C₁₉H₂₂N₂O₃

Molecular weight: 326.40

CAS Registry No.: 3583-64-0, 69365-73-7 (calcium salt hemihydrate)

Merck Index: 1507



SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform: isopropanol:n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 235

CHROMATOGRAM

Retention time: 6.72

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds(all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; tolaxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vindesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loperazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydra-

mine; cyclizine; histapyrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254–262.

Bumetanide


Molecular formula: C₁₇H₂₀N₂O₅S

Molecular weight: 364.42

CAS Registry No.: 28395-03-1

Merck Index: 1508

Lednicer No.: 2 87



SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 30 µL 500 ng/mL methylbumetanide in water + 200 µL 1 M sulfuric acid + 6 mL diethyl ether, shake 20 min, centrifuge at 1000 g for 5 min. Evaporate 4 mL of ether layer at 30° under a stream of nitrogen. Dissolve residue in 150 µL mobile phase, inject a 100 µL aliquot.

HPLC VARIABLES

Guard column: 22 × 3.9 37-50 µm Corasil C18

Column: 75 × 4.6 3 µm Supelco LC-8-DB

Mobile phase: MeCN:30 mM pH 3.0 sodium phosphate buffer 125:200

Flow rate: 1

Injection volume: 100

Detector: F ex 340 em 440

CHROMATOGRAM

Retention time: 4.5

Internal standard: methylbumetanide (6.5)

Limit of detection: 0.1 ng/mL

KEY WORDS

plasma

REFERENCE

Bökens,H.; Bourscheidt,C.; Müller,R.F. Determination of bumetanide in plasma by high-performance liquid chromatography, *J.Chromatogr.*, **1988**, 434, 327-329.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 500 ng naproxen + 1 mL 100 mM HCl + 10 mL dichloromethane, extract. Dry organic layer at 50° under nitrogen, dissolve in 1 mL mobile phase, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 150 × 4.6 3 µm Alltech C8

Mobile phase: MeCN:80 mM phosphoric acid 35:65

Flow rate: 1

Injection volume: 20

Detector: F ex 235 em 405

CHROMATOGRAM

Retention time: 12.3

Internal standard: naproxen (3.5)

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Extracted: furosemide

KEY WORDS

plasma; horse; pharmacokinetics

REFERENCE

Singh,A.K.; McArdle,C.; Gordon,B.; Ashraf,M.; Granley,K. Simultaneous analysis of furosemide and bumetanide in horse plasma using high performance liquid chromatography, *Biomed.Chromatogr.*, **1989**, 3, 262-265.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Plasma + piretanide + 1 mL MTBE + 100 μ L 1 M pH 4 phosphate buffer, extract for 5 min, centrifuge. Remove 850 μ L of the organic layer and evaporate it to dryness, reconstitute the residue in 180 μ L mobile phase, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 50 \times 4.6 LC-8-DB (Supelco)

Mobile phase: MeCN:50 mM pH 3 NaH_2PO_4 30:70

Flow rate: 2

Injection volume: 50

Detector: F (wavelengths not specified)

CHROMATOGRAM

Retention time: 6.2

Internal standard: piretanide (3.9)

Limit of detection: 0.4 ng/mL

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Oberbauer,R.; Krivanek,P.; Turnheim,K. Pharmacokinetics and pharmacodynamics of the diuretic bumetanide in the elderly, *Clin.Pharmacol.Ther.*, **1995**, 57, 42-51.

SAMPLE

Matrix: blood, tissue

Sample preparation: Plasma. 200 μ L Plasma + 50 μ L 250 μ g/mL acetophenone in MeCN: water 50:50, vortex, add 400 μ L MeCN, vortex, sonicate for 2 min, centrifuge for 10 min, inject an aliquot of the supernatant. Tissue. Homogenize (Tissuemizer) rat stomach, kidney, or liver with three volumes cold 250 mM sucrose, centrifuge at 9000 g. Incubate 1 mL supernatant with bumetanide, add 1 mL 1 M NaOH. 100 μ L Mixture + 250 μ L MeCN + 100 μ L 0.3% HCl, vortex, centrifuge, inject an aliquot of the supernatant (from Bio-pharm. Drug Dispos. 1991, 12, 311).

HPLC VARIABLES

Column: 250 \times 4.6 Partisil-10 ODS-3

Mobile phase: MeOH:water:acetic acid 70:30:1

Flow rate: 1.5

Detector: F ex 338 em 433 or UV 254

CHROMATOGRAM

Retention time: 5.5 (F)

Internal standard: acetophenone (4.0, UV)

Limit of detection: 5 ng/mL

OTHER SUBSTANCES

Noninterfering: metabolites

KEY WORDS

plasma; human; rat; pharmacokinetics; stomach; kidney; liver; rabbit (J.Pharm.Sci. 1995; 84; 236)

REFERENCE

Smith,D.E. High-performance liquid chromatographic assay for bumetanide in plasma and urine, *J.Pharm.Sci.*, **1982**, 71, 520–523.

SAMPLE

Matrix: blood, urine

Sample preparation: Serum. 400 μ L Serum + 100 μ L 500 ng/mL IS + 1.0 mL 1.0 M pH 5.0 potassium citrate, vortex, let sit for 1 min. Add 4.0 mL ethyl acetate:cyclohexane 70:30, mix (Rotatorque mixer) at 60 rpm for 10 min, centrifuge at 1150 g for 1 min. Evaporate upper organic layer to dryness under reduced pressure. Redissolve residue in 125 μ L mobile phase, inject a 50 μ L aliquot. Urine. Centrifuge sample at 16000 g for 1 min. Mix one volume urine supernatant with one volume 800 mM pH 3.9 potassium formate buffer, inject an aliquot.

HPLC VARIABLES

Column: 300 \times 4 Zorbax C8

Mobile phase: MeOH:30 mM pH 2.5 potassium phosphate 63:37

Column temperature: 40

Flow rate: 1

Injection volume: 50

Detector: F ex 340 em 440

CHROMATOGRAM

Retention time: 6.8

Internal standard: R021-1825 (6.3)

Limit of quantitation: 3 ng/mL (plasma), 10 ng/mL (urine)

KEY WORDS

serum; pharmacokinetics

REFERENCE

Sullivan,J.E.; Witte,M.K.; Yamashita,T.S.; Myers,C.M.; Blumer,J.L. Pharmacokinetics of bumetanide in critically ill infants, *Clin.Pharmacol.Ther.*, **1996**, 60, 405–413.

SAMPLE

Matrix: blood, urine

Sample preparation: Condition a Bond Elut 1 mL C18 SPE cartridge with 600 μ L MeCN, 2 mL MeOH, and 3 mL water. 200 μ L Plasma or urine + 20 (plasma) or 100 (urine) ng piretanide in methanol + 400 μ L MeCN, vortex 2 min, sonicate 2 min in a water bath, centrifuge 10 min 1000 g. Transfer supernatant to a clean tube and add 4 mL 100 mM pH 5.0 phosphate buffer, add to SPE cartridge, wash with 3 mL water, elute with 300 μ L MeCN. Evaporate to dryness under a stream of nitrogen at 37°, take up in 200 μ L mobile phase, inject a 100 μ L aliquot.

HPLC VARIABLES

Guard column: C18 (unspecified)

Column: 100 \times 8 5 μ m C18 Radial Pak

Mobile phase: MeOH:water:glacial acetic acid 66:34:1

Flow rate: 1.2

Injection volume: 100

Detector: F ex 228 em 418

CHROMATOGRAM

Retention time: 6.5

Internal standard: piretanide (5.1)

Limit of detection: 5 ng/mL

OTHER SUBSTANCES

Simultaneous: dopamine, amikacin, furosemide

Noninterfering: carbamazepine, phenytoin, midazolam, diazepam, lorazepam, ranitidine, gentamicin, erythromycin, ampicillin, cefotaxime, ceftazidime, ceftriaxone, vancomycin, clindamycin, theophylline, digoxin

Interfering: nafcillin

KEY WORDS

plasma; pharmacokinetics; SPE

REFERENCE

Wells, T.G.; Hendry, I.R.; Kearns, G.L. Measurement of bumetanide in plasma and urine by high-performance liquid chromatography and application to bumetanide disposition, *J. Chromatogr.*, **1991**, *570*, 235-242.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A: B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 19.055

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, *763*, 149-163.

SAMPLE

Matrix: bulk, formulations

Sample preparation: Dissolve in mobile phase: acetone 95:5 containing 100 μ g/mL salicylic acid, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 250 × 4.6 10 µm Varian RP-8

Mobile phase: MeOH:water:acetic acid 60:40:0.5

Column temperature: 25

Flow rate: 1

Injection volume: 10

Detector: UV 231

CHROMATOGRAM

Retention time: 9.3

Internal standard: salicylic acid (4.5)

Limit of detection: 2000 ng/mL

KEY WORDS

tablets; ampoules

REFERENCE

Zivanov-Stakic,D.; Solomun,L.J.; Zivanovic,L.J. High-performance liquid chromatographic method for the determination of bumetanide in pharmaceutical preparations, *J.Pharm.Biomed.Anal.*, **1989**, 7, 1889–1892.

SAMPLE

Matrix: urine

Sample preparation: Inject 5 µL urine onto column A and elute to waste with mobile phase A, after 1 min backflush the contents of column A onto column B with mobile phase B. Monitor the effluent from column B.

HPLC VARIABLES

Column: A 20 × 2.1 30 µm Hypersil ODS-C18; B 125 × 4 5 µm LiChrospher 100 RP 18

Mobile phase: A 50 mM pH 3 phosphate buffer; B MeCN:50 mM pH 3 phosphate buffer 60:40 (Prepare buffer as follows. Dissolve 3.45 g NaH₂PO₄ monohydrate in 500 mL water containing 750 µL propylamine hydrochloride, adjust to pH 3 with concentrated phosphoric acid.)

Flow rate: 1

Injection volume: 5

Detector: UV 254, F ex 228 em 418

CHROMATOGRAM

Retention time: 8.9

Limit of detection: 100 pg/mL

OTHER SUBSTANCES

Extracted: amiloride, furosemide, triamterene

KEY WORDS

column-switching; pharmacokinetics

REFERENCE

Campins-Falcó,P.; Herráez-Hernández,R.; Pastor-Navarro,M.D. Analysis of diuretics in urine by column-switching chromatography and fluorescence detection, *J.Liq.Chromatogr.Rel.Technol.*, **1997**, 20, 1867–1885.

SAMPLE

Matrix: urine

Sample preparation: 100 µL Urine + 100 µL water + 50 µL 500 µg/mL acetophenone in MeCN:water 50:50, vortex, inject an aliquot. (For rabbit (*J.Pharm.Sci.* 1995, 84, 236) or

rat (Biopharm. Drug Dispos. 1991, 12, 311) 100 μ L urine + 250 μ L MeCN, vortex, centrifuge, inject an aliquot.)

HPLC VARIABLES

Column: 250 \times 4.6 Partisil-10 ODS-3

Mobile phase: MeCN:15 mM phosphoric acid 50:50

Flow rate: 2

Injection volume: 6.5 (F)

Detector: F ex 338 em 433 or UV 254

CHROMATOGRAM

Internal standard: acetophenone (4.5, UV)

OTHER SUBSTANCES

Noninterfering: metabolites

KEY WORDS

pharmacokinetics; human; rabbit; rat

REFERENCE

Smith,D.E. High-performance liquid chromatographic assay for bumetanide in plasma and urine, *J.Pharm.Sci.*, **1982**, 71, 520-523.

SAMPLE

Matrix: urine

Sample preparation: 2 mL Urine + 2 mL 1 M pH 4.1 NaH_2PO_4 + 4 mL ethyl acetate, vortex for 2 min, centrifuge at 1500 g for 5 min. Remove the organic phase and add it to 5 mL 100 mM pH 7.5 Na_2HPO_4 , vortex for 2 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 60°, reconstitute the residue in 100 μ L MeCN:10 mM pH 3.0 phosphate buffer, inject a 5 μ L aliquot.

HPLC VARIABLES

Column: 125 \times 4 5 μ m LiChrosorb RP-18

Mobile phase: Gradient. MeCN:10 mM pH 3.0 phosphate buffer 10:90 for 1.5 min then to 35:65 over 2 min

Column temperature: 50

Flow rate: 1.5

Injection volume: 5

Detector: UV 271

CHROMATOGRAM

Retention time: 9.3

Limit of quantitation: 1500 ng/mL

OTHER SUBSTANCES

Extracted: chlorothiazide, hydrochlorothiazide, quinethazone, chlorthalidone, clopamide, methyclothiazide, furosemide, metolazone, mefruside, cyclopenthiazide, bendroflumethiazide

Simultaneous: indapamide, clorexolone, ethacrynic acid

Noninterfering: aspirin, albuterol, allopurinol, alprenolol, atenolol, captopril, carbimazole, clonidine, coloxyl, danthron, diazepam, digoxin, doxepin, glibenclamide, hydralazine, indomethacin, labetalol, metformin, methyl dopa, metoprolol, mianserin, minoxidil, nifedipine, nitrazepam, oxazepam, oxprenolol, pindolol, prazosin, propranolol, senokot, theophylline, trifluoperazine

REFERENCE

Fullinaw,R.O.; Bury,R.W.; Moulds,R.F.W. Liquid chromatographic screening of diuretics in urine, *J.Chromatogr.*, **1987**, 415, 347-356.

SAMPLE**Matrix:** urine

Sample preparation: 2 mL Urine + 0.5 g solid buffer I (pH 5-5.5), vortex 15 s, add 4 mL ethyl acetate, agitate for 10 min, centrifuge at 600 g for 5 min. Remove organic layer and vortex it with 2 mL 5% aqueous lead acetate for 10 s, centrifuge at 600 g for 5 min, remove and keep organic phase. 2 mL Urine + 0.5 g solid buffer II (pH 9-9.5), vortex 15 s, add 4 mL ethyl acetate, agitate for 10 min, centrifuge at 600 g for 5 min. Remove organic layer and combine it with previous organic layer. Evaporate to dryness at 50° under a stream of nitrogen, reconstitute in 300 µL 50 µg/mL β-hydroxyethyltheophylline in MeOH, inject 5 µL aliquot. (Solid buffer I was $\text{KH}_2\text{PO}_4\text{:Na}_2\text{HPO}_4$ 99:1, solid buffer II was $\text{NaHCO}_3\text{:K}_2\text{CO}_3$ 3:2.)

HPLC VARIABLES**Column:** 250 × 4.6 5 µm HP Hypersil ODS (A) or HP LiChrosorb RP-18 (B)

Mobile phase: Gradient. MeCN:buffer from 15:85 at 2 min to 80:20 at 20 min (Buffer was 50 mM NaH_2PO_4 containing 16 mM propylamine hydrochloride, adjusted to pH 3 with concentrated phosphoric acid.)

Flow rate: 1**Injection volume:** 5**Detector:** UV 230, UV 275

CHROMATOGRAM**Retention time:** 15.25 (A), 15.8 (B)**Internal standard:** β-hydroxyethyltheophylline (3.7 (A), 4.4 (B))**Limit of detection:** 1000 ng/mL

OTHER SUBSTANCES

Extracted: furosemide, metolazone, amiloride, acetazolamide, chlorothiazide, hydrochlorothiazide, quinethazone, triamterene, hydroflumethiazide, chlorthalidone, dichlorphenamide, trichloromethiazide, methyclothiazide, benzthiazide, cyclothiazide, polythiazide, bendroflumethiazide, probenecid, spironolactone, canrenone, flumethiazide

Noninterfering: acetaminophen, aspirin, caffeine, diflunisal, fenoprofen, ibuprofen, indomethacin, methocarbamol, naproxen, phenylbutazone, sulindac, tetracycline, theobromine, theophylline, tolmetin, trimethoprim, verapamil

Interfering: ethacrynic acid

REFERENCE

Cooper, S.F.; Massé, R.; Dugal, R. Comprehensive screening procedure for diuretics in urine by high-performance liquid chromatography, *J. Chromatogr.*, **1989**, 489, 65–88.

SAMPLE**Matrix:** urine

Sample preparation: 2 mL Urine + 1 mL 10 mM HCl + 2000 ng bendroflumethiazide, extract with 5 mL ethyl acetate, centrifuge at 3000 rpm for 5 min. Remove the organic layer and dry it under a stream of nitrogen at 40°. Reconstitute with 100 µL MeOH, inject a 2 µL aliquot.

HPLC VARIABLES**Column:** 100 × 2.1 5 µm Hypersil ODS

Mobile phase: Gradient. MeOH: 50 mM ammonium acetate from 10:90 to 60:40 over 10 min, maintain at 60:40 for 10 min.

Column temperature: 40**Flow rate:** 0.3**Injection volume:** 2**Detector:** F ex 231 em 426 or UV 230

CHROMATOGRAM**Retention time:** 6.2

Internal standard: Bendroflumethiazide (F ex 223 em 415) (8.6)

Limit of detection: <10 ng/mL

OTHER SUBSTANCES

Extracted: furosemide (UV), piretanide (UV), cyclopenthiiazide (UV), etozolin (UV), canren-one (UV)

REFERENCE

Gradeen,C.Y.; Billay,D.M.; Chan,S.C. Analysis of bumetanide in human urine by high-performance liquid chromatography with fluorescence detection and gas chromatography/mass spectrometry, *J.Anal.Toxicol.*, **1990**, *14*, 123–126.

SAMPLE

Matrix: urine

Sample preparation: Make 5 mL urine alkaline (pH 9-10), add 2 g NaCl, extract twice with 6 mL ethyl acetate. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 200 μ L MeCN/water, inject a 10-20 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 4.5 μ m SGE 100 GL-4 C18P (Scientific Glass Engineering)

Mobile phase: MeCN:MeOH:water:trifluoroacetic acid 15:15:70:0.5

Flow rate: 0.8 or 1

Injection volume: 10-20

Detector: MS, ZAB2-SEQ (VG), PSP source coupled to LC, source 250°, probe 240-260°, scan m/z 200-550 or UV 270

CHROMATOGRAM

Retention time: 3.7

Limit of detection: 50 ng (by MS)

OTHER SUBSTANCES

Extracted: probenecid, ethacrynic acid, spironolactone

REFERENCE

Ventura,R.; Fraisse,D.; Becchi,M.; Paisse,O.; Segura,J. Approach to the analysis of diuretics and masking agents by high-performance liquid chromatography-mass spectrometry in doping control, *J.Chromatogr.*, **1991**, *562*, 723–736.

SAMPLE

Matrix: urine

Sample preparation: Buffer urine to 4.9 by mixing with an equal volume of pH 4.9 200 mM sodium phosphate buffer. Inject a 40 μ L aliquot onto column A with mobile phase A, after 3 min backflush the contents of column A onto column B with mobile phase B and start the gradient. At the end of the run re-equilibrate for 10 min.

HPLC VARIABLES

Column: A 20 \times 4.5 μ m Hypersil octadecylsilica ODS; B 200 \times 4.6 5 μ m Shiseido SG-120 polymer-based C18

Mobile phase: A water; B Gradient. MeCN:buffer from 7:93 to 15:85 over 3.5 min, to 50:50 over 8.5 min, maintain at 50:50 for 11 min (Buffer was 6.9 g NaH₂PO₄·H₂O in 1 L water, pH adjusted to 3.1 with phosphoric acid.)

Flow rate: 1

Injection volume: 40

Detector: UV 230

CHROMATOGRAM

Retention time: 19.5

Limit of detection: 500 ng/mL

OTHER SUBSTANCES

Extracted: acetazolamide, amiloride, bendroflumethiazide, benzthiazide, caffeine, carbamazepine, chlorothiazide, chlorthalidone, clopamide, dichlorfenamide, ethacrynic acid, furosemide, hydrochlorothiazide, metyrapone, probenecid, spironolactone, triamterene, trichlormethiazide

KEY WORDS

column-switching; optimum detection wavelengths vary for each drug

REFERENCE

Saareninen, M.; Sirén, H.; Riekkola, M.-L. A column switching technique for the screening of diuretics in urine by high performance liquid chromatography, *J. Liq. Chromatogr.*, **1993**, *16*, 4063–4078.

SAMPLE

Matrix: urine

Sample preparation: 5 mL Urine + 50 μ L 100 μ g/mL 7-propyltheophylline in MeOH + 200 μ L ammonium chloride buffer + 2 g NaCl, extract with 6 mL ethyl acetate by rocking at 40 movements/min for 20 min and centrifuging at 800 g for 5 min, repeat extraction, combine organic layers, evaporate to dryness at 40° under a stream of nitrogen. Reconstitute in 200 μ L MeCN:water 15:85, inject a 20 μ L aliquot. (Ammonium chloride buffer was 28 g ammonium chloride in 100 mL water with the pH adjusted to 9.5 with concentrated ammonia solution.)

HPLC VARIABLES

Column: 75 \times 4.6 3 μ m Ultrasphere ODS

Mobile phase: Gradient. MeCN:100 mM ammonium acetate adjusted to pH 3 with concentrated phosphoric acid from 10:90 to 15:85 over 2 min to 55:45 over 3 min to 60:40 over 3 min, maintain at 60:40 for 1 min, decrease to 10:90 over 1 min, equilibrate at 10:90 for 2 min.

Flow rate: 1

Injection volume: 20

Detector: UV 270

CHROMATOGRAM

Retention time: 7.2

Internal standard: 7-propyltheophylline (4.5)

Limit of detection: 20 ng/mL

OTHER SUBSTANCES

Extracted: acetazolamide, amiloride, bendroflumethiazide, benzthiazide, buthiazide, caffeine, canrenone, chlorthalidone, clopamide, cyclothiazide, diclofenamide, ethacrynic acid, furosemide, hydrochlorothiazide, mesocarb, morazone, piretanide, polythiazide, probenecid, spironolactone, torsemide, triamterene

Interfering: xipamide

REFERENCE

Ventura, R.; Nadal, T.; Alcalde, P.; Pascual, J.A.; Segura, J. Fast screening method for diuretics, probenecid and other compounds of doping interest, *J. Chromatogr. A*, **1993**, *655*, 233–242.

SAMPLE

Matrix: urine

Sample preparation: Inject 50 μ L urine directly onto column A with mobile phase A and elute to waste, after 1 min backflush the contents of column A onto column B with mobile phase B, elute column B with mobile phase B, monitor the effluent from column B.

HPLC VARIABLES

Column: A 20×2.1 30 μm Hypersil ODS-C18; B 250×4 5 μm Hypersil ODS-C18

Mobile phase: A Water; B Gradient. MeCN:buffer 15:85 for 1.5 min then to 80:20 over 8 min. Keep at 80:20 for 2.5 min then re-equilibrate with 15:85. (Buffer was 50 mM NaH_2PO_4 + 1.4 mL propylamine hydrochloride per liter adjusted to pH 3 with concentrated phosphoric acid.)

Flow rate: 1

Injection volume: 50

Detector: UV 230

CHROMATOGRAM

Retention time: 10

Limit of detection: 4 ng/mL

OTHER SUBSTANCES

Extracted: acetazolamide, amiloride, bendroflumethiazide, chlorthalidone, cyclothiazide, furosemide, hydrochlorothiazide, probenecid, spironolactone, triamterene

Interfering: ethacrynic acid

KEY WORDS

column-switching

REFERENCE

Campíns-Falco,P.; Herráez-Hernández,R.; Sevillano-Cabeza,A. Column-switching techniques for screening of diuretics and probenecid in urine samples, *Anal.Chem.*, **1994**, *66*, 244–248.


Bunazosin

Molecular formula: C₁₉H₂₇N₅O₃

Molecular weight: 373.46

CAS Registry No.: 80755-51-7

Merck Index: 1512



SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 150 × 4.6 5 μm Nucleosil C18

Mobile phase: MeCN:17 mM acetate containing 5 mM sodium laurylsulfate 50:50

Flow rate: 1

Detector: UV 245


KEY WORDS

water; buffer

REFERENCE

Kato,A.; Iwata,S. Studies on improved corneal permeability to bunazosin, *J.Pharmacobiodyn.*, **1988**, *11*, 330-334.

Bupivacaine



Molecular formula: $C_{18}H_{28}N_2O$

Molecular weight: 288.43

CAS Registry No.: 2180-92-9, 14252-80-3 (HCl hydrate), 18010-40-7 (HCl)

Merck Index: 1520

Lednicher No.: 1 17

SAMPLE

Matrix: blood

Sample preparation: Mix 50 μ L water with 500 μ L pH 6.0 buffer and 100 μ L plasma sample, add 10 mL chloroform, shake mechanically for 10 min, centrifuge at 873 g for 10 min, evaporate the organic layer under a stream of nitrogen, reconstitute the residue in 100 μ L mobile phase, inject a 35 μ L aliquot. (Prepare pH 6.0 buffer by mixing 67 mM KH_2PO_4 with 67 mM Na_2HPO_4 in an 87.7:12.3 ratio.)

HPLC VARIABLES

Column: 250 \times 2.1 5 μ m Supelcosil ABZ+plus deactivated reversed-phase

Mobile phase: MeOH:MeCN:50 mM pH 4.5 monobasic ammonium phosphate 5:7:63

Flow rate: 0.4

Injection volume: 35

Detector: UV 235

CHROMATOGRAM

Retention time: 26.90

Internal standard: bupivacaine (26.90)

OTHER SUBSTANCES

Extracted: cocaine, benzoylecgonine, norcocaine, cocaethylene,

Simultaneous: ascorbic acid, morphine, oxymorphone, noroxymorphone, norhydromorphone, norcodeine, codeine, nalorphine, procaine, acetaminophen, oxycodone, hydrocodone, caffeine, ethylmorphine, lidocaine, benzoynorecgonine, ketamine, acepromazine, salicylic acid, benzoic acid, thebaine, cocaine propyl ester, benzocaine, tetracaine, pentobarbital

KEY WORDS

rat; plasma; pharmacokinetics; bupivacaine is IS

REFERENCE

Pan, W.; Hedaya, M.A. Sensitive and specific high-performance liquid chromatographic assay with ultra-violet detection for the determination of cocaine and its metabolites in rat plasma, *J.Chromatogr.B*, **1997**, 703, 129–138.

SAMPLE

Matrix: blood

Sample preparation: Mix 1 mL plasma with 200 μ L 2 μ g/mL IS in MeOH, add 2 mL water and 2 mL MeCN, vortex gently, set aside for 3 min, centrifuge at 2200 g for 20 min. Separate the clear supernatant, add 500 μ L 200 mM NaOH and extract with 6 mL n-hexane by vortexing for 2 min. Centrifuge at 2200 g for 15 min. Evaporate 5 mL of the organic phase to dryness under reduced pressure. Reconstitute the residue in 120 μ L mobile phase. Inject a 100 μ L aliquot.

HPLC VARIABLES

Guard column: 10 \times 3 5 μ m AGP bonded silica (ChromTech, Hagersten, Sweden)

Column: 150 \times 4 5 micro.m AGP bonded silica (ChromTech, Hagersten, Sweden)

Mobile phase: Isopropanol:buffer 4:96 (Prepare mobile phase by adding 4% isopropanol and 0.6% diethylamine to 8 mM sodium dihydrogen phosphate containing 100 mM NaCl, adjust to pH 7.05 with 50% phosphoric acid.)

Flow rate: 0.9

Injection volume: 100

Detector: UV 214

CHROMATOGRAM

Retention time: 29.35 (R-(+)), 38.25 (S-(-))

Internal standard: diazepam(19.21)

Limit of detection: 4 ng/mL

Limit of quantitation: 5 ng/mL

OTHER SUBSTANCES

Extracted: lidocaine

KEY WORDS

plasma; pharmacokinetics; chiral

REFERENCE

Abraham,I.; Fawcett,J.P.; Kennedy,J.; Kumar,A.; Ledger,R. Simultaneous analysis of lignocaine and bupivacaine enantiomers in plasma by high-performance liquid chromatography, *J.Chromatogr.B*, **1997**, 703, 203–208.

SAMPLE

Matrix: blood

Sample preparation: Directly inject 100 μ L serum onto column A and elute to waste with mobile phase A, after 2 min elute to waste with mobile phase B, after 2 min backflush the contents of column A onto column B with mobile phase C, after 2 min remove column A from the circuit, elute column B with mobile phase C, monitor the effluent from column B. Re-equilibrate column A with mobile phase A for 4 min.

HPLC VARIABLES

Column: A 20 \times 4.6 protein-coated Lichrosorb RP-8 (preparation details not given); B 250 \times 4.6 Lichrosorb RP-18

Mobile phase: A pH 7.4 phosphate buffered saline; B MeOH:100 mM pH 5.5 phosphate buffer 20:80; C MeCN:MeOH:pH 6.4 phosphate buffer:ethylamine 30:30:40:0.3

Flow rate: 1

Injection volume: 100

Detector: UV 254

CHROMATOGRAM

Retention time: <10

Limit of detection: 70 ng/mL

OTHER SUBSTANCES

Noninterfering: metabolites, codeine, epinephrine, lignocaine, meperidine, morphine

KEY WORDS

serum; column-switching

REFERENCE

Emara,S.; Khedr,A.; Askal,H. Rapid and specific precolumn extraction high-performance liquid chromatographic assay for bupivacaine in human serum, *Biomed.Chromatogr.*, **1996**, 10, 131–134.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 50 μ L 40 μ g/mL etidocaine hydrochloride in water + 100 μ L 1 M NaOH, vortex for 15 s, add 5 mL diethyl ether, shake on a reciprocating shaker for 10 min, centrifuge. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 100 μ L mobile phase, inject an 80 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 3.2 10 μ m μ Bondapak C18

Mobile phase: MeCN:50 mM pH 5.80 Na₂HPO₄ 25:75

Flow rate: 0.9

Injection volume: 80

Detector: UV 210

CHROMATOGRAM

Retention time: 8.7

Internal standard: etidocaine (12.0)

Limit of detection: 50 ng/mL

OTHER SUBSTANCES

Extracted: 2,6-pipecolylxylidine, mepivacaine, lidocaine

Noninterfering: metabolites, 2,3-chloroprocaine, theophylline, mexiletine, quinidine, disopyramide, verapamil, phenobarbital, phenytoin, carbamazepine, ethosuximide, digoxin, theobromine, caffeine, furosemide, phenprocoumon, aldactone

KEY WORDS

plasma

REFERENCE

Ha,H.-R.; Funk,B.; Gerber,H.R.; Follath,F. Determination of bupivacaine in plasma by high-performance liquid chromatography, *Anesth.Analg.*, **1984**, 63, 448–450.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Serum + 1 μ g diazepam +100 μ L 2 M NaOH + 7 mL n-hexane, extract. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 100 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 4 Enantiopak α 1-acid glycoprotein (LKB)

Mobile phase: Isopropanol:8 mM sodium phosphate buffer containing 100 mM NaCl 9:91

Column temperature: 30

Flow rate: 0.3

Injection volume: 20

Detector: UV 215

CHROMATOGRAM

Retention time: 21.0 (R), 28.3 (S)

Internal standard: diazepam (15.7)

Limit of quantitation: 500 ng/mL

KEY WORDS

serum; chiral

REFERENCE

Lee,E.J.; Ang,S.B.; Lee,T.L. Stereoselective high-performance liquid chromatographic assay for bupivacaine enantiomers, *J.Chromatogr.*, **1987**, 420, 203–206.

SAMPLE**Matrix:** blood**Sample preparation:** 1 mL Plasma + 50 μ L 2 μ g/mL etidocaine in water + 100 μ L 1 M NaOH + 3 mL heptane:ethyl acetate 90:10, shake for 2 min, centrifuge at 1200 g for 10 min. Remove the organic phase and add it to 50 μ L 50 mM sulfuric acid, shake for 2 min, centrifuge at 1200 g for 5 min. Remove the aqueous phase and add it to 820 μ g sodium acetate, inject a 40 μ L aliquot. (The sodium acetate was measured out by adding 50 μ L 200 mM sodium acetate in MeOH to the tube and evaporating the MeOH.)

HPLC VARIABLES**Column:** 250 \times 4 10 μ m μ Bondapak C18**Mobile phase:** MeCN:10 mM NaH_2PO_4 20:80, adjusted to pH 2.1**Column temperature:** 30**Flow rate:** 1**Injection volume:** 40**Detector:** UV 205

CHROMATOGRAM**Retention time:** 12**Internal standard:** etidocaine (10)**Limit of detection:** 2 ng/mL

KEY WORDS

plasma; rabbit

REFERENCE

Le Guévello,P.; Le Corre,P.; Chevanne,P.; Le Verge,R. High-performance liquid chromatographic determination of bupivacaine in plasma samples for biopharmaceutical studies and application to seven other local anaesthetics, *J.Chromatogr.*, **1993**, 622, 284–290.

SAMPLE**Matrix:** blood**Sample preparation:** 1 mL Plasma + 10 μ L 100 μ g/mL lidocaine hydrochloride in water + 100 μ L 2 M NaOH, mix, add 3 mL n-hexane, shake for 1 min, centrifuge at 3500 rpm for 10 min. Remove the organic layer and evaporate it to dryness with nitrogen under vacuum, reconstitute the residue in 200 μ L mobile phase, vortex, inject a 100 μ L aliquot.

HPLC VARIABLES**Column:** 250 \times 4 5 μ m Spherisorb ODS-2**Mobile phase:** MeOH:50 mM pH 5.9 KH_2PO_4 38:62**Flow rate:** 1**Injection volume:** 100**Detector:** UV 254

CHROMATOGRAM**Retention time:** 5.9**Internal standard:** lidocaine (3.8)**Limit of detection:** 25 ng/mL

KEY WORDS

plasma

REFERENCE

Murillo,I.; Costa,J.; Salvá,P. Determination of bupivacaine in human plasma by HPLC, *J.Liq.Chromatogr.*, **1993**, 16, 3509–3514.

SAMPLE**Matrix:** blood

Sample preparation: Condition a Bakerbond cyano SPE cartridge two 1 mL aliquots of eluent. Add 2 mL MeCN to 1 mL plasma slowly while whirlmixing, let stand at room temperature for 5 min, centrifuge at 3000 g for 10 min. Remove the supernatant and dilute it with 15 mL water, add this mixture to the SPE cartridge, wash with three 1 mL aliquots of water, elute with eluent. Evaporate the eluent to about 500 μ L under a stream of nitrogen, add 1 mL water, add 50 μ L 1 M NaOH, whirlmix for 30 s, add 6 mL n-hexane, rotate for 10 min, centrifuge at 3000 g for 5 min. Remove 5 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 80 μ L mobile phase, let stand for 2 h, inject a 50 μ L aliquot. (The eluent was MeOH:50 mM NaH_2PO_4 adjusted to pH 3.0 with 1 M phosphoric acid 50:50.)

HPLC VARIABLES

Guard column: 10 \times 4 Bakerbond chiral α 1-acid glycoprotein

Column: 100 \times 4 Bakerbond chiral α 1-acid glycoprotein

Mobile phase: Isopropanol:10 mM pH 6.8 NaH_2PO_4 6:94

Flow rate: 1

Injection volume: 50

Detector: UV 210

CHROMATOGRAM

Retention time: 20 (R-(+)), 25 (S-(-))

Limit of detection: 10 ng/mL (S), 8 ng/mL (R)

KEY WORDS

plasma; SPE; chiral; pharmacokinetics

REFERENCE

Groen,K.; Zeijlmans,P.W.M.; Burm,A.G.L.; van Kleef,J.W. Improved clean-up procedure for the high-performance liquid chromatographic assay of bupivacaine enantiomers in human plasma and ultrafiltrate in the nanogram per milliliter range, *J.Chromatogr.B*, **1994**, 655, 163–166.

SAMPLE

Matrix: blood

Sample preparation: Condition a 1 mL BondElut C18 SPE cartridge once with 1 M HCl, twice with MeOH, and once with water, remove the liquid completely with suction each time. Add 250 μ L IS solution and 250 μ L serum to the column at 1 mL/min, wash twice with water and once with MeCN draining the column completely after each wash, elute with 250 μ L eluting solution, centrifuge for 20 s to remove last of eluate, inject a 5 μ L aliquot of the eluate. (Prepare IS solution by adding 40 μ L 1 mg/mL N-pentyl-2,6-pipecoloxylidide (1-pentyl-N-(2,6-dimethylphenyl)-2-piperidinecarboxamide, pentyl-PPX) in MeOH to 10 mL 100 mM NaH_2PO_4 . Eluting solution was 2.5 mL 35% perchloric acid in 100 mL MeOH.)

HPLC VARIABLES

Guard column: 15 \times 3.2 7 μ m RP-8 (Applied Biosystems)

Column: 150 \times 4.6 5 μ m Ultrasphere octyl

Mobile phase: MeCN:10 mM KH_2PO_4 25:80, pH 5.2

Flow rate: 1.5

Injection volume: 5

Detector: UV 205

CHROMATOGRAM

Retention time: 7.2

Internal standard: N-pentyl-2,6-pipecoloxylidide (1-pentyl-N-(2,6-dimethylphenyl)-2-piperidinecarboxamide, pentyl-PPX) (14.5)

Limit of detection: 25 ng/mL

OTHER SUBSTANCES

Extracted: metabolites, mepivacaine, meperidine, fentanyl

Noninterfering: acetaminophen, codeine, epinephrine, morphine, diazepam

KEY WORDS

serum; SPE

REFERENCE

Gupta, R.N.; Dauphin, A. Column liquid chromatographic determination of bupivacaine in human serum using solid-phase extraction, *J.Chromatogr.B*, **1994**, 658, 113–119.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Plasma + 100 μ L 2 M NaOH, vortex briefly, add 5 mL anhydrous ethyl ether, vortex for 30 s, rotate for 10 min, centrifuge at 1000 g for 5 min. Remove 4.5 mL ether and add to 250 μ L 12.5 mM sulfuric acid, vortex for 30 s, rotate for 10 min, centrifuge for 5 min, inject a 50 μ L aliquot of the lower aqueous phase.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Octyl 1B (Keystone)

Mobile phase: MeCN:50 mM Na₂HPO₄ 27:73 pH adjusted to 5.8 with 50% phosphoric acid

Flow rate: 1

Injection volume: 50

Detector: UV 210

CHROMATOGRAM

Retention time: 9.18

Internal standard: bupivacaine

OTHER SUBSTANCES

Extracted: lidocaine, prilocaine, o-toluidine

KEY WORDS

plasma; pig; bupivacaine is IS

REFERENCE

Klein, J.; Fernandes, D.; Gazarian, M.; Kent, G.; Koren, G. Simultaneous determination of lidocaine, prilocaine and the prilocaine metabolite o-toluidine in plasma by high-performance liquid chromatography, *J.Chromatogr.B*, **1994**, 655, 83–88.

SAMPLE

Matrix: blood

Sample preparation: 450 μ L Plasma + 50 μ L pentycaine in 10 mM pH 7 ACES buffer, vortex, filter (0.22 μ m nylon syringe filter), inject a 10 μ L aliquot onto column A with mobile phase A, after 3 min elute contents of column A onto column B with mobile phase B, monitor the effluent from column B. After 3 min remove column A from the circuit and re-equilibrate with mobile phase A. (ACES is N-(2-acetamido)-2-aminoethanesulfonic acid.)

HPLC VARIABLES

Column: A 10 \times 4.6 5 μ m Regis semi-permeable surface (SPS) guard cartridge; B 100 \times 3.2 7 μ m Hypercarb pH (Shandon)

Mobile phase: A Isopropanol:buffer 3:97 (Buffer was 10 mM ACES adjusted to pH 7.0 with 5 M ammonia solution. ACES is N-(2-acetamido)-2-aminoethanesulfonic acid.); B 10 mM Acetic acid and 4 mM triethylamine in MeOH

Flow rate: A 0.8; B 0.5

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 5

Internal standard: pentycaine (6)

Limit of detection: 200 ng/mL

OTHER SUBSTANCES

Extracted: ropivacaine

KEY WORDS

plasma; column switching

REFERENCE

Yu,Z.; Abdel-Rehim,M.; Westerlund,D. Determination of amide-type local anaesthetics by direct injection of plasma in a column-switching high-performance liquid chromatographic system using a pre-column with a semipermeable surface, *J.Chromatogr.B*, **1994**, 654, 221–230.

SAMPLE

Matrix: blood

Sample preparation: Adjust pH of plasma to 7.4 with carbon dioxide gas. Filter (Amicon MPS-1 with YMT 30 membrane) while centrifuging at 37° at 500 g for 15 min, inject a 150 µL aliquot of the ultrafiltrate on to column A and elute to waste with mobile phase A, collect the eluate containing the compound in a 1 mL loop and inject it onto column B with mobile phase B, elute column B with mobile phase B and monitor the effluent.

HPLC VARIABLES

Column: A 125 × 4 4 µm Superspher RP-select B (Merck); B 150 × 4.6 5 µm Nucleosil 5 SA

Mobile phase: A MeCN:buffer 30:70 (Buffer was 3.1 mL 1 M phosphoric acid and 20 mL 1 M NaH₂PO₄ made up to 1 L with water, pH 3.); B MeCN:buffer 40:60 (Buffer was 149 mL 2 M ammonium hydroxide and 348 mL 1 M phosphoric acid made up to 1 L with water, pH 2.6.)

Column temperature: 27

Flow rate: 1

Injection volume: 150

Detector: UV 210

CHROMATOGRAM

Retention time: 12

Limit of detection: 10 nM

OTHER SUBSTANCES

Extracted: ropivacaine (elutes in a different fraction from column A)

KEY WORDS

plasma; ultrafiltrate; column-switching; rugged

REFERENCE

Arvidsson,T.; Eklund,E. Determination of free concentration of ropivacaine and bupivacaine in blood plasma by ultrafiltration and coupled-column liquid chromatography, *J.Chromatogr.B*, **1995**, 668, 91–98.

SAMPLE

Matrix: blood

Sample preparation: Condition a 3 mL Clean Screen SPE cartridge (Worldwide Monitoring) with two 2 mL portions of MeOH, with 3 mL water, and with 3 mL 10 mM pH 3.0

phosphate buffer. 1 mL Serum + 500 μ L 10 mM pH 3.0 phosphate buffer, mix, add to the SPE cartridge, air dry for 30 s, wash with 3 mL phosphate buffer, wash with 3 mL 100 mM HCl, wash with 3 mL MeOH, elute with 2 mL chloroform:isopropanol:ammonium hydroxide 22:20.5:2.5. Evaporate the eluate to dryness under a stream of nitrogen at 30°, reconstitute the residue in 100 μ L mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 150 \times 4.5 SemiPermeable Surface (SPS) C8 (Regis)

Mobile phase: THF:2.5 mM potassium phosphate buffer 3.25:96.75 containing 0.0025% triethylamine, final pH adjusted to 2.7-2.8 with 85% orthophosphoric acid

Flow rate: 0.5

Detector: UV 235

CHROMATOGRAM

Retention time: 24

Internal standard: bupivacaine

Limit of detection: 5 ng/mL

OTHER SUBSTANCES

Extracted: acepromazine, atropine, benzoylecgonine, benzoylnorecgonine, cocaine, ketamine, norcocaine

Noninterfering: benzethonium chloride, benzyl alcohol

KEY WORDS

SPE; serum; bupivacaine is IS

REFERENCE

Muztar,J.; Chari,G.; Bhat,R.; Ramaro,S.; Vidyasagar,D. A high-performance liquid chromatographic procedure for the separation of cocaine and some of its metabolites from acepromazine, ketamine, and atropine from serum, *J.Liq.Chromatogr.*, **1995**, *18*, 2635-2645.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 100 μ L 10 μ g/mL lidocaine in 25 mM sulfuric acid + 1 mL 1 M NaOH + 5 mL diethyl ether, shake or rotate for 15 min, centrifuge at 1000 rpm for 5 min, freeze at -20°. Remove the organic layer and add it to 250 μ L 25 mM sulfuric acid, shake for 15 min, centrifuge at 1000 rpm for 5 min, freeze, discard the organic layer. Thaw the aqueous layer, pass air over the aqueous phase at room temperature to remove traces of ether, adjust pH to 5.0-6.5 by adding 10 μ L 1 M NaOH, inject a 50 μ L aliquot.

HPLC VARIABLES

Guard column: 4 \times 4 4 μ m LiChroCART Superspher 60 RP Select B (Merck)

Column: 125 \times 4 4 μ m LiChroCART Superspher 60 RP Select B (Merck)

Mobile phase: MeCN:buffer 30:70 (Buffer was 7.0 g/L K_2HPO_4 in water adjusted to pH 5.8 with 1 M NaOH.)

Column temperature: 40

Flow rate: 1

Injection volume: 50

Detector: UV 202

CHROMATOGRAM

Retention time: 10

Internal standard: lidocaine (5)

Limit of quantitation: 100 ng/mL

KEY WORDS

plasma

REFERENCE

Sattler,A.; Krämer,I.; Jage,J.; Vrana,S.; Kleemann,P.P.; Dick,W. Development of a HPLC-system for quantitative measurement of lidocaine and bupivacaine in patients plasma during postoperative epidural pain therapy, *Pharmazie*, **1995**, *50*, 741–744.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform: isopropanol:n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 264

CHROMATOGRAM

Retention time: 5.43

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds(all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaline; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vandesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phencyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loperazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorphenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrodine; phenylbutazone; demexiptiline; clozapine; proganil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; op-

ipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozone; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254–262.

SAMPLE

Matrix: blood

Sample preparation: Condition a Bond Elut phenyl SPE cartridge with 4 mL MeOH and 4 mL water. 1 mL Serum or 300–500 μ L ultrafiltrate + lidocaine, add to the SPE cartridge, wash with 5 mL water, wash with 2 mL MeOH:water 5:95, wash with 2 mL EtOH:water 2.5:97.5, wash with 2 mL MeCN:water 10:90, elute with 1 mL MeCN:50 mM pH 2.4 phosphate buffer 25:75, inject a 100 μ L aliquot of the eluate.

HPLC VARIABLES

Column: 300 \times 4.6 μ Bondapak

Mobile phase: MeCN:50 mM pH 4.0 KH_2PO_4 25:75

Flow rate: 1.5

Injection volume: 100

Detector: UV 210

CHROMATOGRAM

Internal standard: lidocaine

Limit of detection: 10–15 ng/mL

KEY WORDS

serum; ultrafiltrate; SPE

REFERENCE

Mazoit, J.X.; Cao, L.S.; Samii, K. Binding of bupivacaine to human serum proteins, isolated albumin and isolated α -1-acid glycoprotein. Differences between the two enantiomers are partly due to cooperativity, *J. Pharmacol. Exp. Ther.*, **1996**, *276*, 109–115.

SAMPLE

Matrix: blood

Sample preparation: Filter (0.22 μ m) plasma, add ropivacaine (590 ng per 1 mL plasma), inject a 500 μ L aliquot of the filtrate on to column A and elute to waste with mobile phase A, after 7 min backflush the contents of column A on to column B with mobile phase B, monitor the effluent from column B.

HPLC VARIABLES

Column: A 25 \times 4 C18-alkyl-diol silica (*J. Chromatogr. B* 1995, 666, 315; *J. Pharm. Biomed. Anal.* 1995, 13, 615); B 10 \times 4.6 5 μ m Kromasil C8 (Eka Nobel) + 100 \times 4.6 5 μ m Kromasil C18 (Eka Nobel)

Mobile phase: A Isopropanol:buffer 5:95; B MeOH:buffer 63:37 (Prepare buffer by mixing 36.68 mL 1 M phosphoric acid with water, adjust pH to 7.4 with 1 M NaOH, make up to 2 L with water.)

Column temperature: 30 (column B)

Flow rate: A 1.5; B 1

Injection volume: 500

Detector: UV 210

CHROMATOGRAM

Retention time: 15

Internal standard: ropivacaine (9)

Limit of detection: 10 ng/mL

KEY WORDS

plasma; column-switching

REFERENCE

Yu,Z.; Westerlund,D. Direct injection of large volumes of plasma in a column-switching system for the analysis of local anaesthetics. II. Determination of bupivacaine in human plasma with an alkyl-diol silica precolumn, *J.Chromatogr.A*, **1996**, 725, 149–155.

SAMPLE

Matrix: blood

Sample preparation: Filter (0.22 μ m) plasma, inject a 500 μ L aliquot of the filtrate on to column A and elute to waste with mobile phase A, after 10 min backflush the contents of column A on to column B with mobile phase B, monitor the effluent from column B.

HPLC VARIABLES

Column: A 10 \times 10 5 μ m semi-permeable surface C8 (Regis); B 10 \times 4.6 5 μ m Kromasil C8 (Eka Nobel) + 100 \times 4.6 5 μ m Kromasil C18 (Eka Nobel)

Mobile phase: A Isopropanol:buffer 5:95; B MeOH:buffer 63:37 (Prepare buffer by mixing 36.68 mL 1 M phosphoric acid with 68.34 mL 1 M NaOH and diluting to 2 L with water, pH 7.7.)

Flow rate: A 1.5; B 1

Injection volume: 500

Detector: UV 240

CHROMATOGRAM

Retention time: 16

OTHER SUBSTANCES

Extracted: ropivacaine

KEY WORDS

plasma; column-switching

REFERENCE

Yu,Z.; Westerlund,D. Direct injection of large volumes of plasma in a column-switching system for the analysis of local anaesthetics. I. Optimization of semi-permeable surface precolumns in the system and characterization of some interference peaks, *J.Chromatogr.A*, **1996**, 725, 137–147.

SAMPLE

Matrix: blood, urine

Sample preparation: 500 μ L Serum or urine + 200 μ L 2 μ g/mL etidocaine in water + 1 mL 100 mM pH 9.0 sodium tetraborate + 5 mL diethyl ether, shake for 20 min, centrifuge at 1200 g. Remove the organic layer and add it to 500 μ L 200 mM HCl, shake for 20 min, centrifuge at 1200 g. Remove the aqueous layer and add it to 1 mL 100 mM pH 9.0 sodium tetraborate, add 5 mL diethyl ether, shake for 20 min, centrifuge at 1200 g. Remove the organic layer and evaporate it to dryness under a stream of air at 50°, reconstitute the residue in 100 μ L mobile phase, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 10 μ m Nucleosil C8

Mobile phase: THF:10 mM pH 2.4 potassium phosphate 8:92

Flow rate: 1.6

Injection volume: 50

Detector: UV 210

CHROMATOGRAM

Retention time: 4.15

Internal standard: etidocaine (10.60)

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

serum; pharmacokinetics

REFERENCE

Lindberg,R.L.P.; Kanto,J.H.; Pihlajamäki,K.K. Simultaneous determination of bupivacaine and its two metabolites, desbutyl- and 4'-hydroxybupivacaine, in human serum and urine, *J.Chromatogr.*, **1986**, 383, 357-364.

SAMPLE

Matrix: blood, urine

Sample preparation: 1 mL Plasma or urine + 100 μ L 5 μ g/mL IS in water + 100 μ L 1 M pH 10 sodium carbonate buffer, shake by hand for 5 s, add 6 mL n-hexane:isopropanol 5:1, rotate on a tumble-mixer at 28 rpm for 10 min, centrifuge at 2000 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 120 μ L mobile phase, inject a 100 μ L aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: MeOH:MeCN:20 mM pH 6 sodium phosphate buffer 15:40:45

Flow rate: 1.2

Injection volume: 100

Detector: UV 210

CHROMATOGRAM

Retention time: 9.2

Internal standard: 1-pentyl-2-(2',6'-xylylcarbamoyl)piperidine hydrochloride (13.2)

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

Noninterfering: fentanyl, morphine, oxycodone

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Kastrissios,H.; Hung,M.-F.; Triggs,E.J. High-performance liquid chromatographic method for the quantitation of bupivacaine, 2,6-pipecoloxylidide and 4'-hydroxybupivacaine in plasma and urine, *J.Chromatogr.*, **1992**, 577, 103-107.

SAMPLE

Matrix: blood, urine

Sample preparation: 2 mL Whole blood, plasma, or urine + 1 mL saturated sodium carbonate + 10 μ L 100 μ g/mL etidocaine, add to a 3 mL Extrelut SPE cartridge, elute with

15 mL dichloromethane. Evaporate eluate to dryness under a stream of nitrogen at 40°, reconstitute in 100 µL 10 mM HCl, add 3 mL diethyl ether, vortex for 20 s, centrifuge at 2800 g for 5 min, inject a 40 µL aliquot of the aqueous layer.

HPLC VARIABLES

Guard column: 5 × 6 µBondapak Guard Pak

Column: 300 × 3.9 10 µm µBondapak C18

Mobile phase: MeCN:100 mM ammonium acetate 50:50

Flow rate: 1.5

Injection volume: 40

Detector: UV 230

CHROMATOGRAM

Retention time: 11

Internal standard: etidocaine (14)

Limit of detection: 40 ng/mL

OTHER SUBSTANCES

Extracted: lidocaine, prilocaine, dibucaine

Also analyzed: procaine, butacaine, tetracaine, p-aminobenzoic acid, articaine, o-toluidine, caffeine, amphetamine, ephedrine, epinephrine, morphine, monoacetylmorphine, diamorphine, ethylmorphine, codeine, acetylcodeine

KEY WORDS

whole blood; plasma; SPE

REFERENCE

Rop,P.P.; Grimaldi,F.; Bresson,M.; Fornaris,M.; Viala,A. Liquid chromatographic analysis of cocaine, benzoylecgonine, local anaesthetic agents and some of their metabolites in biological fluids, *J.Liq.Chromatogr.*, **1993**, 16, 2797–2811.

SAMPLE

Matrix: formulations

Sample preparation: Dilute 5% bupivacaine hydrochloride injection and 25 mg/mL morphine sulfate injection with 0.9% NaCl injections to a bupivacaine concentration of 625 µg/mL and a morphine concentration of 100 µg/mL, inject a 10 µL aliquot.

HPLC VARIABLES

Column: 150 × 4.6 5 µm Spherisorb Phenyl

Mobile phase: MeCN:100 mM pH 5 phosphate buffer 40:60

Flow rate: 1.5

Injection volume: 10

Detector: UV 260

CHROMATOGRAM

Retention time: 5.2

OTHER SUBSTANCES

Simultaneous: degradation products

Noninterfering: morphine

KEY WORDS

injections; stability indicating

REFERENCE

Johnson,C.E.; Christen,C.; Perez,M.M.; Ma,M. Compatibility of bupivacaine hydrochloride and morphine sulfate, *Am.J.Health-Syst.Pharm.*, **1997**, 54, 61–64.

SAMPLE**Matrix:** formulations**Sample preparation:** Inject a 20 μ L aliquot

HPLC VARIABLES**Column:** 100 \times 3.9 4 μ m Radial Pak phenyl (Waters)**Mobile phase:** MeOH:buffer 65:35 (Buffer was 5 mM pH 4.8 phosphate buffer containing 1.4 mM tetrabutylammonium hydroxide.)**Flow rate:** 3**Injection volume:** 20**Detector:** UV 210

CHROMATOGRAM**Retention time:** 7.5

OTHER SUBSTANCES**Simultaneous:** degradation products, fentanyl

KEY WORDS

injections; saline; stability-indicating

REFERENCETu, Y.H.; Stiles, M.L.; Allen, L.V., Jr. Stability of fentanyl citrate and bupivacaine hydrochloride in portable pump reservoirs, *Am.J.Hosp.Pharm.*, **1990**, 47, 2037–2040.

SAMPLE**Matrix:** formulations

HPLC VARIABLES**Column:** 75 \times 4.6 3 μ m Ultrasphere XL-ODS C18**Mobile phase:** MeCN:100 mM pH 5 KH_2PO_4 40:60**Flow rate:** 1.5**Detector:** UV 254

KEY WORDS

injections; saline; stability-indicating

REFERENCEJones, J.W.; Davis, A.T. Stability of bupivacaine hydrochloride in polypropylene syringes, *Am.J.Hosp.Pharm.*, **1993**, 50, 2364–2365.

SAMPLE**Matrix:** formulations**Sample preparation:** Inject a 20 μ L aliquot.

HPLC VARIABLES**Column:** 300 \times 3.9 10 μ m μ Bondapak phenyl**Mobile phase:** MeCN:20 mM KH_2PO_4 adjusted to pH 6.0 with 1 M KOH 50:50**Flow rate:** 1**Injection volume:** 20**Detector:** UV 235

CHROMATOGRAM**Retention time:** 14.1**Limit of detection:** 254 ng/mL

OTHER SUBSTANCES

Simultaneous: morphine, hydromorphone

KEY WORDS

saline; injections

REFERENCE

Venkateshwaran,T.G.; Stewart,J.T. HPLC determination of morphine-hydromorphone-bupivacaine and morphine-hydromorphone-tetracaine mixtures in 0.9% sodium chloride injection, *J.Liq.Chromatogr.*, **1995**, *18*, 565–578.

SAMPLE

Matrix: perfusate

HPLC VARIABLES

Column: 100 × 8 4 μm Novapak C18

Mobile phase: MeCN:0.092%phosphoric acid + 0.2% triethylamine 26:74

Flow rate: 2

Detector: UV 214

CHROMATOGRAM

Internal standard: lidocaine

Limit of quantitation: 10 ng/mL

OTHER SUBSTANCES

Also analyzed: diphenhydramine, diltiazem, metabolites

KEY WORDS

rat; liver

REFERENCE

Hussain,M.D.; Tam,Y.K.; Gray,M.R.; Coutts,K.T. Kinetic interactions of lidocaine, diphenhydramine, and verapamil with diltiazem: A study using isolated perfused rat liver, *Drug Metab.Dispos.*, **1994**, *22*, 530–536.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 μg/mL solution in MeOH, inject a 20 μL aliquot.

HPLC VARIABLES

Column: 125 × 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 1.6

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlor-

promazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazepine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipanone, diprenorphine, dipyridamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserin, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypropazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phenotolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tolycaine, translycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripelethnamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J.Chromatogr.*, **1985**, 323, 191-225.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in MeOH:water 1:1 at a concentration of 50 µg/mL, inject a 10 µL aliquot.

HPLC VARIABLES

Column: 300 × 3.9 10 µm µBondapak C18

Mobile phase: MeOH:acetic acid:triethylamine:water 50:1.5:0.5:48

Flow rate: 1.5

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 10

OTHER SUBSTANCES

Simultaneous: butacaine, lidocaine, benzocaine, tetracaine

REFERENCE

Roos,R.W.; Lau-Cam,C.A. General reversed-phase high-performance liquid chromatographic method for the separation of drugs using triethylamine as a competing base, *J.Chromatogr.*, **1986**, 370, 403-418.

SAMPLE

Matrix: solutions

Sample preparation: Inject an aliquot of a 200 μ M solution in MeOH.

HPLC VARIABLES

Column: 100 \times 4.7 7 μ m Hypercarb (Shandon)

Mobile phase: MeOH:isopropanol 90:10 containing 15 mM N-benzoyloxycarbonylglycyl-L-proline and 9 mM NaOH

Column temperature: 17

Injection volume: 20

Detector: UV 270

CHROMATOGRAM

Retention time: k' 3.09 (first enantiomer)

KEY WORDS

chiral; α = 1.08

REFERENCE

Huynh,N.-H.; Karlsson,A.; Pettersson,C. Enantiomeric separation of basic drugs using N-benzoyloxycarbonylglycyl-L-proline as counter ion in methanol, *J.Chromatogr.A*, **1995**, 705, 275-287.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Supelcosil LC-DP (A) or 250 \times 4 5 μ m LiChrospher 100 RP-8 (B)

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 10.60 (A), 5.79 (B)

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordiazepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyridamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenopropfen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, fu-rosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, le-

vorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazindol, mefenamic acid, meperidine, mephenytoin, mepivacaine, mesoridazine, metaproterenol, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl-dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfinpyrazone, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocinide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

KEY WORDS

also details of plasma extraction

REFERENCE

Koves,E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J.Chromatogr.A*, **1995**, 692, 103–119.