


Cabergoline

Molecular formula: C₂₆H₃₇N₅O₂

Molecular weight: 451.61

CAS Registry No.: 81409-90-7, 85329-89-1 (diphosphate)

Merck Index: 1637



SAMPLE

Matrix: blood, urine

Sample preparation: 1 mL Plasma or urine + 0.5 (plasma) or 1 (urine) mL 500 mM pH 9 borate buffer + 2.5 mL dichloromethane:isooctane 40:60, shake on a rotary mixer for 5 min, centrifuge at 1200 g for 15 min, repeat extraction. Combine the organic layers and add them to 100 μ L 100 mM phosphoric acid, vortex for 1 min, centrifuge at 1200 g for 10 min. Remove the aqueous phase and add it to 1 mL n-hexane, vortex, inject a 50 μ L aliquot of the aqueous phase.

HPLC VARIABLES

Column: 250 \times 4.5 μ m Nucleosil C18

Mobile phase: MeCN:75 mM pH 3 phosphate buffer 20:80

Flow rate: 0.6

Injection volume: 50

Detector: E, ESA Model 5100 Coulochem, Model 5020 guard cell +670 mV, Model 5011 analytical cell, detector 1 +350 mV, detector 2 +650 mV

CHROMATOGRAM

Retention time: 13.5

Limit of quantitation: 0.30 ng/mL (urine), 0.25 ng/mL (plasma)

OTHER SUBSTANCES

Noninterfering: levodopa, methyldopa, dopamine, homovanillic acid

KEY WORDS

plasma

REFERENCE

Pianezzola,E.; Bellotti,V.; La Croix,R.; Strolin Benedetti,M. Determination of cabergoline in plasma and urine by high-performance liquid chromatography with electrochemical detection, *J.Chromatogr.*, **1992**, 574, 170-174.

SAMPLE

Matrix: urine

Sample preparation: Add urine to Amberlite XAD-2 resin (resin/urine = 0.35-0.40), wash with water, elute with acetone, elute with acetic acid. Combine aliquots of the eluates, evaporate to dryness under vacuum, dissolve the residue in mobile phase, inject an aliquot.

HPLC VARIABLES

Guard column: Guardpak C18 (Waters)

Column: 300 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: Gradient. MeCN:100 mM pH 2 KCl/HCl buffer 12:88 for 5 min, to 30:70 over 15 min, maintain at 30:70 for 15 min

Flow rate: 1

Injection volume: 200

Detector: UV 280

CHROMATOGRAM**Retention time:** 25.4

OTHER SUBSTANCES**Extracted:** metabolites

KEY WORDSrat; monkey; human; SPE

REFERENCE

Battaglia,R.; Strolin Benedetti,M.; Mantegani,S.; Castelli,M.G.; Cocchiara,G.; Dostert,P. Disposition and urinary metabolite pattern of cabergoline, a potent dopaminergic agonist, in rat, monkey and man, *Xenobiotica*, **1993**, 23, 1377–1389.


Cadralazine

Molecular formula: C₁₂H₂₁N₅O₃

Molecular weight: 283.33

CAS Registry No.: 64241-34-5

Merck Index: 1669



SAMPLE

Matrix: blood

Sample preparation: 200 μ L Plasma + 1.2 mL chloroform:EtOH 95:5, shake horizontally at 300 rpm for 12 min, centrifuge at 1220 g for 5 min. Remove the organic phase and add it to 100 μ L 90 mM KCl in 10 mM HCl (pH 2), shake horizontally at 400 rpm for 12 min, centrifuge at 1220 g for 5 min, inject a 20 μ L aliquot of the aqueous supernatant. (cf *J. Chromatogr.* 1984, 290, 223.)

HPLC VARIABLES

Column: 150 \times 1.5 μ m Nucleosil C18 in a glass-lined stainless steel column

Mobile phase: MeCN:100 mM NaH₂PO₄:1 M NaOH 18:79:3

Flow rate: 0.06

Injection volume: 20

Detector: UV 254 (2.4 μ L flow cell)

CHROMATOGRAM

Retention time: 6

Internal standard: CGP 24 751

Limit of quantitation: 5 ng/mL

KEY WORDS

plasma; microbore

REFERENCE

Rouan, M.C. Microbore liquid chromatographic determination of cadralazine and cephalexin in plasma with large-volume injection, *J. Chromatogr.*, **1988**, 426, 335–344.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Whole blood + 500 μ L MeCN, mix, centrifuge at 2200 g for 2 min. Remove the supernatant and saturate it with anhydrous potassium carbonate (about 2 g per 1 mL of whole blood), centrifuge at 2200 g for 3 min, remove the MeCN layer, inject a 100 μ L aliquot of the MeCN layer.

HPLC VARIABLES

Guard column: 10 μ m octyl (Upchurch)

Column: 150 \times 4.6 3 μ m C8 (Alltech)

Mobile phase: MeCN:5 mM hexanesulfonic acid in 1% acetic acid 30:70

Flow rate: 1

Injection volume: 100

Detector: UV 250

CHROMATOGRAM

Retention time: 11

Limit of quantitation: 10 ng/mL

OTHER SUBSTANCES

Noninterfering: erythromycin, cimetidine, ranitidine, aspirin, acetaminophen, creatinine, theophylline, caffeine

KEY WORDS

whole blood; salting-out

REFERENCE

Rustum, A.M. Determination of cadralazine in human whole blood using reversed-phase high-performance liquid chromatography: utilizing a salting-out extraction procedure, *J. Chromatogr.*, **1989**, *489*, 345–352.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 1 g Plasma + 500 μ L 1 μ g/mL IS in 5 mM sulfuric acid + 1 mL water + 6 mL chloroform:EtOH 95:5, shake horizontally at 300 rpm for 12 min, centrifuge at 1220 g for 5 min. Remove the organic phase and add it to 500 μ L pH 2 buffer, shake horizontally at 400 rpm for 12 min, centrifuge at 1220 g for 5 min, remove the aqueous layer and inject a 100 μ L aliquot. Urine. 1 g Urine + 500 μ L 60 μ g/mL IS in 5 mM sulfuric acid + 1 mL pH 12 phosphate buffer + 6 mL chloroform:EtOH 95:5, shake horizontally at 300 rpm for 12 min, centrifuge at 1220 g for 5 min. Remove the organic phase and add it to 500 μ L pH 2 buffer, shake horizontally at 400 rpm for 12 min, centrifuge at 1220 g for 5 min, remove the aqueous layer and inject a 100 μ L aliquot. (pH 2 Buffer was 90 mmol/L KCl in 1 L 10 mM HCl, pH 2.06. pH 12 Buffer was 3.8 g $\text{Na}_3\text{PO}_4 \cdot 12\text{H}_2\text{O}$ in 1 L water adjusted to pH 12 with about 20 mL 2 M NaOH.)

HPLC VARIABLES

Guard column: 50 \times 3.2 30–40 μ m Perisorb RP-8

Column: 250 \times 4.6 10 μ m LiChrosorb RP-8

Mobile phase: MeCN:100 mM pH 6 phosphate buffer 15:85 (Buffer was 100 mM NaH_2PO_4 adjusted to pH 6 with about 30 mL 1 M NaOH.)

Column temperature: 30

Flow rate: 2.7

Injection volume: 100

Detector: UV 254

CHROMATOGRAM

Retention time: 6

Internal standard: 2-[6-[ethyl(2-hydroxypropyl)amino]-3-pyridazinyl]hydrazine carboxylic acid propyl ester (CGP 24 751) (14)

Limit of quantitation: 10.59 nmol/g (urine), 0.141 nmol/g (plasma)


KEY WORDS

plasma; pharmacokinetics

REFERENCE

Haufler, S.A.; Dubois, J.P. Determination of cadralazine in human plasma and urine by high-performance liquid chromatography, *J. Chromatogr.*, **1984**, *290*, 223–230.

Caffeine



Molecular formula: $C_8H_{10}N_4O_2$

Molecular weight: 194.19

CAS Registry No.: 58-08-2, 5743-12-4 (monohydrate)

Merck Index: 1674

Lednicer No.: 1 111

SAMPLE

Matrix: beverages

Sample preparation: Filter sample.

HPLC VARIABLES

Column: 150×4.5 μm Hiasil C18 (Higgins)

Mobile phase: MeOH:25 mM phosphate buffer 45:55, pH 3.0

Flow rate: 1.0

Injection volume: 20

Detector: UV 218

CHROMATOGRAM

Retention time: 2.0

Limit of detection: 1.2 mg/mL

OTHER SUBSTANCES

Simultaneous: aspartame, benzoic acid

KEY WORDS

comparison with UV spectrophotometry and capillary electrophoresis; soft drinks

REFERENCE

McDevitt,V.L.; Rodriguez,A.; Williams,K.R. Analysis of soft drinks: UV spectrophotometry, liquid chromatography, and capillary electrophoresis, *J.Chem.Educ.*, **1998**, 75, 625–629.

SAMPLE

Matrix: beverages, blood, formulations

Sample preparation: Plasma. 100 μL Plasma + 10 μL 1,3-dimethyl-7-(2-hydroxyethyl)xanthine in water, vortex, add 1 mL acetone, vortex for 1 min, centrifuge at 2500 g for 5 min. Filter (0.45 μm) the supernatant, inject a 20 μL aliquot of the supernatant. Beverages. Dilute beverages 4–25 times with water, inject an aliquot. Tablets. Powder tablets, dissolve in water, filter, inject an aliquot of the filtrate.

HPLC VARIABLES

Column: 150×6 5 μm 3-aminopropylsilyl silica gel with the amino group derivatized with 1,8-naphthalic anhydride (Bunseki Kagaku 1993, 42, 817)

Mobile phase: MeOH:buffer 50:50 (Prepare buffer by dissolving 6.183 g boric acid and 1.461 g NaCl in 500 mL water, adjust pH to 6.4 with sodium borate solution.)

Column temperature: 30

Flow rate: 1

Injection volume: 20

Detector: UV 270

CHROMATOGRAM

Retention time: 24

Internal standard: 1,3-dimethyl-7-(2-hydroxyethyl)xanthine (12)

Limit of detection: 140 ng/mL

OTHER SUBSTANCES

Simultaneous: hypoxanthine, pentoxifylline, propentofylline, theobromine, theophylline, uric acid, xanthine

KEY WORDS

plasma; tablets; pharmacokinetics

REFERENCE

Nakashima,K.; Inoue,K.; Mayahara,K.; Kuroda,N.; Hamachi,Y.; Akiyama,S. Use of 3-(1,8-naphthalimido)propyl-modified silyl silica gel as a stationary phase for the high-performance liquid chromatographic separation of purine derivatives, *J.Chromatogr.A*, **1996**, 722, 107–113.

SAMPLE

Matrix: beverages, syrup

Sample preparation: Dilute syrup ten fold. Filter (0.45 μm) beverages and diluted syrup, inject a 10-20 μL aliquot of the filtrate.

HPLC VARIABLES

Column: 300 \times 3.9 10 μm μ Bondapak C18

Mobile phase: MeOH:acetic acid:water 20:5:75

Flow rate: 2

Injection volume: 10-20

Detector: UV 254

CHROMATOGRAM

Retention time: 4.5

Limit of detection: 20 ng

OTHER SUBSTANCES

Simultaneous: acesulfame, benzoic acid, dulcin, p-hydroxybenzoic acid, saccharin, vanillin

REFERENCE

Veerabhadrarao,M.; Narayan,M.S.; Kapur,O.; Sastry,C.S. Reverse phase liquid chromatographic determination of some food additives, *J.Assoc.Off.Anal.Chem.*, **1987**, 70, 578–582.

SAMPLE

Matrix: blood

Sample preparation: Mix 100 μL serum with 50 μL 1 M sodium hydroxide, extract with 1 mL 1.5 $\mu\text{g/mL}$ IS in isopropanol:chloroform 10:90. Centrifuge at 2000 g for 1 min, decant the organic layer, dry under a gentle air flow at 60°, dissolve the residue in 0.2 mL mobile phase, inject a 50 μL aliquot.

HPLC VARIABLES

Column: 150 \times 3.9 Nova Pak C18

Mobile phase: MeOH:THF:50 mM pH 4.7 ammonium acetate buffer 2:1.5:96.5

Flow rate: 1

Injection volume: 50

Detector: UV 273

CHROMATOGRAM

Internal standard: β -hydroxyethyltheophylline

Limit of quantitation: 200 ng/mL

OTHER SUBSTANCES

Extracted: theobromine, paraxanthine, theophylline

KEY WORDS

pharmakokinetics; serum

REFERENCE

Lee,T.C.; Charles,B.; Steer,P.; Flenady,V.; Shearman,A. Population pharmacokinetics of intravenous caffeine in neonates with apnea of prematurity, *Clin.Pharmacol.Ther.*, **1997**, 61, 628–640.

SAMPLE

Matrix: blood

Sample preparation: Mix 50 μ L plasma with 50 μ L 20 μ g/mL IS in water, vortex for 10 s, add 20 μ L 20% perchloric acid, vortex for 10 s, centrifuge at 2000 g for 5 min, inject a 50 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Spherisorb C18

Mobile phase: MeCN:THF:concentrated acetic acid:water 2:2:0.5:95.5

Column temperature: 35

Flow rate: 1

Injection volume: 50

Detector: UV 273

CHROMATOGRAM

Retention time: 7.5

Internal standard: 7-(β -hydroxypropyl)theophylline (9.2)

Limit of detection: 100 ng/mL

Limit of quantitation: 2 μ g/mL

OTHER SUBSTANCES

Extracted: theophylline

Simultaneous: β -hydroxyethyltheophylline, 8-chlorotheophylline, theobromine

KEY WORDS

plasma

REFERENCE

Schreiber-Deturmeny,E.; Bruguerolle,B. Simultaneous high-performance liquid chromatographic determination of caffeine and theophylline for routine drug monitoring in human plasma, *J.Chromatogr.B*, **1996**, 677, 305–312.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Plasma + 200 μ L 5 μ g/mL IS, extract with 6.0 mL dichloromethane, centrifuge. Aspirate aqueous layer to waste, evaporate organic layer under a stream of nitrogen at 35°. Reconstitute residue in 200 μ L MeOH:water 15:85, inject an aliquot.

HPLC VARIABLES

Column: 100 \times 4.6 C8 Rainin "Short One"

Mobile phase: Gradient. MeCN:MeOH:25 mM pH 3.0 phosphate buffer 3:15:82 to 15:15:70 over 10 min.

Detector: UV 273

CHROMATOGRAM

Retention time: 4.8

Internal standard: 7-(β -chloroethyl)theophylline (9.6)

KEY WORDS

plasma

REFERENCE

Sarich,T.; Kalhorn,T.; Magee,S.; Al-sayegh,F.; Adams,S.; Slattery,J.; Goldstein,J.; Nelson,S.; Wright,J.
The effect of omeprazole pretreatment on acetaminophen metabolism in rapid and slow metabolizers of S-mephenytoin, *Clin.Pharmacol.Ther.*, **1997**, 62, 21–28.

SAMPLE**Matrix:** blood**Sample preparation:** Inject a 5 μ L aliquot of serum directly.

HPLC VARIABLES**Column:** 100 \times 4.6 5-10 μ m Silicalite (by sieving Silicalite, 3M Co.(?))**Mobile phase:** MeCN:20 mM pH 6.9 phosphate buffer 11:89**Flow rate:** 1**Injection volume:** 5**Detector:** UV 254

CHROMATOGRAM**Retention time:** 2.99

OTHER SUBSTANCES**Also analyzed:** metabolites

KEY WORDS

serum

REFERENCE

Ambrose,D.L.; Fntz,J.S. High-performance liquid chromatographic determination of drugs and metabolites in human serum and urine using direct injection and a unique molecular sieve, *J.Chromatogr.B*, **1998**, 709, 89–96.

SAMPLE**Matrix:** blood**Sample preparation:** Mix 250 μ L serum with 250 μ L 800 mM perchloric acid, vortex, centrifuge at 14000 g for 3–4 min. Remove 350 μ L supernatant, adjust to ca. pH 5.0 with 27 μ L 4 M NaOH, inject an aliquot.

HPLC VARIABLES**Guard column:** 10 \times 4.6 Adsorbosphere ODS**Column:** 250 \times 4.6 5 μ m Ultrasphere ODS**Mobile phase:** MeOH:15 mM pH 4.9 potassium phosphate buffer 15:85 (After each elution flush the column with MeOH:water 80:20 for 5 min and then re-equilibrate with the for 5 min.)**Flow rate:** 1.75**Injection volume:** 100**Detector:** UV 274

CHROMATOGRAM**Retention time:** 17.75**Limit of quantitation:** 50 ng/mL

OTHER SUBSTANCES**Extracted:** metabolites**Simultaneous:** acetaminophen, aspirin, hydrocortisone, prednisolone, prednisone, phenylbutazone, phenytoin

KEY WORDS

serum

REFERENCE

Holland,D.T.; Godfredsen,K.A.; Page,T.; Connor,J.D. Simple high-performance liquid chromatography method for the simultaneous determination of serum caffeine and paraxanthine following rapid sample preparation, *J.Chromatogr.B*, **1998**, 707, 105–110.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 100 μ L 0.2 mg/mL β -hydroxyethyltheophylline in buffer + 100 μ L 40% aqueous trichloroacetic acid, vortex for 30 s, let stand for 5 min, centrifuge at 2000 g for 15 min, inject a 25 μ L aliquot of the supernatant. (Buffer was 10 mM sodium acetate adjusted to pH 4.0 with glacial acetic acid.)

HPLC VARIABLES

Column: 10 μ m μ Bondapak C18

Mobile phase: MeCN:buffer 6:94 (Buffer was 10 mM sodium acetate adjusted to pH 4.0 with glacial acetic acid.)

Column temperature: 40

Flow rate: 2

Injection volume: 25

Detector: UV 274

CHROMATOGRAM

Retention time: 13

Internal standard: β -hydroxyethyltheophylline (9)

OTHER SUBSTANCES

Extracted: dyphylline, theophylline, theobromine

KEY WORDS

plasma

REFERENCE

Valia,K.H.; Hartman,C.A.; Kucharczyk,N.; Sofia,R.D. Simultaneous determination of dyphylline and theophylline in human plasma by high-performance liquid chromatography, *J.Chromatogr.*, **1980**, 221, 170–175.

SAMPLE

Matrix: blood

Sample preparation: 50 μ L Serum + 50 μ L 15 μ g/mL β -hydroxyethyltheophylline in MeCN + 2 mL chloroform:isopropanol 95:5, mix for 30 s, centrifuge at 3000 g for 3 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 50 μ L MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: μ Bondapak C18

Mobile phase: MeCN:buffer 9.75:90.25 (Buffer was 100 mM KH_2PO_4 adjusted to pH 4.0 with phosphoric acid.) (At the end of each day clean with water for 20 min and MeOH for 30 min.)

Flow rate: 2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 8.5

Internal standard: β -hydroxyethyltheophylline (5.8)

Limit of detection: 500 ng/mL

OTHER SUBSTANCES

Extracted: dyphylline, theophylline

Simultaneous: acetaminophen, aspirin, salicylic acid, procainamide, N-acetylprocainamide

Noninterfering: benzoic acid

KEY WORDS

serum

REFERENCE

Ou, C.-N.; Frawley, V.L. Theophylline, dyphylline, caffeine, acetaminophen, salicylate, acetylsalicylate, procainamide, and N-acetylprocainamide determined in serum with a single liquid-chromatographic assay, *Clin. Chem.*, **1982**, *28*, 2157–2160.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Serum or plasma + 200 μ L 100 mM pH 7.0 phosphate buffer + 3 mL 0.5 μ g/mL 8-chlorotheophylline in isopropanol, stir at 40000 rpm for 5 s using dental micromotor with a PTFE mixing head, centrifuge at 3500 g for 2 min. Remove the supernatant and evaporate it to dryness under a stream of nitrogen at 60°, reconstitute the residue in 50 μ L MeOH, inject a 10 μ L aliquot.

HPLC VARIABLES

Guard column: 50 \times 3.2 30-38 μ m Co:Pell ODS

Column: 250 \times 4.6 5 μ m Ultrasphere ODS

Mobile phase: MeCN:MeOH:10 mM pH 5.2 sodium acetate buffer 6:3:91

Column temperature: 40

Flow rate: 1.5

Injection volume: 10

Detector: UV 274

CHROMATOGRAM

Retention time: 12.1

Internal standard: 8-chlorotheophylline (8.85)

OTHER SUBSTANCES

Extracted: dyphylline, theophylline, proxyphylline, paraxanthine

Simultaneous: cefoxitin

Noninterfering: carbenicillin, cefoperazone, cephacetril, heparin, penicillin G, phenytoin, phenobarbital

KEY WORDS

plasma; serum; pharmacokinetics

REFERENCE

Wenk, M.; Eggs, B.; Follath, F. Simultaneous determination of diprophylline, proxyphylline and theophylline in serum by reversed-phase high-performance liquid chromatography, *J. Chromatogr.*, **1983**, *276*, 341–348.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Serum + 200 μ L 4 μ g/mL 8-chlorotheophylline in MeCN, mix, centrifuge, evaporate the supernatant to dryness, reconstitute in 400 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 100 × 3.8 µm octadecyl CP-tm-Spher C18 glass column (Chrompack)

Mobile phase: MeCN:20 mM sodium acetate 20:80 adjusted to pH 4.4 with phosphoric acid

Flow rate: 0.8

Injection volume: 20

Detector: UV 273

CHROMATOGRAM

Retention time: 10.0

Internal standard: 8-chlorotheophylline (8.8)

Limit of detection: 1000 ng/mL

OTHER SUBSTANCES

Simultaneous: theophylline

KEY WORDS

serum

REFERENCE

Van Damme,M.; Molle,L.; Abi Khalil,F. Useful sample handlings for reversed phase high performance liquid chromatography in emergency toxicology, *J.Toxicol.Clin.Toxicol.*, **1985**, *23*, 589–614.

SAMPLE

Matrix: blood

Sample preparation: 100 µL Serum + 1 mL reagent, vortex for 10 min, centrifuge for 2–3 min. Remove the lower organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 200 µL mobile phase, inject a 40 µL aliquot. (Reagent was 1.5 mg β-hydroxyethyltheophylline and 25 µL glacial acetic acid in 100 mL chloroform:isopropanol 95:5.)

HPLC VARIABLES

Guard column: 30 mm long 5 µm Ultrasphere ion-pair

Column: 150 × 4.6 5 µm Ultrasphere ion-pair

Mobile phase: MeCN:MeOH:water:1 M tetra-n-butylammonium hydroxide 3.5:3.5:91:2 containing 1.82 g Trizma base (Tris, tris(hydroxymethyl)aminomethane), pH adjusted to 7.50 ± 0.03 with concentrated HCl

Flow rate: 1.2

Injection volume: 40

Detector: UV 280

CHROMATOGRAM

Retention time: k' 6.3

Internal standard: β-hydroxyethyltheophylline (k' 4.2)

OTHER SUBSTANCES

Extracted: theophylline, theobromine, 1,7-dimethyl xanthine

Simultaneous: acetaminophen, acetazolamide, allopurinol, dimethylurea, dyphylline, 3-methylxanthine, oxypurinol, procainamide, sulfadiazine, sulfamethazine, uric acid

Noninterfering: ampicillin, cefazolin, cephalothin, cephapirin, chlorotheophylline, 1,3-dimethyluric acid, gentamicin, lidocaine, methicillin, methylurea, 3-methyluric acid, quinidine, sulfamerazine, 1,3,7-trimethyluric acid

KEY WORDS

serum

REFERENCE

Lauff,J.J. Ion-pair high-performance liquid chromatographic procedure for the quantitative analysis of theophylline in serum samples, *J.Chromatogr.*, **1987**, *417*, 99–109.

SAMPLE**Matrix:** blood**Sample preparation:** 100 μ L Serum + 100 μ L buffer + 1.5 mL IS in 5% isopropanol in chloroform, vortex for 30 s, centrifuge. Remove the organic layer and evaporate it to dryness under a stream of air at room temperature, reconstitute the residue in 100 μ L mobile phase, inject a 6-10 μ L aliquot. (Buffer was 13.6 g KH_2PO_4 in 90 mL water, pH adjusted to 6.8 with about 3 mL 10 M NaOH, made up to 100 mL.)

HPLC VARIABLES**Guard column:** 20 \times 4.6 Supelguard LC-1 (Supelco)**Column:** 250 \times 4.6 5 μ m Supelcosil LC-1 (Supelco)**Mobile phase:** MeOH:MeCN:buffer 17.5:17.5:65 (Buffer was 2.72 g KH_2PO_4 in 1.9 L water, pH adjusted to 6.3 with about 2 mL 1 M NaOH, made up to 2 L.)**Flow rate:** 2**Injection volume:** 6-10**Detector:** UV 273

CHROMATOGRAM**Retention time:** 2.18**Internal standard:** 3-isobutyl-1-methylxanthine (3.15)

OTHER SUBSTANCES**Extracted:** acetaminophen, amobarbital, barbital, carbamazepine, chloramphenicol, ethosuximide, mephobarbital, methsuximide, pentobarbital, phenobarbital, phenytoin, primidone, secobarbital, theophylline, thiopental**Also analyzed:** acetanilide, N-acetylcysteine, N-acetylprocainamide, ampicillin, aspirin, butabarbital, butalbital, chlorpropamide, cimetidine, codeine, cyheptamide, diazoxide, diflunisal, diphylline, disopyramide, ethchlorvynol, gentisic acid, glutethimide, heptabarbital, hexobarbital, ibuprofen, indomethacin, ketoprofen, mefenamic acid, mephentoin, methaqualone, methsuximide, methyl salicylate, methyprylon, morphine, naproxen, nirvanol, oxphenylbutazone, phenacetin, phensuximide, phenylbutazone, procainamide, salicylamide, salicylic acid, sulfamethoxazole, sulindac, tolmetin, trimethoprim, vancomycin**Noninterfering:** amikacin, gentamicin, meprobamate, netilmicin, quinidine, tetracycline, tobramycin, valproic acid

KEY WORDS

serum

REFERENCEMeatherall, R.; Ford, D. Isocratic liquid chromatographic determination of theophylline, acetaminophen, chloramphenicol, caffeine, anticonvulsants, and barbiturates in serum, *Ther. Drug Monit.*, **1988**, *10*, 101-115.

SAMPLE**Matrix:** blood**Sample preparation:** Prepare an SPE cartridge by plugging the end of a 1 mL disposable pipette tip with glass wool and adding about 100 mg Chromosorb P/NAW. Add 50 μ L plasma then 50 μ L 10 μ g/mL tolylphenobarbital in 200 mM HCl to the SPE cartridge, let stand for 2 min, elute with 1 mL chloroform:isopropanol 6:1. Evaporate the eluate to dryness under a stream of nitrogen at 30°, reconstitute the residue in 100 μ L mobile phase, inject a 15 μ L aliquot.

HPLC VARIABLES**Column:** 150 \times 4.6 5 μ m Supelcosil-LC-8**Mobile phase:** MeCN:water 20:80**Flow rate:** 3.3**Injection volume:** 15

Detector: UV 208

CHROMATOGRAM

Retention time: 0.93

Internal standard: tolylphenobarbital (7.57)

Limit of detection: 50-100 ng/mL

OTHER SUBSTANCES

Extracted: theophylline, barbital, ethosuximide, primidone, carbamazepinediol, phenacemide, methypylon, nirvanol, phenobarbital, chloramphenicol, butabarbital, carbamazepine epoxide, mephenytoin, pentobarbital, amobarbital, carbamazepine, glutethimide, phenytoin, secobarbital, methaqualone

Noninterfering: acetaminophen, amikacin, amitriptyline, clonazepam, cyclosporine, desipramine, diazepam, digoxin, disopyramide, gentamicin, imipramine, lidocaine, methotrexate, N-acetylprocainamide, netilmicin, nortriptyline, procainamide, quinidine, salicylic acid, sulfamethoxazole, tobramycin, trimethoprim, valproic acid, p-hydroxyphenobarbital, vancomycin

KEY WORDS

plasma; SPE

REFERENCE

Svinarov,D.A.; Dotchev,D.C. Simultaneous liquid-chromatographic determination of some bronchodilators, anticonvulsants, chloramphenicol, and hypnotic agents, with Chromosorb P columns used for sample preparation, *Clin.Chem.*, **1989**, *35*, 1615-1618.

SAMPLE

Matrix: blood

Sample preparation: Inject 20 μ L serum onto column A with mobile phase A and elute to waste, after 1 min backflush the contents of column A onto column B with mobile phase B, after 1 min remove column A from the circuit, elute column B with mobile phase B, monitor the effluent from column B. Re-equilibrate column A with mobile phase A.

HPLC VARIABLES

Column: A 30 \times 4.6 IRSP silica (for preparation see Anal. Chem. 1989, 61, 2445); B 150 \times 4.6 TSK gel ODS-80TM

Mobile phase: A 20 mM NaH₂PO₄; B MeCN:100 mM NaH₂PO₄ 10:90

Flow rate: A 0.8; B 1

Injection volume: 20

Detector: UV 275

CHROMATOGRAM

Retention time: 12.5

Limit of quantitation: 1 μ g/mL

OTHER SUBSTANCES

Extracted: theobromine, theophylline

KEY WORDS

serum; column-switching

REFERENCE

Haginaka,J.; Wakai,J.; Yasuda,H.; Kimura,Y. Determination of anticonvulsant drugs and methyl xanthine derivatives in serum by liquid chromatography with direct injection: column-switching method using a new internal-surface reversed-phase silica support as a precolumn, *J.Chromatogr.*, **1990**, *529*, 455-461.

SAMPLE**Matrix:** blood**Sample preparation:** 50 μ L Serum + 10 μ L 25 μ g/mL lidocaine in 4 mM HCl, mix, add 100 μ L 1 M pH 9.0 borate buffer, add 1 mL chloroform:EtOH 82.5:17.5, mix, centrifuge. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 50 μ L mobile phase, inject an aliquot.

HPLC VARIABLES**Column:** 300 \times 2 10 μ m μ Bondapak C18**Mobile phase:** MeOH:MeCN:buffer 12:16:72 (Buffer was 31 mM sodium acetate adjusted to pH 5.1 with 40% phosphoric acid containing 0.15 mM tetrabutylammonium phosphate.)**Flow rate:** 0.3**Injection volume:** 20**Detector:** UV 230

CHROMATOGRAM**Retention time:** 5**Internal standard:** lidocaine (10)**Limit of detection:** 12.5 ng/mL

OTHER SUBSTANCES**Extracted:** cocaine**Simultaneous:** barbital, phenobarbital, flumazepil, mazindol, hexobarbital**Noninterfering:** amphetamine, desipramine, tetracaine, methadone, reserpine, buspirone, diazepam, haloperidol, chlordiazepoxide, oxazepam, midazolam, clonazepam, chlorpromazine, pentobarbital**Interfering:** nicotine, procaine, cotinine

KEY WORDSserum; rat

REFERENCELau, C.E.; Ma, F.; Falk, J.L. Simultaneous determination of cocaine and its metabolites with caffeine in rat serum microsomes by high-performance liquid chromatography, *J. Chromatogr.*, **1990**, 532, 95-103.

SAMPLE**Matrix:** blood**Sample preparation:** 100 μ L Plasma or serum + 35 μ L 22 μ g/mL 3-ethylxanthine in 20 mM pH 4.0 acetate buffer, add to a Celute-MX SPE cartridge (Jones Chromatography), let stand for 10 min, elute with two portions of isopropanol:dichloromethane 10:90, evaporate the eluate to dryness under a stream of nitrogen below 37°, reconstitute the residue in 100 μ L mobile phase B, centrifuge at 4400 rpm in a refrigerated centrifuge for 5 min, inject a 20 μ L aliquot.

HPLC VARIABLES**Column:** 250 \times 4.5 3 μ m ODS Apex I (Jones Chromatography)**Mobile phase:** Gradient. A was MeCN:THF:10 mM pH 4.0 acetate buffer 25:2:73. B was THF:10 mM pH 4.0 acetate buffer 0.01:99.99. From A:B 0:100 increasing at 2.1% A/min.**Column temperature:** 50**Flow rate:** 0.8**Injection volume:** 20**Detector:** UV 273

CHROMATOGRAM**Retention time:** 20.17**Internal standard:** 3-ethylxanthine (13.64)

OTHER SUBSTANCES

Simultaneous: theophylline, theobromine, paraxanthine, methylxanthines, dimethyluric acids, trimethyluric acid, acetaminophen

KEY WORDS

plasma; serum; SPE

REFERENCE

Leakey, T.E. Simultaneous analysis of theophylline, caffeine and eight of their metabolic products in human plasma by gradient high-performance liquid chromatography, *J. Chromatogr.*, **1990**, 507, 199–220.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Serum + 500 μ L 220 μ g/mL 4-nitroacetanilide in ethyl acetate, vortex for 30 s, centrifuge at 7000 g for 1 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 100 μ L mobile phase, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 33 \times 4.6 3 μ m C18 (Perkin-Elmer)

Mobile phase: Isopropanol:100 mM pH 5.0 sodium acetate buffer 2:98

Flow rate: 2

Injection volume: 10

Detector: UV 278

CHROMATOGRAM

Retention time: 3.85

Internal standard: 4-nitroacetanilide (8.27)

Limit of quantitation: 10000 ng/mL

OTHER SUBSTANCES

Extracted: chloramphenicol

Simultaneous: theobromine, theophylline, dipylline, chloramphenicol 3-monosuccinate

Noninterfering: acetaminophen, N-acetylprocainamide, amikacin, amitriptyline, carbamazepine, cyclosporine, digoxin, desipramine, disopyramide, ethosuximide, gentamicin, imipramine, lidocaine, lithium, methotrexate, netilmicin, nortriptyline, phenobarbital, phenytoin, primidone, procainamide, quinidine, salicylic acid, theophylline, tobramycin, valproic acid, vancomycin

KEY WORDS

serum

REFERENCE

Markin, R.S.; Wadman, M.C.; Bottjen, P.L.; Haven, M.C.; Huth, J.A. Short-column liquid chromatographic assay for caffeine and chloramphenicol in serum, *J. Chromatogr.*, **1990**, 525, 464–470.

SAMPLE

Matrix: blood

Sample preparation: Centrifuge, filter (0.45 μ m), inject an aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m internal surface reversed phase Pinkerton, silica derivatized with glycine-phenylalanine-phenylalanine (Regis) (periodically reverse the column)

Mobile phase: 100 mM pH 6.8 phosphate buffer

Flow rate: 0.3

Injection volume: 10

Detector: UV 275

CHROMATOGRAM

Retention time: 13.03

Limit of detection: <1000 ng/mL

OTHER SUBSTANCES

Extracted: dyphylline, doxofylline, theophylline

Noninterfering: acetaminophen, amitriptyline, amphetamine, atropine, benzoylecgonine, benzotropine, caffeine, carbamazepine, carisoprodol, chlorpheniramine, chlorpromazine, chlorprothixene, cimetidine, cocaine, codeine, dextromethorphan, diazepam, diphenhydramine, diphenoxilate, disopyramide, doxepin, doxylamine, emetine, erythromycin, flurazepam, glutethimide, hydrocortisone, hydromorphone, hydroxyzine, imipramine, lidocaine, loxapine, meperidine, meprobamate, methadone, methamphetamine, methapyrilene, methaqualone, methocarbamol, methylphenidate, nicotine, nordiazepam, nortriptyline, orphenadrine, papaverine, pentazocine, phenacetin, phencyclidine, phenmetrazine, phenolphthalein, phentermine, phenylpropanolamine, phenytoin, prazepam, procainamide, procaine, propoxyphene, propranolol, protriptyline, pseudoephedrine, pyrilamine, quinine, salicylamide, spironolactone, strychnine, terpin hydrate, thioridazine, thiothixene, triamterene, trifluoperazine, triflupromazine, trihexyphenidyl, trimeprazine, trimethobenzamide, trimethoprim, tripeleennamine

KEY WORDS

plasma; serum; direct injection

REFERENCE

Tagliaro, F.; Dorizzi, R.; Frigerio, A.; Marigo, M. Non-extraction HPLC method for simultaneous measurement of dyphylline and doxofylline in serum, *Clin.Chem.*, **1990**, 36, 113–115.

SAMPLE

Matrix: blood

Sample preparation: 10 μ L Plasma + 300 μ L 100 mM pH 6.0 KH_2PO_4 buffer + 100 μ L 10 μ g/mL diprophylline + 2 mL chloroform:isopropanol 50:50, vortex for 30 s, centrifuge at 2000 rpm for 10 min. Remove the organic layer and evaporate it to dryness under a stream of air at 40°, reconstitute the residue in 100 μ L mobile phase, inject a 25 μ L aliquot.

HPLC VARIABLES

Column: NovaPak C18 radial compression

Mobile phase: MeOH:MeCN:10 mM KH_2PO_4 9:2.5:90

Flow rate: 2

Injection volume: 25

Detector: UV 270

CHROMATOGRAM

Retention time: 17

Internal standard: diprophylline (10)

OTHER SUBSTANCES

Extracted: theophylline, 3-methylxanthine, theobromine

KEY WORDS

plasma

REFERENCE

Augustijns, P.; Verbeke, N. A microassay method for the determination of theophylline in biological samples using HPLC with electrochemical detection, *J.Liq.Chromatogr.*, **1992**, 15, 1303–1313.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Serum + 100 μ L 20 μ g/mL hydroxyethyltheophylline in 2 M perchloric acid, vortex, centrifuge 5 min, inject 50 μ L aliquot of supernatant.

HPLC VARIABLES

Column: 125 \times 4 LiChroSpher RP-8 5 μ m

Mobile phase: MeOH:buffer 15:85 (Buffer was 5 mL 2 M sodium acetate + 845 mL water, pH adjusted to 4.0 with acetic acid.)

Column temperature: 45

Flow rate: 1.5

Injection volume: 50

Detector: UV 282

CHROMATOGRAM

Retention time: 9.6

Internal standard: hydroxyethyltheophylline (5.6)

OTHER SUBSTANCES

Simultaneous: theophylline, chloramphenicol

KEY WORDS

serum

REFERENCE

Hannak,D.; Haux,P.; Scharbert,F.; Kattermann,R. Liquid chromatographic analysis of phenobarbital, phenytoin, and theophylline, *Wien.Klin.Wochenschr.Suppl.*, **1992**, *191*, 27–31.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma + 50 μ L 3 μ g/mL 8-chlorotheophylline in mobile phase + 100 μ L 1 M HCl + 3 mL dichloromethane, vortex for 2 min, centrifuge at 1200 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m TSK gel ODS-80TM (Tosoh)

Mobile phase: MeOH:100 mM NaH₂PO₄ 30:70

Flow rate: 0.8

Injection volume: 20

Detector: UV 274

CHROMATOGRAM

Retention time: 8

Internal standard: 8-chlorotheophylline (11.7)

Limit of quantitation: 5 ng/mL

OTHER SUBSTANCES

Extracted: theophylline, theobromine, paraxanthine, metabolites

Noninterfering: acetaminophen, phenylbutazone, phenacetin, thiamine, salicylic acid, phenobarbital, chlorpheniramine, trimethadione, hydrocortisone, prednisolone, prednisone, phenytoin, aspirin, ethenzamide

KEY WORDS

plasma

REFERENCE

Tanaka,E. Simultaneous determination of caffeine and its primary demethylated metabolites in human plasma by high-performance liquid chromatography, *J.Chromatogr.*, **1992**, 575, 311–314.

SAMPLE

Matrix: blood

Sample preparation: Dilute serum with an equal volume 7.5 µg/mL theobromine, inject a 20 µL aliquot directly.

HPLC VARIABLES

Column: 150 × 4.6 ChromSpher 5 BioMatrix (Chrompack)

Mobile phase: MeCN:water 5:95

Flow rate: 1

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 4.2

Internal standard: theobromine (3.3)

OTHER SUBSTANCES

Simultaneous: theophylline

KEY WORDS

serum

REFERENCE

Helmsing,P.J.; Huisman,R.; van der Weele,A. HPLC determination of caffeine and theophylline by direct serum injection [letter], *Clin.Chem.*, **1993**, 39, 1348–1349.

SAMPLE

Matrix: blood

Sample preparation: Condition an Merck Extrelut-3 glass extraction column with 12 mL dichloromethane the day before. 1.5 mL Serum + 100 µL 3 µg/mL N-ethylnorcotinine in water + 1.4 mL 0.5 M NaOH. Add to column, after 15 min elute under gravity with dichloromethane:isopropanol 9:1. Add 300 µL 25 mM HCl in MeOH to the organic phase and evaporate it to dryness under nitrogen. Redissolve in 100 µL water and inject a 20 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Supelcosil LC8DB

Mobile phase: Gradient. A was MeCN:water 3.6:96.4 containing 2 mL/L triethylamine, 12 mM sodium heptanesulfonate, 12 mM K₂HPO₄, 12 mM citric acid adjusted to pH 4.7 with citric acid. B was MeCN:water 19.7:80.3 containing 2 mL/L triethylamine, 12 mM sodium heptanesulfonate, 12 mM K₂HPO₄, 12 mM citric acid adjusted to pH 5.2 with citric acid. A:B 100:0 for 15 min then to 50:50 over 20 min using a concave gradient. Re-equilibrate for 15 min before next injection.

Flow rate: 1.5

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 27

Internal standard: N-ethylnorcotinine (33)

OTHER SUBSTANCES

Simultaneous: nicotine, cotinine-N-oxide, trans-3'-hydroxycotinine, norcotinine, cotinine

KEY WORDS

serum

REFERENCE

Zuccaro,P.; Altieri,I.; Rosa,M.; Passa,A.R.; Pichini,S.; Ricciarello,G.; Pacifici,R. Determination of nicotine and four metabolites in the serum of smokers by high-performance liquid chromatography with ultraviolet detection, *J.Chromatogr.*, **1993**, 621, 257-261.

SAMPLE

Matrix: blood

Sample preparation: Rock 5 mL whole blood + 10 mL water + 8.5 mL Na_2WO_4 in a 50 mL stoppered tube for 1 min, add 6 mL NiCl_2 , rock for 5 min, add 15 mL dichloromethane: isobutyl alcohol:THF 30:45:25, centrifuge at 2500 g for 15 min. Remove organic phase and repeat the process. Filter all organic phases through a 40-90 μm filter and evaporate to dryness in a 100 mL porcelain dish at a moderate temperature in a sand bath. Take up residue in 500 μL MeCN:water 80:20, inject a 20 μL aliquot. (Na_2WO_4 prepared by mixing 10 g $\text{Na}_2\text{WO}_4 \cdot 2\text{H}_2\text{O}$ in 38 mL of 2 M NaOH and 2.5 g of NaHCO_3 and making up to 100 mL. NiCl_2 was 17% w/v NiCl_2 in water.)

HPLC VARIABLES

Column: 200 \times 4.6 5 μm Hypersil C8

Mobile phase: A = MeCN; B = 20 mM n-propylamine adjusted to pH 5 with 85% phosphoric acid. A:B from 15:85 to 20:80 over 5 min to 45:55 over another 15 min to 65:35 over another 5 min

Injection volume: 20

Detector: UV 230

CHROMATOGRAM

Retention time: 5

Limit of detection: 0.10 ppm

OTHER SUBSTANCES

Extracted: buprenorphine, cocaine, codeine, diamorphine, ethylmorphine, lidocaine, methaqualone, morphine, naloxone, noscapine, papaverine, pentazocine, procaine

Also analyzed: bromazepam, clonazepam, diazepam, flunitrazepam, flurazepam, medazepam, nitrazepam, oxazepam

KEY WORDS

whole blood

REFERENCE

Bernal,J.L.; Del Nozal,M.J.; Rosas,V.; Villarino,A. Extraction of basic drugs from whole blood and determination by high performance liquid chromatography, *Chromatographia*, **1994**, 38, 617-623.

SAMPLE

Matrix: blood

Sample preparation: Wash PCPure cartridge containing 0.4 g hydroxyapatite with 10 mL MeCN and remove MeCN by evaporation. 75 μL Plasma + 25 μL of 50 $\mu\text{g/mL}$ IS in 0.5% aqueous MeCN injected onto PCPure cartridge, elute with MeCN:water 10:90. Use first 600 μL of eluate, inject 20 μL aliquots.

HPLC VARIABLES

Column: 150 \times 4.6 5 μm Inertsil ODS-2

Mobile phase: MeCN:water 5:95

Column temperature: 40

Flow rate: 1

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 15

Internal standard: 7-(2-hydroxyethyl)theophylline (10)

OTHER SUBSTANCES

Simultaneous: theophylline

KEY WORDS

plasma; SPE

REFERENCE

Iwase,H.; Gondo,K.; Koike,T.; Ono,I. Novel precolumn deproteinization method using a hydroxyapatite cartridge for the determination of theophylline and diazepam in human plasma by high-performance liquid chromatography with ultraviolet detection, *J.Chromatogr.B*, **1994**, 655, 73–81.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Plasma + 100 μ L 20 μ g/mL theophylline + 8 mL dichloromethane, shake for 20 min, centrifuge at 2500 rpm for 20 min. Remove 7 mL of the organic layer and evaporate it to dryness under nitrogen or at 60°. Dissolve residue in 200 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 6 Shimpack CLS-ODS (Shimadzu)

Mobile phase: MeCN:0.5 mM phosphoric acid 10:90

Column temperature: 40

Flow rate: 1.5

Injection volume: 20

Detector: UV 273

CHROMATOGRAM

Internal standard: theophylline

KEY WORDS

plasma; rat

REFERENCE

Lee,C.K.; Uchida,T.; Kitagawa,K.; Yagi,A.; Kim,N.-S.; Goto,S. Skin permeability of various drugs with different lipophilicity, *J.Pharm.Sci.*, **1994**, 83, 562–565.

SAMPLE

Matrix: blood

Sample preparation: 250 μ L Plasma + 25 μ L 100 μ g/mL theophylline, extract with chloroform:isopropanol 85:15. Remove the organic layer and evaporate it to dryness, reconstitute the residue in mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 100 \times 4.6 microsphere C18 (Chrompack)

Mobile phase: MeCN:55 mM pH 4.0 sodium acetate buffer 6:94

Flow rate: 1.2

Detector: UV 278

CHROMATOGRAM

Internal standard: theophylline

KEY WORDS

plasma; pig; pharmacokinetics

REFERENCE

Monshouwer,M.; Witkamp,R.F.; Nijmeijer,S.M.; Pijpers,A.; Verheijden,J.H.M.; Van Miert,A.S.J.P.A.M. Selective effects of a bacterial infection (*Actinobacillus pleuropneumoniae*) on the hepatic clearance of caffeine, antipyrine, paracetamol, and indocyanine green in the pig, *Xenobiotica*, **1995**, 25, 491–499.

SAMPLE

Matrix: blood

Sample preparation: Filter (0.22 μm), inject a 20 μL aliquot of the filtrate.

HPLC VARIABLES

Column: 250 \times 4.5 μm LiChrospher 100 Diol

Mobile phase: MeCN:50 mM pH 6.9 phosphate buffer 1.8:98.2

Flow rate: 0.6

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 6

OTHER SUBSTANCES

Extracted: theophylline

KEY WORDS

serum; direct injection

REFERENCE

Nimura,N.; Itoh,H.; Kinoshita,T. Diol-bonded silica gel as a restricted access packing forming a binary-layered phase for direct injection of serum for the determination of drugs, *J.Chromatogr.A*, **1995**, 689, 203–210.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Serum + 100 μL 70% perchloric acid, centrifuge at 2000 g for 10 min. Remove the supernatant and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100 μL water, inject a 20 μL aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 μm Supelcosil LC8DB

Mobile phase: MeCN:MeOH:THF:buffer 4.3:2.3:0.3:93.1 (Buffer was 2 mL/L triethylamine containing 12 mM sodium heptanesulfonate, 12 mM K_2HPO_4 , and 12 mM citric acid, adjusted to pH 4.7 with citric acid.)

Flow rate: 1.4

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 11

Limit of detection: 5 ng/mL

OTHER SUBSTANCES

Extracted: nicotine, cotinine

KEY WORDS

serum

REFERENCE

Pichini, S.; Altieri, I.; Passa, A.R.; Rosa, M.; Zuccaro, P.; Pacifici, R. Use of solvent optimization software for rapid selection of conditions for reversed-phase high-performance liquid chromatography of nicotine and its metabolites, *J. Chromatogr. A*, **1995**, 697, 383–388.

SAMPLE**Matrix:** blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform: isopropanol:n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES**Column:** 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30**Flow rate:** 0.8**Injection volume:** 50**Detector:** UV 273**CHROMATOGRAM****Retention time:** 3.06**Limit of detection:** <120 ng/mL**KEY WORDS**

whole blood; plasma; interferences may occur—compounds(all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; mocllobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vindsine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; tri-

fluperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254–262.

SAMPLE

Matrix: blood, CSF

Sample preparation: 200 μ L Serum, plasma, or CSF + 300 μ L reagent. Flush column A to waste with 500 μ L 500 mM ammonium sulfate, inject sample onto column A, flush column A to waste with 500 μ L 500 mM ammonium sulfate, elute the contents of column A onto column B with mobile phase, monitor the effluent from column B. (Reagent was 8.05 M guanidine hydrochloride and 1.02 M ammonium sulfate in water.)

HPLC VARIABLES

Column: A 30 \times 2.1 40 μ m preparative grade C18 (Analytichem); B 250 \times 4.6 10 μ m Partisil C8

Mobile phase: Gradient. A was 50 mM pH 4.5 KH_2PO_4 . B was MeCN:isopropanol 80:20. A: B 90:10 for 1 min, to 30:70 over 15 min, maintain at 30:70 for 4 min.

Column temperature: 50

Flow rate: 1.5

Detector: UV 280 for 5 min then UV 254

CHROMATOGRAM

Retention time: 6.16

Internal standard: heptanophenone (19.2)

OTHER SUBSTANCES

Extracted: acetazolamide, ampicillin, bromazepam, carbamazepine, chloramphenicol, chlorothiazide, diazepam, droperidol, ethionamide, furosemide, isoniazid, methadone, penicillin G, phenobarbital, phenytoin, prazepam, propoxyphene, pyrazinamide, rifampin, trimepazine, trimethoprim

KEY WORDS

plasma; serum; column-switching

REFERENCE

Seifart,H.I.; Kruger,P.B.; Parkin,D.P.; van Jaarsveld,P.P.; Donald,P.R. Therapeutic monitoring of anti-tuberculosis drugs by direct in-line extraction on a high-performance liquid chromatography system, *J.Chromatogr.*, **1993**, *619*, 285–290.

SAMPLE

Matrix: blood, CSF, gastric contents, urine

Sample preparation: 200 μ L Serum, urine, CSF, or gastric fluid + 300 μ L reagent. Flush column A to waste with 500 μ L 500 mM ammonium sulfate, inject sample onto column A, flush column A to waste with 500 μ L 500 mM ammonium sulfate, backflush the contents of column A onto column B with mobile phase, monitor the effluent from column B. (Reagent was 8.05 M guanidine HCl and 1.02 M ammonium sulfate in water.)

HPLC VARIABLES

Column: A 40 μm preparative grade C18 (Analytichem); B 75 \times 2.1 pellicular C18 (Whatman) + 250 \times 4.6 5 μm C8 end-capped (Whatman)

Mobile phase: Gradient. A was 50 mM pH 4.5 KH_2PO_4 . B was MeCN:isopropanol 80:20. A: B 90:10 for 1 min, to 30:70 over 20 min.

Column temperature: 50

Flow rate: 1.5

Detector: UV 220

CHROMATOGRAM

Retention time: 6.1

Internal standard: heptanophenone (19)

OTHER SUBSTANCES

Extracted: acetaminophen, allobarbitol, azinphos, barbital, brallobarbitone, bromazepam, butethal, carbamazepine, carbaryl, cephaloridine, chloramphenicol, chlordiazepoxide, chlorothiazide, chlorvinphos, clothiapine, cocaine, coomassie blue, desipramine, diazepam, diphenhydramine, dipipanone, ethylbromphos, flufenamic acid, formothion, griseofulvin, indomethacin, lidocaine, lorazepam, malathion, medazepam, midazolam, oxazepam, paraoxon, penicillin G, pentobarbital, prazepam, propoxyphene, prothiophos, quinine, salicylic acid, secobarbital, strychnine, sulfamethoxazole, theophylline, thiopental, thioridazine, trimethoprim

KEY WORDS

serum; column-switching

REFERENCE

Kruger,P.B.; Albrecht,C.F.De V.; Jaarsveld,P.P. Use of guanidine hydrochloride and ammonium sulfate in comprehensive in-line sorption enrichment of xenobiotics in biological fluids by high-performance liquid chromatography, *J.Chromatogr.*, **1993**, 612, 191–198.

SAMPLE

Matrix: blood, gastric contents, pancreatic juice

Sample preparation: 500 μL Serum, pancreatic juice, or gastric juice + 25 μL 200 $\mu\text{g/mL}$ β -hydroxyethyltheophylline in water+ 500 μL MeCN, centrifuge at 6000 rpm for 10 min, remove supernatant and add it to 1.8 mL chloroform, vortex, centrifuge at 6000 rpm for 10 min. Remove the lower organic phase and evaporate it to dryness at 60° under a stream of air, dissolve the residue in 500 μL mobile phase, inject a 50 μL aliquot.

HPLC VARIABLES

Guard column: C18 LichroCART

Column: 250 \times 4.6 7 μm Lichrosorb C18

Mobile phase: THF:MeOH:10 mM pH 3.5 KH_2PO_4 1:20:79

Flow rate: 0.8

Injection volume: 50

Detector: UV 280

CHROMATOGRAM

Retention time: 10.78

Internal standard: β -hydroxyethyltheophylline (7.48)

Limit of quantitation: 250 ng/mL

OTHER SUBSTANCES

Simultaneous: theophylline, 1-methyluric acid, 1-methylxanthine, theobromine, paraxanthine

KEY WORDS

serum; dog

REFERENCE

Casoli,P.; Vérine,H. High performance liquid chromatographic determination of methylxanthines in canine serum, gastric and pancreatic juices, *Biomed.Chromatogr.*, **1990**, 4, 209–213.

SAMPLE

Matrix: blood, milk

Sample preparation: 25 μL serum + 100 μL 250 ng/mL β -hydroxyethyltheophylline in water + 200 μL 200 mM pH 6.0 phosphate buffer + 3 mL dichloromethane, shake at 120 oscillations/min for 20 min, centrifuge at 174 g. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 60°, reconstitute the residue in 230 μL water, heat at 90° for 6 min, vortex, inject a 180 μL aliquot. (No details given for milk extraction.)

HPLC VARIABLES

Guard column: Corasil Bondapak C18

Column: 5 μm radial-compression C18 (Waters)

Mobile phase: MeOH:THF:10 mM KH_2PO_4 9:1:90 adjusted to pH 3.5

Flow rate: 1.2

Injection volume: 180

Detector: UV 214

CHROMATOGRAM

Retention time: 16.5

Internal standard: β -hydroxyethyltheophylline (10.9)

OTHER SUBSTANCES

Extracted: theophylline, paraxanthine, theobromine

KEY WORDS

serum; pharmacokinetics

REFERENCE

Oo,C.Y.; Burgio,D.E.; Kuhn,R.C.; Desai,N.; McNamara,P.J. Pharmacokinetics of caffeine and its metabolites in lactation: Prediction of milk to serum concentration ratios, *Pharm.Res.*, **1995**, 12, 313–316.

SAMPLE

Matrix: blood, saliva

Sample preparation: Add 500 μL plasma, serum, or saliva to 200 mg ammonium sulfate, add 50 μL 15 $\mu\text{g/mL}$ IS in water, add 500 μL 200 mM pH 4.5 sodium acetate buffer, vortex briefly, add 3 mL dichloromethane, shake at 85 cycles/min for 10 min, centrifuge at 2000 g for 10min. Evaporate the organic layer to dryness at 40° under a stream of nitrogen, reconstitute in 250 μL mobile phase, inject a 50 μL aliquot.

HPLC VARIABLES

Guard column: 20 \times 2 pellicular C18

Column: 100 \times 4.6 3 μm Microsorb MV C18

Mobile phase: MeOH:THF:100 mM pH 4.5 sodium acetate:water 6.5:1.4:10:82.1

Flow rate: 0.8

Injection volume: 50

Detector: UV 274

CHROMATOGRAM

Retention time: 9.1

Internal standard: β -hydroxyethyltheophylline (6.5)

Limit of detection: 5 ng/mL

OTHER SUBSTANCES

Extracted: metabolites, paraxanthine, theobromine, theophylline

Noninterfering: chlorzoxazone, dapsone

KEY WORDS

plasma; serum

REFERENCE

Frye,R.F.; Stiff,D.D.; Branch,R.A. A sensitive method for the simultaneous determination of caffeine and its dimethylxanthine metabolites in human plasma: Application to CYP1A2 phenotyping, *J.Liq.Chromatogr.Rel.Technol.*, **1998**, 21, 1161–1171.

SAMPLE

Matrix: blood, saliva

Sample preparation: Saliva. Add 3 mg N-acetylcysteine (mucolytic agent) to 200 μ L saliva, let stand for 10-15 min. 50 μ L Saliva + 50 μ L 1 M NaOH, extract with 1 mL 10% (v/v) isopropanol in chloroform (Caution! Chloroform is a carcinogen!) containing 1.5 mg/L β -hydroxyethyl theophylline. Centrifuge at 2000 g for 1 min. Evaporate organic layer to dryness under a gentle airflow at 60°. Dissolve residue in 200 μ L mobile phase, inject a 50 μ L aliquot. Serum. 50 μ L Serum + 50 μ L 1 M NaOH, extract with 1 mL 10% isopropanol in chloroform containing 1.5 mg/L β -hydroxyethyl theophylline. Centrifuge at 2000 g for 1 min. Evaporate organic layer to dryness under a gentle airflow at 60°. Dissolve the residue in 200 μ L mobile phase, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 3.9 C18 Nova Pak

Mobile phase: MeOH:THF:50 mM pH 4.7 ammonium acetate buffer 2:1.5:96.5

Flow rate: 1

Injection volume: 50

CHROMATOGRAM

Retention time: 10.2

Internal standard: β -hydroxyethyl theophylline (6.6)

Limit of quantitation: 200 μ g/L

OTHER SUBSTANCES

Extracted: paraxanthine, theobromine, theophylline

KEY WORDS

serum

REFERENCE

Lee,T.C.; Charles,B.G.; Steer,P.A.; Flenady,V.J. Saliva as a valid alternative to serum in monitoring intravenous caffeine treatment for apnea of prematurity, *Ther.Drug Monit.*, **1996**, 18, 288–293.

SAMPLE

Matrix: blood, tissue

Sample preparation: Tissue. Homogenize brain in twenty volumes of water, centrifuge at 1500 g, freeze at -20° for 2 h, centrifuge. 200 μ L Supernatant + 100 μ L pH 7.2 phosphate buffer, mix, add 6 mL dichloromethane:propanol 95:5, mix for 2 min, centrifuge at 1500 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 35°, reconstitute the residue in 200 μ L mobile phase, vortex for 1 min, inject a 25 μ L aliquot. Serum. 200 μ L Serum + 100 μ L pH 7.2 phosphate buffer, mix, add 6 mL dichloromethane:propanol 95:5, mix for 2 min, centrifuge at 1500 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 35°, reconstitute the residue in 200 μ L mobile phase, vortex for 1 min, inject a 25 μ L aliquot.

HPLC VARIABLES

Guard column: 10 μ m μ Bondapak C18

Column: 250 × 4.6 10 µm µBondapak C18

Mobile phase: THF:10 mM Na₂HPO₄ 3:97, pH adjusted to 6.5 with phosphoric acid

Flow rate: 2.5

Injection volume: 25

Detector: UV 273

CHROMATOGRAM

Retention time: 7

Limit of detection: 62.5 ng/mL (serum), 62.5 ng/g (tissue)

OTHER SUBSTANCES

Extracted: theophylline, theobromine, paraxanthine

KEY WORDS

serum; rat; brain

REFERENCE

Parra,P.; Limon,A.; Ferre,S.; Guix,T.; Jane,F. High-performance liquid chromatographic separation of caffeine, theophylline, theobromine and paraxanthine in rat brain and serum, *J.Chromatogr.*, **1991**, 570, 185–190.

SAMPLE

Matrix: blood, urine

Sample preparation: Urine. 100 µL Urine + 100 µL 1 M NaOH, let stand at room temperature for 10 min, add 100 µL 1 M HCl, make up to 1 mL with mobile phase, vortex, centrifuge at 13000 g for 2 min, inject a 20 µL aliquot. Plasma. 500 µL Plasma + 100 µL 60% perchloric acid, add 200 µL 1 M NaOH to make pH >9, let stand for 10 min, add 200 µL 1 M HCl, add 200 µL 50 µg/mL β-hydroxyethyltheophylline, add 6 mL chloroform: isopropanol 50:50, agitate for 15 min on a rotating mixer, centrifuge at 4000 g. Remove the organic layer and evaporate it to dryness under vacuum, reconstitute the residue in 100 µL mobile phase, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 100 × 4.5 5 µm ChromSep (Chrompack)

Mobile phase: THF:10 mM sodium acetate 0.8:99.2

Flow rate: 1

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 7.03

Internal standard: β-hydroxyethyltheophylline (4.76)

Limit of quantitation: 1000 ng/mL

OTHER SUBSTANCES

Extracted: 1,3,7-trimethyluric acid, 1,7-dimethylxanthine, 1,3-dimethylxanthine, metabolites

KEY WORDS

plasma

REFERENCE

Dobrocky,P.; Bennett,P.N.; Notarianni,L.J. Rapid method for the routine determination of caffeine and its metabolites by high-performance liquid chromatography, *J.Chromatogr.B*, **1994**, 652, 104–108.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A: B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 205.2

CHROMATOGRAM

Retention time: 6.647

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

SAMPLE

Matrix: bulk

Sample preparation: Prepare a 750 µg/mL solution in 10 mM pH 2.5 orthophosphoric acid, sonicate for 10 min, filter (0.2 µm), inject a 15 µL aliquot.

HPLC VARIABLES

Guard column: 4 × 4 5 µm LiChrospher 100

Column: 125 × 4 3 µm Spherisorb ODS-1

Mobile phase: Gradient. A was water containing 5 mL/L 85% orthophosphoric acid and 0.56 mL/L hexylamine. B was MeCN:water 90:10 containing 5 mL/L 85% orthophosphoric acid and 0.56 mL/L hexylamine. A:B from 91:9 to 86:14 over 4 min, maintain at 86:14 for 13 min, to 55:45 over 11 min, maintain at 55:45 for 8 min, re-equilibrate at initial conditions for 20 min.

Flow rate: 0.7

Injection volume: 15

Detector: UV 210

CHROMATOGRAM

Retention time: 14.3

OTHER SUBSTANCES

Simultaneous: acetaminophen, acetylcodeine, benzocaine, cocaine, codeine, diamorphine, lidocaine, 6-monoacetylmorphine, morphine, noscapine, papaverine, procaine

REFERENCE

Grogg-Sulser,K.; Helmlin,H.-J.; Clerc,J.-T. Qualitative and quantitative determination of illicit heroin street samples by reversed-phase high-performance liquid chromatography: method development by CARTAGO-S, *J.Chromatogr.A*, **1995**, 692, 121–129.

SAMPLE

Matrix: dialysate

Sample preparation: Dialyze blood with Ringer's solution, inject a 0.5 μ L aliquot of the dialysate.

HPLC VARIABLES

Column: 14 \times 1.3 μ m BAS Sep-Stik ODS (Bioanalytical Systems)

Mobile phase: MeCN:50 mM pH 2.5 ammonium phosphate buffer 5:95

Flow rate: 0.2

Injection volume: 0.5

Detector: UV 273

CHROMATOGRAM

Retention time: 0.77

Limit of quantitation: 300 ng/mL

OTHER SUBSTANCES

Extracted: metabolites, theobromine, paraxanthine

KEY WORDS

rat; microbore

REFERENCE

Chen,A.; Lunte,C.E. Microdialysis sampling coupled on-line to fast microbore liquid chromatography, *J.Chromatogr.A*, **1995**, 691, 29–35.

SAMPLE

Matrix: food

Sample preparation: Dilute 5 mL beverage with 5 mL mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 Partisil PXS 5/25 ODS-3 or 150 \times 4.6 5 μ m Ultrasphere ODS

Mobile phase: MeOH:water:33% acetic acid 50:49:1

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 5

Internal standard: aspirin (6.6)

KEY WORDS

drink; beverage; soft drink; tea; coffee

REFERENCE

Galasko,G.T.F.; Furman,K.I.; Alberts,E. The caffeine content of non-alcoholic beverages, *Food Chem.Toxicol.*, **1989**, 27, 49–51.

SAMPLE

Matrix: formulations

Sample preparation: Powder tablets or capsules, weigh 20 mg, add 100 (tablets) or 500 (capsules) mL 50 mM sodium dodecyl sulfate, sonicate for 3 min, filter (sintered glass 4), dilute to a concentration of ca. 6 µg/mL with 50 mM sodium dodecyl sulfate, inject an aliquot.

HPLC VARIABLES

Guard column: 35 × 4.6 5 µm Spherisorb ODS-2 C18

Column: 120 × 4.5 µm Spherisorb ODS-2 C18

Mobile phase: Isopropanol:50 mM sodium dodecyl sulfate 1.5:98.5, pH 7.0

Flow rate: 1

Injection volume: 20

Detector: UV 273

CHROMATOGRAM

Retention time: 11.5

KEY WORDS

tablets; capsules

REFERENCE

Perez-Martinez,I.; Sagrado,S.; Medina-Hernández,M.J. Chromatographic determination of caffeine in pharmaceutical formulations using micellar mobile phases, *Chromatographia*, **1996**, 43, 149–152.

SAMPLE

Matrix: formulations

Sample preparation: Add 50 mL of mobile phase to 0.5 g of sample and swirl to aid dissolution. Dilute to 100 mL with mobile phase. Dilute 1:10, filter (0.22 µm nylon). Inject a 10 µL aliquot.

HPLC VARIABLES

Column: 100 × 2.1 5 µm Hypersil ODS

Mobile phase: MeCN:water:triethylamine:acetic acid 5.5:94.1:0.2:0.2

Flow rate: 1.5

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 2.2

Limit of quantitation: 5 µg/mL

OTHER SUBSTANCES

Simultaneous: acetaminophen, aspirin

KEY WORDS

powder

REFERENCE

Ferguson,G.K. Quantitative HPLC analysis of an analgesic/caffeine formulation: Determination of caffeine, *J.Chem.Educ.*, **1998**, 75, 467–469.

SAMPLE

Matrix: formulations

Sample preparation: Weigh 500 mg homogenized analgesic powder, transfer to 100 mL volumetric flask, add ca. 50 mL mobile phase, swirl and dilute to volume with mobile phase. Dilute an aliquot of this solution 1:10 with mobile phase, filter (0.20 µm Cameo nylon filter, MSI, Westboro, MA) an aliquot, inject an aliquot of the filtrate.

HPLC VARIABLES

Column: 100 × 2.1 5 µm Hypersil ODS

Mobile phase: MeCN:triethylamine:acetic acid:water 5.5:0.2:0.2:94.1 (Prepare mobile phase as follows. Mix 110 mL MeCN, 4 mL triethylamine, 4 mL glacial acetic acid and make up to 2 L with water.)

Flow rate: 1.5

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 2.2

OTHER SUBSTANCES

Extracted: aspirin, acetaminophen

Noninterfering: salicylic acid

KEY WORDS

powder

REFERENCE

Ferguson, G.K. Quantitative HPLC analysis of an analgesic/caffeine formulation: Determination of caffeine, *J.Chem.Educ.*, **1998**, 75, 467–469.

SAMPLE

Matrix: formulations

Sample preparation: Sonicate 75 mg powdered tablets with 25 mL mobile phase for 15 min, filter (paper), inject a 135 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Ultrabase C18 (Scharlau Science, Spain)

Mobile phase: MeOH:20 mM pH 4.0 KH₂PO₄ 30:70 adjusted to pH 4.0 with orthophosphoric acid

Flow rate: 1.5

Injection volume: 135

Detector: UV 274

CHROMATOGRAM

Retention time: 6.6

Limit of quantitation: 2.2 µg/mL

OTHER SUBSTANCES

Simultaneous: aspirin, salicylic acid, thiamine

KEY WORDS

tablets

REFERENCE

Gámiz-Gracia, L.; Luque de Castro, M.D. An HPLC method for the determination of vitamin B1, caffeine, acetylsalicylic acid, and the impurities of salicylic acid in a pharmaceutical preparation, *J.Liq.Chromatogr.Rel.Technol.*, **1997**, 20, 2123–2133.

SAMPLE

Matrix: formulations

Sample preparation: Oils. 1 mL Sample + 25 mL MeOH:water 90:10, shake vigorously for 5 min, centrifuge, inject a 10 µL aliquot of the supernatant. Tablets. Grind a tablet to a fine powder, add 25 mL MeOH, sonicate for 5-10 min, filter (0.45 µm), discard first 5 mL of filtrate, inject a 10 µL aliquot of the remaining filtrate. Suspensions (aqueous).

Make up 5 mL to 50 mL with MeOH, filter (0.45 μm), discard first 5 mL of filtrate, inject a 10 μL aliquot of the remaining filtrate.

HPLC VARIABLES

Column: 250 \times 4.6 5 μm Zorbax ODS

Mobile phase: MeOH:water 75:25

Flow rate: 1.5

Injection volume: 10

Detector: UV 240

CHROMATOGRAM

Retention time: 2.6

Limit of detection: 5 $\mu\text{g/mL}$

OTHER SUBSTANCES

Simultaneous: benzyl alcohol, cortisone, fluoxymesterone, norethindrone, oxandrolone (UV 210), boldenone, ethisterone, methandrostenolone, nandrolone, norgestrel, testosterone, dehydroepiandrosterone (UV 210), mibolerone, methyltestosterone, methandriol (UV 210), norethindrone acetate, calusterone, mesterolone (UV 210), norethandrolone, trenbolone acetate, benzyl benzoate, nandrolone acetate, testosterone acetate, stanozolol, oxymetholone, nandrolone propionate, methenolone acetate, testosterone propionate

Interfering: aspirin, formebolone, testolactone

KEY WORDS

oils; tablets; suspensions

REFERENCE

Walters, M.J.; Ayers, R.J.; Brown, D.J. Analysis of illegally distributed anabolic steroid products by liquid chromatography with identity confirmation by mass spectrometry or infrared spectrophotometry, *J. Assoc. Off. Anal. Chem.*, **1990**, 73, 904–926.

SAMPLE

Matrix: formulations

HPLC VARIABLES

Column: 100 \times 4 3 μm Hypersil BDS-C18

Mobile phase: Gradient. MeCN:water adjusted to pH 2.1 from 0.3:99.7 to 25:75 over 11 min

Flow rate: 0.5

Detector: UV 220

CHROMATOGRAM

Retention time: 9

OTHER SUBSTANCES

Simultaneous: biotin, citric acid, folic acid, niacinamide, niacin, pantothenic acid, riboflavin, saccharin, thiamine, pyridoxine, vitamin B12, ascorbic acid

KEY WORDS

tablets

REFERENCE

Hewlett Packard Leaflet 12-5091-7351 EUS, 1993, **1993**,

SAMPLE

Matrix: formulations

Sample preparation: Powder tablet and add 50 mg to 50 mL MeCN:20 mM pH 3.8 phosphate buffer 3:97, sonicate for 5 min, filter (0.5 μ m), inject a 20 μ L aliquot of the filtrate.

HPLC VARIABLES

Guard column: Supelguard pre-column containing 5 μ m Suplex pKb100 (Supelco)

Column: 150 \times 4.6 5 μ m Suplex pKb100 (Supelco)

Mobile phase: Gradient. MeCN:20 mM pH 3.8 phosphate buffer at 3:97 for 3 min, to 15:85 over 5 min, stay at 15:85 for 4 min, re-equilibrate for 8 min.

Flow rate: 1.5

Injection volume: 20

Detector: UV 220 for 5 min then UV 280

CHROMATOGRAM

Retention time: 10.40

Limit of quantitation: 3 μ g/mL

OTHER SUBSTANCES

Simultaneous: methamphetamine, amphetamine, ephedrine, 3,4-methylenedioxyamphetamine, N-methyl-3,4-methylenedioxyamphetamine, N-ethyl-3,4-methylenedioxyamphetamine

KEY WORDS

tablets

REFERENCE

Longo,M.; Martines,C.; Rolandi,L.; Cavallaro,A. Simple and fast determination of some phenethylamines in illicit tablets by base-activated reversed phase HPLC, *J.Liq.Chromatogr.*, **1994**, 17, 649–658.

SAMPLE

Matrix: formulations

Sample preparation: Finely powder half a tablet, add 9 mL mobile phase, sonicate for 20 min, make up to 10 mL with mobile phase, filter (Whatman type 40 and 0.2 μ m Millipore), inject an aliquot of the filtrate.

HPLC VARIABLES

Column: 250 \times 4 5 μ m LiChrospher 100 CN

Mobile phase: MeCN:THF:buffer 7:6:87 (Buffer was 0.8% acetic acid containing 5 mM sodium hexanesulfonate, 10 mM di-n-butylamine, and 0.12% phosphoric acid, pH 3.3.)

Flow rate: 1

Injection volume: 20

Detector: UV 298

CHROMATOGRAM

Retention time: 3.8

Limit of detection: 1.1 μ g/mL

OTHER SUBSTANCES

Simultaneous: acetaminophen (UV 310), chlorpheniramine (UV 265), guaifenesin (glycerylguaiacolate) (UV 284), phenylpropanolamine (UV 260)

KEY WORDS

tablets

REFERENCE

Indrayanto,G.; Sunarto,A.; Adriani,Y. Simultaneous assay of phenylpropanolamine hydrochloride, caffeine, paracetamol, glycerylguaiacolate and chlorpheniramine in SilabatTM tablet using HPLC with diode array detection, *J.Pharm.Biomed.Anal.*, **1995**, 13, 1555–1559.

SAMPLE**Matrix:** formulations**Sample preparation:** Finely powder tablets, weigh out amount equivalent to 20 mg enalapril maleate, suspend in 100 mL mobile phase, filter, inject a 5 μ L aliquot.

HPLC VARIABLES**Column:** 250 \times 4.6 12 μ m Hypersil C18**Mobile phase:** MeCN:water 20:80 adjusted to pH 3.8 with acetic acid**Flow rate:** 1**Injection volume:** 5**Detector:** UV 215 for 3.5 min, then UV 275

CHROMATOGRAM**Retention time:** 4.8

OTHER SUBSTANCES**Simultaneous:** enalapril, hydrochlorothiazide

KEY WORDS

tablets

REFERENCE

el Walily,A.F.M.; Belal,S.F.; Heaba,E.A.; El Kersh,A. Simultaneous determination of enalapril maleate and hydrochlorothiazide by first-derivative ultraviolet spectrophotometry and high-performance liquid chromatography, *J.Pharm.Biomed.Anal.*, **1995**, 13, 851–856.

SAMPLE**Matrix:** formulations**Sample preparation:** Weigh out powdered sample containing 8 mg caffeine, add 80 mL MeOH, sonicate for 10 min, dilute to 100 mL with MeOH, centrifuge. Remove a 5 mL aliquot of the supernatant and add it to 1 mL 2 mg/mL resorcinol, add 2 mL MeOH, make up to 20 mL with 50 mM pH 3.0 triethylamine phosphate, inject an aliquot.

HPLC VARIABLES**Column:** 150 \times 3.2 5 μ m Hypersil ODS**Mobile phase:** THF:50 mM pH 3.0 triethylamine phosphate 12:88**Flow rate:** 0.6**Injection volume:** 20**Detector:** UV 275

CHROMATOGRAM**Retention time:** 15**Internal standard:** resorcinol (9)

OTHER SUBSTANCES**Simultaneous:** acetaminophen, aspirin (post-column irradiation gives an increase in peak height), propyphenazone

REFERENCE

Di Pietra,A.M.; Gatti,R.; Andrisano,V.; Cavrini,V. Application of high-performance liquid chromatography with diode-array detection and on-line post-column photochemical derivatization to the determination of analgesics, *J.Chromatogr.A*, **1996**, 729, 355–361.

SAMPLE**Matrix:** microsomal incubations**Sample preparation:** 1 mL Microsomal incubation + 1 mL chloroform:isopropanol 85:15 (Caution! Chloroform is a carcinogen!). Add 25 pmole 4-acetoimidophenol. (Metabolites

were extracted with the addition of an additional 3.5 mL chloroform:isopropanol 85:15 and 500 mg ammonium sulfate before centrifugation.). Mix, centrifuge at 3000 rpm for 10 min. Remove the organic phase, evaporate it and reconstitute the residue with 150 μ L mobile phase A. Inject an aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Beckman C18 (Phenomenex, USA)

Mobile phase: Gradient. A was THF:50 mM pH 4.5 sodium acetate 1:99. B was MeCN:50 mM pH 4.5 sodium acetate 60:40. A:B 100:0 for 1 min, to 70:30 over 23 min, to 0:100 over 1 min, maintain at 0:100 for 5 min, to 100:0 over 1 min, maintain at 100:0 for 9 min.

Flow rate: 0.8

Detector: UV 270

CHROMATOGRAM

Retention time: 18

Internal standard: 4-acetoimidophenol (16)

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

liver

REFERENCE

Ring,B.J.; Catlow,J.; Lindsay,T.J.; Gillespie,T.; Roskos,L.K.; Cerimele,B.J.; Swanson,S.P.; Hamman,M.A.; Wrighton,S.A. Identification of the human cytochromes P450 responsible for the in vitro formation of the major oxidative metabolites of the antipsychotic agent olanzapine, *J.Pharmacol.Exp.Ther.*, **1996**, 276, 658–666.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 275 μ L Microsomal incubation + 50 μ L 30% perchloric acid, centrifuge at 2000 g for 10 min, inject a 60 μ L aliquot of the supernatant.

HPLC VARIABLES

Guard column: 20 \times 4.6 40 μ m Supelcosil LC-18

Column: 150 \times 4.6 5 μ m Supelcosil LC-18

Mobile phase: MeCN:THF:20 mM pH 3.5 sodium perchlorate 0.5:0.5:99

Flow rate: 1.5

Injection volume: 60

Detector: UV 280

CHROMATOGRAM

Retention time: 35

OTHER SUBSTANCES

Extracted: paraxanthine, theobromine, theophylline, trimethylurate

KEY WORDS

monkey; liver

REFERENCE

Bullock,P.; Pearce,R.; Draper,A.; Podval,J.; Bracken,W.; Veltman,J.; Thomas,P.; Parkinson,A. Induction of liver microsomal cytochrome P450 in cynomolgus monkeys, *Drug Metab.Dispos.*, **1995**, 23, 736–748.

SAMPLE**Matrix:** milk**Sample preparation:** Centrifuge at 12800 g for 10 min, remove top layer of fat with a spatula. 450 μ L Supernatant + 50 μ L water + 150 μ L 21.65 μ g/mL proxyphylline in 6% perchloric acid, vortex for 5 s, cool on ice for 10-15 min, centrifuge at 12800 g for 10 min, inject a 50 μ L aliquot of the supernatant

HPLC VARIABLES**Guard column:** Co:Pell ODS glass beads (Whatman)**Column:** 250 \times 4.6 5 μ m Ultrasphere ODS**Mobile phase:** Gradient. A was 10 mM sodium acetate + 5 mM tetrabutylammonium hydrogen sulfate, adjusted to pH 4.9 with 2 M NaOH. B was MeOH:water 50:50 containing 10 mM sodium acetate + 5 mM tetrabutylammonium hydrogen sulfate, adjusted to pH 4.8 with glacial acetic acid. A:B 100:0 for 7.5 min, then to 85:15 over 7.5 min, then to 70:30 over 10 min, then to 68:32 over 4 min, then to 100:0 over 3 min**Flow rate:** 1.5**Injection volume:** 50**Detector:** UV 272

CHROMATOGRAM**Retention time:** 30**Internal standard:** proxyphylline (25)

OTHER SUBSTANCES**Simultaneous:** theobromine, paraxanthine, theophylline

REFERENCEBlanchard,J.; Weber,C.W.; Shearer,L.E. HPLC analysis of methylxanthines in human breast milk, *J.Chromatogr.Sci.*, **1990**, 28, 640-642.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 250 \times 4 ODS (Hitachi)**Mobile phase:** MeCN:50 mM phosphoric acid 18:82**Column temperature:** 55**Flow rate:** 0.6**Injection volume:** 20**Detector:** UV 273

REFERENCESugawara,M.; Takekuma,Y; Yamada,H.; Kobayashi,M.; Iseki,K.; Miyazaki,K. A general approach for the prediction of the intestinal absorption of drugs: regression analysis using the physicochemical properties and drug-membrane electrostatic interactions, *J.Pharm.Sci.*, **1998**, 87, 960-966.

SAMPLE**Matrix:** solutions**Sample preparation:** Dissolve in MeOH at a concentration of 1 mg/mL, inject a 20 μ L aliquot.

HPLC VARIABLES**Column:** 250 \times 5 Spherisorb S5W**Mobile phase:** MeOH:buffer 90:10 (Buffer was 94 mL 35% ammonia and 21.5 mL 70% nitric acid in 884 mL water, adjust the pH to 10.1 with ammonia.)**Flow rate:** 2**Injection volume:** 20

Detector: UV 254

CHROMATOGRAM

Retention time: 1.66

OTHER SUBSTANCES

Simultaneous: levallorphan, hydroxypethidine, normethadone, meperidine, dipipanone, diamorphine, pentazocine, acetylcodeine, monoacetylmorphine, thebacon, oxycodone, thebaine, norlevorphanol, methadone, benzylmorphine, ethylmorphine, morphine-N-oxide, codeine, codeine-N-oxide, morphine, ethoheptazine, morphine-3-glucuronide, pholcodeine, norpethidine, hydrocodone, dihydrocodeine, dihydromorphine, levorphanol, norcodeine, normorphine, epinephrine, pipradol, phenylpropanolamine, fencamfamin, chlorphentermine, norpseudoephedrine, phentermine, fenfluramine, methylenedioxyamphetamine, amphetamine, normetanephine, 4-hydroxyamphetamine, bromo-STP, STP, prolintane, 2-phenethylamine, tyramine, trimethoxyamphetamine, phenylephrine, pseudoephedrine, ephedrine, methylephedrine, dimethylamphetamine, methamphetamine, mescaline, mephentermine, buprenorphine, dextromoramide, phenoperidine

Noninterfering: dopamine, levodopa, methyl dopa, methyl dopate, norepinephrine

Interfering: pemoline, benzphetamine, diethylpropion, mazindol, tranlycypromine, fenethyline, phendimetrazine, methylphenidate, phenelzine, fentanyl, etorphine, piritramide, noscapine, papaverine, naloxone, dextropropoxyphene, nalorphine, phenazocine, norpipanone

REFERENCE

Law,B.; Gill,R.; Moffat,A.C. High-performance liquid chromatography retention data for 84 basic drugs of forensic interest on a silica column using an aqueous methanol eluent, *J.Chromatogr.*, **1984**, *301*, 165-172.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve compounds in MeCN:water 80:20, inject a 1 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 1.3 μ m Hitachi-Gel 3057 ODS silica (Hitachi)

Mobile phase: MeCN:water 25:75

Flow rate: 0.03

Injection volume: 1

Detector: UV 254

CHROMATOGRAM

Retention time: 5

OTHER SUBSTANCES

Simultaneous: dipyrone (sulpyrin), guaifenesin (guaiacol glycerol ether), acetaminophen, buccetin (3-hydroxy-p-butyrophenetidine), methyl p-hydroxybenzoate, phenacetin

KEY WORDS

semi-micro

REFERENCE

Matsushima,Y.; Nagata,Y.; Niyomura,M.; Takakusagi,K.; Takai,N. Analysis of antipyretics by semimicro liquid chromatography, *J.Chromatogr.*, **1985**, *332*, 269-273.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in MeOH:water 1:1 at a concentration of 50 μ g/mL, inject a 10 μ L aliquot.

HPLC VARIABLES**Column:** 300 × 3.9 10 µm µBondapak C18**Mobile phase:** MeOH:acetic acid:triethylamine:water 25:1.5:0.5:73**Flow rate:** 1.5**Injection volume:** 10**Detector:** UV 254

CHROMATOGRAM**Retention time:** 11

OTHER SUBSTANCES**Simultaneous:** theobromine, diphylline, theophylline, 8-chlorotheophylline

REFERENCE

Roos,R.W.; Lau-Cam,C.A. General reversed-phase high-performance liquid chromatographic method for the separation of drugs using triethylamine as a competing base, *J.Chromatogr.*, **1986**, 370, 403–418.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 150 × 4.6 Supelcosil LC-ABZ**Mobile phase:** MeCN:25 mM pH 6.9 potassium phosphate buffer 35:65**Flow rate:** 1.5**Injection volume:** 25**Detector:** UV 254

CHROMATOGRAM**Retention time:** 1.329

OTHER SUBSTANCES

Also analyzed: 6-acetylmorphine, amiloride, amphetamine, benzocaine, benzoylecgonine, cocaine, codeine, doxylamine, fluoxetine, glutethimide, hexobarbital, hypoxanthine, levorphanol, LSD, meperidine, mephobarbital, methadone, methylphenidate, methyprylon, N-norcodeine, oxazepam, oxycodone, phenylpropanolamine, prilocaine, procaine, terfenadine

REFERENCE

Ascah,T.L. Improved separations of alkaloid drugs and other substances of abuse using Supelcosil LC-ABZ column, *Supelco Reporter*, **1993**, 12(3), 18–21.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Guard column:** 30 × 3.2 7 µm SI 100 ODS (not commercially available)**Column:** 150 × 3.2 7 µm SI 100 ODS (not commercially available)**Mobile phase:** MeCN:buffer 31.2:68.8 (Buffer was 6.66 g KH₂PO₄ and 4.8 g 85% phosphoric acid in 1 L water, pH 2.3.)**Flow rate:** 0.5-1**Detector:** UV 225, 267

CHROMATOGRAM**Retention time:** 1.5**Internal standard:** 5-(4-methylphenyl)-5-phenylhydantoin (7.3)

OTHER SUBSTANCES

Also analyzed: aspirin, carbamazepine, chlordiazepoxide, chlorprothixene, clonazepam, diazepam, doxylamine, ethosuximide, furosemide, haloperidol, hydrochlorothiazide, methocarbamol, methotrimeprazine, nicotine, oxazepam, procaine, promazine, propafenone, propranolol, salicylamide, temazepam, tetracaine, thiopental, triamterene, verapamil, zolpidem, zopiclone

REFERENCE

Below, E.; Burrmann, M. Application of HPLC equipment with rapid scan detection to the identification of drugs in toxicological analysis, *J. Liq. Chromatogr.*, **1994**, 17, 4131-4144.

SAMPLE

Matrix: solutions

Sample preparation: Dilute a 10 mg/mL solution of caffeine in water 1:1000 with mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 2540 \times 4.6 5 μ m C18 (Supelco)

Mobile phase: MeOH:10 mM ammonium acetate and 2.5 mM sodium heptanesulfonate 20:80, pH adjusted to 5.1 with glacial acetic acid

Flow rate: 2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 10.5

OTHER SUBSTANCES

Simultaneous: benzoic acid

KEY WORDS

stability-indicating; water

REFERENCE

Donnelly, R.F.; Tirona, R.G. Stability of citrated caffeine injectable solution in glass vials, *Am. J. Hosp. Pharm.*, **1994**, 51, 512-514.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, carbamazepine, carbromal, chloramphenicol, chlordiazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheni-

ramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapson, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenoprofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiacol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephen-
termine, mephenytoin, mephesis, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methyprylon, metoprolol, mibolone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscipine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sufadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, **1994**, *18*, 233-242.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4 5 µm LiChrospher 100 RP-8 (B)

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 4.79 (A), 3.84 (B)

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, carbamazepine, cetirizine, chlorcyclizine, chlordiazepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyridamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazindol, mefenamic acid, meperidine, mephénytoin, mepivacaine, mesoridazine, metaproterenol, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl-dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfapyrazone, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocinide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

KEY WORDS

also details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J.Chromatogr.A*, **1995**, 692, 103–119.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 150 mm long 5 μ m Microsorb-MV C18

Mobile phase: MeOH:water 60:40

Flow rate: 1

Detector: UV 270

REFERENCE

Phillips, C.A.; Michniak, B.B. Transdermal delivery of drugs with differing lipophilicities using azone analogs as dermal penetration enhancers, *J.Pharm.Sci.*, **1995**, 84, 1427–1433.

SAMPLE

Matrix: urine

Sample preparation: 40 μ L Urine + 100 μ L 160 μ M IS in 30% EtOH, make up to 400 μ L with 10 mM pH 4.0 acetate buffer. Add 5 mL ethyl acetate:2-propanol 93:7, shake for 10

min, centrifuge at 1000 g for 10 min, freeze the aqueous phase at -30°. Evaporate organic phase to dryness under a stream of nitrogen at 55°. Reconstitute the residue in 300 µL mobile phase, vortex for 10 s, centrifuge at 1000 g for 30 s, inject a 40 µL aliquot. (Keep sample in the autosampler at 5°.)

HPLC VARIABLES

Column: 250 × 4.6 5 µm Ultrasphere ODS

Mobile phase: MeOH:10 mM pH 4.0 acetate buffer 9:91

Column temperature: 30

Flow rate: 1 from 0 to 11 min, 1.5 from 11 to 17 min, 2.5 from 17 to 30 min

Injection volume: 40

Detector: UV 273

CHROMATOGRAM

Retention time: 25

Internal standard: β-hydroxy-ethyl theophylline (19.8)

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

pharmacokinetics

REFERENCE

Rasmussen,B.B.; Brosen,K. Determination of urinary metabolites of caffeine for the assessment of cytochrome P4501A2, xanthine oxidase, and N-acetyltransferase activity in humans, *Ther.Drug Monit.*, **1996**, *18*, 254–262.

SAMPLE

Matrix: urine

Sample preparation: Add 100-120 mg NaCl, 50 µL 100 µg/mL β-hydroxyethyltheophylline, and 100 µL ammonia buffer to 1 mL urine. Extract with 5 mL MeOH:dichloromethane 10:90 for 10 min, centrifuge at 150 g for 5 min, remove the upper aqueous layer and evaporate the organic layer under nitrogen at 40°. Dissolve the residue in 200 µL mobile phase, inject a 20 µL aliquot. (Prepare the pH 9.5 ammonia buffer by adding ammonia to a saturated ammonium chloride solution.)

HPLC VARIABLES

Guard column: 10 × 2 40 µm C18

Column: 100 × 3 5 µm Hypersil 5 ODS (Chrompack)

Mobile phase: THF:water 1:100

Flow rate: 1

Injection volume: 20

Detector: UV 275

CHROMATOGRAM

Retention time: 6.5

Internal standard: β-hydroxyethyltheophylline (4.5)

OTHER SUBSTANCES

Extracted: theophylline, theobromine, paraxanthine

KEY WORDS

horse; human; urine

REFERENCE

Delbeke,F.T.; De Backer,P. Threshold level for theophylline in doping analysis, *J.Chromatogr.B*, **1996**, *687*, 247–252.