

Fentanyl


Molecular formula: $C_{22}H_{28}N_2O$

Molecular weight: 336.48

CAS Registry No.: 437-38-7, 990-73-8 (citrate)

Merck Index: 4043

Lednicer No.: 1 299



SAMPLE

Matrix: blood

Sample preparation: 100 μ L Plasma + 50 μ L 5 M NaOH + 100 μ L MeCN + 600 μ L hexane, mix, centrifuge at 2000 rpm for 5 min. Evaporate the organic phase under nitrogen at 30° for about 10 min. Reconstitute the residue in 100 μ L MeCN:50 mM pH 3.0 phosphate buffer 50:50, inject a 50 μ L aliquot. (Use silanized glassware.)

HPLC VARIABLES

Column: 100 \times 8 4 μ m cyano column (Waters)

Mobile phase: MeCN:50 mM pH 3.2 phosphate buffer 50:50

Flow rate: 2.5

Injection volume: 50

Detector: UV 210

CHROMATOGRAM

Retention time: 6.2

Limit of detection: 150 pg

Limit of quantitation: 2.5 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

Noninterfering: ampicillin, calcium chloride, clafloxan, diazepam, dobutamine, dopamine, furosemide, gentamicin, morphine, midazolam, pavulon, phenytoin, vitamin K

KEY WORDS

plasma

REFERENCE

Bansal,R.; Aranda,J.V. High-performance liquid chromatography microassay for the simultaneous determination of fentanyl and its major metabolites in biological samples, *J.Liq.Chromatogr.Rel.Technol.*, **1996**, *19*, 353–364.

SAMPLE

Matrix: blood

Sample preparation: Condition a 1 mL BondElut C18 SPE cartridge once with 1 M HCl, twice with MeOH, and once with water, remove the liquid completely with suction each time. Add 250 μ L IS solution and 250 μ L serum to the column at 1 mL/min, wash twice with water and once with MeCN draining the column completely after each wash, elute with 250 μ L eluting solution, centrifuge for 20 s to remove last of eluate, inject a 5 μ L aliquot of the eluate. (Prepare IS solution by adding 40 μ L 1 mg/mL N-pentyl-2,6-pipecoloxylidide (1-pentyl-N-(2,6-dimethylphenyl)-2-piperidinecarboxamide, pentyl-PPX) in MeOH to 10 mL 100 mM NaH_2PO_4 . Eluting solution was 2.5 mL 35% perchloric acid in 100 mL MeOH.)

HPLC VARIABLES

Guard column: 15 \times 3.2 7 μ m RP-8 (Applied Biosystems)

Column: 150 \times 4.6 5 μ m Ultrasphere octyl

Mobile phase: MeCN:10 mM KH_2PO_4 25:80, pH 5.2

Flow rate: 1.5

Injection volume: 5

Detector: UV 205

CHROMATOGRAM**Retention time:** 11.4**Internal standard:** N-pentyl-2,6-pipecoloxylidide (1-pentyl-N-(2,6-dimethylphenyl)-2-piperidinecarboxamide, pentyl-PPX) (14.5)

OTHER SUBSTANCES**Extracted:** bupivacaine, mepivacaine, meperidine**Noninterfering:** acetaminophen, codeine, epinephrine, morphine, diazepam

KEY WORDS

serum; SPE

REFERENCEGupta,R.N.; Dauphin,A. Column liquid chromatographic determination of bupivacaine in human serum using solid-phase extraction, *J.Chromatogr.B*, **1994**, 658, 113–119.

SAMPLE**Matrix:** blood**Sample preparation:** Automated SPE by ASPEC system. Condition a C18 Clean-Up SPE cartridge (CEC 18111, Worldwide Monitoring) with 2 mL MeOH then 2 mL water. 1 mL Plasma + 1 mL 400 ng/mL protriptyline in water, vortex, add to column, wash with 3 mL water, wash with 3 mL 750 mL/L methanol. Elute with three aliquots of 300 μ L 0.1 M ammonium acetate in MeOH. Add 0.5 mL 0.5 M NaOH and 4 mL 50 mL/L isopropanol in heptane to eluate, mix thoroughly. Allow 5 min for phase separation. Remove upper heptane phase and add it to 300 μ L 0.1 M phosphoric acid (pH 2.5), mix, separate, inject a 100 μ L aliquot of the aqueous phase.

HPLC VARIABLES**Guard column:** LC-8-DB (Supelco)**Column:** 150 \times 4.6 LC-8-DB (Supelco)**Mobile phase:** MeCN:buffer 35:65 (Buffer was 10 mL/L triethylamine in water adjusted to pH 5.5 with glacial acetic acid.)**Flow rate:** 2**Injection volume:** 100**Detector:** UV 228

CHROMATOGRAM**Retention time:** 2.8**Internal standard:** protriptyline (4)

OTHER SUBSTANCES**Extracted:** acetazolamide, amitriptyline, chlordiazepoxide, chlorimipramine, chlorpromazine, desipramine, dextromethorphan, diazepam, doxepin, encainide, fluoxetine, flurazepam, hydroxyethylflurazepam, ibuprofen, imipramine, lidocaine, maprotiline, methadone, methaqualone, mexiletine, midazolam, norchlorimipramine, nordiazepam, norfluoxetine, nortriptyline, norverapamil, pentazocine, promazine, propafenone, propoxyphene, propranolol, protriptyline, quinidine, temazepam, trimipramine, verapamil**Noninterfering:** acetaminophen, acetylmorphine, amiodarone, amobarbital, amphetamine, benodroflumethiazide, benzocaine, benzoylecgonine, benzthiazide, butalbital, carbamazepine, chlorothiazide, clonazepam, cocaine, codeine, cotinine, cyclosporine, cyclothiazide, desalkylflurazepam, diamorphine, dicumerol, ephedrine, ethacrynic acid, ethanol, ethchlorvynol, ethosuximide, furosemide, glutethimide, hydrochlorothiazide, hydrocodone, hydroflumethiazide, hydromorphone, lorazepam, mephentermine, meprobamate, methamphetamine, metharbital, methoxsalen, methoxyphenteramine, methsuximide, methylcyclothiazide, metoprolol, MHPG, monoacetylmorphine, morphine, normethsuximide, oxazepam, oxycodone, oxymorphone, pentobarbital, phencyclidine, phenteramine, phenylephrine, phenytoin, polythiazide, primidone, prochlorperazine, salicylic acid, sulfanilamide, THC-COOH, theophylline, thiazolam, thiopental, thioridazine, tocainide, trichloromethiazide, trifluoperazine, valproic acid, warfarin**Interfering:** diphenhydramine, flecainide, nordoxepin, haloperidol (reduced), trazodone

KEY WORDS

plasma; SPE

REFERENCE

Nichols, J.H.; Charlson, J.R.; Lawson, G.M. Automated HPLC assay of fluoxetine and norfluoxetine in serum, *Clin. Chem.*, **1994**, *40*, 1312–1316.

SAMPLE

Matrix: blood, urine

Sample preparation: 50 μ L Plasma or urine + 50 μ L 4 M NaOH + 100 μ L MeCN + 500 μ L n-hexane, vortex for 30 s, centrifuge at 2000 rpm for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 30°, reconstitute the residue in 100 μ L mobile phase, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 8 4 μ m Nova pak cyano

Mobile phase: MeCN:5 mM pH 3.2 phosphate buffer 70:30

Flow rate: 2.5

Injection volume: 50

Detector: UV 214

CHROMATOGRAM

Retention time: 6.5

Limit of detection: 2.5 ng/mL

OTHER SUBSTANCES

Extracted: alfentanil, sufentanil

KEY WORDS

plasma

REFERENCE

Bansal, R.; Aranda, J.V. Simultaneous microassay of alfentanil, fentanyl, and sufentanil by high performance liquid chromatography, *J. Liq. Chromatogr.*, **1995**, *18*, 339–348.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10–30

Detector: UV 255.8

CHROMATOGRAM

Retention time: 14.202

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

SAMPLE

Matrix: formulations

Sample preparation: Inject a 20 μL aliquot

HPLC VARIABLES

Column: 100×3.9 4 μm Radial Pak phenyl (Waters)

Mobile phase: MeOH:buffer 65:35 (Buffer was 5 mM pH 4.8 phosphate buffer containing 1.4 mM tetrabutylammonium hydroxide.)

Flow rate: 3

Injection volume: 20

Detector: UV 210

CHROMATOGRAM

Retention time: 11.6

OTHER SUBSTANCES

Simultaneous: degradation products, bupivacaine

KEY WORDS

injections; saline; stability-indicating

REFERENCE

Tu,Y.H.; Stiles,M.L.; Allen,L.V.,Jr. Stability of fentanyl citrate and bupivacaine hydrochloride in portable pump reservoirs, *Am.J.Hosp.Pharm.*, **1990**, 47, 2037–2040.

SAMPLE

Matrix: formulations

Sample preparation: 100 μL Injection solution + 400 μL 2.5 $\mu\text{g/mL}$ haloperidol in water, inject a 100 μL aliquot.

HPLC VARIABLES

Column: 100×4.6 5 μm Brownlee C18

Mobile phase: MeCN:MeOH:10 mM NaH_2PO_4 24:31:45, pH adjusted to 5.0 with 2 M KOH

Flow rate: 1.7

Injection volume: 100

Detector: UV 210

CHROMATOGRAM

Retention time: 7.88

Internal standard: haloperidol (9.74)

OTHER SUBSTANCES

Extracted: sufentanil

KEY WORDS

injections

REFERENCE

Dewell,W.M.,Jr.; Khandaghabadi,M.; D'Souza,M.J.; Solomon,H.M. High-performance liquid chromatographic determination of fentanyl and sufentanil returned from the operating room, *Am.J.Hosp.Pharm.*, **1993**, 50, 2374–2375.

SAMPLE

Matrix: formulations

Sample preparation: Dilute 5-fold with mobile phase, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Spherisorb phenyl

Mobile phase: MeOH:buffer 65:35 (Buffer was 5 mM pH 4.8 KH_2PO_4 containing 1.4 mM tetrabutylammonium hydroxide.)

Flow rate: 1.5

Injection volume: 10

Detector: UV 229

CHROMATOGRAM

Retention time: 5.1

OTHER SUBSTANCES

Noninterfering: midazolam

KEY WORDS

injections; 5% dextrose; stability-indicating

REFERENCE

Bhatt-Mehta,V.; Johnson,C.E.; Leininger,N.; Agarwal,M. Stability of fentanyl citrate and midazolam hydrochloride during simulated intravenous coadministration, *Am.J.Health-Syst.Pharm.*, **1995**, 52, 511–513.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 500 μ L Microsomal incubation + 500 μ L DMSO, vortex, centrifuge, filter, inject an aliquot of the filtrate.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Spherisorb ODS-2

Mobile phase: Gradient. A was 1 M ammonium acetate. B was MeCN:MeOH:THF:1 M ammonium acetate 30:20:40:10. A:B from 100:0 to 35:65 over 40 min.

Flow rate: 1

Detector: UV 230

CHROMATOGRAM

Retention time: 40

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

human; liver; pharmacokinetics

REFERENCE

Tateishi,T.; Wood,A.J.J.; Guengerich,F.P.; Wood,M. Biotransformation of tritiated fentanyl in human liver microsomes. Monitoring metabolism using phenylacetic acid and 2-phenylethanol, *Biochem.Pharmacol.*, **1995**, 50, 1921–1924.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 150 \times 3.9 μ Bondapak C18

Mobile phase: MeCN:water adjusted to pH 3 with phosphoric acid 70:30

Flow rate: 1

Injection volume: 10

Detector: UV 210

CHROMATOGRAM

Internal standard: codeine

REFERENCE

Vanbever,R.; Le Boulengé,E.; Préat,V. Transdermal delivery of fentanyl by electroporation. I. Influence of electrical factors, *Pharm.Res.*, **1996**, *13*, 559–565.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in MeOH at a concentration of 1 mg/mL, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 5 Spherisorb S5W

Mobile phase: MeOH:buffer 90:10 (Buffer was 94 mL 35% ammonia and 21.5 mL 70% nitric acid in 884 mL water, adjust the pH to 10.1 with ammonia.)

Flow rate: 2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 1.46

OTHER SUBSTANCES

Simultaneous: tranlylcypromine, caffeine, fenethyline, phendimetrazine, methylphenidate, phenelzine, epinephrine, pipradol, phenylpropanolamine, fencamfamin, chlorphentermine, norpseudoephedrine, phentermine, fenfluramine, methylenedioxymphetamine, amphetamine, normetanephine, 4-hydroxyamphetamine, bromo-STP, STP, prolintane, 2-phenethylamine, tyramine, trimethoxyamphetamine, phenylephrine, pseudoephedrine, ephedrine, methylephedrine, dimethylamphetamine, methamphetamine, mescaline, mephentermine, nalorphine, phenazocine, norpipanone, levallorphan, hydroxypethidine, normethadone, meperidine, dipipanone, diamorphine, pentazocine, acetylcodeine, monoacetylmorphine, thebacon, oxycodone, thebaine, norlevorphanol, methadone, benzylmorphine, ethylmorphine, morphine-N-oxide, codeine, codeine-N-oxide, morphine, etioheptazine, morphine-3-glucuronide, pholcodine, norpethidine, hydrocodone, dihydrocodeine, dihydromorphine, levorphanol, norcodeine, normorphine

Noninterfering: dopamine, levodopa, methylodpa, methylodopate, norepinephrine

Interfering: pemoline, benzphetamine, diethylpropion, mazindol, buprenorphine, dextromoramide, phenoperidine, etorphine, piritramide, noscapine, papaverine, naloxone, dextropropoxyphene

REFERENCE

Law,B.; Gill,R.; Moffat,A.C. High-performance liquid chromatography retention data for 84 basic drugs of forensic interest on a silica column using an aqueous methanol eluent, *J.Chromatogr.*, **1984**, *301*, 165–172.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 μ g/mL solution in MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 125 \times 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 1.6

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine,

chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpentine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipranone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserin, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylegonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozide, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilcaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolantane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, *323*, 191–225.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in mobile phase.

HPLC VARIABLES

Guard column: 15 × 3.2 7 µm Applied Biosystems pre-column

Column: 100 × 2 10 µm µPorasil

Mobile phase: MeCN:5 mM pH 3.75 sodium acetate 80:20

Flow rate: 1

Injection volume: 200

Detector: UV 214

CHROMATOGRAM

Retention time: 9.39

Limit of detection: 5.5 ng/mL

OTHER SUBSTANCES

Simultaneous: buprenorphine, nalbuphine, ethylmorphine, morphine, codeine, meperidine, tramadol

Noninterfering: thiopentone, succinylcholine, pancuronium, diazepam, atropine, neostigmine

Interfering: butorphanol

REFERENCE

Ho, S.-T.; Wang, J.-J.; Ho, W.; Hu, O. Y.-P. Determination of buprenorphine by high-performance liquid chromatography with fluorescence detection: application to human and rabbit pharmacokinetic studies, *J. Chromatogr.*, **1991**, *570*, 339–350.

SAMPLE**Matrix:** solutions**HPLC VARIABLES****Column:** 250 × 4.6 Zorbax RX**Mobile phase:** Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.**Column temperature:** 30**Flow rate:** 2**Detector:** UV 210**OTHER SUBSTANCES**

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitrityline, amobarbital, amoxapine, amphetamine, amyllocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenoprofen, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaacol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarbostyryl, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebedazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methypyrrol, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, nor-epinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphen-butazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantane, phenacetin, phenazocine, phenazopyridine, phenacyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenyl-butazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopola-mine, scopoletin, secobarbital, strychnine, sulfacetamide, sufadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sul-fasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tol-metin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, tripolidine, tropacocaine, tyramine, verapa-mil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, **1994**, *18*, 233–242.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4 5 µm LiChrospher 100 RP-8 (B)**Mobile phase:** MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)**Flow rate:** 0.6**Injection volume:** 25**Detector:** UV 229**CHROMATOGRAM****Retention time:** 11.40 (A), 6.03 (B)**OTHER SUBSTANCES**

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordi-azepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyridamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenopropfen, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazindol, mefenamic acid, meperidine, mephenytoin, mepivacaine, mesoridazine, metaproterenol, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozone, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfipyrazole, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocinamide, tolbutamide, tometin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone


KEY WORDS

also details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, 692, 103–119.

Fenthion

Molecular formula: C₁₀H₁₅O₃PS₂**Molecular weight:** 278.33**CAS Registry No.:** 55-38-9**Merck Index:** 4044**SAMPLE****Matrix:** blood

Sample preparation: 1.5 mL Serum + 2 mL 200 mM pH 7.0 phosphate buffer, add to an Extrelut No. 3 SPE column, let stand for 10 min, elute with 15 mL n-hexane:diethyl ether 80:20. Evaporate the eluate to dryness under a stream of nitrogen at 40°, reconstitute the residue with 150 μ L MeOH:water 70:30, inject a 100 μ L aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: Gradient. MeOH:water from 70:30 to 90:10.

Flow rate: 1

Injection volume: 200

Detector: MS, Hitachi Model M-2000, APCI non-equilibrium interface, vaporizer 250°, nebulizer 400°, ionization needle electrode current 5 μ A, drift voltage 230 V, vacuum 0.0001 Pa, ion-source slit 500 μ m, collector slit 400 μ m, accelerated electrical potential 4 kV, secondary electronic step-up tube potential 1.3 kV, positive-ion mode

CHROMATOGRAM

Retention time: 10.6

Limit of detection: 20 ng

OTHER SUBSTANCES

Extracted: diazinon, dichlorvos, dimethoate, dimethylvinphos, ediphenphos, IBP, isoxathon, malathion, phenthoate, propaphos, pyridafenthion

KEY WORDS

serum; SPE; m/z 279

REFERENCE

Kawasaki,S.; Ueda,H.; Itoh,H.; Tadano,J. Screening of organophosphorus pesticides using liquid chromatography-atmospheric pressure chemical ionization mass spectrometry, *J.Chromatogr.*, **1992**, 595, 193–202.

SAMPLE

Matrix: bulk

Sample preparation: Dissolve in 1 mL 0.5 M NaOH, heat at 75° for 45 min, cool, add 2 drops methyl isobutyl ketone, add 2 drops 0.1% dansyl chloride in acetone, shake well, heat at 65° for 30 min, cool, acidify with 10% HCl, add 300 μ L benzene (Caution! Benzene is a carcinogen!), extract, inject a 1–10 μ L aliquot.

HPLC VARIABLES

Column: 1000 \times 2.4 Zipax coated with 0.5% β,β' -oxydipropionitrile

Mobile phase: Hexane:MeOH 95:5

Flow rate: 0.78

Injection volume: 1–10

Detector: F ex Turner filter no. 811 em Turner filter no. 817

CHROMATOGRAM

Retention time: 3

KEY WORDS

derivatization; normal phase

REFERENCE

Frei,R.W.; Lawrence,J.F. Fluorogenic labelling in high-speed liquid chromatography, *J.Chromatogr.*, **1973**, 83, 321–330.

SAMPLE

Matrix: solutions

Sample preparation: Condition a 10 \times 2 SPE column packed with 40 μ m octadecylsilica (Spark Holland) with 10 mL MeCN, 10 mL MeOH, and 10 mL water at 2 mL/min. Add nitric acid to a final concentration of 0.5% to water sample, filter (0.45 μ m), add a 150 mL aliquot to the SPE column at 3 mL/min, wash with 3 mL distilled water, elute the contents of the SPE column on to the analytical column with the mobile phase.

HPLC VARIABLES

Column: 250 × 4 μm Superspher 60 RP-8 endcapped C8 (Merck)

Mobile phase: Gradient. A was MeCN:MeOH 80:20. B was water. A:B from 10:90 to 40:60 over 10 min, maintain at 40:60 for 5 min, to 90:10 over 33 min, return to initial conditions over 5 min.

Flow rate: 1

Detector: UV 220

CHROMATOGRAM

Retention time: 39.8

Limit of detection: 65 ng/L

OTHER SUBSTANCES

Simultaneous: azinphos-ethyl, azinphos-methyl, chlorfenvinphos, dichlorvos, fenitrothion, malathion, mevinphos, parathion-ethyl, parathion-methyl

Interfering: diazinon

KEY WORDS

groundwater; wastewater; SPE

REFERENCE

Lacorte,S.; Barceló,D. Improvements in the determination of organophosphorus pesticides in ground- and wastewater samples from interlaboratory studies by automated on-line liquid-solid extraction followed by liquid chromatography-diode array detection, *J.Chromatogr.A*, **1996**, 725, 85–92.

Fentiazac


Molecular formula: C₁₇H₁₂ClNO₂S

Molecular weight: 329.81

CAS Registry No.: 18046-21-4

Merck Index: 4045

Lednicer No.: 4 96

**SAMPLE**

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 μL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 μL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 μm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 247

CHROMATOGRAM

Retention time: 13.49

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoyllecgonine; acetaminophen; diazoxide; dacarbazine; sulfapyrazole; flumazenil; sulpride; morphine; atenolol; tolaxotone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vandesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidamine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254–262.

Fenticonazole


Molecular formula: C₂₄H₂₀Cl₂N₂O₃S

Molecular weight: 455.41

CAS Registry No.: 72479-26-6, 73151-29-8 (nitrate)

Merck Index: 4047

Lednicer No.: 4 93

**SAMPLE**

Matrix: formulations

Sample preparation: Tablets. Powder tablets, weigh out amount equivalent to about 30 mg, add 100 mL MeOH, sonicate for 5 min, filter. Add a 2 mL aliquot of filtrate to 5 mL of 100 µg/mL ketoconazole in MeOH, make up to 25 mL with MeOH, inject 20 µL aliquot. Cream. Condition a 500 mg Bond-Elut diol cartridge with 6 mL dichloromethane. Weigh out cream equivalent to about 5 mg of drug, add 30 mL dichloromethane, sonicate for 3 min, make up to 100 mL with dichloromethane, filter. Add a 2 mL aliquot to the cartridge, wash with 2 mL dichloromethane:methanol 4:1, wash with 2 mL dichloromethane, elute with 3 mL MeOH:buffer 85:15. Add eluate to 0.5 mL 100 µg/mL ketoconazole in MeOH, make up to 5 mL with MeOH,

inject 20 μ L aliquot. (Buffer was 50 mM triethylamine adjusted to pH 7.0 with phosphoric acid.)

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Spherisorb CN

Mobile phase: THF:buffer 30:70 (Buffer was 50 mM triethylamine adjusted to pH 3.0 with phosphoric acid.)

Flow rate: 1

Injection volume: 20

Detector: UV 230

CHROMATOGRAM

Retention time: 32

Internal standard: ketoconazole (7)

OTHER SUBSTANCES

Simultaneous: clotrimazole, ketoconazole, bifonazole, tioconazole, isoconazole, econazole, miconazole

KEY WORDS

tablets; creams

REFERENCE

Di Pietra,A.M.; Cavrini,V.; Andrisano,V.; Gatti,R. HPLC analysis of imidazole antimycotic drugs in pharmaceutical formulations, *J.Pharm.Biomed.Anal.*, **1992**, 10, 873-879.


Feprazone

Molecular formula: C₂₀H₂₀N₂O₂

Molecular weight: 320.39

CAS Registry No.: 30748-29-9

Merck Index: 4053



SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 2 mL Plasma + 250 μ L 1 M HCl + 5 mL diethyl ether, mix, centrifuge at 1100 g for 10 min. Remove the organic phase, repeat the extraction, combine the organic extracts, evaporate under nitrogen at 40°, reconstitute the residue in 200 μ L mobile phase, vortex, inject a 20 μ L aliquot. Urine. 1 mL Urine + 250 μ L 1 M HCl + 5 mL diethyl ether, rotate for 15 min. Remove the organic layer, add 1 mL 1% sodium hydrogen carbonate, vortex for 1 min. Remove the organic layer and evaporate it under nitrogen at 40°, reconstitute the residue in 200 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 3 5 μ m Nucleosil

Mobile phase: MeCN:water:acetic acid 40:59.4:0.6

Flow rate: 0.8

Injection volume: 20

Detector: UV 245

CHROMATOGRAM

Retention time: 8.3

Internal standard: feprazone

OTHER SUBSTANCES

Extracted: tiaprofenic acid (UV 310)

Noninterfering: alclofenac, diclofenac, fenoprofen, flunixin, flurbiprofen, ibuprofen, indomethacin, naproxen, oxyphenbutazone, phenylbutazone, piroxicam

KEY WORDS

plasma; feprazone is IS

REFERENCE

Delbeke, F.T.; Baert, K.; De Backer, P. Disposition of human drug preparations in the horse. VI. Tiaprofenic acid, *J.Chromatogr.B*, **1997**, *704*, 207–214.

Fibrinolysin

Molecular weight: about 90000

CAS Registry No.: 9004-09-5, 9001-90-5

Merck Index: 7678

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 100 × 6 Asahipak GS-520-AHA-ABA (Prepare by suspending 10 g Asahipak-GS gel (Asahi Chemical Industry) in water, sonicate for 5 min, wash with 200 mL water, wash with 200 mL dioxane, suspend in 100 mL dioxane, add 3.24 g 1,1'-carbonyldiimidazole, stir gently for 15 min at room temperature, wash with 200 mL dioxane, suspend in 200 mL 1 M sodium bicarbonate containing 1 M 6-aminohexanoic acid, shake at 4° for 25 h, wash with 200 mL water, wash with 100 mL 1 M NaCl, wash with 200 mL water. Suspend 2 g of the gel in 15 mL 200 mM pH 4.752-(morpholino)ethanesulfonic acid/NaOH buffer, add 288 mg 1-ethyl-3-(3-dimethylaminopropyl)carbodiimide monohydrochloride, stir gently for 30 min, add 28.3 mg p-aminobenzamide monohydrochloride, adjust pH three times to 4.75 with 1 M HCl or 1 M NaOH at 30 min intervals, shake gently at room temperature for 24 h, wash with 150 mL water, wash with 100 mL 50 mM NaOH containing 1 M NaCl, wash with 100 mL 50 mM HCl containing 1 M NaCl, wash with water until washings are neutral. Caution! Dioxane is a carcinogen. It may be possible to use acetone instead.)

Mobile phase: Gradient. Buffer A, after 10 min buffer B, after 25 min buffer B containing 40 mM 6-hexanoic acid, after 40 min buffer B containing 40 mM 6-hexanoic acid and 1 M urea (step gradient). Buffer A was 50 mM pH 6.5 sodium phosphate. Buffer B was 50 mM pH 7.4 sodium phosphate containing 100 mM NaCl.

Flow rate: 1.8

Detector: F ex 285 em 340 (?)

CHROMATOGRAM

Retention time: 50

OTHER SUBSTANCES

Simultaneous: plasminogens

REFERENCE

Ito, N.; Noguchi, K.; Kazama, M.; Shimura, K.; Kasai, K.-I. Separation of human Glu-plasminogen, Lys-plasminogen and plasmin by high-performance affinity chromatography on Asahipak GS gel coupled with p-aminobenzamide, *J.Chromatogr.*, **1985**, *348*, 199–204.

Finasteride


Molecular formula: $C_{23}H_{36}N_2O_2$

Molecular weight: 372.55

CAS Registry No.: 98319-26-7

Merck Index: 4125

Lednicer No.: 5 49



SAMPLE

Matrix: blood

Sample preparation: Condition a 6 mL 300 mg Carbograph SPE cartridge (Alltech) with 2 mL water and 4 mL MeOH. Add 100 μ L 10 μ g/mL IS to 1.0 mL plasma, briefly vortex, add the mixture to the SPE cartridge. Wash with 2 mL water and 4 mL MeOH:water 20:80. Elute with 2 mL MeOH:chloroform 20:80 (Caution! Chloroform is a carcinogen!). Centrifuge eluate at 1000 g for 10 min, filter through a WTP 500 nm filter, evaporate to dryness with a stream of nitrogen under vacuum. Reconstitute the residue in 1.0 mL mobile phase, vortex, inject a 50 μ L aliquot.

HPLC VARIABLES

Guard column: 20 \times 4.6 40 μ m Pelliguard (Supelco)

Column: 250 \times 4.6 5 μ m reversed-phase Symmetry (Waters)

Mobile phase: MeCN:40 mM pH 4.0 orthophosphoric acid 45:55

Flow rate: 1.2

Injection volume: 50

Detector: UV 215

CHROMATOGRAM

Retention time: 14.6

Internal standard: 4-androstene-3,17-dione (9.8)

Limit of detection: 5 ng/mL

Limit of quantitation: 10 ng/mL

KEY WORDS

plasma; SPE

REFERENCE

Carlucci, G.; Mazzeo, P. Finasteride in biological fluids: extraction and separation by a graphitized carbon black cartridge and quantification by high-performance liquid chromatography, *J. Chromatogr. B*, **1997**, 693, 245-248.

SAMPLE

Matrix: blood

Sample preparation: Non-buffered extraction. 1 mL Plasma + 100 μ L 100 ng/mL IS in MeOH. Add 7 mL MTBE, mix and rotate for 15 min, centrifuge, remove the organic layer and evaporate it to dryness under a stream of nitrogen at 50°. Reconstitute the residue in 300 μ L mobile phase, inject an aliquot. Buffered extraction. 1 mL Plasma + 100 μ L 100 ng/mL IS in MeOH. Add 1 mL 200 mM pH 9.8 carbonate buffer and 7 mL MTBE. Extract as described above. Inject an aliquot.

HPLC VARIABLES

Column: 30 \times 4.6 3 μ m BDS Hypersil C18 (A), 50 \times 2 3 μ m BDS Hypersil C18 (B)

Mobile phase: MeCN:water containing 0.1% formic acid 90:10

Column temperature: 60

Flow rate: 1(A), 0.2 (B)

Detector: MS, PE Sciex API IIIplus tandem MS, heated nebulizer 500°, corona discharge needle +4 μ A, APCI positive ion, nebulizing gas air 550 kPa, auxiliary gas at 2 mL/min, curtain gas nitrogen at 0.9 mL/min, orifice +50 V, interface heater 60°, collision gas argon, m/z 373 (A), Turbo ion spray, nebulizing gas air 410 kPa, auxiliary gas at 6 mL/min, curtain gas at 1 mL/min, orifice +30 V, interface heater 60°, m/z 373 (B)

tyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pe-
 cazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone,
 phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide,
 phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine,
 phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimino-
 dine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pi-
 renzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, pri-
 maquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine,
 prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prop-
 thipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, qui-
 nine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldi-
 amine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine,
 thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranlylcypromine, tra-
 zodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, tri-
 methoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, *323*, 191-225.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4 5 µm LiChrospher 100 RP-8 (B)

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 11.30 (A), 5.94 (B)

KEY WORDS


also details of plasma extraction; also acebutolol; acepromazine; acetaminophen; acetazolamide; acetophenazine; albuterol; alprazolam; amitriptyline; amobarbital; amoxapine; antipyrine; atenolol; atropine; azatadine; baclofen; benzocaine; bromocriptine; brompheniramine; brotizolam; bupivacaine; buspirone; butabarbital; butalbital; caffeine; carbamazepine; cetirizine; chlorcyclizine; chlordiazepoxide; chlormezanone; chloroquine; chlorpheniramine; chlorpromazine; chlorpropamide; chlorprothixene; chlorthalidone; chlorzoxazone; cimetidine; cisapride; clomipramine; clonazepam; clonidine; clozapine; cocaine; codeine; colchicine; cyclizine; cyclobenzaprine; dantrolene; desipramine; diazepam; diclofenac; diflunisal; diltiazem; diphenhydramine; diphenidol; diphenoxylate; dipyrindamole; disopyramide; dobutamine; doxapram; doxepin; droperidol; encainide; ethidium bromide; ethopropazine; fenoprofen; fentanyl; fluoxetine; fluphenazine; flurazepam; flurbiprofen; fluvoxamine; furosemide; glutethimide; glyburide; guaifenesin; haloperidol; homatropine; hydralazine; hydrochlorothiazide; hydrocodone; hydromorphone; hydroxychloroquine; hydroxyzine; ibuprofen; imipramine; indomethacin; ketoconazole; ketoprofen; ketorolac; labetalol; levorphanol; lidocaine; loratadine; lorazepam; lovastatin; loxapine; mazindol; mefenamic acid; meperidine; mephénytoin; mepivacaine; mesoridazine; metaproterenol; methadone; methdilazine; methocarbamol; methotrexate; methotrimiprazine; methoxamine; methyl dopa; methylphenidate; metoclopramide; metolazone; metoprolol; metronidazole; midazolam; moclobemide; morphine; nadolol; nalbuphine; naloxone; naphazoline; naproxen; nifedipine; nizatidine; norepinephrine; nortriptyline; oxazepam; oxycodone; oxymetazoline; paroxetine; pemoline; pentazocine; pentobarbital; pentoxifylline; perphenazine; pheniramine; phenobarbital; phenol; phenolphthalein; phentolamine; phenylbutazone; phenyltoloxamine; phenytin; pimozone; pindolol; piroxicam; pramoxine; prazepam; prazosin; probenecid; procainamide; procaine; prochlorperazine; procyclidine; promazine; promethazine; propafenone; propantheline; propiomazine; propofol; propranolol; protriptyline; quazepam; quinidine; quinine; racemethorphan; ranitidine; remoxipride; risperidone; salicylic acid; scopolamine; secobarbital; sertraline; sotalol; spironolactone; sulfapyrazone; sulindac; temazepam; terbutaline; terfenadine

dine; tetracaine; theophylline; thiethylperazine; thiopental; thioridazine; thiothixene; timolol; tocinide; tolbutamide; tolmetin; trazodone; triamterene; triazolam; trifluoperazine; triflupromazine; trimeprazine; trimethoprim; trimipramine; verapamil; warfarin; xylometazoline; yohimbine; zopiclone

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, 692, 103–119.

Flecainide



Molecular formula: C₁₇H₂₀F₆N₂O₃

Molecular weight: 414.35

CAS Registry No.: 54143-55-4, 54143-56-5 (acetate)

Merck Index: 4136

Lednicer No.: 3 59

SAMPLE

Matrix: blood

Sample preparation: Automated SPE by ASPEC system. Condition a C18 Clean-Up SPE cartridge (CEC 18111, Worldwide Monitoring) with 2 mL MeOH then 2 mL water. 1 mL Plasma + 1 mL 400 ng/mL protriptyline in water, vortex, add to column, wash with 3 mL water, wash with 3 mL 750 mL/L methanol. Elute with three aliquots of 300 µL 0.1 M ammonium acetate in MeOH. Add 0.5 mL 0.5 M NaOH and 4 mL 50 mL/L isopropanol in heptane to eluate, mix thoroughly. Allow 5 min for phase separation. Remove upper heptane phase and add it to 300 µL 0.1 M phosphoric acid (pH 2.5), mix, separate, inject a 100 µL aliquot of the aqueous phase.

HPLC VARIABLES

Guard column: LC-8-DB (Supelco)

Column: 150 × 4.6 LC-8-DB (Supelco)

Mobile phase: MeCN:buffer 35:65 (Buffer was 10 mL/L triethylamine in water adjusted to pH 5.5 with glacial acetic acid.)

Flow rate: 2

Injection volume: 100

Detector: UV 228

CHROMATOGRAM

Retention time: 3.1

Internal standard: protriptyline (4)

OTHER SUBSTANCES

Extracted: acetazolamide, amitriptyline, chlordiazepoxide, chlorimipramine, chlorpromazine, desipramine, dextromethorphan, diazepam, encainide, fluoxetine, flurazepam, hydroxyethylflurazepam, ibuprofen, imipramine, lidocaine, maprotiline, methadone, methaqualone, mexiletine, midazolam, norchlorimipramine, nordiazepam, norfluoxetine, nortriptyline, norverapamil, pentazocine, promazine, propafenone, propoxyphene, propranolol, protriptyline, quinidine, temazepam, trimipramine, verapamil

Noninterfering: acetaminophen, acetylmorphine, amiodarone, amobarbital, amphetamine, ben-droflumethiazide, benzocaine, benzoylcegonine, benzthiazide, butalbital, carbamazepine, chlorothiazide, clonazepam, cocaine, codeine, cotinine, cyclosporine, cyclothiazide, desalkylflurazepam, diamorphine, dicumerol, ephedrine, ethacrynic acid, ethanol, ethchlorvynol, ethosuximide, furosemide, glutethimide, hydrochlorothiazide, hydrocodone, hydroflumethiazide, hydromorphone, lorazepam, mephentermine, meprobamate, methamphetamine, metharbital, methoxsalen, methoxyphenteramine, methsuximide, methylcyclothiazide, metoprolol, MHPG, monoacetylmorphine, morphine, normethsuximide, oxazepam, oxycodone, oxymorphone, pentobarbital, phenacyclidine, phenteramine, phenylephrine, phenytoin, polythiazide, primidone, prochlorperazine, salicylic acid, sulfanilamide, THC-COOH, theophylline, thiazolam, thiopental, thioridazine, tocinide, trichloromethiazide, trifluoperazine, valproic acid, warfarin

Interfering: diphenhydramine, doxepin, fentanyl, haloperidol, nordoxepin, trazodone

KEY WORDS

plasma; SPE

REFERENCE

Nichols, J.H.; Charlson, J.R.; Lawson, G.M. Automated HPLC assay of fluoxetine and norfluoxetine in serum, *Clin. Chem.*, **1994**, *40*, 1312–1316.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 299

CHROMATOGRAM

Retention time: 5.22

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoyllecgonine; acetaminophen; diazoxide; dacarbazine; sulfipyrazole; flumazenil; sulpride; morphine; atenolol; tolaxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vindesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loperazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclicine; histapyrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvox-

amine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpi-pramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254-262.

SAMPLE

Matrix: blood, gastric contents, tissue, urine

Sample preparation: Homogenize 15 g liver with water to a final volume of 90 mL, dilute 1 to 10 with water. Homogenize gastric contents, dilute 1 to 100 with water. Dilute urine 1 to 25. 1 mL Whole blood, diluted liver homogenate, diluted gastric contents, or diluted urine + 50 μ L 100 μ g/mL clomipramine, mix, inject an aliquot.

HPLC VARIABLES

Guard column: 20 mm long Hypersil ODS

Column: 150 mm long Hypersil ODS

Mobile phase: MeCN:10 mM Na₂HPO₄:n-nonylamine 40:60:0.12, pH 3

Flow rate: 1

Detector: UV 295

CHROMATOGRAM

Retention time: 3

Internal standard: clomipramine (5.5)

Limit of quantitation: 1 μ g/mL

KEY WORDS

liver; whole blood

REFERENCE

Sadler,D.W.; Quigley,C. Unsuspected self-poisoning with flecainide and alcohol, *J.Forensic Sci.*, **1995**, *40*, 903-905.

SAMPLE

Matrix: blood, urine

Sample preparation: Dilute urine 1:100. 1 mL Plasma or diluted urine + 100 μ L 4 μ g/mL IS in water + 1 mL water + 1 mL 1 M NaOH + 8 mL distilled diethyl ether, shake for 15 min, centrifuge at 1000 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 1 mL freshly prepared 49 mM 1-[(4-nitrophenyl)sulfonyl]-L-prolyl chloride in ethyl acetate, add 100 μ L 0.016% triethylamine in ethyl acetate, vortex for 5 s, heat at 80° for 2 h, cool, wash with 3 mL 600 mM HCl, centrifuge for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 200 μ L mobile phase, inject a 40-180 μ L aliquot. (Synthesis of 1-[(4-nitrophenyl)sulfonyl]-L-prolyl chloride is as follows. Stir 40-45 mmoles L-(-)-proline in 40 mL THF and 200 mL 10% potassium carbonate, add 37-43 mmoles 4-nitrophenylsulfonyl chloride in 40 mL THF dropwise, heat at 50° for 3 h and maintain at pH 8 or above, cool, acidify to pH 2 and extract into chloroform. Extract with 10% potassium carbonate, acidify the aqueous layer, extract into chloroform. Dry the chloroform extract, evaporate, recrystallize from petroleum ether and benzene (Caution! Benzene is a carcinogen!). Stir 15 mmoles of the 4-nitrophenylsulfonylproline in 100 mL benzene under reflux condenser fitted with a calcium sulfate drying tube, add dropwise a five-fold molar excess of thionyl chloride in 50 mL benzene, heat at 35-40° until formation of the acid chloride is complete (about 48 h, monitor by IR spectroscopy), evaporate, recrystallize product from HPLC-grade heptane (mp 110-110.5°).)

HPLC VARIABLES

Guard column: 50 mm long octadecyl Pellicular ODS (Whatman)

Column: 300 \times 3.9 Bondapak C18

Mobile phase: MeCN:water:triethylamine 45:55:0.2

Flow rate: 1

Injection volume: 40-180

Detector: UV 280

CHROMATOGRAM

Retention time: 29.8 (R), 31.9 (S)

Internal standard: (R,S)-N-(2-piperidylmethyl)-2,3-bis-(2,2,2-trifluoroethoxy)benzamide (R, 25.5, S, 27.9)

Limit of quantitation: 50 ng/mL

KEY WORDS

plasma; derivatization; pharmacokinetics; chiral

REFERENCE

Alessi-Severini,S.; Jamali,F.; Pasutto,F.M.; Coutts,R.T.; Gulamhusein,S. High-performance liquid chromatographic determination of the enantiomers of flecainide in human plasma and urine, *J.Pharm.Sci.*, **1990**, 79, 257-260.

SAMPLE

Matrix: formulations

Sample preparation: Extract ground tablets containing 10-50 mg of the compound with 100 mL MeOH, filter. Add 500 μ L 100 μ g/mL IS in MeOH to 1 mL filtrate, make up to 10 mL with MeOH. Inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4 10 μ m LiChrosorb

Mobile phase: MeCN:MeOH:buffer 30:60:10 (Buffer was 67 mM KH_2PO_4 adjusted to pH 2.9 with phosphoric acid.)

Flow rate: 1

Injection volume: 50

Detector: UV 254

CHROMATOGRAM

Retention time: 3.30 (flecainide acetate)

Internal standard: amiodarone (7.46)

Limit of quantitation: 4 μ g/mL

KEY WORDS

tablets

REFERENCE

Paw,B.; Przyborowski,L.; Slawik,T. Determination of flecainide acetate in tablets by HPLC and UV-spectrophotometry, *Pharmazie*, **1998**, 53, 97-98.

SAMPLE

Matrix: solutions

Sample preparation: Mix 1 mL of an aqueous solution with 1 mL 100 mM nickel sulfate in water, 1 mL 20% aqueous ammonia, and 5 mL chloroform:carbon disulfide 98:2, shake vigorously for 1 min, wash the organic layer with three 2 mL portions of water, filter (phase-separation paper). Evaporate the filtrate to dryness under a stream of nitrogen, reconstitute with 1 mL mobile phase, inject a 10 μ L aliquot. (Copper may also be used with electrochemical detection or UV detection at 270 nm.)

HPLC VARIABLES

Guard column: 30 \times 4 40 μ m LiChrosorb RP-18

Column: 250 \times 4 7 μ m LiChrosorb RP-18

Mobile phase: MeOH:20 mM pH 5.8 sodium acetate buffer 80:20 containing 5 mM lithium perchlorate

Flow rate: 1.5

Injection volume: 10

Detector: UV 325, E, Merck-Clevenot E 230, Model LCC 231 thin-layer electrolytic cell with a glassy carbon electrode at +0.7 V, standard calomel reference electrode

CHROMATOGRAM**Retention time:** k' 2.93**Limit of detection:** 1 fmole (E), 1 nmole (UV)

OTHER SUBSTANCES**Also analyzed:** acebutolol, alprenolol, ephedrine, methamphetamine, propranolol

KEY WORDS

derivatization; complexation

REFERENCE

Leroy,P.; Nicolas,A. Determination of secondary amino drugs as their metal dithiocarbamate complexes by reversed-phase high-performance liquid chromatography with electrochemical detection, *J.Chromatogr.*, 1984, 317, 513-521.

SAMPLE**Matrix:** solutions

Sample preparation: Mix a 50 μ L aliquot of a solution in MeOH:triethylamine 99:1 with 20 μ L 0.1% FLOPIC in dry toluene, vortex briefly, let stand at room temperature in the dark for 30 min, add 50 μ L 1% ethanolamine in MeOH, let stand at room temperature for 15 min, evaporate to dryness under reduced pressure, reconstitute with 100 μ L mobile phase, sonicate for 30 s, inject a 20 μ L aliquot. (FLOPIC is (-)-(S)-flunoxaprofen isocyanate; synthesis is as follows. Dissolve 1 g (+)-(S)-flunoxaprofen in 30 mL acetone, cool to 0°, add a solution of 500 μ L triethylamine in 2 mL acetone dropwise, add a solution of 370 μ L ethyl chloroformate in 2 mL acetone dropwise, stir at 0° for 15 min, add a solution of 250 mg sodium azide in 1 mL water dropwise (Caution! Sodium azide is highly toxic!), stir for 1 h, pour into 60 mL ice water, stir for 10 min, filter, wash the solid with two 50 mL aliquots of ice-water, dry under reduced pressure to obtain flunoxaprofen azide. Dissolve 100 mg flunoxaprofen azide in 3 mL dry toluene, reflux for 10-15 min, cool to room temperature, filter. Evaporate the filtrate to dryness under reduced pressure and dry under reduced pressure to obtain FLOPIC as a crystalline solid (mp 93-94°), store in a desiccator under reduced pressure.)

HPLC VARIABLES**Column:** 150 \times 3.9 4 μ m Nova Pak C18 (A) or 200 \times 4.6 5 μ m Nucleosil cyano (B)**Mobile phase:** MeOH:water:THF 62:35:3 (A) or n-hexane:isopropanol:diethylamine 95:5:0.05 (B)**Flow rate:** 1**Injection volume:** 20**Detector:** F ex 296 em 356

CHROMATOGRAM**Retention time:** 29.0 (R (A)), 31.0 (S (A)), 22.2 (R (A)), 26.8 (S (B))

OTHER SUBSTANCES**Simultaneous:** mexiletine (system A only)

KEY WORDS

derivatization; chiral

REFERENCE

Martin,E.; Quinke,K.; Spahn,H.; Mutschler,E. (-)-(S)-Flunoxaprofen and (-)-(S)-naproxen isocyanate: two new fluorescent chiral derivatizing agents for an enantiospecific determination of primary and secondary amines, *Chirality*, 1989, 1, 223-234.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 150 \times 4.6 Spherisorb S5SCX**Mobile phase:** MeOH:MeCN:water 40:40:20 containing 25 mM perchloric acid**Flow rate:** 2**Detector:** F ex 215 no emission filter

CHROMATOGRAM**Retention time:** 4.5**Internal standard:** benzimidazole (7)

OTHER SUBSTANCES**Simultaneous:** propranolol

REFERENCE

Croes, K.; McCarthy, P.T.; Flanagan, R.J. HPLC of basic drugs and quaternary ammonium compounds on micro-particulate strong cation-exchange materials using methanolic or aqueous methanol eluents containing an ionic modifier, *J.Chromatogr.A*, **1995**, 693, 289–306.

SAMPLE**Matrix:** solutions

Sample preparation: Mix 20 μL of a 1 mM solution in MeOH or water with 50 μL pH 8 borate buffer and 50 μL 18 mM 2-(6-methoxy-2-naphthyl)-1-propyl chloroformate in acetone, vortex, let stand at room temperature for 30 min, add 100 μL 10 mM trans-4-hydroxy-L-proline in water, mix, let stand for 2 min, add 2 mL dichloromethane, vortex for 30 s. Remove the organic layer and evaporate it to dryness under reduced pressure, reconstitute the residue in 100 μL mobile phase, inject an aliquot. Prepare 2-(6-methoxy-2-naphthyl)-1-propyl chloroformate as follows. Stir 1.5 mmoles lithium aluminum hydride in THF, slowly add 2 mmoles (S)-naproxen in 20 mL anhydrous THF, reflux for 1 h, evaporate most of the solvent, cautiously add water with stirring, acidify with 6 N HCl, extract three times with diethyl ether. Combine the organic layers and dry them over anhydrous sodium sulfate, evaporate to dryness, chromatograph on silica gel with dichloromethane:MeOH 100:2 (flash chromatography), evaporate eluate to dryness, dry under vacuum over KOH to give 2-(6-methoxy-2-naphthyl)propanol as a white solid (mp 92–3°). Stir 0.5 mmoles 2-(6-methoxy-2-naphthyl)propanol and 0.5 mmoles triethylamine in 10 mL dry toluene at 0°, add 1 mL 20% phosgene in toluene (Caution! Phosgene is highly toxic, perform reaction in a chemical fume hood!) (Fluka), stir for 4 h, filter, evaporate to dryness under reduced pressure, dry under vacuum to give 2-(6-methoxy-2-naphthyl)-1-propyl chloroformate (mp 60°). Store under vacuum over phosphorus pentoxide at room temperature.)

HPLC VARIABLES**Column:** 250 \times 4.5 μm Zorbax-SIL**Mobile phase:** n-Hexane:isopropanol 100:1.5**Flow rate:** 1.5**Injection volume:** 100**Detector:** UV 230, F ex 270 em 365

CHROMATOGRAM**Retention time:** k' 14.4 (R-(-)), k' 15.2 (S-(+))

OTHER SUBSTANCES**Simultaneous:** metoprolol, tocainide**Interfering:** propafenone

KEY WORDS

derivatization; chiral; normal phase

REFERENCE

Büschges, R.; Linde, H.; Mutschler, E.; Spahn-Langguth, H. Chloroformates and isothiocyanates derived from 2-arylpropionic acids as chiral reagents: synthetic routes and chromatographic behaviour of the derivatives, *J.Chromatogr.A*, **1996**, 725, 323–334.

Fleroxacin


Molecular formula: $C_{17}H_{16}F_3N_3O_3$

Molecular weight: 369.34

CAS Registry No.: 79660-72-3

Merck Index: 4137

Lednicer No.: 5 125



SAMPLE

Matrix: bile, blood, urine

Sample preparation: Dilute urine 1:20. Dilute bile 1:10. 500 μ L Serum, diluted urine, or diluted bile + 3.2 mL dichloromethane, vortex, rotate at 20 rpm for 10 min, centrifuge at 1000 g for 10 min. Remove 3 mL of the lower organic phase and add it to 200 μ L 100 mM NaOH, rotate at 20 rpm for 30 min, centrifuge at 1000 g for 10 min, inject a 20 μ L aliquot of the aqueous layer.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Ultrasphere C18

Mobile phase: MeCN:buffer 10:90, pH adjusted to 2 with 14.6 M phosphoric acid (Buffer was 10 mM NaH_2PO_4 containing 5 mM tetrabutylammonium bromide.)

Flow rate: 2

Injection volume: 20

Detector: F ex 277 em 445

CHROMATOGRAM

Retention time: 3

Limit of detection: 25 ng/mL (bile), 50 ng/mL (urine), 2.5 ng/mL (serum)

OTHER SUBSTANCES

Noninterfering: amikacin, aztreonam, carbamazepine, cephalosporins, ciprofloxacin, clavulanic acid, difloxacin, digitoxin, digoxin, fosfomycin, furosemide, gentamycin, imipenem, lidocaine, netilmicin, norfloxacin, pefloxacin, penicillins, phenobarbital, phenytoin, primidone, procainamide, quinidine, rifampin, salicylic acid, teicoplanin, temafloxacin, theophylline, tobramycin, valproic acid, vancomycin

Interfering: ofloxacin

KEY WORDS

serum; human; rabbit

REFERENCE

Koechlin,C.; Jehl,F.; Linger,L.; Monteil,H. High-performance liquid chromatography for the determination of three new fluoroquinolones, fleroxacin, temafloxacin and A-64730, in biological fluids, *J.Chromatogr.*, **1989**, *491*, 379-387.

SAMPLE

Matrix: blood

Sample preparation: Filter 1 mL plasma using a micropartition system (MPS-1, Amicon, MA) while centrifuging at 2000 g for 20 min at 10°, inject an aliquot of the ultrafiltrate.

HPLC VARIABLES

Column: 250 \times 4.6 Spherisorb ODS-2 endcapped

Mobile phase: MeCN:buffer 22:78 containing 5 mM tetrabutylammonium sulfate, adjusted to pH 2.5 with 1 M NaOH (Buffer was 100 mM citric acid containing 200 mM ammonium perchlorate.)

Column temperature: 37

Flow rate: 1

Detector: UV 287

CHROMATOGRAM**Retention time:** 4.55**Internal standard:** difloxacin (9.5)

KEY WORDS

plasma; ultrafiltrate

REFERENCE

Zlotos,G.; B cker,A.; Kinzig-Schippers,M.; Sorgel,F.; Holzgrabe,U. Plasma protein binding of gyrase inhibitors, *J.Pharm.Sci.*, **1998**, 87, 215–220.

SAMPLE**Matrix:** blood**Sample preparation:** 500 μ L Serum + 500 μ L 7% perchloric acid, vortex for 10 s, centrifuge at >700 g for 10 min, inject a 20 μ L aliquot of the supernatant.

HPLC VARIABLES**Column:** 150 \times 3.9 Nova-Pak C18**Mobile phase:** MeOH:18 mM KH_2PO_4 containing 0.13 mM heptanesulfonic acid:concentrated phosphoric acid 30:70:0.1**Flow rate:** 0.8**Injection volume:** 20**Detector:** F ex 288 em 475 bandpass filter

CHROMATOGRAM**Retention time:** 5.5**Limit of detection:** 100 ng/mL

OTHER SUBSTANCES**Extracted:** metabolites

KEY WORDS

serum

REFERENCE

Griggs,D.J.; Wise,R. A simple isocratic high-pressure liquid chromatographic assay of quinolones in serum, *J.Antimicrob.Chemother.*, **1989**, 24, 437–445.

SAMPLE**Matrix:** blood**Sample preparation:** Filter (0.2 μ m cellulose acetate, Schleicher & Schuell) while centrifuging at 4° at 2500 g for 5 min, maintain filtrate at 4°, inject a 5 μ L aliquot of the filtrate onto column A and elute to waste with mobile phase A, collect the effluent containing fleroxacin in a sample loop and backflush the contents of this loop onto column B with mobile phase B, elute column B with mobile phase B and monitor the effluent.

HPLC VARIABLES**Column:** A 40 \times 6 12 μ m TSKgel PWxl methacrylate polymer gel (Toso-Haas); B 250 \times 4.6 5 μ m Zorbax Rx-C8**Mobile phase:** A Isopropanol:30 mM pH 6.2 potassium phosphate buffer 10:90; B MeCN:50 mM pH 2.7 potassium phosphate buffer 18:82**Flow rate:** A 0.6; B 1**Injection volume:** 5**Detector:** UV 287

CHROMATOGRAM**Retention time:** 15**Limit of quantitation:** 270 ng/mL

OTHER SUBSTANCES**Extracted:** Ro 23-9424 (a pro-drug)

KEY WORDS

plasma; column-switching; heart-cut; pharmacokinetics

REFERENCE

Szuna,A.J.; Blain,R.W. Determination of a new antibacterial agent (Ro 23-9424) by multidimensional high-performance liquid chromatography with ultraviolet detection and direct plasma injection, *J.Chromatogr.*, **1993**, 620, 211–216.

SAMPLE

Matrix: blood

Sample preparation: 250 μ L Serum + 20 μ L 30 μ g/mL pipemidic acid + 250 μ L 25% sodium sulfate, vortex, add 3.5 mL chloroform, shake at low speed for 10 min, centrifuge at 1155 g for 10 min. Remove the lower organic layer and add it to 200 μ L 1 M NaOH, shake fast for 20 min, centrifuge at 1155 g for 10 min, inject a 150 μ L aliquot of the aqueous layer.

HPLC VARIABLES

Column: 10 μ m Nucleosil C18

Mobile phase: MeCN:MeOH:10 mM phosphate buffer containing 5 mM tetrabutylammonium hydrogen sulfate 13:6:81

Flow rate: 1

Injection volume: 150

Detector: F ex 274

CHROMATOGRAM

Internal standard: pipemidic acid

Limit of detection: 50 ng/mL

KEY WORDS

serum; pharmacokinetics

REFERENCE

Bertino,J.S.,Jr.; Nafziger,A.N.; Wong,M.; Stragand,L.; Puleo,C. Effect of a fat- and calcium-rich breakfast on pharmacokinetics of fleroxacin administered in single and multiple doses, *Antimicrob.Agents Chemother.*, **1994**, 38, 499–503.

SAMPLE

Matrix: blood

Sample preparation: 250 μ L Serum + 20 μ L (?) 30 μ g/mL pipemidic acid + 250 μ L 25% sodium sulfate, vortex, add 3.5 mL chloroform, shake at low speed for 10 min, centrifuge at 1155 g for 10 min. Remove the organic layer and add it to 200 μ L 1 M NaOH, shake quickly for 20 min, centrifuge at 1155 g for 10 min, inject a 150 μ L aliquot of the aqueous layer.

HPLC VARIABLES

Column: 10 μ m Nucleosil C18

Mobile phase: MeCN:MeOH:10 mM phosphate buffer containing 5 mM tetrabutylammonium sulfate 13:6:81

Flow rate: 1

Injection volume: 150

Detector: F ex 274 (?)

CHROMATOGRAM

Internal standard: pipemidic acid

Limit of quantitation: 50 ng/mL

KEY WORDS

pharmacokinetics; serum

REFERENCE

Bertino,J.S.,Jr.; Nafziger,A.N. Pharmacokinetics of oral fleroxacin in male and premenopausal female volunteers, *Antimicrob.Agents Chemother.*, **1996**, 40, 789–791.

SAMPLE

Matrix: blood, dialysate

Sample preparation: Plasma. 100 μ L Plasma + 25 μ L MeCN:70% perchloric acid 80:20 containing IS, mix vigorously, centrifuge. Inject a 5 or 50 μ L aliquot. Dialysate. Dilute 20 μ L dialysate with 980 μ L 50 mM sodium dihydrogen phosphate buffer containing IS. Mix thoroughly. Inject a 5 or 10 μ L aliquot.

HPLC VARIABLES

Column: Spherisorb ODS II

Mobile phase: MeCN containing 2 mM tetrabutyl ammonium hydrogen sulfate:100 mM citric acid buffer containing 5 mM ammonium perchlorate 87:13, adjusted to pH 2.2

Flow rate: 1.2

Injection volume: 5-50

Detector: F ex 290 em 460

CHROMATOGRAM

Retention time: 6.6

Internal standard: pipemidic acid (3.2)

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Uehlinger,G.E.; Schaedeli,F.; Kinzig,M.; Sörgel,F.; Frey,F.J. Pharmacokinetics of fleroxacin after multiple oral dosing in patients receiving regular hemodialysis, *Antimicrob.Agents Chemother.*, **1996**, 40, 1903-1909.

SAMPLE

Matrix: blood, intestinal efflux

Sample preparation: Intestinal efflux. Freeze intestinal efflux at -80°, lyophilize, reconstitute with 1 mL ofloxacin in MeOH:100 mM phosphoric acid 50:50, centrifuge at 3000 rpm for 10 min, inject a 20 μ L aliquot. Serum. Deproteinize serum with MeOH containing ofloxacin, extract with dichloromethane at pH 7.5.

HPLC VARIABLES

Column: 150 \times 3.9 Novapack C18

Mobile phase: MeCN:buffer 16:84 (Buffer was 10 mM pH 3.0 potassium phosphate buffer containing 25 mM sodium heptanesulfonate (PIC B7) and 20 mM triethylamine.)

Flow rate: 1.5

Injection volume: 20

Detector: F ex 290 em 430

CHROMATOGRAM

Internal standard: ofloxacin

Limit of detection: 10 ng/mL

KEY WORDS

serum; rat

REFERENCE

Rubinstein,E.; Dautrey,S.; Farinoti,R.; St.Julien,L.; Ramon,J.; Carbon,C. Intestinal elimination of sparflloxacin, fleroxacin, and ciprofloxacin in rats, *Antimicrob.Agents Chemother.*, **1995**, 39, 99-102.

SAMPLE

Matrix: blood, urine

Sample preparation: 500 μ L Plasma + 100 μ L 1 μ g/mL IS1 in 10 mM HCl + 1 mL pH 7.5 Sørensen buffer + 7 mL dichloromethane:isopropanol 70:30, extract by turning head over head at 70 rotations/min for 10 min, centrifuge at 2000 g for 10 min. Remove 5 mL of the organic

layer and evaporate it to dryness under a stream of nitrogen at 35°, reconstitute the residue in 500 µL mobile phase, vortex, inject an aliquot. Urine. 500 µL Urine + 100 µL 1 mg/mL IS2 in 10 mM HCl + 500 µL 1 M acetic acid + 500 µL 25 mM sodium dodecyl sulfate in water + 7 mL dichloromethane:isopropanol 70:30, extract by turning head over head at 70 rotations/min for 10 min, centrifuge at 2000 g for 10 min. Remove 5 mL of the organic layer and evaporate it to dryness under a stream of nitrogen at 35°, reconstitute the residue in 500 µL mobile phase, vortex, inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm TSK-Gel ODS 120T (Toyo Soda)

Mobile phase: MeOH:5 mM tetrabutylammonium hydrogen sulfate 21:79 (plasma) or 24.5:74.5 (urine)

Column temperature: ambient (urine), 27° (plasma)

Flow rate: 0.8

Detector: F ex 290 em 450

CHROMATOGRAM

Retention time: 7.5 (plasma), 9 (urine)

Internal standard: IS1 6-fluoro-1-(2-fluoroethyl)-1,4-dihydro-7-(4-methyl-1-piperazinyl)-4-oxo-3-quinolinecarboxylic acid (12), IS2 pipemidic acid (7)

Limit of quantitation: 10-20 ng/mL

KEY WORDS

plasma; protect from light; pharmacokinetics

REFERENCE

Dell,D.; Partos,C.; Portmann,R. The determination of a new trifluorinated quinolone, fleroxacin, its N-demethyl, and N-oxide metabolites in plasma and urine by high performance liquid chromatography with fluorescence detection, *J.Liq.Chromatogr.*, **1988**, *11*, 1299–1312.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 200 µL Plasma + 1 mL 50 ng/mL IS1 in MeCN:water 95:5, vortex, centrifuge at 14° at 1000 g for 10 min. Remove the supernatant and evaporate it to dryness under a stream of nitrogen at 38°, reconstitute the residue in 200-400 µL mobile phase, inject a 20-30 µL aliquot onto column A and column B in series, after 1.5 min remove column A from the circuit and backflush it to waste, elute column B with mobile phase and monitor the effluent. Urine. 100 µL Urine + 100 µL 100 µg/mL IS2 in water + 5 mL mobile phase, inject a 20 µL aliquot onto column A and column B in series, after 1.5 min remove column A from the circuit and backflush it to waste, elute column B with mobile phase and monitor the effluent.

HPLC VARIABLES

Column: A 10 × 4 Nucleosil 5-C18 (replaced every 60-80 samples); B 250 × 4.6 5 µm TSK ODS-120T (Toyo Soda)

Mobile phase: MeOH:buffer 28:72, pH adjusted to 2.6 with 40% phosphoric acid (Buffer was 50 mM KH₂PO₄ containing 10 mM tetrabutylammonium hydrogen sulfate.)

Column temperature: 28

Flow rate: 1

Injection volume: 20-30

Detector: F ex 290 em 450

CHROMATOGRAM

Retention time: 7.7

Internal standard: IS1 6-fluoro-1-(2-fluoroethyl)-1,4-dihydro-7-(4-methyl-1-piperazinyl)-4-oxo-3-quinoline carboxylic acid (AM-735) (13.0), IS2 pipemidic acid (6.0)

Limit of detection: 1 ng/mL (plasma)

Limit of quantitation: 1000 ng/mL (urine), 10 ng/mL (plasma)

OTHER SUBSTANCES

Extracted: metabolites

Noninterfering: ciprofloxacin, norfloxacin, pefloxacin

KEY WORDS

plasma; column-switching; protect from light; pharmacokinetics

REFERENCE

Heizmann,P.; Dell,D.; Eggers,H.; Gora,R. Determination of the new fluoroquinolone fleroxacin and its N-demethyl and N-oxide metabolites in plasma and urine by high-performance liquid chromatography with fluorescence detection, *J.Chromatogr.*, **1990**, 527, 91–101.

SAMPLE

Matrix: cells

Sample preparation: Incubate cells in 2 mL 100 mM pH 3.0 glycine-HCl buffer for 2 h at room temperature, centrifuge at 5600 g for 5 min, inject an aliquot.

HPLC VARIABLES

Column: Bondapak C18

Mobile phase: MeCN:25 mM phosphoric acid adjusted to pH 3.0 with tetrabutylammonium hydroxide 25:75

Flow rate: 1.5

Detector: F ex 340 em 425

OTHER SUBSTANCES

Also analyzed: ciprofloxacin, lomefloxacin, norfloxacin, ofloxacin, temafloxacin

REFERENCE

Pascual,A.; Garcia,I.; Conejo,M.C.; Perea,E.J. Fluorometric and high-performance liquid chromatographic measurement of quinolone uptake by human neutrophils, *Eur.J.Clin.Microbiol.Infect.Dis.*, **1991**, 10, 969–971.

SAMPLE

Matrix: milk

Sample preparation: Milk + pipemidic acid + phosphate buffer, extract with dichloromethane: isopropanol 70:30, centrifuge. Remove an aliquot of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in mobile phase, add n-hexane, mix thoroughly, centrifuge, inject an aliquot of the aqueous phase.

HPLC VARIABLES

Column: 250 × 4.6 5 µm ODS-120T (Toyo Soda)

Mobile phase: MeOH:5 mM tetrabutylammonium hydrogen sulfate 28:72

Flow rate: 0.8

Detector: F ex 290 em 450

CHROMATOGRAM

Internal standard: pipemidic acid

Limit of quantitation: 100 ng/mL

KEY WORDS

protect from light; pharmacokinetics

REFERENCE

Dan,M.; Weidekamm,E.; Sagiv,R.; Portmann,R.; Zakut,H. Penetration of fleroxacin into breast milk and pharmacokinetics in lactating women, *Antimicrob.Agents Chemother.*, **1993**, 37, 293–296.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 125 × 4.6 3 µm ODS-Hypersil

Mobile phase: MeOH:THF:670 mM pH 3.0 phosphate buffer 20:0.8:79.2 plus 2 g/L tetrabutylammonium hydrogen sulfate and 2 mL/L 85% phosphoric acid

Flow rate: 1

Injection volume: 20

Detector: UV 278

CHROMATOGRAM

Retention time: 4.94

OTHER SUBSTANCES

Simultaneous: photodegradation products, ofloxacin

Interfering: ciprofloxacin

REFERENCE

Tiefenbacher,E.M.; Haen,E.; Przybilla,B.; Kurz,H. Photodegradation of some quinolones used as antimicrobial therapeutics, *J.Pharm.Sci.*, **1994**, 83, 463–467.

SAMPLE

Matrix: solutions

Sample preparation: Filter (0.45 μm) a solution in MeCN:water 10:90, inject an aliquot of the filtrate.

HPLC VARIABLES

Column: 250 \times 4.5 μm LiChrospher 100 RP-18

Mobile phase: MeCN:buffer 7:93 (Buffer was 25 mM phosphoric acid adjusted to pH 3.89 with 100 mM tetrabutylammonium hydroxide.)

Flow rate: 1

Injection volume: 10

Detector: UV 280

CHROMATOGRAM

Retention time: 7

OTHER SUBSTANCES

Simultaneous: ciprofloxacin, enoxacin, norfloxacin, ofloxacin (UV 295), pipemidic acid

REFERENCE

Barbosa,J.; Bergés,R.; Sanz-Nebot,V. Solvatochromic parameter values and pH in aqueous-organic mixtures used in liquid chromatography. Prediction of retention of a series of quinolones, *J.Chromatogr.A*, **1996**, 719, 27–36.

SAMPLE

Matrix: urine

Sample preparation: 250 μL Urine + 20 μL 30 $\mu\text{g/mL}$ pipemidic acid + 250 μL 500 mM pH 7.5 phosphate buffer, vortex, add 3.5 mL chloroform, shake at low speed for 10 min, centrifuge at 1155 g for 10 min. Remove the lower organic layer and add it to 200 μL 1/15 M pH 12.5 sodium phosphate, shake fast for 20 min, centrifuge at 1155 g for 10 min, inject a 150 μL aliquot of the aqueous layer.

HPLC VARIABLES

Column: 10 μm Nucleosil C18

Mobile phase: MeCN:MeOH:10 mM phosphate buffer containing 5 mM tetrabutylammonium hydrogen sulfate 13:6:81

Flow rate: 1

Injection volume: 150

Detector: F ex 274

CHROMATOGRAM

Internal standard: pipemidic acid

Limit of detection: 50 ng/mL

KEY WORDS

pharmacokinetics

REFERENCE


Bertino, J.S., Jr.; Nafziger, A.N.; Wong, M.; Stragand, L.; Puleo, C. Effect of a fat- and calcium-rich breakfast on pharmacokinetics of feroxacin administered in single and multiple doses, *Antimicrob. Agents Chemother.*, **1994**, *38*, 499–503.

Flezelastine

Molecular formula: $C_{29}H_{30}F_3O$

Molecular weight: 416.58

CAS Registry No.: 135381-77-0



SAMPLE

Matrix: microsomal incubations

Sample preparation: Adjust pH of 750 μ L microsomal incubation to 10 with 100 mM NaOH, add 3 mL cyclohexane:ethyl acetate 80:20, shake mechanically for 10 min, centrifuge at 2500 g for 10 min, repeat the extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute with a solution of azelastine in MeOH, evaporate to dryness under a stream of nitrogen, reconstitute with 200 μ L MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: $4 \times 4.5 \mu\text{m}$ LiChrospher Si-60

Column: $250 \times 4.5 \mu\text{m}$ LiChrospher Si-60

Mobile phase: MeOH containing 0.033% perchloric acid

Flow rate: 0.5 for 17 min then 0.9

Injection volume: 20

Detector: F ex 210 em 360

CHROMATOGRAM

Retention time: 15

Internal standard: azelastine (27)

Limit of detection: 125 $\mu\text{g/mL}$

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

human; rat; cow; pig; liver

REFERENCE

Paris, S.; Blaschke, G.; Locher, M.; Borbe, H.O.; Engel, J. Investigation of the stereoselective in vitro metabolism of the chiral antiasthmatic/antiallergenic drug flezelastine by high-performance liquid chromatography and capillary zone electrophoresis, *J. Chromatogr. B*, **1997**, *691*, 463–471.

Floctafenine


Molecular formula: $C_{20}H_{17}F_3N_2O_4$

Molecular weight: 406.36

CAS Registry No.: 23779-99-9

Merck Index: 4140

Lednicer No.: 3 184



SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol:n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 353

CHROMATOGRAM

Retention time: 10.73

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfapyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; car-teolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihy-dralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazo-lam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lor-azepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temaze-pam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; almino-profen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acen-ocoumarol; vindesine; mexiletine; dipyrindamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; fle-cainide; phencyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; bupren-orphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazo-lam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclam-ide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; metha-done; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidra-mine; proprietyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvox-amine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpi-pramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254–262.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 209.9

CHROMATOGRAM

Retention time: 17.177

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

Flosequinan


Molecular formula: C₁₁H₁₀FNO₂S

Molecular weight: 239.27

CAS Registry No.: 76568-02-0

Merck Index: 4146

Lednicer No.: 5 127



SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 10 µL 100 µg/mL IS in MeCN:MeOH 50:50 + 5 mL chloroform, shake for 10 min, centrifuge at 1800 g for 10 min. Remove 4 mL of the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100 µL MeOH, inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.6 10 µm Chiralcel OD

Mobile phase: MeOH:EtOH 78:22

Column temperature: 30

Flow rate: 0.7

Detector: UV 320

CHROMATOGRAM

Retention time: 12.1 (R(+)), 7.9 (S(-))

Internal standard: (±)-7-chloro-1-methyl-3-methylsulfinyl-4-quinolone (8.9, 14.6)

Limit of quantitation: 5 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; chiral; pharmacokinetics

REFERENCE

Kashiyama,E.; Odomi,M.; Shimizu,T. Stereospecific and simultaneous high-performance liquid chromatographic assay of flosequinan and its metabolites in human plasma, *J.Chromatogr.B*, **1994**, 652, 179–185.

SAMPLE

Matrix: blood, urine

Sample preparation: Condition a 1 mL Baker-10 octadecyl C18 SPE cartridge (J.T. Baker) with MeOH then water. 100 μ L Plasma or urine + 25 μ L 2 (plasma) or 10 (urine) μ g/mL IS in water, add to the SPE cartridge, wash with 2 column volumes of water, elute with 200 μ L MeOH, inject a 2 (urine) or 5 (plasma) aliquot of the eluate.

HPLC VARIABLES

Column: 150 \times 3.9 5 μ m Nova-Pak C18

Mobile phase: MeCN:MeOH:water 7:20:73

Flow rate: 1.2

Injection volume: 2-5

Detector: UV 254

CHROMATOGRAM

Retention time: 3.4

Internal standard: 7-chloro-1-methyl-3-methylsulfinyl-4-quinolone (7.8)

Limit of quantitation: 50 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

Noninterfering: acetaminophen, amiodarone, diazepam, digoxin, furosemide, lidocaine, nitroglycerin, procainamide, propranolol, quinidine

KEY WORDS

plasma; SPE

REFERENCE

Slegowski,M.B.; Miller,C.; Porter,R.S. Simplified high-performance liquid chromatographic determination of flosequinan and its metabolite in plasma, serum and urine, *J.Chromatogr.*, **1988**, 425, 227–232.

SAMPLE

Matrix: incubations

Sample preparation: 2 mL Incubation + 10 μ g IS + 4 mL chlorobutane:1,2-dichloroethane 80:20, shake for 15 min, centrifuge. Remove the upper organic layer and evaporate it to dryness under a stream of nitrogen at 70°, reconstitute the residue in 800 μ L mobile phase, inject a 100 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 Hypersil 30DS

Mobile phase: MeCN:MeOH:water 7:20:73

Flow rate: 2

Injection volume: 100

Detector: UV 254

CHROMATOGRAM

Internal standard: 7-chloro-1-methyl-3-methylsulfinyl-4-quinolone (BTS 49037)

OTHER SUBSTANCES

Extracted: metabolites

REFERENCE

Lee, S.C.; Renwick, A.G. Sulphoxide reduction by rat intestinal flora and by *Escherichia coli* *in vitro*, *Biochem. Pharmacol.*, **1995**, 49, 1567–1576.


Floxacillin

Molecular formula: C₁₉H₁₇ClFN₃O₅S

Molecular weight: 453.88

CAS Registry No.: 5250-39-5, 1847-24-1 (sodium salt), 34214-51-2 (sodium monohydrate)

Merck Index: 4147



SAMPLE

Matrix: blood

Sample preparation: Condition a 1 mL Bond Elut C18 SPE cartridge with 2 mL MeCN and 1 mL 10 mM pH 2 Na₂HPO₄. 250 µL Plasma + 100 µL 20 µg/mL dicloxacillin sodium in water, add 400 µL MeCN at -15° while vortexing, add 700 µL 10 mM pH 2 Na₂HPO₄, centrifuge at 8000 g for 10 min. Add the supernatant to the SPE cartridge, wash with 1 mL water, elute with two 500 µL portions of MeCN:water 35:65 containing 10 mM Na₂HPO₄ (pH adjusted to 6 with phosphoric acid), inject a 20 µL aliquot of the eluate.

HPLC VARIABLES

Column: 100 × 2.5 µm ODS Hypersil

Mobile phase: MeCN:water 40:60 containing 10 mM Na₂HPO₄, pH adjusted to 2 with orthophosphoric acid

Flow rate: 0.5

Injection volume: 20

Detector: UV 220

CHROMATOGRAM

Retention time: 2.5

Internal standard: dicloxacillin (3.5)

Limit of detection: 100 ng/mL

OTHER SUBSTANCES

Simultaneous: cloxacillin

KEY WORDS

plasma; pharmacokinetics; SPE

REFERENCE

Hung, C.T.; Lim, J.K.C.; Zoest, A.R.; Lam, F.C. Optimization of high-performance liquid chromatographic analysis for isoxazolyl penicillins using factorial design, *J. Chromatogr.*, **1988**, 425, 331–341.

SAMPLE

Matrix: blood

Sample preparation: 100 µL Plasma + 100 µL dicloxacillin in water + 25 µL glacial acetic acid + 2 mL ethyl acetate, vortex for 30 s, centrifuge at 2000 g for 5 min. Remove the supernatant and evaporate it to dryness under a stream of nitrogen at 70°, reconstitute the residue in 250 µL mobile phase, inject a 10–20 µL aliquot.

HPLC VARIABLES

Column: 40 × 3.2 RP18 VeloSep (Brownlee)

Mobile phase: MeCN:10 mM pH 7 phosphate buffer 18:82

Flow rate: 1.2

Injection volume: 10–20

Detector: UV 220