
REFERENCE

Lee, S.C.; Renwick, A.G. Sulphoxide reduction by rat intestinal flora and by *Escherichia coli* *in vitro*, *Biochem. Pharmacol.*, **1995**, 49, 1567–1576.


Floxacillin

Molecular formula: C₁₉H₁₇ClFN₃O₅S

Molecular weight: 453.88

CAS Registry No.: 5250-39-5, 1847-24-1 (sodium salt), 34214-51-2 (sodium monohydrate)

Merck Index: 4147



SAMPLE

Matrix: blood

Sample preparation: Condition a 1 mL Bond Elut C18 SPE cartridge with 2 mL MeCN and 1 mL 10 mM pH 2 Na₂HPO₄. 250 µL Plasma + 100 µL 20 µg/mL dicloxacillin sodium in water, add 400 µL MeCN at -15° while vortexing, add 700 µL 10 mM pH 2 Na₂HPO₄, centrifuge at 8000 g for 10 min. Add the supernatant to the SPE cartridge, wash with 1 mL water, elute with two 500 µL portions of MeCN:water 35:65 containing 10 mM Na₂HPO₄ (pH adjusted to 6 with phosphoric acid), inject a 20 µL aliquot of the eluate.

HPLC VARIABLES

Column: 100 × 2.5 µm ODS Hypersil

Mobile phase: MeCN:water 40:60 containing 10 mM Na₂HPO₄, pH adjusted to 2 with orthophosphoric acid

Flow rate: 0.5

Injection volume: 20

Detector: UV 220

CHROMATOGRAM

Retention time: 2.5

Internal standard: dicloxacillin (3.5)

Limit of detection: 100 ng/mL

OTHER SUBSTANCES

Simultaneous: cloxacillin

KEY WORDS

plasma; pharmacokinetics; SPE

REFERENCE

Hung, C.T.; Lim, J.K.C.; Zoest, A.R.; Lam, F.C. Optimization of high-performance liquid chromatographic analysis for isoxazolyl penicillins using factorial design, *J. Chromatogr.*, **1988**, 425, 331–341.

SAMPLE

Matrix: blood

Sample preparation: 100 µL Plasma + 100 µL dicloxacillin in water + 25 µL glacial acetic acid + 2 mL ethyl acetate, vortex for 30 s, centrifuge at 2000 g for 5 min. Remove the supernatant and evaporate it to dryness under a stream of nitrogen at 70°, reconstitute the residue in 250 µL mobile phase, inject a 10–20 µL aliquot.

HPLC VARIABLES

Column: 40 × 3.2 RP18 VeloSep (Brownlee)

Mobile phase: MeCN:10 mM pH 7 phosphate buffer 18:82

Flow rate: 1.2

Injection volume: 10–20

Detector: UV 220


CHROMATOGRAM**Retention time:** 2.8**Internal standard:** dicloxacillin (4.4)**Limit of detection:** 50 ng/mL**Limit of quantitation:** 300 ng/mL

OTHER SUBSTANCES**Simultaneous:** carbamazepine, phenytoin, phenobarbital**Noninterfering:** acetaminophen, N-acetylprocainamide, amikacin, amitriptyline, caffeine, chloramphenicol, cyclosporine, digoxin, ethosuximide, gentamicin, lidocaine, nortriptyline, methotrexate, primidone, procainamide, quinidine, salicylic acid, theophylline, tobramycin, valproic acid, vancomycin, metabolites

KEY WORDSplasma

REFERENCECharles, B.G.; Foo, C.C.; Gath, J. Rapid column liquid chromatographic analysis of flucloxacillin in plasma on a microparticulate pre-column, *J.Chromatogr.B*, **1994**, 660, 186–190.

Floxuridine

Molecular formula: C₉H₁₁FN₂O₅**Molecular weight:** 246.20**CAS Registry No.:** 50-91-9**Merck Index:** 4148

SAMPLE**Matrix:** blood**Sample preparation:** 1 mL Plasma + 100 µL water, vortex 10 s, add 200 µL water containing 30% w/v trichloroacetic acid and 30% w/v perchloric acid, vortex 10 s, place in an ice bath for 2 min, centrifuge at 3000 g for 20 min, inject a 20 µL aliquot of the supernatant.

HPLC VARIABLES**Column:** 100 × 4.6 5µm Hypersil-ODS**Mobile phase:** 20 mM Na₂HPO₄ adjusted to pH 2.0 with orthophosphoric acid**Flow rate:** 2**Injection volume:** 20**Detector:** UV 254

CHROMATOGRAM**Retention time:** 7**Internal standard:** floxuridine

OTHER SUBSTANCES**Extracted:** allopurinol, oxipurinol

KEY WORDSplasma; floxuridine is IS

REFERENCEHung, C.T.; Zoest, A.R.; Perrier, D.G. Analysis of allopurinol and oxipurinol in plasma by reversed phase HPLC, *J.Liq.Chromatogr.*, **1986**, 9, 2471–2483.

SAMPLE**Matrix:** blood

Sample preparation: 100 μ L Serum + 5 μ L 42 μ M 5-chlorouracil in acetone, mix, let stand at room temperature for 3 min, add 500 μ L 100 mM pH 3.5 potassium phosphate buffer, add 2 mL ethyl acetate, vortex for 2 min, centrifuge at 1000 g for 5 min. Remove a 1.4 mL aliquot of the organic layer and evaporate it to dryness under reduced pressure. Add 20 mg solid potassium bicarbonate:anhydrous sodium sulfate 1:7 to the residue, add 50 μ L 1.3 mM 3-bromomethyl-6,7-dimethoxy-1-methyl-2(1H)-quinoxalinone in acetone, add 50 μ L 1.5 mM 18-crown-6 in acetone, heat at 50° for 20 min, cool, inject a 10 μ L aliquot. (Silanize all glassware. Synthesize 3-bromomethyl-6,7-dimethoxy-1-methyl-2(1H)-quinoxalinone as follows. Stir 483 g veratrole in 1.45 L acetic acid at 15° for 1 h, add 683 g concentrated nitric acid (d 1.05) over 1 h (maintain the temperature below 40° by cooling and regulating the rate of addition of the nitric acid). Continue stirring and add 2.127 L fuming nitric acid (d 1.50) over 1 h while maintaining the temperature below 30°, let stand for 2 h, pour into a large volume of cold water, filter, wash the solid with water until the washings are neutral, recrystallize from EtOH to give 4,5-dinitroveratrole (mp 129.5-130.5°) (J. Am. Chem. Soc. 1946, 68, 1536). Reflux 5 g 4,5-dinitroveratrole in 200 mL benzene (Caution! Benzene is a carcinogen!), add 100 g 60 mesh iron powder and 20 mL concentrated HCl in small portions over 1 h, reflux for 4 h, add 10 mL water, reflux for 2 h, cool, make alkaline with 2.5 M NaOH, extract several times with 200 mL portions of benzene. Combine the organic layers and evaporate them to dryness, add 10 mL concentrated HCl, recrystallize from EtOH to give 1,2-diamino-4,5-dimethoxybenzene monohydrochloride as very slightly pink needles (mp 240°) (Anal. Chim. Acta 1982, 134, 39). Heat 2.5 mmoles 1,2-diamino-4,5-dimethoxybenzene hydrochloride and 2.4 mmoles pyruvic acid in 30 mL 500 mM HCl on a boiling water bath for 2 h, cool with ice-water, filter. Wash the precipitate with water and dry it under vacuum, recrystallize from MeOH:water 90:10 to give 6,7-dimethoxy-3-methyl-2(1H)-quinoxalinone as yellow needles (mp 255°) (Chem. Pharm. Bull. 1985, 33, 3493). Treat 1 g 6,7-dimethoxy-3-methyl-2(1H)-quinoxalinone dissolved in 50 mL anhydrous MeOH with a solution of diazomethane in ether, evaporate to dryness under reduced pressure, dissolve the residue in 5 mL ethyl acetate, chromatograph on a 250 \times 35 column filled with 130 g 70-230 mesh silica gel 60 (Merck) using n-hexane:ethyl acetate 25:75 to give 6,7-dimethoxy-1,3-dimethyl-2(1H)-quinoxalinone as yellow needles (mp 170-171°). Dissolve 350 mg 6,7-dimethoxy-1,3-dimethyl-2(1H)-quinoxalinone in 3 mL acetic acid, add 350 mg anhydrous sodium acetate, add 2 mL 1.5 M bromine in acetic acid, heat at 100° for 15 min, cool, add 10 mL ether, filter, wash the solid 2 or 3 times with small portions of ether. Combine the filtrate and washings and evaporate them to dryness, dissolve the residue in 5 mL ethyl acetate, chromatograph on a 250 \times 35 column filled with 130 g 70-230 mesh silica gel 60 (Merck) using ether, evaporate the main fraction to dryness, recrystallize the residue from n-hexane:ethyl acetate 50:50 to give 3-bromomethyl-6,7-dimethoxy-1-methyl-2(1H)-quinoxalinone as yellow needles (mp 161-163°) (J. Chromatogr. 1985, 346, 227). 3-Bromomethyl-6,7-dimethoxy-1-methyl-2(1H)-quinoxalinone is also available from Dojindo Molecular Technologies, Inc., 3 Bethesda Metro Center, Suite 700, Bethesda MD 20814; (301) 664-8448; www.dojindo.co.jp.)

HPLC VARIABLES

Column: 100 \times 8 10 μ m Radial Pak C18 (Waters) (Wash with MeOH at 2 mL/min for 20 min at the end of each day.)

Mobile phase: Gradient. MeOH:water 35:65 for 15 min, 50:50 for 25 min (step gradient), re-equilibrate at initial conditions for 20 min.

Flow rate: 1.5

Injection volume: 10

Detector: F ex 370 em 455

CHROMATOGRAM

Retention time: 28.1

Internal standard: 5-chlorouracil (32.5)

Limit of detection: 25 ng/mL

OTHER SUBSTANCES

Extracted: 5-fluorouracil

KEY WORDS

derivatization; serum; pharmacokinetics

REFERENCE

Yamaguchi, M.; Nakamura, M.; Kuroda, N.; Ohkura, Y. Determination of 5-fluorouracil and 5-fluoro-2'-deoxyuridine in human serum by high-performance liquid chromatography with fluorescence detection, *Anal. Sci.*, 1987, 3, 75-79.

SAMPLE**Matrix:** blood**Sample preparation:** Add 250 μ L serum to a 20 \times 7 DEAE-Cellulofine AM anion-exchange column (Seikagaku Tokyo), elute with 3.5 mL 1 mM HCl, discard the first 0.5 mL eluate, collect the next 3 mL eluate. Evaporate the eluate to 0.5 mL under reduced pressure, add 15 mL ethyl acetate, shake, centrifuge. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 800 μ L anhydrous acetone, add 100 μ L 750 μ g/mL 4-bromomethyl-6,7-dimethoxycoumarin in acetone, add 100 μ L 250 μ g/mL 18-crown-6 in acetone, add 1.5 mg anhydrous potassium carbonate, heat at 70° for 15 min (protect from atmospheric moisture with a calcium chloride drying tube), cool, inject an aliquot

HPLC VARIABLES**Column:** 200 \times 4.5 μ m Nucleosil 5 C18**Mobile phase:** MeOH:water 60:40**Flow rate:** 0.8**Detector:** F ex 340 em 420

CHROMATOGRAM**Retention time:** 5**Limit of quantitation:** 100 ng/mL

OTHER SUBSTANCES**Extracted:** florasur, 5-fluorouracil

KEY WORDS

derivatization; serum; protect from light; SPE

REFERENCE

Yoshida, S.; Adachi, T.; Hirose, S. 4-Bromomethyl-6,7-dimethoxycoumarin as a fluorescence reagent for pre-column derivatization of 5-fluorouracil compounds in high-performance liquid chromatography, *J. Chromatogr.*, **1988**, 430, 156–162.

SAMPLE**Matrix:** blood**Sample preparation:** Condition a 1 mL LC-SCX Supelclean strong cation-exchange SPE cartridge (Supelco) with 2 mL MeOH, 1 mL 100 mM copper(II) sulfate solution, and 3 mL 50 mM pH 7 phosphate buffer, do not allow to dry. 300 μ L Serum + 5-bromouracil, add to the SPE cartridge, wash with 2 mL 50 mM pH 7 phosphate buffer, wash with 2 mL MeOH, elute with 700 μ L 1.7 M ammonia solution, add 70 μ L glacial acetic acid to the eluate, mix thoroughly, inject a 20 μ L aliquot.

HPLC VARIABLES**Guard column:** 20 \times 4.6 5 μ m Supelguard LC-18-S (Supelco)**Column:** 250 \times 4.6 5 μ m Supelcosil LC-18-S ODS**Mobile phase:** Gradient. A was MeOH:50 mM pH 6.5 phosphate buffer 60:40. B was 50 mM pH 6.5 phosphate buffer.**Flow rate:** 1**Injection volume:** 20**Detector:** UV 269

CHROMATOGRAM**Retention time:** 13**Internal standard:** 5-bromouracil (12)**Limit of detection:** 70 ng/mL

OTHER SUBSTANCES**Extracted:** doxifluridine, 5-fluorouracil, 5-fluorouridine monophosphate, metabolites

KEY WORDS

serum; SPE

REFERENCE

Guerrieri,A.; Palmisano,F.; Zamboni,P.G.; De Lena,M.; Lorusso,V. Solid-phase extraction of fluoropyrimidine derivatives on a copper-modified strong cation exchanger: determination of doxifluridine, 5-fluorouracil and its main metabolites in serum by high-performance liquid chromatography with ultraviolet detection, *J.Chromatogr.*, **1993**, 617, 71-77.

SAMPLE

Matrix: blood, peritoneal fluid

Sample preparation: Plasma. 1 mL Plasma + 10 μ L 100 μ M bromouridine + 70 μ L perchloric acid, mix thoroughly, let stand at 4° for at least 12 h, centrifuge for 5 min. Remove the supernatant and adjust the pH to 7 with 5 M KOH, let stand on ice for 2 h, inject a 20 μ L aliquot. Peritoneal fluid. 1 mL Peritoneal fluid + 10 μ L 100 μ M bromouridine, dilute 1:100 with water, inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: 12.5 \times 4.6 5 μ m Zorbax RX

Column: 250 \times 4.6 5 μ m Zorbax RX

Mobile phase: Gradient. A was 25 mL pH 2.5 ammonium phosphate. B was MeCN:25 mM pH 7.5 ammonium phosphate 7:93. A:B 100:0 for 5 min, to 0:100 over 10 min, maintain at 0:100 for 10 min, return to initial conditions over 1 min, re-equilibrate for 20 min.

Column temperature: 20

Flow rate: 1

Injection volume: 20

Detector: UV 270

CHROMATOGRAM

Retention time: 16-17

Internal standard: bromouridine (18)

Limit of detection: 2.5 nM

OTHER SUBSTANCES

Extracted: 5-fluorouracil

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Smith-Rogers,J.A.; Tong,W.P.; Duafala,M.E.; Markman,M.; Bertino,J.R. High-performance liquid chromatographic method for the simultaneous measurement of floxuridine and fluorouracil in human body fluids, *J.Chromatogr.*, **1991**, 566, 147-154.

SAMPLE

Matrix: blood, tissue

Sample preparation: Add 100 μ L 10% perchloric acid and 20 μ L 200 μ g/mL IS to 100 μ L serum or homogenized tissue. Shake for 2 min and centrifuge at 2000 g for 10 min. Filter (45 μ m) supernatant, inject a 10 μ L aliquot of the filtrate.

HPLC VARIABLES

Guard column: 45 \times 4.6 5 μ m ODS Hypersil (VDS Optilab)

Column: 250 \times 4.6 5 μ m ODS Hypersil (VDS Optilab)

Mobile phase: MeOH:99% acetic acid:water 3:0.5:96.95

Column temperature: 30

Flow rate: 1

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 12.75

Internal standard: 5-bromouracil (10.34)

Limit of quantitation: 200 ng/mL (serum); 600 ng/mL (tissue)

OTHER SUBSTANCES

Extracted: metabolites, fluorouracil

KEY WORDS

serum; rat; liver; tumor; kidney; spleen; peritoneum; gastric mucosa; lung; heart; pancreas

REFERENCE

Jung,M.; Berger,G.; Pohlen,U.; Päuser,S.; Reszka,R.; Buhr,H.J. Simultaneous determination of 5-fluorouracil and its active metabolites in serum and tissue by high-performance liquid chromatography, *J.Chromatogr.B*, **1997**, 702, 193–202.

SAMPLE

Matrix: blood, tissue

Sample preparation: 0.5 g Tissue or 1 mL plasma + 5 μ L 1 M sulfuric acid (lung and heart only) + 2 μ L 1 M sulfuric acid (plasma only) + 500 μ L 200 mg/mL sodium sulfate (liver and kidney only) + 50 μ L 1 M pH 6 sodium acetate (liver only) + 50 μ L 1 M pH 5 sodium acetate (kidney only) + 100 μ L 2% trichloroacetic acid (lung and heart only) + 15 mL n-propanol:ether (liver 16:84, kidney 20:80, lung 88:12, heart 40:60, plasma 88:12), sonicate for 30 s, shake for 15 min, centrifuge for 15 min. Remove the organic layer and evaporate it to dryness under reduced pressure, reconstitute the residue in 1 mL 50 mM ammonium phosphate (liver pH 11, kidney pH 3, lung pH 2.5, heart pH 5, plasma pH 2.5), inject a 20 μ L aliquot. (From *J. Liq.Chromatogr.* 1994, 17, 1621.)

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Spherisorb 5 ODS 2

Mobile phase: MeCN:50 mM phosphate buffer 0.5:99.5 (Liver pH 3, kidney pH 6, lung pH 5, heart pH 5, plasma pH 2.5)

Column temperature: 10 (kidney), 35 (lung), 20 (heart), 25 (plasma), 15 (liver)

Flow rate: 1

Injection volume: 20

Detector: UV 200

CHROMATOGRAM

Retention time: 18 (plasma), 28 (lung), 30 (liver), 32 (kidney), 36 (heart)

Internal standard: flucytosine (for plasma) (4), 4-chlorouracil (for tissue) (17 (lung), 18 (liver), 11 (kidney), 19 (heart))

Limit of quantitation: 670 ng/g plasma, 110 ng/g (heart), 90 ng/g (lung), 210 ng/g (kidney), 500 ng/g (liver)

OTHER SUBSTANCES

Extracted: metabolites, 5-fluorouracil

KEY WORDS

plasma; rabbit; liver; kidney; lung; heart; pharmacokinetics

REFERENCE

Del Nozal,M.J.; Bernal,J.L.; Pampliega,A.; Marinero,P.; Pozuelo,M. Determination of the concentrations of 5-fluorouracil and its metabolites in rabbit plasma and tissues by high-performance liquid chromatography, *J.Chromatogr.B*, **1994**, 656, 397–405.

SAMPLE

Matrix: formulations

Sample preparation: Inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 30 mm long 3 μ m C18 (Perkin-Elmer)

Mobile phase: MeOH:water 0.5:95.5 (sic)

Flow rate: 1

Injection volume: 10

Detector: UV 268

CHROMATOGRAM**Retention time:** 2

KEY WORDS

stability-indicating; saline; injections

REFERENCE

Smith, J.A.; Morris, A.; Duafala, M.E.; Bertino, J.R.; Markman, M.; Kleinberg, M. Stability of floxuridine and leucovorin calcium admixtures for intraperitoneal administration, *Am.J.Hosp.Pharm.*, **1989**, *46*, 985–989.

SAMPLE**Matrix:** formulations**Sample preparation:** Dilute formulation 1:100 with water, inject a 50 µL aliquot.

HPLC VARIABLES**Guard column:** 5 × 4 35-60 µm Perisorb RP18**Column:** 250 × 4 10 µm LiChrosorb RP18**Mobile phase:** MeOH:300 mM sodium acetate 2:98**Injection volume:** 50**Detector:** UV 254

CHROMATOGRAM**Retention time:** 10.5

OTHER SUBSTANCES**Simultaneous:** fluorouracil

KEY WORDS

injections; water

REFERENCE

Sadjak, A.; Wintersteiger, R. Compatibility of morphine, baclofen, floxuridine and fluorouracil in an implantable medication pump, *Arzneimittelforschung*, **1995**, *45*, 93–98.

SAMPLE**Matrix:** solutions**Sample preparation:** Prepare 50 mg/mL or 1 mg/mL solutions in 0.9% sodium chloride injection, dilute 1:1000 or 1:100 with mobile phase, inject an aliquot.

HPLC VARIABLES**Column:** Bakerbond C18**Mobile phase:** MeOH:water 0.5:99.5**Flow rate:** 1**Detector:** UV 268

CHROMATOGRAM**Retention time:** 4.6

OTHER SUBSTANCES**Simultaneous:** degradation products

KEY WORDS

stability-indicating

REFERENCE

Stiles, M.L.; Allen, L.V., Jr.; Prince, S.J. Stability of deferoxamine mesylate, floxuridine, fluorouracil, hydromorphone hydrochloride, lorazepam, and midazolam hydrochloride in polypropylene infusion-pump syringes, *Am.J.Health-Syst.Pharm.*, **1996**, *53*, 1583–1588.

SAMPLE**Matrix:** solutions**Sample preparation:** Add 10 mg of 4-bromomethyl-7-methoxycoumarin and 10 mg potassium carbonate to 5 mL of a solution in DMSO, after 5 min at room temperature add 4-nitrobenzoic acid.

HPLC VARIABLES**Column:** 200 × 4.5 µm Nucleosil 5 C18**Mobile phase:** MeCN:MeOH:water 5:29:66**Flow rate:** 0.6**Detector:** F ex 346 em 395

CHROMATOGRAM**Retention time:** 12.5**Internal standard:** chlorodeoxyuridine (16.5)

OTHER SUBSTANCES**Simultaneous:** inosine, thymine, uracil, uridine


KEY WORDS

derivatization; some details of serum extraction in paper

REFERENCE

Yoshida, S.; Hirose, S.; Iwamoto, M. Use of 4-bromomethyl-7-methoxycoumarin for derivatization of pyrimidine compounds in serum analysed by high-performance liquid chromatography with fluorimetric detection, *J. Chromatogr.*, **1986**, 383, 61–68.

Flubendazole

**Molecular formula:** C₁₆H₁₂FN₃O₃**Molecular weight:** 313.29**Merck Index:** 4154**Lednicer No.:** 2 354

SAMPLE**Matrix:** blood, urine**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES**Guard column:** 20 mm long Symmetry C18**Column:** 250 × 4.6 µm Symmetry C8 (Waters)**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.**Column temperature:** 30**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)**Injection volume:** 10–30**Detector:** UV 211.1

CHROMATOGRAM**Retention time:** 16.745**KEY WORDS**

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149–163.

SAMPLE**Matrix:** solutions**HPLC VARIABLES****Column:** 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water; make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30**Flow rate:** 2**Detector:** UV 210**OTHER SUBSTANCES**

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenopropfen, fenproporex, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiacol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarbostyryl, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephénytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methylpyrrolone, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxyamphetamine, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phenacyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidone, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopalamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfathiazole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole,

thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranylcypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleennamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.


Fluconazole

Molecular formula: C₁₃H₁₂F₂N₆O

Molecular weight: 306.27

CAS Registry No.: 86386-73-4

Merck Index: 4158



SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 100 μ L 700 μ g/mL IS + 1.0 mL 1 M KOH, mix, extract twice with 3 mL portions of ethyl acetate. Combine the organic layers and extract with 2 mL 1 M HCl. Add 1.75 mL 1 M KOH to the aqueous layer and extract twice with 3 mL portions of ethyl acetate. Combine the organic layers and evaporate them to dryness under a stream of dry nitrogen, reconstitute the residue in 100 μ L MeCN:10 mM pH 9.0 ammonium phosphate 50:50, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 3 μ m Microsorb C18

Mobile phase: Gradient. MeCN:10 mM pH 9.0 ammonium phosphate 10:90 for 1 min, to 60:40 over 15 min, maintain at 60:40 for 2 min.

Flow rate: 0.5

Injection volume: 20

Detector: UV 260

CHROMATOGRAM

Retention time: 7

Internal standard: α -(2,4-dichlorophenyl)-1H-imidazole-1-ethanol (9)

Limit of quantitation: 12 μ g/mL

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Black,D.J.; Kunze,K.L.; Wienkers,L.C.; Gidal,B.E.; Seaton,T.L.; McDonnell,N.D.; Evans,J.S.; Bauwens,J.E.; Trager,W.F. Warfarin-fluconazole II. A metabolically based drug interaction: In vivo studies, *Drug Metab.Dispos.*, **1996**, *24*, 422-428.

SAMPLE

Matrix: blood

Sample preparation: Add 500 μ L MeCN containing an excess of sodium carbonate to 500 μ L plasma, vortex for 10 s, centrifuge at 2000 g for 10 min, dilute a 350 μ L aliquot of the supernatant with 700 μ L distilled water, inject an aliquot.

HPLC VARIABLES

Column: 125 \times 4 5 μ m Ultrabase C8 (SFCC, Neuilly Plaisance, France)

Mobile phase: MeCN:water 28:72

Flow rate: 1

Injection volume: 40

Detector: UV 260

CHROMATOGRAM

Retention time: 3.6

Limit of detection: 150 ng/mL

Limit of quantitation: 400 ng/mL

OTHER SUBSTANCES

Noninterfering: amikacin, amoxicillin, carbamazepine, cilastatin, clavulanic acid, epoxycarbamazepine, flunitrazepam, imipenem, methotrexate, midazolam, ofloxacin, phenobarbital, phenytoin, piperacillin, teicoplanin, theophylline, thiopental, tobramycin, valproic acid, vancomycin

KEY WORDS

plasma

REFERENCE

Cociglio,M.; Brandissou,S.; Alric,R.; Bressolle,F. High-performance liquid chromatographic determination of fluconazole in plasma, *J.Chromatogr.B*, **1996**, 686, 11–17.

SAMPLE

Matrix: blood

Sample preparation: Condition a 6 mL 500 mg Bakerbond C18 SPE cartridge with 2 column volumes of MeOH and 2 column volumes of buffer. 1 mL Plasma + 100 μ L 100 μ g/mL IS in buffer + 2 mL buffer, mix, add to the SPE cartridge, wash with 1 mL buffer, wash with 1 mL MeOH:buffer 15:85, dry under vacuum (8-9 inches Hg) for 1 min, elute with 500 μ L MeOH without vacuum, elute with 500 μ L MeOH under 3-4 inches Hg vacuum. Combine the eluates and evaporate them to dryness under a stream of nitrogen at 40°, reconstitute with 200 μ L mobile phase, inject a 15 μ L aliquot. (Buffer was 100 mM pH 6.0 sodium phosphate buffer.)

HPLC VARIABLES

Column: 250 \times 4.5 μ m Nucleosil C18

Mobile phase: MeCN:25 mM pH 7.0 Tris buffer 25:75

Flow rate: 1

Injection volume: 15

Detector: UV 210

CHROMATOGRAM

Retention time: 5

Internal standard: 2-(2,4-difluorophenyl)-1-fluoro-1,3-bis(1H-1,2,4-triazol-1-yl)propan-2-ol (Pfizer UK-54373) (8)

Limit of detection: 100 ng/mL

KEY WORDS

plasma; comparison with capillary electrophoresis; pharmacokinetics; SPE

REFERENCE

von Heeren,F.; Tanner,R.; Theurillat,R.; Thormann,W. Determination of fluconazole in human plasma by micellar electrokinetic capillary chromatography with detection at 190 nm, *J.Chromatogr.A*, **1996**, 745, 165–172.

SAMPLE

Matrix: blood, dialysate

Sample preparation: Plasma. Extract 100 μ L plasma with 5 mL dichloromethane, centrifuge at 1000 g for 15 min. Evaporate the organic to dryness under a stream of nitrogen. Reconstitute the residue with 100 μ L mobile phase, inject a 20 μ L aliquot. Dialysate. Inject a 10 μ L aliquot directly.

HPLC VARIABLES

Column: 200 \times 2.1 5 μ m HP ODS

Mobile phase: MeCN:buffer 13:87 (Buffer was 20 mM ammonium monobasic phosphate, adjusted to pH 7 with 5 M NaOH.)

Flow rate: 0.5

Injection volume: 10-20

Detector: UV 210

CHROMATOGRAM

Internal standard: UK-51,060 (plasma)

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Yang,H.; Wang,Q.; Elmquist,W.F. Fluconazole distribution to the brain: a crossover study in freely-moving rats using in vivo microdialysis, *Pharm.Res.*, **1996**, 13, 1570-1575.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 11.398

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.


Flucytosine

Molecular formula: C₄H₄FN₃O

Molecular weight: 129.09

CAS Registry No.: 2022-85-7

Merck Index: 4161



SAMPLE

Matrix: blood

Sample preparation: 500 μL Plasma + 50 μL 100 $\mu\text{g/mL}$ 5-chlorouracil in MeOH + 30 μL 1 M pH 4.0 sodium acetate buffer + 300 μL 20% sodium sulfate in water + 8 mL diethyl ether:n-propanol 75:25, vortex for 5 min, centrifuge at 3000 rpm for 10 min. Remove the organic layer and evaporate it to dryness under vacuum using a freeze dryer, reconstitute the residue in 1 mL water, vortex, filter (0.22 μm), inject a 20 μL aliquot.

HPLC VARIABLES

Column: 150 \times 6 5 μm YMC A-312 octadecylsilane (Yamamura Chemical)

Mobile phase: MeOH:50 mM pH 3.0 phosphate buffer 1:99

Flow rate: 1.5

Injection volume: 20

Detector: UV 275

CHROMATOGRAM

Retention time: 8

Internal standard: 5-chlorouracil (3)

Limit of quantitation: 200 ng/mL

KEY WORDS

dog; plasma; pharmacokinetics

REFERENCE

Bonny,J.-D.; Kyowa,M. Use of in vitro release tests for the prediction of the in vivo behavior and the development of flucytosine controlled-release capsules, *J.Pharm.Sci.*, **1995**, *84*, 619–623.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10–30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 3.052

KEY WORDS

whole blood

REFERENCE

Gaillard,X.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, *763*, 149–163.

SAMPLE

Matrix: formulations

Sample preparation: Dilute 30 μL of the injection to 100 mL with buffer, add 50 μL 10 mg/mL 5-flucytosine in water, inject a 20 μL aliquot. (Buffer contained 1 g/L $\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$ and 0.363 g/L KH_2PO_4 , pH 7.4.)

HPLC VARIABLES

Column: 150 \times 4.6 5 μm Supelcosil LC 18

Mobile phase: 10 mM KH_2PO_4 adjusted to pH 6.8 with 25% KOH

Flow rate: 1

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: k' 0.9

Internal standard: flucytosine

OTHER SUBSTANCES

Simultaneous: vidarabine phosphate

KEY WORDS

injections; water; flucytosine is IS

REFERENCE

Kwee, M.S.L.; Stolk, L.M.L. Formulation of a stable vidarabine phosphate injection, *Pharm. Weekbl. [Sci.]*, **1984**, 6, 101–104.

SAMPLE

Matrix: formulations

Sample preparation: Tablets. Powder tablets, weigh out amount equivalent to 5 mg flucytosine, add 50 mL water, stir for 15 min, filter. 1 mL Filtrate + 1.5 mL 27 $\mu\text{g}/\text{mL}$ thymine in buffer, make up to 10 mL with buffer, inject an aliquot. Injections. Dilute with water to a drug concentration of 100 $\mu\text{g}/\text{mL}$. 1 mL Sample + 1.5 mL 27 $\mu\text{g}/\text{mL}$ thymine in buffer, make up to 10 mL with buffer, inject an aliquot. (Buffer was 70 mM KH_2PO_4 adjusted to pH 3.0 with HCl.)

HPLC VARIABLES

Column: 250 \times 4.6 5 μm Spherisorb CN

Mobile phase: 9 mM Sodium heptanesulfonate adjusted to pH 2.8 with phosphoric acid

Flow rate: 0.8

Injection volume: 20

Detector: UV 266

CHROMATOGRAM

Retention time: 6.3

Internal standard: thymine (4.85)

OTHER SUBSTANCES

Simultaneous: 5-fluorouracil

KEY WORDS

tablets; injections

REFERENCE

Cavrini, V.; Bonazzi, D.; Di Pietra, A.M. Analysis of flucytosine dosage forms by derivative UV spectroscopy and liquid chromatography, *J. Pharm. Biomed. Anal.*, **1991**, 9, 401–407.

SAMPLE

Matrix: formulations

Sample preparation: Dilute 60 μL to 10 mL with 40 $\mu\text{g}/\text{mL}$ p-aminobenzoic acid in MeCN: water 25:75, inject a 5 μL aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 10 μm $\mu\text{Bondapak C18}$

Mobile phase: MeCN:buffer 25:75 (Buffer was 2 mL glacial acetic acid and 700 mg 1-octanesulfonic acid in 750 mL water.)

Flow rate: 1

Injection volume: 5

Detector: UV 285

CHROMATOGRAM

Retention time: 3.2

Internal standard: p-aminobenzoic acid (4.7)

OTHER SUBSTANCES

Noninterfering: degradation products

KEY WORDS

stability-indicating; oral liquids

REFERENCE

Wintermeyer,S.M.; Nahata,M.C. Stability of flucytosine in an extemporaneously compounded oral liquid, *Antimicrob.Agents Chemother.*, **1996**, 40, 407-409.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Guard column: 10 μ m precolumn (Beckman Instruments Inc.)

Column: 150 \times 4.6 5 μ m C18 Altex Ultrasphere ODS

Mobile phase: Isocratic. MeOH:buffer 15:85 containing 2.5 mM sodium pentanesulfonate and 2.5 mM sodium heptanesulfonate. Gradient. A was MeOH containing 2.5 mM sodium pentanesulfonate and 2.5 mM sodium heptanesulfonate. B was buffer containing 2.5 mM sodium pentanesulfonate and 2.5 mM sodium heptanesulfonate. A:B 0:100 for 2 min, to 30:70 over 1 min, maintain at 30:70. (Buffer was 50 mM phosphoric acid containing 50 mM KH_2PO_4 , pH 2.5.)

Flow rate: 1

Injection volume: 20

Detector: UV 254; UV 285

CHROMATOGRAM

Retention time: 3.76 (isocratic), 9.20 (gradient)

Internal standard: 5-methylcytosine (5.66 (isocratic), 12.62 (gradient))

Limit of quantitation: 12.5 μ M

OTHER SUBSTANCES

Simultaneous: barbituric acid, cytosine, fluorouracil, hydroxycytosine, uracil, urea

REFERENCE

Biondi,L.; Nairn,J.G. High performance liquid chromatographic assay for 5-fluorouracil and 5-fluorocytosine, *J.Liq.Chromatogr.*, **1985**, 8, 1881-1892.

SAMPLE

Matrix: solutions

Sample preparation: Centrifuge and filter cell solutions (0.22 μ m), inject an aliquot.

HPLC VARIABLES

Guard column: Guard-PAK C18 (Waters)

Column: 150 \times 3.9 5 μ m NOVA PAK C18

Mobile phase: MeOH:50 mM pH 6.0 KH_2PO_4 3:97

Flow rate: 1

Detector: UV 280

CHROMATOGRAM

Retention time: 2.5

REFERENCE

Koga,H. High-performance liquid chromatography measurement of antimicrobial concentrations in polymorphonuclear leukocytes, *Antimicrob.Agents Chemother.*, **1987**, 31, 1904–1908.

SAMPLE

Matrix: solutions

Sample preparation: Inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Spherisorb 5 ODS 2

Mobile phase: MeCN:50 mM pH 2.5 phosphate buffer 0.5:99.5

Column temperature: 25

Flow rate: 1

Injection volume: 20

Detector: UV 200

CHROMATOGRAM

Retention time: 4

OTHER SUBSTANCES

Simultaneous: 5-fluorouracil, floxuridine

REFERENCE

Del Nozal,M.J.; Bernal,J.L.; Pampliega,A.; Marinero,P.; Pozuelo,M. Determination of the concentrations of 5-fluorouracil and its metabolites in rabbit plasma and tissues by high-performance liquid chromatography, *J.Chromatogr.B*, **1994**, 656, 397–405.

Fludarabine


Molecular formula: C₁₀H₁₂FN₅O₄

Molecular weight: 285.23

CAS Registry No.: 21679-14-1, 75607-67-9 (phosphate)

Merck Index: 4162

Lednicer No.: 4 167

**SAMPLE**

Matrix: blood

Sample preparation: Condition a 500 mg Bakerbond SPE cartridge packed with LiChrosorb RP-18 with 2 mL MeOH, 2 mL water and 1 mL MeOH. Mix 500 μ L plasma with 100 μ L 30 μ g/mL mercaptopurine in MeOH, make up to 3 mL with MeOH. Centrifuge at 1100 g for 15 min, add 1.5 mL supernatant to the SPE cartridge, elute at 50 μ L/min flow-rate, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 200 \times 4 7 μ m Lichrosorb RP-18

Mobile phase: MeOH: pH 4.15 phosphate buffer 20:80

Flow rate: 1

Injection volume: 20

Detector: UV 262

CHROMATOGRAM

Retention time: 6.30

Internal standard: mercaptopurine (3.30)

Limit of detection: 50 ng/mL

Limit of quantitation: 100 ng/mL

KEY WORDS

plasma; SPE

REFERENCE

Misztal, G.; Paw, B. Determination of fludarabine phosphate in human plasma using reversed phase high-performance liquid chromatography, *Pharmazie*, **1996**, *51*, 733–734.

SAMPLE**Matrix:** blood

Sample preparation: Isolate mononuclear cells from 10 mL blood by a standard step-gradient density centrifugation procedure. Wash once with PBS, resuspend in 500 μ L water, add 500 μ L 800 mM perchloric acid, centrifuge at 400 g for 5 min, wash the pellet with 500 μ L 400 mM perchloric acid, centrifuge at 400 g for 5 min. Combine supernatants, neutralize with 10 M KOH, bring to pH 7 with 1 M KOH (Universal indicator paper), cool in ice, centrifuge at 400 g for 5 min, inject a 50–2000 μ L aliquot of the supernatant. (PBS was 8.1 g NaCl, 0.22 g KCl, 1.14 g NaHPO_4 (sic), 0.27 g KH_2PO_4 in 1 L water, pH 7.4.)

HPLC VARIABLES**Column:** 250 \times 4.6 10 μ m Partisil 10 SAX

Mobile phase: Gradient. A was 5 mM pH 2.8 $(\text{NH}_4)_2\text{HPO}_4$. B was 750 mM pH 3.5 $(\text{NH}_4)_2\text{HPO}_4$. A:B from 70:30 to 0:100 over 30 min (concave gradient, Waters no. 9). (At the start of each day pump through 20 mL 2 M $(\text{NH}_4)_2\text{HPO}_4$ then inject 100 μ L 100 mM disodium EDTA into the initial mobile phase.)

Flow rate: 3**Injection volume:** 50–2000**Detector:** UV 262**CHROMATOGRAM****Retention time:** 29 (as triphosphate, Fara-ATP)**OTHER SUBSTANCES****Extracted:** ara-CTP (cytarabine triphosphate), ATP, CTP, UTP, GTP**KEY WORDS**

mononuclear cells

REFERENCE

Gandhi, V.; Danhauser, L.; Plunkett, W. Separation of 1- β -D-arabinofuranosylcytosine 5'-triphosphate and 9- β -D-arabinofuranosyl-2-fluoroadenine 5'-triphosphate in human leukemia cells by high-performance liquid chromatography, *J. Chromatogr.*, **1987**, *413*, 293–299.

SAMPLE**Matrix:** blood, urine

Sample preparation: Plasma. Filter (Amicon CF) while centrifuging at 1000 g for 20 min, filter (0.45 μ m) the ultrafiltrate, inject an aliquot. Urine. 250 μ L Urine + 500 μ L 150 mM $\text{Ba}(\text{OH})_2$, vortex, add 500 μ L 5% $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$, vortex, centrifuge at 1700 g for 15 min, filter (0.45 μ m), inject an aliquot.

HPLC VARIABLES**Column:** 250 \times 4.6 5 μ m Ultrasphere C18**Mobile phase:** MeOH:10 mM pH 4.15 $(\text{NH}_4)_2\text{HPO}_4$ 6:94**Flow rate:** 1–1.2**Detector:** UV 254**CHROMATOGRAM****Retention time:** 55**Limit of detection:** 12.5 pmole**OTHER SUBSTANCES****Extracted:** fludarabine phosphate

KEY WORDS

plasma; pharmacokinetics; ultrafiltrate

REFERENCE

Hersh,M.R.; Kuhn,J.G.; Phillips,J.L.; Clark,G.; Ludden,T.M.; von Hoff,D.D. Pharmacokinetic study of fludara-bine phosphate (NSC 312887), *Cancer Chemother.Pharmacol.*, **1986**, 17, 277–280.

SAMPLE

Matrix: dialysate

HPLC VARIABLES

Column: 80 × 4.6 C18 (Perkin-Elmer)

Mobile phase: MeOH:10 mM pH 7 KH₂PO₄ 20:80

Flow rate: 1

Detector: UV 265

CHROMATOGRAM

Retention time: 3

OTHER SUBSTANCES

Simultaneous: araA, adenosine, 2-chloro-2'-arabino-fluoro-2'-deoxyadenosine, 2-chloro-2'-deoxy-adenosine, 2-chloroadenosine, 5'-chloro-5'-deoxyadenosine, 2'-deoxyadenosine

REFERENCE

Reichelova,V.; Liliemark,J.; Albertioni,F. Structure-activity relationships of 2-chloro-2'-arabino-fluoro-2'-deoxyadenosine and related analogues: Protein binding, lipophilicity, and retention in reversed-phase LC, *J.Liq.Chromatogr.*, **1995**, 18, 1123–1135.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 80 × 4.6 3 μm (Perkin-Elmer)

Mobile phase: MeOH:10 mM pH 6.8 potassium phosphate buffer 20:80

Flow rate: 1

Detector: UV 265

CHROMATOGRAM

Retention time: 3

OTHER SUBSTANCES

Simultaneous: cladribine, analogs, degradation products

REFERENCE

Reichelova,V.; Liliemark,J.; Albertioni,F. Liquid chromatographic study of acid stability of 2-chloro-2'-arabino-fluoro-2'-deoxyadenosine, 2-chloro-2'-deoxyadenosine and related analogues, *J.Pharm.Biomed.Anal.*, **1995**, 13, 711–714.


Fludeoxyglucose F18

Molecular formula: C₆H₁₁¹⁸FO₅

Molecular weight: 181.15

CAS Registry No.: 105851-17-0

Merck Index: 4163

**SAMPLE**

Matrix: cell incubations

Sample preparation: Wash cells with PBS, freeze in liquid nitrogen, add 3 mL 400 mM perchloric acid, let stand at 4° for 30 min, centrifuge at 12000 g at 4°, extract with 3 mL 500 mM trioctylamine in 1,1,2-trichlorotrifluoroethane, filter (0.22 µm) the inorganic phase, inject an aliquot.

HPLC VARIABLES

Guard column: LiChrosorb RP 18-5

Column: 100 × 8 Radial-PAK Partisil 10 SAX (Waters)

Mobile phase: Gradient. A was MeOH:15 mM pH 3.8 (NH₄)H₂PO₄ 3:97. B was MeOH:750 mM pH 4.8 (NH₄)H₂PO₄ 3:97. A:B from 100:0 to 0:100 over 35 min (concave gradient), re-equilibrate for 15 min.

Flow rate: 2

Injection volume: 150

Detector: Radioactivity (Berthold LB 506 C-1) (with Quick Szint Flow 306 pumped at 4 mL/min)

CHROMATOGRAM

Retention time: 3

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

C14 labeled; chondrocytes

REFERENCE

Fedders,G.; Kock,R.; Van de Leur,E.; Greiling,H. A radiochemical high-performance liquid chromatographic method for the analysis of 2-fluoro-2-deoxy-D-glucose-derived metabolites in human chondrocytes, *Anal.Biochem.*, **1993**, *211*, 81–86.

SAMPLE

Matrix: reaction mixtures

Sample preparation: Evaporate reaction mixture in the presence of active charcoal and pass through a 100 × 10 column of Merck Kieselgel 60 for dry column chromatography using 70 mL MeCN:water 99.7:0.3, evaporate the eluate, take up the residue in 1 mL water, inject an aliquot.

HPLC VARIABLES

Guard column: 30 × 3.6 Aminex HPX-87C (Bio-Rad)

Column: 300 × 7.8 Aminex HPX-87P (Bio-Rad) (The column contains lead which is washed off. Periodically backflush with 100 mM lead nitrate at 0.1 mL/min overnight.)

Mobile phase: Water

Column temperature: 20

Flow rate: 0.1

Detector: RI (at 18°)

CHROMATOGRAM

Retention time: 69.4

KEY WORDS

F18 labeled

REFERENCE

Oberdorfer,F.; Hull,W.E.; Traving,B.C.; Maier-Borst,W. Synthesis and purification of 2-deoxy-2-[18F]fluoro-D-glucose and 2-deoxy-2-[18F]fluoro-D-mannose: characterization of products by 1H- and 19F-NMR spectroscopy, *Int.J.Rad.Appl.Instrum.[A]*, **1986**, *37*, 695–701.

SAMPLE

Matrix: tissue hydrolyzate

Sample preparation: Centrifuge at 12000 g for 10 min, inject an aliquot.

HPLC VARIABLES**Guard column:** 20 × 4 APS Hypersil**Column:** 250 × 4 BioSil Amino 5S (Bio-Rad)**Mobile phase:** MeCN:water 90:10**Flow rate:** 1**Detector:** Radioactivity (with Quick Szint Flow 302)

CHROMATOGRAM**Retention time:** 7


KEY WORDS

C14 labeled

REFERENCE

Fedders,G.; Kock,R.; Van de Leur,E.; Greiling,H. The metabolism of 2-fluoro-2-deoxy-D-glucose in human chondrocytes and its incorporation into keratan sulfate proteoglycans, *Eur.J.Biochem.*, **1994**, 219, 1063–1071.

Fludrocortisone

Molecular formula: C₂₁H₂₉FO₅**Molecular weight:** 380.46**CAS Registry No.:** 127-31-1, 514-36-3 (acetate)**Merck Index:** 4166**Lednicer No.:** 1 192

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 250 × 4.6 5 μm Spherex C18 (Phenomenex USA)**Mobile phase:** MeOH:THF:water 3:25:72**Flow rate:** 1.0**Injection volume:** 60**Detector:** UV 254

CHROMATOGRAM**Retention time:** 15.9

OTHER SUBSTANCES

Simultaneous: 11-deoxycortisol, dexamethasone, hydrocortisone, methylprednisolone, prednisolone

REFERENCE

McWhinney,B.C.; Ward,G.; Hickman,P.E. Improved HPLC method for simultaneous analysis of cortisol, 11-deoxycortisol, prednisolone, methylprednisolone, and dexamethasone in serum and urine, *Clin.Chem.*, **1996**, 42, 979–981.

SAMPLE**Matrix:** solutions**Sample preparation:** Inject 20 μL aliquot of a MeOH solution.

HPLC VARIABLES**Column:** 250 × 4.6 5 μm Hypersil 5-ODS**Mobile phase:** THF:water 23:77**Column temperature:** 30**Flow rate:** 1

Injection volume: 20

Detector: UV 245

CHROMATOGRAM

Retention time: k' 7.87 (fludrocortisone), k' 28.04 (fludrocortisone acetate)

Internal standard: methylprednisolone (k' 11.36)

OTHER SUBSTANCES

Simultaneous: metabolites, betamethasone, corticosterone, cortisone, deflazacort, deoxycorticosterone, dexamethasone, fluorocortisone, fluorocortisone acetate, hydrocortisone, 21-hydroxydeflazacort, 11 α -hydroxyprogesterone, methylprednisolone, prednisolone, prednisone, triamcinolone acetonide, triamcinolone

REFERENCE

Santos-Montes,A.; Gonzalo-Lumbreras,R.; Gasco-Lopez,A.I.; Izquierdo-Hornillos,R. Extraction and high-performance liquid chromatographic separation of deflazacort and its metabolite 21-hydroxydeflazacort. Application to urine samples, *J.Chromatogr.B*, **1994**, 657, 248–253.

SAMPLE

Matrix: urine

Sample preparation: Equilibrate a Sephadex G-25M column with 100 mM pH 7.0 phosphate buffer. Condition a Bond-Elut C18 SPE cartridge with 1 mL MeCN, 4 mL acetone:water 20:80, and 4 mL water. 2 mL Urine + 500 μ L 500 mM pH 5.0 acetate buffer + 50 μ L MeOH + 160 μ L 100000 Fishmann U/mL β -glucuronidase and 800000 Roy U/mL arylsulfatase (from Helix pomatia, Boehringer Mannheim), heat at 37° for 24 h, filter (0.45 μ m), add to the Sephadex column, wash with three 2 mL portions of 100 mM pH 7.0 phosphate buffer, elute with four 2 mL portions of 100 mM pH 7.0 phosphate buffer. Add the eluate to the SPE cartridge, wash with 4 mL water, wash with 4 mL acetone:water 20:80, elute with 1 mL MeCN. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute the residue in 100 μ L MeOH, add 20 μ L cupric acetate solution, let stand at room temperature for 1 h, add 100 μ L reagent, heat at 60° for 40 min, cool, centrifuge briefly at 1000 g, inject a 100 μ L aliquot of the supernatant. (Cupric acetate solution was 0.7 g cupric acetate in water diluted to 100 mL with MeOH. Reagent was 7 mM 1,2-diamino-4,5-methylenedioxybenzene in water containing 200 mM β -mercaptoethanol and 250 mM sodium hydrosulfite, store in the dark at 4°, stable for at least 2 weeks. Prepare 1,2-diamino-4,5-methylenedioxybenzene as follows. Add 5 g 1,2-(methylenedioxy)-4-nitrobenzene to 37.5 mL concentrated nitric acid and 12.5 mL glacial acetic acid, pour the yellow-colored solution into water, recrystallize the 1,2-dinitro-4,5-methylenedioxybenzene from EtOH (Rec.Trav.Chim.Pays-Bas 1930, 49, 45). Dissolve 5 g 1,2-dinitro-4,5-methylenedioxybenzene in 200 mL benzene (Caution! Benzene is a carcinogen!), add 100 g 80 mesh iron powder, add 20 mL concentrated HCl in small portions over 1 h while heating the mixture under reflux. Reflux for 4 h, add 10 mL water, reflux for 2 h, cool, make alkaline with 2.6 M NaOH, extract three times with 200 mL portions of benzene. Combine the extracts, evaporate to dryness to give 1,2-diamino-4,5-methylenedioxybenzene, mix with 10 mL concentrated HCl, recrystallize from EtOH to give 1,2-diamino-4,5-methylenedioxybenzene dihydrochloride, mp 176-9° (Chem.Pharm.Bull. 1987, 35, 687).)

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m L-Column ODS (Chemicals Inspection and Testing Institute, Tokyo)

Mobile phase: MeOH:MeCN:500 mM ammonium acetate 50:10:40 (After each injection wash with MeOH:water 80:20 for 20 min for 20 min, re-equilibrate for 20 min.)

Flow rate: 1

Injection volume: 100

Detector: F ex 350 em 390

CHROMATOGRAM

Retention time: 35.6

Internal standard: fludrocortisone

OTHER SUBSTANCES

Extracted: hydrocortisone, tetrahydroaldosterone, aldosterone

Noninterfering: cortisone, corticosterone, hydroxycorticosteroids


KEY WORDS

SPE; derivatization; fludrocortisone is IS

REFERENCE

Yoshitake,T.; Ishida,J.; Sonezaki,S.; Yamaguchi,M. High performance liquid chromatographic determination of 3 α , 5 β -tetrahydroaldosterone and cortisol in human urine with fluorescence detection, *Biomed. Chromatogr.*, **1992**, 6, 217–221.

Flufenamic acid



Molecular formula: C₁₄H₁₀F₃NO₂

Molecular weight: 281.23

CAS Registry No.: 530-78-9

Merck Index: 4167

Lednicer No.: 1 110

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 40 μ L 100 μ g/mL mefenamic acid in MeOH + 1 mL 1 M hydrochloric acid + 6 mL dichloromethane, shake mechanically for 20 min. Centrifuge at 2000 g for 5 min, evaporate organic phase to dryness under a gentle stream of nitrogen at 40°. Reconstitute the residue with 50 μ L MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 5 μ m Nucleosil C18

Mobile phase: MeOH:water 77:23

Flow rate: 0.8

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 7.2

Internal standard: mefenamic acid (8.6)

Limit of detection: 100 ng/mL

Limit of quantitation: 500 ng/mL

KEY WORDS

plasma; rat; pharmacokinetics

REFERENCE

Cerretani,D.; Micheli,L.; Fiaschi,L.; Giorgi,G. High-performance liquid chromatography of flufenamic acid in rat plasma, *J.Chromatogr.B*, **1996**, 678, 365–368.

SAMPLE

Matrix: formulations

Sample preparation: 100-300 mg Gel ointment + 3 mL MeOH, mix vigorously, filter (0.2 μ m), inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: MeOH

Flow rate: 1

Injection volume: 10

Detector: UV 280

CHROMATOGRAM

Retention time: 4.3

KEY WORDS

ointment

REFERENCE

Yamamura, K.; Yamada, J.-I.; Yotsuyanagi, T. High-performance liquid chromatographic assay of antiinflammatory drugs incorporated in gel ointments. Separation and stability testing, *J.Chromatogr.*, **1985**, *331*, 383–388.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitrityline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenpropofen, fenproporex, fentanyl, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiacol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarboxystiril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazin-dol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephentyoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazepam, methocarbamol, methoxamine, methsuximide, methyl salicylate, methylidopa, methyldopamine, methylphenidate, methylprednisolone, methyltestosterone, methyprylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, nor-epinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphen-butazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, per-santine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenyl-butazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primi-done, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sufadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sul-fasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tol-metin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapa-mil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233–242.

SAMPLE

Matrix: solutions

Sample preparation: Inject a 20 μL aliquot of a 100–500 $\mu\text{g/mL}$ solution in mobile phase.

HPLC VARIABLES

Column: 100 \times 4.6 5 μm Hypersil C8 MOS 100A coated with phosphatidylcholine (95% pure soybean lecithin, Epikuron, Lucas Meyer & Co.) (Coat column by recycling a 1 mM solution of phosphatidylcholine in MeOH:water 80:20 for 24 h.)

Mobile phase: MeCN:35 mM pH 7.4 sodium phosphate buffer 40:60

Flow rate: 0.5–2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: k' 1.35

OTHER SUBSTANCES

Also analyzed: amoxicillin, antipyrine, carbamazepine, chlorpheniramine, chlorpromazine, clonidine, codeine, desipramine, diphenhydramine, dipyridamole, ephedrine, haloperidol, hydroxyzine, imipramine, indomethacin, lidocaine, megestrol acetate, metoprolol, nabumetone, nadolol, phenobarbital, phenol, promazine, propranolol, pyrilamine, quinidine, ropinirole, testosterone, thioridazine, tolfenamic acid, verapamil

Noninterfering: acetaminophen, aspirin, azathioprine, caffeine, carprofen, chlorambucil, cimetidine, fenoterol, flurbiprofen, ibuprofen, ketoprofen, ranitidine, salicylic acid, sulfamethoxazole, theophylline, thioguanine, tiaprofenic acid, trimethoprim, valproic acid


KEY WORDS

comparison with capillary electrophoresis

REFERENCE

Hanna,M.; de Biasi,V.; Bond,B.; Salter,C.; Hutt,A.J.; Camilleri,P. Estimation of the partitioning characteristics of drugs: A comparison of a large and diverse drug series utilizing chromatographic and electrophoretic methodology, *Anal.Chem.*, **1998**, *70*, 2092–2099.

Fluindione



Molecular formula: $\text{C}_{15}\text{H}_9\text{FO}_2$

Molecular weight: 240.23

CAS Registry No.: 957-56-2

Merck Index: 4168

SAMPLE

Matrix: blood

Sample preparation: Add 100 μL plasma to a siliconized tube containing 50 μL 2 $\mu\text{g/mL}$ IS in MeOH and 100 μL MeCN. Vortex for 30 s, centrifuge at 3000 g for 10 min. Mix a 150 μL aliquot of the supernatant with 150 μL 0.9% NaCl, inject a 100 μL aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μm Supelcosil LC-18

Mobile phase: MeCN:67 mM Na_2HPO_4 buffer adjusted to pH 7.2 with orthophosphoric acid 23:77 (At the end of each chromatographic session rinse the HPLC system with 200 mL MeCN: water 50:50.)

Flow rate: 1.5

Injection volume: 100

Detector: UV 280

CHROMATOGRAM

Retention time: 4.2

Internal standard: coumarin (8.5)

Limit of quantitation: 25 ng/mL

OTHER SUBSTANCES

Also analyzed: acenocoumarol, brodifacoum, bromadiolone, chlorphacinone, coumatetralyl, difenacoum, diphenadione, ethyl biscoumacetate, phenidione, warfarin

KEY WORDS

plasma

REFERENCE

Aymard,G.; Legrand,M.; Comets,E.; Mentre,F.; Diquet,B. Rapid and simple micromethod for the quantification of fluindione in human plasma using high-performance liquid chromatography, *J.Chromatogr.B*, **1998**, 707, 169-173.

Flumazenil


Molecular formula: C₁₅H₁₄FN₃O₃

Molecular weight: 303.29

CAS Registry No.: 78755-81-4

Merck Index: 4169

Lednicer No.: 4 220



SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 25 µL EtOH + 25 µL 3 µg/mL IS1 and 7.6 µg/mL IS2 in EtOH + 1 mL 100 mM Na₂HPO₄ adjusted to pH 10.5 with NaOH + 5 mL diethyl ether: dichloromethane 60:40, vortex for 30 s, centrifuge at 4° at 2000 g for 10 min. Remove the organic phase and add it to 1 mL 100 mM Na₂HPO₄ adjusted to pH 10.5 with NaOH, vortex for 30 s, centrifuge at 4° at 2000 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 50 µL mobile phase, inject a 1-15 µL aliquot.

HPLC VARIABLES

Column: 100 × 4.6 3 µm CP-Microspher C18 (Chrompack)

Mobile phase: Gradient. A was MeOH:buffer 1:2. B was MeOH:water 80:20. A:B 93.8:6.2 for 5.5 min, to 60:40 over 0.15 min, maintain at 60:40 for 11.3 min, to 2.5:97.5 over 0.5 min, maintain at 2.5:97.5 for 3.5 min, return to initial conditions over 0.5 min (Buffer was 6 g/L NaH₂PO₄ and 1 mL/L triethylamine adjusted to pH 7.00 with NaOH.)

Column temperature: 40

Flow rate: 1.5

Injection volume: 1-15

Detector: UV 220

CHROMATOGRAM

Retention time: 5.0

Internal standard: IS1 ethyl 7-chloro-5,6-dihydro-5-methyl-6-oxo-4H-imidazo[1,5-a][1,4]benzodiazepine-3-carboxylate (Ro 15-305) (6.4), IS2 clomazepam (17.0)

Limit of detection: 0.3 ng/mL

OTHER SUBSTANCES

Extracted: midazolam

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Vletter, A.A.; Burm, A.G.L.; Breimer, L.T.M.; Spierdijk, J. High-performance liquid chromatographic assay to determine midazolam and flumazenil simultaneously in human plasma, *J. Chromatogr.*, **1990**, 530, 177–185.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Serum or plasma + 1 mL saturated Na_2HPO_4 + 100 μL 10 $\mu\text{g/mL}$ IS + 5 mL diisopropyl ether:isopropanol 95:5 (Caution! Diisopropyl ether readily forms explosive peroxides!), shake for 5 min, centrifuge. Remove the organic layer and evaporate it to dryness at 35°, reconstitute the residue in 50 μL mobile phase, vortex, inject a 20 μL aliquot.

HPLC VARIABLES

Column: 100 \times 3 Chrompack CP spher C8

Mobile phase: MeCN:50 mM NaH_2PO_4 70:30 adjusted to pH 2.2 with 85% orthophosphoric acid

Flow rate: 0.8

Injection volume: 20

Detector: UV 210

CHROMATOGRAM

Retention time: 2.1

Internal standard: N-6-dimethyl-2-(4-methylphenyl)-N-propylimidazo[1,2- α]pyridine-3-acetamide methanesulfonate (Synthelabo France) (6.4)

OTHER SUBSTANCES

Extracted: prothipendyl, zolpidem

KEY WORDS

serum; plasma; pharmacokinetics

REFERENCE

Debailleul, G.; Khalil, F.A.; Lheureux, P. HPLC quantification of zolpidem and prothipendyl in a voluntary intoxication, *J. Anal. Toxicol.*, **1991**, 15, 35–37.

SAMPLE

Matrix: blood

Sample preparation: Condition a Sep-Pak C18 SPE cartridge with water, MeOH, and 100 mM ammonium acetate. Add 200 μL plasma to the SPE cartridge, wash with 100 mM ammonium acetate, elute with MeOH:100 mM ammonium acetate 3:1. Evaporate the eluate to dryness under reduced pressure, dissolve the residue in 200 μL mobile phase, inject a 20 μL aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 Hitachi gel 3056 octadecylsilica

Mobile phase: MeOH:100 mM ammonium acetate 60:40

Flow rate: 1

Injection volume: 20

Detector: MS, Hitachi M1000, APCI, nebulizer 260°, vaporizer 399

CHROMATOGRAM

Retention time: 2.7

Limit of detection: 0.5–2.5 ng/mL

OTHER SUBSTANCES

Simultaneous: atipamezole, atropine, butorphanol, ketamine, medetomidine, midazolam, xylazine

KEY WORDS

plasma; SPE; dog

REFERENCE

Kanazawa,H.; Nagata,Y.; Matsushima,Y.; Takai,N.; Uchiyama,H.; Nishimura,R.; Takeuchi,A. Liquid chromatography-mass spectrometry for the determination of medetomidine and other anaesthetics in plasma, *J.Chromatogr.*, **1993**, *631*, 215–220.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 244

CHROMATOGRAM

Retention time: 3.26

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfapyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vandesine; mexiletine; dipyrindamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrridine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; flotaefenine; alpidem; loratadine; chlorpromazine; clomipramine; caripramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254–262.

SAMPLE

Matrix: blood, urine

Sample preparation: 500 μ L Plasma or urine + 50 μ L 4 μ g/mL flurazepam in MeOH + 1 mL 100 mM pH 9 sodium phosphate buffer + 4 mL dichloromethane:diethyl ether 60:40, shake at 45 rpm for 15 min, centrifuge at 10° at 1870 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 80 μ L MeOH, inject a 30 μ L aliquot. (Deconjugate urine as follows. 250 μ L Urine + 750 μ L pH 5.4 acetate buffer + 500 U β -glucuronidase, heat at 37° for 18 h, add 20 μ L 5 M NaOH, centrifuge, proceed as above using 5 mL dichloromethane:diethyl ether.)

HPLC VARIABLES

Guard column: 30 \times 4.6 30 μ m C8

Column: 100 \times 8 4 μ m Nova Pak C18

Mobile phase: MeCN:40 mM sodium phosphate buffer 32:68 containing 1 mL/L triethylamine, final pH 7.2

Flow rate: 1.5

Injection volume: 30

Detector: UV 220

CHROMATOGRAM

Retention time: 3.8

Internal standard: flurazepam (18.5)

Limit of detection: 4 ng/mL

OTHER SUBSTANCES

Extracted: midazolam, 4-hydroxymidazolam, 1-hydroxymethylmidazolam

Noninterfering: alfentanil, atropine, bupivacaine, lignocaine, neostigmine

KEY WORDS

plasma

REFERENCE

Chan,K.; Jones,R.D.M. Simultaneous determination of flumazenil, midazolam and metabolites in human biological fluids by liquid chromatography, *J.Chromatogr.*, **1993**, *619*, 154–160.

SAMPLE

Matrix: formulations

Sample preparation: Inject a 25 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Spherisorb C18

Mobile phase: MeOH:THF:isopropanol:water 30:3.5:1.5:65

Column temperature: 40

Flow rate: 1

Injection volume: 25

Detector: UV 245

CHROMATOGRAM

Retention time: 6.30

OTHER SUBSTANCES

Simultaneous: methylparaben, propylparaben, phenol, degradation products

Noninterfering: aminophylline, cimetidine, dobutamine, dopamine, famotidine, lidocaine, procainamide, ranitidine

KEY WORDS

injections; 5% dextrose; stability-indicating

REFERENCE

Olsen, K.M.; Gurley, B.J.; Davis, G.A.; Christensen, R.; Monaghan, M.S. Stability of flumazenil with selected drugs in 5% dextrose injection, *Am. J. Hosp. Pharm.*, **1993**, 50, 1907–1912.

Flumequine


Molecular formula: C₁₄H₁₂FNO₃

Molecular weight: 261.25

CAS Registry No.: 42835-25-6

Merck Index: 4172

Lednicer No.: 3 186



SAMPLE

Matrix: blood

Sample preparation: Mix 100 μ L plasma with 100 μ L pH 6.0 phosphate buffer and 1 mL ethyl acetate, shake, centrifuge at 20000 g for 5 min, collect 800 μ L the organic phase, dry under nitrogen at 45°. Add 300 μ L pH 7.8 phosphate buffer and 300 μ L hexane to the residue, shake, centrifuge at 20000 g for 5 min, inject a 100 μ L aliquot of the aqueous phase.

HPLC VARIABLES

Guard column: 4 \times 4 C18 (Merck)

Column: 125 \times 4 5 μ m Lichrospher RP Select B

Mobile phase: MeCN:DMF:water:orthophosphoric acid 18:28:40.5:13.5

Flow rate: 0.8

Injection volume: 100

Detector: UV 324 (A), F ex 320 em 365 (B)

CHROMATOGRAM

Retention time: 6.4

Limit of quantitation: 5 ng/mL (F), 100 ng/mL (UV)

OTHER SUBSTANCES

Extracted: metabolites

Noninterfering: ciprofloxacin, danofloxacin, enrofloxacin, marbofloxacin, nalidixic acid, oxolinic acid

KEY WORDS

sheep; plasma

REFERENCE

Delmas, J.M.; Chapel, A.M.; Sanders, P. Determination of flumequine and 7-hydroxyflumequine in plasma of sheep by high-performance liquid chromatography, *J. Chromatogr. B*, **1998**, 712, 263–268.

SAMPLE

Matrix: tissue

Sample preparation: 2 g Minced tissue + 100 μ L 5 μ g/mL IS in water + 4 g anhydrous sodium sulfate, mix until homogenized. Add 10 mL ethyl acetate, shake mechanically for 10 min, centrifuge at 1500 g for 10 min. Transfer organic layer into another tube and repeat extraction on the tissue pellet with 10 mL ethyl acetate. Evaporate combined organic phases under a stream of nitrogen at 50°. Dissolve residue in 1 mL MeCN:2.7 mM pH 2.5 oxalic acid 50:50, vortex, sonicate for 5 min, filter through 0.45 μ m filter (GHP Acrodisc GF, Gelman Sciences, USA). Inject a 100 μ L aliquot. (1 mg/mL IS stock solution was prepared in 20 mM NaOH, it was diluted to working concentration with water immediately before use.)

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Ultrabase C18 (Shandon, UK)

Mobile phase: Gradient. A was MeCN. B was 2.7 mM pH 2.5 oxalic acid. A:B from 10:90 to 70:30 over 20 min, maintain at 70:30 for 5 min, return to initial conditions over 5 min.

Flow rate: 0.8

Injection volume: 100

Detector: F ex 252 em 356

CHROMATOGRAM

Retention time: 20

Internal standard: ibafloxacin (21.3)

Limit of detection: 15 ng/g

Limit of quantitation: 50 ng/g

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

kidney; pig

REFERENCE

Guyonnet,J.; Pacaud,M.; Richard,M.; Doisi,A.; Spavone,F.; Hellings,P. Routine determination of flumequine in kidney tissue of pig using automated liquid chromatography, *J.Chromatogr.B*, **1996**, 679, 177-184.

SAMPLE

Matrix: tissue

Sample preparation: Homogenize 1 g tissue with 5 mL acetone, centrifuge, decant and save the supernatant, repeat the extraction. Add 2 mL acetone, 3 mL hexane, and 6 mL 3% NaCl to the combined supernatants, extract, centrifuge, discard the hexane layer. Add the extract to 25 mL chloroform (Caution! Chloroform is a carcinogen!), mix, separate the phases, add 2.5 mL 100 mM pH 9.0 Na_3PO_4 buffer and one drop 1 M NaOH to the chloroform phase, mix, separate the phases, discard the chloroform layer. Wash the aqueous layer with 2.5 mL chloroform, centrifuge the aqueous layer. Dialyze a 740 μL aliquot of the supernatant (in 2 portions) against 3.9 mL 20 mM pH 5.0 sodium phosphate buffer pumped at 0.6 mL/min using a Cuprophan 15000 MW cut-off cellulose acetate membrane. After washing column A with 500 μL MeCN:water 50:50 and 500 μL 20 mM pH 5.0 Na_3PO_4 buffer the dialysate flowed through column A to waste. Wash column A with 500 μL 20 mM pH 5.0 sodium phosphate buffer, elute the contents of column A onto column B with mobile phase, elute with mobile phase, monitor the effluent from column B. (Wash the donor channel with 2 mL 0.01% Triton X-100 in 20 mM pH 5.0 sodium phosphate buffer. Wash the acceptor channel with 3 mL 20 mM pH 5.0 sodium phosphate buffer. Regenerate the membrane with 2 mL 100 mM pH 9.0 sodium phosphate buffer (donor channel) and 3 mL 20 mM pH 5.0 sodium phosphate buffer (acceptor channel).)

HPLC VARIABLES

Column: A 5.8×4.6 Hypersil ODS; B 150×4.6 PLRP-S (Polymer Labs., UK)

Mobile phase: MeCN:THF:20 mM pH 5.0 Na_3PO_4 buffer 20:15:65

Flow rate: 0.6

Detector: F ex 318 em 364

CHROMATOGRAM

Retention time: 12

Limit of detection: 5 ng/g

OTHER SUBSTANCES

Extracted: oxolinic acid

KEY WORDS

column switching; chicken; liver; dialysis

REFERENCE

Eng,G.Y.; Maxwell,R.J.; Cohen,E.; Piotrowski,E.G.; Fiddler,W. Determination of flumequine and oxolinic acid in fortified chicken tissue using on-line dialysis and high-performance liquid chromatography with fluorescence detection, *J.Chromatogr.A*, **1998**, 799, 349-354.

Flumethasone


Molecular formula: $C_{22}H_{28}F_2O_5$

Molecular weight: 410.46

CAS Registry No.: 2135-17-3, 2002-29-1 (pivalate)

Merck Index: 4173

Lednicer No.: 1 200



SAMPLE

Matrix: blood

Sample preparation: Add 1 mL serum to a Sep Pak C18 SPE cartridge, wash with 4 mL water, elute with 4 mL MeOH, evaporate to dryness under vacuum, reconstitute in 50 μ L MeCN: water 30:70, inject whole sample.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Ultrasphere ODS

Mobile phase: MeCN:water 30:70

Flow rate: 1

Injection volume: 50

Detector: enzyme immunoassay of fractions

CHROMATOGRAM

Retention time: 18

Limit of detection: 0.3 pg

OTHER SUBSTANCES

Extracted: dexamethasone, betamethasone, triamcinolone

Noninterfering: endogenous steroids

KEY WORDS

serum; SPE; horse

REFERENCE

Friedrich,A.; Schulz,R.; Meyer,H.H. Use of enzyme immunoassay and reverse-phase high-performance liquid chromatography to detect and confirm identity of dexamethasone in equine blood, *Am.J.Vet.Res.*, **1992**, *53*, 2213-2220.

SAMPLE

Matrix: blood

Sample preparation: Condition a Tef Elutor C18 cartridge with two 3 mL portions of MeOH then two 3 mL portions of water. Heat 1 mL plasma at 50° for 10 min, add to cartridge, wash with 2 mL water, 1 mL MeOH:water 10:90, 4 mL acetone:water 20:80, apply suction to cartridge for 10 min to air dry. Elute with 1 mL MeOH, evaporate eluent at 45° under nitrogen, reconstitute with 50 μ L mobile phase, inject 25 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 2 3 μ m C18 Hypersil

Mobile phase: MeCN:THF:water 8:10:82, containing 5 mL/L triethylamine, pH adjusted to 6.5 with citric acid

Flow rate: 0.6

Injection volume: 25

Detector: UV 242

CHROMATOGRAM

Retention time: 11.50

Internal standard: flumethasone

Limit of detection: 300 pg/mL

OTHER SUBSTANCES

Simultaneous: dexamethasone, prednisone, hydrocortisone, adrenosterone, prednisolone, estriol, corticosterone, methylprednisolone, cortisone, hydroxyprogesterone, testosterone, deoxycorticosterone, fluorometholone, spironolactone, equilenin, estrone, estradiol, progesterone, diphenhydramine, propranolol, aspirin, theophylline, imipramine, desipramine, indomethacin, amitriptyline, nortriptyline, nordiazepam, diazepam, chlordiazepoxide, tripeleennamine, carbamazepine, probenecid, phenobarbital

Noninterfering: caffeine, nicotine, cotinine, chlorothiazide, acetazolamide, phenytoin, pheniramine, cephalothin, primidone, acebutolol, hydrochlorothiazide, quinine, acetophenetidin, furosemide, aldosterone, triamcinolone, ephedrine, allopurinol, phenylephrine

KEY WORDS

plasma; SPE; flumethasone is IS

REFERENCE

Hariharan,M.; Naga,S.; VanNoord,T.; Kindt,E.K. Simultaneous assay of corticosterone and cortisol in plasma by reversed-phase liquid chromatography, *Clin.Chem.*, **1992**, 38, 346–352.

SAMPLE

Matrix: blood

Sample preparation: Condition a 2 mL 200 mg Tef Elutor C18 SPE cartridge (Versa Prep) with 3 mL MeOH and two 3 mL portions of water. Heat 1 mL plasma at 50° for 10 min, add to SPE cartridge, wash with 2 mL water, wash with 1 mL MeOH:water 10:90, wash with 4 mL acetone:water 20:80, air-dry for 10 min, elute with 1 mL MeOH. Evaporate the eluate to dryness under a stream of nitrogen at 45°, reconstitute the residue in 50 µL mobile phase, inject a 25 µL aliquot.

HPLC VARIABLES

Column: 100 × 2.3 µm Hypersil

Mobile phase: MeCN:THF:water 8:10:82 containing 5 mL/L triethylamine, pH adjusted to 6.5 with citric acid

Flow rate: 0.6

Injection volume: 25

Detector: UV 242

CHROMATOGRAM

Retention time: 13

Internal standard: flumethasone

OTHER SUBSTANCES

Extracted: cortisone, corticosterone, hydrocortisone

Simultaneous: acebutolol, acetazolamide, acetophenetidin, adrenosterone, aldosterone, amitriptyline, androsten-3,17-dione, aspirin, carbamazepine, cephalothin, chlordiazepoxide, chlorothiazide, dehydrocorticosterone, deoxycorticosterone, deoxycortisol, desipramine, dexamethasone, diazepam, diphenhydramine, equilenin, estradiol, estriol, estrone, fluorometholone, furosemide, hydrochlorothiazide, hydroxycorticosterone, hydroxyprogesterone, hydroxyprogesterone, imipramine, indomethacin, methylhydroxyprogesterone, methylprednisolone, nandrolone, nordiazepam, nortriptyline, pheniramine, phenobarbital, phenytoin, prednisolone, prednisone, primidone, probenecid, progesterone, propranolol, quinine, spironolactone, testosterone, theophylline, triamcinolone, tripeleennamine

Noninterfering: caffeine, nicotine, cotinine, ephedrine, allopurinol, phenylephrine

KEY WORDS

serum; SPE; flumethasone is IS

REFERENCE

Hariharan,M.; Naga,S.; VanNoord,T.; Kindt,E.K. Assay of human plasma cortisone by liquid chromatography: normal plasma concentrations (between 8 and 10 a.m.) of cortisone and corticosterone, *J.Chromatogr.*, **1993**, 613, 195–201.

SAMPLE

Matrix: formulations

Sample preparation: Ointment. Weigh out 1-1.5 g ointment, add 10 mL chloroform:hexane 50:50, warm until the ointment dissolves, add to a silica Sep-Pak SPE cartridge, rinse container with 10 mL hexane:chloroform 50:50, add rinse to SPE cartridge, discard all eluates, rinse container with 5 mL MeOH, add rinse to SPE cartridge, elute with 15 mL MeOH, collect all MeOH eluates, make up to 25 mL with MeOH, inject a 25 μ L aliquot. Cream. Weigh out 1 g cream and add it to 10 mL MeOH, shake for 5 min, sonicate for 5 min, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 C8 (Brownlee)
Mobile phase: MeOH:THF:water 30:30:40
Flow rate: 1.5
Injection volume: 10-25
Detector: UV 254

CHROMATOGRAM

Retention time: 8.4 (flumethasone pivalate)

KEY WORDS

ointment; cream; SPE

REFERENCE

Lodge,B.A.; Vincent,A. Analysis of flumethasone pivalate formulations by high-performance liquid chromatography, *J.Chromatogr.*, **1984**, 301, 477-480.

SAMPLE

Matrix: liposomal preparations
Sample preparation: Dilute 1000-fold with MeOH/water, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: C18
Mobile phase: MeOH:1% acetic acid 70:30
Flow rate: 1
Injection volume: 20
Detector: UV 254

CHROMATOGRAM

Retention time: 5.12

OTHER SUBSTANCES

Interfering: dexamethasone

REFERENCE

Devoisselle,J.-M.; Vion-Dury,J.; Confort-Gouny,S.; Coustaut,D.; Cozzone,P.J. Liposomes containing fluorinated steroids: an analysis based on photon correlation and fluorine-19 nuclear magnetic resonance spectroscopy, *J.Pharm.Sci.*, **1992**, 81, 249-254.

SAMPLE

Matrix: solutions
Sample preparation: Condition a Bond Elut C18 SPE cartridge with 4 mL water then 3 mL MeOH. Add aqueous steroid solution to cartridge, elute with MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.5 μ m Nucleosil C18
Mobile phase: MeCN:water 70:30
Flow rate: 1
Injection volume: 20
Detector: UV 237

CHROMATOGRAM

Limit of detection: 120 ng/mL

OTHER SUBSTANCES

Also analyzed: betamethasone 17-valerate, dexamethasone

KEY WORDS

for flumethasone 21-acetate; SPE

REFERENCE

Valenta,C.; Janout,H. Corticosteroid analysis by HPLC with increased sensitivity by use of precolumn concentration, *J.Liq.Chromatogr.*, **1994**, *17*, 1141–1146.

SAMPLE

Matrix: synovial fluid

Sample preparation: 100 μ L Synovial fluid + 10 μ L 10 μ g/mL flumethasone in MeOH + 1 mL 100 mM NaOH + 10 mL dichloromethane, shake for 10 min, centrifuge at 8400 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100 μ L mobile phase, vortex, inject the whole amount.

HPLC VARIABLES

Column: 100 \times 8 10 μ m Radial Pak B silica (Waters)

Mobile phase: Dichloromethane:MeOH:glacial acetic acid 96.8:2.4:0.8

Flow rate: 1.4

Injection volume: 100

Detector: UV 254

CHROMATOGRAM

Retention time: 6

Internal standard: flumethasone

OTHER SUBSTANCES

Extracted: methylprednisolone, methylprednisolone acetate

KEY WORDS

normal phase; cow; flumethasone is IS

REFERENCE

Alvinerie,M.; Toutain,P.L. Determination of methylprednisolone and methylprednisolone acetate in synovial fluid using high-performance liquid chromatography, *J.Chromatogr.*, **1984**, *309*, 385–390.

Flunarizine


Molecular formula: C₂₆H₂₆F₂N₂

Molecular weight: 404.50

CAS Registry No.: 52468-60-7, 30484-77-6 (2.HCl)

Merck Index: 4179

Lednicer No.: 2 31



SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)


HPLC VARIABLES**Guard column:** 20 mm long Symmetry C18**Column:** 250 × 4.6 5 µm Symmetry C8 (Waters)**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.**Column temperature:** 30**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)**Injection volume:** 10-30**Detector:** UV 200.5**CHROMATOGRAM****Retention time:** 19.317**KEY WORDS**

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, *763*, 149–163.

Flunisolide

Molecular formula: C₂₄H₃₁FO₆**Molecular weight:** 434.50**CAS Registry No.:** 3385-03-3, 77326-96-6 (hemihydrate), 4533-89-5 (acetate)**Merck Index:** 4180**Lednicer No.:** 2 181**SAMPLE****Matrix:** blood**Sample preparation:** Add 30 µL IS to 1 mL serum, add 5 mL diethyl ether/dichloromethane (ratio not given), mix for 30 min. Centrifuge and freeze at -80° for 15 min. Evaporate the supernatant under a stream of nitrogen, reconstitute the residue in 200 µL mobile phase. Inject a 50 µL aliquot.**HPLC VARIABLES****Column:** Lichrospher RP Select B**Mobile phase:** MeCN:20 mM ammonium acetate buffer 80:20**Flow rate:** 1**Injection volume:** 50**Detector:** MS, PE-Sciex API 300, negative ion mode**CHROMATOGRAM****Internal standard:** budesonide**Limit of detection:** 100 pg/mL**KEY WORDS**

serum; pharmacokinetics

REFERENCE

Möllmann,H.; Derendorf,H.; Barth,J.; Meibohm,B.; Wagner,M.; Krieg,M.; Weisser,H.; Knöller,J.; Möllmann,A.; Hochhaus,G. Pharmacokinetic/pharmacodynamic evaluation of systemic effects of flunisolide after inhalation, *J.Clin.Pharmacol.*, **1997**, *37*, 893–903.

SAMPLE**Matrix:** formulations**Sample preparation:** Dissolve in mobile phase, inject an aliquot.

HPLC VARIABLES**Column:** 10 μ m Spherisorb ODS**Mobile phase:** MeCN:water:acetic acid 35:64:1**Injection volume:** 100**Detector:** UV 240

CHROMATOGRAM**Internal standard:** fluocinonide

KEY WORDS

nasal spray

REFERENCE

Yu,C.D.; Jones,R.E.; Henesian,M. Cascade impactor method for the droplet size characterization of a metered-dose nasal spray, *J.Pharm.Sci.*, **1984**, 73, 344–348.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 300 \times 3.9 μ Bondapak C18**Mobile phase:** MeOH:water 55:45**Flow rate:** 2**Detector:** UV 254


CHROMATOGRAM**Retention time:** 10.2

OTHER SUBSTANCES**Simultaneous:** metabolites

REFERENCE

Tőkés,L.; Cho,D.; Maddox,M.L.; Chaplin,M.D.; Chu,N.I. Isolation and identification of an oxidatively defluorinated metabolite of flunitrazepam in man, *Drug Metab.Dispos.*, **1981**, 9, 485–486.

Flunitrazepam

Molecular formula: C₁₆H₁₂FN₃O₃**Molecular weight:** 313.29**CAS Registry No.:** 1622-62-4**Merck Index:** 4181**Lednicer No.:** 2 406

SAMPLE**Matrix:** blood**Sample preparation:** Vortex 1 mL plasma with 3 mL toluene:isoamyl alcohol 95:5 at 1000 rpm for 90 s, centrifuge at 2600 g for 10 min. Evaporate a 2.5 mL aliquot of the upper organic layer to dryness under nitrogen at 40°, reconstitute with 100 μ L mobile phase, inject a 25 μ L aliquot.

HPLC VARIABLES**Guard column:** 10 \times 4.6 5 μ m Nucleosil 120-5 C18