
OTHER SUBSTANCES

Extracted: morphine, codeine, hydromorphone, hydrocodone, caffeine, cotinine, benzoylecgonine, secobarbital, oxazepam, phenobarbital, nordiazepam, diazepam, phenylpropanolamine, phenetermine, amphetamine, phenmetrazine, lidocaine, ephedrine, pentazocine, methamphetamine, desipramine, nortriptyline, diphenhydramine, methadone

Interfering: imipramine, amitriptyline


KEY WORDS

column-switching

REFERENCE

Binder, S.R.; Regalia, M.; Biaggi-McEachern, M.; Mazhar, M. Automated liquid chromatographic analysis of drugs in urine by on-line sample cleanup and isocratic multi-column separation, *J. Chromatogr.*, **1989**, *473*, 325–341.

Flurbiprofen



Molecular formula: C₁₅H₁₃FO₂

Molecular weight: 244.27

CAS Registry No.: 5104-49-4

Merck Index: 4234

Lednicer No.: 1 86

SAMPLE

Matrix: aqueous humor

Sample preparation: 100 μ L Aqueous humor + 500 μ L MeCN + 30 μ L 400 ng/mL (+)-naproxen in MeOH, mix mechanically for 90 s, centrifuge at 3000 g for 20 min. Remove the supernatant and dry it under nitrogen at room temperature, dissolve the residue in 50 μ L mobile phase by swirl-mixing for 1 min, centrifuge at 3000 g for 20 s, reduce volume to 20–30 μ L, inject an aliquot.

HPLC VARIABLES

Column: 150 \times 4.5 5 μ m Ultrasphere octyl

Mobile phase: MeCN:triethylamine:1.65% glacial acetic acid 505:0.65:495, pH 4.35

Column temperature: 30

Flow rate: 1

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 6.04

Internal standard: naproxen (3.89)

Limit of detection: 0.4 ng

OTHER SUBSTANCES

Extracted: diclofenac, indomethacin, meclofenamic acid

Simultaneous: bacitracin, cortisone acetate, diazepam, fluorometholone, hydrocortisone acetate, imipramine, ketoprofen, ketorolac tromethamine, levobunolol, metipranolol, neomycin, prednisolone acetate, proparacaine, propranolol, salicylic acid, sulfacetamide, suprofen

Noninterfering: acebutolol, acetaminophen, acetazolamide, alprenolol, apraclonidine, atenolol, atropine, betamethasone, betaxolol, bupivacaine, caffeine, cyclopentolate, dexamethasone, diphenhydramine, erythromycin, haloperidol, lidocaine, phenylephrine, polymyxin B, procaine, scopolamine, timolol, tropicamide

KEY WORDS

human; rabbit

REFERENCE

Riegel, M.; Ellis, P.P. High-performance liquid chromatography assay for antiinflammatory agents diclofenac and flurbiprofen in ocular fluids, *J. Chromatogr. B*, **1994**, *654*, 140–145.

SAMPLE

Matrix: aqueous humor

Sample preparation: 100 μ L Aqueous humor + 500 μ L MeCN + 30 μ L 400 ng/mL (+)-naproxen, mix mechanically for 90 s, centrifuge at 3000 g for 20 min. Remove supernatant and dry it under nitrogen at room temperature. Dissolve residue in 50 μ L mobile phase by swirl mixing for 1 min, centrifuge at 3000 g for 20 s. For concentrations of < 20 ng/mL, reduce volume to 20–30 μ L under nitrogen.

HPLC VARIABLES

Column: 150 \times 4.5 5 μ m Ultrasphere octyl

Mobile phase: 505 mL MeCN containing 0.65 mL triethylamine + 495 mL 1.65% glacial acetic acid, apparent pH 4.35

Column temperature: 30

Flow rate: 1

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 6.04

Internal standard: naproxen (3.89)

Limit of detection: 4 ng/mL

OTHER SUBSTANCES

Extracted: diclofenac

Simultaneous: bacitracin, cortisone, diazepam, fluorometholone, hydrocortisone, imipramine, indomethacin, ketoprofen, ketorolac, levobunolol, meclofenamic acid, metipranolol, neomycin, prednisolone, proracaine, propranolol, salicylic acid, sulfacetamide, suprofen

Noninterfering: acebutolol, acetaminophen, acetazolamide, alprenolol, apraclonidine, atenolol, atropine, betamethasone, betaxolol, bupivacaine, caffeine, cyclopentolate, dexamethasone, diphenhydramine, erythromycin, haloperidol, lidocaine, phenylephrine, polymyxin B, procaine, scopolamine, timolol, tropicamide.

KEY WORDS

rabbit; human

REFERENCE

Riegel, M.; Ellis, P.P. High-performance liquid chromatography assay for antiinflammatory agents diclofenac and flurbiprofen in ocular fluids, *J. Chromatogr. B*, **1994**, *654*, 140–145.

SAMPLE

Matrix: blood

Sample preparation: Add 125 μ L MeCN to 50 μ L plasma, mix, centrifuge at 9000g for 10 min, inject a 10 μ L aliquot of the supernatant.

HPLC VARIABLES

Guard column: 5 μ m C18 (Brownlee)

Column: 250 \times 4.6 5 μ m RP C18 (Hibar, Merck, Germany)

Mobile phase: MeCN:water:phosphoric acid 60:40:0.05

Flow rate: 1.5

Injection volume: 10

Detector: F ex 250 em 285

CHROMATOGRAM

Retention time: 3.4

Limit of detection: 50 ng/mL

KEY WORDS

rat; plasma; pharmacokinetics

REFERENCE

Park,K.-M.; Gao,Z.-G.; Kim,C.-K. Assay of flurbiprofen in rat plasma using HPLC with fluorescence detection, *J.Liq.Chromatogr.Rel.Technol.*, **1997**, 20, 1849–1855.

SAMPLE

Matrix: blood

Sample preparation: Mix 500 μ L plasma with 50 μ L 200 μ g/mL IS in MeOH. Add 500 μ L 1 M HCl, extract with 10 mL dichloromethane. Separate the organic layer and evaporate it under a stream of nitrogen at 30°. Reconstitute the residue with 500 μ L mobile phase. Inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 3.9 10 μ m C18 μ Bondapak

Mobile phase: MeCN:MeOH:1% pH 5.8 glacial acetic acid 19:19:62

Flow rate: 1.8

Injection volume: 50

Detector: UV 260

CHROMATOGRAM

Retention time: 5.65

Internal standard: 2-(2'-chloro-4-biphenyl)propionic acid (7.99)

Limit of detection: 500 ng/mL

KEY WORDS

plasma

REFERENCE

Pargal,A.; Kelkar,M.G.; Nayak,P.J. The effect of food on the bioavailability of ibuprofen and flurbiprofen from sustained release formulations, *Biopharm.Drug Dispos.*, **1996**, 17, 511–519.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma + 1 mL water + ibuprofen (15 μ g per 100 μ L) + 2 mL 10% trichloroacetic acid, mix, add 5 mL hexane, mix for 15 min, centrifuge at 800 g for 5-10 min, repeat extraction. Combine hexane layers and evaporate them under a stream of nitrogen at 37-40°. Reconstitute with 300 μ L chloroform, add 200 μ L 65 mg/mL 1,1'-carbonyldiimidazole in chloroform, let stand 5-10 min at room temperature, add 10 μ L glacial acetic acid, vortex briefly, let stand 5-10 min at room temperature, add 50 μ L S-(α)-methylbenzylamine, mix briefly, let stand for 30 min at room temperature, add 3 mL 0.5 M ammonium hydroxide, add 5 mL hexane, mix gently for 15 min. Remove hexane and wash it with 3 mL 1 M HCl, 3 mL 0.5 M ammonium hydroxide, and 3 mL 1 M HCl (with 15 min mixing each time). Evaporate hexane under a stream of nitrogen at 37°, dissolve in 150 μ L mobile phase, inject 25 μ L aliquot.

HPLC VARIABLES

Guard column: Brownlee RP18

Column: 250 \times 4.6 5 μ m Ultrasphere ODS

Mobile phase: MeCN:water 62:38

Flow rate: 1

Injection volume: 25

Detector: UV 245

CHROMATOGRAM

Retention time: 12 (S), 14 (R)

Internal standard: ibuprofen (17 S, 19 R)

Limit of detection: 10 ng/mL

Limit of quantitation: 25 ng/mL

OTHER SUBSTANCES

Simultaneous: metabolites

KEY WORDS

plasma; derivatization; chiral

REFERENCE

Knadler,M.P.; Hall,S.D. High-performance liquid chromatographic analysis of the enantiomers of flurbiprofen and its metabolites in plasma and urine, *J.Chromatogr.*, **1989**, 494, 173–182.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma + 1.25 μ g S-Naproxen + 200 μ L 2 M HCl + 6 mL hexane: diethyl ether 8:2 (ice cold), extract, centrifuge at 1500 g for 5 min. Remove 5 mL of organic layer and evaporate it to dryness under a stream of nitrogen, redissolve in 500 μ L 30 mM pH 7.5 phosphate buffer, inject 50-100 μ L aliquot

HPLC VARIABLES

Column: 100 \times 4.5 μ m Grom AGP

Mobile phase: 2-Propanol:20 mM pH 6.5 phosphate buffer 5:95 containing 1 mM dimethyloctylamine

Column temperature: 15

Flow rate: 0.8

Injection volume: 50-100

Detector: UV 246

CHROMATOGRAM

Retention time: 7.6 (R), 15.2 (S)

Internal standard: S-naproxen (5.2)

Limit of quantitation: 50 ng/mL

KEY WORDS

plasma; chiral

REFERENCE

Geisslinger,G.; Menzel-Soglowek,S.; Schuster,O.; Brune,K. Stereoselective high-performance liquid chromatographic determination of flurbiprofen in human plasma, *J.Chromatogr.*, **1992**, 573, 163–167.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma + 150 μ L 1 M phosphoric acid + 5 mL hexane:ether 80:20, extract. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 200 μ L mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 4 μ m Nova Pak C18

Mobile phase: MeCN:water:acetic acid 59:40.5:0.5

Flow rate: 1.3

Detector: UV 233

CHROMATOGRAM

Internal standard: flurbiprofen

OTHER SUBSTANCES

Extracted: ibuprofen

KEY WORDS

plasma; flurbiprofen is IS

REFERENCE

al-Meshal,M.A.; El-Sayed,Y.M.; al-Balla,S.R.; Gouda,M.W. The effect of colestipol and cholestyramine on ibuprofen bioavailability in man, *Biopharm.Drug Dispos.*, **1994**, 15, 463–471.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 100 μ L 1 M HCl + 100 μ L water + 6 mL ether:hexane 20:80, shake for 10 min, centrifuge at 900 g for 5 min. Remove 4 mL of the organic layer and evaporate it to dryness under a stream of nitrogen at 20°, reconstitute the residue in 100 μ L 10 mM NaOH, sonicate for 3 min, add 50 μ L 100 mM pH 7.0 phosphate buffer containing 0.1% dimethyloctylamine, sonicate for 3 min, inject a 5 μ L aliquot.

HPLC VARIABLES

Guard column: 10 \times 3 5 μ m Chiral-AGP (ChromTech)

Column: 100 \times 4 5 μ m Chiral-AGP (ChromTech)

Mobile phase: Gradient. A was 10 mM pH 7.0 phosphate buffer containing 1 mM dimethyloctylamine. B was isopropanol:10 mM pH 7.0 phosphate buffer 50:50 containing 1 mM dimethyloctylamine. A:B 99.2:0.8 for 5 min, to 59:41 over 10 min, re-equilibrate for 10 min.

Flow rate: 0.9

Injection volume: 5

Detector: UV 220 for 7 min, then UV 245

CHROMATOGRAM

Retention time: 9.5, 11.6 (enantiomers)

Internal standard: flurbiprofen

OTHER SUBSTANCES

Extracted: ibuprofen

KEY WORDS

plasma; chiral; flurbiprofen is IS

REFERENCE

de Vries, J.X.; Schmitz-Kummer, E.; Siemon, D. The analysis of ibuprofen enantiomers in human plasma by high-performance liquid chromatography on an α 1-acid glycoprotein chiral stationary phase, *J. Liq. Chromatogr.*, 1994, 17, 2127–2145.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Plasma + + 25 μ L 2 M HCl, vortex for 15 s, add 2 mL isooctane: isopropanol 85:15, rotate for 5 min, centrifuge at 3000 rpm for 10 min. Remove the upper organic layer and evaporate it to dryness under a stream of nitrogen at 45°, reconstitute the residue in 25 μ L 5 mg/mL 5-bromoacetyl acenaphthene in MeCN, add 10 μ L 3% triethylamine in MeCN, vortex for 30 s, heat at 75° for 5 min, evaporate to dryness under reduced pressure, reconstitute with 25 μ L MeCN, inject a 20 μ L aliquot. (Prepare 5-bromoacetyl acenaphthene as follows. Add 43 g bromoacetyl chloride to 43 g acenaphthene dissolved in 200 mL dichloroethane, cool to -5° in an ice/salt bath, stir vigorously and add 38 g aluminum chloride in small portions over 90 min, do not allow temperature to go above 3°, place under reduced pressure for 30 min, add an excess of crushed ice. Separate the dichloroethane layer and wash it with two 100 mL portions of dilute HCl, wash with 100 mL 5% sodium carbonate solution. Dry the organic layer over anhydrous magnesium sulfate, remove the solvent under reduced pressure, allow the oily residue to solidify, remove liquid by blotting with filter paper. Purify the solid by chromatography on a 300 \times 20 column of 60-120 mesh silica gel, elute with toluene, unreacted acenaphthene elutes first followed by 5-bromoacetyl acenaphthene (mp 87-90°).)

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Hypersil C18

Mobile phase: MeCN:water 90:10

Flow rate: 1

Injection volume: 20

Detector: F ex 250 em 450

CHROMATOGRAM

Internal standard: flurbiprofen

OTHER SUBSTANCES

Extracted: ibuprofen

KEY WORDS

rat; plasma; protect from light; derivatization; flurbiprofen is IS

REFERENCE

Gifford, L.A.; Owusu-Daaku, F.T.K.; Stevens, A.J. Acenaphthene fluorescence derivatization reagents for use in high-performance liquid chromatography, *J.Chromatogr.A*, **1995**, 715, 201–212.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma + 125 μ L 40 mM decanoic acid in MeCN, mix. Dialyze a 100 μ L sample against 20 mM pH 7.0 phosphate buffer using a Gilson Cuprophane membrane (molecular mass cut-off 15 kDa). Continuously pump the buffer through the dialysis cell and through column A at 3 mL/min for 9.6 min, backflush the contents of column A onto column B with the mobile phase, monitor the effluent from column B. (After each injection flush plasma channel with 1 mL 0.05% Triton X-100, with 1 mL 1 mM HCl, and with 2 mL water. After each injection flush buffer channel with 3 mL 20 mM pH 7.0 phosphate buffer and condition column A with 1 mL 20 mM pH 7.0 phosphate buffer.)

HPLC VARIABLES

Column: A 10×2.40 μ m Bondesil C18 (Analytichem); B 250×3.15 μ m C18 (RoSil Research Separation Laboratories)

Mobile phase: MeCN:MeOH:20 mM pH 3.2 phosphate buffer 50:10:40

Flow rate: 1

Injection volume: 100

Detector: UV 247

CHROMATOGRAM

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Extracted: fenoprofen (UV 272), ibuprofen (UV 264), ketoprofen (UV 261), naproxen (UV 272)

KEY WORDS

plasma; dialysis; column-switching

REFERENCE

Herráez-Hernández, R.; Van de Merbel, N.C.; Brinkman, U.A.T. Determination of the total concentration of highly protein-bound drugs in plasma by on-line dialysis and column liquid chromatography: application to non-steroidal anti-inflammatory drugs, *J.Chromatogr.B*, **1995**, 666, 127–137.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma + 100 μ L 25 μ g/mL indomethacin + 500 μ L 600 mM sulfuric acid + 15 mL dichloromethane, mix for 20 min, centrifuge at 2000 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at room temperature, reconstitute the residue in 100 μ L 50 mM triethylamine in MeCN + 50 μ L 60 mM ethyl chloroformate in MeCN, vortex for 30 s, add 50 μ L 100 mM L-leucinamide in MeOH:triethylamine 100:14, let stand for 2 min, add 50 μ L water, inject a 10–50 μ L aliquot of the reaction mixture.

HPLC VARIABLES

Column: 250×4.65 μ m Ultrabase C18 (Shandon)

Mobile phase: MeCN:60 mM KH_2PO_4 :triethylamine 49:51:0.1

Flow rate: 1.8

Injection volume: 10–50

Detector: UV 275

CHROMATOGRAM

Retention time: 3.5 (R-(-)), 4.4 (S-(+))

Internal standard: indomethacin (5.3)

Limit of detection: 100 ng/mL

OTHER SUBSTANCES

Extracted: ibuprofen (UV 225), ketoprofen

KEY WORDS

plasma; chiral; derivatization

REFERENCE

Péhourcq,F.; Lagrange,F.; Labat,L.; Bannwarth,B. Simultaneous measurement of flurbiprofen, ibuprofen, and ketoprofen enantiomer concentrations in plasma using L-leucinamide as the chiral coupling component, *J.Liq.Chromatogr.*, **1995**, *18*, 3969–3979.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 247

CHROMATOGRAM

Retention time: 8.01

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfapyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; metoprolol; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vandesine; mexiletine; dipyrindamole; trazodone; pipamperone; pyrimethamine; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrodine; phenylbutazone; demoxiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenopropfen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidra-

mine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpi-pramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254–262.

SAMPLE

Matrix: blood, milk

Sample preparation: 200 μ L Plasma or milk + 20 μ L 30% perchloric acid + 200 μ L 2 μ g/mL IS in MeOH, vortex for 2 min, add 20 μ L 5 M NaOH, centrifuge at 5000 g. Remove a 200 μ L aliquot of supernatant and add it to 200 μ L mobile phase, mix, centrifuge, inject a 50 μ L aliquot.

HPLC VARIABLES

Guard column: Waters C18 Guard Pak

Column: 100 \times 8 10 μ m μ Bondapak C8 Radial Pak

Mobile phase: MeCN:MeOH:1% acetic acid 30:30:40

Flow rate: 2

Injection volume: 50

Detector: UV 254

CHROMATOGRAM

Retention time: 6

Internal standard: 2-(2'-chloro-4-biphenyl)propionic acid (7)

Limit of detection: 50 ng/mL

KEY WORDS

plasma

REFERENCE

Johnson, V.A.; Wilson, J.T. Flurbiprofen analysis in plasma and breast milk by high-performance liquid chromatography, *J. Chromatogr.*, **1986**, *382*, 367–371.

SAMPLE

Matrix: blood, synovial fluid

Sample preparation: 0.5 mL Plasma or synovial fluid + 200 μ L 2 M HCl + 5 mL hexane, tumble 10 min on a rotary mixer, centrifuge at 10 000 g for 5 min. Remove organic layer and evaporate it to dryness under vacuum centrifugation. Reconstitute residue in 150 μ L MeOH + 100 μ L water, vortex mix, inject aliquot.

HPLC VARIABLES

Guard column: 20 \times 2 Perisorb RP18 30-40 μ m pellicular

Column: 125 \times 4.6 5 μ m Spherisorb ODS 1

Mobile phase: MeOH:water 63:37 adjusted to pH 3.3 with phosphoric acid

Flow rate: 1

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 9

Internal standard: flurbiprofen

Limit of detection: <100 ng/mL

OTHER SUBSTANCES

Extracted: diclofenac

KEY WORDS

plasma; flurbiprofen is IS

REFERENCE

Blagbrough, I.S.; Daykin, M.M.; Doherty, M.; Patrick, M.; Shaw, P.N. High-performance liquid chromatographic determination of naproxen, ibuprofen and diclofenac in plasma and synovial fluid in man, *J. Chromatogr.*, **1992**, 578, 251–257.

SAMPLE

Matrix: blood, urine

Sample preparation: 100 μ L Serum or urine + 50 μ L 1 M NaOH, let stand at room temperature for 20 min (to hydrolyse conjugates), add 50 μ L 1 M HCl, add 1 mL 1 μ g/mL IS in MeCN, add 2 mL 50 mM pH 2.6 potassium phosphate buffer, centrifuge at 2000 rpm for 10 min, inject 100 μ L aliquot of supernatant.

HPLC VARIABLES

Guard column: 30 \times 4.6 5 μ m Brownlee RP-8

Column: 300 \times 3.9 μ Bondapak C18

Mobile phase: THF:50 mM pH 2.6 potassium phosphate 45:55

Flow rate: 1.9

Injection volume: 100

Detector: F ex 360 em 320

CHROMATOGRAM

Retention time: 10

Internal standard: (RS)-2-(2-methoxy-4-biphenyl)propionic acid (8)

Limit of quantitation: 100 ng/mL

KEY WORDS

serum

REFERENCE

Adams, W.J.; Bothwell, B.E.; Bothwell, W.M.; VanGiessen, G.J.; Kaiser, D.G. Simultaneous determination of flurbiprofen and its major metabolite in physiological fluids using liquid chromatography with fluorescence detection, *Anal. Chem.*, **1987**, 59, 1504–1509.

SAMPLE

Matrix: blood, urine

Sample preparation: 500 μ L Urine + 100 μ L 1 M NaOH (to hydrolyze conjugates), vortex. 500 μ L Plasma or 600 μ L basified urine + 200 μ L 600 mM sulfuric acid, vortex, add 50 μ L 100 μ g/mL ketoprofen in water + 3 mL isooctane:isopropanol 96:5, mix vigorously for 45 s, centrifuge at 3000 rpm for 5 min. Remove the top layer and add it to 3 mL water, mix vigorously, centrifuge for 5 min, discard the organic layer. Add the aqueous layer to 350 μ L 600 mM sulfuric acid, add 3 mL chloroform, mix, centrifuge for 5 min. Remove the organic layer and evaporate it to dryness under reduced pressure, reconstitute the residue in 100 μ L 50 mM triethylamine in MeCN, after 30 s add 50 μ L 60 mM ethylchloroformate in MeCN, after 30 s add 50 μ L 1 M L-leucinamide, after 2 min add 50 μ L water, inject a 10–50 μ L aliquot.

HPLC VARIABLES

Guard column: 50 mm long 10 μ m Partisil 5 ODS-3

Column: 100 \times 4.6 5 μ m Partisil 5 ODS-3

Mobile phase: MeCN:67 mM KH_2PO_4 :triethylamine 35:65:0.02

Flow rate: 1

Injection volume: 10–50

Detector: UV 250

CHROMATOGRAM

Retention time: 17 (-), 21 (+)

Internal standard: ketoprofen (UV 275) (8 (-), 10 (+))

Limit of quantitation: 250 ng/mL (urine), 100 ng/mL (plasma)

KEY WORDS

plasma; derivatization; pharmacokinetics; chiral

REFERENCE

Berry, B. W.; Jamali, F. Stereospecific high-performance liquid chromatographic (HPLC) assay of flurbiprofen in biological specimens, *Pharm. Res.*, **1988**, *5*, 123-125.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 21.337

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, *763*, 149-163.

SAMPLE

Matrix: blood, urine

Sample preparation: Serum. Activate an SPE cartridge filled with 100 mg 40-63 μ m silica gel (Merck) with 1 mL MeOH and dry in a hot air oven at 100° for 1 h, equilibrate with 1 mL dichloromethane before use. 500 μ L Serum + 100 μ L 1 M HCl, mix, add 1 mL 1 M pH 3.8 sodium phosphate buffer, mix, add 3 mL diethyl ether, rock for 20 min, centrifuge at 1000 g for 2 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100 μ L 1 mg/mL 1-(3-dimethylaminopropyl)-3-ethylcarbodiimide hydrochloride in dichloromethane, add 100 μ L 1 mg/mL 1-hydroxybenzotriazole in dichloromethane, add 100 μ L 1 mg/mL (R)-(+)-1-(1-naphthyl)ethylamine in dichloromethane, vortex briefly, let stand at room temperature for 2 h, add to the SPE cartridge, elute with two 1 mL portions of dichloromethane:MeCN 90:10. Combine all the eluate and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100 μ L mobile phase, inject a 50 μ L aliquot. Urine. Activate an SPE cartridge filled with 100 mg 40-63 μ m silica gel (Merck) with 1 mL MeOH and dry in a hot air oven at 100° for 1 h, equilibrate with 1 mL dichloromethane before use. 500 μ L Urine + 100 μ L 1 M HCl, mix, add 1.5 mL 1 M pH 3.8 sodium phosphate buffer, mix, add 5 mL hexane:isopropanol 90:10, rock for 20 min, centrifuge at 1000 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100 μ L 1 mg/mL 1-(3-dimethylaminopropyl)-3-ethylcarbodiimide hydrochloride in dichloromethane, add 100 μ L 1 mg/mL 1-hydroxybenzotriazole in dichloromethane, add 100 μ L 1 mg/mL (R)-(+)-1-(1-naphthyl)ethylamine in dichloromethane,

vortex briefly, let stand at room temperature for 2 h, add to the SPE cartridge, elute with two 1 mL portions of dichloromethane:MeCN 90:10. Combine all the eluate and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100 µL mobile phase, inject a 50 µL aliquot.

HPLC VARIABLES

Guard column: 10 × 2.1 40-63 µm pellicular C18 (Alltech)

Column: 150 × 3.9 5 µm Resolve C18 (Waters)

Mobile phase: MeCN:10 mM pH 3.5 phosphate buffer 50:50

Flow rate: 1.5

Injection volume: 50

Detector: F ex 290 em 330

CHROMATOGRAM

Retention time: 14.5 (R), 17.8 (S)

Internal standard: flurbiprofen

OTHER SUBSTANCES

Interfering: ibuprofen (a metabolite of ibuprofen interferes with the R enantiomer)

KEY WORDS

flurbiprofen is IS; derivatization; chiral; serum; SPE

REFERENCE

Tan, S.C.; Jackson, S.H.D.; Swift, C.G.; Hutt, A.J. Enantiospecific analysis of ibuprofen by high performance liquid chromatography: Determination of free and total drug enantiomer concentrations in serum and urine, *Chromatographia*, **1997**, *46*, 23-32.

SAMPLE

Matrix: bulk

Sample preparation: 10 mg Compound + 10 mg 1-(3-dimethylaminopropyl)-3-ethylcarbodiimide hydrochloride + 2 drops 3,5-dimethylaniline + 1.5 mL dichloromethane, mix, after 30 min add 1 mL 1 M HCl, shake vigorously. Remove the lower organic layer and dry it over anhydrous magnesium sulfate, inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm D N-(3,5-dinitrobenzoyl)phenylglycine (Regis)

Mobile phase: Hexane:isopropanol 80:20

Flow rate: 2

Injection volume: 20

Detector: UV 254, UV 280

CHROMATOGRAM

Retention time: k' 2.27 (for first enantiomer)

OTHER SUBSTANCES

Also analyzed: carprofen, cicloprofen, etodolac, fenoprofen, ibuprofen, ketoprofen, naproxen, piroprofen, tiaprofenic acid

KEY WORDS

derivatization; $\alpha = 1.26$; chiral

REFERENCE

Pirkle, W.H.; Murray, P.G. The separation of the enantiomers of a variety of non-steroidal anti-inflammatory drugs (NSAIDS) as their anilide derivatives using a chiral stationary phase, *J. Liq. Chromatogr.*, **1990**, *13*, 2123-2134.

SAMPLE

Matrix: cell suspensions

Sample preparation: Centrifuge cell suspension at 2000 g for 4 min. Remove a 2 mL aliquot of the supernatant and add it to 200 µL 200 µg/mL IS in DMF, mix, add 200 µL 5 M HCl, extract

twice with 3 mL portions of toluene. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 1 mL dichloromethane, add 20 μ L 10 mg/mL 1-hydroxybenzotriazole in dichloromethane:pyridine 99:1, add 300 μ L 10 mg/mL 1-(3-dimethylaminopropyl)-3-ethylcarbodiimide hydrochloride in dichloromethane, add 300 μ L 10 mg/mL (-)-S- α -methylbenzylamine in dichloromethane, let stand for 30 min, evaporate to dryness, reconstitute with 500 μ L mobile phase, inject a 10 μ L aliquot.

HPLC VARIABLES

Guard column: 10 mm long Techsphere ODS (HPLC Technology, Macclesfield UK)

Column: 250 \times 5 μ m Techsphere ODS (HPLC Technology, Macclesfield UK)

Mobile phase: MeCN:7.5 mM NaH₂PO₄ 60:40, containing 5 mM sodium pentanesulfonate, pH adjusted to 2.8 with phosphoric acid

Flow rate: 1

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: k' 4.75, 5.35 (enantiomers)

Internal standard: (S)-naproxen (k' 3.15)

Limit of detection: 500 ng/mL

KEY WORDS

derivatization; chiral

REFERENCE

Thomason,M.J.; Hung,Y.-F.; Rhys-Williams,W.; Hanlon,G.W.; Lloyd,A.W. Indirect enantiomeric separation of 2-arylpropionic acids and structurally related compounds by reversed phase HPLC, *J.Pharm.Biomed.Anal.*, 1997, 15, 1765-1774.

SAMPLE

Matrix: dialysate

Sample preparation: Inject a 10 μ L aliquot of dialysate (pH 7.4 isotonic phosphate buffer) containing 250 ng/mL naproxen.

HPLC VARIABLES

Guard column: 37-50 μ m Corasil C18

Column: 100 \times 4 μ m Nucleosil C18

Mobile phase: MeCN:50 mM pH 3.0 phosphate buffer 48:52

Flow rate: 1.1

Injection volume: 10

Detector: F ex 258 em 310

CHROMATOGRAM

Retention time: 4

Internal standard: naproxen (F ex 262 em 356) (2.5)

KEY WORDS

mouse; rat; pharmacokinetics

REFERENCE

Evrard,P.A.; Deridder,G.; Verbeeck,R.K. Intravenous microdialysis in the mouse and the rat: Development and pharmacokinetic application of a new probe, *Pharm.Res.*, 1996, 13, 12-17.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 300 μ L Microsomal incubation + 20 μ L 6 M HCl, centrifuge at 3000 g for 10 min, inject a 100 μ L aliquot of the supernatant.

HPLC VARIABLES

Guard column: μ Bondapak C18

Column: 250 \times 4 μ m Lichrospher RP-18

Mobile phase: MeCN:water:trifluoroacetic acid 65:145:0.08
Column temperature: 37
Flow rate: 1.5
Injection volume: 100
Detector: UV 273

CHROMATOGRAM

Retention time: 70
Internal standard: 1-naphthol- β -D-glucuronide

OTHER SUBSTANCES

Extracted: metabolites glucuronides

KEY WORDS

rat; human; liver

REFERENCE

Hamdoune,M.; Mounie,J.; Magdalou,J.; Masmoudi,T.; Goudonnet,H.; Escousse,A. Characterization of the in vitro glucuronidation of flurbiprofen enantiomers, *Drug Metab.Dispos.*, **1995**, 23, 343–348.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 Spherisorb S10-ODS2
Mobile phase: MeCN:water:acetic acid 60:35:0.5
Flow rate: 1
Injection volume: 20
Detector: UV 280

OTHER SUBSTANCES

Simultaneous: mefenamic acid

REFERENCE

Galia,E.; Nicolaidis,E.; Hörter,D.; Löbenberg,R.; Reppas,C.; Dressman,J.B. Evaluation of various dissolution media for predicting in vivo performance of class I and II drugs, *Pharm.Res.*, **1998**, 15, 698–705.

SAMPLE

Matrix: solutions
Sample preparation: Dilute in an appropriate solvent, inject an aliquot.

HPLC VARIABLES

Guard column: RC18 Guardpak (Waters)
Column: 150 \times 4.5 5 μ m Altex C18
Mobile phase: MeOH:water:acetic acid 67:32.5:0.5
Flow rate: 1.5
Detector: UV 254

CHROMATOGRAM

Retention time: 8

REFERENCE

Richman,J.B.; Tang-Liu,D.D.-S. A corneal perfusion device for estimating ocular bioavailability in vitro, *J.Pharm.Sci.*, **1990**, 79, 153–157.

SAMPLE

Matrix: solutions
Sample preparation: Mix 100 μ L of a 1–200 μ M solution of carboxylic acid in dichloromethane with 100 μ L 800 μ M (-)-APMB in dichloromethane, 100 μ L 1.6 mM 2,2'-dipyridyl disulfide in dichloromethane, and 100 μ L 1.6 mM triphenylphosphine in dichloromethane, let stand at

room temperature for 20 min. Evaporate to dryness under a stream of nitrogen, reconstitute the residue in 400 μ L mobile phase, inject a 10 μ L aliquot. ((-)-APMB is (-)-2-[4-(1-aminoethyl)phenyl]-6-methoxybenzoxazole. Synthesis is as follows. Hydrogenate 5-methoxy-2-nitrophenol in EtOH over platinum oxide to give 2-amino-5-methoxyphenol (J. Org. Chem. 1957, 22, 220). It should be possible to prepare ethyl 4-acetylbenzimidate hydrochloride ($\text{CH}_3\text{COC}_6\text{H}_4\text{C}(=\text{NH})\text{OC}_2\text{H}_5\cdot\text{HCl}$) by passing dry hydrogen chloride into a mixture of 4-acetylbenzonitrile and 1.2-1.5 equivalents EtOH in an inert solvent (e.g., benzene, chloroform, dioxane, ether, nitrobenzene (Caution! Benzene, chloroform, and dioxane are carcinogens!)) at 0-5°, the benzimidate should crystallize from the mixture in 7-10 days (J. Chem. Soc. 1942, 103). Add a solution of 5.5 g 2-amino-5-methoxyphenol in 200 mL MeOH to 9 g ethyl 4-acetylbenzimidate hydrochloride, stir at 60-70° for 4 h, evaporate to dryness under reduced pressure, recrystallize from EtOH to give 4-(6-methoxy-2-benzoxazolyl)acetophenone as fine orange-yellow crystals (mp 167°) (J. Chromatogr. 1990, 532, 65). Add 7.0 g hydroxylamine hydrochloride and 8.2 g sodium acetate to 10.1 g 4-(6-methoxy-2-benzoxazolyl)acetophenone in 500 mL EtOH: water 95:5, reflux for 1 h, pour into ice-water, filter, recrystallize from EtOH:water 90:10 to give 4-(6-methoxy-2-benzoxazolyl)acetophenone oxime as faint reddish needles (mp 212°). Dissolve 4.7 g 4-(6-methoxy-2-benzoxazolyl)acetophenone oxime in 300 mL MeOH, add 3 g 10% palladium on charcoal, add 10.5 g ammonium formate, reflux for 30 min, filter, evaporate the filtrate to dryness under reduced pressure. Take up the residue in 100 mL 5% HCl and wash the aqueous phase with 100 mL ethyl acetate. Adjust the pH of the aqueous layer to 13-14 with 10% NaOH and extract with 200 mL ethyl acetate. Wash the organic layer with 100 mL water and dry it over anhydrous sodium sulfate, evaporate to dryness under reduced pressure to give racemic 2-[4-(1-aminoethyl)phenyl]-6-methoxybenzoxazole. Dissolve 3.6 g racemic 2-[4-(1-aminoethyl)phenyl]-6-methoxybenzoxazole in 50 mL EtOH and add 3.5 g (S)-(-)- α -methoxy- α -trifluoromethylphenylacetic acid, allow to stand overnight at 5°. Collect the precipitate and fractionally crystallize it from EtOH 4 times. Take up the final product in 5% NaOH and extract it with ethyl acetate, wash the organic layer with water, dry over anhydrous sodium sulfate, evaporate to dryness under reduced pressure, recrystallize from EtOH to give (-)-2-[4-(1-aminoethyl)phenyl]-6-methoxybenzoxazole as pale yellow crystals (mp 74°).

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m TSK gel ODS-80TM (Tosoh)

Mobile phase: MeCN:water:acetic acid 60:40:0.1

Flow rate: 1

Injection volume: 10

Detector: F ex 320 em 380

CHROMATOGRAM

Retention time: 16 (S), 18 (R)

Limit of detection: 10 fmole

OTHER SUBSTANCES

Simultaneous: ibuprofen, naproxen

KEY WORDS

derivatization; chiral

REFERENCE

Kondo,J.; Imaoka,T.; Kawasaki,T.; Nakanishi,A.; Kawahara,Y. Fluorescence derivatization reagent for resolution of carboxylic enantiomers by high-performance liquid chromatography, *J.Chromatogr.*, **1993**, 645, 75-81.

SAMPLE

Matrix: solutions

Sample preparation: Mix 100 μ L of a solution of the carboxylic acid in MeCN with 100 μ L 2 mM (S)-1-(4-dansylaminophenyl)ethylamine in MeCN, add 100 μ L 3 mM triphenylphosphine in MeCN, vortex, let stand at room temperature for 3 h, inject a 5 μ L aliquot. (Synthesis of (S)-1-(4-dansylaminophenyl)ethylamine is as follows. Add 2.2 g di-tert-butyl dicarbonate dropwise to a stirred solution of 2 g (S)- α -methyl-4-nitrobenzylamine hydrochloride ((S)-1-(4-nitrophenyl)ethylamine hydrochloride) and 1.1 g triethylamine in 20 mL MeCN at 0°, stir at room temperature for 1 h, evaporate to dryness under reduced pressure. Dissolve the residue in 50 mL ethyl acetate and wash with 10% aqueous citric acid, wash with water, dry over anhydrous sodium sulfate,

evaporate to dryness under reduced pressure to give (S)-N-tert-butoxycarbonyl-1-(4-nitrophenyl)ethylamine as white crystals (mp 86-89°). Add 200 mg 5% Pd/C to a solution of 2 g (S)-N-tert-butoxycarbonyl-1-(4-nitrophenyl)ethylamine in 40 mL MeOH, stir, hydrogenate at room temperature for 3 h, filter, evaporate the filtrate to dryness to obtain (S)-N-tert-butoxycarbonyl-1-(4-aminophenyl)ethylamine. Stir 1.4 g (S)-N-tert-butoxycarbonyl-1-(4-aminophenyl)ethylamine in 10 mL MeCN and 50 mL 100 mM pH 9.0 sodium bicarbonate, add a solution of 1.9 g dansyl chloride in 30 mL MeCN dropwise while maintaining the pH of the solution at 9.0 with 1 M NaOH, stir at room temperature for 20 min, stir at 45° for 2 h, cool to room temperature, extract three times with 30 mL portions of ethyl acetate. Combine the organic layers and evaporate them to dryness under reduced pressure, recrystallize the residue from benzene/hexane to obtain (S)-N-tert-butoxycarbonyl-1-(4-dansylaminophenyl)ethylamine as pale-yellow crystals (mp 97-101°) (Caution! Benzene is a carcinogen!). Add 2 mL concentrated HCl to a solution of 1.5 g (S)-N-tert-butoxycarbonyl-1-(4-dansylaminophenyl)ethylamine in 10 mL MeOH, stir at room temperature for 30 min, evaporate to dryness under reduced pressure. Dissolve the residue in 30 mL water and adjust the pH to 8.0 with sodium bicarbonate, extract three times with 30 mL portions of ethyl acetate. Combine the organic layers and evaporate them to dryness under reduced pressure, recrystallize the residue from EtOH to obtain (S)-1-(4-dansylaminophenyl)ethylamine as pale-yellow crystals (mp 157-160°; $[\alpha]_D^{28} -11.1^\circ$ (c = 0.2 in MeCN)). (The (R)-enantiomer can be prepared in an exactly analogous fashion.)

HPLC VARIABLES

Column: 150 × 4.6 5 μ m ODS-80TM (Tosoh)

Mobile phase: MeCN:50 mM pH 6.5 sodium acetate buffer 65:35

Flow rate: 1

Injection volume: 5

Detector: F ex 338 em 535

CHROMATOGRAM

Retention time: k' 8.71, k' 10.57 (enantiomers)

Limit of detection: 170 fmole

OTHER SUBSTANCES

Also analyzed: ibuprofen, phenoprofen, 2-phenylpropionic acid, pranoprofen, naproxen

KEY WORDS

derivatization; chiral

REFERENCE

Iwaki, K.; Bunrin, T.; Kameda, Y.; Yamazaki, M. Resolution and sensitive detection of carboxylic acid enantiomers using fluorescent chiral derivatization reagents by high-performance liquid chromatography, *J. Chromatogr. A*, 1994, 662, 87-93.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 10 μ m Chiralpak AD (Daicel)

Mobile phase: Carbon dioxide:MeOH 96:4

Column temperature: 30

Flow rate: 2.5

Detector: UV 210

CHROMATOGRAM

Retention time: 10, 17 (enantiomers)

OTHER SUBSTANCES

Simultaneous: fenoprofen, ibuprofen, ketoprofen, naproxen

KEY WORDS

SFC; 250 bar; chiral

REFERENCE

Kot,A.; Sandra,P.; Venema,A. Sub- and supercritical fluid chromatography on packed columns: A versatile tool for the enantioselective separation of basic and acidic drugs, *J.Chromatogr.Sci.*, **1994**, 32, 439–448.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: Ultrasphere ODS

Mobile phase: MeOH:water:acetic acid 67:32.5:0.5

Flow rate: 1.5

Detector: UV 254

CHROMATOGRAM

Retention time: 7

Limit of detection: 2 ng

OTHER SUBSTANCES

Simultaneous: flurbiprofen amide

KEY WORDS

rabbit; buffer

REFERENCE

Tang-Liu,D.D.-S.; Richman,J.B.; Weinkam,R.J.; Takruri,H. Effects of four penetration enhancers on corneal permeability of drugs in vitro, *J.Pharm.Sci.*, **1994**, 83, 85–90.

SAMPLE

Matrix: solutions

Sample preparation: Mix 1 mL 100 µg/mL compound in dichloromethane with 300 µL 100 µg/mL 1-hydroxybenzotriazole in dichloromethane:pyridine 99:1, 300 µL 1.1 mg/mL 1-ethyl-3-dimethylaminopropylcarbodiimide hydrochloride in dichloromethane, and 300 µL 300 µg/mL benzylamine in dichloromethane, vortex, let stand at room temperature for 1.5 h, evaporate to dryness under a stream of nitrogen, reconstitute the residue in 500 µL MeOH, inject an aliquot.

HPLC VARIABLES

Column: 150 × 4.6 10 µm EXP B101 tris(4-methylbenzoate) cellulose on silica (Bio-Rad)

Mobile phase: MeOH:buffer 70:30 (Prepare buffer solution by dissolving 14.05 g sodium perchlorate in water, adjust pH to 2.0, make up to 1 L with water.)

Flow rate: 1

Detector: UV 230

CHROMATOGRAM

Retention time: 27, 35 (enantiomers)

OTHER SUBSTANCES

Also analyzed: benoxaprofen (MeOH:buffer 80:20), carprofen, fenoprofen, ibuprofen, ketoprofen, piroprofen, tiaprofenic acid

KEY WORDS

derivatization; chiral

REFERENCE

Van Overbeke,A.; Baeyens,W.; Van den Bossche,W.; Dewaele,C. Separation of 2-arylpropionic acids on a cellulose based chiral stationary phase by RP-HPLC, *J.Pharm.Biomed.Anal.*, **1994**, 12, 901–909.

SAMPLE

Matrix: solutions

Sample preparation: Mix 1 mL 100 µg/mL compound in dichloromethane with 300 µL 1 mg/mL 1-hydroxybenzotriazole in dichloromethane:pyridine 99:1, 300 µL 1 mg/mL 1-ethyl-3-di-

methylaminopropylcarbodiimide hydrochloride in dichloromethane, and 300 μ L 1-1.5 mg/mL benzylamine in dichloromethane, vortex, let stand at room temperature for 1.5 h, wash with 1 mL 250 mM HCl, wash with 1 mL water. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in MeOH, inject an aliquot.

HPLC VARIABLES

Column: 10 μ m EXP B101 4-methylbenzoate cellulose on silica gel(Bio-Rad)

Mobile phase: MeOH:50 mM pH 1.5 perchlorate buffer 80:20

Flow rate: 1

Detector: UV 230

CHROMATOGRAM

Retention time: 7, 8

OTHER SUBSTANCES

Simultaneous: ketoprofen

KEY WORDS

derivatization; chiral

REFERENCE

Van Overbeke,A.; Baeyens,W.; Van den Bossche,W.; Dewaele,C. Enantiomeric separation of amide derivatives of some 2-arylpropionic acids by HPLC on a cellulose-based chiral stationary phase, *J.Pharm.Biomed.Anal.*, 1994, 12, 911-916.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 125 \times 3 Ecocart LiChrospher 100 RP-18

Mobile phase: Isopropanol:100 mM KH_2PO_4 :formic acid 54:100:0.1

Flow rate: 0.6

Detector: UV 254

CHROMATOGRAM

Retention time: 12.0

Limit of quantitation: 200-500 ng/mL

OTHER SUBSTANCES

Simultaneous: acemetacin; diclofenac; indomethacin; lonazoloc; ketoprofen; naproxen; piroxicam; sulindac; tenoxicam

REFERENCE

Baeyens,W.R.G.; Van Der Weken,G.; Van Overbeke,A.; Zhang,Z.D. Preliminary results on the LC-separation of non-steroidal anti-inflammatory agents in conventional and narrow-bore RP set-ups applying columns with different internal diameters, *Biomed.Chromatogr.*, 1995, 9, 261-262.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Supelcosil LC-DP (A) or 250 \times 4 5 μ m LiChrospher 100 RP-8 (B)

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 8.01 (A), 8.91 (B)

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordi-azepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyridamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazindol, mefenamic acid, meperidine, mephenytoin, mepivacaine, mesoridazine, metaproterenol, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, nor-epinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spirinolactone, sulfinpyrazone, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocanide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

KEY WORDS

details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, 692, 103–119.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 100 μ M solution in buffer, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 4.6 column containing riboflavin binding proteins (Prepare as follows. Add riboflavin to saturate protein of egg yolk, homogenize with 3 volumes buffer, centrifuge, add the supernatant to a 500 \times 30 column of DEAE-cellulose (Whatman) equilibrated with buffer, wash extensively with buffer to remove bound protein, elute riboflavin binding proteins (RFBP) with buffer containing 200 mM NaCl (RFBP has intense yellow color, absorption at 455 nm). Purify RFBP on a Sephadex G-100 column with 50 mM pH 7.5 Tris-HCl buffer as eluent, remove the bound riboflavin by extensive dialysis at pH 3.0. Add 4.5 g N,N-disuccinylimidyl carbonate to 3 g Nucleosil 5NH₂ slurried in MeCN, filter, wash with MeCN, wash with 50 mM pH 7.5 phosphate buffer. Suspend 300 mg RFBP in 50 mM phosphate buffer, add the activated silica, mix gently for 2 h using a rotary evaporator, filter, wash with sterile water, wash with isopropanol:water 1:2, pack in a 100 \times 4.6 column.) (Buffer was 100 mM pH 5.3 sodium acetate.)

Mobile phase: 50 mM pH 5.5 KH₂PO₄

Flow rate: 0.8

Injection volume: 20

Detector: UV

CHROMATOGRAM

Retention time: k' 6.10

OTHER SUBSTANCES

Simultaneous: isradipine, ketoprofen, nimodipine, suprofen

KEY WORDS

chiral; $\alpha = 1.13$

REFERENCE

Massolini, G.; De Lorenzi, E.; Ponci, M.C.; Gandini, C.; Caccialanza, G.; Monaco, H.L. Egg yolk riboflavin binding protein as a new chiral stationary phase in high-performance liquid chromatography, *J. Chromatogr. A*, **1995**, 704, 55–65.

SAMPLE

Matrix: solutions

Sample preparation: 1 mL 5 mM flurbiprofen in dichloromethane + 300 μ L 1 mg/mL hydroxy-benzotriazole in dichloromethane:pyridine 99:1 + 300 μ L 11 mg/mL 1-ethyl-3-dimethylamino-propylcarbodiimide in dichloromethane + 300 μ L 3.47 mg/mL 1-naphthylamine (Caution! 1-Naphthylamine in a carcinogen!) in dichloromethane, vortex, let stand for 1 h, evaporate to dryness under a stream of nitrogen, reconstitute with 5 mL MeOH, inject an aliquot.

HPLC VARIABLES

Column: 150 \times 2.1 Tollycellulose EXP B101 (tris(4-methylbenzoate)cellulose covalently bonded to 10 μ m aminopropylsilica)

Mobile phase: MeOH:buffer 85:15 (Buffer was 14.05 g/L sodium perchlorate adjusted to pH 2.0.)

Flow rate: 0.21

Injection volume: 1

Detector: UV 230, UV 254

CHROMATOGRAM

Retention time: k' 6.46 (first enantiomer)

OTHER SUBSTANCES

Also analyzed: fenoprofen, ibuprofen, ketoprofen, tiaprofenic acid

KEY WORDS

derivatization; narrow-bore; chiral; $\alpha=2.18$; (see Biomed. Chromatogr. 1995; 9; 292)

REFERENCE

Van Overbeke, A.; Baeyens, W.; Van Der Weken, G.; Van de Voorde, I.; Dewaele, C. Comparative chromatographic study on the chiral separation of the 1-naphthylamine derivative of ketoprofen on cellulose-based columns of different sizes, *Biomed. Chromatogr.*, **1995**, 9, 289–290.

SAMPLE

Matrix: urine

Sample preparation: Condition a Sep-Pak C18 SPE cartridge with 10 mL ethyl acetate and dry by aspiration of air. Evaporate an aliquot of 100 μ L 20 μ g/mL IS in MeOH to dryness at 37°. Add 1 mL urine, vortex, add 250 μ L 1 M pH 5.0 acetate buffer, vortex. Add 250 μ L of the mixture to the SPE cartridge, dry by aspiration of air, elute with 3 mL ethyl acetate, evaporate the eluate to dryness under a stream of nitrogen at 37°, reconstitute the residue in 200 μ L mobile phase, inject a 10–30 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Inertsil ODS-2

Mobile phase: MeCN:50 mM pH 5.0 phosphate buffer 42:58

Flow rate: 0.9

Injection volume: 10–30

Detector: UV 230

CHROMATOGRAM

Retention time: 15

Internal standard: indomethacin (18.5)

Limit of quantitation: 50 ng/mL

OTHER SUBSTANCES

Extracted: diclofenac, ibuprofen, felbinac, fenbufen, ketoprofen, loxoprofen, mefenamic acid, naproxen, piroxicam, sulindac

KEY WORDS

SPE

REFERENCE

Hirai,T.; Matsumoto,S.; Kishi,I. Simultaneous analysis of several non-steroidal anti-inflammatory drugs in human urine by high-performance liquid chromatography with normal solid-phase extraction, *J.Chromatogr.B*, **1997**, 692, 375–388.

SAMPLE

Matrix: urine

Sample preparation: 1 mL Urine + 1 mL 6 M HCl, heat at 90° for 30 min, cool to room temperature, add ibuprofen (15 µg per 100 µL), add 5 mL dichloroethane, mix 15 min, centrifuge at 800 g for 5-10 min. Evaporate organic layer under a stream of nitrogen at 37-40°. Reconstitute with 300 µL chloroform, add 200 µL 65 mg/mL 1,1'-carbonyldiimidazole in chloroform, let stand 5-10 min at room temperature, add 10 µL glacial acetic acid, vortex briefly, let stand 5-10 min at room temperature, add 50 µL S-(α)-methylbenzylamine, mix briefly, let stand for 30 min at room temperature, add 3 mL 0.25 M ammonium hydroxide, add 5 mL dichloroethane, mix gently for 15 min. Remove organic layer and wash it with 3 mL 0.25 M ammonium hydroxide and twice with 3 mL 1 M HCl (with 15 min mixing each time). Evaporate organic layer under a stream of nitrogen at 37°, dissolve in 150 µL MeCN:water 45:55, inject 20-30 µL aliquot.

HPLC VARIABLES

Guard column: Brownlee RP18

Column: 250 × 4.6 5 µm Ultrasphere ODS

Mobile phase: MeCN:50 mM acetic acid 55:45

Flow rate: 1

Injection volume: 20-30

Detector: F ex 200 em 320 (cut-off) (flurbiprofen), UV 232 (ibuprofen)

CHROMATOGRAM

Retention time: 17 (S), 19 (R)

Internal standard: ibuprofen (24 S, 27 R)

Limit of detection: 10 ng/mL

Limit of quantitation: 25 ng/mL

OTHER SUBSTANCES

Simultaneous: metabolites

KEY WORDS

derivatization; chiral

REFERENCE

Knadler,M.P.; Hall,S.D. High-performance liquid chromatographic analysis of the enantiomers of flurbiprofen and its metabolites in plasma and urine, *J.Chromatogr.*, **1989**, 494, 173–182.

SAMPLE

Matrix: urine

Sample preparation: Dilute urine 20-fold with 100 mM pH 2.0 phosphate buffer, extract twice with two volumes of ethyl acetate, centrifuge at 5000 g for 5 min. Combine the organic layers and evaporate them to dryness under a stream of nitrogen below 30°. Reconstitute in 0.2-1 mL mobile phase, inject an aliquot.

HPLC VARIABLES

Guard column: 4 × 4 5 µm LiChrospher 100 RP-18

Column: 250 × 4 5 µm LiChrospher CH-18

Mobile phase: MeOH:10 mM pH 6.0 phosphate buffer 80:20 containing 2.5 mM cethexonium bromide (Rinse with 100 mL MeOH:EtOH:water 50:25:25 at the end of the day.)

Flow rate: 1

Injection volume: 10

Detector: UV 273

CHROMATOGRAM

Retention time: k' 4.1


OTHER SUBSTANCES

Extracted: glucuronides, pirprofen

REFERENCE

Liu,H.-F.; Leroy,P.; Nicolas,A.; Magdalou,J.; Siest,G. Evaluation of a versatile reversed-phase high-performance liquid chromatographic system using cethexonium bromide as ion-pairing reagent for the analysis of glucuronic acid conjugates, *J.Chromatogr.*, **1989**, 493, 137-147.

Flutamide



Molecular formula: C₁₁H₁₁F₃N₂O₃

Molecular weight: 276.22

CAS Registry No.: 13311-84-7

Merck Index: 4242

Lednicer No.: 3 57

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 22.248

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

SAMPLE

Matrix: microsomal incubations

Sample preparation: Vigorously mix 500 μ L microsomal incubation with 5 mL dichloromethane. Evaporate the dichloromethane layer under a stream of nitrogen at 30°. Reconstitute the residue in 100 μ L MeOH. Inject an aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 CP-18 μ Bondapak

Mobile phase: Gradient. A was MeOH:water 60:40. B was MeOH. A:B from 100:0 to 60:40 over 40 min, flush column with MeOH for 10 min, re-equilibrate at initial conditions

Flow rate: 1

Detector: UV; Radioactivity (Radiometer Flo-1)

CHROMATOGRAM

Retention time: 18.3

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

radiolabeled

REFERENCE

Shet,M.D.; McPhaul,M.; Fisher,C.W.; Stallings,N.R.; Estabrook,R.W. Metabolism of the antiandrogenic drug (flutamide) by human CYP1A2, *Drug Metab.Dispos.*, **1997**, 25, 1298–1303.


Flutazolam

Molecular formula: C₁₉H₁₈ClFN₂O₃

Molecular weight: 376.81

CAS Registry No.: 27060-91-9

Merck Index: 4243

**SAMPLE**

Matrix: blood

Sample preparation: 500 μ L Serum + 20 μ L 20 μ g/mL IS + 200 μ L 1 M potassium carbonate + 3 mL chloroform, mix for 2 min, centrifuge at 1200 g for 5 min, aspirate aqueous phase. Evaporate the organic phase in a water bath under a stream of nitrogen at 40°. Dissolve residue in 100 μ L mobile phase, inject a 20 μ L aliquot. (Caution! Chloroform is a carcinogen!)

HPLC VARIABLES

Column: 100 \times 4.6 2 μ m TSK gel Super-ODS (A) or 100 \times 4.6 5 μ m Hypersil ODS-C18 (B)

Mobile phase: MeCN:5 mM pH 6 NaH₂PO₄ 45:55

Flow rate: 0.65

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 32.2 (A), 97.4 (B)

Internal standard: diazepam (29.8 (A), 77.5 (B))

Limit of quantitation: 10 ng/mL (A)

OTHER SUBSTANCES

Extracted: bromazepam, chlordiazepoxide, clonazepam, estazolam, etizolam, haloxazolam, lorazepam, nitrazepam, oxazolam, triazolam

Simultaneous: alprazolam

Noninterfering: barbital, carbamazepine, cloxazolam, ethosuximide, hexobarbital, mexazolam, oxazepam, pentobarbital, phenobarbital, phenytoin, primidone, trimethadione

KEY WORDS

serum

REFERENCE

Tanaka,E.; Terada,M.; Misawa,.; Wakasugi,C. Simultaneous determination of twelve benzodiazepines in human serum using a new reversed-phase chromatographic column on a 2- μ m porous microspherical silica gel, *J.Chromatogr.B*, **1996**, 682, 173–178.


Flutoprazepam

Molecular formula: C₁₉H₁₆ClFN₂O

Molecular weight: 342.80

CAS Registry No.: 25967-29-7

Merck Index: 4245



SAMPLE

Matrix: blood

Sample preparation: 200 μ L Plasma + 25 μ L 1 μ g/mL o-chlorodiazepam in MeOH, vortex gently, extract twice with 1 mL portions of benzene (Caution! Benzene is a carcinogen!) with shaking. Combine the extracts evaporate them to dryness under a stream of nitrogen at 30–35°, reconstitute the residue in 100 μ L mobile phase, inject a 10–90 μ L aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: MeCN:10 mM KH₂PO₄ 50:50 adjusted to pH 4.5 with orthophosphoric acid

Flow rate: 1.2

Injection volume: 10–90

Detector: UV 229

CHROMATOGRAM

Retention time: 16

Internal standard: o-chlorodiazepam

Limit of detection: 25 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; rat; pharmacokinetics

REFERENCE

Conti,I.; Sarati,S.; Caccia,S. Propranolol does not alter flutoprazepam kinetics and metabolism in the rat, *Eur.J.Drug Metab.Pharmacokinet.*, **1991**, 16, 53–58.


Flutrimazole

Molecular formula: $C_{22}H_{16}F_2N_2$

Molecular weight: 346.38

CAS Registry No.: 119006-77-8

Merck Index: 4247



SAMPLE

Matrix: formulations, tissue

Sample preparation: Formulations. 500 mg Cream + 30 mL EtOH, heat mixture to 50° to melt it, cool, dilute to 50 mL, filter through a 0.45 μ m filter. Inject a 25 μ L aliquot of the filtrate. Tissue. Add 3 mL MeOH to cleaned and cut skin, sonicate for 15 min, inject an aliquot.

HPLC VARIABLES

Guard column: 10 μ m μ Bondapak C18

Column: 10 μ m μ Bondapak C18

Mobile phase: MeOH:pH 7.0 phosphate buffer 80:20

Flow rate: 1

Injection volume: 25

Detector: UV 225

CHROMATOGRAM

Limit of quantitation: 50 ng/mL

KEY WORDS

cream; skin

REFERENCE

Ramis,J.; Conte,L.; Segado,X.; Forn,J.; Lauroba,J.; Calpena,A.; Escribano,E.; Domenech,J. Influence of formulation on the in vitro transdermal penetration of flutrimazole, *Arzneimittelforschung*, **1997**, 47, 1139–1144.

Fluvastatin


Molecular formula: $C_{24}H_{26}FNO_4$

Molecular weight: 411.47

CAS Registry No.: 93957-54-1, 93957-55-2 (sodium salt)

Merck Index: 4250

Lednicer No.: 5 105



SAMPLE

Matrix: blood

Sample preparation: Mix 500 μ L plasma with 500 μ L MeCN and 500 μ L pH 6.0 phosphate buffer, extract with 5 mL MTBE by shaking for 30 min, centrifuge at 1200 g for 5 min, evaporate the organic phase under nitrogen. Dissolve the residue in 400 μ L MeCN:water 40:60, vortex 3 times for 30 s, inject a 20–300 μ L aliquot.

HPLC VARIABLES

Guard column: 15 \times 3.2 7 μ m Brownlee Si

Column: 250 \times 4.6 10 μ m Chiralcel OD-R (Daicel Chemical Industries, Japan)

Mobile phase: MeCN:buffer 40:60 (Buffer was pH 2.5 phosphate buffer (I = 0.04).)

Column temperature: 15

Flow rate: 0.5

Injection volume: 100
Detector: F ex 305 em 390

CHROMATOGRAM

Retention time: 28.0 (3S, 5R), 30.0 (3R, 5S)
Limit of quantitation: 5 nM

KEY WORDS

chiral; plasma

REFERENCE

Toreson, H.; Eriksson, B.-M. Determination of fluvastatin enantiomers and the racemate in human blood plasma by liquid chromatography and fluorimetric detection, *J. Chromatogr. A*, **1996**, 729, 13–18.

SAMPLE

Matrix: blood

Sample preparation: Mix 500 μ L plasma with 100 μ L 500 nM IS, 500 μ L MeCN, and 500 μ L pH 6.0 phosphate buffer, extract with 5 mL MTBE by shaking for 30 min, centrifuge at 1200 g for 5 min, evaporate the organic phase under nitrogen. Dissolve the residue in 400 μ L mobile phase, vortex 3 times for 30 s, inject an aliquot.

HPLC VARIABLES

Guard column: 15 \times 3.2 7 μ m Brownlee CN

Column: 150 \times 4.6 5 μ m Zorbax Rx-C8

Mobile phase: MeOH:pH 6.0 phosphate buffer:100 mM tetrabutylammonium fluoride 60:25:15

Column temperature: 40

Flow rate: 1

Injection volume: 50

Detector: F ex 305 em 390 following post-column photolysis. The column effluent flowed through a knitted 10 m \times 0.3 mm PTFE coil irradiated at 254 nm to the detector.

CHROMATOGRAM

Retention time: 15.0

Internal standard: Sandoz compound 63-267 ([R,S,-(E)-] (\pm)-7-(3-(4-fluorophenyl)-1-(1-methyl-ethyl)-1H-indol-2-yl)-(3,5-dihydroxy-6-methyl)-6-heptenoic acid, monosodium salt) (20.5)

Limit of quantitation: 500 pM

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; post-column reaction; post-column photochemical derivatization

REFERENCE

Toreson, H.; Eriksson, B.-M. Determination of fluvastatin enantiomers and the racemate in human blood plasma by liquid chromatography and fluorimetric detection, *J. Chromatogr. A*, **1996**, 729, 13–18.


Fluvoxamine

Molecular formula: C₁₅H₂₁F₃N₂O₂

Molecular weight: 318.34

CAS Registry No.: 54739-18-3

Merck Index: 4251



SAMPLE

Matrix: blood

Sample preparation: 2 mL Plasma + 100 μ L 900 ng/mL IS + 4 mL 300 mM trisodium phosphate, extract with 400 μ L diisopropyl ether for 20 min, centrifuge. Evaporate upper organic layer to dryness under a stream of nitrogen, dissolve residue in 100 μ L MeCN, inject a 50 μ L aliquot. (Caution! Diisopropyl ether readily forms explosive peroxides!)

HPLC VARIABLES

Column: 150 \times 4.6 3 μ m Supelcosil LC-SI

Mobile phase: MeCN:MeOH:concentrated ammonia 87.5:12.0:0.5

Flow rate: 2

Injection volume: 50

Detector: UV 254

CHROMATOGRAM

Retention time: 2.1

Internal standard: 5-(pyrrolidinylpropyliden)-10,11-dihydro-5H-dibenzo[a,d]cyclohepten (4.1)

Limit of detection: 2.2 nM

OTHER SUBSTANCES

Noninterfering: amitriptyline, citalopram, clomipramine, desmethylclomipramine, imipramine, nortriptyline, risperidone

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Carrillo, J.A.; Dahl, M.-L.; Svensson, J.-O.; Alm, C.; Rodríguez, I.; Bertilsson, L. Disposition of fluvoxamine in humans is determined by the polymorphic CYP2D6 and also by the CYP1A2 activity, *Clin. Pharmacol. Ther.*, **1996**, 60, 183–190.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 20 μ L 10 μ g/mL clovoxamine + 120 μ L 2 M NaOH + 4 mL heptane:isopropanol 98:2, shake for 30 min, centrifuge at 3000 g for 10 min. Remove the organic layer and add it to 100 μ L 100 mM HCl, shake for 20 min, centrifuge at 3000 g for 10 min, inject a 20 μ L aliquot of the aqueous layer.

HPLC VARIABLES

Column: 120 \times 4.6 5 μ m Nucleosil C8

Mobile phase: MeCN:buffer 36:64 (Buffer was 16 mM KH_2PO_4 adjusted to pH 2.5 with concentrated phosphoric acid.)

Flow rate: 1

Injection volume: 20

Detector: UV 215

CHROMATOGRAM

Retention time: 4.3

Internal standard: clovoxamine (3.3)

Limit of detection: 25 ng/mL

OTHER SUBSTANCES

Simultaneous: amitriptyline, chlorimipramine, desipramine, doxepin, imipramine, nortriptyline, trimipramine

KEY WORDS

plasma

REFERENCE

Foglia, J.P.; Birder, L.A.; Perel, J.M. Determination of fluvoxamine in human plasma by high-performance liquid chromatography with ultraviolet detection, *J. Chromatogr.*, **1989**, 495, 295–302.

SAMPLE**Matrix:** blood**Sample preparation:** 2 mL Plasma + 200 μ L 1 μ g/mL metapramine in MeOH + 2 mL 1 M pH 10.0 phosphate buffer + 6 mL diethyl ether:hexane 50:50, shake for 15 min, centrifuge at 4000 g for 5 min. Remove the organic layer and add it to 2 mL 62.5 mM sulfuric acid, vortex for 5 min, centrifuge at 4000 g for 5 min. Remove the aqueous phase and add it to 1 mL 500 mM NaOH, vortex, add 6 mL hexane:diethyl ether 50:50, shake for 10 min, centrifuge at 4000 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 50°, reconstitute the residue in 100 mM sodium carbonate, add 10 μ L 10 mg/mL dansyl chloride in acetone, vortex for 1 min, heat at 45° for 30 min, evaporate under a stream of nitrogen at 50°. Reconstitute the residue in 200 μ L MeCN:water 45:55, inject a 100 μ L aliquot.

HPLC VARIABLES**Column:** 125 \times 4.6 5 μ m Hypersil ODS**Mobile phase:** Gradient. MeCN:water from 45:55 to 65:35 over 10 min, maintain at 65:35 for 20 min**Column temperature:** 30**Flow rate:** 1.5**Injection volume:** 100**Detector:** F (Fluorichrom ex 7.54 and 7.60 filters, em 3.71 and 4.76 filters)

CHROMATOGRAM**Retention time:** 19**Internal standard:** metapramine (28)**Limit of detection:** 1.5 ng/mL

OTHER SUBSTANCES**Noninterfering:** alimemazine, alprazolam, amineptine, amitriptyline, caffeine, clobazam, clomipramine, clorazepate, cyamemazine, diazepam, demethyldiazepam, flunitrazepam, levomepromazine, loprazolam, lorazepam, meprobamate, nitrazepam, oxazepam, triazolam, viloxazine

KEY WORDS

plasma; protect from light; derivatization

REFERENCEPommery,J.; Lhermitte,M. High performance liquid chromatographic determination of fluvoxamine in human plasma, *Biomed.Chromatogr.*, **1989**, 3, 177–179.

SAMPLE**Matrix:** blood**Sample preparation:** 1 mL Plasma + 500 μ L 2 M sodium bicarbonate + 100 μ L 0.1 or 1 μ g/mL clovoxamine in MeOH + 10 mL hexane, shake horizontally for 20 min, centrifuge at 2000 g for 10 min. Remove the upper organic layer and evaporate it to dryness under a stream of nitrogen at 35°, reconstitute the residue in 200 μ L mobile phase, inject a 50 μ L aliquot.

HPLC VARIABLES**Column:** 150 \times 3.9 5 μ m Resolve spherical silica (Waters)**Mobile phase:** MeOH:MeCN:THF:water:diethylamine 98.59:1:0.2:0.2:0.01**Flow rate:** 1**Injection volume:** 50**Detector:** UV 254

CHROMATOGRAM**Retention time:** 4.7**Internal standard:** clovoxamine (5.2)**Limit of detection:** 0.5 ng/mL**Limit of quantitation:** 2 ng/mL

OTHER SUBSTANCES**Simultaneous:** alimemazine, amitriptyline, chlorpromazine, clomipramine, desipramine, fluoxetine, haloperidol, imipramine, levomepromazine, norfluoxetine, nortriptyline, propericiazine, trimipramine, viloxazine

Noninterfering: metabolites, amineptine, buspirone, clobazam, clonazepam, clorazepate, diazepam, flunitrazepam, lorazepam, nitrazepam, oxazepam

Interfering: imipramine

KEY WORDS

plasma; human; rat; normal phase; pharmacokinetics

REFERENCE

Van Der Meersch-Mougeot,V.; Diquet,B. Sensitive one-step extraction procedure for column liquid chromatographic determination of fluvoxamine in human and rat plasma, *J.Chromatogr.*, **1991**, 567, 441–449.

SAMPLE

Matrix: blood

Sample preparation: Add 10 μL 20 $\mu\text{g/mL}$ oxaprotiline in MeOH to 990 μL plasma or serum. Inject 100 μL plasma or serum onto column A with mobile phase A and elute to waste, after 15 min elute column A onto column B with mobile phase B for 2 min. Remove column A from circuit and re-equilibrate it with mobile phase A for 5 min. Chromatograph on column B with mobile phase B.

HPLC VARIABLES

Column: A 20 \times 4.6 10 μm Hypersil MOS C8; B 20 \times 4.6 5 μm Hypersil CPS CN + 250 \times 4.6 5 μm Nucleosil 100 CN

Mobile phase: A MeOH:water 5:95; B MeOH:MeCN:10 mM pH 6.8 potassium phosphate buffer 188:578:235

Flow rate: 1.5

Injection volume: 100

Detector: UV 214

CHROMATOGRAM

Retention time: 6.2

Internal standard: oxaprotiline (9.5)

Limit of detection: 20 ng/mL

OTHER SUBSTANCES

Simultaneous: metoclopramide, doxepin, amitriptyline, clomipramine, fluoxetine, imipramine, norfluoxetine, nortriptyline, desipramine, maprotiline

Noninterfering: haloperidol, spiroperidol, pimozide, fluspirilene, trifluoperidol, perazine, chlor-diazepoxide, clobazam, diazepam, nordiazepam, flurazepam, lorazepam, nitrazepam, oxazepam, carbamazepine

Interfering: clozapine

KEY WORDS

plasma; serum; column-switching

REFERENCE

Härter,S.; Wetzel,H.; Hiemke,C. Automated determination of fluvoxamine in plasma by column-switching high-performance liquid chromatography, *Clin.Chem.*, **1992**, 38, 2082–2086.

SAMPLE

Matrix: blood

Sample preparation: For each 1 mL plasma or serum add 10 μL 14 $\mu\text{g/mL}$ trimipramine in MeOH. Inject serum or plasma directly onto column A with mobile phase A, elute with mobile phase A to waste. After 15 min elute column A onto column B (foreflush) with mobile phase B. After 2 min remove column A from the circuit, elute column B with mobile phase B, monitor the effluent from column B. Re-equilibrate column A with mobile phase A.

HPLC VARIABLES

Column: A 20 \times 4.6 10 μm Hypersil MOS C8; B 20 \times 4.6 5 μm Hypersil CPS CN + 250 \times 4.6 5 μm Nucleosil 100 CN

Mobile phase: A MeOH:water 5:95; B MeCN:MeOH:buffer 578:188:235 (Buffer was 10 mM K_2HPO_4 adjusted to pH 6.8 with 85% phosphoric acid.)

Flow rate: 1.5
Injection volume: 100
Detector: UV 214

CHROMATOGRAM

Retention time: 6.19

Internal standard: trimipramine (6.5)

Limit of detection: 1 ng/mL (with three injections onto column A before switching), 5-10 ng/mL

OTHER SUBSTANCES

Extracted: metabolites, amitriptyline, clomipramine, desipramine, doxepin, imipramine, maprotiline, nortriptyline

Noninterfering: chlordiazepoxide, clobazam, clozapine, diazepam, flurazepam, fluspirilene, haloperidol, nitrazepam, oxazepam, perazine, pimozone, spiroperidol, trifluoperidol

KEY WORDS

plasma; serum; column-switching

REFERENCE

Härter, S.; Hiemke, C. Column switching and high-performance liquid chromatography in the analysis of amitriptyline, nortriptyline and hydroxylated metabolites in human plasma or serum, *J.Chromatogr.*, **1992**, 578, 273-282.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 10 μ L 100 μ g/mL nortriptyline in water + 100 μ L 1 M pH 7.6 K_2HPO_4 , vortex for 5 s, add 6 mL ethyl acetate, vortex for 1 min, centrifuge at 1900 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 45°, reconstitute the residue in 250 μ L 40 mM sodium bicarbonate, add 20 μ L 10 mg/mL dansyl chloride in acetone, add 750 μ L acetone, vortex for 30 s, let stand at room temperature at 22° for 15 min. Evaporate under a stream of nitrogen at 45°, reconstitute in 1 mL mobile phase, vortex for 1 min, centrifuge at 1900 g for 10 min, inject a 100 μ L aliquot of the supernatant.

HPLC VARIABLES

Guard column: 15 \times 3.2 7 μ m NewGuard RP-8 (Brownlee)

Column: 250 \times 4.6 5 μ m Supelcosil LC-18-DB

Mobile phase: MeCN:10 mM pH 7.2 potassium phosphate 85:15

Flow rate: 1.5

Injection volume: 100

Detector: F (wavelengths not given)

CHROMATOGRAM

Retention time: 5.3

Internal standard: nortriptyline (9.7)

Limit of quantitation: 10 ng/mL

OTHER SUBSTANCES

Simultaneous: desipramine

KEY WORDS

plasma; derivatization

REFERENCE

Pullen, R.H.; Fatmi, A.A. Determination of fluvoxamine in human plasma by high-performance liquid chromatography with fluorescence detection, *J.Chromatogr.*, **1992**, 574, 101-107.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 1 mL 0.6 M pH 9.8 carbonate buffer + 40 μ L 5 μ g/mL maprotiline in 10 mM HCl + 5 mL 200 g/L ethyl acetate in n-heptane, mix by rocking for 10 min, centrifuge at 1500 g for 10 min. Remove organic layer and add it to 150 μ L 100 mM HCl,

mix 10 min, centrifuge at 1500 g for 10 min. Discard organic layer and evaporate aqueous layer at 45° in a vacuum centrifuge for 1 h. Take up residue in 50 µL 1 M pH 10.3 carbonate buffer and 25 µL 10 mg/mL dansyl chloride in MeCN, vortex, allow to react at room temperature for 45 min, evaporate at 45° in a vacuum centrifuge for 20 min, reconstitute in 125 µL MeCN:water 75:25, vortex, centrifuge for 3-5 min, inject a 25-40 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Supelcosil LC-18

Mobile phase: MeCN:25 mM KH₂PO₄ 75:25 + 500 µL/L orthophosphoric acid + 600 µL/L n-butylamine

Flow rate: 2

Injection volume: 25-40

Detector: F ex 235 em 470 (cut-off)

CHROMATOGRAM

Retention time: 6.78

Internal standard: maprotiline (12.8)

OTHER SUBSTANCES

Simultaneous: fluoxetine, propranolol, clovoxamine, fenfluramine, amoxapine, desipramine, protriptyline, nortriptyline, sertraline, norfluoxetine

Noninterfering: amitriptyline, imipramine, clomipramine, trimipramine, mianserin, chlordiazepoxide, trazodone, cyclobenzaprine, nomifensine, bupropion, metoprolol, atenolol, pindolol, tranlycypromine, moclobemide, thioridazine, citalopram, clozapine, carbamazepine, doxepin, loxapine

KEY WORDS

plasma

REFERENCE

Suckow,R.F.; Zhang,M.F.; Cooper,T.B. Sensitive and selective liquid-chromatographic assay of fluoxetine and norfluoxetine in plasma with fluorescence detection after precolumn derivatization, *Clin.Chem.*, **1992**, *38*, 1756-1761.

SAMPLE

Matrix: blood

Sample preparation: Condition a 1 mL BondElut C18 SPE cartridge with 1 mL 1 M HCl, 1 mL MeOH, 1 mL water, and 1 mL 1% potassium carbonate. 700 µL Serum + 50 µL 5 µg/mL trimipramine in 5% potassium bicarbonate + 700 µL MeCN, vortex, centrifuge at 1500 g for 5 min, add supernatant to SPE cartridge (at ca. 1 mL/min). Wash with 2 mL water and 1 mL MeCN, elute with 250 µL MeOH:35% perchloric acid 20:1 by gravity (10 min) then centrifuge for 20 s to remove rest of eluant, inject a 50 µL aliquot of the eluate.

HPLC VARIABLES

Guard column: 15 mm 7 µm Brownlee RP-8

Column: 150 × 4.6 5 µm Ultrasphere Octyl

Mobile phase: MeCN:water 37.5:62.5 containing 0.5 g/L tetramethylammonium perchlorate and 0.5 mL/L 7% perchloric acid

Flow rate: 1.5

Injection volume: 50

Detector: UV 215

CHROMATOGRAM

Retention time: 6.8

Internal standard: trimipramine (9.6)

Limit of quantitation: 5 ng/mL

OTHER SUBSTANCES

Extracted: amitriptyline, clomipramine, doxepin, fluoxetine, maprotiline, nortriptyline

Interfering: desmethylmaprotiline, desipramine, imipramine, protriptyline

KEY WORDS

serum; SPE

REFERENCE

Gupta, R.N. An improved solid phase extraction procedure for the determination of antidepressants in serum by column liquid chromatography, *J. Liq. Chromatogr.*, **1993**, *16*, 2751–2765.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Serum + 5 mL 300 mM sodium phosphate + 400 μ L diisopropyl ether (Caution! Diisopropyl ether readily forms explosive peroxides!) + imipramine, mix for 20 min, centrifuge for 10 min, inject an aliquot of the organic layer.

HPLC VARIABLES

Column: 150 \times 4.6 3 μ m Apex Silica (Jones Chromatography)

Mobile phase: MeCN:MeOH:25% ammonia 345:65:1.7

Flow rate: 1.3

Detector: UV 254

CHROMATOGRAM

Retention time: 3.5

Internal standard: imipramine

Limit of quantitation: 0.5 nM

KEY WORDS

serum; normal phase; pharmacokinetics

REFERENCE

Spigset, O.; Carlborg, L.; Hedenmalm, K.; Dahlqvist, R. Effect of cigarette smoking on fluvoxamine pharmacokinetics in humans, *Clin. Pharmacol. Ther.*, **1995**, *58*, 399–403.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 μ L mobile phase, centrifuge at 2800 g for 5 min, inject a 50 μ L aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 \times 3.9 4 μ m NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH_2PO_4 adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 253

CHROMATOGRAM

Retention time: 8.71

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoyllecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; tolaxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; car-

teolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melfhalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; videsine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phencyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenopropfen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpi-pramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254–262.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Serum + oxprotiline, make alkaline with borate buffer, extract with cyclohexane:dichloromethane 60:40. Remove the organic layer and extract it with acid, inject an aliquot of the acid layer.

HPLC VARIABLES

Column: C18 DB (Supelco)

Mobile phase: MeCN:pH 2.5 phosphate buffer 37:63

Detector: UV 254

CHROMATOGRAM

Internal standard: oxprotiline

OTHER SUBSTANCES

Extracted: bromazepam, clobazam, diazepam, lorazepam, oxazepam

KEY WORDS

serum

REFERENCE

Vandenbergh, H.; MacDonald, J.C. Analysis of fluvoxamine, clobazam and other benzodiazepines on the same HPLC system (Abstract 40), *Ther. Drug Monit.*, **1995**, *17*, 393–393.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the

organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 15.347

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

SAMPLE

Matrix: formulations

Sample preparation: Pulverize a tablet, add 100 mL MeOH, shake mechanically for 5 min, centrifuge an aliquot at 3000 rpm for 5 min. Remove a 50 μ L aliquot of the supernatant and add it to 50 μ L 100 μ g/mL propyl paraben in MeOH, make up to 1 mL with MeCN, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 8 10 μ m μ Bondapak C18

Mobile phase: MeCN:50 mM ammonium acetate 60:40 adjusted to pH 5.2 with glacial acetic acid

Flow rate: 0.8

Injection volume: 20

Detector: UV 273

CHROMATOGRAM

Retention time: 9

Internal standard: propyl paraben (7)

KEY WORDS

tablets

REFERENCE

Foda,N.H. Quantitative analysis of fluvoxamine maleate in tablet formulations by HPLC, *J.Liq.Chromatogr.*, **1995**, 18, 1591–1601.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Supelcosil LC-DP (A) or 250 \times 4 5 μ m LiChrospher 100 RP-8 (B)

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 9.97 (A), 5.94 (B)

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordi-azepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazindol, mefenamic acid, meperidine, mephenytoin, mepivacaine, mesoridazine, metaproterenol, methadone, methdiazine, methocarbital, methotrexate, methotrimeprazine, methoxamine, methylidopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, nor-epinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenidine, quinidine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfipyrazole, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocainide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimetoprim, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

KEY WORDS

details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, 692, 103–119.

SAMPLE

Matrix: urine

Sample preparation: 8 mL Urine + 2 mL 100 mM pH 9.5 borate buffer, mix, filter (0.2 µm). Remove a 500 µL aliquot of the filtrate and add it to 100 µL 10 mM NaCN in 20 mM pH 9.5 borate buffer, add 500 µL 1 mM naphthalene-2,3-dicarboxaldehyde in MeOH, mix, let stand at room temperature for 20 min, add 100 µL 100 mM glycine in 20 mM pH 9.5 borate buffer, mix, let stand for 10 min, add 500 µL hexane:toluene 50:50, extract. Remove a 400 µL aliquot of the organic layer and add it to 800 µL dichloromethane, mix, inject a 35 µL aliquot.

HPLC VARIABLES

Column: 250 × 3.1 5 µm LiChrosorb Si-60

Mobile phase: Dichloromethane:MeOH 99.8:0.2

Flow rate: 0.5

Injection volume: 35

Detector: Chemiluminescence (418 nm cutoff filter) following post-column reaction. The column effluent mixed with 50 mM hydrogen peroxide in MeCN:dichloromethane 50:50 containing 0.5 mM triethylamine pumped at 0.1 mL/min and with 5 mM bis(2,4,6-trichlorophenyl) oxalate in dichloromethane pumped at 0.1 mL/min and the mixture flowed into the detector.

CHROMATOGRAM

Retention time: 7

Limit of detection: 5 nM

KEY WORDS

derivatization; normal phase

REFERENCE

Kwakman,P.J.M.; Koelewijn,H.; Kool,I.; Brinkman,U.A.T.; de Jong,G.J. Naphthalene- and anthracene-2,3-dialdehyde as precolumn labelling reagents for primary amines using reversed- and normal-phase liquid chromatography with peroxyoxalate chemiluminescence detection, *J. Chromatogr.*, **1990**, 511, 155–166.


Folic acid

Molecular formula: C₁₉H₁₉N₇O₆

Molecular weight: 441.40

CAS Registry No.: 59-30-3

Merck Index: 4253



SAMPLE

Matrix: blood, formulations, urine

Sample preparation: Tablets. Powder tablets, dissolve in water, inject a 10 µL aliquot. Injections. Dilute with water, inject a 10 µL aliquot. Plasma, urine. Condition a Lichrolut RP-18 (Merck) SPE cartridge with 3 mL MeOH and 3 mL water. Mix 40 µL plasma or 100 µL urine with twice the volume of MeCN for 2 min, add 100 µL water, centrifuge at 3500 rpm for 15 min, evaporate the supernatant under nitrogen at 45° to remove the organic solvents, add slowly to the SPE cartridge, collect the eluate. Evaporate to dryness under a stream of nitrogen at 45°. Reconstitute the residue with 500 µL MeOH containing 4.2 µg/mL IS. Inject a 10 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.5 µm Lichrosorb RP-18

Mobile phase: Gradient. A was MeOH. B was 50 mM ammonium acetate. A:B from 5:95 to 15:85 over 6 min, to 30:70 over 7 min, maintain at 30:70 over 7 min

Flow rate: 1

Injection volume: 10

Detector: UV 270

CHROMATOGRAM

Retention time: 9.99

Internal standard: xanthine (4.65)

Limit of detection: 3 ng

OTHER SUBSTANCES

Extracted: ascorbic acid, niacin, niacinamide, riboflavin, vitamin B12

KEY WORDS

plasma; SPE; tablets; injections

REFERENCE

Papadoyannis,I.N.; Tsioni,G.K.; Samanidou,V.F. Simultaneous determination of nine water and fat soluble vitamins after SPE separation and RP-HPLC analysis in pharmaceutical preparations and biological fluids, *J. Liq. Chromatogr. Rel. Technol.*, **1997**, 20, 3203–3231.