


**SAMPLE****Matrix:** urine**Sample preparation:** Filter (0.2  $\mu\text{m}$ ), inject an aliquot directly. Hydrolyze conjugates by heating with 6 M HCl at 37° for 18 h, inject an aliquot.**HPLC VARIABLES****Column:** 150  $\times$  3.9 5  $\mu\text{m}$  Resolve C18 (Waters)**Mobile phase:** MeOH:1.5% trifluoroacetic acid on water 10:90**Flow rate:** 0.5**Detector:** UV or radioactivity**CHROMATOGRAM****Retention time:** 1.8**OTHER SUBSTANCES****Extracted:** phenol, phenyl glucuronide, phenyl sulfate**KEY WORDS**

rat

**REFERENCE**Hughes,M.F.; Hall,L.L. Disposition of phenol in rat after oral, dermal, intravenous, and intratracheal administration, *Xenobiotica*, **1995**, 25, 873–883.

# Hydroxychloroquine

**Molecular formula:** C<sub>18</sub>H<sub>26</sub>ClN<sub>3</sub>O**Molecular weight:** 335.88**CAS Registry No.:** 118-42-3, 747-36-4 (sulfate)**Merck Index:** 4863**Lednicer No.:** 1 342**SAMPLE****Matrix:** blood**Sample preparation:** 1 mL Plasma + 500  $\mu\text{L}$  5 M NaOH + 100  $\mu\text{L}$  10  $\mu\text{g/mL}$  chloroquine in MeOH + 5 mL hexane:diethyl ether 50:50, vortex for 1 min, centrifuge at 1000 g for 10 min, freeze in dry ice/acetone, remove the organic layer. Thaw out the aqueous layer and repeat the extraction. Combine the organic layers and evaporate them to dryness under vacuum, reconstitute the residue in 100  $\mu\text{L}$  mobile phase, inject a 50  $\mu\text{L}$  aliquot.**HPLC VARIABLES****Guard column:** 5  $\mu\text{m}$  cyano (Regis)**Column:** 75  $\times$  4.6 3  $\mu\text{m}$  Ultremex cyano (Phenomenex)**Mobile phase:** 20 mM Dimethyloctylamine phosphate:60 mM ammonium acetate 40:60, pH adjusted to 4.5 (Dimethyloctylamine phosphate was prepared by adding phosphoric acid to N,N-dimethyloctylamine to precipitate the salt.)**Flow rate:** 0.6**Injection volume:** 50**Detector:** UV 320**CHROMATOGRAM****Retention time:** 15**Internal standard:** chloroquine (23)**Limit of detection:** 10 ng/mL**OTHER SUBSTANCES****Extracted:** metabolites

**KEY WORDS**

plasma; rabbit; pharmacokinetics

**REFERENCE**

Iredale, J.; Wainer, I.W. Determination of hydroxychloroquine and its major metabolites in plasma using sequential achiral-chiral high-performance liquid chromatography, *J. Chromatogr.*, **1992**, 573, 253–258.

**SAMPLE****Matrix:** blood

**Sample preparation:** Serum. 200  $\mu$ L Serum + 50  $\mu$ L 1  $\mu$ g/mL IS in water + 50  $\mu$ L 4 M NaOH + 200  $\mu$ L MTBE, vortex for 30 s, centrifuge at 9950 g for 2 min, inject 100  $\mu$ L of the organic layer. Whole blood. 100  $\mu$ L Whole blood + 500  $\mu$ L water + 50  $\mu$ L 1  $\mu$ g/mL IS in water + 50  $\mu$ L 4 M NaOH + 200  $\mu$ L MTBE, vortex for 30 s, centrifuge at 9950 g for 2 min, inject 100  $\mu$ L of the organic layer. Dried blood. Spread 100  $\mu$ L whole blood on a 70  $\times$  30 mm piece of filter paper, allow to dry, cut paper into 10  $\times$  5 mm strips, add 100  $\mu$ L 1  $\mu$ g/mL IS in water, add 1.5 mL 0.5 M NaOH, vortex for 30 s, let stand for 30 min at room temperature, add 300  $\mu$ L MTBE, vortex for 30 s, centrifuge at 2000 g for 5 min, inject a 100  $\mu$ L aliquot of the organic layer.

**HPLC VARIABLES****Column:** 150  $\times$  5  $\mu$ m Spherisorb S5SCX sulfophenylpropyl-modified silica**Mobile phase:** MeOH:water 98.5:1.5 containing 9.41 g/L ammonium perchlorate, adjust apparent pH to 8.0 with 220 mL/L 50 mM NaOH in MeOH**Flow rate:** 1.5**Injection volume:** 100**Detector:** F ex 215 em no filter**CHROMATOGRAM****Retention time:** 8**Internal standard:** 6,8-dichloro-4-(1-methyl-4-diethylaminobutylamino)quinoline (5)**Limit of quantitation:** 5 ng/mL (serum), 10 ng/mL (whole blood, dried blood)**OTHER SUBSTANCES****Extracted:** chloroquine, quinine, metabolites**Simultaneous:** acebutolol, N-acetylprocainamide, atenolol, butriptyline, chlorpromazine, desipramine, flecainide, fluoxetine, imipramine, labetalol, maprotiline, mepacrine, metoprolol, mexiletine, norbutriptyline, normaprotiline, procainamide, propranolol, sotalol**Noninterfering:** amitriptyline, amodiaquin, carbamazepine, clomipramine, dapsone, diazepam, dothiepin, doxepin, fluvoxamine, lorazepam, mefloquine, nitrazepam, norclomipramine, nordiazepam, nordothiepin, nordoxepin, nortriptyline, primaquine, proguanil, pyrimethamine**KEY WORDS**

serum; whole blood; dried blood

**REFERENCE**

Croes, K.; McCarthy, P.T.; Flanagan, R.J. Simple and rapid HPLC of quinine, hydroxychloroquine, chloroquine, and desethylchloroquine in serum, whole blood, and filter paper-adsorbed dry blood, *J. Anal. Toxicol.*, **1994**, 18, 255–260.

**SAMPLE****Matrix:** blood

**Sample preparation:** 1 mL Whole blood + 2 mL 20 ng/mL propranolol hydrochloride in water, vortex for 10 s, sonicate for 10 min, centrifuge at 3000 rpm for 20 min. Remove 2 mL of the supernatant and add it to 1 mL 100 mM NaOH and 8 mL ethyl acetate:isopropanol 90:10, vortex for 30 s, centrifuge at 3000 rpm for 15 min. Remove the organic layer and evaporate it to dryness under a stream of air at 35–40°, reconstitute the residue in 250  $\mu$ L MeOH, vortex for 30 s, inject a 10–100  $\mu$ L aliquot.

**HPLC VARIABLES****Guard column:** 10  $\times$  3  $\mu$ m cyano (Applied Biosystems)**Column:** 250  $\times$  4.6  $\mu$ m cyanopropyl (Baxter/Burdick & Jackson)

**Mobile phase:** MeCN:MeOH:buffer 50:30:20 (Buffer was 25 mM K<sub>2</sub>HPO<sub>4</sub> adjusted to pH 6.0 with phosphoric acid.)  
**Column temperature:** 50  
**Flow rate:** 2  
**Injection volume:** 10-100  
**Detector:** F ex 230 em 385 (370 nm cut-off filter)

---

#### CHROMATOGRAM

**Retention time:** 9.72  
**Internal standard:** propranolol (4)  
**Limit of quantitation:** 2 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** metabolites

---

#### KEY WORDS

whole blood

---

#### REFERENCE

Wei, Y.; Nygard, G.A.; Khalil, S.K.W. A HPLC method for the separation and quantification of the enantiomers of hydroxychloroquine and its three major metabolites, *J. Liq. Chromatogr.*, **1994**, *17*, 3479-3490.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol:n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

---

#### HPLC VARIABLES

**Column:** 300 × 3.9 4 µm NovaPack C18

**Mobile phase:** MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH<sub>2</sub>PO<sub>4</sub> adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

**Column temperature:** 30

**Flow rate:** 0.8

**Injection volume:** 50

**Detector:** UV 222

---

#### CHROMATOGRAM

**Retention time:** 4.65

**Limit of detection:** <120 ng/mL

---

#### KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procabazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acen-

ocoumarol; vindesine; mexiletine; dipyrindamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loperazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrridine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

## REFERENCE

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254–262.

## SAMPLE

**Matrix:** blood

**Sample preparation:** 1 mL Serum + 100  $\mu$ L water containing 5  $\mu$ g/mL 2,3-diaminonaphthalene and 3.5  $\mu$ g/mL 18-hydroxy-11-deoxycorticosterone + 1 mL 250 mM NaOH + 7 mL diethyl ether, shake on a rotary shaker for 15 min, repeat extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 30–40°, reconstitute the residue in 70  $\mu$ L MeOH:100 mM perchloric acid 50:50, inject a 20  $\mu$ L aliquot.

## HPLC VARIABLES

**Column:** 150  $\times$  3.9 4  $\mu$ m Nova-Pak C18

**Mobile phase:** Gradient. A was 58 mM NaH<sub>2</sub>PO<sub>4</sub> containing 6 mM sodium heptanesulfonate, adjusted to pH 3.1 with concentrated phosphoric acid. B was MeCN:MeOH 85:15. A:B from 100:0 to 78:22 over 5 min, to 70:30 over 12 min, maintain at 70:30 for 4 min, to 65:35 over 9 min.

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 245, 256, 343

## CHROMATOGRAM

**Retention time:** 11.93

**Internal standard:** 2,3-diaminonaphthalene (10.71), 18-hydroxy-11-deoxycorticosterone (15.85)

**Limit of detection:** 2 ng/mL (343 nm)

## OTHER SUBSTANCES

**Extracted:** betamethasone, chloroquine, corticosterone, cortisone, dexamethasone, fluocinolone acetate, fluendrenolide, fluorometholone, fluprednisolone, hydrocortisone, 17 $\beta$ -hydroxyprogesterone, meprednisolone, methylprednisolone, methylprednisolone acetate, paramethasone, prednisolone, prednisone, progesterone, triamcinolone

**Noninterfering:** aspirin, ibuprofen, indomethacin, phenylbutazone, pregnenolone

## KEY WORDS

serum

## REFERENCE

Volin, P. Simple and specific reversed-phase liquid chromatographic method with diode-array detection for simultaneous determination of serum hydroxychloroquine, chloroquine and some corticosteroids, *J. Chromatogr. B*, **1995**, *666*, 347–353.

## SAMPLE

**Matrix:** blood, erythrocytes, urine

**Sample preparation:** Condition a 3 mL Bond Elut C8 SPE cartridge with 2 mL MeOH and 2 mL buffer. Hemolyze erythrocytes in water 1:3. Dilute urine with water 1:99. Add 1 mL plasma, hemolyzed erythrocytes, or diluted urine to the SPE cartridge, wash with 4 mL buffer, wash with 2 mL MeOH:buffer 50:50, elute with 3 mL MeOH:ammonia 99:1. Evaporate the eluate to dryness under a stream of nitrogen at 30°, reconstitute the residue in the initial mobile phase, vortex, inject a 50  $\mu$ L aliquot. (Prepare buffer by mixing equal volumes of 100 mM ammonium formate and 100 mM ammonia solution, pH 9.2.)

---

**HPLC VARIABLES**

**Guard column:** 10  $\times$  4 Inertsil

**Column:** 250  $\times$  4.5  $\mu$ m Inertsil

**Mobile phase:** Gradient. A was MeCN. B was MeOH:25% ammonia solution 92.5:7.5. A:B 78:22 for 3 min, then to 65:35 over 2 min (Waters curve no. 3), maintain at 65:35 for 20 min, return to 78:22 over 5 min (Waters curve no. 3).

**Flow rate:** 0.85

**Injection volume:** 50

**Detector:** F ex 325 em 375

---

**CHROMATOGRAM**

**Retention time:** 11.5

**Internal standard:** hydroxychloroquine

---

**OTHER SUBSTANCES**

**Extracted:** chloroquine, quinine, monodesethylchloroquine, bidesethylchloroquine

**Simultaneous:** halofantrine, quinidine

**Noninterfering:** proguanil, cycloguanil, 4-chlorophenylbiguanide, amodiaquine, mefloquine, pyrimethamine, sulfadoxine, cinchonine, cinchonidine

---

**KEY WORDS**

SPE; hydroxychloroquine is IS; plasma

---

**REFERENCE**

Chaulet, J.-F.; Robet, Y.; Prevosto, J.-M.; Soares, O.; Brazier, J.-L. Simultaneous determination of chloroquine and quinine in human biological fluids by high-performance liquid chromatography, *J. Chromatogr.*, **1993**, *613*, 303–310.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a solution in mobile phase, inject an aliquot.

---

**HPLC VARIABLES**

**Guard column:** Chiral-AGP

**Column:** 150  $\times$  4.6 Chiral-AGP (Regis)

**Mobile phase:** MeCN:EtOH:30 mM pH 7.0 sodium phosphate buffer 1:20:79 containing 5 mM N,N-dimethyloctylamine

**Flow rate:** 0.9

**Injection volume:** 50

**Detector:** UV 320

---

**CHROMATOGRAM**

**Retention time:** 10 (-), 14 (+)

**Limit of detection:** 5 ng

---

**KEY WORDS**

chiral

---

**REFERENCE**

Iredale, J.; Wainer, I.W. Determination of hydroxychloroquine and its major metabolites in plasma using sequential achiral-chiral high-performance liquid chromatography, *J. Chromatogr.*, **1992**, *573*, 253–258.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a solution in mobile phase, inject a 10–100  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 100 × 4 5 µm Chiral-AGP (ChromTech)

**Mobile phase:** MeCN:isopropanol:buffer 1:5:94 (Buffer was 50 mM (NH<sub>4</sub>)H<sub>2</sub>PO<sub>4</sub> containing 5 mM dihexylamine, pH adjusted to 7.0 with 3 M NaOH.)

**Column temperature:** 40

**Flow rate:** 1

**Injection volume:** 10-100

**Detector:** F ex 230 em 385 (370 nm cut-off filter)

---

**CHROMATOGRAM**

**Retention time:** 10 (S(+)), 13 (R(-))

---

**KEY WORDS**

chiral

---

**REFERENCE**

Wei,Y.; Nygard,G.A.; Khalil,S.K.W. A HPLC method for the separation and quantification of the enantiomers of hydroxychloroquine and its three major metabolites, *J.Liq.Chromatogr.*, **1994**, 17, 3479-3490.

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 150 × 5 Spherisorb S5SCX

**Mobile phase:** MeOH:water 98.5:1.5 containing 80 mM ammonium perchlorate, adjusted to pH 8.0 with 50 mM NaOH in MeOH

**Flow rate:** 1.5

**Detector:** F ex 215 no emission filter

---

**CHROMATOGRAM**

**Retention time:** 8

---

**OTHER SUBSTANCES**

**Simultaneous:** hydroquinine, quinine, chloroquine

---

**REFERENCE**

Croes,K.; McCarthy,P.T.; Flanagan,R.J. HPLC of basic drugs and quaternary ammonium compounds on micro-particulate strong cation-exchange materials using methanolic or aqueous methanol eluents containing an ionic modifier, *J.Chromatogr.A*, **1995**, 693, 289-306.

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4 5 µm LiChrospher 100 RP-8 (B)

**Mobile phase:** MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

**Flow rate:** 0.6

**Injection volume:** 25

**Detector:** UV 229

---

**CHROMATOGRAM**

**Retention time:** 9.60 (A), 3.23 (B)

---

**OTHER SUBSTANCES**

**Also analyzed:** acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordiazepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine,

clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, etidrium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazindol, mefenamic acid, meperidine, mephenytoin, mepivacaine, mesoridazine, metaproterenol, methadone, methdiazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, nor-epinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfipyrazole, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thietilperazine, thiopental, thioridazine, thiothixene, timolol, tocinide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

---

**KEY WORDS**

details of plasma extraction

---

**REFERENCE**

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, 692, 103–119.

---

# Hydroxyprogesterone


**Molecular formula:**  $C_{21}H_{30}O_3$

**Molecular weight:** 330.47

**CAS Registry No.:** 68-96-2

**Merck Index:** 4886

**Lednicer No.:** 1 176



---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 100  $\mu$ L 2  $\mu$ g/mL dexamethasone in EtOH:water 10:90 + 100  $\mu$ L 250 mM NaOH + 7 mL ether:dichloromethane 60:40, vortex for 30 s, centrifuge at 2000 rpm for 5 min. Remove 6 mL of the organic layer and evaporate it to dryness under a stream of air at 40°, reconstitute the residue in 100  $\mu$ L dichloromethane:EtOH:water 95:4:1, inject a 50  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.5 5  $\mu$ m Partisil silica

**Mobile phase:** Dichloromethane:EtOH:water 95:4:1

**Flow rate:** 1.5

**Injection volume:** 50

**Detector:** UV 239

---

**CHROMATOGRAM**

**Retention time:** 2.5 (17-hydroxyprogesterone)

**Internal standard:** dexamethasone (11.5)

**Limit of quantitation:** 25 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** corticosterone, 11-deoxycortisol, hydrocortisone, 6 $\alpha$ -methylprednisolone, prednisolone, prednisone, progesterone

---

**KEY WORDS**

plasma; normal phase

---

**REFERENCE**

Scott,N.R.; Chakraborty,J.; Marks,V. Determination of prednisolone, prednisone, and cortisol in human plasma by high-performance liquid chromatography, *Anal.Biochem.*, **1980**, *108*, 266–268.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Condition a Bond-Elut C18 SPE cartridge with MeOH and water. Add 500  $\mu$ L plasma to the SPE cartridge, wash with 2 mL water, wash with 2 mL MeOH:water 20:80, elute with two 500  $\mu$ L aliquots of MeOH. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute the residue in 50  $\mu$ L MeOH:water 20:80, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 100  $\times$  1.5  $\mu$ m Hypersil ODS

**Mobile phase:** MeCN:MeOH:water 25:25:50

**Flow rate:** 0.1

**Injection volume:** 20

**Detector:** UV (wavelength not given)

---

**CHROMATOGRAM**

**Retention time:** 8 (17 $\alpha$ -hydroxyprogesterone)

**Limit of detection:** 2 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** androstenedione, 20 $\alpha$ -hydroxy-4-pregnen-3-one, norethindrone, progesterone, testosterone

---

**KEY WORDS**

microbore; rat; plasma; SPE

---

**REFERENCE**

Taylor,R.B.; Kendle,K.E.; Reid,R.G.; Hung,C.T. Chromatography of progesterone and its major metabolites in rat plasma using microbore high-performance liquid chromatography columns with conventional injection and detection systems, *J.Chromatogr.*, **1987**, *385*, 383–392.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 100  $\mu$ L Serum + 500  $\mu$ L water + 100  $\mu$ L 10  $\mu$ g/mL 3,7-dimethoxyflavone in EtOH + 8 mL diethyl ether, shake, centrifuge at 4° at 1000 g for 5 min, freeze in acetone/dry ice. Remove the organic layer and dry it over anhydrous sodium sulfate, evaporate to dryness under a stream of nitrogen, reconstitute the residue in 100  $\mu$ L MeOH:water 40:60, inject a 50  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 3  $\mu$ m NS-Gel C18

**Mobile phase:** Gradient. MeOH:water from 40:60 to 55:45, maintain at 55:45 for 24 min, to 80:20 over 25 min

**Column temperature:** 50

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 210, UV 240

---

**CHROMATOGRAM**

**Retention time:** 36.52

**Internal standard:** 3,7-dimethoxyflavone (47)

---



---

**OTHER SUBSTANCES**

**Extracted:** aldosterone, androstenedione, dehydroepiandrosterone, deoxycorticosterone, 11-deoxycortisol, estradiol, estrone, hydrocortisone, pregnenolone, progesterone

---

**KEY WORDS**

serum

---

**REFERENCE**

Ueshiba,H.; Segawa,M.; Hayashi,T.; Miyachi,Y.; Irie,M. Serum profiles of steroid hormones in patients with Cushing's syndrome determined by a new HPLC/RIA method, *Clin.Chem.*, **1991**, 37, 1329-1333.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 1 mL Serum + 100  $\mu$ L water containing 5  $\mu$ g/mL 2,3-diaminonaphthalene and 3.5  $\mu$ g/mL 18-hydroxy-11-deoxycorticosterone + 1 mL 250 mM NaOH + 7 mL diethyl ether, shake on a rotary shaker for 15 min, repeat extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 30-40°, reconstitute the residue in 70  $\mu$ L MeOH:100 mM perchloric acid 50:50, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 150  $\times$  3.9 4  $\mu$ m Nova-Pak C18

**Mobile phase:** Gradient. A was 58 mM NaH<sub>2</sub>PO<sub>4</sub> containing 6 mM sodium heptanesulfonate, adjusted to pH 3.1 with concentrated phosphoric acid. B was MeCN:MeOH 85:15. A:B from 100:0 to 78:22 over 5 min, to 70:30 over 12 min, maintain at 70:30 for 4 min, to 65:35 over 9 min.

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 245, 256, 343

---

**CHROMATOGRAM**

**Retention time:** 29.26 (17 $\beta$ -hydroxyprogesterone)

**Internal standard:** 2,3-diaminonaphthalene (10.71), 18-hydroxy-11-deoxycorticosterone (15.85)

**Limit of detection:** 1-10 ng/mL (245 nm)

---

**OTHER SUBSTANCES**

**Extracted:** betamethasone, chloroquine, corticosterone, cortisone, dexamethasone, fluendrenolide, fluocinolone acetonide, fluorometholone, fluprednisolone, hydrocortisone, hydroxychloroquine, meprednisone, methylprednisolone, methylprednisolone acetate, paramethasone, prednisolone, prednisone, progesterone, triamcinolone

**Noninterfering:** aspirin, ibuprofen, indomethacin, phenylbutazone, pregnenolone

---

**KEY WORDS**

serum

---

**REFERENCE**

Volin,P. Simple and specific reversed-phase liquid chromatographic method with diode-array detection for simultaneous determination of serum hydroxychloroquine, chloroquine and some corticosteroids, *J.Chromatogr.B*, **1995**, 666, 347-353.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Make up 1 mL injection to 100 mL with EtOH, remove a 1.2 mL aliquot and add it to 1 mL 1 mg/mL 17-hydroxyprogesterone in EtOH. Dilute this mixture to 50 mL with EtOH, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 300  $\times$  4  $\mu$ Bondapak CN

**Mobile phase:** MeOH:20 mM KH<sub>2</sub>PO<sub>4</sub> 40:60

**Flow rate:** 2

**Injection volume:** 20

**Detector:** UV 254

---

**CHROMATOGRAM****Retention time:** 7 (hydroxyprogesterone caproate)**Internal standard:** 17-hydroxyprogesterone (3)

---

**OTHER SUBSTANCES****Simultaneous:** benzyl benzoate**Noninterfering:** polyethylene glycol 4000, myristyl-gamma-picolinium chloride, methylcellulose, thimerosal

---

**KEY WORDS**

injections

---

**REFERENCE**

Das Gupta, V. Quantitation of hydroxyprogesterone caproate, medroxyprogesterone acetate, and progesterone by reversed-phase high-pressure liquid chromatography, *J.Pharm.Sci.*, **1982**, *71*, 294–297.

---

**SAMPLE****Matrix:** formulations

**Sample preparation:** Injections. 1 mL Injection (200 µg/mL) + 300 mg NaCl + 200 µL 10% ammonia + 5 mL dichloromethane, shake vigorously for 10 min, let stand for a few min. Remove 4 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 4 mL water, mix an aliquot with an equal volume of 30 µg/mL 17α-hydroxyprogesterone in MeOH, inject a 20 µL aliquot. Tablets. Weigh out amount of powdered tablets equivalent to about 200 µg compound, add 1 mL water, sonicate for 2 min, add 300 mg NaCl, add 200 µL 10% ammonia, add 5 mL dichloromethane, shake vigorously for 10 min, let stand for a few min. Remove 4 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 4 mL water, mix an aliquot with an equal volume of 30 µg/mL 17α-hydroxyprogesterone in MeOH, inject a 20 µL aliquot.

---

**HPLC VARIABLES****Column:** 150 × 4 5 µm LiChrosorb RP-18**Mobile phase:** MeCN:50 mM pH 3.5 acetate buffer 40:60 containing 1.5 mM triethylamine**Column temperature:** 30**Flow rate:** 1**Injection volume:** 20**Detector:** UV 254

---

**CHROMATOGRAM****Retention time:** 12**Internal standard:** 17α-hydroxyprogesterone

---

**OTHER SUBSTANCES****Simultaneous:** benzyl alcohol, ergonovine, methylegonovine**Noninterfering:** ascorbic acid

---

**KEY WORDS**

injections; tablets; 17α-hydroxyprogesterone is IS

---

**REFERENCE**

Tokunaga, H.; Kimura, T.; Kawamura, J. Determination of ergometrine maleate and methylergometrine maleate in pharmaceutical preparations by high-performance liquid chromatography, *Chem.Pharm.Bull.(Tokyo)*, **1983**, *31*, 3988–3993.

---

**SAMPLE****Matrix:** formulations

**Sample preparation:** Grind tablets containing about 1 mg levothyroxine, add 4.5 mL 0.5 mg/mL hydroxyprogesterone caproate in MeOH, add 20.5 mL 10 mM NaOH in MeOH:water 75:25, shake intermittently for 5 min, filter, discard first 5 mL filtrate, inject a 25 µL aliquot.

---

**HPLC VARIABLES****Column:** 300 × 3.9 µBondapak CN

**Mobile phase:** MeCN:0.1% phosphoric acid in water 35:65

**Flow rate:** 3

**Injection volume:** 25

**Detector:** UV 225

---

#### CHROMATOGRAM

**Retention time:** 8 (hydroxyprogesterone caproate)

**Internal standard:** hydroxyprogesterone caproate

---

#### OTHER SUBSTANCES

**Simultaneous:** levothyroxine

---

#### KEY WORDS

tablets; hydroxyprogesterone caproate is IS

---

#### REFERENCE

Das Gupta,V.; Odom,C.; Bethea,C.; Plattenburg,J. Effect of excipients on the stability of levothyroxine sodium tablets, *J.Clin.Pharm.Ther.*, **1990**, *15*, 331–336.

---

#### SAMPLE

**Matrix:** ileostomy effluent

**Sample preparation:** Dilute ileostomy effluent 1:2 by weight with water. Extract 3 g aliquot three times with 10 mL dichloromethane by shaking for 1 min and centrifuging at 2000 rpm for 2 min. Wash combined extracts successively with 2 mL 0.1 M NaOH and 4 mL water by shaking for 30 s and centrifuging for 1 min then dry the organic layer under air at 40°. Take up the extract in 1 mL MeOH, add 1.1 mL water and apply to C18 Bond Elut SPE cartridge. Wash with 10 mL water, wash with 5 mL MeOH:water 45:55, elute with 2 mL MeOH. Add 50  $\mu$ L 20  $\mu$ g/mL progesterone to the eluate, dry at 40°, take up in 100  $\mu$ L MeOH, inject 10  $\mu$ L aliquot.

---

#### HPLC VARIABLES

**Guard column:** Bondapak C18/Corasil

**Column:** 300  $\times$  3.9  $\mu$ Bondapak C18

**Mobile phase:** MeOH:50 mM pH 3.0 sodium phosphate buffer 55:45

**Flow rate:** 3

**Injection volume:** 10

**Detector:** UV 254 and 238

---

#### CHROMATOGRAM

**Retention time:** 6.0

**Internal standard:** progesterone (11.6)

---

#### OTHER SUBSTANCES

**Extracted:** beclomethasone alcohol, beclomethasone 17-monopropionate, beclomethasone dipropionate

---

#### KEY WORDS

SPE; 17-hydroxyprogesterone is IS

---

#### REFERENCE

Levine,D.S.; Raisys,V.A.; Ainardi,V. Coating of oral beclomethasone dipropionate capsules with cellulose acetate phthalate enhances delivery of topically active antiinflammatory drug to the terminal ileum, *Gastroenterology*, **1987**, *92*, 1037–1044.

---

#### SAMPLE

**Matrix:** solutions

---

#### HPLC VARIABLES

**Column:** 250  $\times$  4.6 5  $\mu$ m RP-18 C18 (Brownlee)

**Mobile phase:** MeCN:MeOH:water 30:30:40

**Injection volume:** 20

**Detector:** UV 254

---

## KEY WORDS

for 17 $\alpha$ -hydroxyprogesterone

---

## REFERENCE

Kane,M.P.; Tsuji,K. Radiolytic degradation scheme for <sup>60</sup>Co-irradiated corticosteroids, *J.Pharm.Sci.*, **1983**, *72*, 30–35.

---

## SAMPLE

**Matrix:** solutions

---

## HPLC VARIABLES

**Column:** 150  $\times$  6.5  $\mu$ m Shim-pack CLC-ODS

**Mobile phase:** MeOH:THF:water 26:18:56

**Column temperature:** 48

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 240

---

## CHROMATOGRAM

**Retention time:** 14.6 (17 $\alpha$ )

---

## OTHER SUBSTANCES

**Simultaneous:** cortisone, estriol, cortisol, corticosterone, 11-deoxycortisol, androstenedione, prednisone acetate, 11-deoxycorticosterone, testosterone, dexamethasone acetate, estradiol, estrone, progesterone

---

## REFERENCE

Wei,J.Q.; Wei,J.L.; Zhou,X.T. Optimization of an isocratic reversed phase liquid chromatographic system for the separation of fourteen steroids using factorial design and computer simulation, *Biomed.Chromatogr.*, **1990**, *4*, 3–38.

---

## SAMPLE

**Matrix:** solutions

**Sample preparation:** Prepare a 25  $\mu$ g/mL solution in mobile phase, inject an aliquot.

---

## HPLC VARIABLES

**Column:** 250  $\times$  4.6 Partisil 10 ODS-1

**Mobile phase:** MeOH:water 55:45

**Column temperature:** 40

**Flow rate:** 1.5

**Detector:** UV 240

---

## CHROMATOGRAM

**Retention time:** k' 2.175

---

## OTHER SUBSTANCES

**Also analyzed:** androsterone (UV 210), cortexolone (UV 240), cortisone (UV 240), estradiol (UV 280), estrone (UV 280), ethinyl estradiol (UV 280), ethisterone (UV 240), hydrocortisone (UV 240), lynestrenol (UV 210), medroxyprogesterone acetate (UV 240), medroxyprogesterone (UV 240), methandienone (UV 240), methylandrostenediol (UV 210), methylprednisolone acetate (UV 240), methylprednisolone (UV 240), methyltestosterone (UV 240), nandrolone (UV 240), norethisterone (UV 240), prednisolone acetate (UV 240), prednisolone (UV 240), prednisone (UV 240), pregnenolone (UV 210), progesterone (UV 240), testosterone (UV 240)

---

## REFERENCE

Sadlej-Sosnowska,N. Structure retention relationship for steroid hormones. Functional groups as structural descriptors, *J.Liq.Chromatogr.*, **1994**, *17*, 2319–2330.

---

**SAMPLE****Matrix:** solutions**Sample preparation:** Inject 20  $\mu\text{L}$  aliquot of a MeOH solution.

---

**HPLC VARIABLES****Column:** 250  $\times$  4.6 5  $\mu\text{m}$  Hypersil 5-ODS**Mobile phase:** THF:water 23:77**Column temperature:** 30**Flow rate:** 1**Injection volume:** 20**Detector:** UV 245

---

**CHROMATOGRAM****Retention time:**  $k'$  11.86 (11 $\alpha$ -hydroxyprogesterone)**Internal standard:** methylprednisolone ( $k'$  11.36)

---

**OTHER SUBSTANCES****Simultaneous:** metabolites, betamethasone, corticosterone, cortisone, deflazacort, deoxycorticosterone, dexamethasone, fludrocortisone, fludrocortisone acetate, fluorocortisone, fluorocortisone acetate, hydrocortisone, 21-hydroxydeflazacort, methylprednisolone, prednisolone, prednisone, triamcinolone acetonide, triamcinolone

---

**REFERENCE**

Santos-Montes,A.; Gonzalo-Lumbreras,R.; Gasco-Lopez,A.I.; Izquierdo-Hornillos,R. Extraction and high-performance liquid chromatographic separation of deflazacort and its metabolite 21-hydroxydeflazacort. Application to urine samples, *J.Chromatogr.B*, **1994**, 657, 248–253.

---

**SAMPLE****Matrix:** solutions

---

**HPLC VARIABLES****Column:** 250  $\times$  4.6 5  $\mu\text{m}$  Nucleosil phenyl**Mobile phase:** Gradient. Carbon dioxide:MeOH from 98:2 to 78:22 over 40 min**Column temperature:** 50**Flow rate:** 2**Detector:** UV

---

**CHROMATOGRAM****Retention time:** 9.7

---

**OTHER SUBSTANCES****Simultaneous:** testosterone, estradiol, norethisterone, hydrocortisone, estriol, other steroids

---

**KEY WORDS**

SFC; 200 bar

---

**REFERENCE**

Hanson,M. Aspects of retention behaviour of steroids in packed column supercritical fluid chromatography, *Chromatographia*, **1995**, 40, 58–68.

---

**SAMPLE****Matrix:** tissue**Sample preparation:** Extract 70-125 mg tissue four times with 5 mL portions of ether:chloroform 80:20. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 100  $\mu\text{L}$  MeOH, inject an aliquot.

---

**HPLC VARIABLES****Column:** 80 mm long 10  $\mu\text{m}$  octadecylsilane radial compression (Radial-Pak) (Waters)**Mobile phase:** Gradient. A was MeOH:water 50:50. B was MeOH. A:B from 100:0 to 70:30 over 20 min, to 0:100 over 20 min

**Flow rate:** 2  
**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 25 (17-hydroxyprogesterone)

---

**OTHER SUBSTANCES**

**Extracted:** androstenedione, deoxycortisol, hydrocortisone, testosterone

**Simultaneous:** estriol, estradiol, pregnenolone, progesterone, testosterone enanthate, testosterone propionate

---

**KEY WORDS**

tumor

---

**REFERENCE**

Kessler, M.J. Analysis of steroids from normal and tumor tissue by HPLC, *Clin.Chim.Acta*, **1982**, *125*, 21–30.

---

**SAMPLE**

**Matrix:** urine

**Sample preparation:** 3 mL Urine + 1.5 µg betamethasone + 100 mg K<sub>2</sub>HPO<sub>4</sub> + 500 mg anhydrous sodium sulfate + 5 mL diethyl ether, shake mechanically for 10 min, centrifuge at 2500 g for 5 min. Remove the organic layer and evaporate it to dryness under vacuum, reconstitute the residue in 200 µL MeOH, filter (0.45 µm), inject a 15 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 100 × 4.6 5 µm Hypersil ODS

**Mobile phase:** Gradient. MeCN:water from 4:96 to 30:70 over 10 min, to 45:55 over 5 min, to 50:50 over 3 min

**Column temperature:** 40

**Flow rate:** 1

**Injection volume:** 15

**Detector:** UV 246

---

**CHROMATOGRAM**

**Retention time:** 14.66

**Internal standard:** betamethasone (12.83)

**Limit of detection:** 10 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** corticosterone, cortisone, deoxycorticosterone, hydrocortisone, prednisolone, prednisone, triamcinolone, triamcinolone acetonide

---

**REFERENCE**

Park, S.-J.; Kim, Y.-J.; Pyo, H.-S.; Park, J. Analysis of corticosteroids in urine by HPLC and thermospray LC/MS, *J.Anal.Toxicol.*, **1990**, *14*, 102–108.

---

**SAMPLE**

**Matrix:** urine

**Sample preparation:** 3 mL Urine + 0.25 g NaCl, adjust pH to 9.0 with 0.5 g Na<sub>2</sub>HPO<sub>4</sub>, add 4 mL dichloromethane, vortex 1 min, centrifuge at 3700 g for 3 min. Remove organic phase and dry it over anhydrous sodium sulfate. Evaporate a 3 mL aliquot to dryness under vacuum, reconstitute residue with 200 µL 5 µg/mL IS in MeOH, inject a 20 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 Hypersil ODS

**Mobile phase:** MeCN:water 32:68

**Column temperature:** 30

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 245

---

**CHROMATOGRAM****Retention time:** 19**Internal standard:** methylprednisolone (9)

---

**OTHER SUBSTANCES**

**Simultaneous:** triamcinolone, triamcinolone acetonide, prednisolone, corticosterone, hydroxyprogesterone, fluorocortisone acetate, cortisone, hydrocortisone, fluorocortisone, betamethasone, dexamethasone, prednisone

---

**KEY WORDS**

SPE also discussed

---

**REFERENCE**

Santos-Montes,A.; Gonzalo-Lumbreras,R.; Gasco-Lopez,A.I.; Izquierdo-Hornillos,R. Solvent and solid-phase extraction of natural and synthetic corticoids in human urine, *J.Chromatogr.B*, **1994**, 652, 83–89.

---

# Hydroxypropyl methylcellulose

**CAS Registry No.:** 9004-65-3**Merck Index:** 4889

---

**SAMPLE****Matrix:** formulations

**Sample preparation:** Dilute with an equal volume of MeOH:water 40:60, inject an 80  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 300  $\times$  7.8 Ultrahydrogel 250 250 Å cross-linked methacrylate gel (waters)

**Mobile phase:** MeOH:buffer 20:80 (Buffer was 95 mM boric acid containing 10 mM KCl, 13  $\mu$ M sodium borate, and 1.5 mM dextrose.) (At the end of each set of analyses rinse instrument and column with 0.5% sodium azide. (Caution! Sodium azide is carcinogenic, highly toxic, and can form explosive heavy metal azides. Do not discharge to the plumbing system!))

**Column temperature:** 35**Flow rate:** 1**Injection volume:** 80**Detector:** RI

---

**CHROMATOGRAM****Retention time:** 10.6

---

**OTHER SUBSTANCES****Simultaneous:** PEG 400

---

**KEY WORDS**


ophthalmic solutions; SEC

---

**REFERENCE**

Delker,G.; Chen,C.; Miller,R.B. Size exclusion chromatographic determination of hydroxypropyl methyl cellulose and polyethylene glycol 400 in an ophthalmic solution, *Chromatographia*, **1995**, 41, 263–266.

# Hydroxyquinoline



**Molecular formula:** C<sub>9</sub>H<sub>7</sub>NO

**Molecular weight:** 145.16

**CAS Registry No.:** 148-24-3

**Merck Index:** 4890

## SAMPLE

**Matrix:** solutions

## HPLC VARIABLES

**Column:** 250 × 4.6 Zorbax RX

**Mobile phase:** Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

**Column temperature:** 30

**Flow rate:** 2

**Detector:** UV 210

## OTHER SUBSTANCES

**Also analyzed:** acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitrityline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenpropfen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, keto-profen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepa-crine, meperidine, mephentermine, mephénytoin, mephesin, mephobarbital, mepivacaine, mes-caline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, metha-zolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methylprylon, me-toprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, na-proxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepi-nephine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbi-tal, persantine, phenacetin, phenazocine, phenazopyridine, phenacyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenyl-butazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primi-done, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopola-mine, scopoletin, secobarbital, strychnine, sulfacetamide, sufadiazine, sulfadimethoxine, sul-faethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sul-fasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tol-



metin, tranlylcpromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleennamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine


---

**REFERENCE**

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.

---

# Hydroxyurea



**Molecular formula:** CH<sub>4</sub>N<sub>2</sub>O<sub>2</sub>

**Molecular weight:** 76.06

**CAS Registry No.:** 127-07-1

**Merck Index:** 4896

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Mix 500 µL serum, 20 µL 70% perchloric acid, and 20 µL 15.5 mM methylurea, centrifuge at 12000 g for 5 min. Mix 200 µL supernatant with 700 µL BUN acid reagent (No. 535-3, Sigma) and 600 µL BUN color reagent (No. 535-5, Sigma), place the mixture in a boiling water bath to form colored complexes. After 10 min cool the mixture in ice water, inject a 25-100 µL aliquot of the colored solution.

---

**HPLC VARIABLES**

**Column:** 300 × 3.9 Bondclone 10 C18 (Phenomenex)

**Mobile phase:** MeCN:water 13:87

**Flow rate:** 1.7

**Injection volume:** 25-100

**Detector:** UV 449

---

**CHROMATOGRAM**

**Retention time:** 6.5

**Internal standard:** methylurea (12.2)

---

**KEY WORDS**

serum; derivatization

---

**REFERENCE**

Manouilov,K.K.; McGuire,T.R.; Gwilt,P.R. Colorimetric determination of hydroxyurea in human serum using high-performance liquid chromatography, *J.Chromatogr.B*, **1998**, *708*, 321-324.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Vortex plasma with two volumes of 5% trichloroacetic acid for 30 s, centrifuge at 7000 g for 5 min, inject a 5 µL aliquot of the supernatant.

---

**HPLC VARIABLES**

**Column:** 100 × 3.2 3 µm Phase-II ODS (Bioanalytical Systems)

**Mobile phase:** 50 mM sodium acetate containing 5 mM tetrabutylammonium hydroxide, adjusted to pH 6.75 ± 0.02 with 50 mM acetic acid

**Flow rate:** 0.5

**Injection volume:** 5

**Detector:** E, Bioanalytical Systems LC-4B, glassy carbon electrode + 800 mV, Ag/AgCl reference electrode

---

**CHROMATOGRAM**

**Retention time:** 1.1 (elute for 5 min to remove plasma components)

**Limit of detection:** 20  $\mu$ M

---

**OTHER SUBSTANCES**

**Noninterfering:** acetaminophen, caffeine, carbamazepine, ethosuximide, morphine, phenobarbital, phenytoin, salicylic acid, valproic acid

---

**KEY WORDS**

plasma

---

**REFERENCE**

Havard,J.; Grygiel,J.; Sampson,D. Determination by high-performance liquid chromatography of hydroxyurea in human plasma, *J.Chromatogr.*, **1992**, 584, 270–274.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 100  $\mu$ L Serum + 200  $\mu$ L trichloroacetic acid, vortex for 30 s, centrifuge at 5000 g for 5 min, pass the supernatant through a 1 mL Bond Elut C18 SPE cartridge, inject an aliquot of the eluate.

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Supelguard LC 18 (Supelco)

**Column:** 150  $\times$  4.6 Supelcosil LC 18

**Mobile phase:** 50 mM Sodium acetate adjusted to pH 6.75 with acetic acid

**Flow rate:** 0.8

**Detector:** E, ESA Coulochem II, Model 5020 guard cell 650 mV, Model 5011 analytical cell 600 mV

---

**CHROMATOGRAM**

**Limit of detection:** 600 nM

**Limit of quantitation:** 1.3  $\mu$ M

---

**KEY WORDS**

pharmacokinetics; serum; SPE


---

**REFERENCE**

Villani,P.; Maserati,R.; Regazzi,M.B.; Giacchino,R.; Lori,F. Pharmacokinetics of hydroxyurea in patients infected with human immunodeficiency virus type I, *J.Clin.Pharmacol.*, **1996**, 36, 117–121.

---

# Hydroxyzine



**Molecular formula:**  $C_{21}H_{27}ClN_2O_2$

**Molecular weight:** 374.91

**CAS Registry No.:** 68-88-2, 2192-20-3 (di HCl), 10246-75-0 (pamoate)

**Merck Index:** 4897

**Lednicer No.:** 1 59, 4 118

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 1 mL Serum + 100  $\mu$ L 1  $\mu$ g/mL triprolidine + 250  $\mu$ L 10% KOH + 5 mL ether, vortex, centrifuge. Remove ether layer and add it to 100  $\mu$ L 0.5% phosphoric acid, vortex, centrifuge, remove most of ether layer and discard it, remove traces of ether by nitrogen at room temperature for 2-3 min, inject all of aqueous layer.

---

**HPLC VARIABLES**

**Column:** Waters CN reverse-phase radial compression

**Mobile phase:** MeCN:buffer 27:73 (Buffer was 75 mM pH 3.0 phosphate buffer containing 20 mM dibutylamine and 50 ng/mL triprolidine.)

**Flow rate:** 1

**Injection volume:** 100

**Detector:** UV 229

---

#### CHROMATOGRAM

**Retention time:** 6.9

**Internal standard:** triprolidine (3.6)

**Limit of detection:** 3 ng/mL

---

#### KEY WORDS

serum

---

#### REFERENCE

Simons,F.E.R.; Simons,K.J.; Frith,E.M. The pharmacokinetics and antihistaminic of the H<sub>1</sub> receptor antagonist hydroxyzine, *J.Allerg.Clin.Immunol.*, **1984**, 73, 69–75.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** 1 mL Serum + 50  $\mu$ L 3  $\mu$ g/mL IS + 1 mL 1 M pH 5.0 sodium citrate buffer + 3 mL ethyl acetate, vortex 1 min, centrifuge at 4000 rpm for 15 min, remove organic layer, repeat extraction. Combine organic layers and add 200  $\mu$ L 1.7% phosphoric acid, vortex 1 min, centrifuge 5 min, remove and discard ethyl acetate layer, remove traces of ethyl acetate from aqueous layer using a stream of nitrogen, inject.

---

#### HPLC VARIABLES

**Column:** radial 4  $\mu$ m NovoPak C18 radial compression

**Mobile phase:** MeCN:buffer 46:54 (Buffer was 10 mM pH 2.9 KH<sub>2</sub>PO<sub>4</sub> + 20 mM sodium 1-decanesulfonate.)

**Flow rate:** 1.4

**Detector:** UV 229

---

#### CHROMATOGRAM

**Retention time:** 3.6

**Internal standard:** P-265, an ethoxy derivative of cetirizine (6.8)

**Limit of detection:** 1 ng/mL

---

#### OTHER SUBSTANCES

**Simultaneous:** cetirizine

---

#### KEY WORDS

serum

---

#### REFERENCE

Simons,K.J.; Watson,W.T.A.; Chen,X.Y.; Simons,F.E.R. Pharmacokinetic and pharmacodynamic studies of the H<sub>1</sub>-receptor antagonist hydroxyzine in the elderly, *Clin.Pharmacol.Ther.*, **1989**, 45, 9–14.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol:n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100  $\mu$ L mobile phase, centrifuge at 2800 g for 5 min, inject a 50  $\mu$ L aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

---

#### HPLC VARIABLES

**Column:** 300  $\times$  3.9 4  $\mu$ m NovaPack C18

**Mobile phase:** MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH<sub>2</sub>PO<sub>4</sub> adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

**Column temperature:** 30**Flow rate:** 0.8**Injection volume:** 50**Detector:** UV 230**CHROMATOGRAM****Retention time:** 8.62**Limit of detection:** <120 ng/mL**KEY WORDS**

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfapyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procabazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; ciprolool; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vindesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benzepiril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; caripramine; thioridazine; feniazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

**REFERENCE**

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254–262.

**SAMPLE****Matrix:** blood, CSF

**Sample preparation:** Plasma. Centrifuge blood at 7000 rpm, decant 100  $\mu$ L plasma. Mix 100  $\mu$ L plasma with 200  $\mu$ L acetone, centrifuge at 7000 rpm for 5 min. Evaporate the supernatant under a stream of nitrogen, reconstitute the residue with mobile phase, inject an aliquot. CSF. Add 25  $\mu$ L water to 25  $\mu$ L CSF, mix with 50  $\mu$ L acetone, centrifuge at 7000 rpm for 5 min. Evaporate the supernatant under a stream of nitrogen, reconstitute the residue with mobile phase, inject an aliquot.

**HPLC VARIABLES****Column:** Regis SPS phenyl**Mobile phase:** MeCN: pH 3.0 water 21:79 containing 9 mM decanesulfonic acid**Column temperature:** 60

---

**KEY WORDS**

plasma; rat; pharmacokinetics

---

**REFERENCE**

Chou, K.-J.; Donovan, M.D. Distribution of antihistamines into the CSF following intranasal delivery, *Bio-pharm. Drug Dispos.*, **1997**, *18*, 335–346.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50  $\mu$ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood)  $\mu$ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18

**Column:** 250  $\times$  4.6 5  $\mu$ m Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10–30

**Detector:** UV 200.5

---

**CHROMATOGRAM**

**Retention time:** 15.267

---

**KEY WORDS**

whole blood

---

**REFERENCE**

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, *763*, 149–163.

---

**SAMPLE**

**Matrix:** bulk

**Sample preparation:** Prepare solutions in mobile phase, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 120  $\times$  4.0 5  $\mu$ m Lichrosorb RP-18

**Mobile phase:** MeCN:buffer 30:70 (Buffer was 9.8% triethylamine and 0.35% sodium methane-sulfonate with pH adjusted to 2.85 with sulfuric acid.)

**Flow rate:** 1

**Detector:** UV 230

---

**CHROMATOGRAM**

**Retention time:** 12.3

**Limit of detection:** 180 ng/mL

**Limit of quantitation:** 600 ng/mL

---

**OTHER SUBSTANCES**

**Simultaneous:** impurities

---

**REFERENCE**

Simpson,D.; Skellern,G.G.; Miller,J.H.MB. Method for the control of known impurities in hydroxyzine hydrochloride, *J.Pharm.Biomed.Anal.*, **1996**, 14, 1371-1375.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** 5 mL Formulation + 5 mL IS solution, make up to 50 mL with MeOH, inject 10  $\mu$ L aliquot. (IS solution was 0.2 mg p-nitroacetophenone and 2.5 mg isobutyrophenone per mL of MeOH.)

---

**HPLC VARIABLES**

**Column:** 300  $\times$  4 10  $\mu$ m  $\mu$ Bondapak C18

**Mobile phase:** MeCN:water:reagent 25:60:15, pH 2.6 (Reagent was MeOH containing 0.06% sulfuric acid, 0.5% sodium sulfate, and 0.02% sodium heptanesulfonate.)

**Flow rate:** 2

**Injection volume:** 10

**Detector:** UV 257

---

**CHROMATOGRAM**

**Retention time:** 9

**Internal standard:** p-nitroacetophenone (6) and isobutyrophenone (13)

---

**OTHER SUBSTANCES**

**Simultaneous:** benzyl alcohol, benzoic acid, benzaldehyde, p-chlorobenzoic acid, p-chlorobenzaldehyde, p-chlorobenzophenone

---

**KEY WORDS**

injections; stability-indicating

---

**REFERENCE**

Menon,G.N.; Norris,B.J. Simultaneous determination of hydroxyzine hydrochloride and benzyl alcohol in injection solutions by high-performance liquid chromatography, *J.Pharm.Sci.*, **1981**, 70, 697-698.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a 10  $\mu$ g/mL solution in MeOH, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 125  $\times$  4.9 Spherisorb S5W silica

**Mobile phase:** MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

**Flow rate:** 2

**Injection volume:** 20

**Detector:** E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

---

**CHROMATOGRAM**

**Retention time:** 2.2

---

**OTHER SUBSTANCES**

**Also analyzed:** acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzoctamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazepine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipamnone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, er-

gosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserlin, laudanose, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenamproide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenylglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranylecypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

---

## REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, *323*, 191-225.

---

## SAMPLE

**Matrix:** solutions

**Sample preparation:** Dissolve in MeOH:water 1:1 at a concentration of 50 µg/mL, inject a 10 µL aliquot.

---

## HPLC VARIABLES

**Column:** 300 × 3.9 10 µm µBondapak C18

**Mobile phase:** MeOH:acetic acid:triethylamine:water 58:1.5:0.5:40

**Flow rate:** 1.5

**Injection volume:** 10

**Detector:** UV 230

---

## CHROMATOGRAM

**Retention time:** 12

---

## REFERENCE

Roos, R. W.; Lau-Cam, C. A. General reversed-phase high-performance liquid chromatographic method for the separation of drugs using triethylamine as a competing base, *J. Chromatogr.*, **1986**, *370*, 403-418.

---

## SAMPLE

**Matrix:** solutions

**Sample preparation:** Prepare a 1 mg/mL solution in MeOH, inject a 5 µL aliquot.

---

## HPLC VARIABLES

**Column:** 250 × 4.6 5 µm Lichrosphere cyanopropyl

**Mobile phase:** Carbon dioxide:MeOH:isopropylamine 90:10:0.05

**Column temperature:** 50

**Flow rate:** 3

**Injection volume:** 5

**Detector:** UV 220

---

## CHROMATOGRAM

**Retention time:** 2.42

**OTHER SUBSTANCES**

**Simultaneous:** benactyzine, buclizine, perphenazine, thioridazine, amitriptyline, desipramine, imipramine, nortriptyline, protriptyline

**KEY WORDS**

SFC; pressure 200 bar

**REFERENCE**

Berger, T.A.; Wilson, W.H. Separation of drugs by packed column supercritical fluid chromatography. 2. Anti-depressants, *J.Pharm.Sci.*, **1994**, *83*, 287-290.

**SAMPLE**

**Matrix:** solutions

**HPLC VARIABLES**

**Column:** 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4 5 µm LiChrospher 100 RP-8 (B)

**Mobile phase:** MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

**Flow rate:** 0.6

**Injection volume:** 25

**Detector:** UV 229

**CHROMATOGRAM**

**Retention time:** 11.40 (A), 6.27 (B)

**OTHER SUBSTANCES**

**Also analyzed:** acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlorthalidone, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrizamide, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazindol, mefenamic acid, meperidine, mephentermine, mepivacaine, mesoridazine, metaproterenol, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfapyrazole, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocinamide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, trifluorpromazine, trimetoprim, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

**KEY WORDS**

details of plasma extraction



---

**REFERENCE**

Koves,E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J.Chromatogr.A*, **1995**, 692, 103–119.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Inject a 20  $\mu\text{L}$  aliquot of a 100–500  $\mu\text{g/mL}$  solution in mobile phase.

---

**HPLC VARIABLES**

**Column:** 100  $\times$  4.6 5  $\mu\text{m}$  Hypersil C8 MOS 100A coated with phosphatidylcholine (95% pure soybean lecithin, Epikuron, Lucas Meyer & Co.) (Coat column by recycling a 1 mM solution of phosphatidylcholine in MeOH:water 80:20 for 24 h.)

**Mobile phase:** MeCN:35 mM pH 7.4 sodium phosphate buffer 40:60

**Flow rate:** 0.5–2

**Injection volume:** 20

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** k' 8.32

---

**OTHER SUBSTANCES**

**Also analyzed:** amoxicillin, antipyrine, carbamazepine, chlorpheniramine, chlorpromazine, clonidine, codeine, desipramine, diphenhydramine, dipyridamole, ephedrine, flufenamic acid, haloperidol, imipramine, indomethacin, lidocaine, megestrol acetate, metoprolol, nabumetone, naldolol, phenobarbital, phenol, promazine, propranolol, pyrilamine, quinidine, ropinirole, testosterone, thioridazine, tolfenamic acid, verapamil

**Noninterfering:** acetaminophen, aspirin, azathioprine, caffeine, carprofen, chlorambucil, cimetidine, fenoterol, flurbiprofen, ibuprofen, ketoprofen, ranitidine, salicylic acid, sulfamethoxazole, theophylline, thioguanine, tiaprofenic acid, trimethoprim, valproic acid

---

**KEY WORDS**

comparison with capillary electrophoresis

---

**REFERENCE**

Hanna,M.; de Biasi,V.; Bond,B.; Salter,C.; Hutt,A.J.; Camilleri,P. Estimation of the partitioning characteristics of drugs: A comparison of a large and diverse drug series utilizing chromatographic and electrophoretic methodology, *Anal.Chem.*, **1998**, 70, 2092–2099.

---

---

# Ibafloxacin


**Molecular formula:**  $\text{C}_{15}\text{H}_{14}\text{FNO}_3$

**Molecular weight:** 275.28

**CAS Registry No.:** 91618-36-9

**Merck Index:** 4919

---



---

**SAMPLE**

**Matrix:** tissue

**Sample preparation:** Mix 2 g minced tissue with 4 g anhydrous sodium sulfate until homogenized. Add 10 ml ethyl acetate, shake mechanically for 10 min, centrifuge at 1500 g for 10 min. Transfer organic layer into another tube and repeat extraction on the tissue pellet with 10 mL ethyl acetate. Evaporate combined organic phases under a stream of nitrogen at 50°. Dissolve residue in 1 mL MeCN:2.7 mM pH 2.5 oxalic acid 50:50, vortex, sonicate for 5 min, filter through 0.45  $\mu\text{m}$  filter (GHP Acrodisc GF, Gelman Sciences, USA).. Inject a 100  $\mu\text{L}$  aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu\text{m}$  Ultrabase C18 (Shandon, UK)