

REFERENCE

Jane,I.; McKinnon,A.; Flanagan,R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J.Chromatogr.*, **1985**, 323, 191–225.

Indeloxazine hydrochloride


Molecular formula: $C_{14}H_{18}ClNO_2$

Molecular weight: 267.76

CAS Registry No.: 65043-22-3

Merck Index: 4972

Lednicer No.: 4 59

**SAMPLE**

Matrix: blood

Sample preparation: 1 mL Plasma + 1 mL 300 ng/mL viloxazine in water + 500 μ L 200 mM ammonium hydroxide + 4 mL ether, extract, centrifuge. Remove the organic layer and evaporate it to dryness at 45°, reconstitute the residue in 100 μ L 3 mg/mL sodium bicarbonate and 200 μ L 125 μ g/mL dansyl chloride in acetone, heat at 45° for 20 min, cool, add 4 mL ether. Wash the ether solution twice with 3 mL water and evaporate it to dryness at 45°. Take up the residue in 100 μ L n-heptane, inject a 3–5 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.5 μ m LiChrosorb SI-60

Mobile phase: n-Heptane:ethyl acetate 20:3

Flow rate: 1.5

Injection volume: 3–5

Detector: F ex 365 em 505

CHROMATOGRAM

Retention time: 3

Internal standard: viloxazine (4)

Limit of detection: 5 ng/mL

KEY WORDS

plasma; normal phase; derivatization

REFERENCE

Kamimura,H.; Sasaki,H.; Yokoi,K.; Kawamura,S. Determination of indeloxazine in plasma by liquid chromatography and gas chromatography-mass spectrometry, *J.Pharm.Sci.*, **1985**, 74, 559–561.


Indinavir

Molecular formula: $C_{36}H_{47}N_5O_4$

Molecular weight: 613.80

CAS Registry No.: 150378-17-9, 157810-81-6 (sulfate)

Merck Index: 4979

**SAMPLE**

Matrix: blood

Sample preparation: 300 μ L Plasma + 300 μ L 50 mM pH 9.0 ammonium dihydrogen phosphate buffer + 30 μ L 10.5 μ g/mL IS in water, vortex for 10 s, add 3 mL diethyl ether, vortex for 30

s, keep at -20° for 30 min. Remove the ether layer and evaporate it to dryness under nitrogen. Reconstitute the residue with 150 μ L 10 mM pH 5.5 ammonium dihydrogen phosphate buffer, centrifuge at 750 g for 5 min, inject a 35 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 3.9 5 μ m Delta-pak C4 (Waters)

Mobile phase: MeCN:buffer 35:65 (Buffer was 10 mM ammonium dihydrogen phosphate and 1 mM 1-heptanesulfonic acid sodium salt, pH adjusted to 4.8 with ammonium hydroxide.)

Flow rate: 0.6

Injection volume: 35

Detector: UV 210

CHROMATOGRAM

Retention time: 13.2 \pm 0.5

Internal standard: methylindinavir (methyl at 5 position on pyridine) (14.1 \pm 0.7)

Limit of quantitation: 10 ng/mL

OTHER SUBSTANCES

Simultaneous: nelfinavir, ritonavir, saquinavir

Noninterfering: didanosine, lamivudine, stavudine, zalcitabine, zidovudine

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Iayewardene, A.L.; Zhu, F.; Aweeka, F.T.; Gambertoglio, J.G. Simple high-performance liquid chromatographic determination of the protease inhibitor indinavir in human plasma, *J. Chromatogr. B*, **1998**, 707, 203–211.

SAMPLE

Matrix: blood

Sample preparation: Add 50 μ L 1 μ g/mL IS in MeCN:water 50:50 to 1 mL plasma and vortex. Add 1 mL 100 mM pH 9.5 borate buffer, vortex, add 8 mL isopropanol:chloroform 1:15 (Caution! Chloroform is a carcinogen!), mix on a flat-bed shaker for 15 min, centrifuge at 1500 g for 5 min. Remove the lower organic layer and evaporate it under a stream of nitrogen at 45°. Reconstitute the residue in 1 mL mobile phase, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 50 \times 3.9 5 μ m Symmetry C8

Mobile phase: MeCN:water 40:60 containing 2 mM ammonium acetate

Flow rate: 1

Injection volume: 50

Detector: MS, PE-SCIEX API IIIplus triple quadrupole, heated nebulizer, corona discharge needle (+4.5 μ A), nebulizer 500°, collision gas argon at 260×10^{12} atoms/cm², nebulizing gas nitrogen at 80 psi and 0.6 mL/min, orifice +75 V, electron multiplier -4.0 kV, dwell time 400 ms with a 30 ms pause time between scans, Q1 at m/z 523, Q2 at m/z 512, Q3 at m/z 273

CHROMATOGRAM

Retention time: 1.5

Internal standard: structure given in paper

Limit of quantitation: 5 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma

REFERENCE

Woolf, E.; Haddix, H.M.; Matuszewski, B. Determination of an in vivo metabolite of a human immunodeficiency virus protease inhibitor in human plasma by high-performance liquid chromatography with tandem mass spectrometry, *J. Chromatogr. A*, **1997**, 762, 311–319.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 50 μ L MeCN:water 50:50 + 25 μ L 1 μ g/mL IS in MeCN, vortex, add 1 mL 100 mM pH 9.5 borate buffer, vortex, add 8 mL MTBE, mix on a Glas Col (Terre Haute, IN) RD 350 rotator for 15 min, centrifuge at 1500 g for 5 min, freeze the aqueous phase in a dry ice-acetone bath. Decant the organic layer and evaporate it under a stream of nitrogen at 42°. Reconstitute the residue in 175 μ L mobile phase, inject a 6 μ L aliquot.

HPLC VARIABLES

Column: 50 \times 2.3 μ m BDS Hypersil C8

Mobile phase: MeCN:water 40:60 containing 7 mM pH ammonium acetate, adjusted to pH 4.9 with formic acid

Flow rate: 0.2

Injection volume: 6

Detector: MS, PE-Sciex API IIIplus triple quadrupole, turbo-ion spray, turbo probe at 500°, auxiliary gas nitrogen, nebulizing gas nitrogen at 80 psi, positive ion mode, interface sprayer at +4 kV, sampling orifice at +60 V, m/z 614

CHROMATOGRAM

Retention time: 2.6

Internal standard: indinavir analog (4.2)

Limit of quantitation: 1 ng/mL

OTHER SUBSTANCES

Extracted: hexadeuterated indinavir (m/z 620)

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Woolf, E.J.; Matuszewski, B.K. Simultaneous determination of unlabeled and deuterium-labeled indinavir in human plasma by high-performance liquid chromatography with tandem mass spectrometric detection, *J.Pharm.Sci.*, **1997**, 86, 193–198.

SAMPLE

Matrix: microsomal incubations

Sample preparation: Add four volumes of MeCN to microsomal incubation, centrifuge at 1500 g for 10 min, evaporate supernatant to dryness under nitrogen at 40°, reconstitute the residue in 100–120 μ L mobile phase, inject an 80 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 Hypersil 5C18 (Phenomenex)

Mobile phase: Gradient. A was MeCN. B was 0.3% pH 6 triethylamine in water. A:B from 25:75 to 45:55 over 18 min

Flow rate: 1.5

Injection volume: 80

Detector: UV 240

CHROMATOGRAM

Retention time: 16.1

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

human; liver

REFERENCE

Koudriakova, T.; Iatsimirskaia, E.; Utkin, I.; Gangl, E.; Vouros, P.; Storozhuk, E.; Orza, D.; Marinina, J.; Gerber, N. Metabolism of the human immunodeficiency virus protease inhibitors indinavir and zidovudine by human intestinal microsomes and expressed cytochrome P4503A4/3A5: Mechanism-based inactivation of cytochrome P4503A by zidovudine, *Drug Metab.Dispos.*, **1998**, 26, 552–561.

SAMPLE

Matrix: microsomal incubations

Sample preparation: Add four volumes of 1-chlorobutane containing 50 ng/mL 5-methoxypsoralen to microsomal incubation, mix vigorously, centrifuge at 1500 g for 10 min, evaporate 1-chlorobutane extract to dryness under nitrogen, reconstitute the residue in 100 μ L mobile phase, inject an 80 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 Hypersil 5C18 (Phenomenex)

Mobile phase: Gradient. A was MeCN. B was 0.3% pH 6 triethylamine in water. A:B from 25:75 to 45:55 in 18 min

Flow rate: 1.5

Injection volume: 80

Detector: UV 240

CHROMATOGRAM

Retention time: 16.1

Internal standard: 5-methohypsoralen (14.6)

Limit of detection: 30 nM

OTHER SUBSTANCES

Extracted: metabolites


KEY WORDS

human; liver

REFERENCE

Koudriakova,T.; Iatsimirskaia,E.; Utkin,I.; Gangl,E.; Vouros,P.; Storozhuk,E.; Orza,D.; Marinina,J.; Gerber,N. Metabolism of the human immunodeficiency virus protease inhibitors indinavir and zidovudine by human intestinal microsomes and expressed cytochrome P4503A4/3A5: Mechanism-based inactivation of cytochrome P4503A by zidovudine, *Drug Metab.Dispos.*, **1998**, 26, 552–561.

Indobufen



Molecular formula: C₁₆H₁₇NO₃

Molecular weight: 295.34

CAS Registry No.: 63610-08-2

Merck Index: 4991

SAMPLE

Matrix: blood

Sample preparation: 0.5-1 mL Plasma + 500 μ L water + 100 μ L 600 mM sulfuric acid + 40 mg NaCl + 4 mL diethyl ether, extract on a rotamixer, centrifuge at 1200 g. Remove the organic layer and evaporate it to dryness under a stream of air at 30°, add 5 drops toluene, evaporate to dryness under a stream of air at 30°, reconstitute the residue in 200 μ L 50 mM triethylamine in MeCN, add 100 μ L 60 mM ethyl chloroformate in MeCN, after 30 s add 100 μ L 1 M l-leucinamide hydrochloride in MeOH containing 1 M triethylamine, after 2 min add 500 μ L 250 mM HCl, extract with 4 mL ethyl acetate. Evaporate the organic layer to dryness under a stream of air at 30°, reconstitute the residue with 500 μ L MeOH, inject a 10 μ L aliquot.

HPLC VARIABLES

Guard column: 30 \times 4 Perisorb RP-18 (Merck)

Column: 250 \times 4 7 μ m LiChroCart RP-18 (Merck)

Mobile phase: MeCN:10 mM pH 6.5 phosphate buffer 38:62

Flow rate: 2

Injection volume: 10

Detector: UV 275

CHROMATOGRAM**Retention time:** 6 (-), 8 (+) (assignment tentative)**Internal standard:** indobufen

OTHER SUBSTANCES**Extracted:** indoprofen

KEY WORDSplasma; derivatization; chiral; indobufen is IS

REFERENCEBjörkman, S. Determination of the enantiomers of indoprofen in blood plasma by high-performance liquid chromatography after rapid derivatization by means of ethyl chloroformate, *J. Chromatogr.*, **1985**, 339, 339–346.

SAMPLE**Matrix:** blood**Sample preparation:** 0.5–1 mL Plasma + 100 μ L 600 mM sulfuric acid + 40 mg NaCl + 4 mL diethyl ether, extract by rotation, centrifuge at 1200 g. Remove the organic layer and evaporate it to dryness under a stream of air at 30°, add 5 drops of toluene, evaporate to dryness under a stream of air. Reconstitute the residue in 200 μ L 50 mM triethylamine in MeCN, add 100 μ L 60 mM ethyl chloroformate in MeCN, let stand for 30 s, add 100 μ L 1 M leucinamide hydrochloride in MeOH containing 1 M triethylamine, let stand for 2 min, add 500 μ L 250 mM HCl, extract with 4 mL ethyl acetate. Remove the organic layer and evaporate it to dryness under a stream of air at 30°, reconstitute in 500 μ L MeOH, inject a 10 μ L aliquot (*J. Chromatogr.* 1985, 339, 339).

HPLC VARIABLES**Column:** 300 \times 3.9 10 μ m μ Bondapak C18**Mobile phase:** MeCN:pH 6.4 phosphate buffer 40:60**Flow rate:** 1.5**Injection volume:** 10**Detector:** UV 275

CHROMATOGRAM**Retention time:** 6.977 (R), 8.848 (S)

KEY WORDSchiral; derivatization

REFERENCEPerrone, G.; Farina, M. High-performance liquid chromatographic method for direct resolution of the indobufen enantiomeric components, *J. Chromatogr.*, **1990**, 520, 373–378.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 250 \times 4.6 10 μ m Chiralcel OD**Mobile phase:** Hexane:isopropanol:formic acid 80:20:0.5**Flow rate:** 1.5**Detector:** UV 270

CHROMATOGRAM**Retention time:** 8.20 (R), 10.91 (S)**Limit of detection:** 2 ng

KEY WORDSchiral; $\alpha = 1.33$

REFERENCE

Perrone,G.; Farina,M. High-performance liquid chromatographic method for direct resolution of the indobufen enantiomeric components, *J.Chromatogr.*, **1990**, 520, 373–378.

SAMPLE

Matrix: urine

Sample preparation: 500 μ L Urine + 50 μ L 5000 U/mL β -glucuronidase (Sigma) + 2 mL 100 mM pH 5 acetate buffer, heat at 37° for 16 h, add S-indoprofen, add 1 mL 1 M HCl, extract with diethyl ether. Remove the organic phase and extract it with 500 μ L 1 M NaOH. Remove the aqueous phase and add it to 1 mL 1 M HCl, extract with diethyl ether. Remove the organic layer and evaporate it to dryness under a stream of air at 30°, add 5 drops of toluene, evaporate to dryness under a stream of air. Reconstitute the residue in 200 μ L 50 mM triethylamine in MeCN, add 100 μ L 60 mM ethyl chloroformate in MeCN, let stand for 30 s, add 100 μ L 1 M leucinamide hydrochloride in MeOH containing 1 M triethylamine, let stand for 2 min, add 500 μ L 250 mM HCl, extract with 4 mL ethyl acetate. Remove the organic layer and evaporate it to dryness under a stream of air at 30°, reconstitute in 500 μ L MeOH, inject a 10 μ L aliquot (*J.Chromatogr.* 1985, 339, 339).

HPLC VARIABLES

Column: 125 \times 4 μ m Lichrocart Superspher (Merck)

Mobile phase: MeCN:10 mM pH 6.5 phosphate buffer 35:65

Injection volume: 10

Detector: UV 275

CHROMATOGRAM

Internal standard: indoprofen

Limit of detection: 5 μ g/mL

KEY WORDS

chiral; derivatization; pharmacokinetics

REFERENCE

Strolin Benedetti,M.; Frigerio,E.; Tamassia,V.; Nosedà,G.; Caldwell,J. The dispositional enantioselectivity of indobufen in man, *Biochem.Pharmacol.*, **1992**, 43, 2032–2034.

SAMPLE

Matrix: urine

Sample preparation: Inject an aliquot.

HPLC VARIABLES

Column: 250 \times 4 LiChrosorb C18

Mobile phase: Gradient. MeCN:water 25:75 containing 0.1% trifluoroacetic acid for 30 min then MeCN:water 40:60 containing 0.1% trifluoroacetic acid for 20 min (step gradient).

Flow rate: 2

Detector: F ex 290 em 440

CHROMATOGRAM

Retention time: 36.8

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

rat; mouse

REFERENCE

Grubb,N.; Caldwell,J.; Strolin-Benedetti,M. Excretion balance and urinary metabolism of indobufen in rats and mice, *Biochem.Pharmacol.*, **1993**, 46, 759–761.


Indocyanine green

Molecular formula: $C_{43}H_{47}N_2NaO_6S_2$

Molecular weight: 774.98

CAS Registry No.: 3599-32-4

Merck Index: 4992



SAMPLE

Matrix: blood

Sample preparation: 0.5 mL Plasma + 10 μ L 250 μ g/mL 1-acetamidopyrene in MeOH + 200 μ L 1 M ammonium sulfate + 800 μ L cold MeCN, vortex for 30 s, store at -20° for at least 30 min, vortex, centrifuge at 1500 g for 30 min. Remove 400 μ L of the upper organic layer and evaporate it under a stream of nitrogen. Reconstitute with 100 μ L mobile phase, vortex for 30 s, inject a 75 μ L aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: MeCN:buffer 47:53 (Buffer was 6.805 g potassium monophosphate in 1 L water, adjust pH to 6.00 with 10 M NaOH.)

Flow rate: 1

Injection volume: 75

Detector: UV 214

CHROMATOGRAM

Retention time: 7.2

Internal standard: 1-acetamidopyrene (9.7)

Limit of detection: 39 ng/mL

OTHER SUBSTANCES

Simultaneous: lorazepam, antipyrine

Noninterfering: adenosine, albuterol, alphenal, aspirin, caffeine, carbamazepine, cefazolin, cephalixin, cephalothin, cimetidine, ciprofloxacin, claforan, desipramine, enoxacin, fleroxacin, furosemide, hydralazine, hydrochlorothiazide, minoxidil, norfloxacin, phenytoin, propafenone, sulindac, teicoplanin, theophylline, vancomycin

KEY WORDS

plasma

REFERENCE

Awni, W.M.; Bakker, L.J. Antipyrine, indocyanine green, and lorazepam determined in plasma by high-pressure liquid chromatography, *Clin. Chem.*, **1989**, *35*, 2124–2126.

SAMPLE

Matrix: blood

Sample preparation: 25 μ L Plasma + 25 μ L cold (-20°) 60 μ g/mL IS in acetone, vortex for 30 s, centrifuge at 13000 g for 1 min, inject a 10–20 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m octylsilane (Alltech)

Mobile phase: MeCN:50 mM pH 4.0 phosphate buffer 42:58

Flow rate: 1.5

Injection volume: 10–20

Detector: F ex 214 em 370

CHROMATOGRAM

Retention time: 7.3

Internal standard: 1,1',3,3,3',3'-hexamethyl-4,4',5,5'-dibenzo-2,2'-indotricarbocyanine perchlorate (Eastman Kodak) (8.7)

Limit of detection: 500 ng/mL


KEY WORDS

rat; plasma; pharmacokinetics

REFERENCE

Pollack,G.M.; Brouwer,K.L.R.; Demby,K.B.; Jones,J.A. Determination of hepatic blood flow in the rat using sequential infusions of indocyanine green or galactose, *Drug Metab.Dispos.*, **1990**, 18, 197-202.

Indomethacin

Molecular formula: C₁₉H₁₆ClNO₄**Molecular weight:** 357.79**CAS Registry No.:** 53-86-1, 74252-25-8 (sodium salt trihydrate)**Merck Index:** 4998**Lednicer No.:** 1 318; 2 345; 3 165**SAMPLE****Matrix:** aqueous humor

Sample preparation: 100 µL Aqueous humor + 500 µL MeCN + 30 µL 400 ng/mL (+)-naproxen in MeOH, mix mechanically for 90 s, centrifuge at 3000 g for 20 min. Remove the supernatant and dry it under nitrogen at room temperature, dissolve the residue in 50 µL mobile phase by swirl-mixing for 1 min, centrifuge at 3000 g for 20 s, reduce volume to 20-30 µL, inject an aliquot.

HPLC VARIABLES**Column:** 150 × 4.5 5 µm Ultrasphere octyl**Mobile phase:** MeCN:triethylamine:1.65% glacial acetic acid 505:0.65:495, pH 4.35**Column temperature:** 30**Flow rate:** 1**Injection volume:** 20**Detector:** UV 280**CHROMATOGRAM****Retention time:** 6.64**Internal standard:** naproxen (3.89)**OTHER SUBSTANCES****Extracted:** diclofenac, flurbiprofen, meclofenamic acid

Simultaneous: bacitracin, cortisone acetate, diazepam, fluorometholone, hydrocortisone acetate, imipramine, ketoprofen, ketorolac tromethamine, levobunolol, metipranolol, neomycin, prednisolone acetate, proparacaine, propranolol, salicylic acid, sulfacetamide, suprofen

Noninterfering: acebutolol, acetaminophen, acetazolamide, alprenolol, apraclonidine, atenolol, atropine, betamethasone, betaxolol, bupivacaine, caffeine, cyclopentolate, dexamethasone, diphenhydramine, erythromycin, haloperidol, lidocaine, phenylephrine, polymyxin B, procaine, scopolamine, timolol, tropicamide

KEY WORDS

human; rabbit

REFERENCE

Riegel,M.; Ellis,P.P. High-performance liquid chromatography assay for antiinflammatory agents diclofenac and flurbiprofen in ocular fluids, *J.Chromatogr.B*, **1994**, 654, 140-145.

SAMPLE**Matrix:** aqueous humor

Sample preparation: 200 µL Aqueous humor + 200 µL MeCN, vortex for 15 s, centrifuge at 2000 rpm for 10 min, inject a 20 µL aliquot of the supernatant.

HPLC VARIABLES**Guard column:** 4 × 4.5 µm Lichrospher 100 C18**Column:** 125 × 4.5 µm Lichrospher 100 C18**Mobile phase:** MeOH:50 mM pH 3 sodium phosphate buffer 65:35**Flow rate:** 1**Injection volume:** 20**Detector:** E, Shimadzu L-ECD-6A, glassy carbon working electrode 1.20 V, Ag/AgCl reference electrode or UV 266

CHROMATOGRAM**Retention time:** 7.5**Limit of detection:** 50 ng/mL (UV), 20 ng/mL (E)

KEY WORDS

rabbit

REFERENCE

Baudrit, O.; Fabre, H. Evaluation of electrochemical and fluorescence detection in liquid chromatography for the assay of indomethacin in aqueous humor samples, *J. Liq. Chromatogr.*, **1995**, *18*, 3283–3299.

SAMPLE**Matrix:** bile, blood, gastric contents, urine

Sample preparation: Plasma. 300 µL Plasma + 1 mL MeCN, vortex for 30 s, centrifuge at 3500 rpm. Remove the supernatant and evaporate it to 100 µL under a stream of nitrogen at 50°, inject a 20 µL aliquot. Urine, bile, gastric fluid. 300 µL Urine, bile, or gastric fluid + 100 µL 5 M NaOH, let stand at room temperature for 15 min, adjust the pH with 100 µL 28.3% phosphoric acid, add 1 (urine, bile) or 1.5 (gastric fluid) mL MeCN, vortex for 30 s, centrifuge at 3500 rpm. Remove the supernatant and evaporate it to 100 µL under a stream of nitrogen at 50°, inject a 20 (urine), 35 (bile), or 40 (gastric fluid) µL aliquot. (Borate buffer was 12.4 g boric acid and 10 mL 1 M NaOH made up to 1 L with water, pH 7.2.)

HPLC VARIABLES**Guard column:** Microguard reverse-phase (Bio-Rad)**Column:** 100 × 8 Radial-PAK C18 in a radial compression module**Mobile phase:** MeCN:buffer 70:30 (Buffer was 6.8 g/L KH₂PO₄ adjusted to pH 3.0 with 85% phosphoric acid.)**Flow rate:** 2**Injection volume:** 20–40**Detector:** UV 340

CHROMATOGRAM**Retention time:** 13**Internal standard:** indomethacin

OTHER SUBSTANCES**Extracted:** sulindac

KEY WORDS

plasma; indomethacin is IS

REFERENCE

Musson, D.G.; Vincek, W.C.; Constanzer, M.L.; Detty, T.E. Analytical methods for the determination of sulindac and metabolites in plasma, urine, bile, and gastric fluid by liquid chromatography using ultraviolet detection, *J. Pharm. Sci.*, **1984**, *73*, 1270–1273.

SAMPLE**Matrix:** blood

Sample preparation: Mix 500 µL serum with 100 µL MeOH:50 mM pH 3 sodium phosphate buffer 50:50, vortex for 5 s. Add 1 mL MeCN, vortex for 1 min. Centrifuge the mixture at 14000 rpm for 5 min, decant the clear upper layer. Add 500 µL MeCN to the pellet, mix, centrifuge at 14000 rpm for 5 min. Combine the upper layers and evaporate at 40°. Reconstitute the

residue with 300 μ L MeOH:50 mM pH 3.0 sodium phosphate buffer 40:60 containing 0.5% sodium metabisulfate (sic). Vortex for 30 s and inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Techsil C18 (HPLC Technology, Macclesfield)

Mobile phase: MeCN:50 mM phosphate buffer 46:63, adjusted to pH 3.0 with NaOH

Flow rate: 1.2

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 7.17

Internal standard: indomethacin

OTHER SUBSTANCES

Extracted: sulindac

KEY WORDS

indomethacin is IS; serum

REFERENCE

Kanfer,I.; Brown,C.; Koninigs,M.; Swarbrick,J. Absorption of sulindac from a novel (Pro-SorbTM) liquid formulation, *Biopharm.Drug Dispos.*, **1996**, 17, 209–221.

SAMPLE

Matrix: blood

Sample preparation: Precipitate 100 μ L serum with 200 μ L 2 μ g/mL IS in MeCN, centrifuge at 12 000 g for 5 min. Inject a 50 μ L aliquot of the supernatant.

HPLC VARIABLES

Guard column: 4 \times 4 5 μ m LiChroCART LiChrospher 60 RP Select B

Column: 125 \times 4 5 μ m LiChroCART LiChrospher 60 RP Select B

Mobile phase: MeCN:buffer 60:40 (Buffer was 25 mM pH 3.0 triethylammonium phosphate containing 2% MeCN.)

Flow rate: 1

Injection volume: 50

Detector: UV 280

CHROMATOGRAM

Retention time: 2.28

Internal standard: mefenamin (2.94)

Limit of detection: 110 ng/mL

KEY WORDS

serum

REFERENCE

Hannak,D.; Scharbert,F.; Kattermann,R. Stepwise binary gradient high-performance liquid chromatographic system for routine drug monitoring, *J.Chromatogr.A*, **1996**, 728, 307–310.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma + 250 μ L 1 M sulfuric acid + 5 mL 24 ng/mL p-phenyl-phenol in dichloromethane, vortex for 10 s, centrifuge at 500 g for 5 min. Remove 4 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 100 μ L MeOH, inject a 20–30 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 Zorbax ODS

Mobile phase: Gradient. MeOH:100 mM pH 5 acetate buffer 45:55 for 3 min, to 62:38 over 2 min, maintain at 62:38 for 10 min.

Column temperature: 40

Flow rate: 1.5

Injection volume: 20-30

Detector: UV 254

CHROMATOGRAM

Retention time: 13.83

Internal standard: p-phenylphenol (12.43)

Limit of quantitation: 100 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

Simultaneous: acetaminophen, caffeine, carbamazepine, ethosuximide, fenoprofen, naproxen, phenobarbital, phenytoin, primidone, quinidine, salicylic acid, sulindac, theophylline, tolmetin, valproic acid

Noninterfering: ibuprofen

KEY WORDS

plasma

REFERENCE

Shimek, J.L.; Rao, N.G.S.; Wahba Khalil, S.K. High performance liquid chromatographic analysis of tolmetin, indomethacin and sulindac in plasma, *J.Liq.Chromatogr.*, **1981**, *4*, 1987-2013.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Plasma + 100 μ L pH 2 dilute sulfuric acid + 1 mL MeCN, vortex for 30 s, centrifuge at 2500 rpm for 5 min. Remove the supernatant and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 200 μ L mobile phase, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Zorbax CN

Mobile phase: MeCN:4% aqueous acetic acid 45:55

Flow rate: 1

Injection volume: 10

Detector: UV 340

CHROMATOGRAM

Retention time: 6

Internal standard: indomethacin

OTHER SUBSTANCES

Extracted: sulindac

KEY WORDS

plasma; indomethacin is IS

REFERENCE

Clark, C.R.; McMillian, C.L.; Hoke, J.F.; Campagna, K.D.; Ravis, W.R. Liquid chromatographic determination of sulindac and metabolites in serum, *J.Chromatogr.Sci.*, **1987**, *25*, 247-251.

SAMPLE

Matrix: blood

Sample preparation: Activate a 1 mL Bond-Elut C8 SPE cartridge with 2 mL MeOH then 1 mL 10 mM HCl, do not allow it to dry completely. Sonicate 1 mL whole blood for 20-30 min then apply to cartridge. Wash with 100 μ L water, elute with three 500 μ L portions of MeOH: MeCN:1% aqueous ammonium hydroxide 50:20:30, combine eluents and evaporate to dryness under a stream of nitrogen at 40°. Redissolve in 1 mL MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 250 × 4.5 5 µm Spherisorb ODS

Mobile phase: MeCN:MeOH:buffer 35:13:52 (Buffer was water adjusted to pH 3.2 with orthophosphoric acid)

Flow rate: 1

Injection volume: 20

Detector: UV 250

CHROMATOGRAM

Retention time: 9

OTHER SUBSTANCES

Simultaneous: ketoprofen, acetaminophen, salicylic acid, naproxen, fenoprofen, ibuprofen

KEY WORDS

whole blood; SPE

REFERENCE

Moore,C.M.; Tebbett,I.R. Rapid extraction of anti-inflammatory drugs in whole blood for HPLC analysis, *Forensic Sci.Int.*, **1987**, *34*, 155–158.

SAMPLE

Matrix: blood

Sample preparation: Wash a Sep-Pak C18 cartridge with 2 mL MeOH, 5 mL water, and 1 mL 0.25 mM pH 3.0 ammonium phosphate buffer. 20–200 µL Plasma + 100 µL MeOH + 20 µL 50 µg/mL indomethacin in MeOH + 100 µL 0.25 mM pH 3.0 ammonium phosphate buffer + 100 µL water, vortex for 2 min, centrifuge at 1800 g for 10 min. Add the supernatant to the cartridge, wash with 5 mL water, elute twice with 5 mL portions of MeOH. Evaporate eluate to dryness under vacuum, dissolve the residue in 1 mL mobile phase, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 100 × 4.6 5 µm Brownlee RP18

Mobile phase: MeOH:buffer 75:25 (Buffer prepared by diluting 0.25 mM ammonium phosphate buffer adjusted to pH 3.0 with orthophosphoric acid.)

Injection volume: 20

Detector: E ESA Coulochem Model 5100 A, + 0.9 V

CHROMATOGRAM

Retention time: 14.6

Internal standard: naproxen (10.0)

Limit of detection: 20 ng/mL

OTHER SUBSTANCES

Also analyzed: sulindac, piroxicam, diflunisal

KEY WORDS

plasma

REFERENCE

Kazemifard,A.G.; Moore,D.E. Liquid chromatography with amperometric detection for the determination of non-steroidal anti-inflammatory drugs in plasma, *J.Chromatogr.*, **1990**, *533*, 125–132.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 1 mL 1 M pH 10.0 Tris buffer containing 150 ng/mL indomethacin + 8 mL diethyl ether, vortex 3 min, freeze for 1 h. Remove organic phase and evaporate it to dryness at room temperature. Dissolve residue in 200 µL MeOH, inject 20 µL aliquot.

HPLC VARIABLES

Column: 70 × 4.6 3 µm Ultrasphere XL ODS

Mobile phase: MeOH:20 mM ammonium acetate buffer (pH 5.0) 65:35

Flow rate: 1.5

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 1.64

Internal standard: indomethacin

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Simultaneous: dipyrindamole

KEY WORDS

plasma; indomethacin is IS

REFERENCE

Barberi,M.; Merlin,J.L.; Weber,B. Sensitive determination of free and plasma protein-bound dipyrindamole by high-performance liquid chromatography, *J.Chromatogr.*, **1991**, 565, 511–515.

SAMPLE

Matrix: blood

Sample preparation: Condition a Bond Elut phenyl SPE cartridge with 5 mL MeOH and 5 mL water. Adjust pH of 500 μ L plasma to 3.4 with 345 mM citrate buffer, add to SPE cartridge, wash with water, dry, elute with 5 mL hexane:diethyl ether 50:50. Evaporate the eluate to dryness under a stream of nitrogen at 40°, reconstitute the residue in 200 μ L MeOH, inject a 25 μ L aliquot.

HPLC VARIABLES

Guard column: 6 \times 4 μ m Nova-Pack C18

Column: 150 \times 4.6 μ m Ultrasphere ODS

Mobile phase: MeCN:20 mM ammonium sulfate 55:45

Flow rate: 1.5

Injection volume: 25

Detector: UV 340

CHROMATOGRAM

Retention time: 7.5

Internal standard: indomethacin

OTHER SUBSTANCES

Extracted: phenylbutazone, oxyphenbutazone, suxibuzone

KEY WORDS

plasma; SPE; indomethacin is IS

REFERENCE

Caturia,M.C.; Cusido,E. Solid-phase extraction for the high-performance liquid chromatographic determination of indomethacin, suxibuzone, phenylbutazone and oxyphenbutazone in plasma, avoiding degradation of compounds, *J.Chromatogr.*, **1992**, 581, 101–107.

SAMPLE

Matrix: blood

Sample preparation: 20 μ L serum + 20 μ L MeCN, vortex for a few s, centrifuge at 10000 g for 2 min, inject a 20 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 150 \times 4.6 μ m Inertsil ODS-2

Mobile phase: MeCN:25 mM phosphate buffer 35:65 containing 30 mM hydrogen peroxide, adjusted to pH 7.0 with 1 M NaOH

Flow rate: 1

Injection volume: 20

Detector: F ex 358 em 462 following post-column reaction. The column effluent flowed through a 15 m × 0.5 mm ID stainless steel coil at 180° then a 3 m × 0.5 mm ID stainless steel coil at 15° to the detector.

CHROMATOGRAM

Retention time: 10

Limit of detection: 500 ng/mL

KEY WORDS

post-column reaction; serum; pharmacokinetics

REFERENCE

Kubo,H.; Umiguchi,Y.; Kinoshita,T. Fluorometric determination of indomethacin in serum by high performance liquid chromatography using in-line oxidation with hydrogen peroxide, *J.Liq.Chromatogr.*, **1993**, *16*, 465–474.

SAMPLE

Matrix: blood

Sample preparation: 50 µL Plasma + 250 µL 0.75 µg/mL mefenamic acid in MeCN + 50 µL MeCN, vortex, centrifuge at 9000 g for 3 min. Remove 250 µL of the supernatant and evaporate it to dryness under vacuum, dissolve the residue in 50 µL mobile phase, inject a 20 µL aliquot.

HPLC VARIABLES

Guard column: 10 × 4.6 Alltech 5 µm C18 bonded-phase silica

Column: 250 × 4.6 Vydac column packed with Merck 5 µm C18 bonded-phase silica

Mobile phase: MeCN:10 mM phosphoric acid 60:40, pH 2.6

Flow rate: 0.9

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 6.3

Internal standard: mefenamic acid (9.2)

Limit of detection: 60 ng/mL

KEY WORDS

plasma

REFERENCE

Niopas,I.; Mamzoridi,K. Determination of indomethacin and mefenamic acid in plasma by high-performance liquid chromatography, *J.Chromatogr.B*, **1994**, *656*, 447–450.

SAMPLE

Matrix: blood

Sample preparation: Condition a 1 mL 60 µm Separon SGX C18 SPE cartridge with 5 mL MeOH, 5 mL water, and 5 mL buffer. 250 µL Blood + 500 µL water, shake for 5 min, sonicate for 5 min, let stand at room temperature for 5 min, add 1 mL buffer, shake for 5 min, centrifuge at 1930 g for 10 min, add supernatant to cartridge, wash with 5 mL buffer, wash with 10 mL water, dry with vacuum for 5 min, elute with dichloromethane. Evaporate eluate to dryness under a stream of nitrogen, dissolve in 100 µL mobile phase, inject 10 µL aliquot. (Buffer was 66 mM KH₂PO₄ adjusted to pH 2.0 with phosphoric acid.)

HPLC VARIABLES

Column: 150 × 3.3 5 µm Separon SGX C18 glass column

Mobile phase: MeOH water 220:100, adjusted to pH 3.0 with 5% perchloric acid

Flow rate: 1.3

Injection volume: 10

Detector: UV 222

CHROMATOGRAM**Retention time:** 7.8**Internal standard:** indomethacin

OTHER SUBSTANCES**Simultaneous:** ibuprofen

KEY WORDS**SPE;** indomethacin is IS; rabbit; human

REFERENCE

Sochor,J.; Klimes,J.; Zahradnicek,M.; Sedlacek,J. High-performance liquid chromatographic assay for ibuprofen in whole blood using solid-phase extraction, *J.Chromatogr.B*, **1994**, 654, 282–286.

SAMPLE**Matrix:** blood

Sample preparation: Condition a Bond-Elut C2 SPE cartridge with 1 mL MeOH and 1 mL mobile phase. 1 mL Serum + 1 drop saturated ammonium sulfate solution + 1 drop 1 M HCl, vortex for 3 min, add to the SPE cartridge, wash with six 500 μ L portions of wash solvent, elute with four 250 μ L aliquots of mobile phase, combine the eluates, vortex, inject a 100 μ L aliquot. (Wash solvent was MeCN:water adjusted to pH 3.0 with phosphoric acid 20:80.)

HPLC VARIABLES**Column:** 150 \times 4.6 5 μ m Ultrasphere C8**Mobile phase:** MeCN:68 mM pH 2.5 phosphate buffer 55:45**Flow rate:** 0.5**Injection volume:** 100

Detector: F ex 232 em 335 (filter) following post-column photolysis. The effluent from the column flowed through a 7.9 m \times 0.3 mm i.d. coil of PTFE irradiated by an SC3-9 UV lamp (UVP) (cooled with air) to the detector.

CHROMATOGRAM**Retention time:** 12**Internal standard:** indomethacin

OTHER SUBSTANCES**Extracted:** sulindac

KEY WORDS**serum;** post-column reaction; **SPE;** indomethacin is IS; post-column photochemical derivatization

REFERENCE

Siluveru,M.; Stewart,J.T. Determination of sulindac and its metabolites in human serum by reversed-phase high-performance liquid chromatography using on-line post-column ultraviolet irradiation and fluorescence detection, *J.Chromatogr.B*, **1995**, 673, 91–96.

SAMPLE**Matrix:** blood

Sample preparation: Erythrocytes. 500 μ L Erythrocytes + 900 μ L water, shake for 5 min, sonicate for 5 min, let stand at room temperature for 5 min, add 400 μ L 3 M HCl, shake for 5 min, add 6 mL dichloromethane, shake, centrifuge at 1930 g for 5 min. Remove 5 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 100 μ L mobile phase, inject a 10 μ L aliquot. Plasma. Acidify 500 μ L plasma gradually with 900 μ L 1 M HCl, shake, add 200 μ L 3 M HCl, shake for 5 min, add 6 mL dichloromethane, shake, centrifuge at 1930 g for 5 min. Remove 5 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 100 μ L mobile phase, inject a 10 μ L aliquot.

HPLC VARIABLES**Column:** 150 \times 3.3 5 μ m C18 glass column (Tessek)**Mobile phase:** MeOH:water 66:30 adjusted to pH 3.0 with 5% perchloric acid

Flow rate: 1.3
Injection volume: 10
Detector: UV 222

CHROMATOGRAM

Retention time: 7.8
Internal standard: indomethacin

OTHER SUBSTANCES

Extracted: ibuprofen
Simultaneous: diazepam, phenylanthranilic acid

KEY WORDS

plasma; erythrocytes; rabbit; indomethacin is IS

REFERENCE

Sochor,J.; Klimes,J.; Sedláček,J.; Zahradnicek,M. Determination of ibuprofen in erythrocytes and plasma by high-performance liquid chromatography, *J.Pharm.Biomed.Anal.*, **1995**, *13*, 899–903.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 253

CHROMATOGRAM

Retention time: 8.91

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vindesine; mexiletine; dipyrindamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phencyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; bupren-

orphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loperamide; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; acepromazine; glibenclamide; chlorphenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrodine; phenylbutazone; demoxiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozone; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254–262.

SAMPLE

Matrix: blood, CSF, gastric contents, urine

Sample preparation: 200 μ L Serum, urine, CSF, or gastric fluid + 300 μ L reagent. Flush column A to waste with 500 μ L 500 mM ammonium sulfate, inject sample onto column A, flush column A to waste with 500 μ L 500 mM ammonium sulfate, backflush the contents of column A onto column B with mobile phase, monitor the effluent from column B. (Reagent was 8.05 M guanidine HCl and 1.02 M ammonium sulfate in water.)

HPLC VARIABLES

Column: A 40 μ m preparative grade C18 (Analytichem); B 75 \times 2.1 pellicular C18 (Whatman) + 250 \times 4.6 5 μ m C8 end-capped (Whatman)

Mobile phase: Gradient. A was 50 mM pH 4.5 KH_2PO_4 . B was MeCN:isopropanol 80:20. A:B 90:10 for 1 min, to 30:70 over 20 min.

Column temperature: 50

Flow rate: 1.5

Detector: UV 220

CHROMATOGRAM

Retention time: 15.27

Internal standard: heptanophenone (19)

OTHER SUBSTANCES

Extracted: acetaminophen, allobarbitol, azinphos, barbitol, brallobarbitone, bromazepam, butethal, caffeine, carbamazepine, carbaryl, cephaloridine, chloramphenicol, chlordiazepoxide, chlorothiazide, chlorvinphos, clothiapine, cocaine, coomassie blue, desipramine, diazepam, diphenhydramine, dipipanone, ethylbromphos, flufenamic acid, formothion, griseofulvin, lidocaine, lorazepam, malathion, medazepam, midazolam, oxazepam, paraoxon, penicillin G, pentobarbital, prazepam, propoxyphene, prothiophos, quinine, salicylic acid, secobarbital, strychnine, sulfamethoxazole, theophylline, thiopental, thioridazine, trimethoprim

KEY WORDS

serum; column-switching

REFERENCE

Kruger, P.B.; Albrecht, C.F. De V.; Jaarsveld, P.P. Use of guanidine hydrochloride and ammonium sulfate in comprehensive in-line sorption enrichment of xenobiotics in biological fluids by high-performance liquid chromatography, *J. Chromatogr.*, **1993**, *612*, 191–198.

SAMPLE

Matrix: blood, tissue

Sample preparation: Plasma. 50 μ L Plasma + 150 μ L 8 μ g/mL mefenamic acid in MeOH, mix, centrifuge at 15000 rpm, filter (0.45 μ m), inject an aliquot. Tissue. Homogenize liver in ice-cold 10 mM pH 7.4 phosphate buffer, 1 mL homogenate + 2 mL 8 μ g/mL mefenamic acid in MeOH, mix, centrifuge at 15000 rpm, inject an aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Inertsil ODS-2

Mobile phase: MeCN:MeOH:water:acetic acid 65:10:25:1

Flow rate: 1

Detector: UV 254

CHROMATOGRAM

Internal standard: mefenamic acid

Limit of detection: 30 ng/mL

KEY WORDS

plasma; rat; liver

REFERENCE

Ogiso,T.; Iwaki,M.; Kinoshita,T.; Tanino,T.; Paku,T. Pharmacokinetics of indomethacin octyl ester (prodrug) and indomethacin produced from the prodrug, *J.Pharm.Sci.*, **1994**, *83*, 34–37.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 1 mL Plasma + 50 μ L 500 ng/mL mefenamic acid + 1 mL 100 mM HCl + 10 mL dichloromethane, rotate for 10 min, centrifuge at 1500 g for 15 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 45°. Redissolve the residue in mobile phase, inject a 20 μ L aliquot. Urine. 50 μ L Urine + 1 mL mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 75 \times 4.6 3 μ m Supelcosil LC-8

Mobile phase: MeCN:50 mM phosphoric acid 45:55

Flow rate: 1

Injection volume: 20

Detector: UV 235

CHROMATOGRAM

Retention time: 5

Internal standard: mefenamic acid (8)

Limit of detection: 50-250 ng/mL

OTHER SUBSTANCES

Simultaneous: naproxen, flunixin, thiosalicylic acid, ethacrynic acid, phenylbutazone

KEY WORDS

plasma

REFERENCE

Singh,A.K.; Jang,Y.; Mishra,U.; Granley,K. Simultaneous analysis of flunixin, naproxen, ethacrynic acid, indomethacin, phenylbutazone, mefenamic acid and thiosalicylic acid in plasma and urine by high-performance liquid chromatography and gas chromatography-mass spectrometry, *J.Chromatogr.*, **1991**, *568*, 351–361.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 100 μ L Plasma + 400 μ L 200 mM perchloric acid, centrifuge at 3000 g, inject a 20 μ L aliquot of the supernatant. Urine. Dilute 1:1 with water, inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: 75 × 2.1 pellicular reverse phase (Chrompack no. 28653)

Column: 250 × 4.6 5 µm Cp Spherisorb ODS (Chrompack)

Mobile phase: Gradient; MeCN:5 g/L orthophosphoric acid from 80:20 to 60:40 (sic, probably 20:80 to 40:60) over 30 min then stay there for 5 min, then to initial conditions over 5 min, equilibrate for 2 min before next injection

Flow rate: 1.2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 28.13

Limit of detection: 5.3 ng/mL

Limit of quantitation: 146 ng/mL (urine), 60 ng/mL (plasma)

OTHER SUBSTANCES

Simultaneous: metabolites, glucuronides

KEY WORDS

plasma; pharmacokinetics; also details of preparative procedure

REFERENCE

Vree,T.B.; Van den Biggelaar-Martea,M.; Verwey-van Wissen,C.P.W.G.M. Determination of indomethacin, its metabolites and their glucuronides in human plasma and urine by means of direct gradient high-performance liquid chromatographic analysis. Preliminary pharmacokinetics and effect of probenecid, *J.Chromatogr.*, **1993**, 616, 271–282.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10–30

Detector: UV 201.7

CHROMATOGRAM

Retention time: 21.748

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

SAMPLE

Matrix: bulk

Sample preparation: 500 μ L Plasma + 100 μ L 5 mg/mL indomethacin in MeOH, filter (Sartorius SM 13243 ultrafiltration unit at 4000 g for 30 min). Add filtrate to a dry Chem Elut column (modified diatomaceous earth), leave 3 to 5 min, elute with 6 mL diethyl ether, evaporate eluant under a stream of nitrogen at room temperature, sonicate residue with 100 μ L MeOH for 10 min, vortex, inject 20 μ L aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: MeOH:100 mM pH 4.0 sodium acetate buffer 55:45

Flow rate: 1

Injection volume: 20

Detector: E, ESA Coulochem 5100A dual electrode with an ESA guard cell, + 0.65 V

CHROMATOGRAM

Retention time: 11.74

Internal standard: indomethacin

OTHER SUBSTANCES

Simultaneous: dipyridamole

KEY WORDS

plasma; indomethacin is IS

REFERENCE

Barberi-Heyob,M.; Merlin,J.L.; Pons,L.; Calco,M.; Weber,B. A sensitive isocratic liquid chromatography assay for the determination of dipyridamole in plasma with electrochemical detection, *J.Liq.Chromatogr.*, **1994**, *17*, 1837–1848.

SAMPLE

Matrix: formulations

Sample preparation: Reconstitute 1 mg indomethacin sodium trihydrate injection with 2 mL water, inject a 5 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.2 5 μ m Ultrasphere ODS C18 (Beckman)

Mobile phase: MeCN:50 mM H_3PO_4 60:40

Flow rate: 1

Injection volume: 5

Detector: UV 260

CHROMATOGRAM

Retention time: 6.5-7.0

KEY WORDS

injections

REFERENCE

Walker,S.E.; Gray,S.; Schmidt,B. Stability of reconstituted indomethacin sodium trihydrate in original vials and polypropylene syringes, *Am.J.Health-Syst.Pharm.*, **1998**, *55*, 154–158.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 250 μ L Microsomal incubation, 150 μ L ice-cold MeCN and 2.5 ng keto-profen, centrifuge. Extract the mixture with 4 mL ethyl acetate, centrifuge at 3000 rpm for 10 min, remove the organic fraction and evaporate it under a gentle stream of nitrogen at 40°. Dissolve the residue in 30 μ L MeOH and dilute to 60 μ L with water. Inject a 30 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m RP-18 (Kanto Chemical, Tokyo)

Mobile phase: MeCN:water 40:60 containing 0.6% acetic acid

Column temperature: 35

Flow rate: 1

Injection volume: 30

Detector: UV 260

CHROMATOGRAM

Retention time: 49.0

Internal standard: ketoprofen (18.0)

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

liver; pharmacokinetics

REFERENCE

Nakajima,M.; Inoue,T.; Shimada,N.; Tokudome,S.; Yamamoto,T.; Kuroiwa,Y. Cytochrome P450 2C9 catalyzes indomethacin O-demethylation in human liver microsomes, *Drug Metab.Dispos.*, **1998**, *26*, 261–266.

SAMPLE

Matrix: perfusate

Sample preparation: Mix 100 nL perfusate with 600 nL perfusion fluid containing IS, inject an aliquot. (Perfusion fluid contained 104 mM NaCl, 25 mM sodium bicarbonate, 2.3 mM sodium biphosphate, 10 mM sodium acetate, 1.2 mM calcium chloride, 1 mM magnesium sulfate, 5 mM KCl, 5 mM dextrose, and 5 mM alanine.)

HPLC VARIABLES

Column: 300 × 2 10 µm µBondapak C18

Mobile phase: MeOH:water 52:48

Flow rate: 0.13

Injection volume: 0.2

Detector: F ex 295 em 376 following post-column reaction. The column effluent mixed with 4 M NaOH pumped at 0.0013 mL/min and the mixture flowed through a 130 µL PTFE coil at 64° to the detector.

CHROMATOGRAM

Internal standard: phenylbutazone

Limit of detection: 25 ng/mL

KEY WORDS

post-column reaction; microbore

REFERENCE

De Zeeuw,D.; Leinfelder,J.L.; Brater,D.C. Highly sensitive measurement of indomethacin using a high performance liquid chromatographic technique combined with post column in-line hydrolysis, *J.Chromatogr.*, **1986**, *380*, 157–162.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4 ODS (Hitachi)

Mobile phase: MeCN:50 mM phosphoric acid 40:60 adjusted to pH 5.5 with NaOH

Column temperature: 55

Flow rate: 0.6

Injection volume: 20

Detector: UV 264

OTHER SUBSTANCES

Also analyzed: carbamazepine, fenbufen, ketoprofen, α.-naphthoquinone, naproxen, tolmetin

REFERENCE

Sugawara,M.; Takekuma,Y; Yamada,H.; Kobayashi,M.; Iseki,K.; Miyazaki,K. A general approach for the prediction of the intestinal absorption of drugs: regression analysis using the physicochemical properties and drug-membrane electrostatic interactions, *J.Pharm.Sci.*, **1998**, 87, 960-966.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4 OmniPac PAX-500 (Dionex)

Mobile phase: Gradient. A was MeCN:10 mM sodium carbonate 18:82. B was MeCN:50 mM sodium carbonate 33:67. A:B from 100:0 to 0:100 over 10 min.

Flow rate: 1

Detector: UV 254

CHROMATOGRAM

Retention time: 18

OTHER SUBSTANCES

Simultaneous: tolmetin, aspirin, ibuprofen, fenbufen, naproxen, carprofen, diflunisal

REFERENCE

Slingsby,R.W.; Rey,M. Determination of pharmaceuticals by multi-phase chromatography: Combined reversed phase and ion exchange in one column, *J.Liq.Chromatogr.*, **1990**, 13, 107-134.

SAMPLE

Matrix: solutions

Sample preparation: Sample + 400 μ L 5 mM DBD-PZ + 70 mM diethylphosphorocyanidate in MeCN, react for 6 h, inject a 1 μ L aliquot. (Synthesis of 4-(N,N-dimethylaminosulfonyl)-7-N-piperazino-2,1,3-benzoxadiazole (DBD-PZ) is as follows. Dissolve 0.5 g magnesium sulfate heptahydrate and 6 g NaOH in 60 mL water, throughout the reaction keep the flask at about 20° with cold water cooling, add 15 mL 30% hydrogen peroxide, add 75 mL MeOH, add 12.1 g powdered benzoyl peroxide in one go, stir for 10 min, pour into 150 mL 20% sulfuric acid, extract three times with 50 mL portions of chloroform, determine peroxybenzoic acid concentration by iodometric titration (Tetrahedron 1967, 23, 3327). Slowly add 110 mL 1 M peroxybenzoic acid in chloroform to 7 g 2,6-difluoroaniline dissolved in 100 mL chloroform, stir at room temperature, when reaction is complete (iodometric titration) wash with 2% sodium thiosulfate, wash with 5% sodium carbonate, wash with water, dry over anhydrous sodium sulfate, evaporate to dryness under reduced pressure, recrystallize 2,6-difluoronitrosobenzene form EtOH (mp 108.5-109.5). Stir 8.5 g 2,6-difluoronitrosobenzene in 85 mL DMSO at room temperature and add a solution of 3.91 g sodium azide in 85 mL DMSO dropwise, let stand for about 1 h, add to a large volume of water, extract with ether, dry the extracts over anhydrous sodium sulfate, evaporate to dryness under reduced pressure and distil to give 4-fluoro-2,1,3-benzoxadiazole as a colorless oil (bp 83°/12 mm Hg) (J.Chem.Soc.(C) 1970, 1433). Add 11 mL chlorosulfonic acid dropwise to 3 g 4-fluoro-2,1,3-benzoxadiazole in 10 mL chloroform at 0-10° (use a calcium chloride drying tube), stir at room temperature for 1 h, reflux for 2 h, cool, slowly pour into ice water, remove the organic layer, extract the aqueous layer with chloroform, combine the organic layer, wash, dry over anhydrous magnesium sulfate, evaporate under reduced pressure, take up the residue in 5 mL benzene (Caution! Benzene is a carcinogen!), chromatograph on a 150 × 30 column of silica gel (100-200 mesh Kanto Chemical) with n-hexane:benzene 50:50, evaporate the appropriate fractions to give 4-(chlorosulfonyl)-7-fluoro-2,1,3-benzoxadiazole (CBD-F) as pale yellow needles (mp 64-66°) (Anal. Chem. 1984, 56, 2461). Stir 0.76 g CBD-F in 70 mL MeCN at 0-10° and add 1 g dimethylamine hydrochloride in 10 mL 100 mM pH 10 borax dropwise, adjust pH to 5 with 1 M HCl, concentrate to about 10 mL under reduced pressure, extract three times with 200 mL portions of diethyl ether, wash with water, dry over anhydrous magnesium sulfate, evaporate under reduced pressure, chromatograph on a 500 × 20 column of silica gel with chloroform, isolate the appropriate fraction and re-chromatograph on the same column with ethyl acetate:benzene 1:2 to give 4-(N,N-dimethylaminosulfonyl)-7-fluoro-2,1,3-benzoxadiazole (DBD-F) as white needles (mp 124-125°) (yield = 1% !). On a Merck no. 5714 60F₂₅₄ tlc plate eluted with chloroform DBD-F has R_f 0.32 and lies between two other reaction products (Analyst 1989, 114, 413). It is also reported that DBD-F can be purchased from Tokyo Kasei (TCI America, Portland OR). Add 123 mg 4-(N,N-dimethyl-

aminosulfonyl)-7-fluoro-2,1,3-benzoxadiazole in 20 mL MeCN dropwise to 129 mg piperazine in 20 mL MeCN at room temperature, stir for 30 min, evaporate under reduced pressure, dissolve residue in 50 mL 5% HCl, wash three times with 20 mL ethyl acetate, discard ethyl acetate extracts, adjust pH of aqueous solution to 13-14 with 5% NaOH, extract five times with 50 mL ethyl acetate, combine extracts, wash with 20 mL water, dry over anhydrous sodium sulfate, evaporate under vacuum to give 4-(N,N-dimethylaminosulfonyl)-7-N-piperazino-2,1,3-benzoxadiazole as orange crystals (mp 121-2°).

HPLC VARIABLES

Column: 150 × 4.6 5 µm Inertsil ODS-2

Mobile phase: MeCN:water 65:35

Column temperature: 40

Flow rate: 1

Injection volume: 1

Detector: F ex 437 em 561

CHROMATOGRAM

Retention time: 8

Limit of detection: 3.9 fmol

OTHER SUBSTANCES

Simultaneous: ibuprofen

REFERENCE

Toyooka,T.; Ishibashi,M.; Takeda,Y.; Nakashima,K.; Akiyama,S.; Uzu,S.; Imai,K. Precolumn fluorescence tagging reagent for carboxylic acids in high-performance liquid chromatography: 4-substituted-7-aminoalkyl-amino-2,1,3-benzoxadiazoles, *J.Chromatogr.*, **1991**, 588, 61-71.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, am-triptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenopropfen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminosilbene, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenine acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclufenamic acid, medazepam, mefenamic acid, megestrol,

mepacrine, meperidine, mephentermine, mephentytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methypylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, nor-epinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phenacyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithydione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfathidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlylcypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233–242.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 125 × 3 Ecocart LiChrospher 100 RP-18

Mobile phase: Isopropanol:100 mM KH₂PO₄:formic acid 54:100:0.1

Flow rate: 0.6

Detector: UV 254

CHROMATOGRAM

Retention time: 10.7

Limit of quantitation: 200-500 ng/mL

KEY WORDS

solutions acetaminacin; diclofenac; flurbiprofen; lonazolac; ketoprofen; naproxen; piroxicam; sulindac; tenoxicam

REFERENCE

Baeyens,W.R.G.; Van Der Weken,G.; Van Overbeke,A.; Zhang,Z.D. Preliminary results on the LC-separation of non-steroidal anti-inflammatory agents in conventional and narrow-bore RP set-ups applying columns with different internal diameters, *Biomed.Chromatogr.*, **1995**, *9*, 261–262.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 150 × 3.9 Nova-Pak C18

Mobile phase: MeCN:water 45:55, pH adjusted to 3.5 with acetic acid

Detector: UV 280

CHROMATOGRAM

Internal standard: clomethacin

OTHER SUBSTANCES

Also analyzed: diclofenac, phenylbutazone

REFERENCE

Guterres, S.S.; Fessi, H.; Barratt, G.; Puisieux, F.; Devissaguet, J.-P. Poly(D,L-lactide) nanocapsules containing non-steroidal anti-inflammatory drugs: Gastrointestinal tolerance following intravenous and oral administration, *Pharm.Res.*, **1995**, *12*, 1545–1547.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4.5 µm LiChrospher 100 RP-8 (B)

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 8.45 (A), 9.22 (B)

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordiazepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenopropfen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazindol, mefenamic acid, meperidine, mephentermine, mepivacaine, mesoridazine, metaproterenol, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemioline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phenolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimizide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfapyrazole, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocinamide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, trifluorpromazine, trimetoprim, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

KEY WORDS

details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J.Chromatogr.A*, **1995**, *692*, 103–119.

SAMPLE

Matrix: solutions

Sample preparation: Inject a 20 µL aliquot of a 250 ng/mL solution in pH 2.5 buffer.

HPLC VARIABLES

Column: 150 × 0.32 3 µm Hypersil C18

Mobile phase: MeCN:pH 6.0 acetate/citrate buffer 45:55

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 14

Limit of detection: 0.4 ng/mL

Limit of quantitation: 1.2 ng/mL

KEY WORDS

microcolumn

REFERENCE

Streel,B.; Ceccato,A.; Chiap,P.; Hubert,P.; Crommen,J. Injection-generated solvent and pH gradients for sample enrichment on injection of large volumes in microcolumn liquid chromatography, *Biomed.Chromatogr.*, **1995**, 9, 254–256.

SAMPLE

Matrix: solutions

Sample preparation: Inject a 10 µL aliquot.

HPLC VARIABLES

Guard column: 10 × 4.6 10 µm Spherisorb silica

Column: 120 × 2 5 µm Spherisorb ODS

Mobile phase: MeOH:20 mM pH 7.0 phosphate buffer 58:42

Column temperature: 30

Flow rate: 0.378

Injection volume: 10

Detector: UV 254

REFERENCE

Lough,W.J.; Mills,M.J.; Maltas,J. Analyte adsorption in liquid chromatography valve injectors for samples in non-eluting solvents, *J.Chromatogr.A*, **1996**, 726, 67–75.

SAMPLE

Matrix: solutions

Sample preparation: Filter (0.45 µm), dilute the filtrate with mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 150 × 4.6 5 µm Hypersil ODS

Mobile phase: MeCN:10 mM pH 4 acetate buffer 50:50

Detector: UV 242

REFERENCE

Okimoto,K.; Rajewski,R.A.; Uekama,K.; Jona,J.A.; Stella,V.J. The interaction of charged and uncharged drugs with neutral (HP-β-CD) and anionically charged (SBE7-β-CD) β-cyclodextrins, *Pharm.Res.*, **1996**, 13, 256–264.

SAMPLE

Matrix: solutions

Sample preparation: Inject a 20 µL aliquot of a 100–500 µg/mL solution in mobile phase.

HPLC VARIABLES

Column: 100 × 4.6 5 µm Hypersil C8 MOS 100A coated with phosphatidylcholine (95% pure soybean lecithin, Epikuron, Lucas Meyer & Co.) (Coat column by recycling a 1 mM solution of phosphatidylcholine in MeOH:water 80:20 for 24 h.)

Mobile phase: MeCN:35 mM pH 7.4 sodium phosphate buffer 40:60

Flow rate: 0.5–2

Injection volume: 20**Detector:** UV 254**CHROMATOGRAM****Retention time:** k' 0.58**OTHER SUBSTANCES**

Also analyzed: amoxicillin, antipyrine, carbamazepine, chlorpheniramine, chlorpromazine, clonidine, codeine, desipramine, diphenhydramine, dipyrindamole, ephedrine, flufenamic acid, haloperidol, hydroxyzine, imipramine, lidocaine, megestrol acetate, metoprolol, nabumetone, nadolol, phenobarbital, phenol, promazine, propranolol, pyrilamine, quinidine, ropinirole, testosterone, thioridazine, tolfenamic acid, verapamil

Noninterfering: acetaminophen, aspirin, azathioprine, caffeine, carprofen, chlorambucil, cimetidine, fenoterol, flurbiprofen, ibuprofen, ketoprofen, ranitidine, salicylic acid, sulfamethoxazole, theophylline, thioguanine, tiaprofenic acid, trimethoprim, valproic acid


KEY WORDS

comparison with capillary electrophoresis

REFERENCE

Hanna, M.; de Biasi, V.; Bond, B.; Salter, C.; Hutt, A. J.; Camilleri, P. Estimation of the partitioning characteristics of drugs: A comparison of a large and diverse drug series utilizing chromatographic and electrophoretic methodology, *Anal. Chem.*, **1998**, 70, 2092–2099.

Indoramin

**Molecular formula:** C₂₂H₂₅N₃O**Molecular weight:** 347.46**CAS Registry No.:** 26844-12-2**Merck Index:** 5000**Lednicer No.:** 2 344**SAMPLE****Matrix:** blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES**Guard column:** 20 mm long Symmetry C18**Column:** 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10–30**Detector:** UV 200.5**CHROMATOGRAM****Retention time:** 12.533


KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

Inositol



Molecular formula: $C_6H_{12}O_6$

Molecular weight: 180.16

CAS Registry No.: 87-89-8, 573-35-3 (monophosphate)

Merck Index: 5008

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 300 × 8.7 Aminex HPX-87C Ca^{++} (Bio-rad)

Mobile phase: Water

Column temperature: 50

Flow rate: 0.6

Detector: E, pulsed amperometric detector (Dionex ?), E1 0.1 V, T1 300 ms, E2 0.6 V, T2 120 ms, E3 -0.8 V, T3 300 ms, following post-column reaction. The column effluent mixed with 100 mM NaOH pumped at 0.2 mL/min and the mixture flowed to the detector.

CHROMATOGRAM

Retention time: 13.64 (myo-inositol), 14.16 (chiro-inositol), 11.16 (scyllo-inositol), 19.58 (neo-inositol)

OTHER SUBSTANCES

Also analyzed: 2-deoxygalactitol, 2-deoxyribose, dextrose, fucitol, fucose, galactitol, galactose, mannitol, mannose, perseitol, sorbitol

KEY WORDS

post-column reaction

REFERENCE

Wang,W.T.; Safar,J.; Zopf,D. Analysis of inositol by high-performance liquid chromatography, *Anal.Biochem.*, **1990**, 188, 432–435.

Insulin

Molecular formula: $C_{257}H_{383}N_{65}O_{77}S_6$ (human)

Molecular weight: 5807.69 (human)

CAS Registry No.: 9004-10-8 (injection), 8049-62-5 (zinc suspension), 11061-68-0 (human), 12584-58-6 (pig), 11070-73-8 (cow), 9004-14-2 (neutral insulin), 8049-62-5 (isophane insulin), 9004-17-5 (protamine zinc suspension)

Merck Index: 5011

SAMPLE

Matrix: formulations

HPLC VARIABLES

Column: 250 × 4.5 µm Nucleosil RP18

Mobile phase: MeCN:buffer 24:76 (w/w) (Prepare buffer by dissolving 9.8 g 85% phosphoric acid and 28.4 g sodium sulfate in 1 L water.)

Column temperature: 40

Flow rate: 1

Injection volume: 20

Detector: UV 214

CHROMATOGRAM

Retention time: 10.89

OTHER SUBSTANCES

Simultaneous: degradation products

KEY WORDS

comparison with capillary electrophoresis

REFERENCE

Kunkel,A.; Günter,S.; Dette,C.; Wätzig,H. Quantitation of insulin by capillary electrophoresis and high-performance liquid chromatography. Method comparison and validation, *J.Chromatogr.A*, **1997**, 781, 445–455.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 150 × 4.6 Synchropak C4

Mobile phase: Gradient. A was 0.05% trifluoroacetic acid in water. B was 0.05% trifluoroacetic acid in MeCN. A:B from 74:26 to 38:62 over 15 min.

Flow rate: 1.5

Detector: UV 220

CHROMATOGRAM

Retention time: 4.7

REFERENCE

Ho,H.-O.; Hsiao,C.-C.; Sheu,M.-T. Preparation of microemulsions using polyglycerol fatty acid esters as surfactant for the delivery of protein drugs, *J.Pharm.Sci.*, **1996**, 85, 138–143.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 150 × 4.6 5 µm Zorbax 300 Å SB-C3

Mobile phase: Gradient. A was MeCN:water:trifluoroacetic acid 5:95:0.1. B was MeCN:water:trifluoroacetic acid 5:95:0.085. A:B from 85:15 to 47:53 over 20 min.

Column temperature: 35

Flow rate: 1

Injection volume: 10

Detector: UV 215

CHROMATOGRAM

Retention time: 10

OTHER SUBSTANCES

Simultaneous: angiotensin II, carbonic anhydrase, cytochrome C, leucine enkephalin, lysozyme, myoglobin, RNAase

REFERENCE

Ricker,R.D.; Sandoval,L.A.; Permar,B.J.; Boyes,B.E. Improved reversed-phase high performance liquid chromatography columns for biopharmaceutical analysis, *J.Pharm.Biomed.Anal.*, **1996**, 14, 93–105.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Protein & Peptide C18 (Vydac)

Mobile phase: MeCN:buffer 26:74 (Buffer was 28.4 g sodium sulfate and 2.7 mL phosphoric acid in 1 L water, pH adjusted to 2.3 with ethanolamine (if necessary).)

Column temperature: 40

Flow rate: 0.8

Detector: UV 214

CHROMATOGRAM

Retention time: 15

OTHER SUBSTANCES

Simultaneous: degradation products

REFERENCE

Yomota,C.; Yoshii,Y.; Takahata,T.; Okada,S. Separation of B-3 monodesamidoinulin from human insulin by high-performance liquid chromatography under alkaline conditions, *J.Chromatogr.A*, **1996**, 721, 89–96.

Interferon

Molecular formula: C₈₆₀H₁₃₅₃N₂₂₇O₂₅₅S₉

Molecular weight: 19241.28

CAS Registry No.: 76543-88-9 (αA), 99210-65-8 (α2B), 98059-61-1 (gamma-1B)

Merck Index: 5015

Lednicer No.: 4 1

SAMPLE

Matrix: solutions

Sample preparation: Terminate the enzymatic hydrolysis by adding 10 μL 1% trifluoroacetic acid, 100 μL 8 M guanidine hydrochloride, inject an aliquot of the enzymatic digests.

HPLC VARIABLES

Column: 150 × 4.6 5 μm Vydac C18

Mobile phase: Gradient. A was 0.1% trifluoroacetic acid. B was MeCN:water 90:10 containing 0.1% trifluoroacetic acid. A:B from 52:48 to 45:55 in 45 min, to 20:80 over 2.5 min.

Flow rate: 0.6

Detector: UV 214

CHROMATOGRAM

Retention time: 20

OTHER SUBSTANCES

Simultaneous: impurities

REFERENCE

Gitlin,G.; Tsaropoulos,A.; Patel,S.T.; Sydor,W.; Pramanik,B.N.; Jacobs,S.; Westreich,L.; Mittelman,S.; Bausch,J.N. Isolation and characterization of a monomethioninesulfoxide variant of interferon α-2b, *Pharm.Res.*, **1996**, 13, 762–769.


Iodamide

Molecular formula: $C_{12}H_{11}I_3N_2O_4$

Molecular weight: 627.94

CAS Registry No.: 440-58-4, 18656-21-8 (meglumine salt)

Merck Index: 5031



SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 1 mL Plasma + 100 μ L 500 μ g/mL iodopyracet in water + 5 mL MeOH, centrifuge at 1500 g for 20 min. Remove the supernatant and evaporate it to dryness under a stream of air at 65°, reconstitute the residue in 300 μ L MeOH:water 90:10, let stand for 10 min, inject a 10 μ L aliquot. Urine. 100 μ L Urine + 900 μ L mobile phase + 100 μ L 500 μ g/mL iodopyracet in water, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 4.6 5 μ m LiChrosorb RP-8

Mobile phase: MeOH:water 85:15 containing 10 mM tetrabutylammonium hydrogen sulfate and 10 mM Tris

Flow rate: 1.5

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 8

Internal standard: iodopyracet (11)

Limit of detection: 200 ng/mL

KEY WORDS

plasma

REFERENCE

Hekman,P.; Van Ginneken,C.A. Rapid determination of renal contrast media in biological fluids by means of high-performance liquid chromatography, *J.Chromatogr.*, **1980**, 182, 492-495.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 0.5-1 mL Plasma + 1 mL 1 M HCl + 5 mL ethyl acetate, shake for 10 min, centrifuge, repeat extraction. Combine the organic layers and evaporate them under a stream of nitrogen at 60°, add 3 mL 100 mM NaOH, shake for 10 min, centrifuge, discard the organic layer. Remove the aqueous layer and add it to 500 μ L 1 M HCl, add 5 mL ethyl acetate, shake for 10 min, centrifuge, repeat extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 60°, reconstitute the residue in 100 μ L mobile phase, inject a 10-20 μ L aliquot. Urine. 100 μ L Urine + 1 mL 1 M HCl + add 5 mL ethyl acetate, shake for 10 min, centrifuge, repeat extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 60°, reconstitute the residue in 400 μ L mobile phase, inject a 10-20 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 2.6 10 μ m ODS-HC Sil-X-1 (Perkin-Elmer)

Mobile phase: MeCN:water:85% phosphoric acid 4:96:0.03

Flow rate: 1

Injection volume: 10-20

Detector: UV 235

CHROMATOGRAM

Retention time: 2.5

Internal standard: iodamide

OTHER SUBSTANCES

Extracted: o-iodohippurate, iothalamate

KEY WORDS

plasma; iodamide is IS

REFERENCE

Boschi,S.; Marchesini,B. High-performance liquid chromatographic method for the simultaneous determination of iothalamate and o-iodohippurate, *J.Chromatogr.*, **1981**, 224, 139–143.


Iodinated glycerol

Molecular formula: $C_6H_{11}IO_3$

Molecular weight: 258.06

CAS Registry No.: 5634-39-9

Merck Index: 5033

**SAMPLE**

Matrix: bulk

Sample preparation: Prepare an aqueous solution, inject an aliquot.

HPLC VARIABLES

Guard column: 10 μ m Guard-pak (Waters)

Column: 250 \times 4.6 5 μ m Ultrasphere ODS C18

Mobile phase: MeCN:water 5:95

Flow rate: 1

Detector: RI

KEY WORDS

this procedure determines glycerol, a component of iodinated glycerol

REFERENCE


Cannon,J.M.; Brown,R.D.; Murrill,E.M.; Jameson,C.W. Identification of components in iodinated glycerol, *J.Pharm.Sci.*, **1989**, 78, 48–51.

Iodochlorhydroxyquin

Molecular formula: C_9H_5ClINO

Molecular weight: 305.50

Merck Index: 5052

**SAMPLE**

Matrix: blood

Sample preparation: 1 mL Plasma + 100 μ L concentrated perchloric acid, vortex for 30 s, centrifuge at 15° at 5000 rpm for 15 min. Remove 500 μ L of the supernatant and add it to 5 mL ether, vortex for 10 s, centrifuge at 15° for 10 min, repeat extraction. Add 10 mL ether to the original precipitate, vortex for 1 min, centrifuge at 15° at 5000 rpm for 15 min. Combine the extracts and dry them over anhydrous sodium sulfate, evaporate to dryness under a stream of nitrogen at 38°, reconstitute the residue in 500 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: 40 \times 5 RP-18-MPLC (Brownlee)

Column: 250 \times 2.6 ODS-HC-SIL-X-I (Perkin-Elmer)

Mobile phase: MeOH:50 mM phosphoric acid 80:20 (Flush column with MeOH at the end of each day.)
Column temperature: 40
Flow rate: 1
Injection volume: 20
Detector: UV 256

CHROMATOGRAM

Retention time: 5
Limit of quantitation: 1 µg/mL

OTHER SUBSTANCES

Noninterfering: hydrocortisone

KEY WORDS

plasma

REFERENCE

Ezzedeen,F.W.; Masoud,A.N.; Stohs,S.J.; Lerman,S.J. High-performance liquid chromatographic analysis of iodochlorhydroxyquin in plasma, *J.Pharm.Sci.*, **1981**, *70*, 889–891.

SAMPLE

Matrix: bulk, formulations

Sample preparation: Weigh out 30 mg of bulk drug or an amount of cream equivalent to 30 mg iodochlorhydroxyquin, add 70 mL THF, shake vigorously until the cream has dissolved, make up to 100 mL with THF. Remove a 5 mL aliquot and add it to 1 mL pyridine and 1 mL acetic anhydride, heat at 60° for 15 min, cool, add 15 mL 450 µg/mL testosterone acetate in 94:6 butyl chloride:THF, mix thoroughly. Remove a 3 mL aliquot and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 15 mL mobile phase with gentle warming and vigorous shaking, inject an aliquot.

HPLC VARIABLES

Column: 300 × 4 10 µm µPorasil

Mobile phase: Butyl chloride:water-saturated butyl chloride:THF:glacial acetic acid 55:55:3:2

Flow rate: 2-3

Detector: UV 254

CHROMATOGRAM

Retention time: 6

Internal standard: testosterone acetate (8)

OTHER SUBSTANCES

Simultaneous: impurities

KEY WORDS

cream; normal phase; derivatization; acetylation

REFERENCE

Kubiak,E.J.; Munson,J.W. Analysis of iodochlorhydroxyquin in cream formulations and bulk drugs by high-performance liquid chromatography, *J.Pharm.Sci.*, **1982**, *71*, 872–875.

SAMPLE

Matrix: bulk, formulations

Sample preparation: Weigh out ointment, cream, or bulk drug containing 30 mg of iodochlorhydroxyquin, add 40-50 mL THF, dissolve with heating, make up to 100 mL with THF. Remove a 5 mL aliquot and add it to 1 mL 10 mg/mL nickel chloride in MeOH, add 5 mL 0.4 mg/mL diphenylamine in MeOH, make up to 50 mL with MeOH, filter, inject an aliquot.

HPLC VARIABLES

Column: 300 × 3.9 µBondapak phenyl

Mobile phase: MeCN:MeOH:water 30:20:50 containing 1 mM NiCl₂

Flow rate: 1.2

Detector: UV 273

CHROMATOGRAM

Retention time: 7.5

Internal standard: diphenylamine (11)

OTHER SUBSTANCES

Simultaneous: chloroxine, iodoquinol

KEY WORDS

ointment; creams; separated as Ni chelates

REFERENCE

Wojtowicz,E.J. Reverse-phase high-performance liquid chromatographic determination of halogenated 8-hydroxyquinoline compounds in pharmaceuticals and bulk drugs, *J.Pharm.Sci.*, **1984**, 73, 1430–1433.

SAMPLE

Matrix: formulations

Sample preparation: Ointment. 50 mg Ointment + 10 mL ether, vortex until dissolved. Remove a 200 µL aliquot and add phenyl salicylate in mobile phase, evaporate to dryness under a stream of nitrogen at 40°, reconstitute the residue in 10 mL mobile phase, warm for 1 min on a steam bath, vortex for 1 min, cool. Remove an aliquot, dilute with mobile phase, inject an aliquot. Cream. Suspend 50 mg cream in 10 mL mobile phase by vortexing. Remove an aliquot and add phenyl salicylate in mobile phase, evaporate to dryness under a stream of nitrogen at 40°, suspend the residue in 10 mL mobile phase, warm for 1 min on a steam bath, vortex for 1 min, cool. Remove an aliquot, dilute with mobile phase, inject an aliquot.

HPLC VARIABLES

Guard column: 40 × 5 RP-18-MPLC (Brownlee)

Column: 250 × 2.6 ODS-HC-SIL-X (Perkin-Elmer)

Mobile phase: MeOH:50 mM phosphoric acid 70:30 (Flush column with MeOH at the end of each day.)

Column temperature: 40

Flow rate: 1

Injection volume: 20

Detector: UV 256

CHROMATOGRAM

Retention time: 7.5

Internal standard: phenyl salicylate (5.25)

OTHER SUBSTANCES

Simultaneous: hydrocortisone

KEY WORDS

ointment; cream

REFERENCE

Ezzedeen,F.W.; Stohs,S.J.; Masoud,A.N. High-performance liquid chromatographic analysis of iodochlorhydroxyquin and hydrocortisone in ointments and creams, *J.Pharm.Sci.*, **1983**, 72, 1036–1039.

SAMPLE

Matrix: feces, tissue, urine

Sample preparation: Urine. 1-5 mL Urine + 1-5 mL diethyl ether, vortex for 10 s, centrifuge at 15° at 3000 g for 10 min. Remove the organic layer and dry it over anhydrous sodium sulfate, evaporate to dryness under a stream of nitrogen at 40°, reconstitute the residue in mobile phase, inject a 20 µL aliquot. (Hydrolyze conjugates by adjusting pH to 5 with 1 M acetate buffer, add β-glucuronidase (Sigma) to a final concentration of 200 U/mL, heat at 37° for 2 h, neutralize with 3 M NaOH, add 1-5 mL diethyl ether, vortex for 10 s, centrifuge at 15° at 3000

g for 10 min. Remove the organic layer and dry it over anhydrous sodium sulfate, evaporate to dryness under a stream of nitrogen at 40°, reconstitute in 500 µL benzene (Caution! Benzene is a carcinogen!), add to an alumina column, wash with 2 mL benzene:pyridine 7:1, wash with 2 mL acetone, wash with 2 mL 100 mM acetic acid in MeOH, wash with 2 mL MeOH. Remove the alumina and add it to 1 mL saturated aqueous sodium fluoride, extract twice with 5 mL diethyl ether. Combine the extracts and evaporate them to dryness, reconstitute the residue in mobile phase, inject a 20 µL aliquot.) Tissue. Homogenize (Potter-Elvehjem) 1 g liver and 2 mL mobile phase, place homogenate in another tube, rinse original tube with 1 mL mobile phase, add rinse to homogenate, add 5 mL diethyl ether, vortex for 1 min, centrifuge at 3000 g for 10 min, repeat extraction twice. Combine the organic layers and evaporate them to dryness, reconstitute in 500 µL benzene (Caution! Benzene is a carcinogen!), add to an alumina column, wash with 2 mL benzene:pyridine 7:1, wash with 2 mL acetone, wash with 2 mL 100 mM acetic acid in MeOH, wash with 2 mL MeOH. Remove the alumina and add it to 1 mL saturated aqueous sodium fluoride, extract twice with 5 mL diethyl ether. Combine the extracts and evaporate them to dryness, reconstitute the residue in mobile phase, inject a 20 µL aliquot. Feces. Homogenize (Potter-Elvehjem) 1 g feces and 5 mL mobile phase, adjust pH to 5 with 1 M acetate buffer, add β-glucuronidase (Sigma) to a final concentration of 200 U/mL, heat at 37° for 2 h, neutralize with 3 M NaOH, add 1-5 mL diethyl ether, vortex for 10 s, centrifuge at 15° at 3000 g for 10 min. Remove the organic layer and dry it over anhydrous sodium sulfate, evaporate to dryness under a stream of nitrogen at 40°, reconstitute in 500 µL benzene (Caution! Benzene is a carcinogen!), add to an alumina column, wash with 2 mL benzene:pyridine 7:1, wash with 2 mL acetone, wash with 2 mL 100 mM acetic acid in MeOH, wash with 2 mL MeOH. Remove the alumina and add it to 1 mL saturated aqueous sodium fluoride, extract twice with 5 mL diethyl ether. Combine the extracts and evaporate them to dryness, reconstitute the residue in mobile phase, inject a 20 µL aliquot.

HPLC VARIABLES

Guard column: 40 × 5 RP-18-MPLC (Brownlee)

Column: 250 × 2.6 ODS-HC-SIL-X-I (Perkin-Elmer)

Mobile phase: MeOH:50 mM phosphoric acid 70:30 (Flush column with MeOH at the end of each day.)

Column temperature: 40

Flow rate: 1

Injection volume: 20

Detector: UV 256

CHROMATOGRAM

Retention time: 7.5

Limit of detection: 200 ng

Limit of quantitation: 250 ng/g (feces), 500 ng/g (liver), 200 ng/mL (urine)

KEY WORDS

liver; dog; SPE

REFERENCE

Ezzedeen, F.W.; Stohs, S.J.; Stublar, M. Analysis of iodochlorhydroxyquin in biological materials by high-performance liquid chromatography, *J. Chromatogr.*, **1983**, 276, 121-128.


Iodoquinol

Molecular formula: C₉H₅I₂NO

Molecular weight: 396.95

CAS Registry No.: 83-73-8

Merck Index: 5063



SAMPLE

Matrix: bulk, formulations

Sample preparation: Cream, bulk. Weigh out cream or bulk drug containing 20 mg of iodoquinol, add 40-50 mL THF, dissolve with heating, make up to 100 mL with THF. Remove a 5

mL aliquot and add it to 1 mL 10 mg/mL nickel chloride in MeOH, add 3 mL 0.4 mg/mL diphenylamine in MeOH, make up to 50 mL with MeOH, filter, inject an aliquot. Tablets. Finely powder tablets, weigh out amount equivalent to 200 mg iodoquinol, add 150 mL THF, heat on a steam bath for a few min, shake mechanically for 30 min, make up to 250 mL with THF. Remove a 25 mL aliquot and make up to 100 mL with THF. Remove a 5 mL aliquot and add it to 1 mL 10 mg/mL nickel chloride in MeOH, add 3 mL 0.4 mg/mL diphenylamine in MeOH, make up to 50 mL with MeOH, filter, inject an aliquot.

HPLC VARIABLES

Column: 300 × 3.9 μBondapak phenyl

Mobile phase: MeCN:MeOH:water 30:20:50 containing 1 mM NiCl₂

Flow rate: 1.2

Detector: UV 273

CHROMATOGRAM

Retention time: 9

Internal standard: diphenylamine (11)

OTHER SUBSTANCES

Simultaneous: chloroxine, iodochlorhydroxyquin

KEY WORDS

tablets; creams; separated as Ni chelates

REFERENCE


Wojtowicz, E.J. Reverse-phase high-performance liquid chromatographic determination of halogenated 8-hydroxyquinoline compounds in pharmaceuticals and bulk drugs, *J. Pharm. Sci.*, **1984**, 73, 1430–1433.

Iofratol

Molecular formula: C₃₁H₃₆I₆N₆O₁₃

Molecular weight: 1462.09

CAS Registry No.: 141660-63-1



SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. Add 30 μL 2 mg/mL IS in water to 100 μL plasma. Add 30 μL 35% perchloric acid. Agitate and centrifuge at 3500 g for 10 min. Inject a 10 μL aliquot of the clear supernatant. Urine. Dilute 1 mL urine with 2 mL water, centrifuge at 4500 g for 15 min. Add 100 μL 5 mg/mL IS, 100 μL glacial acetic acid, and an ion-exchange resin mixture (1 g Duolite A-30B + 900 mg Amberlite IR-120). Dilute the suspension to 5 mL with water. Agitate for 30 min and centrifuge at 3500 g for 5 min. Inject a 10 μL aliquot of the supernatant.

HPLC VARIABLES

Guard column: 30 × 4 7 μm LiChrosorb RP-8

Column: 250 × 4.6 5 μm LiChrosorb RP-8

Mobile phase: MeCN:5 mM pH 4.5 potassium dihydrogen phosphate 5:95

Column temperature: 45

Flow rate: 1

Injection volume: 10

Detector: UV 242

CHROMATOGRAM

Retention time: 12.0

Internal standard: iopamidol (4.1)

Limit of detection: 75 ng/mL (plasma), 430 ng/mL (urine)


KEY WORDS

plasma; rat; human

REFERENCE

Arbughi,T.; Bertani,F.; Celeste,R.; Grotti,A.; Tirone,P. High-performance liquid chromatographic determination of the X-ray imaging contrast agent, iofratol, in plasma and urine, *J.Chromatogr.B*, **1997**, 701, 103–113.

Iohexol

Molecular formula: C₁₉H₂₆I₃N₃O₉**Molecular weight:** 821.14**CAS Registry No.:** 66108-95-0**Merck Index:** 5068**SAMPLE****Matrix:** blood**Sample preparation:** Mix serum with an equal volume of MeCN, vortex for 15 s, centrifuge at 14000 g for 30 s, dilute the supernatant 100-fold with mobile phase, inject a 10 µL aliquot.**HPLC VARIABLES****Column:** 40 × 3.2 3 µm Velosep RP-18 (Applied Biosystems)**Mobile phase:** 8 mM pH 2 Phosphoric acid**Flow rate:** 1**Injection volume:** 10**Detector:** UV 254**KEY WORDS**

serum

REFERENCE

Shihabi,Z.K.; Constantinescu,M.S. Iohexol in serum determined by capillary electrophoresis, *Clin.Chem.*, **1992**, 38, 2117–2120.

SAMPLE**Matrix:** blood**Sample preparation:** 50 µL Serum + 50 µL 250 µg/mL acetaminophen in 100 mM HCl, add to SPE cartridge containing 150 mg 80-100 mesh Chromosorb P/NAW, elute with 1 mL ethyl acetate:MeOH 5:1, add the eluate to 50 µL 100 mM HCl, vortex for 15 s, centrifuge at 10000 g for 3 min, inject a 20 µL aliquot of the lower aqueous phase.**HPLC VARIABLES****Column:** 5 µm C8**Mobile phase:** MeCN:20 mM pH 3.3 phosphoric acid 2.5:97.5**Injection volume:** 20**Detector:** UV 254**CHROMATOGRAM****Internal standard:** acetaminophen**Limit of detection:** <1 µg/mL**OTHER SUBSTANCES****Extracted:** aminohippuric acid (PAH)**KEY WORDS**

serum; SPE