


KEY WORDS

plasma; rat; human

REFERENCE

Arbughi,T.; Bertani,F.; Celeste,R.; Grotti,A.; Tirone,P. High-performance liquid chromatographic determination of the X-ray imaging contrast agent, iofratol, in plasma and urine, *J.Chromatogr.B*, **1997**, 701, 103–113.

Iohexol

Molecular formula: $C_{19}H_{26}I_3N_3O_9$ **Molecular weight:** 821.14**CAS Registry No.:** 66108-95-0**Merck Index:** 5068**SAMPLE****Matrix:** blood**Sample preparation:** Mix serum with an equal volume of MeCN, vortex for 15 s, centrifuge at 14000 g for 30 s, dilute the supernatant 100-fold with mobile phase, inject a 10 μ L aliquot.**HPLC VARIABLES****Column:** 40 \times 3.2 3 μ m Velosep RP-18 (Applied Biosystems)**Mobile phase:** 8 mM pH 2 Phosphoric acid**Flow rate:** 1**Injection volume:** 10**Detector:** UV 254**KEY WORDS**

serum

REFERENCE

Shihabi,Z.K.; Constantinescu,M.S. Iohexol in serum determined by capillary electrophoresis, *Clin.Chem.*, **1992**, 38, 2117–2120.

SAMPLE**Matrix:** blood**Sample preparation:** 50 μ L Serum + 50 μ L 250 μ g/mL acetaminophen in 100 mM HCl, add to SPE cartridge containing 150 mg 80-100 mesh Chromosorb P/NAW, elute with 1 mL ethyl acetate:MeOH 5:1, add the eluate to 50 μ L 100 mM HCl, vortex for 15 s, centrifuge at 10000 g for 3 min, inject a 20 μ L aliquot of the lower aqueous phase.**HPLC VARIABLES****Column:** 5 μ m C8**Mobile phase:** MeCN:20 mM pH 3.3 phosphoric acid 2.5:97.5**Injection volume:** 20**Detector:** UV 254**CHROMATOGRAM****Internal standard:** acetaminophen**Limit of detection:** <1 μ g/mL**OTHER SUBSTANCES****Extracted:** aminohippuric acid (PAH)**KEY WORDS**

serum; SPE

REFERENCE

Andreeva,M.; Rapondjieva,A.; Deskova,D.; Tishkov,I.; Svinarov,D. Liquid chromatographic determination of iohexol and PAH with Chromosorb P column used for sample preparation (Abstract 175), *Ther.Drug Monit.*, 1995, 17, 427.

SAMPLE

Matrix: blood, tissue

Sample preparation: Plasma. 50 μ L Plasma + 150 μ L buffer, centrifuge at 1000 g at 0° for 10 min, inject a 10 μ L aliquot. Tissue. Weigh 2 testes, add 5 mL buffer, homogenize (Kinematica type PT 10/35, setting 7.5) for 1 min, vortex, centrifuge at 1000 g at 0° for 10 min. Remove a 50 μ L aliquot of the supernatant and add it to 150 μ L buffer, centrifuge at 1000 g at 0° for 10 min, inject a 10 μ L aliquot. (Buffer was 5 mM metaphosphoric acid and 5 mM disodium EDTA.)

HPLC VARIABLES

Guard column: C18 Bondapak guard column

Column: 300 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: MeCN:buffer 2.5:97.5–5:95 (Buffer was 100 mM NaH_2PO_4 and 0.2 mM Na_2EDTA adjusted to pH 3.1 with orthophosphoric acid.)

Flow rate: 1

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 7.5 (two isomers seen)

Limit of detection: 670 ng/mL

OTHER SUBSTANCES

Also analyzed: iopamidol, diatrizoate

KEY WORDS

plasma; mouse; testes

REFERENCE

Harapanhalli,R.S.; Yaghmai,V.; Patel,Y.D.; Baker,S.R.; Rao,D.V. Assay of radiographic contrast agents in mice plasma and testes by high-performance liquid chromatography, *Anal.Chem.*, 1993, 65, 606–612.


Iopamidol

Molecular formula: $\text{C}_{17}\text{H}_{22}\text{I}_3\text{N}_3\text{O}_8$

Molecular weight: 777.09

CAS Registry No.: 60166-93-0

Merck Index: 5071



SAMPLE

Matrix: blood, tissue

Sample preparation: Plasma. 50 μ L Plasma + 150 μ L buffer, centrifuge at 1000 g at 0° for 10 min, inject a 10 μ L aliquot. Tissue. Weigh 2 testes, add 5 mL buffer, homogenize (Kinematica type PT 10/35, setting 7.5) for 1 min, vortex, centrifuge at 1000 g at 0° for 10 min. Remove a 50 μ L aliquot of the supernatant and add it to 150 μ L buffer, centrifuge at 1000 g at 0° for 10 min, inject a 10 μ L aliquot. (Buffer was 5 mM metaphosphoric acid and 5 mM disodium EDTA.)

HPLC VARIABLES

Guard column: C18 Bondapak guard column

Column: 300 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: MeCN:buffer 2.5:97.5–5:95 (Buffer was 100 mM NaH_2PO_4 and 0.2 mM Na_2EDTA adjusted to pH 3.1 with orthophosphoric acid.)

Flow rate: 1
Injection volume: 10
Detector: UV 254

CHROMATOGRAM

Retention time: 9.5
Limit of detection: 640 ng/mL

OTHER SUBSTANCES

Also analyzed: iohexaol, diatrizoate

KEY WORDS

plasma; mouse; testes

REFERENCE

Harapanhalli,R.S.; Yaghmai,V.; Patel,Y.D.; Baker,S.R.; Rao,D.V. Assay of radiographic contrast agents in mice plasma and testes by high-performance liquid chromatography, *Anal.Chem.*, **1993**, 65, 606–612.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. Add 30 μ L water to 100 μ L plasma. Add 30 μ L 35% perchloric acid. Agitate and centrifuge at 3500 g for 10 min. Inject a 10 μ L aliquot of the clear supernatant. Urine. Dilute 1 mL urine with 2 mL water, centrifuge at 4500 g for 15 min. Add 100 μ L 5 mg/mL IS, 100 μ L glacial acetic acid, and an ion-exchange resin mixture (1 g Duolite A-30B + 900 mg Amberlite IR-120). Dilute the suspension to 5 mL with water. Agitate for 30 min and centrifuge at 3500 g for 5 min. Inject a 10 μ L aliquot of the supernatant.

HPLC VARIABLES

Guard column: 30 \times 4.0 7 μ m LiChrosorb RP-8

Column: 250 \times 4.6 5 μ m LiChrosorb RP-8

Mobile phase: MeCN:5 mM pH 4.5 potassium dihydrogen phosphate 5:95

Column temperature: 45

Flow rate: 1.0

Injection volume: 10

Detector: UV 242

CHROMATOGRAM

Retention time: 4.1

OTHER SUBSTANCES

Extracted: iofratol

KEY WORDS

iopamidol is IS; plasma; rat; human

REFERENCE


Arbuchi,T.; Bertani,F.; Celeste,R.; Grotti,A.; Tirone,P. High-performance liquid chromatographic determination of the X-ray imaging contrast agent, iofratol, in plasma and urine, *J.Chromatogr.B*, **1997**, 701, 103–113.

Iothalamate

Molecular formula: C₁₁H₉I₃N₂O₄.C₇H₇NO₅ (meglumine),
C₁₁H₈I₃N₂NaO₄ (sodium salt)

Molecular weight: 809.13 (meglumine), 635.90 (sodium salt)

CAS Registry No.: 13087-53-1, 2276-90-6 (iothamic acid), 6284-40-8 (meglumine),
17692-74-9 (²⁵¹I radioactive agent), 15845-98-4 (¹³¹I radioactive agent), 1225-20-3 (sodium salt)



SAMPLE

Matrix: blood, urine

Sample preparation: Dilute urine 10-fold with water. Add 200 μL MeCN containing 20 $\mu\text{g/mL}$ p-aminobenzoic acid to 100 μL diluted urine, vortex briefly, centrifuge at 12000 g for 4 min, inject a 20 μL aliquot of the supernatant.

HPLC VARIABLES

Guard column: 5 μm Ultrasphere C18

Column: 250 mm 5 μm Primesphere C18 (Torrance, USA)

Mobile phase: MeOH:buffer 18:82 (Buffer was 50 mM NaH_2PO_4 with 0.5 mM tetrabutyl ammonium hydrogen sulfate with an unadjusted pH of 4.11.)

Flow rate: 0.8

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 13

Internal standard: p-aminobenzoic acid (9.8)

OTHER SUBSTANCES

Extracted: p-aminohippuric acid

KEY WORDS

serum

REFERENCE

Agarwal,R. Chromatographic estimation of iothalamate and p-aminohippuric acid to measure glomerular filtration rate and effective renal plasma flow in humans, *J.Chromatogr.B*, **1998**, 705, 3–9.

SAMPLE

Matrix: formulations

Sample preparation: Mix a dilution of the injectable solution in water with a mixture of 2,4-dinitrobenzenesulfonyl chloride (DNBS-Cl) and sodium carbonate in acetone at ambient temperature for 1 hr. Dilute the reaction mixture with 25% MeOH and inject a 20 μL aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 4 μm Nova-Pak C18

Mobile phase: MeOH:buffer 25:75 (Buffer was 20 mM tetrapropylammonium hydroxide adjusted to pH 3.4 with orthophosphoric acid.)

Flow rate: 1.3

Injection volume: 20

Detector: UV 262

CHROMATOGRAM

Retention time: 2.75

Limit of quantitation: 1 μg

OTHER SUBSTANCES

Simultaneous: diatrizoate, meglumine

KEY WORDS

injections; only meglumine is derivatized under these conditions

REFERENCE

Lau-Cam,C.A.; Roos,R.W. HPLC method with spectrophotometric detection for the simultaneous assay of meglumine and its counterions iothalamic acid or diatrizoic acid in radiographic solutions for injection (Abstract 3372), *Pharm.Res.*, **1997**, 14, S586.

SAMPLE**Matrix:** blood**Sample preparation:** 100 μ L Plasma or urine + 0.5-5 μ g p-aminobenzoic acid + 200 μ L MeCN, vortex for a few s, centrifuge at 800 g for 5 min, inject a 5 μ L aliquot of the supernatant.

HPLC VARIABLES**Column:** 300 \times 3.9 10 μ m μ Bondapak C18**Mobile phase:** MeCN:0.04% pH 2.5 \pm 0.05 phosphoric acid 3.5:96.5**Flow rate:** 1.5**Injection volume:** 5**Detector:** UV 254

CHROMATOGRAM**Retention time:** 6**Internal standard:** p-aminobenzoic acid (8)**Limit of detection:** 500 ng/mL

OTHER SUBSTANCES**Extracted:** p-aminohippuric acid

KEY WORDS

plasma; dog; human; pharmacokinetics

REFERENCE

Prueksaritanont,T.; Chen,M.L.; Chiou,W.L. Simple and micro high-performance liquid chromatographic method for simultaneous determination of p-aminohippuric acid and lothalamate in biological fluids, *J.Chromatogr.*, **1984**, 306, 89-97.

SAMPLE**Matrix:** blood**Sample preparation:** Add barbital to plasma. 100 μ L Plasma + 500 μ L MeOH, vortex for 15 s, centrifuge at 2500 rpm for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in buffer, inject a 15-20 μ L aliquot.

HPLC VARIABLES**Column:** 250 \times 4.6 5 μ m Octyl C8 (Rainin)**Mobile phase:** MeOH:MeCN:buffer 90:10:300 (Buffer was 6.44 g KH_2PO_4 , 7.04 g K_2HPO_4 , and 14 mL 500 mM dodecyltriethylammonium phosphate (Regis) in 4 L water.)**Flow rate:** 1**Injection volume:** 15-20**Detector:** UV 254

CHROMATOGRAM**Retention time:** 18.2**Internal standard:** barbital (15.9)**Limit of quantitation:** 3000 ng/mL

OTHER SUBSTANCES**Extracted:** p-aminohippuric acid

KEY WORDS

plasma

REFERENCE

Jayewardene,A.L.; Seneviratne,A.K.; Gambertoglio,J.G. Paired ion reversed-phase HPLC assay for the simultaneous determination of lothalamate acid and para aminohippuric acid in plasma, *J.Liq.Chromatogr.*, **1994**, 17, 2395-2412.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 0.5-1 mL Plasma + 18 µg iodamide + 2 µg hippuric acid + 1 mL 1 M HCl + 5 mL ethyl acetate, shake for 10 min, centrifuge, repeat extraction. Combine the organic layers and evaporate them under a stream of nitrogen at 60°, add 3 mL 100 mM NaOH, shake for 10 min, centrifuge, discard the organic layer. Remove the aqueous layer and add it to 500 µL 1 M HCl, add 5 mL ethyl acetate, shake for 10 min, centrifuge, repeat extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 60°, reconstitute the residue in 100 µL mobile phase, inject a 10-20 µL aliquot. Urine. 1 mL Urine + 10 µL 10% iodamide in water, mix. Remove a 100 µL aliquot and add it to 1 mL 1 M HCl, add 5 mL ethyl acetate, shake for 10 min, centrifuge, repeat extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 60°, reconstitute the residue in 400 µL mobile phase, inject a 10-20 µL aliquot.

HPLC VARIABLES

Column: 250 × 2.6 10 µm ODS-HC Sil-X-1 (Perkin-Elmer)

Mobile phase: MeCN:water:85% phosphoric acid 4:96:0.03

Flow rate: 1

Injection volume: 10-20

Detector: UV 235

CHROMATOGRAM

Retention time: 2.0

Internal standard: iodamide (2.5), hippuric acid (5.6)

Limit of detection: 500 ng/mL

OTHER SUBSTANCES

Extracted: o-iodohippurate

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Boschi,S.; Marchesini,B. High-performance liquid chromatographic method for the simultaneous determination of iothalamate and o-iodohippurate, *J.Chromatogr.*, **1981**, *224*, 139-143.

SAMPLE

Matrix: blood, urine

Sample preparation: 200 µL Plasma + 200 µL 100 mM NaH₂PO₄ + 400 µL MeCN, mix for 5 s, let stand at 4° for 15 min, centrifuge at 10500 g for 1 min. Remove the supernatant and add it to 2 mL dichloromethane, mix for 5 min, centrifuge at 4800 g for 10 min, inject a 5-50 µL aliquot of the upper aqueous phase. Urine. Centrifuge urine at 4800 g for 10 min, dilute 1:10 with 50 mM NaH₂PO₄, inject an aliquot.

HPLC VARIABLES

Column: 125 × 4 5 µm LiChrosorb RP-18

Mobile phase: MeCN:water 7.5:92.5 containing 5.50 g/L NaH₂PO₄·H₂O, 1.80 g/L Na₂HPO₄·2H₂O, and 20 mg/L tetrabutylammonium bromide, pH 6.4 (plasma) or MeCN:water 5:95 containing 5.50 g/L NaH₂PO₄·H₂O, 1.80 g/L Na₂HPO₄·2H₂O, and 22.5 mg/L tetrabutylammonium bromide, pH 6.4 (urine)

Flow rate: 1

Injection volume: 5-50

Detector: UV 254

CHROMATOGRAM

Retention time: 2.3 (plasma), 2.7 (urine)

Limit of detection: 1000 ng/mL

OTHER SUBSTANCES

Extracted: cefotetan (UV 280)

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Kees,F.; Grobecker,H.; Naber,K.G. High-performance liquid chromatographic analysis of cefotetan epimers in human plasma and urine, *J.Chromatogr.*, **1984**, 305, 363–371.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 500 μ L Plasma + 500 μ L 1 M HCl, vortex for 10 s, add 6 mL ethyl acetate, vortex for 20 s, centrifuge at 4° at 1700 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 400 μ L 25 mM pH 3 KH_2PO_4 , add 500 μ L dichloromethane, shake gently horizontally for 5 min, centrifuge at 1700 g for 5 min, inject a 20 μ L aliquot of the aqueous phase. Urine. Dilute 1:10 with water. Remove a 1 mL aliquot and add it to 1 mL 1 M HCl, vortex for 10 s, add 6 mL ethyl acetate, vortex for 20 s, centrifuge at 4° at 1700 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 400 μ L 25 mM pH 3 KH_2PO_4 , add 500 μ L dichloromethane, shake gently horizontally for 5 min, centrifuge at 1700 g for 5 min, inject a 20 μ L aliquot of the aqueous phase.

HPLC VARIABLES

Guard column: 10 \times 3 anion-exchange guard column (Chrompack)

Column: 250 \times 4.6 Partisil 10 SAX

Mobile phase: MeCN:25 mM pH 3 phosphate buffer 15:85

Flow rate: 1

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 7.05

Limit of detection: 500 ng/mL

OTHER SUBSTANCES

Extracted: p-aminohippuric acid

Noninterfering: acipimox, allopurinol, aspirin, atenolol, captopril, chlorthalidone, clonidine, digoxin, digoxin, diltiazem, dipyridamole, enalapril, furosemide, gemfibrozil, hydralazine, hydrochlorothiazide, ibopamine, insulin, inulin, isosorbide dinitrate, α -methyldopa, nicardipine, nifedipine, prazosin, propranolol, salicylic acid, simvastatin, trinitrin, verapamil

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Gaspari,F.; Mainardi,L.; Ruggerenti,P.; Remuzzi,G. High-performance liquid chromatographic determination of lothalamate acid in human plasma and urine, *J.Chromatogr.*, **1991**, 570, 435–440.

SAMPLE

Matrix: urine

Sample preparation: Dilute urine 1:100 or 1:500. 200 μ L Diluted urine + 50 μ L barbital solution, vortex for 15 s, inject a 20–30 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Ultrasphere C8

Mobile phase: MeCN:MeOH:10 mM pH 7.5 potassium phosphate buffer:0.5 M dodecyl triethylammonium phosphate 6:94:300:0.6 (0.5 M Dodecyl triethylammonium phosphate was Q-12, Ion pair reagent, Regis Chemical Co.)

Flow rate: 1

Injection volume: 20–30

Detector: UV 254

CHROMATOGRAM

Retention time: 16.8

Internal standard: barbital (14.5)

Limit of quantitation: 50000 ng/mL

OTHER SUBSTANCES

Extracted: p-aminohippuric acid

REFERENCE

Seneviratne,A.K.; Jayewardene,A.L.; Gambertoglio,J.G. Paired-ion reversed-phase HPLC assay for the determination of iothalamic acid and para aminohippuric acid in urine, *J.Pharm.Biomed.Anal.*, **1994**, *12*, 1311-1316.


Ipecac

Molecular formula: $C_{28}H_{38}N_2O_4$, $C_{29}H_{40}N_2O_4$

Molecular weight: 466.62, 480.65

CAS Registry No.: 483-17-0, 483-18-1

Merck Index: 5086



SAMPLE

Matrix: blood, emesis, urine

Sample preparation: 2 mL Plasma, whole blood, urine, or emesis (diluted 1:100) + 100 μ L 10 ng/mL N-propylprocainamide in water + 2 mL saturated sodium borate + 7 mL n-butyl chloride, vortex for 30 s, centrifuge. Remove the organic phase and add it to 200 μ L 10 mM HCl, vortex, centrifuge, inject a 30-50 μ L aliquot of the aqueous phase.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Ultrasphere ODS

Mobile phase: MeOH:25 mM pH 8.0 Na_2HPO_4 72:28

Flow rate: 1.7

Injection volume: 30-50

Detector: F ex 285 em 316

CHROMATOGRAM

Retention time: 3 (cephaeline), 4.1 (emetine)

Internal standard: N-propylprocainamide (1.7)

Limit of quantitation: 5 ng/mL

KEY WORDS

plasma; whole blood

REFERENCE

Crouch,D.J.; Moran,D.M.; Finkle,B.S.; Peat,M.A. Quantitative analysis of emetine and cephaeline by reversed-phase high performance liquid chromatography with fluorescence detection, *J.Anal.Toxicol.*, **1984**, *8*, 63-65.

SAMPLE

Matrix: blood, vomit

Sample preparation: Mix 250 μ L vomit or 500 μ L serum with an equal volume of 10% ammonium hydroxide for 5 min, add 4 mL ether, mix. Remove the organic layer and evaporate it

to dryness under vacuum at 25°, reconstitute the residue in 10 (vomit) or 0.2 (serum) mL 0.01% HCl. Mix 1 mL vomit solution with 50 μ L 4 mg/mL acrinol or 200 μ L serum solution with 50 μ L 200 μ g/mL acrinol, inject a 100 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m TSK gel ODS-80TM

Mobile phase: MeOH:buffer 54:46 (Buffer was 10 mM sodium 1-heptanesulfonate adjusted to pH 4 with glacial acetic acid.)

Flow rate: 1

Injection volume: 100

Detector: F ex 285 em 316

CHROMATOGRAM

Internal standard: acrinol

KEY WORDS

dog; serum; pharmacokinetics; assay determines cephaeline, a constituent of ipecac

REFERENCE

Teshima,D.; Suzuki,A.; Otsubo,K.; Higuchi,S.; Aoyama,T.; Shimozone,Y.; Saita,M.; Noda,K. Efficacy of emetic and United State Pharmacopoeia ipecac syrup in prevention of drug absorption, *Chem.Pharm.Bull.(Tokyo)*, **1990**, 38, 2242–2245.

SAMPLE

Matrix: cell cultures

Sample preparation: Dry 1-5 g cell culture at 60°, powder, extract with MeOH in a glass percolator for 72 h. Filter (paper) extract, wash solid with MeOH at 45°, concentrate filtrate under reduced pressure, filter (0.5 μ m), dilute filtrate with MeOH, inject a 2.5 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4 10 μ m silica gel (Spectra Physics)

Mobile phase: Chloroform:MeOH:diethylamine 90:10:0.2

Flow rate: 0.5

Injection volume: 2.5

Detector: UV 280

CHROMATOGRAM

Retention time: 4.5 (emetine), 5 (cephaeline)

KEY WORDS

normal phase

REFERENCE

Jha,S.; Sahu,N.P.; Mahato,S.B. Production of the alkaloids emetine and cephaeline in callus cultures of *Cephaelis ipecacuanha*, *Planta Med.*, **1988**, 54, 504–506.

SAMPLE

Matrix: formulations

Sample preparation: Syrup. Dilute syrup with an equal volume of water. 2 mL Diluted syrup + 2 mL 1% dansyl chloride in acetone + 200 μ L 1.5 M sodium carbonate, heat at 45 \pm 2° in the dark for 20 min, cool, add 3 mL water, add 500 μ L benzene (Caution! Benzene is a carcinogen!), shake. Remove the organic layer and evaporate it to dryness, reconstitute the residue in mobile phase, inject a 10 μ L aliquot. Capsules. Sonicate the contents of a capsule in 20 mL water for 10 min, centrifuge at 2000 g for 10 min. 2 mL Supernatant + 2 mL 1% dansyl chloride in acetone + 200 μ L 1.5 M sodium carbonate, heat at 45 \pm 2° in the dark for 20 min, cool, add 3 mL water, add 500 μ L benzene (Caution! Benzene is a carcinogen!), shake. Remove the organic layer and evaporate it to dryness, reconstitute the residue in mobile phase, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 2.8 10 μ m silica gel SI 100 (Merck)

Mobile phase: Diisopropyl ether:isopropanol:concentrated ammonia 48:2:0.3 (Caution! Diisopropyl ether readily forms explosive peroxides!)

Injection volume: 10

Detector: UV 254 or F ex 358 em 492 (cephaeline) or ex 356 em 481 (emetine)

CHROMATOGRAM

Retention time: 2.5 (cephaeline), 3.5 (emetine)

OTHER SUBSTANCES

Simultaneous: ephedrine, codeine (not derivatized, detect at UV 254 only)

KEY WORDS

derivatization; syrup; capsules; normal phase

REFERENCE

Frei, R.W.; Santi, W.; Thomas, M. Liquid chromatography of dansyl derivatives of some alkaloids and the application to the analysis of pharmaceuticals, *J.Chromatogr.*, **1976**, 116, 365-377.

SAMPLE

Matrix: formulations

Sample preparation: Mix 10 g linctus with 10-20 mL 200 µg/mL ethyl 4-hydroxybenzoate in MeCN:mobile phase 10:90, make up to 100 mL with mobile phase, inject a 10-25 µL aliquot.
Mix 10 pastilles with 50 mL 200 µg/mL ethyl 4-hydroxybenzoate in MeCN:mobile phase 10:90, make up to 200 mL with mobile phase, inject a 10-25 µL aliquot.

HPLC VARIABLES

Column: 150 × 3.9 µBondapak C18

Mobile phase: MeOH:water 40:60 containing 1 g/L sodium 1-heptanesulfonate and 1 mL/L orthophosphoric acid

Column temperature: 35

Flow rate: 2

Injection volume: 10-25

Detector: F ex 276 em 304

CHROMATOGRAM

Retention time: 8 (cephaeline), 15 (emetine)

Internal standard: ethyl 4-hydroxybenzoate (5)

Limit of quantitation: 5 µg/g

OTHER SUBSTANCES

Simultaneous: dihydroemetine (UV 214), emetamine (UV 214), O-methylpsychotrine (UV 214), tetrahydroemetine (UV 214)

KEY WORDS

stability-indicating; linctus; pastilles; rugged

REFERENCE

Elvidge, D.A.; Johnson, G.W.; Harrison, J.R. Selective, stability-indicating assay of the major ipecacuanha alkaloids, emetine and cephaeline, in pharmaceutical preparations by high-performance liquid chromatography using spectrofluorimetric detection, *J.Chromatogr.*, **1989**, 463, 107-118.

SAMPLE

Matrix: formulations

Sample preparation: 2 mL Sample + 1 mL 200 µg/mL emetine hydrochloride in water, make up to 10 mL with mobile phase, filter (0.45 µm), inject a 50-100 µL aliquot.

HPLC VARIABLES

Column: 100 × 4.6 5 µm Technosphere RP C-8 (HPLC Technology)

Mobile phase: MeCN:40 mM tetramethylammonium bromide:1 M acetic acid 80:15:5 (apparent pH 4.5)

Flow rate: 1.5
Injection volume: 50-100
Detector: UV 260

CHROMATOGRAM

Retention time: 1.83
Internal standard: emetine

OTHER SUBSTANCES

Simultaneous: benzalkonium (C12, C14, C16), tetrahydrozoline, naphazoline

KEY WORDS

nasal; ophthalmic; emetine is IS

REFERENCE

Santoni,G.; Medica,A.; Gratteri,P.; Furlanetto,S.; Pinzauti,S. High-performance liquid chromatographic determination of benzalkonium and naphazoline or tetrahydrozoline in nasal and ophthalmic solutions, *Farmaco*, **1994**, 49, 751-754.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 300 × 4 C18 Micropak SP
Mobile phase: MeCN:3 mM pH 2.5 phosphoric acid 20:80
Flow rate: 2
Injection volume: 15
Detector: F ex 285 em 315 or UV 205

CHROMATOGRAM


Retention time: 4.1 (cephaeline), 6.2 (emetine)
Limit of detection: 30 ng/mL

REFERENCE

Lachman,M.F.; Romeo,R.; McComb,R.B. Emetine identified in urine by HPLC, with fluorescence and ultraviolet/diode array detection, in a patient with cardiomyopathy, *Clin.Chem.*, **1989**, 35, 499-502.

Ipriflavone

Molecular formula: C₁₈H₁₆O₃
Molecular weight: 280.32
CAS Registry No.: 35212-22-7
Merck index: 5090



SAMPLE

Matrix: bulk
Sample preparation: Inject a 10 µL aliquot of a solution in MeOH.


HPLC VARIABLES

Column: 30 × 4.6 3 µm ODS C18 (Perkin-Elmer)
Mobile phase: MeCN:buffer 50:50 (Buffer was 10 mM triethylamine adjusted to pH 2.5 with orthophosphoric acid.)
Flow rate: 1.2
Injection volume: 10
Detector: UV 225

CHROMATOGRAM**Retention time:** 3.4**OTHER SUBSTANCES****Simultaneous:** impurities**REFERENCE**

Sustacha,K.; Chacón,M.; Lucero,M.L.; Orjales,A. Determination of ipriflavone and its synthetic impurities by high-performance liquid chromatography using diode-array detection, *J.Chromatogr.A*, **1996**, 719, 245–250.

Iprindole

Molecular formula: C₁₉H₂₈N₂**Molecular weight:** 284.44**CAS Registry No.:** 5560-72-5, 20432-64-8 (HCl)**Merck Index:** 5091**Lednicer No.:** 1 318**SAMPLE****Matrix:** solutions**Sample preparation:** Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.**HPLC VARIABLES****Column:** 125 × 4.9 Spherisorb S5W silica**Mobile phase:** MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7**Flow rate:** 2**Injection volume:** 20**Detector:** E, LeCarbone, V25 glassy carbon electrode, + 1.2 V**CHROMATOGRAM****Retention time:** 4.5**OTHER SUBSTANCES**

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipanonone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, isothipendyl, isoxsuprine, ketanserin, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypropazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide,

phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimino-dine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pi-renzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, pri-maquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, pro-thipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, qui-nine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldi-amine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranlycypromine, tra-zodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, tri-methoprim, trimipramine, tripelethamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, 323, 191–225.

Iproniazid


Molecular formula: C₉H₁₃N₃O

Molecular weight: 179.22

CAS Registry No.: 54-92-2, 305-33-9 (phosphate)

Merck Index: 5094

Lednicer No.: 1 254



SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 264

CHROMATOGRAM

Retention time: 3.00

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order: tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoyllecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; car-teolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihy-dralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazo-lam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam;

tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vandesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phencyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacemone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrridine; phenylbutazone; dexamethylamine; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozone; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpi-pramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254-262.


Irbesartan

Molecular formula: C₂₅H₂₈N₆O

Molecular weight: 428.54

CAS Registry No.: 138402-11-6

Merck Index: 5097



SAMPLE

Matrix: blood, urine

Sample preparation: Condition a 3 mL 100 mg Isolute Cyano SPE cartridge (Jones Chromatography) with 2 mL MeOH and 2 mL 0.85% phosphoric acid. Plasma. Mix 1.0 mL 0.85% phosphoric acid and 100 µL 1 µg/mL IS in 0.85% phosphoric acid with 250 µL plasma. Pass the mixture slowly through the SPE cartridge, wash with 3 mL 0.85% phosphoric acid, wash with 1 mL hexane. Elute with 1 mL of mixture MeOH:0.85% phosphoric acid 50:50. Inject a 20 µL aliquot of the eluate. Urine. Mix 1.0 mL 0.85% phosphoric acid and 100 µL 1 µg/mL IS in 0.85% phosphoric acid with 250 µL urine. Pass the mixture slowly through the SPE cartridge, wash with 4 mL 0.85% phosphoric acid. Elute with 1 mL of mixture MeOH:0.85% phosphoric acid 50:50. Inject a 20 µL aliquot of the eluate.

HPLC VARIABLES

Column: 150 × 4.6 5 µm YMC-ODS-AQ

Mobile phase: MeCN:buffer 50:50 (Prepare buffer by adding 1 mL triethylamine to 1 L water, adjust pH to 3.5 with 85% phosphoric acid.)

Flow rate: 0.8

Injection volume: 20

Detector: F ex 250 em 371


CHROMATOGRAM**Retention time:** 6.6**Internal standard:** BMS-190462 (9.4)**Limit of quantitation:** 1 ng/mL**KEY WORDS**

plasma; urine; SPE; pharmacokinetics

REFERENCE

Chang,S.-Y.; Whigan,D.B.; Vachharajani,N.N.; Patel,R. High-performance liquid chromatographic assay for the quantitation of irbesartan (SR 47436/BMS-186295) in human plasma and urine, *J.Chromatogr.B*, **1997**, 702, 149–155.

Irinotecan

Molecular formula: $C_{33}H_{38}N_4O_6$ **Molecular weight:** 586.69**CAS Registry No.:** 97682-44-5, 100286-90-6 (HCl)**Merck Index:** 5104**SAMPLE****Matrix:** bile**Sample preparation:** Dilute bile, inject an aliquot.**HPLC VARIABLES****Column:** 300 × 7.2 TSKgel-ODS 80Tm (Tosoh)**Mobile phase:** Gradient. MeOH:100 mM pH 4.0 phosphate buffer from 50:50 to 60:40 over 20 min, maintain at 60:40 for 20 min.**Flow rate:** 3**Detector:** UV 254**CHROMATOGRAM****Retention time:** 16**OTHER SUBSTANCES****Extracted:** metabolites**KEY WORDS**

rat; preparative for metabolites

REFERENCE

Atsumi,R.; Suzuki,W.; Hokusui,H. Identification of the metabolites of irinotecan, a new derivative of camptothecin, in rat bile and its biliary excretion, *Xenobiotica*, **1991**, 21, 1159–1169.

SAMPLE**Matrix:** blood

Sample preparation: Mix 10 μ L plasma with 100 μ M diisopropyl fluorophosphate, 125 ng/mL camptothecin and 375 μ L MeOH with an automatic mixer (S-100, Taitec, Saitama, Japan) for 30 s, centrifuge at 3000 rpm for 10 min. Evaporate the supernatant on a Speed Vac Plus SC110A (Savant Instruments. Inc., Farmingdale, NY) and dissolve the residue in 200 μ L THF: 50 mM pH 2.0 KH_2PO_4 containing 5 mM heptanesulfonate 25:75. Centrifuge at 10000 rpm for 3 min, inject 20 or 50 μ L aliquot.

HPLC VARIABLES

Guard column: 15 × 3.5 5 µm TSK-GEL ODS-80TS (Toyo soda, Tokyo)

Column: 150 × 4.6 5 µm TSK-GEL ODS-80TS (Toyo soda, Tokyo)

Mobile phase: THF:50mM pH 4.0 KH₂PO₄ containing 5mM heptanesulfonate 25:75

Flow rate: 0.8

Injection volume: 20 or 50

Detector: F ex 370 em 430

CHROMATOGRAM

Internal standard: camptothecin

KEY WORDS

plasma; mouse; pharmacokinetics

REFERENCE

Ohdo,S.; Makinosumi,T.; Ishizaki,T.; Yukawa,E.; Higuchi,S.; Nakano,S.; Ogawa,N. Cell cycle-dependent chronotoxicity of irinotecan hydrochloride in mice, *J.Pharmacol.Exp.Ther.*, **1997**, 283, 1383–1388.

SAMPLE

Matrix: blood

Sample preparation: Determination of lactone form. 1 mL Plasma + 100 µL 25 ng/mL IS in MeOH:10 mM HCl 40:60 + 800 mg solid NaCl, extract with 7.5 mL MeCN:n-butyl chloride 20:80 for 5 min, centrifuge at 4000 g for 2 min. Rotate quickly by hand to break the gels, centrifuge at 4000 g for 5 min. Mix the organic layer with 50 µL DMSO, dry under a gentle stream of nitrogen at 60° to approximately 50 µL. Reconstitute the residue in 100 µL MeOH and 100 µL perchloric acid:water 1:500, vortex for 5 s, inject a 100 µL aliquot. Total determination. Mix 250 µL plasma with 500 µL cold (-20°) MeOH:perchloric acid:water 20:1:20, centrifuge at 24000 g for 5 min, inject a 100 µL aliquot of the supernatant.

HPLC VARIABLES

Guard column: 4 × 4 5 µm LiChrospher 100 RP-18

Column: 100 × 4.6 5 µm Hypersil ODS

Mobile phase: MeOH:buffer 40:60 adjusted to pH 5.5 with HCl (A, lactone form determination) or MeOH:buffer 35:65 adjusted to pH 5.5 with HCl (B, total determination). (Buffer was 100 mM ammonium acetate containing 10 mM tetrabutylammonium sulfate.)

Column temperature: 50

Flow rate: 1

Injection volume: 100

Detector: F ex 355 em 515

CHROMATOGRAM

Retention time: 4.9 (A), 8.3 (B)

Internal standard: camptothecin (6.5, A)

Limit of quantitation: 200 pg/mL (lactone form), 2.0 ng/mL (total)

OTHER SUBSTANCES

Extracted: active metabolite

Noninterfering: acetaminophen, alizapride, codeine, dexamethasone, domperidone, metoclopramide, morphine, ranitidine

KEY WORDS

plasma; pharmacokinetics

REFERENCE

de Bruijn,P.; Verweij,J.; Loos,W.J.; Nooter,K.; Stoter,G.; Sparreboom,A. Determination of irinotecan (CPT-11) and its active metabolite SN-38 in human plasma by reversed-phase high-performance liquid chromatography with fluorescence detection, *J.Chromatogr.B*, **1997**, 698, 277–285.

SAMPLE

Matrix: blood

Sample preparation: Add 1 mL MeCN:MeOH 50:50 to 500 μ L plasma at -20° , vortex for 10 s. Centrifuge at 10000 g at 4° for 3 min, mix 70 μ L supernatant with 70 μ L mobile phase at 0° , vortex, inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: 10×3 C18 (Chrompack, Netherlands)

Column: 150×4.6 3.5 μ m Zorbax SB C18

Mobile phase: MeCN:100mM pH 6.4 ammonium acetate:triethylamine 15.6:80:0.1 containing 5 mM tetrabutylammonium phosphate

Flow rate: 1.5

Injection volume: 20

Detector: F ex 375 em 460

CHROMATOGRAM

Retention time: 5.7 (carboxylate form), 11.1 (lactone form)

Limit of quantitation: 1 ng/mL

OTHER SUBSTANCES

Extracted: metabolites (F ex 385 em 525)

Noninterfering: lurtotecan, topotecan

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Herben,V.M.M.; Mazee,D.; van Gortel-van Zomeren,D.M.; Zeedijk,S.; Schellens,J.H.M.; ten Bokkel Huinink,W.W.; Beijnen,J.H. Sensitive determination of the carboxylate and lactone forms of the novel antitumor drug irinotecan and its active metabolite in plasma by HPLC, *J.Liq.Chromatogr.Rel.Technol.*, **1998**, *21*, 1541-1558.

SAMPLE

Matrix: blood

Sample preparation: Condition a 100 mg Bond Elut C18 SPE cartridge with 1 mL MeOH and 1 mL water at 0.5 mL/s. 100 μ L Plasma + 50 μ L 1 μ g/mL camptothecin in 10 mM HCl + 850 μ L 10 mM HCl, mix, add to the SPE cartridge at 2.4 mL/min, wash with 1.5 mL water at 6 mL/min, wash with 1 mL 10 mM HCl at 6 mL/min, elute with 1.5 mL acidic MeOH at 2.4 mL/min. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute the residue in 250 μ L mobile phase, inject an aliquot. (Prepare IS solution by sonicating 1 mg camptothecin in 25 mL MeOH at room temperature for 20 min, dilute to 10 μ g/mL with water, dilute to 1 μ g/mL with 10 mM HCl. Prepare acidic MeOH by mixing 100 μ L concentrated HCl with 100 mL MeOH.)

HPLC VARIABLES

Guard column: 22×3.5 10 μ m Nucleosil C18

Column: 300×3.9 10 μ m Nucleosil octadecylsilane

Mobile phase: MeCN:100 mM KH_2PO_4 34:66 containing 3 mM sodium heptanesulfonate, adjusted to pH 4 with 1 M HCl

Flow rate: 1

Detector: F ex 380 em 500

CHROMATOGRAM

Retention time: 5.5

Internal standard: camptothecin (9)

Limit of detection: 1 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; pharmacokinetics; SPE

REFERENCE

Barilero,I.; Gandia,D.; Armand,J.-P.; Mathieu-Boué,A.; Ré,M.; Chabot,G.G. Simultaneous determination of the camptothecin analogue CPT-11 and its active metabolite SN-38 by high-performance liquid chromatography: application to plasma pharmacokinetic studies in cancer patients, *J.Chromatogr.*, **1992**, 575, 275–280.

SAMPLE

Matrix: blood

Sample preparation: Evaporate 50 μL 1 $\mu\text{g/mL}$ camptothecin in acetone into the bottom of a tube under a stream of nitrogen, add 50 μL plasma, add 100 μL ice-cold MeCN:MeOH 50:50, vortex for 5 s, centrifuge at 8000 g briefly. Remove a 100 μL aliquot of the supernatant and add it to 70 μL 75 mM pH 6.4 ammonium acetate buffer, vortex briefly, inject a 5-20 μL aliquot.

HPLC VARIABLES

Guard column: Guard-Pak Nova-Pak C18

Column: 100 \times 5 4 μm Nova-Pak Radial-Pak C18

Mobile phase: MeCN:75 mM pH 6.4 ammonium acetate buffer 22:78 containing 5 mM tetra-butylammonium phosphate (PIC A)

Flow rate: 1.5

Injection volume: 5-20

Detector: F ex 355 em 515

CHROMATOGRAM

Retention time: 4.2 (carboxylate form), 8.2 (lactone form)

Internal standard: camptothecin (6.5 (carboxylate form), 10.5 (lactone form))

Limit of quantitation: 10 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; pharmacokinetics; carboxylate form is the inactive form

REFERENCE

Rivory,L.P.; Robert,J. Reversed-phase high-performance liquid chromatographic method for the simultaneous quantitation of the carboxylate and lactone forms of the camptothecin derivative irinotecan, CPT-11 and its metabolite SN-38 in plasma, *J.Chromatogr.B*, **1994**, 661, 133–141.

SAMPLE

Matrix: blood

Sample preparation: 200 μL Plasma + 50 μL 1.25 $\mu\text{g/mL}$ camptothecin in MeOH:0.5 M HCl 97.5:2.5 + 750 μL MeOH, vortex for 10 s, centrifuge at 20° at 10000 g for 5 min. Remove the supernatant and evaporate it to dryness under reduced pressure, reconstitute with 400 μL mobile phase adjusted to pH 2.0, vortex for 10 s, centrifuge at 20° at 10000 g for 5 min, inject a 100 μL aliquot of the supernatant.

HPLC VARIABLES

Guard column: 15 \times 3.2 TSK guardgel ODS-120T

Column: 150 \times 4.6 5 μm TSK gel ODS-80Ts

Mobile phase: MeCN:50 mM Na_2HPO_4 28:72 containing 5 mM sodium 1-heptanesulfonate, adjusted to pH 3.0 with orthophosphoric acid (Prepared by dissolving 17.9 g $\text{Na}_2\text{HPO}_4 \cdot 12\text{H}_2\text{O}$ and 1.1 g sodium 1-heptanesulfonate hydrate in 1 L water and adding 388 mL MeCN, adjust pH to 3.0 with orthophosphoric acid.)

Column temperature: 30

Flow rate: 1

Injection volume: 100

Detector: F ex 380 em 556 or ex 370 em 430

CHROMATOGRAM

Retention time: 5.4

Internal standard: camptothecin (8.8)

Limit of quantitation: 30 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Sumiyoshi,H.; Fujiwara,Y.; Ohune,T.; Yamaoka,N.; Tamura,K.; Yamakido,M. High-performance liquid chromatographic determination of irinotecan (CPT-11) and its active metabolite (SN-38) in human plasma, *J.Chromatogr.B*, **1995**, 670, 309–316.

SAMPLE

Matrix: blood, feces, urine

Sample preparation: Plasma, urine. Thaw frozen plasma or urine sample in a waterbath, vortex, add 250 μ L plasma or 250 μ L plasma:diluted urine 50:50 to 500 μ L MeOH:5% (w/v) aqueous perchloric acid 50:50. Mix for 5 min, centrifuge at 24000 g for 5 min, dilute the upper aqueous layer from plasma and urine extracts 2-fold and 10-fold, respectively, with mobile phase, inject a 100–200 μ L aliquot. Feces. Homogenize feces with 5 volumes of 5% (w/v) perchloric acid in water using five 1 min bursts from an Ystral X1020 tissue homogenizer at 20500 r.p.m. Centrifuge at 24000 g for 10 min, dilute with one volume of drug-free human plasma, and process further as described for the urine sample.

HPLC VARIABLES

Guard column: 4 \times 4 LiChroCart endcapped RP-18 Merck

Column: 100 \times 4.6 5 μ m Hypersil ODS

Mobile phase: MeOH:100 mM ammonium acetate containing 10 mM tetrabutylammonium sulphate 30:70, pH adjusted to 5.3 with HCl

Column temperature: 50

Flow rate: 1.0

Injection volume: 100–200

Detector: F ex 355 em 515

CHROMATOGRAM

Retention time: 16.1

Limit of quantitation: 10 ng/mL (plasma), 100 ng/mL (feces, urine)

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

pharmacokinetics; plasma

REFERENCE

Sparreboom,A.; de Bruijn,P.; de Jonge,M.J.A.; Loos,W.J.; Stoter,G.; Verweij,J.; Nooter,K. Liquid chromatographic determination of irinotecan and three major metabolites in human plasma, urine and feces, *J.Chromatogr.B*, **1998**, 712, 225–235.

SAMPLE

Matrix: blood, tissue

Sample preparation: Condition an Analytichem C18 SPE cartridge with 1.5 mL MeOH and 1.5 mL water. Homogenize (PTFE homogenizer) tissue with four volumes ice-cold 150 mM KCl, centrifuge at 2° at 9000 g for 20 min. Dilute plasma, serum, or tissue homogenate supernatant 10-fold with 100 mM HCl, add to the SPE cartridge, elute the contents of the SPE cartridge on to the column with the mobile phase.

HPLC VARIABLES

Guard column: RP-18

Column: 250 \times 4 LiChrosorb RP-18

Mobile phase: MeCN:EtOH:0.8% ammonium carbonate 50:25:25

Column temperature: 50

Flow rate: 1

Detector: F ex 373 em 428

CHROMATOGRAM

Retention time: 6.6

Limit of detection: 1 ng

KEY WORDS

mouse; plasma; serum; liver; epithelium; pharmacokinetics; SPE

REFERENCE

Kaneda,N.; Nagata,H.; Furuta,T.; Yokokura,T. Metabolism and pharmacokinetics of the camptothecin analogue CPT-11 in the mouse, *Cancer Res.*, **1990**, 50, 1715-1720.


Isepamicin

Molecular formula: $C_{22}H_{43}N_5O_{12}$

Molecular weight: 569.61

CAS Registry No.: 58152-03-7, 67814-76-0 (sulfate)

Merck Index: 5121



SAMPLE

Matrix: blood

Sample preparation: Mix 1 mL plasma with 0.5-5.0 μ g IS and 2 mL EtOH. Add 7 mL dichloromethane and 1 mL water, mix, centrifuge. Inject an aliquot of the aqueous supernatant onto column A, elute to waste with mobile phase A, after 4 min elute the contents of column A onto column B with mobile phase B, elute column B with mobile phase B. Mix the effluent from column B with o-phthalaldehyde at 0.2 mL/min and monitor.

HPLC VARIABLES

Column: A 3.9×4.0 10 μ m Guard Pak Cyano (Waters); B 150×4.6 5 μ m Shandon Hypersil C18

Mobile phase: A 17 mM acetic acid containing 10 mM hexanesulfonic acid; B 17 mM acetic acid containing 10 mM hexanesulfonic acid, 100 mM sodium acetate, and 3.53 M MeOH

Detector: F ex 338 em 450 (cut-off filter) following post-column reaction. The column effluent mixed with o-phthalaldehyde pumped at 0.2 mL/min and the mixture flowed to the detector.

CHROMATOGRAM

Retention time: 7.4

Internal standard: dibekacin (9.5)

Limit of quantitation: 100 ng/mL

OTHER SUBSTANCES

Simultaneous: aspirin, caffeine, chlorpheniramine, gentamicin, neomycin, netilmicin, sisomicin

KEY WORDS

plasma; pharmacokinetics; comparison with RIA and microbiological assay; post-column reaction; column-switching

REFERENCE

Lin,C.-.; Veals,J.; Korduba,C.; Hilbert,M.J.; Nomeir,A. Analysis of isepamicin in human plasma by radioimmunoassay, microbiologic assay, and high-performance liquid chromatography, *Ther.Drug Monit.*, **1997**, 19, 675-681.

SAMPLE

Matrix: blood

Sample preparation: Dry blood on gauze, soak gauze in 500 μL 500 mM Na_2HPO_4 at 35° for 30 min, inject a 50 μL aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 5 μm $\mu\text{Bondapak C18}$

Mobile phase: 16 mM Sodium sulfate containing 5 mM 1-heptanesulfonic acid (PIC B-7)

Flow rate: 0.6

Injection volume: 50

Detector: F ex 360 em 440 following post-column reaction. The column effluent was mixed with reagent pumped at 0.3 mL/min and flowed through a 5 m \times 0.25 mm i.d. coil of PTFE tubing to the detector. (Prepare reagent by dissolving 300 mg o-phthalaldehyde in 500 mL MeOH, add 1.25 mL β -mercaptoethanol, add 500 mL 400 mM pH 10.4 potassium borate buffer.)

CHROMATOGRAM

Retention time: 12

Limit of detection: 100 ng/mL

KEY WORDS

post-column reaction; dried blood

REFERENCE

Shoshihara,M.; Kase,K.; Yoshizawa,E.; Takao,M.; Fujimoto,T. Column liquid chromatographic determination of isepamicin in nasal cavity using gauze, *J.Chromatogr.*, **1990**, 529, 473–478.

SAMPLE

Matrix: blood

Sample preparation: 50 μL Plasma + 20 μL 10 mg/mL gentamicin C1a in water + 50 μL buffer, vortex for 15 s, add 200 μL MeCN, vortex for 20 s, centrifuge at 2000 g for 5 min. Filter (Millex-HV4) the supernatant. Heat 200 μL filtrate and 20 μL 250 mg/mL 1-fluoro-2,4-dinitrobenzene in MeCN at 80° for 1 h, cool, inject a 50 μL aliquot. (Buffer was 3.81 g disodium tetraborate decahydrate in water, adjust pH to 10 with NaOH, make up to 100 mL with water.)

HPLC VARIABLES

Guard column: 25 \times 4 10 μm LiChroCART RP 18

Column: 250 \times 4 5 μm LiChrosorb RP 18

Mobile phase: MeCN:water 70:30 containing 1 mL/L acetic acid

Flow rate: 2

Injection volume: 50

Detector: UV 365

CHROMATOGRAM

Retention time: 13.0

Internal standard: gentamicin C1a (10.0)

Limit of quantitation: 500 ng/mL

OTHER SUBSTANCES

Noninterfering: ampicillin, aspirin, captopril, cefazolin, cefotaxime, ceftazidime, ceftriaxone, cephalosporins, chlorpromazine, diazepam, heparin, propranolol, sulfamethoxazole, sulpiride, trimethoprim, verapamil

KEY WORDS

plasma; guinea pig; human; derivatization

REFERENCE

Dionisotti,S.; Bamonte,F.; Scaglione,F.; Ongini,E. Simple measurement of isepamicin, a new aminoglycoside antibiotic, in guinea pig and human plasma, using high-performance liquid chromatography with ultraviolet detection, *Ther.Drug Monit.*, **1991**, 13, 73–78.

SAMPLE

Matrix: blood, dialysate, urine

Sample preparation: Plasma. Condition a 3 mL Baker cyanopropylsilane CN SPE cartridge with 2 mL MeOH, 2 mL water, and 2 mL buffer. 1 mL Plasma + 100 μ L 100 μ g/mL dibekacin in water, vortex for 15 s, add 1 mL buffer, vortex for 15 s, centrifuge at 3100 g at 4° for 7 min, add to SPE cartridge, wash with 500 μ L water, wash with 250 μ L mobile phase, elute to dryness. Elute with 250 μ L mobile phase, inject an aliquot of the eluate. Urine, dialysate. Dilute 1:100 with water, add 100 μ L 100 μ g/mL dibekacin per 1 mL of sample, mix well, inject a 100 μ L aliquot. (Buffer was 0.94 g sodium hexanesulfonate in 300 mL water, add 500 μ L glacial acetic acid, dilute to 500 mL with water.)

HPLC VARIABLES

Guard column: 10 \times 4.6 5 μ m Hypersil C18

Column: 150 \times 4.6 5 μ m Hypersil C18

Mobile phase: MeOH:buffer 10:90 (Buffer was 3.76 g sodium hexanesulfonate + 28.4 g sodium sulfate in 2 L water, acidify to pH 3.4 with 2 mL glacial acetic acid.)

Column temperature: 25

Flow rate: 1.1

Injection volume: 100

Detector: F ex 338 em 418 (bandpass filter) following post-column reaction. The column effluent mixed with the reagent pumped at 0.4 mL/min and the mixture flowed through a 3 m \times 0.05 mm i.d. knitted PTFE reaction coil at 25° to the detector (Derivatizing reagent was 0.4 g o-phthalaldehyde in 3 mL MeOH added to 390 mL buffer, add 2 mL β -mercaptoethanol, make up to 500 mL with water, store at 4°. Buffer was 1 M pH 10.4 borate from equal volumes of 1 M KOH and boric acid.)

CHROMATOGRAM

Retention time: 6.7

Internal standard: dibekacin (17)

Limit of quantitation: 100 ng/mL

OTHER SUBSTANCES

Simultaneous: kanamycin, gentamicin, tobramycin, netilmicin

KEY WORDS

post-column reaction; SPE; plasma

REFERENCE

Maloney, J.A.; Awani, W.M. High-performance liquid chromatographic determination of isepamicin in plasma, urine and dialysate, *J. Chromatogr.*, **1990**, 526, 487–496.

SAMPLE

Matrix: blood, urine

Sample preparation: 1 mL Plasma or urine + 0.5–5 μ g dibekacin + 2 mL EtOH, mix, add 7 mL dichloromethane, add 1 mL water, centrifuge, inject an aliquot of the aqueous supernatant on to column A and elute to waste with mobile phase A, after 4 min elute the contents of column A on to column B with mobile phase B, monitor the effluent from column B.

HPLC VARIABLES

Column: A 3.9 \times 4 10 μ m Guard Pak with Cyano insert; B 150 \times 4.6 5 μ m Hypersil C18

Mobile phase: A 17 mM Acetic acid containing 10 mM hexanesulfonate; B MeOH:buffer 14.3:85.7 (Mobile phase contained 100 mM sodium acetate, 17 mM acetic acid, 10 mM hexanesulfonate, and 3.53 M MeOH.)

Detector: F ex 338 (band-pass filter) em 418–700 and 450 cut-off filters following post-column reaction with o-phthalaldehyde derivatizing reagent pumped at 0.2 mL/min.

CHROMATOGRAM

Internal standard: dibekacin

Limit of quantitation: 100 ng/mL

KEY WORDS

plasma; pharmacokinetics; post-column reaction; column-switching

REFERENCE

Lin, C.-C.; Radwanski, E.; Korduba, C.; Cayen, M.; Affrime, M. Pharmacokinetics of intravenously administered isepamicin in men, *Antimicrob. Agents Chemother.*, **1995**, 39, 2774–2778.

Isocarboxazid


Molecular formula: $C_{12}H_{13}N_3O_2$

Molecular weight: 231.25

CAS Registry No.: 59-63-2

Merck Index: 5172

Lednicer No.: 1 233



SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 200 μ L 1.5 M NaOH + 10 μ L 200 μ g/mL IS in 100 mM HCl, vortex, add 5 mL hexane:ethyl acetate 80:20, shake mechanically at low speed for 10 min, centrifuge at 10° at 1100 g for 10 min. Remove the organic layer and add it to 500 μ L 2 M HCl, shake mechanically for 5 min, centrifuge at 10° at 1100 g for 5 min. Remove the aqueous phase and add it to 200 μ L saturated K_2HPO_4 , vortex, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m LC-18 (Supelco)

Mobile phase: MeOH:5 mM octanesulfonic acid 50:50

Flow rate: 1

Injection volume: 50

Detector: UV 230

CHROMATOGRAM

Retention time: 6.4

Internal standard: 5-methyl-3-isoxazolecarboxylic acid 2-(2-propyl-1-phenyl)hydrazide (Ro 5-1226) (15.5)

Limit of quantitation: 100 ng/mL

KEY WORDS

plasma; dog; pharmacokinetics

REFERENCE

Powell, M.L.; Town, C.; Henderson, L.; Buck, C. Determination of marplan in human plasma using high-performance liquid chromatography, *J. Chromatogr.*, **1990**, 529, 237–244.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitrityline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenopropfen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminos-tilbene, imipramine, indomethacin, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazin-dol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyldopa, methyldopamine, methylphenidate, methylprednisolone, methyltestosterone, methyprylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, nor-epinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphen-butazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, per-santine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenyl-butazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primi-done, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopola-mine, scopoletin, secobarbital, strychnine, sulfacetamide, sufadiazine, sulfadimethoxine, sul-faethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sul-fasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tol-metin, tranlycpromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleppamine, triprolidine, tropacocaine, tyramine, verapa-mil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.

Isoetharine


Molecular formula: C₁₃H₂₁NO₃

Molecular weight: 239.31

CAS Registry No.: 530-08-5, 2576-92-3 (HCl), 7279-75-6 (mesylate)

Merck Index: 5185

Lednicer No.: 29



SAMPLE

Matrix: blood

Sample preparation: 25-150 μ L Plasma + 100 μ L 50 ng/mL colterol mesylate in 25 mM sulfuric acid + 2 mL 2% boric acid + 2 g ammonium sulfate + 10 mL 0.5% di(2-ethylhexyl)phosphoric acid in benzene (Caution! Benzene is a carcinogen!), shake, centrifuge. Remove the organic phase and add it to 130 μ L 25 mM sulfuric acid, shake, centrifuge, inject a 100 μ L aliquot of the aqueous phase.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Ultrasphere ODS

Mobile phase: MeOH:buffer 12.5:87.5 (Buffer was 100 mM sodium sulfate adjusted to pH 2.8 phosphoric acid then adjusted to pH 3.0 with NaOH.)

Flow rate: 1.2

Injection volume: 100

Detector: E, Bioanalytical Systems LC-4, TL-5 glassy carbon electrode +0.60 V, Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 9.3

Internal standard: colterol mesylate (7.5)

Limit of quantitation: 0.5 ng/mL

KEY WORDS

plasma; rat; pharmacokinetics

REFERENCE

Park,G.B.; Koss,R.F.; Utter,J.; Mayes,B.A.; Edelson,J. Determination of isoetharine in plasma by reversed-phase chromatography with amperometric detection, *J.Pharm.Sci.*, **1982**, 71, 932-934.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 Chirex 3020 (Phenomenex)

Mobile phase: Hexane:1,2-dichloroethane:EtOH/trifluoroacetic acid 55:35:10 (EtOH/trifluoroacetic acid was premixed 20:1.)

Flow rate: 0.7-1

Injection volume: 20

Detector: UV 282

KEY WORDS

chiral; α = 1.21 for enantiomers

REFERENCE

Cleveland,T. Pirkle-concept chiral stationary phases for the HPLC separation of pharmaceutical racemates, *J.Liq.Chromatogr.*, **1995**, 18, 649-671.


Isofezolac

Molecular formula: $C_{23}H_{19}N_2O_2$

Molecular weight: 354.41

CAS Registry No.: 50270-33-2

Merck Index: 5188



SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 1 mL Plasma + 1 mL 100 mM pH 4.4 citrate buffer + 100 μ L 40 μ g/mL IS in MeOH + 15 mL diethyl ether, agitate, centrifuge at 3000 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 0.5-3 mL mobile phase, inject a 20 μ L aliquot. Urine. 100 μ L Urine + 1 mL 100 mM pH 4.4 citrate buffer + 100 μ L 40 μ g/mL IS in MeOH + 100 μ L enzyme solution containing 100000 U/mL β -glucuronidase and 1000000 U/mL arylsulfatase (I.B.F.), heat at 37° for 16 h, add 15 mL diethyl ether, agitate, centrifuge at 3000 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 1-2 mL mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m LiChrosorb RP 8

Mobile phase: MeCN:water:200 mM pH 3 phosphate buffer 65:15:20

Flow rate: 1.5

Injection volume: 20

Detector: F ex 273 em 335 or UV 265

CHROMATOGRAM

Retention time: 2.7

Internal standard: 1-phenyl-3,4-di-p-chlorophenylpyrazole-5-acetic acid (4.5)

Limit of detection: 10 ng/mL (F)

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Bannier,A.; Brazier,J.L. Determination of isofezolac in biological fluids by reversed-phase liquid column chromatography, *J.Chromatogr.*, **1980**, 182, 369-377.


Isoflupredone

Molecular formula: $C_{21}H_{27}FO_5$

Molecular weight: 378.44

CAS Registry No.: 338-95-4, 338-98-7 (acetate)

Merck Index: 5190



SAMPLE

Matrix: bulk, formulations

Sample preparation: Cream. Weigh out cream containing 1 mg diflorasone diacetate, add 30 mL 40 μ g/mL isoflupredone acetate in water-saturated chloroform, shake for 30 min, centrifuge at 2000 rpm for 15 min, inject a 10 μ L aliquot of the lower chloroform layer. Ointment. Weigh out ointment containing 0.5 mg diflorasone diacetate, add 15 mL 40 μ g/mL isoflupredone acetate in water-saturated chloroform, shake for 30 min, centrifuge at 2000 rpm for 15 min, inject a 10 μ L aliquot of the lower chloroform layer. Bulk. Dissolve 1.5 mg bulk drug in 50 mL 40 μ g/mL isoflupredone acetate in water-saturated chloroform, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 100 × 4.6 3 µm silica gel (Perkin-Elmer part 0258-1500)

Mobile phase: Butyl chloride:dichloromethane:THF:acetic acid 70:25:2:3 (Butyl chloride and dichloromethane were saturated with water.)

Flow rate: 2.5

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 24 (isoflupredone acetate)

Internal standard: isoflupredone acetate

OTHER SUBSTANCES

Simultaneous: related compounds, diflorasone diacetate

KEY WORDS

cream; ointment; normal phase; isoflupredone acetate is IS

REFERENCE

Shaw,M.C.; Vanderwielen,A.J. Liquid chromatographic assay for diflorasone diacetate in cream and ointment formulations, *J.Pharm.Sci.*, **1984**, 73, 1606–1608.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 5 µm SI-100 (Brownlee)

Mobile phase: Butyl chloride:THF:MeOH:glacial acetic acid 88:2.5:2.5:2.5 (Butyl chloride was 50% water saturated.)

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 16 (isoflupredone), 10 (isoflupredone acetate)

OTHER SUBSTANCES

Simultaneous: bromoprednisolone acetate, prednisolone, fluoroprednisone acetate


KEY WORDS

normal phase

REFERENCE

Kane,M.P.; Tsuji,K. Radiolytic degradation scheme for ⁶⁰Co-irradiated corticosteroids, *J.Pharm.Sci.*, **1983**, 72, 30–35.

Isoflurane



Molecular formula: C₃H₂ClF₅O

Molecular weight: 184.49

CAS Registry No.: 26675-46-7

Merck Index: 5191

SAMPLE

Matrix: solutions

Sample preparation: Mix 50 µL phosphate buffer containing isoflurane and 50 µL 0.05 mM toluene in MeOH, inject a 20 µL aliquot.

HPLC VARIABLES

Guard column: RCSS Guard-Pak μ Bondapak C18 precolumn cartridge

Column: 100 \times 8 4 μ m Nova-Pak C18 Radial Compression Module

Mobile phase: MeOH:water 50:50

Flow rate: 3.5

Injection volume: 20

Detector: UV 203

CHROMATOGRAM

Retention time: 5

Internal standard: toluene (12)

Limit of detection: 0.2 mM

OTHER SUBSTANCES

Simultaneous: enflurane, halothane

KEY WORDS

buffer

REFERENCE

Janicki,P.K.; Erskine,W.A.R.; James,M.F.M. High-performance liquid chromatographic method for the direct determination of the volatile anaesthetics halothane, isoflurane and enflurane in water and in physiological buffer solutions, *J.Chromatogr.*, **1990**, 518, 250–253.

Isoniazid


Molecular formula: C₆H₇N₃O

Molecular weight: 137.14

CAS Registry No.: 54-85-3

Merck Index: 5203

Lednicer No.: 1 254



SAMPLE

Matrix: blood

Sample preparation: 2 mL Plasma + 100 μ L 30 μ g/mL IS in water, mix, add 2 g ammonium sulfate, add 40 mL water saturated n-butanol:chloroform 30:70, shake for 10 min, centrifuge at 500 g for 10 min. Transfer the organic layer into a tube, add 1 mL 1 M sulfuric acid, shake for 10 min, centrifuge at 500 g for 10 min, inject a 250 μ L aliquot of the upper aqueous layer. (Caution! Chloroform is a carcinogen !)

HPLC VARIABLES

Guard column: 10 μ m μ Bondapak C18

Column: 115 \times 8 5 μ m μ Bondapak C18 radial compression

Mobile phase: EtOH:1 mM dioctyl sulfosuccinate 45:55, adjusted to pH 2.50

Flow rate: 4

Injection volume: 250

Detector: UV 254

CHROMATOGRAM

Retention time: 10.7

Internal standard: 1-benzoyl-2-isonicotinoylhydrazine (8.9)

Limit of detection: 100 ng/mL

OTHER SUBSTANCES

Extracted: acetylisoniazid

KEY WORDS

plasma

REFERENCE

Holdiness, M.R. High pressure liquid chromatographic determination of isoniazid and acetylisoniazid in human plasma, *J.Liq.Chromatogr.*, **1982**, 5, 707–714.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Plasma + 50 μ L 50 μ g/mL phenelzine + 400 μ L 10% acetic acid + 7 mL diethyl ether:dichloromethane 2:1, shake, centrifuge at 2059 g for 10 min. Remove the aqueous layer and add it to 600 μ L 10% acetic acid, add 300 μ L 0.1% salicaldehyde in EtOH, heat at 60° for 30 min, cool to room temperature, add 1 mL 1 M K_2PO_4 (sic), extract with 7 mL diethyl ether, centrifuge at 2059 g for 10 min, repeat extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 40°, reconstitute the residue in 200 μ L buffer, inject a 20 μ L aliquot. (Buffer was 5 mM heptanesulfonic acid in MeCN:water:triethylamine 70:30:0.4.)

HPLC VARIABLES

Guard column: 50 \times 4.6 30 μ m C8

Column: 250 \times 4.6 Spherisorb S5 ODS2 C18

Mobile phase: Gradient. MeCN:buffer:water 0:75:25 for 5 min, 15:85:0 for 12 min (step gradient). (Buffer was 5 mM heptanesulfonic acid in MeCN:water:triethylamine 70:30:0.4.)

Flow rate: 1

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 4.68

Internal standard: phenelzine (11.09)

Limit of detection: 250 ng/mL

OTHER SUBSTANCES

Extracted: hydrazine, monoacetylhydrazine

KEY WORDS

derivatization

REFERENCE

Walubo, A.; Smith, P.; Folb, P.I. Comprehensive assay for pyrazinamide, rifampicin and isoniazid with its hydrazine metabolites in human plasma by column liquid chromatography, *J.Chromatogr.B*, **1994**, 658, 391–396.

SAMPLE

Matrix: blood

Sample preparation: 250 μ L Plasma + 150 μ L 10% zinc sulfate in water, mix, centrifuge at 1000 g for 1 min. Remove a 250 μ L aliquot of the supernatant and add it to 100 μ L MeOH, mix, centrifuge, inject a 5 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m LC-CN (Supelco)

Mobile phase: Isopropanol:water 5:95 containing 5 g/L ammonium formate

Flow rate: 1

Injection volume: 5

Detector: E, ESA Model 5100A, Model 5011 analytical cell with first detector +0.6 V and second (monitored) detector +0.8 V, Model 5020 guard cell +1.0 V between pump and injector

CHROMATOGRAM

Retention time: 4.4

Limit of detection: 100 ng/mL

KEY WORDS

rat; plasma

REFERENCE

Hansen, E.B., Jr.; Dooley, K.L.; Thompson, H.C., Jr. High-performance liquid chromatographic analysis of the antituberculosis drugs aconiazide and isoniazid, *J. Chromatogr. B*, **1995**, 670, 259–266.

SAMPLE

Matrix: blood

Sample preparation: 250 μ L Plasma + 50 μ L 1.5% p-hydroxybenzaldehyde in MeOH + 40 μ L 20% trichloroacetic acid in water, vortex thoroughly, centrifuge at 8000 g for 10 min, let stand on ice for 30 min, centrifuge at 8000 g for 10 min, inject a 30 μ L aliquot of the supernatant.

HPLC VARIABLES

Guard column: 20 \times 2 30–40 μ m Perisorb RP18

Column: 100 \times 8 4 μ m Radial-Pak Novapak C18

Mobile phase: MeOH:water:20% tetraethylammonium hydroxide:70% perchloric acid 24:76:0.05:0.05, apparent pH 2.3

Flow rate: 2

Injection volume: 30

Detector: UV 350

CHROMATOGRAM

Retention time: 4.4

Limit of detection: <500 ng/mL

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Schall, R.; Müller, F.O.; Duursema, L.; Groenewoud, G.; Hundt, H.K.L.; Middle, M.V.; Mogilnicka, E.M.; Swart, K.J. Relative bioavailability of rifampicin, isoniazid and ethambutol from a combination tablet vs. concomitant administration of a capsule containing rifampicin and a tablet containing isoniazid and ethambutol, *Arzneimittelforschung*, **1995**, 45, 1236–1239.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Serum + 100 μ L 30% trichloroacetic acid, mix, centrifuge at 2000 g for 5 min. Remove a 100 μ L aliquot of the supernatant and add it to 20 μ L 0.1% trans-cinnamaldehyde in MeOH, let stand for 10 min, add 20 μ L 1 M KOH, inject an aliquot.

HPLC VARIABLES

Column: 125 \times 3.9 4 μ m Nova-pak C18

Mobile phase: MeCN:water:triethylamine:acetic acid 40:60:0.2:0.1, pH 5 \pm 1

Flow rate: 1.3

Detector: UV 340

CHROMATOGRAM

Retention time: 1.95

Limit of detection: 20 ng/mL

OTHER SUBSTANCES

Noninterfering: aceprometazine, adrafinil, allopurinol, alprostadiol, altretamine, atenolol (tenormine), baclofen, bendroflumethiazide, benserazide, betamethasone, bisoprolol, bromocriptine, caffeine, captopril, chlorpromazine, clomipramine, clonazepam, cortisone, cyamemazine, difebarbamate, dothiepin (dosulepin), ethosuximide, fenspiride, flumazenil, fluoxetine, fluvoxamine, halofantrine, hydrochlorothiazide, hydroxyzine, ibuprofen, imipramine, levamisole, levodopa, maprotiline, medifoxamine, metopimazine, midazolam, nafronyl (naftidrofuryl), naftazone, naproxen, nicergoline, nitrazepam, nordazepam, nortriptyline, penfluridol, phenobarbital, pimozone, pipamperone, pipotizine, primidone, pyrazinamide, pyridoxine, quinine, rifampin, selegiline (deprenyl), streptomycin, tetrazepam, theophylline, thioproperazide, tiapride, triazolam, trihexyphenidyl, trimeprazine (alimemazine), trimipramine, tropatepine, vigabatrin, zopiclone

KEY WORDSserum; derivatization

REFERENCE

Sadeg,N.; Pertat,N.; Dutertre,H.; Dumontet,M. Rapid, specific and sensitive method for isoniazid determination in serum, *J.Chromatogr.B*, **1996**, 675, 113–117.

SAMPLE**Matrix:** blood, CSF

Sample preparation: 200-500 μL Plasma or CSF + 50 μL 50 $\mu\text{g/mL}$ phenelzine sulfate + 100 μL 10% (?) aqueous acetic acid + 5 mL n-hexane, shake for 30 min, centrifuge at 1870 g for 10 min. Discard the organic layer. Add 300 μL 0.1% salicaldehyde in EtOH and 400 μL 10% aqueous acetic acid to the aqueous layer, heat at 60° for 30 min, cool, add 1 mL 1 M pH 6.5 K_2HPO_4 , shake for 10 s, add 5 mL diethyl ether, shake for 10 min, centrifuge at 1870 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 50 μL mobile phase, inject a 25 μL aliquot.

HPLC VARIABLES**Guard column:** 30 \times 4.6 30 μm C8 (Waters)**Column:** 300 \times 3.9 10 μm $\mu\text{Bondapak}$ C18**Mobile phase:** MeCN:water:triethylamine 70:30:0.4 containing 5 mM heptanesulfonic acid, pH adjusted to 6.0 with acetic acid**Flow rate:** 1**Injection volume:** 25**Detector:** UV 320

CHROMATOGRAM**Retention time:** 1.6**Internal standard:** phenelzine sulfate (3)**Limit of detection:** 200 ng/mL

OTHER SUBSTANCES**Extracted:** hydrazine**Noninterfering:** p-aminosalicylic acid, pyrazinamide, rifampin

KEY WORDSplasma; rabbit; derivatization

REFERENCE

Walubo,A.; Chan,K.; Wong,C.L. Simultaneous assay for isoniazid and hydrazine metabolite in plasma and cerebrospinal fluid in the rabbit, *J.Chromatogr.*, **1991**, 567, 261–266.

SAMPLE**Matrix:** blood, CSF

Sample preparation: 200 μL Serum, plasma, or CSF + 300 μL reagent. Flush column A to waste with 500 μL 500 mM ammonium sulfate, inject sample onto column A, flush column A to waste with 500 μL 500 mM ammonium sulfate, elute the contents of column A onto column B with mobile phase, monitor the effluent from column B. (Reagent was 8.05 M guanidine hydrochloride and 1.02 M ammonium sulfate in water.)

HPLC VARIABLES**Column:** A 30 \times 2.1 40 μm preparative grade C18 (Analytichem); B 250 \times 4.6 10 μm Partisil C8**Mobile phase:** Gradient. A was 50 mM pH 4.5 KH_2PO_4 . B was MeCN:isopropanol 80:20. A:B 90:10 for 1 min, to 30:70 over 15 min, maintain at 30:70 for 4 min.**Column temperature:** 50**Flow rate:** 1.5**Detector:** UV 280 for 5 min then UV 254

CHROMATOGRAM**Retention time:** 3.10

Internal standard: heptanophenone (19.2)

Limit of quantitation: 2500 ng/mL

OTHER SUBSTANCES

Extracted: acetazolamide, ampicillin, bromazepam, caffeine, carbamazepine, chloramphenicol, chlorothiazide, diazepam, droperidol, ethionamide, furosemide, methadone, penicillin G, phenobarbital, phenytoin, prazepam, propoxyphene, pyrazinamide, rifampin, trimeprazine, trimethoprim

KEY WORDS

plasma; serum; column-switching

REFERENCE

Seifart,H.I.; Kruger,P.B.; Parkin,D.P.; van Jaarsveld,P.P.; Donald,P.R. Therapeutic monitoring of antituberculosis drugs by direct in-line extraction on a high-performance liquid chromatography system, *J.Chromatogr.*, **1993**, 619, 285–290.

SAMPLE

Matrix: blood, CSF, urine

Sample preparation: Centrifuge 1.5 mL whole blood, CSF, or urine at 3000 g for 3 min, add 500 μ L of the supernatant to 500 μ L 10% trichloroacetic acid, centrifuge at 10000 g for 1 min. Remove a 200 μ L aliquot and add it to 20 μ L water, add 40 μ L 1% trans-cinnamaldehyde in MeOH, let stand at room temperature for 10 min, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 Partisil 5 C8

Mobile phase: Gradient. A was 50 mM KH_2PO_4 . B was MeCN:isopropanol 80:20. A:B 60:40 for 1 min, to 30:70 over 9 min, maintain at 30:70 for 4.5 min, re-equilibrate at initial conditions for 4 min.

Column temperature: 50

Flow rate: 1

Injection volume: 20

Detector: UV 340

CHROMATOGRAM

Retention time: 7

Limit of detection: 500 ng/mL

OTHER SUBSTANCES

Extracted: hydrazine

KEY WORDS

whole blood; derivatization

REFERENCE

Seifart,H.I.; Gent,W.L.; Parkin,D.P.; van Jaarsveld,P.P.; Donald,P.R. High-performance liquid chromatographic determination of isoniazid, acetylisoniazid and hydrazine in biological fluids, *J.Chromatogr.B*, **1995**, 674, 269–275.

SAMPLE

Matrix: blood, urine

Sample preparation: Condition a 3 mL CN-bonded maxi-clean cartridge (Alltech) with 5 mL MeOH and 2 mL 1% aqueous acetic acid. 2 mL Plasma or urine + 1 mL isopropanol:chloroform 50:50 + 1 μ g rifampin, shake for 30 s, add to the SPE cartridge with rinses, collect the eluate, inject an aliquot of the eluate.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Excalibar C18-CN (Alltech)

Mobile phase: MeOH:5 mM tetra-n-butylammonium hydroxide 80:20 adjusted to pH 3.0 with phosphoric acid

Flow rate: 1.5

Injection volume: 20

Detector: UV 265

CHROMATOGRAM

Retention time: 2.7

Internal standard: rifampin (4.3)

Limit of detection: 250 ng/mL (urine), 200 ng/mL (plasma)

OTHER SUBSTANCES

Extracted: pyrazinamide

KEY WORDS

plasma; SPE

REFERENCE

Gaitonde, C.D.; Pathak, P.V. Rapid liquid chromatographic method for the estimation of isoniazid and pyrazinamide in plasma and urine, *J.Chromatogr.*, **1990**, 532, 418–423.

SAMPLE

Matrix: milk

Sample preparation: Condition a 3 mL 500 mg C18 SPE cartridge (J.T. Baker) with 6 mL MeOH and 6 mL water. 80 mL milk + 20 mL 20% trichloroacetic acid in water, mix, let stand at room temperature for 5 min, centrifuge at 5500 g for 15 min. Filter (0.45 μ m) 75 mL of the liquid, add 1 mL 1% cinnamaldehyde in MeOH to the filtrate, mix for a few s, let stand at room temperature for 15 min, add to the SPE cartridge at 5 mL/min, dry for 10 min, wash with 3 mL mobile phase at 0.2 mL/min, wash with 5 mL n-hexane at 0.3 mL/min, dry for 10 min, elute with 6 mL MeOH. Evaporate the eluate under reduced pressure at 35–40°, reconstitute the residue with 200 μ L mobile phase, inject an aliquot.

HPLC VARIABLES

Guard column: 4 \times 4 5 μ m LiChrospher 100-CN

Column: 250 \times 4 5 μ m LiChrospher 100-CN

Mobile phase: MeOH:water 40:60 containing 0.41 g/L sodium acetate trihydrate and 10 mL/L glacial acetic acid

Flow rate: 1

Injection volume: 100

Detector: UV 330

CHROMATOGRAM

Retention time: 7.4

Limit of detection: 0.1 ng/mL

KEY WORDS

cow; SPE; derivatization

REFERENCE

Defilippi, A.; Piancone, G.; Costa Laia, R.; Balla, S.; Tibaldi, G.P. High-performance liquid chromatography with UV detection and diode-array UV confirmation of isonicotinic acid hydrazide in cattle milk, *J.Chromatogr.B*, **1994**, 656, 466–471.

SAMPLE

Matrix: milk

Sample preparation: Condition a 3 mL 500 mg 40 μ m phenyl SPE cartridge (J.T. Baker) with 6 mL MeOH and 6 mL water. 80 mL Milk + 20 mL 20% trichloroacetic acid in water, let stand at room temperature for 5 min, centrifuge at 5500 g for 15 min, filter (0.45 μ m cellulose acetate) a 75 mL aliquot. Add the filtrate to 1 mL 1% cinnamaldehyde in MeOH, mix for a few s, let stand at room temperature for 15 min. Add 70 mL of the sample to the SPE cartridge at 3 mL/min, wash with 3 mL MeOH:water 40:60 at 3 mL/min, dry under vacuum for 1 min, wash with 3 mL MeOH:water 46:54 at 3 mL/min, dry under vacuum for 1 min, wash with 3 mL n-hexane at 3 mL/min, dry under vacuum for 1 min, elute with 6 mL MeOH. Evaporate the eluate under vacuum, reconstitute in 233 μ L mobile phase, inject an aliquot.

HPLC VARIABLES

Guard column: 4 × 4.5 µm LiChrosphere 100 RP18

Column: 250 × 4.5 µm LiChrosphere 100 RP18

Mobile phase: MeCN:MeOH:buffer 41:13:46 (Buffer was 1% ammonium acetate adjusted to pH 5.5 with acetic acid.)

Flow rate: 1

Injection volume: 100

Detector: UV 330

CHROMATOGRAM

Retention time: 4.02

Limit of detection: 0.05 ng/mL

KEY WORDS

cow; SPE; derivatization

REFERENCE

Defilippi, A.; Piancone, G.; Costa Laia, R.; Tibaldi, G. P. An HPLC screening method for the detection of isonicotinic acid hydrazide in cattle milk, *Chromatographia*, **1995**, *40*, 170–174.

SAMPLE

Matrix: reaction mixtures

Sample preparation: Centrifuge, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Microsorb C8

Mobile phase: MeOH:0.4 g/L (NH₄)H₂PO₄ + 0.1% triethylamine (pH 10.0) 10:90

Flow rate: 1

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 5.3

Limit of detection: 10000 ng/mL

REFERENCE

Lunn, G.; Sansone, E. B. Reductive destruction of dacarbazine, procarbazine hydrochloride, isoniazid, and iproniazid, *Am. J. Hosp. Pharm.*, **1987**, *44*, 2519–2524.

SAMPLE

Matrix: solution

Sample preparation: Mix 1 mL pH 7.4 phosphate buffer or hepatocyte solution with 3 mL ethyl acetate:n-butanol 2:1, 300 µL 500 mM tetra-n-butylammonium hydroxide, and 100 µL 80 µg/mL IS in water. Shake the mixture for 20 min, centrifuge at 850 g for 15 min, remove 2.5 mL of the upper organic layer, mix with 1 mL 0.2% hydrobromic acid, shake for 20 min, centrifuge at 850 g for 15 min, discard the upper organic layer, inject a 5 µL aliquot of the aqueous solution.

HPLC VARIABLES

Column: 250 × 4.6 5 µm TSK-gel ODS-80Ts

Mobile phase: MeOH:67 mM KH₂PO₄ 4:96

Column temperature: 37

Flow rate: 0.8

Injection volume: 5

Detector: UV 265

CHROMATOGRAM

Retention time: 11.42

Internal standard: 6-methylnicotinic acid (8)

Limit of quantitation: 2 µg/mL

OTHER SUBSTANCES

Simultaneous: acetylisoniazid, isonicotinic acid

KEY WORDS

rat; hepatocytes

REFERENCE

Ono,Y.; Noda,A.; Zaima,Y.; Jitsufuchi,N.; Eto,S.; Noda,H. Determination of isonicotinic acid in the presence of isoniazid and acetylisoniazid. Studies on isonicotinic acid formation from isoniazid in isolated rat hepatocytes, *J.Chromatogr.B*, **1996**, 677, 339–343.

SAMPLE

Matrix: solutions

Sample preparation: Centrifuge and filter cell solutions (0.22 μm), inject an aliquot.

HPLC VARIABLES

Guard column: Guard-PAK C18 (Waters)

Column: 300 \times 3.9 5 μm μ Bondapak C18

Mobile phase: MeOH:50 mM pH 6.0 KH_2PO_4 3:97

Flow rate: 3

Detector: UV 254

CHROMATOGRAM

Retention time: 3.5

REFERENCE

Koga,H. High-performance liquid chromatography measurement of antimicrobial concentrations in polymorphonuclear leukocytes, *Antimicrob.Agents Chemother.*, **1987**, 31, 1904–1908.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amyllocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenpropofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiacol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-

stilbene, imipramine, indomethacin, isocarbostyryl, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephénytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methyprylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, nor-epinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233–242.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 12 µm Dynamax C18 (Rainin)

Mobile phase: pH 7.5 sodium phosphate buffer

Flow rate: 2

Detector: UV 266

CHROMATOGRAM

Retention time: 5.4

OTHER SUBSTANCES

Simultaneous: metabolites, acetylisoniazid

REFERENCE

Hickman,D.; Palamanda,J.R.; Unadkat,J.D.; Sim,E. Enzyme kinetic properties of human recombinant arylamine N-acetyltransferase 2 allotypic variants expressed in *Escherichia coli*, *Biochem.Pharmacol.*, **1995**, *50*, 697–703.


Isoproterenol

Molecular formula: C₁₁H₁₇NO₃

Molecular weight: 211.26

CAS Registry No.: 7683-59-2, 51-30-9 (HCl), 6700-39-6 (sulfate dihydrate), 299-95-6 (sulfate)

Merck Index: 5236



SAMPLE

Matrix: blood