

Sample preparation: Condition a 6 mL Bond Elut C18 SPE cartridge with 5 mL MeOH and 5 mL water. Filter 30 mL microsomal incubation, add filtrate to SPE cartridge, elute with MeOH. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute the residue, inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.6 10 μm C18 (Alltech)

Mobile phase: Gradient. MeCN:200 mM pH 6.5 ammonium acetate 10:90 for 5 min, to 20:80 over 5 min, maintain at 20:80 for 10 min, to 50:50 over 10 min.

Flow rate: 1

Detector: UV 254, UV 362

CHROMATOGRAM

Retention time: 35

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

rabbit; liver; SPE

REFERENCE

Woolf, T.F.; Black, A.; Chang, T. In vitro metabolism of isoxicam by horseradish peroxidase, *Xenobiotica*, **1989**, *19*, 1369–1377.

Isoxsuprine


Molecular formula: C₁₈H₂₃NO₃

Molecular weight: 301.39

CAS Registry No.: 395-28-8, 579-56-6 (HCl)

Merck Index: 5259

Lednicer No.: 1 69



SAMPLE

Matrix: blood

Sample preparation: Dilute blood with an equal volume of water. 1 mL Plasma or 900 μL diluted blood + 0.9-1 mL buffer + 5 mL freshly distilled ethyl acetate, vortex for 1 min, centrifuge at 1750 g for 7 min. Remove the organic layer and evaporate it almost to dryness under a stream of nitrogen at 57°, evaporate the final 500 μL at room temperature, reconstitute the residue in 100 μL MeCN, vortex for 15 s, inject the whole amount. (Buffer was 26.5 g sodium carbonate and 21 g sodium bicarbonate in 500 mL water, pH 9.48.)

HPLC VARIABLES

Column: 300 × 3.9 10 μm μBondapak phenyl (plasma) or 200 × 4.6 5 μm Spheri-5 RP-18 (blood)

Mobile phase: MeCN:0.05% orthophosphoric acid 17:83 (plasma) or 63:37 (blood)

Flow rate: 2

Injection volume: 100

Detector: F ex 200 no emission filter or UV 254

CHROMATOGRAM

Retention time: 15.1 (plasma), 16.3 (blood)

Internal standard: isoxsuprine hydrochloride

OTHER SUBSTANCES

Extracted: ritodrine

Simultaneous: fenoterol

Noninterfering: acetaminophen, albuterol, betamethasone, bupivacaine, caffeine, chloral hydrate, dexamethasone, diazepam, lignocaine, meperidine, metoclopramide, morphine, nitrazepam, terbutaline

KEY WORDS

isoxsuprine is IS; plasma; whole blood

REFERENCE

Gross,A.S.; Brown,K.F.; Baird-Lambert,J.A.; Nation,R.L. Determination of ritodrine in blood and plasma by high-performance liquid chromatography with fluorescence detection, *J.Chromatogr.*, **1987**, *416*, 400-408.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma + nylidrin + 500 μ L buffer, mix briefly, add 50-75 mg resin, mix at 10-25 rpm for 30 min, discard the supernatant, wash twice with 1 mL buffer, add 500 μ L 50 mg/mL KOH in MeOH:water 50:50, mix for 30 min, inject a 20 μ L aliquot of the eluate. (Buffer was 100 mM citric acid:200 mM Na_2HPO_4 29:71, pH 6.5 (McIlvaine buffer) Wash 20-50 mesh Dowex HCR-S resin twice with water and allow it to equilibrate in buffer).

HPLC VARIABLES

Guard column: 25 \times 4.6 5 μ m Spherisorb ODS-I

Column: 250 \times 4.6 5 μ m Spherisorb ODS-I

Mobile phase: MeCN:MeOH:buffer 30:18:52 containing 1.8 mM octanesulfonic acid (Buffer was 30 mM KH_2PO_4 adjusted to pH 3.0 with concentrated orthophosphoric acid.)

Flow rate: 1

Injection volume: 20

Detector: E, Gynotek M20, glassy carbon working electrode 950 mV, Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 8.5

Internal standard: nylidrin (11)

Limit of detection: 1 ng/mL (from 2 mL plasma)

Limit of quantitation: 5 ng/mL

KEY WORDS

horse; plasma; pharmacokinetics; SPE

REFERENCE

Hashem,A.; Lubczyk,B. Determination of isoxsuprine in equine plasma by high-performance liquid chromatography with electrochemical detection, *J.Chromatogr.*, **1991**, *563*, 216-223.

SAMPLE

Matrix: blood

Sample preparation: Condition a Baker 500 mg C18 SPE cartridge with 3 mL MeOH and 3 mL water. 1 mL Serum + 714 units β -glucuronidase (E. coli, Sigma), heat at 37° overnight, add 2 mL water, add to the SPE cartridge at 0.15 mL/s, wash with 3 mL water, elute with 2 mL MeOH:triethylamine 99:1. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute the residue in 200 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4 5 μ m LC 8 DB (Supelco)

Mobile phase: Gradient. MeOH:buffer 5:95 for 1 min, to 95:5 over 20 min, maintain at 95:5 for 5 min, re-equilibrate at initial conditions for 10 min. (Buffer was 100 mM sodium acetate containing 0.1% triethylamine adjusted to pH 3.4 with 85% phosphoric acid.)

Flow rate: 1.2

Injection volume: 20

Detector: UV 276

CHROMATOGRAM

Limit of detection: 400 ng/mL

KEY WORDS

serum; horse; SPE; pharmacokinetics

REFERENCE

Pompa,G.; Caloni,F.; Montana,M.; Pasqualucci,C. Prolonged presence of isoxsuprine in equine serum after oral administration, *Xenobiotica*, **1994**, *24*, 339–346.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 125 × 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 1.6

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipiprone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, flupromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, ketanserin, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphane, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, J.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, *323*, 191–225.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bicucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenoprofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminos-tilbene, imipramine, indomethacin, isocarbostyryl, isocarboxazid, isoniazid, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, meth-azolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methylidopa, meth-ylidopamine, methylphenidate, methylprednisolone, methyltestosterone, methypyrrol, meto-prolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, nor-epinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphen-butazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, per-santine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenyl-butazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primi-done, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopola-mine, scopolin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sul-faethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sul-fasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicilic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tol-metin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapa-mil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Chirex 3014 (Phenomenex)

Mobile phase: Hexane:1,2-dichloroethane:EtOH/trifluoroacetic acid 55:35:10 (EtOH/trifluoroacetic acid was premixed 20:1.)

Flow rate: 0.7-1

Injection volume: 20

Detector: UV 276

KEY WORDS

chiral; $\alpha = 1.40$ for enantiomers

REFERENCE

Cleveland,T. Pirkle-concept chiral stationary phases for the HPLC separation of pharmaceutical racemates, *J.Liq.Chromatogr.*, **1995**, *18*, 649-671.

Isradipine


Molecular formula: C₁₉H₂₁N₃O₅

Molecular weight: 371.39

CAS Registry No.: 75695-93-1

Merck Index: 5260

Lednicer No.: 4 107



SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 1 mL 4 ng/mL IS in water + 100 μ L 2 M NaOH + 6 mL dichloromethane, mix on a rotary mixer for 20 min, centrifuge at 1500 g for 5 min (if emulsions form stir with a glass rod and re-centrifuge). Remove the organic layer and evaporate it under a stream of nitrogen at 40°, add 500 μ L dichloromethane, vortex, evaporate under a stream of nitrogen at 40°, reconstitute in 100 μ L mobile phase, vortex for 30 s, allow to stand for 10 min, vortex, inject a 75 μ L aliquot.

HPLC VARIABLES

Column: 150 × 3.9 4 μ m Nova-pak C18

Mobile phase: MeOH:10 mM dibutylamine phosphate (Waters D-4) 50:50, pH 2.8-3.0

Column temperature: 48

Flow rate: 1

Injection volume: 75

Detector: UV 325

CHROMATOGRAM

Retention time: 12

Internal standard: PY 108-068 (diethyl ester of isradipine) (13)

Limit of quantitation: 0.5 ng/mL

OTHER SUBSTANCES

Simultaneous: metabolites

KEY WORDS

plasma

REFERENCE

Boutagy,J.; Rumble,F.; Dunagan,F. Determination of isradipine and the oxidative pyridine metabolite in human plasma by high-performance liquid chromatography, *J.Chromatogr.*, **1989**, 487, 483-488.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 22.352

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

SAMPLE

Matrix: cells, culture medium

Sample preparation: Cells. To 100 μ L containing 40×10^6 cell suspension add 100 μ L of a methanolic solution of p-phenoxyphenol and benzoylated spermine so as to give concentrations of 10.63 μ M and 2815 nM, respectively. Vortex for 1 min, sonicate for 3 min, add 200 μ L MeOH, vortex for 1 min, centrifuge at 2000 g for 5 min, inject a 50 μ L aliquot of the supernatant. Culture medium. To 5 mL culture medium add 200 μ L of a methanolic solution of p-phenoxyphenol and benzoylated spermine so as to give concentrations of 10.63 μ M and 2815 nM, respectively. Adjust the pH to 8.0 with 1 M NaOH, add 5 mL chloroform (Caution! Chloroform is a carcinogen!), mix at 30 rpm for 30 min, centrifuge at 1000 g for 5 min, remove the chloroform layer. Adjust the aqueous phase to pH 3.0 with 1 M HCl, add 1 mL 100 mM pH 3 citrate buffer and 1 mL 37 mM tetrabutylammonium hydrogen sulfate. Extract the aqueous phase with 5 mL chloroform at 30 rpm for 30 min, centrifuge at 1000 g for 5 min. Combine the chloroform layers and evaporate them to dryness under a stream of nitrogen, dissolve the residue in 300 μ L mobile phase, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 4.6 5 μ m Chromspher C8 (Chrompack)

Mobile phase: MeOH:THF:buffer 43.6:1:55.4 (Buffer was water containing 160 mM NaH_2PO_4 and 17 mM tetrabutylammonium hydrogen sulfate, pH 3.0.)

Flow rate: 1

Injection volume: 50

Detector: UV 240

CHROMATOGRAM**Retention time:** 37**Internal standard:** p-phenoxyphenol (19), benzoylated spermine (33)**Limit of detection:** 110 pmol/10⁶ (cells), 210 nM (culture medium)**Limit of quantitation:** 350 pmol/10⁶ (cells), 700 nM (culture medium)

OTHER SUBSTANCES**Extracted:** metabolites

REFERENCE

Bidouil,S.; Dubois,J.; Hanocq,M. Isocratic high-performance liquid chromatographic method for the separation of isradipine and its main metabolites. Application to in vitro metabolization by h3A4/OR cells, *J.Chromatogr.B*, **1997**, 693, 359–366.

SAMPLE**Matrix:** formulations**Sample preparation:** Shake bottle by hand, dilute a 1 mL aliquot with MeOH:95% EtOH 1:1 to an expected isradipine concentration of 100 µg/mL, filter (0.22 µm), inject a 10 µL aliquot.

HPLC VARIABLES**Column:** 250 × 4.6 5 µm Spheri-5 ODS (Applied Biosystems)**Mobile phase:** MeOH:THF:water 42:20:38**Flow rate:** 1**Injection volume:** 10**Detector:** UV 240

CHROMATOGRAM**Retention time:** 6

OTHER SUBSTANCES**Simultaneous:** degradation products

KEY WORDS

suspensions; stability-indicating

REFERENCE

MacDonald,J.L.; Johnson,C.E.; Jacobson,P. Stability of isradipine in an extemporaneously compounded oral liquid, *Am.J.Hosp.Pharm.*, **1994**, 51, 2409–2411.

SAMPLE**Matrix:** solutions**Sample preparation:** Prepare a 100 µM solution in buffer, inject a 20 µL aliquot.

HPLC VARIABLES


Column: 100 × 4.6 column containing riboflavin binding proteins (Prepare as follows. Add riboflavin to saturate protein of egg yolk, homogenize with 3 volumes buffer, centrifuge, add the supernatant to a 500 × 30 column of DEAE-cellulose (Whatman) equilibrated with buffer, wash extensively with buffer to remove bound protein, elute riboflavin binding proteins (RFBP) with buffer containing 200 mM NaCl (RFBP has intense yellow color, absorption at 455 nm). Purify RFBP on a Sephadex G-100 column with 50 mM pH 7.5 Tris-HCl buffer as eluent, remove the bound riboflavin by extensive dialysis at pH 3.0. Add 4.5 g N,N-disuccinylimidyl carbonate to 3 g Nucleosil 5NH₂ slurried in MeCN, filter, wash with MeCN, wash with 50 mM pH 7.5 phosphate buffer. Suspend 300 mg RFBP in 50 mM phosphate buffer, add the activated silica, mix gently for 2 h using a rotary evaporator, filter, wash with sterile water, wash with isopropanol:water 1:2, pack in a 100 × 4.6 column.) (Buffer was 100 mM pH 5.3 sodium acetate.)

Mobile phase: 50 mM pH 5.5 KH₂PO₄**Flow rate:** 0.8**Injection volume:** 20**Detector:** UV

CHROMATOGRAM**Retention time:** k' 8.57**OTHER SUBSTANCES****Simultaneous:** flurbiprofen, ketoprofen, nimodipine, suprofen**KEY WORDS**chiral; $\alpha = 1.28$ **REFERENCE**

Massolini, G.; De Lorenzi, E.; Ponci, M. C.; Gandini, C.; Caccialanza, G.; Monaco, H. L. Egg yolk riboflavin binding protein as a new chiral stationary phase in high-performance liquid chromatography, *J. Chromatogr. A*, **1995**, *704*, 55–65.

Itraconazole

Molecular formula: $C_{35}H_{38}Cl_2N_8O_4$ **Molecular weight:** 705.64**CAS Registry No.:** 84625-61-6**Merck Index:** 5262**SAMPLE****Matrix:** blood**Sample preparation:** 1 mL Plasma + 1 mL 1 mg/mL IS in MeCN, add 500 mg KCl, vortex, centrifuge at 1500 g for 5 min, inject an aliquot of the supernatant.**HPLC VARIABLES****Column:** 250 × 4 5 μ m LiChrospher RP8**Mobile phase:** MeCN:water 55:45**Flow rate:** 1.5**Injection volume:** 40**Detector:** UV 263**CHROMATOGRAM****Retention time:** 10-11**Internal standard:** R51012 (14-15) (Janssen, France)**Limit of detection:** 20 μ g/mL**Limit of quantitation:** 40 μ g/mL**OTHER SUBSTANCES****Extracted:** metabolites**KEY WORDS**

plasma

REFERENCE

Cociglio, M.; Hillaire-Buys, D.; Alric, R. Prevalidation statistical design to assess analytical methods. Example of a quick liquid chromatographic assay of itraconazole in serum, *J. Chromatogr. B*, **1997**, *698*, 225–233.

SAMPLE**Matrix:** blood**Sample preparation:** 250 μ L serum + IS + 50 μ L 0.3 N barium hydroxide + 50 μ L 0.4 N zinc sulfate + 1 mL MeCN, vortex, centrifuge at 3521 g for 15 min, evaporate the supernatant to dryness under a stream of nitrogen, reconstitute with 250 μ L mobile phase, inject an aliquot.

HPLC VARIABLES

Guard column: 7.5 × 4.6 5 µm Alltech Alltima C18

Column: 250 × 4.6 5 µm Alltech Alltima C18

Mobile phase: MeCN:MeOH:50 mM pH 6.7 phosphate buffer 47:8:45

Column temperature: 37

Flow rate: 1

Detector: UV 263

CHROMATOGRAM

Internal standard: saperconazole

Limit of detection: 10 ng/mL

Limit of quantitation: 25 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

serum; pharmacokinetics

REFERENCE

Christensen,K.J.; Gubbins,P.O.; Gurley,B.J.; Bowman,J.L.; Buice,R.G. Relative bioavailability of itraconazole from an extemporaneously prepared suspension and from the marketed capsules,*Am.J.Health-Syst.Pharm.*, **1998**, *55*, 261–265.

SAMPLE

Matrix: blood

Sample preparation: Condition a Bond Elut SPE cartridge (No. 607101) with 1 mL MeOH and 1 mL water. 1 mL Serum + 100 µL 1 µg/mL IS in MeOH, mix, add to the SPE cartridge, wash with 1-2 mL water, wash with 1-2 mL MeOH:water 50:50, elute with 400 µL MeOH:triethylamine:concentrated orthophosphoric acid 99.7:0.3:0.3, inject a 20 µL aliquot of the eluate.

HPLC VARIABLES

Column: 100 × 4.6 MPLC RP-18 Spheri-5

Mobile phase: MeCN:water 60:40 containing 20 mM triethylamine, pH adjusted to 2.3 with phosphoric acid

Flow rate: 1

Injection volume: 20

Detector: F ex 260 em 365

CHROMATOGRAM

Retention time: 3

Internal standard: cis-4-[4-[4-[4-[[2-(2,4-dichlorophenyl)-2-(1H-1,2,4-triazol-1-ylmethyl)-1,3-dioxolan-4-yl]methoxy]phenyl]-1-piperazinyl]phenyl]-2,4-dihydro-5-methyl-2-(3-methylbutyl)-3H-1,2,4-triazol-3-one (R 51 012) (4)

Limit of quantitation: 4 ng/mL

OTHER SUBSTANCES

Noninterfering: amphotericin B, cefaclor, imipenem, flucytosine, gentamycin, ketoconazole, netilmicin, salicylic acid, sulfamethoxazole, tienamycin, tobramycin, trimethoprim

KEY WORDS

serum; SPE

REFERENCE

Allenmark,S.; Edebo,A.; Lindgren,K. Determination of itraconazole in serum with high-performance liquid chromatography and fluorescence detection [letter], *J.Chromatogr.*, **1990**, *532*, 203–206.

SAMPLE

Matrix: blood

Sample preparation: 100 µL Plasma or serum + 300 µL 250 nM IS in MeOH, vortex for 1 min, centrifuge at 2000 g for 2 min, inject a 250 µL aliquot of the supernatant.

HPLC VARIABLES**Guard column:** 50 × 4.6 20 µm Ultrasphere ODS**Column:** 150 × 4.6 5 µm Ultrasphere ODS**Mobile phase:** MeCN:water:diethylamine 60:40:0.05 (At the end of each day flush column with MeOH:DMSO 90:10.)**Flow rate:** 1.5**Injection volume:** 250**Detector:** UV 261

CHROMATOGRAM**Retention time:** 6.1**Internal standard:** cis-4-[4-[4-[[2-(2,4-dichlorophenyl)-2-(1H-1,2,4-triazol-1-ylmethyl)-1,3-dioxolan-4-yl]methoxy]phenyl]-1-piperazinyl]phenyl]-2,4-dihydro-5-methyl-2-(3-methylbutyl)-3H-1,2,4-triazol-3-one (R51012) (8.5)**Limit of detection:** 8 nM

OTHER SUBSTANCES**Noninterfering:** amoxicillin, amphotericin B, ampicillin, cimetidine, diazepam, erythromycin, gentamycin, griseofulvin, ketoconazole, miconazole, nystatin, penicillin G, prednisolone, sulfamethoxazole, trimethoprim, zidovudine

KEY WORDSplasma; serum

REFERENCEBadcock, N.R. Micro-scale method for itraconazole in plasma by reversed-phase high-performance liquid chromatography, *J. Chromatogr.*, **1990**, 525, 478–483.

SAMPLE**Matrix:** blood**Sample preparation:** 1 mL Serum + 1 mL 50 mM sodium borate + 100 µL 10 µg/mL IS in MeOH + 200 µL MeOH, extract twice with 4 mL aliquots of heptane:isoamyl alcohol 95:5 in a rotary mixer. Combine the organic phases and add them to 2 mL 1 M sulfuric acid, extract. Discard the organic phase and add 600 µL concentrated ammonium hydroxide to the aqueous phase. Extract the aqueous phase twice with 2.5 mL portions of heptane:isoamyl alcohol 95:5. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 60°, reconstitute the residue in 100 µL MeCN:water 60:40, inject a 20 µL aliquot.

HPLC VARIABLES**Guard column:** C18**Column:** 100 × 4.5 3 µm Hypersil octadecylsilane**Mobile phase:** MeCN:water 40:60 containing 0.03% diethylamine adjusted to pH 7.8 with dilute orthophosphoric acid**Flow rate:** 1**Injection volume:** 20**Detector:** UV 254

CHROMATOGRAM**Retention time:** k' 7.4**Internal standard:** cis-4-[4-[4-[[2-(2,4-dichlorophenyl)-2-(1H-1,2,4-triazol-1-ylmethyl)-1,3-dioxolan-4-yl]methoxy]phenyl]-1-piperazinyl]phenyl]-2,4-dihydro-5-methyl-2-(3-methylbutyl)-3H-1,2,4-triazol-3-one (R51012) (k' 10.8)**Limit of quantitation:** 20 ng/mL

OTHER SUBSTANCES**Extracted:** metabolites, hydroxyitraconazole

KEY WORDSserum

REFERENCE

Law,D.; Moore,C.B.; Denning,D.W. Bioassay for serum itraconazole concentrations using hydroxyitraconazole standards, *Antimicrob.Agents Chemother.*, **1994**, 38, 1561–1566.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Serum + 1.5 mL MeOH, shake for 5 min, centrifuge. Remove the supernatant and evaporate it to dryness, reconstitute the residue in 100 μ L MeOH, inject a 20 μ L aliquot. Alternatively, condition a 100 mg Bakerbond C-18 SPE cartridge with 2 mL MeOH and 2 mL water. 1 mL Plasma + 1 mL water, add to the SPE cartridge, wash with 2 mL water containing 100 μ L MeCN, elute with 2 mL MeOH or MeCN. Evaporate the eluate to dryness, reconstitute the residue in 100 μ L MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 3.5 μ m Separon SGX C-18 (Tessek, Prague)

Mobile phase: MeOH:water:triethylamine 72:28:0.05

Flow rate: 1.1

Injection volume: 20

Detector: UV (wavelength not specified)

CHROMATOGRAM

Retention time: 9

Internal standard: R 51012 (13)

Limit of detection: 10 ng/mL (?)

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

SPE; plasma; serum

REFERENCE

Brandsteterova,E.; Kubalec,P.; Rády,A.; Krcméry,V. Determination of itraconazole and its metabolites in plasma using SPE-HPLC, *Pharmazie*, **1995**, 50, 597–599.

SAMPLE

Matrix: blood

Sample preparation: Condition a 100 mg Bond-Elut C18 SPE cartridge with 2 mL 95% EtOH and 2 mL MeCN:water 15:85. 200 μ L Plasma or whole blood + 50 μ L 100 μ M testosterone propionate in MeOH + 3 mL MeCN:water 15:85, vortex for 30 s, add to the SPE cartridge, wash with 9 mL MeCN:water 30:70, dry, elute with 200 μ L 95% EtOH, inject a 10 μ L aliquot of the eluate.

HPLC VARIABLES

Column: 50 \times 4.6 μ m Supelcosil LC-8DB

Mobile phase: MeOH:buffer 72.5:27.5 (Buffer was 25 mM K_2HPO_4 adjusted to pH 3 with 670 mM phosphoric acid.)

Flow rate: 1

Injection volume: 10

Detector: UV 210

CHROMATOGRAM

Retention time: 2.55

Internal standard: testosterone propionate (3.60)

OTHER SUBSTANCES

Extracted: doxepin

Noninterfering: acetaminophen, N-acetylprocainamide, amitriptyline, aspirin, barbituric acid, brompheniramine, caffeine, carbamazepine, chloramphenicol, chlorpheniramine, clonazepam, desipramine, desmethyldoxepin, digitoxin, digoxin, disopyramide, ethosuximide, felbamate, gentamicin, ibuprofen, imipramine, lidocaine, maprotiline, mephenytoin, mephobarbital, meth-

arbital, methsuximide, methylsuccinimide, nortriptyline, paramethadione, phenacemide, phenobarbital, phensuximide, phenylpropanolamine, phenytoin, primidone, procainamide, protriptyline, quinidine, theophylline, tobramycin, trimethadione, valproic acid, vancomycin

Interfering: clotrimazole

KEY WORDS

plasma; SPE; whole blood

REFERENCE

Rifai,N.; Sakamoto,M.; Law,T.; Platt,O.; Mikati,M.; Armsby,C.C.; Brugnara,C. HPLC measurement, blood distribution, and pharmacokinetics of oral clotrimazole, potentially useful antisickling agent, *Clin.Chem.*, **1995**, *41*, 387–391.

SAMPLE

Matrix: blood

Sample preparation: Condition a 100 mg Bond-Elut C18 SPE cartridge with 2 mL 95% EtOH and 2 mL MeCN:water 15:85. 100 μ L Plasma + 50 μ L 4 μ g/mL + 3 mL MeCN:water 15:85, vortex for 30 s, add to the SPE cartridge, wash with 9 mL MeCN:water 40:60, dry the SPE cartridge, elute with 500 μ L dichloromethane:MeOH 50:50. Evaporate the eluate to dryness, reconstitute the residue in 100 μ L mobile phase, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 50 \times 4.6 5 μ m Supelcosil LC-1

Mobile phase: MeCN:MeOH:25 mM pH 6.3 K_2HPO_4 30:30:40

Flow rate: 2

Injection volume: 10

Detector: UV 263

CHROMATOGRAM

Retention time: 2.1

Internal standard: R051012 (Jannsen) (2.8)

Limit of detection: 50 ng/mL

OTHER SUBSTANCES

Noninterfering: acetaminophen, N-acetylprocainamide, amikacin, amitriptyline, aspirin, barbituric acid, brompheniramine, caffeine, carbamazepine epoxide, carbamazepine, chloramphenicol, chlorpheniramine, clonazepam, clotrimazole, desipramine, desmethyldoxepin, digitoxin, digoxin, disopyramide, doxepin, ethosuximide, felbamate, gentamicin, ibuprofen, imipramine, lidocaine, maprotiline, mephenytoin, mephobarbital, metharbital, methsuximide, methylsuccinimide, nortriptyline, paramethadione, phenacemide, phenobarbital, phensuximide, phenylpropanolamine, phenytoin, primidone, protriptyline, theophylline, tobramycin, trimethadione, vancomycin

KEY WORDS

SPE; plasma

REFERENCE

Rifai,N.; Sakamoto,M.; Platt,O.; Brugnara,C. A high-performance liquid chromatographic assay for the determination of itraconazole concentration using solid-phase extraction and small sample volume, *Ther.Drug Monit.*, **1995**, *17*, 522–525.

SAMPLE

Matrix: blood, tissue

Sample preparation: Tissue. Homogenize 25 mg tissue with 1 mL MeOH. Add 50 μ L 2.5 μ g/mL IS in MeOH, vortex. Centrifuge the sample at 833 g for 15 min. Evaporate the supernatant to dryness under a stream of nitrogen. Reconstitute the residue in 300 μ L MeOH, inject an aliquot. Plasma. Add 50 μ L 12.5 μ g/mL IS to 100 μ L plasma, vortex. Add 300 μ L MeOH, vortex for 1 min. Centrifuge the sample at 833 g for 5 min. Inject a 190 μ L aliquot.

HPLC VARIABLES

Guard column: Novapak Guard-Pak

Column: 100 × 8 4 μm Novapak C18

Mobile phase: MeCN:diethylamine:water 58:0.05:42, adjusted to pH 2.45 with 85% phosphoric acid

Flow rate: 1.5

Injection volume: 190

Detector: F ex 260 em 365

CHROMATOGRAM

Retention time: 14.10

Internal standard: R51012 (Janssen Research Diagnostics, USA) (18.5)

Limit of detection: 5 ng/mL

Limit of quantitation: 10 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; serum; liver; bird

REFERENCE

Cox,S.K.; Orosz,S.; Burnette,J.; Frazier,D. Microassay for determination of itraconazole and hydroxyitraconazole in plasma and tissue biopsies, *J.Chromatogr.B*, **1997**, 702, 175–180.

SAMPLE

Matrix: blood, tissue

Sample preparation: Homogenize (Ultra-Turrax) tissue in water 1:4. 1-2 mL Plasma or homogenate + 100 μL 2-10 μg/mL IS in MeOH + 50 mM pH 7.8 phosphate buffer, extract twice with 4 mL heptane:isoamyl alcohol 98.5:1.5 for 10 min. Combine the organic layers and add them to 3 mL 50 mM sulfuric acid, extract, centrifuge at 1000 g. Remove the aqueous phase and adjust the pH to 9 with concentrated ammonia, extract twice with 2 mL heptane:isoamyl alcohol 98.5:1.5. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 55°, reconstitute the residue in 100 μL mobile phase, inject a 40 μL aliquot.

HPLC VARIABLES

Column: 150 × 3.1 5 μm RSiL C18HL octadecyl (Alltech)

Mobile phase: MeCN:water 60:40 containing 0.05% diethylamine

Flow rate: 0.5

Injection volume: 40

Detector: UV 263

CHROMATOGRAM

Retention time: 4.3

Internal standard: cis-4-[4-[4-[4-[[2-(2,4-dichlorophenyl)-2-(1H-1,2,4-triazol-1-ylmethyl)-1,3-dioxolan-4-yl]methoxy]phenyl]-1-piperazinyl]phenyl]-2,4-dihydro-5-methyl-2-(3-methylbutyl)-3H-1,2,4-triazol-3-one (R 51 012) (5.8)

Limit of detection: 1 ng/mL

KEY WORDS

plasma; human; rat; pharmacokinetics

REFERENCE

Woestenborghs,R.; Lorreyne,W.; Heykants,J. Determination of itraconazole in plasma and animal tissues by high-performance liquid chromatography, *J.Chromatogr.*, **1987**, 413, 332–337.

SAMPLE

Matrix: blood, tissue

Sample preparation: Tissue. Homogenize (Baxter Scientific disposable tissue grinder) esophageal tissue with 20 volumes of cold phosphate-buffered saline (pH 7.8). 250 μL Homogenate + 1 mL MeCN + 1 μL 1 mg/mL IS, vortex for 1 min, centrifuge at 1000 g for 5 min. Remove the supernatant and dry in air for 30 min, reconstitute in 100 μL mobile phase, inject an aliquot. Plasma. 250 μL Plasma + 1 mL MeCN + 1 μL 1 mg/mL IS, vortex for 1 min, centrifuge at

1000 g for 5 min. Remove the supernatant and dry in air for 30 min, reconstitute in 100 μ L mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 150 \times 3.9 Nova-pak C18

Mobile phase: MeCN:10 mM K_2HPO_4 60:40, final pH adjusted to 7.8 with 85% phosphoric acid

Flow rate: 1.3

Detector: UV 263

CHROMATOGRAM

Retention time: 4.35

Internal standard: R51012 (Research Diagnostics Inc.) (6.00) Remarks

Limit of detection: 10 ng/g, 5 ng/mL

OTHER SUBSTANCES

Noninterfering: ampicillin, cefazolin, cefoperazone, ceftriaxone, cefuroxime, clindamycin, fluconazole, folic acid, minocycline, nafcillin, norfloxacin, rifampin, tetracycline, vancomycin, zalcitabine, zidovudine

KEY WORDS

plasma

REFERENCE

Darouiche,R.O.; Setoodeh,A.; Anaissie,E.J. Potential use of a simplified method for determination of itraconazole levels in plasma and esophageal tissue by using high-performance liquid chromatography, *Antimicrob.Agents Chemother.*, **1995**, 39, 757-759.

SAMPLE

Matrix: formulations

Sample preparation: 250 μ L Syrup + 1.5 mL DMF, make up to 10 mL with mobile phase. 250 μ L Solution + 500 μ L 1 mg/mL IS, make up to 10 mL with mobile phase, inject a 5 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Spheri-5 ODS

Mobile phase: MeCN:water:diethylamine 60:40:0.05

Flow rate: 1

Injection volume: 5

Detector: UV 263

CHROMATOGRAM

Retention time: 7

Internal standard: R51012 (Janssen) (11)

KEY WORDS

syrup; stability-indicating; suspensions

REFERENCE

Jacobson,P.A.; Johnson,C.E.; Walters,J.R. Stability of itraconazole in an extemporaneously compounded oral liquid, *Am.J.Health-Syst.Pharm.*, **1995**, 52, 189-191.


Ivermectin

Molecular formula: C₄₈H₇₄O₁₄ (B1a)

Molecular weight: 875.11 (B1a)

CAS Registry No.: 70288-86-7,
70161-11-4 (B1a), 70209-81-3 (B1b)

Merck Index: 5264



SAMPLE

Matrix: blood

Sample preparation: Mix 5 mL Serum with 5.0 mL MeOH:buffer 20:80. Pass through an immunoaffinity column (80 × 7 mm, 1.0 mL bed volume, IgG + CNBr-activated Sepharose 4B, Pharmacia, Sweden) at 1.2 mL/min. Wash with 20 mL MeOH:buffer 10:90 and 5 mL MeOH:water 10:90. Elute with 5 mL MeOH. Evaporate the eluate to dryness at 55°. Add 1 mL MeOH to the residue, vortex for 15 s, inject an aliquot. (Buffer was prepared by dissolving 200 mg KH₂PO₄, 2.9 g Na₂HPO₄·12 H₂O, 200 mg KCl, and 18.8 g NaCl in 900 mL water, adjusting to pH 7.4 with 2 M NaOH, and making up to 1 L with water.)

HPLC VARIABLES

Column: 220 × 4.6 5 μm Spheri-5 RP-18

Mobile phase: MeOH:water 95:5

Flow rate: 1.0

Injection volume: 10

Detector: UV 245

CHROMATOGRAM

Retention time: 6.5

Limit of detection: 200 pg/mL

Limit of quantitation: 2 ng/mL

KEY WORDS

serum; sheep; immunoaffinity; SPE

REFERENCE

Li, J.; Zhang, S. Immunoaffinity column cleanup and liquid chromatographic method for determining ivermectin in sheep serum, *J. AOAC Int.*, **1996**, 79, 1300–1302.

SAMPLE

Matrix: blood

Sample preparation: Condition a Sep-Pak C18 SPE cartridge with 4 mL MeCN:water 50:50. Condition a Sep-Pak silica SPE cartridge with 4 mL MeCN and 4 mL dichloromethane. 5 mL Serum + 5 mL MeCN:water 50:50, add to the C18 SPE cartridge, wash with 4 mL MeCN:water 50:50, blow out excess solvent, elute the C18 SPE cartridge onto the silica SPE cartridge with 4 mL MeCN:dichloromethane 10:90, discard the C18 cartridge and elute the silica cartridge with 4 mL MeCN. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute the residue in 1 mL mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 100 × 8 4 μm Nova-Pak radial PAK

Mobile phase: MeCN:MeOH:water 45:45:10

Flow rate: 1

Detector: UV 245

CHROMATOGRAM**Retention time:** 4**Limit of detection:** 2 ppb

KEY WORDScow; serum; SPE

REFERENCE

Oehler,D.D.; Miller,J.A. Liquid chromatographic determination of ivermectin in bovine serum, *J.Assoc.Off.Anal.Chem.*, **1989**, 72, 59–59.

SAMPLE**Matrix:** blood

Sample preparation: Condition a 500 mg Sep-Pak C18 SPE cartridge with 4 mL MeCN, 5 mL chloroform, 4 mL MeCN, and 4 mL water. 1 mL Plasma + 500 μ L MeCN + 500 μ L 1 ng/mL IS in MeCN, mix for 15 s, centrifuge at 2500 g for 10 min, add the supernatant to the SPE cartridge, dry under vacuum for 15 min, elute with 5 mL chloroform. Evaporate the eluate to dryness under a stream of nitrogen at <50°, reconstitute with 100 μ L N-methylimidazole:MeCN 1:1, add 150 μ L trifluoroacetic anhydride:MeCN 1:2, let stand for <30 s, inject a 100 μ L aliquot.

HPLC VARIABLES**Column:** 250 \times 4.6 5 μ m Zorbax C8**Mobile phase:** MeCN:THF:water 40:40:20**Column temperature:** 30**Flow rate:** 1**Injection volume:** 100**Detector:** F ex 365 em 475

CHROMATOGRAM**Retention time:** 18 (ivermectin B_{1a})**Internal standard:** avermectin B_{1a} (12.5)**Limit of detection:** 20 pg/mL

KEY WORDSplasma; SPE; cow; derivatization

REFERENCE

de Montigny,P.; Shim,J.S.K.; Pivnichny,J.V. Liquid chromatographic determination of ivermectin in animal plasma with trifluoroacetic anhydride and N-methylimidazole as the derivatization reagent, *J.Pharm.Biomed.Anal.*, **1990**, 8, 507–511.

SAMPLE**Matrix:** blood

Sample preparation: Condition a 3 mL C18 SPE cartridge (J.T. Baker) with 4 mL MeCN, 5 mL chloroform, 4 mL MeCN, and 4 mL water. 1 mL Plasma + 50 μ L 86-285 ng/mL IS in MeCN, vortex for 15 s, add 1 mL MeCN, mix, centrifuge at 1130 g for 15 min, add the supernatant to the SPE cartridge. Reconstitute the residue in 3.5 mL MeCN:water 1:2, vortex for 15 s, centrifuge at 1130 g for 15 min, add the supernatant to the SPE cartridge. Wash the SPE cartridge with 4 mL MeCN:water 1:2, dry under vacuum for 1 h, elute with 5 mL MeCN:chloroform 50:50. Evaporate the eluate to dryness under a stream of nitrogen at 45°, reconstitute with two 100 μ L portions of MeCN, evaporate to dryness under a stream of nitrogen at room temperature. Reconstitute with 100 μ L N-methylimidazole:MeCN 1:1, add 150 μ L trifluoroacetic anhydride:MeCN 1:2, let stand for 1.7 min, inject a 150 μ L aliquot.

HPLC VARIABLES**Column:** 100 \times 2 3 μ m MOS-Hypersil-2**Mobile phase:** Gradient. MeCN:water from 72:28 to 92:8 over 15 min.**Flow rate:** 0.3**Injection volume:** 150**Detector:** F ex 365 em 475

CHROMATOGRAM**Retention time:** 22.5**Internal standard:** ivermectin monosaccharide (20)**Limit of detection:** 10 pg/mL

KEY WORDSplasma; SPE; dog; narrow bore; derivatization

REFERENCE

Rabel, S.R.; Stobaugh, J.F.; Heinig, R.; Bostick, J.M. Improvements in detection sensitivity for the determination of ivermectin in plasma using chromatographic techniques and laser-induced fluorescence detection with automated derivatization, *J.Chromatogr.*, **1993**, 617, 79–86.

SAMPLE**Matrix:** blood, tissue

Sample preparation: Whole blood, serum. Condition a 1 mL Bond Elut C18 SPE cartridge with 2 mL MTBE, 2 mL MeCN, and 2 mL MeCN:water 50:50. 500 μ L Whole blood or serum + 20 μ L IS solution + 50 μ L 200 mM zinc sulfate solution + 500 μ L MeCN, vortex, centrifuge for 5 min, add to the SPE cartridge, wash with 2 mL MeCN:water 50:50, elute with 2 mL MTBE. Evaporate the eluate and dissolve the residue in 150 μ L mobile phase, inject a 50 μ L aliquot. Muscle. Condition a 3 mL Bond Elut C18 SPE cartridge with 4 mL MTBE, 4 mL MeCN, and 4 mL MeCN:water 50:50. Weigh out 1 g muscle, add 50 μ L/g IS solution, add 6 mL MeCN:water 50:50, homogenize (Vertis 45), rinse blades and jar with 2 mL MeCN:water 50:50, add 200 mL (sic) 200 mM zinc sulfate, mix, centrifuge, add the supernatant to the SPE cartridge, wash with 4 mL MeCN:water 50:50, elute with 4 mL MTBE. Evaporate the eluate and dissolve the residue in 150 μ L mobile phase, inject a 50 μ L aliquot.

HPLC VARIABLES**Column:** 70 \times 4.6 3 μ m Ultrasphere XL ODS**Mobile phase:** MeCN:MeOH:water 49:33:18**Column temperature:** 56**Flow rate:** 1**Injection volume:** 50**Detector:** UV 245

CHROMATOGRAM**Retention time:** 5.6

Internal standard: dehydroivermectin (13.4) (Prepare by evaporating 1 mL 278 μ g/mL ivermectin in MeCN into a tube. Add 200 μ L 1-methylimidazole, add 300 μ L acetic anhydride, add 900 μ L DMF, mix well, heat at 60° for 15 min, add 4 mL MeCN, pass through a silica SPE cartridge, Evaporate the eluate to 500 μ L, make up to 20 mL with MeCN, use this solution.)

Limit of detection: 2 ng/g (tissue), 2 ng/mL (whole blood, serum)

KEY WORDSwhole blood; serum; muscle; human; cow; SPE

REFERENCE

Dickinson, C.M. Improved high-performance liquid chromatographic method for quantitation of ivermectin in whole blood, serum or muscle tissue, *J.Chromatogr.*, **1990**, 528, 250–257.

SAMPLE**Matrix:** feces

Sample preparation: Stir 5 g feces with 25 mL MeOH for 25 min, centrifuge at 1500 g for 15 min, concentrate the supernatant to 7 mL under reduced pressure at 80°, centrifuge at 1500 g for 15 min, add procaine to a concentration of 50 ppm, make up to 10 mL with MeOH, filter (0.50 μ m PTFE), inject a 20 μ L aliquot of the filtrate.

HPLC VARIABLES**Column:** 300 \times 9 Bondclone 10C18 (Phenomenex)**Mobile phase:** MeCN:MeOH:water 47:33:20**Flow rate:** 1

Injection volume: 20**Detector:** UV 245

CHROMATOGRAM**Retention time:** 12.6**Internal standard:** procaine (7.3)**Limit of detection:** 20 ng/g

KEY WORDS

cow

REFERENCE

Bernal,J.L.; Del Nozal,M.J.; Salas,M.; Galante,E.; Lumaret,J.P. HPLC determination of ivermectin in cattle dung following subcutaneous injection, *J.Liq.Chromatogr.*, **1994**, 17, 2429–2444.

SAMPLE**Matrix:** milk

Sample preparation: Prepare a SPE cartridge by adding 2 g 40 μ m Bondesil C18 18% load endcapped (Varian) to a 25 mL syringe barrel fitted with a 20 μ m frit, wash with 5 mL petroleum ether, 5 mL acetone, and two 5 mL aliquots of MeOH, aspirate with full vacuum for <5 s (A). Condition a 500 mg Bond Elut LRC silica SPE cartridge with 3 mL hexane:ethyl acetate 60:40 (B). Condition a 500 mg Bond Elut LRC silica SPE cartridge with 4 mL chloroform (C). 25 mL Milk + 200 μ L 500 ng/mL abamectin (avermectins) in MeOH, mix, add 5 mL to the SPE cartridge (A), mix milk with C18 material, let stand for 2 min, wash spatula with water, wash with two 5 mL portions of water, elute with 10 mL ethyl acetate, allow eluate to pass through a 5 cm layer of anhydrous sodium sulfate. Evaporate the eluate to dryness under a stream of nitrogen below 50°, add 2 mL hexane:ethyl acetate 60:40 to the oily residue, vortex, sonicate for 1 min, add mixture to SPE cartridge (B), rinse in with 1 mL hexane:ethyl acetate 60:40, wash with 5 mL hexane:ethyl acetate 60:40, elute with 5 mL MeOH:ethyl acetate 50:50. Evaporate the eluate to dryness under a stream of nitrogen below 60° (this residue should have no moisture in it), reconstitute the residue in 100 μ L reagent, vortex gently for a few s, heat at 95° for 1 h, cool, add 1 mL chloroform, vortex, add to SPE cartridge (C), wash in with three 1 mL portions of chloroform, elute with 2 mL chloroform. Collect all the eluate and evaporate it to dryness under a stream of nitrogen below 60°, reconstitute in 500 μ L MeOH, inject a 50 μ L aliquot. (Prepare reagent by sequentially mixing 900 μ L DMF, 300 μ L acetic anhydride, and 200 μ L N-methylimidazole just before use.)

HPLC VARIABLES**Guard column:** Newguard RP-18 (Brownlee)**Column:** 250 \times 4.6 5 μ m Econosil C18**Mobile phase:** MeOH:THF:water 85:15:5**Flow rate:** 1**Injection volume:** 50**Detector:** F ex 364 em 455

CHROMATOGRAM**Retention time:** 15**Internal standard:** abamectin (avermectins) (10.5)**Limit of detection:** <1 ppb

KEY WORDS

cow; SPE; MSPD; derivatization; silylate glassware

REFERENCE

Schenck,F.J. Isolation and quantification of ivermectin in bovine milk by matrix solid phase dispersion (MSPD) extraction and liquid chromatographic determination, *J.Liq.Chromatogr.*, **1995**, 18, 349–362.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 250 \times 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, ivermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenoprofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminos-tilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclufenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methylidopa, methyldopamine, methylphenidate, methylprednisolone, methyltestosterone, methypyrrolon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, nor-epinephrine, nortriptyline, noscapine, nyldrin, oxazepam, oxycodone, oxymorphone, oxyphen-butazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenyl-butazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopola-mine, scopoletin, secobarbital, strychnine, sulfacetamide, sufadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfa-soxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tol-metin, tranlylcypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapa-mil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233–242.

SAMPLE

Matrix: tissue

Sample preparation: Condition a 3 mL 500 mg Isolute C8 end-capped SPE cartridge (Inter-national Sorbent Technologies) with 5 mL MeCN and 5 mL MeCN:water 1:2 containing 0.1% triethylamine. Condition a 3 mL 500 mg Isolute silica SPE cartridge (International Sorbent Technologies) with 5 mL ethyl acetate:hexane 40:60. Homogenize (Tissuemizer) 5 g tissue in 25 mL MeCN, centrifuge at 3000 rpm for 5 min, decant supernatant into 50 mL water and 75

μ L triethylamine. mix, pass through the C8 SPE cartridge at 2 mL/min, discard eluate, dry column under vacuum for 5 min, elute with 5 mL MeCN. Dry the eluate under a stream of nitrogen at 50-55°, resolubilize dry residues in 5 mL of ethyl acetate:hexane 40:60, vortex briefly, pass through the silica SPE cartridge at 2 mL/min. Rinse the reservoir with 5 mL ethyl acetate:hexane 40:60, add the rinse to the SPE cartridge. Dry the SPE cartridge for 5 min, elute with 5 mL MeOH:ethyl acetate 50:50, dry the eluate under a stream of nitrogen at 50-55°. Reconstitute the extract with 200 μ L of fresh methylimidazole:MeCN 50:50, vortex briefly, add 300 μ L of fresh trifluoroacetic anhydride:MeCN 1:2, vortex briefly, dry under a stream of nitrogen at 50-55° for 15 min. Add 500 μ L MeOH:ammonium acetate:molecular sieves 4:1:1), vortex briefly, dry under a stream of nitrogen at 50-55°. Add 1 mL MeCN, vortex thoroughly, filter (0.45 μ m), inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 200 \times 4.6 5 μ m Hypersil ODS (C18)

Mobile phase: MeCN-water 90:10

Column temperature: 65

Flow rate: 1.0

Injection volume: 50

Detector: F ex 272 em 465

CHROMATOGRAM

Retention time: 11 (B1b), 12.5 (B1a)

Limit of detection: 0.25 ppb

OTHER SUBSTANCES

Extracted: doramectin

KEY WORDS

SPE; derivatization; salmon; muscle

REFERENCE

Rupp,H.S.; Turnipseed,S.B.; Walker,C.C.; Roybal,J.E.; Long,A.R. Determination of ivermectin in salmon muscle tissue by liquid chromatography with fluorescence detection, *J.AOAC Int.*, **1998**, *81*, 549-553.

SAMPLE

Matrix: tissue

Sample preparation: Condition a C8 Bond-Elut SPE cartridge with 5 mL MeCN and 5 mL MeCN:water 30:70 containing 0.1% triethylamine. Mix 5 g muscle with 15 mL MeCN in a high speed blender for 1 min. Centrifuge at 3000 rpm for 5 min. Dilute 15 mL of the supernatant with 35 mL water, add 50 μ L triethylamine, add to the SPE cartridge, elute with 5 mL MeCN, evaporate the eluate to dryness under a stream of nitrogen at 60°. Dissolve the residue in 1 mL MeOH, filter (0.45 μ m Millipore filter). Transfer a 500 μ L aliquot to a silanized 3 mL reaction vial, evaporate to dryness. Add 100 μ L derivatization reagent, heat at 95° for 1 h. Cool to room temperature, add 1 mL chloroform, vortex, add quantitatively to a Sep-Pak silica SPE cartridge with 3-4 mL chloroform. Elute with 9 mL chloroform, evaporate the combined chloroform eluates to dryness. Dissolve the residue in 500 μ L MeOH, filter (0.45 μ m Millipore filter), inject a 50 μ L aliquot. (Caution! Chloroform is a carcinogen! Derivatization reagent was 1-methylimidazole:acetic anhydride:DMF 2:6:9, freshly prepared.)

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Shandon ODS Hypersil C18

Mobile phase: MeOH:water 97:3

Column temperature: 30

Flow rate: 1

Injection volume: 50

Detector: F ex 364 em 470

CHROMATOGRAM

Retention time: 14

Limit of detection: 500-1000 ng/kg

KEY WORDS

pig; muscle; SPE; derivatization

REFERENCE

Nordlander, I.; Johnsson, H. Determination of ivermectin residues in swine tissues—an improved clean-up procedure using solid-phase extraction, *Food Addit. Contam.*, **1990**, 7, 79–82.

SAMPLE

Matrix: tissue

Sample preparation: Add 15 mL MeOH to 5 g ‘homogenized liver, shake thoroughly by hand and again using a shaking apparatus at medium speed for 1 h. Adjust the volume to 20 mL with MeOH, shake, centrifuge at 2000 g for 5 min. Mix 10 mL supernatant with 40 mL phosphate buffer A, add to immunoaffinity column (column preparation described in detail in paper) at 1.2 mL/min, wash with 40 mL MeOH:phosphate buffer B 10:90 and 10 mL MeOH:water 20:80. Elute with 3 mL MeOH, evaporate the eluate to less than 1 mL on a rotary evaporator at 55°, vortex with 5 mL ethyl acetate for 15 s. Collect the organic layer, evaporate to dryness at 55°, redissolve the residue in 1 mL mobile phase by vortexing for 15 s, filter (0.45 µm filter), inject a 100 µL aliquot of the filtrate. (Phosphate buffer A was 200 mg KH₂PO₄, 2900 mg Na₂HPO₄, 200 mg KCl, and 8800 mg NaCl in 900 mL water adjusted to pH 7.4 with 2 M NaOH and diluted to 1 L with water. Phosphate buffer B was 200 mg KH₂PO₄, 2900 mg Na₂HPO₄, 200 mg KCl, and 29.3 g NaCl in 900 mL water adjusted to pH 7.4 with 2 M NaOH and diluted to 1 L with water.)

HPLC VARIABLES

Column: 220 × 4.6 5 µm Brownlee C18

Mobile phase: MeCN:MeOH:water 45:45:10

Flow rate: 1

Injection volume: 100

Detector: UV 245

CHROMATOGRAM

Retention time: 15

Limit of detection: 2 mg/g

Limit of quantitation: 5 mg/g

KEY WORDS

immunoaffinity; liver; pig; SPE

REFERENCE

Li, J.S.; Li, X.W.; Hu, H.B. Immunoaffinity column cleanup procedure for analysis of ivermectin in swine liver, *J. Chromatogr. B*, **1997**, 696, 166–171.

SAMPLE

Matrix: tissue

Sample preparation: Condition a Bond-Elut C8 SPE cartridge (Analytichem International) with 5 mL MeCN and 5 mL MeCN:water 30:70 containing 0.1% triethylamine. Mix 5 g muscle with 15 mL MeCN in a high speed blender for 1 min, centrifuge at 3000 rpm for 5 min. Dilute 15 mL supernatant with 35 mL water, add 50 µL triethylamine, add the entire mixture to the SPE cartridge, elute with 5 mL MeCN, evaporate the eluate to dryness under a stream of nitrogen at 60°. Dissolve the residue in 1 mL MeOH and filter (0.45 µm). Evaporate a 500 µL aliquot of the filtrate to dryness in a silanized vial, add 100 µL derivatization reagent (1-methyl-imidazole:acetic anhydride:dimethylformamide 2:6:9, freshly prepared), heat in an oven at 95° for 1 h. Cool to room temperature, add 1 mL chloroform and vortex briefly. (Caution! Chloroform is a carcinogen!) Transfer quantitatively to a Sep-Pak silica SPE cartridge with 3–4 mL chloroform, elute with 9 mL chloroform, evaporate the combined chloroform eluates to dryness. Dissolve the residue in 500 µL MeOH, filter (0.45 µm), inject a 50 µL aliquot of the filtrate.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Hypersil ODS C18 (Shandon)

Mobile phase: MeOH:water 97:3

Column temperature: 30

Flow rate: 1

Injection volume: 50

Detector: F ex 364 em 470

CHROMATOGRAM**Retention time:** 14 (dihydroavermectin B1A)**Limit of quantitation:** 500-1000 ng/g

KEY WORDS

derivatization; muscle; SPE; pig

REFERENCE

Nordlander, I.; Johnsson, H. Determination of ivermectin residues in swine tissues--an improved clean-up procedure using solid-phase extraction, *Food Addit. Contam.*, **1990**, *7*, 79-82.

SAMPLE**Matrix:** tissue

Sample preparation: Condition a Bond Elut C8 SPE cartridge with 5 mL MeCN and 5 mL MeCN:water:triethylamine 30:70:0.1. Homogenize (Silverson) 5 g frozen minced tissue with 15 mL MeCN at full speed for 1 min, centrifuge at 4° at 2000 g for 10 min. Remove a 13 mL aliquot of the supernatant and add it to 35 mL water and 50 µL triethylamine, elute with 5 mL MeCN. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute with 200 µL 1-methylimidazole:MeCN 1:1, add 300 µL trifluoroacetic anhydride:MeCN 1:2, mix, store cold, inject an aliquot.

HPLC VARIABLES**Column:** 250 × 4.6 Partisil 5 ODS-3**Mobile phase:** MeOH:water 96:4**Flow rate:** 1.8**Detector:** F ex 364 em 470

CHROMATOGRAM**Retention time:** 6 (ivermectin B_{1a})**Limit of detection:** 1 ng/g

KEY WORDS

SPE; salmon; derivatization; brain; gill; kidney; liver; muscle; skin; spleen

REFERENCE

Kennedy, D.G.; Cannavan, A.; Hewitt, S.A.; Rice, D.A.; Blanchflower, W.J. Determination of ivermectin residues in the tissues of Atlantic salmon (*Salmo salar*) using HPLC with fluorescence detection, *Food Addit. Contam.*, **1993**, *10*, 579-584.

SAMPLE**Matrix:** tissue

Sample preparation: Condition a 6 mL 500 mg Bond Elut C18 SPE cartridge with 5 mL MeCN and 5 mL MeCN:water:triethylamine 30:70:0.1. Homogenize (Polytron) 5 g tissue and 15 mL MeCN for 20 s, rinse probe with 5 mL MeCN, shake mechanically at high speed for 5 min, centrifuge at 2000 g for 5 min. Re-extract the solid with 10 mL MeCN. Add the supernatants to the alumina column. Combine the eluates, add 70 mL water, add 100 µL triethylamine, mix, add to the C18 SPE cartridge, pull air through the SPE cartridge for 3 min, elute with 5 mL MeCN. Evaporate the eluate to dryness under a stream of nitrogen at 60°, reconstitute the residue in 100 µL freshly prepared reagent, vortex for 15 s, heat at 95-100° for 45 min, cool, add 1 mL chloroform, vortex, add to a 2.8 mL 500 mg Bond Elut silica SPE cartridge, elute with three 3 mL aliquots of chloroform. Combine the eluates and evaporate them to dryness under a stream of nitrogen at 60°, reconstitute the residue in 1 mL MeOH, filter, inject a 40 µL aliquot. (Prepare alumina column as follows. Shake 94 g Brockman Activity I neutral alumina (Fisher) and 6 mL water for 45 min, add 4.5 g alumina to an 8 mL column with a frit. Reagent was 200 µL 1-methylimidazole, 600 µL acetic anhydride, and 900 µL DMF.)

HPLC VARIABLES**Guard column:** 30 × 4.6 RP-18 (Brownlee)**Column:** 250 × 4.6 RP-18 OD-224 (Brownlee)**Mobile phase:** MeOH:water 97:3**Flow rate:** 1.8

Injection volume: 40
Detector: F ex 365 em 425

CHROMATOGRAM

Retention time: 9.3
Limit of detection: 2 ppb

KEY WORDS

SPE; cow; pig; sheep; fish; liver; muscle; derivatization

REFERENCE

Salisbury, C.D.C. Modified method for the determination of ivermectin residues in animal tissues, *J.AOAC Int.*, **1993**, 76, 1149–1151.

SAMPLE

Matrix: tissue

Sample preparation: Condition a Bond Elut C18 SPE cartridge with 4 mL MeCN and 4 mL MeCN:water 10:90. 4 g Minced meat or liver + 40 mL MeCN + 3.5 mL water, vortex for 2 min, centrifuge at 2000 rpm for 10 min, remove supernatant, repeat extraction with 20 mL MeCN and 3.5 mL water. Combine the supernatants and evaporate them to 6 mL under reduced pressure (all the MeCN should be removed), add 6 mL water to the residue, add to the SPE cartridge, pull air through the cartridge for 10 min, elute with 5 mL MeCN. Evaporate the eluate under a stream of nitrogen, add 150 μ L trifluoroacetic anhydride:MeCN 1:2 to the residue, add 100 μ L N-methylimidazole:MeCN 50:50, shake, store in the dark, inject a 10–50 μ L aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 μ Bondapak C18
Mobile phase: MeOH:water 95:5 (At the end of the day flush column with 30 mL MeOH.)
Flow rate: 1.5
Injection volume: 10–50
Detector: F ex 364 em 470

CHROMATOGRAM

Retention time: 11
Limit of quantitation: 5 ng/g

KEY WORDS

meat; liver; cow; pig; muscle; SPE; derivatization

REFERENCE

Degroodt, J.M.; Wyhowski de Bukanski, B.; Srebrnik, S. Determination of ivermectin residues in meat and liver by HPLC and fluorometric detection, *J.Liq.Chromatogr.*, **1994**, 17, 1419–1426.

SAMPLE

Matrix: tissue

Sample preparation: Condition a 3 mL Bakerbond C8 SPE cartridge with 5 mL MeCN and 5 mL MeCN:water:triethylamine 30:70:0.1. Prepare a 55 mm column of 70–230 mesh Kieselgel 60 (Merck) in a Pasteur pipette, condition with 3 mL hexane:isopropanol 60:40. Homogenize (Polytron) 5 g blended tissue with 12 mL MeCN, rinse homogenizer with 3 mL MeCN, centrifuge the mixture at 4000 rpm for 10 min. Remove the supernatant and make up to 50 mL with water, add 50 μ L triethylamine, shake, add to the SPE cartridge, elute with 5 mL MeCN at 1 drop/s. Evaporate the eluate to 300 μ L under a stream of nitrogen at 40°, transfer to a smaller vial with MeCN, evaporate to dryness under a stream of nitrogen at 60°, reconstitute the residue in 100 μ L 1-methylimidazole:MeCN 50:50, cool in an ice bath, add 150 μ L trifluoroacetic anhydride:MeCN 1:2, shake at room temperature for 1 min, add to the column, rinse the vial with 500 μ L hexane:isopropanol 60:40, add the rinse to the column, elute with 1 mL hexane:isopropanol 60:40. Evaporate the eluate to dryness under a stream of nitrogen with heating, reconstitute with 250 μ L MeCN, inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: 5 μ m LiChrospher 100 RP-18

Column: 125 × 4 5 µm LiChrospher 60 RP-select B
Mobile phase: MeOH:water 95:5
Column temperature: 40
Flow rate: 1
Injection volume: 20
Detector: F ex 365 em 465

CHROMATOGRAM

Retention time: 3.5
Limit of detection: 1 ppb
Limit of quantitation: 2.5 ppb

KEY WORDS

derivatization; muscle; liver; pig; cow; SPE

REFERENCE

Guggisberg,D.; Sievi,M.; Koch,H. Methode zur quantitativen Bestimmung von Ivermectin in Fleisch und Leber mit HPLC und Vorsäulenderivatisierung [Method for the quantitative determination of ivermectin in meat and liver by HPLC and pre-column derivatization], *Mitteilungen aus dem Gebiete der Lebensmitteluntersuchung und Hygiene*, **1994**, 85, 395–405.

SAMPLE

Matrix: tissue

Sample preparation: Condition a 6 mL 500 mg Bakerbond C18 SPE cartridge with three 5 mL portions of MeOH, 5 mL MeCN, and three 5 mL portions of MeCN:water:triethylamine 30:70:0.1. Condition a Waters silica SPE cartridge with 8 mL chloroform. Homogenize (Ultraturrax) 5 g minced tissue with 15 mL MeCN at high speed for 3 min, rinse blade with 2 mL MeCN, sonicate for 15 min, centrifuge at 3000 rpm for 5 min, filter (paper), extract the residue again with 10 mL MeCN, wash the filter with 3 mL MeCN. Combine the organic layers and add 70 mL water and 100 µL triethylamine, stir thoroughly, add to the C18 SPE cartridge, wash with two 5 mL portions of MeCN:water 50:50, elute with 7 mL MTBE at 2 mL/min. Store the eluate overnight at -20°, remove the organic layer and evaporate it to dryness under a stream of nitrogen at 50°, reconstitute the residue in 3 mL MeOH, add 100 µL water, add 3 mL hexane, vortex, remove the hexane layer, repeat the hexane wash. Extract the combined hexane layers with 1 mL MeOH. Combine the MeOH layers and evaporate them to dryness under a stream of nitrogen at 50°, heat in a vacuum oven at 50° for 30 min, reconstitute the residue in 150 µL 1-methylimidazole:acetic anhydride:DMF 2:3:9 (freshly prepared), vortex for 30 s, heat at 100° for 1 h, cool, add 1 mL chloroform, vortex, add to the silica SPE cartridge, elute with three 3 mL portions of chloroform. Evaporate the eluate to dryness under a stream of nitrogen at 50°, reconstitute the residue in 400 µL MeOH, vortex, inject a 20 µL aliquot.

HPLC VARIABLES

Guard column: 20 × 4.6 5 µm Supelcosil LC-18
Column: 150 × 4.6 5 µm Supelcosil LC-18
Mobile phase: MeOH:water 95:5
Flow rate: 1.8
Injection volume: 20
Detector: F ex 360 em 470

CHROMATOGRAM

Retention time: 7
Limit of detection: 2 ng/g

KEY WORDS

derivatization; SPE; liver; muscle; fat; guinea pig; cow; pig; horse; sheep; pharmacokinetics

REFERENCE

Dusi,G.; Curatolo,M.; Fierro,A.; Faggionato,E. Determination of the antiparasitic drug ivermectin in liver, muscle and fat tissue samples from swine, cattle, horses and sheep using HPLC with fluorescence detection, *J.Liq.Chromatogr.Rel.Technol.*, **1996**, 19, 1607–1616.