


Josamycin

Molecular formula: $C_{42}H_{69}NO_{15}$

Molecular weight: 828.01

CAS Registry No.: 16846-24-5,
51016-68-3 (propionate)

Merck Index: 5280



SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 231.1

CHROMATOGRAM

Retention time: 16.763

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

SAMPLE

Matrix: bulk

HPLC VARIABLES

Column: 250 \times 4.6 8 μ m 1000 Å PLRP-S (Polymer Labs., UK)

Mobile phase: MeCN:buffer:water 52:20:28 (Prepare buffer by mixing 200 mM K_3PO_4 with 200 mM K_2HPO_4 to obtain a pH of 10.0.)

Column temperature: 60

Flow rate: 1

Injection volume: 20

Detector: UV 232

CHROMATOGRAM

Retention time: 20 (josamycin propionate)

OTHER SUBSTANCES

Simultaneous: josamycin, leucomycin A4 propionate, josamycin 2',9-dipropionate, josamycin 3'',9-dipropionate, platenomycin A1 propionate

REFERENCE

Roets,E.; Lepoudre,X.; Van Rompaey,V.; Velghe,G.; Liu,L.; Hoogmartens,J. Liquid chromatography of josamycin propionate on poly(styrene-divinylbenzene), *J. Chromatogr.A*, **1998**, 812, 303-308.

SAMPLE

Matrix: tissue

Sample preparation: Condition a 1 mL 100 mg Bond-Elut diol SPE cartridge with 1 mL chloroform (Caution! Chloroform is a carcinogen!). Mix 2 g minced muscle tissue with 800 μ L water. Stir, vortex for 1 min at maximum speed, let stand for 15 min. Add 2 mL pH 8 buffer, mix briefly, add 10 mL chloroform. Stir at 100 rpm for 15 min, centrifuge at 4000 g for 10 min, discard the aqueous layer, filter the chloroform layer through glass wool. Add the filtrate to the SPE cartridge, wash with 500 μ L chloroform, dry under vacuum, elute with three 200 μ L portions of MeOH:100 mM ammonium acetate 50:50, inject a 200 μ L aliquot of the eluate. (Buffer was 33.46 g K_2HPO_4 and 1.046 g KH_2PO_4 in 1 L water.)

HPLC VARIABLES

Guard column: 4 \times 4.5 μ m C18

Column: 125 \times 4.5 μ m Lichrospher RP18

Mobile phase: Gradient. A was MeCN. B was MeOH. C was 0.1% trifluoroacetic acid in water. A:B:C from 20:20:60 to 25:55:20 in 10 (?) min

Flow rate: 0.5

Injection volume: 200

Detector: MS, HP Model 5989 A, desolvation chamber 60°, source 280° and 300° in negative and positive chemical ionization mode, respectively, with methane as reagent, quadrupole 100°, particle beam nebulizer helium 345 kPa, scan m/z 735.4-828.5 in NCI and 828.5-769.4 in PCI

CHROMATOGRAM

Retention time: 7.2

Limit of detection: 50 μ g/kg

OTHER SUBSTANCES

Extracted: erythromycin, spiramycin, tilmicosin, tylosin

KEY WORDS

muscle; cow; SPE

REFERENCE

Delépine,B.; Hurtaud-Pessel,D.; Sanders,P. Multiresidue method for confirmation of macrolide antibiotics in bovine muscle by liquid chromatography/mass spectrometry, *JAOAC Int.*, **1996**, 79, 397-404.

SAMPLE

Matrix: tissue

Sample preparation: Condition a 500 mg Bond Elut SCX SPE cartridge (Varian) with 5 mL MeOH and 10 mL 100 mM pH 4.4 KH_2PO_4 buffer. Homogenize 5 g tissue with 100 mL MeOH: 0.3% metaphosphoric acid 30:70 at high speed for 2 min, filter through 2 mm Hyflo Super-Cel coated on a pressure funnel (when filtering liver or kidney add several grams of Hyflo Super-Cel to the homogenized solution before filtration). Evaporate the filtrate to ca. 20 mL under reduced pressure at 45°, add to the SPE cartridge, wash with 10 mL distilled water and 5 mL 100 mM pH 8.9 K_2HPO_4 buffer, elute with 10 mL MeOH, evaporate the eluate to dryness under reduced pressure at 45°, dissolve the residue in 1 mL MeCN:50 mM pH 4.5 NaH_2PO_4 buffer 30:70, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Puresil 5C18 (Waters)

Mobile phase: Gradient. A:B from 60:40 to 0:100 over 16 min. A was buffer. B was MeCN:buffer 40:60 (Buffer was 2.5 g KH_2PO_4 dihydrate and 0.65 mL 85% phosphoric acid dissolved in 1 L distilled water, pH 2.5.)

Column temperature: 35

Flow rate: 1

Injection volume: 10

Detector: UV 232 for 9 min, UV 287 for 2 min, UV 232 for 4 min

CHROMATOGRAM

Retention time: 13.73

Limit of detection: 50 ng/g

OTHER SUBSTANCES

Extracted: leucomycin (kitasamycin), mirosamicin, spiramycin, tylosin

KEY WORDS

meat; SPE

REFERENCE

Horie,M.; Saito,K.; Ishii,R.; Yoshida,T.; Haramaki,Y.; Nakazawa,H. Simultaneous determination of five macrolide antibiotics in meat by high-performance liquid chromatography, *J.Chromatogr.A*, **1998**, 812, 295–302.

SAMPLE

Matrix: tissue

Sample preparation: Blend (Virtis model 45 with U-shaped blades) 2.5 g tissue with 20 mL MeCN:10 mM pH 6.0 phosphate buffer 65:35 for 10 min, centrifuge at 5° at 8500 g for 5 min. Remove the supernatant and adjust the volume to 25 mL with MeCN:10 mM pH 6.0 phosphate buffer 65:35. Remove a 1.5 mL aliquot and add it to 5 mL isooctane, shake for 10 min, centrifuge at 5° at 3000 g for 5 min, discard the organic layer. Add 500 µL reagent to the aqueous layer, mix, heat at 90° for 2 h, cool, inject a 100 µL aliquot. (Prepare reagent by dissolving 1 g cyclohexa-1,3-dione and 25 g ammonium acetate in 60 mL water and 8 mL concentrated HCl, make up to 100 mL with water. Store at 5°, discard after 1 month.)

HPLC VARIABLES

Guard column: 4 × 4 5 µm LiChrospher 100 RP-18 end capped

Column: 125 × 4 5 µm LiChrospher 100 RP-18 end capped

Mobile phase: MeCN:MeOH:10 mmole pH 6.0 phosphate buffer 45:5:50

Column temperature: 45

Flow rate: 1.5

Injection volume: 100

Detector: F ex 375 em 450

CHROMATOGRAM

Retention time: 10.3

Limit of detection: 25 ng/g

OTHER SUBSTANCES

Simultaneous: spiramycin, tylosin

Noninterfering: acetaldehyde, benzaldehyde, erythromycin, formaldehyde

KEY WORDS

pig; muscle; liver; kidney; fat; derivatization

REFERENCE

Leroy,P.; Decolin,D.; Nicolas,A.; Archimbault,P. Determination of josamycin residues in porcine tissues using high-performance liquid chromatography with pre-column derivatization and spectrofluorometric detection, *Analyst*, **1994**, 119, 2743–2747.