


# Kanamycin

**Molecular formula:**  $C_{18}H_{36}N_4O_{11}$  (A)

**Molecular weight:** 484.51 (A)

**CAS Registry No.:** 8063-07-8, 25389-94-0 (A sulfate), 59-01-8 (A), 4696-78-8 (B)

**Merck Index:** 5293



	R	R'
Kanamycin A	NH <sub>2</sub>	OH
Kanamycin B	NH <sub>2</sub>	NH <sub>2</sub>
Kanamycin C	OH	NH <sub>2</sub>

## SAMPLE

**Matrix:** blood, urine

**Sample preparation:** Plasma. 300  $\mu$ L Plasma + 30  $\mu$ L water + 100  $\mu$ L 2 M perchloric acid, vortex for 2-3 s, centrifuge at 1000 g for 5 min. Remove the supernatant and neutralize it with 1.5 M NaOH, add 300  $\mu$ L buffer, add 400  $\mu$ L DMSO, add 100  $\mu$ L 2% 2,4-dinitrofluorobenzene in EtOH, vortex, heat at 64° for 30 min, add 3 mL toluene, vortex, centrifuge, discard the upper toluene layer, add 3 mL MeCN:toluene 50:50, vortex for 5-10 s. Remove the upper organic layer and evaporate it to dryness under a stream of nitrogen at 30-40°, reconstitute the residue in 1 mL MeCN:water 50:50, inject a 20  $\mu$ L aliquot. Urine. Dilute urine 100-fold with water. 300  $\mu$ L Diluted urine + 30  $\mu$ L water + 300  $\mu$ L buffer + 400  $\mu$ L DMSO + 100  $\mu$ L 2% 2,4-dinitrofluorobenzene in EtOH, vortex, heat at 64° for 30 min, add 3 mL toluene, vortex, centrifuge, discard the upper toluene layer, add 3 mL MeCN:toluene 50:50, vortex for 5-10 s. Remove the upper organic layer and evaporate it to dryness under a stream of nitrogen at 30-40°, reconstitute the residue in 1 mL MeCN:water 50:50, inject a 20  $\mu$ L aliquot. (Prepare buffer by mixing 80 mL 100 mM Na<sub>2</sub>HPO<sub>4</sub> and 20 mL 100 mM NaH<sub>2</sub>PO<sub>4</sub>, adding 1 g Tris HCl, and adjusting the pH to 7.8 with 6 M HCl.)

## HPLC VARIABLES

**Column:** 250  $\times$  4.6 5  $\mu$ m Zorbax SB-C18

**Mobile phase:** MeOH:water 64:36, adjusted to pH 3.0 with phosphoric acid

**Column temperature:** 50

**Flow rate:** 2

**Injection volume:** 10-20

**Detector:** UV 350

## CHROMATOGRAM

**Retention time:** 24.0 (kanamycin B)

**Internal standard:** kanamycin B (24.0)

## OTHER SUBSTANCES

**Extracted:** paromomycin

## KEY WORDS

derivatization; plasma; pharmacokinetics; kanamycin B is IS

## REFERENCE

Lu,J.; Cwik,M.; Kanyok,T. Determination of paromomycin in human plasma and urine by reversed-phase high-performance liquid chromatography using 2,4-dinitrofluorobenzene derivatization, *J.Chromatogr.B*, **1997**, 695, 329-335.

## SAMPLE

**Matrix:** bulk

**Sample preparation:** Prepare a 2 mg/mL solution in 20 mM pH 9.0 borate buffer, remove a 5 mL aliquot and add it to 15 mL 150 mM 2,4-dinitrofluorobenzene in MeOH (prepare fresh daily), heat at 100° for 45 min, cool, make up to 250 mL with mobile phase, discard the upper aqueous phase, inject a 20  $\mu$ L aliquot of the lower organic phase.

## HPLC VARIABLES

**Column:** 250  $\times$  4.6 5  $\mu$ m LiChrosorb SI-100

**Mobile phase:** Chloroform:THF:water 42:56.4:1.6

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 350

---

#### CHROMATOGRAM

**Retention time:** 18

---

#### KEY WORDS

normal phase; derivatization

---

#### REFERENCE

Tsuji,K.; Goetz,J.F.; VanMeter,W.; Gusciora,K.A. Normal-phase high-performance liquid chromatographic determination of neomycin sulfate derivatized with 1-fluoro-2,4-dinitrobenzene, *J.Chromatogr.*, **1979**, 175, 141–152.

---

#### SAMPLE

**Matrix:** fermentation solutions

**Sample preparation:** 5 mL Fermentation broth + 5 mL saturated aqueous solution of Tris + 20 mL MeCN, centrifuge at 3000 rpm for 10 min. Remove a 1 mL aliquot of the supernatant and add it to 3 mL 150 mM 2,4-dinitrofluorobenzene in MeOH, heat at 100° under a reflux condenser for 45 min, make up to 4 mL with mobile phase, inject an aliquot.

---

#### HPLC VARIABLES

**Column:** 200 × 4.6 10 µm LiChrosorb RP-8

**Mobile phase:** MeCN:water:acetic acid 55:45:0.15

**Flow rate:** 1.2

**Injection volume:** 20

**Detector:** UV 350

---

#### CHROMATOGRAM

**Retention time:** 11.28 (kanamycin B)

---

#### OTHER SUBSTANCES

**Extracted:** apramycin, tobramycin

---

#### KEY WORDS

derivatization

---

#### REFERENCE

Harangi,J.; Deák,M.; Nánási,P.; Bacsa,G. Determination of the major factors of fermentation of the nebramycin complex by high performance liquid chromatography, *J.Liq.Chromatogr.*, **1984**, 7, 83–93.

---

#### SAMPLE

**Matrix:** formulations

**Sample preparation:** Dilute 3 mg/mL ophthalmic suspension in water with 10 mM sulfuric acid to a tobramycin concentration of 240 µg/mL. Mix 4 mL diluted suspension with 10 mL 10 mg/mL 2,4-dinitrofluorobenzene in EtOH and 10 mL 15 mg/mL tris(hydroxymethyl)aminomethane in water:dimethylsulfoxide 20:80. Heat at 70 ± 2°. for 20 min, allow to cool slightly for 2 min and add 24 mL MeCN. Allow to cool to room temperature, make up to 50 mL with MeCN, inject a 30 µL aliquot.

---

#### HPLC VARIABLES

**Column:** 150 × 3.9 Nova-Pak C18

**Mobile phase:** MeCN:buffer 55:45 (Prepare mobile phase as follows. Dissolve 2.0 g tris (hydroxymethyl)aminomethane in 960 mL water, add 20 mL 0.5 M sulfuric acid and 1200 mL MeCN.)

**Flow rate:** 1.5

**Injection volume:** 30

**Detector:** UV 365

---

**CHROMATOGRAM****Retention time:** 6

---

**OTHER SUBSTANCES****Simultaneous:** neamine, nebramine, tobramycin

---

**KEY WORDS**

derivatization; ophthalmic suspension

---

**REFERENCE**

Russ,H.; McCleary,D.; Katimy,R.; Montana,J.L.; Miller,R.B.; Krishnamoorthy,R.; Davis,C.W. Development and validation of a stability-indicating HPLC method for the determination of tobramycin and its related substances in an ophthalmic suspension, *J.Liq.Chromatogr.Rel.Technol.*, **1998**, *21*, 2165–2181.

---

**SAMPLE****Matrix:** reaction mixtures

**Sample preparation:** 50  $\mu$ L Buffered reaction mixture + 50  $\mu$ L isopropanol + 50  $\mu$ L reagent, heat at 60° for 10 min, centrifuge at 1000 g for 2 min, immediately inject a 50  $\mu$ L aliquot of the supernatant. (Reagent was 80 mM o-phthalaldehyde and 250 mM thioglycolic acid in 1 M boric acid, pH adjusted to 10.4 with 40% KOH.)

---

**HPLC VARIABLES****Column:** 100  $\times$  5 Hypersil ODS

**Mobile phase:** A was MeOH:water:acetic acid 50:45:5 containing 5 g/L heptanesulfonic acid. B was MeOH:water:acetic acid 75:20:5 containing 5 g/L heptanesulfonic acid. A:B 60:40.

**Flow rate:** 2**Injection volume:** 50**Detector:** UV 330

---

**CHROMATOGRAM****Retention time:** 19 (kanamycin A)

---

**KEY WORDS**

derivatization

---

**REFERENCE**

Lovering,A.M.; White,L.O.; Reeves,D.S. Identification of aminoglycoside-acetylating enzymes by high-pressure liquid chromatographic determination of their reaction products, *Antimicrob.Agents Chemother.*, **1984**, *26*, 10–12.

---

**SAMPLE****Matrix:** solutions

**Sample preparation:** Prepare a solution in MeCN, inject a 50  $\mu$ L aliquot.

---

**HPLC VARIABLES****Column:** 250  $\times$  4.6 5  $\mu$ m Ultrasphere octyl

**Mobile phase:** MeCN:20 mM buffer 52:48 (Buffer was 2.68 g  $\text{KH}_2\text{PO}_4$  in 1 L water adjusted to pH 3.0 with phosphoric acid.)

**Column temperature:** 50**Flow rate:** 2**Injection volume:** 50**Detector:** UV 340

---

**CHROMATOGRAM****Retention time:** 12

---

**OTHER SUBSTANCES****Simultaneous:** amikacin

**Noninterfering:** acetaminophen, acetazolamide, N-acetylprocainamide, amobarbital, ampicillin, amitriptyline, caffeine, cefamandole, cefoxime, cefoxitin, cephalothin, clindamycin, chloram-

phenicol, chlordiazepoxide, diazepam, erythromycin, ethosuximide, gentamicin, nitrofurantoin, penicillin G, pentobarbital, phenobarbital, phenytoin, primidone, procainamide, quinidine, salicylic acid, secobarbital, tetracycline, theophylline, tobramycin, vancomycin

## REFERENCE

Kabra,P.M.; Bhatnager,P.K.; Nelson,M.A. Liquid chromatographic determination of amikacin in serum with spectrophotometric detection, *J.Chromatogr.*, **1984**, 307, 224–229.

## SAMPLE

**Matrix:** solutions

**Sample preparation:** 50  $\mu$ L Buffer solution + 25  $\mu$ L 242 mg/mL pH 10.4 Tris buffer + 100  $\mu$ L MeCN:water 50:50 + 30  $\mu$ L 250 mg/mL 2,4,6-trinitrobenzenesulfonic acid in MeCN:water 80:20, vortex for 10 s, heat at 70° for 15 min, add 2 mL chloroform, shake horizontally at 180 cycles/min for 5 min, centrifuge at 750 g for 5 min. Remove the lower organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 200  $\mu$ L MeCN, vortex, inject a 20  $\mu$ L aliquot.

## HPLC VARIABLES

**Column:** 250  $\times$  4.6 5  $\mu$ m Ultrasphere octyl

**Mobile phase:** MeCN:50 mM  $\text{KH}_2\text{PO}_4$  62:38, pH adjusted to 3.5 with phosphoric acid

**Flow rate:** 2.5

**Injection volume:** 20

**Detector:** UV 340

## CHROMATOGRAM

**Retention time:** 6.4

## OTHER SUBSTANCES

**Simultaneous:** tobramycin

## KEY WORDS

derivatization

## REFERENCE

Dash,A.K.; Suryanarayanan,R. A liquid-chromatographic method for the determination of tobramycin, *J.Pharm.Biomed.Anal.*, **1991**, 9, 237–245.

# Ketamine


**Molecular formula:**  $\text{C}_{13}\text{H}_{16}\text{ClNO}$

**Molecular weight:** 237.73

**CAS Registry No.:** 6740-88-1, 1867-66-9 (HCl)

**Merck Index:** 5306

**Lednicer No.:** 1 57



## SAMPLE

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 20  $\mu$ L IS + 350  $\mu$ L 200 mM pH 13 borate buffer, mix. Extract with 5 mL dichloromethane:ethyl acetate 80:20, vortex at 60 rpm for 10 min, centrifuge at 1500g at 15° for 3 min, remove organic phase and extract again with 3 mL dichloromethane:ethyl acetate 80:20. Combine the organic layers and evaporate them to dryness under a stream of nitrogen. Reconstitute the residue in 500  $\mu$ L dichloromethane:ethyl acetate 80:20, extract with 2 mL 2 M HCl. Remove the aqueous layer and evaporate it to dryness at 45°. Reconstitute the residue in 100  $\mu$ L mobile phase, inject a 60  $\mu$ L aliquot.

## HPLC VARIABLES

**Column:** 125  $\times$  4 5  $\mu$ m Purospher RP-18e (Merck)

**Mobile phase:** MeCN:30 mM pH 7.2 phosphate buffer 23:77

**Column temperature:** 20

**Flow rate:** 1.5

**Injection volume:** 60

**Detector:** UV 210

---

#### CHROMATOGRAM

**Retention time:** 16.26

**Internal standard:** nortilidine (6.29)

**Limit of detection:** 3 ng

**Limit of quantitation:** 5 µg/L

---

#### OTHER SUBSTANCES

**Extracted:** metabolites

**Simultaneous:** atropine, buprenorphine, diazepam, dopamide, furosemide, nalbuphine, omeprazole, phenobarbital, phytomenadione, propofol

---

#### KEY WORDS

plasma; pharmacokinetics

---

#### REFERENCE

Bolze,S.; Bouliou,R. HPLC determination of ketamine, norketamine, and dehydronorketamine in plasma with a high-purity reversed-phase sorbent, *Clin.Chem.*, **1998**, *44*, 560–564.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** Activate two 130 mg Sep-Pak Light C18 SPE cartridges with 1 mL MeOH and 1 mL water. Add 1 mL plasma onto SPE cartridge, wash with 1 mL water, 2 mL 5 mM pH 9.6 ammonium sulfate buffer containing 3% MeCN, and 1 mL 5 mM pH 9.6 ammonium sulfate buffer containing 20% MeCN. Displace washing solution with 200 µL 20 mM pH 2.1 phosphoric acid buffer containing 25% MeCN and elute with 500 µL of the same solution. Mix eluate with 1 mL 40 mM pH 11.5 sodium hydroxide. Add the mixture onto the second SPE cartridge, wash with 2 mL 5 mM pH 9.6 ammonium sulfate buffer containing 3% MeCN and with 1 mL 5 mM pH 9.6 ammonium sulfate buffer containing 20% MeCN. Displace washing solution with 200 µL 20 mM pH 2.1 phosphoric acid buffer containing 25% MeCN and elute with 500 µL of the same solution. Evaporate the eluate in a vacuum centrifuge to about 150 µL, add 12 µL 1 mM sodium hydroxide immediately before injection, inject the whole volume on the column. (Buffers were adjusted with ammonia. SPE flow-rate at all steps was approximately 1.5 mL/min.)

---

#### HPLC VARIABLES

**Column:** 150 × 4.0 5 µm Chiral AGP (Chrom Tech, Sweden)

**Mobile phase:** MeOH:10 mM pH 7.0 KH<sub>2</sub>PO<sub>4</sub> 16:84 (Buffer was adjusted with potassium hydroxide.)

**Column temperature:** 40

**Flow rate:** 1

**Injection volume:** 162

**Detector:** UV 220

---

#### CHROMATOGRAM

**Retention time:** 12 (S), 14 (R)

**Limit of detection:** 1.7 ng/mL (S), 2.0 ng/mL (R)

**Limit of quantitation:** 10 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** metabolite

---

#### KEY WORDS

plasma; SPE; chiral

---

#### REFERENCE

Svensson,J.O.; Gustafsson,L.L. Determination of ketamine and norketamine enantiomers in plasma by solid-phase extraction and high-performance liquid chromatography, *J.Chromatogr.B*, **1996**, *678*, 373–376.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 100  $\mu$ L bisnortilidine in mobile phase + 100  $\mu$ L 3 M NaOH, mix, add 5.75 mL ice-cold cyclohexane, agitate at 4° for 15 min, repeat extraction with 3 mL ice-cold cyclohexane. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 500  $\mu$ L cyclohexane, evaporate to dryness under a stream of nitrogen, reconstitute in 0.2-1 mL mobile phase, mix for 30 s, let stand for 2 min, inject a 50-100  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 100  $\times$  4 AGP (Grom)

**Mobile phase:** Isopropanol:20 mM pH 7 phosphate buffer 2.5:97.5

**Column temperature:** 25

**Flow rate:** 0.5

**Injection volume:** 50-100

**Detector:** UV 215

---

**CHROMATOGRAM**

**Retention time:** 22 (S), 26 (R)

**Internal standard:** bisnortilidine (ethyl trans-2-amino-1-phenyl-3-cyclohexene-1-carboxylate hydrochloride) (one enantiomer only, separated by HPLC) (15)

**Limit of detection:** 20 ng/mL

**Limit of quantitation:** 40 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** metabolites

---

**KEY WORDS**

chiral; plasma; pharmacokinetics

---

**REFERENCE**

Geisslinger, G.; Menzel-Soglowek, S.; Kamp, H.-D.; Brune, K. Stereoselective high-performance liquid chromatographic determination of the enantiomers of ketamine and norketamine in plasma, *J. Chromatogr.*, **1991**, *568*, 165-176.

---

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Condition a Sep-Pak C18 SPE cartridge with water, MeOH, and 100 mM ammonium acetate. Add 200  $\mu$ L plasma to the SPE cartridge, wash with 100 mM ammonium acetate, elute with MeOH:100 mM ammonium acetate 3:1. Evaporate the eluate to dryness under reduced pressure, dissolve the residue in 200  $\mu$ L mobile phase, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 150  $\times$  4.6 Hitachi gel 3056 octadecylsilica

**Mobile phase:** MeOH:100 mM ammonium acetate 60:40

**Flow rate:** 1

**Injection volume:** 20

**Detector:** MS, Hitachi M1000, APCI, nebulizer 260°, vaporizer 399

---

**CHROMATOGRAM**

**Retention time:** 6.2

**Limit of detection:** 0.5-2.5 ng/mL

---

**OTHER SUBSTANCES**

**Simultaneous:** atipamezole, atropine, butorphanol, flumazenil, medetomidine, midazolam, xylazine

---

**KEY WORDS**

plasma; SPE; dog

---

---

**REFERENCE**

Kanazawa,H.; Nagata,Y.; Matsushima,Y.; Takai,N.; Uchiyama,H.; Nishimura,R.; Takeuchi,A. Liquid chromatography-mass spectrometry for the determination of medetomidine and other anaesthetics in plasma, *J.Chromatogr.*, **1993**, 631, 215–220.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 500  $\mu$ L Serum + 500  $\mu$ L MeCN:water 85% phosphoric acid 20:78:2, vortex for 10–15 s, filter (Amicon Centricon-10 microseparation system, 10000 molecular mass cut-off) while centrifuging at 3000 g for 30 min, inject a 20–100  $\mu$ L aliquot of the ultrafiltrate.

---

**HPLC VARIABLES**

**Column:** 100  $\times$  4.6 3  $\mu$ m Spherisorb phenyl

**Mobile phase:** MeCN:MeOH:10 mM  $\text{NaH}_2\text{PO}_4$ :85% phosphoric acid 10:30:59.8:0.2

**Column temperature:** 50

**Injection volume:** 20–100

**Detector:** UV 215

---

**CHROMATOGRAM**

**Retention time:** 5.1

**Limit of detection:** 5 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** metabolites

---

**KEY WORDS**

horse; serum; ultrafiltrate

---

**REFERENCE**

Seay,S.S.; Aucoin,D.P.; Tyczkowska,K.L. Rapid high-performance liquid chromatographic method for the determination of ketamine and its metabolite dehydronorketamine in equine serum, *J.Chromatogr.*, **1993**, 620, 281–287.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100  $\mu$ L mobile phase, centrifuge at 2800 g for 5 min, inject a 50  $\mu$ L aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

---

**HPLC VARIABLES**

**Column:** 300  $\times$  3.9 4  $\mu$ m NovaPack C18

**Mobile phase:** MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic))  $\text{KH}_2\text{PO}_4$  adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

**Column temperature:** 30

**Flow rate:** 0.8

**Injection volume:** 50

**Detector:** UV 269

---

**CHROMATOGRAM**

**Retention time:** 4.19

**Limit of detection:** <120 ng/mL

---

**KEY WORDS**

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfapyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vindesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrodine; phenylbutazone; demexiptiline; clozapine; progauil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thiopropazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

**REFERENCE**

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254–262.

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18

**Column:** 250 × 4.6 5 µm Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10–30

**Detector:** UV 202.8



---

**CHROMATOGRAM****Retention time:** 9.637

---

**KEY WORDS**whole blood

---

**REFERENCE**

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

---

**SAMPLE****Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 5 µm YMC GEL, ODS-AM coated with poly-(R)-1-(α-naphthyl)ethyl methacrylamide (Prepare (R)-1-(α-naphthyl)ethyl methacrylamide by reacting methacryl chloride with (R)-1-(α-naphthyl)ethylamine. Prepare poly-(R)-1-(α-naphthyl)ethyl methacrylamide by polymerizing this compound in anhydrous benzene/THF with 2,2'-azobis(isobutyronitrile)(Caution! Benzene is a carcinogen!). Average molecular weight = 2500. Coat 4 g 5 µm YMC GEL, ODS-AM with 0.8 g of this polymer using dichloromethane as a solvent.)

**Mobile phase:** MeCN:0.5M sodium perchlorate 40:60**Flow rate:** 1

---

**CHROMATOGRAM****Retention time:** k' 3.18 (α = 1.21)

---

**OTHER SUBSTANCES****Also analyzed:** propranolol

---

**KEY WORDS**chiral

---

**REFERENCE**

Oi,N.; Hashimoto,S.; Ishizuka,N.; Ohtake,J. Enantiomer separation with poly-(R)-1-(α-naphthyl)-ethyl-methacrylamide coated on ODS silica gel by reversed phase HPLC, *Biomed.Chromatogr.*, **1997**, 11, 296–297.

---

**SAMPLE****Matrix:** solutions

---

**HPLC VARIABLES****Column:** 250 × 4.6 Zorbax RX

**Mobile phase:** Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

**Column temperature:** 30**Flow rate:** 2**Detector:** UV 210

---

**OTHER SUBSTANCES**

**Also analyzed:** acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepan, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone,

debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenoprofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-stilbene, imipramine, indomethacin, isocarbostyryl, isocarboxazid, isoniazid, isoproterenol, iso-xsuprine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, lox-apine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, meth-aqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methylpyrrolone, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, na-phazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitra-zepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbi-tal, persantine, phenacetin, phenazocine, phenazopyridine, phenacyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenyl-butazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primi-done, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopola-mine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sul-faethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sul-fasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tol-metin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapa-mil, vincamine, warfarin, yohimbine, zoxazolamine

---

## REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233–242.

---

## SAMPLE

**Matrix:** solutions

---

## HPLC VARIABLES

**Column:** 250 × 4.6 Sumchiral CSP 10 (Sumika Chemical Analysis Service)

**Mobile phase:** n-Hexane:EtOH:trifluoroacetic acid 200:40:0.6

**Flow rate:** 1

**Detector:** UV 230-280

---

## CHROMATOGRAM

**Retention time:** k' 7.78 (first enantiomer)

---

## KEY WORDS

chiral;  $\alpha = 1.12$

---

## REFERENCE

Oi,N.; Kitahara,H.; Aoki,F. Direct enantiomer separations by high-performance liquid chromatography with chiral urea derivatives as stationary phases, *J.Chromatogr.A*, **1995**, *694*, 129–134.

---

## SAMPLE

**Matrix:** solutions

**HPLC VARIABLES**

**Column:** 250 × 4.6 CSP-4 (Prepare as follows. Add a solution of 1.07 g L-valyl-L-valyl-L-valine isopropylester (Bunseki Kagaku 1079, 28, 125) in 30 mL dry dioxane (Caution! Dioxane is a carcinogen!) dropwise to a mixture of 2.2 g 2,4,6-trichloro-1,3,5-triazine (cyanuric chloride) in 20 mL dry dioxane stirred at 0°, add 3 g anhydrous sodium carbonate at room temperature, stir, filter, evaporate to give a colorless solid. Dissolve 8.3 g of this solid in 30 mL dry dioxane, add 2 g N-(2-aminoethyl)-3-aminopropyltrimethoxysilane, add 1.5 g anhydrous sodium carbonate, reflux with stirring for 40 h, filter, add 3 g dried 10 µm LiChrosorb Si 100, reflux with slow stirring for 10 h, cool, filter. Wash the solid with dioxane, MeOH, and diethyl ether, dry under reduced pressure (J.Chromatogr. 1984, 292, 427).)

**Mobile phase:** Hexane:EtOH:trifluoroacetic acid 96:4:0.24

**Detector:** UV

**CHROMATOGRAM**

**Retention time:** k' 7.85 (first enantiomer)

**KEY WORDS**

chiral;  $\alpha = 1.09$

**REFERENCE**

Oi,N.; Kitahara,H.; Matsushita,Y.; Kisu,N. Enantiomer separation by gas and high-performance liquid chromatography with tripeptide derivatives as chiral stationary phases, *J.Chromatogr.A*, **1996**, 722, 229–232.

# Ketanserin


**Molecular formula:** C<sub>22</sub>H<sub>22</sub>FN<sub>3</sub>O<sub>3</sub>

**Molecular weight:** 395.43

**CAS Registry No.:** 74050-98-9

**Merck Index:** 5307

**Lednicer No.:** 3 193

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 25 µL Serum + 50 µL MeOH, vortex for 30 s, centrifuge at 13600 g for 5 min, inject a 25 µL aliquot of the supernatant.

**HPLC VARIABLES**

**Column:** 150 × 3.9 4 µm Nova-Pak C18

**Mobile phase:** MeCN:buffer:water 31:50:19 (Buffer was 2% acetic acid adjusted to pH 7.0 with ammonium hydroxide.)

**Flow rate:** 1

**Injection volume:** 25

**Detector:** F ex 225 em no emission filter

**CHROMATOGRAM**

**Retention time:** 7.7

**Limit of quantitation:** 40 ng/mL

**OTHER SUBSTANCES**

**Extracted:** metabolites

**KEY WORDS**

rat; serum; pharmacokinetics

**REFERENCE**

Wong,Y.W.; Skinner,M.H. Rapid method for the determination of ketanserin in rat serum by high-performance liquid chromatography with fluorimetric detection, *J.Chromatogr.*, **1991**, 571, 318–323.

---

**SAMPLE****Matrix:** solutions**Sample preparation:** Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.

---

**HPLC VARIABLES****Column:** 125 × 4.9 Spherisorb S5W silica**Mobile phase:** MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7**Flow rate:** 2**Injection volume:** 20**Detector:** E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

---

**CHROMATOGRAM****Retention time:** 1.4

---

**OTHER SUBSTANCES**

**Also analyzed:** acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazepine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipranone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypropazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindine, pimozide, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

---

**REFERENCE**

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, *323*, 191-225.

# Ketoconazole


**Molecular formula:**  $C_{26}H_{28}Cl_2N_4O_4$

**Molecular weight:** 531.44

**CAS Registry No.:** 65277-42-1

**Merck Index:** 5313

**Lednicer No.:** 3 132



## SAMPLE

**Matrix:** blood

**Sample preparation:** 500  $\mu$ L Plasma + 2  $\mu$ g clotrimazole + hexane:isoamyl alcohol 98.5:1.5, vortex, centrifuge. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 250  $\mu$ L MeCN, inject an aliquot.

## HPLC VARIABLES

**Column:** 150  $\times$  3.9 NovaPak C18

**Mobile phase:** MeCN:MeOH:50 mM phosphate buffer 40:5:55

**Detector:** UV 220

## CHROMATOGRAM

**Limit of detection:** 100-200 ng/mL

## KEY WORDS

plasma; pharmacokinetics

## REFERENCE

von Moltke, L.L.; Greenblatt, D.J.; Harmatz, J.S.; Duan, S.X.; Harrel, L.M.; Cotreau-Bibbo, M.M.; Pritchard, G.A.; Wright, C.E.; Shader, R.I. Triazolam biotransformation by human liver microsomes in vitro: Effects of metabolic inhibitors and clinical confirmation of a predicted interaction with ketoconazole, *J.Pharmacol.Exp.Ther.*, **1996**, 276, 370-379.

## SAMPLE

**Matrix:** blood, microsomal incubations

**Sample preparation:** Vortex 1 mL plasma or microsomal incubation with 200  $\mu$ L 5  $\mu$ g/mL diazepam and 100  $\mu$ L 5 M NaOH solution for 10 s, add 5 mL butan-1-ol:hexane 2:98, vortex for 1 min, centrifuge at 2000 g and 4° for 5 min, evaporate the organic phase to dryness at 40° using a vacuum vortex evaporator, reconstitute the residue in 200  $\mu$ L mobile phase, inject a 50  $\mu$ L aliquot.

## HPLC VARIABLES

**Column:** 5  $\mu$ m Nova-Pak C18

**Mobile phase:** MeCN:buffer 35:65 (Buffer was water containing 1% triethylamine, adjusted to pH 3 with orthophosphoric acid.)

**Flow rate:** 2

**Injection volume:** 50

**Detector:** UV 240

## CHROMATOGRAM

**Retention time:** 6.3

**Internal standard:** diazepam

## OTHER SUBSTANCES

**Extracted:** amitriptyline, nortriptyline

**Noninterfering:** furafylline, hydroxyamitriptyline, hydroxynortriptyline, quinidine, mephentoin, triacetyloleandomycin

## KEY WORDS

human; liver; rat; plasma

---

**REFERENCE**

Ghahramani,P.; Lennard,M.S. Quantitative analysis of amitriptyline and nortriptyline in human plasma and liver microsomal preparations by high-performance liquid chromatography, *J.Chromatogr.B*, **1996**, 685, 307–313.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50  $\mu$ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood)  $\mu$ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18

**Column:** 250  $\times$  4.6 5  $\mu$ m Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10–30

**Detector:** UV 202.8

---

**CHROMATOGRAM**

**Retention time:** 15.738

---

**KEY WORDS**

whole blood

---

**REFERENCE**

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 125  $\times$  4 5  $\mu$ m Hibar RP 18 (Merck)

**Mobile phase:** MeCN:water:diethylamine 48:55:0.02

**Flow rate:** 1.2

**Injection volume:** 20, 100


**Detector:** UV 254

---

**REFERENCE**

Galia,E.; Nicolaidis,E.; Hörter,D.; Löbenberg,R.; Reppas,C.; Dressman,J.B. Evaluation of various dissolution media for predicting in vivo performance of class I and II drugs, *Pharm.Res.*, **1998**, 15, 698–705.

# Ketoprofen



**Molecular formula:** C<sub>16</sub>H<sub>14</sub>O<sub>3</sub>

**Molecular weight:** 254.29

**CAS Registry No.:** 22071-15-4

**Merck Index:** 5316

**Lednicer No.:** 2 64

## SAMPLE

**Matrix:** blood

**Sample preparation:** Mix 50 µL plasma with 150 µL 25 µg/mL ibufenac in MeCN:water 95:5, vortex for 30 s, centrifuge for 1 min using a Beckman Microfuge (Beckman Instruments, Palo Alto, CA). Inject a 20 g/mL aliquot of the supernatant.

## HPLC VARIABLES

**Column:** 33 × 4.6 3 µm C18 (Perkin Elmer, Norwalk, CT)

**Mobile phase:** MeCN:buffer 40:60 (Prepare mobile phase as follows. Dissolve 4 mL concentrated phosphoric acid in 600 mL water and mix with 400 mL MeCN.)

**Flow rate:** 1.5

**Injection volume:** 20

**Detector:** UV 220

## CHROMATOGRAM

**Retention time:** 1.1

**Internal standard:** ibufenac (2.6)

**Limit of detection:** 1 µg/mL

**Limit of quantitation:** 5 µg/mL

## OTHER SUBSTANCES

**Simultaneous:** dicloxacillin, ibuprofen

**Noninterfering:** acetaminophen, N-acetylprocainamide, amikacin, amitriptyline, aspirin, aztreonam, barbituric acid, brompheniramine, cafazolin, caffeine, carbamazepine, carbamazepine epoxide, cephalixin, chlorpheniramine, clonazepam, clotrimazole, desipramine, desmethyldoxepin, digitoxin, digoxin, disopyramide, doxepin, ethosuximide, felbamate, gentamicin, imipenem, imipramine, lidocaine, maprotiline, mephentyoin, mephobarbital, metharbital, methsuximide, methylsuccinimide, nortriptyline, paramethadione, phenacemide, phenobarbital, phensuximide, phenylpropanolamine, phenytoin, primidone, protriptyline, sulfamethoxazole, theophylline, tobramycin, trimethadione, trimethoprim, vancomycin.

## KEY WORDS

plasma

## REFERENCE

Rifai,N.; Lafi,M.; Sakamoto,M.; Law,T. Measurement of plasma ketoprofen by a rapid high-performance liquid chromatography assay, *Ther.Drug Monit.*, **1997**, 19, 175–178.

## SAMPLE

**Matrix:** blood

**Sample preparation:** Condition a 1 mL 100 mg C 18 SPE cartridge (Varian) with 1 mL MeOH and 1 mL 100 mM pH 2.5 phosphate buffer. Mix 100 µL 50 µg/mL 3,4-dimethoxybenzoic acid in MeOH with 1 mL plasma. Vortex with 360 mg ammonium sulfate for 30 min to deproteinate the plasma. Centrifuge at 8000 g at 4°. for 30 min, acidify the supernatant with 3 ml 100 mM sulfuric acid. Add the supernatant to the SPE cartridge, wash twice with 750 µL MeOH:100 mM pH 2.5 phosphate buffer 20:80. Elute with two 500 µL portions of MeOH:50 mM pH 7.4 phosphate buffer 75:25. Inject a 40 µL aliquot.

## HPLC VARIABLES

**Guard column:** Supelco C18

**Column:** 250 × 4.6 Chirex 3005 [(R)-1-naphtylglycine and 3,5-dinitrobenzoic acid] (Phenomenex, Torrance, CA, USA)  
**Mobile phase:** 20 mM ammonium acetate in MeOH  
**Flow rate:** 1.2  
**Injection volume:** 40  
**Detector:** UV 254

---

#### CHROMATOGRAM

**Retention time:** 14.2 (R-(-)), 16.1 (S-(+))  
**Internal standard:** 3,4-dimethoxybenzoic acid (18.1)  
**Limit of quantitation:** 160 ng/mL

---

#### KEY WORDS

SPE; chiral; pharmacokinetics

---

#### REFERENCE

Boisvert, J.; Caillé, G.; McGilveray, I. J.; Qureshi, S. A. Quantification of ketoprofen enantiomers in human plasma based on solid-phase extraction and enantioselective column chromatography, *J. Chromatogr. B*, **1997**, 690, 189–193.

---

#### SAMPLE

**Matrix:** blood  
**Sample preparation:** 500 µL Plasma + 1 mL 1 M pH 7.0 phosphate buffer + IS, extract with 7 mL diethyl ether, evaporate, reconstitute the residue in 250 µL MeCN:5 mM pH 7.0 phosphate buffer 10:90, inject 50 µL aliquot.

---

#### HPLC VARIABLES

**Guard column:** 4 × 4 LiChrocart  
**Column:** 125 × 3.5 µm Ecocart packed with LiChrospher 100 RP-18  
**Mobile phase:** MeCN:50 mM pH 7.0 phosphate buffer 16:84  
**Column temperature:** 35  
**Flow rate:** 0.6  
**Injection volume:** 50  
**Detector:** UV 260

---

#### CHROMATOGRAM

**Retention time:** 14.88  
**Internal standard:** naproxen (10.85)  
**Limit of detection:** 30 pg/mL  
**Limit of quantitation:** 100 pg/mL

---

#### KEY WORDS

horse; plasma

---

#### REFERENCE

Baeyens, W. R. G.; Van Der Weken, G.; Van Overbeke, A.; Corveleyn, S.; Remon, J. P.; Deprez, P. Comparative narrow-bore high-performance liquid chromatographic determination of ketoprofen in horse plasma, *Bio-med. Chromatogr.*, **1998**, 12, 167–169.

---

#### SAMPLE

**Matrix:** blood  
**Sample preparation:** Acidify plasma with an equal volume of 100 mM pH 3.2 phosphoric acid, inject a 20 µL aliquot onto column A and elute to waste with mobile phase A. After 5 min elute the contents of column A onto column B with mobile phase B, monitor the effluent from column B. (Before each run wash column A with MeCN:100 mM pH 3.2 phosphate buffer 20:80 for 3 min, with water for 1 min, with MeOH for 3 min, and with mobile phase A. Equilibrate column B with mobile phase B.)

---

#### HPLC VARIABLES

**Column:** A 10 × 4.5 µm Develosil NH<sub>2</sub>-5 (Nomura Chemical Co., Japan) + 10 × 4 Nucleosil 5CN (Macherey-Nagel, Düren, Germany); B 150 × 4.6 Ultron ES-OVM G bonded silica column (Shinwa Kako Co., Japan)



**Mobile phase:** A 100 mM pH 3.2 phosphate buffer; B MeOH:100 mM pH 3.2 phosphate buffer 20:80  
**Column temperature:** 20  
**Flow rate:** 0.8  
**Injection volume:** 20  
**Detector:** UV 265

---

#### CHROMATOGRAM

**Retention time:** 9.6 (S(+)), 12.0 (R(-))

---

#### KEY WORDS

column-switching; chiral; plasma

---

#### REFERENCE

Tamai,G.; Edani,M.; Imai,H. Determination of ketoprofen enantiomers in plasma by solid phase extraction and column switching high performance liquid chromatography, *Anal.Sci.*, **1991**, 7, 29–32.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 100  $\mu$ L tolmetin solution + 500  $\mu$ L pH 1.8 phosphate buffer, extract with 1-butanol/MTBE. Remove the organic layer and add it to 500  $\mu$ L pH 6.1 ammonium acetate buffer, mix, inject an aliquot of the aqueous layer.

---

#### HPLC VARIABLES

**Column:** 150  $\times$  4.6 5  $\mu$ m Cosmosil C18

**Mobile phase:** MeCN:250 mM pH 5.0 ammonium acetate buffer 20:80

**Flow rate:** 1.8

**Detector:** UV 350

---

#### CHROMATOGRAM

**Internal standard:** tolmetin (UV 258)

**Limit of quantitation:** 5 ng/mL

---

#### KEY WORDS

plasma; pharmacokinetics

---

#### REFERENCE

Shah,A.K.; Wei,G.; Lanman,R.C.; Bhargava,V.O.; Weir,S.J. Percutaneous absorption of ketoprofen from different anatomical sites in man, *Pharm.Res.*, **1996**, 13, 168–172.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** Condition a 1 mL Bond-Elut C18 SPE cartridge with 1 mL MeOH. 1 mL Serum + 1 mL water + 20  $\mu$ L saturated ammonium sulfate solution + 60  $\mu$ L concentrated HCl, vortex for 3 min, add to the SPE cartridge, wash with three 1 mL portions of water, allow to dry for 3 min, elute with five 500  $\mu$ L portions of MeOH. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute the residue in 1 mL mobile phase, inject a 100  $\mu$ L aliquot.

---

#### HPLC VARIABLES

**Column:** 100  $\times$  4.6 5  $\mu$ m Spheri-5 cyano

**Mobile phase:** MeCN:MeOH:water:phosphoric acid 21:22:56.5:0.5

**Flow rate:** 0.5

**Injection volume:** 100

**Detector:** F ex 248 em 335 (filter) following post-column reaction. The column effluent flowed through a knitted 7.9 m  $\times$  0.3 mm ID PTFE coil irradiated with an SC3-9 UV lamp (UVP, San Gabriel CA) and cooled with a fan to the detector.

---

#### CHROMATOGRAM

**Retention time:** 6.5

**Internal standard:** ketoprofen

---

---

**OTHER SUBSTANCES**

Extracted: fenbufen

---

**KEY WORDS**

ketoprofen is IS; post-column reaction; post-column photochemical derivatization; serum; SPE

---

**REFERENCE**

Siluveru,M.; Stewart,J.T. Determination of fenbufen and its metabolites in serum by reversed-phase high-performance liquid chromatography using solid-phase extraction and on-line post-column ultraviolet irradiation and fluorescence detection, *J.Chromatogr.B*, **1996**, 682, 89–94.

---

**SAMPLE**

Matrix: blood

Sample preparation: 100  $\mu$ L Serum + 200  $\mu$ L MeCN, vortex for 30 s, centrifuge at 14000 g for 30 s, inject a 20  $\mu$ L aliquot of the supernatant.

---

**HPLC VARIABLES**

Column: 150  $\times$  4.6 5  $\mu$ m Econosphere CN

Mobile phase: MeCN:water:phosphoric acid 4:100:0.02

Flow rate: 1.5

Injection volume: 20

Detector: UV 254

---

**CHROMATOGRAM**

Retention time: 7

Limit of detection: 100 ng/mL

---

**KEY WORDS**

comparison with capillary electrophoresis; serum

---

**REFERENCE**

Friedberg,M.; Shihabi,Z.K. Ketoprofen analysis in serum by capillary electrophoresis, *J.Chromatogr.B*, **1997**, 695, 193–198.

---

**SAMPLE**

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50  $\mu$ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood)  $\mu$ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

Guard column: 20 mm long Symmetry C18

Column: 250  $\times$  4.6 5  $\mu$ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10–30

Detector: UV 200.5

---

**CHROMATOGRAM**

Retention time: 19.628

---

**KEY WORDS**whole blood

---

**REFERENCE**

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149–163.

---

**SAMPLE**

**Matrix:** cell suspensions

**Sample preparation:** Centrifuge cell suspension at 2000 g for 4 min. Remove a 2 mL aliquot of the supernatant and add it to 200  $\mu$ L 4 mg/mL IS in DMF, mix, add 200  $\mu$ L 5 M HCl, extract twice with 3 mL portions of toluene. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 1 mL dichloromethane, add 20  $\mu$ L 40 mg/mL 1-hydroxybenzotriazole in dichloromethane:pyridine 99:1, add 300  $\mu$ L 40 mg/mL 1-(3-dimethylaminopropyl)-3-ethylcarbodiimide hydrochloride in dichloromethane, add 300  $\mu$ L 40 mg/mL (-)-(*S*)- $\alpha$ -methylbenzylamine in dichloromethane, let stand for 30 min, evaporate to dryness, reconstitute with 500  $\mu$ L mobile phase, inject a 10  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Guard column:** 10 mm long Techsphere ODS (HPLC Technology, Macclesfield UK)

**Column:** 250  $\times$  5  $\mu$ m Techsphere ODS (HPLC Technology, Macclesfield UK)

**Mobile phase:** MeCN:7.5 mM NaH<sub>2</sub>PO<sub>4</sub> 50:50, containing 5 mM sodium pentanesulfonate, pH adjusted to 2.8 with phosphoric acid

**Flow rate:** 1

**Injection volume:** 10

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** *k'* 6.60, 7.55 (enantiomers)

**Internal standard:** phenylacetic acid (*k'* 2.50)

**Limit of detection:** 1  $\mu$ g/mL

---

**KEY WORDS**derivatization; chiral

---

**REFERENCE**

Thomason, M.J.; Hung, Y.-F.; Rhys-Williams, W.; Hanlon, G.W.; Lloyd, A.W. Indirect enantiomeric separation of 2-arylpropionic acids and structurally related compounds by reversed phase HPLC, *J. Pharm. Biomed. Anal.*, **1997**, 15, 1765–1774.

---

**SAMPLE**

**Matrix:** microsomal incubations

**Sample preparation:** 250  $\mu$ L Microsomal incubation, 150  $\mu$ L ice-cold MeCN and 2.5 ng ketoprofen, centrifuge. Extract the mixture with 4 mL ethyl acetate, centrifuge at 3000 rpm for 10 min, remove the organic fraction and evaporate it under a gentle stream of nitrogen at 40°. Dissolve the residue in 30  $\mu$ L MeOH and dilute to 60  $\mu$ L with water. Inject a 30  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6  $\mu$ m RP-18 (Kanto Chemical, Tokyo)

**Mobile phase:** MeCN:water 40:60 containing 0.6% acetic acid

**Column temperature:** 35

**Flow rate:** 1

**Injection volume:** 30

**Detector:** UV 260

---

**CHROMATOGRAM**

**Retention time:** 18.0

---

**OTHER SUBSTANCES**

**Extracted:** metabolites, indomethacin

---

**KEY WORDS**

liver; pharmacokinetics; ketoprofen is IS

---

**REFERENCE**

Nakajima,M.; Inoue,T.; Shimada,N.; Tokudome,S.; Yamamoto,T.; Kuroiwa,Y. Cytochrome P450 2C9 catalyzes indomethacin O-demethylation in human liver microsomes, *Drug Metab.Dispos.*, **1998**, 26, 261–266.

---

**SAMPLE**

**Matrix:** permeate

---

**HPLC VARIABLES**

**Column:** 250 × 4 ODS (Hitachi)

**Mobile phase:** MeCN:50 mM phosphoric acid 40:60 adjusted to pH 5.5 with NaOH

**Column temperature:** 55

**Flow rate:** 0.6

**Injection volume:** 20

**Detector:** UV 260

---

**OTHER SUBSTANCES**

**Also analyzed:** carbamazepine, fenbufen, indomethacin,  $\alpha$ -naphthoquinone, naproxen, tolmetin

---

**REFERENCE**

Sugawara,M.; Takekuma,Y; Yamada,H.; Kobayashi,M.; Iseki,K.; Miyazaki,K. A general approach for the prediction of the intestinal absorption of drugs: regression analysis using the physicochemical properties and drug-membrane electrostatic interactions, *J.Pharm.Sci.*, **1998**, 87, 960–966.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Mix 50  $\mu$ L of a 0.001–5 mM solution in MeCN with 50  $\mu$ L 1 mM DNS-APy in MeCN containing 50 mM 2,2'-dipyridyl disulfide and 50 mM triphenylphosphine, let stand at room temperature for 30 min. Remove a 10  $\mu$ L aliquot and dilute it to 100  $\mu$ L with MeCN, inject a 2  $\mu$ L aliquot. (Synthesis of DNS-APy, 1-(5-dimethylamino-1-naphthalenesulfonyl)-(S)-3-aminopyrrolidine, is as follows. Cool a solution of 16.4 g (R)-(+)-1-benzyl-3-pyrrolidinol in 164 mL pyridine to +5°, add 19.35 g p-toluenesulfonyl chloride, stir at +10° for 48 h, evaporate to dryness, chromatograph using dichloromethane:acetone 95:5 to obtain (3R)-3-[(4-tolylsulfonyl)oxy]-1-(phenylmethyl)pyrrolidine (mp 68°). Heat a solution of (3R)-3-[(4-tolylsulfonyl)oxy]-1-(phenylmethyl)pyrrolidine in 200 mL anhydrous DMF to 65°, add 33.5 g sodium azide (Caution! Sodium azide is highly toxic!), stir at 60° for 7 h, filter, evaporate the filtrate to dryness under reduced pressure, dissolve the residue in ethyl acetate, wash twice with water, dry over anhydrous magnesium sulfate, evaporate to obtain (3S)-3-azido-1-(phenylmethyl)pyrrolidine as an oil. Add 3.5 g 10% palladium on carbon under nitrogen to a solution of 7.05 g (3S)-3-azido-1-(phenylmethyl)pyrrolidine in 34.8 mL 1 M HCl in water and 245 mL EtOH, hydrogenate at atmospheric pressure for 30 min, add 3.5 g catalyst, hydrogenate for 2 h, filter, add 34.8 mL 1 M HCl to the filtrate, evaporate to dryness under reduced pressure, take up the residue in 70 mL EtOH, filter, evaporate the filtrate to dryness under reduced pressure, repeat this operation twice, crystallize with the minimum amount of EtOH to obtain (3S)-(+)-3-aminopyrrolidine dihydrochloride (J. Med. Chem. 1992, 35, 4205). (3S)-(+)-Aminopyrrolidine dihydrochloride is also reported to be available from Tokyo Kasei. Stir 800 mg (3S)-(+)-3-aminopyrrolidine dihydrochloride and 2 mL triethylamine in 800 mL MeCN at 0–10°, add a solution of 440 mg dansyl chloride in 80 mL MeCN dropwise, stir in the dark for 30 min, evaporate to dryness under reduced pressure, dissolve the residue in 200 mL 5% HCl, wash twice with 40 mL portions of dichloromethane. Adjust the pH of the organic layer to 13–14 with 5% NaOH, extract twice with 10 mL portions of dichloromethane. Combine the organic layers and wash them with 80 mL water. Dry the organic layer over anhydrous sodium sulfate, evaporate to dryness under reduced pressure, dissolve the residue in dichloromethane:MeOH 90:10, chromatograph on silica gel with dichloromethane:MeOH 90:10. Collect the greenish-yellow fluorescent band and evaporate it under reduced pressure to obtain DNS-APy as a greenish-yellow oil.)

---

**HPLC VARIABLES**

**Column:** 150 × 4.6 5  $\mu$ m TSK gel ODS-80TM (Tosoh)

**Mobile phase:** MeCN:water 50:50

**Flow rate:** 1  
**Injection volume:** 2  
**Detector:** F ex 340 em 530

---

**CHROMATOGRAM**

**Retention time:** 37 ((S)-(+)), 40 ((R)-(-))  
**Limit of detection:** 0.1 pmole

---

**OTHER SUBSTANCES**

**Simultaneous:** pranoprofen

---

**KEY WORDS**

derivatization; chiral

---

**REFERENCE**

Al-Kindy,S.; Santa,T.; Fukushima,T.; Homma,H.; Imai,K. 1-(5-Dimethylamino-1-naphthalenesulphonyl)-(S)-3-aminopyrrolidine (DNS-Apy) as a fluorescence chiral labelling reagent for carboxylic acid enantiomers, *Biomed.Chromatogr.*, **1997**, 11, 137–142.

---

**SAMPLE**

**Matrix:** urine

**Sample preparation:** Condition a Sep-Pak C18 SPE cartridge with 10 mL ethyl acetate and dry by aspiration of air. Evaporate an aliquot of 100  $\mu$ L 20  $\mu$ g/mL IS in MeOH to dryness at 37°. Add 1 mL urine, vortex, add 250  $\mu$ L 1 M pH 5.0 acetate buffer, vortex. Add 250  $\mu$ L of the mixture to the SPE cartridge, dry by aspiration of air, elute with 3 mL ethyl acetate, evaporate the eluate to dryness under a stream of nitrogen at 37°, reconstitute the residue in 200  $\mu$ L mobile phase, inject a 10–30  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 150  $\times$  4.6 5  $\mu$ m Inertsil ODS-2  
**Mobile phase:** MeCN:50 mM pH 5.0 phosphate buffer 42:5  
**Flow rate:** 0.9  
**Injection volume:** 10–30  
**Detector:** UV 230

---

**CHROMATOGRAM**

**Retention time:** 6  
**Internal standard:** indomethacin (18.5)  
**Limit of quantitation:** 50 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** diclofenac, ibuprofen, felbinac, fenbufen, flurbiprofen, loxoprofen, mefenamic acid, naproxen, piroxicam, sulindac

---

**KEY WORDS**


SPE

---

**REFERENCE**

Hirai,T.; Matsumoto,S.; Kishi,I. Simultaneous analysis of several non-steroidal anti-inflammatory drugs in human urine by high-performance liquid chromatography with normal solid-phase extraction, *J.Chromatogr.B*, **1997**, 692, 375–388.

# Ketorolac



**Molecular formula:**  $C_{15}H_{13}NO_3$

**Molecular weight:** 255.27

**CAS Registry No.:** 74103-06-3, 74103-07-4 (tromethamine)

**Merck Index:** 5318

**Lednicer No.:** 4 81

## SAMPLE

**Matrix:** blood

**Sample preparation:** Dilute 0.1-1 mL plasma to 1 mL with water, add 100  $\mu$ L 2  $\mu$ g/mL IS in MeOH:water 90:10, add 100  $\mu$ L 500 mM pH 3 sodium acetate buffer, add 6 mL hexane:ethyl acetate 70:30, vortex vigorously for 5 min. Centrifuge at 2000 rpm for 2-5 min, place in dry ice/isopropanol or dry ice/MeOH bath to freeze the aqueous layer. Decant the organic layer and evaporate it to dryness under a stream of nitrogen at 38°. Add 500  $\mu$ L MeOH:water 90:10 and 3 mL hexane, sonicate for 15 s, vortex vigorously for 3 min, let stand for at least 5 min, discard the hexane layer. Evaporate the remaining MeOH/water to dryness under a stream of nitrogen at 38°. Add 100  $\mu$ L MeOH and 100  $\mu$ L mobile phase, sonicate for 30 s, vortex for 30 s, inject a 20  $\mu$ L aliquot.

## HPLC VARIABLES

**Column:** 150  $\times$  4.6 4  $\mu$ m Nova Pak C18

**Mobile phase:** MeCN:0.05% phosphoric acid 34:66

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 317

## CHROMATOGRAM

**Retention time:** 5.5

**Internal standard:** RS-37414-000, ( $\pm$ )-5-p-fluorobenzoyl-1,2-dihydro-3H-pyrrolo[1,2a]-pyrrole-1-carboxylic acid (7)

**Limit of quantitation:** 10 ng/mL

## KEY WORDS

plasma; pharmacokinetics

## REFERENCE

Tsina,I.; Chu,F.; Kaloostian,M.; Pettibone,M.; Wu,A. HPLC method for the determination of ketorolac in human plasma, *J.Liq.Chromatogr.Rel.Technol.*, **1996**, 19, 957-967.

## SAMPLE

**Matrix:** blood

**Sample preparation:** Condition a 10  $\times$  3 C18 SPE cartridge (Analytichem) with 2 mL MeOH, 2 mL water, and 2 mL 50 mM pH 3.5 sodium acetate at 2 mL/min. 550  $\mu$ L Plasma + 550  $\mu$ L 0/9% NaCl, vortex vigorously, add 25  $\mu$ L 10  $\mu$ g/mL ketoprofen in MeOH:water 10:90, add 1 mL to the SPE cartridge at 1 mL/min, wash with 1 mL 50 mM pH 3.5 sodium acetate at 1 mL/min, wash with 1.5 mL MeOH:0.1% acetic acid 20:80 at 1.5 mL/min, elute the contents of the cartridge on to the column with mobile phase.

## HPLC VARIABLES

**Guard column:** 15  $\times$  3.2 7  $\mu$ m Newguard RP-18

**Column:** 100  $\times$  8 4  $\mu$ m Nova-pak C18 radial pak

**Mobile phase:** Gradient. MeCN:0.1% acetic acid from 30:70 to 60:40 over 10 min, maintain at 60:40 for 2 min, to 100:0 over 3 min.

**Flow rate:** 2

**Detector:** UV 313 for 7.2 min then UV 258

## CHROMATOGRAM

**Retention time:** 6.7

**Internal standard:** ketoprofen (8.8)**Limit of quantitation:** 5 ng/mL**KEY WORDS**

plasma; SPE

**REFERENCE**

Solà,J.; Pruñonosa,J.; Colom,H.; Peraire,C.; Obach,R. Determination of ketorolac in human plasma by high-performance liquid chromatography after automated on-line solid-phase extraction, *J.Liq.Chromatogr.Rel.Technol.*, **1996**, 19, 89–99.

**SAMPLE****Matrix:** blood

**Sample preparation:** 500  $\mu$ L Plasma + 100  $\mu$ L 600 mM sulfuric acid + 3 mL isooctane:isopropanol 95:5, vortex for 30 s, centrifuge at 1800 g for 5 min. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 100  $\mu$ L 2 mg/mL 4-(dimethylamino)pyridine in MeCN, add 100  $\mu$ L 60 mM trichloroethyl chloroformate in MeCN, add 1 M L-leucinamide in MeCN, let stand for 2 min, add 500  $\mu$ L 250 mM HCl, extract with chloroform. Remove the organic layer and evaporate it to dryness, reconstitute the residue in mobile phase, inject a 10–100  $\mu$ L aliquot. (A 7% conversion of S to R is observed during the derivatization procedure. No racemization is observed using a direct procedure with a chiral column.)


**HPLC VARIABLES****Column:** 100  $\times$  4.6 5  $\mu$ m Partisil 5 ODS-2**Mobile phase:** MeCN:60 mM  $\text{KH}_2\text{PO}_4$ :triethylamine 30:70:0.02**Flow rate:** 1**Injection volume:** 10–100**Detector:** UV 310**CHROMATOGRAM****Retention time:** 10.7 (R), 19.5 (S)**Internal standard:** ketorolac**OTHER SUBSTANCES****Extracted:** tiaprofenic acid**KEY WORDS**

derivatization; chiral; plasma; ketorolac is IS

**REFERENCE**

Vakily,M.; Jamali,F. Pharmacokinetics of tiaprofenic acid in humans: Lack of stereoselectivity in plasma using both direct and precolumn derivatization methods, *J.Pharm.Sci.*, **1996**, 85, 638–642.

# Ketotifen

**Molecular formula:**  $\text{C}_{19}\text{H}_{19}\text{NOS}$ **Molecular weight:** 309.43**CAS Registry No.:** 34580-13-7, 34580-14-8 (fumarate)**Merck Index:** 5319**Lednicer No.:** 3 239**SAMPLE****Matrix:** blood, tissue

**Sample preparation:** Blood. Dilute 1 mL plasma with 100  $\mu$ L 1 M pH 9 phosphate buffer and 100  $\mu$ L water, add 8 mL diethyl ether, extract. Evaporate the organic layer and dissolve the residue in 400  $\mu$ L mobile phase. Inject 200  $\mu$ L aliquot. Tissue. Homogenize the brain with 2

fold the weight of water. Dilute 1500  $\mu$ L brain homogenate with 500  $\mu$ L 1 M pH 9 phosphate buffer, add 8 mL diethyl ether, extract. Evaporate the organic layer, dissolve the residue in 400-500  $\mu$ L mobile phase. Inject a 200  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  6 Intersil ODS-2

**Mobile phase:** MeCN:0.018% TFA 20:80 (plasma), MeCN:0.018% TFA 15:85 (tissue)

**Column temperature:** 40

**Flow rate:** 0.7

**Injection volume:** 200

**Detector:** UV 300

---

**CHROMATOGRAM**

**Limit of quantitation:** 5 ng/mL (plasma), 40 ng/mL (brain)

---

**KEY WORDS**

brain; cat; mouse; pharmacokinetics; plasma; rat

---

**REFERENCE**

Kato,M.; Nishida,A.; Aga,Y.; Kita,J.; Kudo,Y.; Narita,H.; Endo,T. Pharmacokinetic and pharmacodynamic evaluation of central effect of the novel antiallergic agent betotastine besilate, *Arzneimittelforschung*, **1997**, *47*, 1116-1124.

---

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 150  $\times$  4.6 12  $\mu$ m 1-myristoyl-2-[(13-carboxyl)-tridecoyl]-sn-3-glycerophosphocholine chemically bonded to silica (Regis)

**Mobile phase:** MeCN:100 mM pH 7.0 phosphate buffer 20:80

**Flow rate:** 1

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** k' 14.72

---

**OTHER SUBSTANCES**

**Also analyzed:** acebutolol, alprenolol, antazoline, atenolol, betaxolol, bisoprolol, bopindolol, bupranolol, carteolol, celiprolol, chloropyramine, chlorpheniramine, cicloprolol, cimetidine, cinarizine, cirazoline, clonidine, dilevalol, dimethindene, diphenhydramine, doxazosin, esmolol, famotidine, isothipendyl, metiamide, metoprolol, moxonidine, nadolol, naphazoline, nifenalol, nizatidine, oxprenolol, pheniramine, phentolamine, pindolol, pizotyline (pizotifen), practolol, prazosin, promethazine, propranolol, pyrilamine (mepyramine), ranitidine, roxatidine, sotalol, tiamenidine, timolol, tramazoline, tripeleminamine, triprolidine, tymazoline, UK-14,304

---

**REFERENCE**

Kaliszan,R.; Nasal,A.; Turowski,M. Binding site for basic drugs on  $\alpha_1$ -acid glycoprotein as revealed by chemometric analysis of biochromatographic data, *Biomed.Chromatogr.*, **1995**, *9*, 211-215.

---

---

# Labetalol

**Molecular formula:** C<sub>19</sub>H<sub>24</sub>N<sub>2</sub>O<sub>3</sub>

**Molecular weight:** 328.41

**CAS Registry No.:** 36894-69-6, 32780-64-6 (HCl)

**Merck Index:** 5341

**Lednicer No.:** 3 24; 4 20

