


Levallorphan



Molecular formula: C₁₉H₂₅NO

Molecular weight: 283.41

CAS Registry No.: 152-02-3, 71-82-9 (tartrate)

Merck Index: 5485

SAMPLE

Matrix: blood, CSF, urine

Sample preparation: Dilute urine 3:1 or more with water. Mix 1 mL CSF, plasma, or diluted urine with 500 μ L saturated sodium carbonate, add 5 mL hexane containing 0.1% n-octylamine. Vortex for 60 s, centrifuge at 2000 g for 10 min. Re-extract aqueous phase with 5 mL hexane containing 0.1% n-octylamine, evaporate the combined hexane extracts to dryness under a stream of nitrogen in a 50° water bath. Reconstitute residue with 150 μ L 100 mM hydrochloric acid, inject a 100 μ L aliquot.

HPLC VARIABLES

Guard column: 10 \times 4.6 5 μ m CN

Column: 220 \times 4.6 5 μ m Brownlee Spheri-5CN (Applied Biosystems, USA)

Mobile phase: MeCN:n-octylamine:water 19:0.05:80.95 adjusted to pH 2.8 with phosphoric acid

Column temperature: 40

Flow rate: 1

Injection volume: 100

Detector: F ex 230 em 330

CHROMATOGRAM

Retention time: 7.4

Internal standard: levallorphan

OTHER SUBSTANCES

Extracted: dextromethorphan

KEY WORDS

plasma; levallorphan is IS

REFERENCE

Kimiskidis,V.K.; Kazis,A.D.; Niopas,I. Simultaneous determination of dextromethorphan and dextrorphan in human plasma, urine and cerebrospinal fluid by HPLC with fluorescence detection, *J.Liq. Chromatogr.Rel.Technol.*, **1996**, *19*, 1267–1275.

SAMPLE

Matrix: blood, saliva, urine

Sample preparation: Add 100 μ L 28% ammonium hydroxide and 5 mL n-butanol:hexane 10:90 to 5 mL urine, 1 mL plasma or 3 mL saliva. Rotate for 30 min and centrifuge at 4500 rpm for 10 min. Remove the upper organic layer and extract it with 300 μ L 100 mM HCl by vortexing for 20 min, centrifuge for 5 min. Inject a 40 μ L (urine) or 200 μ L (plasma, saliva) aliquot of the acidic layer.

HPLC VARIABLES

Column: Zorbax RP-phenyl

Mobile phase: MeCN:10 mM potassium phosphate 50:50 adjusted to pH 4.0 with 8.5% phosphoric acid

Column temperature: 40

Flow rate: 1.0

Injection volume: 40-200

Detector: F ex 280 em 310

CHROMATOGRAM

Retention time: 8.6

Internal standard: levallorphan

OTHER SUBSTANCES

Extracted: dextromethorphan

KEY WORDS

levallorphan is IS; plasma

REFERENCE

Hu,O.Y.-P.; Tang,H.-S.; Lane,H.-Y.; Chang,W.-H.; Hu,T.-M. Novel single-point plasma or saliva dextromethorphan method for determining CYP2D6 activity, *J.Pharmacol.Exp.Ther.*, **1998**, 285, 955–960.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 500 μ L Microsomal incubation + 500 μ L MeOH:water:zinc sulfate 50:45:5, centrifuge at 10000 g for 3 min, extract the supernatant twice with 3 mL dichloromethane. Combine dichloromethane extracts, dry under a stream of nitrogen, dissolve the residue in 200 μ L MeOH. Inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m polymer C18 (Astec, Whippany, NJ)

Mobile phase: MeCN:50 mM (sic) pH 9.0 ammonium carbonate buffer 60:40 (Buffer was adjusted to pH 9.0 with ammonium hydroxide.)

Column temperature: 30

Flow rate: 0.7

Injection volume: 20

Detector: UV (wavelength not given)

CHROMATOGRAM

Retention time: 11.35

Internal standard: levallorphan

OTHER SUBSTANCES

Extracted: dextromethorphan, dextrorphan


KEY WORDS

liver; levallorphan is IS

REFERENCE

Vielnascher,E.; Spatzenegger,M.; Mayrhofer,A.; Klinger,P.; Jäger,W. Metabolism of dextromethorphan in human liver microsomes: a rapid HPLC assay to monitor cytochrome P450 2D6 activity, *Pharmazie*, **1996**, 51, 586–588.

Levamisole



Molecular formula: C₁₁H₁₂N₂S

Molecular weight: 204.30

CAS Registry No.: 14769-73-4, 16595-80-5 (HCl), 5036-02-2 (racemic), 5086-74-8 (racemic HCl)

Merck Index: 5486

Lednicer No.: 4 217

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the

residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 213.4

CHROMATOGRAM

Retention time: 6.97

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

SAMPLE

Matrix: cell cultures

Sample preparation: Condition a Sep-Pak C18 SPE cartridge with 10 mL MeOH and 20 mL water. Add 2 mL cell culture to the SPE cartridge, wash with 20 mL water, elute with 10 mL MeOH:water 80:20, inject a 10 μ L portion of the eluate.

HPLC VARIABLES

Column: 100 \times 8 4 μ m RCM Nova-Pak C18

Mobile phase: MeCN:100 mM ammonium acetate 46:54

Flow rate: 1.3

Injection volume: 10

Detector: UV 245

CHROMATOGRAM

Retention time: 2.2

Limit of quantitation: 5 μ g/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

SPE

REFERENCE

Shu, Y.-Z.; Kingston, D.G.I.; Van Tassell, R.L.; Wilkins, T.D. Metabolism of levamisole, an anti-colon cancer drug, by human intestinal bacteria, *Xenobiotica*, **1991**, 21, 737-750.

SAMPLE

Matrix: milk

Sample preparation: Condition a Sep-Pak C18 SPE cartridge with 10 mL ethyl acetate, with 10 mL MeOH, and with 10 mL water. Acidify 50 mL milk to pH 4.6 with 6 M HCl, centrifuge at 10° at 4200 g for 15 min, remove the supernatant from the solid and the fat layer. Adjust

the pH of the supernatant to 11.0-11.2 with a few drops of 40% NaOH, add a 25 mL aliquot to the SPE cartridge, elute with 10 mL water-saturated ethyl acetate. Evaporate the eluate to dryness under reduced pressure at 35°, reconstitute the residue in 200 μ L mobile phase.

HPLC VARIABLES

Guard column: normal phase universal guard cartridge kit (Whatman)

Column: 250 \times 2.1 Spherisorb S5W

Mobile phase: Dichloromethane:MeOH 95:5

Flow rate: 0.2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 12

KEY WORDS

normal phase; SPE; detection with GC for higher sensitivity

REFERENCE

Chappell, C.G.; Creaser, C.S.; Shepherd, M.J. Modified on-column interface for coupled high-performance liquid chromatography-gas chromatography and its application to the determination of levamisole in milk, *J. Chromatogr.*, **1992**, 626, 223-230.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbitol, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenopropfen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaacal, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephentyoin, mephesisin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methyprylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone,

naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidone, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233–242.

SAMPLE

Matrix: solutions

Sample preparation: Inject 10 μL of a solution in RPMI-1640.

HPLC VARIABLES

Column: 250 \times 4.6 5 μm Hypersil MS C8

Mobile phase: MeCN:water 75:25 containing 0.02% dimethyloctylamine and 0.02% trifluoroacetic acid

Injection volume: 10

Detector: radioactivity

CHROMATOGRAM

Retention time: 9.9

KEY WORDS

tritium labelled

REFERENCE

Ho,N.F.H.; Sims,S.M.; Vidmar,T.J.; Day,J.S.; Barsuhn,C.L.; Thomas,E.M.; Geary,T.G.; Thompson,D.P. Theoretical perspectives on anthelmintic drug discovery: Interplay of transport kinetics, physicochemical properties, and in vitro activity of anthelmintic drugs, *J.Pharm.Sci.*, **1994**, *83*, 1052–1059.

SAMPLE

Matrix: tissue

Sample preparation: Condition a 3 mL Bakerbond Cyano SPE cartridge with 5 mL ethyl acetate:hexane 50:50, do not allow to dry out. Pulverize frozen tissue samples to a fine powder. 3 g Powdered tissue + 2 g anhydrous sodium sulfate + 9 mL ethyl acetate + 500 μL 50% KOH (w/v), homogenize (Silverson) for 1 min, centrifuge at 2000 rpm for 10 min. Remove 6 mL of the ethyl acetate extract and add it to 6 mL hexane, mix, add to the SPE cartridge, wash with 5 mL chloroform:hexane 50:50, dry under vacuum for 10 min, elute with two 2 mL portions of MeOH. Evaporate the eluate to dryness under a stream of nitrogen at 70°, reconstitute the residue in 200 μL , sonicate for 5 min, inject a 50 μL aliquot.

HPLC VARIABLES

Guard column: 5 μm LiChrospher 60 RP-select B guard column

Column: 125 \times 4 5 μm LiChrospher 60 RP-select B

Mobile phase: MeCN:THF:triethylamine:water 35:5:0.2:59.8 containing 100 mM ammonium acetate

Flow rate: 1

Injection volume: 50

Detector: MS, Hewlett-Packard 5989A, thermospray, vaporizer 15° below take-off temperature, ion source 250°, quadrupole 100°, electron multiplier +200 V with respect to autotune voltage, SIM, m/z 205

CHROMATOGRAM

Retention time: 3.5

Limit of detection: <3 ng/g

KEY WORDS

sheep; liver; kidney; muscle; SPE; LC-MS; thermospray

REFERENCE

Cannavan, A.; Blanchflower, W.J.; Kennedy, D.G. Determination of levamisole in animal tissues using liquid chromatography-thermospray mass spectrometry, *Analyst*, **1995**, *120*, 331-333.

Levobunolol


Molecular formula: C₁₇H₂₅NO₃

Molecular weight: 291.39

CAS Registry No.: 47141-42-4, 27912-14-7 (HCl)

Merck Index: 5488

Lednicer No.: 2 110, 215, 280



SAMPLE

Matrix: blood, urine

Sample preparation: Blood. Add MeCN to blood so that the ratio is 1:5. Remove a 1 mL aliquot and add 10 µL 2 µg/mL metoprolol in water, add 1 mL water, adjust pH to 9.8-10.2 with 7-10 drops 100 mM NaOH, add 2 mL benzene (CAUTION! Benzene is a carcinogen!), shake on an automatic shaker for 30 min, centrifuge, remove the organic phase and extract the aqueous layer again with 2 mL benzene for 20 min. Combine the organic layers and evaporate them under reduced pressure at 30°, dissolve the residue in 50 µL mobile phase, inject a 10-50 µL aliquot and determine amount of dihydrolevobunolol present (DHL-A). Add MeCN to blood so that the ratio was 1:5. Remove a 1 mL aliquot and add 10 µL 2 µg/mL metoprolol in water, add 1 mL water, adjust pH to 9.8-10.2 with 7-10 drops 100 mM NaOH, add 2 mL benzene, shake on an automatic shaker for 30 min, centrifuge, remove the organic phase and extract the aqueous layer again with 2 mL benzene for 20 min. Combine the organic layers and evaporate them under reduced pressure at 30°, dissolve the residue in 200 µL MeOH, add 5 mg sodium borohydride, let stand at room temperature in a closed tube for 30 min, add 1 mL water, add 300 mg NaCl, extract with 3 mL benzene for 20 min, centrifuge. Remove the organic layer and evaporate it under reduced pressure, dissolve the residue in 50 µL mobile phase, inject a 10-50 µL aliquot and determine the amount of dihydrolevobunolol (DHL-B) that is now present which represents the total amount of levobunolol and dihydrolevobunolol originally present. Determine amount of levobunolol originally present by subtracting DHL-A from DHL-B. Urine. 20-1000 µL Urine + 200 ng metoprolol + 1 mL 200 mM pH 10.2 sodium borate buffer (Sørensen), adjust pH to 9.8-10.2 with 100 mM NaOH (if necessary), add 500 mg NaCl, extract with 4 mL benzene for 30 min, proceed as described above for blood samples. (Levobunolol itself is not fluorescent so it must be reduced to the fluorescent derivative dihydrolevobunolol. Addition of MeCN to freshly drawn blood stopped enzymatic conversion of levobunolol to dihydrolevobunolol.)

HPLC VARIABLES

Column: 250 × 4.6 10 µm µBondapak C18

Mobile phase: MeOH:water 48:52 containing 0.4% phosphoric acid and 0.2% heptanesulfonic acid

Flow rate: 2

Injection volume: 10-50

Detector: F ex 225 em 295

CHROMATOGRAM**Retention time:** 3.8 (dihydrolevobunolol)**Internal standard:** metoprolol (4.6)**Limit of detection:** 0.5-1 ng/mL

KEY WORDShuman; dog

REFERENCE

Hengy,H.; K lle,E.-U. Determination of levobunolol and dihydrolevobunolol in blood and urine by high-performance liquid chromatography using fluorescence detection, *J.Chromatogr.*, **1985**, 338, 444-449.

SAMPLE**Matrix:** solutions**Sample preparation:** Dilute in an appropriate solvent, inject an aliquot.


HPLC VARIABLES**Guard column:** RC18 Guardpak (Waters)**Column:** 150 × 4.5 5 µm Altex C18**Mobile phase:** MeOH:water:phosphoric acid 48:52:0.1, containing 2 g/L sodium heptanesulfonate, pH adjusted to 3.5 with 3 M NaOH**Flow rate:** 1.5**Detector:** UV 255

CHROMATOGRAM**Retention time:** 5-5.3

REFERENCE

Richman,J.B.; Tang-Liu,D.D.-S. A corneal perfusion device for estimating ocular bioavailability in vitro, *J.Pharm.Sci.*, **1990**, 79, 153-157.

Levodopa

**Molecular formula:** C₉H₁₁NO₄**Molecular weight:** 197.19**CAS Registry No.:** 59-92-7**Merck Index:** 5490

SAMPLE**Matrix:** amniotic fluid, blood, CSF, urine

Sample preparation: Plasma. Condition a 100 mg Bond Elut SCX (propylbenzenesulfonic acid, H⁺ form) SPE cartridge with 1 mL 50 mM HCl, 1 mL MeOH, 2 mL water, and 1 mL 50 mM HCl. 100 µL Plasma + 100 µL 250 µM norleucine in 100 mM HCl + 10 mg solid sulfosalicylic acid + 800 µL acetone or MeOH, mix, centrifuge, add a 50 µL aliquot to the SPE cartridge, wash with 2 mL water, elute with two 500 µL portions of MeOH:water:triethylamine 40:40:20, dry the eluate under vacuum, add 10 µL MeOH:1 M sodium acetate:triethylamine 40:40:20, dry under vacuum at 70 mTorr, reconstitute with 20 µL MeOH:triethylamine:water:phenylisothiocyanate 70:10:10:10, let stand at room temperature for 20 min, evaporate to dryness under vacuum, reconstitute with 100 µL MeCN:5 mM pH 7.4 sodium phosphate buffer 5:95, inject a 20 µL aliquot. Dried blood. Add 25 µL 250 µM norleucine in 100 mM HCl to a 6 mm filter paper disc containing dried blood, add 100 µL MeCN, let stand for 30 min, centrifuge, remove a 75 µL aliquot of the supernatant, evaporate to dryness under reduced pressure, add 10 µL MeOH:1 M sodium acetate:triethylamine 40:40:20, dry under vacuum at 70 mTorr, reconstitute with 20 µL MeOH:triethylamine:water:phenylisothiocyanate 70:10:10:10, let stand at room temperature for 2 min, evaporate to dryness under vacuum, reconstitute with 50 µL MeCN:5 mM pH 7.4 sodium phosphate buffer 5:95, inject a 20 µL aliquot. Amniotic fluid, CSF. Mix amniotic fluid or CSF with an equal volume of 250 µM norleucine in 100 mM HCl, filter (Centrifree 10000 MW cutoff) while centrifuging at 2200 g. Evaporate a 50 µL aliquot of the

ultrafiltrate to dryness under vacuum, add 10 μL MeOH:1 M sodium acetate:triethylamine 40:40:20, dry under vacuum at 70 mTorr, reconstitute with 20 μL MeOH:triethylamine:water:phenylisothiocyanate 70:10:10:10, let stand at room temperature for 20 min, evaporate to dryness under vacuum, reconstitute with 50 (CSF) or 100 (amniotic fluid) μL MeCN:5 mM pH 7.4 sodium phosphate buffer 5:95, inject a 20 μL aliquot. Urine. Dilute urine with water to a creatinine concentration of 1 mM, mix an aliquot with an equal volume of 250 μM norleucine in 100 mM HCl, filter (Centrifree 10000 MW cutoff) while centrifuging at 2200 g. Evaporate a 50 μL aliquot of the ultrafiltrate to dryness under vacuum, add 10 μL MeOH:1 M sodium acetate:triethylamine 40:40:20, dry under vacuum at 70 mTorr, reconstitute with 20 μL MeOH:triethylamine:water:phenylisothiocyanate 70:10:10:10, let stand at room temperature for 20 min, evaporate to dryness under vacuum, reconstitute with 100 μL MeCN:5 mM pH 7.4 sodium phosphate buffer 5:95, inject a 20 μL aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 Pico-Tag amino acid column (Waters)

Mobile phase: Gradient. A was MeCN:70 mM pH 6.55 sodium acetate 2.5:97.5. B was MeCN:MeOH:water 45:15:40

Column temperature: 46

Flow rate: 1

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 34.14

Internal standard: norleucine (55.07)

OTHER SUBSTANCES

Extracted: β -alanine, alanine, alloisoleucine, α -aminoadipic acid, 4-aminobenzoic acid, gamma-aminobutyric acid, β -amino-n-butyric acid, gamma-amino-n-butyric acid, 4-aminohippuric acid, β -aminoisobutyric acid, 4-aminophenylacetic acid, α -aminophenylacetic acid, 3-amino-3-phenylpropionic acid, δ -amino-n-valeric acid, ammonia, anserine, arginine, asparagine, aspartic acid, aspartylglucosamine, carnosine, citrulline, cystathionine, cysteic acid, cysteine, cysteine-homocysteine (mixed disulfide), cystine, ethanolamine, ethionine, ethylamine, galactosamine, glucosamine, glutamic acid, glutamine, glutathionine (oxidized), glycine, glycyglycine, glycyhistidine, glycyllucine, glycyphenylalanine, glycytyrosine, histidine, homoarginine, homocitrulline, homoserine, homocystine, 3-hydroxyanthranilic acid, 3-hydroxykynurenine, hydroxyproline, isoleucine, kynurenine, leucine, lysine, methionine sulfone, methionine, 3-methylhistidine, 1-methylhistidine, ornithine, phenylalanine, phosphoethanolamine, phosphoserine, proline, sarcosine, serine, serotonin, taurine, threonine, tromethamine, tryptophan, tyrosine, valine

Noninterfering: cadaverine, 2-phenylethylamine

KEY WORDS

derivatization; SPE; ultrafiltrate; plasma; dried blood

REFERENCE

Davey, J.F.; Ersser, R.S. Amino acid analysis of physiological fluids by high-performance liquid chromatography with phenylisothiocyanate derivatization and comparison with ion-exchange chromatography, *J. Chromatogr.*, **1990**, 528, 9–23.

SAMPLE

Matrix: blood

Sample preparation: Mix 450 μL plasma with 50 μL 310 $\mu\text{g/L}$ iso-homovanilic acid in antioxidant solution. Dialyze this solution using a Carnegie Medicin (Stockholm) microdialysis system. Perfuse 10 mm microdialysis probes with Ringer solution at 2 $\mu\text{L/min}$. Collect dialysate for each plasma sample over 20 min in a vial containing 80 μL antioxidant solution, inject a 100 μL aliquot. (Antioxidant solution was 10 mM HCl containing 1 g/L sodium metabisulfite and 0.1 g/L Na_2EDTA .)

HPLC VARIABLES

Guard column: 30 \times 4 Bondapak C18/Corasil 37-50 μm m

Column: 250 \times 4.8 5 μm ODS (Beckman, CA)

Mobile phase: MeOH:buffer 20:80 (Buffer was 70 mM pH 2.55 NaH_2PO_4 containing 2.08 mM sodium octanesulfonate and 80 μM EDTA.)

Flow rate: 1

Injection volume: 100

Detector: E, Waters 460 containing an electrochemical cell fitted with a glassy carbon working electrode and an Ag/AgCl reference electrode; the detector potential is + 0.8 V vs the reference electrode

CHROMATOGRAM

Retention time: 3.5

Internal standard: iso-homovanilic acid (10)

Limit of detection: 1.1 nM/L

OTHER SUBSTANCES

Extracted: dihydroxyphenylacetic acid, dopamine, homovanilic acid

KEY WORDS

plasma; pharmacokinetics; dialysate

REFERENCE

Dethy,S.; Laute,M.A.; Van Blercom,N.; Damhaut,P.; Goldman,S.; Hildebrand,J. Microdialysis-HPLC for plasma levodopa and metabolites monitoring in parkinsonian patients, *Clin.Chem.*, **1997**, 43, 740–744.

SAMPLE

Matrix: blood

Sample preparation: Total levodopa. Add 300 μL 60 mM trichloroacetic acid to 1 mL plasma, place the mixture in an ice bath for 10 min, centrifuge at 5000 g for 10 min, inject a 20 μL aliquot of the supernatant. Unbound levodopa. Add 250 μL plasma to 5 μL 4 M orthophosphoric acid, centrifuge through a membrane Minicent 10 (Bio-Rad) at 8000 g at 25° for 30 min, inject a 20 μL aliquot of the ultrafiltrate.

HPLC VARIABLES

Column: 45 \times 4.6 5 μm C18 (Bio-Rad)

Mobile phase: MeCN:1.4 mM sodium dodecylsulphate:50 mM KH_2PO_4 16:42:42, pH 2.8

Flow rate: 1

Injection volume: 20

Detector: E, ESA Coulochem 5100, analytical cell +50 mV, -300 mV, conditioning cell +300 mV

CHROMATOGRAM

Retention time: 3.7

Limit of detection: 200 pg/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; ultrafiltrate

REFERENCE

Rizzo,V.; Memmi,M.; Moratti,R.; Melzi,E.; Perucca,E. Concentrations of L-dopa in plasma and plasma ultrafiltrates, *J.Pharm.Biomed.Anal.*, **1996**, 14, 1043–1046.

SAMPLE

Matrix: blood

Sample preparation: 0.5 mL Plasma + 50 μL 4 M perchloric acid + 50 μL 1 $\mu\text{g/mL}$ dihydroxybenzylamine in 0.1 M perchloric acid, centrifuge at 1500 g for 10 min. Remove 300 μL supernatant and centrifuge it at 1600 g through a 0.2 μm regenerated cellulose filter, inject a 20 μL aliquot.

HPLC VARIABLES

Guard column: 20 \times 4.6 5 μm Biophase ODS + 50 \times 4.6 Pelliguard LC-18

Column: 250 × 4.6 5 µm Biophase ODS or 250 × 4.6 Phase II ODS (both from Bioanalytical Systems)

Mobile phase: MeOH:buffer 5:95 (Buffer was 20 mM sodium citrate, 100 mM NaH₂PO₄, 0.15 mM, and 1.25 mM heptanesulfonic acid, pH 3.2.)

Column temperature: 28

Flow rate: 1-1.5

Injection volume: 20

Detector: E, Bioanalytical Systems LC-150 in dual-parallel mode, channel 1 700 mV 200 nA f.s. for levodopa and 3-O-methyldopa, channel 2 560 mV 10 nA f.s. for dopamine, carbidopa, and dihydroxyphenylacetic acid, Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 5.1

Internal standard: dihydroxybenzylamine (7)

Limit of detection: 1 ng/mL

OTHER SUBSTANCES

Simultaneous: 3-O-methyldopa, dopamine, carbidopa, dihydroxyphenylacetic acid

KEY WORDS

plasma

REFERENCE

Cedarbaum, J.M.; Williamson, R.; Kutt, H. Simultaneous determination of levodopa, its metabolites and carbidopa in clinical samples, *J. Chromatogr.*, **1987**, *415*, 393-399.

SAMPLE

Matrix: blood

Sample preparation: Prepare a 20 × 5 polypropylene column packed with CM-Sephadex pre-swollen in water, wash with 5 mL 2 M HCl, wash with 10 mL water, wash with 10 mL 100 mM pH 7 phosphate buffer. Add 1 mL plasma to column, elute with 5.5 mL water, discard first 1 mL. Add next 4.5 mL to 0.5 mL 0.5 M perchloric acid, centrifuge, inject 10 µL aliquot of supernatant.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Nucleosil C18

Mobile phase: MeCN:MeOH:25 mM sodium acetate 4:4:92 containing 0.2 mM 1-octanesulfonic acid and 0.3 mM disodium EDTA, pH was adjusted to pH 3 with acetic acid

Flow rate: 0.9

Injection volume: 10

Detector: E, ESA Coulochem 5100 A, 5010 A analytical cell, first electrode +0.25 V, second electrode -0.30 V

CHROMATOGRAM

Retention time: 7

Limit of quantitation: 5 ng/mL

OTHER SUBSTANCES

Simultaneous: O-methyldopa, carbidopa, dihydroxyphenylacetic acid (DOPAC)

KEY WORDS

plasma

REFERENCE

Betto, P.; Ricciarello, G.; Giambenedetti, M.; Lucarelli, C.; Ruggeri, S.; Stocchi, F. Improved high-performance liquid chromatographic analysis with double detection system for L-dopa, its metabolites and carbidopa in plasma of parkinsonian patients under L-dopa therapy, *J. Chromatogr.*, **1988**, *459*, 341-349.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 100 μ L 1.25 μ g/mL α -ethylidopa in 0.1 M HCl + 100 μ L 4 M perchloric acid, vortex, centrifuge at 2000 g for 10 min, inject a 60 μ L aliquot of the supernatant (keep sample tray at $6 \pm 1^\circ$).

HPLC VARIABLES

Guard column: 45 \times 4.6 37-40 μ m Whatman pellicular-ODS followed by 45 \times 4.6 5 μ m Ultra sphere-IP C18

Column: 250 \times 4.6 5 μ m Ultrasphere IP C18

Mobile phase: MeOH:20 mM orthophosphoric acid and 4 mM sodium octanesulfonate 25:75 adjusted to pH 2.8 ± 0.05 with 50% NaOH

Column temperature: 40

Flow rate: 1

Injection volume: 60

Detector: E, BAS LC-4B, 0.75 V vs Ag/AgCl, 5 nA full scale for carbidopa, 20 nA full scale for 3-O-methylidopa and levodopa

CHROMATOGRAM

Retention time: 5.9

Internal standard: α -ethylidopa (15.4)

Limit of quantitation: 25 ng/mL

OTHER SUBSTANCES

Extracted: carbidopa, 3-O-methylidopa

KEY WORDS

plasma; stabilize plasma sample immediately with EDTA and 2 mg/mL sodium metabisulfite

REFERENCE

Titus,D.C.; August,T.F.; Yeh,K.C.; Eisenhandler,R.; Bayne,W.F.; Musson,D.G. Simultaneous high-performance liquid chromatographic analysis of carbidopa, levodopa and 3-O-methylidopa in plasma and carbidopa, levodopa and dopamine in urine using electrochemical detection, *J.Chromatogr.*, **1990**, 534, 87-100.

SAMPLE

Matrix: blood

Sample preparation: 4 mL Plasma + 500 μ L 20 mg/mL ascorbic acid solution, vortex for 30 s. 1 mL Aliquot + 75 mg acid washed alumina + 100 μ L 1 μ g/mL 3,4-dihydroxybenzylamine hydrobromide in buffer, vortex, add 1 mL 1.5 M pH 8.6 TRIS buffer, shake at 230 oscillations/min for 15 min. Allow to settle and discard plasma, wash the alumina twice by shaking with 5 mL water for 10 min. To the washed alumina add 900 μ L buffer, vortex for 20 s, allow to settle, inject a 50 μ L aliquot of the supernatant. (Buffer was 200 mM phosphoric acid containing 3.3 μ M EDTA and 6.7 μ M potassium metabisulfite.)

HPLC VARIABLES

Column: 100 \times 4.6 3 μ m Spherisorb ODS-2

Mobile phase: MeCN:buffer 8:92 (Buffer was pH 2.6 550 mM NaH_2PO_4 containing 1 mM sodium octyl sulfate and 0.7 mM EDTA.)

Flow rate: 1.5

Injection volume: 50

Detector: E, Bioanalytical Systems LC-4B, glassy carbon electrode, Ag/AgCl reference electrode, 0.75 V.

CHROMATOGRAM

Retention time: 3.3

Internal standard: 3,4-dihydroxybenzylamine hydrobromide (4.5)

Limit of quantitation: 5 ng/mL

OTHER SUBSTANCES

Extracted: carbidopa

Noninterfering: caffeine, ibuprofen, aspirin, nicotine, acetaminophen, theophylline

KEY WORDS

plasma; SPE

REFERENCE

Miller,R.B.; Dehelean,L.; Bélanger,L. Determination of carbidopa and levodopa in human plasma by high-performance liquid chromatography with electrochemical detection, *Chromatographia*, **1993**, 35, 607–612.

SAMPLE

Matrix: blood

Sample preparation: 0.5 mL Plasma + 25 μ L 2 μ g/mL 3,4-dihydroxybenzylamine in 0.4 M perchloric acid + 25 μ L 70% perchloric acid, vortex 1 min, keep on ice for 1 min, store at -80°. Allow to thaw at +4°, vortex 1 min, centrifuge at 1200 g at +4° for 15 min. Remove 300 μ L supernatant and add it to 200 μ L 2 M pH 4.5 potassium citrate buffer, centrifuge at 1200 g at +4° for 10 min, inject 50 μ L of the supernatant.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Spherisorb C18

Mobile phase: MeOH:MeCN:buffer 8:4:88 (Buffer was 50 mM pH 3.2 phosphate containing 3.5 mM heptanesulfonic acid and 0.05 mM EDTA.)

Flow rate: 1

Injection volume: 50

Detector: E, ESA Model 5100, Model 5020 guard cell +0.6 V, Model 5010 analytical cell, DET 1 +0.35 V, DET 2 -0.35 V, both DET 1 and DET 2 monitored

CHROMATOGRAM

Retention time: 4.15

Internal standard: 3,4-dihydroxybenzylamine (8.26)

Limit of detection: 19.6 ng/mL

OTHER SUBSTANCES

Simultaneous: 3-O-methyldopa, dopamine, L-DOPA methyl ester

KEY WORDS

plasma; rat; human

REFERENCE

Rondelli,I.; Acerbi,D.; Mariotti,F.; Ventura,P. Simultaneous determination of levodopa methyl ester, levodopa, 3-O-methyldopa and dopamine in plasma by high-performance liquid chromatography with electrochemical detection, *J.Chromatogr.B*, **1994**, 653, 17–23.

SAMPLE

Matrix: blood, tissue

Sample preparation: Plasma. 20 μ L Plasma + 180 μ L isoproterenol in 100 mM perchloric acid containing 10 μ M disodium EDTA, centrifuge at 4° at 4000 rpm for 10 min. Filter (0.45 μ m cellulose acetate) the supernatant and inject a 10 μ L aliquot of the supernatant. Tissue. Homogenize (Polytron PT 10-35) 50 mg brain tissue + 500 μ L isoproterenol in 100 mM perchloric acid containing 10 μ M disodium EDTA at 15000 rpm for 10 s, centrifuge at 4° at 4000 rpm for 10 min. Filter (0.45 μ m cellulose acetate) the supernatant and inject a 10 μ L aliquot of the supernatant.

HPLC VARIABLES

Guard column: Supelcosil LC-18-DB

Column: Supelcosil LC-18-DB

Mobile phase: MeOH:10 mM pH 4.4 citrate buffer 10:90 containing 10 μ M disodium EDTA and 0.5 mM sodium 1-octanesulfonate

Flow rate: 1

Injection volume: 10

Detector: E, EICOM Co. ECD-100, +0.7 V, Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 4

Internal standard: isoproterenol (45)

KEY WORDS

plasma; rat; brain; pharmacokinetics

REFERENCE

Sato,S.; Koitabashi,T.; Koshiro,A. Pharmacokinetic and pharmacodynamic studies of L-dopa in rats. I. Pharmacokinetic analysis of L-dopa in rat plasma and striatum, *Biol.Pharm.Bull.*, **1994**, *17*, 1616–1621.

SAMPLE

Matrix: blood, urine

Sample preparation: Condition a 150 μ L Toyopak IC-SP S (sulfopropyl resin, H⁺ form) SPE cartridge (Tosoh) with 10 mL water. Plasma. 700 μ L Plasma + 50 μ L 700 nM 3,4-dihydroxybenzylamine + 350 μ L 2 M perchloric acid, mix, centrifuge at 4° at 1000 g for 15 min. Remove a 700 μ L aliquot of the supernatant and add it to 30 μ L 2 M potassium carbonate, centrifuge at 4° at 1000 g for 5 min. Add a 500 μ L aliquot of the supernatant to the SPE cartridge, wash with 1 mL water, wash with 500 μ L EtOH:water 50:50, wash with 5 mL water, elute with 500 μ L 2 M sodium perchlorate, filter (0.2 μ m), inject a 50 μ L aliquot of the filtrate. Urine. Acidify urine collected over 24 h with 10 mL 6 M HCl. 500 μ L Urine + 25 μ L 10 μ M 3,4-dihydroxybenzylamine + 25 μ L 40 μ M ferulic acid + 500 μ L 1 M perchloric acid, mix, centrifuge at 4° at 1000 g for 15 min. Remove a 700 μ L aliquot of the supernatant and add it to 30 μ L 2 M potassium carbonate, centrifuge at 4° at 1000 g for 5 min, add a 500 μ L aliquot of the supernatant to the SPE cartridge, wash with 1.5 mL water, wash with 500 μ L EtOH:water 50:50, wash with 5 mL water, elute with 500 μ L 2 M sodium perchlorate, filter (0.2 μ m), inject a 50 μ L aliquot of the filtrate.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m TSK-gel ODS-80TM (Tosoh)

Mobile phase: Gradient. A was buffer. B was MeCN:MeOH:buffer 8:12:80, pH 3.1. A:B 100:0 for 4 min, to 60:40 over 8 min, to 0:100 over 2 min, maintain at 0:100 for 16 min, return to initial conditions (step gradient), re-equilibrate for 20 min. Buffer was 60 mM pH 3.1 citric acid containing 32 mM Na₂HPO₄, 1.7 mM sodium hexanesulfonate, and 0.1 mM disodium EDTA (J. Chromatogr. 1989, 467, 237).

Flow rate: 1

Injection volume: 50

Detector: F ex 345 em 480 following post-column reaction. The column effluent passed through a Hitachi 655A electrochemical detector with carbon cloth electrodes; working electrode at +0.68 V versus reference electrode (200 mM equimolar mixture of potassium hexacyanoferrate(II) and potassium hexacyanoferrate(III) containing 200 mM potassium nitrate and 200 mM KOH). The effluent from the electrochemical detector mixed with 20 mM meso-1,2-diphenylethylenediamine in 50 mM HCl pumped at 0.4 mL/min and with 1 M glycine containing 490 mM KOH and 3 mM potassium hexacyanoferrate(III) pumped at 0.4 mL/min. This mixture flowed through a 10 m \times 0.47 mm ID coil at 80° to the detector (J. Chromatogr. 1989, 467, 237).

CHROMATOGRAM

Retention time: 9.5

Internal standard: 3,4-dihydroxybenzylamine (12.5)

Limit of detection: 12 nM (urine), 10 nM (plasma)

OTHER SUBSTANCES

Extracted: dopamine, epinephrine, metanephrine, 3-methoxytyramine, norepinephrine

KEY WORDS

post-column reaction; plasma; SPE

REFERENCE

Nohta,H.; Yamaguchi,E.; Ohkura,Y.; Watanabe,H. Measurement of catecholamines, their precursor and metabolites in human urine and plasma by solid-phase extraction followed by high-performance liquid chromatography with fluorescence derivatization, *J.Chromatogr.*, **1989**, *493*, 15–26.

SAMPLE

Matrix: blood, urine

Sample preparation: Condition a Toyopak IC-SP S sulfopropyl resin, H⁺ form, SPE cartridge (Tosoh) with 10 mL water and 2 mL 200 mM pH 5.0 sodium phosphate buffer. Plasma. 700 μ L Plasma + 30 μ L 700 nM isoproterenol + 50 μ L 7 μ M 3,4-dihydroxyphenylpropanoic acid + 350 μ L 2 M perchloric acid, mix, centrifuge at 4° at 1000 g for 15 min. Remove a 700 μ L aliquot

of the supernatant and adjust the pH to 1.5-2.0 with about 150 μL 2 M potassium carbonate, centrifuge at 4° at 1000 g for 5 min, add the supernatant to the SPE cartridge, wash with 10 mL water, elute with 300 μL MeOH:2 M sodium perchlorate 7:93, filter (cellulose acetate membrane), inject a 100 μL aliquot of the filtrate. Urine. Collect human urine for 24 h in the presence of 10 mL 6 M HCl. 500 μL Urine + 10 μL 15 μM isoproterenol + 25 μL 800 μM 3,4-dihydroxyphenylpropanoic acid + 500 μL 1 M perchloric acid, mix, centrifuge at 4° at 1000 g for 15 min. Remove a 700 μL aliquot of the supernatant and adjust the pH to 1.5-2.0 with about 130 μL 2 M potassium carbonate, centrifuge at 4° at 1000 g for 5 min, add the supernatant to the SPE cartridge, wash with 1.5 mL water, wash with 500 μL EtOH:water 50:50, wash with 5 mL water, elute with 500 μL 1.5 M KCl in MeOH:100 mM HCl 7:93, filter (cellulose acetate membrane), inject a 100 μL aliquot of the filtrate.

HPLC VARIABLES

Column: 250 \times 4.6 5 μm TSK-gel ODS-80TM (Tosoh)

Mobile phase: MeOH:buffer 7:93 (Buffer was 30 mM pH 2.5 citrate buffer containing 0.4 mM sodium octanesulfonate.)

Flow rate: 0.8

Injection volume: 100

Detector: F ex 350 em 480 following post-column reaction. The column effluent mixed with reagent A pumped at 0.3 mL/min and the mixture flowed through a 3 m \times 0.5 mm ID stainless steel coil at 90°. The effluent from this coil mixed with reagent B pumped at 0.3 mL/min and the mixture flowed through a 10 m \times 0.5 mm ID stainless steel coil at 90° and through a 1 m \times 0.5 mm ID stainless steel cooling coil to the detector (Anal. Sci. 1991, 7, 257). (Reagent A was 10 mM sodium periodate containing 3 mM potassium ferricyanide. Reagent B was 30 mM meso-1,2-diphenylethylenediamine in EtOH:water 70:30 containing 130 mM sodium methylate.)

CHROMATOGRAM

Retention time: 23

Internal standard: isoproterenol (60)

Limit of detection: 7-9 nM

OTHER SUBSTANCES

Extracted: dopamine, epinephrine, metanephrine, 3-methoxytyramine, norepinephrine, normetanephrine

KEY WORDS

post-column reaction; plasma; SPE

REFERENCE

Jeon, H.-K.; Nohta, H.; Ohkura, Y. High-performance liquid chromatographic determination of catecholamines and their precursor and metabolites in human urine and plasma by postcolumn derivatization involving chemical oxidation followed by fluorescence reaction, *Anal. Biochem.*, **1992**, *200*, 332-338.

SAMPLE

Matrix: blood, urine

Sample preparation: 2 mL Plasma or 1 mL urine + dihydroxybenzylamine + 20 mg Sigma WA4 alumina + 200 μL 1 M pH 8.6 Tris-EDTA buffer, mix for 10 min, discard plasma. Wash the alumina three times with 3 mL water and dry it. Add 125 μL 500 mM phosphoric acid, after 1 min inject a 100 μL aliquot. (Ann. Clin. Biochem. 1985, 22, 194-203)

HPLC VARIABLES

Column: 250 \times 4.5 5 μm Ultratechsphere

Mobile phase: Per liter 75 mmol citric acid, 58.5 mmol NaH_2PO_4 , 0.2 mmol disodium EDTA, and 4.4 mmol heptanesulfonic acid, pH adjusted to 3.4, made up to a final volume of 2 L, add 200 mL MeOH

Flow rate: 1

Injection volume: 100

Detector: E, ESA Coulochem conditioning cell +0.35 V, first electrode +0.05 V, second electrode -0.35 V

CHROMATOGRAM**Retention time:** 5.37**Internal standard:** dihydroxybenzylamine (10.53)**Limit of detection:** 50 ng/mL

OTHER SUBSTANCES**Simultaneous:** norepinephrine, metanephrine, epinephrine, 3-methoxytyrosine, normetanephrine, dihydroxyphenylacetic acid, dopamine

KEY WORDSplasma

REFERENCE

Dutton, J.; Copeland, L.G.; Playfer, J.R.; Roberts, N.B. Measuring L-dopa in plasma and urine to monitor therapy of elderly patients with Parkinson disease treated with L-dopa and a dopa decarboxylase inhibitor, *Clin. Chem.*, **1993**, 39, 629–634.

SAMPLE**Matrix:** blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES**Guard column:** 20 mm long Symmetry C18**Column:** 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10–30**Detector:** UV 200.5

CHROMATOGRAM**Retention time:** 3.575

KEY WORDSwhole blood

REFERENCE

Gaillard, X.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149–163.

SAMPLE**Matrix:** formulations

Sample preparation: Powder tablets or contents of capsules, weigh out an amount equivalent to about 100 mg levodopa, add 30 mL 0.1 M HCl, sonicate, make up to 50 mL with 0.1 M HCl, mix, filter (0.45 μ m), discard first 5 mL filtrate. 5 mL Filtrate + 10 mL 2 mg/mL methyl dopa in 0.1 M HCl, make up to 100 mL with mobile phase, mix, inject a 20 μ L aliquot.

HPLC VARIABLES**Column:** 300 \times 3.9 10 μ m μ Bondapak C18**Mobile phase:** 3% aqueous acetic acid

Flow rate: 1.5
Injection volume: 20
Detector: UV 280

CHROMATOGRAM

Retention time: 3
Internal standard: methyl dopa (4.5)

KEY WORDS

capsules; tablets

REFERENCE

Ting,S. Liquid chromatographic determination of levodopa and levodopa-carbidopa in solid dosage forms: collaborative study, *J.Assoc.Off.Anal.Chem.*, **1987**, 70, 987–990.

SAMPLE

Matrix: formulations

Sample preparation: Dissolve in mobile phase, filter, inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.6 10 μm μBondapak C18

Mobile phase: MeOH:50 mM ammonium acetate adjusted to pH 4.1 with 0.6 M acetic acid 1:99

Flow rate: 0.9

Detector: E, Coulochem model 5100A, screen electrode +0.3 V, sample electrode +0.6 V and UV 280

CHROMATOGRAM

Retention time: 6.35

Limit of detection: 1000 ng/mL (UV), 20 ng/mL (E)

OTHER SUBSTANCES

Simultaneous: hydroxydopa, carbidopa, methyl dopa, methoxytyrosine, methylcarbidopa, impurities

KEY WORDS

stability-indicating; tablets

REFERENCE

Kafil,J.B.; Dhingra,B.S. Stability-indicating method for the determination of levodopa, levodopa-carbidopa and related impurities, *J.Chromatogr.A*, **1994**, 667, 175–181.

SAMPLE

Matrix: formulations

Sample preparation: Tablets. Grind tablets, weigh out a portion, dissolve in 50 mL mobile phase, sonicate, filter (No. 4 sintered glass plate), dilute, inject an aliquot. Capsules. Dissolve 10 capsules (without opening) in 100 mL mobile phase, sonicate, inject an aliquot. Injections, ampules, sprays. Dilute, inject an aliquot.

HPLC VARIABLES

Column: 120 × 4.6 Spherisorb C18 ODS-2

Mobile phase: Isopropanol:buffer 5:95 (Buffer was 100 mM sodium dodecyl sulfate containing 25 mM Na₂HPO₄, pH adjusted to 3.0 with HCl.)

Flow rate: 1

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 3

Limit of detection: 4 ng/mL

OTHER SUBSTANCES

Simultaneous: carbidopa, dopamine, epinephrine, hydrochlorothiazide, isoproterenol, methyl-dopa, norepinephrine, phenylephrine

KEY WORDS

tablets; capsules; injections; ampules; sprays

REFERENCE

Villanueva Camañas, R.M.; Sanchis Mallols, J.M.; Torres Lapasió, J.R.; Ramis-Ramos, G. Analysis of pharmaceutical preparations containing catecholamines by micellar liquid chromatography with spectrophotometric detection, *Analyst*, **1995**, *120*, 1767–1772.

SAMPLE

Matrix: solutions

Sample preparation: Inject an aliquot of a 100 µg/mL solution in mobile phase.

HPLC VARIABLES

Column: 150 × 4.5 µm Crownpak CR(+) immobilized crown ether

Mobile phase: MeOH:0.1% pH 1.9 perchloric acid 15:85

Column temperature: 40

Flow rate: 1

Detector: UV 210

CHROMATOGRAM

Retention time: 1.48, 1.70 (enantiomers)

OTHER SUBSTANCES

Simultaneous: baclofen, norephedrine, primaquine

KEY WORDS

chiral; comparison with capillary electrophoresis

REFERENCE

Nishi, H.; Nakamura, K.; Nakai, H.; Sato, T. Separation of enantiomers and isomers of amino compounds by capillary electrophoresis and high-performance liquid chromatography utilizing crown ethers, *J. Chromatogr. A*, **1997**, *757*, 225–235.

SAMPLE

Matrix: tissue

Sample preparation: Prepare a 70 × 5 SPE column of Sephadex G 10 in a Pasteur pipette, wash with 3 mL 20 mM ammonia and 3 mL 10 mM formic acid, let stand for 10 days. Homogenize up to 150 mg rat brain in 1 mL 100 mM perchloric acid, centrifuge at 4000 g at 4° for 15 min, add 500 µL of the supernatant to the SPE column, wash with 2.5 mL 10 mM formic acid, elute with 1 mL 10 mM formic acid followed by 1.5 mL 5 mM Na₂HPO₄, inject an aliquot of the eluate.

HPLC VARIABLES

Column: Nucleosil 5 C18

Mobile phase: pH 5.5 Buffer prepared from 200 mM Na₂HPO₄ and 100 mM citric acid

Flow rate: 0.8

Injection volume: 200

Detector: E, rotating disc electrode, 500 mV

CHROMATOGRAM

Retention time: 4.5

Limit of detection: 0.015 nmole/g

OTHER SUBSTANCES

Extracted: dopamine, uric acid

KEY WORDS

rat; brain; SPE

REFERENCE

Westerink,B.H.C.; Mulder,T.B.A. Determination of picomole amounts of dopamine, noradrenaline, 3,4-dihydroxyphenylacetic acid, homovanillic acid, and 5-hydroxyindolacetic acid in nervous tissue after one-step purification on Sephadex G-10, using high-performance liquid chromatography with a novel type of electrochemical detection, *J.Neurochem.*, **1981**, 36, 1449-1462.

SAMPLE**Matrix:** urine

Sample preparation: 100 μ L Urine + 100 μ L solution containing 55 mM ascorbic acid and 55 mM disodium EDTA + 25 μ L 1.25 μ g/mL α -ethylidopa in 0.1 M HCl + 25 mg alumina + 1 mL 2 M pH 8.6 Tris-HCl buffer in a microfilter tube (Centrex, Schleicher & Schuell), vortex 5 min, allow to stand for 10 min, filter off water, wash with 5 mL water, add 5 mL water, centrifuge at 3000 g, vortex with 400 μ L 0.2 M perchloric acid containing 11 mM disodium EDTA and 0.4 M sodium metabisulfite, centrifuge at 9000 g for 5 min, inject 50 μ L of filtrate. (Stabilize each 10 mL urine sample immediately with 0.5 mL 0.1 M HCl and 1 mL solution containing 55 mM ascorbic acid and 55 mM disodium EDTA.)

HPLC VARIABLES**Guard column:** 40 \times 4.6 Bio-Sil ODS-10 (Bio-Rad)**Column:** 250 \times 4.6 5 μ m Ultrasphere IP C18

Mobile phase: MeOH:water 22.5:77.5 containing 20 mM citric acid, 20 mM Na_2HPO_4 , 4 mM sodium octanesulfonate, and 0.05 mM disodium EDTA, pH adjusted to 2.74 ± 0.01 with 2 M citric acid

Column temperature: 40**Injection volume:** 50**Detector:** E, BAS LC-4B, 0.54 V vs Ag/AgCl, 50 nA full scale

CHROMATOGRAM**Retention time:** 6**Internal standard:** α -ethylidopa (14)**Limit of quantitation:** 250 ng/mL

OTHER SUBSTANCES**Extracted:** carbidopa, dopamine

KEY WORDS

SPE

REFERENCE

Titus,D.C.; August,T.F.; Yeh,K.C.; Eisenhandler,R.; Bayne,W.F.; Musson,D.G. Simultaneous high-performance liquid chromatographic analysis of carbidopa, levodopa and 3-O-methylidopa in plasma and carbidopa, levodopa and dopamine in urine using electrochemical detection, *J.Chromatogr.*, **1990**, 534, 87-100.

Levonordefrin


Molecular formula: C₉H₁₃NO₃

Molecular weight: 183.21

CAS Registry No.: 829-74-3, 138-61-4 (racemic HCl)

Merck Index: 6785

Lednicer No.: 1 68



SAMPLE

Matrix: formulations

Sample preparation: Tablets. Dissolve powdered tablets in 10 mM HCl, filter if necessary, inject an aliquot. Injections, solutions. Dilute with 10 mM HCl, inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 μm Partisil-5 ODS-3

Mobile phase: MeOH:buffer 30:70 (Buffer was 10 mM sodium 1-octanesulfonate in 0.2% acetic acid.)

Flow rate: 1

Injection volume: 20

Detector: UV 220

CHROMATOGRAM

Retention time: 12

Limit of detection: 36 ng

OTHER SUBSTANCES

Simultaneous: norepinephrine, epinephrine, isoproterenol, phenylephrine, metaraminol, impurities

KEY WORDS

tablets; injections; ophthalmic solutions; inhalation solutions

REFERENCE

Smela, M.J., Jr.; Stromberg, R. Liquid chromatographic determination of six sympathomimetic drugs in dosage forms, *J. Assoc. Off. Anal. Chem.*, **1991**, 74, 289–291.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 5 μm Partisil ODS-3

Mobile phase: MeOH:buffer 30:70 (Buffer was 10 mM octanesulfonic acid in 0.2% acetic acid.)

Flow rate: 1

Detector: UV 220

CHROMATOGRAM

Retention time: 12


OTHER SUBSTANCES

Simultaneous: epinephrine, isoproterenol, metaraminol, phenylephrine

REFERENCE

Phenomenex Catalog, **1994**, p. 1.077.

Levorphanol



Molecular formula: $C_{17}H_{23}NO$

Molecular weight: 257.38

CAS Registry No.: 77-07-6, 5985-38-6 (tartrate dihydrate), 125-72-4 (tartrate)

Merck Index: 5496

Lednicer No.: 1 293

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 20 ng levallorphan, mix, add 1 mL 1 M pH 9 borate buffer, add 10 mL hexane:ethyl acetate 90:10, shake for 10 min, centrifuge at 2500 rpm for 10 min. Remove 9 mL of the organic layer and evaporate it to dryness under a stream of nitrogen at 50°, reconstitute the residue in 200 μ L mobile phase, inject a 100 μ L aliquot.

HPLC VARIABLES

Column: 300 \times 4 10 μ m μ Bondapak C18

Mobile phase: MeCN:buffer 30:70 containing 0.1 mM EDTA (Buffer was 10 mM NaCl adjusted to pH 4.8 with 1 M HCl.)

Flow rate: 1

Injection volume: 100

Detector: E, Bioanalytical Systems LC-4B, glassy carbon electrode +1.00 V, Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 6.5

Internal standard: levallorphan (9.5)

Limit of quantitation: 1.25 ng/mL

OTHER SUBSTANCES

Simultaneous: 6-acetylmorphine, diamorphine, oxymorphone, pentazocine

Noninterfering: acetaminophen, 1- α -acetylmethodol, caffeine, codeine, hydrocodone, hydromorphone, meperidine, morphine, oxycodone, propoxyphene

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Lucek,R.; Dixon,R. Quantitation of levorphanol in plasma using high-performance liquid chromatography with electrochemical detection, *J.Chromatogr.*, **1985**, *341*, 239–243.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 100 μ L 2.58 μ g/mL laudanosine in MeCN + 500 μ L saturated sodium carbonate solution, vortex for 10 s, add 5 mL chloroform, vortex for 10 s, mix on a rocking mixer for 40 min, centrifuge at 2000 g for 25 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 400 μ L mobile phase, inject a 300 μ L aliquot. (Hydrolyze conjugates by heating 1 mL plasma with 1 mL 3000 U/mL β -glucuronidase (Helix pomatia type H-1 (Sigma)) in 100 mM pH 5.0 sodium citrate at 37° for 2 h, proceed as above.)

HPLC VARIABLES

Column: 150 \times 4.6 Spherisorb 5-CN

Mobile phase: MeCN:water:triethylamine 10:89:1 adjusted to pH 6 with orthophosphoric acid

Flow rate: 1

Injection volume: 300

Detector: F ex 280 em 315

CHROMATOGRAM**Retention time:** 6.615**Internal standard:** laudanosine (8.603)**Limit of detection:** 30 ng/mL

OTHER SUBSTANCES**Extracted:** guaifenesin

KEY WORDSplasma

REFERENCE

Stavchansky,S.; Demirbas,S.; Reyderman,L.; Chai,C.-K. Simultaneous determination of dextrorphan and guaifenesin in human plasma by liquid chromatography with fluorescence detection, *J.Pharm.Biomed.Anal.*, **1995**, 13, 919-925.

SAMPLE**Matrix:** blood, CSF, urine

Sample preparation: Dilute urine 3:1 or more with water. Vortex 1 mL CSF, plasma, or diluted urine with 100 μ L 100 ng/mL IS, add 500 μ L saturated sodium carbonate, mix, add 5 mL hexane containing 0.1% n-octylamine. Vortex for 60 s, centrifuge at 2000 g for 10 min. Re-extract aqueous phase with 5 mL hexane containing 0.1% n-octylamine, evaporate the combined hexane extracts to dryness under a stream of nitrogen in a 50° water bath. Reconstitute residue with 150 μ L 100 mM HCl, inject a 100 μ L aliquot.

HPLC VARIABLES**Guard column:** 10 \times 4.6 5 μ m CN**Column:** 220 \times 4.6 5 μ m Brownlee Spheri-5CN (Applied Biosystems, USA)**Mobile phase:** MeCN:n-octylamine:water 19:0.05:80.95 adjusted to pH 2.8 with phosphoric acid**Column temperature:** 40**Flow rate:** 1**Injection volume:** 100**Detector:** F ex 230 em 330

CHROMATOGRAM**Retention time:** 5.3**Internal standard:** levallorphan (7.4)**Limit of detection:** 1 ng/mL**Limit of quantitation:** 1 ng/mL (CSF, plasma), 5 ng/mL (urine)

OTHER SUBSTANCES**Extracted:** dextromethorphan

KEY WORDSplasma

REFERENCE

Kimiskidis,V.K.; Kazis,A.D.; Niopas,I. Simultaneous determination of dextromethorphan and dextrorphan in human plasma, urine and cerebrospinal fluid by HPLC with fluorescence detection, *J.Liq. Chromatogr.Rel.Technol.*, **1996**, 19, 1267-1275.

SAMPLE**Matrix:** blood, urine

Sample preparation: Plasma. 1 mL Plasma + 100 μ L 1 μ g/mL pholcodine in water + 500 μ L saturated sodium carbonate, mix, add 4 mL diethyl ether:chloroform:isopropanol 20:9:1, mix on a rotary mixer for 10 min, centrifuge at 2000 g for 10 min. Remove the organic layer and add it to 100 μ L 100 mM HCl, mix on a rotary mixer for 10 min, centrifuge at 2000 g for 5 min, inject a 10-50 μ L aliquot of the aqueous layer. Urine. 500 μ L Urine + 50 μ L 50 μ g/mL pholcodine in water + 500 μ L saturated sodium carbonate, mix, add 4 mL diethyl ether:chloroform:isopropanol 20:9:1, mix on a rotary mixer for 10 min, centrifuge at 2000 g for 10 min. Remove the organic layer and add it to 100 μ L 100 mM HCl, mix on a rotary mixer for 10

min, centrifuge at 2000 g for 5 min, inject a 10-50 μL aliquot of the aqueous layer. (If desired, hydrolyse 500 μL plasma or urine with 500 μL 8000 U/mL β -glucuronidase (*Helix pomatia*, type H-1, Sigma) in 200 mM pH 5 acetate buffer at 37° for 16 h, proceed as above.)

HPLC VARIABLES

Column: 150 \times 4.6 5 μm Spherisorb cyano

Mobile phase: MeCN:water:triethylamine 17:82.94:0.06, adjusted to pH 3.0 with orthophosphoric acid

Flow rate: 1

Injection volume: 10-50

Detector: F ex 230 em 330

CHROMATOGRAM

Retention time: 5.1

Internal standard: pholcodine (6.3)

Limit of detection: 1 ng/mL

OTHER SUBSTANCES

Extracted: metabolites, dextromethorphan

Noninterfering: acetaminophen, clofibrate, codeine, cyclophosphamide, diclofenac, digoxin, doxepin, doxorubicin, estrogens, flucloxacillin, folic acid, furosemide, metformin, metoclopramide, miconazole, minoxidil, morphine, nifedipine, nitroglycerin, norcodeine, norethisterone, oxazepam, oxethazaine, prednisolone, pseudoephedrine, quinine, spironolactone, temazepam, tolbutamide, warfarin

KEY WORDS

plasma

REFERENCE

Chen,Z.R.; Somogyi,A.A.; Bochner,F. Simultaneous determination of dextromethorphan and three metabolites in plasma and urine using high-performance liquid chromatography with application to their disposition in man, *Ther.Drug Monit.*, **1990**, 12, 97-104.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 600 μL Microsomal incubation + 600 μL saturated sodium carbonate, place on ice, add 150 μL 12.5 $\mu\text{g/mL}$ betaxolol, extract with 5 mL ethyl acetate. Remove the organic layer and add it to 300 μL 0.5% orthophosphoric acid, extract, inject an aliquot of the aqueous layer.

HPLC VARIABLES

Column: 5 mm i.d. 4 μm Nova-Pak phenyl radial-Pak

Mobile phase: MeCN:MeOH:0.05% orthophosphoric acid 24:10:66

Flow rate: 1.6

Detector: F ex 261 em 306

CHROMATOGRAM

Retention time: 5.5

Internal standard: betaxolol

Limit of detection: 60 nM

OTHER SUBSTANCES

Extracted: dextromethorphan

KEY WORDS

rat; liver

REFERENCE

Laslett,T.J.; Alvarez,F.; Nation,R.L.; Evans,A.M.; Scott,S.D.; Stupans,I. Effect of cyclophosphamide administration on the activity and relative content of hepatic P4502D1 in rat, *Xenobiotica*, **1995**, 25, 1031-1039.

SAMPLE**Matrix:** solutions**Sample preparation:** Dissolve in MeOH at a concentration of 1 mg/mL, inject a 20 μ L aliquot.

HPLC VARIABLES**Column:** 250 \times 5 Spherisorb S5W**Mobile phase:** MeOH:buffer 90:10 (Buffer was 94 mL 35% ammonia and 21.5 mL 70% nitric acid in 884 mL water, adjust the pH to 10.1 with ammonia.)**Flow rate:** 2**Injection volume:** 20**Detector:** UV 254

CHROMATOGRAM**Retention time:** 5.55

OTHER SUBSTANCES

Simultaneous: buprenorphine, dextromoramide, phenoperidine, fentanyl, etorphine, piritramide, noscapine, papaverine, naloxone, dextropropoxyphene, nalorphine, phenazocine, norpipanone, levallorphan, hydroxypethidine, normethadone, meperidine, dipipanone, diamorphine, pentazocine, acetylcodeine, monoacetylmorphine, thebacon, oxycodone, thebaine, norlevorphanol, methadone, benzylmorphine, ethylmorphine, morphine-N-oxide, codeine, codeine-N-oxide, morphine, ethoheptazine, morphine-3-glucuronide, pholcodeine, norpethidine, hydrocodone, dihydrocodeine, dihydromorphine, norcodeine, normorphine, pemoline, benzphetamine, diethylpropion, mazindol, tranlycypromine, caffeine, fenethyline, phendimetrazine, methylphenidate, phenelzine, epinephrine, pipradol, phenylpropanolamine, fencamfamin, chlorphen-termine, norpseudoephedrine, phentermine, fenfluramine, methylenedioxymphetamine, amphetamine, normetanephine, 4-hydroxyamphetamine, bromo-STP, STP, prolintane, 2-phenethylamine, tyramine, trimethoxyamphetamine, phenylephrine, pseudoephedrine, ephedrine, methylephedrine, dimethylamphetamine, methamphetamine, mescaline, mephentermine

Noninterfering: dopamine, levodopa, methyldopa, methyldopate, norepinephrine

REFERENCE

Law,B.; Gill,R.; Moffat,A.C. High-performance liquid chromatography retention data for 84 basic drugs of forensic interest on a silica column using an aqueous methanol eluent, *J.Chromatogr.*, **1984**, *301*, 165–172.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 150 \times 4.6 Supelcosil LC-ABZ**Mobile phase:** MeCN:25 mM pH 6.9 potassium phosphate buffer 35:65**Flow rate:** 1.5**Injection volume:** 25**Detector:** UV 254

CHROMATOGRAM**Retention time:** 2.755

OTHER SUBSTANCES

Also analyzed: 6-acetylmorphine, amiloride, amphetamine, benzocaine, benzoylecgonine, caffeine, cocaine, codeine, doxylamine, fluoxetine, glutethimide, hexobarbital, hypoxanthine, LSD, meperidine, mephobarbital, methadone, methylphenidate, methyprylon, N-norcodeine, oxazepam, oxycodone, phenylpropanolamine, prilocaine, procaine, terfenadine

REFERENCE

Ascah,T.L. Improved separations of alkaloid drugs and other substances of abuse using Supelcosil LC-ABZ column, *Supelco Reporter*, **1993**, *12*(3), 18–21.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 250 × 4.6 Zorbax RX**Mobile phase:** Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.**Column temperature:** 30**Flow rate:** 2**Detector:** UV 210**OTHER SUBSTANCES**

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisquinone, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenopropfen, fenproporex, fentanyl, flubendazole, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-stilbene, imipramine, indomethacin, isocarboxystiril, isocarboxamide, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclufenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephentoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, meth-azolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, meth-yl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methyprylon, meto-prolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, nor-epinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphen-butazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, per-santine, phenacetin, phenazocine, phenazopyridine, phenacyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenyl-butazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primi-done, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopola-mine, scopoletin, secobarbital, strychnine, sulfacetamide, sufadiazine, sulfadimethoxine, sul-faethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sul-fasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tol-metin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapa-mil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, **1994**, *18*, 233–242.

SAMPLE**Matrix:** solutions**HPLC VARIABLES****Column:** 250 × 4.5 µm LiChrospher 100 RP-8

Mobile phase: MeCN:0.025% phosphoric acid:buffer 60:25:15 (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 4.26

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordiazepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazindol, mefenamic acid, meperidine, mephenytoin, mepivacaine, mesoridazine, metaproterenol, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phenotolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, resperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spirinolactone, sulfinpyrazone, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocainide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

KEY WORDS

details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, 692, 103–119.

SAMPLE

Matrix: urine

Sample preparation: 1 mL Urine + 1 mL 100 mM pH 5 acetate buffer + 20 μ L β -glucuronidase/ β -arylsulfatase (Helix pomatia, Boehringer Mannheim) + 50 μ L 600 mM sodium azide in water, heat at 37° for 12 h. Inject a 100 μ L aliquot onto column A and elute to waste at 0.5 mL/min, after 3 min elute contents of column A onto column B at 1.4 mL/min, monitor the effluent from column B.

HPLC VARIABLES

Column: A 30 \times 4 10 μ m LiChrosorb CN; B 250 \times 4.6 5 μ m Spherisorb phenyl

Mobile phase: MeCN:10 mM KH_2PO_4 50:50, adjusted to pH 4 with phosphoric acid

Flow rate: A 0.5; B 1.4

Injection volume: 100

Detector: UV 200 or F ex 280 em 310

CHROMATOGRAM**Retention time:** 9.5**Limit of detection:** <50 ng/mL (UV)

OTHER SUBSTANCES**Extracted:** metabolites, dextromethorphan

KEY WORDS

column-switching

REFERENCE

Motassim,N.; Decolin,D.; Le Dinh,T.; Nicolas,A.; Siest,G. Direct determination of dextromethorphan and its three metabolites in urine by high-performance liquid chromatography using a precolumn switching system for sample clean-up, *J.Chromatogr.*, **1987**, 422, 340-345.

SAMPLE**Matrix:** urine

Sample preparation: Condition a Bond Elut silica modified with carboxylic acid ion-exchange groups SPE cartridge with 1 mL MeCN:100 mM HCl 40:60 and 1 mL water. Adjust 1 mL urine to pH 5.0-5.5, add β -glucuronidase/arylsulfatase (Helix pomatia (Boehringer Mannheim), heat at 37° for 18 h, add 1 mL to the SPE cartridge, wash with 1 mL water, wash with 500 μ L 100 mM HCl, elute with 1 mL MeCN:100 mM HCl 40:60, inject a 20 μ L aliquot of the eluate.

HPLC VARIABLES**Column:** 250 mm long 5 μ m Zorbax phenyl**Mobile phase:** MeCN:100 mM KH_2PO_4 45:55, adjusted to pH 4**Flow rate:** 1.5**Injection volume:** 20**Detector:** F ex 280 em 310 or UV 280

CHROMATOGRAM**Retention time:** 4.8**Limit of detection:** 20 ng/mL (F)

OTHER SUBSTANCES**Extracted:** metabolites, dextromethorphan

KEY WORDS

SPE

REFERENCE

Jacqz-Aigrain,E.; Menard,Y.; Popon,M.; Mathieu,H. Dextromethorphan phenotypes determined by high-performance liquid chromatography and fluorescence detection, *J.Chromatogr.*, **1989**, 495, 361-363.

SAMPLE**Matrix:** urine

Sample preparation: Condition a 3 mL 200 mg Bond Elut C18 SPE cartridge with 6 mL MeOH, 6 mL water, and 4 mL 100 mM pH 9.2 sodium carbonate buffer. 750 μ L Urine + 750 μ L 100 mM pH 5.0 sodium acetate buffer containing 20 μ L β -glucuronidase-arylsulfatase (Helix pomatia, 100000 Fisherman units/mL, Boehringer Mannheim) + 50 μ L 600 mM sodium azide, heat at 37° for 18 h. 250 μ L Hydrolysed urine + 100 μ L 10 μ g/mL levallorphan tartrate in water + 2 mL 100 mM pH 9.2 sodium carbonate, add to SPE cartridge, wash with 2 mL water, wash with 1 mL MeCN, elute with 3 mL MeOH:MeCN:2% phosphoric acid 50:30:20. Evaporate the eluate to dryness under a stream of nitrogen at 70°, reconstitute the residue in 500 μ L mobile phase, inject a 20 μ L aliquot. (For low concentrations of dextromethorphan: 500 μ L Hydrolysed urine + 100 μ L 1 μ g/mL levallorphan tartrate in water + 2 mL 100 mM pH 9.2 sodium carbonate, add to SPE cartridge, wash with 2 mL water, wash with 1 mL MeCN, elute with 3 mL MeOH:MeCN:2% phosphoric acid 50:30:20. Evaporate the eluate to dryness under a stream of nitrogen at 70°, reconstitute the residue in 250 μ L mobile phase, inject a 120 μ L aliquot.)

HPLC VARIABLES

Guard column: 15 × 3.2 RP-2 (Brownlee)

Column: 250 × 4.6 5 µm Zorbax phenyl

Mobile phase: MeCN:MeOH:10 mM pH 2.5 phosphate buffer containing 2.5 mM 1-octanesulfonic acid 27:13:60

Column temperature: 30

Flow rate: 1

Injection volume: 20-120

Detector: F ex 270 em 312

CHROMATOGRAM

Retention time: 7.0

Internal standard: levallorphan tartrate (9.5)

Limit of detection: 1 ng/mL

OTHER SUBSTANCES

Extracted: metabolites, dextromethorphan

KEY WORDS

SPE

REFERENCE

Wenk,M.; Todesco,L.; Keller,B.; Follath,F. Determination of dextromethorphan and dextrorphan in urine by high-performance liquid chromatography after solid-phase extraction, *J.Pharm.Biomed.Anal.*, **1991**, 9, 341-344.

SAMPLE

Matrix: urine

Sample preparation: 250 µL Urine + 1.25 µg levallorphan + 250 µL 140 mM pH 5 sodium acetate buffer, mix, add 25 µL β-glucuronidase (glucurase, from bovine liver, 5000 U/mL, Sigma), heat at 37° overnight, add 1.5 mL pH 11.3 glycine buffer, add 6 mL hexane:butanol 90:10, shake vigorously for 10 min, centrifuge at 1500 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 37°, reconstitute the residue in 200 µL mobile phase, inject a 10 µL aliquot.

HPLC VARIABLES

Column: Nucleosil 7 C6H5

Mobile phase: MeCN:MeOH:buffer 20:10:70 (Buffer was 10 mM phosphate buffer containing 2.5 mM sodium 1-octanesulfonate, adjusted to pH 2.5 with concentrated phosphoric acid.)

Flow rate: 1.3

Injection volume: 10

Detector: F ex 270 em 312

CHROMATOGRAM

Retention time: 6

Internal standard: levallorphan (9)

OTHER SUBSTANCES

Extracted: dextromethorphan

REFERENCE

Caslavska,J.; Hufschmid,E.; Theurillat,R.; Desiderio,C.; Wolfisberg,H.; Thormann,W. Screening for hydroxylation and acetylation polymorphisms in man via simultaneous analysis of urinary metabolites of mephentanyl, dextromethorphan and caffeine by capillary electrophoretic procedures, *J.Chromatogr.B*, **1994**, 656, 219-231.

SAMPLE

Matrix: urine

Sample preparation: 250 µL Urine + 5000 U β-glucuronidase in 1 M pH 5.0 sodium acetate buffer, heat at 37° for 18 h, add 500 µL saturated sodium carbonate, add 10 mL hexane:

triethylamine 99.9:0.1. Remove the organic layer, dry, reconstitute the residue in 250 μ L mobile phase, inject a 5-50 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 Selectosil (Phenomenex)

Mobile phase: MeCN:10 mM pH 3.0 potassium phosphate buffer 30:70

Flow rate: 1

Injection volume: 5-50

Detector: F ex 280 em 305

CHROMATOGRAM

Limit of detection: 2 μ M

OTHER SUBSTANCES

Extracted: dextromethorphan

REFERENCE

Marinac,J.S.; Foxworth,J.W.; Willisie,S.K. Dextromethorphan polymorphic hepatic oxidation (CYP2D6) in healthy black american adult subjects, *Ther.Drug Monit.*, **1995**, *17*, 120-124.

Levothyroxine sodium


Molecular formula: $C_{15}H_{10}I_4NNaO_4$

Molecular weight: 798.86

CAS Registry No.: 55-03-8, 25416-65-3 (hydrate), 51-48-9 (free acid)

Merck Index: 5497

Lednicer No.: 1 97

**SAMPLE**

Matrix: blood

Sample preparation: Equilibrate a Sep-Pak silica cartridge with 5 mL ethyl acetate. 1 mL Serum + 3 mL 5% trichloroacetic acid + 4 mL ethyl acetate, vortex vigorously, centrifuge at 1500 g for 5 min. Remove organic layer and repeat extraction twice with 3 mL portions of ethyl acetate. Combine extracts, evaporate to about 1.5 mL, add to Sep-Pak cartridge. Wash with 8 mL ethyl acetate, elute with 4 mL MeOH:ammonium hydroxide (90:10). Evaporate the eluent to dryness under nitrogen, reconstitute in 100 μ L MeOH, inject.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Ultrasphere I.P.

Mobile phase: MeCN:buffer 35:65 (Buffer was 13.6 g sodium acetate, 1.0 g cupric sulfate pentahydrate, 0.92 g L-proline, and 0.34 g silver nitrate.)

Flow rate: 1.5

Injection volume: 100

Detector: E, Bioanalytical Systems Inc. TL-5 Kel-F glassy carbon thin-layer cell, LC-4 electronic controller, +0.78 V, 2-5 nA/V

CHROMATOGRAM

Retention time: 8

Limit of detection: 3 ng/mL

OTHER SUBSTANCES

Extracted: dextrothyroxine, triiodothyronine

KEY WORDS

serum; SPE

REFERENCE

Hay, I.D.; Annesley, T.M.; Jiang, N.S.; Gorman, C.A. Simultaneous determination of D- and L-thyroxine in human serum by liquid chromatography with electrochemical detection, *J. Chromatogr.*, **1981**, 226, 383-390.

SAMPLE

Matrix: bulk, formulations

Sample preparation: Weigh out powder equivalent to about 65 mg thyroid, add 5 mL enzyme solution, mix well, incubate at 37° for 28 h, agitate after 4-8 h and after 20-24 h, add 2 mL deactivating solution, mix well, centrifuge at 2000 rpm for 5-10 min, if necessary filter (0.45 µm). (The enzyme solution was about 150 protease units/mL of bacterial protease from *Streptomyces griseus* in 110 mM NaCl + 40 mM Tris buffer + 50 mM methimazole (pH adjusted to 8.4 ± 0.05 with 6 M HCl) reducing buffer. Deactivating solution was 1:100 phosphoric acid: MeCN.)

HPLC VARIABLES

Column: 300 × 4 µm Bondapak C18

Mobile phase: MeCN:0.5% phosphoric acid in water 28:72

Column temperature: 34

Flow rate: 1.5

Injection volume: 200

Detector: UV 225

CHROMATOGRAM

Retention time: 22

OTHER SUBSTANCES

Simultaneous: liothyronine, L-3,3',5'-triiodothyronine

KEY WORDS

tablets

REFERENCE

Richheimer, S.L.; Jensen, C.B. Determination of liothyronine and levothyroxine in thyroid preparations by liquid chromatography, *J. Pharm. Sci.*, **1986**, 75, 215-217.

SAMPLE

Matrix: formulations

Sample preparation: Grind a tablet, add 50 µg 3,3',5'-triiodothyronine, add 20 mL solvent A, stir for 10 min, add 40 mL solvent B, stir for 30 min, filter. Remove the upper layer and wash it six times with 15 mL portions of water saturated with butanol, evaporate under vacuum at 40-42°, reconstitute in 2.5 mL 3% ammonium hydroxide in MeOH, inject an aliquot (Anal. Lett. 1979, 12, 1201). (Prepare the solvents by mixing 1.8 L 1-butanol, 1.35 L water and 450 mL concentrated HCl, shake vigorously for 20 min, allow to separate. The lower layer was solvent A and the upper layer was solvent B.)

HPLC VARIABLES

Guard column: 25 × 2.5 Co: Pel ODS

Column: 300 × 3.9 10 µm µBondapak C18

Mobile phase: MeOH:100 mM pH 5.0 ammonium acetate 50:50

Flow rate: 1.5

Detector: UV 254

CHROMATOGRAM

Retention time: 19.5

Internal standard: 3,3',5'-triiodothyronine

OTHER SUBSTANCES

Simultaneous: liothyronine

KEY WORDS

protect from light; tablets

REFERENCE

Rapaka,R.S.; Knight,P.W.; Prasad,V.K. Reversed-phase high-performance liquid chromatographic analysis of liothyronine sodium and levothyroxine sodium in tablet formulations: preliminary studies on dissolution and content uniformity, *J.Pharm.Sci.*, **1981**, 70, 131–134.

SAMPLE

Matrix: formulations

Sample preparation: Powder tablets, weigh out amount equivalent to about 200 µg sodium levothyroxine, add 10 mL mobile phase, sonicate for 5 min, centrifuge. Filter (0.45 µm, 25 mm Acrodisc CR, Gelman) the supernatant, inject a 200 µL aliquot.

HPLC VARIABLES

Guard column: 40 × 4 40 µm RP 201SC pellicular (Vydac)

Column: 300 × 4 µBondapak C18

Mobile phase: MeCN:buffer 60:40 (Buffer was pH 3.0 containing 5 mM 1-octanesulfonic acid and 5 mM tetramethylammonium chloride.)

Flow rate: 2

Injection volume: 200

Detector: UV 230

CHROMATOGRAM

Retention time: 4.5

OTHER SUBSTANCES

Simultaneous: liothyronine, 3,5-diiodo-L-thyronine

KEY WORDS

tablets; stability-indicating

REFERENCE

Richheimer,S.L.; Amer,T.M. Stability-indicating assay, dissolution, and content uniformity of sodium levothyroxine in tablets, *J.Pharm.Sci.*, **1983**, 72, 1349–1351.

SAMPLE

Matrix: formulations

Sample preparation: Grind tablets containing about 1 mg levothyroxine, add 4.5 mL 0.5 mg/mL hydroxyprogesterone caproate in MeOH, add 20.5 mL 10 mM NaOH in MeOH:water 75:25, shake intermittently for 5 min, filter, discard first 5 mL filtrate, inject a 25 µL aliquot.

HPLC VARIABLES

Column: 300 × 3.9 µBondapak CN

Mobile phase: MeCN:0.1% phosphoric acid in water 35:65

Flow rate: 3

Injection volume: 25

Detector: UV 225

CHROMATOGRAM

Retention time: 5

Internal standard: hydroxyprogesterone caproate (8)

KEY WORDS

tablets

REFERENCE

Das Gupta,V.; Odom,C.; Bethea,C.; Plattenburg,J. Effect of excipients on the stability of levothyroxine sodium tablets, *J.Clin.Pharm.Ther.*, **1990**, 15, 331–336.

SAMPLE

Matrix: formulations

Sample preparation: Condition a 13 mm Empore C18 SPE disk (Baker) with 2.5 mL MeOH and 2.5 mL water at 1.5 mL/min. Dissolve tablet in dissolution medium (?). Pass 40 mL through the SPE disk, wash with 2.5 mL water, dry, add 1 mL MeOH and let it soak in for 3 min, elute at 0.5 mL/min, dilute the eluate to 2 mL with 50 mM pH 2.5 phosphate buffer, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 125 \times 4 5 μ m LiChrospher cyano

Mobile phase: MeCN:water:phosphoric acid 45:55:0.05

Column temperature: 25

Flow rate: 1.6

Injection volume: 20

Detector: UV 225

KEY WORDS

tablets; SPE; comparison with capillary electrophoresis

REFERENCE

Carducci,C.N.; Lucangiolì,S.E.; Rodríguez,V.G.; Fernández Otero,G.C. Application of extraction disks in dissolution tests of clenbuterol and levothyroxine tablets by capillary electrophoresis, *J.Chromatogr.A*, **1996**, 730, 313–319.

SAMPLE

Matrix: solutions

Sample preparation: Take up 1.5 mg levothyroxine in 200 μ L 100 mM sodium bicarbonate and 400 μ L reagent, stir in an ice bath for 30 min, evaporate to dryness below 30°, add 100 μ L trifluoroacetic acid to the dry residue, let stand for 30 min at room temperature, add 2 mL 1 M sodium bicarbonate, centrifuge. Remove the precipitate and dissolve it in 600 μ L MeOH:20 mM NaOH 50:50, inject a 15 μ L aliquot. Reagent was 7 mg/mL BOC-L-Leu-SU (tert-butyloxy-L-leucine-N-hydroxysuccinimide ester) in MeOH, prepared immediately before use.)

HPLC VARIABLES

Column: 150 \times 3.2 7 μ m LiChrosorb RP-18

Mobile phase: MeOH:water 60:40 containing 0.05% methanesulfonic acid

Flow rate: 1

Injection volume: 15

Detector: UV 230

CHROMATOGRAM

Retention time: 9

Limit of detection: 0.05% of the D form

OTHER SUBSTANCES

Simultaneous: dextrothyroxine, impurities

KEY WORDS

derivatization; chiral

REFERENCE

Lankmayr,E.P.; Budna,K.W.; Nachtmann,F. Separation of enantiomeric iodinated thyronines by liquid chromatography of diastereomers, *J.Chromatogr.*, **1980**, 198, 471–479.

SAMPLE

Matrix: solutions

Sample preparation: 20 μ L Solution + 40 μ L 50 mM pH 8.5 borate buffer + 40 μ L 4 mM dabsyl chloride in MeCN, mix, heat at 70° for 15 min, add 100 μ L 25 mM pH 6.5 sodium acetate buffer, inject an aliquot.

HPLC VARIABLES

Guard column: 10 \times 3 30 μ m Chromspher C18 (Chrompack)

Column: 100 \times 3 5 μ m Chromspher C18 (Chrompack)

Mobile phase: Gradient. A was MeOH:25 mM pH 6.5 sodium acetate buffer 70:30. B was MeOH:25 mM pH 6.5 sodium acetate buffer 90:10. A:B from 100:0 to 0:100 over 15 min, maintain at 0:100 for 5 min.

Column temperature: 35

Flow rate: 0.7

Detector: UV 436

CHROMATOGRAM

Retention time: 14.5

Limit of detection: 0.39 pmole

OTHER SUBSTANCES

Simultaneous: 3,5-diiodothyronine, 3,5-diiodotyrosine, liothyronine, 3-monoiodotyrosine, thyronine, 3,3',5-triiodothyronine, tyrosine

KEY WORDS

derivatization; comparison with other derivatization procedures

REFERENCE

Doorn,L.; Jansen,E.H.; Van Leeuwen,F.X. Comparison of high-performance liquid chromatographic detection methods for thyronine and tyrosine residues in toxicological studies of the thyroid, *J.Chromatogr.*, **1991**, 553, 135–142.

SAMPLE

Matrix: solutions

Sample preparation: 100 μ L Solution + 100 μ L 500 mM pH 7.7 borate buffer + 100 μ L 2.5 mM 9-fluorenylmethyl chloroformate in dry acetone, mix, let stand at room temperature for 45 s, add 200 μ L 12 mM 1-adamantamine in MeCN, inject an aliquot.

HPLC VARIABLES

Guard column: 10 \times 3 30 μ m Chrospher C18 (Chrompack)

Column: 100 \times 3 5 μ m Chrospher C18 (Chrompack)

Mobile phase: Gradient. A was MeOH:50 mM pH 4.2 sodium acetate buffer 40:60. B was MeCN:MeOH:50 mM pH 4.2 sodium acetate buffer 20:60:20. A:B from 100:0 to 0:100 over 40 min.

Column temperature: 35

Flow rate: 0.7

Detector: UV 260

CHROMATOGRAM

Retention time: 33

Limit of detection: 1.6 pmole

OTHER SUBSTANCES

Simultaneous: 3,5-diiodothyronine, 3,5-diiodotyrosine, liothyronine, 3-monoiodotyrosine, thyronine, 3,3',5-triiodothyronine

KEY WORDS

derivatization; comparison with other derivatization procedures

REFERENCE

Doorn,L.; Jansen,E.H.; Van Leeuwen,F.X. Comparison of high-performance liquid chromatographic detection methods for thyronine and tyrosine residues in toxicological studies of the thyroid, *J.Chromatogr.*, **1991**, 553, 135–142.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Phenomenex cyano-bonded silica

Mobile phase: MeCN:water:phosphoric acid 400:600:1

Flow rate: 1.5
Detector: UV 225

CHROMATOGRAM

Retention time: 8.6

OTHER SUBSTANCES

Simultaneous: degradation products, liothyronine

REFERENCE

Won,C.M. Kinetics of degradation of levothyroxine in aqueous solution and in solid state, *Pharm.Res.*, **1992**, 9, 131-137.

SAMPLE

Matrix: tissue

Sample preparation: 100 μ L Thyroid tissue + 200 μ L MeCN, mix, centrifuge. Remove a 100 μ L aliquot of the supernatant and add it to 100 μ L 4 nM dabsyl chloride in MeCN, heat at 70° for 10 min, add 400 μ L MeOH:50 mM pH 7.0 phosphate buffer 50:50, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Hypersil ODS

Mobile phase: Gradient. A was MeOH:25 mM pH 6.5 sodium acetate 56:44. B was MeOH. A:B from 80:20 to 35:65 over 15 min, maintain at 35:65 for 3 min, to 0:100 over 1 min, maintain at 0:100 for 2 min.

Flow rate: 1

Injection volume: 20

Detector: UV 436

CHROMATOGRAM

Retention time: 17.5

OTHER SUBSTANCES

Extracted: diiodothyronine (T2), liothyronine (T3)

KEY WORDS

derivatization; thyroid

REFERENCE

Jansen,E.H.J.M.; van den Berg,R.H.; Both-Miedema,R.; Doorn,L. Advantages and limitations of pre-column derivatization of amino acids with dabsyl chloride, *J.Chromatogr.*, **1991**, 553, 123-133.

Lidocaine


Molecular formula: C₁₄H₂₂N₂O

Molecular weight: 234.34

CAS Registry No.: 137-58-6, 6108-05-0 (HCl monohydrate), 73-78-9 (HCl)

Merck Index: 5505

Lednicer No.: 1 16



SAMPLE

Matrix: blood

Sample preparation: Mix 1 mL plasma with 200 μ L 2 μ g/mL IS in MeOH, add 2 mL water and 2 mL MeCN, vortex gently, set aside for 3 min, centrifuge at 2200 g for 20 min. Separate the clear supernatant, add 500 μ L 200 mM NaOH and extract with 6 mL n-hexane by vortexing for 2 min. Centrifuge at 2200 g for 15 min. Evaporate 5 mL of the organic phase to dryness under reduced pressure. Reconstitute the residue in 120 μ L mobile phase. Inject a 100 μ L aliquot.