


Mepivacaine



Molecular formula: $C_{15}H_{22}N_2O$

Molecular weight: 246.35

CAS Registry No.: 96-88-8, 1722-62-9 (HCl)

Merck Index: 5905

Lednicer No.: 1 17

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 25 μ L EtOH + 25 μ L 17.6 μ g/mL ropivacaine in EtOH, mix. Add 5 mL diethyl ether, vortex for 3 min, centrifuge at 2000 g for 10 min at 4°. Add 500 μ L 100 mM hydrochloric acid to organic phase, back extract for 3 min. Discard organic phase, add 5 mL n-pentane containing 100 μ L isoamyl alcohol, wash for 3 min, discard organic phase again. Add 50 μ L 2 M NaOH, extract aqueous phase with 5 mL n-pentane containing 100 μ L isoamyl alcohol, centrifuge at 2000 g for 10 min at 4°. Evaporate organic phase to dryness under a gentle stream of nitrogen at 40°. Reconstitute the residue in 100 μ L mobile phase, inject a 30 μ L aliquot.

HPLC VARIABLES

Guard column: 10 \times 4.0 40 μ m α 1-AGP (J.T.Baker)

Column: 100 \times 4.0 5 μ m α 1-AGP (J.T.Baker)

Mobile phase: 2-Propanol:buffer 6.8:93.2 (Buffer was 3.6 g/L $Na_2HPO_4 \cdot 12 H_2O$ adjusted to pH 6.8 with phosphoric acid.)

Column temperature: 30

Flow rate: 1.1

Injection volume: 30

Detector: UV 210

CHROMATOGRAM

Retention time: 4.6 (R-(-)), 5.8 (S-(+))

Internal standard: S-(-)-ropivacaine (10.5)

Limit of detection: 3 ng/mL

Limit of quantitation: 5 ng/mL

OTHER SUBSTANCES

Noninterfering: bupivacaine

Interfering: lidocaine, prilocaine

KEY WORDS

plasma; chiral; pharmacokinetics

REFERENCE

Vletter, A.A.; Olieman, W.; Burm, A.G.L.; Groen, K.; van Kleef, J.W. High-performance liquid chromatographic assay of mepivacaine enantiomers in human plasma in the nanogram per milliliter range, *J. Chromatogr. B*, 1996, 678, 369–372.

SAMPLE

Matrix: blood

Sample preparation: Condition a 1 mL Bond-Elut C18 SPE cartridge with 1 mL 100 mM HCl, 1 mL MeOH and 1 mL water. Add 350 μ L 100 mM pH 4.5 phosphate buffer and 40 μ L 200 μ g/mL S-bupivacaine in MeOH to 1 mL serum, vortex. Add the mixture to the SPE cartridge and pass it through the cartridge attached to a vacuum manifold, wash 4 times with 250 μ L water, then wash with 500 μ L MeCN and allow to air dry for 1 min between each wash. Elute with four 250 μ L aliquots of MeOH containing 2% HCl, evaporate the eluate to dryness under a stream of nitrogen, reconstitute the residue in 1 mL mobile phase and inject a 100 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4 5 μ m Sumichiral OA-4700 (YMC, Wilmington, NC, USA)

Mobile phase: Hexane:dichloroethane:absolute MeOH 85:10:5

Flow rate: 0.8

Injection volume: 100

Detector: UV 220

CHROMATOGRAM

Retention time: 9.3 (R-+), 11.4 (S-(-))

Internal standard: S-bupivacaine (7.7)

Limit of detection: 100 ng/mL

Limit of quantitation: 150 ng/mL

KEY WORDS

chiral; SPE; serum

REFERENCE

Siluveru,M.; Stewart,J.T. Stereoselective determination of mepivacaine in human serum using a bush-type chiral stationary phase and solid-phase extraction, *J.Chromatogr.B*, **1997**, 690, 359–362.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 50 μ L 40 μ g/mL etidocaine hydrochloride in water + 100 μ L 1 M NaOH, vortex for 15 s, add 5 mL diethyl ether, shake on a reciprocating shaker for 10 min, centrifuge. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 100 μ L mobile phase, inject a 80 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 3.2 10 μ m μ Bondapak C18

Mobile phase: MeCN:50 mM pH 5.80 Na₂HPO₄ 25:75

Flow rate: 0.9

Injection volume: 80

Detector: UV 210

CHROMATOGRAM

Retention time: 3.7

Internal standard: etidocaine (12.0)

OTHER SUBSTANCES

Extracted: 2,6-pipecolylxylidine, bupivacaine, lidocaine

Noninterfering: metabolites, 2,3-chloroprocaine, theophylline, mexiletine, quinidine, disopyramide, verapamil, phenobarbital, phenytoin, carbamazepine, ethosuximide, digoxin, theobromine, caffeine, furosemide, phenprocoumon, aldactone

KEY WORDS

plasma

REFERENCE

Ha,H.-R.; Funk,B.; Gerber,H.R.; Follath,F. Determination of bupivacaine in plasma by high-performance liquid chromatography, *Anesth.Analg.*, **1984**, 63, 448–450.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 100 μ L 1 M NaOH + 3 mL heptane:ethyl acetate 90:10, shake for 2 min, centrifuge at 1200 g for 10 min. Remove the organic phase and add it to 50 μ L 50 mM sulfuric acid, shake for 2 min, centrifuge at 1200 g for 5 min. Remove the aqueous phase and add it to 820 μ g sodium acetate, inject a 40 μ L aliquot. (The sodium acetate was measured out by adding 50 μ L 200 mM sodium acetate in MeOH to the tube and evaporating the MeOH.)

HPLC VARIABLES

Column: 250 \times 4 10 μ m μ Bondapak C18

Mobile phase: MeCN:10 mM NaH₂PO₄ 7:93, adjusted to pH 2.1

Column temperature: 30

Flow rate: 1

Injection volume: 40

Detector: UV 205

CHROMATOGRAM

Retention time: 15

Limit of detection: 5 ng/mL

OTHER SUBSTANCES

Extracted: articaine

KEY WORDS

plasma; rabbit

REFERENCE

Le Guévello,P.; Le Corre,P.; Chevanne,P.; Le Verge,R. High-performance liquid chromatographic determination of bupivacaine in plasma samples for biopharmaceutical studies and application to seven other local anaesthetics, *J.Chromatogr.*, **1993**, 622, 284–290.

SAMPLE

Matrix: blood

Sample preparation: Condition a 1 mL BondElut C18 SPE cartridge once with 1 M HCl, twice with MeOH, and once with water, remove the liquid completely with suction each time. Add 250 μ L bupivacaine in 100 mM NaH₂PO₄ and 250 μ L serum to the column at 1 mL/min, wash twice with water and once with MeCN draining the column completely after each wash, elute with 250 μ L eluting solution, centrifuge for 20 s to remove last of eluate, inject a 5 μ L aliquot of the eluate. (Eluting solution was 2.5 mL 35% perchloric acid in 100 mL MeOH.)

HPLC VARIABLES

Guard column: 15 \times 3.2 7 μ m RP-8 (Applied Biosystems)

Column: 150 \times 4.6 5 μ m Ultrasphere octyl

Mobile phase: MeCN:10 mM KH₂PO₄ 25:80, pH 5.2

Flow rate: 1.5

Injection volume: 5

Detector: UV 205

CHROMATOGRAM

Retention time: 2.6

Internal standard: bupivacaine (7.2)

OTHER SUBSTANCES

Extracted: bupivacaine, meperidine, fentanyl

Noninterfering: acetaminophen, codeine, epinephrine, morphine, diazepam

KEY WORDS

serum; SPE

REFERENCE

Gupta,R.N.; Dauphin,A. Column liquid chromatographic determination of bupivacaine in human serum using solid-phase extraction, *J.Chromatogr.B*, **1994**, 658, 113–119.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol:n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 μ L mobile phase, centrifuge at 2800 g for 5 min, inject a 50 μ L aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 264

CHROMATOGRAM

Retention time: 4.54

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfapyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; ace-nocoumarol; vandesine; mexiletine; dipyrindamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpi-pramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254–262.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 125 × 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2**Injection volume:** 20**Detector:** E, LeCarbone, V25 glassy carbon electrode, + 1.2 V**CHROMATOGRAM****Retention time:** 1.7**OTHER SUBSTANCES**

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzoctamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipiprone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenylglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindine, pimozide, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thiopropazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranylecypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, *323*, 191–225.

SAMPLE**Matrix:** solutions**HPLC VARIABLES****Column:** 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30**Flow rate:** 2**Detector:** UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, meclizine, meclorfenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephentoin, mephesisin, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methyprylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidone, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfafethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, **1994**, *18*, 233–242.

SAMPLE

Matrix: solutions

Sample preparation: Inject an aliquot of a 200 μ M solution in MeOH.

HPLC VARIABLES

Column: 100 \times 4.7 7 μ m Hypercarb (Shandon)

Mobile phase: MeOH containing 5 mM N-benzoyloxycarbonylglycyl-L-proline and 4.5 mM NaOH

Column temperature: 17

Injection volume: 20

Detector: UV 270

CHROMATOGRAM

Retention time: k' 2.32 (first enantiomer)

KEY WORDSchiral; $\alpha = 1.11$ **REFERENCE**

Huynh, N.-H.; Karlsson, A.; Pettersson, C. Enantiomeric separation of basic drugs using N-benzoyloxycarbonyl-glycyl-L-proline as counter ion in methanol, *J. Chromatogr. A*, **1995**, 705, 275–287.

SAMPLE**Matrix:** solutions**HPLC VARIABLES****Column:** 250 × 4.6 5 μ m Supelcosil LC-DP (A) or 250 × 4 5 μ m LiChrospher 100 RP-8 (B)**Mobile phase:** MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)**Flow rate:** 0.6**Injection volume:** 25**Detector:** UV 229**CHROMATOGRAM****Retention time:** 7.71 (A), 4.20 (B)**OTHER SUBSTANCES**

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordiazepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenopropfen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazinol, mefenamic acid, meperidine, mephentoin, mesoridazine, metaproterenol, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimizide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfipyrazole, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocainide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, trifluorpromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone


KEY WORDS

details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, 692, 103–119.

Meprobamate



Molecular formula: C₉H₁₈N₂O₄

Molecular weight: 218.25

CAS Registry No.: 57-53-4

Merck Index: 5908

Lednicer No.: 1 218

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Serum + 100 μ L 200 mM HCl + 200 μ L 100 μ g/mL carisoprodol + 1 mL chloroform, shake for 2 min, sonicate for 1 min, centrifuge for 1 min. Remove the organic layer and add it to 200 mM NaOH, shake for 1 min, sonicate for 1 min, centrifuge for 1 min. Remove 750 μ L of the organic layer and evaporate it to dryness, dissolve the residue in 50 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 3 octyl CP-tm-Spher C8 glass column (Chrompack)

Mobile phase: MeCN:50 mM NaH₂PO₄ 35:65 adjusted to pH 2.2 with phosphoric acid

Flow rate: 0.8

Injection volume: 20

Detector: UV 190

CHROMATOGRAM

Retention time: 1.8

Internal standard: carisoprodol (5.5)

Limit of detection: 2000 ng/mL

KEY WORDS

serum

REFERENCE

Van Damme, M.; Molle, L.; Abi Khalil, F. Useful sample handlings for reversed phase high performance liquid chromatography in emergency toxicology, *J. Toxicol. Clin. Toxicol.*, **1985**, 23, 589–614.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Serum + 100 μ L 200 mM HCl + 200 μ L 100 μ g/mL carisoprodol + 1 mL chloroform, shake for 2 min, sonicate for 1 min, centrifuge for 1 min. Remove the organic layer and add it to 200 mM NaOH, shake for 1 min, sonicate for 1 min, centrifuge for 1 min. Remove 750 μ L of the organic layer and evaporate it to dryness, dissolve the residue in 50 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 3 octyl CP-tm-Spher C8 glass column (Chrompack)

Mobile phase: MeCN:50 mM NaH₂PO₄ 35:65 adjusted to pH 2.2 with phosphoric acid

Flow rate: 0.8

Injection volume: 20

Detector: UV 190

CHROMATOGRAM

Retention time: 1.8

Internal standard: carisoprodol (5.5)

Limit of detection: 2000 ng/mL

KEY WORDS

serum

REFERENCE

Hormazabal,V.; Steffanak,I.; Yndestad,M. Simultaneous extraction and determination of sulfadiazine and trimethoprim in medicated fish feed by high-performance liquid chromatography, *J.Chromatogr.*, **1993**, 648, 183–186.

SAMPLE

Matrix: formulations

Sample preparation: Powder tablets, weigh out amount equivalent to 100 mg meprobamate, dissolve in 50 mL mobile phase, sonicate for 10 min, filter, inject an aliquot.

HPLC VARIABLES

Column: 125 × 4.5 µm Lichrospher RP-18

Mobile phase: MeOH:50 mM pH 1.9 phosphoric acid 30:70 containing 1 mM benzoic acid

Column temperature: 35

Flow rate: 0.9

Injection volume: 50

Detector: UV 273

CHROMATOGRAM

Retention time: 13.5

KEY WORDS

tablets; indirect UV detection

REFERENCE

Bechet,I.; Ceccato,A.; Hubert,P.; Herne,P.; Crommen,J. Determination of meprobamate in pharmaceutical dosage forms also containing carbromal by liquid chromatography and indirect photometric detection, *J.Pharm.Biomed.Anal.*, **1992**, 10, 995–999.


Meptazinol

Molecular formula: C₁₅H₂₃NO

Molecular weight: 233.35

CAS Registry No.: 54340-58-8, 59263-76-2 (HCl)

Merck Index: 5910

**SAMPLE**

Matrix: blood

Sample preparation: 1 mL Plasma + 20 µL 2.5 µg/mL IS in MeOH + 100 µL buffer + 3 mL ethyl acetate, vortex for 2 min, centrifuge at 2000 rpm for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 50°, reconstitute the residue in 100 µL MeOH, add 50 µL water, inject an aliquot. (Buffer was 1 M Na₂HPO₄ adjusted to pH 9.5 with 5 M NaOH.)

HPLC VARIABLES

Column: 250 × 4.6 µm Spherisorb CN

Mobile phase: MeOH:0.5% ammonium acetate 45:55

Flow rate: 2

Detector: F ex 282 em 300

CHROMATOGRAM

Retention time: 5.8

Internal standard: m-(1-cyclopropylmethyl-3-ethylhexahydro-1H-azepin-3-yl)phenol (7.0)

Limit of detection: 3 ng/mL

KEY WORDS

plasma

REFERENCE

Frost, T. Determination of meptazinol in plasma by high-performance liquid chromatography with fluorescence detection, *Analyst*, **1981**, *106*, 999–1000.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Plasma + 100 μ L 40 ng/mL fenethazine in 2 M aqueous Tris + 200 μ L MTBE, vortex for 30 s, centrifuge at 9950 g for 2 min, inject a 100 μ L aliquot of the organic phase.

HPLC VARIABLES

Column: 125 \times 5 μ m Spherisorb S5W silica

Mobile phase: MeOH:glacial acetic acid:ammonia 996:3:1 (Ammonia was 0.88 g/mL.)

Flow rate: 2

Injection volume: 100

Detector: E, EDT Research LCA 15, glassy carbon electrode +1.2 V, Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 3

Internal standard: fenethazine (4)

Limit of detection: 500 ng/mL

KEY WORDS

plasma; cow; human

REFERENCE

Storey, G.C.A.; Schootstra, R.; Henry, J.A. Measurement of meptazinol in plasma by high-performance liquid chromatography with electrochemical detection, *J. Chromatogr.*, **1985**, *341*, 223–227.

SAMPLE

Matrix: microsomal incubations, urine

Sample preparation: Microsomal incubations. 2 mL Microsomal incubation + IS + 1 mL 1 M pH 10.0 ammonia/ammonium chloride buffer + 5 mL n-hexane:ethyl acetate 50:50, extract. Remove 4 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 50 μ L water, inject a 20 μ L aliquot. Urine. Untreated or deconjugated (with 500–2000 IU/mL β -glucuronidase/sulfatase from *Helix pomatia* for 24 h) urine + IS + 1 mL 1 M pH 10.0 ammonia/ammonium chloride buffer + 5 mL n-hexane:ethyl acetate 30:70, stir vigorously for 10 min. Remove 4 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 50 μ L water, inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: 30 \times 4 μ m LiChrospher 60 CN

Column: 250 \times 4 μ m LiChrospher 100 CN

Mobile phase: MeCN:MeOH:1% triethylammonium acetate buffer (pH 5.5) 15:5:80

Flow rate: 0.5

Injection volume: 20

Detector: UV 220

CHROMATOGRAM

Retention time: 31

Internal standard: 3-(3-ethylhexahydro-1-ethyl-1H-azepin-3-yl)phenol (35)

Limit of quantitation: 392 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

rabbit; liver; human; pharmacokinetics; rat

REFERENCE

Rudolphi, C.; Blaschke, G. Determination of the stereoselective aspects in in-vitro and in-vivo metabolism of the analgesic meptazinol by high-performance liquid chromatography, *J. Chromatogr. B*, **1995**, *663*, 315–326.


SAMPLE**Matrix:** solutions**Sample preparation:** Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.**HPLC VARIABLES****Column:** 125 × 4.9 Spherisorb S5W silica**Mobile phase:** MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7**Flow rate:** 2**Injection volume:** 20**Detector:** E, LeCarbone, V25 glassy carbon electrode, + 1.2 V**CHROMATOGRAM****Retention time:** 3.6**OTHER SUBSTANCES**

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazepine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipranone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, mepivacaine, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolitane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, *323*, 191–225.

Mercaptobenzothiazole



Molecular formula: C₇H₅NS₂

Molecular weight: 167.26

CAS Registry No.: 149-30-4

Merck Index: 5916

SAMPLE

Matrix: blood, tissue, urine

Sample preparation: Urine. 1 mL Urine + 4 mL EtOH, centrifuge at 3000 rpm for 10 min. Dilute the supernatant with EtOH, filter (0.2 μ m, Acro LC13), inject a 20 μ L aliquot. (Hydrolyze conjugates by mixing 100 μ L urine with 900 μ L 100 mM pH 5 acetate buffer and 33 units sulfatase or 5000 units β -glucuronidase, heat at 37° for 30 min, proceed as above.) Tissue. Homogenize liver or kidney in 3 volumes of cold 100 mM pH 5 acetate buffer containing 2 mM aminooxyacetic acid, add 4 volumes of EtOH, centrifuge at 3000 rpm for 20 min. Filter (0.2 μ m, Acro LC13) the supernatant and inject a 20 μ L aliquot. Plasma. Add 4 volumes of EtOH to plasma, centrifuge at 3000 rpm for 20 min. Filter (0.2 μ m, Acro LC13) the supernatant and inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: 30 mm long C18

Column: 250 \times 4.6 Partisil 10 ODS 3

Mobile phase: MeCN:water 50:50 adjusted to pH 4.5 with trifluoroacetic acid

Flow rate: 1

Injection volume: 20

Detector: UV 321

CHROMATOGRAM

Retention time: 7.0

Limit of detection: 200 nM

KEY WORDS

rat; hamster; guinea pig; mouse; plasma; liver; kidney

REFERENCE

Elfarra, A.A.; Hwang, I.Y. *In vivo* metabolites of S-(2-benzothiazolyl)-L-cysteine as markers of *in vivo* cysteine conjugate β -lyase and thiol glucuronosyl transferase activities, *Drug Metab. Dispos.*, **1990**, 18, 917–922.

SAMPLE

Matrix: bulk, formulations

Sample preparation: Bulk. Weigh out 100 mg morphine, dissolve in 25 mL MeOH:water:acetic acid 24:72:1, dilute with MeOH:water:acetic acid 24:72:1 to a final concentration of 240 μ g/mL morphine, filter (0.45 μ m), inject a 20 μ L aliquot. Injections. Dilute with MeOH:water:acetic acid 24:72:1 to a final concentration of 240 μ g/mL morphine, filter (0.45 μ m), inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: MeOH:5 mM sodium 1-heptanesulfonate:acetic acid 24:72:1

Flow rate: 1.5

Injection volume: 20

Detector: UV 230

CHROMATOGRAM

Retention time: 16

OTHER SUBSTANCES

Simultaneous: morphine (UV 284), phenol (UV 284), pseudomorphine (UV 230)

KEY WORDSinjections

REFERENCE

Bello,A.C.; Jhangiani,R.K. Liquid chromatographic determination of morphine sulfate and some contaminants in injections and bulk drug material: collaborative study, *J.Assoc.Off.Anal.Chem.*, **1988**, 71, 1046–1048.

SAMPLE**Matrix:** formulations**Sample preparation:** Directly inject a 20 μL aliquot of a 250 $\mu\text{g/mL}$ digoxin injection.

HPLC VARIABLES**Column:** 300 \times 3.9 10 μm $\mu\text{Bondapak C18}$ **Mobile phase:** MeCN:water 29:71**Flow rate:** 2**Injection volume:** 20**Detector:** UV 218

CHROMATOGRAM**Retention time:** 7

OTHER SUBSTANCES**Simultaneous:** digoxin

KEY WORDSinjections

REFERENCE

Reepmeyer,J.C.; Juhl,Y.H. Contamination of injectable solutions with 2-mercaptopbenzothiazole leached from rubber closures, *J.Pharm.Sci.*, **1983**, 72, 1302–1305.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 100 \times 3 5 μm ChromSpher C18 (Chrompack)**Mobile phase:** MeOH:water:acetic acid 45:54:1**Flow rate:** 0.4**Injection volume:** 75**Detector:** UV 321

CHROMATOGRAM**Retention time:** 5.5**Limit of detection:** 5 μM

REFERENCE

Stijntjes,G.J.; te Koppele,J.M.; Vermeulen,N.P.E. High-performance liquid chromatography-fluorescence assay of pyruvic acid to determine cysteine conjugate β -lyase activity: application to S-1,2-dichlorovinyl-L-cysteine and S-2-benzothiazolyl-L-cysteine, *Anal.Biochem.*, **1992**, 206, 334–343.

SAMPLE**Matrix:** solutions**Sample preparation:** Prepare a solution in MeCN, inject a 50 μL aliquot.

HPLC VARIABLES**Column:** 250 \times 4.6 10 μm Spheri-10 RP-18**Mobile phase:** MeCN:water:acetic acid 30:70:0.03**Flow rate:** 1**Injection volume:** 50

Detector: UV 328

CHROMATOGRAM

Retention time: 2

Limit of detection: 500 ng

REFERENCE

Gaind,V.S.; Jedrzejczak,K. HPLC determination of rubber septum contaminants in the iodinated intravenous contrast agent (sodium iothalamate), *J.Anal.Toxicol.*, **1993**, 17, 34–37.

SAMPLE

Matrix: solutions

Sample preparation: Mix 500 μ L of a solution in MeOH:water 90:10 with 500 μ L 50 μ M CY5.4a-IA in MeOH:water 70:30, pass dry nitrogen through the mixture for 1 min, heat at 65° for 1.5 h, inject a 25 μ L aliquot. (Some details for the synthesis of CY5.4a-IA are given in the paper.)

HPLC VARIABLES

Column: 150 \times 3.1 5 μ m LiChrosorb RP-8

Mobile phase: MeOH:10 mM pH 6.8 phosphate buffer 65:35 containing 1 mM triethylamine

Flow rate: 0.75

Injection volume: 25

Detector: F ex 670 (9.5 mW Lasermix LAS200-670-10 diode laser)

CHROMATOGRAM

Retention time: 7

Limit of detection: 1 nM


KEY WORDS

derivatization

REFERENCE

Mank,A.J.G.; Molenaar,E.J.; Lingeman,H.; Gooijer,C.; Brinkman,U.A.T.; Velthorst,N.H. Visible diode laser induced fluorescence detection in liquid chromatography after precolumn derivatization of thiols, *Anal.Chem.*, **1993**, 65, 2197–2203.

Mercaptopurine



Molecular formula: C₅H₄N₄S

Molecular weight: 152.18

CAS Registry No.: 50-44-2, 6112-76-1 (monohydrate)

Merck Index: 5919

SAMPLE

Matrix: blood

Sample preparation: Condition a 500 mg Bakerbond SPE cartridge packed with LiChrosorb RP-18 with 2 mL MeOH, 2 mL water and 1 mL MeOH. Make up 500 μ L plasma to 3 mL with MeOH. Centrifuge at 1100 g for 15 min, add 1.5 mL supernatant to the SPE cartridge, elute at 50 μ L/min flow-rate, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 200 \times 4 7 μ m Lichrosorb RP-18

Mobile phase: MeOH: pH 4.15 phosphate buffer 20:80

Flow rate: 1

Injection volume: 20

Detector: UV 262

CHROMATOGRAM**Retention time:** 3.30**Internal standard:** mercaptopurine

KEY WORDSplasma; SPE; mercaptopurine is IS

REFERENCE

Misztal, G.; Paw, B. Determination of fludarabine phosphate in human plasma using reversed phase high-performance liquid chromatography, *Pharmazie*, **1996**, *51*, 733–734.

SAMPLE**Matrix:** blood

Sample preparation: 2 mL Blood + 120 µg dithiothreitol (to prevent oxidation of 6-mercaptopurine), centrifuge at 2000 g for 5 min. Cool plasma in ice, add a volume of ice-cold 50% trichloroacetic acid equal to 10% of the plasma volume, mix vigorously, keep on ice for 10 min, centrifuge at 2000 g for 10 min. Remove the supernatant and adjust the pH to 6–7 with 4 M KOH, inject a 460 µL aliquot.

HPLC VARIABLES**Column:** Two 250 × 4.6 10 µm Spherisorb 10-ODS columns in series**Mobile phase:** 50 mM pH 6.35 potassium phosphate buffer**Flow rate:** 1.5**Injection volume:** 460**Detector:** UV 312

CHROMATOGRAM**Retention time:** 9**Limit of detection:** 3 ng/mL

OTHER SUBSTANCES**Extracted:** metabolites

KEY WORDSplasma; pharmacokinetics; dog

REFERENCE

De Abreu, R.A.; van Baal, J.M.; Schouten, T.J.; Schretlen, E.D.A.M.; de Bruyn, C.H.M.M. High-performance liquid chromatographic determination of plasma 6-mercaptopurine in clinically relevant concentrations, *J.Chromatogr.*, **1982**, *227*, 526–533.

SAMPLE**Matrix:** blood

Sample preparation: 1 mL Plasma + 25 µL 4 µg/mL 6-thioguanine in water + 100 µL 400 mM NaOH + 1 mL 0.3% phenylmercuric acetate in ethyl acetate + 3 mL diethyl ether, shake on a tumble mixer for 10 min, centrifuge for 5 min. Remove the organic layer and add it to 500 µL 100 mM HCl, whirlmix for 2 min, centrifuge for 5 min, discard the organic layer, evaporate traces of organic solvent under a stream of nitrogen at room temperature for 15 min, add 10 µL 3 mg/mL dithioerythritol in water, inject an aliquot.

HPLC VARIABLES**Column:** 250 × 4.6 LiChrosorb 10 RP-18**Mobile phase:** Isopropanol:water 3:97 containing 13.80 g/L NaH₂PO₄·H₂O, 200 µL/L 85% phosphoric acid, 60 mg/L dithioerythritol, and 500 mg/L sodium octanesulfonate, pH 3.6–3.7**Flow rate:** 1.5

Detector: F ex 295 em 380 following post-column reaction. The column effluent mixed with 8 mM potassium chromate in 500 mM HCl pumped at 0.16 mL/min and with air flowing at 0.32 mL/min and the mixture flowed through a single mixing coil. The effluent from this coil mixed with 1.6% sodium metabisulfite pumped at 0.16 mL/min and this mixture flowed through a single mixing coil. The effluent from this coil mixed with 4 M ammonium hydroxide pumped

at 0.23 mL/min and this mixture flowed through a double mixing coil to a debubbler. The liquid effluent from the debubbler flowed to the detector.

CHROMATOGRAM

Retention time: 6

Internal standard: 6-thioguanine (8)

Limit of detection: <2 ng/mL

KEY WORDS

post-column reaction; plasma; pharmacokinetics

REFERENCE

Jonkers,R.E.; Oosterhuis,B.; ten Berge,R.J.M.; van Bortel,C.J. Analysis of 6-mercaptopurine in human plasma with a high-performance liquid chromatographic method including post-column derivatization and fluorimetric detection, *J.Chromatogr.*, **1982**, 233, 249–255.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 80 μ L 10 μ g/mL 6-thioguanine in water + 10 μ L 1 M dithiothreitol, vortex for 10 s, add 2 mL MeCN, vortex for 30 s, centrifuge at 2000 g for 5 min. Remove the supernatant and add it to 2 mL dichloromethane, shake on a reciprocating shaker for 5 min, centrifuge at 2000 g for 5 min. Remove 750 μ L from the top aqueous layer and evaporate it to dryness under a stream of nitrogen at 37°, reconstitute the residue in 150 μ L distilled water, vortex for 1 min, inject a 15 μ L aliquot.

HPLC VARIABLES

Guard column: 70 \times 2.2 30-38 μ m Co:Pell ODS

Column: 250 \times 4.6 5 μ m Ultrasphere ODS

Mobile phase: MeCN:acetic acid:water 3.5:0.2:96.3

Flow rate: 1.4

Injection volume: 15

Detector: UV 322

CHROMATOGRAM

Retention time: 4.8

Internal standard: 6-thioguanine (6.5)

Limit of detection: 5 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

Noninterfering: caffeine, cytarabine, 5-fluorouracil, prednisone, theophylline, vinblastine, vincristine

KEY WORDS

plasma; protect from light; monkey; pharmacokinetics; human

REFERENCE

Narang,P.K.; Yeager,R.L.; Chatterji,D.C. Quantitation of 6-mercaptopurine in biologic fluids using high-performance liquid chromatography: a selective and novel procedure, *J.Chromatogr.*, **1982**, 230, 373–380.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Serum + 5 μ L 50 μ g/mL 2-ethyl-4-oxoquinazoline in EtOH + 100 μ L reagent, let stand at room temperature for 1 h, add 1.8 mL ethyl acetate, mix, centrifuge at 1800 g for 5 min, remove 1.5 mL of the supernatant, repeat the extraction. Combine the organic layers and evaporate them to dryness under reduced pressure below 30°, reconstitute the residue in 100 μ L initial mobile phase, inject a 90 μ L aliquot. (Reagent was 30 mg N-ethylmaleimide in 2 mL 50 mM pH 7.0 phosphate buffer, prepare fresh daily.)

HPLC VARIABLES

Column: 10 μ m μ Bondapak C18

Mobile phase: Gradient. MeCN:10 mM KH₂PO₄ 9:91 for 26 min, then 50:50 for 1 min (step gradient).

Flow rate: 1.5

Injection volume: 90

Detector: UV 280

CHROMATOGRAM

Retention time: 18.6

Internal standard: 2-ethyl-4-oxoquinazoline (28)

Limit of quantitation: 10 ng/mL

OTHER SUBSTANCES

Extracted: azathioprine

KEY WORDS

serum; derivatization

REFERENCE

Tsutsumi,K.; Otsuki,Y.; Kinoshita,T. Simultaneous determination of azathioprine and 6-mercaptopurine in serum by reversed-phase high-performance liquid chromatography, *J.Chromatogr.*, **1982**, 231, 393–399.

SAMPLE

Matrix: blood

Sample preparation: 250 µL Plasma + 5 µL 2.5 µg/mL 6-thioguanine + 25 µL 1 M aluminum perchlorate in water, let stand at room temperature for 15 min, chill in ice water for 15 min, centrifuge at 15600 g for 15 min, discard the supernatant. Suspend the precipitate in 500 µL 50 mM aluminum perchlorate in water by stirring to break up the precipitate, vortex for 20 s, centrifuge at 15600 g for 15 min, discard the supernatant. Suspend the precipitate in 150 µL 400 mM perchloric acid, add 5 µL freshly prepared 200 mM aqueous sodium hydrosulfite, mix, let stand at room temperature for 30 min, chill in ice water, centrifuge at 15600 g for 15 min, inject a 50 µL aliquot of the supernatant.

HPLC VARIABLES

Guard column: 45 × 2 40 µm ODS

Column: 150 × 4.3 5 µm Ultrasphere ODS

Mobile phase: Water:85% phosphoric acid 99.32:0.68 containing 154.3 mg/L dithiothreitol

Flow rate: 1

Injection volume: 50

Detector: UV 340

CHROMATOGRAM

Retention time: 6.08

Internal standard: 6-thioguanine (4.73)

Limit of quantitation: 3 ng/mL

OTHER SUBSTANCES

Extracted: metabolites, 6-thiouric acid

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Lin,K.T.; Varin,F.; Rivard,G.E.; Leclerc,J.M. Isolation of 6-mercaptopurine in human plasma by aluminum ion complexation for high-performance liquid chromatographic analysis, *J.Chromatogr.*, **1991**, 536, 349–355.

SAMPLE

Matrix: blood

Sample preparation: 200 µL Plasma + 800 µL MeOH, stir for 15 s, let stand for 20 min, centrifuge at 5° at 3000 rpm for 10 min. Remove 800 µL of the supernatant and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 100 µL mobile phase, inject a 40 µL aliquot.

HPLC VARIABLES

Column: 150 × 4.6 ODS-II (Shimadzu)

Mobile phase: MeOH:water 7:100 containing 0.2% glacial acetic acid and 0.1% sodium 1-heptanesulfonate

Flow rate: 1

Injection volume: 40

Detector: UV 325

CHROMATOGRAM

Retention time: 6.8

Limit of detection: 2 ng/mL

KEY WORDS

plasma; pharmacokinetics; rat

REFERENCE

Takeichi, Y.; Kimura, T. Improvement of aqueous solubility and rectal absorption of 6-mercaptopurine by addition of sodium benzoate, *Biol. Pharm. Bull.*, **1994**, *17*, 1391–1394.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 50 µL 2.8 µg/mL 6-mercaptopurine arabinoside + 2–3 mg mercurial cellulose, vortex for 30 s, centrifuge at 550 g for 5 min. Resuspend the pellet in 2 mL phosphate buffered saline, centrifuge, repeat this process twice. Suspend the pellet in 250 µL 20 mM 2-mercaptoethanol, centrifuge, inject a 100 µL aliquot of the supernatant. (Prepare mercurial cellulose as follows. Add 40 g Sigmacell Type 20 cellulose (20 µm) to 175 mL 40% NaOH at 5°, let stand at 0° for 2 h, add 40 mL allyl glycidyl ether dropwise, add 100 mL water, heat at 70–80° for 2 h. Suspend the cellulose in 4 L water, allow to settle, decant the water, wash several more times, filter, wash with water until the pH of the filtrate is about 8, wash with 1 L 10% acetic acid, wash 1 L with EtOH:water 95:5, wash with 1 L MeOH, wash with 1 L diethyl ether, air dry. Add to 400 mL 10% acetic acid containing 2.3 g mercuric acetate, stir at 60° for 1 h, filter, wash with 2 L 10% acetic acid, wash with 1.5 L water, wash with 300 mL EtOH:water 95:5, wash with 700 mL MeOH, wash with 1 L diethyl ether, air dry in the dark, store in the dark (*Anal. Biochem.* 1985, *144*, 514).)

HPLC VARIABLES

Column: 250 × 4.6 5 µm ODS (Beckman)

Mobile phase: MeCN:10 mM pH 3.0 sodium phosphate buffer 2:98

Flow rate: 1.2

Injection volume: 100

Detector: UV 323

CHROMATOGRAM

Retention time: 7.6

Internal standard: 6-mercaptopurine arabinoside (13.6)

Limit of detection: 0.5 ng/mL

OTHER SUBSTANCES

Noninterfering: aspirin, chloroquine, cyclosporin, diltiazem, nifedipine, prednisolone

KEY WORDS

plasma; pharmacokinetics; SPE

REFERENCE

Albertioni, F.; Pettersson, B.; Ohlman, S.; Peterson, C. Analysis of azathioprine and 6-mercaptopurine in plasma in renal transplant recipients after administration with oral azathioprine, *J. Liq. Chromatogr.*, **1995**, *18*, 3991–4005.

SAMPLE

Matrix: blood

Sample preparation: Add 1 volume ice-cold 8 M perchloric acid to 20 volumes plasma, mix, keep on ice for 10 min, centrifuge at 10 000 g for 15 min, remove the supernatant. Adjust the pH of the supernatant to 6-7 with 10 volumes ice-cold 4 M K_2HPO_4 , keep on ice for 10 min, centrifuge at 10 000 g for 5 min, inject a 100 μL aliquot of the supernatant.

HPLC VARIABLES

Column: 250 \times 4.6 5 μm Supelcosil LC-18-DB

Mobile phase: Gradient. A was 25 mM KH_2PO_4 , B was MeOH:50 mM KH_2PO_4 25:75. A:B from 98:2 to 96:4 over 5 min, to 85:15 over 3 min, to 80:20 over 2 min, to 40:60 over 10 min, maintain at 40:60 over 2 min, to 20:80 over 3 min, maintain at 20:80 for 20 min, return to initial conditions over 3 min, re-equilibrate for 12 min.

Flow rate: 1.25

Injection volume: 100

Detector: UV 320

CHROMATOGRAM

Retention time: 11.6

Limit of detection: 20-50 nM

OTHER SUBSTANCES

Extracted: metabolites, thioguanine (UV 342)

KEY WORDS

plasma

REFERENCE

Keuzenkamp-Jansen, S.W.; De Abreu, R.A.; Bökterink, J.P.M.; Trijbels, J.M.F. Determination of extracellular and intracellular thiopurines and methylthiopurines by high-performance liquid chromatography, *J.Chromatogr.B*, **1995**, 672, 53-61.

SAMPLE

Matrix: blood

Sample preparation: Treat whole blood with EDTA, freeze a 2 mL aliquot at -20° , thaw, add 2 mL PBS, centrifuge at 3000 g for 15 min, suspend the pellet in 200 μL PBS, treat with proteinase K at 70° for 10 min, add to a disposable spin column (Diagen QIAamp Blood Kit, Hilden, Germany), centrifuge for 1 min, wash twice with buffer for 1 min, elute with 200 μL 10 mM pH 9.0 Tris-HCl buffer containing 0.1 mM EDTA, heat at 100° for 5 min, cool in ice. Remove a 100 μL aliquot and add it to 10 μL buffer, add 20 μL 25 $\mu\text{g}/\text{mL}$ P_1 nuclease (Boehringer Mannheim) and 12.5 U/mL acid phosphatase (Sigma) in buffer:water 10:90, heat at 42° for 1 h, add 10 μL 400 mM formic acid, add 60 μL MeOH, add 1 μL 5 mM N-[6-(7-amino-4-methylcoumarin-3-acetamido)hexyl]-3'-(2'-pyridyldithio)propionamide (AMCA-HPDP, Pierce) in DMF, inject a 25 μL aliquot. (Buffer was 500 mM pH 4.5 sodium acetate buffer containing 10 mM magnesium chloride.)

HPLC VARIABLES

Guard column: 20 mm long Supelguard (Supelco)

Column: 150 \times 4.6 3 μm Supelcosil LC-8

Mobile phase: MeOH:buffer 37:63 (Between analyses wash column with MeOH:buffer 80:20 for 3 min. Buffer was 200 mM formic acid adjusted to pH 4.0 with 10 M NaOH.)

Column temperature: 45

Flow rate: 1

Injection volume: 25

Detector: F ex 345 em 450

CHROMATOGRAM

Retention time: 11.5 (as 2'-deoxy-6-thioguanosine metabolite)

Limit of detection: 60 pmole/g DNA

KEY WORDS

derivatization; whole blood

REFERENCE

Warren,D.J.; Andersen,A.; Slordal,L. Quantitation of 6-thioguanine residues in peripheral blood leukocyte DNA obtained from patients receiving 6-mercaptopurine-based maintenance therapy, *Cancer Res.*, **1995**, *55*, 1670-1674.

SAMPLE

Matrix: blood, CSF, urine

Sample preparation: Blood, CSF. Collect 2 mL blood in tubes containing heparin and 120 μ g dithiothreitol (DTT). Mix, centrifuge at 2000 g for 5 min, remove plasma. CSF. Collect 0.5 mL CSF in tubes containing 30 μ g DTT. Cool CSF and plasma samples on ice, add freshly prepared ice-cold 50% trichloroacetic acid equal to 10% of sample volume. Urine. Collect 2 mL urine in tubes containing 120 μ g DTT, filter (0.22 μ m), inject an aliquot.

HPLC VARIABLES

Column: two 250 \times 4.6 10 μ m Nucleosil 10 C18 columns in series

Mobile phase: Gradient. A was 25 mM pH 2.75 phosphoric acid. B was MeOH:water 50:50. C was 100 mM pH 6.6 KH₂PO₄. A:B:C from 100:0:0 to 98:2:0 over 5 min, to 30:3.5:66.5 over 5 min, maintain at 30:3.5:66.5 for 10 min, re-equilibrate at initial conditions for 10 min.

Column temperature: 33

Flow rate: 1.7

Injection volume: 195, 500

Detector: UV 312

CHROMATOGRAM

Retention time: 14

Limit of detection: 20 nM

OTHER SUBSTANCES

Extracted: metabolites, 6-mercaptopurine riboside, thioguanine (UV 342), 6-thioguanosine

KEY WORDS

plasma; goat; pharmacokinetics

REFERENCE

van Baal,J.M.; van Leeuwen,M.B.; Schouten,T.J.; De Abreu,R.A. Sensitive high-performance liquid chromatographic determination of 6-mercaptopurine, 6-thioguanine, 6-mercaptopurine riboside and 6-thioguanosine in biological fluids, *J.Chromatogr.*, **1984**, *336*, 422-428.

SAMPLE

Matrix: enzyme incubations

Sample preparation: To 237 μ L enzyme incubation add 850 μ L ice-cold 3.5 mM DL-dithiothreitol immediately followed by 500 μ L 1.5 M sulfuric acid, place tubes on ice. Equilibrate tubes to room temperature, heat at 100° for 2 h. Cool, add 500 μ L 3.4 M NaOH immediately followed by 8 mL toluene:amyl alcohol:phenyl mercury acetate mixture. Shake gently for 10 min, centrifuge at 10° at 900 g for 5 min. Transfer 6 mL toluene to another tube and add 200 μ L 100 mM HCl. Vortex for four 20 s periods, centrifuge at 10° at 900 g for 5 min, inject a 50 μ L aliquot of the aqueous layer. (Prepare toluene:amyl alcohol:phenyl mercury acetate mixture by adding enough phenyl mercury acetate to toluene containing 170 mM amyl alcohol so as to form a solution containing 1.3 mM phenyl mercury acetate, mix gently for 1 h, store in the dark (*J.Chromatogr.* 1992, 583, 83).)

HPLC VARIABLES

Guard column: 5 \times 4 5 μ m Resolve C18

Column: 100 \times 8 5 μ m Resolve C18 radial compression

Mobile phase: MeOH:water 20:80 containing 100 mM triethylamine and 0.5 mM dithiothreitol, pH adjusted to 3.2 with orthophosphoric acid (Add dithiothreitol immediately prior to use.)

Flow rate: 1.5

Injection volume: 50

Detector: UV 303

CHROMATOGRAM

Retention time: 4

OTHER SUBSTANCES

Extracted: methylmercaptopurine

REFERENCE

Lennard,L.; Singleton,H.J. High-performance liquid chromatographic assay of human red blood cell thiopurine methyltransferase activity, *J.Chromatogr.B*, **1994**, 661, 25–33.

SAMPLE

Matrix: formulations

Sample preparation: Dissolve crushed tablets or the freeze-dried compound for injection in 20 mM NaOH. Add a 10 mL aliquot of this solution (or a saline injection) to 10 mL 3 mg/mL theophylline in 20 mM NaOH, inject a 1.5 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 5 μ m ODS-Hypersil

Mobile phase: MeOH:25 mM KH_2PO_4 :glacial acetic acid 20:79:5 adjusted to pH 4.50 (Flush column with MeOH:water 60:40 at the end of each day.)

Flow rate: 1.5

Injection volume: 1.5

Detector: UV 240

CHROMATOGRAM

Retention time: 1.7

Internal standard: theophylline (3.5)

OTHER SUBSTANCES

Simultaneous: azathioprine, impurities, degradation products

KEY WORDS

stability-indicating; injections; tablets

REFERENCE

Fell,A.F.; Plag,S.M.; Neil,J.M. Stability-indicating assay for azathioprine and 6-mercaptopurine by reversed-phase high-performance liquid chromatography, *J.Chromatogr.*, **1979**, 186, 691–704.

SAMPLE

Matrix: solutions

Sample preparation: Inject a 25 μ L aliquot.

HPLC VARIABLES

Guard column: 40 \times 4 C18 Corasil II

Column: 300 \times 4 10 μ m μ Bondapak phenyl

Mobile phase: Propanol:6 mM C_{12} DAPS (Fluka) 3:97 (C_{12} DAPS is 3-(dimethyldodecylammonio)propanesulfonate.)

Injection volume: 25

Detector: UV 273

CHROMATOGRAM

Retention time: 3.3

OTHER SUBSTANCES

Simultaneous: albendazole, aminophylline, antipyrine, caffeine, dipropyline, flubendazole, metronidazole, nimorazole, procaine, β -hydroxytheophylline, theophylline, tinidazole

Interfering: albendazole sulfoxide, amyleine, theobromine

KEY WORDS

micellar chromatography

REFERENCE

Habel,D.; Guermouche,S.; Guermouche,M.H. Direct determination of theophylline in human serum by high-performance liquid chromatography using zwitterionic micellar mobile phase. Comparison with an enzyme multiplied immunoassay technique, *Analyst*, **1993**, *118*, 1511–1513.


Meropenem

Molecular formula: C₁₇H₂₅N₃O₅S

Molecular weight: 383.45

CAS Registry No.: 96036-03-2

Merck Index: 5960



SAMPLE

Matrix: bile, blood

Sample preparation: Condition a 100 mg Bond Elut C8 SPE cartridge with two 1 mL portions of MeOH and two 1 mL portions of 500 mM KH₂PO₄. 100 μ L Plasma or 50 μ L bile + 900 μ L water, mix, add to the SPE cartridge, wash with 1 mL 50 mM KH₂PO₄, elute with 800 (plasma) or 600 (bile) μ L MeOH:50 mM pH 6.0 phosphate buffer 10:90, inject a 100 μ L aliquot of the eluate.

HPLC VARIABLES

Column: 100 \times 4.6 3 μ m Hypersil ODS

Mobile phase: MeOH:buffer 15:85 (Buffer was 50 (plasma) or 57.5 (bile) mM pH 7.0 phosphate buffer.)

Flow rate: 1

Injection volume: 100

Detector: UV 296 (plasma), UV 310 (bile)

CHROMATOGRAM

Retention time: 5-5.5

Limit of detection: 280 ng/mL (bile), 100 ng/mL (plasma)

KEY WORDS

plasma; SPE

REFERENCE

Granai,F.; Smart,H.L.; Triger,D.R. A study of the penetration of meropenem into bile using endoscopic retrograde cholangiography, *J.Antimicrob.Chemother.*, **1992**, *29*, 711–718.

SAMPLE

Matrix: blood

Sample preparation: Dilute 50 μ L plasma with 75 μ L 50mM pH 7.0 phosphate buffer, inject a 100 μ L aliquot onto column A with mobile phase A, elute to waste with mobile phase A, after 4 min backflush the contents of column A onto column B with mobile phase B, monitor the effluent from column B. After 11 min re-equilibrate column A with mobile phase A and column B with mobile phase B for 5 min.

HPLC VARIABLES

Column: 20 \times 3.9 25-40 μ m LiChroprep RP-8; B 4.0 \times 10 Nova -Pak C8 + 150 \times 4.6 5 μ m Inertsil ODS

Mobile phase: A 50 mM pH 7.0 phosphate buffer; B MeCN:50 mM pH 7.0 phosphate buffer 6:94

Flow rate: 1

Injection volume: 100

Detector: UV 300

CHROMATOGRAM

Retention time: 6.5

Limit of detection: 100 ng/mL

KEY WORDS

rat; plasma; column-switching; pharmacokinetics

REFERENCE

Lee,H.S.; Shim,H.O.; Yu,S.R. High-performance liquid chromatographic determination of meropenem in rat plasma using column-switching, *Chromatographia*, **1996**, 42, 405–408.

SAMPLE

Matrix: blood

Sample preparation: Dilute serum with an equal volume of water, inject 20 μ L diluted serum onto column A, elute to waste with mobile phase A, after 1.1 min elute the contents of column A to column B with mobile phase B, elute with mobile phase B, monitor the effluent from column B.

HPLC VARIABLES

Column: A 50 \times 2.1 40 μ m Supelclean LC-NH₂ (Supelco); B 150 \times 4.6 3 μ m C18 (Supelco)

Mobile phase: A MeOH:10 mM pH 7.0 phosphate buffer 5:95; B MeOH:10 mM pH 7.0 phosphate buffer containing 5 mM tetrabutylammonium hydrogen sulfate 30:70

Flow rate: A 0.3; B 1

Injection volume: 20

Detector: UV 298

CHROMATOGRAM

Retention time: 4.2

Limit of detection: 100 ng/mL

OTHER SUBSTANCES

Noninterfering: acetaminophen, acyclovir, digoxin, fluconazole, theophylline, vancomycin

KEY WORDS

column-switching; serum

REFERENCE

Bompadre,S.; Ferrante,L.; de Martinis,M.; Leone,L. Determination of meropenem in serum by high-performance liquid chromatography with column switching, *J.Chromatogr.A*, **1998**, 812, 249–253.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Plasma + 900 μ L water, add to a conditioned Bond Elut C18 SPE cartridge, wash with 1 mL 50 mM KH₂PO₄, elute with 800 μ L mobile phase, inject an aliquot of the eluate.

HPLC VARIABLES

Column: 100 mm long 3 μ m Hypersil ODS

Mobile phase: MeOH:5 mM tetrabutylammonium dihydrogen phosphate 12:88

Flow rate: 1

Detector: UV 296

CHROMATOGRAM

Retention time: 9

Limit of detection: 200 ng/mL

KEY WORDS

SPE; plasma; pharmacokinetics

REFERENCE

Harrison,M.P.; Haworth,S.J.; Moss,S.R.; Wilkinson,D.M.; Featherstone,A. The disposition and metabolic fate of ¹⁴C-meropenem in man, *Xenobiotica*, **1993**, 23, 1311–1323.

SAMPLE

Matrix: blood

Sample preparation: Dilute plasma with water, prepare using a 1 mL C18 SPE cartridge.

HPLC VARIABLES

Column: 100 × 4.6 3 µm Hypersil C18

Mobile phase: MeOH:5 mM tetrabutylammonium dihydrogen phosphate 15:85

Detector: UV 296

CHROMATOGRAM

Limit of quantitation: 500 ng/mL

KEY WORDS

plasma; SPE; pharmacokinetics

REFERENCE

Bedikian,A.; Okamoto,M.P.; Nakahiro,R.K.; Farino,J.; Heseltine,P.N.R.; Appleman,M.D.; Yellin,A.E.; Berne,T.V.; Gill,M.A. Pharmacokinetics of meropenem in patients with intra-abdominal infections, *Antimicrob.Agents Chemother.*, **1994**, 38, 151–154.

SAMPLE

Matrix: blood, formulations, urine

Sample preparation: Dilute plasma 1:4, urine 1:10, and injections 1:100 with water, filter (0.6 µm), inject a 5-50 µL aliquot of the filtrate.

HPLC VARIABLES

Column: 200 × 4 5 µm Nucleosil C18

Mobile phase: MeOH:10 mM pH 7.4 potassium phosphate buffer 18:82 (plasma, injections) or 25:75 (urine)

Flow rate: 1

Injection volume: 5-50

Detector: UV 296

CHROMATOGRAM

Limit of detection: 5 µg/mL (urine), 500 ng/mL (plasma)

KEY WORDS

plasma; injections; saline; pharmacokinetics

REFERENCE

Burman,L.Å.; Nilsson-Ehle,I.; Hutchison,M.; Haworth,S.J.; Norrby,S.R. Pharmacokinetics of meropenem and its metabolite ICI 213,689 in healthy subjects with known renal metabolism of imipenem, *J.Antimicrob.Chemother.*, **1991**, 27, 219–224.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. Condition a Bond Elut C18 SPE cartridge. Add 100 µL plasma to the SPE cartridge, wash with 50 mM KH₂PO₄, elute with 800 µL eluent, inject an aliquot of the eluate. Urine. Dilute urine 1:10 with water, inject an aliquot directly. (Eluent was MeOH: 5 mM tetrabutylammonium dihydrogen phosphate 12:88.)

HPLC VARIABLES

Column: 100 × 4 3 µm Hypersil ODS

Mobile phase: MeOH:5 mM tetrabutylammonium dihydrogen phosphate 12:88 (plasma) or MeCN:10 mM pH 7.4 phosphate buffer 6:100 (urine)

Detector: UV 296

CHROMATOGRAM

Limit of detection: 60 ng/mL (plasma), 10 µg/mL (urine)

KEY WORDS

SPE; plasma; pharmacokinetics

REFERENCE

Bax,R.P.; Bastain,W.; Featherstone,A.; Wilkinson,D.M.; Hutchison,M. Haworth,S.J. The pharmacokinetics of meropenem in volunteers, *J.Antimicrob.Chemother.*, **1989**, 24, 311–320.

SAMPLE

Matrix: enzyme incubations

Sample preparation: Add 2 volumes of MeOH, mix well, centrifuge at 3000 g for 15 min, inject a 10 µL aliquot of the supernatant.

HPLC VARIABLES

Column: Inertsil ODS-2

Mobile phase: MeOH:5 mM tetrabutylammonium phosphate 30:70

Flow rate: 1

Injection volume: 10

Detector: UV 300

CHROMATOGRAM

Limit of quantitation: 1 µg/mL

REFERENCE

Hikida,M.; Kawashima,K.; Yoshida,M.; Mitsuhashi,S. Inactivation of new carbapenem antibiotics by dehydropeptidase-I from porcine and human renal cortex, *J.Antimicrob.Chemother.*, **1992**, 30, 129–134.

SAMPLE

Matrix: urine

Sample preparation: Dilute urine 1:5 with water, inject a 50 µL aliquot.

HPLC VARIABLES

Column: 50 mm long 3 µm Hypersil ODS

Mobile phase: MeCN:10 mM pH 7.4 phosphate buffer 6:100

Flow rate: 1

Injection volume: 50

Detector: UV 296

CHROMATOGRAM

Retention time: 5.5

Limit of detection: 600 ng/mL

KEY WORDS

pharmacokinetics

REFERENCE

Harrison,M.P.; Haworth,S.J.; Moss,S.R.; Wilkinson,D.M.; Featherstone,A. The disposition and metabolic fate of ¹⁴C-meropenem in man, *Xenobiotica*, **1993**, 23, 1311–1323.

SAMPLE

Matrix: urine

Sample preparation: Dilute urine with water, prepare using a 1 mL C18 SPE cartridge.

HPLC VARIABLES

Column: 150 × 4.6 5 µm Partisil silica

Mobile phase: 0.1% Phosphoric acid in water

Detector: UV 313

CHROMATOGRAM

Limit of quantitation: 1 µg/mL


KEY WORDS

SPE; pharmacokinetics

REFERENCE

Bedikian,A.; Okamoto,M.P.; Nakahiro,R.K.; Farino,J.; Heseltine,P.N.R.; Appleman,M.D.; Yellin,A.E.; Berne,T.V.; Gill,M.A. Pharmacokinetics of meropenem in patients with intra-abdominal infections, *Antimicrob.Agents Chemother.*, **1994**, *38*, 151–154.

Mesalamine

**Molecular formula:** C₇H₇NO₃**Molecular weight:** 153.14**CAS Registry No.:** 89-57-6**Merck Index:** 5964**SAMPLE****Matrix:** blood**Sample preparation:** 1 mL Plasma + 4 mL MeOH, let stand for 15 min at 4°, centrifuge at 3000 g for 15 min. Remove 1 mL of the supernatant and add it to 1 mL water, mix, inject an aliquot.**HPLC VARIABLES****Column:** 40 × 4.6 3 μm Hypersil**Mobile phase:** MeOH:water:200 mM pH 6.5 potassium phosphate buffer 40:40:20 containing 4 mM hexadecyltrimethylammonium bromide**Column temperature:** 34**Flow rate:** 1.7**Injection volume:** 20**Detector:** F ex 315 em 470**CHROMATOGRAM****Retention time:** 2.5**KEY WORDS**

plasma

REFERENCE

Tjornelund,J.; Hansen,S.H. Stability of 5-aminosalicylic acid and its metabolites in plasma at -20°C. Formation of N-β-D-glucopyranosyl-5-aminosalicylic acid, *J.Chromatogr.*, **1991**, *570*, 224–228.

SAMPLE**Matrix:** blood, feces, urine**Sample preparation:** Plasma. 1 mL Plasma + 2 mL MeOH, vortex for 20 s, let stand for 15 min at -20° and 10 min at room temperature, mix, centrifuge at 2000 g for 15 min. Remove the supernatant and add it to 5 mL 1,1,1-trichloroethane, shake for 2 min, centrifuge at 2000 g for 15 min, inject a 20 μL aliquot of the supernatant. Urine. 1 mL Urine + 4 mL MeOH, mix, let stand at -20° for 15 min, centrifuge at 2000 g for 15 min. Dilute an aliquot of the supernatant with an equal volume of water, inject a 20 μL aliquot. Feces. Extract with 500-1000 mL MeOH, centrifuge at 10000 g for 4 min. Dilute an aliquot of the supernatant with an equal volume of water, inject a 20 μL aliquot.**HPLC VARIABLES****Column:** 40 × 4.6 3 μm Hypersil**Mobile phase:** MeOH:water:200 mM pH 6.5 potassium phosphate buffer 40:40:20 containing 4 mM hexadecyltrimethylammonium bromide**Column temperature:** 34**Flow rate:** 1.7**Injection volume:** 20

Detector: F ex 315 em 470

CHROMATOGRAM

Retention time: 2.5

Limit of detection: 4 ng/mL

OTHER SUBSTANCES

Simultaneous: metabolites

KEY WORDS

plasma

REFERENCE

Tjornelund, J.; Hansen, S.H. High-performance liquid chromatographic assay of 5-aminosalicylic acid (5-ASA) and its metabolites N- β -D-glucopyranosyl-5-ASA, N-acetyl-5-ASA, N-formyl-5-ASA and N-butyryl-5-ASA in biological fluids, *J. Chromatogr.*, **1991**, 570, 109–117.

SAMPLE

Matrix: blood, urine

Sample preparation: Serum. Add disodium-3,5'-azo-bis-(6-hydroxybenzoate) to serum, treat with proteinase K (0.5 mg/mL protein) for 10 min, add tetrabutylammonium hydrogen sulfate buffered to pH 6.5, add dichloromethane, agitate for 30 min, centrifuge. Remove the organic layer and evaporate it to dryness, reconstitute the residue in mobile phase, inject an aliquot. Urine. Add 2,4-dihydroxybenzoic acid to urine, react with propionic anhydride for 5 min, add perchloric acid, add diethyl ether, shake for 10 min, freeze. Remove the organic layer and add it to pH 7.4 phosphate buffer, extract, inject an aliquot of the aqueous phase.

HPLC VARIABLES

Guard column: 30 \times 4 30-40 μ m Perisorb RP-18

Column: 250 \times 4 10 μ m Nucleosil C18

Mobile phase: MeOH:buffer 30:70 (Buffer was pH 7.4 phosphate buffer containing 20 mM tetrabutylammonium hydrogen sulfate.)

Detector: F ex 312 em 469

CHROMATOGRAM

Retention time: k' 2.4 (propionyl derivative)

Internal standard: disodium-3,5'-azo-bis-(6-hydroxybenzoate), 2,4-dihydroxybenzoic acid

Limit of quantitation: 6.1 μ M urine, 0.4 μ M (serum)

OTHER SUBSTANCES

Extracted: acetylaminosalicylic acid

KEY WORDS

serum; pharmacokinetics; derivatization

REFERENCE

Ryde, E.M.; Ahnfelt, N.-O. The pharmacokinetics of olsalazine sodium in healthy volunteers after a single i.v. dose and after oral doses with and without food, *Eur. J. Clin. Pharmacol.*, **1988**, 34, 481–488.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 500 μ L Plasma + 500 μ L 100 mg/mL (sic) 2,4-dihydroxybenzoic acid in MeOH, mix, centrifuge, inject an aliquot. Urine. 100 μ L Urine + 500 μ L 600 mg/mL (sic) 2,4-dihydroxybenzoic acid in MeOH, mix, centrifuge, inject an aliquot.

HPLC VARIABLES

Guard column: μ Bondapak C18 Guard Pak

Column: 250 \times 4.1 10 μ m Caulab CSC PRP-1

Mobile phase: MeCN:MeOH:50 mM pH 7.9 potassium phosphate buffer 10:10:80

Flow rate: 1

Detector: F ex 315 em 475

CHROMATOGRAM

Retention time: 5.1

Internal standard: 2,4-dihydroxybenzoic acid (3.3)

Limit of quantitation: 1000 ng/mL (urine), 25 ng/mL (plasma)

OTHER SUBSTANCES

Extracted: N-acetyl-5-aminosalicylic acid

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Corey,A.E.; Rose,G.M.; Conklin,J.D. Bioavailability of single and multiple doses of enteric-coated mesalamine and sulphasalazine, *J.Int.Med.Res.*, **1990**, 18, 441-453.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 205.2

CHROMATOGRAM

Retention time: 4.737

KEY WORDS

whole blood

REFERENCE

Gaillard,X.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

SAMPLE

Matrix: bulk

Sample preparation: Dissolve in dilute HCl (pH 2), sonicate if necessary, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Hypersil-BDS C18

Mobile phase: MeOH:THF:buffer 11:4:85 (Buffer was 8.6 g NaH₂PO₄·H₂O, 8.2 g NaCl, 2.2 g sodium 1-heptanesulfonate, and 5.75 mL 85% phosphoric acid in 2 L water, pH 2.0.)

Column temperature: 35

Flow rate: 1.5

Injection volume: 20
Detector: UV 300 or 215

CHROMATOGRAM
Retention time: 3.6

OTHER SUBSTANCES
Simultaneous: impurities, salicylic acid

REFERENCE
Kersten,B.S.; Catalano,T.; Rozenman,Y. Ion-pairing high-performance liquid chromatographic method for the determination of 5-aminosalicylic acid and related impurities in bulk chemical, *J.Chromatogr.*, **1991**, 588, 187–193.

SAMPLE
Matrix: formulations
Sample preparation: 75 μ L Sample + 6 mL mobile phase, vortex for 2 min, add 300 μ L 300 μ g/mL acetaminophen, make up to 10 mL with mobile phase, mix for 2 min, inject a 5 μ L aliquot.

HPLC VARIABLES
Column: 250 \times 4.6 5 μ m Spheri-5 ODS (Applied Biosystems)
Mobile phase: MeOH:buffer 30:70 (Buffer was 900 mL 50 mM Na_2HPO_4 + 18.75 mL tetrabutylammonium phosphate, pH adjusted to 6.8 with 1 N phosphoric acid.)
Flow rate: 1.2
Injection volume: 5
Detector: UV 254

CHROMATOGRAM
Retention time: 6.0
Internal standard: acetaminophen (3.8)

KEY WORDS
stability-indicating; rectal suspension; enema

REFERENCE
Henderson,L.M.; Johnson,C.E.; Berardi,R.R. Stability of mesalamine in rectal suspension diluted with distilled water, *Am.J.Hosp.Pharm.*, **1994**, 51, 2955–2957.

SAMPLE
Matrix: microsomal incubations
Sample preparation: 250 μ L Microsomal incubation + 500 μ L MeOH, centrifuge, inject an aliquot of the supernatant.

HPLC VARIABLES
Column: 120 \times 4.6 5 μ m RP-18
Mobile phase: MeCN:200 mM pH 7.5 phosphate buffer:1.25 mM cetyltrimethylammonium bromide 30:5:65
Flow rate: 1
Detector: F ex 315 em 500

KEY WORDS
rat; for 5-aminosalicylic acid-O-sulfate and 5-aminosalicylic acid

REFERENCE
Herzog,R.; Leuschner,J. Experimental studies on the pharmacokinetics and toxicity of 5-aminosalicylic acid-O-sulfate following local and systemic application, *Arzneimittelforschung*, **1995**, 45, 300–303.

SAMPLE
Matrix: tissue

Sample preparation: Freeze mucosal intestinal biopsy (ca. 5 mg) in liquid nitrogen and crush the sample, allow to warm to room temperature, add 20 μL propionic anhydride, add 100 μL 57.6 $\mu\text{g/mL}$ N-propionyl-4-aminosalicylic acid (purified by HPLC) in mobile phase, add 500 μL 50 mM pH 7.4 phosphate buffer, wash grinding/mixing rod with 500 μL 50 mM pH 7.4 phosphate buffer, sonicate (80 W) by immersing microprobe tip in mixture for 1 min, vortex, let stand at 37° for 1 h, add 500 μL 10% NaCl, add 6 mL MeCN, extract, cool at 4° for 1 h, centrifuge at 1500 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of air, reconstitute the residue in 500 μL mobile phase, filter (0.45 μm), inject a 100 μL aliquot.

HPLC VARIABLES

Guard column: 10 \times 4.6 5 μm Spherisorb ODS-2

Column: 150 \times 4.6 3 μm Spherisorb ODS-2

Mobile phase: MeCN:100 mM acetic acid:triethylamine 8:92:0.2

Flow rate: 1.5

Injection volume: 100

Detector: F ex 315 em 430

CHROMATOGRAM

Retention time: 5.7

Internal standard: N-propionyl-4-aminosalicylic acid (7.3)

Limit of detection: 1 ng

OTHER SUBSTANCES

Extracted: N-acetyl-5-aminosalicylic acid

KEY WORDS

mucosal intestinal biopsy; derivatization

REFERENCE

De Vos, M.; Verdier, H.; Schoonjans, R.; Beke, R.; De Weerd, G. A.; Barbier, F. High-performance liquid chromatographic assay for the determination of 5-aminosalicylic acid and acetyl-5-aminosalicylic acid concentrations in endoscopic intestinal biopsy in humans, *J. Chromatogr.*, **1991**, 564, 296–302.

SAMPLE

Matrix: tissue

Sample preparation: Sonicate tissue in 2 mL 5 ng/mL 3,4-dihydroxybenzylamine in MeOH for two 30 s cycles ($W = 60$), centrifuge at 1800 g for 10 min. Filter (0.5 μm) the supernatant, evaporate the filtrate to dryness under a stream of nitrogen under reduced pressure, reconstitute with 100 μL mobile phase, vortex, inject a 5 μL aliquot.

HPLC VARIABLES

Guard column: 50 \times 4.6 40 μm Pelliguard (Supelco)

Column: 250 \times 4.6 10 μm Erbasil S C18 (Carlo Erba)

Mobile phase: MeOH:buffer 15:85, pH adjusted to 3 with 100 mM NaOH. (Buffer was 10 mM Na_2HPO_4 containing 0.1 mM EDTA, 100 mM citric acid, and 0.1 mM heptanesulfonic acid.)

Flow rate: 1

Injection volume: 5

Detector: E, ESA Model 5100A, Model 5021 conditioning cell +0.35 V (between column and analytical cell), Model 5011 analytical cell, first electrode +0.05 V, second electrode -0.50 V (monitored)

CHROMATOGRAM

Retention time: 3.2

Internal standard: 3,4-dihydroxybenzylamine (4.4)

Limit of detection: 1 ng/mL


OTHER SUBSTANCES

Extracted: metabolites

REFERENCE

Palumbo, G.; Carlucci, G.; Mazzeo, P.; Frieri, G.; Pimpo, M.T.; Fanini, D. Simultaneous determination of 5-amino-salicylic acid, acetyl-5-aminosalicylic acid and 2,5-dihydroxybenzoic acid in endoscopic intestinal biopsy samples in humans by high-performance liquid chromatography with electrochemical detection, *J.Pharm.Biomed.Anal.*, **1996**, 14, 175-180.

Mescaline



Molecular formula: $C_{11}H_{17}NO_3$

Molecular weight: 211.26

CAS Registry No.: 54-04-6

Merck Index: 5965

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES


Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitrityline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenpropofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiacol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-stilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephénytoin, mephesin, mephobarbital, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methyprylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nyldrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phenacyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primi-

done, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfathiazole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranylcypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.

Mesna



Molecular formula: $C_2H_5NaO_3S_2$

Molecular weight: 164.18

CAS Registry No.: 19767-45-4, 3375-50-6 (free acid)

Merck Index: 5969

SAMPLE

Matrix: blood

Sample preparation: 50 μ L Plasma (within 3 min of collection) + penicillamine + 10 μ L 25 mM monobromobimane in MeCN, let stand for 5 min at room temperature, add 20 μ L 20% perchloric acid. (For total mesna add 100 μ L 20 mM dithiothreitol in 200 mM pH 8.5 Tris/HCl buffer to 50 μ L plasma, let stand for 40 min at room temperature, add 50 μ L 15% sulfosalicylic acid, centrifuge. Remove 200 μ L of the supernatant and wash it three times with ethyl acetate. Remove 60 μ L of the aqueous phase and add it to 300 μ L 200 mM pH 8.5 Tris/HCl buffer and 10 μ L 15 mM monobromobimane in MeCN, let stand at room temperature for 5 min, add 20 μ L 20% perchloric acid.)

HPLC VARIABLES

Column: 150 \times 4.6 7 μ m Nucleosil RP-18

Mobile phase: Gradient. A was MeCN. B was 1% aqueous acetic acid containing 1 g/L octane-sulfonic acid. A:B from 5:95 to 8:92 over 2 min, to 10:90 over 13 min (Waters convex), to 30:70 over 20 min (Waters convex), maintain at 30:70 for 6 min, re-equilibrate at initial conditions for 9 min.

Flow rate: 1.4

Detector: F (wavelengths not given)

CHROMATOGRAM

Internal standard: penicillamine

Limit of detection: 10 μ M

KEY WORDS

plasma; derivatization; pharmacokinetics

REFERENCE

Stofer-Vogel,B.; Cerny,T.; Borner,M.; Lauterburg,B.H. Oral bioavailability of mesna tablets, *Cancer Chemother.Pharmacol.*, **1993**, *32*, 78-81.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 400 μ L Plasma + 400 μ L 333 mM sulfuric acid + 400 μ L 5% sodium hexametaphosphate, vortex for 5 s, centrifuge at 15000 g for 2 min, inject a 20 μ L aliquot of the supernatant. (Process plasma as rapidly as possible after collection.) (For total mesna mix 100 μ L plasma, 100 μ L EDTA solution, and 100 μ L 4% sodium borohydride, heat to 50° for 30 min, cool, add 200 μ L 9.7% acetic acid, store in the dark for 3-4 days, inject an aliquot. (EDTA solution was 1% disodium EDTA in a mixture of 50 mL 500 mM Na_2HPO_4 and 22 mL 1 M NaOH.)) Urine. Dilute 1:10 to 1:400 with mobile phase, inject an aliquot. (For total mesna dilute 1:10 to 1:1000 then proceed as above.)

HPLC VARIABLES

Column: 100 \times 4.6 3 μ m Hypersil ODS

Mobile phase: MeCN:buffer 3:97 (Buffer was 100 mM pH 7.0 phosphate buffer containing 5 mM n-tetrabutylammonium hydroxide.)

Flow rate: 1

Injection volume: 20

Detector: E, Coulochem Model 5100A, Model 5011 analytical cell, detector 1 not used, detector 2 -200 mV

CHROMATOGRAM

Retention time: 4.5

Limit of quantitation: 6.1 μ M

KEY WORDS

plasma

REFERENCE

James,C.A.; Rogers,H.J. Estimation of mesna and dimesna in plasma and urine by high-performance liquid chromatography with electrochemical detection, *J.Chromatogr.*, **1986**, 382, 394-398.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 100 μ L Plasma + 25 μ L 330 mM sulfuric acid + 25 μ L 5% sodium hexametaphosphate + 10 μ L 40 μ g/mL p-aminobenzoic acid, make up to 200 μ L, vortex for 2 s, centrifuge at 3400 g for 6 min, inject a 50-100 μ L aliquot of the supernatant. Urine. Dilute 1:50 with water. 100 μ L Diluted urine + 100 μ L 150 μ g/mL p-aminobenzoic acid in 1.25% sodium hexametaphosphate, vortex for 2 s, inject a 50 μ L aliquot. (Reduce dimesna to mesna as follows. 100 μ L Plasma or diluted urine + 100 μ L 1% EDTA in buffer + 100 μ L 1.06 M sodium borohydride in water, mix thoroughly, heat at 50° for 30 min, cool to room temperature, add 200 μ L 1.74 M acetic acid, vortex for 20 s, centrifuge for 10 min, inject a 50-100 μ L aliquot of the supernatant. Buffer was 50 mL 500 mM Na_2HPO_4 + 22 mL NaOH.)

HPLC VARIABLES

Guard column: Guard PAK C18 (Waters)

Column: 100 \times 8 10 μ m Resolve C18 (Waters)

Mobile phase: 100 mM Sodium citrate containing 1 mM tetrabutylammonium phosphate and 0.71 mM triethylamine, pH adjusted to 5 with 85% phosphoric acid

Flow rate: 2

Injection volume: 50-100

Detector: E, ESA Coulochem II model 5100, model 5011 analytical cell +450 mV, model 5021 conditioning cell +500 mV

CHROMATOGRAM

Retention time: 6.73

Internal standard: p-aminobenzoic acid (8.87)

Limit of quantitation: 1 μ g/mL

OTHER SUBSTANCES

Simultaneous: acetaminophen, 5-fluorouracil, procarbazine, prochlorperazine

Noninterfering: aspirin, bleomycin, carboplatin, carmustine, chlorambucil, cyclophosphamide, cyclosporine A, cytarabine, doxorubicin, etoposide, hydrocortisone, ifosfamide, lomustine, methotrexate, mitomycin, mitoxantrone, teniposide, thiotepa, vincristine

KEY WORDS

plasma; pharmacokinetics

REFERENCE

el-Yazigi,A.; Yusuf,A.; Al-Rawithi,S. Liquid chromatographic analysis of mesna and dimesna in plasma and urine of patients treated with mesna, *Ther.Drug Monit.*, **1995**, 17, 153-158.

SAMPLE

Matrix: formulations

Sample preparation: Dilute formulation with 15 µg/mL disodium EDTA solution to a mesna concentration of 3.2 µg/mL. Remove a 1 mL aliquot and add it to 300 µL reagent solution, let stand at room temperature for 20 min, add 500 µL 300 mM phosphoric acid solution, add 3 mL 4 µg/mL IS solution, make up to 10 mL with water, inject a 50 µL aliquot. (Prepare the reagent solution by dissolving 3.5 mg methyl 4-(6-methoxynaphthalen-2-yl)-4-oxo-2-butenate in 10 mL THF, make up to 25 mL with pH 7.5 borate buffer. Prepare methyl 4-(6-methoxynaphthalen-2-yl)-4-oxo-2-butenate as follows. Dissolve 5 g 6'-methoxy-2'-acetone naphthone in warm glacial acetic acid and add 2.5 g glyoxylic acid, reflux for 24 h, evaporate to dryness under reduced pressure. Take up the residue in chloroform and extract it three times with 5% sodium carbonate solution. Combine the aqueous layers and acidify them with concentrated HCl, collect the product by filtration, recrystallize from MeOH/water or acetic acid to give 4-(6-methoxynaphthalen-2-yl)-4-oxo-2-butenic acid (mp 167-9°) (Farmaco, Ed. Sci. 1982, 37, 171). Reflux 0.5 g 4-(6-methoxynaphthalen-2-yl)-4-oxo-2-butenic acid, 2.5 mL MeOH, and 2-3 drops sulfuric acid in 25 mL anhydrous benzene (Caution! Benzene is a carcinogen!) for 1 h, add 20 mL water, wash the organic layer with 10 mL 5% sodium bicarbonate solution, wash the organic layer with 20 mL water. Dry the organic layer over anhydrous sodium sulfate, evaporate to dryness under reduced pressure, purify by flash chromatography on silica gel using ethyl acetate:light petroleum (bp 40-70°) 40:60 to give methyl 4-(6-methoxynaphthalen-2-yl)-4-oxo-2-butenate as a pale yellow compound (mp 116-120°).)

HPLC VARIABLES

Column: 150 × 4.5 µm Spherisorb RP-8

Mobile phase: MeOH:50 mM pH 3.0 triethylammonium phosphate 53:47

Flow rate: 1

Injection volume: 50

Detector: F ex 310 em 450

CHROMATOGRAM

Retention time: 6

Internal standard: 4-(6-methoxynaphthalen-2-yl)-4-oxobutanoic acid (8)

OTHER SUBSTANCES

Simultaneous: acetylcysteine, cysteamine, cysteine, glutathione, homocysteine

Noninterfering: bacitracin, biotin, calcium pantothenate, cystine, glycine, magnesium oxide, neomycin, starch, threonine, vitamin E, pyridoxine, riboflavin phosphate

KEY WORDS

solutions; derivatization

REFERENCE

Gatti,R.; Cavrini,V.; Roveri,P.; Pinzauti,S. High-performance liquid chromatographic determination of aliphatic thiols with acryloyl acids as fluorogenic precolumn derivatization reagents, *J.Chromatogr.*, **1990**, 507, 451-458.

SAMPLE

Matrix: formulations

Sample preparation: Dilute with mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.6 10 µm cyano

Mobile phase: MeCN:100 mM NaH₂PO₄ 20:80 adjusted to pH 4.2 with phosphoric acid

Flow rate: 2.5

Injection volume: 20

Detector: UV 200

CHROMATOGRAM

Retention time: 1.16

OTHER SUBSTANCES

Simultaneous: diphenhydramine, granisetron (UV 300)

KEY WORDS

stability-indicating; injections; saline

REFERENCE

Mayron,D.; Gennaro,A.R. Stability and compatibility of granisetron hydrochloride in i.v. solutions and oral liquids and during simulated Y-site injection with selected drugs, *Am.J.Health-Syst.Pharm.*, **1996**, *53*, 294–304.

SAMPLE

Matrix: urine

Sample preparation: Dilute urine 1:3 to 1:39, inject a 50 µL aliquot.

HPLC VARIABLES

Column: 250 × 4 5 µm Hypersil ODS

Mobile phase: MeOH:250 mM pH 7.4 phosphate buffer 5:95 containing 5 mM tetrabutylammonium phosphate

Flow rate: 1

Injection volume: 50

Detector: UV 412 following post-column reaction. The column effluent mixed with the reagent pumped at 0.5 mL/min and the mixture flowed through a 200 × 4 column packed with 100–120 mesh glass beads (dichlorodimethylsilane treated) to the detector. (Prepare reagent by diluting 0.2% 5,5'-dithiobis(2-nitrobenzoic acid) in 250 mM pH 7.4 phosphate buffer containing 10% tripotassium citrate 1:10 with water.)

CHROMATOGRAM

Retention time: 7.5

Limit of detection: 75 ng

KEY WORDS

post-column reaction; comparison with electrochemical detection without derivatization

REFERENCE

Sidau,B.; Shaw,I.C. Determination of sodium 2-mercaptoethanesulphonate by high-performance liquid chromatography using post-column reaction colorimetry or electrochemical detection, *J.Chromatogr.*, **1984**, *311*, 234–238.

Mesoridazine


Molecular formula: C₂₁H₂₆N₂OS₂

Molecular weight: 386.58

CAS Registry No.: 5588-33-0, 32672-69-8 (besylate)

Merck Index: 5970

Lednicer No.: 1 389



SAMPLE

Matrix: blood

Sample preparation: 2 mL Plasma + 100 µL 1 M HCl, vortex for 30 s, add 4 mL isopropanol, mix for 5 min, centrifuge at 5000 rpm at 0° for 20 min. Remove the supernatant and adjust