

KEY WORDS

chiral

REFERENCE

Kristensen,K.; Angelo,H.R. A stereoselective HPLC method for the determination of methadone and its main metabolite in urine (using an AGP column) (Abstract 39), *Ther.Drug Monit.*, **1995**, 17, 393–393.

Methamphetamine

**Molecular formula:** C₁₀H₁₅N**Molecular weight:** 149.24**CAS Registry No.:** 537-46-2**Merck Index:** 6015**Lednicer No.:** 1 37

SAMPLE**Matrix:** blood**Sample preparation:** Inject a 5 µL aliquot of serum directly.

HPLC VARIABLES**Column:** 100 × 4.6 5-10 µm Silicalite (by sieving Silicalite, 3M Co.(?))**Mobile phase:** MeNC:20 mM pH 6.9 phosphate buffer 10:90**Flow rate:** 1**Injection volume:** 5**Detector:** UV 254

CHROMATOGRAM**Retention time:** 3.79

OTHER SUBSTANCES**Extracted:** ethosuximide, sulfamethoxazole, primidone

KEY WORDS

serum

REFERENCE

Ambrose,D.L.; Fntz,J.S. High-performance liquid chromatographic determination of drugs and metabolites in human serum and urine using direct injection and a unique molecular sieve, *J.Chromatogr.B*, **1998**, 709, 89–96.

SAMPLE**Matrix:** blood

Sample preparation: 1 mL Plasma + 1 mL 100 mM NaOH + 3 mL n-hexane, shake for 20 min, centrifuge for 10 min. Remove 2 mL of the organic layer and evaporate it to dryness using a vacuum centrifuge, reconstitute the residue in 500 µL 100 µg/mL (S)-(+)-benoxaprofen chloride in dried dichloromethane, let stand at room temperature for 30 min, inject a 10 µL aliquot. (Synthesis of benoxaprofen chloride is as follows. Dissolve 600 mg benoxaprofen in 50 mL toluene, slowly add 5 mL freshly-distilled thionyl chloride, reflux for 30 min, evaporate to dryness, recrystallize benoxaprofen chloride from dichloromethane.)

HPLC VARIABLES**Column:** 250 × 4.6 7 µm Zorbax-Sil**Mobile phase:** Cyclohexane:dichloromethane:THF 50:10:10**Flow rate:** 1**Injection volume:** 10**Detector:** F ex 312 em 365

CHROMATOGRAM

Retention time: 10.5 (R-(-)), 11.5 (S-(+))

OTHER SUBSTANCES

Extracted: amphetamine

Interfering: tranlycypromine

KEY WORDS

plasma; derivatization; normal phase; chiral

REFERENCE

Weber,H.; Spahn,H.; Mutschler,E.; Möhrke,W. Activated α -alkyl- α -arylacetic acid enantiomers for stereoselective thin-layer chromatographic and high-performance liquid chromatographic determination of chiral amines, *J.Chromatogr.*, **1984**, 307, 145–153.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Serum + 1 mL acetone, vortex for 2 min, centrifuge at 2000 g for 10 min. Remove a 500 μ L aliquot of the supernatant and add it to 500 μ L reagent, add 20 μ L 0.5% triethylamine in MeOH, vortex, heat at 80° for 30 min, cool to room temperature, inject a 20 μ L aliquot. (Prepare reagent by mixing 2 mL 2.5 mM N-(4-aminobutyl)-N-ethylisoluminol in MeOH with 2 mL 2.5 mM N,N'-disuccinimidyl carbonate in MeCN, let stand for 2 h.)

HPLC VARIABLES

Guard column: 30 \times 4.6 TSKm Guardgel ODS-80TM (Toyo Soda)

Column: 150 \times 6 5 μ m Shimpack CLC C18 (Shimadzu)

Mobile phase: MeOH:water 54:46 containing 30 mM sodium 1-octanesulfonate

Flow rate: 1

Injection volume: 20

Detector: Chemiluminescence following post-column reaction. The column effluent mixed with 15 mM potassium ferricyanide in 2.5 M NaOH pumped at 1 mL/min and this mixture flowed through a 200 mm \times 0.5 mm ID stainless steel coil. The effluent from this coil mixed with 300 mM hydrogen peroxide containing 10 mM β -cyclodextrin pumped at 1 mL/min and this mixture flowed through a 100 mm \times 0.5 mm ID stainless steel coil to the detector.

CHROMATOGRAM

Retention time: 30

Limit of detection: 20 pM

KEY WORDS

derivatization; serum

REFERENCE

Nakashima,K.; Suetsugu,K.; Akiyama,S.; Yoshida,M. High-performance liquid chromatography-chemiluminescence determination of methamphetamine in human serum using N-(4-aminobutyl)-N-ethylisoluminol as a chemiluminogen, *J.Chromatogr.*, **1990**, 530, 154–159.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Serum + 50 μ L 100 ng/mL aniline sulfate in water + 200 μ L 20 mM pH 10.6 carbonate buffer + 2 mL ethyl acetate, shake for 15 min, centrifuge at 1200 g for 5 min. Remove the organic layer and add it to 200 μ L 50 mM HCl, shake for 15 min, centrifuge at 1200 g for 5 min. Remove the aqueous layer and add it to 40 μ L 250 mM NaOH, add 50 μ L 330 mM pH 7.8 phosphate buffer, add 250 μ L MeCN, add 25 μ L 1 mM (-)-1-(9-fluorenyl)ethyl chloroformate in acetone, let stand overnight at room temperature, add 30 μ L 100 mM glycine in water, add 750 μ L n-pentane, vortex for 2 min, centrifuge at 1200 g for 5 min. Remove the organic layer and evaporate it to dryness under vacuum, reconstitute the residue in MeCN: water 50:50, inject a 100 μ L aliquot.

HPLC VARIABLES

Guard column: Direct-Connect column prefilter (Alltech)

Column: 150 × 4.6 3 µm Adsorbosphere HS C18 (Alltech)
Mobile phase: MeCN:THF:20 mM pH 3.6 acetate buffer 39:15:46
Flow rate: 1
Injection volume: 100
Detector: F ex 265 em 330

CHROMATOGRAM

Retention time: 27.7 (R), 29.0 (S)
Internal standard: aniline (21.0)
Limit of quantitation: 5 ng/mL

OTHER SUBSTANCES

Extracted: amphetamine

KEY WORDS

serum; rat; chiral; derivatization

REFERENCE

Hutchaleelaha,A.; Walters,A.; Chow,H.-H.; Mayersohn,M. Sensitive enantiomer-specific high-performance liquid chromatographic analysis of methamphetamine and amphetamine from serum using precolumn fluorescent derivatization, *J.Chromatogr.B*, **1994**, 658, 103–112.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 1 mL 500 mM pH 11 borate buffer, mix, add 2.5 mL diethyl ether, vortex for 5 min, centrifuge at 1200 g for 5 min, remove organic layer, repeat extraction. Combine the organic layers and add them to 200 µL 100 mM HCl, vortex for 2 min, centrifuge at 1200 g for 5 min. Remove the aqueous phase and add it to 150 µL 1 M pH 8 borate buffer and 100 µL 4 mM 9-fluorenylmethyl chloroformate in MeCN, shake, allow to react at 50° for 5 min, add 20 µL 20 mM proline in water, allow to react at 50° for 2 min, inject a 200 µL aliquot.

HPLC VARIABLES

Column: 150 × 3.9 4 µm Nova-Pak phenyl
Mobile phase: MeCN:50 mM pH 6.0 sodium phosphate buffer 50:50
Flow rate: 1
Injection volume: 200
Detector: F ex 260 em 315

CHROMATOGRAM

Retention time: 15
Limit of quantitation: 0.5 ng/mL

OTHER SUBSTANCES

Extracted: amphetamine, desmethyldeprenyl

KEY WORDS

plasma

REFERENCE

La Croix,R.; Pianezzola,E.; Strolin Benedetti,M. Sensitive high-performance liquid chromatographic method for the determination of the three main metabolites of selegiline (L-deprenyl) in human plasma, *J.Chromatogr.B*, **1994**, 656, 251–258.

SAMPLE

Matrix: blood, urine

Sample preparation: Adjust pH of 10 mL plasma or 20 mL urine to 11.4 with 5 M NaOH, add to a column containing 1.5 g Amberlite XAD-2, wash with 10 mL water, elute with 20 (plasma) or 40 (urine) mL chloroform:isopropanol 75:25, add 100 µL 6 M HCl in EtOH to the eluate, evaporate to dryness under reduced pressure, reconstitute with 1 mL 8% sodium bicarbonate,

add 1 mL 0.5% sodium 1,2-naphthoquinone-4-sulfonate, heat at 70° for 20 min, add an equal volume of chloroform, vortex for 1 min, inject a 50 µL aliquot of the organic layer.

HPLC VARIABLES

Column: 150 × 5 Partisil 5

Mobile phase: Hexane:chloroform:ethyl acetate:EtOH 50:25:35:1

Column temperature: 20

Flow rate: 2.5

Injection volume: 50

Detector: UV 248

CHROMATOGRAM

Retention time: 3

Internal standard: phenylethylamine (6)

Limit of detection: 2 ng

OTHER SUBSTANCES

Extracted: amphetamine, hydroxyamphetamine, norephedrine

KEY WORDS

derivatization; plasma; normal phase; SPE; comparison with other derivatization reagents and with ion-pair chromatography

REFERENCE

Farrell,B.M.; Jefferies,T.M. An investigation of high-performance liquid chromatographic methods for the analysis of amphetamines, *J.Chromatogr.*, **1983**, 272, 111–128.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 206.4

CHROMATOGRAM

Retention time: 8.433

KEY WORDS

whole blood

REFERENCE

Gaillard,X.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

SAMPLE**Matrix:** bulk**Sample preparation:** Mix a 1 mg/mL solution in 1 M sodium carbonate with 2 mL 5 mg/mL 8-quinolinesulfonyl chloride in acetone, heat at 65° for 20 min, cool, extract twice with 30 mL portions of chloroform. Combine the extracts and dry them over anhydrous magnesium sulfate, evaporate to dryness under a stream of air, reconstitute, inject an aliquot.

HPLC VARIABLES**Guard column:** 70 × 2.1 Co:Pell ODS**Column:** 300 × 3.9 µBondapak C18**Mobile phase:** MeCN:water:acetic acid 40:59:1**Flow rate:** 1.5**Detector:** UV 254, UV 280

CHROMATOGRAM**Retention time:** 25

OTHER SUBSTANCES**Simultaneous:** amphetamine, ephedrine, phenmetrazine, phentermine, phenylpropanolamine, pseudoephedrine

KEY WORDSderivatization

REFERENCENoggle,F.T.,Jr.; Clark,C.R. Liquid chromatographic determination of primary and secondary amines as 8-quinolinesulfonyl chloride derivatives, *J.Assoc.Off.Anal.Chem.*, **1984**, 67, 687–691.

SAMPLE**Matrix:** bulk**Sample preparation:** Prepare an aqueous solution. Adjust pH of 500 µL aqueous solution to 12 with 1 M NaOH, add 200 µL (+)-1-(9-fluorenyl)ethyl chloroformate:dichloromethane 1:100, let stand for 20 min, add 500 µL dichloromethane, shake for 30 min. Remove the organic phase and wash it twice with 500 µL aliquots of water, evaporate the organic phase to dryness under a stream of nitrogen, reconstitute the residue in 3 mL MeOH, filter (0.45 µm), inject a 25 µL aliquot.

HPLC VARIABLES**Column:** 250 × 3.9 5 µm 5C18-AR (Waters)**Mobile phase:** MeCN:50 mM pH 6.0 phosphate buffer 65:35**Flow rate:** 1**Injection volume:** 25**Detector:** F ex 295 em 315

CHROMATOGRAM**Retention time:** 44 ((S)-(+)), 47 ((R)-(-))**Limit of quantitation:** 16.7 ng/mL

OTHER SUBSTANCES**Extracted:** ephedrine, pseudoephedrine

KEY WORDSderivatization; chiral

REFERENCEChen,Y.-P.; Hsu,M.-C.; Chien,C.S. Analysis of forensic samples using precolumn derivatization with (+)-1-(9-fluorenyl)ethyl chloroformate and liquid chromatography with fluorimetric detection, *J.Chromatogr.A*, **1994**, 672, 135–140.

SAMPLE**Matrix:** bulk

Sample preparation: Dissolve 10 μmole compound (as free base or hydrochloride) in 500 μL MeCN, add 250 μL 5% sodium carbonate (for hydrochlorides only), add 500 μL 100 mM reagent in MeCN, vortex for 1 min, heat at 60° for 2 h, add 100 μmole L-proline, heat at 60° for 30 min. Remove a 100 μL aliquot and dilute it with mobile phase, neutralize with acetic acid, inject a 10 μL aliquot. Prepare the reagent ((R,R)-N-(3,5-dinitrobenzoyl)-2-aminocyclohexylisothiocyanate) as follows. Add 0.7 mL carbon disulfide to 6 mL (1R,2R)-(-)-1,2-diaminocyclohexane, 12 mL water, and 12 mL EtOH, heat the oil bath to 80°, add 2.8 mL carbon disulfide dropwise (making sure that the product does not start to precipitate), when addition is complete reflux for 1 h, acidify with 500 μL 5 M HCl, reflux for 12 h, cool, filter, wash the solid with a little cold EtOH to give trans-4,5-tetramethyleneimidazolidine-2-thione as a white fluffy solid (mp 148-150°) (Tetrahedron 1993, 49, 4419). Stir 7.97 g 3,5-dinitrobenzoyl chloride in 30 mL dichloroethane at 50°, add a solution of 6 g trans-4,5-tetramethyleneimidazolidine-2-thione in 120 mL dichloroethane containing a catalytic amount of 4-(dimethylamino)pyridine over 15 min, reflux for 2 h, remove the crystals of (R,R)-N-(3,5-dinitrobenzoyl)-2-aminocyclohexylisothiocyanate by filtration, evaporate the filtrate to dryness and dissolve the residue in 60 mL dichloroethane, reflux for 16 h to obtain more (R,R)-N-(3,5-dinitrobenzoyl)-2-aminocyclohexylisothiocyanate (mp >250°, $[\alpha]_{\text{D}}^{25} = -133^\circ$ (c = 1) in MeCN).

HPLC VARIABLES

Column: 125 \times 4.5 μm Lichrospher 60 RP Select B
Mobile phase: MeCN:20 mM ammonium acetate 55:45
Flow rate: 1
Injection volume: 10
Detector: UV 254

CHROMATOGRAM

Retention time: k' 14.33, k' 15.90 (enantiomers)

OTHER SUBSTANCES

Also analyzed: acebutolol, alprenolol, atenolol, carazolol, carvedilol, formoterol, metipranolol, metoprolol, nifedipine, nifedipine, nifedipine, nifedipine, nifedipine, nifedipine, nifedipine, nifedipine, nifedipine, nifedipine

KEY WORDS

derivatization; chiral

REFERENCE

Kleidernigg, O.P.; Posch, K.; Lindner, W. Synthesis and application of a new isothiocyanate as a chiral derivatizing agent for the indirect resolution of chiral amino alcohols and amines, *J. Chromatogr. A*, **1996**, 729, 33-42.

SAMPLE

Matrix: formulations

Sample preparation: Powder tablet and add 50 mg to 50 mL MeCN:20 mM pH 3.8 phosphate buffer 3:97, sonicate for 5 min, filter (0.5 μm), inject a 20 μL aliquot of the filtrate.

HPLC VARIABLES

Guard column: Supelguard pre-column containing 5 μm Suplex pKb100 (Supelco)
Column: 150 \times 4.6 5 μm Suplex pKb100 (Supelco)
Mobile phase: Gradient. MeCN:20 mM pH 3.8 phosphate buffer at 3:97 for 3 min, to 15:85 over 5 min, stay at 15:85 for 4 min, re-equilibrate for 8 min.
Flow rate: 1.5
Injection volume: 20
Detector: UV 220 for 5 min then UV 280

CHROMATOGRAM

Retention time: 4.3
Limit of quantitation: 10 $\mu\text{g/mL}$

OTHER SUBSTANCES

Simultaneous: ephedrine, amphetamine, caffeine, 3,4-methylenedioxyamphetamine, N-methyl-3,4-methylenedioxyamphetamine, N-ethyl-3,4-methylenedioxyamphetamine

KEY WORDS

tablets

REFERENCE

Longo,M.; Martinez,C.; Rolandi,L.; Cavallaro,A. Simple and fast determination of some phenethylamines in illicit tablets by base-activated reversed phase HPLC, *J.Liq.Chromatogr.*, **1994**, *17*, 649-658.

SAMPLE

Matrix: solutions

Sample preparation: Mix 1 mL 20-300 µg/mL amine solution in water with 2 mL 50 mg/mL 4-nitrobenzoyl chloride in THF (freshly prepared) and 1 mL 1 M NaOH, heat at 65° for 1 h, cool, adjust pH to 12 with 1 M NaOH, extract with two 10 mL portions of chloroform. Combine the extracts and wash them with two 20 mL portions of 10% potassium carbonate, wash with water, dry over anhydrous magnesium sulfate. Evaporate to dryness under a stream of air, reconstitute the residue in MeOH, inject a 5 µL aliquot.

HPLC VARIABLES

Column: 300 × 3.9 µBondapak C18

Mobile phase: MeCN:water 35:65

Flow rate: 1.5

Injection volume: 5

Detector: UV 254

CHROMATOGRAM

Retention time: 15

OTHER SUBSTANCES

Simultaneous: amphetamine, benzylamine, α-methylbenzylamine, n-propylamphetamine

KEY WORDS

derivatization

REFERENCE

Clark,R.C.; Teague,J.D.; Wells,M.M.; Ellis,J.H. Gas and high-pressure liquid chromatographic properties of some 4-nitrobenzamides of amphetamines and related arylalkylamines, *Anal.Chem.*, **1977**, *49*, 912-915.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in MeOH at a concentration of 1 mg/mL, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 250 × 5 Spherisorb S5W

Mobile phase: MeOH:buffer 90:10 (Buffer was 94 mL 35% ammonia and 21.5 mL 70% nitric acid in 884 mL water, adjust the pH to 10.1 with ammonia.)

Flow rate: 2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 4.05

OTHER SUBSTANCES

Simultaneous: dihydrocodeine, dihydromorphine, levorphanol, norcodeine, normorphine, pemo-line, benzphetamine, diethylpropion, mazindol, tranlycypromine, caffeine, fenethyline, phen-dimetrazine, methylphenidate, phenelzine, epinephrine, pipradol, phenylpropanolamine, fencamfamin, chlorphentermine, norpseudoephedrine, fenfluramine, methylenedioxyamphet-amine, amphetamine, normetanephine, 4-hydroxyamphetamine, bromo-STP, STP, prolintane, 2-phenethylamine, tyramine, trimethoxyamphetamine, phenylephrine, pseudoephedrine, ephedrine, methylephedrine, mephentermine, buprenorphine, dextromoramide, phenoperidine, fentanyl, etorphine, piritramide, noscapine, papaverine, naloxone, dextropropoxyphene, nalor-phine, phenazocine, norpipanone, levallorphan, hydroxypethidine, normethadone, meperidine, dipipanone, diamorphine, pentazocine, acetylcodeine, monoacetylmorphine, thebacon, oxy-codone, thebaine, norlevorphanol, methadone, benzylmorphine, ethylmorphine, morphine-N-oxide, codeine, codeine-N-oxide, morphine, ethoheptazine, morphine-3-glucuronide, pholcodine

Noninterfering: dopamine, levodopa, methyl dopa, methyl dopate, norepinephrine

Interfering: dimethylamphetamine, mescaline, norpethidine, hydrocodone

REFERENCE

Law, B.; Gill, R.; Moffat, A.C. High-performance liquid chromatography retention data for 84 basic drugs of forensic interest on a silica column using an aqueous methanol eluent, *J. Chromatogr.*, **1984**, *301*, 165–172.

SAMPLE

Matrix: solutions

Sample preparation: Mix 1 mL of an aqueous solution with 1 mL 100 mM nickel sulfate in water, 1 mL 20% aqueous ammonia, and 5 mL chloroform:carbon disulfide 98:2, shake vigorously for 1 min, wash the organic layer with three 2 mL portions of water, filter (phase-separation paper). Evaporate the filtrate to dryness under a stream of nitrogen, reconstitute with 1 mL mobile phase, inject a 10 μ L aliquot. (Copper may also be used with electrochemical detection or UV detection at 270 nm.)

HPLC VARIABLES

Guard column: 30 \times 4 40 μ m LiChrosorb RP-18

Column: 250 \times 4 7 μ m LiChrosorb RP-18

Mobile phase: MeOH:20 mM pH 5.8 sodium acetate buffer 80:20 containing 5 mM lithium perchlorate

Flow rate: 1.5

Injection volume: 10

Detector: UV 325, E, Merck-Clevenot E 230, Model LCC 231 thin-layer electrolytic cell with a glassy carbon electrode at +0.7 V, standard calomel reference electrode

CHROMATOGRAM

Retention time: 8.49

Limit of detection: 1 fmole (E), 1 nmole (UV)

OTHER SUBSTANCES

Simultaneous: ephedrine

Also analyzed: acebutolol, alprenolol, flecainide, propranolol

KEY WORDS

derivatization; complexation

REFERENCE

Leroy, P.; Nicolas, A. Determination of secondary amino drugs as their metal dithiocarbamate complexes by reversed-phase high-performance liquid chromatography with electrochemical detection, *J. Chromatogr.*, **1984**, *317*, 513–521.

SAMPLE

Matrix: solutions

Sample preparation: 2 mL THF + 1 mL 33.5 mM reagent in THF (freshly prepared) + 1 mL 1 mg/mL amphetamine in water + 700 μ L 10% sodium bicarbonate in water, heat at 65° for 1 h, cool, extract three times with 10 mL aliquots of chloroform. Combine the extracts and wash them with 10 mL water, dry over anhydrous magnesium sulfate, evaporate to dryness, reconstitute with 2.5 mL mobile phase, inject a 5 μ L aliquot. (Prepare reagent 1-[(4-nitrophenyl)sulfonyl]propyl chloride) as follows. Mix 40–45 mmoles L-(-)-proline, 40 mL THF, and 200 mL 10% potassium carbonate, add 37–43 mmoles 4-nitrobenzenesulfonyl chloride in 40 mL THF dropwise, heat at 50° and maintain at pH 8 or above for 3 h, cool, acidify to pH 2, extract with chloroform. Extract the organic layers with potassium carbonate in water. Acidify the aqueous layer and extract it with chloroform. Dry the chloroform layer and evaporate it to dryness, recrystallize the resulting 1-[(4-nitrophenyl)sulfonyl]proline from petroleum ether and benzene (Caution! Benzene is a carcinogen!). Stir 15 mmoles 1-[(4-nitrophenyl)sulfonyl]proline in 100 mL benzene and add 75 mmoles thionyl chloride in 50 mL benzene dropwise, heat at 35–40° until the reaction is complete (about 48 h; monitor by IR), evaporate to dryness, recrystallize from n-heptane to give 1-[(4-nitrophenyl)sulfonyl]propyl chloride (Anal. Chem. **1984**, *56*, 958) (mp 110–110.5°).)

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Zorbax ODS

Mobile phase: MeOH:water 60:40**Flow rate:** 1.5**Injection volume:** 5**Detector:** UV 254**CHROMATOGRAM****Retention time:** 12 (S), 14 (R)**KEY WORDS**

derivatization; chiral

REFERENCE

Barksdale, J.M.; Clark, C.R. Liquid chromatographic determination of the enantiomeric composition of amphetamine and related drugs by diastereomeric derivatization, *J.Chromatogr.Sci.*, **1985**, *23*, 176–180.

SAMPLE**Matrix:** solutions**Sample preparation:** Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.**HPLC VARIABLES****Column:** 125 × 4.9 Spherisorb S5W silica**Mobile phase:** MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7**Flow rate:** 2**Injection volume:** 20**Detector:** E, LeCarbone, V25 glassy carbon electrode, + 1.2 V**CHROMATOGRAM****Retention time:** 2.7**OTHER SUBSTANCES**

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzoctamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazepine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipranone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methapyrilene, methdilazene, methotrimprazine, methoxamine, methoxyphenamine, methoxypropazine, methylephedrine, methyl-ergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampro-mide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, pheninda-mine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxy-benzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozide, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenze-pine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, prima-quine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolin-tane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, qui-nine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldi-

amine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleppamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J.Chromatogr.*, **1985**, 323, 191–225.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in MeOH:water 1:1 at a concentration of 50 µg/mL, inject a 10 µL aliquot.

HPLC VARIABLES

Column: 300 × 3.9 10 µm µBondapak C18

Mobile phase: MeOH:acetic acid:triethylamine:water 15:1.5:0.5:83

Flow rate: 1.5

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 13

OTHER SUBSTANCES

Simultaneous: phenylpropanolamine, ephedrine, hydroxyamphetamine, amphetamine, phentermine, mephentermine

REFERENCE

Roos, R.W.; Lau-Cam, C.A. General reversed-phase high-performance liquid chromatographic method for the separation of drugs using triethylamine as a competing base, *J.Chromatogr.*, **1986**, 370, 403–418.

SAMPLE

Matrix: solutions

Sample preparation: 100 µL 100 µM Solution + 300 µL buffer + 500 µL 1 mM dansyl chloride in acetone, mix, heat at 45° in the dark for 1 h, dilute with MeCN:water 50:50, inject a 20 µL aliquot. (Prepare buffer by adjusting 10 mM sodium bicarbonate to pH 9.0 with NaOH.)

HPLC VARIABLES

Guard column: 30 × 4.6 Spheri-5 RP-18

Column: 250 × 4.6 Inertsil ODS-2

Mobile phase: MeCN:water 70:30 containing 1 mM imidazole, pH adjusted to 7.0 with nitric acid

Flow rate: 1

Injection volume: 20

Detector: Chemiluminescence following post-column reaction. The column effluent mixed with the reagent pumped at 1 mL/min and the mixture flowed through a 300 mm × 0.25 mm ID coil to the detector. (Prepare the reagent by dissolving 112 mg bis(2,4,6-trichlorophenyl) oxalate in 500 mL MeCN, add 8.6 mL 30% hydrogen peroxide, sonicate.), F ex 343 em 530

CHROMATOGRAM

Retention time: 18.5

Limit of detection: 4 fmole (chemiluminescence), 50 fmole (F)

OTHER SUBSTANCES

Simultaneous: benzylamine, ephedrine, N-isopropylbenzylamine, N-methylphenethylamine, phenylbutylamine, phenylethylamine, phenylpropanolamine, phenylpropylamine

KEY WORDS

derivatization; post-column reaction; comparison with other derivatization reagents

REFERENCE

Hayakawa,K.; Hasegawa,K.; Imaizumi,N.; Wong,O.S.; Miyazaki,M. Determination of amphetamine-related compounds by high-performance liquid chromatography with chemiluminescence and fluorescence detections, *J.Chromatogr.*, **1989**, *464*, 343–352.

SAMPLE

Matrix: solutions

Sample preparation: 100 μ L 100 μ M Solution + 400 μ L buffer + 500 μ L 80 mM 4-fluoro-7-nitrobenzoxadiazole in EtOH, mix, heat at 60° in the dark for 1 min, dilute with MeCN:water 50:50, inject a 20 μ L aliquot. (Prepare buffer by adjusting 100 mM boric acid to pH 8.0 with NaOH.)

HPLC VARIABLES

Guard column: 30 \times 4.6 Spheri-5 RP-18

Column: 250 \times 4.6 Inertsil ODS-2

Mobile phase: MeCN:water 60:40 containing 1 mM imidazole, pH adjusted to 7.0 with nitric acid

Flow rate: 1

Injection volume: 20

Detector: Chemiluminescence following post-column reaction. The column effluent mixed with the reagent pumped at 1 mL/min and the mixture flowed through a 300 mm \times 0.25 mm ID coil to the detector. (Prepare the reagent by dissolving 112 mg bis(2,4,6-trichlorophenyl) oxalate in 500 mL MeCN, add 8.6 mL 30% hydrogen peroxide, sonicate.), F ex 470 em 530

CHROMATOGRAM

Retention time: 17

Limit of detection: 20 nmole (chemiluminescence), 30 nmole (F)

OTHER SUBSTANCES

Simultaneous: benzylamine, ephedrine, N-isopropylbenzylamine, N-methylphenethylamine, phenylbutylamine, phenylethylamine, phenylpropanolamine, phenylpropylamine

KEY WORDS

derivatization; post-column reaction; comparison with other derivatization reagents

REFERENCE

Hayakawa,K.; Hasegawa,K.; Imaizumi,N.; Wong,O.S.; Miyazaki,M. Determination of amphetamine-related compounds by high-performance liquid chromatography with chemiluminescence and fluorescence detections, *J.Chromatogr.*, **1989**, *464*, 343–352.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 500 μ g/mL solution in MeOH:water 50:50, inject a 5 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 Zorbax C8

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L MeCN:water 20:80. A:B from 100:0 to 0:100 over 30 min. (Purify triethylamine as follows. Wash neutral alumina (Merck) 3 times with 2 bed volumes of pentane, 3 times with 2 bed volumes of dichloromethane, and 3 times with 2 bed volumes of MeOH, allow solvent to evaporate in a fume hood overnight, heat alumina at 130° for 2 h. Prepare a 14 cm column of the washed alumina in a 290 \times 22 tube, pass through a head volume of MeOH, pass through triethylamine. When triethylamine starts to elute discard the first 20 mL, use the next 20 mL, discard the column.)

Flow rate: 2

Injection volume: 5

Detector: UV 210

CHROMATOGRAM

Retention time: 10

OTHER SUBSTANCES

Simultaneous: acetophenone, amphetamine, desipramine, ethylmorphine, imipramine, mefenamic acid, morphine, phenylbutazone, salicylic acid

KEY WORDS

also details of isocratic elution

REFERENCE

Hill,D.W. Evaluation of alkyl bonded silica and solvent phase modifiers for the efficient elution of basic drugs on HPLC, *J.Liq.Chromatogr.*, **1990**, *13*, 3147–3175.

SAMPLE

Matrix: solutions

Sample preparation: Mix 50 μ L of a 200 ppm solution in MeCN:500 mM pH 9.0 borate buffer 50:50 with 25 mg reagent, after 1 min elute with 1 mL hexane:THF 75:25, inject a 5 μ L aliquot. (Reagent was dinitrophenyl carbamate benzotriazole polymeric reagent, synthesized as follows. (Caution! Chloroform, dichloromethane, dioxane, and hydrazine are carcinogenic in experimental animals! DMF may be carcinogenic! 3,5-dinitrobenzoyl chloride and aluminum chloride are corrosive! Nitrobenzene is toxic!) 10 g Dried macroporous polystyrene (Xe-305, Rohm and Haas) + 10 g 3-nitro-4-chlorobenzyl alcohol + 10 g anhydrous aluminum chloride + 50 mL nitrobenzene, heat at 65-70° for 3 days, cool, filter, wash polymer with three 50 mL portions of 1 M HCl in dioxane, with three 50 mL portions of DMF, with three 50 mL portions of MeOH, and with three 50 mL portions of dichloromethane, dry under vacuum at 100°. Reflux 19 g of this polymer in 60 mL hydrazine hydrate:ethylene glycol monoethyl ether 40:60 for 20 h, cool to room temperature, filter off the polymer and wash it thoroughly with water. Suspend the polymer in 100 mL concentrated HCl:dioxane 50:50, reflux for 20 h, filter the polymer and wash it with five 100 mL portions of water, with three 100 mL portions of MeOH, and with three 50 mL portions of ether, dry under vacuum at 80°. Functionalization was 1.17 mmoles/g (Eur.J.Biochem. 1975, 59, 55). Dissolve 3,5-dinitrobenzoyl chloride in the minimum amount of glacial acetic acid, add an equimolar amount of sodium azide, stir for 30 min, dilute with water, filter to obtain 3,5-dinitrobenzoyl azide (Caution! Azides are toxic and potentially explosive!) (J. Liq. Chromatogr. 1986,9, 443). Heat 71 mg 3,5-dinitrobenzoyl azide in 15 mL toluene (dried over calcium hydride) at ??? for 30 min, cool using an ice bath, add 200 mg polymer, allow to warm to room temperature with stirring for 1 h, filter, wash the polymer with four 10 mL portions of warm (40°) dichloromethane, dry under high vacuum for 1 h.)

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m LC-(R)-naphthylurea (Supelco)

Mobile phase: Hexane:EtOH:MeCN 93:7:0.5

Flow rate: 2

Injection volume: 5

Detector: UV 254

CHROMATOGRAM

Retention time: 8.7, 10.7 (enantiomers)

OTHER SUBSTANCES

Simultaneous: amphetamine

KEY WORDS

derivatization; chiral

REFERENCE

Bourque,A.J.; Krull,I.S. Immobilized isocyanates for derivatization of amines for chiral recognition in liquid chromatography with UV detection, *J.Pharm.Biomed.Anal.*, **1993**, *11*, 495–503.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Guard column: 30 \times 2.1 Spheri-5 RP-8

Column: 220 × 2.1 Spheri-5 RP-8

Mobile phase: Gradient. A was 0.08% diethylamine and 0.09% phosphoric acid in water, pH 2.3. B was MeCN:water 90:10 containing 0.08% diethylamine and 0.09% phosphoric acid. A:B 95:5 for 2 min, to 0:100 over 15 min (?), maintain at 0:100 for 5 min.

Column temperature: 50

Flow rate: 0.5

Detector: UV 200

CHROMATOGRAM

Retention time: 7

OTHER SUBSTANCES

Simultaneous: diethylpropion, phenylpropanolamine, ephedrine, amphetamine, phentermine, fenfluramine

Also analyzed: amitriptyline, chlordiazepoxide, chlorpromazine, desalkylflurazepam, desipramine, desmethyldoxepin, diazepam, doxepin, flurazepam, imipramine, mesoridazine, norchlor-diazepoxide, nordiazepam, nortriptyline, oxazepam, prazepam, promazine, thioridazine, thiothixene, trifluoperazine

REFERENCE

Rainin Catalog, C1-94, 1994, p. 7.24.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlordiazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenopropfen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiacol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-stilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclufenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephentyoin, mephesisin, mephobarbital, mepivacaine, mescaline, mesoridazine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methypyrrolon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepi-

nephrine, nortriptyline, noscapine, nyldrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phenacyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopalamine, scopolin, secobarbital, strychnine, sulfacetamide, sufadiazine, sulfadimethoxine, sulfathiazole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233–242.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.5 µm Chiradex (immobilized β-cyclodextrin) (Merck)

Mobile phase: MeOH:100 mM pH 7 ammonium acetate 5:95

Column temperature: 20

Flow rate: 0.5

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 30.28 (S(+)), 31.93 (R(-))

KEY WORDS

chiral

REFERENCE

Rizzi,A.M.; Hirz,R.; Cladrowa-Runge,S.; Jonsson,H. Enantiomeric separation of amphetamine, methamphetamine and ring substituted amphetamines by means of a β-cyclodextrin-chiral stationary phase, *Chromatographia*, **1994**, *39*, 131–137.

SAMPLE

Matrix: solutions

Sample preparation: 1 mL Solution + 500 µL 0.5% sodium 1,2-naphthoquinone-4-sulfonate in water + 500 µL buffer, let stand for 10 min, extract three times with 2 mL aliquots of n-hexane:ethyl acetate 50:50. Combine the organic layers and evaporate them to dryness at 80°, reconstitute with 2 mL MeCN:water 50:50, inject a 50 µL aliquot. (Buffer was 4% sodium bicarbonate adjusted to pH 10 with 10% NaOH.)

HPLC VARIABLES

Column: 250 × 4.5 µm Hypersil ODS C18

Mobile phase: Gradient. MeCN:0.5% propylamine hydrochloride in water from 40:60 to 50:50 over 2.5 min, to 70:30 over 1 min, maintain at 70:30 for 4.5 min.

Flow rate: 1

Injection volume: 50

Detector: UV 280

CHROMATOGRAM

Retention time: 7

OTHER SUBSTANCES**Extracted:** amphetamine

KEY WORDS

derivatization

REFERENCE

Herráez-Hernández,R.; Campíns-Falcó,P.; Sevillano-Cabeza,A. On-line derivatization into precolumns for the determination of drugs by liquid chromatography and column switching: Determination of amphetamines in urine, *Anal.Chem.*, **1996**, 68, 734–739.

SAMPLE**Matrix:** solutions

Sample preparation: 500 μ L Solution + 250 μ L buffer + 250 μ L 20 mM 9-fluorenylmethyl chloroformate in MeCN, mix, add 1 mL MeCN, inject a 50 μ L aliquot. (Prepare buffer by adjusting the pH of 4% sodium bicarbonate to 10 with 10% NaOH.)

HPLC VARIABLES**Column:** 125 \times 4.5 μ m LiChrospher 100 RP 18

Mobile phase: Gradient. MeCN:water from 40:60 to 50:50 over 2.5 min, to 70:30 over 2.5 min, to 100:0 over 5 min.

Flow rate: 1.5**Injection volume:** 50**Detector:** F ex 264 em 313

CHROMATOGRAM**Retention time:** 8.8

OTHER SUBSTANCES**Extracted:** amphetamine

KEY WORDS

derivatization

REFERENCE

Herráez-Hernández,R.; Campíns-Falcó,P.; Sevillano-Cabeza,A. On-line derivatization into precolumns for the determination of drugs by liquid chromatography and column switching: Determination of amphetamines in urine, *Anal.Chem.*, **1996**, 68, 734–739.

SAMPLE**Matrix:** urine

Sample preparation: Condition a Sep-Pak C18 SPE cartridge with 5 mL MeOH, 3 mL MeCN: 10 mM ammonium acetate 40:60 adjusted to pH 3 with acetic acid, and 5 mL water. 5 mL Urine + 5 mL 500 mM ammonium acetate, adjusted to pH 9.5 with ammonia, mix, add to the SPE cartridge, wash with 20 mL 5 mM pH 9.5 ammonium acetate, wash with 0.5 mL water. Elute with 2 mL MeCN:10 mM ammonium acetate 40:60 adjusted to pH 3 with acetic acid, inject a 50 μ L aliquot of the eluate.

HPLC VARIABLES**Column:** 150 \times 4.6 L-column ODS (Chemical Inspection & Testing Institute, Tokyo)

Mobile phase: Gradient. MeCN:100 mM ammonium acetate 0:100 for 1 min, to 40:60 over 20 min.

Flow rate: 1**Injection volume:** 50

Detector: UV 210; MS Shimadzu model QP-1100EX thermospray, vaporizer temperature from 170 to 150° over 20 min. SIM, m/z 150

CHROMATOGRAM**Retention time:** 15**Limit of detection:** 2–40 ng/mL

OTHER SUBSTANCES

Extracted: 6-acetylmorphine, amphetamine, benzoylecgonine, cocaine, ephedrine, methylephedrine, morphine, morphine-3-glucuronide, morphine-6-glucuronide

KEY WORDS

SPE

REFERENCE

Tatsuno, M.; Nishikawa, M.; Katagi, M.; Tsuchihashi, H. Simultaneous determination of illicit drugs in human urine by liquid chromatography-mass spectrometry, *J. Anal. Toxicol.*, **1996**, 20, 281-286.

SAMPLE

Matrix: urine

Sample preparation: 500 μ L Urine + N-ethylnordiazepam + chlorpheniramine + 100 μ L buffer, centrifuge at 11000 g for 30 s, inject a 500 μ L aliquot onto column A with mobile phase A, after 0.6 min backflush column A with mobile phase A to waste for 1.6 min, elute column A with 250 μ L mobile phase B, with 200 μ L mobile phase C, and with 1.15 mL mobile phase D. Elute column A to waste until drugs start to emerge then elute onto column B. Elute column B to waste until drugs started to emerge, then elute onto column C. When all the drugs have emerged from column B remove it from the circuit, elute column C with mobile phase D, monitor the effluent from column C. Flush column A with 7 mL mobile phase E, with mobile phase D, and mobile phase A. Flush column B with 5 mL mobile phase E then with mobile phase D. (Buffer was 6 M ammonium acetate adjusted to pH 8.0 with 2 M KOH.)

HPLC VARIABLES

Column: A 10×2.1 12-20 μ m PRP-1 spherical poly(styrene-divinylbenzene) (Hamilton); B 10×3.2 11 μ m Aminex A-28 (Bio-Rad); C 25×3.2 5 μ m C8 (Phenomenex) + 150×4.6 5 μ m silica (Macherey-Nagel)

Mobile phase: A 0.1% pH 8.0 potassium borate buffer; B 6 mM KH_2PO_4 containing 5 mM tetramethylammonium hydroxide, and 2 mM dimethyloctylamine, pH adjusted to 6.50 with phosphoric acid; C MeCN:buffer 40:60 (Buffer was 6 mM KH_2PO_4 containing 5 mM tetramethylammonium hydroxide, and 2 mM dimethyloctylamine, pH adjusted to 6.50 with phosphoric acid.); D MeCN:buffer 33:67 (Buffer was 6 mM KH_2PO_4 containing 5 mM tetramethylammonium hydroxide, and 2 mM dimethyloctylamine, pH adjusted to 6.50 with phosphoric acid.); E MeCN:buffer 70:30 (Buffer was 6 mM KH_2PO_4 containing 5 mM tetramethylammonium hydroxide, and 2 mM dimethyloctylamine, pH adjusted to 6.50 with phosphoric acid.)

Column temperature: ambient (column A), 40 (columns B and C)

Flow rate: A 5; B-E 1

Injection volume: 500

Detector: UV 210, UV 235

CHROMATOGRAM

Retention time: k' 3.1

Internal standard: N-ethylnordiazepam (k' 2.1), chlorpheniramine (k' 5.9)

Limit of detection: 300 ng/mL

OTHER SUBSTANCES

Extracted: methadone, imipramine, flurazepam, amitriptyline, morphine, codeine, hydromorphone, hydrocodone, caffeine, cotinine, benzoylecgonine, secobarbital, oxazepam, phenobarbital, nordiazepam, diazepam, phenylpropanolamine, phentermine, amphetamine, phenmetrazine, lidocaine, ephedrine

Interfering: pentazocine, desipramine, nortriptyline, diphenhydramine

KEY WORDS

column-switching

REFERENCE

Binder, S.R.; Regalia, M.; Biaggi-McEachern, M.; Mazhar, M. Automated liquid chromatographic analysis of drugs in urine by on-line sample cleanup and isocratic multi-column separation, *J. Chromatogr.*, **1989**, 473, 325-341.

SAMPLE**Matrix:** urine

Sample preparation: 500 μ L Urine + 100 μ L 25 μ g/mL N-n-propylaniline + 6 mL pH 10.0 carbonate buffer + 15 mL water, add mixture to an Extrelut SPE cartridge, let stand for 20 min, elute with 40 mL hexane:ethyl acetate 90:10. Add the eluate to 3 mL 100 mM sulfuric acid and 500 mg NaCl, stir for 20 min, centrifuge at 1000 g for 5 min. Remove the lower layer and add it to 3 mL 2.5 M NaOH and 20 μ L benzoyl chloride, stir vigorously for 30 min. Extract the mixture with 1.5 mL chloroform. Wash the chloroform layer twice with 5 mL water and evaporate it to dryness at 40°, reconstitute the residue in 200 μ L hexane:isopropanol 90:10, inject an aliquot.

HPLC VARIABLES**Column:** 250 \times 4.6 Chiralcel OB + 250 \times 4.6 Chiralcel OJ**Mobile phase:** Hexane:isopropanol 90:10**Column temperature:** 48**Flow rate:** 1-1.4**Detector:** UV 220

CHROMATOGRAM**Retention time:** 13 (d), 15 (l)**Internal standard:** N-n-propylaniline (11)**Limit of detection:** 25 ng

OTHER SUBSTANCES**Extracted:** amphetamine

KEY WORDSrat; SPE; derivatization; chiral

REFERENCE

Nagai, T.; Kamiyama, S. Assay of the optical isomers of methamphetamine and amphetamine in rat urine using high-performance liquid chromatography with chiral cellulose-based columns, *J. Chromatogr.*, **1990**, 525, 203-209.

SAMPLE**Matrix:** urine

Sample preparation: 200-500 μ L Rat urine + 200-500 μ L pH 3.8 acetate buffer + 25 μ L 40 μ g/mL β -glucuronidase and 20 μ g/mL arylsulfatase (Merck), heat at 37° for 24 h, add 100 μ L 25 μ g/mL 3-methoxytyramine in water, add 100 μ L water, adjust pH to 9.0 with 1.9 M sodium carbonate, add to an Extrelut SPE cartridge, let stand for 20 min, elute with 6 mL ethyl acetate. Add the eluate to 1 mL 100 mM sulfuric acid, extract. Add the aqueous layer to 3 mL 2.5 M NaOH, add 25 μ L benzoyl chloride, extract with 5 mL ethyl acetate. Wash the ethyl acetate layer with water, evaporate to dryness under a stream of nitrogen at 40°, reconstitute the residue in 50 μ L EtOH, 3.5 mL 50 mM pH 8.0 Tris/HCl buffer, and 35 μ L esterase (Type 1 porcine liver, Sigma). Heat at 25° for 45 min, add to an activated Sep-Pak C18 SPE cartridge, wash with 5 mL water, elute with 5 mL acetone. Evaporate the eluate to dryness under a stream of nitrogen at 40°, reconstitute the residue in 200 μ L hexane:EtOH 89:11, inject a 20 μ L aliquot.

HPLC VARIABLES**Column:** 250 \times 4.6 Chiralcel OB + 250 \times 4.6 Chiralcel OJ**Mobile phase:** n-Hexane:EtOH 89:11**Column temperature:** 48**Flow rate:** 1.4**Injection volume:** 20**Detector:** UV 220

CHROMATOGRAM**Retention time:** 10 (D), 11.5 (L)**Internal standard:** 3-methoxytyramine (54)

OTHER SUBSTANCES

Extracted: metabolites, amphetamine

KEY WORDS

rat; SPE; derivatization; chiral

REFERENCE

Nagai,T.; Kamiyama,S. Simultaneous HPLC analysis of optical isomers of methamphetamine and its metabolites, and stereoselective metabolism of racemic methamphetamine in rat urine, *J.Anal.Toxicol.*, **1991**, *15*, 299–304.

SAMPLE

Matrix: urine

Sample preparation: Adjust pH of 3 mL urine to 11 with 10 M KOH, add to an Extrelut 3 column, let stand for 10 min, elute with 15 mL n-hexane into a tube containing one drop of 3 M HCl. Evaporate the eluate to dryness under a stream of nitrogen at 35°. Add 1.5 mL 8% sodium bicarbonate in water and 1 mL 0.5% sodium naphthoquinone-4-sulfonate in water to the residue, heat at 70° for 20 min, cool, extract with 5 mL carbon tetrachloride. Evaporate the eluate to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100 µL mobile phase, inject a 20 µL aliquot.

HPLC VARIABLES

Guard column: 4 × 4.5 µm Lichrospher 100 RP8

Column: 250 × 4.5 µm Lichrospher 100 RP8

Mobile phase: MeCN:buffer 55:45 (Buffer was 1.361 g KH₂PO₄ in 950 mL, add 1.3 mL methanesulfonic acid, adjust pH to 3 with 5 M KOH, make up to 1 L with water.)

Flow rate: 1

Injection volume: 20

Detector: UV 460

CHROMATOGRAM

Retention time: 7.8

Limit of detection: 60 ng/mL

OTHER SUBSTANCES

Extracted: amphetamine

Noninterfering: acetaminophen, aspirin, amitriptyline, buprenorphine, caffeine, carbamazepine, chlorpromazine, desipramine, dextromethorphan, doxepin, ephedrine, fenfluramine, imipramine, lidocaine, loxapine, meperidine, methadone, methaqualone, naloxone, naltrexone, nicotine, orphenadrine, oxycodone, papaverine, pentazocine, phendimetrazine, phenmetrazine, phentermine, phenylpropanolamine, phenytoin, primidone, procaine, promethazine, propoxyphene, propyphenazone, theobromine, theophylline, trazodone, triflupromazine, trimethoprim, trimipramine

KEY WORDS

SPE; derivatization

REFERENCE

Ferrara,S.D.; Tedeschi,L.; Frison,G.; Castagna,F. Solid-phase extraction and HPLC-UV confirmation of drugs of abuse in urine, *J.Anal.Toxicol.*, **1992**, *16*, 217–222.

SAMPLE

Matrix: urine

Sample preparation: Condition a 100 mg Adsorbex SCX cation-exchange SPE cartridge (Merck) with 2 mL MeOH, 1 mL water, and 1 mL 17 mM KH₂PO₄, do not allow to dry. Centrifuge urine at 2000 g for 5 min. 1 mL Urine + 500 µL 50 mM KH₂PO₄, sonicate for 1 min, add to the SPE cartridge, rinse vial with 50 µL 50 mM KH₂PO₄ and add to cartridge, dry cartridge for 1 min, wash with three 500 µL portions of 17 mM KH₂PO₄, wash with 1 mL MeOH, dry under vacuum for 1 min, elute with four 500 µL portions of MeOH:7.3% HCl (97.5:2.5) at a flow rate of 0.5 mL/min, inject a 10 µL aliquot.

HPLC VARIABLES

Column: 125 × 4 3 µm Spherisorb ODS-1

Mobile phase: Gradient. A was water containing 5 mL (8.5 g) 85% orthophosphoric acid and 280 µL (0.22 g) hexylamine per liter. B was MeCN containing 100 mL water, 5 mL (8.5 g) 85% orthophosphoric acid, and 280 µL (0.22 g) hexylamine per liter. A:B 94.5:5.5 for 10.6 min, then to 61:39 over 11 min.

Column temperature: 40

Flow rate: 0.8

Injection volume: 10

Detector: UV 198

CHROMATOGRAM

Retention time: 9

Limit of detection: 30 ng/mL

OTHER SUBSTANCES

Extracted: 3,4-methylenedioxyamphetamine, amphetamine, 4-methoxyamphetamine, phentermine, 3,4-methylenedioxymethamphetamine, 5-methoxy-3,4-methylenedioxyamphetamine, 3,4,5-trimethoxyamphetamine, 3,4-methylenedioxyethylamphetamine, 2,5-dimethoxyamphetamine, 4-bromo-2,5-dimethoxyphenylethylamine, 2,5-dimethoxy-4-methylamphetamine, 4-bromo-2,5-dimethoxyamphetamine, 2,5-dimethoxy-4-ethylamphetamine, mescaline, methoxamine

KEY WORDS

SPE

REFERENCE

Helmlin, H.-J.; Brenneisen, R. Determination of psychotropic phenylalkylamine derivatives in biological matrices by high-performance liquid chromatography with photodiode-array detection, *J. Chromatogr.*, **1992**, 593, 87-94.

SAMPLE

Matrix: urine

Sample preparation: 5 mL Urine + 4.5 mL MeCN + 500 µL 1 M KOH, centrifuge at 2500 rpm for 10 min, filter (45 µm) the supernatant. Inject on to column A at 180 µL/min 25 µL 50 mM KOH in MeCN:water 20:80, 50 µL filtrate, 25 µL 50 mM KOH in MeCN:water 20:80, and 200 µL MeCN:water 20:80, backflush the contents of column A on to column B with mobile phase, after 18 s remove column A from the circuit, elute column B with mobile phase, monitor the effluent from column B. Wash column A with 400 µL MeCN.

HPLC VARIABLES

Column: A 20 × 2 polymeric reagent; B 250 × 4.6 5 µm Supelcosil LC-18-DB (with a guard column) (Prepare polymeric reagent as follows. Prepare a porous rigid resin using a divinylbenzene:ethylstyrene:styrene 24:6:70 mixture with trimethylsilyl modified silica (102 Å average pore size, 1.08 mL/g pore volume, 366 m²/g surface area, 16-20 µm irregular particle shape, IMPAQ RG 1020 Si silica, PQ Co., Valley Forge PA). Further preparation details are not given but a typical procedure given in the cited reference is as follows. Aerate a mixture of 10 g modified silica in 100 mL water with nitrogen for 15 min, add 10 mL styrene:80% divinylbenzene:t-butyl peroxybenzoate 49:49:2 (remove preservative by passing through a butylcatechol remover (Scientific Polymer, Ontario NY), shake vigorously at room temperature for 4 h, add 150 mL 0.75% polyvinyl alcohol, shake for 4 h, heat at 120° for 24 h while shaking on a Parr instrument, cool to room temperature, filter, wash with 100 mL water, wash with 50 mL MeOH. Add the solid to 500 mL 3 M NaOH in MeOH:water 40:60, shake at room temperature for 14 h (to dissolve the silica), filter, wash with water until the washings are neutral, wash with 100 mL MeOH, dry at 60°. The polymer has similar properties to the template silica (US Pat. 4 933 372 (1990)). Soxhlet extract the resin with dioxane for 8 h (Caution! Dioxane is a carcinogen!). Add 25 g aluminum trichloride in 300 mL dry nitrobenzene to 50 g resin and 100 g 4-chloro-3-nitrobenzoyl chloride, stir mechanically at 60° for 5 h, pour into a mixture of 150 mL DMF, 100 mL concentrated HCl, and 150 g ice, filter. Wash the solid with 300 mL portions of DMF:water 75:25 until the washings are colorless, wash with warm (60°) DMF, wash with six 300 mL portions of dichloromethane:MeOH 2:1. Stir the product in 130 mL 40% benzyltrimethylammonium hydroxide in water, 130 mL water, and 260 mL dioxane at 90° for 8 h, filter, repeat the process. Wash the product with four portions of warm (60°) dioxane. Stir the solid

with 30 mL acetic acid for 15 min, filter. Wash the solid with dioxane until the washings are neutral, wash with six 300 mL portions of dichloromethane:MeOH 2:1 to give a nitrobenzophenol-substituted polymer (J. Org. Chem. 1984, 49, 924). Heat 4 g 9-fluoreneacetic acid, 3.9 mL oxalyl chloride, 30 mL benzene (dried over anhydrous sodium sulfate, Caution! Benzene is a carcinogen!), and 3 drops of triethylamine at 55° for 1 h, evaporate under reduced pressure to remove oxalyl chloride, dissolve the product in 35 mL dichloromethane to give a 120 mg/mL solution of 9-fluoreneacetyl chloride, dilute to obtain a 2 mM solution. Stir 1.3 g nitrobenzophenol-substituted polymer, 4.2 mL 2 mM 9-fluoreneacetyl chloride solution, 300 μ L triethylamine, and 20 mL dichloromethane at room temperature for 1 h, filter, wash with three 20 mL portions of MeCN to obtain the reagent, polymer-bound nitrobenzophenol 9-fluoreneacetate (J. Chromatogr. 1992, 609, 103).

Mobile phase: Gradient. MeCN:water 50:50 for 3.5 min, to 70:30 over 12 min, maintain at 70:30 for 2.5 min, return to initial conditions over 1 min. (Place a 100 \times 4.6 30-40 μ m silica column before the injector.)

Column temperature: 60 (column A only)

Injection volume: 25-50

Detector: F ex 254 em 305-395

CHROMATOGRAM

Retention time: 13.6

Limit of quantitation: 25 ng/mL

KEY WORDS

derivatization; column-switching

REFERENCE

Bourque, A.J.; Krull, I.S.; Feibush, B. Automated HPLC analyses of drugs of abuse via direct injection of biological fluids followed by simultaneous solid-phase extraction and derivatization with fluorescence detection, *Biomed. Chromatogr.*, **1994**, 8, 53-62.

SAMPLE

Matrix: urine

Sample preparation: Condition a 100 mg Bond-Elut C18 SPE cartridge with 500 μ L MeOH and 500 μ L water. Adjust pH of urine to 10, centrifuge at 1500 g. 2 mL Supernatant + 100 μ L 75 μ g/mL β -phenylethylamine hydrochloride in water, add to the SPE cartridge, wash with 2.5 mL water, elute with 2 mL MeOH, evaporate the eluate to dryness. Reconstitute in water, add 500 μ L 8% sodium bicarbonate, add 500 μ L 0.5% 1,2-naphthoquinone-4-sulfonic acid sodium salt, make up to 1.5 mL with water, heat at 70° for 20 min, cool, add an equal volume of chloroform, shake for 2 min, centrifuge at 1500 g for 5 min. Remove the organic layer and dry it over anhydrous sodium sulfate, filter (0.45 μ m), inject a 25 μ L aliquot of the filtrate.

HPLC VARIABLES

Column: 125 \times 4.5 μ m LiChrospher Si-60

Mobile phase: EtOH:chloroform:ethyl acetate:n-hexane 1:22:32:45

Flow rate: 2

Injection volume: 25

Detector: UV 280

CHROMATOGRAM

Retention time: 2.6

Internal standard: β -phenylethylamine hydrochloride (4.9)

OTHER SUBSTANCES

Extracted: amphetamine

KEY WORDS

SPE; normal phase; derivatization

REFERENCE

Campins Falcó, P.; Molins Legua, C.; Herráez Hernández, R.; Sevillano Cabeza, A. Improved amphetamine and methamphetamine determination in urine by normal-phase high-performance liquid chromatography with sodium 1,2-naphthoquinone 4-sulphonate as derivatizing agent and solid-phase extraction for sample clean-up, *J. Chromatogr. B*, **1995**, 663, 235-245.

SAMPLE**Matrix:** urine

Sample preparation: Condition a Bond Elut C18 SPE cartridge with 1 mL MeOH and 1 mL 50 mM pH 11 phosphate buffer. 500 μ L Urine + 500 μ L 2000 U/mL β -glucuronidase with sulfatase activity (Type H-1, Sigma) in 100 mM pH 5 acetate buffer, heat at 37° overnight, add 500 mg NaCl, add 500 μ L 50 mM pH 11 potassium phosphate buffer, adjust pH to 11 with ammonium hydroxide, mix. Add the mixture to the SPE cartridge, wash with 1 mL 50 mM pH 11 potassium phosphate buffer, wash with 1 mL freshly prepared MeOH:water 30:70, wash with 1 mL MeCN, elute with 1 mL freshly prepared MeCN:acetic acid 98:2, elute with 1 mL MeCN:HCl 98:2. Combine the eluates and evaporate them to dryness under a stream of air at room temperature, reconstitute the residue in mobile phase (?), inject a 10 μ L aliquot.

HPLC VARIABLES**Guard column:** phenyl**Column:** Microsorb phenyl**Mobile phase:** MeCN:MeOH:50 mM pH 3 potassium phosphate 5:10:85**Flow rate:** 1**Injection volume:** 10**Detector:** UV 215

CHROMATOGRAM**Retention time:** 15.6**Internal standard:** methamphetamine

OTHER SUBSTANCES**Extracted:** amphetamine

KEY WORDSSPE; rat; methamphetamine in IS

REFERENCE

Law,M.Y.L.; Moody,D.E. Simultaneous quantitation of amphetamine and 4'-hydroxyamphetamine by high performance liquid chromatography, *J.Liq.Chromatogr.*, **1995**, 18, 2029–2043.

SAMPLE**Matrix:** urine

Sample preparation: Condition an Extra-Sep C18 SPE cartridge (Teknokroma) with 1 mL MeOH and 1 mL buffer. Adjust pH of 2 mL urine to ca. 10 with 100 μ L concentrated ammonium hydroxide, add 5 μ g β -phenylethylamine, add to the SPE cartridge, wash with 5 mL water, wash with 1 mL MeCN, elute with 2 mL MeOH. Add 100 μ L EtOH:concentrated HCl 6:1 to the eluate, evaporate to dryness. Reconstitute with 1 mL buffer and 1 mL 0.5% 1,2-naphthoquinone-4-sulfonic acid sodium salt, let stand at room temperature for 10 min, add 2 mL n-hexane:ethyl acetate 50:50, shake for 2 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 500 μ L MeCN:water 50:50, inject a 50 μ L aliquot. (Buffer was 1% aqueous sodium bicarbonate adjusted to pH 10 with 5 M NaOH.)

HPLC VARIABLES**Column:** 250 \times 4 5 μ m Hypersil ODS-C18**Mobile phase:** Gradient. MeCN:0.5% propylamine in water from 40:60 to 50:50 over 2.5 min, to 70:30 over 1 min, maintain at 70:30.**Flow rate:** 1**Injection volume:** 50**Detector:** UV 280

CHROMATOGRAM**Retention time:** 6.3**Internal standard:** β -phenylethylamine (4.1)**Limit of detection:** 2 ng/mL**Limit of quantitation:** 10 ng/mL

OTHER SUBSTANCES

Extracted: amphetamine

KEY WORDS

SPE; derivatization

REFERENCE

Molins Legua,C.; Campíns Falcó,P.; Sevillano Cabeza,A. Amphetamine and methamphetamine determination in urine by reversed-phase high-performance liquid chromatography with sodium 1,2-naphthoquinone 4-sulfonate as derivatizing agent and solid-phase extraction for sample clean-up, *J.Chromatogr.B*, **1995**, 672, 81–88.

SAMPLE

Matrix: urine

Sample preparation: 100-300 μ L Urine + 100 μ L 1.5 M NaOH + 5 μ g IS, make up to 1 mL with water, add to an Extrelut 1 SPE cartridge, let stand for 20 min, elute with 6 mL benzene (Caution! Benzene is a carcinogen!). Add the eluate to 1 mL 100 mM sulfuric acid, extract. Remove the aqueous layer and add it to 3 mL 1.5 M NaOH, add 5 μ L benzoyl chloride, vortex vigorously, extract twice with 3 mL portions of n-hexane. Combine the organic layers and wash them twice with 3 mL portions of water, evaporate the organic to dryness under a stream of nitrogen at 40°, reconstitute the residue in 200 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: Chiralcel OB-H

Mobile phase: n-Hexane:isopropanol 90:10

Column temperature: 55

Flow rate: 1

Injection volume: 20

Detector: UV 220

CHROMATOGRAM

Retention time: 6.5 (D), 10 (L)

Internal standard: l-p-methoxyamphetamine (12)

Limit of detection: 30 ng

OTHER SUBSTANCES

Extracted: amphetamine, ethylamphetamine

KEY WORDS

rat; derivatization; SPE; chiral

REFERENCE

Nagai,T.; Kamiyama,S.; Matsushima,K. Analysis of time-lapse changes of d- and l-enantiomers of racemic ethylamphetamine and stereoselective metabolism in rat urine by HPLC determination, *J.Anal.Toxicol.*, **1995**, 19, 225–228.

SAMPLE

Matrix: urine

Sample preparation: Dilute urine 20-fold or more. 500 μ L Diluted urine + 50 μ L 500 ng/mL (+)-2,5-dimethoxyamphetamine hydrochloride + 100 μ L 500 mM NaOH + 2 mL benzene (Caution! Benzene is a carcinogen!), shake for 15 min, centrifuge at 1200 g for 5 min. Remove 1.8 mL of the organic phase and add it to 220 μ L 50 mM HCl, shake for 15 min, centrifuge at 1200 g for 5 min. Remove 200 μ L of the aqueous phase and add it to 40 μ L 250 mM NaOH, add 50 μ L 330 mM pH 7.8 phosphate buffer, add 250 μ L MeCN, add 25 μ L 3 mM (-)-1-(9-fluorenyl)ethyl chloroformate, let stand at room temperature for 24 h, add 30 μ L 100 mM glycine in water, let stand for 30 min, add 750 μ L pentane, vortex for 2 min, centrifuge at 1200 g for 5 min. Remove the organic layer and evaporate it to dryness in a centrifugal evaporator at room temperature, reconstitute the residue in 300 μ L MeCN:water 50:50, inject a 100 μ L aliout.

HPLC VARIABLES

Column: 150 \times 4.6 3 μ m Adsorbosphere HS C18

Mobile phase: MeCN:THF:20 mM pH 3.6 sodium acetate buffer 25:21:54

Flow rate: 1

Injection volume: 100

Detector: F ex 265 em 330

CHROMATOGRAM

Retention time: 42 (L), 44 (D)

Internal standard: (+)-2,5-dimethoxyamphetamine (33)

Limit of detection: 5 ng/mL

OTHER SUBSTANCES

Extracted: amphetamine

KEY WORDS

rat; derivatization; chiral

REFERENCE

Sukbuntherng,J.; Hutchaleelaha,A.; Chow,H.-H.; Mayersohn,M. Separation and quantitation of the enantiomers of methamphetamine and its metabolites in urine by HPLC: Precolumn derivatization and fluorescence detection, *J.Anal.Toxicol.*, **1995**, 19, 139–147.

SAMPLE

Matrix: urine

Sample preparation: Condition a 200 mg Bond Elut C18 SPE cartridge with 1 mL MeOH and 1 mL 1% pH 10 sodium bicarbonate buffer. 2 mL Urine + 400 μ L 8% pH 10 sodium bicarbonate buffer, mix, centrifuge at 1500 g for 2 min, add a 2 mL aliquot of the supernatant to the SPE cartridge, wash with 3 mL water, pass 500 μ L 2% sodium 1,2-naphthoquinone 4-sulfonate through the cartridge, pass 500 μ L 1% pH 10 sodium bicarbonate buffer through the cartridge, let stand at room temperature for 15 min, wash with 3 mL water, elute with 1 mL MeCN: water 50:50, inject a 20 μ L aliquot of the eluate.

HPLC VARIABLES

Column: 250 \times 4.5 μ m Hypersil ODS

Mobile phase: Gradient. MeCN:buffer from 40:60 to 50:50 over 2.5 min, to 70:30 over 0.5 min, maintain at 70:30 for 1.5 min, to 85:15 over 1 min, maintain at 85:15 for 1.5 min. (Buffer was 5 mL/L propylamine in water.)

Flow rate: 1

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 5.6

Internal standard: β -phenylethylamine (3.6)

Limit of detection: 400 ng/mL

OTHER SUBSTANCES

Extracted: amphetamine

KEY WORDS

derivatization; SPE

REFERENCE

Campins-Falcó,P.; Sevillano-Cabeza,A.; Molins-Lagua,C.; Kohlmann,M. Amphetamine and methamphetamine determination in urine by reversed-phase high-performance liquid chromatography with simultaneous sample clean-up and derivatization with 1,2-naphthoquinone 4-sulphonate on solid-phase cartridges, *J.Chromatogr.B*, **1996**, 687, 239–246.

SAMPLE

Matrix: urine

Sample preparation: Inject 15 μ L urine, inject a mixture of 5 μ L 20 mM 9-fluorenylmethyl chloroformate in MeCN and 45 μ L water, and inject 10 μ L buffer on to column A and elute to

waste with mobile phase A. After 2.8 min backflush the contents of column A on to column B with mobile phase B and start the gradient, monitor the effluent from column B. At the end of the run condition column A with 1 mL mobile phase A. (Buffer was 4% sodium bicarbonate adjusted to pH 10 with 10% NaOH.)

HPLC VARIABLES

Column: A 20 × 2.1 30 µm Hypersil ODS-C18; B 125 × 4 5 µm LiChrospher 100 PR-C18

Mobile phase: A water; B Gradient. MeCN:water from 40:60 to 70:30 over 15 min. to 100:0 over 5 min.

Flow rate: A 0.35; B 1

Injection volume: 15

Detector: F ex 264 em 313

CHROMATOGRAM

Retention time: 17.8

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Extracted: amphetamine, ephedrine, norephedrine, 3-phenylpropylamine, pseudoephedrine

KEY WORDS

column-switching; derivatization; on-column derivatization

REFERENCE

Herráez-Hernández,R.; Campíns-Falcó,P.; Sevillano-Cabeza,A. Determination of amphetamine and related compounds in urine using on-line derivatization in octadecyl silica columns with 9-fluorenylmethyl chloroformate and liquid chromatography, *J.Chromatogr.B*, **1996**, 679, 69–78.

SAMPLE

Matrix: urine

Sample preparation: Inject 50 µL urine on to column A and elute to waste with mobile phase A, after 2 min inject a mixture of 25 µL 0.5% sodium 1,2-naphthoquinone-4-sulfonate in water and 25 µL buffer on to column A, stop the flow of mobile phase A, after 10 min start pump A, after 5 min backflush the contents of column A on to column B with mobile phase B and start the gradient, monitor the effluent from column B. After each run flush column A with ethyl acetate for 1 min, n-hexane for 1 min, and ethyl acetate for 1 min, re-equilibrate with mobile phase A. (Buffer was 4% sodium bicarbonate adjusted to pH 10 with 10% NaOH.)

HPLC VARIABLES

Column: A 20 × 2.1 30 µm Hypersil ODS-C18; B 250 × 4 5 µm Hypersil ODS C18

Mobile phase: A water; B Gradient. MeCN:0.5% propylamine hydrochloride in water from 40:60 to 50:50 over 2.5 min, to 70:30 over 1 min, maintain at 70:30 for 4.5 min.

Flow rate: 1

Injection volume: 50

Detector: UV 280

CHROMATOGRAM

Retention time: 7

Limit of detection: 25 ng/mL

OTHER SUBSTANCES

Extracted: amphetamine

KEY WORDS

column-switching; derivatization; on-column derivatization

REFERENCE

Herráez-Hernández,R.; Campíns-Falcó,P.; Sevillano-Cabeza,A. On-line derivatization into precolumns for the determination of drugs by liquid chromatography and column switching: Determination of amphetamines in urine, *Anal.Chem.*, **1996**, 68, 734–739.

SAMPLE**Matrix:** urine**Sample preparation:** Condition a Bond Elut SCX SPE cartridge with 10 mL MeOH:aqueous ammonia 98:2, 10 mL MeOH, and 30 mL water. 10 mL Urine + 10 mL 50 mM pH 7.8 potassium phosphate buffer, mix, add to the SPE cartridge, wash with 2 mL water, wash with 10 mL MeOH, elute with 3 mL MeOH:2% aqueous ammonia 98:2. Add 30 μ L glacial acetic acid and 30 μ L 1 mg/mL N-ethylaniline in MeOH to the eluate, inject a 5-50 μ L aliquot.

HPLC VARIABLES**Column:** 150 \times 6 ULTRON ES-PhCD β -cyclodextrin phenylcarbamate-bonded silica (Shinwa Chemical Industries, Kyoto)**Mobile phase:** MeCN:MeOH:buffer 10:30:60 (Buffer was 50 mM pH 6.0 potassium phosphate for UV detection or 100 mM pH 6.0 ammonium acetate for MS detection.)**Column temperature:** 25**Flow rate:** 1**Injection volume:** 5-50**Detector:** UV 220, MS, Shimadzu LCMS-QP1100EX, thermospray, positive ion mode, filament off, vaporizer 230°, ion source 300°, m/z 150

CHROMATOGRAM**Retention time:** 12 (D, UV), 15 (L, UV), 15 (D, MS), 17.5 (L, MS)**Internal standard:** N-ethylaniline (17, UV only)**Limit of detection:** 50 ng/mL (D, UV detection), 100 ng/mL (L, UV detection), 10 ng/mL (D, scan MS), 20 ng/mL (L, scan MS), 0.8 ng/mL (D, SIM MS), 1.0 ng/mL (L, SIM MS)

OTHER SUBSTANCES**Extracted:** metabolites, amphetamine, p-hydroxymethamphetamine

KEY WORDSSPE; chiral

REFERENCEKatagi,M.; Nishioka,H.; Nakajima,K.; Tsuchihashi,H.; Fujima,H.; Wada,H.; Nakamura,K.; Makino,K. Direct high-performance liquid chromatographic and high-performance liquid chromatographic-thermospray-mass spectrometric determination of enantiomers of methamphetamine and its main metabolites amphetamine and p-hydroxymethamphetamine in human urine, *J.Chromatogr.B*, **1996**, 676, 35-43.

SAMPLE**Matrix:** urine**Sample preparation:** 2 mL Urine + 20 μ L 100 μ M 1-phenylethylamine in water + 400 μ L concentrated HCl, heat at 80° for 1 h, cool, neutralize with 600 μ L 25% ammonia, add 5 mL 10% sodium carbonate solution, add 2 mL 500 mM pH 10.5 sodium borate buffer, add 2 mL chloroform:isopropanol 75:25, vortex for 1 min, centrifuge at 12.5° at 1500 g for 10 min, repeat the extraction. Combine the organic layers and remove a 100 μ L aliquot, add 10 μ L acetic acid to the aliquot, evaporate it to dryness under a stream of nitrogen at room temperature, reconstitute the residue in 50 μ L carbonate buffer, add 50 μ L 10 mM fluorescein-4-isothiocyanate in EtOH, mix, heat in the dark at 80° for 15 min, inject a 20 μ L aliquot. (Prepare carbonate buffer by adjusting the pH of 200 mM sodium bicarbonate to 9.0 with 200 mM sodium carbonate.)

HPLC VARIABLES**Column:** 250 \times 4.6 5 μ m Daisopak SP-120-5-ODS (Daiso, Osaka)**Mobile phase:** Gradient. MeCN:20 mM pH 7.9 sodium phosphate buffer 20:80 for 16 min then 24:76 (step-gradient).**Flow rate:** 0.8**Injection volume:** 20**Detector:** F ex 496 em 518

CHROMATOGRAM**Retention time:** 35.2**Internal standard:** 1-phenylethylamine (26.6)**Limit of detection:** 5.5 nM

OTHER SUBSTANCES

Extracted: metabolites, amphetamine, norepinephrine

KEY WORDS

derivatization

REFERENCE

Al-Dirbashi,O.; Kuroda,N.; Akiyama,S.; Nakashima,K. High-performance liquid chromatography of methamphetamine and its related compounds in human urine following derivatization with fluorescein isothiocyanate, *J.Chromatogr.B*, **1997**, 695, 251-258.

Methapyrilene


Molecular formula: C₁₄H₁₉N₃S

Molecular weight: 261.39

CAS Registry No.: 91-80-5, 33032-12-1 (fumarate), 135-23-9 (HCl)

Merck Index: 6027

Lednicer No.: 1 54

**SAMPLE**

Matrix: solutions

Sample preparation: Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 125 × 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 4.6

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipiprone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethiopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, flupromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypropazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone,

phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimino-dine, pimizole, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldi-amine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tolycaine, tranlycypromine, tra-zodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, tri-methoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J.Chromatogr.*, **1985**, *323*, 191–225.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, ami-triptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspi-rin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benz-phetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clen-buterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapson, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamor-phine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, dil-tiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, dox-apram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenopropfen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-stilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, loraze-pam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medaze-pam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, me-phesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methylglopamine, methylphenidate, methylprednisolone, methyltestosterone, methylpyrrol, me-toprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, na-proxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepi-nephine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbi-tal, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenyl-

butazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopalamine, scopolin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfathiazole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamine, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, **1994**, *18*, 233–242.

Methaqualone


Molecular formula: C₁₆H₁₄N₂O

Molecular weight: 250.30

CAS Registry No.: 72-44-6, 340-56-7 (HCl)

Merck Index: 6028

Lednicer No.: 1 353



SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 226

CHROMATOGRAM

Retention time: 4.19

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine;

prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; ace-nocoumarol; vandesine; mexiletine; dipyrindamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; fle-cainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; bupren-orphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazo-lam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclam-ide; chlorphenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; metha-done; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidra-mine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvox-amine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpi-pramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254-262.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES


Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, ami-triptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspi-rin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benz-phetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clen-buterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamor-phine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, dil-tiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, dox-apram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenoprofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiacol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-stilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, loraze-pam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medaze-pam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, me-

phesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methypylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nyldrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phenacyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidone, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopalamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfathiazole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripelethnamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233–242.

Methazolamide



Molecular formula: $C_5H_8N_4O_3S_2$

Molecular weight: 236.28

CAS Registry No.: 554-57-4

Merck Index: 6031

Lednicer No.: 1 250

SAMPLE

Matrix: blood, urine

Sample preparation: To 500 μ L whole blood, plasma, or urine add 20 μ L 1 mg/mL acetazolamide solution and 2.5 mL 50% ammonium sulfamate, vortex for 30 s. (Place the tube containing whole blood in boiling water for 30 s and then quickly in cold water.) Add 5 mL ethyl acetate, vortex, centrifuge at 3000 g for 10 min, transfer the organic layer to 5 mL pH 8.0 phosphate buffer, vortex, centrifuge at 3000 g for 10 min, transfer the organic layer to 500 μ L pH 10.0 glycine buffer, vortex for 30 s, centrifuge at 3000 g for 5 min. Aspirate and discard the organic layer, add 500 μ L ether to the remaining glycine buffer layer, vortex for 1 min, discard the ether phase. Vent the aqueous layer for about 30 min and inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Altima C18 (Alltech)

Mobile phase: MeCN:100 mM pH 4.0 sodium acetate 20:80

Flow rate: 1

Injection volume: 20

Detector: UV 285

CHROMATOGRAM

Retention time: 7.18

Internal standard: acetazolamide (4.55)

KEY WORDS

plasma; whole blood