

**Mobile phase:** MeOH:buffer 85:15 (Buffer was 90.7 mL 66.7 mM Na<sub>2</sub>HPO<sub>4</sub> and 9.3 mL 66.7 mM KH<sub>2</sub>PO<sub>4</sub> made up to 1 L with water, pH 7.6.)

**Flow rate:** 5 (sic)

**Injection volume:** 20

**Detector:** UV (wavelength not given)

---

#### CHROMATOGRAM

**Retention time:** 4.32

**Limit of detection:** 100 nM

---

#### OTHER SUBSTANCES

**Simultaneous:** chlordiazepoxide, diazepam, flurazepam

---

#### KEY WORDS

comparison with capillary electrophoresis; capillary GC; and polarography

---

#### REFERENCE

McGrath,G.; McClean,S.; O'Kane,E.; Smyth,W.F.; Tagliaro,F. Study of the capillary zone electrophoretic behaviour of selected drugs, and its comparison with other analytical techniques for their formulation assay, *J.Chromatogr.A*, **1996**, 735, 237–247.

---


## Nitrendipine

**Molecular formula:** C<sub>18</sub>H<sub>20</sub>N<sub>2</sub>O<sub>6</sub>

**Molecular weight:** 360.37

**CAS Registry No.:** 39562-70-4

**Merck Index:** 6669



---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** 2 mL Serum + 2 mL 20 ng/mL nimodipine in water, mix, add 1 mL 1 M NaOH, add 12 mL hexane:ethyl ether 50:50, extract, centrifuge at 300 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 250 µL mobile phase, add 250 µL water, inject a 50 µL aliquot.

---

#### HPLC VARIABLES

**Column:** 150 × 4.6 3 µm Spherisorb ODS

**Mobile phase:** MeOH:THF:water 47:15:38

**Flow rate:** 0.5

**Injection volume:** 50

**Detector:** UV 238

---

#### CHROMATOGRAM

**Retention time:** 11.87

**Internal standard:** nimodipine (13.08)

**Limit of detection:** 1 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** metabolites

---

#### KEY WORDS

serum

---

#### REFERENCE

Janis,R.A.; Krol,G.J.; Noe,A.J.; Pan,M. Radioreceptor and high-performance liquid chromatographic assays for the calcium channel antagonist nitrendipine in serum, *J.Clin.Pharmacol.*, **1983**, 23, 266–273.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 100  $\mu$ L Plasma + 1 mL ethyl acetate, shake horizontally for 3 min, centrifuge at 10000 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40-50°, reconstitute the residue in 200  $\mu$ L mobile phase, inject a 100  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Guard column:** 20  $\times$  2 30-40  $\mu$ m Perisorb RP-18 (Upchurch)

**Column:** 150  $\times$  4.6 5  $\mu$ m Hypersil ODS

**Mobile phase:** MeOH:water 65:35

**Flow rate:** 1

**Injection volume:** 100

**Detector:** UV 238

---

**CHROMATOGRAM**

**Retention time:** 7.5

**Internal standard:** nitrendipine

---

**OTHER SUBSTANCES**

**Extracted:** nimodipine

---

**KEY WORDS**

monkey; plasma; protect from light; nitrendipine is IS

---

**REFERENCE**

Qian,M.; Gallo,J.M. High-performance liquid chromatographic determination of the calcium channel blocker nimodipine in monkey plasma, *J.Chromatogr.*, **1992**, 578, 316-320.

---

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100  $\mu$ L mobile phase, centrifuge at 2800 g for 5 min, inject a 50  $\mu$ L aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

---

**HPLC VARIABLES**

**Column:** 300  $\times$  3.9 4  $\mu$ m NovaPack C18

**Mobile phase:** MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic))  $\text{KH}_2\text{PO}_4$  adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

**Column temperature:** 30

**Flow rate:** 0.8

**Injection volume:** 50

**Detector:** UV 236

---

**CHROMATOGRAM**

**Retention time:** 6.00

**Limit of detection:** <120 ng/mL

---

**KEY WORDS**

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazo-

lam; prazosin; flunitrazepam; clonazepam; metoclopramide; melfhalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vandesine; mexiletine; dipyrindamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cycizine; histapyrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

## REFERENCE

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254-262.

## SAMPLE

**Matrix:** blood, gastric contents, tissue, urine

**Sample preparation:** 1 mL Blood, urine, or gastric contents or 1 g tissue homogenate + 500  $\mu$ L buffer + 8 mL n-hexane:ethyl acetate 70:30, mix on a rotary mixer for 10 min, centrifuge at 3000 g for 8 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 100  $\mu$ L 12.5 mM NaOH in MeOH:water 50:50, inject a 50  $\mu$ L aliquot. (Buffer was 13.8 g potassium carbonate in 100 mL water, pH adjusted to 9.5 with concentrated HCl.)

## HPLC VARIABLES

**Guard column:** 4  $\times$  4 30  $\mu$ m LiChrocart Aluspher RP-select B (Merck)

**Column:** 125  $\times$  4 5  $\mu$ m Aluspher RP-select B (Merck)

**Mobile phase:** Gradient. A was 12.5 mM NaOH in MeOH. B was 12.5 mM NaOH in water. A:B 10:90 for 5 min, to 90:10 over 15 min, maintain at 90:10 for 5 min, return to initial conditions over 1 min, re-equilibrate for 5 min.

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 230, 254

## CHROMATOGRAM

**Retention time:** 17

## OTHER SUBSTANCES

**Extracted:** alprenolol, amitriptyline, bromazepam, carbamazepine, chlordiazepoxide, chlorpromazine, clonazepam, desipramine, diazepam, flunitrazepam, haloperidol, nordiazepam, nortriptyline, pindolol, zolpidem

**Also analyzed:** acebutolol, acetaminophen, alprazolam, amphetamine, atenolol, betaxolol, brotizolam, caffeine, camazepam, captopril, chloroquine, clobazam, clomipramine, clothiapine, clotiazepam, cloxazolam, cocaine, codeine, diclofenac, dihydralazine, dihydrocodeine, dihydroergotamine, diphenhydramine, domperidone, doxepin, droperidol, ergotamine, ethyl loflazepate, fenethylamine, fluoxetine, flupentixol, flurazepam, furosemide, glyclazide, hydrochlorothiazide,

hydroxyzine, ibuprofen, imipramine, ketazolam, loprazolam, lorazepam, lormetazepam, maprotiline, medazepam, mepyramine, methadone, methaqualone, methyl dopa, methylphenidate, metoclopramide, metoprolol, mexiletine, mianserin, midazolam, minoxidil, morphine, nadolol, nitrazepam, oxprenolol, papaverine, pentazocine, phenprocoumon, phenylbutazone, pipamperone, piritramide, practolol, prazepam, prazosin, promazine, promethazine, propoxyphene, propranolol, prothipendyl, quinine, sotalol, sulpride, thioridazine, trazodone, triazolam, trimipramine, tripeleminamine, tyramine, verapamil, yohimbine

---

## REFERENCE

Lambert, W.E.; Meyer, E.; De Leenheer, A.P. Systematic toxicological analysis of basic drugs by gradient elution of an alumina-based HPLC packing material under alkaline conditions, *J. Anal. Toxicol.*, **1995**, *19*, 73-78.

---

## SAMPLE

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

## HPLC VARIABLES

**Guard column:** 20 mm long Symmetry C18

**Column:** 250 × 4.6 5 µm Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 236.9

---

## CHROMATOGRAM

**Retention time:** 22.087

---

## KEY WORDS

whole blood

---

## REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, *763*, 149-163.

---

## SAMPLE

**Matrix:** cytosol incubations

**Sample preparation:** 2 mL Incubation + 5 mL chlorobutane:1,2-dichloroethane 80:20, shake for 15 min, centrifuge at 4000 g for 10 min. Remove the upper organic layer and evaporate it to dryness under a stream of nitrogen at 50°, reconstitute the residue in 400 µL mobile phase, inject an 80 µL aliquot.

---

## HPLC VARIABLES

**Column:** 300 × 3.9 µm Bondapak C18

**Mobile phase:** MeOH:10 mM ammonium phosphate buffer 66:34, pH 5.8

**Flow rate:** 1.5

**Injection volume:** 80

**Detector:** UV 230

---

## OTHER SUBSTANCES

**Extracted:** diphenyl sulfoxide

---

**KEY WORDS**

rat; rabbit; nitrendipine is IS

---

**REFERENCE**

Lee,S.C.; Renwick,A.G. Sulphoxide reduction by rat and rabbit tissues *in vitro*, *Biochem.Pharmacol.*, **1995**, 49, 1557–1565.

---

**SAMPLE**

**Matrix:** microsomal incubations

**Sample preparation:** Condition a 100 mg Bond Elut C18 SPE cartridge with 2 mL MeOH and 1 mL 1% aqueous formic acid. Add the microsomal incubation to the SPE cartridge, wash with 1 mL 1% aqueous formic acid, elute with 1 mL MeCN. Evaporate the eluate to dryness under a stream of nitrogen in the dark at room temperature, reconstitute the residue in 1 mL MeCN: 0.5% phosphoric acid 10:90, inject a 50  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4 10  $\mu$ m Nucleosil RP C8

**Mobile phase:** Gradient. MeCN:0.5% phosphoric acid from 40:60 to 60:40 over 12 min.

**Flow rate:** 1.4

**Injection volume:** 50

**Detector:** UV 234, UV 345

---

**CHROMATOGRAM**

**Retention time:** 8

---

**OTHER SUBSTANCES**

**Extracted:** metabolites

---

**KEY WORDS**

rat; liver; SPE

---

**REFERENCE**

Böcker,R.H.; Preuss,E.; Peter,R. High-performance liquid chromatography of the metabolites of nitrendipine and investigation into the metabolic pathways of this dihydropyridine, *J.Chromatogr.*, **1990**, 530, 206–211.

---

**SAMPLE**

**Matrix:** microsomal incubations, perfusate

**Sample preparation:** Extract using a 100 mg 1 mL Bond Elut C2 SPE cartridge, elute with MeOH, evaporate eluate to dryness under a stream of nitrogen, reconstitute the residue in 100  $\mu$ L mobile phase, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  5 Hichrom HiRPB deactivated reverse-phase

**Mobile phase:** MeOH:25 mM pH 5.0 N,N,N',N'-tetramethylethylenediamine phosphate buffer 75:25

**Flow rate:** 1

**Detector:** UV 245

---

**CHROMATOGRAM**

**Internal standard:** felodipine

**Limit of detection:** 5 ng/mL

---

**KEY WORDS**

rat; liver; SPE

---

**REFERENCE**

Walker,D.K.; Humphrey,M.J.; Smith,D.A. Importance of metabolic stability and hepatic distribution to the pharmacokinetic profile of amlodipine, *Xenobiotica*, **1994**, 24, 243–250.

---

**SAMPLE**

**Matrix:** perfusate

---

**HPLC VARIABLES**

**Column:** 100 × 8 5 μm Novapak C18 radial compression

**Mobile phase:** MeCN:10 mM pH 4.5 phosphate buffer 70:30

**Flow rate:** 2

**Detector:** UV 237

---

**OTHER SUBSTANCES**

**Also analyzed:** felodipine, nicardipine, nifedipine, nimodipine

---

**REFERENCE**

Diez,I.; Colom,H.; Moreno,J.; Obach,R.; Peraire,C.; Domenech,J. A comparative in vitro study of transdermal absorption of a series of calcium channel antagonists, *J.Pharm.Sci.*, **1991**, *80*, 931–934.

---

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a solution in mobile phase, inject a 20 μL aliquot.

---

**HPLC VARIABLES**

**Guard column:** present but not specified

**Column:** 75 × 4.6 3 μm XL octyl (Beckman)

**Mobile phase:** MeCN:water 44:56 containing 5 mM heptanesulfonic acid

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 6.2

---

**REFERENCE**

Greiner,P.O.; Angignard,D.; Cahn,J. High performance liquid chromatography of a new 1,4-dihydropyridine: applications to pharmacokinetic study in dogs, *J.Pharm.Sci.*, **1988**, *77*, 387–389.

---

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Inject a 10 μL aliquot of a solution in MeCN.

---

**HPLC VARIABLES**

**Column:** 100 × 4.6 CS-MP Spheri-5 cyano

**Mobile phase:** Gradient. MeCN:buffer from 10:90 to 40:60 over 10 min, re-equilibrate for 5 min. (Buffer was 50 mM KH<sub>2</sub>PO<sub>4</sub> adjusted to pH 3 with phosphoric acid.)

**Flow rate:** 1.5

**Injection volume:** 10

**Detector:** UV 200, 272, 276, 280, 314

---

**CHROMATOGRAM**

**Retention time:** 11

---

**OTHER SUBSTANCES**

**Simultaneous:** nifedipine

---

**REFERENCE**

Logan,B.K.; Patrick,K.S. Photodegradation of nifedipine relative to nitrendipine evaluated by liquid and gas chromatography, *J.Chromatogr.*, **1990**, *529*, 175–181.

---

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Centrifuge at 2000 g at 37° for 15 min.

---

**HPLC VARIABLES**

**Column:** 100 × 8 5 µm C18 Novapak

**Mobile phase:** MeCN:10 mM pH 4.5 phosphate buffer 70:30

**Flow rate:** 2

**Detector:** UV 237

---

**OTHER SUBSTANCES**

**Also analyzed:** nicardipine, nifedipine, nimodipine, felodipine

---

**KEY WORDS**

buffers

---

**REFERENCE**

Diez,I.; Colom,H.; Moreno,J.; Obach,R.; Peraire,C.; Domenech,J. A comparative in vitro study of transdermal absorption of a series of calcium channel antagonists, *J.Pharm.Sci.*, **1991**, 80, 931–934.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a solution in mobile phase, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 Nucleosil 5C18

**Mobile phase:** MeCN:water 60:40 containing 5.25 mL/L 1% phosphoric acid, 2.5 g/L ammonium phosphate, and 10 mL/L 10% tetra-n-butylammonium hydroxide

**Detector:** UV 254

---

**OTHER SUBSTANCES**

**Simultaneous:** nilvadipine

---

**KEY WORDS**

protect from light

---

**REFERENCE**

Kobayashi,D.; Matsuzawa,T.; Sugibayashi,K.; Morimoto,Y.; Kobayashi,M.; Kimura,M. Feasibility of use of several cardiovascular agents in transdermal therapeutic systems with *l*-menthol-ethanol system on hairless rat and human skin, *Biol.Pharm.Bull.*, **1993**, 16, 254–258.

---

# Nitrofurantoin


**Molecular formula:** C<sub>8</sub>H<sub>6</sub>N<sub>4</sub>O<sub>6</sub>

**Molecular weight:** 238.16

**CAS Registry No.:** 67-20-9

**Merck Index:** 6696

**Lednicer No.:** 1 230



---

**SAMPLE**

**Matrix:** cell cultures

**Sample preparation:** Centrifuge cell culture at 4000 g for 30 min, filter (0.2 µm Acrodisk), inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 5 µm Prodigy ODS-3

**Mobile phase:** Gradient. A was MeCN containing 2 g/L ammonium acetate. B was water. A:B 5:95 for 2 min, to 80:20 over 30 min, 80:20 for 5 min

**Flow rate:** 1  
**Injection volume:** 20  
**Detector:** UV 376

---

**CHROMATOGRAM**

**Retention time:** 24.5

---

**REFERENCE**

Rafii,F.; Hansen,E.B.,Jr. Isolation of nitrofurantoin-resistant mutants of nitroreductase-producing *Clostridium* sp. Strains from the human intestinal tract, *Antimicrob.Agents Chemother.*, **1998**, 42, 1121–1126.

---

**SAMPLE**

**Matrix:** milk

**Sample preparation:** Mix 50 mL cow milk with 25 mL 20% trichloroacetic acid, let stand for 15 min. Filter the samples and wash with water. Adjust the pH to 4.5-5 with NaOH, make up to 100 mL with water. Take a 25 mL aliquot and add it to a Sep-Pak Plus C18 SPE cartridge, elute with 2.5 mL mobile phase, pass nitrogen through eluate for at least 2 min (to remove oxygen), inject an aliquot.

---

**HPLC VARIABLES**

**Guard column:** Symmetry C18

**Column:** 150 × 3.9 4 μm Nova Pak C18

**Mobile phase:** MeCN:100mM aqueous sodium perchlorate:glacial acetic acid 28:72:0.5

**Flow rate:** 1

**Injection volume:** 20

**Detector:** E, ESA Coulochem II, Model 5011 analytical cell, porous carbon electrode -600 V, Model 5021 conditioning cell

---

**CHROMATOGRAM**

**Retention time:** 2.2

**Limit of detection:** 4 ppb

---

**OTHER SUBSTANCES**

**Extracted:** furaltadone, furazolidone

---

**KEY WORDS**

cow; SPE


---

**REFERENCE**

Galeano Diaz,T.; Guiberteau Cabanillas,A.; Acedo Valenzuela,M.I.; Correa,C.A.; Salinas,F. Determination of nitrofurantoin, furazolidone and furaltadone in milk by high-performance liquid chromatography with electrochemical detection, *J.Chromatogr.A*, **1997**, 764, 243–248.

---

# Nitrofurazone



**Molecular formula:** C<sub>6</sub>H<sub>6</sub>N<sub>4</sub>O<sub>4</sub>

**Molecular weight:** 198.14

**CAS Registry No.:** 59-87-0

**Merck Index:** 6697

**Lednicer No.:** 1 229

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject



a 10 (urine) or 30 (blood)  $\mu\text{L}$  aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18

**Column:** 250  $\times$  4.6 5  $\mu\text{m}$  Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 260.6

---

**CHROMATOGRAM**

**Retention time:** 10.323

---

**KEY WORDS**

whole blood

---

**REFERENCE**

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

---

**SAMPLE**

**Matrix:** eggs, tissue

**Sample preparation:** Add 30 mL MeCN to 10 g homogenized shelled eggs, liver, or muscle, blend at low speed for 2 min, centrifuge at 1000 g for 5 min, add 10 mL 10% NaCl solution and 50 mL dichloromethane to the supernatant, shake for a few min. Filter the lower organic layer through 5 g anhydrous sodium sulfate, evaporate the filtrate to dryness using a rotary vacuum evaporator at 45°, redissolve the residue in 1 mL MeCN:MeOH:20 mM pH 4.6 sodium acetate 10:50:40, inject an aliquot. (Protect from light. Wash the 1 mL of MeCN:MeOH:20 mM pH 4.6 sodium acetate 10:50:40 three times with 1 mL n-hexane before use.)

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu\text{m}$  Supelcosil L C18-DB

**Mobile phase:** MeCN:water 50:50 containing 1 mM ammonium acetate and 0.025% acetic acid

**Flow rate:** 0.6

**Injection volume:** 20

**Detector:** MS, PESCIEX API I, ionspray interface 5500 V, OR 60 V, m/z 199, split the column effluent so that 0.03 mL/min enters the MS

---

**CHROMATOGRAM**

**Retention time:** 5.2

**Limit of detection:** 3.2 ng/g

---

**KEY WORDS**

chicken; liver; muscle

---

**REFERENCE**

Draisci,R.; Giannetti,L.; Lucentini,L.; Palleschi,L.; Brambilla,G.; Serpe,L.; Gallo,P. Determination of nitrofurans residues in avian eggs by liquid chromatography-UV photodiode array detection and confirmation by liquid chromatography-ionspray mass spectrometry, *J.Chromatogr.A*, **1997**, 777, 201-211.

---

**SAMPLE**

**Matrix:** eggs, tissue

**Sample preparation:** Add 30 mL MeCN to 10 g homogenized shelled eggs, liver, or muscle, blend at low speed for 2 min, centrifuge at 1000 g for 5 min, add 10 mL 10% NaCl solution and 50 mL dichloromethane to the supernatant, shake for a few min. Filter the lower organic layer through 5 g anhydrous sodium sulfate, evaporate the filtrate to dryness using a rotary vacuum evaporator at 45°, redissolve the residue in 1 mL MeCN:MeOH:20 mM pH 4.6 sodium acetate 10:50:40, inject an aliquot. (Protect from light. Wash the 1 mL of MeCN:MeOH:20 mM pH 4.6 sodium acetate 10:50:40 three times with 1 mL n-hexane before use.)

---

**HPLC VARIABLES**

**Guard column:** 10 × 4.6 µm Bondapak C18

**Column:** 150 × 4.6 5 µm Spherisorb ODS2 S5

**Mobile phase:** MeCN:20 mM pH 4.6 sodium acetate 21:79

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 362

---

**CHROMATOGRAM**

**Retention time:** 5.3

**Limit of detection:** 2.5 ng/g

---

**OTHER SUBSTANCES**

**Extracted:** furazolidone, furaltadone

---

**KEY WORDS**

chicken; liver; muscle

---

**REFERENCE**

Draisci,R.; Giannetti,L.; Lucentini,L.; Palleschi,L.; Brambilla,G.; Serpe,L.; Gallo,P. Determination of nitrofurans residues in avian eggs by liquid chromatography-UV photodiode array detection and confirmation by liquid chromatography-ion spray mass spectrometry, *J.Chromatogr.A*, **1997**, 777, 201-211.

---

**SAMPLE**

**Matrix:** eggs, tissue

**Sample preparation:** Blend (Stomacher) 10 g homogenized tissue with 30 mL saline for 3 min, centrifuge at 2000 g, mix 20 mL of the supernatant with 2 mL 1% sodium azide. Dilute 10 mL homogenized egg with 10 mL saline, add 3 mL 10% sodium azide solution. Dialyze sample using a Cuprophane membrane (10000-15000 dalton cut-off) against water pumped at 0.36-1.44 mL/min for 3-9 min, pass the water through column A, flush the column with pure water for 8 min, backflush the contents of column A onto column B with mobile phase, after 5 min remove column A from the circuit, elute column B with mobile phase, monitor the effluent from column B. To increase sensitivity a number of sample batches can be dialyzed before the contents of column A are analyzed. (Caution! Sodium azide is carcinogenic, mutagenic, and highly toxic! Do not discharge to the sink!)

---

**HPLC VARIABLES**

**Column:** A 60 × 4.6 37-50 µm Bondapak C18/Corasil; B 250 × 4.6 5 µm Hypersil ODS

**Mobile phase:** MeCN:100 mM pH 5 acetate buffer 20:80

**Flow rate:** 1

**Detector:** UV 365

---

**CHROMATOGRAM**

**Retention time:** 6

**Limit of detection:** 2 ng/g (tissue), 1 ng/g (eggs)

---

**OTHER SUBSTANCES**

**Extracted:** furaltadone, furazolidone, nitrofurantoin

---

**KEY WORDS**

protect from light; cow; muscle; dialysis; column-switching

---

**REFERENCE**

Aerts, M.M.; Beek, W.M.; Brinkman, U.A. On-line combination of dialysis and column-switching liquid chromatography as a fully automated sample preparation technique for biological samples. Determination of nitrofurans residues in edible products, *J. Chromatogr.*, **1990**, *500*, 453–468.

---

**SAMPLE**

**Matrix:** feed

**Sample preparation:** Grind feed to pass 20 mesh. 10 g Feed + 5 mL water, swirl, let stand for 5 min, add 50 mL DMF:water 95:5, shake vigorously for 15 s, let stand in the dark at room temperature overnight, filter (paper). Add 15 mL of the filtrate to 5 g alumina (Alcoa F-20, 80–200 mesh) in a 300 × 10 glass column, discard first several mL of eluate, collect remaining eluate, inject an aliquot

---

**HPLC VARIABLES**

**Guard column:** 100 × 2 µBondapak C18/Corasil

**Column:** 300 × 4 µBondapak C18

**Mobile phase:** MeCN:1% acetic acid 20:80

**Detector:** UV 280, UV 365

---

**CHROMATOGRAM**

**Retention time:** 5

---

**OTHER SUBSTANCES**

**Extracted:** furazolidone, carbadox

---

**KEY WORDS**

protect from light

---

**REFERENCE**

Thorpe, V.A. Sample preparation of carbadox, furazolidone, nitrofurazone, and ethopabate in medicated feeds for high pressure liquid chromatography, *J. Assoc. Off. Anal. Chem.*, **1980**, *63*, 981–984.

---

**SAMPLE**

**Matrix:** food

**Sample preparation:** Condition a 6 mL Bond Elut C18 SPE cartridge with 5 mL MeOH and 5 mL water. Homogenize (Tissuemizer) 5 g finely-chopped (by hand) shrimp and 20 mL MeCN at medium speed for 45 s, centrifuge at 3000 rpm for 5 min, decant the supernatant. Add 30 mL hexane saturated with MeCN to the supernatant, shake for 30 s, discard the hexane layer. Add 10 mL EtOH to the MeCN layer, evaporate under reduced pressure at 45° to 2–5 mL (until liquid looks milky), add 2 mL EtOH, continue evaporation until there is 2 mL of a thick liquid, add 2 mL EtOH, evaporate to dryness. Add 2 mL water to the residue, sonicate for 5 min, add to the SPE cartridge, wash with 4 mL water, elute at ≤3 mL/min with 5 mL MeCN. Evaporate the eluate to dryness under a stream of nitrogen at ≤45° (remove promptly when dry), reconstitute the residue in 1 mL mobile phase, filter (0.45 µm), inject a 50 µL aliquot.

---

**HPLC VARIABLES**

**Guard column:** 20 × 4 5 µm ODS Hypersil C18

**Column:** 200 × 4.6 5 µm ODS Hypersil C18

**Mobile phase:** MeCN:1% aqueous acetic acid 25:75

**Column temperature:** 40

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 375

---

**CHROMATOGRAM**

**Retention time:** 3.7

**Limit of quantitation:** 4 ng/g

---

**OTHER SUBSTANCES**

**Extracted:** furazolidone

---

**KEY WORDS**

SPE; shrimp

---

**REFERENCE**

Rupp,H.S.; Munns,R.K.; Long,A.R. Simultaneous determination of nitrofurazone and furazolidone in shrimp (*Penaeus vannamei*) muscle tissue by liquid chromatography with UV detection, *JAOAC Int.*, **1993**, 76, 1235-1239.

---

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a solution in MeOH:water 30:70, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 150  $\times$  3.9 5  $\mu$ m Resolve spherical C18 (Waters)

**Mobile phase:** MeOH:water 35:65

**Column temperature:** 40

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 305

---

**CHROMATOGRAM**

**Retention time:** 3

---

**OTHER SUBSTANCES**

**Simultaneous:** carbadox, furazolidone

**Noninterfering:** pyrantel

---

**KEY WORDS**

protect from light

---

**REFERENCE**

Roybal,J.E.; Munns,R.K.; Shimoda,W. Liquid chromatographic determination of carbadox residues in animal feed, *J.Assoc. Off. Anal. Chem.*, **1985**, 68, 653-657.

---

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Dissolve in chloroform at a concentration of 1  $\mu$ g/mL, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4 5  $\mu$ m Lichrospher RP-18

**Mobile phase:** MeCN:10 mM sodium acetate 20:80, pH 5

**Column temperature:** 30

**Flow rate:** 1.6

**Injection volume:** 20

**Detector:** UV 365

---

**CHROMATOGRAM**

**Retention time:** 7.5

**Limit of detection:** 15 ng/mL

---

**OTHER SUBSTANCES**

**Simultaneous:** degradation products, carbadox, nitrofurantoin, furazolidone, furaltadone

---

**REFERENCE**

Kaniou,I.; Zachariadis,G.; Kalligas,G.; Tsoukali,H.; Stratis,J. Separation and determination of carbadox, nitrofurazone, nitrofurantoin, furazolidone, and furaltadone in their mixtures by thin layer and high performance liquid chromatography, *J.Liq.Chromatogr.*, **1994**, 17, 1385-1398.

---

---

**SAMPLE**

**Matrix:** tissue

**Sample preparation:** Homogenize (Waring blender) 10 g muscle, liver, or kidney in 100 mL EtOH for 5 min, let stand for 5 min, filter through 10 g Celite 545 on top of a sintered glass filter, rinse blender with 100 mL EtOH and filter rinse. Add 25 mL 3.6% aqueous metaphosphoric acid to the combined filtrates, evaporate to 25 mL under reduced pressure at 45°. Remove residue, rinse out flask with 5 mL hexane and 3 mL water, combine, centrifuge at 0° at 27000 g for 30 min, discard hexane, rinse surface with 5 mL hexane, discard hexane. Remove aqueous layer, rinse out tube twice with 3 mL portions of water, combine, add 10 mL 1 M  $\text{KH}_2\text{PO}_4$ , make up to 100 mL with water, extract three times for 5 min with 50 mL ethyl acetate. Combine the extracts and dry them over 15 g anhydrous sodium sulfate, filter through glass wool, evaporate to dryness under reduced pressure at 45°. Take up residue in 3 mL ethyl acetate and add to alumina column, rinse flask with 2 mL ethyl acetate and add rinse to column. Elute with 20 mL EtOH:MeOH:ethyl acetate 10:10:80 and combine all the eluate. Evaporate to dryness under reduced pressure at 45°, reconstitute in 500  $\mu\text{L}$  mobile phase, inject a 100  $\mu\text{L}$  aliquot. (Prepare alumina column by slurrying 1 g aluminum oxide (Baker) in 20 mL ethyl acetate and adding to a 200  $\times$  6 glass chromatographic column.)

---

#### HPLC VARIABLES

**Guard column:** Brownlee 10  $\mu\text{m}$  RP-GU MPLC C-8

**Column:** 250  $\times$  4.6 Brownlee RP-10A C-8

**Mobile phase:** MeCN:EtOH:10 mM ammonium acetate 25:5:70, pH 6.8

**Flow rate:** 1

**Injection volume:** 100

**Detector:** UV 350

---

#### CHROMATOGRAM

**Retention time:** 6.6

**Limit of detection:** 2 ng

**Limit of quantitation:** 10 ng

---

#### OTHER SUBSTANCES

**Extracted:** quinoxaline-2-carboxylic acid, furazolidone, carbadox, desoxycarbadox

---

#### KEY WORDS

protect from light; pig; muscle; liver; kidney

---

#### REFERENCE

MacIntosh, A.I.; Neville, G.A. Liquid chromatographic determination of carbadox, desoxycarbadox, and nitrofurazones in pork tissues, *J. Assoc. Off. Anal. Chem.*, **1984**, 67, 958-962.

---

#### SAMPLE

**Matrix:** tissue

**Sample preparation:** Homogenize (Ultra-Turrax TP 18/2) 3 g ground muscle + 6.8 mL MeCN for 6 s, centrifuge at 5000 rpm for 5 min. Remove 6.5 mL of the supernatant and add it to 2 mL 5 M NaCl, shake vigorously for 10 s, centrifuge at 3000 rpm for 2 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 43°, reconstitute the residue in 250  $\mu\text{L}$  MeCN:buffer 20:80, add 1 mL hexane, mix (Whirlimixer), centrifuge for 4 min, discard the hexane layer, filter (Costar Spin-X 0.22  $\mu\text{m}$  cellulose acetate) while centrifuging at 5600 g for 4 min, inject a 20  $\mu\text{L}$  aliquot of the filtrate. (Buffer was 20 mM sodium 1-heptanesulfonate and 10 mM  $\text{Na}_2\text{HPO}_4$ , pH adjusted to 6.0 with phosphoric acid.)

---

#### HPLC VARIABLES

**Guard column:** 20  $\times$  4.6 5  $\mu\text{m}$  Supelcosil LC-ABZ

**Column:** 250  $\times$  4.6 5  $\mu\text{m}$  Supelcosil LC-ABZ

**Mobile phase:** MeCN:buffer 25:75 (Buffer was 4.45 g sodium 1-heptanesulfonate and 9.5 g  $\text{Na}_3\text{PO}_4 \cdot 12\text{H}_2\text{O}$  in 750 mL water, adjust pH to 2.5 with 5 M phosphoric acid, make up to 1 L with water.)

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 365

---

#### CHROMATOGRAM

**Retention time:** 5.5

**Internal standard:** nitrofurazone

## OTHER SUBSTANCES

**Extracted:** furazolidone


## KEY WORDS

cow; muscle; nitrofurazone is IS

## REFERENCE

Hormazábal,V.; Yndestad,M. Simple and rapid method of analysis for furazolidone in meat tissues by high-performance liquid chromatography, *J.Liq.Chromatogr.*, **1995**, *18*, 1871–1877.

# Nizatidine



**Molecular formula:** C<sub>12</sub>H<sub>21</sub>N<sub>5</sub>O<sub>2</sub>S<sub>2</sub>

**Molecular weight:** 331.46

**CAS Registry No.:** 76963-41-2

**Merck Index:** 6758

**Lednicer No.:** 4 95

## SAMPLE

**Matrix:** blood, tissue

**Sample preparation:** Plasma. 100  $\mu$ L Plasma + 100  $\mu$ L 75  $\mu$ M IS + 100  $\mu$ L 5 M NaOH + 5 mL dichloromethane, shake for 10 min, centrifuge at 2000 rpm for 10 min. Evaporate 4 mL of the organic phase to dryness under a stream of nitrogen. Dissolve residue in 100  $\mu$ L mobile phase, inject a 50  $\mu$ L aliquot. Tissue. Homogenize 500 mg liver with saline on ice for 1 min. Add 100  $\mu$ L 75  $\mu$ M IS, 100  $\mu$ L 0.5 M NaOH, and 5 mL dichloromethane, shake for 10 min, centrifuge at 3000 rpm for 10 min. Evaporate 3 mL of the organic phase to dryness under a stream of nitrogen. Dissolve residue in 100  $\mu$ L mobile phase, pass through a Ministar-RC 15 cartridge (Sartorius, Germany), inject a 50  $\mu$ L aliquot.

## HPLC VARIABLES

**Guard column:** 30  $\times$  4.6 Senshu Pak ODS-1031 (Senshu Sciences, Japan)

**Column:** 250  $\times$  4.6 Senshu Pak ODS -1251 (Senshu Sciences, Japan)

**Mobile phase:** MeCN:water 5:95 containing 5 mM NaH<sub>2</sub>PO<sub>4</sub> and 5 mM tetramethylammonium chloride

**Flow rate:** 1.5

**Injection volume:** 50

**Detector:** UV 228

## CHROMATOGRAM

**Internal standard:** cimetidine

**Limit of detection:** 50-100  $\mu$ g/mL (sic)

## KEY WORDS

plasma; rat; pharmacokinetics

## REFERENCE

Takedomi,S.; Matsuo,H.; Yamano,K.; Yamamoto,K.; Iga,T.; Sawada,Y. Quantitative prediction of the interaction of midazolam and histamine H<sub>2</sub> receptor antagonists in rats, *Drug Metab.Dispos.*, **1998**, *26*, 318–323.

## SAMPLE

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50  $\mu$ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject

a 10 (urine) or 30 (blood)  $\mu\text{L}$  aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18

**Column:**  $250 \times 4.6$  5  $\mu\text{m}$  Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 316.4

---

**CHROMATOGRAM**

**Retention time:** 3.302

---

**KEY WORDS**

whole blood

---

**REFERENCE**

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

---

# Nomifensine

**Molecular formula:**  $\text{C}_{16}\text{H}_{18}\text{N}_2$


**Molecular weight:** 238.33

**CAS Registry No.:** 24526-64-5, 32795-47-4 (maleate)

**Merck Index:** 6768

**Lednicer No.:** 4 146

---



---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol:n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100  $\mu\text{L}$  mobile phase, centrifuge at 2800 g for 5 min, inject a 50  $\mu\text{L}$  aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

---

**HPLC VARIABLES**

**Column:**  $300 \times 3.9$  4  $\mu\text{m}$  NovaPack C18

**Mobile phase:** MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic))  $\text{KH}_2\text{PO}_4$  adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

**Column temperature:** 30

**Flow rate:** 0.8

**Injection volume:** 50

**Detector:** UV 240

---

**CHROMATOGRAM**

**Retention time:** 4.87

**Limit of detection:** <120 ng/mL

---

**KEY WORDS**

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfapyrazole; flumazenil; sulpride; morphine; atenolol; tolaxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vandesine; mexiletine; dipyrindamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phencyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluperidol; madaazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thiopropazine; methadone; amoxapine; quinupramine; pipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

---

**REFERENCE**

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254–262.

---

**SAMPLE**

**Matrix:** bulk

**Sample preparation:** Mix 8 mg nomifensine maleate with 200  $\mu$ L trifluoroacetic anhydride, stir at room temperature for 10 min, evaporate to dryness under a stream of nitrogen at 40°, reconstitute the residue in 2 mL mobile phase, inject a 5  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  3.6 5  $\mu$ m Spherisorb silica

**Mobile phase:** Chloroform:MeOH:water 100:3:0.15 containing 2 mM (+)-camphor-10-sulfonic acid

**Flow rate:** 1

**Injection volume:** 5

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 16 (D), 18 (L)

**Limit of quantitation:** 25 ng

---

**KEY WORDS**

derivatization; chiral; normal phase



---

**REFERENCE**

Tsujiyama, T.; Tsuchiya, M.; Hamachi, Y.; Kuriki, T.; Fukunaga, T.; Suzuki, N. Ion-pair chromatographic separation of nomifensine maleate enantiomers, *Anal. Sci.*, **1989**, 5, 285–288.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 125 × 4.9 Spherisorb S5W silica

**Mobile phase:** MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

**Flow rate:** 2

**Injection volume:** 20

**Detector:** E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

---

**CHROMATOGRAM**

**Retention time:** 1.7

---

**OTHER SUBSTANCES**


**Also analyzed:** acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclicizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazepine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipranone, diprenorphine, dipyrizamide, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, flupromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserin, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimetopazine, methoxamine, methoxyphenazine, methoxypropazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindine, pimozide, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluorperazine, trifluoridol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

---

**REFERENCE**

Jane, I.; McKinnon, A.; Flanagan, R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, 323, 191–225.

# Nonoxynol-9



**Molecular formula:** C<sub>33</sub>H<sub>60</sub>O<sub>10</sub>

**Molecular weight:** 616.83

**CAS Registry No.:** 26027-38-3

**Merck Index:** 6772

## SAMPLE

**Matrix:** blood, urine, vaginal fluid

**Sample preparation:** Urine. Centrifuge at 4000 rpm for 5 min, inject a 20  $\mu$ L aliquot of the supernatant. Serum. Mix serum with an equal volume of THF, centrifuge at 4000 rpm for 5 min, remove an aliquot of the supernatant and mix it with an equal volume of THF, centrifuge at 4000 rpm for 5 min, inject a 20  $\mu$ L aliquot of the supernatant. Vaginal fluid. Vortex swab with 1 mL MeCN, let stand at room temperature for 15 min, weigh, vortex, remove swab, centrifuge at 4000 rpm for 5 min, inject a 20  $\mu$ L aliquot of the supernatant.

### HPLC VARIABLES

**Column:** 250 × 4.6 10 μm R-SIL-amine (Alltech)

**Mobile phase:** THF:MeCN 95:5

**Flow rate: 1**

**Injection volume: 20**

**Detector:** F ex 275 em 575

## CHROMATOGRAM

**Retention time: 4**

**Limit of detection:** 460 ng (vaginal fluid), 230 ng/mL (urine), 1.01 µg/mL (serum)

## KEY WORDS

place a 2  $\mu\text{m}$  filter at the column inlet

## REFERENCE

Beck, G.J.; Kossak, D.; Saxena, S.J. A simple, sensitive assay for the spermicide nonoxonyl-9 in biological fluids by high-performance liquid chromatography, *J.Pharm.Sci.*, **1990**, 79, 1029-1031.

## SAMPLE

**Matrix:** solutions

**Sample preparation:** Prepare a 0.2% solution, inject a 10  $\mu$ L aliquot.

## HPLC VARIABLES

**Column:** 300 × 75 PLgel 5 mixed-C (polymer Laboratories) in series with 300 × 78 Ultrastyrigel 100 Å styrene-divinylbenzene (Waters)

**Mobile phase:** THF

**Column temperature: 40**

**Flow rate: 1**

**Injection volume: 10**

**Detector:** RI

## KEY WORDS

SEC

## REFERENCE

Ysambertt, F.; Cabrera, W.; Marquez, N.; Salager, J.L. Analysis of ethoxylated nonylphenol surfactants by high-performance size-exclusion chromatography (HPSEC). *J. Liq. Chromatogr.* **1995**, *18*, 1157-1171.

## SAMPLE

**Matrix:** vaginal lavage fluid

**Sample preparation:** 500  $\mu$ L Vaginal lavage fluid (water) + 500  $\mu$ L 25  $\mu$ g/mL 4-octylphenol in mobile phase, vortex briefly, inject a 20  $\mu$ L aliquot.

**HPLC VARIABLES**

**Column:** 250 × 4.6 10 μm R-SIL-Amine (NH<sub>2</sub>) (Alltech) (Change 2 μm filter in front of column before each run.)

**Mobile phase:** MeCN:THF 2:98

**Flow rate:** 1

**Injection volume:** 20

**Detector:** F ex 227 em 612

**CHROMATOGRAM**

**Retention time:** 5.7

**Internal standard:** 4-octylphenol (4.0)


**Limit of detection:** 480 ng/mL

**Limit of quantitation:** 3.125 μg/mL

**REFERENCE**

McPherson, J.L.; Nichols, J.H.; Barditch-Crovo, P.; Hamzeh, F.M. Determination of the spermicide nonoxonyl-9 in vaginal lavage by high-performance liquid chromatography, *J. Chromatogr. B*, **1996**, 677, 204–208.

# Norepinephrine



**Molecular formula:** C<sub>8</sub>H<sub>11</sub>NO<sub>3</sub>

**Molecular weight:** 169.18

**CAS Registry No.:** 51-41-2, 69815-49-2 (bitartrate monohydrate), 51-40-1 (bitartrate)

**Merck Index:** 6788

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 20 mL Whole blood + 1 mL 20 mg/mL EDTA solution containing 10 mg/mL sodium metabisulfite, mix, centrifuge at 4° at 4000 g for 10 min. Remove the plasma and add concentrated perchloric acid until the concentration of perchloric acid is 400 mM, mix, let stand in the cold for 15 min, centrifuge at 4° at 20000 g for 20 min. Adjust pH of 2 mL supernatant to 7.0 ± 0.2 with 500 mM KOH, add 400 μL 875 μg/mL o-phthalaldehyde in pH 10.40 ± 0.02 buffer (containing 2-mercaptoethanol ?), add 2 g NaCl, add 2 mL ethyl acetate, shake for 1 min, centrifuge at 3400 g, repeat the extraction. Combine the organic layers and add them to 2 mL 35 mM pH 10.0 ± 0.1 Na<sub>2</sub>HPO<sub>4</sub> buffer, shake for 1 min, centrifuge at 3400 g, discard the aqueous layer, wash the ethyl acetate layer again with phosphate buffer. Reduce the ethyl acetate volume to 100 μL under a stream of nitrogen, inject a 10-50 μL aliquot.

**HPLC VARIABLES**

**Guard column:** Co:Pell ODS

**Column:** 300 × 4 10 μm μBondapak phenyl

**Mobile phase:** Gradient. MeCN:25 mM pH 5.10 NaH<sub>2</sub>PO<sub>4</sub> buffer 25:75 for 15 min then MeOH: 25 mM pH 5.10 NaH<sub>2</sub>PO<sub>4</sub> buffer 45:55 (step gradient).

**Column temperature:** 26

**Flow rate:** 1.5

**Injection volume:** 10-50

**Detector:** F ex 340 em 480

**CHROMATOGRAM**

**Retention time:** 11

**Internal standard:** tyramine (44)

**Limit of detection:** 0.5 ng/mL

**OTHER SUBSTANCES**

**Extracted:** dopamine, serotonin