

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 240.5

CHROMATOGRAM

Retention time: 24.038

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

Norethynodrel


Molecular formula: C₂₀H₂₆O₂

Molecular weight: 298.43

CAS Registry No.: 68-23-5

Merck Index: 6791

Lednicer No.: 1 168



SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 50 × 4.6 5 μm Supelcosil LC-18

Mobile phase: MeOH:THF:water 10:20:70

Flow rate: 2

Injection volume: 20

Detector: UV 220

CHROMATOGRAM

Retention time: 6.2 (norethynodrel acetate)

OTHER SUBSTANCES

Simultaneous: ethinyl estradiol, norethindrone, norethindrone acetate, norgestrel

REFERENCE

Supelco Catalog, **1994**, p. 779.

Norfloxacin


Molecular formula: C₁₆H₁₈FN₃O₃

Molecular weight: 319.34

CAS Registry No.: 70458-96-7, 100587-52-8 (norfloxacin succinil)

Merck Index: 6793

Lednicer No.: 4 141, 143



SAMPLE

Matrix: aqueous humor, blood, tissue

Sample preparation: Homogenize sample in 1 M HCl, add MeCN, centrifuge, add an equal volume of dichloromethane, centrifuge, inject a 50-100 μ L aliquot of the aqueous supernatant.

HPLC VARIABLES

Guard column: Corasil C18

Column: μ Bondapak C18

Mobile phase: MeCN:25 mM phosphoric acid adjusted to pH 3.0 with tetrabutylammonium hydroxide 9.5:90.5

Flow rate: 1.5

Injection volume: 50-100

Detector: F ex 272 em 370 (cut-off)

CHROMATOGRAM

Retention time: 6.5

Limit of detection: 0.1 ng/g

KEY WORDS

serum; rabbit; human; cornea

REFERENCE

Bron,A.M.; Péchinot,A.; Garcher,C.; Guyonnet,G.; Kazmierczak,A. Ocular penetration of topically applied norfloxacin 0.3% in the rabbits and in humans, *J.Ocul.Pharmacol.*, **1992**, 8, 241-246.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Serum + 50 μ L 400 μ g/mL IS in water, vortex for 30 s. Add 500 μ L MeCN, vortex for 1 min. Centrifuge at 6000 rpm for 10 min. Evaporate the supernatant to 200 μ L at 40° under a stream of nitrogen, vortex for 30 s. Inject a 30-80 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 8.0 4 μ m Radial-pak Novapak C18

Mobile phase: MeCN:buffer 14:86 (Buffer was 2 g citric acid, 2 g sodium acetate, and 1 mL triethylamine in 1 L water.)

Flow rate: 2.5

Injection volume: 30-80

Detector: F ex 330 em 440

CHROMATOGRAM

Retention time: 4.0

Internal standard: acebutolol (7.4)

Limit of quantitation: 50 ng/mL

OTHER SUBSTANCES

Extracted: pefloxacin

Simultaneous: ciprofloxacin, lomefloxacin, ofloxacin

KEY WORDS

serum; pharmacokinetics

REFERENCE

Abanmi,N.; Zaghloud,I.; El Sayed,N.; al-Khamis,K.I. Determination of pefloxacin and its main active metabolite in human serum by high-performance liquid chromatography, *Ther.Drug Monit.*, **1996**, 18, 158-163.

SAMPLE

Matrix: blood

Sample preparation: Add 20 μ L MeOH:0.1% trifluoroacetic acid 15:85 and 5 μ L (sic) MeCN to 300 μ L plasma. Centrifuge at 600 g for 10 min. Evaporate the supernatant under nitrogen at 40° for 30 min. Reconstitute the residue in 200 μ L MeOH:0.1% trifluoroacetic acid 15:85. Inject a 50 μ L aliquot.

HPLC VARIABLES**Guard column:** 12.5 × 4 Zorbax RX-C18**Column:** 150 × 4.6 5 µm Zorbax SB-C8**Mobile phase:** MeCN:water:trifluoroacetic acid 19:81:0.02**Flow rate:** 1**Injection volume:** 50**Detector:** UV 279

CHROMATOGRAM**Retention time:** 3.8-3.9**Internal standard:** norfloxacin

OTHER SUBSTANCES**Extracted:** ciprofloxacin, enrofloxacin

KEY WORDScat; plasma; norfloxacin is IS

REFERENCE

Kordick,D.L.; Papich,M.G.; Breitschwerdt,E.B. Efficacy of enrofloxacin or doxycycline for treatment of Bartonella henselae or Bartonella clarridgeiae infection in cats, *Antimicrob.Agents Chemother.*, **1997**, 41, 2448-2455.

SAMPLE**Matrix:** blood**Sample preparation:** Filter 1 mL plasma using a micropartition system (MPS-1, Amicon, MA) while centrifuging at 2000 g for 20 min at 10°, inject an aliquot of the ultrafiltrate.

HPLC VARIABLES**Column:** 250 × 4.6 Spherisorb ODS-2 endcapped**Mobile phase:** MeCN:buffer 11:89 containing 5 mM tetrabutylammonium sulfate, adjusted to pH 2.5 with 1 M NaOH (Buffer was 100 mM citric acid containing 200 mM ammonium perchlorate.)**Column temperature:** 37**Flow rate:** 1.2**Detector:** UV 272

CHROMATOGRAM**Retention time:** 8.07**Internal standard:** pipemidic acid (4.14)

OTHER SUBSTANCES**Simultaneous:** pefloxacin

KEY WORDSplasma; ultrafiltrate

REFERENCE

Zlotos,G.; Bucker,A.; Kinzig-Schippers,M.; Sorgel,F.; Holzgrave,U. Plasma protein binding of gyrase inhibitors, *J.Pharm.Sci.*, **1998**, 87, 215-220.

SAMPLE**Matrix:** blood**Sample preparation:** 50 µL Plasma + 1 mL 100 mM pH 7.0 K₂HPO₄ adjusted to pH 7.0 with 85% orthophosphoric acid + 100 µL 300 µg/mL nalidixic acid in water + 3 mL dichloromethane: isoamyl alcohol 9:1, shake vigorously for 10 min, centrifuge at 2270 g for 10 min. Remove 2 mL of the organic phase and evaporate it to dryness under a stream of nitrogen at 40°. Reconstitute residue in 100 µL MeOH:50 mM NaOH 2:1, vortex, inject a 10 µL aliquot.

HPLC VARIABLES**Column:** 150 × 4.6 5 µm Chemcosorb 5-ODS-H

Mobile phase: MeOH:5 mM sodium lauryl sulfate 2:1, adjusted to pH 2.5 with 85% phosphoric acid (Better separation obtained at pH 2.35, *J.Chromatogr.* 1990, 530, 186.)

Column temperature: 40

Flow rate: 0.6

Injection volume: 10

Detector: UV 300

CHROMATOGRAM

Retention time: 6.5

Internal standard: nalidixic acid (5.0)

OTHER SUBSTANCES

Extracted: fenbufen, felbinac

Interfering: ofloxacin, enoxacin

KEY WORDS

plasma; rat

REFERENCE

Katagiri,Y.; Naora,K.; Ichikawa,N.; Hayashibara,M.; Iwamoto,K. Simultaneous determination of ofloxacin, fenbufen and felbinac in rat plasma by high-performance liquid chromatography, *J.Chromatogr.*, **1988**, 431, 135-142.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Serum + 250 μ L 10% trichloroacetic acid, vortex for 10 s, centrifuge at >700 g for 10 min, inject a 20 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 100 \times 8 μ Bondapak C18 Radial-PAK

Mobile phase: MeOH:18 mM KH_2PO_4 containing 0.13 mM heptanesulfonic acid:concentrated phosphoric acid 30:70:0.1

Injection volume: 20

Detector: F ex 278 em 475

CHROMATOGRAM

Retention time: 5.6

KEY WORDS

serum

REFERENCE

Griggs,D.J.; Wise,R. A simple isocratic high-pressure liquid chromatographic assay of quinolones in serum, *J.Antimicrob.Chemother.*, **1989**, 24, 437-445.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Plasma + 100 μ L pH 7.4 phosphate buffer + 50 μ L 10 μ g/mL β -hydroxypropyltheophylline in pH 7.4 phosphate buffer + 5 mL chloroform:isopropanol 80:20, shake on a rotary mixer for 15 min, centrifuge at 800 g for 5 min. Evaporate organic layer under nitrogen at 45°, sonicate residue with 100 μ L mobile phase, inject 25 μ L aliquot.

HPLC VARIABLES

Guard column: 10 \times 4.9 Spherisorb ODS

Column: 250 \times 4.9 Sherisorb S5 ODS2

Mobile phase: MeCN:buffer 15:85 adjusted to pH 3.0 with 85% phosphoric acid immediately before use (Buffer was 4.54 g KH_2PO_4 + 5.94 g $\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$ + 1.49 g tetrabutylammonium hydrogen sulfate per L.)

Flow rate: 1.3

Injection volume: 25

Detector: UV 280

CHROMATOGRAM**Retention time:** 5.5**Internal standard:** β -hydroxypropyltheophylline**Limit of detection:** 500 ng/mL

OTHER SUBSTANCES**Simultaneous:** theophylline, enoxacin, ciprofloxacin

KEY WORDSplasma; rat

REFERENCE

Davis, J.D.; Aarons, L.; Houston, J.B. Simultaneous assay of fluoroquinolones and theophylline in plasma by high-performance liquid chromatography, *J. Chromatogr.*, **1993**, 621, 105–109.

SAMPLE**Matrix:** blood

Sample preparation: 150 μ L Plasma + 10 μ L 100 μ g/mL enoxacin in water + 75 μ L 10% tri-chloroacetic acid + 600 μ L chloroform, vortex for 5 min, centrifuge at 13800 g for 10 min. Remove 500 μ L of the organic layer and evaporate it to dryness under reduced pressure, re-constitute the residue in 150 μ L mobile phase, inject a 50 μ L aliquot.

HPLC VARIABLES**Column:** 80 \times 4.6 5 μ m Zorbax C8**Mobile phase:** MeOH:0.01% trifluoroacetic acid 25:75**Flow rate:** 1.2**Injection volume:** 50**Detector:** F ex 280 em 418

CHROMATOGRAM**Retention time:** 6.9**Internal standard:** enoxacin (5.8)**Limit of quantitation:** 25 ng/mL

KEY WORDSplasma; rat; pharmacokinetics

REFERENCE

Hussain, M.S.; Chukwumaeze-Obiajunwa, V.; Micetich, R.G. Sensitive high-performance liquid chromatographic assay for norfloxacin utilizing fluorescence detection, *J. Chromatogr. B*, **1995**, 663, 379–384.

SAMPLE**Matrix:** blood

Sample preparation: 100 μ L Serum + 1 mL 2.6 μ g/mL N-ethylnorfloxacin in chloroform, vortex, centrifuge at 11000–12000 g for 1 min. Remove the organic layer and evaporate it to dryness under a stream of air at 60°, reconstitute the residue in 200 μ L mobile phase, inject a 25–50 μ L aliquot.

HPLC VARIABLES**Column:** 40 \times 3.2 3 μ m RP-18 Spheri-3**Mobile phase:** MeCN:10 mM pH 2.5 NaH_2PO_4 containing 1 mM triethylamine 11:89**Column temperature:** 30**Flow rate:** 1**Injection volume:** 25–50**Detector:** UV 279

CHROMATOGRAM**Retention time:** 1.9

Internal standard: N-ethylnorfloxacin (Heat 3.2 g norfloxacin, 1.5 g triethylamine, and 12–20 mmoles ethyl iodide in 40 mL DMF at 80–90° with stirring for 2 h, concentrate to dryness,

recrystallize from chloroform/benzene to give N-ethylnorfloxacin (mp 251-3°) (J.Med.Chem. 1980, 23, 1358)) (2.9)

Limit of detection: 20 ng/mL

KEY WORDS

serum

REFERENCE

Wallis,S.C.; Charles,B.G.; Gahan,L.R. Rapid and economical high-performance liquid chromatographic method for the determination of norfloxacin in serum using a microparticulate C18 guard cartridge, *J.Chromatogr.B*, 1995, 674, 306-309.

SAMPLE

Matrix: blood, dialysate

Sample preparation: 100 μ L Plasma or dialysate + 400 μ L MeOH, vortex, centrifuge, inject an aliquot of the supernatant.

HPLC VARIABLES

Column: strong cation exchange

Mobile phase: MeCN:100 mM pH 3 citrate buffer 20:80

Detector: F ex 278 em 440

CHROMATOGRAM

Limit of quantitation: 10 ng/mL

OTHER SUBSTANCES

Extracted: pefloxacin

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Rose,T.F.; Bremner,D.A.; Collins,J.; Ellis-Pegler,R.; Isaacs,R.; Richardson,R.; Small,M. Plasma and dialysate levels of pefloxacin and its metabolites in CAPD patients with peritonitis, *J.Antimicrob.Chemother.*, 1990, 25, 657-664.

SAMPLE

Matrix: blood, formulations

Sample preparation: Blood. Centrifuge 200 μ L fresh blood at 3000 rpm for 10 min. Inject an aliquot of the plasma. Formulations. Completely dissolve 50 mg sample in 20 mL MeOH, sonicate. Filter insoluble material and adjust filtrate to 50 mL with MeOH. Inject an aliquot of the filtrate.

HPLC VARIABLES

Column: 150 \times 4.6 Cosmosil 5C18 (Nacalai Tesque, Japan)

Mobile phase: MeCN:solution 1:6.5 (Solution was 1 g sodium acetate trihydrate, 2 g citric acid monohydrate and 1 mL triethylamine in 1 L water (?).)

Column temperature: 40

Flow rate: 1.5

Injection volume: 100

Detector: UV 277

CHROMATOGRAM

Internal standard: p-nitrophenylacetic acid

KEY WORDS

freeze-dried formulations; plasma; rat; egg albumin; olive oil

REFERENCE

Tsuji,Y.; Kakegawa,H.; Miyataka,H.; Nishiki,M.; Matsumoto,H.; Satoh,T. Pharmaceutical properties of freeze-dried formulations of egg albumin, several drugs and olive oil, *Biol.Pharm.Bull.*, 1996, 19, 636-640.

SAMPLE

Matrix: blood, tissue

Sample preparation: Plasma. Mix 500 μ L plasma with 5 μ g IS, add 4 mL dichloromethane and 100 μ L pH 7.4 phosphate buffer, agitate for 10 min, centrifuge at 5300 g for 10 min. Collect 3.5 mL organic phase. Add 4 mL dichloromethane to the aqueous phase again, agitate, centrifuge. Combine the organic phases, evaporate at 60°, reconstitute the residue in 100 μ L mobile phase, inject a 20 μ L aliquot. Tissue. Mix 500 μ L epiploic-fat and 4 mL dichloromethane and keep at 4°, add 5 μ g IS, mix by using an automatic grinder (Ultra Turrax, Ika-Werk, Stauffen, Germany). Collect the mixture, centrifuge at 5300 g for 10 min. Add 4 mL 100 mM NaOH to the dichloromethane, agitate for 10 min, centrifuge at 5300 g for 5 min. Eliminate the organic phase, adjust the aqueous phase to pH 7.4 with concentrated trichloroacetic acid, add 4 mL dichloromethane, agitate for 10 min and centrifuge at 5300 g for 5 min. Evaporate the organic phase at 60°. Reconstitute the residue in 100 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: 25 \times 4 5 μ m 100 RP-18 Lichrosphere

Column: 125 \times 4 5 μ m 100 RP-18 endcapped Lichrosphere

Mobile phase: MeCN:pH 4.8 citrate buffer 85:15

Flow rate: 1

Injection volume: 20

Detector: F ex 330 em 418

CHROMATOGRAM

Internal standard: 4844P (pefloxacin analog)

Limit of quantitation: 25 ng/mL

OTHER SUBSTANCES

Extracted: pefloxacin

KEY WORDS

epiploic-fat; plasma; pharmacokinetics; fat

REFERENCE

Jacoberger,B.; Ubeaud,G.; Freys,G.; Pottecher,T.; Jung,L.; Koffel,J.C. Concentrations of pefloxacin in plasma and tissue after administration as surgical prophylaxis, *Antimicrob.Agents Chemother.*, **1998**, 42, 425-427.

SAMPLE

Matrix: blood, tissue

Sample preparation: Tissue. Homogenize 100-250 mg tissue with 5 mL 500 mM pH 7.0 phosphate buffer, remove a 1 mL aliquot, add 100 μ L 10 μ g/mL IS, mix, add 10 mL chloroform:isopentanol 90:10, shake for 15 min, centrifuge, repeat extraction. Combine the organic layers and evaporate them to dryness under a stream of air at 60°, reconstitute the residue in 100 μ L 1% ammonia, inject a 25 μ L aliquot. Plasma. 250-500 μ L Plasma + 100 μ L 10 μ g/mL IS + 1 mL 500 mM pH 7.0 sodium phosphate buffer, mix, add 10 mL chloroform:isopentanol 90:10, shake for 15 min, centrifuge, repeat extraction. Combine the organic layers and evaporate them to dryness under a stream of air at 60°, reconstitute the residue in 100 μ L 1% ammonia, inject a 25 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 5 10 μ m Nucleosil C18

Mobile phase: MeCN:water 15:85 containing 2 g/L sodium acetate trihydrate, 2 g/L citric acid monohydrate, and 1 mL/L triethylamine

Flow rate: 2

Injection volume: 25

Detector: F ex 330 em 440

CHROMATOGRAM

Retention time: 2.8

Internal standard: 1-allyl-6-fluoro-1,4-dihydro-4-oxo-7-(4-methyl-1-piperazinyl)quinoline-3-carboxylic acid (Roger Bellon Laboratories) (6.6)

Limit of detection: 30 ng/mL

OTHER SUBSTANCES

Extracted: metabolites, pefloxacin

KEY WORDS

plasma; prostate

REFERENCE

Montay,G.; Tassel,J.P. Improved high-performance liquid chromatographic determination of pefloxacin and its metabolite norfloxacin in human plasma and tissue, *J.Chromatogr.*, **1985**, 339, 214–218.

SAMPLE

Matrix: blood, tissue, urine

Sample preparation: Serum, plasma. Dilute serum or plasma 1:2 to 1:10 with 30 mM phosphoric acid, centrifuge, inject a 20 μ L aliquot of supernatant. Urine. Dilute urine 1:10 to 1:100 with 30 mM phosphoric acid, centrifuge, inject a 20 μ L aliquot of supernatant. Tissue (lung, gut). Cut tissue with a scalpel, homogenize with 1-3 mL buffer, centrifuge at 9600 g for 5 min three times, inject a 20 μ L aliquot. Tissue (chondral). Cut tissue with a scalpel, homogenize with 3-6 mL buffer in an ice bath for 2-3 min, centrifuge at 9600 g for 5 min four or five times, inject a 100 μ L aliquot. Dilute human pleural samples with buffer, centrifuge, inject a 20 μ L aliquot. (Buffer was 66.6 mM K_2HPO_4 adjusted to pH 7.40 with KH_2PO_4 .)

HPLC VARIABLES

Column: 200 \times 4.5 μ m Nucleosil C18

Mobile phase: MeOH:MeCN:buffer 13:7:80, adjusted to pH 3.0 with phosphoric acid (Buffer was 15 mM phosphoric acid adjusted to pH 3.0 with tetrabutylammonium hydroxide.)

Flow rate: 1

Injection volume: 20-100

Detector: F ex 278 em 446

CHROMATOGRAM

Retention time: 4

Limit of detection: 2.5 ng/mL

OTHER SUBSTANCES

Simultaneous: ofloxacin, ciprofloxacin

KEY WORDS

serum; plasma; lung; gut; pleural; chondral

REFERENCE

Knöller,J.; König,W.; Schönfeld,W.; Bremm,K.D.; Köller,M. Application of high-performance liquid chromatography of some antibiotics in clinical microbiology, *J.Chromatogr.*, **1988**, 427, 257–267.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 500 μ L MeCN to 500 μ L plasma or urine diluted with water, vortex vigorously for 20 s, centrifuge at 4000 rpm for 2 min, mix a 250 μ L aliquot of the supernatant with 500 μ L 100 mM perchloric acid containing 20 mM triethylamine, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 80 \times 4.0 3 μ m Nucleosil 3C18

Mobile phase: MeOH:100 mM perchloric acid containing 20 mM triethylamine 30:70

Column temperature: 40

Flow rate: 1

Injection volume: 20

Detector: F ex 300 em 450

CHROMATOGRAM

Retention time: 2.09

OTHER SUBSTANCES

Noninterfering: acyclovir, ampicillin, amoxicillin, clavulanic acid, doxycycline, erythromycin, lansoprazole, metronidazole, minocycline, omeprazole, penicillin V, trimethoprim

KEY WORDS

plasma

REFERENCE

Mascher,H.J.; Kikuta,C. Determination of norfloxacin in human plasma and urine by high-performance liquid chromatography and fluorescence detection, *J.Chromatogr.A*, **1998**, 812, 381–385.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. Wash a C18 Sep-Pak cartridge with 1 mL 4% methanolic phosphoric acid and 10 mL water. Add 1 mL plasma to the cartridge, wash with 6 mL water, elute with 1 mL 4% methanolic phosphoric acid then with 1 mL water. Combine the eluates, make up to 2 mL with water, inject a 50 µL aliquot. Urine. Dilute 100 µL urine with 900 µL mobile phase and directly inject a 10 µL aliquot.

HPLC VARIABLES

Guard column: µBondapak C18-Corasil

Column: 300 × 3.9 µBondapak C18

Mobile phase: MeOH:MeCN:water 300:50:700 + 1.74 g K₂HPO₄ + 20 mg sodium heptanesulfonate, pH adjusted to 3 with phosphoric acid

Flow rate: 2

Injection volume: 10 (urine), 50 (plasma)

Detector: F ex 285 em 440

CHROMATOGRAM

Retention time: 4

Limit of detection: 500 ng/mL (urine), 20 ng/mL (plasma)

KEY WORDS

plasma; SPE

REFERENCE

Gutzler,F.; de Vries,J.X. Bestimmung von Norfloxacin in Plasma und Urin durch Hochdruckflüssigkeitschromatographie, *Fortschr.Antimikr.Antineoplast.Chemother.*, **1984**, 3, 673–677.

SAMPLE

Matrix: blood, urine

Sample preparation: Dilute with one or more volumes of water, filter (0.6 µm)

HPLC VARIABLES

Column: 200 × 4 5 µm Nucleosil C18

Mobile phase: MeCN:25 mM orthophosphoric acid adjusted to pH 3.0 with tetrabutylammonium hydroxide 11:89

Flow rate: 1.5

Injection volume: 10-20

Detector: F ex 278 em 445

CHROMATOGRAM

Retention time: 2.8

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Noninterfering: trimethoprim, sulfamethoxazole, netilmicin, metronidazole, penicillin G, cloxacillin, doxycycline, cefuroxime, erythromycin, salicylic acid, digoxin, furosemide, acetaminophen, prednisolone, warfarin, dextropropoxyphene

Interfering: ciprofloxacin

KEY WORDS

serum

REFERENCE

Nilsson-Ehle,I. Assay of ciprofloxacin and norfloxacin in serum and urine by high-performance liquid chromatography, *J.Chromatogr.*, **1987**, 416, 207–211.

SAMPLE**Matrix:** blood, urine

Sample preparation: Plasma. 250 μ L Plasma + 250 μ L 10 μ g/mL IS in water + 750 μ L MeOH, stir, centrifuge at 2000 rpm for 5 min, inject a 50 μ L aliquot of the supernatant. Urine. 500 μ L Urine + 3.5 mL 8 μ g/mL IS in water, inject a 50 μ L aliquot.

HPLC VARIABLES**Column:** 200 \times 4.6 Nucleosil C8

Mobile phase: MeCN:water:triethylamine:formic acid 11:87.5:0.1:1 containing 0.2% sodium acetate and 0.1% formic acid

Flow rate: 1**Injection volume:** 50**Detector:** F ex 280 em 450

CHROMATOGRAM

Internal standard: 1-ethyl-6-chloro-1,4-dihydro-7-(4-methyl-1-piperazinyl)-4-oxo-3-quinoline-carboxylic acid (RP 41983)

Limit of quantitation: 500 ng/mL (urine), 100 ng/mL (plasma)

OTHER SUBSTANCES**Extracted:** pefloxacin

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Humbert,G.; Brumpt,I.; Montay,G.; Le Liboux,A.; Frydman,A.; Borsa-Lebas,F.; Moore,N. Influence of rifampin on the pharmacokinetics of pefloxacin, *Clin.Pharmacol.Ther.*, **1991**, 50, 682–687.

SAMPLE**Matrix:** cells

Sample preparation: Incubate cells in 2 mL 100 mM pH 3.0 glycine-HCl buffer for 2 h at room temperature, centrifuge at 5600 g for 5 min, inject an aliquot.

HPLC VARIABLES**Column:** Bondapak C18

Mobile phase: MeCN:25 mM phosphoric acid adjusted to pH 3.0 with tetrabutylammonium hydroxide 25:75

Flow rate: 1.5**Detector:** F ex 340 em 425

OTHER SUBSTANCES

Also analyzed: ciprofloxacin, fleroxacin, lomefloxacin, ofloxacin, temafloxacin

REFERENCE

Pascual,A.; Garcia,I.; Conejo,M.C.; Perea,E.J. Fluorometric and high-performance liquid chromatographic measurement of quinolone uptake by human neutrophils, *Eur.J.Clin.Microbiol.Infect.Dis.*, **1991**, 10, 969–971.

SAMPLE**Matrix:** hair

Sample preparation: Wash hair successively with 0.1% sodium dodecyl sulfate and water for 30 min, repeat twice, blot between 2 sheets of paper towel, allow to dry at room temperature. Take a 1 cm fragment of hair, add 500 μ L 1 M NaOH, heat at 80° for 30 min, cool, add 500

μL 1 M HCl, add 1 mL 100 mM pH 4.6 potassium hydrogen citrate buffer, add 50 μL 1 $\mu\text{g/mL}$ IS in water. Add the mixture to a Bond-Elut C8 cartridge, elute with 2 mL THF:25 mM orthophosphoric acid 20:80, evaporate eluate to dryness in vacuum, dissolve residue in 150 μL mobile phase, vortex, inject a 60 μL aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 Tosoh 5 μm TSKgel ODS-80Ts

Mobile phase: MeCN:25 mM orthophosphoric acid adjusted to pH 3.0 with 0.5 M tetra-n-butylamine hydroxide 5:95

Column temperature: 40

Flow rate: 1

Injection volume: 60

Detector: F ex 280 em 445

CHROMATOGRAM

Retention time: 11.7

Internal standard: (R)-9-fluoro-2,3-dihydro-3-methyl-10-(4-ethyl-1-piperazinyl)-7-oxo-7H-pyrido[1,2,3-de][1,4]benzoxazine-6-carboxylic acid (DS-4632) (10.2)

Limit of detection: 0.2 ng/mL

OTHER SUBSTANCES

Simultaneous: ciprofloxacin, ofloxacin (determine at F ex 295 em 490)

KEY WORDS

SPE

REFERENCE

Mizuno,A.; Uematsu,T.; Nakashima,M. Simultaneous determination of ofloxacin, norfloxacin and ciprofloxacin in human hair by high-performance liquid chromatography and fluorescence detection, *J.Chromatogr.B*, **1994**, 653, 187–193.

SAMPLE

Matrix: perfusate

HPLC VARIABLES

Column: 250 \times 4.6 Spheris C18 (Phase Separations)

Mobile phase: MeCN:15 mM tetrabutylammonium iodide 5:95

Flow rate: 1

Detector: UV 275

CHROMATOGRAM

Internal standard: pipemidic acid

REFERENCE

Lin,H.-H.; Hsu,L.-R.; Wu,P.-C.; Tsai,Y.-H. Increased norfloxacin skin permeability for fatty alcohol propylene glycol (FAPG) ointment by optimized process of preparation: Behavior of stearic acid in stratum corneum lipids, *Biol.Pharm.Bull.*, **1995**, 18, 1560–1565.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 450 $\mu\text{g/mL}$ solution in MeCN:water 50:50. 5 mL Solution + 5 mL THF + 200 molar excess of acetic anhydride + 3 molar excess of 1 M NaOH, sonicate for 15 min, add 15 mL mobile phase, sonicate for 15 min, cool to room temperature, make up to 50 mL with mobile phase, inject a 50 μL aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μm Nucleosil C18

Mobile phase: MeCN:buffer 35:65 (Buffer was prepared by mixing equal volumes of 20 mM citric acid and 20 mM sodium citrate, pH adjusted to 2.4 with perchloric acid.)

Flow rate: 1

Injection volume: 50

Detector: UV 280

CHROMATOGRAM

Retention time: 5.9

OTHER SUBSTANCES

Simultaneous: ciprofloxacin, sarafloxacin, temafloxacin

KEY WORDS

derivatization

REFERENCE

Morley, J.A.; Elrod, L., Jr. Determination of fluoroquinolone antibacterials as N-Acyl derivatives, *Chromatographia*, **1993**, *37*, 295–299.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 20 µg/mL solution in MeCN:water 10:90, filter (0.45 µm), inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.5 µm LiChrospher 100 RP-18

Mobile phase: MeCN:25 mM phosphoric acid 7:93, adjusted to pH 3.09 with 100 mM tetrabutylammonium hydroxide

Flow rate: 1

Injection volume: 10

Detector: UV 280

CHROMATOGRAM

Retention time: 8.6

OTHER SUBSTANCES

Simultaneous: ciprofloxacin, ofloxacin (UV 295), pipemidic acid

REFERENCE

Barbosa, J.; Bergés, R.; Sanz-Nebot, V. Linear solvation energy relationships in reversed-phase liquid chromatography. Prediction of retention of several quinolones, *J. Liq. Chromatogr.*, **1995**, *18*, 3445–3463.

SAMPLE

Matrix: solutions

Sample preparation: Filter (0.45 µm) a solution in MeCN:water 10:90, inject an aliquot of the filtrate.

HPLC VARIABLES

Column: 250 × 4.5 µm LiChrospher 100 RP-18

Mobile phase: MeCN:buffer 7:93 (Buffer was 25 mM phosphoric acid adjusted to pH 3.89 with 100 mM tetrabutylammonium hydroxide.)

Flow rate: 1

Injection volume: 10

Detector: UV 280

CHROMATOGRAM

Retention time: 8.8

OTHER SUBSTANCES

Simultaneous: ciprofloxacin, enoxacin, fleroxacin, ofloxacin (UV 295), pipemidic acid

REFERENCE

Barbosa, J.; Bergés, R.; Sanz-Nebot, V. Solvatochromic parameter values and pH in aqueous-organic mixtures used in liquid chromatography. Prediction of retention of a series of quinolones, *J. Chromatogr. A*, **1996**, *719*, 27–36.

SAMPLE**Matrix:** tissue**Sample preparation:** Homogenize (Ultra-Turrax) eye tissue with 3 mL 50 mM pH 5.8 sodium phosphate-citrate buffer and IS, centrifuge. Add the supernatant to 7 mL chloroform, agitate, centrifuge at 1000 g for 10 min. Remove the lower organic layer and evaporate it to dryness under a stream of nitrogen at 37°, reconstitute the residue in 100 µL mobile phase, inject a 5-20 µL aliquot.

HPLC VARIABLES**Column:** 100 × 4.6 3 µm Nucleosil C8**Mobile phase:** MeCN:water 26:74 containing 2 g/L sodium acetate trihydrate, 2 g/L citric acid monohydrate, 4 mL/L triethylamine, and 2 mL/L formic acid, pH 4.8**Flow rate:** 1**Injection volume:** 5-20**Detector:** UV 280

CHROMATOGRAM**Retention time:** 2.19**Internal standard:** 1-allyl-6-fluoro-1,4-dihydro-4-oxo-7-(4-methyl-1-piperazinyl)quinoline-3-carboxylic acid (?) 4662 P (Roger-Bellon Laboratories) (2.95)**Limit of detection:** 5 ng

OTHER SUBSTANCES**Extracted:** pefloxacin

KEY WORDSrabbit; eye; pharmacokinetics

REFERENCE

Cochereau-Massin,I.; Bauchet,J.; Faurisson,F.; Vallois,J.M.; Lacombe,P.; Pocidalò,J.J. Ocular kinetics of pefloxacin after intramuscular administration in albino and pigmented rabbits, *Antimicrob.Agents Chemother.*, **1991**, 35, 1112-1115.

SAMPLE**Matrix:** urine**Sample preparation:** Dilute with water, inject an aliquot.

HPLC VARIABLES**Column:** µBondapak C18**Mobile phase:** MeOH:MeCN:100 mM pH 5.75 phosphate buffer 24.1:2.6:73.3**Flow rate:** 1**Detector:** F ex 275 em 415

CHROMATOGRAM**Limit of quantitation:** 3.13 µg/mL

OTHER SUBSTANCES**Extracted:** pefloxacin

KEY WORDSpharmacokinetics

REFERENCE

Jaehde,U.; Sörgel,F.; Stephan,U.; Schunack,W. Effect of an antacid containing magnesium and aluminum on absorption, metabolism, and mechanism of renal elimination of pefloxacin in humans, *Antimicrob.Agents Chemother.*, **1994**, 38, 1129-1133.


Norgestimate

Molecular formula: $C_{23}H_{31}NO_3$

Molecular weight: 369.50

CAS Registry No.: 35189-28-7

Merck Index: 6796



SAMPLE

Matrix: blood

Sample preparation: 1 mL Serum + 3 mL MTBE, vortex for 1 min, centrifuge at 1500 rpm for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at room temperature, reconstitute the residue in 100 μ L MeOH, inject an aliquot.

HPLC VARIABLES

Guard column: 18 mm long (Brownlee)

Column: 300 \times 3.6 10 μ m μ Bondapak C18

Mobile phase: MeOH:water 80:20

Flow rate: 1

Detector: UV 254 or RIA

CHROMATOGRAM

Retention time: 9.3, 9.8 (syn and anti)

OTHER SUBSTANCES

Extracted: metabolites, norgestrel

KEY WORDS

serum

REFERENCE

Wong,F.A.; Juzwin,S.J.; Tischio,N.S.; Flor,S.C. Determination of norgestimate in serum by automated high-performance liquid chromatography and subsequent radioimmunoassay, *J.Liq. Chromatogr.*, **1995**, *18*, 1851-1861.

SAMPLE

Matrix: culture media

Sample preparation: Extract culture medium twice with 2 volumes of ether, combine the extracts and evaporate them to dryness, reconstitute with MeOH, inject an aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 Techopak 10 C18 (HPLC Technology)

Mobile phase: MeOH:water 70:30

Flow rate: 1.5

Detector: UV 240, radioactivity

CHROMATOGRAM

Retention time: 13.5, 15.5 (racemate)

OTHER SUBSTANCES

Extracted: metabolites, norgestrel

KEY WORDS

tritium labeled

REFERENCE

Wild,M.J.; Rudland,P.S.; Back,D.J. Metabolism of the oral contraceptive steroids ethynylestradiol and norgestimate by normal (Huma 7) and malignant (MCF-7 and ZR-75-1) human breast cells in culture, *J.Steroid Biochem.Mol.Biol.*, **1991**, *39*, 535-543.

SAMPLE**Matrix:** formulations**Sample preparation:** 5 Tablets + 2 glass beads + 25 mL 50 µg/mL dibutyl phthalate in MeOH, vortex 15 min or until tablets have completely disintegrated, sonicate 5 min, filter (2 µm), inject 25 µL aliquot.

HPLC VARIABLES**Column:** 50 × 4.5 5µm IBM C18**Mobile phase:** MeOH:THF:water 10:25:65**Flow rate:** 2.1**Injection volume:** 25**Detector:** UV 230

CHROMATOGRAM**Retention time:** 3.5**Internal standard:** dibutyl phthalate

OTHER SUBSTANCES**Simultaneous:** ethinylestradiol, degradation products

KEY WORDStablets; stability-indicating

REFERENCELane,P.A.; Mayberry,D.O.; Young,R.W. Determination of norgestimate and ethinyl estradiol in tablets by high-performance liquid chromatography, *J.Pharm.Sci.*, **1987**, 76, 44–47.

SAMPLE**Matrix:** microsomal incubations, mucosal fluid**Sample preparation:** Mucosal fluid. Extract 1 mL mucosal fluid twice with 5 mL diethyl ether, evaporate extracts to dryness, resuspend residue in 100 µL MeOH, inject an aliquot. Microsomal incubations. Extract 2.5 mL microsomal incubation with 5 mL diethyl ether, proceed as before.

HPLC VARIABLES**Guard column:** on-line guard column**Column:** 100 × 8 µBondapak radial compression module**Mobile phase:** MeOH:water 70:30**Flow rate:** 1.5**Injection volume:** 100**Detector:** UV 240

CHROMATOGRAM**Retention time:** 17.5

OTHER SUBSTANCES**Simultaneous:** 3-ketonorgestimate, 17-deacetylnorgestimate, norgestrel

REFERENCEMadden,S.; Back,D.J. Metabolism of norgestimate by human gastrointestinal mucosa and liver microsomes in vitro, *J.Steroid Biochem.Mol.Biol.*, **1991**, 38, 497–503.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 250 × 4.6 µBondapak C18**Mobile phase:** Dioxane:water 50:50 (CAUTION! Dioxane is a carcinogen!)**Flow rate:** 1.4**Detector:** UV 254

OTHER SUBSTANCES

Simultaneous: levonorgestrel

REFERENCE

Killinger, J.; Hahn, D.W.; Phillips, A.; Heteyi, N.S.; McGuire, J.L. The affinity of norgestimate for uterine progesterone receptors and its direct action on the uterus, *Contraception*, **1985**, 32, 311–319.

SAMPLE

Matrix: tissue

Sample preparation: Incubate endometrial tissue with buffer, remove tissue, extract medium twice with 2 volumes of diethyl ether, evaporate to dryness, reconstitute in a small volume of MeOH, inject an aliquot.

HPLC VARIABLES

Column: 300 × 3.9 Technopak 10 C18

Mobile phase: MeOH:water 70:30

Flow rate: 1.5

Detector: UV 240

CHROMATOGRAM

Retention time: 15, 18 (syn and anti)

OTHER SUBSTANCES

Simultaneous: norgestrel, metabolites

KEY WORDS

endometrial tissue

REFERENCE

Wild, M.J.; Rudland, P.S.; Back, D.J. Metabolism of the oral contraceptive steroids ethynylestradiol, norgestimate and 3-ketodesogestrel by a human endometrial cancer cell line (HEC-1A) and endometrial tissue *in vitro*, *J.Steroid Biochem.Mol.Biol.*, **1993**, 45, 407–420.

Norgestrel


Molecular formula: C₂₁H₂₈O₂

Molecular weight: 312.45

CAS Registry No.: 797-63-7, 797-64-8 ((-) form), 6533-00-2

Merck Index: 6797

Lednicer No.: 1 167, 2 151, 3 84

**SAMPLE**

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 241.7

CHROMATOGRAM

Retention time: 21.565

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

SAMPLE

Matrix: formulations

Sample preparation: Centrifuge oil formulation at 30° at 2000 rpm for 30 min, filter (Whatman No. 1 paper), collect the last 4 mL of the filtrate. Dilute a 10 µL aliquot to 10 mL with MeCN: water 60:40 containing 0.3% Tween 80, add a 2 mL aliquot and add it to 1 mL 3.33 µg/mL progesterone, vortex for 10 s, inject a 50 µL aliquot.

HPLC VARIABLES

Column: 300 × 3.9 Novapak C18

Mobile phase: MeCN:water 60:40

Flow rate: 2

Injection volume: 50

Detector: UV 248

CHROMATOGRAM

Internal standard: progesterone

KEY WORDS

oils

REFERENCE

Gao,Z.-H.; Shukla,A.J.; Johnson,J.R.; Crowley,W.R. Controlled release of a contraceptive steroid from biodegradable and injectable gel formulations: In vitro evaluation, *Pharm.Res.*, **1995**, 12, 857–863.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 200 × 4.6 5 µm Hypersil ODS

Mobile phase: MeCN:water 50:50

Column temperature: 37

Flow rate: 2.0

Detector: UV 243

CHROMATOGRAM

Retention time: 4.07

REFERENCE

Kim,D.-D.; Kim,J.L.; Chien,Y.W. Mutual hairless rat skin permeation-enhancing effect of ethanol/water system and oleic acid, *J.Pharm.Sci.*, **1996**, 85, 1191–1195.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 μBondapak C18

Mobile phase: Dioxane:water 50:50 (CAUTION! Dioxane is a carcinogen!)

Flow rate: 1.4

Detector: UV 254

OTHER SUBSTANCES

Simultaneous: norgestimate

REFERENCE

Killinger,J.; Hahn,D.W.; Phillips,A.; Heteyi,N.S.; McGuire,J.L. The affinity of norgestimate for uterine progesterogen receptors and its direct action on the uterus, *Contraception*, **1985**, 32, 311–319.

SAMPLE

Matrix: solutions

Sample preparation: Direct injection

HPLC VARIABLES

Column: 250 × 4.6 10 μm Partisil C-18 ODS-3

Mobile phase: MeCN:water 50:50

Flow rate: 2

Detector: UV 243

CHROMATOGRAM

Retention time: 6.0

KEY WORDS

see also J.Pharm.Sci. 1989; 78; 477

REFERENCE

Catz,P.; Friend,D.R. In vitro evaluations of transdermal levonorgestrel, *Drug Des.Deliv.*, **1990**, 6, 49–60.

Nortriptyline


Molecular formula: C₁₉H₂₁N

Molecular weight: 263.38

CAS Registry No.: 72-69-5, 894-71-3 (HCl)

Merck Index: 6812

Lednicer No.: 1 151



SAMPLE

Matrix: bile, blood, gastric contents, tissue, urine

Sample preparation: Chop 5-g tissue and homogenize (Ultra Turrax T25) at 8500, 9500, 13500, 20500, and 24000 rpm for 1 min each. Add homogenate to 20 mL water. Dilute blood, urine, gastric contents, and bile four times with water. Mix 4 mL sample with 100 μL 400 μg/mL IS and 2 mL 500 mM NaOH, vortex briefly, add 4 mL heptane:isoamyl alcohol 98.5:1.5 and mix for 15 min (Spiramix 10, Denley, UK). Separate the organic layer, add 4 mL heptane:isoamyl alcohol 98.5:1.5 to extraction sample, mix. Combine the organic layers and extract them with 2 mL 50 mM sulfuric acid. Make the acid layer alkaline with 1 mL 1.0 M pH 9.0 carbonate/bicarbonate buffer and mix with 2 mL toluene:isoamyl alcohol 85:15 for 15 min. Evaporate the organic layer to dryness, reconstitute the residue in 100 μL MeOH and inject a 20 μL aliquot.

HPLC VARIABLES

Guard column: 20 × 4.6 5 μm Apex II ODS

Column: 150 × 4.6 5 µm Apex II OD

Mobile phase: MeCN:pH 3 phosphate buffer:n-nonylamine 40-50:60:0.12

Flow rate: 1

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 3.70

Internal standard: doxepin (2.99)

OTHER SUBSTANCES

Extracted: amitriptyline

KEY WORDS

liver; lung; muscle; urine; pericardial fluid

REFERENCE

Pounder,D.J.; Adams,E.; Fuke,C.; Langford,A.M. Site to site variability of postmortem drug concentrations in liver and lung, *J.Forensic Sci.*, **1996**, *41*, 927–932.

SAMPLE

Matrix: blood

Sample preparation: Add 250 µL 2 M sodium carbonate to 500 µL plasma. Add 100 µL 1 µg/mL IS in MeOH, extract with 10 mL n-hexane. Shake for 30 min and centrifuge at 3000 g for 10 min. Cool in a dry ice-acetone bath. Add 200 µL 0.3% phosphoric acid to upper organic layer. Shake for 10 min and centrifuge at 3000 g for 10 min. Separate the organic layer. Inject a 100 µL aliquot of the acidic aqueous layer.

HPLC VARIABLES

Column: 250 × 4.6 5 µm C18 Symmetry (Waters Millipore, USA)

Mobile phase: MeCN:67 mM potassium phosphate buffer adjusted to pH 3.0 with phosphoric acid 35:65 (After each chromatographic session wash the column with 200 mL MeCN:water 50:50.)

Flow rate: 1.2

Injection volume: 100

Detector: UV 226, UV 254, UV 400

CHROMATOGRAM

Retention time: 10.45

Internal standard: clovoxamine (6.5)

Limit of quantitation: 5 ng/mL (UV 226, UV 400 nm); 7 ng/mL (UV 254)

OTHER SUBSTANCES

Extracted: metabolites, amitriptyline, clomipramine, desipramine, fluoxetine imipramine, maprotiline

Simultaneous: amineptine, carbamazepine, chlordiazepoxide, chlorpromazine, clonazepam, clorazepate, clozapine, cyamemazine, desmethylmaprotiline, desmethylvenlafaxine, doxepin, flunitrazepam, fluvoxamine, haloperidol, levomepromazine, lorazepam, loxapine, mianserine, sulpiride, trimipramine, venlafaxine, viloxazine, zolpidem, zopiclone

Noninterfering: diazepam, valproic acid

KEY WORDS

plasma

REFERENCE

Aymard,G.; Livi,P.; Pham,Y.T.; Diquet,B. Sensitive and rapid method for the simultaneous quantification of five antidepressants with their respective metabolites in plasma using high-performance liquid chromatography with diode-array detection, *J.Chromatogr.B*, **1997**, *700*, 183–189.

SAMPLE

Matrix: blood

Sample preparation: Condition a 1 mL 30 mg Oasis HLB SPE cartridge with 1 mL MeOH and 1 mL water. Acidify (?) mL serum with 20 μ L phosphoric acid, vortex for 5 s, add to the SPE cartridge, wash with 1 mL MeOH :water 5:95, elute with 1 mL MeOH. Evaporate the eluate to dryness at 40° under a stream of nitrogen. Reconstitute the residue with 200 μ L MeOH:20 mM pH 7 phosphate buffer 20:80, inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: 20 \times 3.9 Sentry

Column: 150 \times 3.9 5 μ m Symmetry C18 (Waters)

Mobile phase: MeOH:20 mM pH 7 potassium phosphate 70:30

Column temperature: 35

Flow rate: 1

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 8

Internal standard: nordoxepin (4.9)

OTHER SUBSTANCES

Simultaneous: amitriptyline, doxepin

KEY WORDS

pig; serum; SPE

REFERENCE

Cheng,Y.-F.; Phillips,D.J.; Neue,U.; Bean,L. Solid-phase extraction for the determination of tricyclic antidepressants in serum using a novel polymeric extraction sorbent, *J.Liq.Chromatogr.Rel.Technol.*, **1997**, 20, 2461-2473.

SAMPLE

Matrix: blood, microsomal incubations

Sample preparation: Vortex 1 mL plasma or microsomal incubation with 200 μ L 1 μ g/mL desipramine and 100 μ L 5 M NaOH for 10 s, add 5 mL butan-1-ol:hexane 2:98, vortex for 1 min, centrifuge at 2000 g and 4° for 5 min, evaporate the organic phase to dryness at 40° using a vacuum vortex evaporator, reconstitute the residue in 200 μ L mobile phase, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 5 μ m Nova-Pak C18

Mobile phase: MeCN:buffer 30:70 (Buffer was water containing 1% triethylamine, adjusted to pH 3 with orthophosphoric acid.)

Flow rate: 2

Injection volume: 50

Detector: UV 240

CHROMATOGRAM

Retention time: 8.5

Internal standard: desipramine (6.3)

Limit of quantitation: 2 ng/mL

OTHER SUBSTANCES

Extracted: amitriptyline

Noninterfering: diazepam, furafylline, hydroxyamitriptyline, hydroxynortriptyline, quinidine, mephenytoin, triacetyloleandomycin

Interfering: ketoconazole

KEY WORDS

human; liver; rat; plasma

REFERENCE

Ghahramani,P.; Lennard,M.S. Quantitative analysis of amitriptyline and nortriptyline in human plasma and liver microsomal preparations by high-performance liquid chromatography, *J.Chromatogr.B*, **1996**, 685, 307-313.

SAMPLE

Matrix: blood, tissue, urine

Sample preparation: Serum, urine. 500 μ L Serum or urine + 100 μ L 2 μ g/mL diazepam + 200 μ L 20% sodium carbonate + 500 μ L water + 3 mL n-hexane:isoamyl alcohol 98.5:1.5, mix for 2 min, centrifuge at 1200 g for 5 min. Remove the organic phase and evaporate it under a gentle stream of nitrogen at about 40°. Dissolve the residue in 100 μ L mobile phase, inject a 10 μ L aliquot. Tissue. Homogenize 1 g sample with 9 mL 100 mM HCl and 100 μ L 20 μ g/mL diazepam, centrifuge at 15000 g for 10 min. Add 500 μ L 20% sodium carbonate and 4 mL n-hexane:isoamyl alcohol 98.5:1.5 to 1 mL of the supernatant, mix for 5 min. Remove the organic phase and evaporate it under a gentle stream of nitrogen at about 40°. Dissolve the residue in 100 μ L mobile phase, filter by microconcentrator (Microcon-30, Grace). Inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 4.6 2 μ m TSK gel Super-Octyl (A) or 100 \times 4.6 5 μ m Hypersil MOS-C8 (B), (Yokogawa, Japan)

Mobile phase: MeOH:20 mM pH 7 KH_2PO_4 60:40

Flow rate: 0.6

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 7.0 (A), 9.2 (B)

Internal standard: diazepam (4.4, A)

Limit of quantitation: 50 ng/mL (serum, urine), 500 ng/mL (tissue)

OTHER SUBSTANCES

Extracted: amitriptyline, amoxapine, clomipramine, desipramine, dothiepin, doxepin, imipramine, maprotiline, melitracen, mianserin

Noninterfering: barbital, carbamazepine, ethosuximide, hexobarbital, lofepramine, pentobarbital, phenobarbital, phenytoin, primidone, sulpiride, trimethadione, trimipramine

KEY WORDS

serum; brain; liver

REFERENCE

Tanaka,E.; Terada,M.; Nakamura,T.; Misawa,S.; Wakasugi,C. Forensic analysis of eleven cyclic antidepressants in human biological samples using a new reversed-phase chromatographic column of 2 μ m porous microspherical silica gel, *J.Chromatogr.B*, **1997**, 692, 405-412.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 500 μ L 3 M ammonia solution and 7 mL n-pentane:isopropanol 95:5 to 2 mL plasma or urine, shake in an overhead shaker for 20 min, let stand for 10 min. Transfer the upper organic layer to a tube containing 1 mL 100 mM HCl, shake for 20 min, let stand for 5 min. Aspirate the organic phase to waste, wash the remaining aqueous layer with 3 mL pentane by shaking for 10 min. Add 500 μ L 3 M ammonia solution and 6 mL a-pentane: isopropanol 95:5 to the washed aqueous layer, shake for 20 min, evaporate the organic layer under a stream of nitrogen at 65°, reconstitute the residue with 160 μ L mobile phase, inject 60 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.5 3 μ m Spherisorb silica

Mobile phase: MeOH:hexane:nonylamine 5:95:0.3

Flow rate: 1

Injection volume: 60

Detector: UV 254

CHROMATOGRAM**Retention time:** 10**Internal standard:** nortriptyline

OTHER SUBSTANCES**Extracted:** N-desmethyldoxepin, doxepin

KEY WORDSdog; human; plasma; normal phase; nortriptyline is IS

REFERENCE

Yan, J.; Hubbard, J.W.; McKay, G.; Midha, K.K. Stereoselective and simultaneous measurement of cis- and trans-isomers of doxepin and N-desmethyldoxepin in plasma or urine by high-performance liquid chromatography, *J.Chromatogr.B*, **1997**, 691, 131–138.

SAMPLE**Matrix:** blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES**Guard column:** 20 mm long Symmetry C18**Column:** 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10–30**Detector:** UV 206.4

CHROMATOGRAM**Retention time:** 15.603

KEY WORDSwhole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

SAMPLE**Matrix:** serum

Sample preparation: 1 mL Serum + 500 µL 750 mM pH 10 sodium bicarbonate/carbonate buffer + 50 µL IS in EtOH:water 50:50 + 8 mL heptane:isoamyl alcohol 98:2, shake at 250 cycles/min for 5 min, centrifuge at 1500 g for 10 min, freeze in dry ice/EtOH. Remove the organic layer and add it to 150 µL 22 mM pH 2.5 KH₂PO₄/phosphoric acid buffer, shake at 250 cycles/min for 5 min, centrifuge at 1500 g for 10 min, freeze in dry ice/EtOH. Discard the organic layer, inject a 65 µL aliquot of the aqueous layer.

HPLC VARIABLES**Column:** 250 × 4.6 Supelco C18

Mobile phase: MeCN:buffer 45:55 (Buffer was 44 mM KH_2PO_4 containing 1.5 mL/L triethylamine, adjusted to pH 2.5 with phosphoric acid.)

Flow rate: 1.5

Injection volume: 65

Detector: UV 240

CHROMATOGRAM

Retention time: 9.73

Internal standard: 1-(3-(dimethylamino)propyl)-1-(p-chlorophenyl)-1,3-dihydroisobenzofuran-5-carbonitrile (LU 10-202) (Lundbeck, Copenhagen) (8.33)

OTHER SUBSTANCES

Extracted: metabolites, citalopram, amitriptyline

Simultaneous: chlorprothixene, clomipramine, clozapine, flupenthixol, haloperidol, levomepromazine, perphenazine, zuclopenthixol

Noninterfering: benzodiazepines

Interfering: desmethyllevomepromazine

KEY WORDS

serum

REFERENCE

Olesen, O.V.; Linnet, K. Simplified high-performance liquid chromatographic method for the determination of citalopram and desmethylcitalopram in serum without interference from commonly used psychotropic drugs and their metabolites, *J. Chromatogr. B*, **1996**, 675, 83–88.


Noscapine

Molecular formula: $\text{C}_{22}\text{H}_{23}\text{NO}_7$

Molecular weight: 413.43

CAS Registry No.: 128-62-1, 6035-40-1 (dl-form), 912-60-7 (HCl)

Merck Index: 6815



SAMPLE

Matrix: blood

Sample preparation: Condition a C18 SPE cartridge with 2 mL MeOH, 2 mL water, and 1 mL 50 mM pH 7.4 sodium phosphate buffer at 1.5 mL/min. Centrifuge rapidly thawed plasma samples. Combine 1 mL plasma with 40 μL 903 ng/mL papaverine free base and 2 mL 50 mM pH 7.4 sodium phosphate buffer, vortex. Add 2.75 mL diluted plasma at a rate of 0.36 mL/min to the SPE cartridge, wash with 1 mL 50 mM pH 7.4 sodium phosphate buffer and 2 mL water. Dry by passing 5 mL of air through the SPE cartridge at 5 mL/min. Elute with 1.5 mL MeCN, evaporate the eluate to dryness under nitrogen at 30°. Reconstitute with 200 μL MeCN, inject a 50 μL aliquot.

HPLC VARIABLES

Column: 125 \times 4.0 5 μm Nucleosil 120-5-C18

Mobile phase: MeCN:25 mM pH 4.5 sodium phosphate buffer 40:60

Column temperature: 25

Flow rate: 1.0

Injection volume: 50

Detector: UV 211

CHROMATOGRAM

Retention time: 4.99

Internal standard: papaverine (3.98)

Limit of detection: 0.9 ng/mL

Limit of quantitation: 7.2 ng/mL

KEY WORDS

plasma; SPE

REFERENCE

Chollet,D.F.; Ruols,C.; Arnera,V. Determination of noscapine in human plasma using solid-phase extraction and high-performance liquid chromatography, *J.Chromatogr.B*, **1997**, 701, 81–85.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 213.4

CHROMATOGRAM

Retention time: 12.827

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 μ g/mL solution in MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 125 \times 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 1.2

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquin-

amide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipiprone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenazpromide, phenazocine, phenbutazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindone, pimizole, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, propriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinamine, sotalol, tacrine, terazosin, terbitaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, 323, 191–225.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapson, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, dil-

tiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenpropfen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxygenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephentoin, mephesis, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methypylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nyldrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, rcinamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopalamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, thebromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycpromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.


Novobiocin

Molecular formula: $C_{31}H_{36}N_2O_{11}$

Molecular weight: 612.63

CAS Registry No.: 303-81-1,
1476-53-5 (Na salt), 4309-70-0 (Ca salt)

Merck Index: 6818



SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 1 mL buffer + 4.5 mL MTBE:isopropanol 97.5:2.5, shake at 70 strokes/min for 5 min, centrifuge at 10-15° at 1100 g for 10 min, repeat extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 20-25°, reconstitute the residue in 0.5 mL mobile phase, vortex, sonicate for 5 min, inject a 100 µL aliquot. (Buffer was 400 mL 1 M KH_2PO_4 , 400 mL K_2HPO_4 solution, and 250 g KCl, adjust pH to 6.5.)

HPLC VARIABLES

Column: 100 × 8 5 µm Nova-Pak C18 Radial-Pak

Mobile phase: MeOH:2-methoxyethanol:buffer 80:5:15 (Buffer was 4.33 g sodium lauryl sulfate and 2 mL 1 M phosphoric acid in 150 mL water, pH 2.8.)

Flow rate: 1.8

Injection volume: 100

Detector: UV 330

CHROMATOGRAM

Retention time: 3.0

Internal standard: novobiocin

OTHER SUBSTANCES

Extracted: coumermycin A1

KEY WORDS

dog; plasma; novobiocin is IS

REFERENCE

Strojny,N.; Conzentino,P.; de Silva,J.A. Determination of coumermycin A1 in plasma by reversed-phase high-performance liquid chromatographic analysis, *J.Chromatogr.*, **1985**, 342, 145–158.

SAMPLE

Matrix: blood

Sample preparation: 0.5-1 mL Plasma + 10 μ L 5 mM prednisone in MeOH + 5 mL MeOH, vortex for 10 s, centrifuge at 2000 g for 10 min, inject a 15 μ L aliquot of the supernatant.

HPLC VARIABLES

Guard column: CN Guard-Pak

Column: 150 \times 4.6 5 μ m APEX octyl EC C8 (Jones Chromatography)

Mobile phase: Gradient. MeOH:buffer from 55:45 to 20:80 over 20 min. (Buffer was water acidified to pH 3.0 with trifluoroacetic acid.)

Flow rate: 1.5

Injection volume: 15

Detector: UV 254

CHROMATOGRAM

Retention time: 18.1

Internal standard: prednisone (5.5)

Limit of detection: 5 μ M

KEY WORDS

plasma

REFERENCE

Chen,T.-L.; Kennedy,M.J.; Dunlap,V.M.; Colvin,O.M. Determination of plasma novobiocin levels by a reversed-phase high-performance liquid chromatographic assay, *J.Chromatogr.B*, **1994**, 652, 109–113.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Serum + 5 μ L 4 mM mitomycin C in MeCN + 400 μ L MeCN, vortex, centrifuge at 12000 g for 2 min, inject a 50 μ L aliquot of the supernatant.

HPLC VARIABLES

Guard column: 15 \times 3.2 NewGuard RP18

Column: 220 \times 4.6 5 μ m Spheri-5 RP18

Mobile phase: MeCN:10 mM phosphoric acid 80:20

Flow rate: 0.8

Injection volume: 50

Detector: UV 340

CHROMATOGRAM

Retention time: 4.53

Internal standard: mitomycin C (6.66)

Limit of quantitation: 1 μ g/mL

KEY WORDS

serum

REFERENCE

Zuhowski,E.G.; Gutheil,J.C.; Egorin,M.J. Rapid and sensitive high-performance liquid chromatographic assay for novobiocin in human serum, *J.Chromatogr.B*, **1994**, 655, 147–152.

SAMPLE

Matrix: blood, milk, tissue

Sample preparation: Blend 10 g tissue with 30 mL 200 mM $(\text{NH}_4)_2\text{H}_2\text{PO}_4$. Dilute each 1 mL of milk or serum with 3 mL 200 mM $(\text{NH}_4)_2\text{H}_2\text{PO}_4$. Add 10 mL Tissue homogenate, diluted milk, or diluted serum to 10 (muscle, milk, serum) or 20 (liver, kidney) mL MeOH, swirl vigorously, let stand for 5 min, filter (paper), refilter if not clear, to each 3 mL of filtrate from liver or kidney add 1 mL water, inject a 200 μL aliquot.

HPLC VARIABLES

Guard column: Supelcosil LC-18-DB

Column: 150 \times 4.6 Supelcosil LC-18-DB

Mobile phase: Gradient. MeCN:MeOH:10 mM phosphoric acid 0:50:50 for 1 min, to 80:0:20 over 19 min, return to initial conditions. (At the end of each day flush the column with the final mobile phase, store in this mobile phase.)

Flow rate: 1 (at end of run set flow rate to 2 mL/min for 5 min then 1 mL/min for 5 min)

Injection volume: 200

Detector: UV 340

CHROMATOGRAM

Retention time: 21

Limit of detection: 10 ppb

KEY WORDS

serum; cow; muscle; liver; kidney

REFERENCE

Moats,W.A.; Leskinen,L. Determination of novobiocin residues in milk, blood, and tissues by liquid chromatography, *J.Assoc.Off.Anal.Chem.*, **1988**, 71, 776–778.

SAMPLE

Matrix: formulations

Sample preparation: Add 10 mL THF to the peanut oil formulation, make up to 100 mL with mobile phase, sonicate, shake at high speed for 5 min, centrifuge at 2000 g for 5 min, inject a 20 μL aliquot of the clear supernatant.

HPLC VARIABLES

Column: 250 \times 4.6 5 μm LiChrosorb SI-100

Mobile phase: Butyl chloride:THF:MeOH:acetic acid 88:5:4:3 (Butyl chloride was 50% water saturated prepared by mixing equal volumes of water-saturated butyl chloride and anhydrous butyl chloride.)

Flow rate: 1

Injection volume: 20

Detector: UV 340

CHROMATOGRAM

Retention time: 16

Internal standard: prednisone (UV 254) (10)

Limit of detection: 10 ng

OTHER SUBSTANCES

Simultaneous: impurities, degradation products

KEY WORDS

stability-indicating; normal phase; oils

REFERENCE

Tsuji,K.; Rahn,P.D.; Kane,M.P. High-performance liquid chromatographic method for the determination of novobiocin, *J.Chromatogr.*, **1982**, 235, 205-214.

SAMPLE

Matrix: milk

Sample preparation: 10 g Milk + 30 mL 200 mM ammonium phosphate, shake. Remove a 10 mL aliquot and add it to 10 mL MeOH, swirl vigorously, swirl occasionally for 5 min, filter (paper), refilter if necessary, inject a 1 mL aliquot.

HPLC VARIABLES

Guard column: 20 × 4.6 5 µm Supelcosil LC-18-DB

Column: 150 × 4.6 5 µm Supelcosil LC-18-DB

Mobile phase: Gradient. MeOH:5 mM phosphoric acid 50:50 for 1 min, to MeCN:MeOH:5 mM phosphoric acid 80:0:20 over 20 min, maintain at MeCN:MeOH:5 mM phosphoric acid 80:0:20 for 5 min, re-equilibrate at initial conditions for 6 min.

Flow rate: 1 (re-equilibrate at 2 mL/min for 5 min and 1 mL/min for 1 min)

Injection volume: 1000

Detector: UV 340

CHROMATOGRAM

Retention time: 23-25

Limit of detection: <0.05 ppm

KEY WORDS

cow

REFERENCE

Reeves,V.B. Liquid chromatographic procedure for the determination of novobiocin residues in bovine milk: Interlaboratory study, *JAOAC Int.*, **1995**, 78, 55-58.

Nylidrin

Molecular formula: C₁₉H₂₅NO₂

Molecular weight: 299.41

CAS Registry No.: 447-41-6, 849-55-8 (HCl)

Merck Index: 6830

Lednicer No.: 1 69

SAMPLE

Matrix: blood

Sample preparation: 500 µL Plasma + 500 µL buffer, mix briefly, add 50-75 mg resin, mix at 10-25 rpm for 30 min, discard the supernatant, wash twice with 1 mL buffer, add 500 µL 50 mg/mL KOH in MeOH:water 50:50, mix for 30 min, inject a 20 µL aliquot of the eluate. (Buffer was 100 mM citric acid:200 mM Na₂HPO₄ 29:71, pH 6.5 (McIlvaine buffer) Wash 20-50 mesh Dowex HCR-S resin twice with water and allow it to equilibrate in buffer).

HPLC VARIABLES

Guard column: 25 × 4.6 5 µm Spherisorb ODS-I

Column: 250 × 4.6 5 µm Spherisorb ODS-I

Mobile phase: MeCN:MeOH:buffer 30:18:52 containing 1.8 mM octanesulfonic acid (Buffer was 30 mM KH₂PO₄ adjusted to pH 3.0 with concentrated orthophosphoric acid.)


Flow rate: 1

Injection volume: 20

Detector: E, Gynotek M20, glassy carbon working electrode 950 mV, Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 11



Internal standard: nylidrin

OTHER SUBSTANCES

Extracted: isoxsuprine

KEY WORDS

horse; plasma; nylidrin is IS; SPE

REFERENCE

Hashem, A.; Lubczyk, B. Determination of isoxsuprine in equine plasma by high-performance liquid chromatography with electrochemical detection, *J. Chromatogr.*, **1991**, 563, 216–223.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, am-triptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenpropofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephentyol, mephesisin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methypylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phenacyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, pumycin, pyrilamine, pyrrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, re-cinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sufadiazine, sulfadimethoxine, sulfaethi-dole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasox-

izole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripelennamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.


Nystatin

Molecular formula: C₄₇H₇₅NO₁₇

Molecular weight: 926.11

CAS Registry No.: 1400-61-9

Merck Index: 6834



SAMPLE

Matrix: CSF

Sample preparation: Condition a BakerBond C18 SPE cartridge with 3 mL MeOH and 3 mL 100 mM pH 9 carbonate buffer. Add 1 mL CSF to the SPE cartridge, wash with 2 mL 100 mM pH 9 carbonate buffer, air dry for 2 min, elute with two 500 µL aliquots of MeOH. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute with 200 µL MeOH, inject a 100 µL aliquot.

HPLC VARIABLES

Column: 150 × 3.9 4 µm Nova-Pak C18

Mobile phase: MeCN:10 mM pH 5 EDTA 35:65

Flow rate: 0.5

Injection volume: 100

Detector: UV 410

CHROMATOGRAM

Retention time: 8.5

Internal standard: nystatin

OTHER SUBSTANCES

Extracted: amphotericin B

KEY WORDS

dog; human; SPE; nystatin is IS

REFERENCE

Liu,H.; Davoudi,H.; Last,T. Determination of Amphotericin B in cerebrospinal fluid by solid-phase extraction and liquid chromatography, *J.Pharm.Biomed.Anal.*, **1995**, *13*, 1395-1400.

SAMPLE

Matrix: bulk

HPLC VARIABLES

Column: 500 × 1 7-8 µm Zorbax B.P. Sil

Mobile phase: MeOH:DMF:water:acetic acid 72:25:3:0.4

Flow rate: 0.125
Injection volume: 1
Detector: UV 308

CHROMATOGRAM

Retention time: 13

KEY WORDS

microbore; normal phase

REFERENCE

Milhaud,J.; Gareil,P.; Rosset,R. Separation of filipin and nystatin complexes by semi-preparative and microbore high-performance liquid chromatography, *J.Chromatogr.*, **1986**, 358, 284–287.

SAMPLE

Matrix: bulk

Sample preparation: Make up a 10% solution in DMF then dilute to 20 mg/mL with MeOH: water 66:34, inject a 200 μ L aliquot.

HPLC VARIABLES

Column: 500 \times 9 10 μ m Partisil ODS 2

Mobile phase: MeOH:water:acetic acid 66:33:1

Flow rate: 8

Injection volume: 200

Detector: UV 308

CHROMATOGRAM

Retention time: 22

KEY WORDS

semi-preparative

REFERENCE

Milhaud,J.; Gareil,P.; Rosset,R. Separation of filipin and nystatin complexes by semi-preparative and microbore high-performance liquid chromatography, *J.Chromatogr.*, **1986**, 358, 284–287.

SAMPLE

Matrix: bulk

HPLC VARIABLES

Column: 300 \times 4.5 μ Bondapak C18

Mobile phase: MeOH:water 70:30

Flow rate: 2

Detector: UV 280

CHROMATOGRAM

Retention time: 10

REFERENCE

Mehta,R.T.; Hopfer,R.L.; Gunner,L.A.; Juliano,R.L.; Lopez-Berestein,G. Formulation, toxicity, and antifungal activity in vitro of liposome-encapsulated nystatin as therapeutic agent for systemic candidiasis, *Antimicrob.Agents Chemother.*, **1987**, 31, 1897–1900.

SAMPLE

Matrix: bulk

Sample preparation: Dissolve in mobile phase to a concentration of 1 mg/mL.

HPLC VARIABLES

Column: 250 \times 4 7 μ m LiChrosorb RP-18

Mobile phase: MeOH:water 90:10 + 1% formic acid, pH 3.5

Flow rate: 1

Injection volume: 10

Detector: UV 340

CHROMATOGRAM

Retention time: 6.3

OTHER SUBSTANCES

Simultaneous: impurities

REFERENCE

Sauer,B.; Matusch,R. High-performance liquid chromatographic separations of nystatin and their influence on the antifungal activity, *J.Chromatogr.A*, **1994**, 672, 247-253.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4 5 µm ODS-Hypersil

Mobile phase: MeOH:DMF:10 mM pH 7.0 Tris in water 56:9.6:34.4

Detector: UV 305

OTHER SUBSTANCES

Simultaneous: degradation products

REFERENCE

Egodage,K.L.; Haslam,J.S.; Rajewski,R.A.; Stella,V.J. Correlation and validation to the USP bioassay of a RP-HPLC assay for nystatin (Abstract 3373), *Pharm.Res.*, **1997**, 14, S587.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in DMSO to 10 mg/mL, dilute 1:20 with MeOH.

HPLC VARIABLES

Column: 250 × 4.6 10 µm µBondapak C18

Mobile phase: MeCN:50 mM phosphate buffer (pH 3.5-8.1) 30:70 to 35:65

Flow rate: 0.4-2

Detector: UV 313

OTHER SUBSTANCES

Simultaneous: amphotericin A

REFERENCE

Aszalos,A.; Bax,A.; Burlinson,N.; Roller,P.; McNeal,C. Physico-chemical and microbiological comparison of nystatin, amphotericin A and amphotericin B, and structure of amphotericin A, *J.Antibiot.(Tokyo)*, **1985**, 38, 1699-1713.


Octhilinone

Molecular formula: C₁₁H₁₉NOS

Molecular weight: 213.34

CAS Registry No.: 26530-20-1

Merck Index: 6853



SAMPLE

Matrix: blood, urine