

Mobile phase: MeOH:water 90:10 + 1% formic acid, pH 3.5

Flow rate: 1

Injection volume: 10

Detector: UV 340

CHROMATOGRAM

Retention time: 6.3

OTHER SUBSTANCES

Simultaneous: impurities

REFERENCE

Sauer,B.; Matusch,R. High-performance liquid chromatographic separations of nystatin and their influence on the antifungal activity, *J.Chromatogr.A*, **1994**, 672, 247-253.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4 5 µm ODS-Hypersil

Mobile phase: MeOH:DMF:10 mM pH 7.0 Tris in water 56:9.6:34.4

Detector: UV 305

OTHER SUBSTANCES

Simultaneous: degradation products

REFERENCE

Egodage,K.L.; Haslam,J.S.; Rajewski,R.A.; Stella,V.J. Correlation and validation to the USP bioassay of a RP-HPLC assay for nystatin (Abstract 3373), *Pharm.Res.*, **1997**, 14, S587.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in DMSO to 10 mg/mL, dilute 1:20 with MeOH.

HPLC VARIABLES

Column: 250 × 4.6 10 µm µBondapak C18

Mobile phase: MeCN:50 mM phosphate buffer (pH 3.5-8.1) 30:70 to 35:65

Flow rate: 0.4-2

Detector: UV 313

OTHER SUBSTANCES

Simultaneous: amphotericin A

REFERENCE

Aszalos,A.; Bax,A.; Burlinson,N.; Roller,P.; McNeal,C. Physico-chemical and microbiological comparison of nystatin, amphotericin A and amphotericin B, and structure of amphotericin A, *J.Antibiot.(Tokyo)*, **1985**, 38, 1699-1713.


Octhilinone

Molecular formula: C₁₁H₁₉NOS

Molecular weight: 213.34

CAS Registry No.: 26530-20-1

Merck Index: 6853



SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 3.112

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

Octopamine


Molecular formula: C₈H₁₁NO₂

Molecular weight: 153.18

CAS Registry No.: 104-14-3, 876-04-0 (D-(-)), 770-05-8 (DL HCl)

Merck Index: 6856

Lednicer No.: 5 23



SAMPLE

Matrix: solutions

Sample preparation: Inject an aliquot of a 100 µg/mL solution in mobile phase.

HPLC VARIABLES

Column: 150 × 4 5 µm Crownpak CR(+) immobilized crown ether

Mobile phase: 0.1% pH 1.9 Perchloric acid

Column temperature: 25

Flow rate: 1

Detector: UV 210

CHROMATOGRAM

Retention time: 4.64, 4.95

OTHER SUBSTANCES

Simultaneous: norepinephrine


KEY WORDS

chiral; comparison with capillary electrophoresis

REFERENCE

Nishi,H.; Nakamura,K.; Nakai,H.; Sato,T. Separation of enantiomers and isomers of amino compounds by capillary electrophoresis and high-performance liquid chromatography utilizing crown ethers, *J.Chromatogr.A*, **1997**, 757, 225-235.

Octreotide

Molecular formula: $C_{49}H_{66}N_{10}O_{10}S_2$ **Molecular weight:** 1019.26**CAS Registry No.:** 83150-76-9**Merck Index:** 6859**SAMPLE****Matrix:** bile, blood, feces, urine

Sample preparation: Hydrolyze pooled plasma with 5 mg Subtilisin A at 50° for 1 h. Lyophilize feces and extract twice with 4 volumes MeOH. Lyophilize urine and bile and reconstitute with water.

HPLC VARIABLES**Guard column:** 30 × 4.6 5 μm RP-18 (Brownlee)**Column:** 100 × 4.6 5 μm RP-18 (Brownlee)

Mobile phase: Gradient. MeCN:water:trifluoroacetic acid from 0:99.8:0.2 to 99.8:0:0.2 over 40 min.

Flow rate: 1.5**Detector:** UV 210 or radioactivity**CHROMATOGRAM****Retention time:** 45**KEY WORDS**

rat; plasma

REFERENCE

Lemaire,M.; Azria,M.; Dannecker,R.; Marbach,P.; Schweitzer,A.; Maurer,G. Disposition of sandostatin, a new synthetic somatostatin analogue, in rats, *Drug Metab.Dispos.*, **1989**, 17, 699-703.

SAMPLE**Matrix:** formulations**Sample preparation:** Inject a 20 μL aliquot.**HPLC VARIABLES****Column:** 250 × 2.6 5 μm Bakerbond C18

Mobile phase: Gradient. A was MeCN:water:1 M tetramethylammonium hydroxide 10:88:2 adjusted to pH 4.5 with concentrated orthophosphoric acid. B was MeCN:water:1 M tetramethylammonium hydroxide 60:38:2 adjusted to pH 4.5 with concentrated orthophosphoric acid. A:B from 100:0 to 0:100 over 14 min.

Flow rate: 1.3**Injection volume:** 20**Detector:** UV 210

CHROMATOGRAM**Retention time:** 11.5

OTHER SUBSTANCES**Simultaneous:** degradation products, des-threninol

KEY WORDSprotect from light; injections

REFERENCE

Stiles,M.L.; Allen,L.V.,Jr.; Resztak,K.E.; Prince,S.J. Stability of octreotide acetate in polypropylene syringes, *Am.J.Hosp.Pharm.*, **1993**, 50, 2356–2358.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 250 × 4.6 5 µm Spherisorb RP18**Mobile phase:** MeCN:buffer 31:69, pH 4.7 (Buffer was 20 mM tetramethylammonium hydroxide adjusted to pH 4.7 with phosphoric acid.)**Column temperature:** 40**Flow rate:** 1.2**Injection volume:** 20**Detector:** UV 210

CHROMATOGRAM**Retention time:** 18

OTHER SUBSTANCES**Simultaneous:** derivatives

REFERENCE

Kuhn,R.; Morin,C.; Erni,F. A simple model describing the retention behavior of octreotide and its glycosylated derivatives in reversed phase HPLC, *Chromatographia*, **1995**, 41, 516–520.

SAMPLE**Matrix:** solutions**Sample preparation:** Inject a 30 µL aliquot.

HPLC VARIABLES**Column:** 250 × 4.6 Spheri-5 C18**Mobile phase:** MeCN:water:1 M tetramethylammonium hydroxide pentahydrate 33:65:2, adjusted to pH 4.5 with phosphoric acid**Flow rate:** 0.8**Injection volume:** 30**Detector:** UV 280

CHROMATOGRAM**Retention time:** 6.2

OTHER SUBSTANCES**Simultaneous:** degradation products (UV 210)

KEY WORDSinjections

REFERENCE

Ripley,R.G.; Ritchie,D.J.; Holstad,S.G. Stability of octreotide acetate in polypropylene syringes at 5 and -20°C, *Am.J.Health-Syst.Pharm.*, **1995**, 52, 1910–1911.


Octyl methoxycinnamate

Molecular formula: $C_{18}H_{26}O_3$

Molecular weight: 290.40

CAS Registry No.: 5466-77-3

Merck Index: 6864



SAMPLE

Matrix: formulations

Sample preparation: 1-1.5 g sun-screen lotion + 50 mL isopropanol, dissolve. Remove a 5 mL aliquot and make up to 50 mL with mobile phase, filter (0.45 μ m) inject an aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 10 μ m C8 Hypersil

Mobile phase: Isopropanol:buffer 10:90 (Buffer was 100 mM sodium dodecyl sulfate (electrophoresis grade) containing 0.3% triethylamine adjusted to pH 3.0 with phosphoric acid.)

Flow rate: 1.5

Injection volume: 20

Detector: UV 254, UV 300

CHROMATOGRAM

Retention time: 19.40

OTHER SUBSTANCES

Simultaneous: 2-ethylhexyl p-dimethylaminobenzoate, oxybenzone, methyl paraben, propyl paraben

KEY WORDS

lotion

REFERENCE

Tomasella,F.P.; Zuting,P.; Love,L.J. Determination of sun-screen agents in cosmetic products by micellar liquid chromatography, *J.Chromatogr.*, **1991**, 587, 325-328.

SAMPLE

Matrix: formulations

Sample preparation: Weigh out lotion, add 20 mL EtOH, heat to 60°, stir for 30 min at room temperature, make up to 25 mL, stir for 5 min, centrifuge at 14500 g for 5 min, inject an aliquot of the supernatant.

HPLC VARIABLES

Column: 200 \times 4.6 5 μ m Ultrasphere C8

Mobile phase: Gradient. MeOH:1% aqueous acetic acid from 80:20 to 100:0 over 10 min, maintain at 100:0 for 2 min, re-equilibrate for 4 min.

Column temperature: 25

Flow rate: 1

Detector: UV 325

CHROMATOGRAM

Retention time: 8.7

OTHER SUBSTANCES

Simultaneous: benzophenone-3, butylmethoxydibenzoylmethane, impurities

KEY WORDS

lotion; sunscreen; an isomer forms on exposure to light

REFERENCE


Meijer, J.; Lodén, M. Stability analysis of three UV-filters using HPLC, *J. Liq. Chromatogr.*, **1995**, *18*, 1821–1832.

Odapipam

Molecular formula: C₁₉H₂₀ClNO₂

Molecular weight: 329.83

CAS Registry No.: 131796-63-9



SAMPLE

Matrix: microsomal incubations

Sample preparation: Condition a 500 mg C18 Bond-Elut SPE cartridge with MeOH and water. Mix 2 mL microsomal incubation with 4.5 mL ice-cold MeOH, centrifuge at 4000 g for 10 min, dilute the supernatant with 15 mL water. Add a 2 mL aliquot to the SPE cartridge, wash with 2 mL water, elute with 2 mL MeOH:25% aqueous ammonia 96:4, evaporate the eluate to dryness under reduced pressure, reconstitute with mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 100 × 4.6 3 µm ChromSpher Si (Chrompack)

Mobile phase: n-Heptane:2-propanol:water:25% ammonia 90:10:0.2:0.1

Detector: UV 280; MS, VG TRIO 1000, particle beam interface at 50°, helium at 25–30 psi, ion source 200°, positive ionization mode, electron current 150 µA, electron energy 70 eV

CHROMATOGRAM

Retention time: 2.3

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

normal phase; rat; liver; SPE

REFERENCE

Andersen, J.V.; Hansen, K.T. Normal-phase liquid chromatography-particle-beam mass spectrometry in drug metabolism studies of the dopamine receptor antagonist Odapipam and the muscarine M1 receptor agonist Xanomeline, *Xenobiotica*, **1997**, *27*, 901–912.

Ofloxacin


Molecular formula: C₁₈H₂₀FN₃O₄

Molecular weight: 361.37

CAS Registry No.: 82419-36-1

Merck Index: 6865

Lednicer No.: 4 141-145



SAMPLE

Matrix: aqueous humor, blood

Sample preparation: Aqueous humor. Inject a 10 µL aliquot directly. Plasma. Condition a 3 mL C18 SPE cartridge (Varian) with two 3 mL portions of MeCN and 3 mL buffer. Add 2 mL 625 ng/mL ciprofloxacin in buffer to 500 µL of plasma, mix, add to the SPE cartridge. Wash with 3 mL buffer. Remove moisture with vacuum (200 mbar) for 10 min. Elute with two 500 µL

portions of MeCN:buffer 40:60. Vortex the eluate, inject a 10 μ L aliquot. (Buffer was 100 mM Tris adjusted to pH 5.0 with HCl).

HPLC VARIABLES

Column: 300 \times 4.6 5 μ m endcapped ODS-Hypersil

Mobile phase: MeCN:DMF:10 mM NaH₂PO₄ 15:6:79, adjusted to pH 3.0 with 85% phosphoric acid

Flow rate: 1

Injection volume: 10

Detector: UV 285

CHROMATOGRAM

Retention time: 10.2

Internal standard: ciprofloxacin (12.0)

Limit of detection: 80 ng/mL (aqueous humor), 310 ng/mL (plasma)

OTHER SUBSTANCES

Extracted: cefotaxime

KEY WORDS

plasma; SPE

REFERENCE

Kraemer,H.-J.; Gehrke,R.; Breithaupt,A.; Breithaupt,H. Simultaneous quantification of cefotaxime, desacetyl-cefotaxime, ofloxacin and ciprofloxacin in ocular aqueous humor and in plasma by high-performance liquid chromatography, *J.Chromatogr.B*, **1997**, *700*, 147–153.

SAMPLE

Matrix: bile, blood, perfusate

Sample preparation: Intestinal perfusate. Centrifuge before analysis. Plasma. Mix 200 μ L plasma with 200 μ L 100 mM pH 6.8 phosphate buffer, add 4 mL dichloromethane, shake at 100 cycles/min for 10 min. Remove the organic layer and dry it under nitrogen at 40°, reconstitute with 200 μ L mobile phase, inject a 20 μ L aliquot (*J. Chromatogr.* 1988, 434, 320). Bile. Mix 100 μ L bile with 900 μ L pH 7.0 phosphate buffer, add 500 μ L MeOH. Shake for 1 min and centrifuge at 1000 g for 5 min. Inject a 20 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: Bovine Serum Albumin (Macherey Nagel)

Mobile phase: 200 mM Potassium dihydrogen phosphate containing 5 mM N,N-dimethyloctylamine, adjusted to pH 8.0 with KOH

Flow rate: 1

Injection volume: 20

Detector: F ex 298 em458

CHROMATOGRAM

Retention time: 6 (S(-)), 7.8 (R(+))

Limit of detection: 100 ng/mL

OTHER SUBSTANCES

Also analyzed: cefoperazone, ciprofloxacin, ofloxacin, quinidine, verapamil

KEY WORDS

chiral; intestinal efflux; pharmacokinetics; plasma

REFERENCE

Rabbaa,L.; Dautrey,S.; Colas-Linhart,N.; Carbon,C.; Farinotti,R. Intestinal elimination of ofloxacin enantiomers in the rat: Evidence of a carrier-mediated process, *Antimicrob.Agents Chemother.*, **1996**, *40*, 2126–2130.

SAMPLE

Matrix: blood

Sample preparation: Mix 1 mL blood with 1 mL 0.25 mM Triton, vortex for 30 s, add 4 mL 6% trichloroacetic acid. Vortex for 30 s, centrifuge at 2000 g for 10 min, inject a 100 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 10 μ m C18

Mobile phase: MeCN:buffer 10:90 (Buffer was 1 L 25 mM phosphoric acid and 15 mL 40% tetrabutyl ammonium hydrogen sulfate, adjusted to pH 3.0 with 66.6 mM phosphate buffer.)

Flow rate: 2

Injection volume: 100

Detector: F ex 330 em 450

CHROMATOGRAM

Retention time: 5

Limit of detection: 30 ng/mL (plasma), 50 ng/mL (blood)

KEY WORDS

plasma; pharmacokinetics; rabbit

REFERENCE

Colino,C.I.; Garc a Turi o,A.; Sanchez Navarro,A.; Lanao,J.M. A comparative study of ofloxacin and ciprofloxacin erythrocyte distribution, *Biopharm.Drug Dispos.*, **1998**, 19, 71–77.

SAMPLE

Matrix: blood

Sample preparation: Mix equal volumes serum and MeOH. Centrifuge sample, inject a 20 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 250 \times 4 Spherisorb 5 ODS II

Mobile phase: MeCN:water:phosphoric acid 5:94.6:0.32 adjusted to pH 3.0 with tetrabutylammonium hydroxide

Column temperature: 50

Flow rate: 2

Injection volume: 20

Detector: F ex 310 em 489

CHROMATOGRAM

Limit of detection: 80 ng/mL

KEY WORDS

serum; pharmacokinetics

REFERENCE

Bethell,D.B.; Day,N.P.; Dung,N.M.; McMullin,C.; Loan,H.T.; Tam,D.T.; Minh,L.T.; Linh,N.T.; Dung,N.Q.; Vinh,H.; MacGowan,A.P.; White,L.O.; White,N.J. Pharmacokinetics of oral and intravenous ofloxacin in children with multidrug-resistant typhoid fever, *Antimicrob.Agents Chemother.*, **1996**, 40, 2167–2172.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Whole blood or 500 μ L plasma, vortex for 30 s with 1 mL 250 μ M Triton, add 4 mL 6% trichloroacetic acid. Vortex for 30 s, centrifuge at 2000 g for 10 min. Inject a 100 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 250 \times 4.6 10 μ m C18

Mobile phase: MeCN:buffer 10:90 (Buffer was mixture of 1 L 25 mM phosphoric acid and 15 mL 40% tetrabutylammonium hydrogen sulfate, adjusted to pH 3 with 66.6 mM phosphate buffer.)

Flow rate: 2

Injection volume: 100

Detector: F ex 330 em 450

CHROMATOGRAM

Retention time: 5

Limit of detection: 30 ng/mL (plasma), 50 ng/mL (blood)

KEY WORDS

plasma; whole blood; pharmacokinetics

REFERENCE

Colino,C.I.; García Turiño,A.; Sanchez Navarro,A.; Lanao,J.M. A comparative study of ofloxacin and ciprofloxacin erythrocyte distribution, *Biopharm.Drug Dispos.*, **1998**, 19, 71–77.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Blood or 500 μ L plasma, vortex for 30 s with 1 mL 250 μ M Triton, add 4 mL 6% trichloroacetic acid. Vortex for 30 s, centrifuge at 2000 g for 10 min. Inject a 100 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 250 \times 4.6 10 μ m C18

Mobile phase: MeCN:buffer 10:90 (Buffer was mixture of 1 L 25 mM phosphoric acid and 15 mL 40% tetrabutylammonium hydrogen sulfate, adjusted to pH 3 with 66.6 mM phosphate buffer.)

Flow rate: 2

Injection volume: 100

Detector: F ex 277 em 445

CHROMATOGRAM

Retention time: 9

Limit of detection: 30 ng/mL (plasma), 50 ng/mL (blood)

KEY WORDS

plasma; whole blood; pharmacokinetics

REFERENCE

Colino,C.I.; García Turiño,A.; Sanchez Navarro,A.; Lanao,J.M. A comparative study of ofloxacin and ciprofloxacin erythrocyte distribution, *Biopharm.Drug Dispos.*, **1998**, 19, 71–77.

SAMPLE

Matrix: blood

Sample preparation: Filter 1 mL plasma using a micropartition system (MPS-1, Amicon, MA) while centrifuging at 2000 g for 20 min at 10°, inject an aliquot of the ultrafiltrate.

HPLC VARIABLES

Column: 250 \times 4.6 Spherisorb ODS-2 endcapped

Mobile phase: MeCN:buffer 13:87 containing 5 mM tetrabutylammonium sulfate, adjusted to pH 2.5 with 1 M NaOH (Buffer was 100 mM citric acid containing 200 mM ammonium perchlorate.)

Column temperature: 37

Flow rate: 1

Detector: UV 295

CHROMATOGRAM

Retention time: 8.54

Internal standard: pipemic acid (4.67)

KEY WORDS

plasma; ultrafiltrate

REFERENCE

Zlotos,G.; Bucker,A.; Kinzig-Schippers,M.; Sorgel,F.; Holzgrabe,U. Plasma protein binding of gyrase inhibitors, *J.Pharm.Sci.*, **1998**, 87, 215–220.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 295

CHROMATOGRAM

Retention time: 8.648

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

SAMPLE

Matrix: formulations

Sample preparation: Dilute a 25 mg/mL levofloxacin injection with an infusion fluid to a concentration of 100 μ g/mL. Inject a 20 μ L aliquot. (Infusion fluids were 0.9% NaCl, 5% dextrose, 5% dextrose and 0.9% NaCl, 5% dextrose and lactated Ringer's injections, 5% sodium bicarbonate, Plasma-Lyte 56 and 5% dextrose, 5% dextrose, 0.45% NaCl, and 0.15% KCl, 1/6 M sodium lactate, sterile water, and 20% mannitol.)

HPLC VARIABLES

Column: 250 \times 4.6 Zorbax SB-Phenyl

Mobile phase: MeCN:MeOH:94 mM KH_2PO_4 :trifluoroacetic acid 15:5:80:0.3

Column temperature: 45

Flow rate: 1

Injection volume: 25

Detector: UV 294

CHROMATOGRAM

Retention time: 15

OTHER SUBSTANCES

Simultaneous: degradation products

KEY WORDS

injections; stability indicating; for levofloxacin

REFERENCE

Williams,N.A.; Bornstein,M.; Johnson,K. Stability of levofloxacin in intravenous solutions in polyvinyl chloride bags, *Am.J.Health-Syst.Pharm.*, **1996**, 53, 2309-2313.

SAMPLE

Matrix: growth medium

Sample preparation: 500 μ L Sample + 500 μ L 100 μ g/mL IS in cold (4°) MeCN, vortex, centrifuge at 3000 g for 5 min. Remove a 500 μ L aliquot of the supernatant, filter (0.45 μ m Acrodisc syringe filter), inject a 30 μ L aliquot. (Protect all specimens from light.)

HPLC VARIABLES

Guard column: C18 5U (Alltech)

Column: 150 \times 4.6 7 μ m Adsorbosphere HS C18 7U

Mobile phase: MeCN:20 mM pH 3.0 phosphate buffer 35:65 containing 0.2% triethylamine and 0.2% sodium dodecyl sulfate, adjusted to pH 3.0 with 85% phosphoric acid

Flow rate: 1.75

Injection volume: 30

Detector: UV 280

CHROMATOGRAM

Retention time: 4.15

Internal standard: sparfloxacin (7.09)

Limit of quantitation: 100 ng/mL

OTHER SUBSTANCES

Also analyzed: ciprofloxacin, clinafloxacin, levofloxacin, sparfloxacin, temafloxacin, trovafloxacin

KEY WORDS

Mueller-Hinton broth

REFERENCE

Wright,D.H.; Herman,V.K.; Konstantinides,F.N.; Rotschafer,J.C. Determination of quinolone antibiotics in growth media by reversed-phase high-performance liquid chromatography, *J.Chromatogr.B*, **1998**, 709, 97–104.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 4 μ m NovaPak C18

Mobile phase: MeCN:MeOH:buffer:acetic acid 2.5:10:86.5:1 containing 20 mM triethylamine (The pH 2.7 buffer was 0.4% diammonium hydrogen phosphate in water containing 0.4% (?) tetrabutylammonium hydrogen sulfate.)

Flow rate: 1

Detector: UV 279

CHROMATOGRAM

Retention time: 11.6

OTHER SUBSTANCES

Extracted: ciprofloxacin, enrofloxacin

REFERENCE

Cester,C.C.; Toutain,P.L. A comprehensive model for enrofloxacin to ciprofloxacin transformation and disposition in dog, *J.Pharm.Sci.*, **1997**, 86, 1148–1155.

SAMPLE

Matrix: solutions

Sample preparation: Filter (0.45 μ m) a solution in MeCN:water 10:90, inject an aliquot of the filtrate.

HPLC VARIABLES

Column: 250 \times 4 5 μ m LiChrospher 100 RP-18

Mobile phase: MeCN:buffer 7:93 (Buffer was 25 mM phosphoric acid adjusted to pH 3.89 with 100 mM tetrabutylammonium hydroxide.)

Flow rate: 1

Injection volume: 10

Detector: UV 295

CHROMATOGRAM

Retention time: 8.8

OTHER SUBSTANCES

Simultaneous: ciprofloxacin (UV 280), enoxacin (UV 280), fleroxacin (UV 280), norfloxacin (UV 280), pipemidic acid (UV 280)

REFERENCE

Barbosa,J.; Bergés,R.; Sanz-Nebot,V. Solvatochromic parameter values and pH in aqueous-organic mixtures used in liquid chromatography. Prediction of retention of a series of quinolones, *J.Chromatogr.A*, **1996**, 719, 27–36.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Guard column: 20 mm long Supelguard LC-18S (Supelco)

Column: 250 × 4.6 Suplecasil LC-18S

Mobile phase: MeCN:buffer:water 10:3.5:86.5 (Buffer was 400 mM tetrabutylammonium hydroxide adjusted to pH 2.85.)

Flow rate: 1.8

Detector: UV 280

REFERENCE

Sinko,P.J.; Hu,P. Determining intestinal metabolism and permeability for several compounds in rats. Implications on regional bioavailability in humans, *Pharm.Res.*, **1996**, 13, 108–113.


Olanzapine

Molecular formula: C₁₇H₂₀N₄S

Molecular weight: 312.44

CAS Registry No.: 132539-06-1

Merck Index: 6959



SAMPLE

Matrix: blood

Sample preparation: Condition a Certify mixed bed (RP and ion exchange) SPE cartridge with MeOH and phosphate buffer. Add serum, wash with 1 M acetic acid, wash with MeOH, elute with 3% ammonium hydroxide in ethyl acetate. Evaporate the eluate, reconstitute, inject an aliquot.

HPLC VARIABLES

Column: 100 × 4.6 Microsorb CN

Mobile phase: MeCN:MeOH:50 mM pH 6.5 sodium phosphate buffer 5:28:67

Column temperature: 37

Flow rate: 1.5

CHROMATOGRAM

Retention time: 9.6

Internal standard: clozapine (12.5)

Limit of detection: 1 ng/mL

OTHER SUBSTANCES

Simultaneous: imipramine

Also analyzed: haloperidol, risperidone

Interfering: paroxetine

KEY WORDS

SPE; serum

REFERENCE

Prieto,I.V.; Hoffman,D.W. HPLC monitoring of olanzapine (Abstract 131), *Ther.Drug Monit.*, **1997**, *19*, 580.

SAMPLE

Matrix: blood

Sample preparation: Mix 1 mL Plasma with 100 μ L 100 ng/mL IS in MeCN, add 500 μ L saturated sodium carbonate, mix. Add 7 mL pentane:dichloromethane 85:15, shake for 10 min, centrifuge at 18° for 10 min. Evaporate the supernatant to dryness under a slow stream of nitrogen at 60°. Dissolve the residue in 150 μ L MeCN, inject a 30-120 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Ultrasphere CN

Mobile phase: MeCN:MeOH:130 mM pH 6.8 ammonium acetate 86:6:8

Column temperature: 40

Flow rate: 0.8

Injection volume: 30-120

Detector: E, ESA Coulochem model 5100A, model 5020 guard cell +1 V, model 5011 dual electrode analytical cell, electrode 1 +300 mV, electrode 2 +930 mV

CHROMATOGRAM

Retention time: 14.6

Internal standard: 2-ethyl analog of olanzapine LY170222 (Lilly Research Labs) (13)

Limit of detection: 250 pg/mL

OTHER SUBSTANCES

Simultaneous: acetaminophen, benztropine, chlorpromazine, clonazepam, clozapine, fluphenazine, ibuprofen, lorazepam, perphenazine, pseudoephedrine, risperidone, spiperone, sulpiride, trifluoperazine, trihexyphenidyl

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Aravagiri,M.; Ames,D.; Wirshing,W.C.; Marder,S.R. Plasma level monitoring of olanzapine in patients with schizophrenia: Determination by high-performance liquid chromatography with electrochemical detection, *Ther.Drug Monit.*, **1997**, *19*, 307-313.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 200 μ L Microsomal incubation + 200 μ L cold MeCN, mix, centrifuge in a microfuge at maximum speed for 5 min, inject an aliquot of the supernatant.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m YMC basic column (YMC, USA)

Mobile phase: MeCN:MeOH:75 mM pH 7 sodium phosphate buffer 20:30:50

Column temperature: 40

Flow rate: 1

Detector: E, ESA Coulochem Dual Electrode Detector, guard cell model 5020 + 300 mV, analytical cell model 5014, detector 1 + 0.0 V, detector 2 + 300 mV

CHROMATOGRAM**Retention time:** 12**Internal standard:** LY170222 (2-ethylolanzapine)

OTHER SUBSTANCES**Extracted:** metabolites


KEY WORDS

liver; human

REFERENCE

Ring,B.J.; Catlow,J.; Lindsay,T.J.; Gillespie,T.; Roskos,L.K.; Cerimele,B.J.; Swanson,S.P.; Hamman,M.A.; Wrighton,S.A. Identification of the human cytochromes P450 responsible for the in vitro formation of the major oxidative metabolites of the antipsychotic agent olanzapine, *J.Pharmacol.Exp.Ther.*, **1996**, 276, 658–666.

Olaquinox

Molecular formula: C₁₂H₁₃N₃O₄**Molecular weight:** 263.25**CAS Registry No.:** 23696-28-8**Merck Index:** 6960

SAMPLE**Matrix:** feed

Sample preparation: Grind and sieve feed with a Moulinex blender. Weigh out 10 g and add it to 20 mL DMF and 60 mL carbon tetrachloride, stir magnetically at 500 rpm at 60° for 30 min, cool, filter (100 µm glass), wash the residue with a little carbon tetrachloride. Remove 25 mL of the filtrate and add it to 45 mL water, stir vigorously for 2 min, centrifuge at 320 g for 5 min, inject an aliquot of the aqueous layer.

HPLC VARIABLES**Column:** 250 × 4.1 10 µm Versapack C18 (Alltech)**Mobile phase:** Gradient. MeOH:water 15:85 for 4 min, to 50:50 over 2 min, maintain at 50:50 for 4 min, return to initial conditions over 2 min.**Flow rate:** 1.5**Injection volume:** 20**Detector:** UV 262

CHROMATOGRAM**Retention time:** 4

OTHER SUBSTANCES**Extracted:** carbadox (UV 305 nm)

KEY WORDS

protect from light

REFERENCE

dos Ramos,F.J.; da Silveira,I.N.; de Graaf,G. Column liquid chromatographic determination of carbadox and olaquinox in feeds, *J.Chromatogr.*, **1991**, 558, 125–130.


Oleandomycin

Molecular formula: $C_{35}H_{61}NO_{12}$

Molecular weight: 687.87

CAS Registry No.: 3922-90-5, 6696-47-5 (HCl), 7060-74-4 (phosphate)

Merck Index: 6962



SAMPLE

Matrix: blood, tissue

Sample preparation: Homogenize (Phycotron) liver with 4 volumes of ice-cold saline. Put 200 μ L plasma or liver homogenate into the tube. Add 2 mL MTBE and 5 μ L 1 M NaOH and shake mechanically for 5 min. Centrifuge at 1500 g for 10 min, transfer the upper layer into a glass tube and evaporate it to dryness under dry nitrogen. Rinse the inner wall of the tube with 200 μ L MeOH and evaporate to dryness. Dissolve the residue in 30 μ L MeOH and inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Cosmosil 5-C18 (Nacalai Tesque)

Mobile phase: MeCN:100 mM pH 6.6 sodium acetate buffer 50:50

Flow rate: 0.6

Injection volume: 10

Detector: E, BAS LC-4C, 1.1 V, Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 8.5

OTHER SUBSTANCES

Extracted: erythromycin


KEY WORDS

rat; plasma; liver; pharmacokinetics; oleandomycin is IS

REFERENCE

Hanada,E.; Ohtani,H.; Kotaki,H.; Sawada,Y.; Iga,T. Determination of erythromycin concentrations in rat plasma and liver by high-performance liquid chromatography with amperometric detection, *J.Chromatogr.B*, **1997**, 692, 478–482.

Olsalazine



Molecular formula: $C_{14}H_{10}N_2O_6$

Molecular weight: 302.24

CAS Registry No.: 15722-48-2, 6054-98-4 (Na salt)

Merck Index: 6976

Lednicer No.: 4 42

SAMPLE

Matrix: blood, urine

Sample preparation: Serum. Add disodium-3,5'-azo-bis-(6-hydroxybenzoate) to serum, treat with proteinase K (0.5 mg/mL protein) for 10 min, add tetrabutylammonium hydrogen sulfate buffered to pH 6.5, add dichloromethane, agitate for 30 min, centrifuge. Remove the organic layer and evaporate it to dryness, reconstitute the residue in mobile phase, inject an aliquot.

Urine. Add 2,4-dihydroxybenzoic acid to urine, add perchloric acid, add diethyl ether, shake for 10 min, freeze. Remove the organic layer and add it to pH 7.4 phosphate buffer, extract, inject an aliquot of the aqueous phase.

HPLC VARIABLES

Guard column: 30 × 4 30-40 µm Perisorb RP-18

Column: 250 × 4 10 µm Nucleosil C18

Mobile phase: MeOH:buffer 52:48 (Buffer was pH 7.4 phosphate buffer containing 20 mM tetrabutylammonium hydrogen sulfate.)

Detector: UV 365

CHROMATOGRAM

Retention time: k' 5.2

Internal standard: disodium-3,5'-azo-bis-(6-hydroxybenzoate), 2,4-dihydroxybenzoic acid

Limit of quantitation: 1 µM urine, 0.5 µM (serum)

OTHER SUBSTANCES

Extracted: olsalazine sulfate

KEY WORDS

serum; pharmacokinetics

REFERENCE

Ryde, E.M.; Ahnfelt, N.-O. The pharmacokinetics of olsalazine sodium in healthy volunteers after a single i.v. dose and after oral doses with and without food, *Eur.J.Clin.Pharmacol.*, **1988**, 34, 481-488.

Omeprazole


Molecular formula: C₁₇H₁₉N₃O₃S

Molecular weight: 345.42

CAS Registry No.: 73590-58-6

Merck Index: 6977

Lednicer No.: 4 133



SAMPLE

Matrix: blood

Sample preparation: Extract 500 µL plasma with dichloromethane containing 50 µL 1.0 M dibasic sodium phosphate and 100 µL 25 µg/mL IS. Centrifuge and aspirate aqueous layer to waste. Evaporate organic layer under a stream of nitrogen, reconstitute residue in 200 µL MeCN:20 mM pH 8.0 phosphate buffer 25:75. Inject a 10 µL aliquot.

HPLC VARIABLES

Column: 150 × 4.6 Zorbax SB C18

Mobile phase: Gradient. MeCN:20 mM pH 7.5 phosphate buffer from 30:70 to 40:60 over 10 min. (Mobile phase was contained 0.1% triethylamine.)

Injection volume: 10

Detector: UV 302, to UV 280 at 6.4 min

CHROMATOGRAM

Retention time: 5.7

Internal standard: carbamazepine (7.3)

KEY WORDS

plasma

REFERENCE

Sarich,T.; Kalhorn,T.; Magee,S.; Al-sayegh,F.; Adams,S.; Slattery,J.; Goldstein,J.; Nelson,S.; Wright,J. The effect of omeprazole pretreatment on acetaminophen metabolism in rapid and slow metabolizers of S-mephenytoin, *Clin.Pharmacol.Ther.*, **1997**, 62, 21-28.

SAMPLE

Matrix: blood

Sample preparation: Add 10 μ L 72 μ g/mL flunitrazepam in MeOH to 1 mL plasma, shake briefly, add 3 mL toluene:isoamyl alcohol 95:5, vortex at 1000 rpm for 90 s, centrifuge at 2600 g for 10 min. Evaporate a 2.5 mL aliquot of the upper organic layer to dryness under nitrogen at 40°, reconstitute with 100 μ L mobile phase, inject a 25 μ L aliquot.

HPLC VARIABLES

Guard column: 10 \times 4.6 5 μ m Nucleosil 120-5 C18

Column: 250 \times 4 5 μ m Nucleosil 120-5 C 18

Mobile phase: MeOH:buffer 47:53 (Buffer was 100 mM Na₂HPO₄ adjusted to pH 7.8 with orthophosphoric acid.)

Column temperature: 37

Flow rate: 1.2

Injection volume: 25

Detector: UV 302

CHROMATOGRAM

Retention time: 10.1

Internal standard: flunitrazepam (11.4)

Limit of quantitation: 9.7 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma

REFERENCE

Macek,J.; Ptáček,P.; Klíma,J. Determination of omeprazole in human plasma by high-performance liquid chromatography, *J.Chromatogr.B*, **1997**, 689, 239-243.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 14.065


KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

Ondansetron

Molecular formula: $C_{18}H_{19}N_3O$ **Molecular weight:** 293.37**CAS Registry No.:** 99614-02-5, 116002-70-1, 99614-01-4
(HCl dihydrate), 103639-04-9 (HCl dihydrate)**Merck Index:** 6979**Lednicer No.:** 5 164**SAMPLE****Matrix:** formulations**HPLC VARIABLES****Column:** 250 × 4.6 5 μ m Spherisorb ODS-2**Mobile phase:** MeCN:20 mM KH_2PO_4 43:57**Flow rate:** 1.8**Injection volume:** 20**Detector:** UV 254**CHROMATOGRAM****Retention time:** 9.6**OTHER SUBSTANCES****Simultaneous:** idarubicin**KEY WORDS**

0.9% NaCl; injections

REFERENCE

Zhang,H.; Ye,L.; Stewart,J.T. HPLC determination of idarubicin-etoposide and idarubicin-ondansetron mixtures in 0.9% sodium chloride injection USP, *J.Liq.Chromatogr.Rel.Technol.*, **1998**, 21, 979–988.

SAMPLE**Matrix:** formulations**Sample preparation:** If necessary, dilute injection 1:9 with mobile phase (for 50 mL admixtures) and 1:4 (for 100 mL admixtures), inject a 20 μ L aliquot.**HPLC VARIABLES****Column:** 125 × 4 5 μ m LiChrospher 60 RP-Select B**Mobile phase:** MeCN:buffer 25:75 (Buffer was 20 mM KH_2PO_4 adjusted to pH 6.0 with NaOH solution.)**Injection volume:** 20**Detector:** UV 241**CHROMATOGRAM****Retention time:** 12.5**OTHER SUBSTANCES****Simultaneous:** dexamethasone

KEY WORDS

5% dextrose; 0.9% sodium chloride; injections; stability-indicating

REFERENCE

Evrard,B.; Ceccato,A.; Gaspard,O.; Delattre,L.; Delporte,J.-P. Stability of ondansetron hydrochloride and dexamethasone sodium phosphate in 0.9% sodium chloride injection and in 5% dextrose injection, *Am.J.Health-Syst.Pharm.*, **1997**, *54*, 1065–1068.

SAMPLE

Matrix: solutions

Sample preparation: Add 380 µg propofol and 83 µg ondansetron hydrochloride to 0.9% sodium chloride, shake vigorously for 2 min, make up to 10 mL with 0.9% sodium chloride, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 300 × 4.6 10 µm T-Bondapak phenyl

Mobile phase: MeCN:buffer 50:50 (Buffer was 10 mM KH₂PO₄ adjusted to pH 4.0 with 10% phosphoric acid.)

Flow rate: 1

Injection volume: 20

Detector: UV 268

CHROMATOGRAM

Retention time: 8.5

Limit of detection: 61 ng/mL

OTHER SUBSTANCES

Simultaneous: propofol

REFERENCE

King,D.T.; Stewart,J.T.; Venkateshwaran,T.G. HPLC determination of propofol-thiopental sodium and propofol-ondansetron mixtures, *J.Liq.Chromatogr.Rel.Technol.*, **1996**, *19*, 2285–2294.

SAMPLE

Matrix: solutions

Sample preparation: Inject a 20 µL aliquot of a solution in 0.9% sodium chloride.

HPLC VARIABLES

Column: 220 × 4.6 5 µm underivatized silica (Brownlee Silica Applied Biosystems, Inc., San Jose)

Mobile phase: MeOH:buffer 40:60 (Buffer was 10 mM KH₂PO₄ adjusted to pH 4.0 with 10% phosphoric acid.)

Flow rate: 1

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 6.7

Limit of detection: 47 ng/mL

OTHER SUBSTANCES

Simultaneous: meperidine, morphine (UV 233)

REFERENCE

Venkateshwaran,T.G.; Stewart,J.T.; King,D.T. HPLC determination of morphine-ondansetron and meperidine-ondansetron mixtures in 0.9% sodium chloride injection, *J.Liq.Chromatogr.Rel.Technol.*, **1996**, *19*, 1329–1338.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 10 µm Partisil ODS1

Mobile phase: MeOH:50 mM pH 3.0 phosphoric acid 40:60

Column temperature: 30

Flow rate: 1.5

Detector: radioactivity detection

KEY WORDS

¹⁴C labeled

REFERENCE


Collett,A.; Sims,E.; Walker,D.; He,Y.-L.; Ayrton,J.; Rowland,M.; Warhurst,G. Comparison of HT29-18-C₁ and Caco-2 cell lines as models for studying intestinal paracellular drug absorption, *Pharm.Res.*, **1996**, *13*, 216–221.

Ontazolast

Molecular formula: C₂₁H₂₅N₃O

Molecular weight: 335.45

CAS Registry No.: 147432-77-7

**SAMPLE**

Matrix: solutions

HPLC VARIABLES

Column: 150 × 4.6 5 µm Zorbax LC-8

Mobile phase: MeCN:MeOH:50 mM pH 6 phosphate buffer 1:3:8

Flow rate: 2

Detector: UV 245

CHROMATOGRAM

Retention time: 7

Internal standard: BIRM 390 BS (10.5)

Limit of quantitation: 5 ng/mL

REFERENCE

Hauss,D.J.; Fogal,S.E.; Ficorilli,J.V.; Price,C.A.; Roy,T.; Jayaraj,A.A.; Keirns,J.J. Lipid-based delivery systems for improving the bioavailability and lymphatic transport of a poorly water-soluble LTB₄ inhibitor, *J.Pharm.Sci.*, **1998**, *87*, 164–169.


Opipramol

Molecular formula: C₂₃H₂₉N₃O

Molecular weight: 363.50

CAS Registry No.: 315-72-0, 909-39-7 (2.HCl)

Merck Index: 6985

**SAMPLE**

Matrix: blood

Sample preparation: Deproteinize 100 µL plasma with 100 µL MeCN containing 2.0 mg/L IS. Inject a 50 µL aliquot of the supernatant.

HPLC VARIABLES

Guard column: 4 × 4 Superspher

Column: 125 × 4 Superspher 60 select B
Mobile phase: MeCN:0.07% orthophosphoric acid 20:80
Column temperature: 30
Flow rate: 1.2
Injection volume: 50
Detector: UV 210, UV 255

CHROMATOGRAM

Retention time: 4.1
Internal standard: methylphenyl-phenylhydantion (7.8)
Limit of detection: 50 µg/L

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma

REFERENCE

Lapenberg-Pelzer, M. Identification and determination of opipramol metabolites in plasma and urine, *J. Anal. Toxicol.*, **1998**, 22, 215–219.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 256

CHROMATOGRAM

Retention time: 7.62

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melfalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acen-

ocoumarol; vindesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrridine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254-262.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 255.8

CHROMATOGRAM

Retention time: 14.163

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, *763*, 149-163.


Ormetoprim

Molecular formula: C₁₄H₁₈N₄O₂

Molecular weight: 274.32

CAS Registry No.: 6981-18-6

Lednicer No.: 2 302



SAMPLE

Matrix: blood, tissue

Sample preparation: Plasma. 200 μ L Plasma + 90 μ L 24% trichloroacetic acid in MeOH + 10 μ L 10 μ g/mL sulfamethoxazole in MeOH, vortex for 30 s, centrifuge at 14000 g for 5 min, inject a 50 μ L aliquot of the supernatant. Muscle. Homogenize 1 g muscle in 1.5 mL MeOH:buffer 20:80, add 50 μ L 10 μ g/mL sulfamethoxazole in MeOH, mix thoroughly for 1 min, centrifuge at 14000 g for 5 min, inject a 50 μ L aliquot of the supernatant. (Buffer was 25 mM NaH₂PO₄ containing 15 mM sodium 1-heptanesulfonate adjusted to pH 2.8 with 5 M phosphoric acid.)

HPLC VARIABLES

Guard column: 20 \times 4.6 40 μ m ODS-Hypersil

Column: 150 \times 4.6 3 μ m ODS-Hypersil C18

Mobile phase: MeCN:buffer:triethylamine 20:80:0.02 (Buffer was 25 mM NaH₂PO₄ containing 15 mM sodium 1-heptanesulfonate adjusted to pH 2.8 with 5 M phosphoric acid.)

Flow rate: 1

Injection volume: 50

Detector: UV 270

CHROMATOGRAM

Retention time: 12.5

Internal standard: sulfamethoxazole (8)

Limit of detection: 50 ng/g (muscle), 30 ng/mL (plasma)

OTHER SUBSTANCES

Extracted: sulfadimethoxine

KEY WORDS

plasma; muscle; fish; salmon

REFERENCE

Samuelsen, O.B. Simultaneous determination of ormetoprim and sulphadimethoxine in plasma and muscle of Atlantic salmon (*Salmo salar*), *J.Chromatogr.B*, **1994**, 660, 412–417.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a solution in MeCN:MeOH:100 mM pH 4.5 acetate buffer 12.5:12.5:75, inject an aliquot.

HPLC VARIABLES

Column: C18 (Rainin)

Mobile phase: MeCN:MeOH:50 mM phosphate buffer 12.5:12.5:75, pH 3.0

Flow rate: 1.5

Detector: UV 280

CHROMATOGRAM

Retention time: 8

OTHER SUBSTANCES

Simultaneous: sulfamethoxazole, trimethoprim

REFERENCE

Brown,M.P.; Gronwall,R.; Castro,L. Pharmacokinetics and body fluid and endometrial concentrations of trimethoprim-sulfamethoxazole in mares, *Am.J.Vet.Res.*, **1988**, 49, 918-922.

SAMPLE

Matrix: tissue

Sample preparation: Condition a Sep-Pak C18 SPE cartridge with 5 mL MeOH and 5 mL water. 5 g Fish + 1 mL 100 µg/mL carbamazepine diol in MeOH + 15 mL MeCN + 500 µL 50% trichloroacetic acid, homogenize (Brinkmann Polytron PT 10/35) at medium speed for 30 s, centrifuge at 4° at 7800 g for 25 min. Remove the supernatant and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 5 mL water, vortex for 30 s, filter (13 mm dia. 8 µm Membra-Fil (Nucleopore)), add the filtrate to the SPE cartridge, elute with 5 mL MeOH. Evaporate the eluate to dryness under a stream of nitrogen at 40°, reconstitute the residue in 1 mL mobile phase, vortex for 30 s, inject a 20 µL aliquot. (Flush injection valve with 1 mL mobile phase between analyses.)

HPLC VARIABLES

Guard column: 15 × 3.2 NewGuard RP-18

Column: 250 × 4.6 5 µm Ultrasphere

Mobile phase: MeCN:MeOH:100 mM pH 4.0 phosphate buffer 17:10:73

Flow rate: 1

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 5.5

Internal standard: carbamazepine diol (10.5)

Limit of quantitation: 0.2 ppm

OTHER SUBSTANCES

Extracted: sulfadimethoxine

Simultaneous: sulfacetamide, sulfadiazine, sulfamerazine, sulfamethazine, sulfoxazole

KEY WORDS

fish; salmon; SPE; pharmacokinetics

REFERENCE

Walisser,J.A.; Burt,H.M.; Valg,T.A.; Kitts,D.D.; McErlane,K.M. High-performance liquid chromatographic analysis of Romet-30 in salmon following administration of medicated feed, *J.Chromatogr.*, **1990**, 518, 179-188.

Orphenadrine

Molecular formula: C₁₆H₂₃NO

Molecular weight: 269.39

CAS Registry No.: 83-98-7, 4682-36-4 (citrate)

Merck Index: 7007

Lednicer No.: 1 42


SAMPLE

Matrix: blood

Sample preparation: Add 300 µL MeCN to 100 µL plasma, vortex, centrifuge for 2 min. Remove the supernatant and add it to 300 µL pH 5.9 sodium phosphate buffer and 3 mL hexane and vortex for 45 s. Centrifuge at 2500 g for 3 min. Evaporate the supernatant under a stream of nitrogen and reconstitute with mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 100 × 4.6 ODS



Mobile phase: MeCN:MeOH:25 mM potassium phosphate 25:25:50 containing 0.75 mL/L 2 M sulfuric acid and 0.25 mL/L triethylamine

Flow rate: 1.0

Detector: UV 250

CHROMATOGRAM

Retention time: 5.7

Internal standard: orphenadrine

OTHER SUBSTANCES

Extracted: ethopropazine

KEY WORDS

plasma; rat; pharmacokinetics; orphenadrine is IS

REFERENCE

Padovani,P.K.; Timby,D.M.; Wright,M.R.; Kapil,R.P. Quantitative analysis of DMP 851 in rat and dog plasma by liquid-liquid extraction and reverse-phase high performance liquid chromatography with ultraviolet detection (Abstract 3318), *Pharm.Res.*, **1997**, *14*, S568.

SAMPLE

Matrix: blood

Sample preparation: 10 mL Plasma or whole blood + 1 mL 1 M NaOH, extract twice with 10 mL hexane for 30 min. Remove the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 1 mL 100 mM HCl, add 5 mL chloroform, vortex for 1 min, centrifuge. Remove a 4.5 mL aliquot of the organic layer and evaporate it to dryness, reconstitute the residue in 100 μ L mobile phase, inject a 50 μ L aliquot. (It is implied, but not explicitly stated in the paper, that this extraction procedure works for this compound.)

HPLC VARIABLES

Column: 10 μ m Micropak CN (Varian)

Mobile phase: MeCN:20 mM ammonium acetate 90:10

Flow rate: 2.5

Injection volume: 50

Detector: UV 254

CHROMATOGRAM

Retention time: 8.2

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Simultaneous: acetophenazine, amitriptyline, benztrapine, butaperazine, carphenazine, fluphenazine, haloperidol, imipramine, mesoridazine, nortriptyline, piperacetazine, promazine, promethazine, thioridazine, thiothixene, trifluoperazine, triflupromazine, trihexyphenidyl, trimeprazine

Interfering: chlorpromazine

KEY WORDS

plasma; whole blood

REFERENCE

Curry,S.H.; Brown,E.A.; Hu,O.Y.-P.; Perrin,J.H. Liquid chromatographic assay of phenothiazine, thioxanthene and butyrophenone neuroleptics and antihistamines in blood and plasma with conventional and radial compression columns and UV and electrochemical detection, *J.Chromatogr.*, **1982**, *231*, 361-376.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 μ g/mL solution in MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 125 \times 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 3.6

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipamnone, diprenorphine, dipyridamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenylglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindine, pimozide, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, 323, 191–225.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 cellulose tris(3,5-dimethylphenylcarbamate)

Mobile phase: Hexane:isopropanol 90:10

Flow rate: 0.5

Detector: UV

CHROMATOGRAM

Retention time: k' 0.56 (of first (+) enantiomer)


KEY WORDSchiral; α 1.89**REFERENCE**

Okamoto,Y.; Aburatani,R.; Hatano,K.; Hatada,K. Optical resolution of racemic drugs by chiral HPLC on cellulose and amylose tris(phenylcarbamate) derivatives, *J.Liq.Chromatogr.*, **1988**, *11*, 2147–2163.

SAMPLE**Matrix:** solutions**HPLC VARIABLES****Column:** 150 × 3.9 5 μ m Spherisorb C8**Mobile phase:** MeCN:buffer 60:40 (Buffer was 1.5 mL triethylamine in 1 L water adjusted to pH 3.0 with 85% phosphoric acid.)**Flow rate:** 1.5**Detector:** UV 199**CHROMATOGRAM****Retention time:** k' 3.90**OTHER SUBSTANCES****Simultaneous:** hyoscyamine, bromocriptine, benztropine, biperiden**Noninterfering:** amantadine, carbidopa, levodopa**REFERENCE**

Selinger,K.; Lebel,G.; Hill,H.M.; Discenza,C. High-performance liquid chromatographic method for the analysis of benztropine in human plasma, *J.Chromatogr.*, **1989**, *491*, 248–252.

Oryzalin

Molecular formula: $C_{12}H_{18}N_4O_6S$ **Molecular weight:** 346.36**CAS Registry No.:** 19044-88-3**Merck Index:** 7015**SAMPLE****Matrix:** blood, microsomal incubations

Sample preparation: Condition a Bond Elut C18 SPE cartridge with 3 mL MeCN and 3 mL water. Plasma. Add 500 μ L plasma to the SPE cartridge, dry with vacuum. Elute with three 200 μ L portions of MeCN. Inject a 10-50 μ L aliquot of the eluate. Microsomal incubations. Centrifuge microsomal incubation at 12000 g. Add 900 μ L portion of the supernatant to the SPE cartridge, dry with vacuum. Elute with three 200 μ L portions of MeCN. Inject a 10-50 μ L aliquot of the eluate.

HPLC VARIABLES**Column:** 150 × 4.6 5 μ m Adsorbosphere HS C18**Mobile phase:** MeCN:water 53:47**Flow rate:** 2**Injection volume:** 10-50**Detector:** UV 254; MS, Finnigan TSQ tandem mass, APCI, vaporizer 450°, corona discharge current 5 μ A, CID, neutral argon -31 eV**CHROMATOGRAM****Retention time:** 4.7**Limit of detection:** 100 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

SPE; plasma; mouse; rat; liver; pharmacokinetics

REFERENCE

Dvorakova,K.; Dorr,R.T.; Gallegos,A.; McClure,T.; Powis,G. Pharmacokinetic studies of the herbicide and antitumor compound oryzalin in mice, *J.Chromatogr.B*, **1997**, 696, 275–281.

Oxacillin


Molecular formula: C₁₉H₁₉N₃O₅S

Molecular weight: 401.44

CAS Registry No.: 66-79-5, 7240-38-2 (Na salt monohydrate), 1173-88-2 (Na salt)

Merck Index: 7036

Lednicer No.: 1 413



SAMPLE

Matrix: blood

Sample preparation: 400 µL Serum + 400 µL MeCN, vortex for 10 s, shake slowly for 15 min, centrifuge at 3000 g for 10 min. Remove the supernatant and add it to 4 mL dichloromethane, vortex for 10 s, shake for 15 min, centrifuge at 3000 g for 10 min, inject a 50 µL aliquot of the upper aqueous layer.

HPLC VARIABLES

Column: µBondapak C18

Mobile phase: MeCN:water:200 mM ammonium acetate 28:62:10, pH 5.6

Flow rate: 1

Injection volume: 50

Detector: UV 254

CHROMATOGRAM

Retention time: 5.7

Limit of detection: 500 ng/mL

OTHER SUBSTANCES

Extracted: cloxacillin, dicloxacillin, methicillin, nafcillin

Noninterfering: amdinocillin (mecillinam), amikacin, amoxicillin, ampicillin, carbenicillin, cefachmandole, cefazolin, ceforanide, cefatoxamine, cefoxitin, cephalixin, cephaloridine, cephalothin, cephradine, cepharin, chloramphenicol, clindamycin, co-trimoxazole, fluorocytosine, gentamicin, metronidazole, moxalactam, penicillin, piperacillin, sulfamethoxazole, theophylline, ticarcillin, tobramycin, trimethoprim, vancomycin

KEY WORDS

serum

REFERENCE

Rudrik,J.T.; Bawdon,R.E. Determination of penicillinase-resistant penicillins in serum using high-pressure liquid chromatography, *J.Liq.Chromatogr.*, **1981**, 4, 1525–1545.

SAMPLE

Matrix: blood, milk

Sample preparation: Milk. Adjust 5 mL milk to pH 6.3 with 100 mM HCl, deproteinize with 10 mL MeCN. Centrifuge at 1932 g for 20 min and extract the aqueous phase with two 5 mL

portions of chloroform for 20 min (Caution! Chloroform is a carcinogen!). Centrifuge at 1932 g for 20 min. Evaporate the organic phase to dryness, reconstitute the residue in 200 μ L mobile phase, inject a 100 μ L aliquot. Serum. Adjust 2.5 mL serum to pH 6.3 with 100 mM HCl, deproteinize with 10 mL MeCN. Centrifuge at 1932 g for 20 min and extract the aqueous phase with two 5 mL portions of dichloromethane for 20 min. Centrifuge at 1932 g for 20 min. Evaporate the organic phase to dryness, reconstitute the residue in 200 μ L mobile phase, inject a 100 μ L aliquot

HPLC VARIABLES

Column: 15 \times 3.9 4 μ m Nova Pack C18
Mobile phase: MeCN:20 mM KH_2PO_4 21:79, pH 5
Flow rate: 1.2
Injection volume: 100
Detector: UV 225

CHROMATOGRAM

Retention time: 3.8
Internal standard: oxacillin

OTHER SUBSTANCES

Extracted: cloxacillin

KEY WORDS

serum; cow; oxacillin is IS

REFERENCE

Pérez,B.; Prats,C.; Castells,E.; Arboix,M. Determination of cloxacillin in milk and blood of dairy cows by high-performance liquid chromatography, *J.Chromatogr.B*, **1997**, 698, 155–160.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma, serum. 100 μ L Plasma or serum + cloxacillin + 100 μ L 500 mM pH 2.2 citric acid buffer + 20 μ L 500 mM HCl + 2.5 mL dichloromethane, extract. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 35°, reconstitute the residue in mobile phase, inject an aliquot. Urine. Dilute urine with water, inject an aliquot.

HPLC VARIABLES

Guard column: 50 \times 2.1 ODS pellicular
Column: 250 \times 4.6 5 μ m Lichrosorb RP-8
Mobile phase: MeCN:20 mM pH 6.6 sodium acetate 34:100
Flow rate: 1
Detector: UV 220

CHROMATOGRAM

Retention time: 7
Internal standard: cloxacillin (9)
Limit of detection: 400 ng/mL

OTHER SUBSTANCES

Extracted: metabolites, dicloxacillin, flucloxacillin

KEY WORDS

plasma; serum

REFERENCE

Thijssen,H.H.W. Analysis of isoxazolyl penicillins and their metabolites in body fluids by high-performance liquid chromatography, *J.Chromatogr.*, **1980**, 183, 339–345.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 1 mL Plasma + 20 μ L 4% aqueous sodium dodecyl hydrogen sulfate solution, shake for 30 min, filter (Amicon MPS-1 micropartition system, YMT membrane) while centrifuging, adjust the pH of the ultrafiltrate to 6.3-6.5 with pH 4 citrate buffer, inject a 500 μ L aliquot onto column A with mobile phase A and elute to waste, after 10 min elute the contents of column A onto column B with mobile phase B, elute with mobile phase B, monitor the effluent from column B. Urine. Make up 5-100 μ L urine to 500 μ L with water, inject onto column A with mobile phase A and elute to waste, after 10 min elute the contents of column A onto column B with mobile phase B, elute with mobile phase B, monitor the effluent from column B.

HPLC VARIABLES

Column: A 50 \times 4 Nucleosil 5-C18; B 250 \times 5 Nucleosil 5-C18

Mobile phase: A MeCN:33 mM NaH_2PO_4 5:95; B MeCN:33 mM NaH_2PO_4 20:80

Injection volume: 500

Detector: UV 210

CHROMATOGRAM

Internal standard: oxacillin

OTHER SUBSTANCES

Extracted: penicillin V

KEY WORDS

plasma; column-switching; ultrafiltrate; oxacillin is IS

REFERENCE

Lintz,W.; Hirsch,I.; Osterloh,G.; Schmidt-Böthelt,E.; Sous,H. Bioverfügbarkeit von Penicillin V in einer wäßrigen Zubereitungsform [Bioavailability of penicillin V in aqueous dosage forms], *Arzneimittelforschung*, 1984, 34, 66-71.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 14.76

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, 1997, 763, 149-163.

SAMPLE**Matrix:** formulations**Sample preparation:** Blend tablets and capsules with water in a high-speed blender for 5 min, filter, dilute with mobile phase, inject a 20 μ L aliquot. Dilute oral suspensions and injections with mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES**Guard column:** 70 mm long Co:Pell ODS**Column:** 300 \times 4.6 10 μ m Chromegabond C18 (E.S. Industries)**Mobile phase:** MeCN:MeOH:10 mM KH_2PO_4 19:11:70**Flow rate:** 1**Injection volume:** 20**Detector:** UV 225

CHROMATOGRAM**Retention time:** 10.0**Limit of detection:** 200 ng/mL

OTHER SUBSTANCES**Simultaneous:** amoxicillin, ampicillin, cloxacillin, dicloxacillin, methicillin, nafcillin, penicillin G, penicillin V

KEY WORDStablets; capsules; oral suspensions; injections

REFERENCEBriguglio, G.T.; Lau-Cam, C.A. Separation and identification of nine penicillins by reverse phase liquid chromatography, *J. Assoc. Off. Anal. Chem.*, **1984**, 67, 228-231.

SAMPLE**Matrix:** milk**Sample preparation:** Mix 10 mL milk with 2 mL 100 mM tetraethylammonium chloride, add 40 mL MeCN slowly with continual stirring, let stand for 10 min, decant the supernatant through a plug of glass wool. Collect 40 mL filtrate, add 2 mL buffer, evaporate to 1-2 mL under reduced pressure at 40-50°, dilute to 4 mL with water, filter (0.45 μ m PVDF). Inject a 2 mL aliquot onto a 150 \times 4.6 5 μ m Supelcosil LC-18 column, elute with MeCN:10 mM KH_2PO_4 0:100 for 3 min, to 60:40 over 37 min at 1 mL/min, collect a 1.5-2 mL aliquot containing the compound (ca. 27.0 min), evaporate to <1 mL under reduced pressure, make up to 1 mL with water, inject an aliquot. (Prepare the buffer by mixing 10 mM KH_2PO_4 and 10 mM Na_2HPO_4 in a 5:1 ratio, pH 6.)

HPLC VARIABLES**Column:** 150 \times 4.6 5 μ m Supelcosil LC-18-DB**Mobile phase:** MeCN:buffer 38:62 (Buffer was 2 mM phosphoric acid containing 8 mM potassium dihydrogen phosphate.)**Flow rate:** 1**Injection volume:** 200**Detector:** UV 215

REFERENCEMoats, W.A.; Romanowski, R.D. Multiresidue determination of β -lactam antibiotics in milk and tissues with the aid of high-performance liquid chromatographic fractionation for clean up, *J. Chromatogr. A*, **1998**, 812, 237-247.

SAMPLE**Matrix:** milk**Sample preparation:** Condition a 6 mL 500 mg Bond Elut C18 SPE cartridge with 10 mL MeOH, 10 mL water, 5 mL 2% NaCl, and 5 mL 100 mM pH 8 phosphate extraction buffer. Add 30 mL 100 mM pH 8 phosphate extraction buffer to 5 mL milk, add 1.65 mL 1 M sulfuric acid to reach pH 4.0-4.5, vortex for 30 s, centrifuge at 2400 g for 10 min, add 600 μ L 5 M NaOH to the supernatant to reach pH 8, vortex, centrifuge at 2400 g for 5 min. Add the supernatant

to a reservoir attached to the SPE cartridge, pull through the SPE cartridge at 3 mL/min, remove the reservoir and elute with 1 mL MeCN:water 40:60. Add 500 μ L derivatizing reagent to the eluate, vortex, heat at 65° for 10 min, cool to room temperature (protect from light), inject a 100 μ L aliquot of the derivatized sample. (Prepare the 100 mM pH 8 phosphate extraction buffer as follows. Dissolve 15.6 g K_2HPO_4 dihydrate in 800 mL water, adjust pH to 8 with 10 M NaOH, make up to 1 L. Prepare the derivatizing reagent as follows. Weigh out 13.78 g 1,2,4-triazole, add 60 mL water, stir, add 10 mL 100 mM mercuric chloride solution, mix, adjust pH to 9.0 ± 0.5 with 5 M NaOH, dilute to 100 mL with water.)

HPLC VARIABLES

Column: 150 \times 3.9 5 μ m Symmetry C8 (Waters)

Mobile phase: MeCN:MeOH:buffer 37:5:58 (Prepare the 100 mM pH 6.5 phosphate buffer containing 15 mM thiosulfate and 30 mM tetrabutylammonium hydrogen sulfate as follows. Weigh 4.969 g anhydrous NaH_2PO_4 , 10.139 g Na_2HPO_4 dihydrate, 3.894 g sodium thiosulfate pentahydrate, and 10.186 g tetrabutylammonium hydrogen sulfate, dissolve in 800 mL water, adjust pH to 6.5 with 5 M NaOH, dilute to 1 L with water, mix thoroughly, filter under vacuum (0.45 μ m).)

Flow rate: 1

Injection volume: 100

Detector: UV 340

CHROMATOGRAM

Retention time: 10

Limit of detection: 2 ng/mL

OTHER SUBSTANCES

Extracted: cloxacillin, dicloxacillin

KEY WORDS

derivatization; SPE

REFERENCE

Verdon,E.; Couedor,P. Determination of isoxazolympenicillins residues in milk by ion-pair reversed-phase high-performance liquid chromatography after precolumn derivatization, *J.Chromatogr.B*, **1998**, 705, 71–78.

SAMPLE

Matrix: milk

Sample preparation: 50 g Milk + 2 drops penicillinase (Difco Laboratories), let stand 1 h at 37°, add 50 mL MeCN, shake vigorously for 1 min, centrifuge at 9000 g for 10 min, decant, add 5 g NaCl, swirl to dissolve, add 100 mL dichloromethane, shake for 1 min, centrifuge at 1000 g for 10 min. Remove top aqueous layer and extract organic layer with 25 mL 10% NaCl by shaking and centrifuging as before. Combine aqueous layers, add 1 mL 0.3% mercuric chloride in water, let stand 30 min, add 1 mL 2 M HCl, extract with three 50 mL portions of dichloromethane by shaking each portion for 1 min and centrifuging at 1000 g for 10 min, filter dichloromethane extracts through 30 g anhydrous sodium sulfate, evaporate to dryness under reduced pressure at 35°, if water remains add 5–10 mL MeOH to flask and complete evaporation. Dissolve residue in 1 mL 10% acetic acid, add 0.5 mL 0.08% dansyl hydrazine in 10% acetic acid, let stand 90 min to overnight in the dark, transfer reaction mixture to a separatory funnel with three 25 mL portions of dichloromethane, add 5 mL 2 M HCl, shake for 1 min, wash organic layer with 5 mL 5% $NaHCO_3$ solution, filter through 10–20 g anhydrous sodium sulfate. Extract acid aqueous layer again with 25 mL dichloromethane. Combine dichloromethane layers and evaporate to dryness at 35° under reduced pressure. Dissolve residue in 2 mL IS solution, inject a 20 μ L aliquot. (Prepare IS solution by dissolving 10 μ L benzaldehyde in 100 mL dichloromethane, evaporate 1 mL to dryness under reduced pressure, dissolve residue in 1 mL 10% acetic acid, add 0.5 mL 0.08% dansyl hydrazine in 10% acetic acid, let stand 90 min to overnight in the dark, transfer reaction mixture to a separatory funnel with three 25 mL portions of dichloromethane, add 5 mL 2 M HCl, shake for 1 min, wash organic layer with 5 mL 5% $NaHCO_3$ solution, filter through 10–20 g anhydrous sodium sulfate. Extract acid aqueous layer again with 25 mL dichloromethane. Combine dichloromethane layers and evaporate to dryness at 35° under reduced pressure. Dissolve residue in 100 mL MeCN then dilute an aliquot 1:4 with MeCN.)

HPLC VARIABLES**Column:** 250 × 4 10 µm Lichrosorb RP-18**Mobile phase:** MeCN:water 58:42**Flow rate:** 1**Injection volume:** 20**Detector:** F ex 254 em 500 filter

CHROMATOGRAM**Retention time:** 7.21**Internal standard:** benzaldehyde (derivatized) (12.18)**Limit of detection:** 5 ng/g

OTHER SUBSTANCES**Extracted:** penicillin G, phenethicillin, methicillin, cloxacillin, dicloxacillin, nafcillin**Interfering:** penicillin V, phenethicillin

KEY WORDS

derivatization

REFERENCE

Munns,R.K.; Shimoda,W.; Roybal,J.E.; Vieira,C. Multiresidue method for determination of eight neutral β -lactam penicillins in milk by fluorescence-liquid chromatography, *J.Assoc.Off.Anal.Chem.*, **1985**, *68*, 968–971.

SAMPLE**Matrix:** milk

Sample preparation: Add 2 volumes MeCN to milk, stand 5 min, decant aqueous portion, suction filter, extract with an equal volume of 1:1 methylene chloride:hexane, centrifuge aqueous phase at 3000 rpm for 10 min. Dilute 3:1 with 20 mM sodium acetate buffer and filter (0.2 µm nylon). Inject 50 µL onto column with mobile phase A, run mobile phase A for 30 min and elute to waste. After 30 min switch to mobile phase B and elute through detector.

HPLC VARIABLES**Column:** 100 × 8 Radial-Pak 10 µm µBondapak C18**Mobile phase:** A 20 mM sodium acetate buffer; B Gradient. MeCN:MeOH:20 mM sodium acetate buffer from 15:10:75 to 30:0:70 over 15 min and hold at 30:0:70**Flow rate:** A 3; B 2**Injection volume:** 50

Detector: E, Waters 464 pulsed electrochemical detector using a thin layer cell with a Ag/AgCl reference electrode. E1 = 1300 mV for 0.166 s, E2 = 1500 mV for 0.166 s, E3 = -200 mV for 0.333 s.

CHROMATOGRAM**Retention time:** 13.4**Limit of detection:** 0.2 ppm

OTHER SUBSTANCES**Extracted:** penicillin V, ampicillin, methicillin, penicillin G, cloxacillin, nafcillin, dicloxacillin

KEY WORDS

column-switching

REFERENCE

Kirchmann,E.; Earley,R.L.; Welch,L.E. The electrochemical detection of penicillins in milk, *J.Liq.Chromatogr.*, **1994**, *17*, 1755–1772.

SAMPLE**Matrix:** milk

Sample preparation: Condition a Bond Elut C8 SPE cartridge with 5 mL MeOH and 5 mL water. 20 mL Milk + 20 mL buffer, heat at 60° for 20 min or until milk curdles, centrifuge for 10 min, add the supernatant to the SPE cartridge, wash with two 2.5 mL portions of water,

elute with 2.5 mL MeOH. Evaporate the eluate to dryness under a stream of nitrogen, extract the residue with three 100 μ L portions of 50 mM pH 6.0 potassium phosphate buffer, filter (0.2 μ m), inject an aliquot of the filtrate. (Buffer was 545 mL 100 mM citric acid, 455 mL 200 mM Na₂HPO₄, and 74.4 g EDTA, adjust to pH 4.5 with ammonium hydroxide, make up to 2 L with water.)

HPLC VARIABLES

Column: 250 \times 4.6 10 μ m Lichrosorb RP-8

Mobile phase: MeOH:50 mM pH 6.0 potassium phosphate buffer 35:65

Flow rate: 1

Injection volume: 200

Detector: UV 210 or Charm II assay

CHROMATOGRAM

Retention time: 33.10

OTHER SUBSTANCES

Extracted: ampicillin, ceftiofur, cephapirin, cloxacillin, dicloxacillin, nafcillin, penicillin G

Simultaneous: amoxicillin

KEY WORDS

SPE

REFERENCE

Al-Obaidy,S.S.; Po,A.L.W.; McKiernan,P.J.; Glasgow,J.F.T.; Millership,J. Assay of paracetamol and its metabolites in urine, plasma and saliva of children with chronic liver disease, *J.Pharm.Biomed.Anal.*, **1995**, 13, 1033–1039.

SAMPLE

Matrix: milk

Sample preparation: Condition a Bond Elut C8 SPE cartridge with 5 mL MeOH and 5 mL water. 20 mL Milk + 20 mL buffer, heat at 60° for 20 min or until milk curdles, centrifuge for 10 min, add the supernatant to the SPE cartridge, wash with two 2.5 mL portions of water, elute with 2.5 mL MeOH. Evaporate the eluate to dryness under a stream of nitrogen, extract the residue with three 100 μ L portions of 50 mM pH 6.0 potassium phosphate buffer, filter (0.2 μ m), inject an aliquot of the filtrate. (Buffer was 545 mL 100 mM citric acid, 455 mL 200 mM Na₂HPO₄, and 74.4 g EDTA, adjust to pH 4.5 with ammonium hydroxide, make up to 2 L with water.)

HPLC VARIABLES

Column: 250 \times 4.6 10 μ m Lichrosorb RP-8

Mobile phase: MeOH:50 mM pH 6.0 potassium phosphate buffer 35:65

Flow rate: 1

Injection volume: 200

Detector: UV 210 or Charm II assay

CHROMATOGRAM

Retention time: 33.10

OTHER SUBSTANCES

Extracted: ampicillin, ceftiofur, cephapirin, cloxacillin, dicloxacillin, nafcillin, penicillin G

Simultaneous: amoxicillin

KEY WORDS

SPE

REFERENCE

Zomer,E.; Quintana,J.; Saul,S.; Charm,S.E. LC-Receptograms: A method for identification and quantitation of β -lactams in milk by liquid chromatography with microbial receptor assay, *J.AOAC Int.*, **1995**, 78, 1165–1172.

SAMPLE**Matrix:** milk

Sample preparation: Condition a 500 mg tC18 SPE cartridge (Waters) with 20 mL MeOH, 20 mL water, and 10 mL 2% NaCl. Centrifuge 30 mL milk at 1500 g for 10 min. Dilute a 10 mL portion of the defatted milk with 20 mL water, add 200 μ L 2 μ g/mL penicillin V in pH 9.0 buffer, add 6 mL 170 mM sulfuric acid, add 5.6 mL 5% sodium tungstate, shake vigorously for 1 min, allow to stand for 5 min, check that the pH is in the range 4.6-4.8 (if it is outside this range start again using a different volume of sodium tungstate solution), centrifuge at 1500 g for 10 min, adjust the pH of the supernatant to 8.1-8.2 with 5 M and 0.1 M NaOH, filter (glass fiber) the clear liquid phase. Pass the filtrate through the SPE cartridge at 2 mL/min, wash with 2 mL water, dry by pulling air through the cartridge for 1 min, elute with 2 mL MeCN. Add 150 μ L pH 9.0 buffer to the eluate and evaporate to about 100 μ L under a stream of nitrogen at 45-50°, add 400 μ L pH 9.0 buffer, add 75 μ L reagent I, vortex for 30 s, let stand at room temperature for 10 min, use 500 μ L water to transfer the mixture to a separatory funnel, add 20 mL dichloromethane, add 5 mL pH 2.45 buffer, shake for 1 min, let stand for no more than 5 min. Remove the organic layer and evaporate it to dryness under reduced pressure at 35-40°, dissolve the residue in 500 μ L pH 9.0 buffer, add 75 μ L reagent I, vortex for 30 s, let stand at room temperature for 10 min, add 450 μ L reagent II, vortex for 1 min, heat at $55 \pm 1^\circ$ for 30 min, cool, filter (0.45 μ m), inject a 150 μ L aliquot. (Prepare pH 9.0 buffer by dissolving 0.34 g KH_2PO_4 in water, adjusting the pH to 9.0 with NaOH, and making up to 100 mL with water. Prepare pH 2.45 buffer by dissolving 2.72 g KH_2PO_4 in water, adjusting the pH to 2.45 with phosphoric acid, and making up to 100 mL with water. Prepare reagent I by dissolving 1.13 g benzoic anhydride in MeCN, make up to 25 mL with MeCN. Prepare reagent II by dissolving 6.905 g 1,2,4-triazole in 30 mL water and adding 5 mL 26 mM mercuric chloride in water, adjust pH to 9.0 ± 0.05 with 5 M NaOH, make up to 50 mL. Prepare reagents I and II 1-4 h before use. Silanize glassware with dichlorodimethylsilane.)

HPLC VARIABLES**Column:** 150 \times 3.9 μ m Nova-Pak C18

Mobile phase: Gradient. A as MeCN:buffer 10:90. B was MeCN:buffer 30:70. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 13 min, return to initial conditions over 2 min, re-equilibrate at initial conditions for 5 min. (Prepare buffer by dissolving 9.938 g Na_2HPO_4 , 17.938 g $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$, and 4.964 g sodium thiosulfate in water, make up to 2 L with water, pH 6.5.)

Column temperature: 30**Flow rate:** 1**Injection volume:** 150**Detector:** UV 323

CHROMATOGRAM**Retention time:** 34**Internal standard:** penicillin V (28.5)**Limit of detection:** 1.4 ng/mL**Limit of quantitation:** 1.9 ng/mL

OTHER SUBSTANCES**Extracted:** amoxicillin, ampicillin, cloxacillin, dicloxacillin, penicillin G

KEY WORDS

derivatization; cow; SPE

REFERENCE

Sorensen, L.K.; Rasmussen, B.M.; Boison, J.O.; Keng, L. Simultaneous determination of six penicillins in cows' raw milk by a multiresidue high-performance liquid chromatographic method, *J. Chromatogr. B*, **1997**, 694, 383-391.

SAMPLE**Matrix:** perfusate

Sample preparation: 200 μ L Perfusate + 300 μ L MeCN, mix, centrifuge for 5 min, inject a 15 μ L aliquot.

HPLC VARIABLES**Guard column:** 50 \times 3.6 LiChrosorb RP-2

Column: Chemcosorb 5-ODS-H (Chemco, Japan)
Mobile phase: MeOH:pH 5.2 acetate buffer 50:50
Column temperature: 40
Flow rate: 1
Injection volume: 15
Detector: UV 230

CHROMATOGRAM

Limit of quantitation: 2 µg/mL

KEY WORDS

rat; liver; pharmacokinetics

REFERENCE

Yasui,H.; Yamaoka,K.; Fukuyama,T.; Nakagawa,T. Effect of liver intoxication by carbon tetrachloride in hepatic local disposition of oxacillin using moment characteristics as index, *Drug Metab.Dispos.*, **1995**, 23, 779-785.

SAMPLE

Matrix: perfusate
Sample preparation: Inject a 15 µL aliquot.

HPLC VARIABLES

Column: 150 × 4.6 Chemcosorb 5-ODS-H (Chemco, Osaka)
Mobile phase: MeOH:acetate buffer 50:50, pH 5.2
Flow rate: 1
Injection volume: 15
Detector: UV 220

CHROMATOGRAM

Limit of quantitation: 3 µg/mL

KEY WORDS

rat

REFERENCE

Ohata,Y.; Yamaoka,K.; Yasui,H.; Nakagawa,T. Consideration on moments of outflow profile in liver perfusion system with change in perfusate flow rate using oxacillin as model drug, *Biol.Pharm.Bull.*, **1996**, 19, 83-87.

SAMPLE

Matrix: solutions
Sample preparation: React the antibiotic, triethylamine, and 1-(2,5-dihydroxyphenyl)-2-bromoethanone in a 1:2:4 molar ratio in DMF at 45° for 2 h (use dibenzo-18-crown-6 to make the sodium salt soluble), inject a 10 µL aliquot. (Preparation of 1-(2,5-dihydroxyphenyl)-2-bromoethanone is as follows. Stir 27.6 g 1,4-dimethoxybenzene and 28 mL bromoacetyl bromide at 0°, add 53.4 g aluminum bromide over 10 min (an exothermic reactions ensues), let stand at room temperature for 12 h, add 100 mL 48% HBr, add 100 g ice, stir for 1 h, extract twice with 200 mL portions of diethyl ether. Combine the extracts and wash them 3 times with 200 mL portions of water, dry over 40 g anhydrous magnesium sulfate, evaporate to dryness, recrystallize the product 3 times from EtOH to yield 1-(2,5-dihydroxyphenyl)-2-bromoethanone monobromoacetate (mp 105-107°). Dissolve 11 g 1-(2,5-dihydroxyphenyl)-2-bromoethanone monobromoacetate in 200 mL warm dry MeOH saturated with HBr, stir for 18 h, add 200 mL water, cool to -10°. Collect the yellow solid and dry it under vacuum at 50° for 48 h, recrystallize from toluene:heptane 50:50 then toluene to obtain 1-(2,5-dihydroxyphenyl)-2-bromoethanone as yellow needles (mp 117-119°).)

HPLC VARIABLES

Column: 250 × 4 7 µm RP-18 LiChrocart (Merck)
Mobile phase: MeOH:100 mM pH 6.5 sodium acetate 58:42
Flow rate: 1

Injection volume: 10

Detector: E, Bioanalytical Systems Model LC4B, glassy carbon electrode 0.8 V, Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 19

OTHER SUBSTANCES

Extracted: carbenicillin, cephalirin, cloxacillin, dicloxacillin, hetacillin, methicillin, nafcillin, penicillin G

KEY WORDS

derivatization

REFERENCE

Munns,R.K.; Roybal,J.E.; Shimoda,W.; Hurlbut,J.A. 1-(4-Hydroxyphenyl)-, 1-(2,4-dihydroxyphenyl)- and 1-(2,5-dihydroxyphenyl)-2-bromoethanones: new labels for determination of carboxylic acids by high-performance liquid chromatography with electrochemical and ultraviolet detection, *J.Chromatogr.*, **1988**, *442*, 209-218.

SAMPLE

Matrix: solutions

Sample preparation: Separate buffer containing drug from human serum albumin by centrifuging at 37° at 700 g for 3 min using a Micropartition System MPS-1 (Amicon) unit, inject a 10-20 µL aliquot of the ultrafiltrate.

HPLC VARIABLES

Guard column: C18/Corasil (Waters)

Column: 300 × 3.9 µBondapak C18

Mobile phase: MeCN:10 mM ammonium acetate 25:75

Flow rate: 1.5

Injection volume: 10-20

Detector: UV 220

OTHER SUBSTANCES

Also analyzed: penicillin V

REFERENCE

Terasaki,T.; Nouda,H.; Tsuji,A. Relationship between lipophilicity and binding affinity with human serum albumin for penicillin and cephem antibiotics, *J.Pharmacobiodyn.*, **1992**, *15*, 99-106.

SAMPLE

Matrix: tissue

Sample preparation: Homogenize (Ultra-Turrax) 25 g tissue with 25 mL MeCN for 1 min, add 5 mL 500 mM pH 2.2 phosphate buffer while the homogenizer is still running, add 65 mL MeCN, homogenize for 1 min, centrifuge at 4000 g for 10 min. Remove the supernatant and add it to 7 g NaCl and 50 mL dichloromethane, shake for 2 min, allow to stand for 30 min. Remove the upper organic layer and add it to 5 g anhydrous sodium sulfate, shake for 30 s, filter through a cotton-wool plug, evaporate to about 4 mL under reduced pressure at 30°, add 3 mL dichloromethane, evaporate to about 4 mL, add 3 mL light petroleum, evaporate to about 0.5 mL. Suspend this residue with sonication in three 3 mL portions of light petroleum and place these fractions in a separate tube, rinse the original tube with 2 mL pH 7 phosphate buffer. Add the phosphate buffer rinse to the light petroleum extracts, vortex for 30 s, centrifuge, remove the aqueous layer. Extract the light petroleum layer with 2 mL pH 7 phosphate buffer and with two 1.5 mL portions of pH 7 phosphate buffer, combine all the aqueous phase, centrifuge, inject a 200 µL aliquot on to column A and elute to waste with mobile phase B, after 15 min elute to waste with mobile phase C at 2 mL/min, after 10 min elute the contents of column A on to column B with mobile phase D, after 2 min remove column A from the circuit, elute column B with mobile phase D, monitor the effluent from column B. (Wash column A with mobile phase A at 2 mL/min for 7 min, with mobile phase A at 1 mL/min for 5 min, with mobile phase B at 2 mL/min for 8 min, and with mobile phase B at 1 mL/min for 6 min.)

HPLC VARIABLES

Column: A 4 × 4.5 µm LiChrospher 100 RP-18e; B 250 × 4.5 µm LiChrospher 100 RP-18e

Mobile phase: A MeCN:water 50:50; B 20 mM pH 7 phosphate buffer; C MeCN:20 mM pH 3 phosphate buffer 10:90; D MeCN:200 mM pH 3.0 phosphate buffer 35:65 containing 2 mM disodium EDTA

Column temperature: 35

Flow rate: 1 (except where indicated)

Injection volume: 200

Detector: E, Merck Model L3500, glassy carbon working electrode +0.65 V, stainless-steel auxiliary electrode, Ag/AgCl reference electrode following post-column reaction. The column effluent flowed through a 10 m × 0.3 mm ID woven PTFE coil illuminated by a UV 254 low-pressure mercury lamp to the detector.

CHROMATOGRAM

Retention time: 7.1

Limit of detection: 2.7 ng

OTHER SUBSTANCES

Extracted: cloxacillin, dicloxacillin, penicillin V, penicillin G

KEY WORDS

post-column reaction; post-column photochemical derivatization; cow; muscle; column-switching

REFERENCE

Lihl,S.; Rehorek,A.; Petz,M. High-performance liquid chromatographic determination of penicillins by means of automated solid-phase extraction and photochemical degradation with electrochemical detection, *J.Chromatogr.A*, **1996**, 729, 229–235.

SAMPLE

Matrix: urine

Sample preparation: Filter (0.45 µm), inject a 5 µL aliquot.

HPLC VARIABLES

Guard column: 50 × 5 LiChrosorb RP-2

Column: 250 × 4.6 LiChrosorb RP-18

Mobile phase: MeOH:30 mM pH 5.6 acetate buffer 33:66

Flow rate: 1.5

Injection volume: 5

Detector: UV 254

CHROMATOGRAM

Retention time: 28

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

rat; human; pharmacokinetics

REFERENCE

Murai,Y.; Nakagawa,T.; Yamaoka,K.; Uno,T. High performance liquid chromatographic analysis and pharmacokinetic investigation of oxacillin and its metabolites in man, *Chem.Pharm.Bull.(Tokyo)*, **1981**, 29, 3290–3297.