


**HPLC VARIABLES****Guard column:** 50 × 6.4 25-37 μm Whatman Co-Pell ODS**Column:** 250 × 4.6 5 μm Ultrasphere ODS**Mobile phase:** MeOH:200 mM ammonium acetate buffer:water 55:10:35**Flow rate:** 1**Injection volume:** 20**Detector:** UV 254**CHROMATOGRAM****Retention time:** 15.2**Internal standard:** nitroglycerin (8)**OTHER SUBSTANCES****Simultaneous:** isosorbide mononitrate, saccharin, isosorbide dinitrate**KEY WORDS**

tablets; capsules

**REFERENCE**

Carlson,M.; Thompson,R.D.; Snell,R.P. Determination of isosorbide dinitrate in pharmaceutical products by HPLC, *J.Chromatogr.Sci.*, **1988**, *26*, 574–578.

# Pentamidine

**Molecular formula:** C<sub>19</sub>H<sub>24</sub>N<sub>4</sub>O<sub>2</sub>**Molecular weight:** 340.43**CAS Registry No.:** 100-33-4, 6823-79-6 (dimethanesulfonate), 140-64-7 (isethionate)**Merck Index:** 7254**SAMPLE****Matrix:** blood

**Sample preparation:** 1 mL Serum + 50 μL 4.92 μg/mL hexamidine, vortex briefly, add 500 μL 2 M NaOH, vortex, add 500 μL 2 M HCl, vortex, add 1 mL pH 10 carbonate buffer, vortex, add 4 mL 20 mM di(2-ethylhexyl)phosphoric acid in chloroform, vortex for 1 min, centrifuge at 700 g for 15 min. Remove the chloroform layer and add it to 1 mL 20 mM HCl, vortex for 1 min. Remove the aqueous layer and adjust the pH to 12 with 4 drops 2 M NaOH, add 2 mL dichloromethane, vortex for 1 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 35°, reconstitute the residue in 100 μL solvent, vortex for 30 s, inject a 20 μL aliquot. (Solvent was MeOH:buffer 50:50. Buffer was 50 mM sodium heptanesulfonate containing 0.4% triethylamine, pH adjusted to 3.0 with phosphoric acid.)

**HPLC VARIABLES****Column:** 150 × 2.1 5 μm solvent miser C18 (Alltech)**Mobile phase:** MeOH:buffer 60:40 (Buffer was 50 mM sodium heptanesulfonate containing 14 mM triethylamine, pH adjusted to 3.0 with phosphoric acid.)**Flow rate:** 0.3**Injection volume:** 20**Detector:** UV 280**CHROMATOGRAM****Retention time:** 5.7**Internal standard:** hexamidine (8.2)**Limit of detection:** 5 ng/mL**KEY WORDS**

serum; dog; human; pharmacokinetics

---

**REFERENCE**

Dickinson, C.M.; Navin, T.R.; Churchill, F.C. High-performance liquid chromatographic method for quantitation of pentamidine in blood serum, *J. Chromatogr.*, **1985**, *345*, 91–97.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Condition an SPE-C8 SPE cartridge (Jones Chromatography) with two 1 mL aliquots of MeOH, two 1 mL aliquots of MeCN:1 M pH 3 ammonium acetate 75:25, and five 1 mL aliquots of water, do not allow to dry. 100  $\mu$ L Plasma + 100  $\mu$ L water, mix, add to the SPE cartridge, wash with three 1 mL aliquots of MeOH, add 80  $\mu$ L 500 ng/mL melphalan in MeOH to the SPE cartridge, elute with 1 mL MeCN:1 M pH 3 ammonium acetate 75:25. Evaporate the eluate to dryness under reduced pressure, reconstitute in 200  $\mu$ L mobile phase, inject a 50  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 150  $\times$  3.9 5  $\mu$ m Resolve C18 (Waters)

**Mobile phase:** MeCN:MeOH:triethylamine:200 mM ammonium acetate 18:2:0.5:79.5, pH adjusted to 3.8 with acetic acid

**Flow rate:** 1.4

**Injection volume:** 50

**Detector:** F ex 265 em 345

---

**CHROMATOGRAM**

**Retention time:** 5.5

**Internal standard:** melphalan (10)

**Limit of detection:** 8.6 ng/mL

---

**KEY WORDS**

plasma; rat; SPE; pharmacokinetics

---

**REFERENCE**

Yeh, T.-K.; Dalton, J.T.; Au, J.L.-S. High-performance liquid chromatographic determination of pentamidine in plasma, *J. Chromatogr.*, **1993**, *622*, 255–261.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Let blood stand at 4° for 4 h, centrifuge at 15000 g for 10 min, add sodium azide to a concentration of 0.01%, inject a 25  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu$ m Zorbax RX-C8

**Mobile phase:** Gradient. A was 4.2 mM phosphoric acid containing 10 mM sodium heptanesulfonate and 10 mM tetramethylammonium chloride. B was MeCN:water 75:25 containing 4.2 mM phosphoric acid, 10 mM sodium heptanesulfonate, and 10 mM tetramethylammonium chloride. A:B from 90:10 to 10:90 over 30 min, return to initial conditions over 7 min, re-equilibrate for 3 min.

**Flow rate:** 1.5

**Injection volume:** 25

**Detector:** UV 265

---

**CHROMATOGRAM**

**Retention time:** 18

**Internal standard:** pentamidine

---

**OTHER SUBSTANCES**

**Extracted:** guanyldiazones

---

**KEY WORDS**

serum; mouse; pentamidine is IS; rat

---

**REFERENCE**

Cerami,C.; Zhang,X.; Ulrich,P.; Bianchi,M.; Tracey,K.J.; Berger,B.J. High-performance liquid chromatographic method for guanyldiazone compounds, *J.Chromatogr.B*, **1996**, 675, 71–75.

---

**SAMPLE**

**Matrix:** blood, broncho-alveolar lavage fluid, cells, urine

**Sample preparation:** Plasma. 0.5 mL Plasma + 1 mL 30 ng/mL hexamidine in MeCN, vortex, centrifuge at 1000 g for 10 min, add the supernatant to a Bond Elut C8 SPE cartridge, wash with 1 mL water, wash with 1 mL MeOH:water 50:50, wash with 1 mL MeOH, elute with MeOH:water:sodium heptanesulfonate:tetramethylammonium chloride:phosphoric acid 97.2:2.2:0.5:0.02:0.1, evaporate the eluate to 200  $\mu$ L under a stream of nitrogen, inject a 50  $\mu$ L aliquot. Bronchoalveolar lavage fluid. Centrifuge bronchoalveolar lavage fluid, add 1 mL supernatant to 1 mL Sorensen's buffer and 1 mL 30 ng/mL hexamidine in MeCN, vortex, centrifuge, add the supernatant to a Bond Elut C8 SPE cartridge, wash with water, wash with MeOH:water 50:50, wash with MeOH, elute with MeOH:water:sodium heptanesulfonate:tetramethylammonium chloride:phosphoric acid 97.2:2.2:0.5:0.02:0.1, evaporate the eluate to 200  $\mu$ L under a stream of nitrogen, inject a 50  $\mu$ L aliquot. Alveolar cells. Wash alveolar cells (from 10 mL bronchoalveolar lavage fluid) twice with phosphate buffered saline, resuspend in 1 L Sorensen's buffer, vortex for 1 min. 250  $\mu$ L Suspension + 250  $\mu$ L Sorensen's buffer + 1 mL 30 ng/mL hexamidine in MeCN, vortex, centrifuge, add the supernatant to a Bond Elut C8 SPE cartridge, wash with water, wash with MeOH:water 50:50, wash with MeOH, elute with MeOH:water:sodium heptanesulfonate:tetramethylammonium chloride:phosphoric acid 97.2:2.2:0.5:0.02:0.1, evaporate the eluate to 200  $\mu$ L under a stream of nitrogen, inject a 50  $\mu$ L aliquot. Urine. 200  $\mu$ L Urine + 1 mL 750 ng/mL hexamidine in MeCN, vortex for 1 min, centrifuge at 1000 g for 5 min, inject a 20  $\mu$ L aliquot (*J.Infect.Dis.* 1986, 154, 823).

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu$ m Ultrasphere Octyl 5

**Mobile phase:** MeCN:water 21:79 containing 0.02% tetramethylammonium chloride and 0.1% phosphoric acid

**Flow rate:** 1

**Injection volume:** 50

**Detector:** F ex 275 em 340

---

**CHROMATOGRAM**

**Internal standard:** hexamidine

**Limit of detection:** 2.29 ng/mL

---

**KEY WORDS**

plasma; SPE; pharmacokinetics

---

**REFERENCE**

Conte,J.E.,Jr.; Golden,J.A. Intrapulmonary and systemic pharmacokinetics of aerosolized pentamidine used for prophylaxis of *Pneumocystis carinii* pneumonia in patients infected with the human immunodeficiency virus, *J.Clin.Pharmacol.*, **1995**, 35, 1166–1173.

---

**SAMPLE**

**Matrix:** cells

**Sample preparation:** 400  $\mu$ L Cells in phosphate saline glucose (99:1) buffer + 10  $\mu$ L 10% orthophosphoric acid, vortex, filter (Ultrafree-MC polysulfone membrane, 100000 molecular mass cut-off) while centrifuging at 5000 g for 5 min, inject a 20  $\mu$ L aliquot of the filtrate.

---

**HPLC VARIABLES**

**Guard column:** 20  $\times$  2 30-40  $\mu$ m pellicular Spherisorb RP-18

**Column:** 200  $\times$  2 5  $\mu$ m Nucleosil C18

**Mobile phase:** MeOH:5 mM citric acid 50:50 containing 5 mM sodium pentanesulfonate, pH adjusted to 4.0 with NaOH

**Flow rate:** 0.5

**Injection volume:** 20

**Detector:** UV 261

---

**CHROMATOGRAM****Retention time:** 6**Limit of detection:** 5.7 ng/mL

---

**REFERENCE**

Rabanal,B.; De Arriba,R.G.; Garzón,M.J.; Reguera,R.M.; Balaña-Fouce,R.; Negro,A. Determination of pentamidine in *Leishmania* infantum promastigotes by ion-paired liquid chromatography, *J.Liq.Chromatogr.*, **1994**, 17, 2017–2029.

---

**SAMPLE****Matrix:** microsomal incubations

**Sample preparation:** 6 mL microsomal incubation + 500  $\mu$ L MeCN + hexamidine, add to a Prep-Sep C18 SPE cartridge (Fisher), wash with water, wash with MeCN, elute with 1 mL MeCN: buffer 95:5. Evaporate the eluate to near dryness under a stream of air, reconstitute the residue in 200  $\mu$ L mobile phase, inject a 5  $\mu$ L aliquot. (Buffer was 4.2 mM phosphoric acid containing 10 mM heptanesulfonate and 10 mM tetramethylammonium chloride.)

---

**HPLC VARIABLES****Column:** 250 X 4.6 Zorbax 5  $\mu$ m RX diisopropyl C8

**Mobile phase:** Gradient. MeCN:buffer from 22.5:77.5 to 45:55 over 25 min (Buffer was 4.2 mM phosphoric acid containing 10 mM heptanesulfonate and 10 mM tetramethylammonium chloride.)

**Column temperature:** 40**Injection volume:** 5**Detector:** UV 265

---

**CHROMATOGRAM****Retention time:** 16**Internal standard:** hexamidine

---

**OTHER SUBSTANCES****Extracted:** metabolites

---

**KEY WORDS**rat; liver; SPE

---

**REFERENCE**

Berger,B.J.; Reddy,V.V.; Le,S.T.; Lombardy,R.J.; Hall,J.E.; Tidwell,R.R. Hydroxylation of pentamidine by rat liver microsomes, *J.Pharmacol.Exp.Ther.*, **1991**, 256, 883–889.

---

**SAMPLE****Matrix:** microsomal incubations

**Sample preparation:** Lyophilize microsomal incubation. Dissolve lyophilizate in 150  $\mu$ L mobile phase, mix thoroughly for 5 min, centrifuge at 6000 g for 5 min, inject a 15  $\mu$ L aliquot of the supernatant.

---

**HPLC VARIABLES****Guard column:** 4  $\times$  4.5  $\mu$ m Lichrospher 60 RP-select B**Column:** 125  $\times$  4.5  $\mu$ m Lichrospher 60 RP-select B

**Mobile phase:** Gradient. A was 10 mM sodium octylsulfonate, 10 mM tetramethylammonium chloride, and 20 mM phosphoric acid adjusted to pH 3.0 with ammonia. B was MeOH. A:B from 60:40 to 30:70 over 30 min.

**Flow rate:** 1**Injection volume:** 15**Detector:** UV 260

---

**CHROMATOGRAM****Retention time:** 22**Limit of detection:** 0.02 nmole

---

**OTHER SUBSTANCES****Extracted:** metabolites**KEY WORDS**

human; rabbit; liver

**REFERENCE**

Clement,B.; Jung,F. *N*-Hydroxylation of the antiprotozoal drug pentamidine catalyzed by rabbit liver cytochrome P-450 2C3 or human liver microsomes, microsomal retroreduction, and further oxidative transformation of the formed amidoximes Possible relationship to the biological oxidation of arginine to NG-hydroxyarginine, citrulline, and nitric oxide, *Drug Metab.Dispos.*, **1994**, 22, 486–497.


**SAMPLE****Matrix:** urine**Sample preparation:** Directly inject an aliquot of urine. Alternatively add urine to a Sep-Pak SPE cartridge, elute with 2 mL MeCN:water 50:50, inject an aliquot of the eluate.**HPLC VARIABLES****Column:** 150 × 4.6 5 µm Ultrasphere ODS**Mobile phase:** Gradient. A was MeCN:50 mM pH 6 ammonium acetate buffer:triethylamine 3:96.8:0.2. B was MeCN:50 mM pH 6 ammonium acetate buffer:triethylamine 50:49.8:0.2. A:B from 100:0 to 0:100 in 60 min (?)**Injection volume:** 100-150**Detector:** MS, Finnigan TSQ 700, chemical ionization APCI interface, vaporizer 400°, heated capillary 200°, current 5 µA, m/z 357 (The effluent from the column was directed to waste for 2 min, then to the MS.)**OTHER SUBSTANCES****Extracted:** metabolites**KEY WORDS**

human; rat; SPE

**REFERENCE**

Nordin,J.; Wilkström,I.; Bronner,U.; Gustafsson,L.L.; Ericsson,. Liquid chromatography-tandem mass spectrometry applied to a study of the metabolism of pentamidine. Discussion of possibilities and problems, *J.Chromatogr.A*, **1997**, 777, 73–79.

# Pentazocine

**Molecular formula:** C<sub>19</sub>H<sub>27</sub>NO**Molecular weight:** 285.45**CAS Registry No.:** 359-83-1, 64024-15-3 (HCl), 17146-95-1 (lactate)**Merck Index:** 7261**Lednicer No.:** 1 297; 2 325**SAMPLE****Matrix:** blood**Sample preparation:** 1 mL Serum + 200 ng doxepin or desipramine + 100 µL 1 M NaOH + 9 mL freshly prepared hexane:isoamyl alcohol 99:1, shake vigorously for 5 min, centrifuge. Remove 8.5 mL of the organic phase and add it to 200 µL 50 mM HCl, shake well for 1 min, centrifuge, inject a 50 µL aliquot of the aqueous phase.**HPLC VARIABLES****Column:** 300 × 4 µBondapak phenyl**Mobile phase:** MeCN:0.01% phosphoric acid containing 0.01% NaCl 35:65, final pH 2.8

**Flow rate:** 1.5  
**Injection volume:** 50  
**Detector:** UV 210

---

#### CHROMATOGRAM

**Retention time:** 7.4  
**Internal standard:** doxepin (12.2), desipramine (14.2)  
**Limit of detection:** 10 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** cocaine, dextromoramide, meperidine, methadone, normeperidine, norpropoxyphene, propoxyphene  
**Simultaneous:** amitriptyline, buprenorphine, chlorpromazine, codeine, desmethyldoxepin, diphenhydramine, ephedrine, imipramine, nortriptyline, oxazepam, oxycodone, pericyazine, pheniramine, propranolol, quinine, thiopropazate, thioridazine

---

#### KEY WORDS

serum

---

#### REFERENCE

Hackett,L.P.; Dusci,L.J.; Ilett,K.F. The analysis of several nonopiate narcotic analgesics and cocaine in serum using high-performance liquid chromatography, *J.Anal.Toxicol.*, **1987**, 11, 269-271.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 200  $\mu$ L 1 M NaOH + 6 mL 20 ng/mL levallorphan tartrate in diethyl ether, agitate at 4° for 15 min. Remove 5 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, dissolve the residue in 250-500  $\mu$ L mobile phase, inject a 50  $\mu$ L aliquot.

---

#### HPLC VARIABLES

**Column:** 250  $\times$  4.5 5  $\mu$ m Nucleosil RP 18  
**Mobile phase:** MeCN:5 mM phosphoric acid 33:67  
**Flow rate:** 1  
**Injection volume:** 50  
**Detector:** F ex 278 em 324

---

#### CHROMATOGRAM

**Retention time:** 6.2  
**Internal standard:** levallorphan tartrate (4.9)  
**Limit of detection:** 1 ng/mL  
**Limit of quantitation:** 4 ng/mL

---

#### KEY WORDS

plasma

---

#### REFERENCE

Moeller,N.; Dietzel,K.; Nuernberg,B.; Geisslinger,G.; Brune,K. High-performance liquid chromatographic determination of pentazocine in plasma, *J.Chromatogr.*, **1990**, 530, 200-205.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** Rock 5 mL whole blood + 10 mL water + 8.5 mL Na<sub>2</sub>WO<sub>4</sub> in a 50 mL stoppered tube for 1 min, add 6 mL NiCl<sub>2</sub>, rock for 5 min, add 15 mL dichloromethane:isobutyl alcohol:THF 30:45:25, centrifuge at 2500 g for 15 min. Remove organic phase and repeat the process. Filter all organic phases through a 40-90  $\mu$ m filter and evaporate to dryness in a 100 mL porcelain dish at a moderate temperature in a sand bath. Take up residue in 500  $\mu$ L MeCN: water 80:20, inject a 20  $\mu$ L aliquot. (Na<sub>2</sub>WO<sub>4</sub> prepared by mixing 10 g Na<sub>2</sub>WO<sub>4</sub>·2H<sub>2</sub>O in 38 mL of 2 M NaOH and 2.5 g of NaHCO<sub>3</sub> and making up to 100 mL. NiCl<sub>2</sub> was 17% w/v NiCl<sub>2</sub> in water.)

---

**HPLC VARIABLES**

**Column:** 200 × 4.6 5 µm Hypersil C8

**Mobile phase:** Gradient A = MeCN; B = 20 mM n-propylamine adjusted to pH 5 with 85% phosphoric acid. A:B from 15:85 to 20:80 over 5 min to 45:55 over another 15 min to 65:35 over another 5 min

**Injection volume:** 20

**Detector:** UV 230

---

**CHROMATOGRAM**

**Retention time:** 20

**Limit of detection:** 0.50 ppm

---

**OTHER SUBSTANCES**

**Extracted:** buprenorphine, caffeine, cocaine, codeine, diamorphine, ethylmorphine, lidocaine, methaqualone, morphine, naloxone, noscapine, papaverine, procaine

**Also analyzed:** bromazepam, clonazepam, diazepam, flunitrazepam, flurazepam, medazepam, nitrazepam, oxazepam

---

**KEY WORDS**

whole blood

---

**REFERENCE**

Bernal, J.L.; Del Nozal, M.J.; Rosas, V.; Villarino, A. Extraction of basic drugs from whole blood and determination by high performance liquid chromatography, *Chromatographia*, **1994**, 38, 617–623.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Automated SPE by ASPEC system. Condition a C18 Clean-Up SPE cartridge (CEC 18111, Worldwide Monitoring) with 2 mL MeOH then 2 mL water. 1 mL Plasma + 1 mL 400 ng/mL protriptyline in water, vortex, add to column, wash with 3 mL water, wash with 3 mL 750 mL/L methanol. Elute with three aliquots of 300 µL 0.1 M ammonium acetate in MeOH. Add 0.5 mL 0.5 M NaOH and 4 mL 50 mL/L isopropanol in heptane to eluate, mix thoroughly. Allow 5 min for phase separation. Remove upper heptane phase and add it to 300 µL 0.1 M phosphoric acid (pH 2.5), mix, separate, inject a 100 µL aliquot of the aqueous phase.

---

**HPLC VARIABLES**

**Guard column:** LC-8-DB (Supelco)

**Column:** 150 × 4.6 LC-8-DB (Supelco)

**Mobile phase:** MeCN:buffer 35:65 (Buffer was 10 mL/L triethylamine in water adjusted to pH 5.5 with glacial acetic acid.)

**Flow rate:** 2

**Injection volume:** 100

**Detector:** UV 228

---

**CHROMATOGRAM**

**Retention time:** 1.9

**Internal standard:** protriptyline (4)

---

**OTHER SUBSTANCES**

**Extracted:** acetazolamide, amitriptyline, chlordiazepoxide, chlorimipramine, chlorpromazine, desipramine, dextromethorphan, diazepam, diphenhydramine, doxepin, fentanyl, flecainide, fluoxetine, flurazepam, haloperidol, hydroxyethylflurazepam, ibuprofen, imipramine, maprotiline, methadone, methaqualone, mexiletine, midazolam, norchlorimipramine, nordiazepam, nordoxepin, norfluoxetine, nortriptyline, norverapamil, promazine, propafenone, propoxyphene, protriptyline, quinidine, temazepam, trazodone, trimipramine, verapamil

**Noninterfering:** acetaminophen, acetylmorphine, amiodarone, amobarbital, amphetamine, ben-droflumethiazide, benzocaine, benzoylecgonine, benzthiazide, butalbital, carbamazepine, chlorothiazide, clonazepam, cocaine, codeine, cotinine, cyclosporine, cyclothiazide, desalkylflurazepam, diamorphine, dicumerol, ephedrine, ethacrynic acid, ethanol, ethchlorvynol, ethosuximide, furosemide, glutethimide, hydrochlorothiazide, hydrocodone, hydroflumethiazide, hydromorphone, lorazepam, mephentermine, meprobamate, methamphetamine, metharbital, methoxsalen, methoxyphenteramine, methsuximide, methylcyclothiazide, metoprolol,

MHPG, monoacetylmorphine, morphine, normethsuximide, oxazepam, oxycodone, oxymorphone, pentobarbital, phencyclidine, phenteramine, phenylephrine, phenytoin, polythiazide, primidone, prochlorperazine, salicylic acid, sulfanilamide, THC-COOH, theophylline, thiazolam, thiopental, thioridazine, tocanide, trichloromethiazide, trifluoperazine, valproic acid, warfarin

**Interfering:** encainide, lidocaine, propranolol

---

## KEY WORDS

plasma; SPE

---

## REFERENCE

Nichols, J.H.; Charlson, J.R.; Lawson, G.M. Automated HPLC assay of fluoxetine and norfluoxetine in serum, *Clin. Chem.*, **1994**, *40*, 1312-1316.

---

## SAMPLE

**Matrix:** blood

**Sample preparation:** 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol:n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

---

## HPLC VARIABLES

**Column:** 300 × 3.9 4 µm NovaPack C18

**Mobile phase:** MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH<sub>2</sub>PO<sub>4</sub> adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

**Column temperature:** 30

**Flow rate:** 0.8

**Injection volume:** 50

**Detector:** UV 220

---

## CHROMATOGRAM

**Retention time:** 4.33

**Limit of detection:** <120 ng/mL

---

## KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylcegonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melfalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vandesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phencyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetrizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine;



bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozone; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; caripramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

## REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254-262.

## SAMPLE

**Matrix:** blood, gastric contents

**Sample preparation:** 1 mL Whole blood or gastric contents + 50  $\mu$ L 400  $\mu$ g/mL IS in MeOH + 1 mL EtOH + 5 drops 1 M pH 9 potassium carbonate + 2 mL water + 8 mL n-hexane:MTBE 25:75, rotate for 15 min, centrifuge at 2500 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 100  $\mu$ L mobile phase, inject a 50  $\mu$ L aliquot.

## HPLC VARIABLES

**Column:** 150  $\times$  4.6 5  $\mu$ m Hypersil BDS C18

**Mobile phase:** Gradient. A was MeCN:MeOH:1.5 M ammonium acetate:water 10:10:3:77. B was MeCN:MeOH:1.5 M ammonium acetate:water 40:40:3:17. A:B from 95:5 to 50:50 over 20 min.

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 230

## CHROMATOGRAM

**Retention time:** 15

**Internal standard:** N-allylnormetazocine (Sigma) (9)

## OTHER SUBSTANCES

**Extracted:** zopiclone

## KEY WORDS

whole blood

## REFERENCE

Van Bocxlaer,J.; Meyer,E.; Clauwaert,K.; Lambert,W.; Piette,M.; De Leenheer,A. Analysis of zopiclone (Imovane) in postmortem specimens by GC-MS and HPLC with diode-array detection, *J.Anal.Toxicol.*, **1996**, *20*, 52-54.

## SAMPLE

**Matrix:** blood, saliva, tissue, urine

**Sample preparation:** Homogenize (Polytron) tissue with 4 (whole brain) or 8 (brain striata) volumes of 100 mM pH 4.5  $\text{NaH}_2\text{PO}_4$  containing 0.5% NaF. Add 500  $\mu$ L brain homogenate or 500  $\mu$ L plasma, saliva, or urine containing 15  $\mu$ L saturated NaF solution to 75  $\mu$ L 150  $\mu$ g/mL IS, add 50  $\mu$ L 50% perchloric acid, mix vigorously for 10 s, let stand at room temperature for 10 min, add 1 mL water, mix briefly, centrifuge at 10° at 2500 (?) for 30 min. Remove the supernatant and add it to 750  $\mu$ L saturated sodium carbonate solution, mix briefly, add 7.5 mL pentane:chloroform 95:5, rock gently for 10 min, centrifuge in a desk-top centrifuge for 2 min, freeze in dry ice/acetone for 2 min. Remove the organic layer and add it to 250  $\mu$ L 100 mM HCl, mix vigorously for 10 s, centrifuge in a desk-top centrifuge for 1-2 min, freeze in dry ice/acetone for 3-5 min, discard the organic layer. Allow the aqueous layer to thaw, remove any trace of organic solvent with a stream of nitrogen, inject a 75  $\mu$ L aliquot of the aqueous layer.

## HPLC VARIABLES

**Guard column:** 15  $\times$  3.2 7  $\mu$ m Brownlee RP-8

**Column:** 250 × 4.6 5 µm Zorbax RX-C18

**Mobile phase:** MeCN:buffer 18:82 (Buffer was 100 mM K<sub>2</sub>HPO<sub>4</sub> containing 0.5% triethylamine, adjusted to pH 2.7 with phosphoric acid.)

**Flow rate:** 2

**Injection volume:** 75

**Detector:** UV 235

---

#### CHROMATOGRAM

**Retention time:** 11.5

**Internal standard:** 2β-carbomethoxy-3β-(4-chlorophenyl)tropane (RTI-31) (Research Biochemical International, Natick MA) (11.4)

---

#### OTHER SUBSTANCES

**Extracted:** chlordiazepoxide, clozapine, cocaine, gepirone, methylphenidate, pseudococaine

**Simultaneous:** acetaminophen, acetophenazine, amoxapine, amphetamine, atropine, benperidol, buspirone, caffeine, carbamazepine, chlorpheniramine, codeine, dextromethorphan, diazepam, diphenhydramine, flupenthixol, flurazepam, haloperidol, hydergine, hydrocodone, hydromorphone, lidocaine, loxapine, mepazine, meperidine, mesoridazine, methaqualone, 3,4-methylenedioxymethamphetamine, 3,4-methylenedioxyethylamphetamine, 3,4-methylenedioxymethamphetamine, morphine, norcocaine, oxazepam, pentobarbital, phenylpropanolamine, procainamide, procaine, propyl benzoylecgonine, quinidine, quinine, salicylic acid, secobarbital, theophylline, trazodone, 3-tropanyl-3,5-dichlorobenzoate, vancomycin, WIN 35428

**Noninterfering:** amitriptyline, benztrapine methanesulfonate, butaperazine, butriptyline, carphenazine, chlorpromazine, clomipramine, cyclobenzaprine, dextropropoxyphene, dronabinol, ephedrine, ethchlorvynol, fluoxetine, fluphenazine, imipramine, meprobamate, methadone, methamphetamine, nicotine, norfluoxetine, nortriptyline, PCP, phenothiazine, pseudoephedrine

---

#### KEY WORDS

rat; cow; plasma; brain

---

#### REFERENCE

Bonate,P.L.; Davis,C.M.; Silverman,P.B.; Swann,A. Determination of cocaine in biological matrices using reversed phase HPLC: Application to plasma and brain tissue, *J.Liq.Chromatogr.*, **1995**, 18, 3473–3494.

---

#### SAMPLE

**Matrix:** blood, tissue, vitreous humor

**Sample preparation:** Blood, vitreous humor. Mix 1 mL sample with 500 µL 2% pH 9.5 sodium tetraborate and 8 mL n-butanol:hexane 5:95. Extract for 30 min, centrifuge. Remove the organic layer and add it to 200 µL 0.2% orthophosphoric acid. Extract for 30 min, inject a 30 µL aliquot of the aqueous layer. Tissue. Mix 500 µL liver homogenate with 500 µL 2% pH 9.5 sodium tetraborate and 8 mL n-butanol:hexane 5:95. Extract for 30 min, centrifuge. Remove the organic layer and add it to 400 µL 0.2% orthophosphoric acid. Extract for 30 min, inject a 30 µL aliquot of the aqueous layer.

---

#### HPLC VARIABLES

**Column:** 150 × 3.9 5 µm NovaPak-Phenyl

**Mobile phase:** MeCN:10 mM KH<sub>2</sub>PO<sub>4</sub> 55:45, adjusted to pH 3.0

**Flow rate:** 1.5

**Injection volume:** 30

**Detector:** UV 214

---

#### OTHER SUBSTANCES

**Extracted:** pimozone, sertraline

---

#### KEY WORDS

liver; pentazocine is IS

---

#### REFERENCE

McIntyre,I.M.; King,C.V.; Staikos,V.; Gall,J.; Drummer,O.H. A fatality involving moclobemide, sertraline, and pimozone, *J.Forensic Sci.*, **1997**, 42, 951–953.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50  $\mu$ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood)  $\mu$ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18

**Column:** 250  $\times$  4.6 5  $\mu$ m Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 200.5

---

**CHROMATOGRAM**

**Retention time:** 12.522

---

**KEY WORDS**

whole blood

---

**REFERENCE**

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

---

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Tablets. Powder tablets, weigh out amount equivalent to about 10 mg, add 75 mL mobile phase, sonicate for 20 min, dilute to 100 mL with mobile phase, mix, filter (0.45  $\mu$ m) (discard first 10 mL of filtrate), inject a 20  $\mu$ L aliquot of the filtrate. Syrups, elixirs, injectables. Measure out amount equivalent to about 10 mg, add 75 mL mobile phase, sonicate for 20 min, dilute to 100 mL with mobile phase, mix, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 300  $\times$  3.9 10  $\mu$ m  $\mu$ Bondapak CN

**Mobile phase:** MeOH:3 mM ammonium acetate 90:10

**Flow rate:** 1.3

**Injection volume:** 20

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 5.3

---

**OTHER SUBSTANCES**

**Also analyzed:** chlorpheniramine, cyclizine, doxylamine, mesoridazine, promethazine, protriptyline, pyrilamine, pyrimethamine, tripeleminamine

---

**KEY WORDS**

tablets; syrups; elixirs; injections

---

**REFERENCE**

Walker, S.T. Liquid chromatographic determination of organic nitrogenous bases in dosage forms: a progress report, *J. Assoc. Off. Anal. Chem.*, **1985**, 68, 539-542.

---

**SAMPLE****Matrix:** solutions**Sample preparation:** Dissolve in MeOH at a concentration of 1 mg/mL, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES****Column:** 250  $\times$  5 Spherisorb S5W**Mobile phase:** MeOH:buffer 90:10 (Buffer was 94 mL 35% ammonia and 21.5 mL 70% nitric acid in 884 mL water, adjust the pH to 10.1 with ammonia.)**Flow rate:** 2**Injection volume:** 20**Detector:** UV 254

---

**CHROMATOGRAM****Retention time:** 2.21

---

**OTHER SUBSTANCES**

**Simultaneous:** thebacon, oxycodone, thebaine, norlevorphanol, methadone, benzylmorphine, ethylmorphine, morphine-N-oxide, codeine, codeine-N-oxide, morphine, ethoheptazine, morphine-3-glucuronide, pholcodeine, norpethidine, hydrocodone, dihydrocodeine, dihydromorphine, levorphanol, norcodeine, normorphine, pemoline, benzphetamine, diethylpropion, mazindol, tranlylcypromine, caffeine, fenethyline, phendimetrazine, methylphenidate, phenelzine, norpseudoephedrine, phentermine, fenfluramine, methylenedioxyamphetamine, amphetamine, normetanephine, 4-hydroxyamphetamine, bromo-STP, STP, prolintane, 2-phenethylamine, tyramine, trimethoxyamphetamine, phenylephrine, pseudoephedrine, ephedrine, methylephedrine, dimethylamphetamine, methamphetamine, mescaline, mephentermine, buprenorphine, dextromoramide, phenoperidine, fentanyl, etorphine, piritramide, noscapine, papaverine, naloxone, dextropropoxyphene, nalorphine, phenazocine, norpipanone, levallorphan, hydroxypethidine

**Noninterfering:** dopamine, levodopa, methylodpa, methylodopate, norepinephrine

**Interfering:** epinephrine, pipradol, phenylpropanolamine, fencamfamin, chlorphentermine, normethadone, meperidine, dipipanone, diamorphine, acetylcodeine, monoacetylmorphine

---

**REFERENCE**

Law,B.; Gill,R.; Moffat,A.C. High-performance liquid chromatography retention data for 84 basic drugs of forensic interest on a silica column using an aqueous methanol eluent, *J.Chromatogr.*, **1984**, *301*, 165-172.

---

**SAMPLE****Matrix:** solutions**Sample preparation:** Prepare a 10  $\mu$ g/mL solution in MeOH, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES****Column:** 125  $\times$  4.9 Spherisorb S5W silica**Mobile phase:** MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7**Flow rate:** 2**Injection volume:** 20**Detector:** E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

---

**CHROMATOGRAM****Retention time:** 2.4

---

**OTHER SUBSTANCES**

**Also analyzed:** acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzoctamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipa-

none, diprenorphine, dipyridamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserin, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypropazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxymbenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, propidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

## REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, *323*, 191-225.

## SAMPLE

**Matrix:** solutions

## HPLC VARIABLES

**Column:** 250 × 4.6 Zorbax RX

**Mobile phase:** Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

**Column temperature:** 30

**Flow rate:** 2

**Detector:** UV 210

## OTHER SUBSTANCES

**Also analyzed:** acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbomal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenopropfen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxinigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiacol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-

stillbene, imipramine, indomethacin, isocarboxystiril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephénytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methylpyrrol, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidone, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethiodole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

---

## REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, **1994**, *18*, 233–242.

---

## SAMPLE

**Matrix:** solutions

---

## HPLC VARIABLES

**Column:** 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4.5 µm LiChrospher 100 RP-8 (B)

**Mobile phase:** MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

**Flow rate:** 0.6

**Injection volume:** 25

**Detector:** UV 229

---

## CHROMATOGRAM

**Retention time:** 9.91 (A), 5.47 (B)

---

## OTHER SUBSTANCES

**Also analyzed:** acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordi-azepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxy-chloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazindol, mefenamic acid, meperidine, mephénytoin, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, met-

ronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfapyrazole, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocainide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

## KEY WORDS

details of plasma extraction

## REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, *692*, 103–119.

## SAMPLE

**Matrix:** urine

**Sample preparation:** 500  $\mu$ L Urine + N-ethylnordiazepam + chlorpheniramine + 100  $\mu$ L buffer, centrifuge at 11000 g for 30 s, inject a 500  $\mu$ L aliquot onto column A with mobile phase A, after 0.6 min backflush column A with mobile phase A to waste for 1.6 min, elute column A with 250  $\mu$ L mobile phase B, with 200  $\mu$ L mobile phase C, and with 1.15 mL mobile phase D. Elute column A to waste until drugs start to emerge then elute onto column B. Elute column B to waste until drugs started to emerge, then elute onto column C. When all the drugs have emerged from column B remove it from the circuit, elute column C with mobile phase D, monitor the effluent from column C. Flush column A with 7 mL mobile phase E, with mobile phase D, and mobile phase A. Flush column B with 5 mL mobile phase E then with mobile phase D. (Buffer was 6 M ammonium acetate adjusted to pH 8.0 with 2 M KOH.)

## HPLC VARIABLES

**Column:** A  $10 \times 2.1$  12-20  $\mu$ m PRP-1 spherical poly(styrene-divinylbenzene) (Hamilton); B  $10 \times 3.2$  11  $\mu$ m Aminex A-28 (Bio-Rad); C  $25 \times 3.2$  5  $\mu$ m C8 (Phenomenex) +  $150 \times 4.6$  5  $\mu$ m silica (Macherey-Nagel)

**Mobile phase:** A 0.1% pH 8.0 potassium borate buffer; B 6 mM  $\text{KH}_2\text{PO}_4$  containing 5 mM tetramethylammonium hydroxide, and 2 mM dimethyloctylamine, pH adjusted to 6.50 with phosphoric acid; C MeCN:buffer 40:60 (Buffer was 6 mM  $\text{KH}_2\text{PO}_4$  containing 5 mM tetramethylammonium hydroxide, and 2 mM dimethyloctylamine, pH adjusted to 6.50 with phosphoric acid.); D MeCN:buffer 33:67 (Buffer was 6 mM  $\text{KH}_2\text{PO}_4$  containing 5 mM tetramethylammonium hydroxide, and 2 mM dimethyloctylamine, pH adjusted to 6.50 with phosphoric acid.); E MeCN:buffer 70:30 (Buffer was 6 mM  $\text{KH}_2\text{PO}_4$  containing 5 mM tetramethylammonium hydroxide, and 2 mM dimethyloctylamine, pH adjusted to 6.50 with phosphoric acid.)

**Column temperature:** ambient (column A), 40 (columns B and C)

**Flow rate:** A 5; B-E 1

**Injection volume:** 500

**Detector:** UV 210, UV 235

## CHROMATOGRAM

**Retention time:** k' 2.8

**Internal standard:** N-ethylnordiazepam (k' 2.1), chlorpheniramine (k' 5.9)

**Limit of detection:** 300 ng/mL

## OTHER SUBSTANCES

**Extracted:** nortriptyline, diphenhydramine, methadone, imipramine, flurazepam, amitriptyline, morphine, codeine, hydromorphone, hydrocodone, caffeine, cotinine, benzoylecgonine, secobarbital, oxazepam, phenobarbital, nordiazepam, diazepam, phenylpropanolamine, phentermine

**Interfering:** amphetamine, phenmetrazine, lidocaine, ephedrine, methamphetamine, desipramine

**KEY WORDS**

column-switching

**REFERENCE**

Binder, S.R.; Regalia, M.; Biaggi-McEachern, M.; Mazhar, M. Automated liquid chromatographic analysis of drugs in urine by on-line sample cleanup and isocratic multi-column separation, *J. Chromatogr.*, **1989**, 473, 325-341.


# Penthienate bromide

**Molecular formula:** C<sub>18</sub>H<sub>30</sub>BrNO<sub>3</sub>S

**Molecular weight:** 420.41

**CAS Registry No.:** 60-44-6

**Merck Index:** 7268

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.

**HPLC VARIABLES**

**Column:** 125 × 4.9 Spherisorb S5W silica

**Mobile phase:** MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

**Flow rate:** 2

**Injection volume:** 20

**Detector:** E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

**CHROMATOGRAM**

**Retention time:** 3.9

**OTHER SUBSTANCES**

**Also analyzed:** acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazepine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipranone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxylbenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, pri-



maquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleppamine, triprolidine, tryptamine, verapamil, xylometazoline

## REFERENCE

Jane,I.; McKinnon,A.; Flanagan,R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J.Chromatogr.*, **1985**, 323, 191-225.

# Pentobarbital


**Molecular formula:**  $C_{11}H_{18}N_2O_3$

**Molecular weight:** 226.28

**CAS Registry No.:** 76-74-4, 57-33-0 (Na salt)

**Merck Index:** 7272

**Lednicer No.:** 1 269



## SAMPLE

**Matrix:** blood

**Sample preparation:** 200  $\mu$ L Serum + 200  $\mu$ L 50  $\mu$ g/mL hexobarbital in MeCN + 25  $\mu$ L glacial acetic acid, vortex for 10 s, centrifuge for 1 min, inject a 30-100  $\mu$ L aliquot of the supernatant.

## HPLC VARIABLES

**Column:**  $\mu$ Bondapak C18

**Mobile phase:** Gradient. MeCN:7.5 g/L  $NaH_2PO_4$  adjusted to pH 3.2 with phosphoric acid 5:95 to 22:78 over 24 min, to 45:55 over 10 min, maintain at 45:55 for 5 min. Re-equilibrate with 5:95 for 5 min.

**Column temperature:** 50

**Flow rate:** 3

**Injection volume:** 30-100

**Detector:** UV 210

## CHROMATOGRAM

**Retention time:** 22.4

**Internal standard:** hexobarbital (20.6)

**Limit of detection:** 200-2000 ng/mL

## OTHER SUBSTANCES

**Extracted:** acetaminophen, amobarbital, butabarbital, butalbital, chlordiazepoxide, diazepam, ethchlorvynol, flurazepam, glutethimide, methaqualone, methypyrrol, nitrazepam, phenobarbital, phenytoin, primidone, salicylic acid, secobarbital, theophylline

**Simultaneous:** amitriptyline, caffeine, clomipramine, codeine, desipramine, ethotoin, imipramine, lidocaine, mesantoin, methsuximide, nirvanol, nortriptyline, oxazepam, procainamide, phenylpropanolamine, propranolol, quinidine

## KEY WORDS

serum

## REFERENCE

Kabra,P.M.; Stafford,B.E.; Marton,L.J. Rapid method for screening toxic drugs in serum with liquid chromatography, *J.Anal.Toxicol.*, **1981**, 5, 177-182.

---

**SAMPLE****Matrix:** blood**Sample preparation:** 2 mL Serum + 1 mL buffer, vortex, add 10 mL n-butyl chloride containing 10 µg/mL barbital and 4 µg/mL thiamylal, extract vigorously for 3 min, centrifuge at 3000 g for 5 min. Remove the upper organic layer and add it to 100 µL 450 mM NaOH, extract vigorously for 3 min, centrifuge for 10 min or until lower aqueous phase is clear, inject a 15 µL aliquot of the lower aqueous phase. (Soak glassware in 1 M HCl overnight, rinse with water, dry. Buffer was prepared from equal volumes of 22 g/L KH<sub>2</sub>PO<sub>4</sub> and 18 g/L Na<sub>2</sub>HPO<sub>4</sub>, pH 6.6 ± 0.2.)

---

**HPLC VARIABLES****Column:** 125 × 4.6 5 µm C-18 (Perkin-Elmer)**Mobile phase:** MeOH:THF:buffer 50:7:43 (Buffer was prepared from equal volumes of 22 g/L KH<sub>2</sub>PO<sub>4</sub> and 18 g/L Na<sub>2</sub>HPO<sub>4</sub>, pH 6.6 ± 0.2.)**Flow rate:** 2**Injection volume:** 6**Detector:** UV 240

---

**CHROMATOGRAM****Retention time:** 2.4**Internal standard:** barbital (0.8), thiamylal (5.2)**Limit of detection:** 2 µg/mL

---

**OTHER SUBSTANCES****Extracted:** thiopental**Simultaneous:** acetaminophen, acetazolamide, aspirin, butabarbital, cefazolin, chlorothiazide, chlorzoxazone, dicloxacillin, ethosuximide, furosemide, hydrochlorothiazide, ibuprofen, oxacillin, phenobarbital, phenylbutazone, phenytoin, salicylic acid, secobarbital, sulfamethoxazole, theophylline, ascorbic acid**Noninterfering:** ampicillin, penicillin G, valproic acid**Interfering:** amobarbital

---

**KEY WORDS**

serum

---

**REFERENCE**Kelner,M.; Bailey,D.N. Reversed-phase liquid-chromatographic simultaneous analysis for thiopental and pentobarbital in serum, *Clin.Chem.*, **1983**, 29, 1097–1100.

---

**SAMPLE****Matrix:** blood**Sample preparation:** 50 µL 25 µg/mL 5-Ethyl-5-p-tolylbarbituric acid in MeOH added to 150 × 10 mm glass centrifuge tube and blow dry under a stream of nitrogen, add 500 µL plasma, add 5 mL dichloromethane, mix on rotary mixer for 5 min, centrifuge. Remove organic layer, evaporate to dryness under nitrogen, take up in 500 µL mobile phase, inject a 40 µL aliquot.

---

**HPLC VARIABLES****Guard column:** 50 × 4 µm Bondapak C18 Corasil B**Column:** 300 × 4 5 µm µBondapak C18**Mobile phase:** MeOH:10 mM potassium phosphate adjusted to pH 4.40 ± 0.05 with 150 mM phosphoric acid 50:50**Flow rate:** 1.7**Injection volume:** 40**Detector:** UV 212

---

**CHROMATOGRAM****Retention time:** 8.90**Internal standard:** 5-ethyl-5-p-tolylbarbituric acid (4.62)**Limit of detection:** 500 ng/mL

---

**OTHER SUBSTANCES****Extracted:** thiopental (UV 284)

**Simultaneous:** acetaminophen, amobarbital, barbital, butalbital, butabarbital, caffeine, carbamazepine, phenacetin, phenobarbital, phenytoin, secobarbital, theobromine, theophylline, vinbarbital

---

**KEY WORDS**

plasma

---

**REFERENCE**

Houdret,N.; Lhermitte,M.; Lalau,G.; Izydorczak,J.; Roussel,P. Determination of thiopental and pentobarbital in plasma using high-performance liquid chromatography, *J.Chromatogr.*, **1985**, 343, 437-442.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Condition a 100 mg Bond-Elut C8 SPE cartridge with 2 volumes of MeOH, 2 volumes of water, and 1 volume of 100 mM pH 5.59 Sørensen's phosphate buffer. 500  $\mu$ L Plasma + 10  $\mu$ L 1 mg/mL sodium secobarbital in EtOH, add to the SPE cartridge, wash with 2 volumes of 100 mM pH 5.59 Sørensen's phosphate buffer, wash with 1 volume of water, elute with 500  $\mu$ L MeOH. Evaporate the eluate to dryness under vacuum, reconstitute in 50  $\mu$ L MeOH, inject an aliquot.

---

**HPLC VARIABLES**

**Guard column:** 10  $\mu$ m Guard-Pak C18 (Waters)

**Column:** 100  $\times$  8 10  $\mu$ m Radial-Pak C8 (Waters)

**Mobile phase:** MeOH:THF:100 mM pH 7.72 Sørensen's phosphate buffer 28:16:52

**Flow rate:** 2.5

**Injection volume:** 50

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 5.21

**Internal standard:** secobarbital (6.39)

**Limit of quantitation:** 500 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** methohexital, thiopental

**Noninterfering:** ketamine

---

**KEY WORDS**

plasma; dog; pharmacokinetics; SPE

---

**REFERENCE**

Avram,M.J.; Krejcie,T.C. Determination of sodium pentobarbital and either sodium methohexital or sodium thiopental in plasma by high-performance liquid chromatography with ultraviolet detection, *J.Chromatogr.*, **1987**, 414, 484-491.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 100  $\mu$ L Serum + 100  $\mu$ L buffer + 1.5 mL IS in 5% isopropanol in chloroform, vortex for 30 s, centrifuge. Remove the organic layer and evaporate it to dryness under a stream of air at room temperature, reconstitute the residue in 100  $\mu$ L mobile phase, inject a 6-10  $\mu$ L aliquot. (Buffer was 13.6 g  $\text{KH}_2\text{PO}_4$  in 90 mL water, pH adjusted to 6.8 with about 3 mL 10 M NaOH, made up to 100 mL.)

---

**HPLC VARIABLES**

**Guard column:** 20  $\times$  4.6 Supelguard LC-1 (Supelco)

**Column:** 250  $\times$  4.6 5  $\mu$ m Supelcosil LC-1 (Supelco)

**Mobile phase:** MeOH:MeCN:buffer 17.5:17.5:65 (Buffer was 2.72 g  $\text{KH}_2\text{PO}_4$  in 1.9 L water, pH adjusted to 6.3 with about 2 mL 1 M NaOH, made up to 2 L.)

**Flow rate:** 2

**Injection volume:** 6-10

**Detector:** UV 204

---

**CHROMATOGRAM****Retention time:** 5.86**Internal standard:** 5-ethyl-5-p-tolybarbituric acid (tolybarb) (4.80)

---

**OTHER SUBSTANCES****Extracted:** acetaminophen, amobarbital, barbital, caffeine, carbamazepine, chloramphenicol, ethosuximide, mephobarbital, methsuximide, phenobarbital, phenytoin, primidone, secobarbital, theophylline, thiopental**Simultaneous:** acetanilide, N-acetylcysteine, N-acetylprocainamide, ampicillin, aspirin, butabarbital, butalbital, chlorpropamide, cimetidine, codeine, cyheptamide, diazoxide, diflunisal, diphylline, disopyramide, ethchlorvynol, gentisic acid, glutethimide, heptabarbital, hexobarbital, ibuprofen, indomethacin, ketoprofen, mefenamic acid, mephenytoin, methaqualone, methsuximide, methyl salicylate, methypylon, morphine, naproxen, nirvanol, oxphenylbutazone, phenacetin, phensuximide, phenylbutazone, procainamide, salicylamide, salicylic acid, sulfamethoxazole, sulindac, tolmetin, trimethoprim, vancomycin**Noninterfering:** amikacin, gentamicin, meprobamate, netilmicin, quinidine, tetracycline, tobramycin, valproic acid

---

**KEY WORDS**

serum

---

**REFERENCE**Meatherall,R.; Ford,D. Isocratic liquid chromatographic determination of theophylline, acetaminophen, chloramphenicol, caffeine, anticonvulsants, and barbiturates in serum, *Ther.Drug Monit.*, **1988**, *10*, 101-115.

---

**SAMPLE****Matrix:** blood**Sample preparation:** Vigorously shake equal volumes of plasma and MeCN, centrifuge at 10000 g for 3 min, inject a 20  $\mu$ L aliquot of the supernatant.

---

**HPLC VARIABLES****Column:** 110  $\times$  4.6 PartiSphere C8 (Whatman)**Mobile phase:** MeCN:120 mM pH 6.2 phosphate buffer 50:50**Flow rate:** 1**Injection volume:** 20**Detector:** UV 270 following post-column reaction. The column effluent flowed through a 6 m  $\times$  0.25 mm ID crocheted PTFE coil irradiated with a Sylvania G8-T5 lamp at 254 nm to the detector.

---

**CHROMATOGRAM****Retention time:** 3.28**Limit of detection:** 200 ng/mL

---

**OTHER SUBSTANCES****Extracted:** thiopental

---

**KEY WORDS**

post-column reaction; post-column photochemical derivatization; plasma

---

**REFERENCE**Schmid,R.W.; Wolf,C. Simultaneous determination of thiopental and its metabolite, pentobarbital, in blood by high-performance liquid chromatography and post-column photochemical reaction, *J.Pharm.Biomed.Anal.*, **1989**, *7*, 1749-1755.

---

**SAMPLE****Matrix:** blood**Sample preparation:** Prepare an SPE cartridge by plugging the end of a 1 mL disposable pipette tip with glass wool and adding about 100 mg Chromosorb P/NAW. Add 50  $\mu$ L plasma then 50  $\mu$ L 10  $\mu$ g/mL tolylphenobarbital in 200 mM HCl to the SPE cartridge, let stand for 2 min, elute with 1 mL chloroform:isopropanol 6:1. Evaporate the eluate to dryness under a stream of nitrogen at 30°, reconstitute the residue in 100  $\mu$ L mobile phase, inject a 15  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 150 × 4.6 5 µm Supelcosil-LC-8

**Mobile phase:** MeCN:water 20:80

**Flow rate:** 3.3

**Injection volume:** 15

**Detector:** UV 208

---

**CHROMATOGRAM**

**Retention time:** 9.18

**Internal standard:** tolylphenobarbital (7.57)

**Limit of detection:** 50-100 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** theophylline, caffeine, barbital, ethosuximide, primidone, carbamazepinediol, phen-acemide, methypyrrolon, nirvanol, phenobarbital, chloramphenicol, butabarbital, carbamazepine epoxide, mephentyoin, amobarbital, carbamazepine, glutethimide, phenytoin, secobarbital, methaqualone

**Noninterfering:** acetaminophen, amikacin, amitriptyline, clonazepam, cyclosporine, desipramine, diazepam, digoxin, disopyramide, gentamicin, imipramine, lidocaine, methotrexate, N-acetylprocainamide, netilmicin, nortriptyline, procainamide, quinidine, salicylic acid, sulfamethoxazole, tobramycin, trimethoprim, valproic acid, p-hydroxyphenobarbital, vancomycin

---

**KEY WORDS**

plasma; SPE

---

**REFERENCE**

Svinarov,D.A.; Dotchev,D.C. Simultaneous liquid-chromatographic determination of some bronchodilators, anticonvulsants, chloramphenicol, and hypnotic agents, with Chromosorb P columns used for sample preparation, *Clin.Chem.*, **1989**, *35*, 1615-1618.

---

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 200-500 µL Whole blood + 1 mL 100 mM pH 7.5 phosphate buffer, vortex for 1 min, add 7 mL n-hexane:diethyl ether 50:50, add 50 µL 100 µg/mL secobarbital in EtOH: water 75:25, shake for 15 min, centrifuge at 4° at 2000 g for 15 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 37°, reconstitute the residue in 100 µL mobile phase, vortex, inject a 5-20 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 100 × 3 5 µm Nucleosil C18

**Mobile phase:** MeCN:water 32:68

**Flow rate:** 0.3

**Injection volume:** 5-20

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 7

**Internal standard:** secobarbital (9)

**Limit of detection:** 10 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** thiopental

---

**KEY WORDS**

whole blood

---

**REFERENCE**

Celardo,A.; Bonati,M. Determination of thiopental measured in human blood by reversed-phase high-performance liquid chromatography, *J.Chromatogr.*, **1990**, *527*, 220-225.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Mix plasma with an equal volume of MeCN, centrifuge at 10000 g, dilute supernatant with an equal volume of water, inject a 50  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 110  $\times$  4.7 5  $\mu$ m PartiSphere C18 (Whatman)

**Mobile phase:** MeCN:15 mM pH 7.0 phosphate buffer 30:70

**Flow rate:** 0.8

**Injection volume:** 50

**Detector:** UV 270 following post-column reaction. The column effluent flowed through a 6 m  $\times$  0.25 mm ID crocheted coil of PTFE tubing irradiated by an 8 W low-pressure mercury lamp to the detector.

---

**CHROMATOGRAM**

**Retention time:** 7

---

**OTHER SUBSTANCES**

**Extracted:** aprobarbital, butethal, secobarbital

---

**KEY WORDS**

plasma; post-column reaction; post-column photochemical derivatization

---

**REFERENCE**

Wolf,C.; Schmid,R.W. Enhanced UV-detection of barbiturates in HPLC analysis by on-line photochemical reaction, *J.Liq.Chromatogr.*, **1990**, 13, 2207–2216.

---

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 300  $\mu$ L Plasma + 20  $\mu$ L 500  $\mu$ g/mL phenylbutazone in 2 mM NaOH + 1.1 mL ether:n-hexane 20:80 + 20  $\mu$ L 3 M phosphoric acid, vortex at 1200 rpm for 1 min, centrifuge at 2000 g for 5 min, freeze in dry ice for 5 min. Remove the organic layer and evaporate it to dryness under reduced pressure (16 mbar) at 40° for 10 min, reconstitute the residue in 50  $\mu$ L 400  $\mu$ M NaOH, inject a 2-40  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Guard column:** 10  $\times$  3 AGP (ChromTech)

**Column:** 100  $\times$  4 AGP-CSP (ChromTech)

**Mobile phase:** Isopropanol:100 mM pH 6.2 phosphate buffer 4.5:95.5

**Flow rate:** 0.9

**Injection volume:** 2-40

**Detector:** UV 220

---

**CHROMATOGRAM**

**Retention time:** 3.1 (R(+)), 4.1 (S(-))

**Internal standard:** phenylbutazone (15.7)

---

**OTHER SUBSTANCES**

**Extracted:** thiopental (UV 287)

---

**KEY WORDS**

sheep; plasma; chiral

---

**REFERENCE**

Huang,J.L.; Mather,L.E.; Duke,C.C. High-performance liquid chromatographic determination of thiopentone enantiomers in sheep plasma, *J.Chromatogr.B*, **1995**, 673, 245–250.

---

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Mix serum with an equal volume of 1 M pH 5.0 phosphate buffer, add IS, add 3 mL ethyl acetate, rotate for 10 min, centrifuge at 3500 rpm for 5 min. Remove 2 mL of the organic layer and evaporate it to dryness, reconstitute the residue in 200  $\mu$ L mobile phase, inject a 10-20  $\mu$ L aliquot.

---

#### HPLC VARIABLES

**Column:** 150  $\times$  6 Shim-pack CLC-ODS (Shimadzu)

**Mobile phase:** MeOH:10 mM pH 5.0 sodium phosphate buffer 55:45

**Flow rate:** 1

**Injection volume:** 10-20

**Detector:** UV

---

#### CHROMATOGRAM

**Internal standard:** 5-(p-methylphenyl)-5-phenylhydantoin

**Limit of quantitation:** 1  $\mu$ g/mL

---

#### OTHER SUBSTANCES

**Also analyzed:** thiopental

---

#### KEY WORDS

serum; rat; pharmacokinetics

---

#### REFERENCE

Nakashima,E.; Matsushita,R.; Ohshima,T.; Tsuji,A.; Ichimura,F. Quantitative relationship between structure and peritoneal membrane transport based on physiological pharmacokinetic concepts for acidic drugs, *Drug Metab.Dispos.*, **1995**, 23, 1220-1224.

---

#### SAMPLE

**Matrix:** blood, CSF, gastric contents, urine

**Sample preparation:** 200  $\mu$ L Serum, urine, CSF, or gastric fluid + 300  $\mu$ L reagent. Flush column A to waste with 500  $\mu$ L 500 mM ammonium sulfate, inject sample onto column A, flush column A to waste with 500  $\mu$ L 500 mM ammonium sulfate, backflush the contents of column A onto column B with mobile phase, monitor the effluent from column B. (Reagent was 8.05 M guanidine HCl and 1.02 M ammonium sulfate in water.)

---

#### HPLC VARIABLES

**Column:** A 40  $\mu$ m preparative grade C18 (Analytichem); B 75  $\times$  2.1 pellicular C18 (Whatman) + 250  $\times$  4.6 5  $\mu$ m C8 end-capped (Whatman)

**Mobile phase:** Gradient. A was 50 mM pH 4.5  $\text{KH}_2\text{PO}_4$ . B was MeCN:isopropanol 80:20. A:B 90:10 for 1 min, to 30:70 over 20 min.

**Column temperature:** 50

**Flow rate:** 1.5

**Detector:** UV 220

---

#### CHROMATOGRAM

**Retention time:** 12.12

**Internal standard:** heptanophenone (19)

---

#### OTHER SUBSTANCES

**Extracted:** acetaminophen, allobarbitol, azinphos, barbital, brallobarbitone, bromazepam, but-ethal, caffeine, carbamazepine, carbaryl, cephaloridine, chloramphenicol, chlordiazepoxide, chlorothiazide, chlorvinphos, clothiapine, cocaine, coomassie blue, desipramine, diazepam, diphenhydramine, dipipanone, ethylbromphos, flufenamic acid, formothion, griseofulvin, indomethacin, lidocaine, lorazepam, malathion, medazepam, midazolam, oxazepam, paraoxon, penicillin G, prazepam, propoxyphene, prothiophos, quinine, salicylic acid, secobarbital, strychnine, sulfamethoxazole, theophylline, thiopental, thioridazine, trimethoprim

---

#### KEY WORDS

serum; column-switching

---

**REFERENCE**

Kruger,P.B.; Albrecht,C.F.De V.; Jaarsveld,P.P. Use of guanidine hydrochloride and ammonium sulfate in comprehensive in-line sorption enrichment of xenobiotics in biological fluids by high-performance liquid chromatography, *J.Chromatogr.*, **1993**, 612, 191-198.

---

**SAMPLE**

**Matrix:** blood, tissue

**Sample preparation:** Blood or serum. 1 mL Blood or serum + 1 µg cianopramine + 1 mL water, vortex, add 1 mL 200 mM sodium carbonate, vortex, add 6 mL hexane:1-butanol 95:5, gently agitate for 30 min, centrifuge at 2500 g for 5 min. Remove the organic layer and add it to 100 µL 0.2% phosphoric acid, agitate gently for 30 min, centrifuge for 5 min. Remove the organic layer and inject a 30 µL aliquot of the aqueous layer. Liver homogenate. 0.5 mL Liver homogenate + 10 µg cianopramine + 500 µL 2% sodium tetraborate + 8 mL hexane:1-butanol 95:5, gently agitate for 30 min, centrifuge at 2500 g for 5 min. Remove the organic layer and add it to 400 µL 0.2% phosphoric acid, agitate gently for 30 min, centrifuge for 5 min. Remove the organic layer and inject a 30 µL aliquot of the aqueous layer.

---

**HPLC VARIABLES**

**Guard column:** 15 × 3.2 7 µm RP-18 Newguard (Applied Biosystems)

**Column:** 100 × 4.6 5 µm Brownlee Spheri-5 RP-18

**Mobile phase:** MeCN:100 mM NaH<sub>2</sub>PO<sub>4</sub>:diethylamine 40:57.5:2.5

**Flow rate:** 2

**Injection volume:** 30

**Detector:** UV 220

---

**CHROMATOGRAM**

**Retention time:** 1.50

**Internal standard:** cianopramine (8.93)

---

**OTHER SUBSTANCES**

**Simultaneous:** amitriptyline, amoxapine, benztropine, brompheniramine, chlorpheniramine, chlorpromazine, clomipramine, cyproheptadine, desipramine, diphenhydramine, dothiepin, doxepin, fluoxetine, haloperidol, imipramine, loxapine, maprotiline, meperidine, mesoridazine, methadone, mianserin, nomifensine, nordoxepin, norfluoxetine, norpropoxyphene, northiaden, nortriptyline, pheniramine, promethazine, propoxyphene, propranolol, protriptyline, quinidine, quinine, sulforidazine, thioridazine, thiothixene, trazodone, trihexyphenidyl, trimipramine, triprolidine

**Noninterfering:** dextromethorphan, norphethidine, phenoxybenzamine, prochlorperazine, trifluoperazine

**Interfering:** moclobemide, tranlycypromine, metoclopramide

---

**KEY WORDS**

serum; whole blood; liver

---

**REFERENCE**

McIntyre,I.M.; King,C.V.; Skafidis,S.; Drummer,O.H. Dual ultraviolet wavelength high-performance liquid chromatographic method for the forensic or clinical analysis of seventeen antidepressants and some selected metabolites, *J.Chromatogr.*, **1993**, 621, 215-223.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18



**Column:** 250 × 4.6 5 µm Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 200.5

---

#### CHROMATOGRAM

**Retention time:** 16.437

---

#### KEY WORDS

whole blood

---

#### REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

---

#### SAMPLE

**Matrix:** formulations

**Sample preparation:** Dissolve injection in mobile phase to give a pentobarbital sodium concentration of 1 mg/mL, inject a 20 µL aliquot.

---

#### HPLC VARIABLES

**Column:** 250 × 4 10 µm Partisil ODS-3 or 300 × 4 10 µm µBondapak C18

**Mobile phase:** MeOH:buffer:propylene glycol 55:45:4 (Buffer was 4.1 g anhydrous sodium acetate and 15 mL acetic acid in 1 L water.)

**Flow rate:** 2

**Injection volume:** 20

**Detector:** UV 230

---

#### CHROMATOGRAM

**Retention time:** 3

---

#### OTHER SUBSTANCES

**Simultaneous:** degradation products

---

#### KEY WORDS

rugged; injections

---

#### REFERENCE

Reif,V.D.; Kaufmann,K.L.; DeAngelis,N.J.; Frankhouser,M.C. Liquid chromatographic assays for barbiturate injections, *J.Pharm.Sci.*, **1986**, 75, 714-716.

---

#### SAMPLE

**Matrix:** solutions

**Sample preparation:** Mix 50 µL of a 20-200 µg/mL solution in acetone with 50 µL of a 0.4-1.6 mg/mL solution of 2-bromo-2'-acetoneaphthone in acetone, add 5-10 mg cesium carbonate, heat at 30° for 30 min, add 50 µL glacial acetic acid, mix, inject an aliquot.

---

#### HPLC VARIABLES

**Column:** 300 × 4 µBondapak C18

**Mobile phase:** MeOH:water 80:20

**Flow rate:** 2

**Detector:** UV 249

---

**CHROMATOGRAM****Retention time:** 8.75**Limit of detection:** 1 ng

---

**OTHER SUBSTANCES****Simultaneous:** amobarbital, barbital, butobarbital, heptobarbital, hexobarbital, mephobarbital, phenobarbital, secobarbital

---

**KEY WORDS**

derivatization

---

**REFERENCE**Hulshoff,A.; Roseboom,H.; Renema,J. Improved detectability of barbiturates in high-performance liquid chromatography by pre-column labelling and ultraviolet detection, *J.Chromatogr.*, **1979**, 186, 535–541.

---

**SAMPLE****Matrix:** solutions

---

**HPLC VARIABLES****Column:** 300 × 3.9 µBondapak C18**Mobile phase:** MeCN:10 mm KH<sub>2</sub>PO<sub>4</sub> + 5 mM 1-decanesulfonic acid 30:70 adjusted to pH 3.2 with 85% phosphoric acid**Flow rate:** 1**Injection volume:** 10**Detector:** UV 214

---

**CHROMATOGRAM****Retention time:** 11.5**Internal standard:** methyl paraben (7.0)**Limit of detection:** 100 ng/mL

---

**OTHER SUBSTANCES****Simultaneous:** allobarbital, barbital, butalbital, aprobarbital, mephobarbital, phenobarbital, secobarbital, talbutal, vinbarbital

---

**KEY WORDS**

stability-indicating

---

**REFERENCE**Ibrahim,F.B. Simultaneous determination and separation of several barbiturates and analgesic products by ion-pair high-performance liquid chromatography, *J.Liq.Chromatogr.*, **1993**, 16, 2835–2851.

---

**SAMPLE****Matrix:** solutions**Sample preparation:** Dissolve in mobile phase to a concentration of 50 µg/mL.

---

**HPLC VARIABLES****Column:** 250 × 4 β-cyclodextrin polymer-coated silica (Chromatographia 1993, 36, 373)**Mobile phase:** MeOH:water 50:50**Flow rate:** 0.6**Injection volume:** 20**Detector:** UV 240

---

**CHROMATOGRAM****Retention time:** k' 2.09

---

**OTHER SUBSTANCES****Simultaneous:** aprobarbital, amobarbital, butabarbital, butalbital, secobarbital, thiopental, phenobarbital

---

**REFERENCE**

Forgács, E.; Cserhádi, T. Retention behaviour of barbituric acid derivatives on a  $\beta$ -cyclodextrin polymer-coated silicon column, *J. Chromatogr. A*, **1994**, *668*, 395–402.

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 Zorbax RX

**Mobile phase:** Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

**Column temperature:** 30

**Flow rate:** 2

**Detector:** UV 210

---

**OTHER SUBSTANCES**

**Also analyzed:** acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenopropfen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiacol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarboxystiril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methylpyrrolone, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nyldrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrrolamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopalamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

---

**REFERENCE**

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Dissolve in mobile phase at a concentration of 100 µg/mL, inject a 5 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 300 × 2 µBondapak C18

**Mobile phase:** MeCN:water 30:70 adjusted to pH 3.0 with formic acid

**Flow rate:** 0.27

**Injection volume:** 5

**Detector:** MS, VG TRIO 2000 single quadrupole MS with EI or CI or UV 270

---

**CHROMATOGRAM**

**Retention time:** 15

---

**OTHER SUBSTANCES**

**Extracted:** butethal, butabarbital, talbutal, butalbital, amobarbital

---

**KEY WORDS**

mass spectra given

---

**REFERENCE**

Ryan,T.W. Identification of barbiturates using high performance liquid chromatography-particle beam EI/CI mass spectrometry, *J.Liq.Chromatogr.*, **1994**, *17*, 867-881.

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4.5 µm LiChrospher 100 RP-8 (B)

**Mobile phase:** MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

**Flow rate:** 0.6

**Injection volume:** 25

**Detector:** UV 229

---

**CHROMATOGRAM**

**Retention time:** 5.92 (A), 5.58 (B)

---

**OTHER SUBSTANCES**

**Also analyzed:** acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordi-azepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenopropfen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxy-chloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazin-dol, mefenamic acid, meperidine, mephentoin, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine,

methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfinpyrazone, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocainide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

## KEY WORDS

details of plasma extraction

## REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J.Chromatogr.A*, **1995**, 692, 103–119.

## SAMPLE

**Matrix:** urine

**Sample preparation:** 2 mL Urine +1 mL 500 mM pH 5.5 phosphate buffer, add to an Extrelut 3 SPE cartridge, let stand for 10 min, elute with 15 mL dichloromethane:isopropanol 95:5. Evaporate the eluate to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100 µL mobile phase, inject a 20 µL aliquot.

## HPLC VARIABLES

**Guard column:** 4 × 4 5 µm Lichrospher 100 RP8

**Column:** 250 × 4 5 µm Lichrospher 100 RP8

**Mobile phase:** Gradient. MeCN:10 mM pH 4.4 phosphate buffer from 30:70 to 40:60 over 8 min, maintain at 40:60 for 6 min, to 30:70 over 1 min

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 212

## CHROMATOGRAM

**Retention time:** 11.1

**Limit of detection:** 300 ng/mL

## OTHER SUBSTANCES

**Extracted:** barbital, allobarbital, butabarbital, phenobarbital, secobarbital

**Noninterfering:** acetaminophen, aspirin, amitriptyline, buprenorphine, caffeine, carbamazepine, chlorpromazine, desipramine, dextromethorphan, doxepin, ephedrine, fenfluramine, imipramine, lidocaine, loxapine, meperidine, methadone, methaqualone, naloxone, naltrexone, nicotine, orphenadrine, oxycodone, papaverine, pentazocine, phendimetrazine, phenmetrazine, phentermine, phenylpropanolamine, phenytoin, primidone, procaine, promethazine, propoxyphene, propylphenazone, theobromine, theophylline, trazodone, triflupromazine, trimethoprim, trimipramine

## KEY WORDS

SPE

## REFERENCE

Ferrara, S.D.; Tedeschi, L.; Frison, G.; Castagna, F. Solid-phase extraction and HPLC-UV confirmation of drugs of abuse in urine, *J.Anal.Toxicol.*, **1992**, 16, 217–222.


# Pentostatin

**Molecular formula:**  $C_{11}H_{16}N_4O_4$

**Molecular weight:** 268.27

**CAS Registry No.:** 53910-25-1

**Merck Index:** 7277



## SAMPLE

**Matrix:** blood

**Sample preparation:** 200  $\mu$ L Serum + 20  $\mu$ L isoamyl alcohol + 50  $\mu$ L chloroform, vortex for 30 s, centrifuge at 20931 g for 10 min. Remove the aqueous layer and add it to 1 mL cold acetone (0°), vortex for 10 s, centrifuge for 5 min. Remove the supernatant and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 50  $\mu$ L water, inject an aliquot.

## HPLC VARIABLES

**Column:** 250  $\times$  4.6 5  $\mu$ m LiChrosorb RP-8

**Mobile phase:** MeCN:buffer 4:96 (Buffer was 5 mM sodium pentanesulfonate, pH 7.2.)

**Column temperature:** 40

**Flow rate:** 1

**Detector:** UV 250

## CHROMATOGRAM

**Retention time:** 5.82

## OTHER SUBSTANCES

**Extracted:** vidarabine

**Noninterfering:** chlorothiazide, cytosine arabinoside, guanine arabinoside, hydroxyzine, kanamycin, metaproterenol, nystatin, penicillin G, phenobarbital, prednisone, sulfamethoxazole, theophylline, trimethoprim, uracil arabinoside

## REFERENCE

Bowman,D.B.; Kauffman,R.E. Reversed-phase high-performance liquid chromatographic method to determine vidarabine and hypoxanthine arabinoside in biological fluids, *J.Chromatogr.*, **1982**, 229, 487-491.

## SAMPLE

**Matrix:** fermentation solutions

## HPLC VARIABLES

**Column:** 250  $\times$  4.6 5  $\mu$ m Phenomenex C8

**Mobile phase:** MeCN:MeOH:50 mM  $(\text{HN}_4)_2\text{HPO}_4$  2.5:2.5:95 adjusted to pH 7.4 with phosphoric acid

**Flow rate:** 1.5

**Detector:** UV 258

## CHROMATOGRAM

**Retention time:** 6.9


## OTHER SUBSTANCES

**Extracted:** cytosine, coformycin, (8S)-pentostatin, 2'-deoxyguanosine, Ara-A

## REFERENCE

Showalter,H.D.H.; Bunge,R.H.; French,J.C.; Hurley,T.R.; Leeds,R.L.; Leja,B.; McDonnell,P.D.; Edmunds,C.R. Improved production of pentostatin and identification of fermentation cometabolites, *J.Antibiot.(Tokyo)*, **1992**, 45, 1914-1918.

# Pentoxifylline



**Molecular formula:**  $C_{13}H_{18}N_4O_3$

**Molecular weight:** 278.31

**CAS Registry No.:** 6493-05-6

**Merck Index:** 7278

**Lednicer No.:** 2 466

## SAMPLE

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50  $\mu$ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood)  $\mu$ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

## HPLC VARIABLES

**Guard column:** 20 mm long Symmetry C18

**Column:** 250  $\times$  4.6 5  $\mu$ m Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 206.4

## CHROMATOGRAM

**Retention time:** 11.477

## KEY WORDS

whole blood

## REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

## SAMPLE

**Matrix:** microsomal incubations

**Sample preparation:** Add 6 mL ice-cold dichloromethane to 500  $\mu$ L microsomal incubation, add 100  $\mu$ L 10  $\mu$ g/mL CT-2410 R, shake on a reciprocal shaker for 10 min, centrifuge at 3000 g for 10 min, evaporate the organic layer to dryness under a stream of nitrogen at 50°, reconstitute the residue in 100  $\mu$ L mobile phase, inject 30-60  $\mu$ L aliquot.

## HPLC VARIABLES

**Guard column:** Opti-Guard (Optimize Technologies, Inc.)

**Column:** 100  $\times$  4.6 3  $\mu$ m Microsorb-MV C18

**Mobile phase:** MeOH:buffer 35:65 (Buffer was 25 mM ammonium phosphate containing 0.25% acetic acid, pH adjusted to 4.5 with ammonium hydroxide.)

**Flow rate:** 0.7

**Injection volume:** 30-60

**Detector:** UV 273

## CHROMATOGRAM

**Retention time:** 9.3

**Internal standard:** CT-2410 R (20.0)

**Limit of detection:** 10 nM

## OTHER SUBSTANCES

**Extracted:** lisofylline

## KEY WORDS

human; liver

## REFERENCE

Lee, S.H.; Slattery, J.T. Cytochrome P450 isozymes involved in lisofylline metabolism to pentoxifylline in human liver microsomes, *Drug Metab. Dispos.*, **1997**, 25, 1354–1358.

## SAMPLE

**Matrix:** solutions

**Sample preparation:** Inject a 20  $\mu$ L aliquot.

## HPLC VARIABLES

**Column:** 150  $\times$  6.5  $\mu$ m 3-aminopropylsilyl silica gel with the amino group derivatized with 1,8-naphthalic anhydride (Bunseki Kagaku 1993, 42, 817)

**Mobile phase:** MeOH:buffer 50:50 (Prepare buffer by dissolving 6.183 g boric acid and 1.461 g NaCl in 500 mL water, adjust pH to 6.4 with sodium borate solution.)

**Column temperature:** 30

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 270

## CHROMATOGRAM

**Retention time:** 20

**Internal standard:** 1,3-dimethyl-7-(2-hydroxyethyl)xanthine (12)

## OTHER SUBSTANCES

**Simultaneous:** caffeine, hypoxanthine, propentofylline, theobromine, theophylline, uric acid, xanthine

## REFERENCE

Nakashima, K.; Inoue, K.; Mayahara, K.; Kuroda, N.; Hamachi, Y.; Akiyama, S. Use of 3-(1,8-naphthalimido)propyl-modified silyl silica gel as a stationary phase for the high-performance liquid chromatographic separation of purine derivatives, *J. Chromatogr. A*, **1996**, 722, 107–113.

# Pergolide


**Molecular formula:**  $C_{19}H_{26}N_2S$

**Molecular weight:** 314.49

**CAS Registry No.:** 66104-22-1, 66104-23-2 (mesylate)

**Merck Index:** 7304

**Lednicer No.:** 3 249



## SAMPLE

**Matrix:** cell cultures


**Sample preparation:** Evaporate 100  $\mu$ L 1 mg/mL lergotril in MeOH in the bottom of a glass tube under a stream of nitrogen, add 2 mL culture homogenate (Polytron), add 2 mL 100 mM pH 8.5 sodium carbonate/sodium bicarbonate buffer, add 4 mL isoamyl alcohol, shake at 30 oscillations/min for 30 min, centrifuge at 1230 g for 10 min. Remove 2 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 250–500  $\mu$ L mobile phase, filter, inject an aliquot.



**HPLC VARIABLES****Column:** 300 × 3.9 µBondapak C18**Mobile phase:** MeCN:10 mM pH 8.4 ammonium carbonate buffer 65.2:34.8**Flow rate:** 2**Injection volume:** 50**Detector:** UV 290**CHROMATOGRAM****Retention time:** 6.74**Internal standard:** lergotril (2.42)**Limit of quantitation:** 10 µg/mL**OTHER SUBSTANCES****Extracted:** metabolites**REFERENCE**

Kerr,K.M.; Smith,R.V.; Davis,P.J. High-performance liquid chromatographic determination of pergolide and its metabolite, pergolide sulfoxide, in microbial extracts, *J.Chromatogr.*, **1981**, 219, 317–320.

# Perhexiline

**Molecular formula:** C<sub>19</sub>H<sub>35</sub>N**Molecular weight:** 277.49**CAS Registry No.:** 6621-47-2, 6724-53-4 (maleate)**Merck Index:** 7305**SAMPLE****Matrix:** blood

**Sample preparation:** 500 µL Plasma + 200 µL 1.4 µg/mL hexadiline hydrochloride in 100 mM HCl + 50 µL 2 M pH 8.75 Tris-HCl buffer + 4 mL n-hexane, shake horizontally at 100 oscillations/min for 15 min, centrifuge at 2500 rpm for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 60°, reconstitute the residue in 100 µL 100 mM pH 10 sodium bicarbonate buffer, add 100 µL 5 mM dansyl chloride in acetone (freshly prepared), vortex, heat at 37° for 20 min, add 1.5 mL n-hexane, vortex, centrifuge at 10° at 2500 rpm for 3 min, freeze in dry ice/EtOH. Remove the organic layer and evaporate it to dryness at 60°, reconstitute the residue in 100 µL mobile phase, inject a 25 µL aliquot.

**HPLC VARIABLES****Column:** 100 × 3.2 3 µm Velosep (Brownlee)**Mobile phase:** MeOH:water 86:14**Column temperature:** 45**Flow rate:** 0.5**Injection volume:** 25**Detector:** F ex 366 em 470**CHROMATOGRAM****Retention time:** 15.9, 16.8 (isomers)**Internal standard:** hexadiline (19.5)**Limit of quantitation:** 150 ng/mL**OTHER SUBSTANCES****Noninterfering:** metabolites**KEY WORDS**

derivatization; plasma; pharmacokinetics

---

**REFERENCE**

Morris,R.G.; Sallustio,B.C.; Saccoia,N.C.; Mangas,S.; Fergusson,L.K.; Kassapidis,C. Application of an improved HPLC perhexiline assay to human plasma specimens, *J.Liq.Chromatogr.*, **1992**, 15, 3219–3232.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 400  $\mu$ L Plasma + 400  $\mu$ L 1 M diammonium hydrogen phosphate + 50  $\mu$ L 2.5  $\mu$ g/mL di-n-hexylamine in water + 100  $\mu$ L 1 mg/mL diisopropylethylamine + 5 mL isopentane:dichloromethane 60:40, vortex for 1.5 min, centrifuge at 1500 g for 5 min, freeze in dry ice/acetone. Remove the organic layer and evaporate it to dryness under reduced pressure at room temperature, reconstitute the residue in 100  $\mu$ L 10 mg/mL trans-4-nitrocinnamoyl chloride in dry MeCN, vortex for 15 s, let stand at room temperature for 30 min, add 100  $\mu$ L 25 mM sodium carbonate, vortex for 15 s, let stand at room temperature for 5 min, add 100  $\mu$ L MeCN:50 mM ammonium acetate 50:50, vortex for 30 s, centrifuge at 1500 g for 5 min, inject a 130  $\mu$ L aliquot of the supernatant.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu$ m Spherisorb phenyl

**Mobile phase:** MeCN:water:glacial acetic acid 46:53.5:0.5

**Flow rate:** 2

**Injection volume:** 130

**Detector:** UV 340

---

**CHROMATOGRAM**

**Retention time:** 18.4

**Internal standard:** di-n-hexylamine (10.2)

**Limit of detection:** 30 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** metabolites

**Simultaneous:** dealkyldisopyramide, desethylamiodarone, desipramine, flecainide, mexiletine, nortriptyline, protriptyline, sotalol

---

**KEY WORDS**

derivatization; plasma

---

**REFERENCE**

Grgurinovich,N. Method for the analysis of perhexiline and its hydroxy metabolite in plasma using high-performance liquid chromatography with precolumn derivatization, *J.Chromatogr.B*, **1997**, 696, 75–80.

---

---

# Pericyazine


**Molecular formula:** C<sub>21</sub>H<sub>23</sub>N<sub>3</sub>OS

**Molecular weight:** 368.47

**CAS Registry No.:** 2622-26-6

**Merck Index:** 7306

---



---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a 10  $\mu$ g/mL solution in MeOH, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 125  $\times$  4.9 Spherisorb S5W silica

**Mobile phase:** MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

**Flow rate:** 2

**Injection volume:** 20