

REFERENCE

Ferrara,S.D.; Tedeschi,L.; Frison,G.; Castagna,F. Solid-phase extraction and HPLC-UV confirmation of drugs of abuse in urine, *J.Anal.Toxicol.*, **1992**, 16, 217-222.

Phenol



Molecular formula: C₆H₆O

Molecular weight: 94.11

CAS Registry No.: 108-95-2

Merck Index: 7390

SAMPLE

Matrix: air

Sample preparation: Pull 5-150 L air through 10 mL 0.06% NaOH in water at 1-2 L/min. Add 1 mL buffer and 3 mL 0.1% 4-nitrobenzenediazonium tetrafluoroborate in water, make up to 20 mL with water, let stand for 15 min, inject a 2-40 µL aliquot. (Buffer was 0.51% sodium bicarbonate containing 0.21% sodium carbonate, pH 11.5.)

HPLC VARIABLES

Column: 200 × 4.6 Polygosil 60-5 C18 (Macherey-Nagel)

Mobile phase: MeOH:water 85:15

Flow rate: 1.1

Injection volume: 2-40

Detector: UV 365

CHROMATOGRAM

Retention time: 4

Limit of detection: 50 pg

OTHER SUBSTANCES

Also analyzed: o-cresol, m-cresol, p-cresol, 1-naphthol, 2-naphthol, 2,3-xyleneol, 2,4-xyleneol, 2,5-xyleneol, 2,6-xyleneol, 3,4-xyleneol, 3,5-xyleneol

KEY WORDS

derivatization

REFERENCE

Kuwata,K.; Uebori,M.; Yamazaki,Y. Determination of phenol in polluted air as *p*-nitrobenzeneazophenol derivative by reversed phase high performance liquid chromatography, *Anal.Chem.*, **1980**, 52, 857-860.

SAMPLE

Matrix: air

Sample preparation: Pull 5-150 L air through 10 mL 0.12% NaOH in water at 1-2 L/min. Remove a 5 mL aliquot, add 1 mL buffer, add 3 mL 0.1% 4-nitrobenzenediazonium tetrafluoroborate in water, make up to 20 mL with water, let stand for 30 min, add 1 mL 1% NaOH, add 1 mL carbon tetrachloride, shake, centrifuge, inject a 2-40 µL aliquot of the aqueous (p-unsubstituted phenols) layer or a 2-10 µL aliquot of the organic (p-substituted phenols) layer. (Buffer was 0.51% sodium bicarbonate containing 0.21% sodium carbonate, pH 11.5.)

HPLC VARIABLES

Column: 200 × 4.6 5 µm LiChrosorb RP-18

Mobile phase: MeOH:water 85:15

Flow rate: 1.1

Injection volume: 2-40

Detector: UV 365

CHROMATOGRAM

Retention time: 3.89

Limit of detection: 0.05 ng

OTHER SUBSTANCES

Simultaneous: o-chlorophenol, m-chlorophenol, p-chlorophenol, o-cresol, m-cresol, p-cresol, o-ethylphenol, m-ethylphenol, p-ethylphenol, α -naphthol, β -naphthol, 2,3-xyleneol, 2,4-xyleneol, 2,5-xyleneol, 2,6-xyleneol, 3,4-xyleneol, 3,5-xyleneol

KEY WORDS

derivatization

REFERENCE

Kuwata,K.; Uebori,M.; Yamazaki,Y. Reversed-phase liquid chromatographic determination of phenols in auto exhaust and tobacco smoke as p-nitrobenzeneazophenol derivatives, *Anal.Chem.*, **1981**, 53, 1531-1534.

SAMPLE

Matrix: air

Sample preparation: Condition a Sep Pak silica SPE cartridge with 10 mL dichloromethane and dry with helium at 5 L/min. Pull air through a 0.80 μ m cellulose ester membrane filter and the SPE cartridge at 2 L/min for 1 h, desorb the filter with 5 mL 1% acetic acid with sonication for 10 min, elute the SPE cartridge with 5 mL 1% acetic acid, inject aliquots of the eluates.

HPLC VARIABLES

Guard column: 30 \times 4.6 Spheri-5 RP-18

Column: 250 \times 4.6 5 μ m Ultrasphere ODS

Mobile phase: Gradient. A was 1% acetic acid. B was MeCN:acetic acid 99:1. A:B from 0:100 to 90:10 over 10.5 min, to 78:22 to 24.5 min, to 0:100 (step gradient), maintain at 0:100 for 5 min, re-equilibrate for 12 min.

Flow rate: 2

Injection volume: 200

Detector: F ex 304 em 338 for 6.3 min, ex 280 em 325 for 7.7 min, ex 257 em 330 for 5.3 min, ex 342 em 464 for 4.7 min, ex 285 em 310 for 11 min

CHROMATOGRAM

Retention time: 15.5

Limit of detection: 0.29 μ g/cu.m.

OTHER SUBSTANCES

Simultaneous: catechol, cresol, hydroquinone, 3-methylcatechol, scopoletin

KEY WORDS

SPE

REFERENCE

Risner,C.H. The quantification of hydroquinone, catechol, phenol, 3-methylcatechol, scopoletin, m+p-cresol and o-cresol in indoor air samples by high-performance liquid chromatography, *J.Liq.Chromatogr.*, **1993**, 16, 4117-4140.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Serum + 10 μ L 100 μ g/mL p-ethylphenol in water, acidify to pH 1.0 with 1 M HCl, saturate with 100 mg NaCl, add 300 μ L ethyl acetate, shake for 10 min, centrifuge at 1300 g for 5 min, inject a 10 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 150 \times 4.6 Finepak Sil C18S (Jasco)

Mobile phase: MeCN:water 30:70

Flow rate: 1

Injection volume: 10

Detector: F ex 260 em 305 or MS, Hitachi M-1000S quadrupole, APCI, desolvation 399°, vaporization 300°, drift voltage 40 V, negative-ion mode, m/z = 93

CHROMATOGRAM**Retention time:** 5.8**Internal standard:** p-ethylphenol ($m/z = 121$) (16)**Limit of detection:** $<1 \mu\text{M}$

OTHER SUBSTANCES**Extracted:** p-cresol

KEY WORDS

serum

REFERENCE

Niwa, T. Phenol and p-cresol accumulated in uremic serum measured by HPLC with fluorescence detection, *Clin. Chem.*, 1993, 39, 108–111.

SAMPLE**Matrix:** blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 μL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 μL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES**Column:** 300 \times 3.9 4 μm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH_2PO_4 adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30**Flow rate:** 0.8**Injection volume:** 50**Detector:** UV 272

CHROMATOGRAM**Retention time:** 3.40**Limit of detection:** $<120 \text{ ng/mL}$

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vandesine; mexiletine; dipyrindamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclam-

ide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozone; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254–262.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10–30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 13.422

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, *763*, 149–163.

SAMPLE

Matrix: bulk, formulations

Sample preparation: Bulk. Weigh out 100 mg morphine, dissolve in 25 mL MeOH:water:acetic acid 24:72:1, dilute with MeOH:water:acetic acid 24:72:1 to a final concentration of 240 µg/mL morphine, filter (0.45 µm), inject a 20 µL aliquot. Injections. Dilute with MeOH:water:acetic acid 24:72:1 to a final concentration of 240 µg/mL morphine, filter (0.45 µm), inject a 20 µL aliquot.

HPLC VARIABLES

Column: 300 × 3.9 10 µm µBondapak C18

Mobile phase: MeOH:5 mM sodium 1-heptanesulfonate:acetic acid 24:72:1

Flow rate: 1.5
Injection volume: 20
Detector: UV 284

CHROMATOGRAM

Retention time: 4

OTHER SUBSTANCES

Simultaneous: morphine, 2-mercaptobenzothiazole (UV 230), pseudomorphine (UV 230)

KEY WORDS

injections

REFERENCE

Bello, A.C.; Jhangiani, R.K. Liquid chromatographic determination of morphine sulfate and some contaminants in injections and bulk drug material: collaborative study, *J.Assoc. Off. Anal. Chem.*, **1988**, *71*, 1046–1048.

SAMPLE

Matrix: feces, urine

Sample preparation: Urine. 5 mL Urine + 2 mL concentrated HCl + 1 mL 250 µg/mL p-chlorophenol in water, boil for 1 h, cool, add 4 mL diethyl ether, extract by repeated inversion for 1 min, centrifuge at 500 g for 10 min. Remove the organic layer and add it to 3 mL 50 mM NaOH in MeOH, evaporate to dryness under a stream of nitrogen, reconstitute the residue in 500 µL water, filter (0.45 µm), inject a 50 µL aliquot of the filtrate. Feces. 500 mg Feces + 5 mL 100 mM pH 5.5 phosphate buffer + 50 µL 250 µg/mL p-chlorophenol in water, vortex, centrifuge at 500 g for 10 min. Remove the top layer and add it to 2 mL concentrated HCl, boil for 1 h, cool, add 4 mL diethyl ether, extract by repeated inversion for 1 min, centrifuge at 500 g for 10 min. Remove the organic layer and add it to 3 mL 50 mM NaOH in MeOH, evaporate to dryness under a stream of nitrogen, reconstitute the residue in 250 µL water, filter (0.45 µm), inject a 50 µL aliquot of the filtrate.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Econosil RP-18

Mobile phase: MeOH:20 mM pH 4.0 phosphate buffer 48:52

Flow rate: 0.7

Injection volume: 50

Detector: UV 270

CHROMATOGRAM

Retention time: 14

Internal standard: p-chlorophenol (37)

Limit of detection: 200 ng

OTHER SUBSTANCES

Extracted: p-cresol

REFERENCE

Birkett, A.M.; Jones, G.P.; Muir, J.G. Simple high-performance liquid chromatographic analysis of phenol and p-cresol in urine and feces, *J.Chromatogr.B*, **1995**, *674*, 187–191.

SAMPLE

Matrix: formulations

Sample preparation: Dilute a 1 mL aliquot to a theoretical concentration of 30 µg/mL with MeCN:water 2:50, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 150 mm long 5 µm Waters Resolve C18

Mobile phase: MeOH:100 mM ammonium acetate:diethylamine 1280:2720:4, adjusted to a pH of 7.5

Column temperature: 30

Flow rate: 1

Injection volume: 20

Detector: UV 230

CHROMATOGRAM

Retention time: 4

OTHER SUBSTANCES

Simultaneous: nizatidine sulfoxide, nizatidine amide, nizatidine

KEY WORDS

injections

REFERENCE

Raineri,D.L.; Cwik,M.J.; Rodvold,K.A.; Deyo,K.L.; Scaros,L.P.; Fischer,J.H. Stability of nizatidine in commonly used intravenous fluids and containers, *Am.J.Hosp.Pharm.*, **1988**, *45*, 1523–1529.

SAMPLE

Matrix: formulations

Sample preparation: Inject a 25 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Spherisorb C18

Mobile phase: MeOH:THF:isopropanol:water 30:3.5:1.5:65

Column temperature: 40

Flow rate: 1

Injection volume: 25

Detector: UV 245

CHROMATOGRAM

Retention time: 17

OTHER SUBSTANCES

Simultaneous: methylparaben, propylparaben, flumazenil

Noninterfering: aminophylline, cimetidine, dobutamine, dopamine, famotidine, lidocaine, procainamide, ranitidine

KEY WORDS

injections; 5% dextrose

REFERENCE

Olsen,K.M.; Gurley,B.J.; Davis,G.A.; Christensen,R.; Monaghan,M.S. Stability of flumazenil with selected drugs in 5% dextrose injection, *Am.J.Hosp.Pharm.*, **1993**, *50*, 1907–1912.

SAMPLE

Matrix: formulations

Sample preparation: Dilute with water, inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: 12.5 \times 4 5 μ m Zorbax RX-C18

Column: 250 \times 4.6 5 μ m Zorbax RX-C18

Mobile phase: MeCN:buffer 20:80 (Buffer was 50 mM NaH_2PO_4 + 1 mM tetramethylammonium chloride + 0.5 mM 1-octanesulfonic acid adjusted to pH 3.5 with concentrated orthophosphoric acid.)

Column temperature: 25

Flow rate: 1

Injection volume: 20

Detector: UV 203

CHROMATOGRAM

Retention time: 12

OTHER SUBSTANCES

Simultaneous: atropine, tropic acid, obidoxime, HI-6

Also analyzed: pralidoxime chloride

KEY WORDS

nerve agent antidote mixtures

REFERENCE

Paddle, B.M.; Dowling, M.H. Simple high-performance liquid chromatographic method for assessing the deterioration of atropine-oxime mixtures employed as antidotes in the treatment of nerve agent poisoning, *J. Chromatogr.*, **1993**, 648, 373–380.

SAMPLE

Matrix: formulations

Sample preparation: Injections and ophthalmic solutions. Dilute with water to an atropine concentration of 80 µg/mL, inject a 20 µL aliquot. Ointment. Weigh out ointment equivalent to about 4 mg atropine sulfate, add 10 mL THF:water 80:20, sonicate and swirl until the ointment is completely dispersed, make up to 50 mL with water, filter (0.45 µm), inject a 20 µL aliquot

HPLC VARIABLES

Column: 250 × 4.6 5 µm Spherisorb CN

Mobile phase: MeCN:50 mM NaH₂PO₄ 10:90, pH adjusted to 4.0 with 10% phosphoric acid

Flow rate: 1

Injection volume: 20

Detector: UV 220

CHROMATOGRAM

Retention time: 5

OTHER SUBSTANCES

Simultaneous: atropine, tropic acid

Noninterfering: benzyl alcohol, methylparaben, benzalkonium chloride, chlorobutanol

KEY WORDS

injections; ophthalmic solutions; ointments

REFERENCE

Lehr, G.J.; Yuen, S.M.; Lawrence, G.D. Liquid chromatographic determination of atropine in nerve gas antidotes and other dosage forms, *J. AOAC Int.*, **1995**, 78, 339–343.

SAMPLE

Matrix: honey

Sample preparation: Condition a 3 mL Baker-10 C18 SPE cartridge with 3 mL MeOH and 3 mL water. Dissolve 10 g honey in 90 mL water, add 30 g NaCl, add 2 mL 10% phosphoric acid, distil about 30 mL at about 5 mL/min. Add a 15 mL aliquot of the distillate to 6 g NaCl and 5 mL 5% sodium bicarbonate, extract with 5 mL benzene (Caution! Benzene is a carcinogen!). Wash the organic layer with 3 mL 1% sodium bicarbonate, extract with 3 mL 100 mM NaOH, extract with 2 mL 100 mM NaOH. Combine the aqueous extracts and adjust the pH to 3.0 with 0.33 M phosphoric acid, add 4.5 g NaCl, add to the SPE cartridge, let cartridge dry under vacuum for 3 min, elute with 1 mL MeOH, inject a 10 µL aliquot of the eluate.

HPLC VARIABLES

Guard column: 70 × 2 30 µm Co-Pell ODS

Column: 150 × 4.6 5 µm Inertsil ODS

Mobile phase: MeCN:buffer 20:80 (Buffer was 10 mM NaH₂PO₄ containing 2 mM EDTA adjusted to pH 5.0.)

Flow rate: 1

Injection volume: 10

Detector: E, Irica Model E-520, glassy carbon electrode 0.7 V, Ag/AgCl reference electrode or UV 270

CHROMATOGRAM**Retention time:** 11**Limit of detection:** 2 ppb

KEY WORDSSPE

REFERENCE

Takeba,K.; Matsumoto,M.; Shida,Y.; Nakazawa,H. Determination of phenol in honey by liquid chromatography with amperometric detection, *J.Assoc.Off.Anal.Chem.*, **1990**, 73, 602–604.

SAMPLE**Matrix:** perfusate

Sample preparation: 100 μ L Perfusate (Kreb's Henselit buffer containing 1% bovine serum albumin) + 200 μ L 5 μ g/mL p-cresol in MeOH, vortex for 20 s, centrifuge at 10000 g for 10 min, inject a 20 μ L aliquot of the supernatant.

HPLC VARIABLES**Column:** 300 \times 4.6 5 μ m Spherisorb C18**Mobile phase:** MeOH:water:orthophosphoric acid 40:60:0.1, pH 2.7 \pm 0.02**Flow rate:** 1**Injection volume:** 20**Detector:** UV 254

CHROMATOGRAM**Retention time:** 6.5**Internal standard:** p-cresol (10.3)**Limit of quantitation:** 500 ng/mL

OTHER SUBSTANCES**Extracted:** phenyl- β -D-glucuronide

REFERENCE

Thompson,M.J.; Ballinger,L.N.; Cross,S.E.; Roberts,M.S. High-performance liquid chromatographic determination of phenol, 4-nitrophenol, β -naphthol and a number of their glucuronide and sulphate conjugates in organ perfusate, *J.Chromatogr.B*, **1996**, 677, 117–122.

SAMPLE**Matrix:** solutions

Sample preparation: Mix solution with 5 mL 200 g/L sodium acetate trihydrate and 650 μ L reagent, let stand for 1 min, add 5 mL 160 g/L sodium carbonate monohydrate, add 150 mg tetrabutylammonium bromide, extract with 2 mL n-butanol, centrifuge, inject an aliquot of the organic layer. (Prepare reagent by mixing 2.566 g p-aminobenzonitrile and 108 mL concentrated HCl in 1 L water. Cool a 25 mL aliquot in an ice bath, slowly add with stirring 3 mL of a 25 g/L sodium nitrite solution. Use reagent within 1 h.)

HPLC VARIABLES**Column:** 250 \times 2.6 HC ODS/SIL-X C18 (Perkin-Elmer)**Mobile phase:** MeOH:water 64:36**Flow rate:** 1**Injection volume:** 10**Detector:** UV 370

CHROMATOGRAM**Retention time:** 3.20**Limit of detection:** 10 ppb

OTHER SUBSTANCES

Simultaneous: 2-sec-butylphenol, 2-tert-butylphenol, 3-tert-butylphenol, catechol, o-cresol, m-cresol, p-cresol, 2,3-dimethylphenol, 2,5-dimethylphenol, 2,6-dimethylphenol, 3,5-dimethylphenol, 3-hydroxybenzoic acid, 3-nitrophenol, salicylic acid, 2,3,5,6-tetramethylphenol, 3-trifluoromethylphenol

KEY WORDSderivatization

REFERENCE

Baiocchi,C.; Campi,E.; Gennaro,M.; Mentasti,E.; Mirti,P. Reversed phase liquid chromatographic separation of phenolic compounds with a new derivatizing reaction, *Chromatographia*, **1982**, 15, 660–664.

SAMPLE

Matrix: solutions

Sample preparation: Mix 10 mL of a 50 µg/mL solution with 2 mL 100 mg/mL NaOH and 5 mL reagent, mix, let stand for 15 min, add 3.4 mL 15 mg/mL tetrabutylammonium bromide in butanol (saturated with water), extract, inject a 10 µL aliquot of the organic layer. (Prepare reagent by mixing 5 volumes 7.6 mg/mL sulfanilic acid with 1 volume 470 mg/mL sulfuric acid, cool in an ice bath, slowly add 5 volumes 3.4 mg/mL sodium nitrite. Discard the reagent after 10 min.)

HPLC VARIABLES

Column: 300 × 3.9 10 µm µBondapak RP phenyl

Mobile phase: MeOH:water 48:52 containing 3 mM tetrabutylammonium bromide

Flow rate: 2

Injection volume: 10

Detector: UV 370

CHROMATOGRAM

Retention time: 5.12

Limit of quantitation: 100 ppb

OTHER SUBSTANCES

Simultaneous: 2-sec-butylphenol, 3-t-butylphenol, 3-chlorophenol, 2-methylphenol, 3-methylphenol, 4-nitrophenol, 2,3-xyleneol, 2,5-xyleneol, 2,6-xyleneol, 3,5-xyleneol

Noninterfering: 2-t-butylphenol, 4-methylphenol

KEY WORDSderivatization

REFERENCE

Baiocchi,C.; Gennaro,M.C.; Campi,E.; Mentasti,E.; Aruga,R. HPLC identification and separation of phenolic compounds derivatized with diazotized sulfanilic acid. Structural effects on retention times, *Anal.Lett.*, **1982**, 15, 1539–1548.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 10 µm LiChrosorb RP 18

Mobile phase: MeOH:10 mM pH 5.5 potassium phosphate buffer 3.5:96.5

Flow rate: 2-3

Detector: UV 254

CHROMATOGRAM

Retention time: 27

OTHER SUBSTANCES

Simultaneous: catechol, hydroquinone (quinol), phenyl glucuronide, phenyl glucoside, phenyl galactopyranoside, phenyl sulfate, resorcinol

REFERENCE

Beyer,J.; Frank,G. Hydroxylation and conjugation of phenol by the frog *Rana temporaria*, *Xenobiotica*, **1985**, 15, 277–280.

SAMPLE**Matrix:** solutions

Sample preparation: Mix 3.5 mL of a 40 nM-50 μ M solution in benzene with 500 μ L 3 mM reagent in benzene and 100 μ L 3.7 mM pyridine in benzene (Caution! Benzene is a carcinogen!), heat in the dark at 100° for 40 min (primary and secondary alcohols) or at 140° for 50 min (tertiary alcohols), cool, dilute 100-fold with mobile phase, inject a 20 μ L aliquot. (The reagent is 2-(5-chlorocarbonyl-2-oxazolyl)-5,6-methylenedioxybenzofuran (Dojindo Laboratories, Kumamoto, Japan). Synthesis is as follows. Add ethyl oxalyl chloride in ether to a solution of diazomethane in ether at 0° to give ethyl diazopyruvate (Caution! Diazo compounds are explosive and toxic!) (cf. Buehler, C.A.; Pearson, D.E. *Survey of Organic Syntheses*, Wiley, New York, 1970, p. 179). Heat 100 mg ethyl diazopyruvate, a few mg copper(II) acetylacetonate, and 400 μ L chloroacetonitrile in benzene at 60° overnight (Caution! Benzene is a carcinogen!), cool, add to sodium bicarbonate solution, extract with ether, dry the organic layer, evaporate, chromatograph on silica with petroleum ether:ethyl acetate 90:10, distil the product at 90°/12 mm Hg to give ethyl 2-chloromethyl-5-oxazolecarboxylate as an oil in 18% yield (US Patent 4 603 209 (July 29, 1986)). Add 2 mL phosphorus oxychloride dropwise to a solution of 2 g sesamol in 3 mL DMF at 0°, heat on a steam bath with frequent shaking for 1 h, cool in ice, add 50 mL saturated sodium acetate solution, heat on a steam bath for 30 min, cool, filter, recrystallize the solid from EtOH to give 2-hydroxy-4,5-methylenedioxybenzaldehyde as colorless needles (mp 125-126°) (Bull. Chem. Soc. Jpn. 1962, 35, 1321). Stir 1.4 g ethyl 2-chloromethyl-5-oxazolecarboxylate, 1.5 g 2-hydroxy-4,5-methylenedioxybenzaldehyde, 2 g potassium carbonate, and 50 mL anhydrous DMF at 120° overnight, cool, filter. Evaporate the filtrate to dryness under reduced pressure to give 2-(5-ethoxycarbonyl-2-oxazolyl)-5,6-methylenedioxybenzofuran as a colorless crystalline powder (mp 186°) (yield 39%). Reflux 260 mg 2-(5-ethoxycarbonyl-2-oxazolyl)-5,6-methylenedioxybenzofuran, 100 mg KOH, 20 mL EtOH, and 30 mL water for 2 h, concentrate under reduced pressure, dissolve the residue in 100 mL water, wash with ethyl acetate, treat the aqueous layer with activated carbon, acidify the aqueous layer to pH 2 with 2 M HCl. Filter the precipitate and recrystallize it from EtOH to give 2-(2-oxazole-5-carboxylic acid)-5,6-methylenedioxybenzofuran as a colorless crystalline powder (mp 294-295°). Reflux 150 mg 2-(2-oxazole-5-carboxylic acid)-5,6-methylenedioxybenzofuran and 5 mL thionyl chloride for 2 h, pour the reaction mixture into 300 mL petroleum ether. Filter the precipitate and dry it over KOH to give 2-(5-chlorocarbonyl-2-oxazolyl)-5,6-methylenedioxybenzofuran (mp 290°).)

HPLC VARIABLES**Column:** 150 \times 4.6 5 μ m Cosmosil 5C18 (Nacalai Tesque)**Mobile phase:** MeCN:water 70:30**Column temperature:** 40**Flow rate:** 1**Injection volume:** 20**Detector:** F ex 360 em 460

CHROMATOGRAM**Retention time:** 4.3**Limit of detection:** 750 fmole

OTHER SUBSTANCES**Simultaneous:** benzyl alcohol, cyclohexanol, 1-butanol, 1-hexanol, 3-methyl-1-butanol, 2-methyl-2-butanol, 2-methyl-1-propanol, 2-methyl-2-propanol, 1-propanol**Noninterfering:** aldehydes, amino acids, aromatic amines, carboxylic acids, ketones, sulfhydryl compounds**Interfering:** 2-propanol

KEY WORDS

derivatization

REFERENCE

Nagaoka, H.; Nohta, H.; Kaetsu, Y.; Saito, M.; Ohkura, Y. 2-(5-Chlorocarbonyl-2-oxazolyl)-5,6-methylenedioxybenzofuran as fluorescence derivatization reagent for alcohols in high performance liquid chromatography, *Anal. Sci.*, **1989**, 5, 525-530.

SAMPLE**Matrix:** solutions

Sample preparation: Inject a 20 μ L aliquot of a solution in 1% acetic acid.

HPLC VARIABLES

Guard column: 30 \times 4.6 Spheri-5 RP-18

Column: 250 \times 4.6 5 μ m Ultrasphere-ODS C18

Mobile phase: Gradient. A was MeCN:acetic acid 99:1. B was 1% acetic acid in water. A:B from 0:100 to 10:90 over 10 min, to 20:80 over 25 min, wash with A for 6 min, re-equilibrate for 14 min.

Flow rate: 2

Injection volume: 20

Detector: F ex 274 em 298

CHROMATOGRAM

Retention time: 6

OTHER SUBSTANCES

Simultaneous: catechol (F ex 280 em 325), hydroquinone (F ex 304 em 338), resorcinol (F ex 284 em 313)

REFERENCE

Risner, C.H.; Cash, S.L. A high-performance liquid chromatographic determination of major phenolic compounds in tobacco smoke, *J. Chromatogr. Sci.*, **1990**, 28, 239–244.

SAMPLE

Matrix: solutions

Sample preparation: Inject 50 μ L onto column A in series with column B, after 1.1 min switch so that column A comes after column B, continue to elute.

HPLC VARIABLES

Column: A 15 \times 3.2 7 μ m New Guard RP-18 (Applied Biosystems); B 100 \times 4.6 3 μ m Econosphere C18 (Alltech)

Mobile phase: MeCN:water:phosphoric acid 25:75:0.2

Flow rate: 1

Injection volume: 50

Detector: UV 200

CHROMATOGRAM

Retention time: 5

Limit of detection: 0.05 ppm

OTHER SUBSTANCES

Extracted: toluene, cresol, benzoic acid

KEY WORDS

groundwater; water; column-switching

REFERENCE

Chamkasem, N.; Hill, K.D.; Sewell, G.W. High-performance liquid chromatographic column-switching technique for the determination of intermediates of anaerobic degradation of toluene in ground water microcosm, *J. Chromatogr.*, **1991**, 587, 185–191.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Supelcosil LC-DP (A) or 250 \times 4 5 μ m LiChrospher 100 RP-8 (B)

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 5.85 (A), 5.29 (B)

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordi-azepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenopropfen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazin-
dol, mefenamic acid, meperidine, mephentermine, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymet-
azoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, phen-
iramine, phenobarbital, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimo-
zide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procain-
amide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, pro-
pantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfapyrazole, sulindac, temazepam, terbutaline, terfena-
dine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocinide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, trifluopro-
mazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yo-
himbine, zopiclone

KEY WORDS

details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J.Chromatogr.A*, **1995**, 692, 103–119.

SAMPLE**Matrix:** solutions

Sample preparation: Filter (0.2 μm) water, remove a 100 μL aliquot and add it to 500 μL 2.59 μM 2-(9-anthryl)ethyl chloroformate in MeCN and 200 μL 25 mM pH 9.6 borate buffer, heat at 43° for 35 min, inject a 10 μL aliquot. (Prepare 2-(9-anthryl)ethyl chloroformate as follows. Stir a solution of 3 g of 9-bromoanthracene in 100 mL ether at 0° under argon or nitrogen, add 9 mL 1.6 M n-butyllithium over 5 min, stir for 30 min, add an ice-cold solution of 3 g ethylene oxide (Caution! Ethylene oxide is a carcinogen!) in 16 mL ether, stir for 1 h, add 70 mL water, add 50 mL ether, remove the organic layer, extract the aqueous layer with 100 mL dichloro-
methane. Combine the organic layers and wash them with water, dry over anhydrous sodium
sulfate, evaporate to dryness, chromatograph on silica gel with dichloromethane to give 2-(9-
anthryl)ethanol as pale yellow crystals (mp 106–8°) (*J.Org.Chem.* 1986, 51, 2956). Stir a solu-
tion of 2-(9-anthryl)ethanol in ether in the presence of pyridine (as an HCl scavenger) at 0°,
add a solution of phosgene in toluene. 2-(9-anthryl)ethyl chloroformate is obtained as colorless
crystals (mp 86–87° from pentane). Protect stock solutions from light and store them in the
refrigerator (*Anal.Chem.* 1991, 63, 292).)

HPLC VARIABLES**Column:** 125 \times 4.5 μm LiChrospher 100 RP-18

Mobile phase: Gradient. MeCN:water from 70:30 to 100:0 over 10 min, maintain at 100:0 for 10 min

Flow rate: 0.75

Injection volume: 10

Detector: F ex 256 em 418 (cut-off filter)

CHROMATOGRAM

Retention time: 8.57

Limit of detection: 7 nM

OTHER SUBSTANCES

Simultaneous: 4-tert-butylphenol, 3,4-dimethylphenol, 4-methylphenol

KEY WORDS

derivatization; wastewater

REFERENCE

Landzettel, W.J.; Hargis, K.J.; Caboot, J.B.; Adkins, K.L.; Strein, T.G.; Veening, H.; Becker, H.-D. High-performance liquid chromatographic separation and detection of phenols using 2-(9-anthrylethyl) chloroformate as a fluorophoric derivatizing agent, *J. Chromatogr. A*, **1995**, 718, 45–51.

SAMPLE

Matrix: solutions

Sample preparation: Wash column A with MeOH at 2 mL/min for 1 min, wash column A with 5 mM tetrabutylammonium bromide at 2 mL/min for 1 min, pump sample through column A at 2 mL/min for 2.5 min and elute to waste, backflush the contents of column A on to column B with MeOH for 1 min, remove column A from the circuit, elute column B with mobile phase, monitor the effluent from column B.

HPLC VARIABLES

Column: A 10 × 2 15-25 µm PLRP-S (Spark Holland); B 250 × 4 ODS-2

Mobile phase: Gradient. MeOH:1% acetic acid from 25:75 to 60:40 over 25 min, to 100:0 over 5 min, maintain at 100:0 for 2 min, return to initial conditions over 2 min.

Column temperature: 65

Flow rate: 1

Injection volume: 5000

Detector: UV 280

CHROMATOGRAM

Retention time: 7

Limit of detection: 0.1-2 ng/mL

OTHER SUBSTANCES

Simultaneous: 4-chloro-3-methylphenol, 2-chlorophenol, 2,4-dichlorophenol, 2,4-dimethylphenol, 2,6-dimethylphenol, 2,4-dinitrophenol, 2-methyl-4,6-dinitrophenol, 2-nitrophenol, 4-nitrophenol, pentachlorophenol (UV 302), 2,4,6-trichlorophenol, 2,4,6-trimethylphenol

KEY WORDS

waste water; river water; column-switching

REFERENCE

Pocurull, E.; Marcé, R.M.; Borrull, F. Improvement of on-line solid-phase extraction for determining phenolic compounds in water, *Chromatographia*, **1995**, 41, 521–526.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 1-10 µg/mL solution in water, inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Hypersil SCX/C18

Mobile phase: MeCN:25 mM pH 3 Na₂HPO₄ 50:50

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: k' 0.79

OTHER SUBSTANCES

Also analyzed: amitriptyline, barbitol, benzoic acid, butabarbital, clomipramine, clonazepam, desipramine, diazepam, flurazepam, furosemide, imipramine, nitrazepam, phenobarbital, phenolphthalein, pindolol, propranolol, resorcinol, salicylic acid, secobarbital, terbutaline, xylazine

KEY WORDS

effect of mobile phase pH on capacity factor is discussed

REFERENCE

Walshe,M.; Kelly,M.T.; Smyth,M.R.; Ritchie,H. Retention studies on mixed-mode columns in high-performance liquid chromatography, *J.Chromatogr.A*, **1995**, 708, 31–40.

SAMPLE

Matrix: solutions

Sample preparation: Inject a 20 µL aliquot of a 100–500 µg/mL solution in mobile phase.

HPLC VARIABLES

Column: 100 × 4.6 5 µm Hypersil C8 MOS 100A coated with phosphatidylcholine (95% pure soybean lecithin, Epikuron, Lucas Meyer & Co.) (Coat column by recycling a 1 mM solution of phosphatidylcholine in MeOH:water 80:20 for 24 h.)

Mobile phase: MeCN:35 mM pH 7.4 sodium phosphate buffer 40:60

Flow rate: 0.5–2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: k' 1.15

OTHER SUBSTANCES

Also analyzed: amoxicillin, antipyrine, carbamazepine, chlorpheniramine, chlorpromazine, clonidine, codeine, desipramine, diphenhydramine, dipyrindamole, ephedrine, flufenamic acid, haloperidol, hydroxyzine, imipramine, indomethacin, lidocaine, megestrol acetate, metoprolol, nabumetone, nadolol, phenobarbital, promazine, propranolol, pyrilamine, quinidine, ropinirole, testosterone, thioridazine, tolafenamic acid, verapamil

Noninterfering: acetaminophen, aspirin, azathioprine, caffeine, carprofen, chlorambucil, cimetidine, fenoterol, flurbiprofen, ibuprofen, ketoprofen, ranitidine, salicylic acid, sulfamethoxazole, theophylline, thioguanine, tiaprofenic acid, trimethoprim, valproic acid

KEY WORDS

comparison with capillary electrophoresis

REFERENCE

Hanna,M.; de Biasi,V.; Bond,B.; Salter,C.; Hutt,A.J.; Camilleri,P. Estimation of the partitioning characteristics of drugs: A comparison of a large and diverse drug series utilizing chromatographic and electrophoretic methodology, *Anal.Chem.*, **1998**, 70, 2092–2099.

SAMPLE

Matrix: urine

Sample preparation: Add 1 mL 12 M HCl to 5 mL urine. Heat at 100° for 30 min. Centrifuge at 2500 rpm for 5 min and dilute 40 fold with water. Extract a 1 mL aliquot of this solution with 1 mL isoamyl alcohol saturated with 6 M HCl. Mix for 2 min. Remove a 500 µL aliquot of the organic layer and add it to 500 µL 500 mM NaOH, vortex for 2 min. Inject a 10 µL aliquot of the aqueous layer.

HPLC VARIABLES

Column: 300 × 4.6 10 µm MicroPak RP18

Mobile phase: MeCN:10 mM HCl 20:80

Flow rate: 1

Injection volume: 10

Detector: UV 220

CHROMATOGRAM

Retention time: 11.93

Limit of detection: 50 ng/mL

REFERENCE

Menezes,M.L.; Demarchi,A.C.C.O. Off line extraction of phenol from human urine sample with isoamyl alcohol and determination by HPLC, *J.Liq.Chromatogr.Rel.Technol.*, **1998**, 21, 2355–2363.

SAMPLE

Matrix: urine

Sample preparation: Adjust pH to 5, hydrolyze with β-glucuronidase/arylsulfatase at 37° for 12 h, add p-chlorophenol at a final concentration of 50 µg/mL, adjust pH to 2 with HCl. 2 mL Sample + 4 mL dichloromethane, vortex, centrifuge at 2400 g for 15 min. Remove 2.5 mL of the organic layer and add it to 500 µL 200 mM NaOH, vortex. Remove 300 µL of the aqueous layer, adjust pH to 7.0 with HCl, inject a 10 µL aliquot.

HPLC VARIABLES

Column: 33 × 4.6 3 µm Pecosphere 3 × 3 C18

Mobile phase: MeOH:water:orthophosphoric acid 30:70:0.1

Flow rate: 1

Injection volume: 10

Detector: UV 215

CHROMATOGRAM

Retention time: 2.02

Internal standard: p-chlorophenol (8)

Limit of detection: 500 ng/mL

OTHER SUBSTANCES

Extracted: p-nitrophenol, cresols, p-methoxyphenol

REFERENCE

Brega,A.; Prandini,P.; Amaglio,C.; Pafumi,E. Determination of phenol, m-, o- and p-cresol, p-aminophenol and p-nitrophenol in urine by high-performance liquid chromatography, *J.Chromatogr.*, **1990**, 535, 311–316.

SAMPLE

Matrix: urine

Sample preparation: Filter (0.22 µm) urine. Dilute a 100 µL aliquot with 900 µL mobile phase, inject a 200 µL aliquot. Hydrolyze conjugates as follows. 1 mL Urine + 1 mL 200 mM pH 4.8 acetate buffer + 200 µL Helix pomatia juice (containing 100000 Fishman Units of β-glucuronidase and 1000000 Roy Units of sulfatase, IBF, France), heat at 37° overnight, dilute 10-fold with mobile phase, centrifuge at 10000 g for 5 min, inject a 200 µL aliquot.

HPLC VARIABLES

Guard column: 15 × 4.6 10 µm ODS 2

Column: 250 × 4.6 5 µm Sup-Rs Classic ODS 2 (Prolabo)

Mobile phase: MeCN:1% phosphoric acid 10:90

Flow rate: 1

Injection volume: 200

Detector: UV 280

CHROMATOGRAM

Retention time: 20.1

Limit of detection: 1 µg/mL

OTHER SUBSTANCES

Extracted: hippuric acid, mandelic acid, 3-methylhippuric acid, phenylglyoxylic acid

REFERENCE

Astier, A. Simultaneous high-performance liquid chromatographic determination of urinary metabolites of benzene, nitrobenzene, toluene, xylene and styrene, *J. Chromatogr.*, **1992**, 573, 318–322.

SAMPLE

Matrix: urine

Sample preparation: Condition a 500 mg Bond Elut SAX SPE cartridge with 3 mL MeOH and 3 mL water. Dilute 125 μ L urine to 4 mL with water, adjust to pH 4.5 with ascorbic acid, add 12.5 μ L enzyme solution, heat at 37° for 48 h, add to the SPE cartridge, wash with 3 mL 5 mM pH 7 phosphate buffer. Acidify the eluate to pH <3 with concentrated HCl, add 5 mL ether, vortex, repeat the extraction twice. Combine the organic layers and evaporate them to dryness under reduced pressure at 30°, reconstitute the residue in 1 mL 1% aqueous phosphoric acid, inject a 20 μ L aliquot. (The enzyme solutions used to deconjugate glucuronides and sulfate esters were β -glucuronidase/arylsulfatase (Merck, 4114), arylsulfatase (Sigma, S 1629), and β -glucuronidase diluted 1:6 with water (Boehringer, 127051).)

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Nucleosil ODS

Mobile phase: MeOH:5 mM pH 3.4 phosphate buffer 30:70

Flow rate: 1

Injection volume: 20

Detector: UV 270

CHROMATOGRAM

Retention time: 9.3

Limit of detection: 36 μ g/mL

OTHER SUBSTANCES

Extracted: catechol, hydroquinone

KEY WORDS

mouse; SPE

REFERENCE

Schad, H.; Schäfer, F.; Weber, L.; Seidel, H. J. Determination of benzene metabolites in urine of mice by solid-phase extraction and high-performance liquid chromatography, *J. Chromatogr.*, **1992**, 593, 147–151.

SAMPLE

Matrix: urine

Sample preparation: Freeze urine, thaw, centrifuge at 6000 g for 10 min, filter (0.45 μ m PVDF), dilute with 3 volumes of 40 mM pH 6.8, inject an aliquot onto column A with mobile phase A (phenol conjugates are hydrolyzed to phenol on this column) and elute onto column B, after 10 min remove column A from the circuit and elute column B with mobile phase B, monitor the effluent from column B.

HPLC VARIABLES

Column: A 50 \times 4.6 immobilized enzyme reactor (Prepare as follows. Stir 1 g Supelclean LC-NH₂ (Supelco) in 25 mL 1% glutaraldehyde for 1 h, rinse solid with water and 40 mM pH 6.8 phosphate buffer, stir solid in 25 mL 1 mg/mL β -glucuronidase/sulfatase (type H-2 from Helix pomatia, Sigma) in 40 mM pH 6.8 buffer at 4° for 24 h, wash the solid with water, wash with 100 mM NaCl. Make a slurry of 0.5 g of the immobilized enzyme in 40 mM pH 6.8 buffer and pack in a 50 \times 4.6 column using high pressure nitrogen, rinse with 40 mM pH 6.8 phosphate buffer. When not in use store in 100 mM NaCl at 4°); B 150 \times 4.6 5 μ m ODS (J & W)

Mobile phase: A MeOH:40 mM pH 6.8 phosphate buffer 30:70; B MeOH:40 mM pH 6.8 phosphate buffer 30:70

Flow rate: 1

Injection volume: 1000

Detector: F ex 270 em 300

CHROMATOGRAM**Retention time:** 20**Limit of detection:** 10 ppb**Limit of quantitation:** 50 ppb**KEY WORDS**

column-switching; immobilized enzyme reactor

REFERENCE

Jen, J.-F.; Tsai, M.-Y. Determination of phenol in urine by high-performance liquid chromatography with on-line precolumn enzymatic hydrolysis of the conjugates, *J. Chromatogr. B*, **1994**, 658, 87–92.

SAMPLE**Matrix:** urine**Sample preparation:** Filter (0.2 μm), inject an aliquot directly. Hydrolyze conjugates by heating with 6 M HCl at 37° for 18 h, inject an aliquot.**HPLC VARIABLES****Column:** 150 \times 3.9 5 μm Resolve C18 (Waters)**Mobile phase:** MeOH:1.5% trifluoroacetic acid on water 10:90**Flow rate:** 0.5**Detector:** UV or radioactivity**CHROMATOGRAM****Retention time:** 5.9**OTHER SUBSTANCES****Extracted:** metabolites, hydroquinone, phenyl glucuronide, phenyl sulfate**KEY WORDS**

rat

REFERENCE

Hughes, M.F.; Hall, L.L. Disposition of phenol in rat after oral, dermal, intravenous, and intratracheal administration, *Xenobiotica*, **1995**, 25, 873–883.

SAMPLE**Matrix:** water**Sample preparation:** Add disodium EDTA to sample. 100 mL Water + 10 mL pH 8-9 Britton-Robinson buffer ($\mu = 0.09$) + 1 mL 1.5% 4-aminoantipyrine in water + 5 mL 2% potassium ferricyanide in water + 10 mL chloroform, stir for 10 min, inject a 10 μL aliquot of the organic layer.**HPLC VARIABLES****Column:** 300 \times 3.9 μm Bondapak phenyl**Mobile phase:** MeOH:water 60:40**Flow rate:** 1**Injection volume:** 10**Detector:** UV 480**CHROMATOGRAM****Retention time:** k' 1.36**Limit of detection:** 20 ng/mL**OTHER SUBSTANCES****Simultaneous:** 4-chloro-3-methylphenol, 2-chlorophenol, 3-chlorophenol, 4-chlorophenol, 2,3-dichlorophenol, 2,4-dichlorophenol, 2,6-dichlorophenol, 2,3-dimethylphenol, 2,5-dimethylphenol, 2,6-dimethylphenol, 3,5-dimethylphenol, 2-ethylphenol, 3-ethylphenol, 2-methylphenol, 3-methylphenol, 2-nitrophenol, 3-nitrophenol, 2,3,5,6-tetramethylphenol, 2,4,6-trichlorophenol, 2,3,5-trimethylphenol, 2,3,6-trimethylphenol

KEY WORDSderivatization

REFERENCE

Blo,G.; Dondi,F.; Betti,A.; Bigli,C. Determination of phenols in water samples as 4-aminoantipyrine derivatives by high-performance liquid chromatography, *J.Chromatogr.*, **1983**, 257, 69-79.

SAMPLE**Matrix:** water

Sample preparation: 100 mL Water + 2 mL 50 mM KCl in water + 1 mL 0.6% dextrin in water + 1 mL 80 mM silver nitrate in water + 2 mL pH 9 borax buffer + 500 μ L 2% 4-aminoantipyrine in water + 10 mL chloroform, stir at 40° for 40 min, inject an aliquot of the organic layer.

HPLC VARIABLES**Column:** 300 \times 3.9 μ Bondapak phenyl**Mobile phase:** MeOH:water 60:40**Flow rate:** 1**Injection volume:** 10**Detector:** UV 254

CHROMATOGRAM**Retention time:** 7**Limit of quantitation:** 0.5 ppm

KEY WORDSderivatization

REFERENCE

Blo,G.; Dondi,F.; Bigli,C. High-performance liquid chromatographic determination of phenols as 4-aminoantipyrine derivatives; silver chloride as oxidizing agent in the derivatization reaction, *J.Chromatogr.*, **1984**, 295, 231-235.

SAMPLE**Matrix:** water

Sample preparation: Condition a Baker amino SPE cartridge with dichloromethane, dry with nitrogen. Adjust pH to 12 with 1 M NaOH. Remove a 500 μ L aliquot and add it to 100 μ L 30 mg/mL pH 12 tetrabutylammonium bromide in water, add 500 μ L pH 12 water, add 600 μ L 100 μ g/mL dansyl chloride in dichloromethane, vortex vigorously for 2 min, add a 500 μ L aliquot of the dichloromethane layer to the amino SPE cartridge, let stand for 10 min, elute with 3 mL dichloromethane. Evaporate the eluate to dryness, reconstitute the residue in 500 μ L MeOH:water 50:50, inject a 100 μ L aliquot. (Excess dansyl chloride reacts with the amino groups in the SPE cartridge.)

HPLC VARIABLES**Column:** 200 \times 3.1 3 μ m LiChrosorb RP-18

Mobile phase: Gradient. A was MeOH:100 mM pH 7.0 imidazole buffer 97.5:2.5. B was MeOH: 2.5 mM pH 7.0 imidazole buffer 2.5:97.5. A:B 75:25 for 9.5 min, to 85:15 over 0.5 min, maintain at 85:15 for 4.5 min, to 95:5 over 0.5 min, maintain at 95:5 for 4.5 min, to 100:0 over 0.5 min, maintain at 100:0 for 20 min, return to initial conditions over 1 min, re-equilibrate for 20 min.

Flow rate: 0.5**Injection volume:** 100

Detector: Chemiluminescence (470 nm cut-off filter) following post-column reaction. The column effluent flowed through a 130 \times 0.3 mm ID PTFE coil irradiated with a fan-cooled 90 w Philips Model 93110E mercury lamp. The effluent from this coil mixed with 50 mM hydrogen peroxide in MeCN containing 5 mM bis(2-nitrophenyl) oxalate pumped at 0.3 mL/min and this mixture flowed to the detector.

CHROMATOGRAM**Retention time:** 14**Limit of detection:** 0.1 ng/mL

OTHER SUBSTANCES

Simultaneous: 4-chloro-3-methylphenol, 2-chlorophenol, 2,4-dichlorophenol, 2,4-dimethylphenol, 2-nitrophenol, 4-nitrophenol, pentachlorophenol, 2,4,6-trichlorophenol

KEY WORDS

derivatization; post-column reaction; SPE; post-column photochemical derivatization

REFERENCE

Kwakman,P.J.M.; Kamminga,D.A.; Brinkman,U.A.T.; de Jong,G.J. Sensitive liquid chromatographic determination of alkyl-, nitro- and chlorophenols by precolumn derivatization with dansyl chloride, postcolumn photolysis and peroxyoxalate chemiluminescence detection, *J.Chromatogr.*, **1991**, 553, 345-356.

SAMPLE

Matrix: water

Sample preparation: Prepare an SPE cartridge by adding 500 mg 120-400 mesh CarboGraph 4 graphitized carbon black (210 m²/g, Carbochimica Romana, Rome) to a 65 × 13 polypropylene tube using polyethylene frits. Condition with 10 mL 10 mM tetrabutylammonium chloride in dichloromethane:MeOH 80:20, 2 mL MeOH, and 14 mL water acidified to pH 2 with HCl. Filter (Whatman GF/C 1.5 µm glass fiber) river water, pass 4 L through the SPE cartridge at 100 mL/min, wash with 7 mL water at 5-7 mL/min, pull air through the SPE cartridge for 1 min, wash with 800 µL MeOH, dry under vacuum for 1 min, elute in a reverse fashion with 6 mL 10 mM tetrabutylammonium chloride in dichloromethane:MeOH 80:20 at 6 mL/min. Remove a 3 mL aliquot of the eluate and evaporate it to dryness under a stream of nitrogen at 27°, reconstitute the residue in 150 µL 100 mM sodium carbonate in MeCN:water 20:80, add 40 µL acetic anhydride, heat at 50° for 6 min, inject a 50 µL aliquot. (Phenols can also be determined without derivatization. Derivatization provides confirmation of peak identity.)

HPLC VARIABLES

Column: 250 × 4.6 5 µm Alltima LC-18 (Alltech)

Mobile phase: Gradient. A was 0.025% trifluoroacetic acid in water. B was 0.0125% trifluoroacetic acid in MeCN. A:B from 78:22 to 10:90 over 27 min.

Flow rate: 1

Injection volume: 50

Detector: UV 280 for 18 min then UV 220

CHROMATOGRAM

Retention time: 12.4

Limit of detection: <50 ng/mL

OTHER SUBSTANCES

Simultaneous: 4-chloro-3-methylphenol, 2-chlorophenol, 2,4-dichlorophenol, 2,4-dimethylphenol, 4,6-dinitro-2-methylphenol, 2,4-dinitrophenol, 2-nitrophenol, 4-nitrophenol, pentachlorophenol, 2,4,6-trichlorophenol

KEY WORDS

derivatization; SPE; river water

REFERENCE

Di Corcia,A.; Bellioni,A.; Madbouly,M.D.; Marchese,S. Trace determination of phenols in natural waters. Extraction by a new graphitized carbon black cartridge followed by liquid chromatography and re-analysis after phenol derivatization, *J.Chromatogr.A*, **1996**, 733, 383-393.

SAMPLE

Matrix: water

Sample preparation: Condition a 500 mg ENVI Chrom P highly cross-linked styrene-divinylbenzene SPE cartridge (Supelco) with 10 mL MeOH, 10 mL water, and 2 mL 5 mM tetrabutylammonium bromide, dry. Filter (0.45 µm) sample, add 3 mL 10% sodium sulfite solution to each 1 L of water, adjust pH to 9 with 1 M NaOH, add tetrabutylammonium bromide to a final concentration of 5 mM, pass a 500 mL aliquot through the SPE cartridge, elute with 5 mL MeOH, acidify with 1% acetic acid, evaporate to 1 mL under reduced pressure, inject a 20 µL aliquot.


HPLC VARIABLES**Column:** 250 × 4.5 µm Spherisorb ODS-2**Mobile phase:** Gradient. MeOH:1% pH 2.8 acetic acid from 25:75 to 40:60 over 25 min, to 100:0 over 5 min, maintain at 200:0 for 2 min, return to initial conditions over 2 min**Column temperature:** 65**Flow rate:** 1**Injection volume:** 20**Detector:** UV 280**CHROMATOGRAM****Retention time:** 6**Limit of detection:** 0.1 ng/mL**OTHER SUBSTANCES****Simultaneous:** other phenolic compounds**KEY WORDS**

river water; SPE

REFERENCE

Pocurull, E.; Calull, M.; Marcé, R. M.; Borrull, F. Determination of phenolic compounds at low µg l⁻¹ levels by various solid-phase extractions followed by liquid chromatography and diode-array detection, *J. Chromatogr. A*, **1996**, 719, 105–112.

Phenolphthalein

Molecular formula: C₂₀H₁₄O₄**Molecular weight:** 318.33**CAS Registry No.:** 77-09-8**Merck Index:** 7392**SAMPLE****Matrix:** bile, blood, urine

Sample preparation: Condition a 1 mL Bond-Elut C18 SPE cartridge with 2 column volumes of MeOH and 2 column volumes of water, do not allow to go dry. Serum. 1 mL Serum + 25 µL 400 µg/mL bromocresol purple in water + 3 mL acidified acetone, vortex for 5 min, centrifuge at 1200 g for 10 min. Remove 3 mL of the supernatant and evaporate it to dryness under reduced pressure, reconstitute the residue in 1 mL 100 mM pH 7.4 sodium phosphate buffer containing 0.9% NaCl, sonicate for 10 min, add to the SPE cartridge, wash with 2 column volumes of water, allow to dry. Elute with two 800 µL aliquots of MeOH, evaporate the eluate to dryness under reduced pressure, reconstitute in 1 mL MeCN:50 mM NaH₂PO₄ 10:90, centrifuge at 1200 g for 10 min, inject a 100 µL aliquot. Urine. 100 µL Urine + 25 µL 400 µg/mL bromocresol purple in water + 875 µL 100 mM pH 7.4 sodium phosphate buffer containing 0.9% NaCl, mix, add 3 mL acetone:water 88:12, vortex for 5 min, centrifuge at 1200 g for 10 min. Remove 3 mL of the supernatant and evaporate it to dryness under reduced pressure, reconstitute the residue in 1 mL 100 mM pH 7.4 sodium phosphate buffer containing 0.9% NaCl, sonicate for 10 min, add to the SPE cartridge, wash with 2 column volumes of water, allow to dry. Elute with two 800 µL aliquots of MeOH, evaporate the eluate to dryness under reduced pressure, reconstitute in 1 mL MeCN:50 mM NaH₂PO₄ 10:90, centrifuge at 1200 g for 10 min, inject a 100 µL aliquot. Bile. 50 µL Bile + 950 µL MeCN:50 mM NaH₂PO₄ 10:90, centrifuge at 1200 g for 10 min, inject a 50 µL aliquot of the supernatant. (Acidified acetone was 0.5 mL glacial acetic acid in 880 mL acetone, made up to 1 L with water.)

HPLC VARIABLES**Guard column:** 5.3 × 4.10 µm Bondapak C18 guard-pak**Column:** 150 × 4.6 5 µm LC-18DB (Supelco)

Mobile phase: Gradient. A was MeCN:50 mM NaH₂PO₄ 10:90. B was MeCN:50 mM NaH₂PO₄. A:B 65:35 for 5 min, to 40:60 over 0.1 min, maintain at 40:60 for 8 min, re-equilibrate at initial conditions for 5 min.

Flow rate: 1

Injection volume: 100

Detector: UV 230

CHROMATOGRAM

Retention time: 11.4

Internal standard: bromocresol purple (9.5)

Limit of detection: 10 µg/mL (bile), 500 ng/mL (urine), 100 ng/mL (serum)

OTHER SUBSTANCES

Extracted: metabolites, glucuronides

KEY WORDS

serum; dog; SPE; pharmacokinetics

REFERENCE

Wilhelm,J.A.; Bailey,L.C.; Shepard,T.A.; Venturella,V.S. Simultaneous determination of phenolphthalein and phenolphthalein glucuronide from dog serum, urine and bile by high-performance liquid chromatography, *J.Chromatogr.*, **1992**, 578, 231–238.

SAMPLE

Matrix: blood, tissue

Sample preparation: Tissue. Homogenize tissue with 2 volumes of water. 1 mL Homogenate + 300 µL MeOH + 1 mL buffer + 10 mL dichloromethane, rotate, centrifuge. Remove the organic layer and add it to 3 mL 100 mM NaOH, rotate, centrifuge. Remove the aqueous layer and acidify it with 1 mL 1 M HCl, extract with 10 mL dichloromethane. Remove the organic layer and evaporate it to dryness at 40°, reconstitute the residue in 300 µL MeOH, inject a 20 µL aliquot. Blood. 500 µL Whole blood + 500 µL buffer + 200 µL MeOH + 5 mL dichloromethane, rotate, centrifuge. Remove the organic layer and evaporate it to dryness at 40°, reconstitute the residue in 300 µL MeOH, inject a 20 µL aliquot. (Buffer was prepared by mixing 500 mM Na₂HPO₄ with 500 mM KH₂PO₄ to pH 5.5.)

HPLC VARIABLES

Column: 250 × 4.6 10 µm RP-8 (Hewlett-Packard) or 260 × 4.6 10 µm Spherisorb C18

Mobile phase: MeOH:water 60:40

Flow rate: 2

Injection volume: 20

Detector: UV 290

CHROMATOGRAM

Retention time: 2.96

Internal standard: phenolphthalein

OTHER SUBSTANCES

Extracted: thiopental

Simultaneous: carbamazepine, pentobarbital (UV 210)

Noninterfering: amobarbital, butabarbital, glutethimide, meprobamate, methaqualone, meth-
pyrlyon, phenobarbital, phenytoin, secobarbital

KEY WORDS

phenolphthalein is IS; whole blood

REFERENCE

Levine,B.; Blanke,R.; Valentour,J. Liquid chromatographic analysis of thiopental in blood and tissues, *J.Anal.Toxicol.*, **1983**, 7, 207–208.

SAMPLE

Matrix: formulations

Sample preparation: Weigh out formulation containing 500 mg phenolphthalein, make up to 100 mL with MeOH, sonicate for 5 min, filter (0.45 μm). Dilute 5 mL of the filtrate to 100 mL with MeOH, inject a 20 μL aliquot.

HPLC VARIABLES

Column: 100 \times 4.7 μm Nucleosil C18

Mobile phase: MeOH:water:acetic acid 50:50:1

Flow rate: 1.5

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 3.9

Limit of detection: 4 $\mu\text{g/mL}$

OTHER SUBSTANCES

Simultaneous: aloin

KEY WORDS

granules; tablets; dragees; emulsions; comparison with spectrophotometric method

REFERENCE

Torrado,S.; Fraile,S.; Torrado,J.J.; Selles,V.E. Comparison of reversed-phase liquid chromatography with colorimetry for analysis of phenolphthalein preparations, *J.Pharm.Biomed.Anal.*, **1995**, 13, 1167–1172.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Guard column: 4 \times 4.5 μm LiChrospher 100 RP-18

Column: 250 \times 4.5 μm LiChrospher CH-18

Mobile phase: MeOH:10 mM pH 6.0 phosphate buffer 75:25 containing 2.5 mM cethexonium bromide (Rinse with 100 mL MeOH:EtOH:water 50:25:25 at the end of the day.)

Flow rate: 1

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: k' 0.3

OTHER SUBSTANCES

Extracted: glucuronides

REFERENCE

Liu,H.-F.; Leroy,P.; Nicolas,A.; Magdalou,J.; Siest,G. Evaluation of a versatile reversed-phase high-performance liquid chromatographic system using cethexonium bromide as ion-pairing reagent for the analysis of glucuronic acid conjugates, *J.Chromatogr.*, **1989**, 493, 137–147.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a solution in MeOH, inject a 20 μL aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 μm ODS (Altex)

Mobile phase: MeOH:water 5:1

Flow rate: 2

Injection volume: 20

Detector: UV 290

OTHER SUBSTANCES

Simultaneous: thiamylal

REFERENCE

Costantino, A.G.; Caplan, Y.H.; Levine, B.S.; Dixon, A.M.; Smialek, J.E. Thiamylal: review of the literature and report of a suicide, *J. Forensic Sci.*, **1990**, *35*, 89–96.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4.5 µm LiChrospher 100 RP-8 (B)

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 6.05 (A), 5.66 (B)

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordiazepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazinol, mefenamic acid, meperidine, mephénytoin, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfinpyrazone, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocainide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimetoprim, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

KEY WORDS

details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, *692*, 103–119.

SAMPLE

Matrix: solutions


Sample preparation: Prepare a 1–10 µg/mL solution in water, inject an aliquot.

HPLC VARIABLES**Column:** 250 × 4.6 5 µm Hypersil SCX/C18**Mobile phase:** MeCN:25 mM pH 3 Na₂HPO₄ 50:50**Injection volume:** 20**Detector:** UV 254**CHROMATOGRAM****Retention time:** k' 0.96**OTHER SUBSTANCES****Also analyzed:** amitriptyline, barbital, benzoic acid, butabarbital, clomipramine, clonazepam, desipramine, diazepam, flurazepam, furosemide, imipramine, nitrazepam, phenobarbital, phenol, pindolol, propranolol, resorcinol, salicylic acid, secobarbital, terbutaline, xylazine**KEY WORDS**

effect of mobile phase pH on capacity factor is discussed

REFERENCEWalshe,M.; Kelly,M.T.; Smyth,M.R.; Ritchie,H. Retention studies on mixed-mode columns in high-performance liquid chromatography, *J.Chromatogr.A*, **1995**, 708, 31–40.

Phenomorphan

Molecular formula: C₂₄H₂₉NO**Molecular weight:** 347.50**CAS Registry No.:** 468-07-5**Merck Index:** 7399**Lednicer No.:** 1 294**SAMPLE****Matrix:** solutions**Sample preparation:** Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.**HPLC VARIABLES****Column:** 125 × 4.9 Spherisorb S5W silica**Mobile phase:** MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7**Flow rate:** 2**Injection volume:** 20**Detector:** E, LeCarbone, V25 glassy carbon electrode, + 1.2 V**CHROMATOGRAM****Retention time:** 2.3**OTHER SUBSTANCES****Also analyzed:** acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzoctamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipranone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, er-

gosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylegonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, 323, 191–225.

Phenoperidine


Molecular formula: $C_{23}H_{29}NO_3$

Molecular weight: 367.49

CAS Registry No.: 562-26-5, 3627-49-4 (HCl)

Merck Index: 7400

Lednicer No.: 1 302



SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 125 × 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 1.7

OTHER SUBSTANCES


Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine,

chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipiparone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, flupromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserin, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazine, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleppamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, *323*, 191-225.

Phenothiazine



Molecular formula: C₁₂H₉NS

Molecular weight: 199.28

CAS Registry No.: 92-84-2

Merck Index: 7404

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 204

CHROMATOGRAM

Retention time: 14.122, 14.265

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

SAMPLE

Matrix: reaction mixtures

Sample preparation: Evaporate reaction mixture (if necessary, chromatograph residue on deactivated alumina with MeOH), dissolve residue in MeOH with IS, inject an aliquot.

HPLC VARIABLES

Column: 150 × 5 Spherisorb A5Y alumina

Mobile phase: Hexane:ethyl acetate containing 0.6% water 95:5

Flow rate: 2

Detector: UV 254

CHROMATOGRAM

Retention time: k' 2.32

Internal standard: 4-methoxy-2-nitroaniline (k' 8.60)

KEY WORDS

normal phase

REFERENCE

Lunn,G. *Nitrogen-containing Reactive Intermediates in Heterocyclic Synthesis*, Ph.D. Thesis, University of Edinburgh, 1975.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 125 × 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 1.0

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquin-

amide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipiprone, diprenorphine, dipyrindimole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserin, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methidiazene, methotrimprazine, methoxamine, methoxyphenamine, methoxypropazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindine, pimizide, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, 323, 191–225.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in MeOH:water 1:1 at a concentration of 50 µg/mL, inject a 10 µL aliquot.

HPLC VARIABLES

Column: 300 × 3.9 10 µm µBondapak C18

Mobile phase: MeOH:acetic acid:triethylamine:water 70:1.5:0.5:28

Flow rate: 1.5

Injection volume: 10

Detector: UV

CHROMATOGRAM

Retention time: k' 1.56

REFERENCE

Roos, R. W.; Lau-Cam, C. A. General reversed-phase high-performance liquid chromatographic method for the separation of drugs using triethylamine as a competing base, *J. Chromatogr.*, **1986**, 370, 403–418.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitrityline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenpropfen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-stilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isox-suprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, mepidine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, meth-apyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methylodopa, methylodopamine, methylphenidate, methylprednisolone, meth-yltestosterone, methypyrrol, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, ox-ymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendi-metrazine, phenelzine, pheniramine, phensuximide, phentermine, phenylbutazone, phenyl-ephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyril-amine, pyriithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sufadiazine, sulfadimethoxine, sulfaethidole, sulfa-merazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulin-dac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycy-promine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleennamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.

SAMPLE

Matrix: solutions

Sample preparation: Mix 100 μ L of a 0.05-1 μ g/mL solution in MeOH with 20 μ L 5% sodium carbonate and about 35 mg of a suspension of W-2 Raney nickel in EtOH (80 μ L), stir at 70° for 15 min, cool in ice-water, add 50 μ L 400 ng/mL 1-naphthol in MeOH, mix well, filter (0.45 μ m), inject a 5 μ L aliquot of the filtrate. (Raney nickel can be purchased from Aldrich or

prepared from aluminum-nickel alloy (Org. Syn. 1955, Coll. Vol. 3, 181). Phenothiazine is desulfurized to diphenylamine.)

HPLC VARIABLES

Column: 150 × 4.6 5 µm Shim-pack CLC-ODS (Shimadzu)

Mobile phase: MeOH:water 75:25 containing 50 mM sodium perchlorate and 5 mL/L glacial acetic acid

Column temperature: 35

Flow rate: 0.6

Injection volume: 5

Detector: E, Showa Denko Shodex EC-1, glassy carbon working electrode 0.8 V, Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 12.5

Internal standard: 1-naphthol (8.5)

Limit of detection: 10 pg

KEY WORDS

derivatization; desulfurization

REFERENCE

Shimada,K.; Mino,T.; Nakajima,M.; Wakabayashi,H.; Yamato,S. Application of the desulfurization of phenothiazines for a sensitive detection method by high-performance liquid chromatography, *J.Chromatogr.B*, **1994**, 661, 85–91.

SAMPLE

Matrix: tissue

Sample preparation: Homogenize (Polytron) 1 g tissue + 5 mL MeOH, centrifuge at 1600 g for 5 min. Remove 1 mL of the supernatant and evaporate it to dryness, reconstitute the residue in 1 mL cyclohexane, centrifuge, inject a 50 µL aliquot of the supernatant.

HPLC VARIABLES

Guard column: 30 × 2.8 30 µm pellicular beads

Column: 250 × 4.6 Partisil-10

Mobile phase: Cyclohexane:n-propanol 99.7:0.3

Flow rate: 1

Injection volume: 50

Detector: UV 254

CHROMATOGRAM

Retention time: 9

Limit of detection: 50 ppb


KEY WORDS

sheep; kidney; muscle; liver; fat

REFERENCE

Blackman,G.L.; Ho,A.C.; Jozsa,A.; Kelly,J.D. High performance liquid chromatographic determination of phenothiazine residues in sheep tissues, *J.Assoc.Off.Anal.Chem.*, **1980**, 63, 988–991.

Phenoxybenzamine



Molecular formula: C₁₈H₂₂ClNO

Molecular weight: 303.83

CAS Registry No.: 59-96-1, 63-92-3 (HCl)

Merck Index: 7409

Lednicer No.: 1 55

SAMPLE

Matrix: blood, tissue

Sample preparation: Homogenize tissue with 3 parts of 50 mM pH 7.4 phosphate buffer, sonicate adipose tissue homogenates for 30 s. 1 mL Plasma or tissue homogenate + 500 μ L 500 mM sodium carbonate + 6 mL ethyl acetate, extract for 1 h, centrifuge. Remove the organic layer and dry it by mixing with anhydrous sodium sulfate for 20 s, centrifuge, evaporate 5 mL to dryness under a stream of nitrogen at 50°, reconstitute the residue in 500 μ L n-heptane, add 200 μ L MeOH:1 M HCl 90:10, mix for 4 min, centrifuge, discard the heptane layer, inject a 10 μ L aliquot of the MeOH/HCl layer.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Spherisorb ODS2

Mobile phase: MeCN:MeOH:10 mM pH 8.0 phosphate buffer 20:66:14

Flow rate: 1

Injection volume: 10

Detector: UV 220

CHROMATOGRAM

Limit of detection: 10 ng

KEY WORDS

rat; plasma; liver; kidney; lung; brain; heart; muscle; adipose tissue

REFERENCE

Moor, M.J.; Bickel, M.H. Tissue distribution of phenoxybenzamine in the rat. Lack of adipose tissue storage, *Life Sci.*, **1987**, *41*, 2041–2046.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 μ g/mL solution in MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 125 \times 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 1.0

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzoctamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextro-

propoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipiprone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylethylgonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozide, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripelethamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, J.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, 323, 191–225.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in MeOH:water 1:1 at a concentration of 50 µg/mL, inject a 10 µL aliquot.

HPLC VARIABLES

Column: 300 × 3.9 10 µm µBondapak C18

Mobile phase: MeOH:acetic acid:triethylamine:water 60:1.5:0.5:38

Flow rate: 1.5

Injection volume: 10

Detector: UV

CHROMATOGRAM

Retention time: k' 1.78

REFERENCE

Roos, R. W.; Lau-Cam, C. A. General reversed-phase high-performance liquid chromatographic method for the separation of drugs using triethylamine as a competing base, *J. Chromatogr.*, **1986**, 370, 403–418.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 cellulose tris(3,5-dimethylphenylcarbamate)

Mobile phase: Hexane:isopropanol 90:10

Flow rate: 0.5

Detector: UV

CHROMATOGRAM

Retention time: k' 0.70 (of first (-) enantiomer)


KEY WORDS

chiral; α 1.13

REFERENCE

Okamoto,Y.; Aburatani,R.; Hatano,K.; Hatada,K. Optical resolution of racemic drugs by chiral HPLC on cellulose and amylose tris(phenylcarbamate) derivatives, *J.Liq.Chromatogr.*, **1988**, *11*, 2147–2163.

Phensuximide



Molecular formula: $C_{11}H_{11}NO_2$

Molecular weight: 189.21

CAS Registry No.: 86-34-0

Merck Index: 7414

Lednicer No.: 1 226

SAMPLE

Matrix: solutions

Sample preparation: Inject a 6-10 μ L aliquot.

HPLC VARIABLES

Guard column: 20×4.6 Supelguard LC-1 (Supelco)

Column: 250×4.6 5 μ m Supelcosil LC-1 (Supelco)

Mobile phase: MeOH:MeCN:buffer 17.5:17.5:65 (Buffer was 2.72 g KH_2PO_4 in 1.9 L water, pH adjusted to 6.3 with about 2 mL 1 M NaOH, made up to 2 L.)

Flow rate: 2

Injection volume: 6-10

Detector: UV 204

CHROMATOGRAM

Retention time: 4.22

Internal standard: 5-ethyl-5-p-tolybarbituric acid (tolybarb) (4.80)

OTHER SUBSTANCES

Simultaneous: acetaminophen, acetanilide, N-acetylcysteine, N-acetylprocainamide, amobarbital, ampicillin, aspirin, barbital, butalbital, caffeine, carbamazepine, chloramphenicol, chlorpropamide, cimetidine, codeine, cyheptamide, diazoxide, diflunisal, diphyllyne, disopyramide, ethchlorvynol, ethosuximide, gentisic acid, glutethimide, heptabarbital, hexobarbital, ibuprofen, indomethacin, ketoprofen, mefenamic acid, mephentyoin, mephobarbital, methaqualone, methsuximide, methyl salicylate, methyprylon, morphine, naproxen, nirvanol, oxphenylbutazone, pentobarbital, phenacetin, phenobarbital, phenylbutazone, phenytoin, primidone, procainamide, salicylamide, salicylic acid, secobarbital, sulfamethoxazole, sulindac, theophylline, thiopental, tolmetin, vancomycin

Noninterfering: amikacin, gentamicin, meprobamate, netilmicin, quinidine, tetracycline, tobramycin, valproic acid

Interfering: butabarbital, trimethoprim

REFERENCE

Meatherall,R.; Ford,D. Isocratic liquid chromatographic determination of theophylline, acetaminophen, chloramphenicol, caffeine, anticonvulsants, and barbiturates in serum, *Ther.Drug Monit.*, **1988**, *10*, 101–115.

SAMPLE

Matrix: solutions

HPLC VARIABLES**Column:** 100 × 4.6 µm 208HS3410 (Vydac)**Mobile phase:** Gradient. MeCN:water from 15:85 to 60:40 over 10 min.**Flow rate:** 1.5**Detector:** UV 210 (?)

CHROMATOGRAM**Retention time:** 5.2

OTHER SUBSTANCES**Simultaneous:** barbitol, carbamazepine, diazepam, ethotoin, mephentyoin, methsuximide, phenacemide, phenobarbital

REFERENCE*Vydac HPLC Catalog, 1994-5, 1994, p. 26.*

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 250 × 4.6 Zorbax RX**Mobile phase:** Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.**Column temperature:** 30**Flow rate:** 2**Detector:** UV 210


OTHER SUBSTANCES**Also analyzed:** acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitrityline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbitol, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenopropfen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiacol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarbostyryl, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclufenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephentyoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methypylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nyldrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phenacyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrrol-

amine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233–242.

Phentermine



Molecular formula: C₁₀H₁₅N

Molecular weight: 149.24

CAS Registry No.: 122-09-8, 1197-21-3 (HCl)

Merck Index: 7415

Lednicer No.: 1 72

SAMPLE

Matrix: bulk

Sample preparation: Mix a 1 mg/mL solution in 1 M sodium carbonate with 2 mL 5 mg/mL 8-quinolinesulfonyl chloride in acetone, heat at 65° for 20 min, cool, extract twice with 30 mL portions of chloroform. Combine the extracts and dry them over anhydrous magnesium sulfate, evaporate to dryness under a stream of air, reconstitute, inject an aliquot.

HPLC VARIABLES

Guard column: 70 × 2.1 Co:Pell ODS

Column: 300 × 3.9 μBondapak C18

Mobile phase: MeCN:water:acetic acid 40:59:1

Flow rate: 1.5

Detector: UV 254, UV 280

CHROMATOGRAM

Retention time: 30

OTHER SUBSTANCES

Simultaneous: amphetamine, ephedrine, methamphetamine, phenmetrazine, phenylpropanolamine, pseudoephedrine

KEY WORDS

derivatization

REFERENCE

Noggle,F.T.,Jr.; Clark,C.R. Liquid chromatographic determination of primary and secondary amines as 8-quinolinesulfonyl chloride derivatives, *J.Assoc.Off.Anal.Chem.*, **1984**, *67*, 687–691.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a solution in mobile phase, inject 75-100 μL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 μm Supelco