
SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Chirex 3014 (Phenomenex)

Mobile phase: Hexane:1,2-dichloroethane:EtOH/trifluoroacetic acid 55:35:10 (EtOH/trifluoroacetic acid was premixed 20:1.)

Flow rate: 1

Injection volume: 20

Detector: UV 268

CHROMATOGRAM

Retention time: 10, 12 (enantiomers)

KEY WORDS

chiral

REFERENCE

Cleveland,T. Pirkle-concept chiral stationary phases for the HPLC separation of pharmaceutical racemates, *J.Liq.Chromatogr.*, **1995**, 18, 649–671.

SAMPLE

Matrix: solutions

Sample preparation: Inject an aliquot of a 100 µg/mL solution in mobile phase.

HPLC VARIABLES

Column: 150 × 4 5 µm Crownpak CR(+) immobilized crown ether

Mobile phase: MeOH:0.1% pH 1.9 perchloric acid 15:85

Column temperature: 40

Flow rate: 1

Detector: UV 210

CHROMATOGRAM

Retention time: 35.17, 39.67

OTHER SUBSTANCES

Simultaneous: baclofen, levodopa, norephedrine

KEY WORDS

chiral; comparison with capillary electrophoresis

REFERENCE

Nishi,H.; Nakamura,K.; Nakai,H.; Sato,T. Separation of enantiomers and isomers of amino compounds by capillary electrophoresis and high-performance liquid chromatography utilizing crown ethers, *J.Chromatogr.A*, **1997**, 757, 225–235.

Primidone

Molecular formula: C₁₂H₁₄N₂O₂

Molecular weight: 218.26


CAS Registry No.: 125-33-7

Merck Index: 7927

Lednicer No.: 1 276

SAMPLE

Matrix: blood



Sample preparation: Mix 500 μL plasma with 500 μL MeCN and 2 μg IS for 30 s, centrifuge at 2700 g for 5 min, inject an aliquot of the supernatant.

HPLC VARIABLES

Column: 150 \times 4.6 5 μm Ultrasphere C18

Mobile phase: MeCN:MeOH:10 mM pH 7.4 phosphate buffer 15:35:50

Column temperature: 25

Flow rate: 1

Detector: UV 219

CHROMATOGRAM

Internal standard: 2-hydroxy-2-ethyl-2-phenylacetamide

Limit of detection: 50 ng/mL

OTHER SUBSTANCES

Extracted: carbamazepine, clonazepam, ethosuximide, D,L-2-hydroxy-2-ethyl-2-phenylpropionamide (HEPP), phenobarbital, phenytoin

KEY WORDS

rat; plasma

REFERENCE

Martínez de Muñoz,D.; Arenas,R.; Chávez González,O. Liquid chromatographic assay in plasma of one of the members of a new series of anticonvulsants: D,L-3-hydroxy-3-ethyl-3-phenylpropionamide, *J.Chromatogr.B*, **1996**, 678, 377–383.

SAMPLE

Matrix: blood

Sample preparation: Add 200 μL 2 $\mu\text{g/mL}$ thymol in MeCN to 200 μL serum, vortex for 10 s, centrifuge at 7000 g for 5 min, inject 20 μL aliquot.

HPLC VARIABLES

Column: 150 \times 3.9 Resolve C18-5 (Waters)

Mobile phase: MeCN:isopropanol:50 mM pH 3.0 phosphate buffer 25:15:60

Column temperature: 30

Flow rate: 0.7

Injection volume: 20

Detector: UV 220

CHROMATOGRAM

Retention time: 3.0

Internal standard: thymol (18.5)

OTHER SUBSTANCES

Extracted: ethosuximide, phenobarbital, phenytoin, carbamazepine, valproic acid

KEY WORDS

human; plasma

REFERENCE

Kondo,K.; Nakamura,M.; Nishioka,R.; Kawai,S. Direct method of determination of valproic acid in serum by high performance liquid chromatography, *Anal.Sci.*, **1985**, 1, 385–387.

SAMPLE

Matrix: blood

Sample preparation: 500 μL Plasma + 100 μL heptabarbital in MeOH + 500 μL 400 mM pH 7.0 sodium phosphate buffer + 10 mL ethyl acetate, extract. Evaporate the extract to dryness at 50°, reconstitute the residue in 20 μL MeOH, inject a 3 μL aliquot.

HPLC VARIABLES

Guard column: 50 \times 2.1 Whatman Co:Pell ODS

Column: 125 × 4.5 5 µm SAS Hypersil

Mobile phase: MeCN:buffer 20:80 (Buffer was 5 mM tetrabutylammonium hydroxide adjusted to pH 7.5 with phosphoric acid.)

Flow rate: 1.6

Injection volume: 3

Detector: UV 200

CHROMATOGRAM

Retention time: 4.2

Internal standard: heptabarbital (9.8)

Limit of quantitation: 2.5 µM

OTHER SUBSTANCES

Extracted: ethosuximide, phenobarbital, pheneturide, carbamazepine, phenytoin

Simultaneous: phenylethylmalonamide, sulthiame, sulfamethoxazole, ethotoin, butabarbital, pentobarbital, methsuximide, cyclobarbital, ethylphenacemide, amobarbital, glutethimide, secobarbital, barbital

KEY WORDS

plasma; horse

REFERENCE

Christofides, J.A.; Fry, D.E. Measurement of anticonvulsants in serum by reversed-phase ion-pair liquid chromatography, *Clin. Chem.*, **1980**, *26*, 499–501.

SAMPLE

Matrix: blood

Sample preparation: 200 µL Serum or plasma + 200 µL 20 µg/mL IS in MeOH:water 10:90 + 75 µL glacial acetic acid, vortex for 30 s, add 5 mL chloroform, shake for 5 min, centrifuge at 2000 rpm for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 200 µL mobile phase, inject a 40 µL aliquot.

HPLC VARIABLES

Guard column: 30 × 2.1 Permaphase ETH (DuPont)

Column: 250 × 4.6 CLC 1 C8 (DuPont)

Mobile phase: MeCN:buffer 35:65 (Buffer was 20 mM KH₂PO₄ and 1 mM K₂HPO₄ adjusted to pH 5.6.)

Column temperature: 25

Flow rate: 2

Injection volume: 40

Detector: UV 220

CHROMATOGRAM

Retention time: 2.0

Internal standard: alphenal (5-allyl-5-phenylbarbituric acid) (4.4)

Limit of quantitation: 500 ng/mL

OTHER SUBSTANCES

Extracted: ethosuximide, phenytoin, carbamazepine, phenobarbital

Simultaneous: amobarbital, barbital, chlordiazepoxide, codeine, cortisol, ethotoin, glutethimide, hexobarbital, mephentyoin, mephobarbital, metharbital, methsuximide, nitrazepam, pentobarbital, phenacetin, phensuximide, secobarbital

Noninterfering: acetaminophen, acetazolamide, amphetamine, bilirubin, caffeine, diazepam, dimenhydrinate, meperidine, meprobamate, methamphetamine, methaqualone, methylphenidate, nicotine, propoxyphene, theophylline, valproate

KEY WORDS

plasma; serum

REFERENCE

Ryzdewski, R.S.; Gadsden, R.H.; Phelps, C.A. Simultaneous rapid HPLC determination of anticonvulsant drugs in plasma and correlation with EMIT, *Ann. Clin. Lab. Sci.*, **1980**, *10*, 89–94.

SAMPLE**Matrix:** blood**Sample preparation:** 200 μ L Serum + 200 μ L 50 μ g/mL hexobarbital in MeCN + 25 μ L glacial acetic acid, vortex for 10 s, centrifuge for 1 min, inject a 30-100 μ L aliquot of the supernatant.

HPLC VARIABLES**Column:** μ Bondapak C18**Mobile phase:** Gradient. MeCN:7.5 g/L NaH_2PO_4 adjusted to pH 3.2 with phosphoric acid 5:95 to 22:78 over 24 min, to 45:55 over 10 min, maintain at 45:55 for 5 min. Re-equilibrate with 5:95 for 5 min.**Column temperature:** 50**Flow rate:** 3**Injection volume:** 30-100**Detector:** UV 210

CHROMATOGRAM**Retention time:** 10.0**Internal standard:** hexobarbital (20.6)**Limit of detection:** 200-2000 ng/mL

OTHER SUBSTANCES**Extracted:** acetaminophen, amobarbital, butabarbital, butalbital, chlordiazepoxide, diazepam, ethchlorvynol, flurazepam, glutethimide, methaqualone, methypylon, nitrazepam, pentobarbital, phenobarbital, phenytoin, salicylic acid, secobarbital, theophylline**Simultaneous:** amitriptyline, caffeine, clomipramine, codeine, desipramine, ethotoin, imipramine, lidocaine, mesantoin, methsuximide, nirvanol, nortriptyline, oxazepam, procainamide, phenylpropanolamine, propranolol, quinidine

KEY WORDS

serum

REFERENCEKabra,P.M.; Stafford,B.E.; Marton,L.J. Rapid method for screening toxic drugs in serum with liquid chromatography, *J.Anal.Toxicol.*, **1981**, 5, 177-182.

SAMPLE**Matrix:** blood**Sample preparation:** 400 μ L Serum or plasma + 400 μ L 10 μ g/mL IS in acetone, vortex for 10 s, centrifuge at 4500-5000 g for 1 min, remove the supernatant to another tube, centrifuge for 30 s, inject a 5-7.5 μ L aliquot.

HPLC VARIABLES**Column:** 300 \times 3.9 μ Bondapak C18**Mobile phase:** MeCN:MeOH:buffer 17:28:55, final pH 6.8-7.0 (Buffer was 400 μ L 1 M KH_2PO_4 in 1 L water, pH adjusted to 6.0 with 900 mM phosphoric acid.)**Column temperature:** 30**Flow rate:** 0.7**Injection volume:** 5-7.5**Detector:** UV 195

CHROMATOGRAM**Retention time:** 7.4**Internal standard:** tolybarb (5-ethyl-5-(p-methylphenyl)barbituric acid) (13.8)

OTHER SUBSTANCES**Extracted:** carbamazepine, N-desmethylnmethsuximide, ethosuximide, phenobarbital, phenytoin**Simultaneous:** acetaminophen, butalbital, caffeine, hexobarbital, methsuximide, phenacetin, phenylethylmalonamide, salicylic acid

KEY WORDS

plasma; serum

REFERENCE

Szabo, G.K.; Browne, T.R. Improved isocratic liquid-chromatographic simultaneous measurement of phenytoin, phenobarbital, primidone, carbamazepine, ethosuximide, and N-desmethylnmethsuximide in serum, *Clin. Chem.*, **1982**, 28, 100-104.

SAMPLE

Matrix: blood

Sample preparation: 50 μ L Serum + 50 μ L 10 μ g/mL IS in MeCN, vortex for 10 s, centrifuge at 3000 g for 1 min, remove the supernatant and place it in another tube, centrifuge for 1 min, inject a 20 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 100 \times 8.5 μ m Nova Pak C18 Radial pak

Mobile phase: MeCN:MeOH:acetone:buffer 8:21:10:61 adjusted to pH 7.95 \pm 0.02 with NaOH (Buffer was 1.36 g/L KH_2PO_4)

Flow rate: 2.8

Injection volume: 20

Detector: UV 200

CHROMATOGRAM

Retention time: 2.03

Internal standard: tolybarb (5-ethyl-5-(p-methylphenyl)barbituric acid) (4.89)

Limit of detection: 500 ng/mL

OTHER SUBSTANCES

Extracted: ethosuximide, phenobarbital, carbamazepine, phenytoin, metabolites

Simultaneous: acetaminophen, N-acetylprocainamide, aspirin, ampicillin, caffeine, cephalirin, chloramphenicol, digoxin, disopyramide, hexobarbital, indomethacin, lidocaine, mephobarbital, methsuximide, nafcillin, pentobarbital, phenylethylmalonamide, procainamide, quinidine, salicylic acid, secobarbital, sulfamerazine, sulfamethazine, terbutaline, tetracycline, theobromine, theophylline

Noninterfering: acetazolamide, amikacin, cephalosporin C, gentamicin, propranolol, sulfadiazine, sulfamethoxazole, sulfisoxazole, tobramycin, valproic acid, verapamil

KEY WORDS

serum

REFERENCE

Ou, C.-N.; Rognerud, C.L. Simultaneous measurement of ethosuximide, primidone, phenobarbital, phenytoin, carbamazepine, and their bioactive metabolites by liquid chromatography, *Clin. Chem.*, **1984**, 30, 1667-1670.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Plasma + 200 μ L 1 M HCl saturated with ammonium sulfate, vortex for 20 s, add 60 μ L 10 μ g/mL 4-methylprimidone in MeCN, vortex for 20 s, centrifuge at 2700 g for 5 min, inject a 5-10 μ L aliquot of the MeCN layer.

HPLC VARIABLES

Column: 250 \times 4.5 μ m LiChrosorb RP-18

Mobile phase: MeOH:THF:50 mM pH 5.9 phosphate buffer 44:1:55

Column temperature: 50

Flow rate: 1.1

Injection volume: 5-10

Detector: UV 210

CHROMATOGRAM

Retention time: 3

Internal standard: 4-methylprimidone (5)

Limit of detection: 50 ng/mL

OTHER SUBSTANCES

Extracted: carbamazepine, phenobarbital, phenytoin, valproic acid

Simultaneous: acetaminophen, salicylic acid, ethylphenylmalonamide, theophylline, caffeine, ethosuximide, chloramphenicol, methylphenobarbital, glutethimide, pentobarbital, lidocaine, diazepam

KEY WORDS

plasma

REFERENCE

Kushida,K.; Ishizaki,T. Concurrent determination of valproic acid with other antiepileptic drugs by high-performance liquid chromatography, *J.Chromatogr.*, **1985**, 338, 131–139.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Serum + 1 mL 5 μ g/mL 5-ethyl-5-tolylhydantoin in MeCN, agitate for 3 min. Remove the supernatant and evaporate it to dryness, dissolve the residue in 1 mL mobile phase, inject a 50 μ L aliquot.

HPLC VARIABLES

Guard column: pellicular reversed phase (Chrompack 28653)

Column: 100 \times 3 octyl CP-tm-Spher C8 glass column (Chrompack)

Mobile phase: MeCN:50 mM NaH₂PO₄ 25:75 adjusted to pH 2.2 with phosphoric acid

Flow rate: 0.9

Injection volume: 50

Detector: UV 210

CHROMATOGRAM

Retention time: 1.7

Internal standard: 5-ethyl-5-tolylhydantoin (7.8)

Limit of detection: 1000 ng/mL

OTHER SUBSTANCES

Simultaneous: phenobarbital, ethylphenylmalonamide

KEY WORDS

serum

REFERENCE

Van Damme,M.; Molle,L.; Abi Khalil,F. Useful sample handlings for reversed phase high performance liquid chromatography in emergency toxicology, *J.Toxicol.Clin.Toxicol.*, **1985**, 23, 589–614.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Serum + 100 μ L buffer + 1.5 mL IS in 5% isopropanol in chloroform, vortex for 30 s, centrifuge. Remove the organic layer and evaporate it to dryness under a stream of air at room temperature, reconstitute the residue in 100 μ L mobile phase, inject a 6-10 μ L aliquot. (Buffer was 13.6 g KH₂PO₄ in 90 mL water, pH adjusted to 6.8 with about 3 mL 10 M NaOH, made up to 100 mL.)

HPLC VARIABLES

Guard column: 20 \times 4.6 Supelguard LC-1 (Supelco)

Column: 250 \times 4.6 5 μ m Supelcosil LC-1 (Supelco)

Mobile phase: MeOH:MeCN:buffer 17.5:17.5:65 (Buffer was 2.72 g KH₂PO₄ in 1.9 L water, pH adjusted to 6.3 with about 2 mL 1 M NaOH, made up to 2 L.)

Flow rate: 2

Injection volume: 6-10

Detector: UV 204

CHROMATOGRAM

Retention time: 2.71

Internal standard: 5-ethyl-5-p-tolybarbituric acid (tolybarb) (4.80)

OTHER SUBSTANCES

Extracted: acetaminophen, amobarbital, barbital, caffeine, carbamazepine, chloramphenicol, ethosuximide, mephobarbital, methsuximide, pentobarbital, phenobarbital, phenytoin, secobarbital, theophylline, thiopental

Also analyzed: acetanilide, N-acetylcysteine, N-acetylprocainamide, ampicillin, aspirin, butabarbital, butalbital, chlorpropamide, cimetidine, codeine, cyheptamide, diazoxide, diflunisal, diphylline, disopyramide, ethchlorvynol, gentisic acid, glutethimide, heptabarbital, hexobarbital, ibuprofen, indomethacin, ketoprofen, mefenamic acid, mephentoin, methaqualone, methsuximide, methyl salicylate, methypylon, morphine, naproxen, nirvanol, oxphenylbutazone, phenacetin, phensuximide, procainamide, salicylamide, salicylic acid, sulfamethoxazole, tolmetin, trimethoprim, vancomycin

Noninterfering: amikacin, gentamicin, meprobamate, netilmicin, quinidine, tetracycline, tobramycin, valproic acid

Interfering: sulindac, phenylbutazone

KEY WORDS

serum

REFERENCE

Meatherall,R.; Ford,D. Isocratic liquid chromatographic determination of theophylline, acetaminophen, chloramphenicol, caffeine, anticonvulsants, and barbiturates in serum, *Ther.Drug Monit.*, **1988**, *10*, 101-115.

SAMPLE

Matrix: blood

Sample preparation: Prepare an SPE cartridge by plugging the end of a 1 mL disposable pipette tip with glass wool and adding about 100 mg Chromosorb P/NAW. Add 50 μ L plasma then 50 μ L 10 μ g/mL tolylphenobarbital in 200 mM HCl to the SPE cartridge, let stand for 2 min, elute with 1 mL chloroform:isopropanol 6:1. Evaporate the eluate to dryness under a stream of nitrogen at 30°, reconstitute the residue in 100 μ L mobile phase, inject a 15 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Supelcosil-LC-8

Mobile phase: MeCN:water 20:80

Flow rate: 3.3

Injection volume: 15

Detector: UV 208

CHROMATOGRAM

Retention time: 1.61

Internal standard: tolylphenobarbital (7.57)

Limit of detection: 50-100 ng/mL

OTHER SUBSTANCES

Extracted: theophylline, caffeine, barbital, ethosuximide, carbamazepinediol, phenacetamide, methypylon, nirvanol, phenobarbital, chloramphenicol, butabarbital, carbamazepine epoxide, mephentoin, pentobarbital, amobarbital, carbamazepine, glutethimide, phenytoin, secobarbital, methaqualone

Noninterfering: acetaminophen, amikacin, amitriptyline, clonazepam, cyclosporine, desipramine, diazepam, digoxin, disopyramide, gentamicin, imipramine, lidocaine, methotrexate, N-acetylprocainamide, netilmicin, nortriptyline, procainamide, quinidine, salicylic acid, sulfamethoxazole, tobramycin, trimethoprim, valproic acid, p-hydroxyphenobarbital, vancomycin

KEY WORDS

plasma; SPE

REFERENCE

Svinarov,D.A.; Dotchev,D.C. Simultaneous liquid-chromatographic determination of some bronchodilators, anticonvulsants, chloramphenicol, and hypnotic agents, with Chromosorb P columns used for sample preparation, *Clin.Chem.*, **1989**, *35*, 1615-1618.

SAMPLE

Matrix: blood

Sample preparation: 250 μ L Plasma + 2 μ g 10-methoxycarbamazepine + 25 μ L 1 M NaOH + 1.2 mL dichloromethane, mix for 15 min, centrifuge. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 20 μ L MeCN, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.9 10 μ m LiChrosorb RP8

Mobile phase: MeCN:water 32:68

Flow rate: 1.8

Injection volume: 10

Detector: UV 215

CHROMATOGRAM

Retention time: 2.7

Internal standard: 10-methoxycarbamazepine (9.3)

OTHER SUBSTANCES

Extracted: carbamazepine, oxcarbazepine, phenobarbital

Noninterfering: clobazam, clonazepam, diazepam, ethosuximide, phenytoin, valproic acid

KEY WORDS

plasma

REFERENCE

Elyas, A.A.; Goldberg, V.D.; Patsalos, P.N. Simple and rapid micro-analytical high-performance liquid chromatographic technique for the assay of oxcarbazepine and its primary active metabolite 10-hydroxycarbamazepine, *J. Chromatogr.*, **1990**, 528, 473–479.

SAMPLE

Matrix: blood

Sample preparation: Inject 20 μ L serum onto column A with mobile phase A and elute to waste, after 1.5 min backflush the contents of column A onto column B with mobile phase B, after 1 min remove column A from the circuit, elute column B with mobile phase B, monitor the effluent from column B. Re-equilibrate column A with mobile phase A.

HPLC VARIABLES

Column: A 30 \times 4.6 ISRP silica (for preparation see Anal. Chem. 1989, 61, 2445); B 150 \times 4.6 5 μ m Nucleosil C18

Mobile phase: A 14 mM NaH₂PO₄ containing 6 mM Na₂HPO₄; B MeCN:MeOH:14 mM NaH₂PO₄ containing 6 mM Na₂HPO₄ 15:20:65

Flow rate: 0.8

Injection volume: 20

Detector: UV 230

CHROMATOGRAM

Retention time: 8

Limit of quantitation: 1 μ g/mL

OTHER SUBSTANCES

Extracted: carbamazepine, phenobarbital, phenytoin

KEY WORDS

serum; column-switching

REFERENCE

Haginaka, J.; Wakai, J.; Yasuda, H.; Kimura, Y. Determination of anticonvulsant drugs and methyl xanthine derivatives in serum by liquid chromatography with direct injection: column-switching method using a new internal-surface reversed-phase silica support as a precolumn, *J. Chromatogr.*, **1990**, 529, 455–461.

SAMPLE

Matrix: blood

Sample preparation: Add two volumes of MeCN to the mouse serum, mix, centrifuge at 1500 g for 5 min, inject a 5 μ L aliquot of the supernatant.

HPLC VARIABLES

Guard column: Sentry (Waters)

Column: 150 \times 4.6 Nova-Pak C18

Mobile phase: MeCN:MeOH:10 mM pH 7.0 phosphate buffer 50:30:110

Column temperature: 40

Flow rate: 0.5

Injection volume: 5

Detector: UV 214

CHROMATOGRAM

Retention time: 4.5

OTHER SUBSTANCES

Extracted: phenylethyl malonamide, phenobarbital, carbamazepine, phenytoin, carbamazepine-10,11-epoxide

KEY WORDS

serum; mouse

REFERENCE

Capparella,M.; Foster,W.,III; Larrousse,M.; Phillips,D.J.; Pomfret,A.; Tuvim,Y. Characteristics and applications of a new high-performance liquid chromatography guard column, *J.Chromatogr.A*, **1995**, 691, 141–150.

SAMPLE

Matrix: blood, saliva, urine

Sample preparation: Serum. 100 μ L Serum + 200 μ L MeCN, vortex for 10 s, centrifuge at 1500 g for 5 min, inject a 2 μ L aliquot of the supernatant. Saliva. 250 μ L Saliva + 50 μ L MeCN, centrifuge at 1500 g for 5 min, inject a 2 μ L aliquot of the supernatant. Urine. Condition a Sep-Pak SPE cartridge with 5 mL MeCN then 20 mL water. Add 2 mL urine to the cartridge, wash with 20 mL water, elute with 500 μ L MeCN, inject 2 μ L of the eluent.

HPLC VARIABLES

Guard column: 20 \times 2 3 μ m ODS-Hypersil

Column: 250 \times 2 3 μ m ODS-Hypersil

Mobile phase: MeCN:MeOH:10 mM pH 7.0 phosphate buffer 50:30:110

Column temperature: 40

Flow rate: 0.2

Injection volume: 2

Detector: UV 200

CHROMATOGRAM

Retention time: 5.0

Limit of quantitation: 780 ng/mL

OTHER SUBSTANCES

Simultaneous: p-hydroxyphenobarbital, phenylethylmaleimide, phenobarbital, dihydrodihydroxycarbamazepine, 5-(p-hydroxyphenyl)-5-phenylhydantoin, carbamazepine, 5-(m-hydroxyphenyl)-5-phenylhydantoin, phenytoin, carbamazepine-10,11-epoxide, hexobarbital, nitrazepam, clonazepam

Noninterfering: oxazepam, nordiazepam, cyheptamide, diazepam, prezepam, temazepam, lorazepam, chlordiazepoxide

KEY WORDS

serum; SPE

REFERENCE

Liu,H.; Delgado,M.; Forman,L.J.; Eggers,C.M.; Montoya,J.L. Simultaneous determination of carbamazepine, phenytoin, phenobarbital, primidone and their principal metabolites by high-performance liquid chromatography with photodiode-array detection, *J.Chromatogr.*, **1993**, 616, 105-115.

SAMPLE

Matrix: blood, tissue

Sample preparation: Tissue. Homogenize 20-200 mg brain tissue with 1 mL 1.5 µg/mL IS in 1% ammonium acetate + 1% sodium azide buffer:MeCN 99:1, flush apparatus with 1 mL extraction buffer, add 1 mL acetone, shake for 5 min, centrifuge for 10 min. Add sample to an Extrelut-3 SPE cartridge (Kieselguhr), add 1 mL extraction buffer, wait for 10 min, elute with 15 mL extraction solvent. Evaporate the eluate, take up residue in 50 µL MeOH, add 50 µL water, inject a 10-25 µL aliquot. Serum. 100 µL Serum + 1 mL 1.5 µg/mL IS in 1% ammonium acetate + 1% sodium azide buffer:MeCN 99:1,mix, add 1 mL extraction buffer, mix, add 1 mL acetone, shake for 5 min, centrifuge for 10 min. Add sample to an Extrelut-3 SPE cartridge (Kieselguhr), add 1 mL extraction buffer, wait for 10 min, elute with 15 mL extraction solvent. Evaporate the eluate, take up residue in 50 µL MeOH, add 50 µL water, inject a 10-25 µL aliquot. (Extraction buffer was 20 g NaH₂PO₄·2H₂O + 4.5 g Na₂HPO₄·2H₂O + 1.5 NaN₃ in 1 L water, pH 6. Extraction solvent was dichloromethane:isopropanol 97:3.)

HPLC VARIABLES

Column: 200 × 2.1 5 µm Hypersil ODS

Mobile phase: Gradient. A was MeCN:50 mM (NH₄)H₂PO₄ (pH 4.4) 10:90. B was MeCN:50 mM (NH₄)H₂PO₄ (pH 4.4) 60:40. A:B from 85:15 to 55:45 over 9.5 min, keep at 55:45 for 0.5 min, return to 85:15 over 0.5 min.

Column temperature: 65

Flow rate: 0.3

Injection volume: 10-25

Detector: UV 207

CHROMATOGRAM

Retention time: 3.73

Internal standard: 5-ethyl-5-(p-tolyl)barbituric acid (9.07)

Limit of quantitation: 500 ng/mL

OTHER SUBSTANCES

Simultaneous: phenobarbital, N-desmethylnethsuximide, carbamazepine-10,11-epoxide, phenytoin, carbamazepine

KEY WORDS

serum; SPE; brain

REFERENCE

Juergens,U.; Rambeck,B. Sensitive analysis of antiepileptic drugs in very small portions of human brain by microbore HPLC, *J.Liq.Chromatogr.*, **1987**, 10, 1847-1863.

SAMPLE

Matrix: blood, urine

Sample preparation: Serum. Inject a 5 or 20 µL aliquot directly onto the column with mobile phase A or C. Urine. Inject a 5 µL aliquot directly onto the column with mobile phase C.

HPLC VARIABLES

Column: 100 × 4.6 5-10 µm Silicalite (by sieving Silicalite, 3M Co.(?))

Mobile phase: MeCN:20 mM pH 6.9 phosphate buffer 10:90 (A) or Gradient. MeCN:20 mM pH 6.9 phosphate buffer from 5:95 to 20:80 over 2 min, to 25:75 over 2 min, to 30:70 over 4 min, to 50:50 over 2 min, maintain at 50:50 for 10 min (B) or Gradient. MeCN:20 mM pH 6.9 phosphate buffer 14:86 for 5 min, to 25:75 over 1 min, to 30:70 over 2 min, to 50:50 over 3 min, maintain at 50:50 for 6 min (C)

Flow rate: 1

Injection volume: 5 (A, C), 20 (B)

Detector: UV 254 (serum); UV 230 (urine)

CHROMATOGRAM

Retention time: 2.16 (A, serum), 9.7 (B, serum), 4.5 (C, urine)

Limit of detection: 1 ng (urine)

OTHER SUBSTANCES

Extracted: acetaminophen (B), barbitol (B), carbamazepine (B,C), ethosuximide (A), methamphetamine (A), phenobarbital (B,C), phenytoin (B,C), sulfamethoxazole (A), sulfapyridine (B)

Also analyzed: metabolites

KEY WORDS

serum

REFERENCE

Ambrose,D.L.; Fntz,J.S. High-performance liquid chromatographic determination of drugs and metabolites in human serum and urine using direct injection and a unique molecular sieve, *J.Chromatogr.B*, **1998**, 709, 89–96.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 11.13

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitrityline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenpropofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-stilbene, imipramine, indomethacin, isocarboxystyrene, isocarboxazid, isoniazid, isoproterenol, isox-suprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megesterol, mepacrine, meperidine, mephentermine, mephenytoin, mephesisin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, meth-apyrilene, methaqualone, methazalamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methylodopa, methylodopamine, methylphenidate, methylprednisolone, methyl-testosterone, methypyrrolon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, ox-ymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phenacyclidine, phendi-metrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phenter-mine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, proben-ecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinna-mine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, sco-poletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, thebaine, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thiorida-zine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tran-lycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphen-idyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, **1994**, *18*, 233-242.

Pristinamycin

Merck Index: 7933

SAMPLE**Matrix:** blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the

residue with 50 μL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 227.5

CHROMATOGRAM

Retention time: 17.235

KEY WORDS

whole blood

REFERENCE


Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

Proadifen

Molecular formula: $\text{C}_{23}\text{H}_{31}\text{NO}_2$

Molecular weight: 353.50

CAS Registry No.: 302-33-0, 62-68-0 (HCl)



SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 $\mu\text{g/mL}$ solution in MeOH, inject a 20 μL aliquot.

HPLC VARIABLES

Column: 125 \times 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 2.4

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cy-

clizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipiprone, diprenorphine, dipyridamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypropazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozide, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, 323, 191–225.

Probenecid


Molecular formula: C₁₃H₁₉NO₄S

Molecular weight: 285.36

CAS Registry No.: 57-66-9

Merck Index: 7934

Lednicer No.: 1 135



SAMPLE

Matrix: blood

Sample preparation: Filter plasma (0.22 μm), inject a 10 μL aliquot.

HPLC VARIABLES

Column: 150 × 4.6 5 μm GFF-S5-80 internal-surface reversed phase "Pinkerton" (Regis)

Mobile phase: THF:100 mM potassium phosphate 5:95, pH 7.0

Flow rate: 1

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 7.8

KEY WORDS

plasma; direct injection

REFERENCE

Nakagawa,T.; Shibukawa,A.; Shimono,N.; Kawashima,T.; Tanaka,H.; Haginaka,J. Retention properties of internal-surface reversed-phase silica packing and recovery of drugs from human plasma, *J.Chromatogr.*, **1987**, *420*, 297–311.

SAMPLE

Matrix: blood

Sample preparation: 50 μ L Plasma + 350 μ L MeOH, mix thoroughly, centrifuge at 10000 g for 5 min. Remove the supernatant and evaporate it to dryness under a stream of nitrogen, re-constitute the residue in 200 μ L mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 150 \times 4 5 μ m Cosmosil 5C18 (Nacalai Tesque)

Mobile phase: MeCN:20 mM pH 7.5 phosphate buffer 25:75

Flow rate: 1

Detector: UV 254

CHROMATOGRAM

Retention time: 10

Limit of quantitation: 1000 ng/mL

KEY WORDS

plasma

REFERENCE

Terashita,S.; Sawamoto,T.; Deguchi,S.; Tokuma,Y.; Hata,T. Sex-dependent and independent renal excretion of nilvadipine metabolites in rat: evidence for a sex-dependent active secretion in kidney, *Xenobiotica*, **1995**, *25*, 37–47.

SAMPLE

Matrix: blood

Sample preparation: 150 μ L Plasma + 150 μ L MeCN, vortex, rotate at 20 rpm for 10 min; centrifuge at 1000 g for 10 min. Transfer supernatant to another tube and add 7 volumes dichloromethane, equilibrate for 10 min, rotate at 20 rpm for 10 min; centrifuge at 1000 g for 10 min, inject an aliquot of the upper aqueous layer (*J.Chromatogr.* 1987, 413, 109).

HPLC VARIABLES

Guard column: C18

Column: 150 \times 1.6 Spherisorb S5-ODS2 C18

Mobile phase: MeOH:100 mM pH 3 acetate buffer 50:50

Flow rate: 1

Detector: UV 254

CHROMATOGRAM

Limit of detection: 300 ng/mL

KEY WORDS

plasma; rat

REFERENCE

Gimeno,M.J.; Martínez,M.; Granero,L.; Torres-Molina,F.; Peris,J.-E. Influence of probenecid on the renal excretion mechanisms of cefadroxil, *Drug Metab.Dispos.*, **1996**, *24*, 270–272.

SAMPLE

Matrix: blood, CSF

Sample preparation: Mix 100 μ L plasma or 50 μ L CSF with an equal volume of 5 μ g/mL n-butyl p-hydroxybenzoate in MeCN, centrifuge at 14000 rpm, inject a 20 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 250 \times 4 LiChrospher RP-18e

Mobile phase: MeCN:water:acetic acid 50:49.9:0.1

Flow rate: 1

Injection volume: 20

Detector: UV 244

CHROMATOGRAM

Internal standard: n-butyl p-hydroxybenzoate

KEY WORDS

rat; plasma; pharmacokinetics

REFERENCE

Seki,T.; Sato,N.; Hasegawa,T.; Kawaguchi,T.; Juni,K. Nasal absorption of zidovudine and its transport to cerebrospinal fluid in rats, *Biol.Pharm.Bull.*, **1994**, 17, 1135–1137.

SAMPLE

Matrix: blood, CSF, urine

Sample preparation: Plasma. Add 500 μ L MeCN to 500 μ L plasma while mixing on a Whirl-mixer, rotate for 10 min, centrifuge at 1000 g for 5 min. Remove a 700 μ L aliquot of the supernatant and add it to 3.5 mL dichloromethane, mix for 30 s, centrifuge at 1000 g for 1 min, inject a 20 μ L aliquot of the aqueous layer. Urine. Dilute with Sørensen buffer, inject an aliquot. CSF. Inject an aliquot directly.

HPLC VARIABLES

Column: 100 \times 3 5 μ m MOS-Hypersil C8

Mobile phase: MeCN:MeOH:buffer 12:26:62 containing 3 mM tetrabutylammonium bromide (Buffer was 5 mM pH 5.0 sodium acetate.)

Column temperature: 22

Flow rate: 1

Injection volume: 20

Detector: UV 231

CHROMATOGRAM

Retention time: 11

Limit of detection: 250 ng/mL

OTHER SUBSTANCES

Extracted: penicillin G

KEY WORDS

plasma

REFERENCE

van Gulpen,C.; Brokerhof,A.W.; van der Kaay,M.; Tjaden,U.R.; Mattie,H. Determination of benzylpenicillin and probenecid in human body fluids by high-performance liquid chromatography, *J.Chromatogr.*, **1986**, 381, 365–372.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 350 μ L 2 μ g/mL Naproxen in 10 mM pH 6.0 phosphate buffer containing 0.05% MeOH + 650 μ L pH 6 phosphate buffer + 100 μ L plasma + 0.5 mL 1 M pH 2 phosphate buffer + 10 mL diethyl ether, vortex 1 min, centrifuge at 2000 g for 3 min. Remove organic phase and evaporate it to dryness under a stream of nitrogen. Dissolve residue in 250 μ L mobile phase, vortex for 15 s, inject aliquot. Urine. 350 μ L 20 μ g/mL Naproxen in 10 mM pH 6.0 phosphate buffer containing 0.5% MeOH + 650 μ L pH 6 phosphate buffer + 100 μ L urine + 1 mL 0.5 M pH 7 phosphate buffer + 10 mL diethyl ether, vortex 1 min, centrifuge at 2000 g for 3 min. Remove organic phase and evaporate it to dryness under a stream of nitrogen. Dissolve residue in 250 μ L mobile phase, vortex for 15 s, inject aliquot.

HPLC VARIABLES

Guard column: 40 \times 3.2 30-44 μ m Vydac reverse-phase

Column: 40 × 4.6 5 µm Spherisorb ODS

Mobile phase: MeCN:50 mM pH 7.0 phosphate buffer 6:94 to 8:92

Flow rate: 2

Injection volume: 5-200

Detector: UV 262

CHROMATOGRAM

Retention time: 18

Internal standard: naproxen (10)

Limit of quantitation: 200 ng/mL

OTHER SUBSTANCES

Simultaneous: ketoprofen, fenoprofen, salicylic acid

KEY WORDS

plasma

REFERENCE

Upton,R.A.; Buskin,J.N.; Guentert,T.W.; Williams,R.L.; Riegelman,S. Convenient and sensitive high-performance liquid chromatography assay for ketoprofen, naproxen and other allied drugs in plasma or urine, *J.Chromatogr.*, **1980**, *190*, 119-128.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 100 µL Plasma + 100 µL 100 µg/mL indoprofen in water + 100 µL 600 mM sulfuric acid + 5 mL isooctane:isopropanol 95:5, vortex 30 s, centrifuge at 1800 g for 5 min. Remove organic layer and add 5 mL water to it. Vortex for 30 s, centrifuge for 3 min. Remove organic layer and evaporate it to dryness on a Speed Vac concentrator. Reconstitute residue in 100 µL 50 mM triethylamine in MeCN, vortex 30 s, add 50 µL 60 mM ethyl chloroformate in MeCN, let stand 30 s, add 50 µL 1 M L-leucinamide hydrochloride and 1 M triethylamine in MeOH, let stand 2 min, add 50 µL water, inject 10-60 µL aliquots. Urine. 100 µL Urine + 25 µL 1 M NaOH, add 125 µL 600 mM sulfuric acid, proceed as for plasma.

HPLC VARIABLES

Guard column: 50 × 5 37-53 µm C18 material

Column: 100 × 4.6 5 µm Partisil 5 ODS-3

Mobile phase: MeCN:60 mM KH₂PO₄:triethylamine 35:65:0.1

Flow rate: 1

Injection volume: 10-60

Detector: UV 275

CHROMATOGRAM

Retention time: 18

Internal standard: indoprofen (6(R), 7(S))

Limit of quantitation: 500 ng/mL

OTHER SUBSTANCES

Simultaneous: ketoprofen

Also analyzed: fenoprofen, cicloprofen, pirprofen, flurbiprofen, indoprofen, carprofen

KEY WORDS

plasma; rat; derivatization

REFERENCE

Palylyk,E.L.; Jamali,F. Simultaneous determination of ketoprofen enantiomers and probenecid in plasma and urine by high-performance liquid chromatography, *J.Chromatogr.*, **1991**, *568*, 187-196.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 100 µL Plasma + 100 µL 600 mM sulfuric acid + 4 mL isooctane:isopropanol 95:5, vortex 30 s, centrifuge at 1800 g for 5 min. Remove organic layer and add 4

mL water to it. Vortex for 30 s, centrifuge for 3 min. Remove organic layer and evaporate it to dryness on a Speed Vac concentrator. Reconstitute residue in 200 μ L MeOH, vortex 30 s, add 100 μ L 100 μ g/mL indoprofen in water, 20 μ L aliquots. Urine. 100 μ L Urine + 25 μ L 1 M NaOH, add 125 μ L 600 mM sulfuric acid, proceed as for plasma.

HPLC VARIABLES

Guard column: 50 \times 5 37-53 μ m C18 material

Column: 100 \times 4.6 5 μ m Partisil 5 ODS-3

Mobile phase: MeCN:60 mM KH_2PO_4 :triethylamine 25:75:0.1

Flow rate: 1

Injection volume: 20

Detector: UV 275

CHROMATOGRAM

Retention time: 6.5

Internal standard: indoprofen (3.4)

OTHER SUBSTANCES

Simultaneous: ketoprofen

KEY WORDS

plasma; rat

REFERENCE

Palylyk, E.L.; Jamali, F. Simultaneous determination of ketoprofen enantiomers and probenecid in plasma and urine by high-performance liquid chromatography, *J. Chromatogr.*, **1991**, 568, 187-196.

SAMPLE

Matrix: solutions

Sample preparation: Inject a 50 μ L aliquot of a solution in mobile phase.

HPLC VARIABLES

Column: 100 \times 4.6 5 μ m Spheri-5 RP-8

Mobile phase: MeOH:buffer 30:70 (Prepare buffer by mixing 4 mM Na_2HPO_4 and 7 mM KH_2PO_4 to achieve pH 7.)

Flow rate: 1

Injection volume: 50

Detector: F ex 355 em 460 (408 nm cutoff filter) following post-column extraction. The column effluent mixed with 50 μ g/mL reagent in mobile phase pumped at 0.5 mL/min and then with chloroform pumped at 1 mL/min and the mixture flowed through a 1.8 m \times 0.3 mm ID knitted PTFE coil to a 50 μ L membrane phase separator using a polyethylene-backed 0.5 μ m Fluoropore membrane filter (design in paper). The organic phase flowed to the detector. (Synthesize the reagent, α -(3,4-dimethoxyphenyl)-4'-trimethylammoniummethylcinnamonitrile methosulfate, as follows. Stir 20 mmoles 3,4-dimethoxyphenylacetone and 20 mmoles p-toluidine in 50 mL EtOH at 50°, add 5 mL 50% aqueous KOH slowly, stir at 50° for 5 min, cool to room temperature, filter, dry the precipitate of α -(3,4-dimethoxyphenyl)-4'-methylcinnamonitrile. Dissolve 20 mmoles α -(3,4-dimethoxyphenyl)-4'-methylcinnamonitrile, 20 mmoles N-bromosuccinimide, and 20 mg benzoyl peroxide in 100 mL carbon tetrachloride (Caution! Carbon tetrachloride is a carcinogen!), reflux with stirring for 1.5 h, cool, filter, evaporate to dryness under reduced pressure, recrystallize from MeOH to give α -(3,4-dimethoxyphenyl)-4'-bromomethylcinnamonitrile. Vigorously stir 30 mmoles anhydrous dimethylamine in 100 mL dry benzene (Caution! Benzene is a carcinogen!), very slowly add 10 mmoles α -(3,4-dimethoxyphenyl)-4'-bromomethylcinnamonitrile while stirring at 0°, stir at room temperature overnight, add 150 mL water, remove the organic phase, extract the aqueous phase twice with 100 mL portions of diethyl ether, wash the organic layers with saturated NaCl solution, dry over anhydrous magnesium sulfate, evaporate under reduced pressure to give α -(3,4-dimethoxyphenyl)-4'-dimethylaminomethylcinnamonitrile (*J. Chem. Eng. Data* 1987, 32, 387). Reflux 10 mmoles α -(3,4-dimethoxyphenyl)-4'-dimethylaminomethylcinnamonitrile, 20 mmoles dimethyl sulfate (Caution! Dimethyl sulfate is a carcinogen and acutely toxic!), and 5 g potassium carbonate in 50 mL acetone for 1 h, cool to room temperature, filter, dry the precipitate under vacuum at room temperature overnight, recrystallize from chloroform containing 2-3 drops of 95% EtOH to give α -(3,4-dimethoxyphenyl)-4'-trimethylammoniummethylcinnamonitrile methosulfate (mp 212-215°). Protect solutions from light.)

CHROMATOGRAM**Retention time:** k' 1.2791**Limit of detection:** 10 ng/mL

OTHER SUBSTANCES**Simultaneous:** ibuprofen, ketoprofen, mefenamic acid, naproxen, salicylic acid, valproic acid

KEY WORDS

post-column extraction; post-column reaction

REFERENCE

Kim, M.; Stewart, J. T. HPLC post-column ion-pair extraction of acidic drugs using a substituted α -phenylcinamonitrile quaternary ammonium salt as a new fluorescent ion-pair reagent, *J. Liq. Chromatogr.*, **1990**, *13*, 213-237.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 250 \times 4.6 Zorbax RX**Mobile phase:** Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.**Column temperature:** 30**Flow rate:** 2**Detector:** UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, flurazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenpropofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarboxystiril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methylpyrrolone, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrrolamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinna-

mine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, transylcypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleonnamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, **1994**, *18*, 233-242.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 mm Supelcosil LC-DP

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 7.49

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordiazepoxide, chlorfentanyl, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenopropfen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazinol, mefenamic acid, meperidine, mephentermine, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methylodopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfinpyrazole, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocainide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, trifluopromazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

KEY WORDS

details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J.Chromatogr.A*, **1995**, 692, 103-119.

SAMPLE

Matrix: urine

Sample preparation: 2 mL Urine + 0.5 g solid buffer I (pH 5-5.5), vortex 15 s, add 4 mL ethyl acetate, agitate for 10 min, centrifuge at 600 g for 5 min. Remove organic layer and vortex it with 2 mL 5% aqueous lead acetate for 10 s, centrifuge at 600 g for 5 min, remove and keep organic phase. 2 mL Urine + 0.5 g solid buffer II (pH 9-9.5), vortex 15 s, add 4 mL ethyl acetate, agitate for 10 min, centrifuge at 600 g for 5 min. Remove organic layer and combine it with previous organic layer. Evaporate to dryness at 50° under a stream of nitrogen, reconstitute in 300 μ L 50 μ g/mL β -hydroxyethyltheophylline in MeOH, inject 5 μ L aliquot. (Solid buffer I was KH_2PO_4 : Na_2HPO_4 99:1, solid buffer II was NaHCO_3 : K_2CO_3 3:2.)

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m HP Hypersil ODS (A) or HP LiChrosorb RP-18 (B)

Mobile phase: Gradient. MeCN:buffer from 15:85 at 2 min to 80:20 at 20 min (Buffer was 50 mM NaH_2PO_4 containing 16 mM propylamine hydrochloride, adjusted to pH 3 with concentrated phosphoric acid.)

Flow rate: 1

Injection volume: 5

Detector: UV 230, UV 275

CHROMATOGRAM

Retention time: 15.78 (A), 16.1 (B)

Internal standard: β -hydroxyethyltheophylline (3.7 (A), 4.4 (B))

Limit of detection: 5000 ng/mL

OTHER SUBSTANCES

Extracted: furosemide, metolazone, amiloride, acetazolamide, chlorothiazide, hydrochlorothiazide, quinethazone, triamterene, hydroflumethiazide, chlorthalidone, dichlorphenamide, trichloromethiazide, methyclothiazide, benzthiazide, cyclothiazide, polythiazide, bendroflumethiazide, spironolactone, canrenone, flumethiazide, bumetanide, ethacrynic acid

Noninterfering: acetaminophen, aspirin, caffeine, diflunisal, fenoprofen, ibuprofen, indomethacin, methocarbamol, naproxen, phenylbutazone, sulindac, tetracycline, theobromine, theophylline, tolmetin, trimethoprim, verapamil

REFERENCE

Cooper, S.F.; Massé, R.; Dugal, R. Comprehensive screening procedure for diuretics in urine by high-performance liquid chromatography, *J.Chromatogr.*, **1989**, 489, 65-88.

SAMPLE

Matrix: urine

Sample preparation: Make 5 mL urine alkaline (pH 9-10), add 2 g NaCl, extract twice with 6 mL ethyl acetate. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 200 μ L MeCN/water, inject a 10-20 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 4 5 μ m SGE 100 GL-4 C18P (Scientific Glass Engineering)

Mobile phase: MeCN:MeOH:water:trifluoroacetic acid 15:15:70:0.5

Flow rate: 0.8 or 1

Injection volume: 10-20

Detector: MS, ZAB2-SEQ (VG), PSP source coupled to LC, source 250°, probe 240-260°, scan m/z 200-550 or UV 270

CHROMATOGRAM

Retention time: 3.2

Limit of detection: 50 ng (by MS)

OTHER SUBSTANCES

Extracted: bumetanide, ethacrynic acid, spironolactone

REFERENCE

Ventura,R.; Fraisse,D.; Becchi,M.; Paise,O.; Segura,J. Approach to the analysis of diuretics and masking agents by high-performance liquid chromatography-mass spectrometry in doping control, *J.Chromatogr.*, **1991**, 562, 723-736.

SAMPLE

Matrix: urine

Sample preparation: Buffer urine to 4.9 by mixing with an equal volume of pH 4.9 200 mM sodium phosphate buffer. Inject a 40 μ L aliquot onto column A with mobile phase A, after 3 min backflush the contents of column A onto column B with mobile phase B and start the gradient. At the end of the run re-equilibrate for 10 min.

HPLC VARIABLES

Column: A 20 \times 4.5 μ m Hypersil octadecylsilica ODS; B 200 \times 4.6 5 μ m Shiseido SG-120 polymer-based C18

Mobile phase: A water; B Gradient. MeCN:buffer from 7:93 to 15:85 over 3.5 min, to 50:50 over 8.5 min, maintain at 50:50 for 11 min (Buffer was 6.9 g $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$ in 1 L water, pH adjusted to 3.1 with phosphoric acid.)

Flow rate: 1

Injection volume: 40

Detector: UV 230

CHROMATOGRAM

Retention time: 19.9

Limit of detection: 1000 ng/mL

OTHER SUBSTANCES

Extracted: acetazolamide, amiloride, bendroflumethiazide, benzthiazide, bumetanide, caffeine, carbamazepine, chlorothiazide, chlorthalidone, clopamide, dichlorfenamide, ethacrynic acid, furosemide, hydrochlorothiazide, metyrapone, spironolactone, triamterene, trichlormethiazide

KEY WORDS

column-switching; optimum detection wavelengths vary for each drug

REFERENCE

Saarinén,M.; Sirén,H.; Riekkola,M.-L. A column switching technique for the screening of diuretics in urine by high performance liquid chromatography, *J.Liq.Chromatogr.*, **1993**, 16, 4063-4078.

SAMPLE

Matrix: urine

Sample preparation: 5 mL Urine + 50 μ L 100 μ g/mL 7-propyltheophylline in MeOH + 200 μ L ammonium chloride buffer + 2 g NaCl, extract with 6 mL ethyl acetate by rocking at 40 movements/min for 20 min and centrifuging at 800 g for 5 min, repeat extraction, combine organic layers, evaporate to dryness at 40° under a stream of nitrogen. Reconstitute in 200 μ L MeCN: water 15:85 and inject 20 μ L aliquots. (Ammonium chloride buffer was 28 g ammonium chloride in 100 mL water with the pH adjusted to 9.5 with concentrated ammonia solution.)

HPLC VARIABLES

Column: 75 \times 4.6 3 μ m Ultrasphere ODS

Mobile phase: Gradient. MeCN:100 mM ammonium acetate adjusted to pH 3 with concentrated phosphoric acid. From 10:90 to 15:85 over 2 min to 55:45 over 3 min to 60:40 over 3 min. Kept at 60:40 for 1 min, decreased to 10:90 over 1 min and equilibrated at 10:90 for 2 min.

Flow rate: 1

Injection volume: 20

Detector: UV 270

CHROMATOGRAM

Retention time: 7.3

Internal standard: 7-propyltheophylline (4.5)

Limit of detection: 100 ng/mL

OTHER SUBSTANCES

Simultaneous: xipamide, bumetanide, acetazolamide, amiloride, bendroflumethiazide, buthiazide, benzthiazide, canrenone, caffeine, clopamide, chlorthalidone, cyclothiazide, diclofenamide, ethacrynic acid, furosemide, hydrochlorothiazide, mesocarb, morazone, piretanide, polythiazide, spironolactone, torsemide, triamterene

REFERENCE

Ventura,R.; Nadal,T.; Alcalde,P.; Pascual,J.A.; Segura,J. Fast screening method for diuretics, probenecid and other compounds of doping interest, *J.Chromatogr.A*, **1993**, 655, 233–242.

SAMPLE

Matrix: urine

Sample preparation: Direct injection into column A with mobile phase A for 1 min then back flush onto column B with mobile phase B.

HPLC VARIABLES

Column: A 20 × 2.1 30 μm Hypersil ODS-C18; B 250 × 4 Hypersil ODS-C18

Mobile phase: A Water; B Gradient. MeCN:buffer 15:85 for 1.5 min then to 80:20 over 8 min.

Keep at 80:20 for 2.5 min then re-equilibrate with 15:85. (Buffer was 50 mM NaH₂PO₄ + 1.4 mL propylamine hydrochloride per liter adjusted to pH 3 with concentrated phosphoric acid.)

Flow rate: 1

Injection volume: 50

Detector: UV 230

CHROMATOGRAM

Retention time: 10.3

Limit of detection: 200 ng/mL.

OTHER SUBSTANCES

Simultaneous: bumetanide, ethacrynic acid, acetazolamide, amiloride, bendroflumethiazide, chlorthalidone, cyclothiazide, furosemide, hydrochlorothiazide, spironolactone, triamterene


KEY WORDS

column-switching

REFERENCE

Campíns-Falco,P.; Herráez-Hernández,R.; Sevillano-Cabeza,A. Column-switching techniques for screening of diuretics and probenecid in urine samples, *Anal.Chem.*, **1994**, 66, 244–248.

Probucol



Molecular formula: C₃₁H₄₈O₂S₂

Molecular weight: 516.85

CAS Registry No.: 23288-49-5

Merck Index: 7935

Lednicer No.: 2 126

SAMPLE

Matrix: blood

Sample preparation: Dilute 1 mL serum 0.5-5 times with saline. Add 1 mL 10 μg/mL IS in EtOH to 1 mL diluted serum dropwise while vortexing, add 1.5 mL n-heptane, vortex for 1 min, centrifuge at 3000 rpm (Labofuge) for 15 min. Remove 1.3 mL of the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 40 μL MeCN:THF 50:50, inject a 5 μL aliquot.

HPLC VARIABLES

Column: 200 × 2.1 5 μm ODS Hypersil

Mobile phase: MeCN:water:THF 81.3:5.7:13

Column temperature: 40

Flow rate: 0.4

Injection volume: 5

Detector: UV 244

CHROMATOGRAM

Retention time: 2.763

Internal standard: 2-pentanone bis(3,5-di-tert)mercaptole (3.666)

Limit of detection: 2 ng

OTHER SUBSTANCES

Extracted: vitamin E (α -tocopherol), gamma-tocopherol, vitamin A (retinol), lycopene, α -carotene, β -carotene, metabolites

KEY WORDS

serum

REFERENCE

Schäfer Elinder,L.; Walldius,G. Simultaneous measurement of serum probucol and lipid-soluble antioxidants, *J.Lipid Res.*, **1992**, 33, 131–137.

SAMPLE

Matrix: blood

Sample preparation: 50 μ L Plasma + 400 μ L 20 μ g/mL IS in EtOH + 800 μ L isooctane + 500 μ L water, vortex for 20 s, centrifuge at 960 g for 5 min. Remove the upper organic layer and evaporate it to dryness under a stream of nitrogen at room temperature, reconstitute the residue in 50 μ L MeCN, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 200 \times 3 5 μ m 5 μ m Hypersil ODS

Mobile phase: MeCN:water 96:4

Flow rate: 1

Injection volume: 10

Detector: UV 241

CHROMATOGRAM

Retention time: 5

Internal standard: 4,4'-[1-methylbutylidene-bis(thio)]-bis[2,6-bis(1,1-dimethylethyl)]phenol (MDL 27272) (7)

Limit of detection: 500 ng/mL

KEY WORDS

plasma; rabbit

REFERENCE

Nourooz-Zadeh,J.; Gopaul,N.K.; Forster,L.A.; Ferns,G.A.; Ånggård,E.E. Measurement of plasma probucol levels by high-performance liquid chromatography, *J.Chromatogr.B*, **1994**, 654, 55–60.

Procainamide


Molecular formula: C₁₃H₂₁N₃O

Molecular weight: 235.33

CAS Registry No.: 51-06-9, 614-39-1 (HCl)

Merck Index: 7936

Lednicer No.: 1 14



SAMPLE

Matrix: activated neutrophils

Sample preparation: Centrifuge, inject a 100 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 150 mm long 5 μ m Spherisorb ODS2

Mobile phase: MeCN:water:acetic acid:triethylamine 20:80:1:0.05

Flow rate: 1

Injection volume: 100

Detector: UV

CHROMATOGRAM

Retention time: 3.6

OTHER SUBSTANCES

Extracted: metabolites

REFERENCE

Uetrecht, J.P.; Zahid, N.; Whitfield, D. Metabolism of vesnarinone by activated neutrophils: Implications for vesnarinone-induced agranulocytosis, *J. Pharmacol. Exp. Ther.*, **1994**, 270, 865–872.

SAMPLE

Matrix: blood

Sample preparation: 50 μ L Serum + 50 μ L buffer + 50 μ L 10 μ g/mL N-propionylprocainamide in MeCN + 3 mL chloroform:isopropanol 95:5, mix for 30 s, centrifuge at 3000 g for 3 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 50 μ L MeOH, inject a 20 μ L aliquot. (Buffer was 2.9 g/L sodium bicarbonate in water, adjusted to pH 11.0 with NaOH.)

HPLC VARIABLES

Column: μ Bondapak C18

Mobile phase: MeCN:buffer 9.75:90.25 (Buffer was 100 mM KH_2PO_4 adjusted to pH 4.0 with phosphoric acid.) (At the end of each day clean with water for 20 min and MeOH for 30 min.)

Flow rate: 2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 3.8

Internal standard: N-propionylprocainamide (6)

Limit of detection: 500 ng/mL

OTHER SUBSTANCES

Extracted: N-acetylprocainamide

Simultaneous: dyphylline, theophylline, caffeine, aspirin, salicylic acid, acetaminophen

Noninterfering: benzoic acid

KEY WORDS

serum

REFERENCE

Ou, C.-N.; Frawley, V.L. Theophylline, dyphylline, caffeine, acetaminophen, salicylate, acetylsalicylate, procainamide, and N-acetylprocainamide determined in serum with a single liquid-chromatographic assay, *Clin. Chem.*, **1982**, 28, 2157–2160.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma + 5 mL MTBE, place on a reciprocating shaker at low speed for 15 min, centrifuge at 1250 g for 10 min. Remove organic layer and evaporate under a stream of nitrogen. Reconstitute in 200 μ L mobile phase, inject 100 μ L aliquot.

HPLC VARIABLES

Guard column: 20 × 4.6 Supelco Pelliguard C18

Column: 150 × 4.6 3 μm Hypersil ODS

Mobile phase: MeCN:7 mM sodium heptanesulfonate 18:82, containing 1% glacial acetic acid and 0.035% triethylamine

Flow rate: 0.8

Injection volume: 100

Detector: UV 272

CHROMATOGRAM

Retention time: 10.1

Internal standard: procainamide.HCl

OTHER SUBSTANCES

Simultaneous: hydrochlorothiazide

Noninterfering: aspirin, acetaminophen, ibuprofen

KEY WORDS

plasma; procainamide is IS

REFERENCE

Azumaya,C.T. Sensitive liquid chromatographic method for the determination of hydrochlorothiazide in human plasma, *J.Chromatogr.*, **1990**, 532, 168–174.

SAMPLE

Matrix: blood

Sample preparation: 100 μL Serum or plasma + 100 μL 500 mM sodium carbonate + 100 μL 15 μg/mL N-propionylprocainamide in water, vortex for 5 s, add 0.5 (procainamide) or 1 (tocainide) mL dichloromethane, vortex for 30 (procainamide) or 60 (tocainide) s, centrifuge at 9500 g for 1 min. Remove the lower organic layer and add it to 200 μL 10 mM HCl, vortex for 15 s, centrifuge, inject a 20 μL aliquot of the aqueous layer.

HPLC VARIABLES

Column: 100 × 5 NovaPak cyano HP radial compression

Mobile phase: MeCN:buffer 10:90, final pH adjusted to 6.0 (Buffer was 5 mM acetate buffer containing 0.05% triethylamine.)

Flow rate: 1.5

Injection volume: 20

Detector: UV 280

CHROMATOGRAM

Retention time: 3.7

Internal standard: N-propionylprocainamide (6.1)

Limit of detection: 1 μg/mL

OTHER SUBSTANCES

Extracted: N-acetylprocainamide, tocainide

Simultaneous: disopyramide, lidocaine, mexiletine, quinidine

Noninterfering: carbamazepine, desmethyldoxepin, digoxin, doxepin, ethosuximide, lithium, phenobarbital, phenytoin, primidone, propranolol, theophylline, valproic acid

KEY WORDS

serum; plasma

REFERENCE

vasBinder,E., Annesley,T. Liquid chromatographic analysis of mexiletine in serum, with alternate application to tocainide, procainamide, and N-acetylprocainamide, *Biomed.Chromatogr.*, **1991**, 5, 19–22.

SAMPLE

Matrix: blood

Sample preparation: Condition a 100 mg LRC Bond Elut unendcapped cyanopropyl SPE cartridge with 1 mL MeOH and 1 mL water. 1 mL Serum or plasma + 250 μ L buffer, vortex, centrifuge at 1300 rpm for 10 min (plasma only), add 1 mL sample to SPE cartridge, wash twice with 1 mL water, dry with nitrogen for 30 s, elute with two 250 μ L aliquots of MeOH:buffer 60:40, inject 100 μ L aliquot of buffer. (Buffer was 50 mM KH_2PO_4 and Na_2HPO_4 50:50 v/v adjusted to pH 6.0 with phosphoric acid.)

HPLC VARIABLES

Guard column: 15 \times 4.6 7 μ m Brownlee RP-18

Column: 250 \times 4.6 5 μ m Spherisorb ODS-1

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 50 mM KH_2PO_4 and Na_2HPO_4 50:50 v/v adjusted to pH 6.0 with phosphoric acid.)

Column temperature: 48

Flow rate: 1.25

Injection volume: 100

Detector: UV 320

CHROMATOGRAM

Retention time: 5.4

Internal standard: procainamide

OTHER SUBSTANCES

Simultaneous: ranitidine

KEY WORDS

serum; plasma; robotic sample preparation; procainamide is IS; SPE

REFERENCE

Lloyd,T.L.; Perschy,T.B.; Gooding,A.E.; Tomlinson,J.J. Robotic solid phase extraction and high performance liquid chromatographic analysis of ranitidine in serum or plasma, *Biomed.Chromatogr.*, **1992**, 6, 311–316.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 μ L mobile phase, centrifuge at 2800 g for 5 min, inject a 50 μ L aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 \times 3.9 4 μ m NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH_2PO_4 adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 287

CHROMATOGRAM

Retention time: 3.65

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfapyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihy-

dralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; mclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vindesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loperazole; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacine; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozone; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254-262.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 25 μ L MeOH + 15 μ L 6 M NaOH + 2 mL ethyl acetate: isopropanol 96:4, shake mechanically for 20 min, centrifuge at 3000 rpm for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100 μ L MeOH, inject a 75 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: MeCN:10 mM pH 5.1 KH_2PO_4 8:92

Flow rate: 2.5

Injection volume: 75

Detector: UV 330

CHROMATOGRAM

Retention time: 3

Internal standard: procainamide

OTHER SUBSTANCES

Extracted: ranitidine

Simultaneous: lidocaine

Noninterfering: brompheniramine, chlorpheniramine, cimetidine, diazepam, diclofenac, glyburide, ibuprofen, ketoprofen, metoclopramide, naproxen, phenylbutazone, verapamil

KEY WORDS

plasma; procainamide is IS

REFERENCE

al-Khamis, K.I.; El-Sayed, Y.M.; Al-Rashood, K.A.; Bawazir, S.A. High-performance liquid chromatographic determination of ranitidine in human plasma, *J. Liq. Chromatogr.*, **1995**, *18*, 277-286.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. 2 mL Plasma + 100 μ L water + 660 μ L 2 M perchloric acid, shake briefly, centrifuge at 3000 rpm for 5 min. Remove 1.5 mL of the supernatant and adjust the pH to 9 with 150 μ L 4 M NaOH and 4 mL 500 mM boric acid/KCl buffer, add 12 mL water-saturated n-pentanol:chloroform 20:60, shake for 20 min, centrifuge at 3000 rpm for 10 min. Remove the organic phase and dry it over anhydrous sodium sulfate, add 8 mL to 300 μ L 100 mM HCl, shake for 10 min, centrifuge at 3000 rpm for 10 min, inject a 25-100 μ L aliquot of the aqueous phase. Urine. Dilute urine with pH 9 borate buffer, add 12 mL water-saturated n-pentanol:chloroform 20:60, shake for 20 min, centrifuge at 3000 rpm for 10 min. Remove the organic phase and dry it over anhydrous sodium sulfate, add 8 mL to 300 μ L 100 mM HCl, shake for 10 min, centrifuge at 3000 rpm for 10 min, inject a 25-100 μ L aliquot of the aqueous phase.

HPLC VARIABLES

Column: μ Bondapak C18

Mobile phase: MeOH:water:acetic acid 38:61:1 containing 1-heptanesulfonic acid (PIC B7)

Flow rate: 1.5

Injection volume: 25-100

Detector: F ex 235 em no filter

CHROMATOGRAM

Retention time: 7

Internal standard: procainamide

OTHER SUBSTANCES

Extracted: sotalol

KEY WORDS

plasma; procainamide is IS

REFERENCE

Lefebvre, M.A.; Girault, J.; Saux, M.C.; Fourtillan, J.B. Fluorometric high-performance liquid chromatographic determination of sotalol in biological fluids, *J.Pharm.Sci.*, **1980**, 69, 1216-1217.

SAMPLE

Matrix: formulations

Sample preparation: 1 mL Sample + 9 mL mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Spherisorb phenyl

Mobile phase: MeCN:500 mM KH₂PO₄:water 22:10:68 adjusted to pH 7.1 with 10 M NaOH

Flow rate: 2

Injection volume: 20

Detector: UV 268

CHROMATOGRAM

Retention time: 6

OTHER SUBSTANCES

Simultaneous: milrinone, degradation products

KEY WORDS

stability-indicating; 5% dextrose; injections

REFERENCE

Riley, C.M. Stability of milrinone and digoxin, furosemide, procainamide hydrochloride, propranolol hydrochloride, quinidine gluconate, or verapamil hydrochloride in 5% dextrose injection, *Am.J.Hosp.Pharm.*, **1988**, 45, 2079-2091.

SAMPLE**Matrix:** formulations**Sample preparation:** Dilute a 1 mL sample to 10 mL with mobile phase, inject an aliquot.

HPLC VARIABLES**Column:** 150 × 4.6 5 µm Spherisorb Phenyl**Mobile phase:** MeCN:water:500 mM KH₂PO₄ 15:75:10, pH 6.8**Flow rate:** 2**Injection volume:** 20**Detector:** UV 268

CHROMATOGRAM**Retention time:** 6

OTHER SUBSTANCES**Simultaneous:** amrinone

KEY WORDS

injections; stability-indicating; 5% dextrose; 0.45% NaCl

REFERENCE

Riley,C.M.; Junkin,P. Stability of amrinone and digoxin, procainamide hydrochloride, propranolol hydrochloride, sodium bicarbonate, potassium chloride, or verapamil hydrochloride in intravenous admixtures, *Am.J.Hosp.Pharm.*, **1991**, 48, 1245–1252.

SAMPLE**Matrix:** formulations**Sample preparation:** Make up 1 mL syrup to 50 mL with water. Remove a 2 mL aliquot and add it to 2 mL 0.1% procaine hydrochloride in water, make up to 100 mL with water, filter (0.45 µm), inject a 50 µL aliquot.

HPLC VARIABLES**Column:** 250 × 4.6 Partisil ODS-3 C18**Mobile phase:** MeCN:buffer:triethanolamine:water 15:80:0.2:4.8, pH adjusted to 4.5 with glacial acetic acid (Buffer was 4.72 g sodium acetate and 1.8 mL acetic acid in 1 L water.)**Flow rate:** 1.2**Injection volume:** 50**Detector:** UV 254

CHROMATOGRAM**Retention time:** 4.6**Internal standard:** procaine (7.5)

KEY WORDS

syrup; stability-indicating

REFERENCE

Alexander,K.S.; Pudipeddi,M.; Parker,G.A. Stability of procainamide hydrochloride syrups compounded from capsules, *Am.J.Hosp.Pharm.*, **1993**, 50, 693–698.

SAMPLE**Matrix:** solutions**Sample preparation:** Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.

HPLC VARIABLES**Column:** 125 × 4.9 Spherisorb S5W silica**Mobile phase:** MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7**Flow rate:** 2**Injection volume:** 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 3.6

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipranone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethiopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypropazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenylglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procaine, prochlorperazine, procyclidine, proheptazine, prolantane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, *323*, 191–225.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 12 µm Dynamax C18 (Rainin)

Mobile phase: MeOH:21 mM pH 4.4 NaH₂PO₄ 13.5:86.5

Flow rate: 1.5

Detector: UV 266

CHROMATOGRAM

Retention time: 2.3

OTHER SUBSTANCES

Simultaneous: metabolites, N-acetylprocainamide

REFERENCE

Hickman,D.; Palamanda,J.R.; Unadkat,J.D.; Sim,E. Enzyme kinetic properties of human recombinant arylamine N-acetyltransferase 2 allotypic variants expressed in *Escherichia coli*, *Biochem.Pharmacol.*, **1995**, 50, 697–703.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4 5 µm LiChrospher 100 RP-8 (B)

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 6.32 (A), 3.52 (B)

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordi-azepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxy-chloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazi-ndol, mefenamic acid, meperidine, mephentoin, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymet-azoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, phen-iramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltolox-amine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propanthe-line, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, race-methorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, ser-traline, sotalol, spirinolactone, sulfipyrazone, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocin-ide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

KEY WORDS

details of plasma extraction

REFERENCE

Koves,E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J.Chromatogr.A*, **1995**, 692, 103–119.