

---

**CHROMATOGRAM****Retention time:**  $k'$  0.38

---

**OTHER SUBSTANCES**

**Also analyzed:** amitriptyline, barbital, benzoic acid, butabarbital, clomipramine, clonazepam, desipramine, diazepam, flurazepam, furosemide, imipramine, nitrazepam, phenobarbital, phenol, phenolphthalein, pindolol, propranolol, salicylic acid, secobarbital, terbutaline, xylazine

---

**KEY WORDS**

effect of mobile phase pH on capacity factor is discussed

---


**REFERENCE**

Walshe,M.; Kelly,M.T.; Smyth,M.R.; Ritchie,H. Retention studies on mixed-mode columns in high-performance liquid chromatography, *J.Chromatogr.A*, **1995**, 708, 31–40.

---

---

# Retinoic acid

**Molecular formula:**  $C_{20}H_{28}O_2$ **Molecular weight:** 300.44**CAS Registry No.:** 302-79-4 (tretinoin (all-trans)), 4759-48-2 (isotretinoin (13-cis)), 5300-03-8 (alitretinoin (9-cis))**Merck Index:** 8333**Lednicer No.:** 3 12 (isotretinoin)

---

**SAMPLE****Matrix:** blood

**Sample preparation:** Evaporate 25-100  $\mu$ L IS in MeOH to near dryness under a stream of nitrogen. Add 0.2-1 mL plasma and 100  $\mu$ L buffer, extract with 2 mL diethyl ether:ethyl acetate 50:50 for 5 min. Centrifuge at 2000 g for 10 min at 4°, evaporate the organic phase to dryness. Dissolve the residue in 30-100  $\mu$ L MeOH, inject an aliquot. (Solution was prepared in yellow amber glass and all handling was performed in a room with dim yellow light! Buffer was 25 mM  $KH_2PO_4$  containing 40 mM  $Na_2HPO_4$ , pH 7.)

---

**HPLC VARIABLES****Column:** 250  $\times$  4.6 4  $\mu$ m Nova-Pak C18

**Mobile phase:** Gradient. A was MeCN:MeOH:THF 33.25:61.75:5. B was 2% acetic acid. A:B from 75:25 to 88:12 over 11 min, maintain at 88:12 for 19 min, return to initial condition at 30 min, equilibrate for 10 min.

**Flow rate:** 1**Injection volume:** 25**Detector:** UV 350

---

**CHROMATOGRAM****Retention time:** 24.9 (isotretinoin), 26.7 (tretinoin), 26.2 (alitretinoin (9-cis))**Internal standard:** acitretin (21), 13-cis-acitretin (19.8)**Limit of quantitation:** 2 ng/mL

---

**OTHER SUBSTANCES****Extracted:** metabolites

**Noninterfering:** acetaminophen, acyclovir, alprazolam, amikacin, amitriptyline, amphotericin B, aspirin, atenolol, bromazepam, caffeine, carbamazepine, ceftriaxone, chlorpromazine, cimetidine, clonazepam, dextromethorphan, diazepam, erythromycin, flunitrazepam, haloperidol, ketoconazole, lorazepam, meprobamate, metronidazole, methylprednisolone, miconazole, midazolam, nifedipine, nitrazepam, netilmicin, nordiazepam, nystatin, oxazepam, phenytoin, prednisolone, prednisone, sulconazole, theophylline, thiopental, zidovudine

---

**KEY WORDS**

plasma; rabbit; rat

---

**REFERENCE**

Disdier,B.; Bun,H.; Catalin,J.; Durand,A. Simultaneous determination of all-trans-, 13-cis, 9-cis-retinoic acid and their 4-oxometabolites in plasma by high-performance liquid chromatography, *J.Chromatogr.B*, **1996**, 683, 143-154.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 400  $\mu$ L Plasma + 1.5 mL EtOH containing IS, freeze at -20° for 30 min, centrifuge. Inject a 1.4 mL aliquot of the supernatant onto column A and elute to waste with mobile phase A (time not given). Elute the contents of column A onto column B with mobile phase B, monitor the effluent from column B.

---

**HPLC VARIABLES**

**Column:** A C18-Corasil or C18-Lichrospher; B two 250  $\times$  4 Supersphere 100 RP-18 endcapped columns in series

**Mobile phase:** A. MeCN:acetic acid (ratio not given) containing 1% ammonium acetate; B. Gradient. MeCN:acetic acid:10% ammonium acetate:water 60:1:6:30, 95:2:0.5:2, and 99:0.5:0.0.5 (times not given).

**Detector:** UV 360

---

**CHROMATOGRAM**

**Internal standard:** acitretin

**Limit of quantitation:** 300 pg/mL

---

**OTHER SUBSTANCES**

**Extracted:** metabolites, 4-oxo-isotretinoin, 4-oxo-tretinoin

---

**KEY WORDS**

plasma; column-switching; pharmacokinetics; for isotretinoin and tretinoin

---

**REFERENCE**

Chen,C.; Mistry,G.; Jensen,B.; Heizmann,P.; Timm,U.; van Brummelen,P.; Rakhit,A.K. Pharmacokinetics of retinoids in women after meal consumption or vitamin A supplementation, *J.Clin.Pharmacol.*, **1996**, 36, 799-808.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Mix 100  $\mu$ L serum with 200  $\mu$ L MeCN and 10  $\mu$ L 20 mM ascorbic acid, centrifuge at 16000 g for 5 min. Mix the supernatant with 200  $\mu$ L water, inject a 2  $\mu$ L aliquot. (Protect all solutions from light.)

---

**HPLC VARIABLES**

**Column:** 50  $\times$  0.18 3  $\mu$ m ODS-AQ (YMC, Wilmington, NC) (The separation capillary column was formed from fused-silica capillaries (Polymicro Technologies, Phoenix) by inserting a small 50  $\mu$ m I.D. capillary ca. 15 mm into a larger 180  $\mu$ m I.D. capillary and fixed by applying epoxy (No. 353ND, Epoxy Technology, Billerica MA). A glass filter paper frit (Whatman GF/A) was inserted into the larger capillary and forced against the smaller capillary with a stream of isopropanol. The stationary phase was suspended in 3 mL isopropanol and pumped into the larger capillary until a 50 mm bed was formed. The larger and smaller diameter capillaries extended no more than 100 and 16 mm from the frit, respectively.)

**Mobile phase:** MeCN:MeOH:water 65:2.5:32.5 containing 1% tetrabutylammonium perchlorate, adjusted to pH 5.0 with acetic acid and 174 mM sodium acetate

**Flow rate:** 0.004

**Injection volume:** 2

**Detector:** E, carbon-fiber working electrode +900 mV, Ag/AgCl reference electrode (details of preparation in paper)

---

**CHROMATOGRAM**

**Retention time:** 9.0 (isotretinoin), 9.5 (tretinoin)

**Limit of detection:** 410 pg/mL (isotretinoin), 64 pg/mL (tretinoin)

**Limit of quantitation:** 49.6 fmol (isotretinoin), 87.6 fmol (tretinoin)

---

---

**OTHER SUBSTANCES**

**Extracted:** retinaldehyde, retinol

---

**KEY WORDS**

cow; serum; capillary HPLC

---

**REFERENCE**

Hagen, J.J.; Washco, K.A.; Monnig, C.A. Determination of retinoids by reversed-phase capillary liquid chromatography with amperometric electrochemical detection, *J. Chromatogr. B*, **1996**, 677, 225–231.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 500  $\mu$ L Plasma + 750  $\mu$ L 100 ng/mL acitretin in MeCN:9 mM NaOH 20:80, centrifuge at 1500 g for 3 min, inject a 500  $\mu$ L aliquot onto column A with mobile phase A and elute for 7 min, elute column A in backflush mode with mobile phase A for 3 min, backflush contents of column A onto column B with mobile phase B and start the gradient for mobile phase B. At the end of the process flush the lines with component B of mobile phase B, re-equilibrate columns for 4 min. (Keep sample at 10° in the autosampler.)

---

**HPLC VARIABLES**

**Column:** A 14  $\times$  4.6 37-50  $\mu$ m Bondapak C18 Corasil (column fitted with 3  $\mu$ m sieves not glass fiber filters); B 30  $\times$  4.5  $\mu$ m Spherisorb ODS 1 + 125  $\times$  4.5  $\mu$ m Spherisorb ODS 1 + 125  $\times$  4.5  $\mu$ m Spherisorb ODS 1

**Mobile phase:** A MeCN:1% ammonium acetate 10:90; B Gradient. A was MeCN:water:10% ammonium acetate:acetic acid 600:400:4:30. B was MeCN:water:10% ammonium acetate:acetic acid 850:146:4:10. A:B 100:0 to 70:30 over 6 min, then to 0:100 over 5 min, stay at 0:100 for 11 min.

**Flow rate:** A 1.5; B 1

**Injection volume:** 500

**Detector:** UV 360

---

**CHROMATOGRAM**

**Retention time:** 25 (isotretinoin), 27 (tretinoin)

**Internal standard:** acitretin (23)

**Limit of detection:** 0.5-1 ng/mL

**Limit of quantitation:** 2 ng/mL

---

**OTHER SUBSTANCES**

**Simultaneous:** tretinoin, 4-oxoisotretinoin, 4-oxotretinoin, metabolites

---

**KEY WORDS**

plasma; column-switching

---

**REFERENCE**

Wyss, R. Determination of retinoids in plasma by high-performance liquid chromatography and automated column switching, *Methods Enzymol.*, **1990**, 189, 146–155.

---

**SAMPLE**

**Matrix:** culture media

**Sample preparation:** 100  $\mu$ L Culture media + 200  $\mu$ L ice-cold EtOH, mix thoroughly, let stand for 15 min, centrifuge at 12000 g for 15 min, inject an aliquot of the supernatant.

---

**HPLC VARIABLES**

**Guard column:** Whatman CO:PELL ODS guard column

**Column:** 100  $\times$  8.5  $\mu$ m Nova-Pak C18 (radial-packed)

**Mobile phase:** MeOH:100 mM pH 7.0 ammonium acetate 90:10

**Flow rate:** 1

**Detector:** UV 340

---

**CHROMATOGRAM**

**Retention time:** 9.44 (tretinoin)

---

---

**OTHER SUBSTANCES**

**Extracted:** isotretin, motretinid, acitretin, Vitamin A (retinol), retinal, etretinate

---

**REFERENCE**

Kochhar,D.M.; Penner,J.D.; Minutella,L.M. Biotransformation of etretinate and developmental toxicity of etretin and other aromatic retinoids in teratogenesis bioassays, *Drug Metab.Dispos.*, **1989**, 17, 618–624.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Weigh out 1.25 g 0.1% isotretinoin cream or 0.05% isotretinoin gel, make up to 10 mL with MeOH, vortex for 15 min. Cut open two 10 mg isotretinoin capsules, add about 40 mL MeOH, sonicate for 10 min, filter, make up to 100 mL with MeOH. Grind 10% isotretinoin beadlets to powder, weigh out a 1.25 mg aliquot. Add 100  $\mu$ L cream, gel, or capsule solution to 0.5 mL stainless steel extraction cartridge partially filled with Celite, complete filling with Celite, load the cartridge into the extraction chamber. Alternatively, directly add 1.25 mg beadlets powder to the extraction chamber (Supercritical Fluid Extractor, Isco, model SFX 2-10). Extract using the following conditions: chamber and restrictor temperature 45°, pressure 325 atm, static extraction time 2.5 min, dynamic extraction time 5 min, mobile phase MeOH:CO<sub>2</sub> 5:95, solvent trap 17 mL MeOH. After extraction cool the MeOH extract to room temperature for 10 min, make up to 25 mL with MeOH, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu$ m Spherisorb ODS-2

**Mobile phase:** MeCN:MeOH:0.05% glacial acid 42.5:32.5:25

**Flow rate:** 1.2

**Injection volume:** 20

**Detector:** UV 360

---

**CHROMATOGRAM**

**Retention time:** 25.0 (isotretinoin), 30.3 (alitretinoin (9-*cis*)), 33.2 (tretinoin)

---

**OTHER SUBSTANCES**

**Extracted:** 11,13-di-*cis* retinoic acid, 9,13-di-*cis* retinoic acid, 9-*cis* retinoic acid, tretinoin

---

**KEY WORDS**

beadlets; capsule; cream; gel; SFE

---

**REFERENCE**

Simmons,B.R.; Chukwumerije,O.; Stewart,J.T. Supercritical fluid extraction of 13-*cis* retinoic acid and its photoisomers from selected pharmaceutical dosage forms, *J.Pharm.Biomed.Anal.*, **1997**, 16, 395–403.

---

**SAMPLE**

**Matrix:** milk

**Sample preparation:** Mix 50 mL milk with 30 mL ethanolic KOH (10:30). Saponify the mixture at 80° for 20 min. Extract twice with 10 mL n-hexane. Evaporate the extracts to dryness, reconstitute the residue with 1 mL mobile phase, inject a 5  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5 $\mu$ m C18 (Alltech)

**Mobile phase:** MeOH:EtOH 80:20

**Flow rate:** 1

**Injection volume:** 5

**Detector:** UV 250

---

**CHROMATOGRAM**

**Retention time:** 2.97 (isotretinoin), 3.38 (tretinoin)

---

**OTHER SUBSTANCES**

**Extracted:** retinal, vitamin A, vitamin D2, vitamin D3, vitamin E, vitamin K1, vitamin K2

---

**REFERENCE**

Gong,B.Y.; Ho,J.W. Simultaneous separation and detection of ten common fat-soluble vitamins in milk, *J.Liq.Chromatogr.Rel.Technol.*, **1997**, 20, 2389–2397.

---

**SAMPLE**

**Matrix:** silicone oils

**Sample preparation:** Condition a 1 g Si Bond-Elut SPE cartridge with 5 mL n-hexane. Mix 1 g silicone oil with 2 mL dichloromethane, vortex for 2 min, centrifuge at 3000 g. Withdrawn the supernatant, repeat this procedure twice, filter (0.45  $\mu$ m), heat the filtrate at 50°, expose to a stream of helium for 30 min. Add 2.5  $\mu$ g retinol acetate, 2.5  $\mu$ g  $\alpha$ -tocopherol acetate, and 25  $\mu$ g BHT. Add the mixture to the SPE cartridge, elute with 500  $\mu$ L MeOH, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 150  $\times$  4.6 5  $\mu$ m Zorbax C8

**Mobile phase:** Gradient. A was MeCN:200 mM ammonium acetate 72:25. B was MeOH:water 95:5. A:B 100:0 for 10 min, to 0:100 over 1 min, maintain at 0:100 for 14 min

**Flow rate:** 2 for 10 min then 1.5

**Injection volume:** 20

**Detector:** UV 350

---

**CHROMATOGRAM**

**Retention time:** 2

**Internal standard:** retinol acetate (9.5)

**Limit of detection:** 79.2 ng/mL

**Limit of quantitation:** 264.1 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** cholesterol (UV 210), retinal (UV 350),  $\alpha$ -tocopherol acetate (UV 210), vitamin A (UV 350), vitamin E (UV 210)

---

**KEY WORDS**

ophthalmic silicone oils; SPE

---

**REFERENCE**

Del Nozal,M.J.; Bernal,J.L.; Marinero,P. Simultaneous HPLC determination of cholesterol,  $\alpha$ -tocopherol, retinol, retinal and retinoic acid in silicone oils used as vitreous substitutes in eye surgery, *J.Liq.Chromatogr.Rel.Technol.*, **1996**, 19, 1151–1167.

---

---


# Ribavirin

**Molecular formula:** C<sub>8</sub>H<sub>12</sub>N<sub>4</sub>O<sub>5</sub>

**Molecular weight:** 244.21

**CAS Registry No.:** 36791-04-5

**Merck Index:** 8365



---

**SAMPLE**

**Matrix:** blood, CSF

**Sample preparation:** Prepare a boronate affinity gel SPE column by packing 100 mg Affi-Gel 601 (Bio-Rad) into a 65  $\times$  6 1.5 mL polypropylene column, condition with 10 mL buffer, store at 4°, immediately before use condition further with two 1 mL aliquots of buffer. 500  $\mu$ L Serum, plasma, or CSF + 500  $\mu$ L buffer + 25  $\mu$ L 100  $\mu$ g/mL 3-methylcytidine methosulfate in water, mix, add to the SPE column, wash with five 1 mL portions of buffer, elute with two 1 mL portions of 100 mM formic acid. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute the residue in 100  $\mu$ L mobile phase, inject the whole amount. Buffer was 250 mM pH 8.8 ammonium acetate.)

---

**HPLC VARIABLES**

**Column:** two 300 × 4 µm Bondapak C18 columns in series

**Mobile phase:** 10 mM Ammonium phosphate adjusted to pH 2.5 with 85% phosphoric acid

**Flow rate:** 1.5

**Injection volume:** 100

**Detector:** UV 235

---

**CHROMATOGRAM**

**Retention time:** 8.0

**Internal standard:** 3-methylcytidine methosulfate (11.0)

**Limit of detection:** 100 nM

---

**OTHER SUBSTANCES**

**Noninterfering:** acetaminophen, N-acetylprocainamide, acyclovir, amikacin, 5-aminoimidazole-4-carboxamide-1-β-D-ribofuranoside, amitriptyline, azidothymidine, caffeine, carbamazepine, chloramphenicol, cyclosporin A, cytidine, desipramine, diazepam, digoxin, disopyramide, ethosuximide, gentamicin, imipramine, kanamycin, lidocaine, lithium, methotrexate, 1-methyladenosine, 5-methylcytidine, 7-methylguanosine, 7-methylinosine, netilmicin, nortriptyline, phenobarbital, phenytoin, primidone, procainamide, quinidine, salicylic acid, theophylline, 2-thiocytidine, tobramycin, uridine, valproic acid, vancomycin, zalcitabine

---

**KEY WORDS**

serum; plasma; SPE

---

**REFERENCE**

Granich, G.G.; Krogstad, D.J.; Connor, J.D.; Desrochers, K.L.; Sherwood, C. High-performance liquid chromatography (HPLC) assay for ribavirin and comparison of the HPLC assay with radioimmunoassay, *Antimicrob. Agents Chemother.*, **1989**, 33, 311–315.

---

**SAMPLE**

**Matrix:** blood, tissue, tracheal aspirate

**Sample preparation:** Dilute tracheal aspirates 10–100 fold with 250 mM pH 8.8 ammonium acetate buffer. Homogenize (glass/PTFE homogenizer) mouse lung with 1 mL water, centrifuge at 13000 g. Filter (Amicon CF25) 1 mL plasma, serum, diluted tracheal aspirate, or lung homogenate while centrifuging at 500 relative centrifugal force (RCF) for 30 min, dilute ultrafiltrate with a volume of 2.5 M pH 8.8 ammonium acetate buffer equal to 10% of the volume of ultrafiltrate, add the mixture to the SPE cartridge, wash with 7 mL 250 mM pH 8.8 ammonium acetate buffer, elute with 6 mL 100 mM formic acid. Lyophilize the eluate, reconstitute with a volume of mobile phase equal to the original volume of ultrafiltrate, inject a 50 µL aliquot. (Soak the CF25 filters in water for 1 h before use. Prepare SPE cartridges by slurrying Matrex PBA-60 gel (88.87 µmole boron per mL, Amicon) in 250 mM pH 8.8 ammonium acetate and packing it into a 40 × 7 column to a bed volume of 1 mL, pass 50 mL 250 mM pH 8.8 ammonium acetate through the column.)

---

**HPLC VARIABLES**

**Guard column:** Guard Pak C18 (Waters)

**Column:** 250 × 4.6 5 µm Microsorb C18

**Mobile phase:** MeOH:buffer 1:99 adjusted to pH 5.10 with 10% ammonium hydroxide or 8% phosphoric acid (Buffer was 20 mM ammonium phosphate. Following each analysis flush column with a gradient of MeCN:water from 0:100 to 80:20.)

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 207

---

**CHROMATOGRAM**

**Retention time:** 5.3

**Limit of detection:** 100 ng/mL

---

**KEY WORDS**

mouse; lung; ultrafiltrate; SPE

---

**REFERENCE**

Smith, R.H.A.; Gilbert, B.E. Quantification of ribavirin in biological fluids and tissues by high-performance liquid chromatography, *J. Chromatogr.*, **1987**, *414*, 202–210.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Serum. 1 mL Serum + 5 µg uridine, vortex, add 1 mL water, vortex, filter (Centricon 30 membrane, 30 000 exclusion limit) while centrifuging at 2750 g for 2 h, wash the filtrate with 1.5 mL dichloromethane, pass the aqueous layer through a 10 mm i.d. column containing 2.3 mL wet 50-100 mesh Dowex 1X4-100 (chloride form) (Sigma) on top of 0.8 mL wet 60-150 mesh Lewatit S-1080 (hydrogen form) (Merck), elute with 3 mL water. Collect all the eluate and evaporate it to dryness under vacuum at 30°. reconstitute in 600 µL water, inject a 20-100 µL aliquot. Urine. 200 µL Urine + 5 µg uridine, vortex for a few s, add the urine mixture dropwise with stirring to 3 mL MeOH, centrifuge at 10° at 10000 g for 20 min. Remove the supernatant and add it to 2 mL n-hexane, mix, discard the hexane layer. Evaporate the MeOH layer to dryness under a nitrogen and vacuum, reconstitute the residue in 3 mL water, inject a 100 µL aliquot (*J. Chromatogr.* 1978, 160 169).

---

**HPLC VARIABLES**

**Guard column:** Guard Pak RP-18 (Waters)

**Column:** 250 × 5 7 µm LiChrosorb C18

**Mobile phase:** Water

**Flow rate:** 1

**Injection volume:** 20-100

**Detector:** UV 207

---

**CHROMATOGRAM**

**Retention time:** 6.9

**Internal standard:** uridine (13)

**Limit of detection:** 100 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** metabolites

---

**KEY WORDS**

serum; ultrafiltrate

---

**REFERENCE**

Paroni, R.; Sirtori, C.R.; Borghi, C.; Kienle, M.G. High-performance liquid chromatographic determination of ribavirin in serum and urine and of its urinary metabolite 1,2,4-triazole-3-carboxamide, *J. Chromatogr.*, **1987**, *420*, 189–196.

---

**SAMPLE**

**Matrix:** tissue

**Sample preparation:** Homogenize hamster brain with phosphate-buffered saline, centrifuge at 1600 g for 10 min. Remove the supernatant and add it to EtOH, heat at 90° to remove EtOH, suspend the protein-free sample in mobile phase, centrifuge at 8300 g for 5 min, inject a 10 µL aliquot of the supernatant.

---

**HPLC VARIABLES**

**Column:** TSKgel ODS-120T (Tosoh)

**Mobile phase:** MeCN:25 mM pH 2.5 buffer 2:98

**Flow rate:** 1

**Injection volume:** 10

**Detector:** UV 226

---

**CHROMATOGRAM**

**Retention time:** 4.7

**Limit of detection:** 10 µg/g

---

**KEY WORDS**

hamster; brain; some interference from endogenous substances; pharmacokinetics

---

**REFERENCE**

Ishii,T.; Hosoya,M.; Mori,S.; Shigeta,S.; Suzuki,H. Effective ribavirin concentration in hamster brains for antiviral chemotherapy for subacute sclerosing panencephalitis, *Antimicrob.Agents Chemother.*, **1996**, *40*, 241–243.

---

**SAMPLE**

**Matrix:** urine

**Sample preparation:** Dilute urine ten-fold with water, filter (0.2  $\mu\text{m}$ ), inject a 50  $\mu\text{L}$  aliquot of the filtrate.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 Lichrosorb RP18

**Mobile phase:** Water

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 220

---

**CHROMATOGRAM**

**Limit of quantitation:** 5  $\mu\text{g/mL}$

---

**KEY WORDS**

mouse

---

**REFERENCE**

Ryan,D.M.; Ticehurst,J.; Dempsey,M.H.; Penn,C.R. Inhibition of influenza virus replication in mice by GG167 (4-guanidino-2,4-dideoxy-2,3-dehydro-*N*-acetylneuramic acid) is consistent with extracellular activity of viral neuraminidase (sialidase), *Antimicrob.Agents Chemother.*, **1994**, *38*, 2270–2275.

---


# Riboflavin

**Molecular formula:**  $\text{C}_{17}\text{H}_{20}\text{N}_4\text{O}_6$

**Molecular weight:** 376.37

**CAS Registry No.:** 83-88-5, 130-40-5 (phosphate sodium)

**Merck Index:** 8367



---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Filter (Amicon 25000 molecular-weight cut-off) serum, inject an aliquot.

---

**HPLC VARIABLES**

**Guard column:** 50  $\times$  4.6 25  $\mu\text{m}$  pellicular reversed-phase (Whatman)

**Column:** 300  $\times$  3.9  $\mu\text{m}$  Bondapak

**Mobile phase:** Gradient. A was 20 mM pH 5.6  $\text{KH}_2\text{PO}_4$ . B was MeOH:water 60:40. A:B from 100:0 to 60:40 over 35 min.

**Flow rate:** 1.5

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 6

---

**OTHER SUBSTANCES**

**Extracted:** creatinine, cytidine, hypoxanthine, inosine, kynurenine, 5-methylcytidine, tryptophan, tyrosine, uric acid, xanthine



---

**KEY WORDS**

dog; human; serum

---

**REFERENCE**

Assenza, S.P.; Brown, P.R. Comparison of high-performance liquid chromatographic serum profiles of humans and dogs, *J. Chromatogr.*, **1980**, *181*, 169–176.

---

**SAMPLE**

**Matrix:** blood, formulations, urine

**Sample preparation:** Tablets. Powder tablets, dissolve in water, inject a 10  $\mu$ L aliquot. Injections. Dilute with water, inject a 10  $\mu$ L aliquot. Plasma, urine. Condition a Lichrolut RP-18 (Merck) SPE cartridge with 3 mL MeOH and 3 mL water. Mix 40  $\mu$ L plasma or 100  $\mu$ L urine with twice the volume of MeCN for 2 min, add 100  $\mu$ L water, centrifuge at 3500 rpm for 15 min, evaporate the supernatant under nitrogen at 45° to remove the organic solvents, add slowly to the SPE cartridge, collect the eluate. Evaporate to dryness under a stream of nitrogen at 45°. Reconstitute the residue with 500  $\mu$ L MeOH containing 4.2  $\mu$ g/mL IS. Inject a 10  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.5  $\mu$ m Lichrosorb RP-18

**Mobile phase:** Gradient. A was MeOH. B was 50 mM ammonium acetate. A:B from 5:95 to 15:85 over 6 min, to 30:70 over 7 min, maintain at 30:70 over 7 min

**Flow rate:** 1

**Injection volume:** 10

**Detector:** UV 270

---

**CHROMATOGRAM**

**Retention time:** 19.67

**Internal standard:** xanthine (4.65)

**Limit of detection:** 2.5 ng

---

**OTHER SUBSTANCES**

**Extracted:** ascorbic acid, folic acid, niacin, niacinamide, vitamin B12

---

**KEY WORDS**

plasma; SPE; tablets; injections

---

**REFERENCE**

Papadoyannis, I.N.; Tsioni, G.K.; Samanidou, V.F. Simultaneous determination of nine water and fat soluble vitamins after SPE separation and RP-HPLC analysis in pharmaceutical preparations and biological fluids, *J. Liq. Chromatogr. Rel. Technol.*, **1997**, *20*, 3203–3231.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50  $\mu$ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood)  $\mu$ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18

**Column:** 250  $\times$  4.6  $\mu$ m Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 267.7

---

#### CHROMATOGRAM

**Retention time:** 7.182

---

#### KEY WORDS

whole blood

---

#### REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

---

#### SAMPLE

**Matrix:** formula, milk

**Sample preparation:** Mix 8.0 g powdered infant milk with 10 mL water to it. Mix the diluted powder or 10.5 g liquid infant milk with 1 g solid trichloroacetic acid, shake thoroughly with magnetic stirring for 10 min, centrifuge at 1250 g for 10 min, add 3 mL 4% trichloroacetic acid to the solid residue, mix thoroughly for 10 min, centrifuge, discard the solid phase. Combine the two acid extracts and make up to 10 mL with 4% trichloroacetic acid, filter (0.45  $\mu$ m), inject an aliquot of the filtrate.

---

#### HPLC VARIABLES

**Guard column:** 5  $\mu$ m Tracer Spherisorb ODS 2 C18 (Teknokroma, Spain)

**Column:** 250  $\times$  4.6 5  $\mu$ m Tracer Spherisorb ODS 2 C18 (Teknokroma, Spain)

**Mobile phase:** MeOH:buffer 15:85 (Buffer was 5 mM octanesulfonic acid and 0.5% triethylamine, pH 3.6.)

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 261 for 6 min, UV 287 for 2 min, UV 290 for 5 min, UV 282 for 3 min, UV 268 for 3.5 min, UV 361 for 20.5 min, UV 246 for 20 min

---

#### CHROMATOGRAM

**Retention time:** 17

**Limit of quantitation:**  $\leq$ 50 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** thiamine, pyridoxine, vitamin B12, folic acid, niacinamide, pyridoxal, pyridoxamine

---

#### REFERENCE

Albalá-Hurtado,S.; Veciana-Nogués,M.; Izquierdo-Pulido,M.; Mariné-Font,A. Determination of water-soluble vitamins in infant milk by high-performance liquid chromatography, *J.Chromatogr.A*, **1997**, 778, 247-253.

---

#### SAMPLE

**Matrix:** formulations

**Sample preparation:** Pulverize tablets if necessary. Add tablets to 100 mL 5 mM pH 4.5 potassium phosphate buffer, sonicate at 75 W for 2 min, cool to room temperature, make up to 200 mL with buffer, filter (0.45  $\mu$ m), inject a 10  $\mu$ L aliquot of the filtrate.

---

#### HPLC VARIABLES

**Column:** 250  $\times$  4.6 10  $\mu$ m LiChrosorb NH2 aminopropyl

**Mobile phase:** MeCN:5 mM  $\text{KH}_2\text{PO}_4$  87:13 (Wash column with MeCN:water 10:90 at the end of the day.)

**Column temperature:** 25

**Flow rate:** 2

**Injection volume:** 10

**Detector:** UV 210

---

**CHROMATOGRAM**

**Retention time:** 2.5

---

**OTHER SUBSTANCES**

**Simultaneous:** pantothenic acid, thiamine, niacinamide, pyridoxine

---

**KEY WORDS**

tablets

---

**REFERENCE**

Hudson,T.J.; Allen,R.J. Determination of pantothenic acid in multivitamin pharmaceutical preparations by reverse-phase high-performance liquid chromatography, *J.Pharm.Sci.*, **1984**, 73, 113–115.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Tablets without iron. Grind 5 tablets to a fine powder, add 10 mL mono-thioglycerol and 800 mL buffer, sonicate for 30 min, add 150 mL MeOH, make up to 1 L with buffer, filter (GF/C paper), discard first few mL, remove a 10 mL aliquot, make up to 25 mL with mobile phase, inject an aliquot. Tablets with dioctyl sodium sulfosuccinate. Grind 5 tablets to a fine powder, add 10 mL 2-monothioglycerol and 1 g barium chloride, make up to 1 L with buffer, stir vigorously for 30 min, filter (GF/C paper), discard first few mL, inject an aliquot. Capsules with iron. Contents of one capsule + 5 mL 2-monothioglycerol + 2 mL glacial acetic acid + 75 mL buffer, sonicate for 5 min, make up to 100 mL with buffer, stir vigorously for 30 min, filter (GF/C paper), add 300 mg cupferron, stir for 10 min, let stand for 1 h at room temperature, filter (GF/C paper), let stand for 30 min, filter again (if necessary), discard first few mL, inject an aliquot. (Buffer was 48 mL glacial acetic acid and 10 mL triethylamine in 1 L water, adjust pH to  $3.6 \pm 0.05$  with acetic acid or triethylamine, make up to 1.7 L with water.)

---

**HPLC VARIABLES**

**Column:** 100 × 8 Radial Pak A C18 (Waters)

**Mobile phase:** MeOH:buffer 15:85 (Buffer was 2.20 g sodium heptanesulfonate, 100 mg EDTA, 48 mL glacial acetic acid, and 10 mL triethylamine made up to 1.7 L with water, adjust pH to  $3.6 \pm 0.05$  with acetic acid or triethylamine.)

**Flow rate:** 2

**Injection volume:** 10

**Detector:** UV 280

---

**CHROMATOGRAM**

**Retention time:** 15

---

**OTHER SUBSTANCES**

**Simultaneous:** niacinamide, thiamine, pyridoxine, ascorbic acid (UV 254)

---

**KEY WORDS**

multi-vitamin; protect from light; tablets; capsules

---

**REFERENCE**

Lam,F.-L.; Holcomb,I.J.; Fusari,S.A. Liquid chromatographic assay of ascorbic acid, niacinamide, pyridoxine, thiamine, and riboflavin in multivitamin-mineral preparations, *J.Assoc.Off.Anal.Chem.*, **1984**, 67, 1007–1011.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Dilute injections with water, inject a 50  $\mu$ L aliquot. Dissolve tablets or capsule contents in water (warm if necessary), filter (0.5  $\mu$ m PTFE), inject a 50  $\mu$ L aliquot of the filtrate. (Dissolve tablets or other formulations containing proteinaceous material in water at 60°, add 5% trichloroacetic acid (to pH 4.4), filter, inject a 50  $\mu$ L aliquot.)

---

**HPLC VARIABLES**

**Guard column:** pellicular Corasil

**Column:** 10  $\mu$ m  $\mu$ Bondapak C18

**Mobile phase:** Gradient. A was prepared by dissolving 1 g sodium dioctylsulfosuccinate in 170 mL MeOH, add 10 mL concentrated formic acid, make up to 800 mL with water, adjust pH to 2.5 with 1 M KOH, make up to 1 L. B was prepared by dissolving 1 g sodium dioctylsulfosuccinate in 450 mL MeOH, add 10 mL concentrated formic acid, make up to 800 mL with water, adjust pH to 4.6, make up to 1 L. A:B 100:0 for 19 min then 0:100 (step gradient) or A:B from 100:0 to 0:100 over 25 min (concave curve 9), maintain at 0:100 for 3 min, return to initial conditions over 2 min.

**Flow rate:** 1.5

**Injection volume:** 50

**Detector:** UV 254

---

#### CHROMATOGRAM

**Retention time:** 4 (step gradient), 8 (curve gradient), 2 (or riboflavin-5'-phosphate (either system))

---

#### OTHER SUBSTANCES

**Simultaneous:** folic acid (UV 280), niacin, niacinamide, pyridoxamine (UV 280), thiamine, pyridoxine (UV 280), ascorbic acid (UV 280)

---

#### KEY WORDS

injections; capsules; tablets

---

#### REFERENCE

Woollard, D.C. New ion-pair reagent for the high-performance liquid chromatographic separation of B-group vitamins in pharmaceuticals, *J.Chromatogr.*, **1984**, 301, 470-476.

---

#### SAMPLE

**Matrix:** formulations

**Sample preparation:** Weigh out 500 mg ground tablets, extract with water, make up to 50 or 100 mL with water, filter, inject an aliquot.

---

#### HPLC VARIABLES

**Column:** 250  $\times$  4.6 Nucleosil 10 C18

**Mobile phase:** MeOH:1% acetic acid 25:75

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 270

---

#### CHROMATOGRAM

**Retention time:** 20

---

#### OTHER SUBSTANCES

**Simultaneous:** menadione hydrogen sulfite, niacinamide, pyridoxine, thiamine, ascorbic acid

---

#### KEY WORDS

tablets; multi-vitamin

---

#### REFERENCE

Sadlej-Sosnowska, N.; Blitek, D.; Wilczynska-Wojtulewicz, I. Determination of menadione sodium hydrogen sulfite and nicotinamide in multivitamin formulations by high-performance liquid chromatography, *J.Chromatogr.*, **1986**, 357, 227-232.

---

#### SAMPLE

**Matrix:** formulations

---

#### HPLC VARIABLES

**Column:** 100  $\times$  4.3  $\mu$ m Hypersil BDS-C18

**Mobile phase:** Gradient. MeCN:water adjusted to pH 2.1 from 0.3:99.7 to 25:75 over 11 min

**Flow rate:** 0.5

**Detector:** UV 220

---

**CHROMATOGRAM**

**Retention time:** 10

---

**OTHER SUBSTANCES**

**Simultaneous:** biotin, caffeine, citric acid, folic acid, niacinamide, niacin, pantothenic acid, saccharin, thiamine, pyridoxine, vitamin B12, ascorbic acid

---

**KEY WORDS**

tablets

---

**REFERENCE**

Hewlett Packard Leaflet 12-5091-7351 EUS, 1993.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Dilute liquid multivitamin formulations, filter (0.45  $\mu\text{m}$ ), inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 100  $\times$  2.1 3  $\mu\text{m}$  Spherisorb ODS-2

**Mobile phase:** MeOH:buffer 20:80 containing 0.1% triethylamine (Buffer was 10 mM  $\text{KH}_2\text{PO}_4$  containing 5 mM sodium hexanesulfonate adjusted to pH 2.8 with phosphoric acid.)

**Flow rate:** 0.2 for 5 min, to 0.3 over 0.5 min, 0.3 for 12.5 min

**Injection volume:** 5

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 16

**Limit of detection:** 1.80 ng

---

**OTHER SUBSTANCES**

**Simultaneous:** folic acid (UV 280), pyridoxine (UV 280), niacin, niacinamide, thiamine

---

**KEY WORDS**

liquid multivitamins; degas solutions with helium; protect from light; narrow bore

---

**REFERENCE**

Blanco,D.; Sánchez,L.A.; Gutiérrez,M.D. Determination of water soluble vitamins by liquid chromatography with ordinary and narrow-bore columns, *J.Liq.Chromatogr.*, 1994, 17, 1525–1539.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Dilute liquid multivitamin formulations, filter (0.45  $\mu\text{m}$ ), inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4 5  $\mu\text{m}$  Lichrosorb RP-8

**Mobile phase:** Gradient. A was 10 mM  $\text{KH}_2\text{PO}_4$  containing 5 mM sodium hexanesulfonate adjusted to pH 2.8 with phosphoric acid. B was MeOH. A:B from 90:10 to 71.8:28.2 over 4 min, maintain at 71.8:28.2 for 1.5 min, to 50:50 over 6.5 min, maintain at 50:50 for 5 min, return to initial conditions over 5 min

**Flow rate:** 1

**Injection volume:** 5

**Detector:** UV 272

---

**CHROMATOGRAM**

**Retention time:** 11.12

**Internal standard:** theobromine (8)

**Limit of detection:** 0.465 ng

---

**OTHER SUBSTANCES**

**Simultaneous:** folic acid, niacin, niacinamide, thiamine, pyridoxine (UV 290)

---

**KEY WORDS**

liquid multivitamins; degas solutions with helium; protect from light

---

**REFERENCE**

Blanco,D.; Sánchez,L.A.; Gutiérrez,M.D. Determination of water soluble vitamins by liquid chromatography with ordinary and narrow-bore columns, *J.Liq.Chromatogr.*, **1994**, 17, 1525–1539.

---

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 5 µm Accubond Amino (J & W)

**Mobile phase:** MeCN:buffer 10 :90 (Buffer was 20 mM phosphoric acid adjusted to pH 3.0 with 20 mM NaOH.)

**Flow rate:** 1

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 2.3

---

**OTHER SUBSTANCES**

**Simultaneous:** p-aminobenzoic acid, niacinamide, pyridoxal, pyridoxamine, thiamine, pyridoxine, vitamin B12

---

**REFERENCE**

*J & W Catalog*, 1992-3, p. 277.

---

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 33 × 4.6 3 µm Supelcosil LC-8-DB

**Mobile phase:** MeOH:buffer 15:85 (Buffer was 4.3 mM sodium hexanesulfonate containing 0.1% triethylamine, adjusted to pH 2.8 with phosphoric acid.)

**Column temperature:** 35

**Flow rate:** 1

**Detector:** UV 200

---

**CHROMATOGRAM**

**Retention time:** 3.7

---

**OTHER SUBSTANCES**

**Simultaneous:** niacin, pantothenic acid, pyridoxine, thiamine, niacinamide, ascorbic acid

---

**REFERENCE**

*Rainin Catalog*, C1-94, **1994**, p. 780.

---

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 100 × 4.6 Spheri-5 RP-8

**Mobile phase:** Gradient. A was 100 mM pH 4.7 acetate buffer. B was MeCN:100 mM pH 4.7 acetate buffer 25:75.

**Column temperature:** 26

**Flow rate:** 4

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 4.2

---

**OTHER SUBSTANCES**

**Simultaneous:** niacin, pyridoxine, thiamine, niacinamide, ascorbic acid

---

**REFERENCE**

*Rainin Catalog, C1-94, 1994, p. 7.21.*

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 150 × 4.6 5 µm Inertsil ODS-2

**Mobile phase:** MeCN:50 mM KH<sub>2</sub>PO<sub>4</sub> 90:10

**Flow rate:** 1

**Detector:** UV 210

---

**CHROMATOGRAM**

**Retention time:** 10.5

---

**OTHER SUBSTANCES**

**Simultaneous:** biotin, folic acid, niacin, pantothenic acid, niacinamide

---

**REFERENCE**

*MetaChem Catalog, 1995, p. 21.*

---

**SAMPLE**

**Matrix:** tissue

**Sample preparation:** Condition a 500 mg Nucleosil C18 SPE cartridge with 2 mL MeOH, 2 mL MeOH containing 5 mM sodium heptanesulfonate, and two 2 mL portions of water. Suspend 5 g homogenized tissue with 35 mL 10 mM HCl, autoclave at 121° for 30 min, add 2 mL 25 mg/mL taka-distase (Fluka) in 2.5 M sodium acetate, add 2 mL 10 (muscle) or 20 (liver) mg/mL clara-distase (Fluka) in water, add 2 mL 50 mg/mL papain (Merck) in water, adjust pH to 4.5, heat at 37° for 16-18 h, filter (paper), adjust pH to 6.5, filter again, make up to 50 mL with water, add 4 mL to the SPE cartridge, wash with 2 mL MeOH:water 20:80 containing 5 mM sodium heptanesulfonate, elute with 2 mL MeOH:water 50:50 containing 5 mM sodium heptanesulfonate, inject a 50 µL aliquot.

---

**HPLC VARIABLES**

**Guard column:** 20 × 4.6 10 µm Nucleosil C18

**Column:** 150 × 4.6 3 µm Nucleosil C18

**Mobile phase:** MeCN:10 mM pH 3.0 KH<sub>2</sub>PO<sub>4</sub> 16:84 (muscle) or 15:85 (liver) containing 5 mM sodium heptane sulfonate (Wash with MeCN:water 20:80 at the end of the day, store column in MeCN.)

**Column temperature:** 45

**Injection volume:** 50

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 4.5 (muscle), 5 (liver)

**Limit of detection:** 160 ng/g

---

**OTHER SUBSTANCES**

**Extracted:** thiamine

---

**KEY WORDS**

pig; muscle; liver; protect from light; SPE

---

**REFERENCE**

Barna,I.; Dworschák,E. Determination of thiamine (vitamin B1) and riboflavin (vitamin B2) in meat and liver by high-performance liquid chromatography, *J.Chromatogr.A*, **1994**, 668, 359–363.

---

**SAMPLE**

**Matrix:** urine

**Sample preparation:** Centrifuge urine at 1400 g for 10 min, inject a 50 µL aliquot.

---

**HPLC VARIABLES**

**Guard column:** 70 × 2.1 Co:Pell ODS

**Column:** 300 × 4 10 µm µBondapak C18

**Mobile phase:** MeOH:buffer 35:65 (Buffer was 10 mM KH<sub>2</sub>PO<sub>4</sub> adjusted to pH 5.0 with 1 M NaOH.)

**Flow rate:** 2

**Injection volume:** 50

**Detector:** F ex 320-400 em 400-700 (filter)

---

**CHROMATOGRAM**

**Retention time:** 4.5 (riboflavin), 3 (riboflavin-5-phosphate)

**Limit of detection:** 50 ng/mL

---

**KEY WORDS**

protect from light

---

**REFERENCE**

Smith,M.D. Rapid method for determination of riboflavin in urine by high-performance liquid chromatography, *J.Chromatogr.*, **1980**, 182, 285–291.

---

**SAMPLE**

**Matrix:** urine

**Sample preparation:** Preserve urine with oxalic acid, centrifuge at 2000 g for 10 min, inject a 25-100 µL aliquot of the supernatant.

---

**HPLC VARIABLES**

**Guard column:** Bondapak C18/Corasil

**Column:** 300 × 3.9 µBondapak C18

**Mobile phase:** MeOH:water 34:66

**Flow rate:** 1

**Injection volume:** 25-100

**Detector:** F ex 450 em 530 (long-pass filter)

---

**CHROMATOGRAM**

**Retention time:** 6.5

**Limit of detection:** 10 ng/mL

---

**OTHER SUBSTANCES**

**Simultaneous:** lumichrome, lumiflavin

**Noninterfering:** acriflavine, apresoline, corticosterone, dipyridamole, estradiol, estriol, estrone, ethylenediamine, hydrocortisone, prednisolone, succinic acid, tetrahydrocortisone, testosterone, urobilin

---

**KEY WORDS**

protect from light


---

**REFERENCE**

Gatautis,V.J.; Naito,H.K. Liquid-chromatographic determination of urinary riboflavin, *Clin.Chem.*, **1981**, 27, 1672–1675.



# Ricinoleic acid



**Molecular formula:**  $C_{18}H_{34}O_3$

**Molecular weight:** 298.47

**CAS Registry No.:** 141-22-0

**Merck Index:** 8378

## SAMPLE

**Matrix:** solutions

**Sample preparation:** Prepare a 0.17-1.7 mg/mL solution in MeOH, inject a 30  $\mu$ L aliquot.

## HPLC VARIABLES

**Column:** 250  $\times$  4.6 5  $\mu$ m Ultrasphere ODS

**Mobile phase:** Gradient. A was MeOH containing 0.05% acetic acid. B was water containing 0.05% acetic acid. A:B 85:15 to 100:0 over 40 min, maintain at 100:0.

**Flow rate:** 1

**Injection volume:** 30

**Detector:** evaporative light scattering (ELSD, MK VIII, Varex), drift tube 75°, nitrogen flow 1 L/min, nitrogen pressure 22 psi or UV 205

## CHROMATOGRAM

**Retention time:** 8.95

## OTHER SUBSTANCES

**Simultaneous:** linoleic acid, fatty acids

## REFERENCE

Lin, J.-T.; McKeon, T.A.; Stafford, A.E. Gradient reversed-phase high-performance liquid chromatography of saturated, unsaturated and oxygenated free fatty acids and their methyl esters, *J.Chromatogr.A*, **1995**, 699, 85-91.


# Rifabutin

**Molecular formula:**  $C_{46}H_{62}N_4O_{11}$

**Molecular weight:** 847.02

**CAS Registry No.:** 72559-06-9

**Merck Index:** 8380



## SAMPLE

**Matrix:** blood

**Sample preparation:** Condition a 1 mL Bond Elut C8 SPE cartridge with 1 mL MeOH and 1 mL water. 1 mL Plasma + 25  $\mu$ L 2  $\mu$ g/mL sulindac in water, vortex briefly, centrifuge at 630 rpm for 5 min, add to the SPE cartridge, wash with 1 mL water, elute with 1 mL MeOH. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute the residue in 250  $\mu$ L mobile phase, inject a 100  $\mu$ L aliquot.

## HPLC VARIABLES

**Guard column:** 15  $\times$  3.2 RP-8 (Brownlee, ABI)

**Column:** 250 × 4.6 5 µm Zorbax RX C8

**Mobile phase:** MeCN:buffer 47:53 (Buffer was 50 mM KH<sub>2</sub>PO<sub>4</sub> containing 50 mM sodium acetate, pH adjusted to 4.0 with acetic acid.)

**Flow rate:** 1

**Injection volume:** 100

**Detector:** UV 275

---

#### CHROMATOGRAM

**Retention time:** 10.8

**Internal standard:** sulindac (6.9)

**Limit of quantitation:** 5 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** metabolites

**Simultaneous:** atevirdine, delavirdine, U-89,255, U-96,183

---

#### KEY WORDS

plasma; protect from light; pharmacokinetics; rugged; SPE

---

#### REFERENCE

Lau,Y.Y.; Hanson,G.D.; Carel,B.J. Determination of rifabutin in human plasma by high-performance liquid chromatography with ultraviolet detection, *J.Chromatogr.B*, **1996**, 676, 125–130.

---

#### SAMPLE

**Matrix:** blood, urine

**Sample preparation:** Plasma. 1 mL Plasma + 500 µL 50 mM pH 7.4 KH<sub>2</sub>PO<sub>4</sub>, mix, add 2 mL dichloromethane:isooctane 40:60, vortex for 10 min, centrifuge at 1200 g for 10 min, repeat extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 300 µL mobile phase, add 1 mL n-hexane, wash, centrifuge at 1200 g for 5 min, discard hexane layer, repeat wash, inject a 200 µL aliquot of the aqueous phase. Urine. Inject an aliquot directly.

---

#### HPLC VARIABLES

**Guard column:** 30-38 µm CO:PELL ODS

**Column:** 300 × 3.9 10 µm µBondapak C18

**Mobile phase:** MeCN:50 mM KH<sub>2</sub>PO<sub>4</sub> 40:60

**Flow rate:** 1

**Injection volume:** 200

**Detector:** UV 275

---

#### CHROMATOGRAM

**Retention time:** 24.0

**Limit of detection:** 15 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** metabolites

---

#### KEY WORDS

plasma; rat; human; rabbit; monkey

---

#### REFERENCE

Battaglia,R.; Pianezzola,E.; Salgarollo,G.; Zini,G.; Strolin Benedetti,M. Absorption, disposition and preliminary metabolic pathway of <sup>14</sup>C-rifabutin in animals and man, *J.Antimicrob.Chemother.*, **1990**, 26, 813–822.

---

#### SAMPLE

**Matrix:** blood, urine

**Sample preparation:** Plasma. 1 mL Plasma + 20 µL 5 µg/mL medazepam in MeOH + 1 mL buffer + 7 mL hexane:ethyl acetate 80:20, vortex for 10 s, centrifuge at 1200 rpm for 10 min, repeat extraction. Combine the organic layers and evaporate them to dryness under reduced pressure at 40°, reconstitute with 250 µL hexane:ethyl acetate 80:20, add 200 µL 50 mM phosphoric acid, vortex for 10 s, centrifuge at 1200 rpm for 1-2 min, freeze in dry ice/isopro-

panol, discard the organic layer, wash the frozen aqueous layer with 1 mL hexane, thaw, place under reduced pressure at room temperature for 10 min, add 100  $\mu$ L 250 mM ammonium acetate in MeOH, mix, inject a 100  $\mu$ L aliquot. Urine. 1 mL Urine + 20  $\mu$ L 50  $\mu$ g/mL medazepam in MeOH + 50  $\mu$ L 3.6 M aqueous sulfuric acid, mix. Remove a 200  $\mu$ L aliquot and add it to 800  $\mu$ L 100 mM ammonium acetate in MeOH:water 40:60, mix, inject a 100  $\mu$ L aliquot. (Buffer was 250 mM  $\text{KH}_2\text{PO}_4$  containing 50 mM sodium 1-heptanesulfonate adjusted to pH 7.4.)

---

**HPLC VARIABLES**

**Guard column:** 15  $\times$  3.2 NewGuard RP-18 (Brownlee)

**Column:** 250  $\times$  4.6 5  $\mu$ m Ultrasphere ODS

**Mobile phase:** MeCN:buffer 38:62 (Buffer was 50 mM  $\text{KH}_2\text{PO}_4$ :triethylamine 61.5:0.5 adjusted to pH 4.2 with phosphoric acid.)

**Column temperature:** 40

**Flow rate:** 1

**Injection volume:** 100

**Detector:** UV 275

---

**CHROMATOGRAM**

**Retention time:** 22

**Internal standard:** medazepam (19)

**Limit of quantitation:** 100 ng/mL (urine), 5 ng/mL (plasma)

---

**OTHER SUBSTANCES**

**Extracted:** metabolites

**Simultaneous:** clofazimine

**Noninterfering:** amikacin, isoniazid, streptomycin, zidovudine

---

**KEY WORDS**

plasma; pharmacokinetics

---

**REFERENCE**

Lewis, R.C.; Hatfield, N.Z.; Narang, P.K. A sensitive method for quantitation of rifabutin and its desacetyl metabolite in human biological fluids by high-performance liquid chromatography (HPLC), *Pharm.Res.*, **1991**, 8, 1434-1440.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50  $\mu$ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood)  $\mu$ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18

**Column:** 250  $\times$  4.6 5  $\mu$ m Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 208.7

---

**CHROMATOGRAM**

**Retention time:** 17.583

---

**KEY WORDS**whole blood

---

**REFERENCE**

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

---

**SAMPLE**

**Matrix:** microsomal incubations, tissue

**Sample preparation:** Tissue. Homogenize liver slices with 200  $\mu$ L water and about 250 mg 1 mm glass beads, using a Mini Bead Beater (Biospec Products), extract with 1 mL MeCN. Extract a 250  $\mu$ L aliquot of the slice incubation medium with 1 mL MeCN. Mix, centrifuge at 11000 g for 4 min at 4°. Dry the supernatant under vacuum, reconstitute the residue with 150  $\mu$ L mobile phase, inject an aliquot. Microsomal incubations. Add 1 mL MeCN to 1 mL microsomal incubation, mix, centrifuge at 1900 g for 15 min at 4°. Extract a 200  $\mu$ L aliquot of the supernatant with 1 mL MeCN, centrifuge at 11000 g for 4 min at 4°. Dry the supernatant under vacuum, reconstitute the residue with 150  $\mu$ L mobile phase, inject an aliquot.

---

**HPLC VARIABLES**

**Guard column:** Symmetry C18 Sentry Guard (Waters)

**Column:** 250  $\times$  4.6 5  $\mu$ m Symmetry C18 (Waters)

**Mobile phase:** MeCN:10 mM pH 4.0 ammonium phosphate triethylamine50:49.9:0.1

**Flow rate:** 1

**Detector:** UV 278

---

**CHROMATOGRAM**

**Retention time:** 31.6

---

**OTHER SUBSTANCES**

**Extracted:** metabolites

---

**KEY WORDS**liver

---

**REFERENCE**

Jamis-Dow,C.A.; Katki,A.G.; Collins,J.M.; Klecker,R.W. Rifampin and rifabutin and their metabolism by human liver esterases, *Xenobiotica*, **1997**, 27, 1015–1024.

---

**SAMPLE**

**Matrix:** tissue

**Sample preparation:** Vortex 200  $\mu$ L incubation mixture and 1 mL MeOH containing 0.1% acetic acid, immerse in ice for 15 min, centrifuge at 6500 g for 5 min, evaporate the supernatant to dryness under nitrogen, reconstitute the residue in 50  $\mu$ L mobile phase, combine two aliquots, inject a 75  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu$ m Hypersil C18

**Mobile phase:** Gradient. A was 3% triethylamine containing 0.3% trifluoroacetic acid adjusted to pH 2 with phosphoric acid. B was MeCN. A:B from 73:27 to 58:42 over 45 min

**Flow rate:** 1.5

**Injection volume:** 75

**Detector:** Radioactivity, Radiomatic A-500 (Flow Scintillation Analyzer)

---

**CHROMATOGRAM**

**Retention time:** 35

---

**OTHER SUBSTANCES**

**Extracted:** metabolites

**KEY WORDS**

liver; enterocytes; intestine

**REFERENCE**

Iatsumirskaja,E.; Tulebaev,S.; Storozhuk,E.; Utkin,I.; Smith,D.; Gerber,N.; Koudriakova,T. Metabolism of rifabutin in human enterocyte and liver microsomes: kinetic parameters, identification of enzyme systems, and drug interactions with macrolides and antifungal agents, *Clin.Pharmacol.Ther.*, **1997**, *61*, 554–562.

**SAMPLE**

**Matrix:** urine

**Sample preparation:** Inject a 50-200  $\mu\text{L}$  aliquot of urine directly.

**HPLC VARIABLES**

**Guard column:** 37-50  $\mu\text{m}$  Corasil C18

**Column:** 300  $\times$  3.9 10  $\mu\text{m}$   $\mu\text{Bondapak C18}$

**Mobile phase:** Gradient. MeCN:50 mM pH 4.5  $\text{KH}_2\text{PO}_4$  40:60 for 10 min, to 60:40 over 5 min, maintain at 60:40 for 5 min.

**Flow rate:** 1 for 10 min then 1.5

**Injection volume:** 50-200

**Detector:** UV 275 or radioactivity

**CHROMATOGRAM**

**Retention time:** 16.5

**OTHER SUBSTANCES**

**Extracted:** metabolites

**KEY WORDS**

rat

**REFERENCE**

Battaglia,R.; Salgarollo,G.; Zini,G.; Montesanti,L.; Strolin Benedetti,M. Absorption, disposition, and urinary metabolism of  $^{14}\text{C}$ -rifabutin in rats, *Antimicrob.Agents Chemother.*, **1991**, *35*, 1391–1396.


# Rifampin

**Molecular formula:**  $\text{C}_{43}\text{H}_{58}\text{N}_4\text{O}_{12}$

**Molecular weight:** 822.95

**CAS Registry No.:** 13292-46-1

**Merck Index:** 8382

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Add 400  $\mu\text{L}$  1 M pH 4.0 phosphate buffer and 40  $\mu\text{L}$  10 mg/mL ascorbic acid solution to 1 mL plasma. Mix for 10 s, add 4 mL ethyl acetate. Vortex for 5 min, centrifuge at 1800 g for 10 min. Remove a 3 mL aliquot of the aqueous layer and evaporate it to dryness under nitrogen at 100°. Reconstitute the residue with 500  $\mu\text{L}$  MeOH, vortex for 1 min. Inject a 20  $\mu\text{L}$  aliquot.

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu\text{m}$  Spherisorb

**Mobile phase:** MeOH:10 mM pH 5.5 phosphate buffer 68:32

**Column temperature:** 28

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 336

---

**CHROMATOGRAM**

**Retention time:** 4.33

**Internal standard:** rifampin

---

**OTHER SUBSTANCES**

**Extracted:** rifapentine

---

**KEY WORDS**

serum; rifampin is IS

---

**REFERENCE**

He,X.; Wang,J.; Liu,X.; Chen,X. High-performance liquid chromatography assay of rifapentine in human serum, *J.Chromatogr.B*, **1996**, 681, 412–415.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Add 100  $\mu$ L MeCN to 100  $\mu$ L plasma, vortex for 10 s, centrifuge for 10 min, inject a 20  $\mu$ L aliquot of the supernatant.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 4  $\mu$ m Nova-Pak C8

**Mobile phase:** MeCN:MeOH:10 mM potassium dihydrogen phosphate 35:5:60

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 334

---

**CHROMATOGRAM**

**Retention time:** 6.5

**Limit of quantitation:** 200 ng/mL

---

**KEY WORDS**

plasma; pharmacokinetics

---

**REFERENCE**

Le Guellec,C.; Gaudet,M.-L.; Lamanetre,S.; Breteau,M. Stability of rifampin in plasma: Consequences for therapeutic monitoring and pharmacokinetic studies, *Ther.Drug Monit.*, **1997**, 19, 669–674.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Condition a 100 mg Bond Elut C2 SPE cartridge with 1 mL MeOH and 1 mL 100 mM HCl just before use. 500  $\mu$ L Plasma + 500  $\mu$ L 100 mM HCl + 50  $\mu$ L 19  $\mu$ g/mL sulindac in MeOH, vortex briefly, add to SPE cartridge, wash with 1 mL 100 mM HCl, elute with 400  $\mu$ L MeOH:MeCN 60:40, mix, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Guard column:**  $\mu$ Bondapak C18 RCSS Guard Pak

**Column:** 100  $\times$  8 4  $\mu$ m Nova Pak C18 Radial Pak

**Mobile phase:** MeCN:buffer 42:58 (Buffer was 50 mM sodium citrate adjusted to pH 4.3 with 50 mM HCl.)

**Flow rate:** 2.3

**Injection volume:** 20

**Detector:** UV 342

---

**CHROMATOGRAM**

**Retention time:** 4.65

**Internal standard:** sulindac (3.07)

**Limit of quantitation:** 160 ng/mL

---

**OTHER SUBSTANCES**

Extracted: metabolites

---

**KEY WORDS**

plasma; pharmacokinetics; SPE

---

**REFERENCE**

Swart,K.J.; Paggis,M. Automated high-performance liquid chromatographic method for the determination of rifampicin in plasma, *J.Chromatogr.*, **1992**, 593, 21–24.

---

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 2 mL Plasma + 50  $\mu$ L 50  $\mu$ g/mL N-butarylaminophenol + 400  $\mu$ L 10% acetic acid + 7 mL diethyl ether:dichloromethane 2:1, shake, centrifuge at 2059 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 200  $\mu$ L MeOH, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Guard column:** 50  $\times$  4.6 30  $\mu$ m C8

**Column:** 250  $\times$  4.6 5  $\mu$ m Spherisorb C8

**Mobile phase:** Gradient. MeCN:5 mM pH 3.5 phosphate buffer from 6:94 to 90:10 over 5 min, maintain at 90:10 for 12 min.

**Flow rate:** 2

**Injection volume:** 20

**Detector:** UV 248

---

**CHROMATOGRAM**

**Retention time:** 11.9

**Internal standard:** N-butarylaminophenol (4.22)

**Limit of detection:** 200 ng/mL

---

**OTHER SUBSTANCES**

Extracted: pyrazinamide

---

**KEY WORDS**

plasma

---

**REFERENCE**

Walubo,A.; Smith,P.; Folb,P.I. Comprehensive assay for pyrazinamide, rifampicin and isoniazid with its hydrazine metabolites in human plasma by column liquid chromatography, *J.Chromatogr.B*, **1994**, 658, 391–396.

---

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Condition a 3 mL 15 mg SPEC.C8 SPE disc (ANSYS) or a Bond-Elut C8 SPE cartridge with two 1 mL portions of MeCN and two 1 mL portions of 100 mM HCl. 500  $\mu$ L Serum + 200  $\mu$ L 5% ascorbic acid in 100 mM HCl + 100  $\mu$ L 50  $\mu$ g/mL papaverine hydrochloride in water, mix well, add to the SPE disc or cartridge, apply a vacuum for 20 min, elute with five 100  $\mu$ L portions of MeCN:MeOH 60:40. Evaporate the eluate to dryness under a stream of nitrogen at room temperature, reconstitute the residue in 300  $\mu$ L 100 mM HCl:water 67:33, inject a 50  $\mu$ L aliquot. (Disc gave higher recoveries than cartridge.)

---

**HPLC VARIABLES**

**Guard column:** ODS

**Column:** 250  $\times$  4.6 Spherisorb ODS-1

**Mobile phase:** MeCN:100 mM pH 4.7 KH<sub>2</sub>PO<sub>4</sub> 35:65

**Column temperature:** 40

**Flow rate:** 1.2

**Injection volume:** 50

**Detector:** UV 340

---

---

**CHROMATOGRAM****Retention time:** 11.5**Internal standard:** papaverine (14.9)

---

**OTHER SUBSTANCES****Extracted:** metabolites

---

**KEY WORDS**

serum; SPE

---

**REFERENCE**

Ye,L.; Stewart,J.T.; Zhang,H. A comparison of disc and cartridge solid-phase extraction for the LC determination of rifampin and 25-desacetyl rifampin in human serum, *J.Pharm.Biomed.Anal.*, **1995**, 13, 1185–1188.

---

**SAMPLE****Matrix:** blood

**Sample preparation:** Condition a 1 mL Bond Elut C2 SPE cartridge with 1 mL MeOH and 1 mL water. 500  $\mu$ L Plasma + 25  $\mu$ L 2  $\mu$ g/mL sulindac in water + 100  $\mu$ L 100 mM HCl, vortex briefly, centrifuge at 630 rpm for 5 min, add to the SPE cartridge, wash with 1 mL water, elute with 1 mL MeOH, add the eluate to 500  $\mu$ L 3 mg/mL ascorbic acid, inject a 150  $\mu$ L aliquot.

---

**HPLC VARIABLES****Guard column:** 15  $\times$  3.2 7  $\mu$ m RP-8 (Brownlee, ABI)**Column:** 250  $\times$  4.6 5  $\mu$ m Zorbax RX C8**Mobile phase:** MeCN:50 mM  $\text{KH}_2\text{PO}_4$  45:55**Flow rate:** 1**Injection volume:** 150**Detector:** UV 340

---

**CHROMATOGRAM****Retention time:** 4.4**Internal standard:** sulindac (7.8)**Limit of quantitation:** 50 ng/mL

---

**OTHER SUBSTANCES****Simultaneous:** atevirdine, delavirdine, U-89,255, U-96,183

---

**KEY WORDS**

plasma; protect from light; pharmacokinetics; rugged; SPE

---

**REFERENCE**

Lau,Y.Y.; Hanson,G.D.; Carel,B.J. Determination of rifampin in human plasma by high-performance liquid chromatography with ultraviolet detection, *J.Chromatogr.B*, **1996**, 676, 147–152.

---

**SAMPLE****Matrix:** blood, CSF

**Sample preparation:** 200  $\mu$ L CSF or plasma + dimethylaminobenzoic acid + 200  $\mu$ L MeOH, mix, add 1 mL buffer, add 7 mL dichloromethane:diethyl ether 40:60, shake for 15 min, centrifuge at 2500 g for 10 min. Remove the organic layer and evaporate it to dryness AT 45°, reconstitute the residue in 50  $\mu$ L MeOH, vortex for 10 s, inject a 25  $\mu$ L aliquot. (Buffer was 1 M  $\text{KH}_2\text{PO}_4$  containing 0.2% ascorbic acid, pH adjusted to 4.2 with 1 M HCl.)

---

**HPLC VARIABLES****Guard column:** 50  $\times$  4.6 30  $\mu$ m C8**Column:** 250  $\times$  4.6 5  $\mu$ m Lichrosorb RP-8**Mobile phase:** MeCN:10 mM pH 3.5 phosphate buffer 48:52**Flow rate:** 1.5**Injection volume:** 25**Detector:** UV 215



---

**CHROMATOGRAM****Retention time:** 5.8**Internal standard:** dimethylaminobenzoic acid (3.8)**Limit of detection:** 500 ng/mL

---

**OTHER SUBSTANCES****Noninterfering:** p-aminosalicylic acid, pyrazinamide, isoniazid

---

**KEY WORDS**plasma; rabbit

---

**REFERENCE**Chan, K. Rifampicin concentrations in cerebrospinal fluid and plasma of the rabbit by high performance liquid chromatography, *Methods Find. Exp. Clin. Pharmacol.*, **1986**, 8, 721-726.

---

**SAMPLE****Matrix:** blood, CSF**Sample preparation:** 200  $\mu$ L Serum, plasma, or CSF + 300  $\mu$ L reagent. Flush column A to waste with 500  $\mu$ L 500 mM ammonium sulfate, inject sample onto column A, flush column A to waste with 500  $\mu$ L 500 mM ammonium sulfate, elute the contents of column A onto column B with mobile phase, monitor the effluent from column B. (Reagent was 8.05 M guanidine hydrochloride and 1.02 M ammonium sulfate in water.)

---

**HPLC VARIABLES****Column:** A 30  $\times$  2.1 40  $\mu$ m preparative grade C18 (Analytichem); B 250  $\times$  4.6 10  $\mu$ m Partisil C8**Mobile phase:** Gradient. A was 50 mM pH 4.5  $\text{KH}_2\text{PO}_4$ . B was MeCN:isopropanol 80:20. A:B 90:10 for 1 min, to 30:70 over 15 min, maintain at 30:70 for 4 min.**Column temperature:** 50**Flow rate:** 1.5**Detector:** UV 280 for 5 min then UV 254

---

**CHROMATOGRAM****Retention time:** 12.93, 16.28, 17.16 (compound undergoes decomposition)**Internal standard:** heptanophenone (19.2)

---

**OTHER SUBSTANCES****Extracted:** acetazolamide, ampicillin, bromazepam, caffeine, carbamazepine, chloramphenicol, chlorothiazide, diazepam, droperidol, ethionamide, furosemide, isoniazid, methadone, penicillin G, phenobarbital, phenytoin, prazepam, propoxyphene, pyrazinamide, trimeprazine, trimethoprim

---

**KEY WORDS**plasma; serum; column-switching

---

**REFERENCE**Seifart, H.I.; Kruger, P.B.; Parkin, D.P.; van Jaarsveld, P.P.; Donald, P.R. Therapeutic monitoring of antituberculosis drugs by direct in-line extraction on a high-performance liquid chromatography system, *J. Chromatogr.*, **1993**, 619, 285-290.

---

**SAMPLE****Matrix:** blood, saliva, urine**Sample preparation:** 1 mL Plasma, urine, or saliva + 50  $\mu$ L MeOH + 1 mL buffer + 1 mL isooctane:dichloromethane 60:40, shake mechanically at 350 rpm for 10 min, centrifuge at 2000 g for 5 min, inject a 5-200  $\mu$ L aliquot of the supernatant. (Buffer was 2 g ascorbic acid and 10 g anhydrous sodium sulfate made up to 50 mL with concentrated pH 6 buffer solution (Merck Cat. No. 9886), stir continuously, prepare fresh each day.)

---

**HPLC VARIABLES****Column:** 100  $\times$  7.5 5  $\mu$ m LiChrosorb Si 60**Mobile phase:** Isooctane:dichloromethane:EtOH:water:acetic acid 45:36.6:16.8:1.65:0.002

---

**Flow rate:** 3  
**Injection volume:** 5-200  
**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 3  
**Limit of detection:** 20 ng/mL (Clin.Pharmacol.Ther. 1991, 50, 682)  
**Limit of quantitation:** 100 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** metabolites

---

**KEY WORDS**

plasma; pharmacokinetics

---

**REFERENCE**

Lecaillon, J.B.; Febvre, N.; Metayer, J.P.; Souppart, C. Quantitative assay of rifampicin and three of its metabolites in human plasma, urine and saliva by high-performance liquid chromatography, *J.Chromatogr.*, **1978**, *145*, 319-324.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50  $\mu$ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood)  $\mu$ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18

**Column:** 250  $\times$  4.6 5  $\mu$ m Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 236.9

---

**CHROMATOGRAM**

**Retention time:** 16.167

---

**KEY WORDS**

whole blood

---

**REFERENCE**

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, *763*, 149-163.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50  $\mu$ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood)  $\mu$ L aliquot. (The detector wavelength shown is the wavelength of

maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18

**Column:** 250 × 4.6 5 µm Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 225.2

---

**CHROMATOGRAM**

**Retention time:** 20.85

---

**KEY WORDS**

whole blood

---

**REFERENCE**

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

---

**SAMPLE**

**Matrix:** cells

**Sample preparation:** 100 µL Cell suspension + 100 µL cefoperazone solution + 100 µL Hanks balanced salt solution, sonicate 30 min, add 800 µL MeCN, centrifuge at 13000 g for 5 min, remove supernatant. Dry supernatant under air, dissolve in 100 µL mobile phase, inject 75 µL.

---

**HPLC VARIABLES**

**Column:** µBondapak C18

**Mobile phase:** MeCN:50 mM pH 4.7 KH<sub>2</sub>PO<sub>4</sub> 40:60

**Flow rate:** 1

**Injection volume:** 75

**Detector:** UV 340

---

**CHROMATOGRAM**

**Retention time:** 6.2

**Internal standard:** papaverine

**Limit of detection:** 100-1000 ng/mL

---

**REFERENCE**

Darouiche, R.O.; Hamill, R.J. Antibiotic penetration of and bactericidal activity within endothelial cells, *Antimicrob. Agents Chemother.*, **1994**, 38, 1059-1064.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 5 µm Vydac 208TP54 C8

**Mobile phase:** MeCN:buffer 32.8:67.2 (Buffer was 50 mM K<sub>2</sub>HPO<sub>4</sub> adjusted to pH 6.5 with phosphoric acid.)

**Flow rate:** 1

**Detector:** UV 254

---

**CHROMATOGRAM****Retention time:** 6.8

---

**OTHER SUBSTANCES****Simultaneous:** minocycline, degradation products

---

**KEY WORDS**injections; saline; 5% dextrose; stability-indicating

---

**REFERENCE**Pearson, S.D.; Trissel, L.A. Stability and compatibility of minocycline hydrochloride and rifampin in intravenous solutions at various temperatures, *Am.J.Hosp.Pharm.*, **1993**, 50, 698–702.

---

**SAMPLE****Matrix:** microsomal incubations, tissue**Sample preparation:** Tissue. Homogenize liver slices with 200  $\mu$ L water and about 250 mg 1 mm glass beads, using a Mini Bead Beater (Biospec Products), extract with 1 mL MeCN. Extract a 250  $\mu$ L aliquot of the slice incubation medium with 1 mL MeCN. Mix, centrifuge at 11000 g for 4 min at 4°. Dry the supernatant under vacuum, reconstitute the residue with 150  $\mu$ L mobile phase, inject an aliquot. Microsomal incubations. Add 1 mL MeCN to 1 mL microsomal incubation, mix, centrifuge at 1900 g for 15 min at 4°. Extract a 200  $\mu$ L aliquot of the supernatant with 1 mL MeCN, centrifuge at 11000 g for 4 min at 4°. Dry the supernatant under vacuum, reconstitute the residue with 150  $\mu$ L mobile phase, inject an aliquot.

---

**HPLC VARIABLES****Guard column:** Symmetry C18 Sentry Guard (Waters)**Column:** 250  $\times$  4.6 5  $\mu$ m Symmetry C18 (Waters)**Mobile phase:** MeCN:10 mM pH 4.0 ammonium phosphate:triethylamine 30:69.9:0.1**Flow rate:** 1**Detector:** UV 240

---

**CHROMATOGRAM****Retention time:** 34.4

---

**OTHER SUBSTANCES****Extracted:** metabolites

---

**KEY WORDS**liver

---

**REFERENCE**Jamis-Dow, C.A.; Katki, A.G.; Collins, J.M.; Klecker, R.W. Rifampin and rifabutin and their metabolism by human liver esterases, *Xenobiotica*, **1997**, 27, 1015–1024.

---

**SAMPLE****Matrix:** solutions**Sample preparation:** Centrifuge and filter cell solutions (0.22  $\mu$ m), inject an aliquot.

---

**HPLC VARIABLES****Guard column:** Guard-PAK C18 (Waters)**Column:** 150  $\times$  3.9 5  $\mu$ m NOVA PAK C18**Mobile phase:** MeCN:50 mM pH 6.0  $\text{KH}_2\text{PO}_4$  40:60**Flow rate:** 1**Detector:** UV 254

---

**CHROMATOGRAM****Retention time:** 3.3

---

**REFERENCE**Koga, H. High-performance liquid chromatography measurement of antimicrobial concentrations in polymorphonuclear leukocytes, *Antimicrob.Agents Chemother.*, **1987**, 31, 1904–1908.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a 5 mg/mL solution in MeOH, inject a 10  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Guard column:** present but not specified

**Column:** 100  $\times$  4.6 Hypersil C18

**Mobile phase:** MeOH

**Flow rate:** 1

**Injection volume:** 10

**Detector:** MS, HP 5985B quadrupole MS, direct liquid introduction interface, 5  $\mu$ m orifice, water-cooled stainless steel diaphragm, ion source pressure 1 Torr, ion source 250°, LC eluent used as reagent gas

---

**OTHER SUBSTANCES**

**Also analyzed:** rifapentine, rifamycin SV

---

**REFERENCE**

Vékey, K.; Edwards, D.M.F.; Zerilli, L.F. Liquid chromatographic-mass spectrometric studies on rifamycin antibiotics, *J. Chromatogr.*, **1989**, 474, 317–327.

---

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu$ m Excalibar C18-CN (Alltech)

**Mobile phase:** MeOH:5 mM tetra-n-butylammonium hydroxide 80:20 adjusted to pH 3.0 with phosphoric acid

**Flow rate:** 1.5

**Injection volume:** 20

**Detector:** UV 265

---

**CHROMATOGRAM**

**Retention time:** 4.3

**Internal standard:** rifampin

---

**OTHER SUBSTANCES**

**Simultaneous:** isoniazid, pyrazinamide

---

**KEY WORDS**

rifampin is IS

---

**REFERENCE**

Gaitonde, C.D.; Pathak, P.V. Rapid liquid chromatographic method for the estimation of isoniazid and pyrazinamide in plasma and urine, *J. Chromatogr.*, **1990**, 532, 418–423.

---

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a 50  $\mu$ g/mL solution in MeOH:water 50:50 (electrospray) or a 400  $\mu$ g/mL solution in MeOH:200 mM ammonium acetate 50:50 (thermospray), inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 150  $\times$  4.6 5  $\mu$ m Spherisorb ODS-2

**Mobile phase:** MeOH:200 mM ammonium acetate 60:40

**Flow rate:** 1.2

**Injection volume:** 20

**Detector:** MS, Delsi/Nermag R3010 triple quadrupole, Delsi/Nermag electrospray interface, flow rate 1  $\mu$ L/min into source or Finnigan MAT TSQ 70 triple quadrupole, Finnigan TSPI interface and source, positive ion, discharge off, filament off, vaporizer 75°, block (jet) 330

**CHROMATOGRAM****Retention time:** 13.5**OTHER SUBSTANCES****Simultaneous:** rifamycin SV, rifamycin B**KEY WORDS**

electrospray; thermospray


**REFERENCE**

Korfmacher, W.A.; Bloom, J.; Churchwell, M.I.; Getek, T.A.; Hansen, E.B., Jr.; Holder, C.L.; McManus, K.T. Characterization of three rifamycins via electrospray mass spectrometry and HPLC-thermospray mass spectrometry, *J. Chromatogr. Sci.*, **1993**, *31*, 498–501.

**SAMPLE****Matrix:** solutions**Sample preparation:** Prepare a solution in mobile phase, inject a 10  $\mu$ L aliquot.**HPLC VARIABLES****Column:** 160  $\times$  4.5  $\mu$ m Zorbax Rx-C8**Mobile phase:** MeOH:buffer 70:30 (Buffer was 5 mM tetrabutylammonium hydroxide adjusted to pH 3.0 with phosphoric acid.)**Flow rate:** 0.7**Injection volume:** 10**Detector:** UV 265**CHROMATOGRAM****Retention time:** 3.8**OTHER SUBSTANCES****Simultaneous:** pyrazinamide**REFERENCE**

Nahata, M.C.; Morosco, R.S.; Peritore, S.P. Stability of pyrazinamide in two suspensions, *Am. J. Health-Syst. Pharm.*, **1995**, *52*, 1558–1560.

# Rifapentine

**Molecular formula:** C<sub>47</sub>H<sub>64</sub>N<sub>4</sub>O<sub>12</sub>**Molecular weight:** 877.04**CAS Registry No.:** 61379-65-5**Merck Index:** 8385**SAMPLE****Matrix:** blood

**Sample preparation:** Add 400  $\mu$ L 1 M pH 4.0 phosphate buffer, 40  $\mu$ L 10 mg/mL ascorbic acid solution, and 50  $\mu$ L 100  $\mu$ g/mL rifampin solution to 1 mL plasma. Mix for 10 s, add 4 mL ethyl acetate. Vortex for 5 min, centrifuge at 1800 g for 10 min. Remove a 3 mL aliquot of the aqueous layer and evaporate it to dryness under nitrogen at 100°. Reconstitute the residue with 500  $\mu$ L MeOH, vortex for 1 min. Inject a 20  $\mu$ L aliquot.

**HPLC VARIABLES****Column:** 250  $\times$  4.6 5  $\mu$ m Spherisorb C18

**Mobile phase:** MeOH:10 mM pH 5.5 phosphate buffer 68:32  
**Column temperature:** 28  
**Flow rate:** 1  
**Injection volume:** 20  
**Detector:** UV 336

---

#### CHROMATOGRAM

**Retention time:** 6.83  
**Internal standard:** rifampin (4.33)  
**Limit of quantitation:** 5 ng

---

#### KEY WORDS

serum

---

#### REFERENCE

He,X.; Wang,J.; Liu,X.; Chen,X. High-performance liquid chromatography assay of rifapentine in human serum, *J.Chromatogr.B*, **1996**, 681, 412–415.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** Filter (0.45  $\mu$ m) plasma, inject a 50  $\mu$ L aliquot.

---

#### HPLC VARIABLES

**Guard column:** 20  $\times$  4.6 Supelcosil LC HISEP

**Column:** 150  $\times$  4.6 Supelcosil LC HISEP (Invert column between runs to prevent clogging.)

**Mobile phase:** Gradient. A was MeCN:2.6 mM 2-(N-morpholino)ethanesulfonic acid 5:95 containing 1% trifluoroacetic acid, adjusted to pH 6.5 with ethanolamine. B was MeCN:THF:2.6 mM 2-(N-morpholino)ethanesulfonic acid 20:10:70 containing 1% trifluoroacetic acid, adjusted to pH 6.5 with ethanolamine. A:B 70:30 for 4 min, to 30:70 over 24 min, re-equilibrate for 6 min.

**Flow rate:** 1.4

**Injection volume:** 50

**Detector:** UV 254

---

#### CHROMATOGRAM

**Retention time:** 21

**Limit of detection:** 200 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** metabolites

---

#### KEY WORDS

plasma; direct injection

---

#### REFERENCE

Riva,E.; Merati,R.; Cavenaghi,L. High-performance liquid chromatographic determination of rifapentine and its metabolite in human plasma by direct injection into a shielded hydrophobic phase column, *J.Chromatogr.*, **1991**, 553, 35–40.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** 100  $\mu$ L Serum + 300  $\mu$ L 10  $\mu$ g/mL rifampin in 1% ascorbic acid solution, mix, keep at 4°, inject a 100  $\mu$ L aliquot on column A with mobile phase A and elute to waste, after 6 min backflush the contents of column A onto column B with mobile phase B, after 4 min remove column A from the circuit, elute column B with mobile phase B and monitor the effluent from column B. Re-equilibrate column A with mobile phase A for 15 min.

---

#### HPLC VARIABLES

**Column:** A 40  $\times$  20 37-50  $\mu$ m Corasil RP C18; B 20  $\times$  4.6 25-40  $\mu$ m LiChrosorb RP-8 + 300  $\times$  3.9 10  $\mu$ m  $\mu$ Bondapak C18

**Mobile phase:** A 50 mM pH 7.0 phosphate buffer; B MeCN:THF:50 mM pH 7.0 phosphate buffer 42:5:53

**Flow rate:** 1

**Injection volume:** 100

**Detector:** UV 332

## CHROMATOGRAM

**Retention time:** 12.6

**Internal standard:** rifampin (7.2)

**Limit of detection:** 100 ng/mL

## OTHER SUBSTANCES

**Noninterfering:** acetaminophen, p-aminosalicylic acid, aspirin, caffeine, cefuroxime, ciprofloxacin, ibuprofen, isonicotinic acid, theophylline, vitamin A, thiamine, riboflavin, pyridoxine, ascorbic acid

## KEY WORDS

serum; column-switching; dog; pharmacokinetics

## REFERENCE

Lee,H.S.; Shin,H.C.; Han,S.S.; Roh,J.K. High-performance liquid chromatographic determination of rifapentine in serum using column switching, *J.Chromatogr.*, **1992**, 574, 175–178.

## SAMPLE

**Matrix:** solutions

**Sample preparation:** Prepare a 5 mg/mL solution in MeOH, inject a 10  $\mu$ L aliquot.

## HPLC VARIABLES

**Guard column:** present but not specified

**Column:** 100  $\times$  4.6 Hypersil C18

**Mobile phase:** MeOH

**Flow rate:** 1

**Injection volume:** 10

**Detector:** MS, HP 5985B quadrupole MS, direct liquid introduction interface, 5  $\mu$ m orifice, water-cooled stainless steel diaphragm, ion source pressure 1 Torr, ion source 250°, LC eluent used as reagent gas

## OTHER SUBSTANCES

**Also analyzed:** rifampin, rifamycin SV

## REFERENCE

Vékey,K.; Edwards,D.M.F.; Zerilli,L.F. Liquid chromatographic-mass spectrometric studies on rifamycin antibiotics, *J.Chromatogr.*, **1989**, 474, 317–327.


# Rifaximin

**Molecular formula:** C<sub>43</sub>H<sub>51</sub>N<sub>3</sub>O<sub>11</sub>

**Molecular weight:** 785.89

**CAS Registry No.:** 80621-81-4

**Merck Index:** 8386



## SAMPLE

**Matrix:** bile



**Sample preparation:** 1 mL Bile + 4 mL pH 4.5 buffer + 200  $\mu$ L 50  $\mu$ g/mL IS in MeOH, extract twice with chloroform. Combine the organic layers and dry them over anhydrous sodium sulfate, evaporate to dryness, reconstitute the residue in 100  $\mu$ L MeOH, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** HSC18

**Mobile phase:** MeOH:1 mM pH 3 phosphate buffer 65:35

**Column temperature:** 35

**Flow rate:** 1

**Detector:** UV 276

---

**CHROMATOGRAM**

**Internal standard:** M 302

---

**REFERENCE**

Verardi,S.; Verardi,V. Bile rifaximin concentration after oral administration in patients undergoing cholecystectomy, *Farmaco*, **1990**, *45*, 131–135.

---


# Rilmazafone

**Molecular formula:**  $C_{21}H_{20}Cl_2N_6O_3$

**Molecular weight:** 475.33

**CAS Registry No.:** 99593-25-6

**Merck Index:** 8387



---

**SAMPLE**

**Matrix:** feed

**Sample preparation:** 1 g Feed + 10 mL MeOH, shake for 45 min, centrifuge at 2000 rpm for 10 min. Remove 5 mL of the supernatant and evaporate it to dryness under a stream of nitrogen at 30°, reconstitute the residue in 1 mL MeCN, add 100  $\mu$ L 1 M tartaric acid containing 50 mM sodium 1-dodecanesulfonate, vortex for 10 min, add 5 mL n-hexane, vortex for 10 min, filter (0.45  $\mu$ m) the aqueous layer, inject a 10  $\mu$ L aliquot of the filtrate.

---

**HPLC VARIABLES**

**Column:** two 300  $\times$  4 10  $\mu$ m Nucleosil 10C18 columns in series

**Mobile phase:** MeCN:EtOH:buffer 35:12.5:52.5 (Buffer was 100 mM tartaric acid containing 5 mM sodium 1-dodecanesulfonate and 0.29 mM triethylamine.)

**Flow rate:** 1.2

**Injection volume:** 10

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 30

---

**OTHER SUBSTANCES**

**Simultaneous:** degradation products

---

**REFERENCE**

Ikenishi,R.; Kitagawa,T. Analytical studies on 1-(2-o-chlorobenzoyl-4-chlorophenyl)-3-dimethylcarbamoyl-5-glycyl-aminomethyl-1H-1,2,4-triazole hydrochloride dihydrate. III. High-performance liquid chromatographic method applicable to animal feed, *Chem.Pharm.Bull.*, **1987**, *35*, 4544–4551.