


Rilmenidine



Molecular formula: $C_{10}H_{16}N_2O$

Molecular weight: 180.25

CAS Registry No.: 54187-04-1

Merck Index: 8388

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 9.8


KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

Riluzole



Molecular formula: $C_8H_5F_3N_2OS$

Molecular weight: 234.20

CAS Registry No.: 1744-22-5

Merck Index: 8389

SAMPLE

Matrix: blood

Sample preparation: Add 100-200 ng IS to 1 mL plasma adjusted to pH 6.0. Vortex, centrifuge at 2000 g for 5 min. Add to a 1 mL C2 Bond Elut SPE cartridge or CBA Bond Elut ion exchange SPE cartridge. Wash with 1 mL 10 mM dipotassium hydrogen phosphate, 500 μ L water, and MeOH:water 10:90. Elute with 500 μ L mobile phase, dry under a stream of nitrogen. Inject a 300 μ L aliquot.

HPLC VARIABLES

Guard column: 39 × 4.6 μBondapak C18

Column: 300 × 4.6 μBondapak C18

Mobile phase: MeCN:MeOH:acetic acid:dipotassium hydrogen phosphate 10:55:1:35

Flow rate: 1

Injection volume: 300

Detector: UV 265

CHROMATOGRAM

Retention time: 9.2

Internal standard: RP 61307 (riluzole N-methyl derivative) (12.4)

Limit of detection: 5 ng/mL

KEY WORDS

plasma; pharmacokinetics; SPE

REFERENCE

Le Liboux,A.; Lefebvre,P.; Le Roux,Y.; Truffinet,P.; Aubeneau,M.; Kirkesseli,S.; Montay,G. Single- and multiple-dose pharmacokinetics of riluzole in white subjects, *J.Clin.Pharmacol.*, **1997**, 37, 820-827.

SAMPLE

Matrix: urine

Sample preparation: Add 20 µg IS to 1 mL urine adjusted to pH 6.0. Vortex, centrifuge at 2000 g for 5 min. Add to a 1 mL C2 Bond Elut SPE cartridge or CBA Bond Elut ion exchange cartridge. Wash with 1 mL 10 mM dipotassium hydrogen phosphate, 500 µL water, MeOH: water 10:90. Elute with 500 µL mobile phase, inject a 300 µL aliquot.

HPLC VARIABLES

Guard column: 39 × 4.6 μBondapak C18

Column: 300 × 4.6 μBondapak C18

Mobile phase: MeCN:MeOH:acetic acid:dipotassium hydrogen phosphate 10:55:1:35

Flow rate: 0.8-1

Injection volume: 300

Detector: UV 265

CHROMATOGRAM

Retention time: 9.2

Internal standard: RP 61 307 (riluzole N-methyl derivative) (12.4)

Limit of detection: 10 ng/mL

KEY WORDS

pharmacokinetics; SPE

REFERENCE

Le Liboux,A.; Lefebvre,P.; Le Roux,Y.; Truffinet,P.; Aubeneau,M.; Kirkesseli,S.; Montay,G. Single- and multiple-dose pharmacokinetics of riluzole in white subjects, *J.Clin.Pharmacol.*, **1997**, 37, 820–827.

Rimantadine


Molecular formula: C₁₂H₂₁N

Molecular weight: 179.31

CAS Registry No.: 13392-28-4, 1501-84-4 (HCl)

Merck Index: 8390

Lednicher No.: 2 19



SAMPLE

Matrix: solutions

Sample preparation: 50 μ L 5 mg/mL Rimantadine in 100 mM HCl + 50 μ L buffer + 100 μ L reagent, swirl for 1 min, place on ice for 5 min, add 2 mL mobile phase, inject a 5 μ L aliquot. (Buffer was 100 mM sodium borate adjusted to pH 9.50 with 2 M NaOH. Reagent was 13.40 g o-phthalaldehyde and 36.4 mg 2,3,4,6-tetra-O-acetyl-1-thio- β -D-glucopyranoside in 1 mL MeOH, protect from light, keep on ice.)

HPLC VARIABLES

Column: 150 \times 3.9 4 μ m Nova-Pak C18

Mobile phase: MeOH:buffer 85:15 (Buffer was 3 mL/L glacial acetic acid in water, pH adjusted to 7.20 with 2 M NaOH.)

Flow rate: 1

Injection volume: 5

Detector: F ex 338 em 420 or UV 254

CHROMATOGRAM

Retention time: 5.88, 7.11 (enantiomers)

Limit of detection: 6 ng (UV)

KEY WORDS

derivatization; protect from light; chiral

REFERENCE

Desai,D.M.; Gal,J. Enantiospecific drug analysis via the *ortho*-phthalaldehyde/homochiral thiol derivatization method, *J.Chromatogr.*, **1993**, 629, 215–228.

SAMPLE

Matrix: urine

Sample preparation: Adjust pH of 8 mL urine to 5.5 with 2 M HCl, add 80 μ L Glusulase (DuPont), heat at 37° for 18 h. Remove a 2 mL aliquot and add it to 500 μ L 5 M NaOH, add 8 mL cyclohexane saturated with triethanolamine:chloroform 2:1, add 100 μ L 2% pentafluorobenzoyl chloride in cyclohexane, shake at 30 strokes/min for 20 min, centrifuge at 1500 g. Remove 7.5 mL of the upper organic layer and evaporate it to dryness under a stream of nitrogen at 50°, reconstitute the residue in 160 μ L heptane:isopropanol 95:5, inject an aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 Zorbax silica

Mobile phase: Gradient. Hexane:isopropanol 95:5 for 10 min, to 90:10 over 20 min, re-equilibrate for 5 min.

Flow rate: 1.5

Detector: radioactivity

CHROMATOGRAM

Retention time: 6.4

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

14C-labeled; derivatization; normal phase

REFERENCE

Rubio,F.R.; Fukuda,E.K.; Garland,W.A. Urinary metabolites of rimantadine in humans, *Drug Metab.Dispos.*, **1988**, 16, 773–777.

Risperidone


Molecular formula: C₂₃H₂₇FN₄O₂

Molecular weight: 410.49

CAS Registry No.: 106266-06-2

Merck Index: 8397

Lednicer No.: 5 150



SAMPLE

Matrix: blood

Sample preparation: 1 mL Serum + 500 μ L 600 mM pH 10 sodium carbonate/bicarbonate buffer + 50 μ L 3.76 mg/mL IS in MeOH + 8 mL heptane:isoamyl alcohol 98:2, shake at 250 cycles/min for 5 min, centrifuge at 1500 g for 10 min. Freeze the aqueous layer, evaporate the heptane layer to dryness under a gentle stream of nitrogen at 60°. Dissolve the residue in 75 μ L mobile phase, inject a 65 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 LiChroCart

Mobile phase: MeOH:40 mM pH 7.0 ammonium acetate buffer 90:10

Flow rate: 1

Injection volume: 65

Detector: UV 280

CHROMATOGRAM

Retention time: 4.42

Internal standard: haloperidol (5.10)

Limit of quantitation: 1.2 ng/mL

OTHER SUBSTANCES

Extracted: metabolites IN clozapine, diltiazem, fluphenazine, hydroxyzine, mianserine, perphenazine, zuclopenthixol

Simultaneous: amitriptyline, citalopram, chlorprothixene, clomipramine, desipramine, desmethylcitalopram, desmethylclomipramine, desmethylsertraline, fluoxetine, 10-hydroxyamitriptyline, 8-hydroxyclopmipramine, 8-hydroxydesmethylclomipramine, 10-hydroxynortriptyline, imipramine, methotrimeprazine sulfoxide, norfluoxetine, nortriptyline, paroxetine, sertraline

Noninterfering: carbamazepine, clonazepam, flunitrazepam, nitrazepam, oxazepam, oxcarbazepine

KEY WORDS

serum

REFERENCE

Olesen, O.V.; Linnet, K. Simplified high-performance liquid chromatographic method for determination of risperidone and 9-hydroxyrisperidone in serum from patients comedicated with other psychotropic drugs, *J.Chromatogr.B*, **1997**, 698, 209–216.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 100 μ L 25 ng/mL remoxipride in MeCN, mix for 5 s, add 1 mL saturated sodium carbonate, mix for 5 s, add 7 mL pentane:dichloromethane 75:25, shake gently for 10 min, centrifuge at 18° at 1735 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 60°, reconstitute the residue in 140 μ L MeCN, inject an aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Ultrasphere cyano

Mobile phase: MeCN:MeOH:40 mM pH 6.8 ammonium acetate 82:8:10

Column temperature: 40

Flow rate: 1.5

Injection volume: 150

Detector: E, ESA Coulochem model 5100A, model 5011 analytical cell, screening electrode 0.6 V, detection electrode 0.92 V, model 5020 guard cell 1 V

CHROMATOGRAM

Retention time: 14

Internal standard: remoxipride (17)

Limit of detection: 0.1 ng/mL

OTHER SUBSTANCES

Simultaneous: pseudoephedrine

Noninterfering: metabolites, acetaminophen, benztropine, clonazepam, clozapine, fluphenazine, haloperidol, ibuprofen, lorazepam, trihexyphenidyl

KEY WORDS

plasma

REFERENCE

Aravagiri, M.; Marder, S.R.; Van Putten, T.; Midha, K.K. Determination of risperidone in plasma by high-performance liquid chromatography with electrochemical detection: application to therapeutic drug monitoring in schizophrenic patients, *J. Pharm. Sci.*, **1993**, *82*, 447-449.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 100 μ L 200 ng/mL IS in MeOH + 1 mL 50 mM pH 10 borate buffer, vortex briefly, add to an Extrelut 3 SPE cartridge, let stand for 5 min, elute with 15 mL hexane:dichloromethane 50:50. Add the eluate to 3 mL 50 mM sulfuric acid, mix for 10 min, centrifuge at 3000 g for 10 min. Remove the aqueous layer and add it to 6 mL hexane:dichloromethane 50:50, wash for 5 min, centrifuge. Make the aqueous layer basic with 150 μ L 28% ammonia, extract twice with 3 mL hexane:dichloromethane 50:50. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 60°, reconstitute the residue in 100 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: 30 \times 4.6 5 μ m Spherisorb cyano

Column: 250 \times 4.6 5 μ m Ultrasphere cyano

Mobile phase: MeCN:buffer 60:40 (Buffer was 50 mM KH_2PO_4 adjusted to pH 6.5 with 28% ammonia.)

Flow rate: 1

Injection volume: 20

Detector: E, 5100 A Coulochem, 5020 guard cell 1.00 V, 5011 analytical cell, detector 1 0.55 V, detector 2 0.80 V, output of detector 2 is monitored

CHROMATOGRAM

Retention time: 7.8

Internal standard: methylrisperidone (R68808) (14.3)

Limit of quantitation: 2 ng/mL

OTHER SUBSTANCES

Extracted: metabolites, chlorpromazine, clomipramine, cyamemazine, desipramine, droperidol, flunitrazepam, haloperidol, imipramine, trihexyphenidyl

Noninterfering: alprazolam, bromazepam, carbamazepine, chlorazepate, diazepam, diphenylhydantoin, estazolam, ethylbenzatropine, oxazepam, phenobarbital, triazolam, valproic acid

Interfering: pipamperone

KEY WORDS

plasma; SPE

REFERENCE

Le Moing, J.P.; Edouard, S.; Levron, J.C. Determination of risperidone and 9-hydroxyrisperidone in human plasma by high-performance liquid chromatography with electrochemical detection, *J. Chromatogr.*, **1993**, *614*, 333-339.

SAMPLE

Matrix: blood, feces, urine

Sample preparation: Urine. Inject up to 1.95 mL urine directly. Feces. Homogenize (Ultra-Turrax TP-25) feces with MeOH, centrifuge, repeat extraction twice more. Evaporate a 10 mL aliquot of the extracts to dryness under a stream of nitrogen, reconstitute in DMSO, inject a 200 μ L aliquot. Plasma. Add an equal volume of MeCN to plasma, centrifuge. Remove the supernatant and reduce the volume under a stream of nitrogen at 40°, inject a 2 mL aliquot.

HPLC VARIABLES

Column: 300 \times 4.6 5 μ m Hypersil C18

Mobile phase: Gradient. A was 100 mM ammonium acetate containing 0.2% diethylamine, pH 6.0. B was MeCN:MeOH:1 M ammonium acetate containing 2% diethylamine, pH 6.0 80:10:10. A:B 100:0 to 50:50 over 1 h, maintain at 50:50 for 5 min, to 0:100 over 5 min.

Injection volume: 200-2000

Detector: UV 280 or radioactivity

CHROMATOGRAM

Retention time: 49

OTHER SUBSTANCES

Extracted: metabolites

REFERENCE

Mannens,G.; Huang,M.-L.; Meuldermans,W.; Hendrickx,J.; Woestenborghs,R.; Heykants,J. Absorption, metabolism, and excretion of risperidone in humans, *Drug Metab.Dispos.*, **1993**, 21, 1134-1141.

SAMPLE

Matrix: blood, tissue, urine

Sample preparation: Plasma, urine. 1 mL Plasma or urine + 200 μ L MeOH + 100 μ L 2 μ g/mL IS in MeOH + 1 mL 50 mM sodium borate, add 4 mL ethyl acetate, rotate at 10 rpm for 10 min, centrifuge at 1000 g for 10 min, repeat extraction. Combine the organic layers and add them to 3 mL 50 mM sulfuric acid, extract. Add the aqueous layer to 150 μ L concentrated ammonia, extract twice with 2.5 mL portions of heptane:isoamyl alcohol 90:10. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 60°, reconstitute the residue in 80 μ L mobile phase, inject a 30 μ L aliquot. Tissue. Grind (Waring blender) tissue, homogenize (Ultra-Turrax) with 4 volumes of water. 1 mL Homogenate + 100 μ L 2 μ g/mL IS in MeOH + 1 mL 100 mM sodium borate, add 4 mL ethyl acetate, rotate at 10 rpm for 10 min, centrifuge at 1000 g for 10 min, repeat extraction. Combine the organic layers and add them to 3 mL 50 mM sulfuric acid, extract. Add the aqueous layer to 150 μ L concentrated ammonia, extract twice with 2.5 mL portions of heptane:isoamyl alcohol 90:10. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 60°, reconstitute the residue in 80 μ L mobile phase, inject a 30 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 2.1 5 μ m ODS-Hypersil

Mobile phase: MeCN:water:diethylamine 35:65:0.02

Flow rate: 0.8

Injection volume: 30

Detector: UV 280

CHROMATOGRAM

Retention time: 2.9

Internal standard: 3-[2-[4-([6-fluoro-1,2-benzisoxazol-3-yl]methyl)-1-piperidinyl]ethyl]-2,7-dimethyl-4H-pyrido[1,2-a]pyrimidin-4-one (Janssen R 68808) (5.4)

Limit of quantitation: 10 ng/g (tissue), 2 ng/mL (plasma, tissue)

OTHER SUBSTANCES

Extracted: metabolites

Simultaneous: diazepam, oxazepam

Interfering: desmethyldiazepam

KEY WORDS

plasma; dog; human; muscle; pharmacokinetics

REFERENCE

Woestenborghs,R.; Lorreyne,W.; Van Rompaey,F.; Heykants,J. Determination of risperidone and 9-hydroxyrisperidone in plasma, urine and animal tissues by high-performance liquid chromatography, *J.Chromatogr.*, **1992**, 583, 223-230.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4 5 µm LiChrospher 100 RP-8 (B)

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 9.12 (A), 4.63 (B)

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordi-azepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyridamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenopropfen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazin-
dol, mefenamic acid, meperidine, mephenytoin, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl-
dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymet-
azoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, phen-
iramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltolox-
amine, phenytoin, pimo-
zide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quin-
ine, racemethorphan, ranitidine, remoxipride, salicylic acid, scopolamine, secobarbital, ser-
traline, sotalol, spironolactone, sulfapyrazone, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocin-
ide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

KEY WORDS

details of plasma extraction

REFERENCE

Koves,E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J.Chromatogr.A*, **1995**, 692, 103-119.

Ritodrine


Molecular formula: $C_{17}H_{21}NO_3$

Molecular weight: 287.36

CAS Registry No.: 26652-09-5, 23239-51-2 (HCl)

Merck Index: 8401

Lednicer No.: 2 39



SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 1 mL water + 1 mL 600 mM potassium carbonate + 8 mL ethyl acetate, shake mechanically for 5 min, centrifuge at 2000 rpm for 15 min. Remove the organic layer and add it to 1.2 mL 100 mM HCl, shake for 5 min, centrifuge at 2000 rpm for 10 min. Remove the aqueous layer and add it to 500 μ L 600 mM potassium carbonate and 1 mL ethyl acetate, shake for 5 min, centrifuge at 2000 rpm for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at room temperature, reconstitute the residue in 300 μ L mobile phase, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 5 μ m Biophase ODS

Mobile phase: MeCN:buffer:water 20:10:70 containing 3 mM sodium 1-heptanesulfonate, pH 3.7 (Buffer was 2.1 M acetic acid containing 400 mM ammonium acetate.)

Flow rate: 1

Injection volume: 50

Detector: E, Bioanalytical Systems Model 4B, glassy carbon electrode 0.95 V, Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 9.33

Internal standard: nalbuphine (13.58)

Limit of detection: 0.6 ng/mL

KEY WORDS

plasma

REFERENCE

Kuhnert,P.; Erhard,P.; Dixon,A.; Kuhnert,B.; Gross,T. Determination of ritodrine in plasma using HPLC, *J.Liq.Chromatogr.*, **1983**, 6, 2775-2783.

SAMPLE

Matrix: blood

Sample preparation: Adjust pH of 1 mL serum to 9.4 with 2 M sodium carbonate buffer, add 200 ng IS, vortex, add 6 mL ethyl acetate, shake for 15 min, centrifuge. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 35°, reconstitute the residue in 400 μ L mobile phase, vortex vigorously, inject a 40 μ L aliquot.

HPLC VARIABLES

Guard column: 71 \times 1.2 HC Pellosil (Whatman)

Column: 260 \times 4.6 μ Bondapak C18

Mobile phase: MeOH:buffer 25:75 (Buffer was 10 mM KH_2PO_4 containing 0.3 mM sodium octanesulfonate and 0.1 mM EDTA, pH 4.5.)

Flow rate: 1.5

Injection volume: 40

Detector: E, Bioanalytical Systems Model LC-4B, glassy carbon electrode +0.90 V, Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 7.8

Internal standard: 1-(3,5-dihydroxyphenyl)-2-(1,1-dimethylbutylamino)ethanol (12.0)

Limit of detection: 0.2 ng

KEY WORDS

serum

REFERENCE

Lin,L.S.; Caritis,S.N.; Wong,L.K. Analysis of ritodrine in serum by high-performance liquid chromatography with electrochemical detection, *J.Pharm.Sci.*, **1984**, 73, 131-133.

SAMPLE

Matrix: blood

Sample preparation: Dilute blood with an equal volume of water. 1 mL Plasma or 900 μ L diluted blood + 50 μ L 1.6 μ g/mL isoxsuprine hydrochloride + 0.9-1 mL buffer + 5 mL freshly distilled ethyl acetate, vortex for 1 min, centrifuge at 1750 g for 7 min. Remove the organic layer and evaporate it almost to dryness under a stream of nitrogen at 57°, evaporate the final 500 μ L at room temperature, reconstitute the residue in 100 μ L MeCN, vortex for 15 s, inject the whole amount. (Buffer was 26.5 g sodium carbonate and 21 g sodium bicarbonate in 500 mL water, pH 9.48.)

HPLC VARIABLES

Column: 300 \times 3.9 10 μ m μ Bondapak phenyl (plasma) or 200 \times 4.6 5 μ m Spheri-5 RP-18 (blood)

Mobile phase: MeCN:0.05% orthophosphoric acid 17:83 (plasma) or 63:37 (blood)

Flow rate: 2

Injection volume: 100

Detector: F ex 200 no emission filter or UV 254

CHROMATOGRAM

Retention time: 3.7 (plasma), 9.3 (blood)

Internal standard: isoxsuprine hydrochloride (15.1 (plasma), 16.3 (blood))

Limit of detection: 1 ng/mL (F)

OTHER SUBSTANCES

Simultaneous: fenoterol

Noninterfering: acetaminophen, albuterol, betamethasone, bupivacaine, caffeine, chloral hydrate, dexamethasone, diazepam, lignocaine, meperidine, metoclopramide, morphine, nitrazepam, terbutaline

KEY WORDS

plasma; whole blood; pharmacokinetics

REFERENCE

Gross,A.S.; Brown,K.F.; Baird-Lambert,J.A.; Nation,R.L. Determination of ritodrine in blood and plasma by high-performance liquid chromatography with fluorescence detection, *J.Chromatogr.*, **1987**, 416, 400-408.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol:n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 μ L mobile phase, centrifuge at 2800 g for 5 min, inject a 50 μ L aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 \times 3.9 4 μ m NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH_2PO_4 adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 225

CHROMATOGRAM

Retention time: 3.43

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfapyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procabazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melfalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vindsesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpi-pramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254–262.


Ritonavir

Molecular formula: C₃₇H₄₈N₆O₅S₂

Molecular weight: 720.96

CAS Registry No.: 155213-67-5

Merck Index: 8402



SAMPLE

Matrix: bile, blood, feces, urine

Sample preparation: Rat, dog plasma. Add 3 volumes MeCN to plasma, vortex, centrifuge, evaporate the supernatant to dryness under a stream of nitrogen, reconstitute the residue in mobile phase, filter (0.45 μ m) while centrifuging, inject an aliquot. Human plasma. Add 3 volumes MeCN to plasma, vortex, centrifuge, rinse the protein pellet with MeCN and MeOH. Evaporate the combined supernatants to dryness under a stream of nitrogen, reconstitute the residue with 200 μ L MeCN:MeOH:25 mM ammonium acetate 15:15:70 adjusted to pH 4.8 with trifluoroacetic acid, inject an aliquot. Feces. Centrifuge fecal homogenate, rinse the pellet twice with EtOH, filter the supernatants (0.45 μ m), inject an aliquot. Bile. Filter bile and inject an aliquot. Urine. Concentrate urine under a stream of nitrogen or via lyophilization, reconstitute in mobile phase, filter while centrifuging, inject an aliquot.

HPLC VARIABLES

Guard column: 5 μ m Nucleosil C18

Column: 250 \times 4.6 5 μ m Beckman Ultrasphere C18

Mobile phase: Gradient. A was 0.1% trifluoroacetic acid adjusted to pH 4.8 with ammonium acetate. B was MeCN. A:B from 75:25 to 40:60 over 50 min

Flow rate: 1

Detector: A UV 220; B Radioactivity, Flo-One/Beta Model A-500 radioactivity with 0.5 mL flow cell

CHROMATOGRAM

Retention time: 43

Limit of detection: 150 ng/mL (plasma)

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; pharmacokinetics; rat; dog; human; radiolabeled

REFERENCE

Denissen,J.F.; Grabowski,B.A.; Johnson,M.K.; Buko,A.M.; Kempf,D.J.; Thomas,S.B.; Surber,B.W. Metabolism and disposition of the HIV-1 protease inhibitor ritonavir (ABT-538) in rats, dogs, and humans, *Drug Metab. Dispos.*, **1997**, *25*, 489–501.

SAMPLE

Matrix: blood

Sample preparation: Add 1 mL of 500 mM sodium carbonate and 100 μ L IS to 1 mL of plasma, extract with 6 mL of MTBE. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at room temperature, dissolve in 300 μ L MeCN:10 mM perchloric acid 45:55, add 4 mL of hexane, vortex, centrifuge to separate the layers, freeze in a dry ice/acetone bath, discard the organic layer, inject a 100 μ L aliquot of the aqueous layer.

HPLC VARIABLES

Guard column: 23 \times 4 5 μ m YMC ODS-AQ

Column: 50 \times 4 3 μ m YMC ODS-AQ

Mobile phase: MeCN:MeOH:10 mM tetramethylammonium perchlorate in water:trifluoroacetic acid 40:5:55:0.1

Flow rate: 1.5

Injection volume: 100

Detector: UV 205

CHROMATOGRAM

Limit of quantitation: 5.6 ng/mL

OTHER SUBSTANCES

Extracted: ABT-378

KEY WORDS

plasma

REFERENCE

Bryan,P.; el-Shourbagy,T.; Emry,M.; Marsh,K.; McDonald,E.; Brooks,R.; Sapochak,L.; Hsu-Beischer,R.; Chu,S.
A sensitive and specific HPLC method for the simultaneous determination of ABT-378 and ritonavir in human plasma using uv detection (Abstract 2646), *Pharm.Res.*, **1997**, *14*, S427.

SAMPLE

Matrix: blood, CSF, saliva

Sample preparation: Add 400 μ L MeCN to 100 μ L plasma, CSF, or saliva, vortex for 30 s, centrifuge at 10500 g for 3 min, evaporate 400 μ L supernatant to dryness under a gentle stream of nitrogen at 40°, redissolve the residue in 150 μ L mobile phase, vortex for 60 s, centrifuge at 10500 g for 3 min, inject an aliquot.

HPLC VARIABLES

Guard column: 10 \times 3 Chromguard C18 (Chrompack, Netherlands)

Column: 75 \times 4.6 3.5 μ m Zorbax SB-C18

Mobile phase: MeCN:buffer 44:56 (Buffer was 25 mM sodium acetate with 25 mM hexane-1-sulfonic acid adjusted to pH 4.0 with 37% HCl.)

Flow rate: 1

Injection volume: 100

Detector: UV 239

CHROMATOGRAM

Retention time: 9

Limit of detection: 20 ng/mL

Limit of quantitation: 50 ng/mL

OTHER SUBSTANCES

Noninterfering: didanosine, fluconazole, folinic acid, ganciclovir, indinavir, lamivudine, methadone, methotrexate, nelfinavir, oxazepam, pyrazinamide, ranitidine, rifampin, saquinavir, stavudine, sulfamethoxazole, trimethoprim, zalcitabine, zidovudine, zidovudine glucuronide

KEY WORDS

plasma

REFERENCE

Hoetelmans,R.M.W.; van Essenberg,M.; Profijt,M.; Meenhorst,P.L.; Mulder,J.W.; Beijnen,J.H. High-performance liquid chromatographic determination of ritonavir in human plasma, cerebrospinal fluid and saliva, *J.Chromatogr.B*, **1998**, *705*, 119–126.

SAMPLE

Matrix: microsomal incubations

Sample preparation: Add four volumes of MeCN to microsomal incubation, centrifuge at 1500 g for 10 min, evaporate supernatant to dryness under nitrogen at 40°, resuspend the residue in 100-120 μ L mobile phase, inject an 80 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 Hypersil 5C18 (Phenomenex)

Mobile phase: Gradient. A was MeCN. B was 0.1% triethylamine in water adjusted to pH 4.8 with ammonium acetate. A:B from 25:75 to 75:25 over 15 min

Flow rate: 1.5

Injection volume: 80

Detector: UV 210

CHROMATOGRAM

Retention time: 12.9

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

human; liver

REFERENCE

Koudriakova,T.; Iatsimirskaia,E.; Utkin,I.; Gangl,E.; Vouros,P.; Storozhuk,E.; Orza,D.; Marinina,J.; Gerber,N. Metabolism of the human immunodeficiency virus protease inhibitors indinavir and ritonavir by human intestinal microsomes and expressed cytochrome P4503A4/3A5: Mechanism-based inactivation of cytochrome P4503A by ritonavir, *Drug Metab.Dispos.*, **1998**, 26, 552–561.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 150 × 3.9 5 μm Delta-pak C4 (Waters)

Mobile phase: MeCN:buffer 35:65 (Buffer was 10 mM ammonium dihydrogen phosphate and 1 mM 1-heptanesulfonic acid sodium salt, pH adjusted to 4.8 with ammonium hydroxide.)

Flow rate: 0.6

Injection volume: 35

Detector: UV 210

CHROMATOGRAM

Retention time: 22-27

OTHER SUBSTANCES


Simultaneous: indinavir, nelfinavir, saquinavir

Noninterfering: didanosine, lamivudine, stavudine, zalcitabine, zidovudine

REFERENCE

Iayewardene,A.L.; Zhu,F.; Aweeka,F.T.; Gambertoglio,J.G. Simple high-performance liquid chromatographic determination of the protease inhibitor indinavir in human plasma, *J.Chromatogr.B*, **1998**, 707, 203–211.

Rizatriptan



Molecular formula: C₁₅H₁₉N₅

Molecular weight: 269.35

CAS Registry No.: 144034-80-0, 145202-66-0 (benzoate), 59776-67-7 (sulfate)

SAMPLE

Matrix: bulk

Sample preparation: Prepare an 1 mg/mL solution of the free base in MeCN:water 10:90. Inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 μm Zorbax SB-phenyl

Mobile phase: MeCN:0.1% trifluoroacetic acid in water 16: 84 (A) or MeCN:0.1% phosphoric acid in water 12:88 (B)

Flow rate: 1.5

Detector: UV 280

CHROMATOGRAM

Retention time: 6.3 (A), 6.7 (B)

OTHER SUBSTANCES

Extracted: impurities


REFERENCE

Antonucci,V.; Wright,L.; Toma,P. The reversed-phase liquid chromatographic behavior of the new 5-HT_{1D} receptor agonist rizatriptan benzoate and its potential process impurities, *J.Liq.Chromatogr.Rel.Technol.*, **1998**, 21, 1649–1670.

SAMPLE**Matrix:** bulk**Sample preparation:** Prepare an 1 mg/mL solution of the free base in MeCN:water 10:90. Inject an aliquot.**HPLC VARIABLES****Column:** 150 × 3.9 5 µm Symmetry C8**Mobile phase:** MeCN:0.1% trifluoroacetic acid in water 16: 84 (A) or MeCN:0.1% phosphoric acid in water 12:88 (B)**Flow rate:** 1.5**Detector:** UV 280**CHROMATOGRAM****Retention time:** 5.5 (A), 6.7 (B)**OTHER SUBSTANCES****Extracted:** impurities**REFERENCE**

Antonucci,V.; Wright,L.; Toma,P. The reversed-phase liquid chromatographic behavior of the new 5-HT_{1D} receptor agonist rizatriptan benzoate and its potential process impurities, *J.Liq.Chromatogr.Rel.Technol.*, **1998**, *21*, 1649–1670.

Robenidine

Molecular formula: C₁₅H₁₃Cl₂N₅**Molecular weight:** 334.21**CAS Registry No.:** 25875-51-8,
25875-50-7 (HCl)**Merck Index:** 8403**SAMPLE****Matrix:** bulk**Sample preparation:** Prepare a 200 µg/mL solution in dichloromethane:MeOH 90:10, evaporate an aliquot to dryness under a stream of nitrogen at 60°, reconstitute with 1 mL 4 mg/mL 4-dimethylaminopyridine in dichloromethane, add 1 mL 2 mg/mL dansyl chloride in dichloromethane, heat in a capped tube at 80° for 1 h, cool on ice, add to a 3 mL Sep-Pak silica SPE cartridge, rinse the tube with 5 mL dichloromethane, add the rinse to the SPE cartridge, elute with 5 mL ethyl acetate. Evaporate the eluate to dryness under a stream of nitrogen at 60°, reconstitute the residue in 2 mL mobile phase, inject a 20 µL aliquot.**HPLC VARIABLES****Guard column:** 150 × 3.2 7 µm silica (Brownlee)**Column:** 250 × 4.6 5 µm Zorbax Sil**Mobile phase:** Hexane:chloroform:THF:MeOH 50:50:2:1**Flow rate:** 2**Injection volume:** 20**Detector:** F ex 320 em 485**CHROMATOGRAM****Retention time:** 4**Limit of detection:** 400 ng/mL**KEY WORDS**

derivatization; normal phase; protect from light; SPE

REFERENCE

Cohen, H.; Armstrong, F.; Campbell, H. Sensitive fluorescence detection of robenidine by derivatization with dansyl chloride and high-performance liquid chromatography, *J. Chromatogr. A*, **1995**, 694, 407-413.


Rocuronium bromide

Molecular formula: $C_{32}H_{53}BrN_2O_4$

Molecular weight: 609.69

CAS Registry No.: 119302-91-9

Merck Index: 8407



SAMPLE

Matrix: bile, blood, stoma fluid, tissue, urine

Sample preparation: Homogenize (Ultra-Turrax) 1 g tissue with 9 mL 1 M NaH_2PO_4 for 10 min. Acidify 1 mL plasma, urine, or bile with 200 μ L 1 M NaH_2PO_4 . Homogenize 1 mL stoma fluid with 200 μ L 1 M NaH_2PO_4 . Make up 50-1000 μ L plasma, 200-1000 μ L urine, 5-200 μ L bile, 1000 μ L stoma fluid, or 100-1000 μ L tissue homogenate to 2 mL with water, add 1 mL buffer, add 150 ng IS, add 7 mL dichloromethane, vortex for 15 s, centrifuge at 740 g for 5 min. Remove the lower organic layer and evaporate it to dryness under a stream of air at 37°, reconstitute the residue in 200 μ L mobile phase, inject a 100 μ L aliquot (or less). (Buffer was prepared by mixing 6 mL of an aqueous solution containing 7.505 mg/mL glycine and 5.85 mg/mL NaCl, 4 mL 100 mM NaOH and 6.2 g KI.)

HPLC VARIABLES

Guard column: 4 × 6 μ Bondapak C18

Column: 150 × 3.9 5 μ m Lichrospher 100-RP18

Mobile phase: Dioxane:buffer 16:84 (Caution! Dioxane is a carcinogen!) (Buffer was 100 mM NaH_2PO_4 containing 0.11 mM 9,10-dimethoxyanthracene-2-sulfonate and 0.11 mM 1-heptanesulfonic acid, pH adjusted to 3.0 with orthophosphoric acid. After each series of analyses flush column with 15 mL water and 75 mL MeOH.)

Flow rate: 1

Injection volume: 100

Detector: F ex 385 em 452 following post-column extraction. The column effluent mixed with dichloroethane pumped at 1 mL/min and the mixture flowed through a 1 m × 0.25 mm i.d. stainless steel coil to a phase separator (Organon International) then the organic phase flowed through the detector (*J. Chromatogr.* 1987, 421, 327; *Anal. Chim. Acta* 1987, 192, 267).

CHROMATOGRAM

Retention time: 9

Internal standard: 1-(3 α ,17 β -dihydroxy-2 β -morpholino-5 α -androstan-16 β -yl)-1-methylpiperidinium bromide (Org 7402, Organon) (21)

Limit of detection: 5 ng tissue, 4 ng (urine, bile), 3 ng (plasma)

Limit of quantitation: 20 ng/mL (stoma fluid), 100 ng/mL (bile), 25 ng/mL (urine), 10 ng/mL (plasma)

OTHER SUBSTANCES

Extracted: metabolites

Simultaneous: cefotaxime, ceftazidime, ceftriaxone, cefuroxime, cephamandole, cephradine, dixyrazine, meperidine, metocurine, metoprolol, sulfamethoxazole, trimethoprim

Noninterfering: alfentanil, aprotinin, atropine, bupivacaine, chlorpromazine, daltaparin, dexamethasone, diazepam, dopamine, droperidol, etomidate, fentanyl, furosemide, gallamine, haloperidol, midazolam, morphine, neostigmine, nitroglycerin, nitroprusside, oxytocin, pancuron-

ium, pentobarbital, phenylephrine, phenytoin, pipecuronium, piperacillin, promethazine, propofol, ranitidine, succinylcholine, sufentanil, terbutaline, thiopental, vecuronium, verapamil
Interfering: alizapride, atracurium, ketamine, ketogan, lidocaine, metoclopramide, nimodipine, prochlorperazine, tubocurarine

KEY WORDS

human; dog; plasma; liver; lung

REFERENCE

Kleef, U.W.; Proost, J.H.; Roggevel, J.; Wierda, J.M.K.H. Determination of rocuronium and its putative metabolites in body fluids and tissue homogenates, *J. Chromatogr.*, **1993**, 621, 65–76.

SAMPLE

Matrix: blood

Sample preparation: Condition a 1 mL Bond-Elut C18 SPE cartridge with 1 mL MeOH:MeCN 2:1 and 1 mL water. Acidify 1 mL plasma with 200 μ L 1 M NaH_2PO_4 . Add 20–200 ng IS to 1 mL acidified plasma, add to the SPE cartridge, wash with 1 mL water, wash with 1 mL 100 mM pH 3 NaH_2PO_4 , elute with 400 μ L mobile phase, discard first 100 μ L eluate, inject a 200 μ L aliquot of the remaining eluate (from *J. Chromatogr.* 1987, 421, 327; modifications may be necessary).

HPLC VARIABLES

Column: 150 \times 3.9 4 μ m Nova-Pak C18

Mobile phase: Dioxane:water 20:80 containing 100 mM NaH_2PO_4 and 0.44 mM 9,10-dimethoxyanthracene-2-sulfonate, pH adjusted to 3 with phosphoric acid. (Caution! Dioxane is a carcinogen!) (After each series of analyses flush column with 200 mL MeOH then re-equilibrate with 120 mL mobile phase.)

Flow rate: 1

Injection volume: 200

Detector: F ex 380 em 452 following post-column extraction. The column effluent mixed with dichloroethane pumped at 1.6 mL/min and the mixture flowed through a 1 m \times 0.25 mm i.d. stainless steel coil to a phase separator (*Anal. Chim. Acta* 1987, 192, 267) then the organic phase flowed through the detector.

CHROMATOGRAM

Internal standard: 3,17-didesacetyl vecuronium (Org 7402)

Limit of quantitation: 10 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

SPE; plasma; pharmacokinetics; post-column reaction; post-column extraction

REFERENCE

Cooper, R.A.; Maddineni, V.R.; Mirakhur, R.K.; Wierda, J.M.K.H.; Brady, M.; Fitzpatrick, K.T.J. Time course of neuromuscular effects and pharmacokinetics of rocuronium bromide (Org 9426) during isoflurane anaesthesia in patients with and without renal failure, *Br. J. Anaesth.*, **1993**, 71, 222–226.


Rokitamycin

Molecular formula: C₄₂H₆₉NO₁₅

Molecular weight: 828.01

CAS Registry No.: 74014-51-0

Merck Index: 8408



SAMPLE

Matrix: blood

Sample preparation: Collect blood in tubes containing neostigmine methyl sulfate, final concentration 0.2 mM. 1 mL Plasma + 100 μ L 2 μ g/mL josamycin in MeOH, add 50 μ L 100 mM pH 4.65 acetate buffer + 100 μ L saturated NaCl + 100 μ L 10 mM sodium lauryl sulfate + 5 mL hexane:isoamyl alcohol 90:10, shake for 15 min, centrifuge at 1000 g for 10 min, repeat extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen with mild heating, reconstitute the residue in 200 μ L 0.002% dansylhydrazine in toluene:MeOH:acetic acid 90:10:1.13 (freshly prepared), heat at 60° for 20 min, evaporate to dryness under a stream of nitrogen at 40°, reconstitute with 20 μ L ethyl acetate, add 60 μ L 100 mM HCl, inject a 20 μ L aliquot of the lower aqueous phase. (Silanize glassware with 1% trimethylchlorosilane in toluene for 1 h at 70°, wash twice with MeOH.)

HPLC VARIABLES

Column: 150 × 4.6 3 μm Nucleosil C18

Mobile phase: MeCN:50 mM ammonium acetate 72:28

Column temperature: 32.5

Flow rate: 0.8

Injection volume: 20

Detector: F ex 352 em 537

CHROMATOGRAM

Retention time: 11.5

Internal standard: josamycin (9.5)

Limit of detection: 20 ng/mL

Limit of quantitation: 50 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; derivatization

REFERENCE

Tod,M.; Biarez,O.; Nicolas,P.; Petitjean,O. Sensitive determination of josamycin and rokitamycin in plasma by high-performance liquid chromatography with fluorescence detection. *J.Chromatogr.* **1992**, 575, 171-176.


Romurtide

Molecular formula: $C_{43}H_{78}N_6O_{13}$

Molecular weight: 887.12

CAS Registry No.: 78113-36-7

Merck Index: 8412



SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 5 μ m octadecyl silica

Mobile phase: MeCN:MeOH:20 mM ammonium acetate 66:4:40

Column temperature: 25

Flow rate: 1.5

Detector: UV 220


CHROMATOGRAM

Retention time: 8.6

REFERENCE

Moroi,R.; Yamazaki,K.; Hirota,T.; Watanabe,S.; Kataoka,K.; Ichinose,M. Physico-chemical properties of murec-tasin, *Arzneimittelforschung*, **1988**, 38, 953-959.

Ronidazole



Molecular formula: $C_6H_8N_4O_4$

Molecular weight: 200.15

CAS Registry No.: 7681-76-7

Merck Index: 8413

Lednicer No.: 2 245

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 308

CHROMATOGRAM

Retention time: 8.265

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.


Ropinirole

Molecular formula: $C_{16}H_{24}N_2O$

Molecular weight: 260.38

CAS Registry No.: 91374-21-9, 91374-20-8 (HCl)

Merck Index: 8416



SAMPLE

Matrix: microsomal incubations

Sample preparation: 250 μ L Microsomal incubation + 50 μ L 5% trichloroacetic acid, centrifuge, inject an aliquot of the supernatant.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Supelcosil LC-ABZ (Supelco)

Mobile phase: Gradient. A was MeCN. B was 100 mM pH 4 ammonium acetate. A:B from 0:100 to 15:85 over 10 min, to 100:0 over 12 min, maintain at 100:0 for 15 min

Column temperature: 40

Flow rate: 1

Detector: UV 250; Radioactivity, Ramona-5 (Lablogic, Inc., UK)

CHROMATOGRAM

Retention time: 15

OTHER SUBSTANCES

Simultaneous: metabolites

KEY WORDS

liver

REFERENCE

Bloomer, J.C.; Clarke, S.E.; Chenery, R.J. In vitro identification of the P450 enzymes responsible for the metabolism of ropinirole, *Drug Metab. Dispos.*, **1997**, 25, 840-844.

SAMPLE

Matrix: solutions

Sample preparation: Inject a 20 μ L aliquot of a 100-500 μ g/mL solution in mobile phase.

HPLC VARIABLES

Column: 100 \times 4.6 5 μ m Hypersil C8 MOS 100A coated with phosphatidylcholine (95% pure soybean lecithin, Epikuron, Lucas Meyer & Co.) (Coat column by recycling a 1 mM solution of phosphatidylcholine in MeOH:water 80:20 for 24 h.)

Mobile phase: MeCN:35 mM pH 7.4 sodium phosphate buffer 40:60

Flow rate: 0.5–2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: k' 0.95

OTHER SUBSTANCES

Also analyzed: amoxicillin, antipyrine, carbamazepine, chlorpheniramine, chlorpromazine, clonidine, codeine, desipramine, diphenhydramine, dipyrindamole, ephedrine, flufenamic acid, haloperidol, hydroxyzine, imipramine, indomethacin, lidocaine, megestrol acetate, metoprolol, nabumetone, nadolol, phenobarbital, phenol, promazine, propranolol, pyrilamine, quinidine, testosterone, thioridazine, tolfenamic acid, verapamil

Noninterfering: acetaminophen, aspirin, azathioprine, caffeine, carprofen, chlorambucil, cimetidine, fenoterol, flurbiprofen, ibuprofen, ketoprofen, ranitidine, salicylic acid, sulfamethoxazole, theophylline, thioguanine, tiaprofenic acid, trimethoprim, valproic acid


KEY WORDS

comparison with capillary electrophoresis

REFERENCE

Hanna,M.; de Biasi,V.; Bond,B.; Salter,C.; Hutt,A.J.; Camilleri,P. Estimation of the partitioning characteristics of drugs: A comparison of a large and diverse drug series utilizing chromatographic and electrophoretic methodology, *Anal.Chem.*, **1998**, *70*, 2092–2099.

Ropivacaine



Molecular formula: $C_{17}H_{26}N_2O$

Molecular weight: 274.41

CAS Registry No.: 84057-95-4, 98717-15-8 (HCl),
132112-35-7 (HCl monohydrate)

Merck Index: 8417

SAMPLE

Matrix: blood, urine

Sample preparation: Condition a Bond Elut SCX SPE cartridge. Hydrolyze 1 mL urine with 6 M HCl in a water-bath at 95° for 1 hr, dilute with water 1:5. Acidify plasma with phosphoric acid. Add the sample to the SPE cartridge, wash with 200 mM pH 4.5 acetate buffer, wash with MeOH:buffer 50:50, elute with MeOH:2 M ammonia 80:20, evaporate to dryness, reconstitute with 200 μ L mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 125 \times 2 Superspher RP-Select B (plasma) or 125 \times 4 Superspher RP-Select B (urine)

Mobile phase: MeCN:pH 2 phosphate buffer 15-25:85-75 containing 5-15 mM octanesulfonic acid

Detector: UV 210

CHROMATOGRAM

Limit of quantitation: 10 ng/mL (plasma), 300 nM (urine)

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

SPE; plasma; pharmacokinetics

REFERENCE

Halldin,M.M.; Bredberg,E.; Angelin,B.; Arvidsson,T.; Askemark,Y.; Elofsson,S.; Widman,M. Metabolism and excretion of ropivacaine in humans, *Drug Metab.Dispos.*, **1996**, *24*, 962–968.

SAMPLE

Matrix: formulations

Sample preparation: Dilute with mobile phase to a ropivacaine hydrochloride concentration of 75 µg/mL, inject an aliquot.

HPLC VARIABLES

Guard column: 10 × 3 chiral-AGP

Column: 100 × 4.5 µm chiral-AGP (Chromtech)

Mobile phase: Isopropanol:pH 7.2 phosphate buffer (µ = 0.05) 7:93

Flow rate: 1

Injection volume: 20

Detector: UV 220

CHROMATOGRAM

Retention time: 9 (R), 12 (S)

KEY WORDS

chiral; injections; comparison with capillary electrophoresis

REFERENCE

Sänger-van de Griend,C.E.; Wahlström,H.; Gröningson,K.; Widahl-Näsman,M. A chiral capillary electrophoresis method for ropivacaine hydrochloride in pharmaceutical formulations: validation and comparison with chiral liquid chromatography, *J.Pharm.Biomed.Anal.*, **1997**, *15*, 1051–1061.

Rosoxacin


Molecular formula: C₁₇H₁₄N₂O₃

Molecular weight: 294.31

CAS Registry No.: 40034-42-2

Merck Index: 8426

Lednicer No.: 3 185

**SAMPLE**

Matrix: blood

Sample preparation: Filter 1 mL plasma using a micropartition system (MPS-1, Amicon, MA) while centrifuging at 2000 g for 20 min at 10°, inject an aliquot of the ultrafiltrate.

HPLC VARIABLES

Column: 250 × 4.6 Spherisorb ODS-2 endcapped

Mobile phase: MeCN:buffer 20:80 containing 5 mM tetrabutylammonium sulfate, adjusted to pH 2.5 with 1 M NaOH (Buffer was 100 mM citric acid containing 200 mM ammonium perchlorate.)

Column temperature: 37

Flow rate: 1

Detector: UV 280

CHROMATOGRAM

Retention time: 5.13

Internal standard: sparfloxacin (8.3)

KEY WORDS

plasma; ultrafiltrate

REFERENCE

Zlotos,G.; B cker,A.; Kinzig-Schippers,M.; Sorgel,F.; Holzgrabe,U. Plasma protein binding of gyrase inhibitors, *J.Pharm.Sci.*, **1998**, *87*, 215–220.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10–30

Detector: UV 271.2

CHROMATOGRAM

Retention time: 13.467

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; P pin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, *763*, 149–163.


Rotenone

Molecular formula: C₂₃H₂₂O₆

Molecular weight: 394.42

CAS Registry No.: 83-79-4

Merck Index: 8427



SAMPLE

Matrix: feed

Sample preparation: 10 g Ground feed + 100 mL MeCN:glacial acetic acid 99:1, shake for 1 h, centrifuge, filter (0.45 µm) an aliquot of the supernatant, inject a 20 µL aliquot of the filtrate.

HPLC VARIABLES

Column: 300 × 3.9 µm Bondapak C18

Mobile phase: MeCN:water 60:40

Flow rate: 1

Injection volume: 20

Detector: UV 280

CHROMATOGRAM**Retention time:** 8.2

OTHER SUBSTANCES**Simultaneous:** rotenonone

REFERENCE

Kline, D.A.; Hanna, G.R.; Honaker, C.B.; Kuhn, G.O.; Jameson, C.W. Preparation and stability of animal feed mixtures dosed with rotenone, *J. Assoc. Off. Anal. Chem.*, **1986**, 69, 660–663.

SAMPLE**Matrix:** sediment, tissue

Sample preparation: Fish, crayfish, mussels. Homogenize (Waring blender) with dry ice, place in freezer overnight. 10 g Sample + 60 (fish), 70 (fish offal, crayfish) or 100 (mussels) g anhydrous sodium sulfate, mix (Sorval Omni-mixer), place sample on top of 50 mm anhydrous sodium sulfate in a 22 mm diameter glass column, place a 50 mM layer of anhydrous sodium sulfate on top of the column, pass 100 mL diethyl ether through the column at 1.5 mL/min. Evaporate the eluate to dryness under reduced pressure at 30°, reconstitute the residue in 10 mL dichloromethane:cyclohexane 50:50, add a 5 mL aliquot to a 400 × 28 column packed with SX-3 BioBeads (Analytical BioChemistry Labs) in dichloromethane:cyclohexane 50:50, elute with dichloromethane:cyclohexane 50:50 at 4 mL/min, discard the first 112 mL of eluate, collect the second 112 mL and evaporate it to dryness under reduced pressure at 30°. Reconstitute with 5 mL benzene (Caution! Benzene is a carcinogen!) and add to the silica column, rinse flask with five 5 mL portions of benzene, add rinses to the column (do not allow column to go dry), elute with 70 mL benzene:acetone 97:3. Evaporate the eluate to dryness under reduced pressure at 30°, reconstitute in 5 mL MeOH, inject a 100 µL aliquot. Sediment. 10 g Sediment + 20 mL MeOH, mix (Sorval Omni-mixer), centrifuge at 1000 g for 5 min, repeat extraction 3 more times with 10 mL MeOH. Filter (Gelman type A/E glass fiber) the combined supernatants, evaporate the filtrate under reduced pressure at 30° to about 25 mL. Add the residue to 500 mL 100 mM HCl, extract three times with 20 mL hexane. Combine the hexane extracts and evaporate them to dryness under reduced pressure at 30°, reconstitute with 5 mL benzene (Caution! Benzene is a carcinogen!) and add to the silica column, rinse flask with five 5 mL portions of benzene, add rinses to the column (do not allow column to go dry), elute with 70 mL benzene:acetone 97:3. Evaporate the eluate to dryness under reduced pressure at 30°, reconstitute in 5 mL MeOH, inject a 100 µL aliquot. (Column was 5 g anhydrous sodium sulfate, 5 g 3% deactivated silica gel, and 5 g anhydrous sodium sulfate, all slurry packed in benzene. Silica gel was 40-140 mesh from J.T. Baker activated at 130° for 24 h.)

HPLC VARIABLES**Column:** 150 × 3.9 4 µm Nova-Pak C18**Mobile phase:** MeOH:water 70:30 (60:40 for fish offal)**Flow rate:** 1**Injection volume:** 100**Detector:** UV 295

CHROMATOGRAM**Retention time:** 5, 17 (for fish offal)**Limit of detection:** 5 ng/g (tissue), 25 ng/g (sediment)

KEY WORDSfish; crayfish; mussels; shellfish; sediment

REFERENCE

Dawson, V.K.; Allen, J.L. Liquid chromatographic determination of rotenone in fish, crayfish, mussels, and sediments, *J. Assoc. Off. Anal. Chem.*, **1988**, 71, 1094–1096.

SAMPLE**Matrix:** solutions**Sample preparation:** Filter (0.45 µm) water, inject a 200 µL aliquot.

HPLC VARIABLES**Column:** 250 × 5 Zorbax ODS

Mobile phase: MeCN:water 70:30
Flow rate: 1.3
Injection volume: 200
Detector: UV 210

CHROMATOGRAM

Retention time: 6.5
Limit of quantitation: 7.5 ppb

OTHER SUBSTANCES

Simultaneous: rotenonone


KEY WORDS

drinking water; stream water; pond water; water

REFERENCE

Bushway, R.J. High-performance liquid chromatographic analysis of rotenone and rotenonone in water by direct injection, *J. Chromatogr.*, **1984**, *303*, 263–266.

Roxarsone



Molecular formula: C₆H₆AsNO₆

Molecular weight: 263.04

CAS Registry No.: 121-19-7

Merck Index: 8430

SAMPLE

Matrix: feed

Sample preparation: Grind feed to pass through No. 20 sieve. 15 g Ground feed + 50 mL 20 g/L K₂HPO₄, shake for 5 min, centrifuge at 2385 rpm for 10 min. REmove a 30 mL aliquot of the supernatant and add it to 1 mL 2.5 M HCl, let stand for 15 min, centrifuge for 10 min. Remove a 25 mL aliquot and add it to 1 mL 6 M NaOH, add 2 g activated charcoal (Darco G-60 or equivalent), swirl, let stand for 30 min with periodic swirling, filter (0.45 μm) an aliquot, add the filtrate to an activated carbon SPE cartridge (Analtech No. 01-97), discard the first 2 mL, inject a 15 μL aliquot of the next 500 μL of eluate.

HPLC VARIABLES

Guard column: Bondapak C18

Column: 8NVC18 Radial Pak radial compression (Waters)

Mobile phase: MeOH:water 25:75 containing 3% PIC-A (tetrabutylammonium sulfate)

Flow rate: 1.2

Injection volume: 15

Detector: UV 243

CHROMATOGRAM

Retention time: 8

Limit of detection: <0.3 ng

OTHER SUBSTANCES

Noninterfering: bacitracin, bambemycins, BHA, BHT, chlortetracycline, erythromycin, ethoxyquin, furazolidone, hygromycin B, monensin, niacin, nicarbazin, ormetoprim, oxytetracycline, riboflavin, sulfadimethoxine, tetracycline, thiamine, vitamin D3, vitamin E, pyridoxine, vitamin A


KEY WORDS

SPE

REFERENCE

Sapp,R.E.; Davidson,S. Determination of Roxarsone in feeds using solid phase extraction and liquid chromatography with ultraviolet detection, *JAOAC Int.*, **1993**, *76*, 956–961.

Roxatidine acetate



Molecular formula: $C_{19}H_{28}N_2O_4$

Molecular weight: 348.44

CAS Registry No.: 78628-28-1, 93793-83-0 (HCl)

Merck Index: 8431

Lednicer No.: 5 26

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 150×4.6 12 μm 1-myristoyl-2-[(13-carboxyl)-tridecoyl]-sn-3-glycerophosphocholine chemically bonded to silica (Regis)

Mobile phase: MeCN:100 mM pH 7.0 phosphate buffer 20:80

Flow rate: 1

Detector: UV 254

CHROMATOGRAM

Retention time: k' 2.29

OTHER SUBSTANCES

Also analyzed: acebutolol, alprenolol, antazoline, atenolol, betaxolol, bisoprolol, bopindolol, bupranolol, carteolol, celiprolol, chlorpyramine, chlorpheniramine, cicloprolol, cimetidine, cinarizine, cirazoline, clonidine, dilevalol, dimethindene, diphenhydramine, doxazosin, esmolol, famotidine, isothipendyl, ketotifen, metiamide, metoprolol, moxonidine, nadolol, naphazoline, nifenalol, nizatidine, oxprenolol, pheniramine, phentolamine, pindolol, pizotyline (pizotifen), practolol, prazosin, promethazine, propranolol, pyrilamine (mepyramine), ranitidine, sotalol, tiamenidine, timolol, tramazoline, tripeleminamine, triprolidine, tymazoline, UK-14,304

REFERENCE

Kaliszan,R.; Nasal,A.; Turowski,M. Binding site for basic drugs on α_1 -acid glycoprotein as revealed by chemometric analysis of biochromatographic data, *Biomed.Chromatogr.*, **1995**, *9*, 211–215.


Roxithromycin

Molecular formula: $C_{41}H_{76}N_2O_{15}$

Molecular weight: 837.06

CAS Registry No.: 80214-83-1

Merck Index: 8433



SAMPLE

Matrix: blood

Sample preparation: 500 μ L Serum + 50 μ L 20 μ g/mL IS in MeOH + 2.5 mL hexane:isoamyl alcohol 95:5, vortex for 20 s, centrifuge at 1500 g for 5 min. Remove 2 mL of the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 200 μ L mobile phase, vortex for 20 s, sonicate for 1 min, centrifuge at 1500 g for 5 min, inject a 50 μ L aliquot of the supernatant.

HPLC VARIABLES

Guard column: 20 mm long Supelguard C18 (Supelco)

Column: 50 \times 4.6 5 μ m Supelcosil C18

Mobile phase: MeCN:MeOH:water 50:20:30 containing 4 g/L (?) ammonium acetate

Flow rate: 1

Injection volume: 50

Detector: E, ESA Coulochem 5100A, 5010 analytical cell, +0.89 V

CHROMATOGRAM

Retention time: 3.5

Internal standard: RU 29767 (7.5)

Limit of detection: 50 ng/mL

KEY WORDS

serum; pharmacokinetics

REFERENCE

Nilsen, O.G.; Aamo, T.; Zahlsen, K.; Svarva, P. Macrolide pharmacokinetics and dose scheduling of roxithromycin, *Diagn. Microbiol. Infect. Dis.*, **1992**, 15, 71S–76S.

SAMPLE

Matrix: blood

Sample preparation: Condition a 10 \times 2 20 mg 30–40 μ m Baker CN SPE cartridge with 2 mL MeOH, 2 mL MeOH:water 10:90, and 4 mL MeCN:pH 10.5 phosphate buffer (I = 0.10) 10:90 at 2 mL/min. Centrifuge plasma at 1300 g for 5 min, 100 μ L plasma + 100 μ L clarithromycin in MeCN:pH 10.5 phosphate buffer (I = 0.10) 10:90, mix, add a 20–100 μ L aliquot to the SPE cartridge, wash SPE cartridge with MeCN:pH 10.5 phosphate buffer (I = 0.10) 10:90 at 0.5 mL/min, after 5 min backflush the contents of the SPE cartridge onto the column with the mobile phase, elute the column with the mobile phase and monitor the effluent.

HPLC VARIABLES

Column: 100 \times 4.6 3 μ m Hypersil BDS C18

Mobile phase: MeCN:water 54:46 containing 4.5 mM NaH_2PO_4 and 6.8 mM Na_2HPO_4 , pH 7

Column temperature: 55

Flow rate: 1

Detector: E, ESA Coulochem II, Model 5011 dual analytical cell, upstream +0.65 V, downstream +0.85 V (monitored), analytical cell protected by an ESA carbon in-line filter

CHROMATOGRAM

Retention time: 7

Internal standard: clarithromycin (6)

Limit of quantitation: 500 nM

KEY WORDS

plasma; SPE

REFERENCE

Hedenmo, M.; Eriksson, B.-M. Liquid chromatographic determination of the macrolide antibiotics roxithromycin and clarithromycin in plasma by automated solid-phase extraction and electrochemical detection, *J. Chromatogr. A*, **1995**, 692, 161–166.

SAMPLE

Matrix: blood, gastric juice, gastric mucosa, saliva, vitreous humor

Sample preparation: Homogenize 5–20 mg gastric mucosa in 300 μ L 10 mM pH 7.4 sodium phosphate buffer with sonication. Add 500 ng clarithromycin in MeOH:water 50:50 to 500 μ L plasma, serum, saliva, gastric juice, leucocytes lysate, vitreous humor or 300 μ L gastric mucosa homogenate, vortex, add 200 μ L 100 mM sodium carbonate and 3 mL MTBE, shake thoroughly (5×2 s in an SMI Multi-tube vortexer), centrifuge at 1000 g for 5 min, freeze the aqueous layer in liquid nitrogen or in a freezer at -70° for 15 min. Evaporate the upper organic layer to dryness in a centrifugal vacuum evaporator (Jouan RC 10.22), reconstitute the residue in 250 μ L MeOH:water 50:50, inject a 20–50 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Zorbax SB CN

Mobile phase: MeCN:MeOH:50 mM Na_2HPO_4 and NaH_2PO_4 buffer 53:41:6 (pH 7.0) (The mobile phase was a mixture of 350 mL MeCN, 50 mL MeOH and 450 mL 50 mM Na_2HPO_4 and NaH_2PO_4 buffer.)

Column temperature: 30

Flow rate: 1

Injection volume: 20–50

Detector: E, ESA Coulochem II, guard cell +1.0 V, screening cell E1 +0.50 V, analytical cell E2 +0.80 V

CHROMATOGRAM

Retention time: 11.5

Internal standard: clarithromycin (10)

Limit of detection: 50 ng/mL

OTHER SUBSTANCES

Extracted: azithromycin

KEY WORDS

pharmacokinetics; plasma; saliva; serum; leucocytes

REFERENCE

Kees, F.; Spangler, S.; Wellenhofer, M. Determination of macrolides in biological matrices by high-performance liquid chromatography with electrochemical detection, *J. Chromatogr. A*, **1998**, 812, 287–293.

SAMPLE

Matrix: blood, saliva, urine

Sample preparation: Plasma. 2 mL Plasma + 20 μ L 750 μ g/mL roxithromycin in MeCN + 5 mL diethyl ether, shake vigorously for 3 min, centrifuge at 900 g at 4° for 5 min. Remove upper layer and evaporate it to dryness under a stream of nitrogen at 45° . Reconstitute residue with 100 μ L MeCN, vortex 5 s, inject 40 μ L aliquot. Urine. 1.5 mL Urine + 100 μ L 750 μ g/mL roxithromycin in saturated K_2HPO_4 + 4 mL diethyl ether, shake vigorously for 3 min, centrifuge at 900 g at 4° for 5 min. Remove upper layer and evaporate it to dryness under a stream of nitrogen at 45° . Reconstitute residue with 100 μ L MeCN, vortex 5 s, inject 40 μ L aliquot. Saliva. 1.5 mL Saliva + 100 μ L 750 μ g/mL roxithromycin in saturated K_2HPO_4 + 4 mL diethyl ether, shake vigorously for 3 min, centrifuge at 900 g at 4° for 15 min. Remove upper layer and evaporate it to dryness under a stream of nitrogen at 45° . Reconstitute residue with 100 μ L MeCN, vortex 5 s, inject 40 μ L aliquot.

HPLC VARIABLES

Column: Nova-Pak C18

Mobile phase: MeCN:MeOH:56 mM sodium acetate buffer 50:4:56, final pH adjusted to 7.0 with glacial acetic acid

Flow rate: 1.1

Injection volume: 40

Detector: E, Waters M460, +0.9 V versus Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 14.7

Internal standard: roxithromycin

Limit of detection: 12.5 ng/mL

OTHER SUBSTANCES

Simultaneous: 4'-acetylerythromycin, 6-O-methylerythromycin, erythromycin base, erythromycin B, erythromycin estolate, erythromycin ethylsuccinate

KEY WORDS

plasma; roxithromycin is IS

REFERENCE

Croteau,D.; Vallée,F.; Bergeron,M.G.; LeBel,M. High-performance liquid chromatographic assay of erythromycin and its esters using electrochemical detection, *J.Chromatogr.*, **1987**, 419, 205–212.

SAMPLE

Matrix: blood, urine

Sample preparation: Dilute urine 1:2 with isotonic NaCl. 200 μ L Plasma or diluted urine + 100 μ L 10 μ g/mL erythromycin in water + 600 μ L pH 9 phosphate buffer + 3 mL dichloromethane, shake for 10 min, centrifuge at 2000 g for 5 min. Remove 2.5 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 50 μ L MeOH, vortex for 10 s, inject a 15 μ L aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: MeCN:MeOH:83 mM ammonium acetate 55:22:23, pH adjusted to 7.5 with acetic acid

Flow rate: 1

Injection volume: 15

Detector: E, ESA Coulochem Model 5100A, Model 5020 guard cell 1.0 V (before injector), Model 5010 dual-electrode cell, screen electrode E1 + 0.7 V, sample electrode E2 +0.9 V, 0.5 μ m ESA carbon filters placed before guard and analytical cells

CHROMATOGRAM

Retention time: 9.8

Internal standard: erythromycin (7.0)

Limit of detection: 100 ng/mL

OTHER SUBSTANCES

Simultaneous: amitriptyline, clomipramine, disopyramide, erythromycin estolate, erythromycin ethylsuccinate, erythromycin stearate, imipramine, josamycin, lidocaine, spiramycin

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Demotes-Mainaird,F.M.; Vinçon,G.A.; Jarry,C.H.; Albin,H.C. Micro-method for the determination of roxithromycin in human plasma and urine by high-performance liquid chromatography using electrochemical detection, *J.Chromatogr.*, **1989**, 490, 115–123.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 15.833

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.


Rifaxacin

Molecular formula: C₁₇H₁₈FN₃O₃S

Molecular weight: 363.41

CAS Registry No.: 101363-10-4

Merck Index: 8448



SAMPLE

Matrix: bile, blood, urine

Sample preparation: Plasma. 500 µL Plasma + 25 µL 350 µg/mL pipemidic acid in water + 50 µL 70% perchloric acid, vortex for 20 s, centrifuge at 14926 g for 5 min, inject a 10 µL aliquot of the supernatant. Urine. 500 µL Urine + 500 µL water + 1 mL 35 µg/mL pipemidic acid in water + 1.5 g Na₂HPO₄:NaH₂PO₄ 50:50, mix, add 8 mL dichloromethane, agitate horizontally at low speed for 40 min, centrifuge at 6590 g for 5 min. Remove 5 mL of the upper organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 1 mL mobile phase, inject a 10 µL aliquot. Alternatively, dilute 1 mL urine to 50 mL with water, inject a 10 µL aliquot. Bile. Centrifuge 500 µL bile at 14926 g for 15 min, filter (0.2 µm) the supernatant, inject a 5 µL aliquot of the supernatant.

HPLC VARIABLES

Guard column: 15 × 3.2 7 µm Polymer RP (Brownlee)

Column: 50 × 4.1 10 µm PRP1 poly(styrene-divinylbenzene) (Hamilton)

Mobile phase: MeCN:0.17% phosphoric acid 12:88, adjust to pH 5.6 with triethylamine then add 5 mL/L THF

Flow rate: 1

Injection volume: 5-10

Detector: F ex 350 em 510 or UV 300

CHROMATOGRAM**Retention time:** 8.9**Internal standard:** pipemidic acid (3.5)**Limit of detection:** 10 ng/mL (F)**Limit of quantitation:** 50 µg/mL (F)

OTHER SUBSTANCES**Extracted:** metabolites

KEY WORDS

plasma; pharmacokinetics; protect from light

REFERENCE

Lombardi,F.; Ardemagni,R.; Colzani,V.; Visconti,M. High-performance liquid chromatographic determination of rufloxacin and its main active metabolite in biological fluids, *J.Chromatogr.*, **1992**, 576, 129–134.

SAMPLE**Matrix:** blood

Sample preparation: Mix 1 mL plasma with 20 µL 100 µg/mL IS in 20 mM NaOH. Add 500 µL 50 mM pH 7.0 phosphate buffer, vortex for 1 min, add 2 mL dichloromethane, vortex for 1 min, shake for 5 min. Centrifuge at 100 g for 10 min, separate the organic layer, repeat the extraction procedure twice, evaporate the combined organic lowers to dryness under reduced pressure. Reconstitute with 200 µL 20 mM NaOH, inject a 20 µL aliquot.

HPLC VARIABLES**Guard column:** 20 × 4.6 10 µm Vydac AXGU**Column:** 250 × 4.6 5 µm Supelcosil LC-SAX**Mobile phase:** MeCN:50 mM pH 7.0 phosphate buffer 10:90**Flow rate:** 1.5**Injection volume:** 20**Detector:** UV 280

CHROMATOGRAM**Retention time:** 7.3**Internal standard:** fenbufen (3.5)**Limit of detection:** 50 ng/mL

OTHER SUBSTANCES**Extracted:** furprofen

KEY WORDS

SPE; plasma

REFERENCE

Carlucci,G.; Mazzeo,P. Simultaneous determination of furprofen and rufloxacin in human plasma by high-performance liquid chromatography, *J.Chromatogr.Sci.*, **1996**, 34, 182–184.

SAMPLE**Matrix:** blood

Sample preparation: Filter 1 mL plasma using a micropartition system (MPS-1, Amicon, MA) while centrifuging at 2000 g for 20 min at 10°, inject an aliquot of the ultrafiltrate..

HPLC VARIABLES**Column:** 250 × 4.6 Spherisorb ODS-2 endcapped

Mobile phase: MeCN:buffer 19:81 containing 5 mM tetrabutylammonium sulfate, adjusted to pH 2.5 with 1 M NaOH (Buffer was 100 mM citric acid containing 200 mM ammonium perchlorate.)

Column temperature: 37**Flow rate:** 1**Detector:** UV 248

CHROMATOGRAM**Retention time:** 4.68**Internal standard:** enrofloxacin (7.53)

KEY WORDS

plasma; ultrafiltrate

REFERENCE

Zlotos,G.; Bucker,A.; Kinzig-Schippers,M.; Sorgel,F.; Holzgrabe,U. Plasma protein binding of gyrase inhibitors, *J.Pharm.Sci.*, **1998**, *87*, 215–220.

SAMPLE**Matrix:** blood

Sample preparation: 1 mL Plasma + 10 μ L 2 mg/mL IS in MeOH, mix, add 2 mL dichloromethane:diethyl ether 80:20, vortex for 15 s, centrifuge at 1500 g for 5 min, remove a 1.7 mL aliquot of the organic phase, repeat the extraction twice more with 2 mL portions of dichloromethane:diethyl ether 80:20. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at room temperature, reconstitute the residue in 200 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES**Guard column:** 20 \times 4.6 40 μ m Pelliguard C18 (Supelco)**Column:** 250 \times 4.6 5 μ m Viosfer C18 (Violet, Rome)**Mobile phase:** MeCN:buffer 15:85 (Buffer was 100 mM phosphate buffer containing 5 mM tetrabutylammonium hydrogen sulfate, adjusted to pH 3 with orthophosphoric acid.)**Flow rate:** 1.3**Injection volume:** 20**Detector:** UV 280

CHROMATOGRAM**Retention time:** 6.2**Internal standard:** 2-[4-(2'-furoyl)phenyl]propionic acid (7.3)**Limit of detection:** 30 ng/mL

OTHER SUBSTANCES**Extracted:** theophylline

KEY WORDS

plasma

REFERENCE

Carlucci,G.; Mazzeo,P.; Palumbo,G. Simultaneous determination of rufloxacin and theophylline by high-performance liquid chromatography in human plasma, *Analyst*, **1995**, *120*, 2493–2495.

SAMPLE**Matrix:** blood, urine

Sample preparation: 1 mL Serum + 1 mL 100 mM pH 7 phosphate buffer + 100 (serum) or 150 (urine) μ L 5 μ g/mL ofloxacin in 500 mM NaOH, mix, add 2.5 mL dichloromethane, shake for 10 min, centrifuge at 1500 g for 10 min, repeat the extraction twice. Combine the organic layers and evaporate them to dryness with nitrogen under vacuum, reconstitute the residue in 200 μ L mobile phase, vortex, inject a 20 μ L aliquot.

HPLC VARIABLES**Guard column:** 20 \times 4.6 10 μ m AXGU anion-exchange (Rainin)**Column:** 250 \times 4.6 10 μ m anion-exchange (Vydac)**Mobile phase:** MeCN:50 mM pH 7 phosphate buffer 20:80**Flow rate:** 1.8**Injection volume:** 20**Detector:** UV 296

CHROMATOGRAM**Retention time:** 7.5**Internal standard:** ofloxacin (5.6)**Limit of detection:** 50 ng/mL (urine), 100 ng/mL (serum)**KEY WORDS**

serum; pharmacokinetics

REFERENCE

Carlucci,G.; Palumbo,G. Analytical procedure for the determination of rufloxacin, a new pyridobenzothiazine, in human serum and urine by high-performance liquid chromatography, *J.Chromatogr.*, **1991**, 564, 346–351.

SAMPLE**Matrix:** urine

Sample preparation: Dilute ten times with water, inject a 20 μ L aliquot. Deconjugate by heating 100 μ L urine, 100 μ L 50000 U/mL β -glucuronidase (Type II from limpets *Patella vulgata*, Sigma), and 800 μ L pH 3.8 KH_2PO_4 buffer at 37° for 16 h.

HPLC VARIABLES**Guard column:** 75 \times 2.1 10 μ m pellicular reversed phase (Chrompack)**Column:** 250 \times 4.6 5 μ m CP Spher 5-ODS (Chrompack)

Mobile phase: Gradient. MeCN:buffer from 4:96 to 26:74 over 37 min, return to initial conditions over 5 min. (Buffer was 6.75 mL orthophosphoric acid and 2 mL diethylamine in 1 L water, pH 3.1.)


Flow rate: 1.5**Injection volume:** 20**Detector:** UV 246**CHROMATOGRAM****Retention time:** 22**Limit of detection:** 500 ng/mL**OTHER SUBSTANCES****Extracted:** metabolites**KEY WORDS**

monkey; pharmacokinetics

REFERENCE

Vree,T.B.; Van den Biggelaar-Martea,M.; Peeters,A.; Imbimbo,B.P. High-performance liquid chromatography and preliminary pharmacokinetics of rufloxacin and its metabolites, N-desmethylrufloxacin and rufloxacin-sulfoxide, in urine of rhesus monkey *Macaca mulatta*, *J.Chromatogr.*, **1992**, 573, 168–172.

Rutin

Molecular formula: $\text{C}_{27}\text{H}_{30}\text{O}_{16}$ **Molecular weight:** 610.53**CAS Registry No.:** 153-18-4**Merck Index:** 8456**SAMPLE****Matrix:** blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 204

CHROMATOGRAM

Retention time: 10.1


KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

Saccharin



Molecular formula: C₇H₅NO₃S

Molecular weight: 183.19

CAS Registry No.: 81-07-2, 6485-34-3 (Ca salt), 6381-91-5 (Ca salt hydrate), 128-44-9 (Na salt), 6155-57-3 (Na salt dihydrate)

Merck Index: 8463

SAMPLE

Matrix: beverage

Sample preparation: Sonicate 25 mL beverage for 15-20 min, filter (0.45 µm) if necessary, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 300 × 3.9 10 µm µBondapak C18

Mobile phase: MeCN:MeOH:water:acetic acid 10:20:70:1

Flow rate: 1.5

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 3.5

OTHER SUBSTANCES

Simultaneous: benzoic acid, hydroquinine, quinine