

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18

**Column:** 250 × 4.6 5 µm Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 204

---

**CHROMATOGRAM**

**Retention time:** 10.1

---

**KEY WORDS**

whole blood


---

**REFERENCE**

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

---

# Saccharin



**Molecular formula:** C<sub>7</sub>H<sub>5</sub>NO<sub>3</sub>S

**Molecular weight:** 183.19

**CAS Registry No.:** 81-07-2, 6485-34-3 (Ca salt), 6381-91-5 (Ca salt hydrate), 128-44-9 (Na salt), 6155-57-3 (Na salt dihydrate)

**Merck Index:** 8463

---

**SAMPLE**

**Matrix:** beverage

**Sample preparation:** Sonicate 25 mL beverage for 15-20 min, filter (0.45 µm) if necessary, inject a 20 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 300 × 3.9 10 µm µBondapak C18

**Mobile phase:** MeCN:MeOH:water:acetic acid 10:20:70:1

**Flow rate:** 1.5

**Injection volume:** 20

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 3.5

---

**OTHER SUBSTANCES**

**Simultaneous:** benzoic acid, hydroquinine, quinine

---

**KEY WORDS**

tonic water; soft drinks

---

**REFERENCE**

Valenti, L.P. Liquid chromatographic determination of quinine, hydroquinine, saccharin, and sodium benzoate in quinine beverages, *J.Assoc.Off.Anal.Chem.*, **1985**, 68, 782–784.

---

**SAMPLE**

**Matrix:** beverages, sweetener

**Sample preparation:** Sweetener. Dissolve 30 mg powdered tabletop sweetener in water and dilute to 25 mL, filter (0.2  $\mu\text{m}$  PTFE). Beverages. Dilute fruit juice 1:10 with water. Degas carbonated beverages in a ultrasonic bath for 5 min, dilute 1:10 with water, filter. Inject a 50  $\mu\text{L}$  aliquot.

---

**HPLC VARIABLES**

**Guard column:** 50  $\times$  4 Dionex IonPak AG4A-SC

**Column:** 250  $\times$  4 Dionex IonPak AS4A-SC

**Mobile phase:** Gradient. A was 1 mM sodium carbonate. B was 12.5 mM sodium carbonate. A: B 100:0 for 4.5 min, from 100:0 to 0:100 over 1 min, maintain at 0:100 for 22.5 min, from 0:100 to 100:0 over 0.1 min

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 190 for 6 min, UV 206 22 min, then UV 190; Conductivity, Dionex ED40 in conductivity mode preceded by a Dionex ASRS-I suppressor (external water mode, 300 mA)

---

**CHROMATOGRAM**

**Retention time:** 18

**Limit of detection:** 19 ng/mL (UV), 260 ng/mL (conductivity)

---

**OTHER SUBSTANCES**

**Simultaneous:** acesulfame, aspartame

---

**REFERENCE**

Chen, Q.-C.; Mou, S.; Liu, K.; Yang, Z.; Ni, Z. Separation and determination of four artificial sweeteners and citric acid by high-performance liquid chromatography, *J.Chromatogr.A*, **1997**, 771, 135–143.

---

**SAMPLE**

**Matrix:** beverages, syrup

**Sample preparation:** Dilute syrup ten fold. Filter (0.45  $\mu\text{m}$ ) beverages and diluted syrup, inject a 10–20  $\mu\text{L}$  aliquot of the filtrate.

---

**HPLC VARIABLES**

**Column:** 300  $\times$  3.9 10  $\mu\text{m}$   $\mu$ Bondapak C18

**Mobile phase:** MeOH:acetic acid:water 20:5:75

**Flow rate:** 2

**Injection volume:** 10–20

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 2.5

**Limit of detection:** 30 ng

---

**OTHER SUBSTANCES**

**Simultaneous:** acesulfame, benzoic acid, caffeine, dulcin, p-hydroxybenzoic acid, vanillin

---

**REFERENCE**

Veerabhadra Rao, M.; Narayan, M.S.; Kapur, O.; Sastry, C.S. Reverse phase liquid chromatographic determination of some food additives, *J.Assoc.Off.Anal.Chem.*, **1987**, 70, 578–582.

---

**SAMPLE**

**Matrix:** blood, feces, feed, milk, urine

**Sample preparation:** Condition a 3 mL quaternary amine SPE cartridge (J.T. Baker) with 3 mL MeOH, 3 mL concentrated ammonium hydroxide, and 3 mL 200 mM HCl. Feed, feces. Air dry feces at 50° for 17 h. 350 mg Feed or dried feces + 350  $\mu$ L 100 mM NaOH, mix, add 2.5 mL 100 mM NaOH, mix, centrifuge at 500 g for 10 min, repeat extraction twice with 2.5 mL aliquots and once with a 2 mL aliquot of 100 mM NaOH. Combine the supernatants and make up to 10 mL with water. Remove a 500  $\mu$ L aliquot and add it to the SPE cartridge, wash with two 1 mL aliquots of MeOH, wash with 2 mL water (add 3  $\mu$ L 2-octanol to top of water as an anti-foaming agent), elute with two 2 mL aliquots of 200 mM pH 8.8  $K_2HPO_4$ , add 80  $\mu$ L phosphoric acid to the eluate, make up to 5 mL with water, inject an aliquot. Urine. 1 mL Urine + 8.5 mL 100 mM NaOH, make up to 10 mL with water. Remove a 500  $\mu$ L aliquot and add it to the SPE cartridge, wash with two 1 mL aliquots of MeOH, wash with 2 mL water (add 3  $\mu$ L 2-octanol to top of water as an anti-foaming agent), elute with two 2 mL aliquots of 200 mM pH 8.8  $K_2HPO_4$ , add 80  $\mu$ L phosphoric acid to the eluate, make up to 5 mL with water, inject an aliquot. Serum. 500  $\mu$ L Serum + 4.25 mL 100 mM NaOH, make up to 5 mL with water. Remove a 500  $\mu$ L aliquot and add it to the SPE cartridge, wash with two 1 mL aliquots of MeOH, wash with 2 mL water (add 3  $\mu$ L 2-octanol to top of water as an anti-foaming agent), elute with two 2 mL aliquots of 200 mM pH 8.8  $K_2HPO_4$ , add 80  $\mu$ L phosphoric acid to the eluate, make up to 5 mL with water, inject an aliquot. Milk. Add 200  $\mu$ L milk to the SPE cartridge, wash with two 1 mL aliquots of MeOH, wash with 2 mL water (add 3  $\mu$ L 2-octanol to top of water as an anti-foaming agent), elute with two 2 mL aliquots of 200 mM pH 8.8  $K_2HPO_4$ , add 80  $\mu$ L phosphoric acid to the eluate, make up to 5 mL with water, inject an aliquot.

---

#### HPLC VARIABLES

**Column:** 150  $\times$  3.9 4  $\mu$ m Novapak C18

**Mobile phase:** MeCN:10 mM pH 2.2  $KH_2PO_4$ , 16:84

**Column temperature:** 30

**Flow rate:** 1

**Detector:** UV 254 or F ex 232 em 432

---

#### CHROMATOGRAM

**Retention time:** 1.9

**Limit of quantitation:** 2.45 pmole (F), 245 pmole (UV)

---

#### KEY WORDS

rat; SPE

---

#### REFERENCE

Tibbels,T.S.; Smith,R.A.; Cohen,S.M. Determination of saccharin in diet and biological materials, *J.Chromatogr.*, 1988, 441, 448-453.

---

#### SAMPLE

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50  $\mu$ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood)  $\mu$ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

#### HPLC VARIABLES

**Guard column:** 20 mm long Symmetry C18

**Column:** 250  $\times$  4.6 5  $\mu$ m Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 200.5

---

**CHROMATOGRAM****Retention time:** 5.907

---

**KEY WORDS**whole blood

---

**REFERENCE**

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

---

---

**SAMPLE****Matrix:** food

**Sample preparation:** 2 g Soy sauce or 1 g sugared fruit or roast beef + 10 g NaCl, make up to 50 mL with acetone, swirl vigorously, let stand for 30 min, filter, wash the solid with acetone, make up filtrate to 50 mL with acetone, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES****Column:** 250  $\times$  4.6 5  $\mu$ m monomeric C18 (Shoko, Kyoto)

**Mobile phase:** MeCN:50 mM pH 4.5  $\alpha$ -hydroxyisobutyric acid in water 22:34 containing 2.5 mM hexadecyltrimethylammonium bromide

**Flow rate:** 1**Injection volume:** 20**Detector:** UV 233

---

**CHROMATOGRAM****Retention time:** 33

---

**OTHER SUBSTANCES**

**Simultaneous:** acesulfame-K, benzoic acid, 3-t-butyl-4-hydroxyanisole, butyl p-hydroxybenzoate, t-butylhydroxyquinone, dulcin, ethyl p-hydroxybenzoate, isobutyl p-hydroxybenzoate, isopropyl p-hydroxybenzoate, methyl p-hydroxybenzoate, sodium dehydroacetate, sorbic acid

---

**KEY WORDS**soy sauce; roast beef; sugared fruit

---

**REFERENCE**

Chen,B.H.; Fu,S.C. Comparison of extraction methods and column types for the determination of additives by liquid chromatography, *J.Liq.Chromatogr.Rel.Technol.*, **1996**, 19, 625–643.

---

---

**SAMPLE****Matrix:** formulations

**Sample preparation:** Pulverize tablets and weigh out 1 g, add 1 mL formic acid, add 25 mL MeOH, shake mechanically for 10 min, make up to 50 mL with methanol. Remove 10 mL and centrifuge. 5 mL Supernatant + 5 mL 0.0025% p-hydroxybenzoic acid in MeOH:water 20:80, make up to 25 mL with water, inject an aliquot. (Analyze within 1 h.)

---

**HPLC VARIABLES****Column:** 250  $\times$  4.6 LiChrosorb RP8

**Mobile phase:** MeOH:200 mM pH 3.5 phosphate buffer:water 20:10:70

**Flow rate:** 1**Injection volume:** 10**Detector:** UV 254

---

**CHROMATOGRAM****Retention time:** 5**Internal standard:** p-hydroxybenzoic acid (18)

---

**OTHER SUBSTANCES**

**Simultaneous:** aspirin, p-aminophenol, 3-O-acetylascorbic acid, 2-O-acetylascorbic acid, Ascorbic acid, acetaminophen, O-acetyl-p-aminophenol, salicylic acid (UV 280), diacetyl-p-aminophenol (UV 280)

---

**KEY WORDS**

tablets

---

**REFERENCE**

Thomis,R.; Roets,E.; Hoogmartens,J. Analysis of tablets containing aspirin, acetaminophen, and ascorbic acid by high-performance liquid chromatography, *J.Pharm.Sci.*, **1984**, 73, 1830–1833.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Grind tablet, add 80 mL water, shake at 240 oscillations/min for 30 min, make up to 100 mL with water, mix well, centrifuge at 1500 rpm, dilute supernatant with water, filter (0.45  $\mu$ m), inject an aliquot of the filtrate.

---

**HPLC VARIABLES**

**Column:** 150  $\times$  4.6 5  $\mu$ m Hypersil ODS

**Mobile phase:** MeCN:buffer 3:97 (Buffer was 250 mM NaH<sub>2</sub>PO<sub>4</sub> adjusted to pH 2.5 with phosphoric acid.) (Flush column with MeCN:water 5:95 at the end of each day.)

**Flow rate:** 2

**Injection volume:** 20

**Detector:** UV 205

---

**CHROMATOGRAM**

**Retention time:** 7

---

**OTHER SUBSTANCES**

**Simultaneous:** pantothenic acid, pantoyllactone, panthenol

---

**KEY WORDS**

tablets

---

**REFERENCE**

Timmons,J.A.; Meyer,J.C.; Steible,D.J.; Assenza,S.P. Reverse phase liquid chromatographic assay for calcium pantothenate in multivitamin preparations and raw materials, *J.Assoc.Off.Anal.Chem.*, **1987**, 70, 510–513.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Dilute 1 mL syrup to 50 mL with mobile phase, filter (0.45  $\mu$ m), inject 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 Zorbax CN

**Mobile phase:** MeCN:water:formic acid:methanesulfonic acid 500:500:1:1, pH adjusted to 3.5 with 10% NaOH

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 290

---

**CHROMATOGRAM**

**Retention time:** 3

---

**OTHER SUBSTANCES**

**Simultaneous:** guaifenesin, dextromethorphan, benzoic acid

---

**KEY WORDS**

syrup

---

**REFERENCE**

Chen, T.M.; Pacifico, J.R.; Daly, R.E. High-pressure liquid chromatographic assay of dextromethorphan hydrobromide, guaifenesin, and sodium benzoate in an expectorant syrup, *J. Chromatogr. Sci.*, **1988**, *26*, 636-639.

---

**SAMPLE**

**Matrix:** formulations

---

**HPLC VARIABLES**

**Column:** 100 × 4.3 µm Hypersil BDS-C18

**Mobile phase:** Gradient. MeCN:water adjusted to pH 2.1 from 0.3:99.7 to 25:75 over 11 min

**Flow rate:** 0.5

**Detector:** UV 220

---

**CHROMATOGRAM**

**Retention time:** 7.5

---

**OTHER SUBSTANCES**

**Simultaneous:** biotin, caffeine, citric acid, folic acid, niacinamide, niacin, pantothenic acid, riboflavin, thiamine, pyridoxine, vitamin B12, ascorbic acid

---

**KEY WORDS**

tablets

---

**REFERENCE**

*Hewlett Packard Leaflet 12-5091-7351 EUS, 1993.*

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 Zorbax RX

**Mobile phase:** Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

**Column temperature:** 30

**Flow rate:** 2

**Detector:** UV 210

---

**OTHER SUBSTANCES**

**Also analyzed:** acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitrityline, amobarbital, amoxapine, amphetamine, amyllocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clonbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenopropfen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiacol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarbostyryl, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclufenamic acid, medazepam,

mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methylidopa, methylidopamine, methylphenidate, methylprednisolone, methyltestosterone, methypyrrolon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidone, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethiodole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlylcypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

## REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233–242.

## SAMPLE

**Matrix:** yogurt

**Sample preparation:** 500 mg Homogenized yogurt + 7 mL buffer, sonicate for 2 min, shake mechanically for 20 min, centrifuge at 2000 rpm for 10 min, repeat extraction twice more. Combine the supernatants and make up to 25 mL with buffer. Remove a 5 mL extract and add it to 5 mL 10 mM tri-n-octylamine in chloroform, shake for 20 min, centrifuge at 2000 rpm for 10 min. Remove 2.5 mL of the organic phase and add it to 2.5 mL 100 mM sodium perchlorate in water, extract, centrifuge, inject an aliquot of the aqueous phase. (Buffer was 24.650 g  $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$  and 1.260 g  $\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$  made up to 2 L with water, pH 5.5.)

## HPLC VARIABLES

**Column:** 250 × 4 10  $\mu\text{m}$  RP-18 (Merck)

**Mobile phase:** MeOH:buffer 40:60 (Buffer was 900  $\mu\text{L}$  1 M phosphoric acid and 27.598 g  $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$  made up to 2 L with water, pH 4.5.)

**Injection volume:** 100

**Detector:** UV 270 for 4 min then UV 240

## CHROMATOGRAM

**Retention time:** 2

**Limit of detection:** 20  $\mu\text{g/g}$


## OTHER SUBSTANCES

**Simultaneous:** benzoic acid, sorbic acid

## REFERENCE

Puttemans,M.L.; Branders,C.; Dryon,L.; Massart,D.L. Extraction of organic acids by ion-pair formation with tri-n-octylamine. Part 6. Determination of sorbic acid, benzoic acid, and saccharin in yogurt, *J.Assoc.Off.Anal.Chem.*, **1985**, *68*, 80–82.

# Salicylamide



**Molecular formula:** C<sub>7</sub>H<sub>7</sub>NO<sub>2</sub>

**Molecular weight:** 137.14

**CAS Registry No.:** 65-45-2

**Merck Index:** 8480

**Lednicer No.:** 1 109

## SAMPLE

**Matrix:** bile, blood, urine

**Sample preparation:** Plasma. 200  $\mu$ L Plasma + 100  $\mu$ L 1 M perchloric acid, mix, centrifuge at 12000 g for 3 min, inject a 100  $\mu$ L aliquot of the supernatant. Bile. Dilute bile 1:10 with water, inject a 20  $\mu$ L aliquot. Urine. Filter (0.45  $\mu$ m) urine, dilute with water 1:4, inject a 30  $\mu$ L aliquot.

## HPLC VARIABLES

**Guard column:** 40  $\mu$ m  $\mu$ Bondapak C18

**Column:** 250  $\times$  4.6 10  $\mu$ m  $\mu$ Bondapak C18

**Mobile phase:** 85 mM KH<sub>2</sub>PO<sub>4</sub>, pH adjusted to 3.35 with glacial acetic acid

**Flow rate:** 0.8 for 34 min, to 1.8 over 5 min, maintain at 1.8 for 11 min

**Injection volume:** 20-100

**Detector:** UV 313

## CHROMATOGRAM

**Retention time:** 48

**Limit of detection:** 1.09 nmole

## OTHER SUBSTANCES

**Extracted:** metabolites, conjugates

## KEY WORDS

rat; plasma; whole blood

## REFERENCE

Xu,X.; Pang,K.S. High-performance liquid chromatographic method for the quantitation of salicylamide and its metabolites in biological fluids, *J.Chromatogr.*, **1987**, *420*, 313-327.

## SAMPLE

**Matrix:** blood

**Sample preparation:** 3 mL Whole blood + 6 mL ethyl acetate + 100  $\mu$ L 10  $\mu$ g/mL p-nitrophenol, shake, centrifuge at 3000 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 100  $\mu$ L MeOH, inject a 20  $\mu$ L aliquot.

## HPLC VARIABLES

**Guard column:** 40  $\mu$ m  $\mu$ Bondapak C18

**Column:** 250  $\times$  4.6 10  $\mu$ m  $\mu$ Bondapak C18

**Mobile phase:** MeOH:0.5% acetic acid 40:60

**Flow rate:** 0.6

**Injection volume:** 20

**Detector:** UV 313

## CHROMATOGRAM

**Internal standard:** p-nitrophenol

**Limit of detection:** 83 ng/mL

## OTHER SUBSTANCES

**Extracted:** metabolites



---

**KEY WORDS**

rat; whole blood

---

**REFERENCE**

Xu,X.; Pang,K.S. High-performance liquid chromatographic method for the quantitation of salicylamide and its metabolites in biological fluids, *J.Chromatogr.*, **1987**, 420, 313–327.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Whole blood. Mix whole blood with 2 volumes cold MeOH, centrifuge, inject an aliquot of the supernatant. Urine. Dilute 24 h urine with 5 mL water, centrifuge, inject an aliquot of the supernatant.

---

**HPLC VARIABLES**

**Guard column:** 25 × 3.9 37-50 µm Bondapak C18/Corasil

**Column:** 300 × 3.9 10 µm µBondapak C18

**Mobile phase:** THF:1% acetic acid 5:95 containing 1.5 mM tetrabutylammonium hydroxide

**Flow rate:** 1.6

**Injection volume:** 20

**Detector:** UV 240

---

**CHROMATOGRAM**

**Retention time:** 20

---

**OTHER SUBSTANCES**

**Extracted:** metabolites, conjugates

---

**KEY WORDS**

mouse; whole blood

---

**REFERENCE**

Howell,S.R.; Kotkoskie,L.A.; Dills,R.L.; Klaassen,C.D. 3-Hydroxylation of salicylamide in mice, *J.Pharm.Sci.*, **1988**, 77, 309–313.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Dissolve compounds in MeOH, inject a 1 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 150 × 1 3 µm Hitachi-Gel 3011 porous polymer (Hitachi)

**Mobile phase:** MeOH:ammonia 99:1

**Flow rate:** 0.03

**Injection volume:** 1

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 3.16

---

**OTHER SUBSTANCES**

**Also analyzed:** acetaminophen, caffeine, bucetin (3-hydroxy-p-butyrophenetidine), phenacetin, dipyrone (sulpyrin), mefenamic acid, aspirin, salicylic acid, ethenzamide (o-ethoxybenzamide), theobromine, theophylline

---

**KEY WORDS**

semi-micro; porous polymer

---

**REFERENCE**

Matsushima,Y.; Nagata,Y.; Niyomura,M.; Takakusagi,K.; Takai,N. Analysis of antipyretics by semimicro liquid chromatography, *J.Chromatogr.*, **1985**, 332, 269–273.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Dissolve in MeOH:water 1:1 at a concentration of 50 µg/mL, inject a 10 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 300 × 3.9 10 µm µBondapak C18

**Mobile phase:** MeOH:acetic acid:triethylamine:water 30:1.5:0.5:68

**Flow rate:** 1.5

**Injection volume:** 10

**Detector:** UV

---

**CHROMATOGRAM**

**Retention time:** k' 2.19

---

**REFERENCE**

Roos,R.W.; Lau-Cam,C.A. General reversed-phase high-performance liquid chromatographic method for the separation of drugs using triethylamine as a competing base, *J.Chromatogr.*, **1986**, 370, 403–418.

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Guard column:** 30 × 3.2 7 µm SI 100 ODS (not commercially available)

**Column:** 150 × 3.2 7 µm SI 100 ODS (not commercially available)

**Mobile phase:** MeCN:buffer 31.2:68.8 (Buffer was 6.66 g KH<sub>2</sub>PO<sub>4</sub> and 4.8 g 85% phosphoric acid in 1 L water, pH 2.3.)

**Flow rate:** 0.5-1

**Detector:** UV 230, 296

---

**CHROMATOGRAM**

**Retention time:** 2.4

**Internal standard:** 5-(4-methylphenyl)-5-phenylhydantoin (7.3)

---

**OTHER SUBSTANCES**

**Also analyzed:** aspirin, caffeine, carbamazepine, chlordiazepoxide, chlorprothixene, clonazepam, diazepam, doxylamine, ethosuximide, furosemide, haloperidol, hydrochlorothiazide, methocarbamol, methotrimeprazine, nicotine, oxazepam, procaine, promazine, propafenone, propranolol, temazepam, tetracaine, thiopental, triamterene, verapamil, zolpidem, zopiclone

---

**REFERENCE**

Below,E.; Burrmann,M. Application of HPLC equipment with rapid scan detection to the identification of drugs in toxicological analysis, *J.Liq.Chromatogr.*, **1994**, 17, 4131–4144.

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 Zorbax RX

**Mobile phase:** Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

**Column temperature:** 30

**Flow rate:** 2

**Detector:** UV 210

---

**OTHER SUBSTANCES**

**Also analyzed:** acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbitol, benzocaine, benzoic acid, benzotropine, benzphet-

amine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clen-buterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamor-phine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, dil-tiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, dox-apram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenopropen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-stilbene, imipramine, indomethacin, isocarboxystyl, isocarboxazid, isoniazid, isoproterenol, isox-suprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, meth-apyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, meth-yltestosterone, methypylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nyldrin, oxazepam, oxycodone, ox-ymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendi-metrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phenter-mine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, predni-solone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidone, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethi-dole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasox-izole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, the-baine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tol-metin, tranilcypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapa-mil, vincamine, warfarin, yohimbine, zoxazolamine

## REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, **1994**, *18*, 233-242.

## SAMPLE

**Matrix:** solutions

**Sample preparation:** Prepare a 50 µg/mL solution in the mobile phase, inject a 10 µL aliquot.

## HPLC VARIABLES

**Column:** 250 × 4.6 7 µm Lichrosorb RP 18

**Mobile phase:** MeOH:water 45:55 containing 1% acetic acid

**Flow rate:** 1

**Injection volume:** 10

**Detector:** UV 230

## CHROMATOGRAM

**Retention time:** 4.20

## OTHER SUBSTANCES

**Simultaneous:** acetaminophen, aspirin, phenacetin, salicylic acid

---

**REFERENCE**

Nivaud-Guernet,E.; Guernet,M.; Ivanovic,D.; Medenica,M. Effect of eluent pH on the ionic and molecular forms of the non-steroidal anti-inflammatory agents in reversed-phase high-performance liquid chromatography, *J.Liq.Chromatogr.*, **1994**, *17*, 2343–2357.

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 8 μm Unisphere-PBD (polybutadiene on alumina) (Biotage, Charlottesville, VA)

**Mobile phase:** MeCN:water 20:80

**Flow rate:** 1

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 3.5

---

**OTHER SUBSTANCES**

**Simultaneous:** acetaminophen, aspirin, phenacetin

---

**REFERENCE**

Jedrejewski,P.T.; Taylor,L.T. Comparison of silica-, alumina-, and polymer-based stationary phases for reversed-phase liquid chromatography, *J.Chromatogr.Sci.*, **1995**, *33*, 438–445.

---

**SAMPLE**

**Matrix:** tissue

**Sample preparation:** Homogenize (Cole-Parmer Model 4420) one hemisphere with 4 mL 100 mM pH 5.0 sodium phosphate buffer and p-nitrophenol at 1000 rpm for 2 min, add 5 mL ethyl acetate, vortex for 2 min, centrifuge. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 200 μL MeOH, inject a 10 μL aliquot.

---

**HPLC VARIABLES**

**Column:** 5 μm Partisil-5 ODS-3 RAC

**Mobile phase:** MeOH:25 mM KH<sub>2</sub>PO<sub>4</sub>:acetic acid 30:69:1

**Flow rate:** 2

**Injection volume:** 10

**Detector:** UV 313

---

**CHROMATOGRAM**

**Retention time:** 7.5

**Internal standard:** p-nitrophenol (12.5)

---

**OTHER SUBSTANCES**

**Extracted:** gentisamide

---

**KEY WORDS**


rat; brain

---

**REFERENCE**

Morris,M.E.; Levy,G. Pharmacodynamics of the hypnotic effect of salicylamide in rats, *J.Pharm.Sci.*, **1985**, *74*, 599–602.

# Salicylic acid



**Molecular formula:**  $C_7H_6O_3$

**Molecular weight:** 138.12

**CAS Registry No.:** 69-72-7, 54-21-7 (Na salt)

**Merck Index:** 8484

## SAMPLE

**Matrix:** blood

**Sample preparation:** Mix 100  $\mu$ L plasma with 100  $\mu$ L 2.5  $\mu$ M IS, add 200  $\mu$ L EtOH, vortex for 2 min, centrifuge at 1600 g for 15 min, dilute 50  $\mu$ L of the supernatant with 950  $\mu$ L mobile phase, filter (0.45  $\mu$ m), inject a 50  $\mu$ L aliquot of the filtrate.

## HPLC VARIABLES

**Guard column:** 30  $\times$  4.6 10  $\mu$ m Spheri-10 RP18 (Alltech)

**Column:** 150  $\times$  4.6 5  $\mu$ m Ultrasphere ODS

**Mobile phase:** MeOH:buffer 15:85 (Buffer was 30 mM sodium citrate and 30 mM sodium acetate, pH 5.45.)

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 295

## CHROMATOGRAM

**Retention time:** 7.2

**Internal standard:** 2,6-dihydroxybenzoic acid (11.4)

**Limit of detection:** 3.89  $\mu$ M

## OTHER SUBSTANCES

**Extracted:** metabolites

## KEY WORDS

plasma

## REFERENCE

Coudray,C.; Mangournet,C.; Bouhadjeb,S.; Faure,H.; Favier,A. Rapid high-performance liquid chromatographic assay for salicylic acid in plasma without solvent extraction, *J.Chromatogr.Sci.*, **1996**, 34, 166–173.

## SAMPLE

**Matrix:** blood

**Sample preparation:** Add 250  $\mu$ L 200 mM orthophosphoric acid to 250  $\mu$ L chilled plasma within 10 min of centrifuging (if fresh plasma) or within 10 min of thawing (if frozen plasma), vortex for 20 s, centrifuge at 5800 g for 3 min. Inject a 200  $\mu$ L aliquot onto column A and elute to waste with mobile phase A, after 2 min backflush the contents of column A onto column B with mobile phase B, elute with mobile phase B, monitor the effluent from column B.

## HPLC VARIABLES

**Column:** A 10  $\times$  4.3 30  $\mu$ m Hypersil C18 PEEK cartridge (Shandon, England); B 10  $\times$  4 30  $\mu$ m Hypersil C8 + 250  $\times$  4.6 5  $\mu$ m Nucleosil C8

**Mobile phase:** A Water:orthophosphoric acid 1000:1, pH 2.5; B MeCN:MeOH:water:orthophosphoric acid 150:200:650:1 (pH 2.6)

**Flow rate:** 1

**Injection volume:** 200

**Detector:** UV 225

## CHROMATOGRAM

**Retention time:** 60

**Limit of detection:** 40 ng/mL

**Limit of quantitation:** 100 ng/mL

---

**OTHER SUBSTANCES****Extracted:** aspirin**Noninterfering:** barbital, butobarbital, caffeine, 8-chlorotheophylline, clonazepam, cocaine, diazepam, flurazepam, furosemide, hydralazine, imipramine, nitrazepam, phenytoin, pindolol, propranolol, quinidine, theophylline**Interfering:** xylazine, prazosin

---

**KEY WORDS**column-switching; plasma

---

**REFERENCE**McMahon, G.P.; Kelly, M.T. Determination of aspirin and salicylic acid in human plasma by column-switching liquid chromatography using on-line solid-phase extraction, *Anal. Chem.*, **1998**, *70*, 409–414.

---

**SAMPLE****Matrix:** blood**Sample preparation:** Mix plasma with an equal volume of MeCN, vortex for 30 s, centrifuge at 900 g for 5 min, inject a 100  $\mu$ L aliquot of the supernatant.

---

**HPLC VARIABLES****Column:** 300  $\times$  3.9 10  $\mu$ m  $\mu$ Bondapak C18**Mobile phase:** MeOH:acetic acid:water 22:5:73**Flow rate:** 2.6**Injection volume:** 100**Detector:** UV 280

---

**CHROMATOGRAM****Retention time:** 7.0**Limit of quantitation:** 5000 ng/mL

---

**OTHER SUBSTANCES****Extracted:** aspirin**Noninterfering:** acetaminophen, albuterol, aminophylline, amitriptyline, amoxicillin, ampicillin, amobarbital, beclomethasone, carbamazepine, carbenicillin, chlordiazepoxide, cimetidine, clonazepam, cyproheptadine, debrisoquine, dextropropoxyphene, diazepam, digoxin, dihydroxyanthraquinone, ergotamine, ethosuximide, fluphenazine, furosemide, gentamicin, gentisic acid, guaifenesin, haloperidol, heparin, hydrocortisone, indomethacin, methdilazine, methyclothiazide, methylphenobarbitone, methysergide, metoclopramide, naproxen, nitrazepam, nystatin, penicillin, pentobarbitone, phenytoin, phenytoin, pizotifen, prazosin, prednisone, prochlorperazine, propranolol, spironolactone, sulfamethoxazole, theophylline, trifluoperazine, trimethoprim, valproic acid

---

**KEY WORDS**plasma

---

**REFERENCE**Cham, B.E.; Ross-Lee, L.; Bochner, F.; Imhoff, D.M. Measurement and pharmacokinetics of acetylsalicylic acid by a novel high performance liquid chromatographic assay, *Ther. Drug Monit.*, **1980**, *2*, 365–372.

---

**SAMPLE****Matrix:** blood**Sample preparation:** 200  $\mu$ L Serum + 200  $\mu$ L 50  $\mu$ g/mL hexobarbital in MeCN + 25  $\mu$ L glacial acetic acid, vortex for 10 s, centrifuge for 1 min, inject a 30–100  $\mu$ L aliquot of the supernatant.

---

**HPLC VARIABLES****Column:**  $\mu$ Bondapak C18**Mobile phase:** Gradient. MeCN:7.5 g/L  $\text{NaH}_2\text{PO}_4$  adjusted to pH 3.2 with phosphoric acid 5:95 to 22:78 over 24 min, to 45:55 over 10 min, maintain at 45:55 for 5 min. Re-equilibrate with 5:95 for 5 min.**Column temperature:** 50**Flow rate:** 3

**Injection volume:** 30-100

**Detector:** UV 210

---

#### CHROMATOGRAM

**Retention time:** 8.0

**Internal standard:** hexobarbital (20.6)

**Limit of detection:** 200-2000 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** acetaminophen, amobarbital, butabarbital, butalbital, chlordiazepoxide, diazepam, ethchlorvynol, flurazepam, glutethimide, methaqualone, methypylon, nitrazepam, pentobarbital, phenobarbital, phenytoin, primidone, secobarbital, theophylline

**Simultaneous:** amitriptyline, caffeine, clomipramine, codeine, desipramine, ethotoin, imipramine, lidocaine, mesantoin, methsuximide, nirvanol, nortriptyline, oxazepam, procainamide, phenylpropanolamine, propranolol, quinidine

---

#### KEY WORDS

serum

---

#### REFERENCE

Kabra,P.M.; Stafford,B.E.; Marton,L.J. Rapid method for screening toxic drugs in serum with liquid chromatography, *J.Anal.Toxicol.*, **1981**, 5, 177-182.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** 50  $\mu$ L Serum + 50  $\mu$ L 15  $\mu$ g/mL  $\beta$ -hydroxyethyltheophylline in MeCN, mix for 30 s, centrifuge at 13000 g for 5 min. Inject the supernatant (about 20  $\mu$ L).

---

#### HPLC VARIABLES

**Column:**  $\mu$ Bondapak C18

**Mobile phase:** MeCN:buffer 9.75:90.25 (Buffer was 100 mM  $\text{KH}_2\text{PO}_4$  adjusted to pH 4.0 with phosphoric acid.) (At the end of each day clean with water for 20 min and MeOH for 30 min.)

**Flow rate:** 2

**Injection volume:** 20

**Detector:** UV 254

---

#### CHROMATOGRAM

**Retention time:** 6.5

**Internal standard:**  $\beta$ -hydroxyethyltheophylline (5.8)

**Limit of detection:** 500 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** aspirin

**Simultaneous:** dyphylline, theophylline, caffeine, acetaminophen, procainamide, N-acetylprocainamide

**Noninterfering:** benzoic acid

---

#### KEY WORDS

serum

---

#### REFERENCE

Ou,C.-N.; Frawley,V.L. Theophylline, dyphylline, caffeine, acetaminophen, salicylate, acetylsalicylate, procainamide, and N-acetylprocainamide determined in serum with a single liquid-chromatographic assay, *Clin.Chem.*, **1982**, 28, 2157-2160.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** Plasma. 1-2 mL Plasma + 600-1000 mg potassium carbonate plasma + 1 mL MeOH, shake for 1 min, centrifuge. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 500  $\mu$ L mobile phase, inject a 25-50  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 610 × 2 25-50 µm Bondapak AX/Corasil anion exchange

**Mobile phase:** MeOH:pH 5.6 phosphate buffer 15:85

**Flow rate:** 0.5

**Injection volume:** 10-50

**Detector:** UV 254, 280

---

**CHROMATOGRAM**

**Retention time:** 11.2 (at UV 280)

**Internal standard:** sodium salicylate

---

**OTHER SUBSTANCES**

**Extracted:** dipyrone (at UV 254)

---

**KEY WORDS**

plasma; sodium salicylate is IS

---

**REFERENCE**

Asmardi,G.; Jamali,F. High-performance liquid chromatography of dipyrone and its active metabolite in biological fluids, *J.Chromatogr.*, **1983**, 277, 183-189.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 500 µL Whole blood + 1 mL 1 M HCl + 200 µL 1 mg/mL β-hydroxyethyl-theophylline in MeOH + 10 mL dichloromethane, rotate for 10 min, centrifuge for 5 min. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 200 µL MeOH, inject a 20 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 150 × 4.6 5 µm ODS (Altex)

**Mobile phase:** MeOH:water:glacial acetic acid 40:66:1

**Flow rate:** 1.5

**Injection volume:** 20

**Detector:** UV 250

---

**CHROMATOGRAM**

**Retention time:** 6.4

**Internal standard:** β-hydroxyethyltheophylline (2)

**Limit of quantitation:** 10 µg/mL

---

**OTHER SUBSTANCES**

**Simultaneous:** methyl salicylate, theophylline

**Noninterfering:** acetaminophen, amobarbital, butalbital, caffeine, carbamazepine, glutethimide, ibuprofen, indomethacin, meprobamate, methaqualone, pentobarbital, phenobarbital, phenylbutazone, phenytoin, primidone, secobarbital

---

**KEY WORDS**

whole blood

---

**REFERENCE**

Levine,B.; Caplan,Y.H. Liquid chromatographic determination of salicylate and methyl salicylate in blood and application to a postmortem case, *J.Anal.Toxicol.*, **1984**, 8, 239-241.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 0.2 mL 1 M HCl + 10 mL diethyl ether, gently mix for 10 min, centrifuge at 1500 rpm for 4 min. Remove the organic phase, evaporate it to dryness at 0° under a stream of nitrogen, add 200 µL mobile phase, vortex 90 s, inject a 5-100 µL aliquot.

---

**HPLC VARIABLES**

**Guard column:** 23 × 3.9 µBondapak C18/Porasil B



**Column:** 300 × 3.9 10 µm µBondapak C18

**Mobile phase:** MeOH:water:1-butanol:orthophosphoric acid 270:720:10:0.13

**Column temperature:** 27

**Flow rate:** 1.8

**Injection volume:** 5-100

**Detector:** UV 234

---

#### CHROMATOGRAM

**Retention time:** 8.0

**Internal standard:** m-anisic acid (9.6)

**Limit of quantitation:** 500 ng/mL

---

#### OTHER SUBSTANCES

**Simultaneous:** aspirin

**Noninterfering:** acetaminophen, albuterol, aminophylline, amitriptyline, atenolol, beclomethasone, bromazepam, caffeine, carbamazepine, chloral hydrate, chlordiazepoxide, cimetidine, clonazepam, codeine, desipramine, dexamethasone, dextropropoxyphene, diazepam, dicyclomine, digoxin, disopyramide, doxycycline, ergotamine, ethosuximide, furosemide, gentisic acid, haloperidol, hydrocortisone, imipramine, indomethacin, levodopa, lignocaine, lithium carbonate, meperidine, methdilazine, methylphenobarbitone, methylprednisolone, methysergide, metoclopramide, metoprolol, mexiletine, midazolam, naphthoxyacetic acid, nitrazepam, nitroglycerin, nortriptyline, oxazepam, oxpranolol, pentobarbitone, pethidine, phenytoin, prednisolone, prednisone, primidone, procainamide, prochlorperazine, propranolol, quinidine, salicylic acid, spironolactone, sulfamethoxazole, theophylline, trimethoprim, valproic acid, verapamil, warfarin

**Interfering:** methclothiazide

---

#### KEY WORDS

plasma; pharmacokinetics

---

#### REFERENCE

Brandon,R.A.; Eadie,M.J.; Smith,M.T. A sensitive liquid chromatographic assay for plasma aspirin and salicylate concentrations after low doses of aspirin, *Ther.Drug Monit.*, **1985**, 7, 216-221.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** Activate a 1 mL Bond-Elut C8 SPE cartridge with 2 mL MeOH then 1 mL 10 mM HCl, do not allow it to dry completely. Sonicate 1 mL whole blood for 20-30 min then apply to cartridge. Wash with 100 µL water, elute with three 500 µL portions of MeOH: MeCN:1% aqueous ammonium hydroxide 50:20:30, combine eluents and evaporate to dryness under a stream of nitrogen at 40°. Redissolve in 1 mL MeOH, inject a 20 µL aliquot.

---

#### HPLC VARIABLES

**Column:** 250 × 4.5 5 µm Spherisorb ODS

**Mobile phase:** MeCN:MeOH:buffer 35:13:52 (Buffer was water adjusted to pH 3.2 with orthophosphoric acid)

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 250

---

#### CHROMATOGRAM

**Retention time:** 3.5

---

#### OTHER SUBSTANCES

**Simultaneous:** ketoprofen, acetaminophen, naproxen, fenoprofen, ibuprofen, indomethacin

---

#### KEY WORDS

whole blood; SPE

---

#### REFERENCE

Moore,C.M.; Tebbett,I.R. Rapid extraction of anti-inflammatory drugs in whole blood for HPLC analysis, *Forensic Sci.Int.*, **1987**, 34, 155-158.

---

**SAMPLE****Matrix:** blood**Sample preparation:** 1 mL Blood + 1 mL water + 50  $\mu$ L 76 mg/L allobarbitol in EtOH:water 10:90 + 5 mL ethyl acetate, shake by hand, add 2 mL of 0.1 M HCl. Mix by inversion with a mechanical shaker for 5 min. Centrifuge at 2700 rpm for 5-10 min. Remove ethyl acetate and evaporate to dryness under a stream of nitrogen at room temperature. Take up in 200  $\mu$ L MeOH, filter (0.45  $\mu$ m), inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES****Guard column:** Guard-Pak precolumn insert**Column:** 200  $\times$  4.6 5  $\mu$ m Hypersil octadecylsilane**Mobile phase:** Gradient. MeCN:1 mM pH 3.2 phosphate buffer from 20:80 to 40:60 over 10 min, stay at 40:60 for 6 min, to 20:80 over 4 min**Column temperature:** 60**Flow rate:** 3**Injection volume:** 20**Detector:** UV 202

---

**CHROMATOGRAM****Retention time:** 3.6**Internal standard:** allobarbitol (2.9)**Limit of quantitation:** 10000 ng/mL pharmacokinetics

---

**OTHER SUBSTANCES****Extracted:** butalbital**Simultaneous:** caffeine, aspirin

---

**REFERENCE**

Drost, M.L.; Walter, L. Blood and plasma concentrations of butalbital following single oral doses in man, *J. Anal. Toxicol.*, **1988**, *12*, 322-324.

---

**SAMPLE****Matrix:** blood**Sample preparation:** 200  $\mu$ L Serum + 400  $\mu$ L 100  $\mu$ g/mL 8-chlorotheophylline in MeCN, vortex for 10 s, centrifuge at 15000 rpm for 3 min, inject a 20  $\mu$ L aliquot of the supernatant.

---

**HPLC VARIABLES****Column:** 250  $\times$  4.6 5  $\mu$ m Ultrasphere ODS**Mobile phase:** MeCN:water:glacial acetic acid 4.84:12 containing 4.84 g/L Trizma, pH 2.3**Flow rate:** 2**Injection volume:** 20**Detector:** UV 254

---

**CHROMATOGRAM****Retention time:** 8.26**Internal standard:** 8-chlorotheophylline (5.29)**Limit of quantitation:** 8  $\mu$ g/mL

---

**OTHER SUBSTANCES****Extracted:** acetaminophen, theophylline**Simultaneous:** caffeine, cefazolin, cimetidine, ergotamine, glutethimide, heparin, methamphetamine, propranolol, sulfamethoxazole, theobromine, tobutamide, trimethoprim**Noninterfering:** amitriptyline, amobarbital, ampicillin, butabarbital, butalbital, celbenine, chlordiazepoxide, chlorpromazine, clorazepate, desipramine, diazepam, doxepin, ethchlorvynol, fluphenazine, hydroxyzine, ibuprofen, imipramine, isoniazid, lidocaine, mephobarbital, mesoridazine, methaqualone, methyluric acid, naprotyline, nordiazepam, nortriptyline, oxazepam, pentobarbital, perphenazine, phenelzine, phenmetrazine, phenobarbital, phenylbutazone, phenytoin, prednisolone, prednisone, procainamide, prochlorperazine, promazine, promethazine, propoxyphene, protriptyline, pyrilamine, secobarbital, thioridazine, thiothixene, timolol, trazodone, triazolam, trifluoperazine

---

**KEY WORDS**

serum

---

**REFERENCE**

Osterloh,J.; Yu,S. Simultaneous ion-pair and partition liquid chromatography of acetaminophen, theophylline and salicylate with application to 500 toxicologic specimens, *Clin.Chim.Acta*, **1988**, 175, 239–248.

---

**SAMPLE****Matrix:** blood

**Sample preparation:** 500  $\mu$ L Plasma + 2.5 mL 10  $\mu$ L/mL salicylamide in MeOH, stir, centrifuge. Remove a 1.8 mL aliquot of the supernatant and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 300  $\mu$ L mobile phase, inject an aliquot.

---

**HPLC VARIABLES****Guard column:** 10  $\mu$ m  $\mu$ Bondapak C18**Column:** 300  $\times$  3.9 10  $\mu$ m  $\mu$ Bondapak C18**Mobile phase:** MeCN:50 mM acetic acid 40:60**Flow rate:** 1**Detector:** UV 237

---

**CHROMATOGRAM****Retention time:** 5.6**Internal standard:** salicylamide (7.7)**Limit of detection:** 500 ng/mL

---

**KEY WORDS**

plasma; rat; dog; pharmacokinetics; also for humans (see *Int.J.Clin.Pharmacol.Ther.Toxicol.* 1988, 26, 421)

---

**REFERENCE**

Ramis,J.; Mis,R.; Forn,J. Pharmacokinetics of fosfosal in rats and dogs, *Arzneimittelforschung*, **1989**, 39, 74–77.

---

**SAMPLE****Matrix:** blood

**Sample preparation:** 100  $\mu$ L Plasma + 100  $\mu$ L 20  $\mu$ g/mL anisic acid + 8 mL dichloromethane, shake for 20 min, centrifuge at 2500 rpm for 20 min. Remove 7 mL of the organic layer and evaporate it to dryness under nitrogen or at 60°. Dissolve residue in 200  $\mu$ L mobile phase, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES****Column:** 150  $\times$  6 Shimpack CLS-ODS (Shimadzu)**Mobile phase:** MeCN:0.5 mM phosphoric acid 30:70**Column temperature:** 40**Flow rate:** 1.5**Injection volume:** 20**Detector:** UV 300

---

**CHROMATOGRAM****Internal standard:** anisic acid

---

**KEY WORDS**

plasma; rat

---

**REFERENCE**

Lee,C.K.; Uchida,T.; Kitagawa,K.; Yagi,A.; Kim,N.-S.; Goto,S. Skin permeability of various drugs with different lipophilicity, *J.Pharm.Sci.*, **1994**, 83, 562–565.

---

**SAMPLE****Matrix:** blood

**Sample preparation:** Add o-anisic acid to 1 mL plasma, acidify with HCl, extract with diethyl ether. Remove the organic layer and evaporate it to dryness, reconstitute the residue, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 150 × 4.6 5 µm Supelcosil LC-8

**Mobile phase:** MeCN:20 mM phosphoric acid 15:85

**Flow rate:** 1.6

**Detector:** UV 237

---

**CHROMATOGRAM**

**Internal standard:** o-anisic acid

**Limit of quantitation:** 100 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** aspirin

---

**KEY WORDS**

pharmacokinetics; plasma

---

**REFERENCE**

Benedek, I.H.; Joshi, A.S.; Pieniaszek, H.J.; King, S.-Y.P.; Kornhauser, D.M. Variability in the pharmacokinetics and pharmacodynamics of low dose aspirin in healthy male volunteers, *J.Clin.Pharmacol.*, **1995**, *35*, 1181–1186.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Filter (0.22 µm), inject a 2 µL aliquot of the filtrate.

---

**HPLC VARIABLES**

**Column:** 250 × 4 5 µm LiChrospher 100 Diol

**Mobile phase:** MeCN:50 mM pH 3.0 phosphate buffer 1.5:98.5

**Flow rate:** 0.5

**Injection volume:** 2

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 9

---

**OTHER SUBSTANCES**

**Extracted:** aspirin

---

**KEY WORDS**

serum; direct injection

---

**REFERENCE**

Nimura, N.; Itoh, H.; Kinoshita, T. Diol-bonded silica gel as a restricted access packing forming a binary-layered phase for direct injection of serum for the determination of drugs, *J.Chromatogr.A*, **1995**, *689*, 203–210.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Plasma. 200 µL Plasma + 200 µL 5 µg/mL IS in 200 mM HCL:200 mM orthophosphoric acid 50:50, vortex for 1-2 s, add 400 µL MeCN, vortex, let stand at 4° for 15 min, centrifuge at 10500 g for 1 min. Remove the supernatant and add it to 100-120 mg NaCl, vortex briefly, let stand at 4° for 10 min, vortex, centrifuge at 10500 g for 1 min, inject a 10 µL aliquot of the upper organic layer. Whole blood. 200 µL Lysed whole blood + 400 µL 5 µg/mL IS in 200 mM HCL:200 mM orthophosphoric acid 50:50, vortex for 1-2 s, add 600 µL MeCN, vortex, let stand at 4° for 15 min, centrifuge at 10500 g for 1 min. Remove the supernatant and add it to 200 mg NaCl, vortex briefly, let stand at 4° for 10 min, vortex, centrifuge at 10500 g for 1 min, inject a 10 µL aliquot of the upper organic layer.

---

**HPLC VARIABLES**

**Column:** 150 × 3.9 4 µm Novapak C18

**Mobile phase:** MeCN:water:85% orthophosphoric acid 18:74:0.09 (Before use prime column by recycling 200 mL mobile phase + 400 µL di-n-butylamine overnight at 0.3 mL/min.)

**Column temperature:** 30

**Flow rate:** 1

**Injection volume:** 10

**Detector:** UV 237

---

**CHROMATOGRAM**

**Retention time:** 6.8

**Internal standard:** 2-methylbenzoic acid (8.9)

**Limit of quantitation:** 100 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** metabolites, aspirin, gentisic acid, salicyluric acid

---

**KEY WORDS**

plasma; whole blood; pharmacokinetics

---

**REFERENCE**

Kees,F.; Jehnich,D.; Grobecker,H. Simultaneous determination of acetylsalicylic acid and salicylic acid in human plasma by high-performance liquid chromatography, *J.Chromatogr.B*, **1996**, 677, 172–177.

---

**SAMPLE**

**Matrix:** blood, CSF, gastric contents, urine

**Sample preparation:** 200 µL Serum, urine, CSF, or gastric fluid + 300 µL reagent. Flush column A to waste with 500 µL 500 mM ammonium sulfate, inject sample onto column A, flush column A to waste with 500 µL 500 mM ammonium sulfate, backflush the contents of column A onto column B with mobile phase, monitor the effluent from column B. (Reagent was 8.05 M guanidine HCl and 1.02 M ammonium sulfate in water.)

---

**HPLC VARIABLES**

**Column:** A 40 µm preparative grade C18 (Analytichem); B 75 × 2.1 pellicular C18 (Whatman) + 250 × 4.6 5 µm C8 end-capped (Whatman)

**Mobile phase:** Gradient. A was 50 mM pH 4.5 KH<sub>2</sub>PO<sub>4</sub>. B was MeCN:isopropanol 80:20. A:B 90:10 for 1 min, to 30:70 over 20 min.

**Column temperature:** 50

**Flow rate:** 1.5

**Detector:** UV 220

---

**CHROMATOGRAM**

**Retention time:** 5.2

**Internal standard:** heptanophenone (19)

---

**OTHER SUBSTANCES**

**Extracted:** acetaminophen, allobarbitol, azinphos, barbital, brallobarbitone, bromazepam, butethal, caffeine, carbamazepine, carbaryl, cephaloridine, chloramphenicol, chlordiazepoxide, chlorothiazide, chlorvinphos, clothiapine, cocaine, coomassie blue, desipramine, diazepam, diphenhydramine, dipipanone, ethylbromphos, flufenamic acid, formothion, griseofulvin, indomethacin, lidocaine, lorazepam, malathion, medazepam, midazolam, oxazepam, paraoxon, penicillin G, pentobarbital, prazepam, propoxyphene, prothiophos, quinine, secobarbital, strychnine, sulfamethoxazole, theophylline, thiopental, thioridazine, trimethoprim

---

**KEY WORDS**

serum; column-switching

---

**REFERENCE**

Kruger,P.B.; Albrecht,C.F.De V.; Jaarsveld,P.P. Use of guanidine hydrochloride and ammonium sulfate in comprehensive in-line sorption enrichment of xenobiotics in biological fluids by high-performance liquid chromatography, *J.Chromatogr.*, **1993**, 612, 191–198.

---

**SAMPLE**

**Matrix:** blood, luminal contents, tissue

**Sample preparation:** Homogenize intestinal tissue with water, make up to 25 mL with water, sonicate for 20 min. 100  $\mu$ L Whole blood, luminal contents or tissue homogenate + 200  $\mu$ L MeCN + 100  $\mu$ L 5 mM HCl, vortex for 30 s, centrifuge at 3000 rpm for 5 min, inject a 30  $\mu$ L aliquot of the clear supernatant.

---

**HPLC VARIABLES**

**Column:** 300  $\times$  3.9 10  $\mu$ m  $\mu$ Bondapak C18

**Mobile phase:** MeOH:water:glacial acetic acid 35:65:4

**Flow rate:** 1.5

**Injection volume:** 30

**Detector:** UV 313

---

**CHROMATOGRAM**

**Limit of detection:** <1  $\mu$ g/mL

---

**KEY WORDS**

rat; whole blood; intestine; pharmacokinetics

---

**REFERENCE**

Choi,Y.M.; Chung,S.M.; Chiou,W.L. First-pass accumulation of salicylic acid in gut tissue after absorption in anesthetized rat, *Pharm.Res.*, **1995**, *12*, 1323–1327.

---

---

**SAMPLE**

**Matrix:** blood, plants, urine, water

**Sample preparation:** Urine, plasma, water. Acidify 2 mL urine, plasma, or drinking water to pH 1 with 2 M HCl, shake with two 5 mL aliquots of diethyl ether, centrifuge. Remove the organic layer and add it to a few drops of saturated sodium bicarbonate solution, evaporate to dryness under a stream of nitrogen at 40°, reconstitute the residue in 2 mL MeOH, inject a 25  $\mu$ L aliquot. Hay, grain. Homogenize 30 g hay or grain with 10 mL 10% NaOH and 290 mL water, heat on a steam bath for 1 h, filter. Acidify a 10 mL aliquot of the filtrate to pH 1 with 2 M HCl, shake with two 5 mL aliquots of diethyl ether, centrifuge. Remove the organic layer and add it to a few drops of saturated sodium bicarbonate solution, evaporate to dryness under a stream of nitrogen at 40°, reconstitute the residue in 2 mL MeOH, inject a 25  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:**  $\mu$ Bondapak C18 radial compression

**Mobile phase:** MeCN:0.05% phosphoric acid 30:70 (plasma) or MeCN:0.2% acetic acid 20:80 (urine)

**Flow rate:** 1.5

**Injection volume:** 25

**Detector:** F ex 313 em 425

---

**CHROMATOGRAM**

**Limit of detection:** 100 ng/mL (plasma), 5  $\mu$ g/mL (urine)

---

**KEY WORDS**

horse; plasma; hay; grain; pharmacokinetics

---

**REFERENCE**

Beaumier,P.M.; Fenwick,J.D.; Stevenson,A.J.; Weber,M.P.; Young,L.M. Presence of salicylic acid in standardbred horse urine and plasma after various feed and drug administrations, *Equine.Vet.J.*, **1987**, *19*, 207–213.

---

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Plasma. 500  $\mu$ L Plasma + 900  $\mu$ L 270 mM HCl + 100  $\mu$ L 100  $\mu$ g/mL  $\alpha$ -phenylcinnamic acid in MeOH + 10 mL dichloromethane, shake at 125 cycles/min for 15 min, centrifuge at 750 g for 5 min. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 500  $\mu$ L MeOH, inject a 25  $\mu$ L aliquot. Urine. 2 mL Urine + 900  $\mu$ L 270 mM HCl + 100  $\mu$ L 100  $\mu$ g/mL  $\alpha$ -phenylcinnamic acid in MeOH + 10 mL hexane, shake at 125 cycles/

min for 15 min, centrifuge at 750 g for 5 min. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 500  $\mu\text{L}$  MeOH, inject a 25  $\mu\text{L}$  aliquot.

---

**HPLC VARIABLES**

**Column:** 300  $\times$  4  $\mu\text{m}$  Bondapak C18  
**Mobile phase:** MeOH:1% acetic acid 60:40  
**Flow rate:** 2  
**Injection volume:** 25  
**Detector:** UV 300

---

**CHROMATOGRAM**

**Retention time:** 3.0  
**Internal standard:**  $\alpha$ -phenylcinnamic acid (8.0)  
**Limit of detection:** 1  $\mu\text{g/mL}$

---

**OTHER SUBSTANCES**

**Extracted:** aspirin (UV 280), salsalate

---

**KEY WORDS**

plasma; pharmacokinetics

---

**REFERENCE**

Harrison, L.I.; Funk, M.L.; Ober, R.E. High-pressure liquid chromatographic determination of salicylsalicylic acid, aspirin, and salicylic acid in human plasma and urine, *J.Pharm.Sci.*, **1980**, 69, 1268–1271.

---

**SAMPLE**

**Matrix:** blood, urine  
**Sample preparation:** Plasma. 200  $\mu\text{L}$  Plasma + 200  $\mu\text{L}$  330 mM perchloric acid, mix thoroughly, centrifuge at 2600 g for 10 min, inject an aliquot of the supernatant. Urine. 100  $\mu\text{L}$  Urine + 900  $\mu\text{L}$  200 mM pH 5 phosphate buffer, mix, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 150  $\times$  4.6 8  $\mu\text{m}$  Cp Spher C8 (Chrompack) + 250  $\times$  4.6 5  $\mu\text{m}$  Spherisorb 5 ODS  
**Mobile phase:** Gradient. MeCN:buffer from 1:99 to 35:65 over 35 min, return to initial conditions over 3 min, re-equilibrate for 2 min. (Buffer was 6 g orthophosphoric acid and 1 mL glacial acetic acid in 1 L distilled water.)  
**Flow rate:** 1.5  
**Injection volume:** 20  
**Detector:** UV 236

---

**CHROMATOGRAM**

**Retention time:** 32.09  
**Limit of quantitation:** 5  $\mu\text{g/mL}$  (urine), 200 ng/mL (plasma)

---

**OTHER SUBSTANCES**

**Extracted:** metabolites, gentisic acid, hippuric acid, salicylic acid acyl glucuronide, salicylic acid phenolic glucuronide, salicyluric acid, salicyluric acid phenolic glucuronide

---

**KEY WORDS**

pharmacokinetics

---

**REFERENCE**

Vree, T.B.; van Ewijk-Beneken Kolmer, E.W.J.; Verwey-van Wissen, C.P.W.G.M.; Hekster, Y.A. Direct gradient reversed-phase high-performance liquid chromatographic determination of salicylic acid, with the corresponding glycine and glucuronide conjugates in human plasma and urine, *J.Chromatogr.B*, **1994**, 652, 161–170.

---

**SAMPLE**

**Matrix:** blood, urine  
**Sample preparation:** Plasma. Adjust pH of plasma to 3 with 43% phosphoric acid. 500  $\mu\text{L}$  Acidified plasma + 300  $\mu\text{L}$  10  $\mu\text{g/mL}$  m-hydroxybenzoic acid in MeCN + 1.2 mL MeCN, vortex

for 1-2 min, centrifuge at 2000 g for 10 min. Remove the supernatant and evaporate it to dryness under a stream of nitrogen at 37°, reconstitute the residue in 500 µL mobile phase, centrifuge at 1500 g for 2 min, inject a 50 µL aliquot of the supernatant. Urine. Dilute 25 µL urine to 500 µL with 40 µg/mL m-hydroxybenzoic acid in MeCN:water 10:90, inject a 25 µL aliquot.

---

**HPLC VARIABLES**

**Guard column:** 15 × 3.2 7 µm RP-18 (Brownlee)

**Column:** 150 × 4.6 5 µm Axxiom ODS (Springfield VA)

**Mobile phase:** MeCN:MeOH:25 mM acetic acid 8.5:8.5:83

**Flow rate:** 1

**Injection volume:** 25-50

**Detector:** UV 310

---

**CHROMATOGRAM**

**Retention time:** 16.5

**Internal standard:** m-hydroxybenzoic acid (8.2)

**Limit of quantitation:** 5 µg/mL (urine), 100 ng/mL (plasma)

---

**OTHER SUBSTANCES**

**Extracted:** metabolites

**Noninterfering:** etodolac, ibuprofen, ketorolac

---

**KEY WORDS**

plasma; paper also contains details of preparative HPLC; pharmacokinetics

---

**REFERENCE**

Liu, J.-H.; Smith, P.C. Direct analysis of salicylic acid, salicyl acyl glucuronide, salicyluric acid and gentisic acid in human plasma and urine by high-performance liquid chromatography, *J.Chromatogr.B*, **1996**, 675, 61-70.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18

**Column:** 250 × 4.6 5 µm Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 202.8

---

**CHROMATOGRAM**

**Retention time:** 12.12

---

**KEY WORDS**

whole blood



---

**REFERENCE**

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

---

**SAMPLE**

**Matrix:** bulk

**Sample preparation:** Dissolve in dilute HCl (pH 2), sonicate if necessary, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu$ m Hypersil-BDS C18

**Mobile phase:** MeOH:THF:buffer 11:4:85 (Buffer was 8.6 g  $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$ , 8.2 g NaCl, 2.2 g sodium 1-heptanesulfonate, and 5.75 mL 85% phosphoric acid in 2 L water, pH 2.0.)

**Column temperature:** 35

**Flow rate:** 1.5

**Injection volume:** 20

**Detector:** UV 215

---

**CHROMATOGRAM**

**Retention time:** 38

---

**OTHER SUBSTANCES**

**Simultaneous:** impurities, 5-aminosalicylic acid

---

**KEY WORDS**

impurity in 5-aminosalicylic acid

---

**REFERENCE**

Kersten,B.S.; Catalano,T.; Rozenman,Y. Ion-pairing high-performance liquid chromatographic method for the determination of 5-aminosalicylic acid and related impurities in bulk chemical, *J.Chromatogr.*, **1991**, 588, 187–193.

---

**SAMPLE**

**Matrix:** bulk, formulations

**Sample preparation:** Bulk. Prepare a 20 mg/mL solution of bulk aspirin in dichloromethane, inject a 10  $\mu$ L aliquot as soon as dissolution is complete. Tablets. Prepare a 20 mg/mL solution of ground aspirin tablets in dichloromethane, filter (0.45  $\mu$ m) immediately, immediately inject a 10  $\mu$ L aliquot of the filtrate.

---

**HPLC VARIABLES**

**Column:** 150  $\times$  4.6 6  $\mu$ m Zorbax SIL

**Mobile phase:** Hexane:chloroform:acetic acid 80:19:3 (Before first use pump 10 column volumes of dichloromethane:acetic acid:2,3-dimethoxypropane 96:2:2 through column at 3 mL/min.)

**Flow rate:** 3

**Injection volume:** 10

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 1.6

**Limit of detection:** 14 ppm

---

**OTHER SUBSTANCES**

**Simultaneous:** aspirin, salsalate

---

**KEY WORDS**

normal phase; tablets

---

**REFERENCE**

Pfeiffer,C.D.; Pankey,J.W. Determination of related compounds in aspirin by liquid chromatography, *J.Pharm.Sci.*, **1982**, 71, 511–514.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Vaccine. Centrifuge vaccine at 3400 g for 15 min, inject a 25  $\mu$ L aliquot of the supernatant.

---

**HPLC VARIABLES**

**Guard column:** 5  $\times$  4.5  $\mu$ m Hypersil C18

**Column:** 210  $\times$  4.6 5  $\mu$ m Hypersil C18

**Mobile phase:** MeOH:water:orthophosphoric acid 35:35:0.9, pH 2.5

**Flow rate:** 0.6

**Injection volume:** 25

**Detector:** UV 222

---

**CHROMATOGRAM**

**Retention time:** 6.5

**Internal standard:** salicylic acid

---

**OTHER SUBSTANCES**

**Simultaneous:** dithiosalicylic acid, thimerosal, thiosalicylic acid,

---

**KEY WORDS**

vaccine; salicylic acid is IS

---

**REFERENCE**

Tleugabulova,D.; Gonzalez Perez,I. Reversed-phase high-performance liquid chromatographic study of thimerosal stability in Cuban recombinant hepatitis B vaccine, *J.Chromatogr.A*, **1996**, 729, 219–227.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Sonicate 75 mg powdered tablets with 25 mL mobile phase for 15 min, filter (paper), inject a 135  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu$ m Ultrabase C18 (Scharlau Science, Spain)

**Mobile phase:** MeOH:20 mM pH 4.0  $\text{KH}_2\text{PO}_4$  30:70 adjusted to pH 4.0 with orthophosphoric acid

**Flow rate:** 1.5

**Injection volume:** 135

**Detector:** UV 224

---

**CHROMATOGRAM**

**Retention time:** 9.7

**Limit of quantitation:** 1.5 mg/L

---

**OTHER SUBSTANCES**

**Simultaneous:** aspirin, caffeine, thiamine

---

**KEY WORDS**

tablets

---

**REFERENCE**

Gámiz-Gracia,L.; Luque de Castro,M.D. An HPLC method for the determination of vitamin B1, caffeine, acetylsalicylic acid, and the impurities of salicylic acid in a pharmaceutical preparation, *J.Liq.Chromatogr.Rel.Technol.*, **1997**, 20, 2123–2133.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Place 5 tablets in MeCN:MeOH:85% phosphoric acid 92:8:0.5, sonicate 15 min, shake 15 min, dilute to 250 mL. Centrifuge an aliquot in 50 mL tube at 2000 rpm for 15 min and filter supernatant (0.45  $\mu$ m), inject an aliquot.

---

**HPLC VARIABLES****Column:** 150 × 3.9 Resolve (Waters)**Mobile phase:** MeCN:water:85% phosphoric acid 24:76:0.5**Column temperature:** 35**Flow rate:** 2**Injection volume:** 10**Detector:** UV 295

---

**CHROMATOGRAM****Retention time:** 3.3

---

**OTHER SUBSTANCES****Simultaneous:** aspirin, caffeine, acetaminophen

---

**KEY WORDS**film coated tablets; tablets

---

**REFERENCE**Fogel,J.; Epstein,P.; Chen,P. Simultaneous high-performance liquid chromatography assay of acetylsalicylic acid and salicylic acid in film-coated aspirin tablets, *J.Chromatogr.*, **1984**, 317, 507–511.

---

---

**SAMPLE****Matrix:** formulations**Sample preparation:** Grind tablets to a fine powder and add about 250 mg aspirin to 100 mL chloroform saturated with citric acid containing 500 µL formic acid, add 500 mg solid citric acid, sonicate for 2 min, centrifuge or filter, inject an aliquot. (If buffers or antacid are present add ground tablets equivalent to about 500 mg aspirin to 3 g acid-washed siliceous earth, mix, add 2 mL 6 M HCl, mix, add to a 200 × 25 column, dry wash container with siliceous earth, add to column, elute column with chloroform saturated with citric acid at 10 mL/min. Collect 150 mL eluent, add 1 mL formic acid, make up to 200 mL with chloroform saturated with citric acid, add 500 mg citric acid, shake, inject an aliquot.)

---

**HPLC VARIABLES****Column:** 250 × 4.6 5µm Zorbax-Sil**Mobile phase:** Chloroform:methylene chloride:acetonitrile:formic acid 700:300:30:4 (At the end of the day wash the column with 200 mL MeOH.)**Flow rate:** 2**Injection volume:** 20**Detector:** UV 300

---

**CHROMATOGRAM****Retention time:** 4.5

---

**OTHER SUBSTANCES****Simultaneous:** aspirin, acetylsalicylsalicylic acid, acetylsalicylic acid anhydride, excipients

---

**KEY WORDS**normal phase; tablets; SPE

---

**REFERENCE**Galante,R.N.; Visalli,A.J.; Grim,W.M. Stabilized normal-phase high-performance liquid chromatographic analysis of aspirin and salicylic acid in solid pharmaceutical dosage forms, *J.Pharm.Sci.*, **1984**, 73, 195–197.

---

---

**SAMPLE****Matrix:** formulations**Sample preparation:** Pulverize tablets and weigh out 1 g, add 1 mL formic acid, add 25 mL MeOH, shake mechanically for 10 min, make up to 50 mL with methanol. Remove 10 mL and centrifuge. 5 mL Supernatant + 5 mL 0.0025% p-hydroxybenzoic acid in MeOH:water 20:80, make up to 25 mL with water, inject an aliquot. (Analyze within 1 h.)

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 LiChrosorb RP8

**Mobile phase:** MeOH:200 mM pH 3.5 phosphate buffer:water 20:10:70

**Flow rate:** 1

**Injection volume:** 10

**Detector:** UV 280

---

**CHROMATOGRAM**

**Retention time:** 3

**Internal standard:** p-hydroxybenzoic acid (18)

---

**OTHER SUBSTANCES**

**Simultaneous:** aspirin (UV 254), p-aminophenol (UV 254), 3-O-acetylascorbic acid (UV 254), 2-O-acetylascorbic acid (UV 254), saccharin (UV 254), acetaminophen (UV 254), O-acetyl-p-aminophenol (UV 254), Ascorbic acid (UV 254), diacetyl-p-aminophenol

---

**KEY WORDS**

tablets

---

**REFERENCE**

Thomis,R.; Roets,E.; Hoogmartens,J. Analysis of tablets containing aspirin, acetaminophen, and ascorbic acid by high-performance liquid chromatography, *J.Pharm.Sci.*, **1984**, 73, 1830–1833.

---

**SAMPLE**

**Matrix:** formulations

---

**HPLC VARIABLES**

**Column:** phenyl

**Mobile phase:** MeCN:buffer 30:70 (Buffer was 10 mM KH<sub>2</sub>PO<sub>4</sub> containing 0.52% tetrabutylammonium phosphate, pH adjusted to 7.1 with KOH.)

**Flow rate:** 1

**Detector:** UV 220

---

**CHROMATOGRAM**

**Retention time:** 6

**Internal standard:** salicylic acid

---

**OTHER SUBSTANCES**

**Simultaneous:** nitroprusside

---

**KEY WORDS**

injections; 5% dextrose; salicylic acid is IS

---

**REFERENCE**

Pramar,Y.; Das Gupta,V.; Gardner,S.N.; Yau,B. Stabilities of dobutamine, dopamine, nitroglycerin and sodium nitroprusside in disposable plastic syringes, *J.Clin.Pharm.Ther.*, **1991**, 16, 203–207.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Condition a 500 mg Extract Clean silica SPE cartridge (Alltech stock no. 209250) with 2 mL hexane. Allow a solution of aspirin in 10 mM sorbitan trioleate in CFC-11 to evaporate, dissolve the residue in 5 mL hexane. Add 1 mL to the SPE cartridge, elute with two 2 mL portions of mobile phase, make up eluate to 5 mL with mobile phase, inject a 20 µL aliquot.

---

**HPLC VARIABLES**

**Guard column:** 20 × 2 3040 µm Perisorb RP-8 Pellicular (Upchurch)

**Column:** 250 × 4.6 5 µm Econosphere C8

**Mobile phase:** MeOH:THF:1 M phosphoric acid:water 44:5:5:46

**Flow rate:** 1

---

**Injection volume:** 20

**Detector:** UV 275

---

**CHROMATOGRAM**

**Retention time:** 10.6

---

**OTHER SUBSTANCES**

**Extracted:** aspirin, degradation products

---

**KEY WORDS**

aerosols; SPE

---

**REFERENCE**

Blondino,F.E.; Byron,P.R. The quantitative determination of aspirin and its degradation products in a model solution aerosol, *J.Pharm.Biomed.Anal.*, **1995**, *13*, 111–119.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Weigh out amount containing 100 mg salicylic acid, make up to 50 mL with mobile phase. Remove a 2 mL aliquot and dilute it to 10 mL with mobile phase, filter, inject a 20  $\mu$ L aliquot of the filtrate.

---

**HPLC VARIABLES**

**Column:** 125  $\times$  4.5  $\mu$ m LiChrospher 100 RP-18

**Mobile phase:** MeOH:water:96% acetic acid 55:44:1, pH 3.0

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 2.44

---

**OTHER SUBSTANCES**

**Simultaneous:** triamcinolone acetonide

---

**KEY WORDS**

topical solution

---

**REFERENCE**

Kedor-Hackmann,E.R.M.; Gianotto,E.A.S.; Santoro,M.I.R.M. Determination of triamcinolone acetonide and salicylic acid in pharmaceutical formulations by high performance liquid chromatography, *Pharmazie*, **1996**, *51*, 63–63.

---

**SAMPLE**

**Matrix:** perfusate

---

**HPLC VARIABLES**

**Column:** Lichrospher 100 RP-18

**Mobile phase:** MeOH:water containing 10 mM citric acid 60:40

**Column temperature:** 45

**Flow rate:** 1

**Detector:** F ex 300 em 408

---

**KEY WORDS**

rat; intestine; Caco-2 monolayer

---

**REFERENCE**

Takagi,M.; Taki,Y.; Sakane,T.; Nadai,T.; Sezaki,H.; Oku,N.; Yamashita,S. A new interpretation of salicylic acid transport across the lipid bilayer: Implication of pH-dependent but not carrier-mediated absorption from the gastrointestinal tract, *J.Pharmacol.Exp.Ther.*, **1998**, *285*, 1175–1180.

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 5 $\mu$ m Brownlee C18 (Applied Biosystems)

**Mobile phase:** MeCN:0.03% phosphoric acid:triethylamine 64:35:1, pH 2.0

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 237

---

**OTHER SUBSTANCES**

**Simultaneous:** diflunisal

---

**REFERENCE**

Hung,D.Y.; Mellick,G.D.; Anissimov,Y.G.; Weiss,M.; Roberts,M.S. Hepatic disposition and metabolite kinetics of a homologous series of diflunisal esters, *J.Pharm.Sci.*, **1998**, 87, 943–951.

---

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 250 × 4 ODS (Hitachi)

**Mobile phase:** MeOH:50 mM phosphoric acid 25:75

**Column temperature:** 55

**Flow rate:** 0.6

**Injection volume:** 20

**Detector:** UV 230

---

**OTHER SUBSTANCES**

**Also analyzed:** indoleacetic acid, phenylacetic acid

---

**REFERENCE**

Sugawara,M.; Takekuma,Y; Yamada,H.; Kobayashi,M.; Iseki,K.; Miyazaki,K. A general approach for the prediction of the intestinal absorption of drugs: regression analysis using the physicochemical properties and drug-membrane electrostatic interactions, *J.Pharm.Sci.*, **1998**, 87, 960–966.

---

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:**  $\mu$ Bondapak C18

**Mobile phase:** MeOH:8 mM pH 7.0 phosphate buffer 5:95

**Flow rate:** 1.2

**Injection volume:** 25

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 10

---

**OTHER SUBSTANCES**

**Simultaneous:** lodoxamide

---

**KEY WORDS**

aerosol

---

**REFERENCE**

Havel,H.A.; Beaubien,L.J.; Haaland,P.D. Analysis of the variance components in a pharmaceutical aerosol product: lodoxamide tromethamine, *J.Pharm.Sci.*, **1985**, 74, 978–982.

---

**SAMPLE****Matrix:** solutions**Sample preparation:** Dissolve compounds in MeOH, inject a 1  $\mu$ L aliquot.

---

**HPLC VARIABLES****Column:** 150  $\times$  1.3  $\mu$ m Hitachi-Gel 3011 porous polymer (Hitachi)**Mobile phase:** MeOH:ammonia 99:1**Flow rate:** 0.03**Injection volume:** 1**Detector:** UV 254

---

**CHROMATOGRAM****Retention time:** 3.56

---

**OTHER SUBSTANCES****Also analyzed:** acetaminophen, caffeine, buccetin (3-hydroxy-p-butyrophenetidine), phenacetin, dipyrone (sulpyrin), mefenamic acid, aspirin, salicylamide, ethenzamide (o-ethoxybenzamide), theobromine, theophylline

---

**KEY WORDS**

semi-micro; porous polymer

---

**REFERENCE**Matsushima, Y.; Nagata, Y.; Niyomura, M.; Takakusagi, K.; Takai, N. Analysis of antipyretics by semimicro liquid chromatography, *J. Chromatogr.*, **1985**, 332, 269–273.

---

**SAMPLE****Matrix:** solutions**Sample preparation:** Add 500  $\mu$ L of a solution in MeCN to 100 mg finely powdered potassium carbonate, add 250  $\mu$ L 3.8 mM 18-crown-6 in MeCN, add 250  $\mu$ L 0.8 mM reagent in MeCN, heat at 80° in the dark for 20 min, cool, inject a 5  $\mu$ L aliquot. (Synthesize the reagent, 3-bromomethyl-6,7-dimethoxy-1-methyl-2(1H)-quinoxalinone, as follows. Stir 483 g veratrole in 1.45 L acetic acid at 15° for 1 h, add 683 g concentrated nitric acid (d 1.05) over 1 h (maintain the temperature below 40° by cooling and regulating the rate of addition of the nitric acid). Continue stirring and add 2.127 L fuming nitric acid (d 1.50) over 1 h while maintaining the temperature below 30°, let stand for 2 h, pour into a large volume of cold water, filter, wash the solid with water until the washings are neutral, recrystallize from EtOH to give 4,5-dinitroveratrole (mp 129.5–130.5°) (*J. Am. Chem. Soc.* 1946, 68, 1536). Reflux 5 g 4,5-dinitroveratrole in 200 mL benzene (Caution! Benzene is a carcinogen!), add 100 g 60 mesh iron powder and 20 mL concentrated HCl in small portions over 1 h, reflux for 4 h, add 10 mL water, reflux for 2 h, cool, make alkaline with 2.5 M NaOH, extract several times with 200 mL portions of benzene. Combine the organic layers and evaporate them to dryness, add 10 mL concentrated HCl, recrystallize from EtOH to give 1,2-diamino-4,5-dimethoxybenzene monohydrochloride as very slightly pink needles (mp 240°) (*Anal. Chim. Acta* 1982, 134, 39). Heat 2.5 mmoles 1,2-diamino-4,5-dimethoxybenzene hydrochloride and 2.4 mmoles pyruvic acid in 30 mL 500 mM HCl on a boiling water bath for 2 h, cool with ice-water, filter. Wash the precipitate with water and dry it under vacuum, recrystallize from MeOH:water 90:10 to give 6,7-dimethoxy-3-methyl-2(1H)-quinoxalinone as yellow needles (mp 255°) (*Chem. Pharm. Bull.* 1985, 33, 3493). Treat 1 g 6,7-dimethoxy-3-methyl-2(1H)-quinoxalinone dissolved in 50 mL anhydrous MeOH with a solution of diazomethane in ether, evaporate to dryness under reduced pressure, dissolve the residue in 5 mL ethyl acetate, chromatograph on a 250  $\times$  35 column filled with 130 g 70-230 mesh silica gel 60 (Merck) using n-hexane:ethyl acetate 25:75 to give 6,7-dimethoxy-1,3-dimethyl-2(1H)-quinoxalinone as yellow needles (mp 170–171°). Dissolve 350 mg 6,7-dimethoxy-1,3-dimethyl-2(1H)-quinoxalinone in 3 mL acetic acid, add 350 mg anhydrous sodium acetate, add 2 mL 1.5 M bromine in acetic acid, heat at 100° for 15 min, cool, add 10 mL ether, filter, wash the solid 2 or 3 times with small portions of ether. Combine the filtrate and washings and evaporate them to dryness, dissolve the residue in 5 mL ethyl acetate, chromatograph on a 250  $\times$  35 column filled with 130 g 70-230 mesh silica gel 60 (Merck) using ether, evaporate the main fraction to dryness, recrystallize the residue from n-hexane:ethyl acetate 50:50 to give 3-bromomethyl-6,7-dimethoxy-1-methyl-2(1H)-quinoxalinone as yellow needles (mp 161–163°).)

---

**HPLC VARIABLES**

**Column:** 100 × 4 10 μm Radial-Pak C18 (Waters)

**Mobile phase:** Gradient. MeOH:water from 57:43 to 100:0 over 20 min, maintain at 100:0 for 12 min

**Flow rate:** 2

**Injection volume:** 5

**Detector:** F ex 370 em 450

---

**CHROMATOGRAM**

**Retention time:** 7.1

**Limit of detection:** 0.3-1 fmole

---

**OTHER SUBSTANCES**

**Simultaneous:** p-aminobenzoic acid, arachidic acid, arachidonic acid, benzoic acid, butyric acid, capric acid, caproic acid, caprylic acid, deoxyuridine, glucuronic acid, imidazole-4-acetic acid, lauric acid, linoleic acid, linolenic acid, margaric acid, 1-methyl-4-imidazoleacetic acid, myristic acid, myristoleic acid, oleic acid, palmitic acid, palmitoleic acid, propionic acid, stearic acid, thymidine, uridine, valeric acid

---

**KEY WORDS**

derivatization

---

**REFERENCE**

Yamaguchi, M.; Hara, S.; Matsunaga, R.; Nakamura, M.; Ohkura, Y. 3-Bromomethyl-6,7-dimethoxy-1-methyl-2(1H)-quinoxalinone as a new fluorescence derivatization reagent for carboxylic acids in high-performance liquid chromatography, *J. Chromatogr.*, **1985**, 346, 227-236.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Dissolve in MeOH:water 1:1 at a concentration of 50 μg/mL, inject a 10 μL aliquot.

---

**HPLC VARIABLES**

**Column:** 300 × 3.9 10 μm μBondapak C18

**Mobile phase:** MeOH:acetic acid:triethylamine:water 35:1.5:0.5:63

**Flow rate:** 1.5

**Injection volume:** 10

**Detector:** UV 261

---

**CHROMATOGRAM**

**Retention time:** 6

---

**OTHER SUBSTANCES**

**Simultaneous:** physostigmine

---

**REFERENCE**

Roos, R.W.; Lau-Cam, C.A. General reversed-phase high-performance liquid chromatographic method for the separation of drugs using triethylamine as a competing base, *J. Chromatogr.*, **1986**, 370, 403-418.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** React the carboxylic acid, triethylamine, and 1-(2,5-dihydroxyphenyl)-2-bromoethanone in a 1:2:4 molar ratio in MeCN at 45° for 2 h, inject a 10 μL aliquot. (Preparation of 1-(2,5-dihydroxyphenyl)-2-bromoethanone is as follows. Stir 27.6 g 1,4-dimethoxybenzene and 28 mL bromoacetyl bromide at 0°, add 53.4 g aluminum bromide over 10 min (an exothermic reaction ensues), let stand at room temperature for 12 h, add 100 mL 48% HBr, add 100 g ice, stir for 1 h, extract twice with 200 mL portions of diethyl ether. Combine the extracts and wash them 3 times with 200 mL portions of water, dry over 40 g anhydrous magnesium sulfate, evaporate to dryness, recrystallize the product 3 times from EtOH to yield 1-(2,5-dihydroxyphenyl)-2-bromoethanone monobromoacetate (mp 105-107°). Dissolve 11 g 1-(2,5-dihydroxyphenyl)-2-bromoethanone monobromoacetate in 200 mL warm dry MeOH satu-



rated with HBr, stir for 18 h, add 200 mL water, cool to  $-10^{\circ}$ . Collect the yellow solid and dry it under vacuum at  $50^{\circ}$  for 48 h, recrystallize from toluene:heptane 50:50 then toluene to obtain 1-(2,5-dihydroxyphenyl)-2-bromoethanone as yellow needles (mp  $117-119^{\circ}$ ).)

---

**HPLC VARIABLES**

**Column:**  $250 \times 4.7 \mu\text{m}$  RP-18 LiChrocart (Merck)

**Mobile phase:** MeOH:100 mM pH 6.5 sodium acetate 58:42

**Flow rate:** 1

**Injection volume:** 10

**Detector:** E, Bioanalytical Systems Model LC4B, glassy carbon electrode 0.6 V, Ag/AgCl reference electrode

---

**CHROMATOGRAM**

**Retention time:** 6

**Limit of detection:** 1 pmole

---

**OTHER SUBSTANCES**

**Simultaneous:** benzoic acid, quinoxaline-2-carboxylic acid

---

**KEY WORDS**

derivatization

---

**REFERENCE**

Munns, R.K.; Roybal, J.E.; Shimoda, W.; Hurlbut, J.A. 1-(4-Hydroxyphenyl)-, 1-(2,4-dihydroxyphenyl)- and 1-(2,5-dihydroxyphenyl)-2-bromoethanones: new labels for determination of carboxylic acids by high-performance liquid chromatography with electrochemical and ultraviolet detection, *J. Chromatogr.*, **1988**, *442*, 209–218.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** 10  $\mu\text{L}$  Solution + 500  $\mu\text{L}$  micelle solution + 25  $\mu\text{L}$  28 mg/mL 9-bromomethylacridine in acetone, mix, heat at  $60^{\circ}$  for 6 min, inject a 20  $\mu\text{L}$  aliquot. (Micelle solution was 25 mM Arkopal N-130 (a polyoxyethylene(13)nonylphenol, Hoechst Holland, Amsterdam) in 10 mM pH 7.0 phosphate buffer containing 6 mM tetrakis(decyl)ammonium bromide. Synthesize 9-bromomethylacridine as follows. Heat 10 g diphenylamine, 10 mL glacial acetic acid, and 40 g anhydrous zinc chloride to  $220^{\circ}$ , evaporate excess acetic acid with stirring, heat at  $220-230^{\circ}$  for 6 h, digest with hot 10% sulfuric acid, make strongly alkaline with 25% ammonia to dissolve the zinc chloride. Extract the insoluble residue with toluene. Extract the organic layer with 10% sulfuric acid, make the aqueous layer alkaline with aqueous ammonia. Collect the yellow precipitate that separates and recrystallize it twice from petroleum ether to give 9-methylacridine as pale yellow needles. Reflux 560 mg 9-methylacridine, 445 mg N-bromosuccinimide, and 10 mg benzoyl peroxide in 30 mL carbon tetrachloride for more than 2 h, cool, chromatograph on silica gel with benzene:ethyl acetate 30:1 (Caution! Benzene is a carcinogen!) to obtain 9-bromomethylacridine as yellow crystals (mp  $147-151^{\circ}$ ) (Anal. Lett. 1987, 20, 1581).)

---

**HPLC VARIABLES**

**Guard column:**  $10 \times 2.1$  40  $\mu\text{m}$  Chromsep C18 (Chrompack)

**Column:**  $100 \times 3.5 \mu\text{m}$  Chromspher C18 (Chrompack)

**Mobile phase:** Gradient. MeOH:water from 20:80 to 100:0 over 13 min.

**Injection volume:** 20

**Detector:** UV 254, F ex 362 em 418

---

**CHROMATOGRAM**

**Retention time:** 6

---

**OTHER SUBSTANCES**

**Simultaneous:** cholic acid, ibuprofen, valproic acid

---

**KEY WORDS**

derivatization

---

**REFERENCE**

van der Horst, F.A.L.; Post, M.H.; Holthuis, J.J.M.; Brinkman, U.A.T. Derivatization of carboxylic acids with 9-bromomethylacridine using micellar phase-transfer catalysis, *Chromatographia*, **1989**, *28*, 267–273.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a 500 µg/mL solution in MeOH:water 50:50, inject a 5 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 Zorbax C8

**Mobile phase:** Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L MeCN:water 20:80. A:B from 100:0 to 0:100 over 30 min. (Purify triethylamine as follows. Wash neutral alumina (Merck) 3 times with 2 bed volumes of pentane, 3 times with 2 bed volumes of dichloromethane, and 3 times with 2 bed volumes of MeOH, allow solvent to evaporate in a fume hood overnight, heat alumina at 130° for 2 h. Prepare a 14 cm column of the washed alumina in a 290 × 22 tube, pass through a head volume of MeOH, pass through triethylamine. When triethylamine starts to elute discard the first 20 mL, use the next 20 mL, discard the column.)

**Flow rate:** 2

**Injection volume:** 5

**Detector:** UV 210

---

**CHROMATOGRAM**

**Retention time:** 16.5

---

**OTHER SUBSTANCES**

**Simultaneous:** acetophenone, amphetamine, desipramine, ethylmorphine, imipramine, mefenamic acid, methamphetamine, morphine, phenylbutazone

---

**KEY WORDS**

also details of isocratic elution

---

**REFERENCE**

Hill, D.W. Evaluation of alkyl bonded silica and solvent phase modifiers for the efficient elution of basic drugs on HPLC, *J. Liq. Chromatogr.*, **1990**, *13*, 3147–3175.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Inject a 50 µL aliquot of a solution in mobile phase.

---

**HPLC VARIABLES**

**Column:** 100 × 4.6 5 µm Spheri-5 RP-8

**Mobile phase:** MeOH:buffer 30:70 (Prepare buffer by mixing 4 mM Na<sub>2</sub>HPO<sub>4</sub> and 7 mM KH<sub>2</sub>PO<sub>4</sub> to achieve pH 7.)

**Flow rate:** 1

**Injection volume:** 50

**Detector:** F ex 355 em 460 (408 nm cutoff filter) following post-column extraction. The column effluent mixed with 50 µg/mL reagent in mobile phase pumped at 0.5 mL/min and then with chloroform pumped at 1 mL/min and the mixture flowed through a 1.8 m × 0.3 mm ID knitted PTFE coil to a 50 µL membrane phase separator using a polyethylene-backed 0.5 µm Fluoropore membrane filter (design in paper). The organic phase flowed to the detector. (Synthesize the reagent, α-(3,4-dimethoxyphenyl)-4'-trimethylammoniummethylcinnamonitrile methosulfate, as follows. Stir 20 mmoles 3,4-dimethoxyphenylacetonitrile and 20 mmoles p-toluamide in 50 mL EtOH at 50°, add 5 mL 50% aqueous KOH slowly, stir at 50° for 5 min, cool to room temperature, filter, dry the precipitate of α-(3,4-dimethoxyphenyl)-4'-methylcinnamonitrile. Dissolve 20 mmoles α-(3,4-dimethoxyphenyl)-4'-methylcinnamonitrile, 20 mmoles N-bromosuccinimide, and 20 mg benzoyl peroxide in 100 mL carbon tetrachloride (Caution! Carbon tetrachloride is a carcinogen!), reflux with stirring for 1.5 h, cool, filter, evaporate to dryness under reduced pressure, recrystallize from MeOH to give α-(3,4-dimethoxyphenyl)-4'-bromomethylcinnamonitrile. Vigorously stir 30 mmoles anhydrous dimethylamine in 100 mL dry benzene (Caution! Benzene is a carcinogen!) at 0°, very slowly add 10 mmoles α-(3,4-dime-

thoxyphenyl)-4'-bromomethylcinnamionitrile while stirring at 0°, stir at room temperature overnight, add 150 mL water, remove the organic phase, extract the aqueous phase twice with 100 mL portions of diethyl ether, wash the organic layers with saturated NaCl solution, dry over anhydrous magnesium sulfate, evaporate under reduced pressure to give  $\alpha$ -(3,4-dimethoxyphenyl)-4'-dimethylaminomethylcinnamionitrile (J.Chem.Eng.Data 1987, 32, 387). Reflux 10 mmoles  $\alpha$ -(3,4-dimethoxyphenyl)-4'-dimethylaminomethylcinnamionitrile, 20 mmoles dimethyl sulfate (Caution! Dimethyl sulfate is a carcinogen and acutely toxic!), and 5 g potassium carbonate in 50 mL acetone for 1 h, cool to room temperature, filter, dry the precipitate under vacuum at room temperature overnight, recrystallize from chloroform containing 2-3 drops of 95% EtOH to give  $\alpha$ -(3,4-dimethoxyphenyl)-4'-trimethylammoniummethylcinnamionitrilemethosulfate (mp 212-215°). Protect solutions from light.)

---

#### CHROMATOGRAM

**Retention time:**  $k'$  0.2791

**Limit of detection:** 1  $\mu$ g/mL

---

#### OTHER SUBSTANCES

**Simultaneous:** ibuprofen, ketoprofen, mefenamic acid, naproxen, probenecid, valproic acid

---

#### KEY WORDS

post-column extraction; post-column reaction

---

#### REFERENCE

Kim, M.; Stewart, J. T. HPLC post-column ion-pair extraction of acidic drugs using a substituted  $\alpha$ -phenylcinnamionitrile quaternary ammonium salt as a new fluorescent ion-pair reagent, *J. Liq. Chromatogr.*, **1990**, *13*, 213-237.

---

#### SAMPLE

**Matrix:** solutions

---

#### HPLC VARIABLES

**Column:** 250  $\times$  4.6 Zorbax RX

**Mobile phase:** Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

**Column temperature:** 30

**Flow rate:** 2

**Detector:** UV 210

---

#### OTHER SUBSTANCES

**Also analyzed:** acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenpropofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiacol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminosilbene, imipramine, indomethacin, isocarboxystiril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin,

mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methypyrrolon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, scopalamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleannamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

---

## REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233–242.

---

## SAMPLE

**Matrix:** solutions

**Sample preparation:** Prepare a 50 µg/mL solution in the mobile phase, inject a 10 µL aliquot.

---

## HPLC VARIABLES

**Column:** 250 × 4.6 7 µm Lichrosorb RP 18

**Mobile phase:** MeOH:water 45:55 containing 1% acetic acid

**Flow rate:** 1

**Injection volume:** 10

**Detector:** UV 230

---

## CHROMATOGRAM

**Retention time:** 7.63

---

## OTHER SUBSTANCES

**Simultaneous:** acetaminophen, aspirin, phenacetin, salicylamide

---

## REFERENCE

Nivaud-Guernet,E.; Guernet,M.; Ivanovic,D.; Medenica,M. Effect of eluent pH on the ionic and molecular forms of the non-steroidal anti-inflammatory agents in reversed-phase high-performance liquid chromatography, *J.Liq.Chromatogr.*, **1994**, *17*, 2343–2357.

---

## SAMPLE

**Matrix:** solutions

**Sample preparation:** 50 µL Solution + 50 µL pH 7.4 PBS + 100 µL MeOH, centrifuge at 12000 g for 10 min, inject a 50 µL aliquot.

---

## HPLC VARIABLES

**Column:** 150 × 4.6 Cosmosil 5C18-P (Nacalai Tesque)

**Mobile phase:** MeOH:50 mM NaH<sub>2</sub>PO<sub>4</sub> 20:80

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 220

---

## CHROMATOGRAM

**Internal standard:** salicylic acid

---

**OTHER SUBSTANCES****Simultaneous:** atenolol

---

**KEY WORDS**buffer; Earle's balanced salt solution; salicylic acid is IS

---

**REFERENCE**Sasaki,H.; Igarishi,Y.; Nishida,K.; Nakamura,J. Intestinal permeability of ophthalmic  $\beta$ -blockers for predicting ocular permeability, *J.Pharm.Sci.*, **1994**, 83, 1335-1338.

---

**SAMPLE****Matrix:** solutions

---

**HPLC VARIABLES****Column:** 250  $\times$  4.6 5  $\mu$ m Supelcosil LC-DP (A) or 250  $\times$  4.5  $\mu$ m LiChrospher 100 RP-8 (B)**Mobile phase:** MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)**Flow rate:** 0.6**Injection volume:** 25**Detector:** UV 229

---

**CHROMATOGRAM****Retention time:** 5.20 (A), 4.35 (B)

---

**OTHER SUBSTANCES**

**Also analyzed:** acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordiazepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrizidine, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenopropfen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazinol, mefenamic acid, meperidine, mephentoin, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methylodopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfapyrazone, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocainide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimetoprim, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

---

**KEY WORDS**details of plasma extraction

---

**REFERENCE**Koves,E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J.Chromatogr.A*, **1995**, 692, 103-119.

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 150 mm long 5  $\mu$ m Microsorb-MV C18

**Mobile phase:** MeOH:water 30:70

**Flow rate:** 1

**Detector:** UV 300

---

**REFERENCE**

Phillips,C.A.; Michniak,B.B. Transdermal delivery of drugs with differing lipophilicities using azone analogs as dermal penetration enhancers, *J.Pharm.Sci.*, **1995**, *84*, 1427–1433.

---

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a 1-10  $\mu$ g/mL solution in water, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu$ m Hypersil SCX/C18

**Mobile phase:** MeCN:25 mM pH 3  $\text{Na}_2\text{HPO}_4$  50:50

**Injection volume:** 20

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** k' 0.50

---

**OTHER SUBSTANCES**

**Also analyzed:** amitriptyline, barbital, benzoic acid, butabarbital, clomipramine, clonazepam, desipramine, diazepam, flurazepam, furosemide, imipramine, nitrazepam, phenobarbital, phenol, phenolphthalein, pindolol, propranolol, resorcinol, secobarbital, terbutaline, xylazine

---

**KEY WORDS**

effect of mobile phase pH on capacity factor is discussed

---

**REFERENCE**

Walshe,M.; Kelly,M.T.; Smyth,M.R.; Ritchie,H. Retention studies on mixed-mode columns in high-performance liquid chromatography, *J.Chromatogr.A*, **1995**, *708*, 31–40.

---

---

**SAMPLE**

**Matrix:** wine

**Sample preparation:** Adjust pH of wine to 7-8 with potassium bicarbonate. Remove a 1 mL aliquot and add it to 1 mL 170 mM phenacyl bromide in acetone, add 1 mL 17 mM 18-crown-6 in acetone, add 1 mL acetone, heat in a boiling water bath for 75 min, cool, inject a 10  $\mu$ L aliquot. (Recrystallize phenacyl bromide from n-heptane.)

---

**HPLC VARIABLES**

**Guard column:** 37-50  $\mu$ m Bondapak C18/Corasil

**Column:** 250  $\times$  4 7  $\mu$ m RP-18 (Merck)

**Mobile phase:** Gradient. MeOH:water from 35:65 to 85:15 over 20 min.

**Flow rate:** 2

**Injection volume:** 10

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 11.7

---

---

**OTHER SUBSTANCES**

**Extracted:** acetic acid, anisic acid, benzoic acid, benzoic acid, butyric acid, caprylic acid, cinnamic acid, citramalic acid, citric acid, enanthic acid, fumaric acid, galacturonic acid, gallic acid, glutaric acid, glycolic acid, glyoxylic acid, p-hydroxybenzoic acid, isocitric acid,  $\alpha$ -ketoglutaric acid, lactic acid, malic acid, mandelic acid, phenylacetic acid, propionic acid, protocatechuic acid, pyruvic acid, sorbic acid, succinic acid, tartaric acid, valeric acid, vanillic acid, ascorbic acid

---

**KEY WORDS**

derivatization

---


**REFERENCE**

Mentasti,E.; Gennaro,M.C.; Sarzanini,C.; Baiocchi,C.; Savigliano,M. Derivatization, identification and separation of carboxylic acids in wines and beverages by high-performance liquid chromatography, *J.Chromatogr.*, **1985**, 322, 177–189.

---

---

# Salinomycin



**Molecular formula:**  $C_{42}H_{70}O_{11}$

**Molecular weight:** 751.01

**CAS Registry No.:** 53003-10-4, 55721-31-8 (sodium salt)

**Merck Index:** 8488

---

**SAMPLE**

**Matrix:** albumen, eggs, feed, premix, tissue

**Sample preparation:** Eggs, feed, premix, tissue. 1 g Sample + 15 mL acetone, vortex for 5 min, centrifuge, decant the acetone layer. Re-extract residue (2x), evaporate the combined acetone layers to dryness. Partition the residue between 25 mL aqueous EtOH (water:EtOH 80:20) and 25 mL petroleum ether (50–110°). Evaporate EtOH fraction to dryness, dissolve residue in MeOH, adjust volume to 500  $\mu$ L. Inject a 50 or 100  $\mu$ L aliquot. Albumen. 1 g Sample + 15 mL MeOH, extract for 15 min, centrifuge at  $1290 \times g$  at  $-5^\circ$ . Evaporate to dryness, dissolve residue in MeOH, adjust to a final volume of 500  $\mu$ L. Inject a 50 or 100  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Guard column:**  $10 \times 4.6$  5  $\mu$ m Inertsil ODS-2 (Mandel Scientific)

**Column:**  $250 \times 3.2$  5  $\mu$ m CSC-Inertsil 150A/ODS-2 (Mandel Scientific)

**Mobile phase:** MeOH:water:acetic acid 95:5:0.1

**Flow rate:** 0.5

**Injection volume:** 50–100

**Detector:** UV 520 nm following post-column reaction. The column effluent mixed with vanillin reagent pumped at 1 mL/min. The mixture flowed through a Teflon knitted reactor coil ( $10 \text{ m} \times 0.5 \text{ mm i.d.}$ , total volume 2 mL) at  $95^\circ$  to the detector. (Prepare the vanillin reagent as follows. Slowly and carefully add 20 mL concentrated sulphuric acid (95–98%) to 950 mL chilled MeOH, mix well, allow to cool to room temperature. Add 30 g vanillin with constant stirring. Filter solution, degas using vacuum, and store in amber-colored bottle.)

---

**CHROMATOGRAM**

**Retention time:** <10

**Limit of detection:** 5 ng/g

**Limit of quantitation:** 10 ng/g

---

**KEY WORDS**

post-column reaction