


# Tazarotene

**Molecular formula:** C<sub>21</sub>H<sub>21</sub>NO<sub>2</sub>S

**Molecular weight:** 351.47

**CAS Registry No.:** 118292-40-3

**Merck Index:** 9249



## SAMPLE

**Matrix:** blood, microsomal incubations

**Sample preparation:** Blood. 500 µL Whole blood + 2.5 mL MeCN:1-butanol 50:50 + 200 µL 1 µg/mL IS, vortex for 30 s, centrifuge at 1500 g for 5 min. Evaporate supernatant to dryness with nitrogen. Microsomal incubations. 1 mL Microsomal incubation + 2.5 mL MeCN:1-butanol 50:50, centrifuge at 1500 g for 5 min. Decant the supernatant, dry it with nitrogen and reconstitute in 300 µL mobile phase. Inject a 100 µL aliquot.

## HPLC VARIABLES

**Guard column:** dry packed C18

**Column:** 150 × 4.6 5 µm Ultrasphere C18

**Mobile phase:** MeCN:water:acetic acid 65:34.5:0.5

**Flow rate:** 1.2

**Injection volume:** 100

**Detector:** UV 345

## CHROMATOGRAM

**Retention time:** 17

**Internal standard:** AGN 190252 (7)

## OTHER SUBSTANCES

**Extracted:** metabolites

## KEY WORDS

liver; human; rat; whole blood

## REFERENCE

Madhu,C.; Duff,S.; Baumgarten,V.; Rix,P.; Small,D.; Tang-Liu,D. Metabolic deesterification of tazarotene in human blood and rat and human liver microsomes, *J.Pharm.Sci.*, **1997**, *86*, 972–974.

# Tazobactam


**Molecular formula:** C<sub>10</sub>H<sub>12</sub>N<sub>4</sub>O<sub>5</sub>S

**Molecular weight:** 300.30

**CAS Registry No.:** 89786-04-9

**Merck Index:** 9251

**Lednicer No.:** 5 156



## SAMPLE

**Matrix:** bile, blood, urine

**Sample preparation:** Plasma, serum. 200 µL Plasma or serum + 200 µL 25 µg/mL penicillin G in 50 mM pH 6.0 sodium phosphate buffer + 800 µL MeCN, vortex for 30 s, centrifuge at 3000 g for 10 min. Remove the clear supernatant and add it to 2 mL dichloromethane, mix for 30 s, centrifuge at 3000 g for 10 min, inject a 25 µL aliquot of the upper aqueous layer. Bile. 200 µL Bile + 400 µL 50 mM pH 7.0 sodium phosphate buffer + 2 mL MeCN, vortex for 30 s, centrifuge at 3000 g for 10 min. Remove the clear supernatant and add it to 2 mL dichloromethane, mix for 30 s, centrifuge at 3000 g for 10 min, inject a 25 µL aliquot of the upper

aqueous layer. Urine. 100  $\mu$ L Urine + 50  $\mu$ L 5 mg/mL penicillin G in water, vortex for 30 s, make up to 10 mL with 50 mM pH 6.0 sodium phosphate buffer, inject a 25  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Guard column:** 15  $\times$  3.2 7  $\mu$ m Brownlee C 18 guard column

**Column:** 250  $\times$  4.6 5  $\mu$ m Hypersil ODS (Keystone)

**Mobile phase:** Gradient. A was MeCN:10 mM sodium phosphate adjusted to pH 2.7 with concentrated phosphoric acid 3:97. B was MeCN:10 mM sodium phosphate adjusted to pH 2.7 with concentrated phosphoric acid 90:10. A:B from 95:5 to 50:50 over 9 min and then to 95:5 over 1 min.

**Flow rate:** 1.5

**Injection volume:** 25

**Detector:** UV 220

---

**CHROMATOGRAM**

**Retention time:** 5.7

**Internal standard:** penicillin G (12.5)

**Limit of quantitation:** 50000 ng/mL (urine), 1000 ng/mL (plasma, serum, bile)

---

**OTHER SUBSTANCES**

**Extracted:** piperacillin

**Simultaneous:** amoxicillin, ampicillin, cefoperazone, cefometazole, cefotaxime, cefotetan, cefuroxime, mezlocillin

---

**KEY WORDS**

plasma; serum

---

**REFERENCE**

Ocampo,A.P.; Hoyt,K.D.; Wadgaonkar,N.; Carver,A.H.; Puglisi,C.V. Determination of tazobactam and piperacillin in human plasma, serum, bile and urine by gradient elution reversed-phase high-performance liquid chromatography, *J.Chromatogr.*, **1989**, 496, 167–179.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Homogenize tissue in water at 4° for 30 s, centrifuge at 15000 rpm (Sorvall centrifuge) for 10 min. 250  $\mu$ L Plasma or tissue homogenate + 500  $\mu$ L 10  $\mu$ g/mL cefpodoxime in MeCN, mix, centrifuge at 15000 rpm (Sorvall centrifuge), extract with 1 mL dichloromethane, inject 120  $\mu$ L of the aqueous phase onto column A with mobile phase A, elute the contents of column A onto column B with mobile phase B, monitor the effluent from column B.

---

**HPLC VARIABLES**

**Column:** A 40  $\times$  4.6 10  $\mu$ m Lichrosorb RP 2; B 250  $\times$  4.6 5  $\mu$ m Spherisorb ODS II

**Mobile phase:** A MeCN:100 mM NaH<sub>2</sub>PO<sub>4</sub> and 5 mM tetrabutylammonium hydrogen sulfate 5:95, pH adjusted to 6.5; B MeCN:100 mM NaH<sub>2</sub>PO<sub>4</sub> and 5 mM tetrabutylammonium hydrogen sulfate 10:90, pH adjusted to 6.5

**Column temperature:** 25

**Flow rate:** A 1; B 1.5

**Injection volume:** 120

**Detector:** UV 210 (UV 300 for IS)

---

**CHROMATOGRAM**

**Retention time:** 18.6

**Internal standard:** cefpodoxime (24.9)

**Limit of detection:** 96 ng/mL

---

**KEY WORDS**

plasma; column-switching

---

**REFERENCE**

Kinzig,M.; Sörgel,F.; Brismar,B.; Nord,C.E. Pharmacokinetics and tissue penetration of tazobactam and piperacillin in patients undergoing colorectal surgery, *Antimicrob.Agents Chemother.*, **1992**, 36, 1997–2004.

**SAMPLE****Matrix:** blood**Sample preparation:** 1 mL Serum + 400  $\mu$ L 10 mg/mL zinc sulfate containing 350  $\mu$ g/mL benzoic acid, vortex for 30 s, centrifuge at 5500 g for 5 min, inject a 20  $\mu$ L aliquot of the supernatant.

---

**HPLC VARIABLES****Column:** 150  $\times$  4.3 5  $\mu$ m Nova Pak**Mobile phase:** MeOH:pH 6.30 phosphate buffer 5:95**Column temperature:** 45**Flow rate:** 2**Injection volume:** 20**Detector:** UV 225

---

**CHROMATOGRAM****Retention time:** 2.4**Internal standard:** benzoic acid (1.9)**Limit of detection:** 5  $\mu$ g/mL

---

**OTHER SUBSTANCES****Extracted:** sulbactam

---

**KEY WORDS**serum

---

**REFERENCE**Guillaume,Y.; Peyrin,E.; Guinchard,C. Rapid determination of sulbactam and tazobactam in human serum by high-performance liquid chromatography, *J.Chromatogr.B*, **1995**, 665, 363–371.

---

**SAMPLE****Matrix:** blood**Sample preparation:** 1 mL Serum or hemofiltrate + 100  $\mu$ L 3 M sulfuric acid + 4 mL diethyl ether, vortex for 1 min, centrifuge at 5500 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 37°, reconstitute the residue in 100  $\mu$ L reagent and 10  $\mu$ L 1% mercuric chloride in water, heat at 37° for 30 min, inject a 20  $\mu$ L aliquot. (Prepare reagent by dissolving 4 g 1,2,4-triazole in water, adjusting the pH to 9.0 with 10 M NaOH, and making up to 20 mL with water.)

---

**HPLC VARIABLES****Column:** 125  $\times$  4 5  $\mu$ m Lichrospher RP (18)**Mobile phase:** MeOH:(NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub> buffer 0.5:99.5, adjusted to pH 6.00 with phosphoric acid (Buffer concentration not given.)**Flow rate:** 2.2**Injection volume:** 20**Detector:** UV 325

---

**CHROMATOGRAM****Retention time:** 8**Limit of detection:** 50 ng/mL

---

**KEY WORDS**serum; hemofiltrate; derivatization

---

**REFERENCE**Peyrin,E.; Guillaume,Y.; Guinchard,C. High-performance liquid chromatographic determination of tazobactam by precolumn derivatization, *J.Chromatogr.B*, **1995**, 672, 160–164.

---

**SAMPLE****Matrix:** blood, tissue

**Sample preparation:** Plasma. 250  $\mu$ L Plasma + 500  $\mu$ L 10  $\mu$ g/mL cefpodoxime in MeCN, mix, centrifuge at 15000 rpm (Sorvall centrifuge). Remove the supernatant and add it to 1 mL dichloromethane, extract, inject a 120  $\mu$ L aliquot of the aqueous phase. Tissue. Blot tissue, homogenize with two volumes of water (IKA-Ultra-Turrax) at 4° for 30 s, centrifuge at 15000 rpm (Sorvall centrifuge) for 10 min, remove 250  $\mu$ L of the supernatant, add 500  $\mu$ L 10  $\mu$ g/mL cefpodoxime in MeCN, mix, centrifuge at 15000 rpm (Sorvall centrifuge). Remove the supernatant and add it to 1 mL dichloromethane, extract, inject a 120  $\mu$ L aliquot of the aqueous phase. (A column-switching technique is used but no details are given in the paper.)

---

#### HPLC VARIABLES

**Column:** A 40  $\times$  4.6 10  $\mu$ m RP-2 (Bischoff Chromatography)); B 250  $\times$  4.6 5  $\mu$ m Spherisorb ODS II

**Mobile phase:** A MeCN:100 mM NaH<sub>2</sub>PO<sub>4</sub> + 5 mM tetrabutylammonium hydrogen sulfate 5:95, pH 6.5; B MeCN:100 mM NaH<sub>2</sub>PO<sub>4</sub> + 5 mM tetrabutylammonium hydrogen sulfate 10:90, pH 6.5

**Column temperature:** 25

**Flow rate:** A 1; B 1.5

**Injection volume:** 120

**Detector:** UV 210 for tazobactam, UV 300 for cefpodoxime

---

#### CHROMATOGRAM

**Retention time:** 18.6

**Internal standard:** cefpodoxime (24.9)

**Limit of quantitation:** 76 ng/g (tissue), 96 ng/mL (plasma)

---

#### KEY WORDS

plasma; column-switching; pharmacokinetics; fat; muscle; skin; intestinal mucosa; appendix

---

#### REFERENCE

Kinzig,M.; Sörgel,F.; Brismar,B.; Nord,C.E. Pharmacokinetics and tissue penetration of tazobactam and piperacillin in patients undergoing colorectal surgery, *Antimicrob.Agents Chemother.*, **1992**, 36, 1997–2004.

---

#### SAMPLE

**Matrix:** formulations

**Sample preparation:** Dilute 18:1 with mobile phase, inject a 20  $\mu$ L aliquot.

---

#### HPLC VARIABLES

**Column:** 250  $\times$  4.6 5  $\mu$ m Adsorbosphere

**Mobile phase:** MeCN:10 mM NaH<sub>2</sub>PO<sub>4</sub> 7:93, pH adjusted to 2.7 with 85% phosphoric acid

**Flow rate:** 1.2

**Injection volume:** 20

**Detector:** UV 215

---

#### CHROMATOGRAM

**Retention time:** 5.0

---

#### KEY WORDS

injections; saline; stability-indicating

---

#### REFERENCE

Choi,J.-S.; Burm,J.-P.; Jhee,S.S.; Chin,A.; Ulrich,R.W.; Gill,M.A. Stability of piperacillin sodium-tazobactam sodium and ranitidine hydrochloride in 0.9% sodium chloride injection during simulated Y-site administration, *Am.J.Hosp.Pharm.*, **1994**, 51, 2273–2276.

---

#### SAMPLE

**Matrix:** formulations

**Sample preparation:** Dilute 8-fold to 20-fold with saline, filter (0.2  $\mu$ m), inject a 20  $\mu$ L aliquot of the filtrate.

---

#### HPLC VARIABLES

**Column:** 250  $\times$  4.6 5  $\mu$ m Adsorbosphere C18

**Mobile phase:** MeCN:100 mM sodium phosphate buffer 7:93 adjusted to pH 2.7 with 85% phosphoric acid  
**Flow rate:** 1.2  
**Injection volume:** 20  
**Detector:** UV 215

---

**CHROMATOGRAM**

**Retention time:** 7-8

---

**KEY WORDS**

saline; injections; stability-indicating

---

**REFERENCE**

Chung,K.C.; Moon,Y.S.K.; Chin,A.; Ulrich,R.W.; Gill,M.A. Compatibility of ondansetron hydrochloride and piperacillin sodium tazobactam sodium during simulated Y-site administration, *Am.J.Health-Syst.Pharm.*, **1995**, 52, 1554-1556.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Dilute a 25-50  $\mu$ L sample with 2 mL mobile phase, filter (0.2  $\mu$ m), inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu$ m Adsorbosphere C18

**Mobile phase:** MeCN:10 mM sodium phosphate 7:93 adjusted to pH 2.7 with 85% phosphoric acid

**Flow rate:** 1.2

**Injection volume:** 20

**Detector:** UV 215

---

**CHROMATOGRAM**

**Retention time:** 6.40

---

**OTHER SUBSTANCES**

**Simultaneous:** degradation products

---

**KEY WORDS**

injections; saline; 5% dextrose

---

**REFERENCE**

Moon,Y.S.K.; Chung,K.C.; Chin,A.; Gill,M.A. Stability of piperacillin sodium-tazobactam sodium in polypropylene syringes and polyvinyl chloride minibags, *Am.J.Health-Syst.Pharm.*, **1995**, 52, 999-1001.


---

# Tazofelone

**Molecular formula:** C<sub>18</sub>H<sub>27</sub>NO<sub>2</sub>S

**Molecular weight:** 321.49

**CAS Registry No.:** 136433-51-7



---

**SAMPLE**

**Matrix:** microsomal incubations

**Sample preparation:** Add 1 mL cold MeCN to 250  $\mu$ L microsomal incubation, centrifuge at 13000 rpm for 3 min, evaporate the supernatant under nitrogen, reconstitute the residue in 100  $\mu$ L MeCN:water 40:60, inject 60  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 Prodigy C 18 (Phenomenex)

**Mobile phase:** MeCN:water 60:40  
**Column temperature:** 40  
**Flow rate:** 1  
**Injection volume:** 60  
**Detector:** UV 214

---

**CHROMATOGRAM**

**Retention time:** 15

---

**OTHER SUBSTANCES**

**Extracted:** metabolites

---

**KEY WORDS**

human; liver

---

**REFERENCE**

Surapaneni,S.S.; Clay,M.P.; Spangle,L.A.; Paschal,J.W.; Lindstrom,T.D. In vitro biotransformation and identification of human cytochrome P450 isozyme-dependent metabolism of tazofelone, *Drug Metab.Dispos.*, **1997**, *25*, 1383–1388.

---


# Tegafur

**Molecular formula:** C<sub>8</sub>H<sub>9</sub>FN<sub>2</sub>O<sub>3</sub>

**Molecular weight:** 200.17

**CAS Registry No.:** 17902-23-7

**Merck Index:** 9267



---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Mix 20  $\mu$ L plasma with 50  $\mu$ L 500 mM pH 8.0 phosphate buffer, 100  $\mu$ L 500 ng/mL 5-chlorouracil in water, and 7 mL ethyl acetate. Shake for 30 min and centrifuge at 2200 g for 10 min. Remove the ethyl acetate layer and evaporate it. Reconstitute the residue in 500  $\mu$ L n-hexane:mobile phase 40:60 with sonication for 10 min. Inject the whole amount.

---

**HPLC VARIABLES**

**Column:** 100  $\times$  4.6 Develosil 60-3 (Nomura Chemical, Japan)

**Mobile phase:** n-Hexane:ethyl acetate:formic acid:water 50:0.5:0.5:0.3

**Flow rate:** 0.9

**Injection volume:** 500

**Detector:** UV 264

---

**CHROMATOGRAM**

**Internal standard:** 5-chlorouracil

**Limit of detection:** 50 ng/mL

---

**KEY WORDS**

plasma; rat; normal phase

---

**REFERENCE**

Fuse,E.; Takai,K.; Okuno,K.; Kobayashi,S. Hepatic extraction ratio of 5-fluorouracil in rats. Dose dependence and effect of uracil and interleukin-2, *Biochem.Pharmacol.*, **1996**, *52*, 561–568.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Dilute 250  $\mu$ L plasma or urine 1:10 with water, add 100  $\mu$ L 200 mM pH 7.0 phosphate buffer, add 100  $\mu$ L 10  $\mu$ g/mL  $\beta$ -hydroxyethyltheophylline, add 4 mL dichloromethane, shake for 10 min, centrifuge at 2000 g for 5 min, transfer the organic layer to another

tube, repeat the extraction. Evaporate the combined organic layers to dryness under a stream of nitrogen at 40°, dissolve the residue in 200 µL mobile phase, inject a 40 µL aliquot.

#### HPLC VARIABLES

**Column:** 150 × 4.6 5 µm Inertsil ODS-2 (GL Sciences)

**Mobile phase:** MeOH:10 mM pH 5.5 phosphate buffer 15:85

**Flow rate:** 1

**Injection volume:** 40

**Detector:** UV 270

#### CHROMATOGRAM

**Retention time:** 6.95

**Internal standard:** β-hydroxyethyltheophylline (10.53)

**Limit of quantitation:** 10 ng/mL (plasma), 100 ng/mL (urine)

#### KEY WORDS

plasma

#### REFERENCE

Matsushima,E.; Yoshida,K.; Kitamura,R.; Yoshida,K.-. Determination of S-1 (combined drug of tegafur, 5-chloro-2,4-dihydroxypyridine and potassium oxonate) and 5-fluorouracil in human plasma and urine using high-performance liquid chromatography and gas-chromatography-negative ion chemical ionization mass spectrometry, *J.Chromatogr.B*, **1997**, 691, 95–104.


# Teicoplanin

**Molecular formula:** C<sub>88</sub>H<sub>95</sub>Cl<sub>2</sub>N<sub>9</sub>O<sub>33</sub>(A<sub>2</sub>-1)

**Molecular weight:** 1877.68 (A<sub>2</sub>-1)

**CAS Registry No.:** 61036-64-4,  
61036-62-2

**Merck Index:** 9269



Teicoplanin A<sub>2</sub>-1 (Z)-4-decanoic acid  
A<sub>2</sub>-2 8-methylnonanoic acid  
A<sub>2</sub>-3 n-decanoic acid  
A<sub>2</sub>-4 8-methyldcanoic acid  
A<sub>2</sub>-5 9-methyldcanoic acid

#### SAMPLE

**Matrix:** blood

**Sample preparation:** Add 1 mL MeCN to 500 µL plasma, vortex briefly, centrifuge at 2000 g for 3 min, add 2 mL chloroform (Caution! Chloroform is a carcinogen!) to the supernatant, vortex, inject an aliquot of the aqueous supernatant layer.

#### HPLC VARIABLES

**Column:** 250 × 4 5 µm LiChrocart RP8

**Mobile phase:** MeCN:20 mM pH 4.4 ammonium acetate buffer 26:74

**Flow rate:** 1.3

**Injection volume:** 30

**Detector:** UV 220

---

#### CHROMATOGRAM

**Retention time:** 6.5

**Limit of detection:** 30 ng/mL

**Limit of quantitation:** 90 ng/mL

---

#### KEY WORDS

plasma

---

#### REFERENCE

Cociglio,M.; Peyrière,H.; Hillaire-Buys,D.; Alric,R. Application of a standardized coextractive cleanup procedure to routine high-performance liquid chromatography assays of teicoplanin and ganciclovir in plasma, *J.Chromatogr.B*, **1998**, 705, 79–85.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** Condition a Bond Elut SAX SPE cartridge with 4 mL MeOH, 1 mL water, and 1 mL 10 mM n-heptanesulfonic acid. 1 mL Plasma + 1 mL 10 mM n-heptanesulfonic acid, vortex, centrifuge at 4000 g for 2 min, add the supernatant to the SPE cartridge, wash with 3 mL 10 mM heptanesulfonic acid, elute with 1 mL MeOH, inject a 20 µL aliquot of the eluate.

---

#### HPLC VARIABLES

**Guard column:** 10 × 4.2 10 µm LiChrosorb RP-8

**Column:** 150 × 4.2 5 µm LiChrosorb RP-8

**Mobile phase:** MeOH:water 5:95 containing 10 mM disodium n-heptanesulfonate, adjusted to pH 4.0 with 2 g/L sodium acetate and glacial acetic acid

**Column temperature:** 30

**Flow rate:** 0.8

**Injection volume:** 20

**Detector:** UV 240

---

#### CHROMATOGRAM

**Retention time:** 6.5

**Limit of detection:** 200-400 ng/mL

---

#### KEY WORDS

SPE; plasma

---

#### REFERENCE

Georgopoulos,A.; Czejka,M.J.; Starzengruber,N.; Jäger,W.; Lackner,H. High-performance liquid chromatographic determination of teicoplanin in plasma: comparison with a microbiological assay, *J.Chromatogr.*, **1989**, 494, 340–346.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** Add plasma to a Bond Elut C8 SPE cartridge, wash with water, elute with MeCN:pH 6 buffer 20:80, inject a 50 µL aliquot of the eluate into a 500 µL sample loop filled with water, inject 200 µL water into the sample loop, operate the valve.

---

#### HPLC VARIABLES

**Column:** 100 × 1 5 µm Hypersil ODS

**Mobile phase:** MeCN:25 mM pH 6.0 NaH<sub>2</sub>PO<sub>4</sub> 20:80

**Flow rate:** 0.05

**Injection volume:** 50

**Detector:** UV 210

---

#### CHROMATOGRAM

**Retention time:** 3 (TA3-1), 10 (TA2-1), 13 (TA2-2), 16 (TA2-3), 30 (TA2-4), 33 (TA2-5)



**Limit of detection:** 50 ng/mL

---

**KEY WORDS**

SPE; plasma

---

**REFERENCE**

Taylor,R.B.; Reid,R.G.; Gould,I.M. Determination of teicoplanin in plasma using microbore high-performance liquid chromatography and injection-generated gradients, *J.Chromatogr.*, **1991**, 563, 451-457.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Add 1-10 mL plasma or urine to the affinity column, wash with 1 mL buffer, elute with three 500  $\mu$ L portions of 1.5% ammonium hydroxide. Add 100  $\mu$ L 18% HCl to the eluate, inject a 170  $\mu$ L aliquot. (Affinity columns consisted of 0.5 mL D-alanyl-D-alanine- $\epsilon$ -aminocaproyl Sepharose CL-6B stored in 100 mM pH 8.5 Tris buffer containing 0.004% sodium merthiolate, equilibrate with buffer before use. After use wash with 5 mL buffer containing 0.004% sodium merthiolate, store at 4°. Buffer was 50 mM  $\text{NaH}_2\text{PO}_4$  containing 200 mM NaCl adjusted to pH 7.5 with 1 M NaOH. Prepare the column material by coupling D-alanyl-D-alanine to activated CH-Sepharose 4B using  $\epsilon$ -aminocaproic acid according to the manufacturer's instructions. 30 mg D-Alanyl-D-alanine in 5 mL 100 mM sodium bicarbonate buffer containing 500 mM NaCl (pH 8) was coupled to 3 mL gel in 1 h, block unreacted ester groups with 1 M ethanolamine hydrochloride at pH 9 for 1 h, wash repeatedly with 100 mM pH 4 sodium acetate buffer containing 500 mM NaCl and 100 mM pH 8 tris-HCl buffer containing 500 mM NaCl (Appl. Biochem. Biotechnol. 1985, 11, 101.))

---

**HPLC VARIABLES**

**Guard column:** 30  $\times$  4.6 5  $\mu$ m Nucleosil C18

**Column:** 100  $\times$  4.6 5  $\mu$ m Nucleosil C18

**Mobile phase:** Gradient. A was 10 mM pH 4.9  $\text{NaH}_2\text{PO}_4$ . B was MeCN:10 mM pH 4.9  $\text{NaH}_2\text{PO}_4$  50:50. A:B 78:22 of 0.5 min, to 47:53 over 44 min, maintain at 47:53 for 2 min, to 0:100 over 1.5 min, maintain at 0:100 for 5 min, return to initial conditions over 2 min, re-equilibrate for 5 min.

**Flow rate:** 1.3

**Injection volume:** 170

**Detector:** UV 240

---

**CHROMATOGRAM**

**Retention time:** 9.1 (A3), 28.8 (A2-1), 31.6 (A2-2), 32.8 (A2-3), 37.1 (A2-4), 38.0 (A2-5)

**Limit of detection:** 100 ng/mL (A2-2)

---

**OTHER SUBSTANCES**

**Noninterfering:** acebutolol, acetaminophen, aspirin, atenolol, benfluorex, chlorazepam, clofibrate, diazepam, furosemide, hydrochlorothiazide, nitrazepam, oxprenolol

---

**KEY WORDS**

plasma; SPE; pharmacokinetics

---

**REFERENCE**

Riva,E.; Ferry,N.; Cometti,A.; Cuisinaud,G.; Gallo,G.G.; Sassard,J. Determination of teicoplanin in human plasma and urine by affinity and reversed-phase high-performance liquid chromatography, *J.Chromatogr.*, **1987**, 421, 99-110.

---

**SAMPLE**

**Matrix:** urine

**Sample preparation:** Adjust pH of 1 mL urine to 7.4, add 1 mL of a homogeneous suspension of D-alanyl-D-alanine-epsilon-aminocaproyl-Sepharose resin:50 mM pH 7.4 phosphate buffer containing 200 mM NaCl 50:50, let stand at 4° overnight, place on top of a column containing 6 mL of a homogeneous suspension of D-alanyl-D-alanine-epsilon-aminocaproyl-Sepharose resin:50 mM pH 7.4 phosphate buffer containing 200 mM NaCl 50:50, wash with 6 mL 50 mM pH 7.4 phosphate buffer containing 200 mM NaCl, elute with four 3 mL portions of 1.5% aqueous ammonia, immediately neutralize the eluate with 100 mM HCl, inject an aliquot.

**HPLC VARIABLES**

**Column:** 250 × 4.6 5 µm Ultrasphere ODS

**Mobile phase:** Gradient. A was 25 mM pH 6.0 sodium phosphate buffer. B was MeCN:water 90:10. A:B from 84:16 to 55:45 over 2 h, wash with 49:51 for 5 min.

**Column temperature:** 20

**Flow rate:** 1.5

**Detector:** UV 254

**CHROMATOGRAM**

**Retention time:** 8 (A3), 19 (A2-1), 20 (A2-2), 22 (A2-3,3a), 27 (A2-4), 28 (A2-5,5a)

**Limit of detection:** 0.5-5 µg/mL

**OTHER SUBSTANCES**

**Extracted:** metabolites

**KEY WORDS**

rat; SPE

**REFERENCE**

Zerilli, L.F.; Cavenaghi, L.; Bernareggi, A.; Assandri, A. Teicoplanin metabolism in rats, *Antimicrob. Agents Chemother.*, **1989**, 33, 1791-1794.

# Temafloracin


**Molecular formula:** C<sub>21</sub>H<sub>18</sub>F<sub>3</sub>N<sub>3</sub>O<sub>3</sub>

**Molecular weight:** 417.39

**CAS Registry No.:** 108319-06-8, 105784-61-0 (HCl)

**Merck Index:** 9284

**Lednicer No.:** 5 125

**SAMPLE**

**Matrix:** bile, blood, urine

**Sample preparation:** Dilute urine 1:20. Dilute bile 1:10. 500 µL Serum, diluted urine, or diluted bile + 3.2 mL dichloromethane, vortex, rotate at 20 rpm for 10 min, centrifuge at 1000 g for 10 min. Remove 3 mL of the lower organic phase and add it to 200 µL 100 mM NaOH, rotate at 20 rpm for 10 min, centrifuge at 1000 g for 10 min, inject a 20 µL aliquot of the aqueous layer.

**HPLC VARIABLES**

**Column:** 150 × 4.6 5 µm Ultrasphere C18

**Mobile phase:** MeCN:buffer 19:81, pH adjusted to 2 with 14.6 M phosphoric acid (Buffer was 10 mM NaH<sub>2</sub>PO<sub>4</sub> containing 5 mM tetrabutylammonium bromide.)

**Flow rate:** 2

**Injection volume:** 20

**Detector:** F ex 275 em 450

**CHROMATOGRAM**

**Retention time:** 3

**Limit of detection:** 100 ng/mL (bile), 200 ng/mL (urine), 10 ng/mL (serum)

**OTHER SUBSTANCES**

**Noninterfering:** amikacin, aztreonam, carbamazepine, cephalosporins, ciprofloxacin, clavulanic acid, difloxacin, digitoxin, digoxin, fleroxacin, fosfomycin, furosemide, gentamycin, imipenem, lidocaine, netilmicin, norfloxacin, ofloxacin, pefloxacin, penicillins, phenobarbital, phenytoin, primidone, procainamide, quinidine, rifampin, salicylic acid, teicoplanin, theophylline, tobramycin, valproic acid, vancomycin

---

**KEY WORDS**

serum; human; rabbit

---

**REFERENCE**

Koechlin,C.; Jehl,F.; Linger,L.; Monteil,H. High-performance liquid chromatography for the determination of three new fluoroquinolones, fleroxacin, temafloxacin and A-64730, in biological fluids, *J.Chromatogr.*, **1989**, *491*, 379-387.

---

**SAMPLE**

**Matrix:** blood, tissue

**Sample preparation:** Homogenize (Ultra-Turrax T25) mouse lung in 1-3 mL pH 6.8 Soerensen phosphate buffer, centrifuge. Add 1 µg sparfloxacin to serum or lung homogenate supernatant, extract using a Bond Elut C2 SPE cartridge, inject a 100 µL aliquot of the extract.

---

**HPLC VARIABLES**

**Column:** 150 × 4.6 5 µm Ultrabase C8 (SFCC, Neuilly Plaisance, France)

**Mobile phase:** MeCN:MeOH:5% acetic acid 15:10:75

**Flow rate:** 1

**Injection volume:** 100

**Detector:** UV 280

---

**CHROMATOGRAM**

**Retention time:** 6

**Internal standard:** sparfloxacin

**Limit of detection:** 15 ng

---

**KEY WORDS**

serum; lung; mouse; pharmacokinetics; SPE

---

**REFERENCE**

Vallée,E.; Azoulay-Dupuis,E.; Bauchet,J.; Pocidal,J.-J. Kinetic disposition of temafloxacin and ciprofloxacin in a murine model of pneumococcal pneumonia. Relevance for drug efficacy, *J.Pharmacol.Exp.Ther.*, **1992**, *262*, 1203-1208.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Plasma. Mix equal volumes of plasma and 100-500 ng/mL IS in MeCN: 75 mM pH 7.4 phosphate buffer containing 0.5% sodium dodecyl sulfate 30:70, filter (Amicon Centrifree with YMT membrane) while centrifuging at less than 1000 g for 20 min, inject an aliquot of the ultrafiltrate. Alternatively, mix 400 µL plasma and 400 µL 500 mM pH 7 phosphate buffer, add 6 mL dichloromethane:EtOH 90:10, shake slowly horizontally for 10 min, centrifuge in a refrigerated centrifuge at 900 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of air at 50-60°, reconstitute the residue in mobile phase, inject an aliquot. Urine. Dilute urine 20-100 fold with mobile phase containing IS, inject an aliquot. Hydrolyze conjugates by heating urine in 1 M NaOH at 60° for 30 min, neutralize, dilute 20-100 fold with mobile phase containing IS, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 7 µm Adsorbosphere HS C18

**Mobile phase:** MeCN:water 53:47 containing 40 mM phosphoric acid, 10 mM NaH<sub>2</sub>PO<sub>4</sub>, 0.2% sodium dodecyl sulfate, and 5 mM N-acetylhydroxamic acid

**Flow rate:** 1.5

**Injection volume:** 50

**Detector:** F ex 280 em 389 (or UV 280, Antimicrob. Agents Chemother. 1991, 35, 436)

---

**CHROMATOGRAM**

**Retention time:** 7.5

**Internal standard:** 1-(4-bromophenyl)-6-fluoro-1,4-dihydro-7-(4-methyl-1-piperazinyl)-oxo-3-quinoline carboxylic acid (A57084) (8.5)

**Limit of quantitation:** 1 ng/mL (extraction), 10 ng/mL (ultrafiltrate)

---

**OTHER SUBSTANCES**

Extracted: metabolites

---

**KEY WORDS**

plasma; ultrafiltrate; pharmacokinetics

---

**REFERENCE**

Granneman, G.R.; Varga, L.L. High-performance liquid chromatographic procedures for the determination of temafloracin in biological matrices, *J. Chromatogr.*, **1991**, 568, 197–206.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** 1 mL Serum or 500  $\mu$ L urine + 500  $\mu$ L 1 M ammonium acetate + 5 mL dichloromethane, shake vigorously for 10 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 500  $\mu$ L MeOH:thionyl chloride 90:10, heat at 60° for 30 min, evaporate to dryness under a stream of nitrogen at 40°, add 500  $\mu$ L 500 mM sulfuric acid, add 1 mL hexane, shake vigorously for 10 min, centrifuge for 5 min. Remove the aqueous layer and add it to 1 mL 2 M sodium carbonate and 3 mL ether, shake vigorously for 10 min, centrifuge for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 20  $\mu$ L 2 mg/mL (S)-(-)-N-1-(2-naphthylsulfonyl)-2- pyrrolidinecarbonyl chloride in dichloromethane, add 5  $\mu$ L triethylamine, let stand at room temperature for 10 min. Evaporate to dryness under a stream of nitrogen, reconstitute the residue in 500  $\mu$ L dichloromethane, inject a 50  $\mu$ L aliquot. (No details of synthesis of chiral reagent given in paper.)

---

**HPLC VARIABLES**

**Column:** 150  $\times$  4.6 Zorbax SIL

**Mobile phase:** Hexane:methyl acetate:MeOH:aqueous ammonia 150:100:10:1

**Flow rate:** 0.8

**Injection volume:** 50

**Detector:** UV 280

---

**CHROMATOGRAM**

**Retention time:** 18.5 (S), 19.5 (R)

**Limit of detection:** 5 ng/mL

---

**KEY WORDS**

serum; chiral; pharmacokinetics; derivatization; normal phase

---

**REFERENCE**

Matsuoka, M.; Banno, K.; Sato, T. Analytical chiral separation of a new quinolone compound in biological fluids by high-performance liquid chromatography, *J. Chromatogr. B*, **1996**, 676, 117–124.

---

**SAMPLE**

**Matrix:** bulk

**Sample preparation:** Dissolve 100 mg bulk drug in 25 mL MeCN:water 50:50 containing 1 mL 1 M NaOH, add 15 mL 20 mM  $\text{KH}_2\text{PO}_4$  adjusted to pH 2.4 with orthophosphoric acid, make up to 50 mL with MeCN:water 50:50, inject a 50  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 150  $\times$  4.6 5  $\mu$ m Nucleosil C18

**Mobile phase:** Gradient. MeCN:THF:buffer from 5:5:90 to 5:60:35 over 50 min, maintain at 5:60:35 for 10 min, return to initial conditions over 5 min, re-equilibrate for 15 min. (Buffer was 20 mM  $\text{KH}_2\text{PO}_4$  adjusted to pH 2.4 with orthophosphoric acid.)

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 325

---

**CHROMATOGRAM**

**Retention time:** 11

**Limit of detection:** 0.05% (of temafloracin)

---

---

**OTHER SUBSTANCES**

**Simultaneous:** degradation products, impurities

---

**REFERENCE**

Elrod,L.,Jr.; Linton,C.L.; Shelat,B.P.; Wong,C.F. Determination of minor impurities in temafloxacin hydrochloride by high-performance liquid chromatography, *J.Chromatogr.*, **1990**, 519, 125–136.

---

**SAMPLE**

**Matrix:** cells

**Sample preparation:** Incubate cells in 2 mL 100 mM pH 3.0 glycine-HCl buffer for 2 h at room temperature, centrifuge at 5600 g for 5 min, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** Bondapak C18

**Mobile phase:** MeCN:25 mM phosphoric acid adjusted to pH 3.0 with tetrabutylammonium hydroxide 25:75

**Flow rate:** 1.5

**Detector:** F ex 340 em 425

---

**OTHER SUBSTANCES**

**Also analyzed:** ciprofloxacin, fleroxacin, lomefloxacin, norfloxacin, ofloxacin

---

**REFERENCE**

Pascual,A.; Garcia,I.; Conejo,M.C.; Perea,E.J. Fluorometric and high-performance liquid chromatographic measurement of quinolone uptake by human neutrophils, *Eur.J.Clin.Microbiol.Infect.Dis.*, **1991**, 10, 969–971.

---

**SAMPLE**

**Matrix:** hair

**Sample preparation:** Hair. Cut 10 hairs into 2 mm sections, weigh, add 1 mL 1 M NaOH, heat at 60° for 1 h, cool, add 500  $\mu$ L 2 M HCl, add 1 mL 200 mM pH 7.0 phosphate buffer, add 500  $\mu$ L 4–20 ng/mL IS in MeOH, add 5 mL chloroform, shake for 30 min, centrifuge at 2000 rpm for 10 min. Remove 4 mL of the organic layer and evaporate it to dryness by aspiration at 40° for 30 min, reconstitute the residue in 200–500  $\mu$ L mobile phase, inject a 100  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 150  $\times$  4.6 TSKgel ODS-80TM (Tosoh)

**Mobile phase:** MeCN:50 mM citric acid:1 M ammonium acetate 22:78:1

**Flow rate:** 1

**Injection volume:** 100

**Detector:** F ex 280 em 406

---

**CHROMATOGRAM**

**Retention time:** 12

**Internal standard:** 1-(4-bromophenyl)-6-fluoro-1,4-dihydro-7-(4-methyl-1-piperazinyl)-oxo-3-quinoline carboxylic acid (A57084) (14)

**Limit of detection:** 0.5 ng/mL

---

**REFERENCE**

Uematsu,T.; Kondo,K.; Yano,S.; Yamaguchi,T.; Umemura,K.; Nakashima,M. Measurement of temafloxacin in human scalp hair as an index of drug exposure, *J.Pharm.Sci.*, **1994**, 83, 42–45.

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a 450  $\mu$ g/mL solution in MeCN:water 50:50. 5 mL Solution + 5 mL THF + 200 molar excess of acetic anhydride + 3 molar excess of 1 M NaOH, sonicate for 15 min, add 15 mL mobile phase, sonicate for 15 min, cool to room temperature, make up to 50 mL with mobile phase, inject a 50  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 150  $\times$  4.6 5  $\mu$ m Nucleosil C18

**Mobile phase:** MeCN:buffer 35:65 (Buffer was prepared by mixing equal volumes of 20 mM citric acid and 20 mM sodium citrate, pH adjusted to 2.4 with perchloric acid.)

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 280

---

#### CHROMATOGRAM

**Retention time:** 21.4

---

#### OTHER SUBSTANCES

**Simultaneous:** ciprofloxacin, norfloxacin, sarafloxacin

---

#### KEY WORDS

derivatization

---

#### REFERENCE

Morley, J.A.; Elrod, L., Jr. Determination of fluoroquinolone antibacterials as N-Acyl derivatives, *Chromatographia*, 1993, 37, 295–299.

---

# Temazepam


**Molecular formula:**  $C_{16}H_{13}ClN_2O_2$

**Molecular weight:** 300.74

**CAS Registry No.:** 846-50-4

**Merck Index:** 9285

**Lednicer No.:** 2 402



---

#### SAMPLE

**Matrix:** bile, blood, gastric contents, tissue, urine

**Sample preparation:** Chop 5-g tissue and homogenize (Ultra Turrax T25) at 8500, 9500, 13500, 20500, and 24000 rpm for 1 min each. Add homogenate to 20 mL water. Dilute blood, urine, gastric contents, and bile four times with water. Mix 4 mL sample with 10  $\mu$ L 1 mg/mL prazepam and 1 mL pH 7.4 phosphate buffer, vortex briefly, add 4 mL diethyl ether and mix for 15 min (Spiramix 10, Denley, UK). Separate the organic layer, add 4 mL diethyl ether to extraction sample, mix. Evaporate combined organic layers to dryness under a stream of dry air at 50°. Purify extracts by partition between 1 mL MeCN and 2 mL heptane, separate MeCN layer, evaporate it to dryness, reconstitute the residue in 100  $\mu$ L MeOH and inject a 20  $\mu$ L aliquot of the solution.

---

#### HPLC VARIABLES

**Guard column:** 20  $\times$  4.6 5  $\mu$ m Apex II ODS

**Column:** 150  $\times$  4.6 5  $\mu$ m Apex II ODS

**Mobile phase:** MeCN:MeOH:10 mM phosphoric acid:10 mM Na<sub>2</sub>HPO<sub>4</sub> 40:20:36:4

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 240

---

#### CHROMATOGRAM

**Retention time:** 4.8

**Internal standard:** prazepam (14.5)

**Limit of quantitation:** 100 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** diazepam, nitrazepam, oxazepam

---

#### KEY WORDS

liver; lung; muscle; urine; pericardial fluid

---

**REFERENCE**

Pounder,D.J.; Adams,E.; Fuke,C.; Langford,A.M. Site to site variability of postmortem drug concentrations in liver and lung, *J.Forensic Sci.*, **1996**, *41*, 927-932.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 200 mg Extrelut + 400  $\mu$ L blood + 100  $\mu$ L 50  $\mu$ g/mL IS in MeOH, mix, let dry at room temperature for 1-2 h. Add to a 30  $\times$  4.6 stainless steel extraction column, extract with carbon dioxide:ethyl acetate 95:5 at 2 mL/min, 65°, and 300 psi. for 10 min, collect by expansion into MeOH. Dry the collected extract at 65° under nitrogen. Reconstitute the residue in 50  $\mu$ L mobile phase. Inject a 20  $\mu$ L aliquot. Condition ca. 10 g Extrelut in a 10 mL plastic syringe with dichloromethane. Add 250  $\mu$ L 5% ammonia to the top. Mix 900  $\mu$ L blood with 100  $\mu$ L 50  $\mu$ g/mL IS in MeOH. Add 1 mL pH 4 phosphate buffer and 250  $\mu$ L 5% ammonia solution, mix thoroughly, add to the extraction column. After 5 min elute with diethyl ether under the influence of gravity. Collect 8 mL eluate, evaporate to dryness at 65° under nitrogen. Reconstitute the residue in 180  $\mu$ L mobile phase, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Guard column:** 20  $\times$  4.6 5  $\mu$ m Hypersil ODS

**Column:** 250  $\times$  4.6 5  $\mu$ m Hypersil ODS

**Mobile phase:** MeOH:Na<sub>2</sub>HPO<sub>4</sub> 70:30

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 3

**Internal standard:** prazepam (11.5)

---

**OTHER SUBSTANCES**

**Also analyzed:** diazepam, chlordiazepoxide, nordiazepam, oxazepam

---

**KEY WORDS**

SFE; SPE; whole blood

---

**REFERENCE**

Scott,K.S.; Oliver,J.S. Development of a supercritical fluid extraction method for the determination of temazepam in whole blood, *J.Anal.Toxicol.*, **1997**, *21*, 297-300.

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50  $\mu$ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood)  $\mu$ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18

**Column:** 250  $\times$  4.6 5  $\mu$ m Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 200.5


**CHROMATOGRAM****Retention time:** 18.562**KEY WORDS**

whole blood

**REFERENCE**

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

# Temocillin

**Molecular formula:**  $C_{16}H_{18}N_2O_7S_2$ **Molecular weight:** 414.46**CAS Registry No.:** 66148-78-5**Merck Index:** 9288**Lednicer No.:** 4 178**SAMPLE****Matrix:** blood

**Sample preparation:** 350  $\mu$ L Serum + 150  $\mu$ L water + 250  $\mu$ L 400 mM HCl + 3.5 mL chloroform: n-amyl alcohol (3:1), mix for 5 min, centrifuge for 5 min. Remove the organic layer and add it to 350  $\mu$ L 10 mM pH 7.0 phosphate buffer, mix for 5 min, centrifuge for 5 min, inject a 20  $\mu$ L aliquot of the upper aqueous layer.

**HPLC VARIABLES****Column:** 300  $\times$  3.9 10  $\mu$ m  $\mu$ Bondapak C18

**Mobile phase:** MeOH:buffer 15:85 (Buffer was 100 mM ammonium acetate adjusted to pH 4.0 with glacial acetic acid.)

**Flow rate:** 1.8**Injection volume:** 20**Detector:** UV 242**CHROMATOGRAM****Retention time:** 5.4**Internal standard:** temocillin**OTHER SUBSTANCES****Extracted:** cefoxitin, cefuroxime, cephalothin, ticarcillin

**Noninterfering:** acetaminophen, acetazolamide, allopurinol, amikacin, ampicillin, azlocillin, caffeine, cefamandole, cefoperazone, cefotaxime, cefsulodin, ceftazidime, ceftizoxime, chloramphenicol, chlorpromazine, clindamycin, dicloxacillin, 5-fluorocytosine, flurazepam, gentamicin, methicillin, metronidazole, mezlocillin, moxalactam, nafcillin, penicillin, phenobarbital, piperacillin, procainamide, rifampin, sulfamethoxazole, theophylline, thienamycin, tobramycin, trimethoprim, vancomycin

**KEY WORDS**

serum; temocillin is IS

**REFERENCE**

Shull,V.H.; Dick,J.D. Determination of ticarcillin levels in serum by high-pressure liquid chromatography, *Antimicrob.Agents Chemother.*, **1985**, 28, 597–600.

**SAMPLE****Matrix:** bulk



**Sample preparation:** Prepare a 300 µg/mL solution in 100 mM pH 7.0 phosphate buffer, inject a 20 µL aliquot.

#### HPLC VARIABLES

**Column:** 300 × 3.9 µBondapak C18

**Mobile phase:** MeOH:buffer 10:90 (Buffer was 15.6 g NaH<sub>2</sub>PO<sub>4</sub>·2H<sub>2</sub>O in 900 mL water, adjust pH to 7.0 with NaOH, make up to 1 L with water.)

**Flow rate:** 2

**Injection volume:** 20

**Detector:** UV 230

#### CHROMATOGRAM

**Retention time:** 6.0 (R), 7.5 (S)

#### OTHER SUBSTANCES

**Simultaneous:** isomers, diastereomers, impurities, degradation products, ticarcillin

#### REFERENCE

Bird, A.E.; Charsley, C.H.; Jennings, K.R.; Marshall, A.C. High-performance liquid chromatographic assay of temocillin and epimerisation of its diastereoisomers, *Analyst*, **1984**, 109, 1209–1212.

#### SAMPLE

**Matrix:** solutions

#### HPLC VARIABLES

**Column:** 300 × 3.9 µBondapak C18

**Mobile phase:** MeOH:100 mM pH 6.5 potassium phosphate 10:90

**Flow rate:** 2

**Injection volume:** 50

**Detector:** UV 254

#### KEY WORDS

radiolabelled compounds

#### REFERENCE

Morecombe, D.J. High-efficiency preparative-scale reversed-phase high-performance liquid chromatographic purification of <sup>14</sup>C-labelled antibiotics, *J.Chromatogr.*, **1987**, 389, 389–395.


# Teniposide

**Molecular formula:** C<sub>32</sub>H<sub>32</sub>O<sub>13</sub>S

**Molecular weight:** 656.66

**CAS Registry No.:** 29767-20-2

**Merck index:** 9291



#### SAMPLE

**Matrix:** blood

**Sample preparation:** Filter 1 mL plasma (Centrifree micropartition device, molecular mass cut-off 30000, Amicon, USA) using a 33° fixed angle centrifuge (Beckman Model G56R) at 2000 g and 25° for 30 min. Add 1 mL chloroform to the ultrafiltrate (Caution! Chloroform is a carcinogen!), agitate slowly for 20 min, centrifuge at 1000 g for 5 min. Evaporate the organic phase to dryness under vacuum at 40°, dissolve the dry extract in 50 µL MeOH, inject a 25 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 300 × 3.9 10 μm 125 Å μBondapak Phenyl (Waters)

**Mobile phase:** MeCN:water:glacial acid 35:64:1

**Flow rate:** 1

**Injection volume:** 25

**Detector:** F ex 288 em 328

---

**CHROMATOGRAM**

**Retention time:** 18

**Internal standard:** teniposide

---

**OTHER SUBSTANCES**

**Extracted:** etoposide

**Noninterfering:** alizapride, doxorubicin, furosemide, idarubicin, ranitidine, vinblastine, vinorelbine

---

**KEY WORDS**

plasma; ultrafiltrate; teniposide is IS

---

**REFERENCE**

Robieux,I.; Aita,P.; Sorio,R.; Toffoli,G.; Boiocchi,M. Determination of unbound etoposide concentration in ultrafiltered plasma by high-performance liquid chromatography with fluorimetric detection, *J.Chromatogr.B*, **1996**, 686, 35–41.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 10 μL 1 mg/mL etoposide in MeOH + 5 mL chloroform, rock gently for 15 min, centrifuge. Remove 4.5 mL of the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, dissolve residue in 50 μL MeOH, vortex, centrifuge for 5-10 min, inject a 20 μL aliquot.

---

**HPLC VARIABLES**

**Column:** 300 × 3.9 10 μm μBondapak C18

**Mobile phase:** MeOH:water 60:40

**Flow rate:** 1-1.2

**Injection volume:** 20

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 7.5

**Internal standard:** etoposide (5)

**Limit of quantitation:** 500 ng/mL

---

**KEY WORDS**

plasma

---

**REFERENCE**

Strife,R.J.; Jardine,I.; Colvin,M. Analysis of the anticancer drugs VP 16-213 and VM 26 and their metabolites by high-performance liquid chromatography, *J.Chromatogr.*, **1980**, 182, 211–220.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 10 μL 100 μg/mL etoposide in MeOH + 5 mL chloroform, rock gently for 15 min, centrifuge. Remove 4.5 mL of the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, dissolve residue in 50 μL MeOH, vortex, centrifuge for 5-10 min, inject a 25 μL aliquot.

---

**HPLC VARIABLES**

**Guard column:** 70 × 2.1 30 μm Co:Pell (Whatman)

**Column:** 300 × 3.9 10 μm μBondapak C18

**Mobile phase:** MeOH:water 60:40

---

Flow rate: 1  
Injection volume: 25  
Detector: F ex 215 em 328

---

**CHROMATOGRAM**

Retention time: 8  
Internal standard: etoposide (5.5)  
Limit of detection: 25 ng/mL  
Limit of quantitation: 50 ng/mL

---

**KEY WORDS**

plasma

---

**REFERENCE**

Strife, R.J.; Jardine, I.; Colvin, M. Analysis of the anticancer drugs etoposide (VP 16-213) and teniposide (VM 26) by high-performance liquid chromatography with fluorescence detection, *J. Chromatogr.*, **1981**, 224, 168-174.

---

**SAMPLE**

Matrix: blood

**Sample preparation:** Mix plasma or serum with an equal volume of proteinase K, let stand for 10 min. (Alternatively, heat serum or plasma with an equal volume of 1 mg/mL subtilisin A at 50° for 15 min.) Inject 1.6 mL hydrolyzed blood or filtered serum on to column A with mobile phase A at 1 mL/min, backflush column A with mobile phase A to waste for 2 min at 2 mL/min, backflush the contents of column A on to column B with mobile phase B, after 30 s remove column A from the circuit, elute column B with mobile phase B, monitor the effluent from column B. (Clean column A by backflushing with MeOH at 2 mL/min for 3 min then forward flush with water at 1 mL/min for 2 min.)

---

**HPLC VARIABLES**

**Column:** A 2 × 4.6 10 µm PRP-1 divinylbenzene-styrene copolymer (Hamilton); B 125 × 4 10 µm LiChrosorb RP-18

**Mobile phase:** A water; B MeOH:water 55:45

**Flow rate:** A 1-2; B 1

**Injection volume:** 1600

**Detector:** F ex 230 em 328 following post-column extraction. The column effluent mixed with dichloroethane pumped at 0.6 mL/min, the mixture flowed through a 2 mm i.d. glass reactor (Technicon) to a phase separator (Technicon (*J. Chromatogr.* 1979, 185, 473)) with a PTFE insert and 0.3 mL/min of the organic phase flowed through the detector.

---

**CHROMATOGRAM**

Retention time: 8

---

**OTHER SUBSTANCES**

Extracted: etoposide

---

**KEY WORDS**

column-switching; post-column extraction; plasma; serum

---

**REFERENCE**

Werkhoven-Goewie, C.E.; Brinkman, U.A.T.; Frei, R.W.; de Ruiter, C.; de Vries, J. Automated liquid chromatographic analysis of the anti-tumorigenic drugs etoposide (VP 16-213) and teniposide (VM 26), *J. Chromatogr.*, **1983**, 276, 349-357.

---

**SAMPLE**

Matrix: blood

**Sample preparation:** 1 mL Plasma + 5 µg etoposide + 3 mL chloroform, shake for 10 min, centrifuge at 1500 g for 10 min, repeat extraction. Combine the organic layers and evaporate a 5 mL aliquot to dryness under a stream of nitrogen at 50°, reconstitute the residue in 200 µL mobile phase, inject a 20 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 300 × 3.9 10 µm µBondapak C18

**Mobile phase:** MeOH:250 mM ammonium acetate:acetic acid 54:45:1

**Flow rate:** 1.5

**Injection volume:** 20

**Detector:** E, Bioanalytical Systems LC4, TL5 glassy carbon electrode, + 900 mV, Ag/AgCl reference electrode

---

**CHROMATOGRAM**

**Retention time:** 8.3

**Internal standard:** etoposide (4.1)

**Limit of detection:** 10 ng/mL

**Limit of quantitation:** 20 ng/mL

---

**KEY WORDS**

plasma; pharmacokinetics

---

**REFERENCE**

Canal,P.; Michel,C.; Bugat,R.; Soula,G.; Carton,M. Quantification of teniposide in human serum by high-performance liquid chromatography with electrochemical detection, *J.Chromatogr.*, **1986**, 375, 451–456.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 450 µL Plasma + 50 µL 380 mM sodium dodecyl sulfate in 59 mM pH 7 sodium phosphate buffer + 5 µL 500 ng/mL etoposide in MeOH, sonicate for 5 min. Inject a 100 µL aliquot onto column A with mobile phase A, elute with mobile phase A for 7.5 min, elute the contents of column A onto column B with mobile phase B for 1 min. After 1 min remove column A from the circuit and re-equilibrate it with mobile phase A for 1.5 min, monitor the effluent from column B.

---

**HPLC VARIABLES**

**Column:** A 10 × 2.1 40 µm Chromsep C18 (Chrompack); B 300 × 4.6 10 µm µBondapak phenyl

**Mobile phase:** A 10 mM pH 7.0 sodium phosphate; B MeOH:10 mM pH 7.0 sodium phosphate buffer 55:45

**Flow rate:** A 0.4; B 1

**Injection volume:** 100

**Detector:** E, +500 mV vs Ag/AgCl or UV 254

---

**CHROMATOGRAM**

**Retention time:** 6.4

**Internal standard:** etoposide (4.2)

**Limit of detection:** 20 ng/mL (electrochemical), 150 ng/mL (UV)

---

**KEY WORDS**

plasma; column-switching

---

**REFERENCE**

van Opstal,M.A.J.; van der Horst,F.A.L.; Holthuis,J.J.M.; Van Bennekom,W.P.; Bult,A. Automated reversed-phase chromatographic analysis of etoposide and teniposide in plasma by using on-line surfactant-mediated sample clean-up and column-switching, *J.Chromatogr.*, **1989**, 495, 139–151.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 50 µL 100 µg/mL etoposide in MeOH, vortex, add 2 mL dichloroethane, shake thoroughly for 1 min, centrifuge at 3000 g for 5 min. Remove 1.5 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 150 µL MeOH:water 70:30, sonicate for 6 min, inject a 10 µL aliquot.

---

**HPLC VARIABLES**

**Guard column:** 10 × 4.6 10 µm LiChrosorb C18

**Column:** 100 × 4.6 10 µm Novapak phenyl

---

**Mobile phase:** MeOH:10 mM pH 7.0 phosphate buffer 55:45

**Flow rate:** 0.7

**Injection volume:** 10

**Detector:** E, Metrohm Model 641 VA, EA 286/1 glassy carbon electrode + 500 mV, stainless-steel auxiliary electrode, Ag/AgCl reference electrode

---

#### CHROMATOGRAM

**Retention time:** 7

**Internal standard:** etoposide (4.5)

**Limit of detection:** 10 ng/mL

---

#### KEY WORDS

plasma

---

#### REFERENCE

van der Horst, F.A.L.; van Opstal, M.A.J.; Teeuwssen, J.; Post, M.H.; Holthuis, J.J.M.; Brinkman, U.A.T. Comparative study on the determination of the anti-neoplastic drug teniposide in plasma using micellar liquid chromatography and surfactant-mediated plasma clean-up, *J. Chromatogr.*, **1991**, 567, 161-174.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** Sonicate 50 million leukemic cells in 1 mL phosphate buffered saline. 500  $\mu$ L Plasma or 1 mL sonicated cells + 2 mL chloroform, mix. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 200  $\mu$ L MeOH: water 50:50, sonicate for 5 min, inject a 100  $\mu$ L aliquot. To measure non-protein-bound etoposide filter (Amicon Centrifree) while centrifuging at 20°, inject a 100-200  $\mu$ L aliquot of the ultrafiltrate.

---

#### HPLC VARIABLES

**Column:** 250  $\times$  4.6 5  $\mu$ m Spherisorb Phenyl

**Mobile phase:** MeOH:water:acetic acid 45:54:1

**Flow rate:** 1

**Injection volume:** 100-200

**Detector:** F ex 220 em 330

---

#### CHROMATOGRAM

**Retention time:** 9.3

**Internal standard:** teniposide

---

#### OTHER SUBSTANCES

**Extracted:** etoposide, cis-etoposide

---

#### KEY WORDS

plasma; cells; ultrafiltrate; teniposide is IS

---

#### REFERENCE

Liliemark, E.; Petterson, B.; Peterson, C.; Liliemark, J. High-performance liquid chromatography with fluorometric detection for monitoring of etoposide and its *cis*-isomer in plasma and leukaemic cells, *J. Chromatogr. B*, **1995**, 669, 311-317.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** 400  $\mu$ L Plasma + 1.6 mL diethyl ether, vortex, centrifuge at 1000 g for 2 min. Remove the organic layer and evaporate it to dryness under a stream of air at 30°, reconstitute the residue in 200  $\mu$ L mobile phase, inject a 100  $\mu$ L aliquot.

---

#### HPLC VARIABLES

**Column:** 250  $\times$  4 5  $\mu$ m LiChrospher 100 RP-18

**Mobile phase:** MeCN:MeOH:buffer 30:20:50 (Buffer was 20 g/L NaH<sub>2</sub>PO<sub>4</sub> containing 0.8 g/L heptanesulfonic acid, pH adjusted to 3.0 with orthophosphoric acid.)

**Flow rate:** 1

**Injection volume:** 100

**Detector:** E, Environmental Sciences Coulochem 5100 A, guard cell +0.90 V (before injector), clean-up cell +0.40 V, detection cell +0.90 V

---

#### CHROMATOGRAM

**Retention time:** 10.6

**Internal standard:** teniposide

---

#### OTHER SUBSTANCES

**Extracted:** vinorelbine

---

#### KEY WORDS

plasma; rabbit; teniposide is IS

---

#### REFERENCE

Mouchard-Delmas,C.; Gourdier,B.; Vistelle,R. Determination of vinorelbine in rabbit plasma by high-performance liquid chromatography with coulometric detection, *J.Chromatogr.B*, **1995**, 663, 390–394.

---

#### SAMPLE

**Matrix:** blood, CSF, urine

**Sample preparation:** 500  $\mu$ L Plasma, urine, or CSF + 500  $\mu$ L saturated ammonium sulfate + 4 mL ethyl acetate + 10  $\mu$ L 100  $\mu$ g/mL etoposide, vortex for 5 min, centrifuge at 3000 rpm for 15 min, repeat the extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 200  $\mu$ L MeOH, inject an aliquot.

---

#### HPLC VARIABLES

**Guard column:** 10  $\mu$ m  $\mu$ Bondapak phenyl

**Column:** 250  $\times$  4.6 10  $\mu$ m  $\mu$ Bondapak phenyl

**Mobile phase:** MeCN:water:acetic acid 30:68:2

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 284 or E, Bioanalytical Systems LC-4A, 0.75 V

---

#### CHROMATOGRAM

**Retention time:** 18

**Internal standard:** etoposide (7)

**Limit of detection:** 20 ng/mL

**Limit of quantitation:** 50 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** metabolites

---

#### KEY WORDS

plasma; pharmacokinetics

---

#### REFERENCE

Sinkule,J.A.; Evans,W.E. High-performance liquid chromatographic analysis of the semisynthetic epipodophyllotoxins teniposide and etoposide using electrochemical detection, *J.Pharm.Sci.*, **1984**, 73, 164–168.

---

#### SAMPLE

**Matrix:** solutions

**Sample preparation:** Prepare a solution in MeCN:water:650 mM pH 4.0 sodium citrate 40:60:0.4, inject a 30  $\mu$ L aliquot.

---

#### HPLC VARIABLES


**Guard column:** 10  $\times$  3.9 10  $\mu$ m  $\mu$ Bondapak phenyl

**Column:** 150  $\times$  3.9 10  $\mu$ m  $\mu$ Bondapak phenyl

**Mobile phase:** MeCN:water:650 mM pH 4.0 sodium citrate 35:57.3:7.7. When run is over wash with MeCN:water:650 mM pH 4.0 sodium citrate 70:27.3:7.7 for 4 min, re-equilibrate with initial mobile phase for 9 min.

**Flow rate:** 2**Injection volume:** 30**Detector:** E, ESA 5100A detector, Model 5020 guard cell between pump and autosampler +0.7 V, Model 5011 dual electrode analytical cell, upstream (screening) electrode +0.2 V, downstream electrode +0.45 V against Ag/AgCl**CHROMATOGRAM****Retention time:** 5.8**OTHER SUBSTANCES****Simultaneous:** etoposide**REFERENCE**Eisenberg, E.J.; Eickhoff, W.M. Determination of etoposide in blood by liquid chromatography with electrochemical detection, *J.Chromatogr.*, **1993**, 621, 110–114.

# Tenoxicam

**Molecular formula:** C<sub>13</sub>H<sub>11</sub>N<sub>3</sub>O<sub>4</sub>S<sub>2</sub>**Molecular weight:** 337.38**CAS Registry No.:** 59804-37-4**Merck Index:** 9293**Lednicer No.:** 4 173**SAMPLE****Matrix:** blood**Sample preparation:** 1 mL Plasma + 1 mL pH 2 phosphate buffer + 10 mL diethyl ether, vortex for 1 min, centrifuge at 1300 g for 5 min. Evaporate the organic layer to dryness under a stream of nitrogen. Reconstitute the residue in 100 µL 10 mM HCl in MeOH, inject a 40 µL aliquot.**HPLC VARIABLES****Column:** 300 × 3.9 10 µm µBondapak ODS**Mobile phase:** MeOH:10 mM pH 2 phosphate buffer 45:55**Flow rate:** 1.5**Injection volume:** 40**Detector:** UV 361**CHROMATOGRAM****Retention time:** 5.81**Internal standard:** tenoxicam**KEY WORDS**

plasma; rat; tenoxicam is IS

**REFERENCE**Amanlou, M.; Dehpour, A.R. Rapid method for the determination of piroxicam in rat plasma using high-performance liquid chromatography, *J.Chromatogr.B*, **1997**, 696, 317–319.**SAMPLE****Matrix:** blood**Sample preparation:** Mix plasma with 800 µL 5 mM pH 4 phosphate buffer. Centrifuge at 2500 rpm for 5 min. Add the supernatant to Extrelut-1 cartridge (Merck, Darmstadt, Germany). Elute with 10 mL dichloromethane, dry eluate under a stream of nitrogen at 35°, dissolve the residue in 300 µL mobile phase, centrifuge and inject 100 µL aliquot of the supernatant.

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 5 µm Hypersil ODS

**Mobile phase:** MeOH:50 mM pH 6 phosphate buffer 50:50

**Column temperature:** 35

**Flow rate:** 1.3

**Injection volume:** 100

**Detector:** UV 371

---

**CHROMATOGRAM**

**Retention time:** 2.2

**Internal standard:** tenoxicam

---

**KEY WORDS**

plasma; tenoxicam is IS

---

**REFERENCE**

Bareggi,S.R.; Gambaro,V.; Valenti,M.; Benvenuti,C. Absorption of oral lornoxicam in healthy volunteers using a granular formulation in comparison with standard tablets, *Arzneimittelforschung*, **1997**, *47*, 755–757.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 1 mL 1 M HCl + 1 mL water + 100 µL 2 µg/mL piroxicam + 100 µL MeOH + 5 mL dichloromethane, mix for 5 min, centrifuge at 1800 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 30°, reconstitute the residue in 100 µL mobile phase, inject a 40 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 125 × 4 5 µm LiChrospher 100 RP-18

**Mobile phase:** MeOH:100 mM pH 7.4 phosphate buffer 40:60

**Flow rate:** 1.1

**Injection volume:** 40

**Detector:** UV 355

---

**CHROMATOGRAM**

**Retention time:** 3.4

**Internal standard:** piroxicam (4.5)

**Limit of detection:** 5 ng/mL

**Limit of quantitation:** 10 ng/mL

---

**KEY WORDS**

plasma

---

**REFERENCE**

Múnera-Jaramillo,M.I.; Botero-Garcés,S. Determination of tenoxicam in plasma by high-performance liquid chromatography, *J.Chromatogr.*, **1993**, *616*, 349–352.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 10 µg ketorolac in water + 100 µL 5% zinc sulfate in water, vortex for 2 min, add 440 µL buffer, vortex for 1 min, centrifuge at 2000 g for 10 min, inject a 100 µL aliquot of the supernatant. (Buffer was 100 mM NaH<sub>2</sub>PO<sub>4</sub> and 10 mM sodium lauryl sulfate, pH adjusted to 2.8 with phosphoric acid.)

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 5 µm Nucleosil C18

**Mobile phase:** MeCN:water 35:65 containing 10 mM NaH<sub>2</sub>PO<sub>4</sub> and 1 mM sodium lauryl sulfate, pH adjusted to 2.8 with phosphoric acid

**Flow rate:** 1.5

**Injection volume:** 100

**Detector:** UV 355



**CHROMATOGRAM****Retention time:** 4.6**Internal standard:** ketorolac (10.3)**Limit of detection:** 40 ng/mL

---

**KEY WORDS**plasma; protect from light; pharmacokinetics

---

**REFERENCE**

Mason, J.L.; Hobbs, G.J. Simple method for the analysis of tenoxicam in human plasma using high-performance liquid chromatography, *J. Chromatogr. B*, **1995**, 665, 410-415.

---

**SAMPLE****Matrix:** blood

**Sample preparation:** 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

---

**HPLC VARIABLES****Column:** 300 × 3.9 4 µm NovaPack C18

**Mobile phase:** MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH<sub>2</sub>PO<sub>4</sub> adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

**Column temperature:** 30**Flow rate:** 0.8**Injection volume:** 50**Detector:** UV 376

---

**CHROMATOGRAM****Retention time:** 2.99**Limit of detection:** <120 ng/mL

---

**KEY WORDS**

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vandesine; mexiletine; dipyrindamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; metha-

done; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpi-pramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

---

**REFERENCE**

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254–262.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 150  $\mu$ L Plasma + 300  $\mu$ L 1 mg/mL piroxicam in MeCN, mix, centrifuge at 3500 rpm for 20 min, inject a 10  $\mu$ L aliquot of the supernatant.

---

**HPLC VARIABLES**

**Column:** 125  $\times$  4 Nucleosil C18

**Mobile phase:** MeCN:water:acetic acid 58:38:4

**Flow rate:** 1

**Injection volume:** 10

**Detector:** UV 365

---

**CHROMATOGRAM**

**Internal standard:** piroxicam

**Limit of detection:** 200 ng/mL

---

**KEY WORDS**

rat; pharmacokinetics; plasma

---

**REFERENCE**

Troconiz,J.F.; Lopez-Bustamante,L.G.; Fos,D. Tenoxicam pharmacokinetics in rats: A population model, *J.Pharm.Sci.*, **1995**, *84*, 1482–1487.

---

**SAMPLE**

**Matrix:** blood, dialysate

**Sample preparation:** Blood, plasma. 150  $\mu$ L Blood or plasma + 300  $\mu$ L 1 mg/mL piroxicam in MeCN, centrifuge at 3500 rpm for 20 min, inject a 10  $\mu$ L aliquot of the supernatant. Dialysate. 100  $\mu$ L Dialysate + 200  $\mu$ L 1 mg/mL piroxicam in MeCN, centrifuge at 3500 rpm for 20 min, inject a 10  $\mu$ L aliquot of the supernatant.

---

**HPLC VARIABLES**

**Column:** 125  $\times$  4 Nucleosil C18

**Mobile phase:** MeCN:water:acetic acid 55:38:4

**Flow rate:** 1

**Injection volume:** 10

**Detector:** UV 365

---

**CHROMATOGRAM**

**Retention time:** 1.5

**Internal standard:** piroxicam (2.2)

**Limit of detection:** 200 ng/mL

---

**KEY WORDS**

rat; whole blood; plasma; pharmacokinetics

---

**REFERENCE**

Lopez-Bustamante,L.G.; Troconiz,J.I.; Fos,D. Tenoxicam: acute dose-dependent disposition studies in rats, *J.Pharm.Sci.*, **1993**, *82*, 851–853.

**SAMPLE****Matrix:** blood, synovial fluid**Sample preparation:** Plasma. 500  $\mu$ L Plasma + 500  $\mu$ L 500 mM pH 4 phosphate buffer, mix, homogenize, centrifuge at 900 g. Add 1 mL of the mixture to an Extrelut-1 SPE cartridge. Elute with two 5 mL portions of dichloromethane. Evaporate the eluate under a stream of nitrogen at 35°. Dissolve the residue in 120  $\mu$ L mobile phase by vortexing, inject an aliquot. Plasma, synovial fluid. Condition a 100 mg C18 SPE cartridge (Phenomenex) with 1 mL MeOH, 1 mL water, and 500  $\mu$ L 500 mM pH 2 phosphate buffer. 500  $\mu$ L Plasma or synovial fluid + 500  $\mu$ L 500 mM pH 2 phosphate buffer. Mix, homogenize, centrifuge at 2500 rpm. Add 1 mL of the clean supernatant to the SPE cartridge, wash with 1 mL water, dry with 2 mL air, elute with 1.25 mL MeCN:25% ammonia 90:10. Evaporate under reduced pressure at 35°. Reconstitute the residue in 100  $\mu$ L MeOH:100 mM pH 8 phosphate buffer 50:50, inject an aliquot.

---

**HPLC VARIABLES****Column:** 250  $\times$  4.6 5  $\mu$ m Hypersil ODS**Mobile phase:** MeOH:100 mM pH 6 sodium dihydrogenphosphate buffer 50:50**Flow rate:** 1.5**Injection volume:** 50-100**Detector:** UV 372

---

**CHROMATOGRAM****Retention time:** 2.6-2.8**Internal standard:** tenoxicam

---

**OTHER SUBSTANCES****Extracted:** lornoxicam

---

**KEY WORDS**plasma; mouse; rat; rabbit; dog; monkey; human; SPE; lornoxicam is IS

---

**REFERENCE**Radhofer-Welte,S.; Dittrich,P. Determination of the novel non-steroidal anti-inflammatory drug lornoxicam and its main metabolite in plasma and synovial fluid, *J.Chromatogr.B*, **1998**, *707*, 151-159.

---

**SAMPLE****Matrix:** blood, urine**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50  $\mu$ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood)  $\mu$ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES****Guard column:** 20 mm long Symmetry C18**Column:** 250  $\times$  4.6 5  $\mu$ m Symmetry C8 (Waters)**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.**Column temperature:** 30**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)**Injection volume:** 10-30**Detector:** UV 200.5

---

**CHROMATOGRAM****Retention time:** 12.733

---

**KEY WORDS**

whole blood

---

**REFERENCE**

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:** 125 × 3 Ecocart LiChrospher 100 RP-18

**Mobile phase:** Isopropanol:100 mM KH<sub>2</sub>PO<sub>4</sub>:formic acid 54:100:0.1

**Flow rate:** 0.6

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 1.5

**Limit of quantitation:** 200-500 ng/mL

---

**OTHER SUBSTANCES**

**Simultaneous:** acemetacin; diclofenac; flurbiprofen; indomethacin; lonazolac; ketoprofen; naproxen; piroxicam; sulindac

---

**REFERENCE**

Baeyens,W.R.G.; Van Der Weken,G.; Van Overbeke,A.; Zhang,Z.D. Preliminary results on the LC-separation of non-steroidal anti-inflammatory agents in conventional and narrow-bore RP set-ups applying columns with different internal diameters, *Biomed.Chromatogr.*, **1995**, 9, 261–262.

---

# Terazosin


**Molecular formula:** C<sub>19</sub>H<sub>25</sub>N<sub>5</sub>O<sub>4</sub>

**Molecular weight:** 387.44

**CAS Registry No.:** 65390-64-7, 70024-40-7  
(HCl dihydrate), 63074-08-8 (HCl)

**Merck Index:** 9297

**Lednicer No.:** 3 194



---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 100 µL 100 ng/mL IS in water, vortex briefly, add 1 mL 0.9% NaCl and 100 µL 2 M NaOH, vortex briefly. Add 5 mL pentane:dichloromethane 50:50, mix at 40 rpm for 20 min, centrifuge at 500 g for 10 min, freeze the lower aqueous layer in a acetone-dry ice bath, evaporate the organic layer to dryness under a stream of nitrogen. Re-constitute the residue in 100 µL mobile phase, vortex vigorously for 15 s, inject a 65 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 10 µm Chiralpak AD (Chiral Technologies, USA)

**Mobile phase:** Hexane:2-propanol:diethylamine 70:30:0.1

**Column temperature:** 30

**Flow rate:** 1

**Injection volume:** 65

**Detector:** F em 238 ex 370

---

**CHROMATOGRAM**

**Internal standard:** (+)-glaucine

**Limit of quantitation:** 500 pg/mL

**OTHER SUBSTANCES****Interfering:** prazosin**KEY WORDS**

plasma; chiral

**REFERENCE**

Zavitsanos,A.P.; Alebic-Kolbah,T. Enantioselective determination of terazosin in human plasma by normal phase high-performance liquid chromatography--electrospray mass spectrometry, *J.Chromatogr.A*, **1998**, 794, 45-56.

**SAMPLE****Matrix:** blood

**Sample preparation:** 1 mL Plasma + 100  $\mu$ L 200 ng/mL IS in water, vortex briefly, add 1 mL 0.9% NaCl and 100  $\mu$ L 2 M NaOH, vortex briefly. Add 5 mL pentane:dichloromethane 50:50, mix at 40 rpm for 20 min, centrifuge at 500 g for 10 min, freeze the lower aqueous layer in a acetone-dry ice bath, evaporate the organic layer to dryness under a stream of nitrogen. Re-constitute the residue in 70  $\mu$ L hexane:2-propanol 90:10, vortex vigorously for 15 s, inject a 20  $\mu$ L aliquot.

**HPLC VARIABLES****Column:** 100  $\times$  2.1 10  $\mu$ m Chiralpak AD (Chiral Technologies, USA)**Mobile phase:** Hexane:2-propanol containing 0.05% diethylamine 65:35**Flow rate:** 0.15**Injection volume:** 20

**Detector:** MS, HP 1100 electrospray, positive ion mode, drying gas nitrogen 310°, 10 L/min, nebulizer pressure 60 p.s.i.g., quadropole temperature 100°, capillary voltage 4500 V, SIM, m/z 388.2, post-column solvent addition of isopropanol:5 mM ammonium acetate 75:25


**CHROMATOGRAM****Retention time:** 4.8, 5.8 (enantiomers)**Internal standard:** prazosin (4.8)**Limit of quantitation:** 62.5 pg/mL**KEY WORDS**

plasma; chiral

**REFERENCE**

Zavitsanos,A.P.; Alebic-Kolbah,T. Enantioselective determination of terazosin in human plasma by normal phase high-performance liquid chromatography--electrospray mass spectrometry, *J.Chromatogr.A*, **1998**, 794, 45-56.

# Terbinafine

**Molecular formula:** C<sub>21</sub>H<sub>25</sub>N**Molecular weight:** 291.44**CAS Registry No.:** 91161-71-6**Merck Index:** 9299**Lednicer No.:** 4 55**SAMPLE****Matrix:** blood, tissue

**Sample preparation:** Adjust pH of 100  $\mu$ L plasma or 500  $\mu$ g tissue to 9.0, add IS, extract with hexane. Extract the organic layer with sulfuric acid:isopropanol, inject an aliquot of this solution.

---

**HPLC VARIABLES****Column:** 150 × 2 Ultrasphere C18**Mobile phase:** MeCN:pH 4.0 phosphate buffer 40:60**Detector:** UV 224

---

**CHROMATOGRAM****Limit of detection:** 1.1 ng/g (fetal tissue), 2.23 ng/mL (plasma)

---

**KEY WORDS**rabbit; plasma; placenta; fetus

---


**REFERENCE**

Gries, W.J.; Wan, W.; Matos, F.J.; de Meireles, J.C.; Pimplaskar, H.K.; Sileno, A.P.; Romeo, V.D.; Xia, W.J.; Behl, C.R. A specific and sensitive method for quantitating buprenorphine hydrochloride in a nasal solution (Abstract 2517), *Pharm. Res.*, **1997**, *14*, S381.

---

---

# Terbutaline

**Molecular formula:** C<sub>12</sub>H<sub>19</sub>NO<sub>3</sub>**Molecular weight:** 225.29**CAS Registry No.:** 23031-25-6, 23031-32-5 (sulfate)**Merck Index:** 9302

---

**SAMPLE****Matrix:** blood

**Sample preparation:** Condition a Bond-Elut C18 SPE cartridge with MeCN and water. 1 mL Plasma + 1 mL 20 mM pH 9.0 Na<sub>2</sub>HPO<sub>4</sub>, mix, add to the SPE cartridge, wash with 5 column volumes of water, dry the SPE cartridge for 5 min, wash with 3 mL dichloromethane:n-butanol 97:3, elute with two 1 mL portions of 0.09% HCl in MeCN, evaporate the eluate to dryness under a stream of nitrogen at 37°, dissolve the residue in 300 µL mobile phase, inject a 200 µL aliquot. (To hydrolyze 3-O-metaproterenol sulfate mix 1 mL plasma and 200 µL 2 µg/mL terbutaline sulfate, add 1 mL 6% trichloroacetic acid, vortex, centrifuge at 1000 g for 15 min. Remove the supernatant and add it to 200 µL 2 M HCl, heat at 65° for 90 min, cool, adjust pH to 10 with 400 µL 2 M carbonate buffer, proceed as above.)

---

**HPLC VARIABLES****Column:** 250 × 4.9 5 µm Spherisorb C8**Mobile phase:** MeCN:buffer:water 4:1.5:94.5**Flow rate:** 1.8**Injection volume:** 200**Detector:** F ex 200 em 300 (cut-off filter)

---

**CHROMATOGRAM****Retention time:** 14.8**Internal standard:** terbutaline

---

**OTHER SUBSTANCES****Extracted:** metaproterenol

---

**KEY WORDS**SPE; plasma; terbutaline is IS

---

**REFERENCE**

Selinger, K.; Hill, H.M.; Matheou, D.; Dehelean, L. Determination of free and total metaproterenol in human plasma by high-performance liquid chromatography with fluorimetric detection, *J. Chromatogr.*, **1989**, *493*, 230–238.

**SAMPLE****Matrix:** blood**Sample preparation:** Inject 1 mL plasma onto column A with mobile phase A, wash with mobile phase A for 1.7 min, backflush contents of column A onto column B with mobile phase B, monitor the effluent from column B.

---

**HPLC VARIABLES****Column:** A  $10 \times 1.5$  37  $\mu\text{m}$  Corasil C18; B C18 guard column +  $250 \times 4.6$  10  $\mu\text{m}$  Spherisorb ODS**Mobile phase:** A water; B MeOH:67 mM pH 5 phosphate buffer:40 g/L sodium dodecyl sulfate: diethylamine 55:45:0.5:0.02 (before use condition column with 55:45:5:0.05)**Flow rate:** 1**Injection volume:** 1000**Detector:** E, EG & G Princeton Applied Research Model 400 EC, carbon fiber working electrode + 1.3 V, silver phosphate reference electrode (construction details given)

---

**CHROMATOGRAM****Retention time:** 7**Limit of detection:** 0.8 ng/mL

---

**KEY WORDS**plasma; column-switching

---

**REFERENCE**Sagar,K.A.; Kelly,M.T.; Smyth,M.R. Analysis of terbutaline in human plasma by high-performance liquid chromatography with electrochemical detection using a micro-electrochemical flow cell, *J.Chromatogr.*, **1992**, 577, 109–116.

---

**SAMPLE****Matrix:** blood**Sample preparation:** Condition a Bond Elut Si SPE cartridge by washing twice with 1 mL MeOH, twice with 1 mL water, and once with 1 mL 100 mM pH 9.2  $\text{K}_2\text{HPO}_4$ . Add 1 mL plasma + 100  $\mu\text{L}$  500 ng/mL atenolol in water, wash twice with 1 mL water, centrifuge at 1000 g for 5 min, elute with 1 mL MeOH. Evaporate MeOH to dryness at 40° under a stream of air and dissolve residue in 200  $\mu\text{L}$  mobile phase, inject an aliquot.

---

**HPLC VARIABLES****Column:**  $250 \times 4.6$  Spherisorb S5 SCX**Mobile phase:** MeOH:MeCN:water 40:40:20 containing 0.2% perchloric acid (apparent pH 1.7)**Flow rate:** 1.5**Injection volume:** 100**Detector:** F ex 200 no emission filter

---

**CHROMATOGRAM****Retention time:** 7**Internal standard:** atenolol (13)**Limit of detection:** 2500 ng/mL

---

**OTHER SUBSTANCES****Extracted:** albuterol**Noninterfering:** aminophylline, beclomethasone, cloprednol, dexamethasone, fenoterol, ipratropium bromide, methylprednisolone, orciprenaline, prednisolone, reproterol, rimiterol, salmeterol, sodium cromoglycate, theophylline

---

**KEY WORDS**plasma; SPE

---

**REFERENCE**McCarthy,P.T.; Atwal,S.; Sykes,A.P.; Ayres,J.G. Measurement of terbutaline and salbutamol in plasma by high performance liquid chromatography with fluorescence detection, *Biomed.Chromatogr.*, **1993**, 7, 25–28.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

---

**HPLC VARIABLES**

**Column:** 300 × 3.9 4 µm NovaPack C18

**Mobile phase:** MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH<sub>2</sub>PO<sub>4</sub> adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

**Column temperature:** 30

**Flow rate:** 0.8

**Injection volume:** 50

**Detector:** UV 225

---

**CHROMATOGRAM**

**Retention time:** 3.36

**Limit of detection:** <120 ng/mL

---

**KEY WORDS**

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylcegonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melfalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vandesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetraacaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cycizine; histapyrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenopropfen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thiopropazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

---

**REFERENCE**

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254–262.



**SAMPLE****Matrix:** blood, urine**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES****Guard column:** 20 mm long Symmetry C18**Column:** 250 × 4.6 5 µm Symmetry C8 (Waters)**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.**Column temperature:** 30**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)**Injection volume:** 10-30**Detector:** UV 200.5

---

**CHROMATOGRAM****Retention time:** 3.683

---

**KEY WORDS**whole blood

---

**REFERENCE**Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

---

**SAMPLE****Matrix:** formulations**Sample preparation:** Tablets, capsules. Mix tablets or capsules with 10 mL water, sonicate 30 min, centrifuge, inject an aliquot. Liquid formulations. Dilute liquid formulations with water, inject an aliquot.

---

**HPLC VARIABLES****Column:** 125 × 4 5 µm LiChrospher 100 RP-18 endcapped**Mobile phase:** MeOH:water 40:60 containing 2 mM KOH + 10 mM hexanoic acid**Flow rate:** 0.4**Injection volume:** 20**Detector:** UV 214

---

**CHROMATOGRAM****Limit of detection:** 1030 ng/mL

---

**OTHER SUBSTANCES****Also analyzed:** albuterol, fenoterol

---

**KEY WORDS**tablets; capsules; liquid formulations

---

**REFERENCE**Ackermans,M.T.; Beckers,J.L.; Everaerts,F.M.; Seelen,I.G. Comparison of isotachopheresis, capillary zone electrophoresis and high-performance liquid chromatography for the determination of salbutamol, terbutaline sulphate and fenoterol hydrobromide in pharmaceutical dosage forms, *J.Chromatogr.*, **1992**, 590, 341-353.

---

**SAMPLE**

**Matrix:** microsomal incubations

**Sample preparation:** Condition a C18 Sep-Pak SPE cartridge with three 3 mL portions of EtOH, two 3 mL portions of water, and with 3 mL 10 mM pH 7.5 phosphate buffer. Add the incubation mixture to the SPE cartridge, wash with two 3 mL portions of water, elute with two 1 mL portions of EtOH:50 mM pH 8.5 ammonium chloride buffer 95:5. Evaporate the eluate to dryness under a stream of nitrogen at 60°, reconstitute the residue in 500 µL of the initial mobile phase, inject a 200 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 150 × 5 Nucleosil 10SA

**Mobile phase:** Gradient. A was 250 mM pH 4.6 ammonium acetate buffer. B was MeCN:500 mM pH 4.6 ammonium acetate buffer 50:50. From A:B 90:10 to 10:90 over 20 min.

**Flow rate:** 1

**Injection volume:** 200

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 9.5

---

**OTHER SUBSTANCES**

**Extracted:** bambuterol

---

**KEY WORDS**

rat; SPE

---

**REFERENCE**

Lindberg,C.; Roos,C.; Tunek,A.; Svensson,L.Å. Metabolism of bambuterol in rat liver microsomes: identification of hydroxylated and demethylated products by liquid chromatography mass spectrometry, *Drug Metab.Dispos.*, **1989**, 17, 311–322.

---

---

**SAMPLE**

**Matrix:** perfusate, tissue

**Sample preparation:** Perfusate. 400 µL Lung perfusate + 400 µL 5% perchloric acid, mix, centrifuge, inject a 200 µL aliquot of the supernatant. Tissue. Homogenize lung in 2 volumes of water (Polytron Homogenizer), mix 400 µL homogenate with 400 µL 5% perchloric acid, mix, centrifuge, inject a 200 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 150 × 5 Nucleosil 10SA

**Mobile phase:** Gradient. A was 19.3 g ammonium acetate and 14.4 mL acetic acid in 1 L water. B was 19.3 g ammonium acetate and 14.4 mL acetic acid in 1 L MeCN:water 50:50. A:B from 90:10 to 10:90 over 20 min, stay at 10:90 for 3 min, return to initial conditions over 3 min, re-equilibrate for 7 min.

**Flow rate:** 1

**Injection volume:** 200

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:** 8

---

**OTHER SUBSTANCES**

**Extracted:** metabolites, bambuterol

---

**KEY WORDS**

lung; guinea pig

---

**REFERENCE**

Ryrfeldt,Å.; Nilsson,E.; Tunek,A.; Svensson,L.-Å. Bambuterol: uptake and metabolism in guinea pig isolated lungs, *Pharm.Res.*, **1988**, 5, 151–155.

**SAMPLE****Matrix:** solutions**Sample preparation:** Dilute with 5% dextrose, inject a 40  $\mu$ L aliquot.**HPLC VARIABLES****Column:** Waters microparticulate C18**Mobile phase:** MeOH:350 mM acetic acid and 5 mM sodium heptanesulfonate 35:65**Flow rate:** 1.6-2.0**Injection volume:** 40**Detector:** F ex 280 em 310**CHROMATOGRAM****Retention time:** 5.33**OTHER SUBSTANCES****Simultaneous:** theophylline, methyl dopate, isoproterenol**REFERENCE**Williams, D.A.; Fung, E.Y.Y.; Newton, D.W. Ion-pair high-performance liquid chromatography of terbutaline and catecholamines with aminophylline in intravenous solutions, *J.Pharm.Sci.*, **1982**, *71*, 956-958.**SAMPLE****Matrix:** solutions**Sample preparation:** Prepare a 10  $\mu$ g/mL solution in MeOH, inject a 20  $\mu$ L aliquot.**HPLC VARIABLES****Column:** 125  $\times$  4.9 Spherisorb S5W silica**Mobile phase:** MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7**Flow rate:** 2**Injection volume:** 20**Detector:** E, LeCarbone, V25 glassy carbon electrode, + 1.2 V**CHROMATOGRAM****Retention time:** 1.7**OTHER SUBSTANCES**

**Also analyzed:** acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipiparone, diprenorphine, dipyridamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine,

phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocinamide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

---

## REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, 323, 191–225.

---

## SAMPLE

**Matrix:** solutions

---

## HPLC VARIABLES

**Column:** 250 × 4.6 Chirex 3020 (Phenomenex)

**Mobile phase:** Hexane:1,2-dichloroethane:EtOH/trifluoroacetic acid 55:35:10 (EtOH/trifluoroacetic acid was premixed 20:1.)

**Flow rate:** 0.7-1

**Injection volume:** 20

**Detector:** UV 278

---

## KEY WORDS

chiral;  $\alpha = 1.28$  for enantiomers

---

## REFERENCE

Cleveland, T. Pirkle-concept chiral stationary phases for the HPLC separation of pharmaceutical racemates, *J. Liq. Chromatogr.*, **1995**, 18, 649–671.

---

## SAMPLE

**Matrix:** solutions

**Sample preparation:** Inject an aliquot of a 200  $\mu$ M solution in MeOH.

---

## HPLC VARIABLES

**Column:** 100 × 4.7 7  $\mu$ m Hypercarb (Shandon)

**Mobile phase:** MeOH containing 5 mM N-benzoyloxycarbonylglycyl-L-proline and 4.5 mM NaOH

**Column temperature:** 17

**Injection volume:** 20

**Detector:** UV 270

---

## CHROMATOGRAM

**Retention time:** k' 3.6 (first enantiomer)

---

## KEY WORDS

chiral;  $\alpha = 1.17$

---

## REFERENCE

Huynh, N.-H.; Karlsson, A.; Pettersson, C. Enantiomeric separation of basic drugs using N-benzoyloxycarbonylglycyl-L-proline as counter ion in methanol, *J. Chromatogr. A*, **1995**, 705, 275–287.

---

## SAMPLE

**Matrix:** solutions

---

## HPLC VARIABLES

**Column:** 250 × 4.6 5  $\mu$ m Supelcosil LC-DP (A) or 250 × 4 5  $\mu$ m LiChrospher 100 RP-8 (B)

**Mobile phase:** MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

**Flow rate:** 0.6

**Injection volume:** 25

**Detector:** UV 229

---

## CHROMATOGRAM

**Retention time:** 5.62 (A), 3.27 (B)

---

## OTHER SUBSTANCES

**Also analyzed:** acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordiazepoxide, chlorhexanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamine, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazinol, mefenamic acid, meperidine, mephénytoin, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyldopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfapyrazone, sulindac, temazepam, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocainide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone

---

## KEY WORDS

details of plasma extraction

---

## REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, 692, 103–119.

---

## SAMPLE

**Matrix:** solutions

**Sample preparation:** Inject a 20  $\mu$ L aliquot of an 8  $\mu$ g/mL solution.

---

## HPLC VARIABLES

**Column:** 250  $\times$  4.4  $\mu$ m Superspher 100 RP-18

**Mobile phase:** Buffer containing 10 mM  $\beta$ -cyclodextrin substituted with 2-hydroxy-3-trimethylammoniumpropyl groups (Roquette Frères, Lestrem, France) (Buffer was 0.8% triethylamine adjusted to pH 5.9 with acetic acid.)

**Column temperature:** 22.5

**Flow rate:** 0.8

**Injection volume:** 20

**Detector:** UV 275

---

**CHROMATOGRAM**

**Retention time:** 18.02, 19.88 (enantiomers)

---

**KEY WORDS**

chiral

---

**REFERENCE**

Roussel, C.; Favrou, A. Cationic  $\beta$ -cyclodextrin: a new versatile chiral additive for separation of drug enantiomers by high-performance liquid chromatography, *J.Chromatogr.A*, **1995**, 704, 67–74.

---

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Prepare a 1–10  $\mu\text{g/mL}$  solution in water, inject an aliquot.

---

**HPLC VARIABLES**

**Column:**  $250 \times 4.6$  5  $\mu\text{m}$  Hypersil SCX/C18

**Mobile phase:** MeCN:25 mM pH 3  $\text{Na}_2\text{HPO}_4$  50:50

**Injection volume:** 20

**Detector:** UV 254

---

**CHROMATOGRAM**

**Retention time:**  $k'$  3.39

---

**OTHER SUBSTANCES**

**Also analyzed:** amitriptyline, barbital, benzoic acid, butabarbital, clomipramine, clonazepam, desipramine, diazepam, flurazepam, furosemide, imipramine, nitrazepam, phenobarbital, phenol, phenolphthalein, pindolol, propranolol, resorcinol, salicylic acid, secobarbital, xylazine

---

**KEY WORDS**

effect of mobile phase pH on capacity factor is discussed

---

**REFERENCE**

Walshe, M.; Kelly, M.T.; Smyth, M.R.; Ritchie, H. Retention studies on mixed-mode columns in high-performance liquid chromatography, *J.Chromatogr.A*, **1995**, 708, 31–40.

---

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Column:**  $150 \times 4.6$  Zorbax phenyl

**Mobile phase:** MeCN:water 15:85 containing 10 mM  $\text{KH}_2\text{PO}_4$ , adjusted to pH 3.1 with phosphoric acid

**Flow rate:** 1.5

**Detector:** UV 212

---

**CHROMATOGRAM**

**Retention time:** 2.2

---

**REFERENCE**

Liu, P.; Bergstrom, T.K. Quantitative evaluation of aqueous isopropyl alcohol enhancement on skin flux of terbutaline (sulfate). 2. Permeability contributions of equilibrated drug species across human skin in vitro, *J.Pharm.Sci.*, **1996**, 85, 320–325.

---

---

**SAMPLE**

**Matrix:** tissue

**Sample preparation:** 300 mg Skin + propranolol, homogenize with 5 mL MeOH four times, combine the homogenates, filter. Evaporate the filtrate to dryness, reconstitute in mobile phase, inject a 5  $\mu\text{L}$  aliquot.

**HPLC VARIABLES**

**Column:** 250 × 4.6 5 µm Spherisorb cyano

**Mobile phase:** MeCN:pH 5.6 buffer 30:70

**Flow rate:** 1.4

**Injection volume:** 5

**Detector:** UV 225

**CHROMATOGRAM**

**Retention time:** 3.3

**Internal standard:** propranolol hydrochloride (7.9)

**Limit of detection:** 100 ng/mL

**KEY WORDS**

stability-indicating; skin

**REFERENCE**

Tenjarla,S.N.; Allen,R.; Mitchell,B. High-performance liquid chromatographic assay of terbutaline for preformulation studies, *J.Liq.Chromatogr.*, **1995**, *18*, 1603–1615.

**SAMPLE**

**Matrix:** urine

**Sample preparation:** Condition a Bond Elut C18 SPE cartridge with 3 volumes of MeOH and 2 volumes of water, dry under vacuum. Add 500 µL urine to the SPE cartridge, wash with 5 volumes of water, elute with 200 µL MeOH:50 mM pH 6 potassium phosphate buffer 50:50, add 50 µL 50 mM Na<sub>3</sub>PO<sub>4</sub> to the eluate, pass argon through the mixture, inject a 25 µL aliquot.

**HPLC VARIABLES**

**Column:** 300 mm long µBondapak phenyl

**Mobile phase:** MeCN:50 mM pH 5 phosphate buffer 6:94

**Flow rate:** 2.8

**Injection volume:** 25

**Detector:** F ex 280 em 310

**CHROMATOGRAM**

**Retention time:** 4.1

**Internal standard:** terbutaline

**OTHER SUBSTANCES**

**Extracted:** metaproterenol

**KEY WORDS**

SPE; protect from light; terbutaline is IS

**REFERENCE**

MacGregor,T.R.; Nastasi,L.; Farina,P.R.; Keirns,J.J. Isolation and characterization of metaproterenol-3-O-sulfate: a conjugate of metaproterenol in human urine, *Drug Metab.Dispos.*, **1983**, *11*, 568–573.

# Terfenadine


**Molecular formula:** C<sub>32</sub>H<sub>41</sub>NO<sub>2</sub>

**Molecular weight:** 471.68

**CAS Registry No.:** 50679-08-8

**Merck Index:** 9307

**Lednicer No.:** 4 48, 104

**SAMPLE**

**Matrix:** blood