

HPLC VARIABLES

Column: 250 × 4.6 5 µm Spherisorb cyano

Mobile phase: MeCN:pH 5.6 buffer 30:70

Flow rate: 1.4

Injection volume: 5

Detector: UV 225

CHROMATOGRAM

Retention time: 3.3

Internal standard: propranolol hydrochloride (7.9)

Limit of detection: 100 ng/mL

KEY WORDS

stability-indicating; skin

REFERENCE

Tenjarla,S.N.; Allen,R.; Mitchell,B. High-performance liquid chromatographic assay of terbutaline for preformulation studies, *J.Liq.Chromatogr.*, **1995**, *18*, 1603–1615.

SAMPLE

Matrix: urine

Sample preparation: Condition a Bond Elut C18 SPE cartridge with 3 volumes of MeOH and 2 volumes of water, dry under vacuum. Add 500 µL urine to the SPE cartridge, wash with 5 volumes of water, elute with 200 µL MeOH:50 mM pH 6 potassium phosphate buffer 50:50, add 50 µL 50 mM Na₃PO₄ to the eluate, pass argon through the mixture, inject a 25 µL aliquot.

HPLC VARIABLES

Column: 300 mm long µBondapak phenyl

Mobile phase: MeCN:50 mM pH 5 phosphate buffer 6:94

Flow rate: 2.8

Injection volume: 25

Detector: F ex 280 em 310

CHROMATOGRAM

Retention time: 4.1

Internal standard: terbutaline

OTHER SUBSTANCES

Extracted: metaproterenol

KEY WORDS

SPE; protect from light; terbutaline is IS

REFERENCE

MacGregor,T.R.; Nastasi,L.; Farina,P.R.; Keirns,J.J. Isolation and characterization of metaproterenol-3-O-sulfate: a conjugate of metaproterenol in human urine, *Drug Metab.Dispos.*, **1983**, *11*, 568–573.

Terfenadine


Molecular formula: C₃₂H₄₁NO₂

Molecular weight: 471.68

CAS Registry No.: 50679-08-8

Merck Index: 9307

Lednicer No.: 4 48, 104

**SAMPLE**

Matrix: blood

Sample preparation: 500 μ L Plasma + 50 μ L 250 ng/mL IS + 5 mL MTBE:dichloromethane:n-butyl chloride 3:2:1, vortex for 1 min, centrifuge at 2500 rpm for 1 min. Froze in an acetone/dry ice bath, transfer organic layer to another tube and extract plasma layer with 5 mL MTBE:dichloromethane:n-butyl chloride 3:2:1 for a second time. Evaporate combined organic layers to dryness under a stream of nitrogen at 40°. Reconstitute in 100 μ L MeCN:20 mM pH 3.5 ammonium acetate 50:50, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 2.0 5 μ m TSK gel ODS-80TS (Tosho, Japan)

Mobile phase: MeCN:1% formic acid:10 mM pH 4.0 ammonium acetate 85:13:2

Flow rate: 0.2

Injection volume: 10

Detector: MS, Fisons, VG Quattro I triple quadrupole, electrospray, source 200°, cone voltage 35 V, collision energy 30 eV, parent/daughter ions 472.2/436 for terfenadine, 470.2/203.2 for IS

CHROMATOGRAM

Retention time: 2.5

Internal standard: α -[4-(1,1-dimethylethyl)phenyl]-4-(hydroxydiphenylmethyl)-1-piperidinebutanone (2.8)

Limit of quantitation: 200 pg/mL

KEY WORDS

plasma

REFERENCE

Lau, Y.Y.; Anderson, P.H.; Talaat, R. High sensitivity high performance liquid chromatography electrospray tandem mass spectrometry determination of terfenadine in human plasma, *J. Liq. Chromatogr. Rel. Technol.*, **1996**, *19*, 2669–2679.

SAMPLE

Matrix: blood

Sample preparation: Add 100 μ L 20 ng/mL terfenadine- d_{10} to 1 mL heparinized plasma, vortex briefly, add 50 μ L 1 M ammonia solution and 4 mL hexane, shake for 5 min on a horizontal shaker, centrifuge at 2500 rpm for 5 min, freeze the aqueous layer in a dry ice-acetone bath, evaporate the organic phase to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100 μ L mobile phase, inject a 30 μ L aliquot.

HPLC VARIABLES

Guard column: 10 \times 2 5 μ m BDS Hypersil C18

Column: 50 \times 3 3 μ m BDS Hypersil C18

Mobile phase: MeCN:MeOH:10 mM ammonium acetate 54.4:32.6:13

Flow rate: 0.8

Injection volume: 30

Detector: MS, Perkin Elmer Sciex API-III, APCI interface, nebulizer 480°, nitrogen flow 1.2 L/min, interface heater 55°, multiplier 4 kV, m/z 472.7, 437

CHROMATOGRAM

Retention time: 1.10-1.15

Internal standard: terfenadine- d_{10} (1.1)

Limit of quantitation: 100 pg/mL

KEY WORDS

plasma

REFERENCE

Xu, A.; Linderholm, K.; Peng, L.; Hulse, J. Development and validation of an LC-MS-MS method for the determination of terfenadine in human plasma, *J. Pharm. Biomed. Anal.*, **1996**, *14*, 1675–1680.

SAMPLE

Matrix: blood, tissue

Sample preparation: Blood. Dilute 1 mL plasma with 100 μ L 1 M pH 7 phosphate buffer, add 100 μ L MeCN and 8 mL diethyl ether, extract. Evaporate the organic layer, dissolve the residue in 400 μ L mobile phase. Inject 200 μ L aliquot. Tissue. Homogenize the brain with 2 fold the weight of 1 M glycine-5 M NaOH pH 10 buffer. Extract 1500 μ L brain homogenate with 22 mL n-hexane. Evaporate the organic layer and dissolve the residue in 400-500 μ L mobile phase. Inject a 200 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 Intersil PH

Mobile phase: MeCN:0.018% TFA 32:68 (plasma), MeCN:0.018% TFA 37.5:62.5 (tissue)

Column temperature: 40

Flow rate: 0.7

Injection volume: 200

Detector: UV 220

CHROMATOGRAM

Limit of quantitation: 10 ng/mL (plasma), 25 ng/mL (brain)

KEY WORDS

brain; cat; mouse; pharmacokinetics; plasma; rat

REFERENCE

Kato,M.; Nishida,A.; Aga,Y.; Kita,J.; Kudo,Y.; Narita,H.; Endo,T. Pharmacokinetic and pharmacodynamic evaluation of central effect of the novel antiallergic agent betotastine besilate, *Arzneimittelforschung*, **1997**, 47, 1116-1124.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 19.118

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 500 μ L Microsomal incubation + 3 mL dichloromethane, add 100 μ L buffer containing 1 M sodium carbonate, 10 mM EDTA, and 2 M NaCl, vortex for 5 min, centrifuge at 2000 g for 10 min. Evaporate the organic phase to dryness under a stream of nitrogen. Reconstitute the residue with mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Supelco cyano

Mobile phase: MeCN:MeOH:12 mM pH 4.3 ammonium acetate buffer 30:30:40

Flow rate: 1.3

Detector: F ex 230 em 280

CHROMATOGRAM

Retention time: 21

OTHER SUBSTANCES

Simultaneous: metabolites

KEY WORDS

Caco-2-TC7 cells; intestine; liver; pharmacokinetics

REFERENCE

Raeissi,S.D.; Guo,Z.; Dobson,G.L.; Artursson,P.; Hidalgo,I.J. Comparison of CYP3A activities in a subclone of Caco-2 cells (TC7) and human intestine, *Pharm.Res.*, **1997**, *14*, 1019–1025.


Teroxirone

Molecular formula: C₁₂H₁₅N₂O₆

Molecular weight: 297.27

CAS Registry No.: 59653-73-5

Lednicer No.: 4 122



SAMPLE

Matrix: blood, urine

Sample preparation: 1 mL Whole blood, plasma, or urine + 500 μ L 10 mM pH 7.4 phosphate buffer + 5 mL chloroform, shake mechanically for 15 min, centrifuge at 3000 g. Remove 4 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 500 μ L 10 mM pH 7.4 phosphate buffer, add 500 μ L 5% diethyldithiocarbamate in water (freshly prepared), let stand at room temperature for 1 h, add 5 mL chloroform, shake, centrifuge. Remove 4 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 50–100 μ L toluene, inject an aliquot.

HPLC VARIABLES

Column: 10 μ m PAC-CN (Whatman)

Mobile phase: Gradient. Heptane:EtOH 60:40 for 10 min, to 10:90 over 6 min.

Flow rate: 2

Detector: UV 254

CHROMATOGRAM

Retention time: 7.7

Limit of detection: 15 ng/mL

KEY WORDS

derivatization; whole blood; plasma; normal phase; human; rabbit; pharmacokinetics

REFERENCE

Ames,M.M.; Kovach,J.S.; Rubin,J. Pharmacological characterization of teroxirone, a triepoxide antitumor agent, in rats, rabbits, and humans, *Cancer Res.*, **1984**, *44*, 4151–4156.


Terpin

Molecular formula: $C_{10}H_{20}O_2$

Molecular weight: 172.27

CAS Registry No.: 80-53-5

Merck Index: 9314



SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 312.8

CHROMATOGRAM

Retention time: 13.868

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.


Tertatolol

Molecular formula: $C_{16}H_{25}NO_2S$

Molecular weight: 295.45

CAS Registry No.: 34784-64-0

Merck Index: 9318



SAMPLE

Matrix: blood, urine

Sample preparation: Evaporate 25 μ L 1 μ g/mL (-)-alprenolol hydrochloride in EtOH into the bottom of a tube, add 1 mL plasma or urine, add 100 μ L 1 M NaOH, vortex for 1 min, add to an Extrelut SPE cartridge, elute with two 4 mL portions of diethyl ether. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute the residue in 100 μ L dichloromethane, add 10 μ L 0.1% S-(+)-naphthylethylisocyanate in dichloromethane, shake for 1 min, let stand

at room temperature for 12 h, add 10 μ L tert-butylamine, evaporate to dryness under a stream of nitrogen, reconstitute the residue in 20 μ L MeCN, vortex for 1 min, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 70 \times 4.6 3 μ m Ultrasphere XLODS

Mobile phase: MeCN:water 40:60

Flow rate: 2

Injection volume: 20

Detector: F ex 220 em 320

CHROMATOGRAM

Retention time: 15 (-), 17 (+)

Internal standard: (-)-alprenolol (13)

Limit of quantitation: 6 ng/mL

KEY WORDS

SPE; derivatization; chiral; plasma; pharmacokinetics

REFERENCE

Lave,T.; Efthymiopoulos,C.; Koffel,J.C.; Jung,L. Determination of tertatolol enantiomers in biological fluids by high-performance liquid chromatography, *J.Chromatogr.*, **1991**, 572, 203–210.

SAMPLE

Matrix: saliva

Sample preparation: Condition a 100 mg 1 mL Bond-Elut C2 SPE cartridge with 1 mL MeOH, 1 mL water, and 1 mL pH 9.0 borate buffer. Centrifuge a cotton roll soaked with saliva at 1000 g for 5 min, remove the liquid supernatant. Add 1 mL supernatant to the SPE cartridge, wash with 500 μ L water, wash with 500 μ L MeCN, elute with two 500 μ L portions of acidified MeOH. Evaporate the eluate to dryness under a stream of nitrogen at 60°, reconstitute the residue in 50 μ L mobile phase, mix for 15 s, inject a 40 μ L aliquot. (Acidified MeOH was 50 mL MeOH + 300 μ L 96% acetic acid.)

HPLC VARIABLES

Guard column: RCSS silica guard-pack (Waters)

Column: 250 \times 4.6 Chiralcel OD-H

Mobile phase: n-Hexane:EtOH:diethylamine 50:50:1

Flow rate: 1

Injection volume: 40

Detector: F ex 225 em 290 cut-off filter

CHROMATOGRAM

Internal standard: (R,S)-tertatolol

OTHER SUBSTANCES

Extracted: atenolol, pindolol

KEY WORDS

SPE; chiral; tertatolol is IS

REFERENCE

Hödl,K.M.; de Boer,D.; Zuidema,J.; Maes,R.A.A. Evaluation of the Salivette as sampling device for monitoring β -adrenoceptor blocking drugs in saliva, *J.Chromatogr.B*, **1995**, 663, 103–110.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 300 \times 3.9 5 μ m Nova-Pak C18

Mobile phase: MeOH:buffer 50:50 (Buffer was pH 4.0 phosphate buffer (ionic strength = 0.1) containing 4 mM N,N-dimethyloctylamine, pH readjusted to 4.00 with 85% phosphoric acid.)

Column temperature: 30

Flow rate: 1

Injection volume: 100

Detector: UV 220

CHROMATOGRAM

Retention time: k' 3.3

OTHER SUBSTANCES

Also analyzed: alprenolol, betaxolol, bopindolol, propranolol

REFERENCE

Hamoir,T.; Verlinden,Y.; Massart,D.L. Reversed-phase liquid chromatography of β -adrenergic blocking drugs in the presence of a tailing suppressor, *J.Chromatogr.Sci.*, **1994**, 32, 14–20.

SAMPLE

Matrix: solutions

Sample preparation: Inject a 20 μ L aliquot of a 1 mg/mL solution.

HPLC VARIABLES

Column: 250 \times 4.6 10 μ m Chiralcel OD

Mobile phase: Hexane:isopropanol:diethylamine 40:60:0.1

Flow rate: 0.5

Injection volume: 20

Detector: UV 290

CHROMATOGRAM

Retention time: k' 0.32, 1.69 (enantiomers)

KEY WORDS

chiral

REFERENCE

Ekelund,J.; van Arkens,A.; Bronnum-Hansen,K.; Fich,K.; Olsen,L.; Petersen,P.V. Chiral separations of β -blocking drug substances using chiral stationary phases, *J.Chromatogr.A*, **1995**, 708, 253–261.

Testolactone


Molecular formula: C₁₉H₂₄O₃

Molecular weight: 300.40

CAS Registry No.: 968-93-4

Merck Index: 9321

Lednicer No.: 1 160



SAMPLE

Matrix: blood, urine

Sample preparation: 500 μ L Plasma or urine + 25 μ L 10 μ g/mL testosterone in MeOH + 4 mL dichloromethane, shake for 15 min, centrifuge at 2000 g for 5 min. Remove the lower organic layer and evaporate it to dryness under a stream of air at 40°, reconstitute the residue in 150–300 (plasma) or 300–500 (urine) μ L mobile phase, inject a 100 μ L aliquot.

HPLC VARIABLES

Guard column: 70 \times 3 30–38 μ m Co:Pell ODS

Column: 250 \times 4.6 5 μ m Zorbax C8

Mobile phase: MeCN:0.1% acetic acid 40:60

Flow rate: 1.2

Injection volume: 100

Detector: UV 242

CHROMATOGRAM

Retention time: 9.2

Internal standard: testosterone (13.6)

Limit of detection: 20 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Pascucci, V.L.; Yeager, R.L.; Sherins, R.J.; Clark, R.V.; Gallelli, J.F.; Chatterji, D.C. Quantitation of testolactone and 4,5-dihydrotestolactone in plasma and urine using high-performance liquid chromatography, *J.Chromatogr.*, **1983**, 277, 79–85.

SAMPLE

Matrix: formulations

Sample preparation: Tablets. Weigh out finely-powdered tablets equivalent to about 50 mg testolactone, add 30 mL MeOH, shake vigorously for 30 min, make up to 50 mL with water, mix well, filter. Remove a 20 mL aliquot of the filtrate and mix it with 20 mL 1 mg/mL propylparaben in MeOH, make up to 100 mL with MeOH, mix, inject a 10 μ L aliquot. Suspensions. Measure out an amount of suspension equivalent to about 100 mg of testolactone, add 60 mL MeOH, make up to 100 mL with water, mix well, filter (if necessary). Remove a 20 mL aliquot and mix it with 20 mL 1 mg/mL propylparaben in MeOH, make up to 100 mL with MeOH, mix, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 10 μ m RP-8 (Brownlee)

Mobile phase: MeOH:water 55:45

Flow rate: 1.75

Injection volume: 10

Detector: UV 240

CHROMATOGRAM

Retention time: k' 1.89

Internal standard: propylparaben (k' 4.78)

Limit of detection: 0.1% (of testolactone)

OTHER SUBSTANCES

Simultaneous: impurities, degradation products, benzyl alcohol, 11-deoxycortisol, progesterone, testosterone

KEY WORDS

tablets; suspensions

REFERENCE

Carignan, G.; Lodge, B.A.; Skakum, W. Analysis of testolactone and its formulations by high-performance liquid chromatography, *J.Chromatogr.*, **1981**, 206, 174–176.

SAMPLE

Matrix: formulations

Sample preparation: Oils. 1 mL Sample + 25 mL MeOH:water 90:10, shake vigorously for 5 min, centrifuge, inject a 10 μ L aliquot of the supernatant. Tablets. Grind a tablet to a fine powder, add 25 mL MeOH, sonicate for 5–10 min, filter (0.45 μ m), discard first 5 mL of filtrate, inject a 10 μ L aliquot of the remaining filtrate. Suspensions (aqueous). Make up 5 mL to 50 mL with MeOH, filter (0.45 μ m), discard first 5 mL of filtrate, inject a 10 μ L aliquot of the remaining filtrate.

HPLC VARIABLES**Column:** 250 × 4.6 5 µm Zorbax ODS**Mobile phase:** MeOH:water 75:25**Flow rate:** 1.5**Injection volume:** 10**Detector:** UV 240**CHROMATOGRAM****Retention time:** 2.9**Limit of detection:** 5 µg/mL**OTHER SUBSTANCES**

Simultaneous: fluoxymesterone, norethindrone, oxandrolone (UV 210), boldenone, ethisterone, methandrostenolone, nandrolone, norgestrel, testosterone, dehydroepiandrosterone (UV 210), mibolerone, methyltestosterone, methandriol (UV 210), norethindrone acetate, calusterone, mesterolone (UV 210), norethandrolone, trenbolone acetate, benzyl benzoate, nandrolone acetate, testosterone acetate, stanozolol, oxymetholone, nandrolone propionate, methenolone acetate, testosterone propionate

Interfering: aspirin, caffeine, formebolone, benzyl alcohol, cortisone

KEY WORDS

oils; tablets; suspensions


REFERENCE

Walters,M.J.; Ayers,R.J.; Brown,D.J. Analysis of illegally distributed anabolic steroid products by liquid chromatography with identity confirmation by mass spectrometry or infrared spectrophotometry, *J.Assoc.Off.Anal.Chem.*, **1990**, 73, 904–926.

Testosterone

Molecular formula: C₁₉H₂₈O₂**Molecular weight:** 288.43

CAS Registry No.: 58-22-0, 53608-96-1 (17-chloral hemiacetal), 58-20-8 (17β-cypionate), 315-37-7 (enantate), 668-56-4 (nicotinate), 5704-03-0 (phenylacetate), 57-85-2 (propionate), 5874-98-6 (ketolaurate)

Merck Index: 9322**Lednicer No.:** 1 172**SAMPLE****Matrix:** blood

Sample preparation: Condition a Bond-Elut C18 SPE cartridge with MeOH and water. Add 500 µL plasma to the SPE cartridge, wash with 2 mL water, wash with 2 mL MeOH:water 20:80, elute with two 500 µL aliquots of MeOH. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute the residue in 50 µL MeOH:water 20:80, inject a 20 µL aliquot.

HPLC VARIABLES**Column:** 100 × 1 5 µm Hypersil ODS**Mobile phase:** MeCN:MeOH:water 25:25:50**Flow rate:** 0.1**Injection volume:** 20**Detector:** UV (wavelength not given)**CHROMATOGRAM****Retention time:** 6**Limit of detection:** 2 ng/mL

OTHER SUBSTANCES

Extracted: androstenedione, 20 α -hydroxy-4-pregnen-3-one, 17 α -hydroxyprogesterone, norethindrone, progesterone

KEY WORDS

microbore; rat; plasma; SPE

REFERENCE

Taylor,R.B.; Kendle,K.E.; Reid,R.G.; Hung,C.T. Chromatography of progesterone and its major metabolites in rat plasma using microbore high-performance liquid chromatography columns with conventional injection and detection systems, *J.Chromatogr.*, **1987**, 385, 383–392.

SAMPLE

Matrix: blood

Sample preparation: Extract 1 mL serum twice with 5 volumes ether by vortexing for 2 min, evaporate extracts to dryness under a stream of nitrogen at 35°, reconstitute in 100 μ L MeOH.

HPLC VARIABLES

Column: 240 \times 4.5 Bio-Rad ODS-5S

Mobile phase: Gradient. MeOH:MeCN:water at 20:60:20 for 3 min then to 5:85:10 over 26 min

Flow rate: 1

Injection volume: 50

Detector: UV 230

OTHER SUBSTANCES

Simultaneous: estradiol, androstenedione, progesterone

KEY WORDS

serum

REFERENCE

Yu,F.H.; Yun,Y.W.; Yuen,B.H.; Moon,Y.S. Effects of hydroxyflutamide on rats treated with a superovulatory dose of pregnant mare serum gonadotropin, *Can.J.Physiol.Pharmacol.*, **1991**, 69, 185–190.

SAMPLE

Matrix: blood

Sample preparation: Condition a 100 mg Bond-Elut C18 SPE cartridge with 2 mL 95% EtOH and 2 mL MeCN:water 15:85. 200 μ L Plasma or whole blood + 50 μ L MeOH + 3 mL MeCN:water 15:85, vortex for 30 s, add to the SPE cartridge, wash with 9 mL MeCN:water 30:70, dry, elute with 200 μ L 95% EtOH, inject a 10 μ L aliquot of the eluate.

HPLC VARIABLES

Column: 50 \times 4.6 5 μ m Supelcosil LC-8DB

Mobile phase: MeOH:buffer 72.5:27.5 (Buffer was 25 mM K₂HPO₄ adjusted to pH 3 with 670 mM phosphoric acid.)

Flow rate: 1

Injection volume: 10

Detector: UV 210

CHROMATOGRAM

Retention time: 3.60 (testosterone propionate)

Internal standard: testosterone propionate

OTHER SUBSTANCES

Extracted: clotrimazole, doxepin, itraconazole

Noninterfering: acetaminophen, N-acetylprocainamide, amitriptyline, aspirin, barbituric acid, brompheniramine, caffeine, carbamazepine, chloramphenicol, chlorpheniramine, clonazepam, desipramine, desmethyldoxepin, digitoxin, digoxin, disopyramide, ethosuximide, felbamate, gentamicin, ibuprofen, imipramine, lidocaine, maprotiline, mephenytoin, mephobarbital, metharbital, methsuximide, methylsuccinimide, nortriptyline, paramethadione, phenacemide, phen-

obarbital, phenisuximide, phenylpropanolamine, phenytoin, primidone, procainamide, protriptyline, quinidine, theophylline, tobramycin, trimethadione, valproic acid, vancomycin

KEY WORDS

testosterone propionate is IS; plasma; SPE; whole blood

REFERENCE

Rifai,N.; Sakamoto,M.; Law,T.; Platt,O.; Mikati,M.; Armsby,C.C.; Brugnara,C. HPLC measurement, blood distribution, and pharmacokinetics of oral clotrimazole, potentially useful antisickling agent, *Clin.Chem.*, **1995**, *41*, 387-391.

SAMPLE

Matrix: blood, urine

Sample preparation: 500 μ L Plasma or urine + 4 mL dichloromethane, shake for 15 min, centrifuge at 2000 g for 5 min. Remove the lower organic layer and evaporate it to dryness under a stream of air at 40°, reconstitute the residue in 150-300 (plasma) or 300-500 (urine) μ L mobile phase, inject a 100 μ L aliquot.

HPLC VARIABLES

Guard column: 70 \times 3 30-38 μ m Co:Pell ODS

Column: 250 \times 4.6 5 μ m Zorbax C8

Mobile phase: MeCN:0.1% acetic acid 40:60

Flow rate: 1.2

Injection volume: 100

Detector: UV 242

CHROMATOGRAM

Retention time: 13.6

Internal standard: testosterone

OTHER SUBSTANCES

Extracted: testolactone

KEY WORDS

plasma; testosterone is IS

REFERENCE

Pascucci,V.L.; Yeager,R.L.; Sherins,R.J.; Clark,R.V.; Gallelli,J.F.; Chatterji,D.C. Quantitation of testolactone and 4,5-dihydrotestolactone in plasma and urine using high-performance liquid chromatography, *J.Chromatogr.*, **1983**, *277*, 79-85.

SAMPLE

Matrix: blood, urine

Sample preparation: Condition a Bond-Elut C18 SPE cartridge with two 500 μ L portions of MeOH and two 500 μ L portions of water. 1 mL Plasma or urine + 1 mL water, add to the SPE cartridge, let stand for 2 min, wash with two 1 mL portions of water, wash with two 1 mL portions of MeOH:water 10:90, elute with 1 mL MeOH. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute with 50 μ L dry benzene (Caution! Benzene is a carcinogen!), add 5 mg potassium carbonate, add 50 μ L 200 mM salicylic acid chloride in dry benzene, add 50 μ L 250 mM 18-crown-6 in dry benzene, shake, heat at 70° for 1 h, cool, centrifuge, evaporate to dryness under a stream of nitrogen, reconstitute with 100 μ L mobile phase, inject a 20 μ L aliquot on to column A and elute to waste with mobile phase A, after 8 min elute the contents of column A on to column B with mobile phase B, elute with mobile phase B, monitor the effluent from column B. (Prepare salicylic acid chloride by stirring 27.5 g freshly distilled thionyl chloride in 30 mL dry benzene at 0° (Caution! Benzene is a carcinogen!), protect the reaction with a calcium chloride drying tube and a nitrogen atmosphere, add 25 g sodium salicylate, stir at 0° for 1 h, remove solvent by vacuum distillation, take up the residue in 50 mL dry petroleum ether, stir for 15 min, centrifuge. Remove the petroleum ether layer and evaporate it to give salicylic acid chloride.)

HPLC VARIABLES

Column: A 5 \times 4 35-40 μ m RP8 Perisorb (Merck); B 100 \times 4.6 Spheri 5 RP8

Mobile phase: A MeOH:water 30:70; B MeOH:water 70:30 containing 2 g/L lithium perchlorate trihydrate and 2 mL/L glacial acetic acid

Flow rate: A 0.8; B 1

Injection volume: 20

Detector: E, LKB (Bromma) 2143, glassy carbon electrode +1.0 V, palladium reference electrode

CHROMATOGRAM

Retention time: 10

Limit of quantitation: 12.5 ng/mL

OTHER SUBSTANCES

Extracted: androsterone

KEY WORDS

derivatization; plasma; SPE; column-switching

REFERENCE

Wintersteiger,R.; Sepulveda,M.J. Electrochemical detection of anabolics in human plasma and urine, *Anal.Chim.Acta*, **1993**, 273, 383-390.

SAMPLE

Matrix: formulations

Sample preparation: 1 mL Injection + 30 mL MeOH:water 90:10, shake for 15 min, centrifuge, remove the MeOH layer, repeat the extraction three times. Combine the extracts, make up to 200 mL with MeOH:water 90:10, cool to -8° for 1 h, filter an aliquot immediately. Remove an aliquot equivalent to 5 mg testosterone, add 5 mL 2 mg/mL 1,2,4,5-tetrachlorobenzene in MeOH, make up to 50 mL with MeOH, inject a 10 µL aliquot.

HPLC VARIABLES

Column: 250 × 3.2 RP8 Express Series (Altex)

Mobile phase: MeOH:THF:water 57:11:32

Flow rate: 1

Injection volume: 10

Detector: UV 240

CHROMATOGRAM

Retention time: k' 0.60 (testosterone), k' 1.80 (testosterone acetate), k' 2.68 (testosterone propionate), k' 5.00 (testosterone benzoate), k' 6.20 (testosterone phenylpropionate), k' 12.81 (testosterone enanthate), k' 13.67 (testosterone cypionate)

Internal standard: 1,2,4,5-tetrachlorobenzene (k' 4.20)

OTHER SUBSTANCES

Simultaneous: benzyl alcohol, sesamin (from sesame oil), sesamol (from sesame oil)

KEY WORDS

oils; injections

REFERENCE

Carignan,G.; Lodge,B.A.; Skakum,W. High-performance liquid chromatographic analysis of testosterone esters in oily solution, *J.Pharm.Sci.*, **1980**, 69, 1214-1217.

SAMPLE

Matrix: formulations

Sample preparation: Oils. 1 mL Sample + 25 mL MeOH:water 90:10, shake vigorously for 5 min, centrifuge, inject a 10 µL aliquot of the supernatant. Tablets. Grind a tablet to a fine powder, add 25 mL MeOH, sonicate for 5-10 min, filter (0.45 µm), discard first 5 mL of filtrate, inject a 10 µL aliquot of the remaining filtrate. Suspensions (aqueous). Make up 5 mL to 50 mL with MeOH, filter (0.45 µm), discard first 5 mL of filtrate, inject a 10 µL aliquot of the remaining filtrate.

HPLC VARIABLES**Column:** 250 × 4.6 5 µm Zorbax ODS**Mobile phase:** MeOH:water 75:25**Flow rate:** 1.5**Injection volume:** 10**Detector:** UV 240

CHROMATOGRAM**Retention time:** 6.3 (testosterone), 10.8 (testosterone acetate), 25.6 (testosterone propionate)**Limit of detection:** 5 µg/mL

OTHER SUBSTANCES**Simultaneous:** dehydroepiandrosterone (UV 210), mibolerone, methyltestosterone, methandriol (UV 210), norethindrone acetate, calusterone, mesterolone (UV 210), norethandrolone, trenbolone acetate, benzyl benzoate, nandrolone acetate, stanozolol, oxymetholone, nandrolone propionate, methenolone acetate, aspirin, caffeine, formebolone, benzyl alcohol, testolactone, cortisone, fluoxymesterone, norethindrone, oxandrolone (UV 210), boldenone, ethisterone, methandrostenolone**Interfering:** nandrolone, norgestrel

KEY WORDSoils; tablets; suspensions

REFERENCE

Walters,M.J.; Ayers,R.J.; Brown,D.J. Analysis of illegally distributed anabolic steroid products by liquid chromatography with identity confirmation by mass spectrometry or infrared spectrophotometry, *J.Assoc.Off.Anal.Chem.*, **1990**, 73, 904–926.

SAMPLE**Matrix:** formulations**Sample preparation:** Crush tablets, weigh out amount equivalent to 10 mg steroid, dissolve in 10 mL MeOH, sonicate for 15 min, filter. 1 mL Filtrate + 5 mL MeOH + 4 mL water, inject a 25 µL aliquot.

HPLC VARIABLES**Column:** 250 × 4.6 5 µm Zorbax ODS**Mobile phase:** Gradient. MeOH:water from 70:30 to 100:0 over 15 min, maintain at 100:0 for 15 min.**Flow rate:** 1**Injection volume:** 25**Detector:** UV 240

CHROMATOGRAM**Retention time:** 9.2 (testosterone), 16.2 (testosterone acetate), 24.2 (testosterone cypionate), 23.9 (testosterone enanthate), 20.0 (testosterone isobutyrate), 17.3 (testosterone propionate), 29.3 (testosterone undecanoate)

OTHER SUBSTANCES**Simultaneous:** boldenone, boldenone acetate, boldenone undecylenate, clostebol acetate, danazol (UV 280), fluoxymesterone, methandriol, methandriol-3-acetate, methandriol dipropionate, methandrostenolone, methyltestosterone, nandrolone, nandrolone decanoate, nandrolone phenylpropionate, nandrolone propionate, stanolone, stanozolol**Noninterfering:** oxandrolone, oxymetholone, testosterone decanoate, testosterone isocaproate

KEY WORDStablets

REFERENCE

Lurie,I.S.rling,A.R.; Meyers,R.P. The determination of anabolic steroids by MECC, gradient HPLC, and capillary GC, *J.Forensic Sci.*, **1994**, 39, 74–85.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 1 mL Microsomal incubation + 1 mL dichloromethane, extract, centrifuge, remove organic layer and evaporate it under vacuum, dissolve residue in 30 μ L MeCN:water 50:50, centrifuge for 3 min, inject supernatant. After each run wash column with MeCN for 1 min, re-equilibrate for 1 min.

HPLC VARIABLES

Column: 50 \times 4.6 3 μ m Spherisorb ODS-2

Mobile phase: MeCN:water 50:50

Column temperature: 60

Flow rate: 2

Injection volume: 30

Detector: UV 200

CHROMATOGRAM

Retention time: 0.9

Limit of detection: <0.1 nmol/mL

OTHER SUBSTANCES

Simultaneous: estradiol

Interfering: estrone, androstenedione

KEY WORDS

human; placenta

REFERENCE

Taniguchi,H.; Feldmann,H.R.; Kaufmann,M.; Pyerin,W. Fast liquid chromatographic assay of androgen aromatase activity, *Anal.Biochem.*, **1989**, 181, 167–171.

SAMPLE

Matrix: microsomal incubations

Sample preparation: Add 11 β -hydroxytestosterone to microsomal incubation, extract with 5 volumes dichloromethane, evaporate the organic layer to dryness under a stream of nitrogen, reconstitute the residue in 250 μ L mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 Nucleosil 5 C18

Mobile phase: MeOH:water 50:50

Flow rate: 0.8

Detector: UV 254

CHROMATOGRAM

Internal standard: 11 β -hydroxytestosterone

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

rat; liver

REFERENCE

Anderson,C.D.; Wang,J.; Kumar,G.N.; McMillan,J.M.; Walle,U.K.; Walle,T. Dexamethasone induction of taxol metabolism in the rat, *Drug Metab.Dispos.*, **1995**, 23, 1286–1290.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 4 mL Microsomal incubation + 100 μ L 10 μ g/mL 11 β -hydroxytestosterone + 6 mL dichloromethane, rotate for 30 min, centrifuge at 2000 g for 5 min. Remove the organic

layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 300 μL MeOH:water 25:75, vortex, inject a 100 μL aliquot.

HPLC VARIABLES

Guard column: 5 μm Hypersil C18

Column: 150 \times 3.9 4 μm Novapack C18

Mobile phase: Gradient. A was MeOH:water 25:75. B was MeCN:MeOH:water 1.5:63.5:35. A:B 75:25 for 10 min, to 49.5:50.5 over 8.5 min, to 30:70 over 6.5 min, to 0:100 over 10 min, maintain at 0:100 for 15 min, return to initial conditions over 10 min, re-equilibrate for 10 min.

Flow rate: 1

Injection volume: 100

Detector: UV 254

CHROMATOGRAM

Retention time: 29.5

Internal standard: 11 β -hydroxytestosterone (20.8)

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

rat; liver

REFERENCE

Ekins,S.; Murray,G.I.; Burke,M.D.; Williams,J.A.; Marchant,N.C.; Hawksorth,G.M. Quantitative differences in phase I and II metabolism between rat precision-cut liver slices and isolated hepatocytes, *Drug Metab Dispos.*, **1995**, 23, 1274–1279.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 1 mL Microsomal incubation + 2.5 mL cold ethyl acetate (4°), mix, centrifuge. Remove the supernatant and evaporate it to dryness under a stream of nitrogen, reconstitute the residue, inject an aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μm C18 (Supelco)

Mobile phase: MeCN:MeOH:water 1.3:37.8:60.9

Flow rate: 1.3

Detector: UV 240

OTHER SUBSTANCES

Extracted: metabolites, phenacetin

KEY WORDS

rat; liver

REFERENCE

Lin,J.H.; Chiba,M.; Chen,I.-W.; Vastag,K.J.; Nishime,J.A.; Dorsey,B.D.; Michelson,S.R.; McDaniel,S.L. Time- and dose-dependent pharmacokinetics of L-754,394, an HIV protease inhibitor, in rats, dogs and monkeys, *J.Pharmacol.Exp.Ther.*, **1995**, 274, 264–269.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 200 μL Microsomal incubation + 200 μL ice-cold MeOH containing 3 nmoles corticosterone, centrifuge at 1500 g for 5 min, inject a 200 μL aliquot of the supernatant.

HPLC VARIABLES

Guard column: 6 \times 4 10 μm μ Bondapak C18 Guard-Pak

Column: 250 \times 4.6 5 μm Ultrasphere IP

Mobile phase: MeOH:THF:water 35:10:55 adjusted to pH 4.0 with glacial acetic acid

Flow rate: 1
Injection volume: 200
Detector: UV 245

CHROMATOGRAM

Retention time: 26.1
Internal standard: corticosterone (15.1)

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

rat; liver

REFERENCE

Sanwald,P.; Blankson,E.A.; Duléry,B.D.; Schoun,J.; Huebert,N.D.; Dow,J. Isocratic high-performance liquid chromatographic method for the separation of testosterone metabolites, *J.Chromatogr.B*, **1995**, 672, 207–215.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Spherisorb S5-ODS2 (A), 125 × 4.5 µm LiChrospher ODS-3 (B), 250 × 4.6 10 µm Partisil 10 ODS-3 (C)
Mobile phase: MeCN:water 65:35 (A), MeCN:MeOH:water 40:30:30 (B and C)
Flow rate: 1(A), 1.5 (B), 2 (C)
Injection volume: 20 (A), 100 (B and C)
Detector: UV 270

CHROMATOGRAM

Internal standard: testosterone propionate

OTHER SUBSTANCES

Simultaneous: danazol

REFERENCE

Galia,E.; Nicolaidis,E.; Hörter,D.; Löbenberg,R.; Reppas,C.; Dressman,J.B. Evaluation of various dissolution media for predicting in vivo performance of class I and II drugs, *Pharm.Res.*, **1998**, 15, 698–705.

SAMPLE

Matrix: solutions
Sample preparation: Prepare solutions in MeCN:water 50:50, inject a 30 µL aliquot.

HPLC VARIABLES

Guard column: 20 × 3.9 4 µm Sentry guard column (Waters)
Column: 150 × 3.9 4 µm mean pore diameter 60 Å Nova-Pak Phenyl
Mobile phase: THF:water 21.5:78.5
Flow rate: 1
Injection volume: 30
Detector: UV 250

CHROMATOGRAM

Retention time: 22.5

OTHER SUBSTANCES

Simultaneous: spironolactone

REFERENCE

Kaukonen,A.M.; Vuorela,P.; Vuorela,H.; Mannermaa,J.-P. High-performance liquid chromatography methods for the separation and quantitation of spironolactone and its degradation products in aqueous formulations and of its metabolites in rat serum, *J.Chromatogr.A*, **1998**, 797, 271–281.

SAMPLE

Matrix: solutions

Sample preparation: Prepare solutions in MeCN, dilute to an appropriate concentration with mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 120 \times 4.6 5 μ m octadecyl Bakerbond

Mobile phase: MeCN:water 30:70 containing 16 mM β -cyclodextrin

Column temperature: 5

Flow rate: 1

Injection volume: 20

Detector: UV 240

CHROMATOGRAM

Retention time: 0.5

OTHER SUBSTANCES

Simultaneous: hydrocortisone, prednisone, cortisone, 17 α -methyltestosterone, 17 α -hydroxyprogesterone

REFERENCE

Zarzycki,P.K.; Wierzbowska,M.; Lamparczyk,H. The influence of temperature on the high performance liquid chromatographic separation of steroids using mobile phases modified with β -cyclodextrin, *J.Pharm.Biomed.Anal.*, **1996**, 14, 1305–1311.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a solution in mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 300 \times 4 10 μ m μ Porasil

Mobile phase: Butyl chloride:water-saturated butyl chloride:THF:glacial acetic acid 55:55:3:2

Flow rate: 2-3

Detector: UV 254

CHROMATOGRAM

Retention time: 8

Internal standard: testosterone acetate

OTHER SUBSTANCES

Simultaneous: iodochlorhydroxyquin

KEY WORDS

normal phase; testosterone acetate is IS

REFERENCE

Kubiak,E.J.; Munson,J.W. Analysis of iodochlorhydroxyquin in cream formulations and bulk drugs by high-performance liquid chromatography, *J.Pharm.Sci.*, **1982**, 71, 872–875.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in MeOH:water 1:1 at a concentration of 50 μ g/mL, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 300 × 3.9 10 μm μBondapak C18

Mobile phase: MeOH:acetic acid:triethylamine:water 80:1.5:0.5:18

Flow rate: 1.5

Injection volume: 10

Detector: UV

CHROMATOGRAM

Retention time: k' 0.52 (testosterone), k' 8.07 (testosterone cypionate), k' 7.12 (testosterone enanthate), k' 1.95 (testosterone propionate)

REFERENCE

Roos,R.W.; Lau-Cam,C.A. General reversed-phase high-performance liquid chromatographic method for the separation of drugs using triethylamine as a competing base, *J.Chromatogr.*, **1986**, 370, 403–418.

SAMPLE

Matrix: solutions

Sample preparation: Inject an aliquot of a solution in MeOH.

HPLC VARIABLES

Column: Radial-PAK μBondapak C18

Mobile phase: MeCN:water 50:50

Flow rate: 2

Injection volume: 100

Detector: UV 254 or 214

CHROMATOGRAM

Retention time: 5.5

OTHER SUBSTANCES

Simultaneous: estrone, estriol, progesterone

Interfering: estradiol

REFERENCE

Erkoc,F.U.; Özsar,S.; Güven,B.; Kalkandelen,G.; Ugrar,E. High-performance liquid chromatographic analysis of steroid hormones, *J.Chromatogr.Sci.*, **1989**, 27, 86–90.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Guard column: 4 × 4 5 μm LiChrospher 100 RP-18

Column: 250 × 4 5 μm LiChrospher CH-18

Mobile phase: MeOH:10 mM pH 6.0 phosphate buffer 75:25 containing 2.5 mM cethexonium bromide (Rinse with 100 mL MeOH:EtOH:water 50:25:25 at the end of the day.)

Flow rate: 1

Injection volume: 10

Detector: UV 220

CHROMATOGRAM

Retention time: k' 1.3

OTHER SUBSTANCES

Extracted: glucuronides

REFERENCE

Liu,H.-F.; Leroy,P.; Nicolas,A.; Magdalou,J.; Siest,G. Evaluation of a versatile reversed-phase high-performance liquid chromatographic system using cethexonium bromide as ion-pairing reagent for the analysis of glucuronic acid conjugates, *J.Chromatogr.*, **1989**, 493, 137–147.

SAMPLE**Matrix:** solutions**Sample preparation:** Dissolve in MeOH at a concentration of 100 µg/mL, inject a 5 µL aliquot.

HPLC VARIABLES**Guard column:** 70 × 2.1 CO:PELL ODS**Column:** 300 × 3.9 Bondex C18 (Phenomenex)**Mobile phase:** MeOH:water 85:15**Flow rate:** 1**Injection volume:** 5**Detector:** UV 254

CHROMATOGRAM**Retention time:** 4.5 (testosterone), 6 (testosterone acetate), 8 (testosterone propionate), 13 (testosterone benzoate), 19 (testosterone enanthate), 21 (testosterone cypionate)

OTHER SUBSTANCES**Simultaneous:** boldenone acetate, nandrolone propionate, boldenone benzoate, nandrolone phenylpropionate

REFERENCENoggle, F.T., Jr.; Clark, C.R.; DeRuiter, J. Liquid chromatographic and mass spectral analysis of the anabolic 17-hydroxy steroid esters, *J.Chromatogr.Sci.*, **1990**, 28, 263–268.

SAMPLE**Matrix:** solutions**Sample preparation:** Inject an aliquot of a 100 µg/mL solution in MeOH.

HPLC VARIABLES**Guard column:** 70 × 2.1 Whatman CO:PELL ODS**Column:** 300 × 3.9 Bondex C18**Mobile phase:** MeOH:water 70:30**Flow rate:** 1**Injection volume:** 5**Detector:** UV 254

CHROMATOGRAM**Retention time:** 8

OTHER SUBSTANCES**Simultaneous:** methyltestosterone, nandrolone, methandrostenolone, boldenone, danazol, fluoxymesterone

REFERENCENoggle, F.T., Jr.; Clark, C.R.; DeRuiter, J. Liquid chromatographic and spectral analysis of the 17-hydroxy anabolic steroids, *J.Chromatogr.Sci.*, **1990**, 28, 162–166.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 150 × 6.5 µm Shim-pack CLC-ODS**Mobile phase:** MeOH:THF:water 26:18:56**Column temperature:** 48**Flow rate:** 1**Injection volume:** 20**Detector:** UV 240

CHROMATOGRAM**Retention time:** 11.5

OTHER SUBSTANCES

Simultaneous: cortisone, estriol, cortisol, corticosterone, 11-deoxycortisol, androstenedione, prednisone acetate, 11-deoxycorticosterone, 17 α -hydroxyprogesterone, dexamethasone acetate, estradiol, estrone, progesterone

REFERENCE

Wei,J.Q.; Wei,J.L.; Zhou,X.T. Optimization of an isocratic reversed phase liquid chromatographic system for the separation of fourteen steroids using factorial design and computer simulation, *Biomed.Chromatogr.*, 1990, 4, 34–38.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 0.5 mg/mL solution in MeOH, inject a 5 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 Zorbax RX

Mobile phase: Gradient. A was 150 mM phosphoric acid and 50 mM triethylamine. B was MeCN: water 80:20 containing 150 mM phosphoric acid and 50 mM triethylamine. A:B 100:0 for 2.2 min then to 0:100 over 30 min.

Column temperature: 30

Flow rate: 2

Injection volume: 5

Detector: UV 210

CHROMATOGRAM

Retention time: 22.9 (testosterone), 30.7 (testosterone propionate)

OTHER SUBSTANCES

Simultaneous: acetaminophen, aprobarbital, butabarbital, chlordiazepoxide, chloroxylenol, chlorpromazine, clenbuterol, cortisone, danazol, doxapram, fluoxymesterone, mefenamic acid, methyltestosterone, nicotine, oxazepam, phentermine, phenylpropanolamine, progesterone, sulfamethazine, sulfanilamide, tranlycypromine, tripeleennamine

Interfering: (with testosterone) diflunisal, estrone

KEY WORDS

details for purification of triethylamine in paper

REFERENCE

Hill,D.W.; Kind,A.J. The effects of type B silica and triethylamine on the retention of drugs in silica based reverse phase high performance chromatography, *J.Liq.Chromatogr.*, 1993, 16, 3941–3964.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepan, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlordiazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine,

chlorthalidone, chlorpromazine, chlorpropamide, chlorthalidone, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarol, danazol, danthron, dapson, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fenclonamine, fenpropofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarboxystyrene, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephentermine, mephentermine, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methyprylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidone, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfathiazole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, transylpromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 25 µg/mL solution in mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.6 Partisil 10 ODS-1

Mobile phase: MeOH:water 55:45

Column temperature: 40

Flow rate: 1.5

Detector: UV 240

CHROMATOGRAM

Retention time: k' 4.146

OTHER SUBSTANCES

Also analyzed: androsterone (UV 210), cortexolone (UV 240), cortisone (UV 240), estradiol (UV 280), estrone (UV 280), ethinyl estradiol (UV 280), ethisterone (UV 240), hydrocortisone (UV 240), hydroxyprogesterone (UV 240), lynestrenol (UV 210), medroxyprogesterone acetate (UV 240), medroxyprogesterone (UV 240), methandienone (UV 240), methylandrostenediol (UV 210), methylprednisolone acetate (UV 240), methylprednisolone (UV 240), methyltestosterone (UV 240), nandrolone (UV 240), norethisterone (UV 240), prednisolone acetate (UV 240), prednisolone (UV 240), prednisone (UV 240), pregnenolone (UV 210), progesterone (UV 240)

REFERENCE

Sadlej-Sosnowska, N. Structure retention relationship for steroid hormones. Functional groups as structural descriptors, *J.Liq.Chromatogr.*, **1994**, *17*, 2319–2330.

SAMPLE

Matrix: solutions

Sample preparation: Inject a 20 μ L aliquot of a solution in MeOH:water 50:50.

HPLC VARIABLES

Column: 250 \times 4 7 μ m LichroCART RP-8 (Merck)

Mobile phase: MeCN:MeOH:water 32:37:31

Flow rate: 1

Injection volume: 20

Detector: UV 230

CHROMATOGRAM

Retention time: 16 (testosterone propionate)

OTHER SUBSTANCES

Simultaneous: fluoxymesterone, medrogestone, mestranol, norethindrone, progesterone

REFERENCE

Gau, Y.S.; Sun, S.W.; Chen, R.R.-L. Optimization of high-performance liquid chromatographic separation for progestogenic, estrogenic, and androgenic steroids using factorial design, *J.Liq.Chromatogr.*, **1995**, *18*, 2373–2382.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Nucleosil phenyl

Mobile phase: Gradient. Carbon dioxide:MeOH from 98:2 to 78:22 over 40 min

Column temperature: 50

Flow rate: 2

Detector: UV

CHROMATOGRAM

Retention time: 9.2 (testosterone), 5.6 (testosterone enanthate)

OTHER SUBSTANCES

Simultaneous: estradiol, hydrocortisone, norethisterone, hydroxyprogesterone, estriol, other steroids

KEY WORDS

SFC; 200 bar

REFERENCE

Hanson, M. Aspects of retention behaviour of steroids in packed column supercritical fluid chromatography, *Chromatographia*, **1995**, *40*, 58–68.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a solution in n-propanol:water 80:20 or DMF:water 80:20, inject an aliquot.

HPLC VARIABLES

Column: 250 \times 4 5 μ m LiChrospher 100 Diol

Mobile phase: Gradient. A was hexane. B was ethyl acetate. C was 0.1% formic acid in MeCN.

D was 0.1% formic acid in water. A:B:C:D 100:0:0:0 for 5 min, to 0:100:0:0 over 15 min, maintain

at 0:100:0:0 for 5 min, to 0:0:100:0 over 5 min, maintain at 0:0:100:0 for 5 min; to 0:0:0:100 over 25 min, maintain at 0:0:0:100 for 5 min.

Flow rate: 0.9

Detector: Evaporative light scattering (Sédex 55, Sédéré)

CHROMATOGRAM

Retention time: 17.71

OTHER SUBSTANCES

Simultaneous: acetylcholine, cholesterol, choline, cortisone, dextrose, estradiol, glycine, phenylalanine, sodium

REFERENCE

Treiber, L.R. Normal-phase high-performance liquid chromatography with relay gradient elution. I. Description of the method, *J. Chromatogr. A*, **1995**, 696, 193–199.

SAMPLE

Matrix: solutions

Sample preparation: Inject a 20 µL aliquot of a 100–500 µg/mL solution in mobile phase.

HPLC VARIABLES

Column: 100 × 4.6 5 µm Hypersil C8 MOS 100A coated with phosphatidylcholine (95% pure soybean lecithin, Epikuron, Lucas Meyer & Co.) (Coat column by recycling a 1 mM solution of phosphatidylcholine in MeOH:water 80:20 for 24 h.)

Mobile phase: MeCN:35 mM pH 7.4 sodium phosphate buffer 40:60

Flow rate: 0.5–2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: k' 3.31

OTHER SUBSTANCES

Also analyzed: amoxicillin, antipyrine, carbamazepine, chlorpheniramine, chlorpromazine, clonidine, codeine, desipramine, diphenhydramine, dipyrindamole, ephedrine, flufenamic acid, haloperidol, hydroxyzine, imipramine, indomethacin, lidocaine, megestrol acetate, metoprolol, nabumetone, nadolol, phenobarbital, phenol, promazine, propranolol, pyrilamine, quinidine, ropinirole, thioridazine, tolfenamic acid, verapamil

Noninterfering: acetaminophen, aspirin, azathioprine, caffeine, carprofen, chlorambucil, cimetidine, fenoterol, flurbiprofen, ibuprofen, ketoprofen, ranitidine, salicylic acid, sulfamethoxazole, theophylline, thioguanine, tiaprofenic acid, trimethoprim, valproic acid

KEY WORDS

comparison with capillary electrophoresis

REFERENCE

Hanna, M.; de Biasi, V.; Bond, B.; Salter, C.; Hutt, A.J.; Camilleri, P. Estimation of the partitioning characteristics of drugs: A comparison of a large and diverse drug series utilizing chromatographic and electrophoretic methodology, *Anal. Chem.*, **1998**, 70, 2092–2099.

SAMPLE

Matrix: tissue

Sample preparation: Extract 70–125 mg tissue four times with 5 mL portions of ether:chloroform 80:20. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 100 µL MeOH, inject an aliquot.

HPLC VARIABLES

Column: 80 mm long 10 µm octadecylsilane radial compression (Radial-Pak) (Waters)

Mobile phase: Gradient. A was MeOH:water 50:50. B was MeOH. A:B from 100:0 to 70:30 over 20 min, to 0:100 over 20 min

Flow rate: 2

Detector: UV 254

CHROMATOGRAM

Retention time: 28

OTHER SUBSTANCES

Extracted: androstenedione, deoxycortisol, hydrocortisone, 17-hydroxyprogesterone

Simultaneous: estriol, estradiol, pregnenolone, progesterone, testosterone enanthate, testosterone propionate

KEY WORDS

tumor

REFERENCE

Kessler, M.J. Analysis of steroids from normal and tumor tissue by HPLC, *Clin. Chim. Acta*, **1982**, *125*, 21–30.

SAMPLE

Matrix: tissue

Sample preparation: 1 g Tissue + 10 mL Chloroform:MeOH 2:1, homogenize for 1 min (Polytron setting 5), filter, rinse tube with an additional 10 mL chloroform:MeOH, filter, combine filtrates, add 4 mL water, vortex for 1 min, centrifuge at 600 g for 10 min. Remove organic layer and dry it under air at 40°. Reconstitute with 200 μ L MeOH, add to an activated Sep-Pak C18 cartridge, wash tube onto cartridge with 200 μ L MeOH, elute with 5 mL each 0, 25, 50, 75, 100% MeOH, collect 5 mL fractions, inject a 100 μ L aliquot of each fraction. (Elutes in 75% MeOH fraction.)

HPLC VARIABLES

Guard column: in line guard column

Column: 80 mm long 10 μ m μ Bondapak C18 radially compressed

Mobile phase: MeCN:water 45:55

Flow rate: 3

Injection volume: 100

Detector: UV 254

CHROMATOGRAM

Retention time: 5.05

Limit of quantitation: 1000 ng/g

OTHER SUBSTANCES

Simultaneous: methyltestosterone

KEY WORDS

fish; muscle; tilapia aurea; SPE

REFERENCE

Goudie, C.A. Extraction of a synthetic androgen from fish muscle and quantitation by high performance liquid chromatography, *Steroids*, **1984**, *44*, 241–252.

SAMPLE

Matrix: tissue

Sample preparation: Dry pack 60 \times 8 mm glass columns with 250 mg Carbowpack B (200–400 mesh) and 60 \times 4 mm glass columns with 50 mg Amberlite CG-400 I (100–200 mesh). Wash Carbowpack column with 5 mL MeOH, 15 mL dichloromethane:MeOH 70:30, and MeOH:water 85:15. Wash Amberlite column with 3 mL 0.5 M NaOH, 8 mL dichloromethane:MeOH 70:30, 1 mL water, and 3 mL 1 M HCl. Repeat this cycle 4 times. Finally pass through 20 mL 50 mM NaOH then 1 mL water. Keep column in water. (Process converts Amberlite to OH form.) Homogenize 1 g of tissue in 5 mL MeOH, sonicate 5 min, centrifuge at 6000 rpm for 10 min. Add another 5 mL MeOH to pellet and repeat. Combine supernatants, make up to 6.8 mL with MeOH, add 1.2 mL water. Pass through Carbowpack column, wash column with 2 mL MeOH:water 85:15 then 2 mL MeOH, elute column with 8 mL dichloromethane:MeOH 70:30. Pass eluate onto Amberlite column, add 1 mL MeOH to column, collect all eluates from column,

evaporate to dryness under nitrogen at 40°, take up in 200 µL MeOH:water 50:50, add 25 µL 10 µg/mL p-chlorophenol, inject 50 µL aliquot.

HPLC VARIABLES

Guard column: 20 × 4.6 5 µm Supelguard LC-18

Column: 250 × 4.6 5 µm Supelco C18

Mobile phase: Gradient. MeCN:water from 40:60 to 65:35 in 30 min

Flow rate: 1.2

Injection volume: 50

Detector: UV 242

CHROMATOGRAM

Retention time: 11

Internal standard: p-chlorophenol (7)

Limit of detection: 1 ng/g

OTHER SUBSTANCES

Simultaneous: trenbolone, progesterone

KEY WORDS

muscle; liver; chicken; ox; SPE

REFERENCE

Laganà,A.; Marino,A. General and selective isolation procedure for high-performance liquid chromatographic determination of anabolic steroids in tissues, *J.Chromatogr.*, **1991**, 588, 89–98.

SAMPLE

Matrix: urine

Sample preparation: Add 10 mL urine to a Supelclean LC-18 SPE tube at a flow rate of 2 mL/min, wash with 4 mL 25 mM sodium borate buffer, wash with 4 mL 40% MeOH, wash with 4 mL 20% acetone, elute with two 500 µL aliquots of 73% MeOH, evaporate under nitrogen at 40°, reconstitute with 1 mL mobile phase, inject a 200 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Microsorb silica

Mobile phase: Cyclohexane:ethyl acetate 40:60

Injection volume: 200

Detector: F ex 247 em 547, after post-column reaction with 30 mM Tb(NO₃)₃ in ethyl acetate using a 50 cm tightly coiled capillary tube to ensure mixing

CHROMATOGRAM

Retention time: 10 (testosterone acetate), 17 (testosterone)

Limit of detection: 130 pg/mL (testosterone acetate), 85 pg/mL (testosterone)

OTHER SUBSTANCES

Extracted: progesterone, bolasterone

Simultaneous: 17-methyltestosterone

KEY WORDS

SPE; normal phase; post-column reaction

REFERENCE

Amin,M.; Harrington,K.; von Wandruszka,R. Determination of steroids in urine by micellar HPLC with detection by sensitized terbium fluorescence, *Anal.Chem.*, **1993**, 65, 2346–2351.

SAMPLE

Matrix: urine

Sample preparation: Inject 200 µL directly.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Rainin C18

Mobile phase: MeCN:100 mM sodium dodecyl sulfate + 10 mM Tb(NO₃)₃ 20:80 (Sodium dodecyl sulfate was 99.99%.)

Column temperature: 40

Injection volume: 200

Detector: F ex 247 em 547

CHROMATOGRAM

Retention time: 8 (testosterone), 27 (testosterone acetate)

Limit of detection: 50 ng/mL (testosterone), 10 ng/mL (testosterone acetate)

OTHER SUBSTANCES

Extracted: progesterone, 17-methyltestosterone, bolasterone

KEY WORDS

micellar chromatography

REFERENCE

Amin,M.; Harrington,K.; von Wandruszka,R. Determination of steroids in urine by micellar HPLC with detection by sensitized terbium fluorescence, *Anal.Chem.*, **1993**, *65*, 2346–2351.

SAMPLE

Matrix: urine

Sample preparation: 5 mL Urine + 10 µL 100 µg/mL methyltestosterone + 1 mL 200 mM pH 7.0 sodium phosphate buffer + 50 µL β-glucuronidase (E. coli K12, Boehringer Mannheim), heat at 55° for 1 h, cool to room temperature, add 1 g sodium bicarbonate:potassium carbonate 1:2, add 1 g sodium sulfate, add 5 mL n-pentane, shake for 20 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 50 µL MeOH, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.5 µm Hypersil BDS-C18

Mobile phase: Gradient. MeCN:1 mM phosphoric acid from 25:75 to 30:70 over 10 min, to 35:65 over 6.5 min, maintain at 35:65 for 11.5 min.

Column temperature: 40

Flow rate: 1.2

Injection volume: 20

Detector: UV 240

CHROMATOGRAM

Retention time: 20

Internal standard: methyltestosterone (25)

Limit of quantitation: 20 ng/mL

OTHER SUBSTANCES

Extracted: epitestosterone

REFERENCE

Navajas,R.; Imaz,C.; Carreras,D.; García,M.; Pérez,M.; Rodríguez,C.; Rodríguez,A.F.; Cortés,R. Determination of epitestosterone and testosterone in urine by high-performance liquid chromatography, *J.Chromatogr.B*, **1995**, *673*, 159–164.

Tetracaine


Molecular formula: $C_{15}H_{24}N_2O_2$

Molecular weight: 264.37

CAS Registry No.: 94-24-6, 136-47-0 (HCl)

Merck Index: 9330

Lednicer No.: 1 110



SAMPLE

Matrix: blood

Sample preparation: Collect 10 mL blood in a bottle containing 200 μ L 10% sodium metabisulfite and 200 μ L 2 M NaOH, centrifuge to separate plasma. 2 mL Plasma + salicylic acid + propiophenone + 8 mL diethyl ether:dichloromethane 10:8, shake for 10 min, centrifuge at 3000 g for 5 min. Remove the organic layer and evaporate it to dryness at 50°, reconstitute the residue in 200 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: μ Bondapak C18

Mobile phase: MeCN:MeOH:water 20:20:60 containing 600 μ L/L sulfuric acid, 5 g/L sodium sulfate, and 200 mg/L sodium heptanesulfonate, pH 2.6

Flow rate: 2

Injection volume: 20

Detector: UV 220

CHROMATOGRAM

Retention time: 5

Internal standard: salicylic acid (6), propiophenone (10)

Limit of detection: 20 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma

REFERENCE

Mazumdar,B.; Tomlinson,A.A.; Faulder,G.C. Preliminary study to assay plasma amethocaine concentrations after topical application of a new local anaesthetic cream containing amethocaine, *Br.J.Anaesth.*, **1991**, 67, 432-436.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 100 μ L 1 M NaOH + 3 mL heptane:ethyl acetate 90:10, shake for 2 min, centrifuge at 1200 g for 10 min. Remove the organic phase and add it to 50 μ L 50 mM sulfuric acid, shake for 2 min, centrifuge at 1200 g for 5 min. Remove the aqueous phase and add it to 820 μ g sodium acetate, inject a 40 μ L aliquot. (The sodium acetate was measured out by adding 50 μ L 200 mM sodium acetate in MeOH to the tube and evaporating the MeOH.)

HPLC VARIABLES

Column: 250 \times 4 10 μ m μ Bondapak C18

Mobile phase: MeCN:10 mM NaH₂PO₄ 35:65, adjusted to pH 2.1

Column temperature: 30

Flow rate: 1

Injection volume: 40

Detector: UV 205

CHROMATOGRAM

Retention time: 8

Limit of detection: 5 ng/mL

OTHER SUBSTANCES

Extracted: pramocaine

KEY WORDS

plasma; rabbit

REFERENCE

Le Guévello,P.; Le Corre,P.; Chevanne,P.; Le Verge,R. High-performance liquid chromatographic determination of bupivacaine in plasma samples for biopharmaceutical studies and application to seven other local anaesthetics, *J.Chromatogr.*, **1993**, 622, 284–290.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Plasma + 100 μ L 20 μ g/mL propentofylline + 8 mL dichloromethane, shake for 20 min, centrifuge at 2500 rpm for 20 min. Remove 7 mL of the organic layer and evaporate it to dryness under nitrogen or at 60°. Dissolve residue in 200 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 6 Shimpack CLS-ODS (Shimadzu)

Mobile phase: MeCN:MeOH:0.5 mM phosphoric acid 23:20:57

Column temperature: 40

Flow rate: 1.5

Injection volume: 20

Detector: UV 290

CHROMATOGRAM

Internal standard: propentofylline

KEY WORDS

plasma; rat

REFERENCE

Lee,C.K.; Uchida,T.; Kitagawa,K.; Yagi,A.; Kim,N.-S.; Goto,S. Skin permeability of various drugs with different lipophilicity, *J.Pharm.Sci.*, **1994**, 83, 562–565.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol:n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 μ L mobile phase, centrifuge at 2800 g for 5 min, inject a 50 μ L aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 \times 3.9 μ m NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH_2PO_4 adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 313

CHROMATOGRAM

Retention time: 6.16

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; tolaxotone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanin; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vandesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254–262.

SAMPLE

Matrix: blood, urine

Sample preparation: Urine. Adjust pH of urine to 5 before freezing. Adjust pH of 5 mL urine to 9.5 with borax buffer, extract twice with 7 mL portions of dichloromethane. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 60°, reconstitute the residue in 1 mL MeOH, inject a 20–100 µL aliquot. Blood. Adjust pH of 4 mL plasma or whole blood to 9.5 with borax buffer, extract twice with 7 mL portions of dichloromethane. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 60°, reconstitute the residue in 0.2–2 mL MeOH, inject a 125 µL aliquot.

HPLC VARIABLES

Column: 10 µm Radial Pak C18

Mobile phase: MeCN:16.5 mM triethylamine 85:15, pH adjusted to 3 with concentrated phosphoric acid

Flow rate: 2

Injection volume: 20–125

Detector: UV 288

CHROMATOGRAM

Internal standard: tetracaine

Limit of detection: 10 ng/mL (urine), 1 ng/mL (plasma)

OTHER SUBSTANCES

Extracted: procaine

KEY WORDS

horse; plasma; whole blood; pharmacokinetics; tetracaine is IS

REFERENCE

Stevenson,A.J.; Weber,M.P.; Todi,F.; Mendonca,M.; Fenwick,J.D.; Young,L.; Kwong,E.; Chen,F.; Beaumier,P.; Timmings,S.; Woodard,W.; Kacew,S. Determination of procaine in equine plasma and urine by high-performance liquid chromatography, *J.Anal.Toxicol.*, **1992**, 16, 93-96.

SAMPLE

Matrix: blood, urine

Sample preparation: 2 mL Whole blood, plasma, or urine + 1 mL saturated sodium carbonate + 50 μ L 10 μ g/mL butacaine, add to a 3 mL Extrelut SPE cartridge, elute with 15 mL dichloromethane. Evaporate eluate to dryness under a stream of nitrogen at 40°, reconstitute in 100 μ L mobile phase, inject a 40 μ L aliquot.

HPLC VARIABLES

Guard column: 5 \times 6 μ Bondapak Guard Pak

Column: 300 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: MeCN:100 mM ammonium acetate 50:50

Flow rate: 1.5

Injection volume: 40

Detector: UV 280

CHROMATOGRAM

Retention time: 14

Internal standard: butacaine (10)

Limit of detection: 40 ng/mL

OTHER SUBSTANCES

Extracted: procaine, p-aminobenzoic acid

Also analyzed: articaine, prilocaine, o-toluidine, lidocaine, bupivacaine, etidocaine, dibucaine, caffeine, amphetamine, ephedrine, epinephrine, morphine, monoacetylmorphine, diamorphine, ethylmorphine, codeine, acetylcodeine

KEY WORDS

whole blood; plasma; SPE

REFERENCE

Rop,P.P.; Grimaldi,F.; Bresson,M.; Fornaris,M.; Viala,A. Liquid chromatographic analysis of cocaine, benzoyl-ecgonine, local anaesthetic agents and some of their metabolites in biological fluids, *J.Liq.Chromatogr.*, **1993**, 16, 2797-2811.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 312.8

CHROMATOGRAM

Retention time: 13.69

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

SAMPLE

Matrix: formulations

Sample preparation: 5 mL Injection + 20 mL 10 mg/mL salicylic acid in MeOH:water 50:50 + 5 mL 10 mg/mL propiophenone in MeOH:water 50:50, make up to 50 mL with MeOH:water 50:50, homogenize (if necessary), inject a 5 μ L aliquot.

HPLC VARIABLES

Column: 300 \times 4 10 μ m μ Bondapak C18

Mobile phase: MeCN:MeOH:water 20:20:60 containing 0.06% sulfuric acid, 0.5% sodium sulfate, and 0.02% sodium heptanesulfonate, pH 2.6

Flow rate: 2

Injection volume: 5

Detector: UV 305

CHROMATOGRAM

Retention time: 6

Internal standard: salicylic acid (4), propiophenone (8)

OTHER SUBSTANCES

Simultaneous: degradation products

KEY WORDS

injections; saline

REFERENCE

Menon, G.N.; Norris, B.J. Simultaneous determination of tetracaine and its degradation product, p-n-butylaminobenzoic acid, by high-performance liquid chromatography, *J. Pharm. Sci.*, **1981**, 70, 569-570.

SAMPLE

Matrix: formulations

Sample preparation: Weigh out 1.25 g ground pastilles, add 50 mL mobile phase, stir mechanically until dissolved, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: MeCN:100 mM pH 5 potassium phosphate buffer 25:75, containing 5.9 g/L NaCl and 30 mM tetrabutylammonium hydrogen sulfate

Flow rate: 2

Injection volume: 50

Detector: UV 294

CHROMATOGRAM**Retention time:** 5

OTHER SUBSTANCES**Simultaneous:** degradation products, n-butyl p-aminobenzoic acid, p-chloroaniline, chlorhexidine

KEY WORDS

pastilles

REFERENCEBauer,M.; Degude,C.; Mailhe,L. Simultaneous determination of chlorhexidine, tetracaine and their degradation products by ion-pair liquid chromatography, *J.Chromatogr.*, **1984**, 315, 457–464.

SAMPLE**Matrix:** formulations**Sample preparation:** Inject a 20 μ L aliquot.

HPLC VARIABLES**Column:** 220 \times 4.6 5 μ m Brownlee silica (Applied Biosystems)**Mobile phase:** MeOH:10 mM KH_2PO_4 adjusted to pH 4.0 with 10% phosphoric acid 25:75**Flow rate:** 1**Injection volume:** 20**Detector:** UV 235

CHROMATOGRAM**Retention time:** 13.3**Limit of detection:** 334 ng/mL

OTHER SUBSTANCES**Simultaneous:** morphine, hydromorphone

KEY WORDS

saline; injections

REFERENCEVenkateshwaran,T.G.; Stewart,J.T. HPLC determination of morphine-hydromorphone-bupivacaine and morphine-hydromorphone-tetracaine mixtures in 0.9% sodium chloride injection, *J.Liq.Chromatogr.*, **1995**, 18, 565–578.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Guard column:** 4 \times 4 5 μ m LiChrospher100RP-18**Column:** 125 \times 4 5 μ m Spherisorb ODS 2**Mobile phase:** MeCN:buffer 35:65 (Buffer was 20 mM sodium acetate containing 0.28% triethylamine, adjusted to pH 4.5 with acetic acid.)**Flow rate:** 1.5**Detector:** UV 280

CHROMATOGRAM**Retention time:** k' 2.85

OTHER SUBSTANCES**Simultaneous:** 4-butylaminobenzoic acid

REFERENCEYang,H.; Thyron,F.C. Determination of six pharmaceuticals and their degradation products in reversed-phase high performance liquid chromatography by using amine additives, *J.Liq.Chromatogr.Rel.Technol.*, **1998**, 21, 1347–1357.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Guard column:** 45 × 4.6 5 µm Ultrasphere ODS**Column:** 150 × 4.6 5 µm Ultrasphere ODS**Mobile phase:** MeCN:MeOH:water containing 2.5 mM hexanesulfonic acid 35:40:25, adjusted to pH 6.0 with 100 mM acetic acid**Flow rate:** 2**Detector:** UV 310

CHROMATOGRAM**Retention time:** 4.5**Internal standard:** procaine (2.6)**Limit of detection:** 800 pg**Limit of quantitation:** 5 ng

OTHER SUBSTANCES**Simultaneous:** metabolites

KEY WORDScomparison with capillary electrophoresis

REFERENCE

Asavapichayont,P.; Hu,J.; Foldvari,M. Development of an HPLC method for simultaneous analysis of tetracaine and its metabolite in dosage forms and biological fluids, with comparison to capillary electrophoresis method (Abstract 3307), *Pharm.Res.*, **1997**, *14*, S565.

SAMPLE**Matrix:** solutions**Sample preparation:** Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.

HPLC VARIABLES**Column:** 125 × 4.9 Spherisorb S5W silica**Mobile phase:** MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7**Flow rate:** 2**Injection volume:** 20**Detector:** E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM**Retention time:** 2.7

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazepine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipanone, diprenorphine, dipyrriamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, flupromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine,

mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J.Chromatogr.*, **1985**, 323, 191-225.

SAMPLE

Matrix: solutions

Sample preparation: Dissolve in MeOH:water 1:1 at a concentration of 50 µg/mL, inject a 10 µL aliquot.

HPLC VARIABLES

Column: 300 × 3.9 10 µm µBondapak C18

Mobile phase: MeOH:acetic acid:triethylamine:water 50:1.5:0.5:48

Flow rate: 1.5

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 14

REFERENCE

Roos, R.W.; Lau-Cam, C.A. General reversed-phase high-performance liquid chromatographic method for the separation of drugs using triethylamine as a competing base, *J.Chromatogr.*, **1986**, 370, 403-418.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a solution in EtOH, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 300 × 4.1 10 µm Versapak C18

Mobile phase: MeCN:250 mM pH 2.7 potassium phosphate buffer 30:70

Flow rate: 2

Injection volume: 20

Detector: UV 235

CHROMATOGRAM

Retention time: 4.9

OTHER SUBSTANCES

Simultaneous: desmethylsertraline, sertraline

REFERENCE

Wiener,H.L.; Kramer,H.K.; Reith,M.E.A. Separation and determination of sertraline and its metabolite, des-methylsertraline, in mouse cerebral cortex by reversed-phase high-performance liquid chromatography, *J.Chromatogr.*, **1990**, 527, 467-472.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Guard column: 30 × 3.2 7 µm SI 100 ODS (not commercially available)

Column: 150 × 3.2 7 µm SI 100 ODS (not commercially available)

Mobile phase: MeCN:buffer 31.2:68.8 (Buffer was 6.66 g KH₂PO₄ and 4.8 g 85% phosphoric acid in 1 L water, pH 2.3.)

Flow rate: 0.5-1

Detector: UV 222, UV309

CHROMATOGRAM

Retention time: 2.4

Internal standard: 5-(4-methylphenyl)-5-phenylhydantoin (7.3)

OTHER SUBSTANCES

Also analyzed: aspirin, caffeine, carbamazepine, chlordiazepoxide, chlorprothixene, clonazepam, diazepam, doxylamine, ethosuximide, furosemide, haloperidol, hydrochlorothiazide, methocarbamol, methotrimeprazine, nicotine, oxazepam, procaine, promazine, propafenone, propranolol, salicylamide, temazepam, thiopental, triamterene, verapamil, zolpidem, zopiclone

REFERENCE

Below,E.; Burrmann,M. Application of HPLC equipment with rapid scan detection to the identification of drugs in toxicological analysis, *J.Liq.Chromatogr.*, **1994**, 17, 4131-4144.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a solution in MeOH, inject a 25 µL aliquot.

HPLC VARIABLES

Guard column: 4 × 4 Lichrospher 100 RP-18

Column: 125 × 4 Lichrospher 100 RP-18

Mobile phase: MeOH:buffer 75:25 (Buffer was 320 mL 20 mM K₂HPO₄ + 680 mL 20 mM KH₂PO₄, pH 7.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 235

CHROMATOGRAM

Retention time: 9.40

OTHER SUBSTANCES

Simultaneous: benzoylecgonine, cocaine

REFERENCE

Fernández,P.; Rodríguez,P.; Bermejo,A.M.; López-Rivadulla,M.; Cruz,A. Simultaneous determination of cocaine and benzoylecgonine in vitreous humor by HPLC, *J.Liq.Chromatogr.*, **1994**, 17, 883-890.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200

mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, am-triptyline, amobarbital, amoxapine, amphetamine, amyllocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenoprofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-stilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isox-suprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, meth-apyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyl-testosterone, methylpyrrolon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, ox-ymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendi-metrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phenter-mine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, predni-solone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidone, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, sal-icylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sufadiazine, sul-fadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, tetracycline, tetramisole, theba-ine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thior-idazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexy-phenidyl, trimethoprim, triptelennamine, triprolidine, tropacocaine, tyramine, verapamil, vin-camine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, **1994**, *18*, 233-242.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4.5 µm LiChrospher 100 RP-8 (B)

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 10.81 (A), 5.44 (B)

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordi-azepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxy-chloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazin-dol, mefenamic acid, meperidine, mephénytoin, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, met-ronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymet-azoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, phen-iramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltolox-amine, phenytoin, pimizole, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quin-ine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, seco-barbital, sertraline, sotalol, spironolactone, sulfinpyrazone, sulindac, temazepam, terbutaline, terfenadine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, to-caicaine, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, trifluop-razine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yo-himbine, zopiclone


KEY WORDS

details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, 692, 103–119.

Tetracycline

Molecular formula: $C_{22}H_{24}N_2O_8$ **Molecular weight:** 444.44**CAS Registry No.:** 60-54-8, 6416-04-2 (trihydrate), 64-75-5 (HCl), 1336-20-5 (phosphate)**Merck Index:** 9337**Lednicer No.:** 1 212**SAMPLE****Matrix:** blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the