


CHROMATOGRAM**Retention time:** 2.6**OTHER SUBSTANCES**

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazacine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipanone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserin, laudanose, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclophenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypropazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozide, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripelennamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, 323, 191–225.

Thiamine

**Molecular formula:** C₁₂H₁₇ClN₄OS**Molecular weight:** 300.81**CAS Registry No.:** 59-43-8, 67-03-8 (HCl), 532-43-4 (mononitrate)**Merck Index:** 9430**SAMPLE****Matrix:** blood

Sample preparation: 1 mL Hemolysate + 30 µL 4 µM IS in 100mM HCl, shake thoroughly. Slowly add 2 mL MeOH, mix, let stand for 30 min. Centrifuge at 2000 g for 10 min. Add 50 µL freshly prepared 30.4 mM potassium ferricyanide and 50 µL 0.8 mM NaOH to 1 mL supernatant. Filter (0.45 µm) and inject a 50 µL aliquot.

HPLC VARIABLES**Guard column:** 50 × 4.0 Spherisorb NH2**Column:** 125 × 4.0 5 µm Spherisorb NH2**Mobile phase:** MeOH:100 mM pH 7.5 potassium phosphate buffer 45:55**Flow rate:** 1**Injection volume:** 50**Detector:** F ex 375 em 430

CHROMATOGRAM**Retention time:** 4.0**Internal standard:** acetylcholine (3.0)**Limit of detection:** 2 nM

KEY WORDSerythrocytes

REFERENCE

Lynch, P.L.M.; Trimble, E.R.; Young, I.S. High-performance liquid chromatographic determination of thiamine diphosphate in erythrocytes using internal standard methodology, *J. Chromatogr. B*, **1997**, 701, 120–123.

SAMPLE**Matrix:** blood**Sample preparation:** 200 µL Plasma, whole blood, or erythrocytes + 200 µL 100 mg/mL trichloroacetic acid, vortex vigorously, centrifuge at 35000 g for 5 min, inject a 100 µL aliquot of the supernatant.

HPLC VARIABLES**Column:** 250 × 4 µm Bondapak C18**Mobile phase:** MeCN:buffer 3.8:96.2 (Mobile phase was 200 mM NaH₂PO₄ in 3 g/L MeCN in water.)**Flow rate:** 1**Injection volume:** 100**Detector:** F ex 375 em 450 following post-column reaction. The column effluent mixed with the reagent pumped at 0.5 mL/min and flowed to the detector. (Reagent was 100 µg/mL potassium ferricyanide in 150 g/L NaOH.)

CHROMATOGRAM**Retention time:** 8.0 (thiamine), 3.1 (thiamine triphosphate), 3.8 (thiamine pyrophosphate), 5.0 (thiamine monophosphate)**Limit of detection:** 30 fmole

KEY WORDSpost-column reaction; pharmacokinetics; plasma; whole blood; erythrocytes

REFERENCE

Kimura, M.; Itokawa, Y. Determination of thiamin and thiamin phosphate esters in blood by liquid chromatography with post-column derivatization, *Clin. Chem.*, **1983**, 29, 2073–2075.

SAMPLE**Matrix:** blood**Sample preparation:** Hemolyze whole blood by freezing at -20° for 20 min, thaw, homogenize. Add a 200 µL aliquot of hemolyzed blood or serum to 200 µL chilled 10% perchloric acid, let stand below 4° for 15 min, centrifuge at 10000 g for 1 min. Remove a 200 µL aliquot of the supernatant and add it to 1.8 M sodium acetate containing 600 mM NaOH, mix, filter (Costar filter unit) while centrifuging for 30 s, inject a 20 µL aliquot of the filtrate.

HPLC VARIABLES**Guard column:** present but not specified**Column:** 110 × 4.7 Partisphere 5 C18 (Whatman)**Mobile phase:** Gradient. A was 15 mM citric acid adjusted to pH 4.2 with 50% ammonium hydroxide, prepare fresh each day. B was 100 mM formic acid containing 4% diethylamine, pH

3.2. A:B from 90:10 to 50:50 over 2.5 min, to 5:95 over 0.5 min, maintain at 5:95 over 3.5 min, to 95:5 over 0.5 min, maintain at 95:5 for 3 min.

Flow rate: 1

Injection volume: 20

Detector: F ex 365 em 435 following post-column reaction. The column effluent mixed with the reagent pumped at 0.2 mL/min and the mixture flowed through a 90 cm \times 0.5 mm ID PTFE coil to the detector. (Prepare reagent by dissolving 100 mg potassium ferricyanide in 120 mL 3 M NaOH.)

CHROMATOGRAM

Retention time: 7.6

Limit of quantitation: 2 nM

OTHER SUBSTANCES

Extracted: thiamine monophosphate, thiamine pyrophosphate, thiamine triphosphate

KEY WORDS

post-column reaction; whole blood; serum

REFERENCE

Lee,B.L.; Ong,H.Y.; Ong,C.N. Determination of thiamine and its phosphate esters by gradient-elution high-performance liquid chromatography, *J.Chromatogr.*, **1991**, 567, 71–80.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 150 μ L 3 M perchloric acid, vortex, centrifuge at >1500 g. Remove 500 μ L of the supernatant and add it to 300 μ L 1 M pH 4.6 acetate buffer, add 100 μ L 10 mg/mL acidic phosphatase (2 U/mg, grade II, Boehringer Mannheim) in water, heat at 40° for 16 h, add 150 μ L 3 M perchloric acid, vortex, centrifuge at >1500 g, inject a 20 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 125 \times 4 Nucleosil 120 5 C18

Mobile phase: MeCN:buffer 25:75 (Buffer was 10 mM perchloric acid containing 10 mM octanesulfonic acid.)

Flow rate: 2

Injection volume: 20

Detector: F ex 365 em 435 following post-column reaction. The column effluent mixed with 0.8 g/L potassium ferricyanide in 3 M NaOH pumped at 1 mL/min, the mixture flowed through a 10 m \times 0.3 mm i.d. PTFE coil at 30° to the detector.

CHROMATOGRAM

Retention time: 2.3

Limit of detection: 2 ng/mL

KEY WORDS

post-column reaction; plasma; pharmacokinetics

REFERENCE

Mascher,H.; Kikuta,C. High-performance liquid chromatographic determination of total thiamine in human plasma for oral bioavailability studies, *J.Pharm.Sci.*, **1993**, 82, 56–59.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation.

Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 2.597

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

SAMPLE

Matrix: feed

Sample preparation: 1 g Ground feed + 20 mL 100 mM HCl, shake vigorously, heat on a boiling water bath for 30 min (shake every 5 min), cool in an ice bath for 5 min, centrifuge at 1000 rpm for 10 min. Remove a 5 mL aliquot and make it up to 50 mL with buffer, centrifuge at 1000 rpm for 5 min, inject a 10 µL aliquot. (Buffer was water adjusted to pH 4.0 with acetic acid.)

HPLC VARIABLES

Column: 250 × 4.6 SynChropack SCD-100 (SynChrom Inc.)

Mobile phase: MeOH:water 40:60 containing 50 mM sodium pentanesulfonate, pH adjusted to 4.0 with acetic acid

Flow rate: 1.5

Injection volume: 10

Detector: F ex 370 em 430 following post-column derivatization. The column effluent was mixed with 200 mM KOH and 0.01% potassium ferricyanide, each pumped at 0.5 mL/min. The mixture flowed in the dark through a 3 m × 0.8 mm i.d. knotted coil of PTFE tubing to the detector.

CHROMATOGRAM

Retention time: 5

Limit of detection: 5 pg

KEY WORDS

post-column reaction; derivatization

REFERENCE

Gehring, T.A.; Cooper, W.M.; Holder, C.L.; Thompson, H.C., Jr. Liquid chromatographic determination of thiamine in rodent feed by postcolumn derivatization and fluorescence detection, *JAOAC Int.*, **1995**, 78, 307-309.

SAMPLE

Matrix: food

Sample preparation: Condition a CBA (methylcarboxylate in acid form) SPE cartridge with 1 mL MeOH and 1 mL 10 mM pH 4.0 phosphate buffer. Add a sample of finely-ground food to 20 mL 100 mM HCl and heat at 100° for 30 min. Cool and adjust pH to 4.4-4.5 with sodium acetate. Add a 6 mg/mL solution of takadiastase (Fluka). Heat at 47° for 3 h. Filter through a

cellulose acetate filter and dilute with water to 50 mL. Add a 2 mL aliquot to the SPE cartridge. Wash twice with 500 μ L portions of pH 4.0 phosphate buffer. Elute with three 200 μ L portions of 100 mM barium chloride. Inject an aliquot.

HPLC VARIABLES

Column: 125 \times 4.5 μ m LiChrospher 100RP-18

Mobile phase: MeOH:10 mM pH 2.8 phosphate buffer:triethylamine 25:85:0.1 containing 5 mM hexanesulfonic acid

Flow rate: 1

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 6.25

Limit of detection: 408 ng/g

KEY WORDS

baby food; cereal; dietetic cookies; SPE

REFERENCE

Blanco,D.; Llanaez,M.B.; Gutierrez,M.D. A paired ion liquid chromatographic method for thiamine determination in selected foods, *J.Liq.Chromatogr.Rel.Technol.*, **1996**, 19, 2155–2164.

SAMPLE

Matrix: formula, milk

Sample preparation: Mix 8.0 g powdered infant milk with 10 mL water to it. Mix the diluted powder or 10.5 g liquid infant milk with 1 g solid trichloroacetic acid, shake thoroughly with magnetic stirring for 10 min, centrifuge at 1250 g for 10 min, add 3 mL 4% trichloroacetic acid to the solid residue, mix thoroughly for 10 min, centrifuge, discard the solid phase. Combine the two acid extracts and make up to 10 mL with 4% trichloroacetic acid, filter (0.45 μ m), inject an aliquot of the filtrate.

HPLC VARIABLES

Guard column: 5 μ m Tracer Spherisorb ODS 2 C18 (Teknokroma, Spain)

Column: 250 \times 4.6 5 μ m Tracer Spherisorb ODS 2 C18 (Teknokroma, Spain)

Mobile phase: MeOH:buffer 15:85 (Buffer was 5 mM octanesulfonic acid and 0.5% triethylamine, pH 3.6.)

Flow rate: 1

Injection volume: 20

Detector: UV 261 for 6 min, UV 287 for 2 min, UV 290 for 5 min, UV 282 for 3 min, UV 268 for 3.5 min, UV 361 for 20.5 min, UV 246 for 20 min

CHROMATOGRAM

Retention time: 48

Limit of quantitation: \leq 100 ng/mL

OTHER SUBSTANCES

Extracted: riboflavin, pyridoxine, vitamin B12, folic acid, niacinamide, pyridoxal, pyridoxamine

REFERENCE

Albal-Hurtado,S.; Veciana-Nogués,M.; Izquierdo-Pulido,M.; Mariné-Font,A. Determination of water-soluble vitamins in infant milk by high-performance liquid chromatography, *J.Chromatogr.A*, **1997**, 778, 247–253.

SAMPLE

Matrix: formulations

Sample preparation: Sonicate 75 mg powdered tablets with 25 mL mobile phase for 15 min, filter (paper), inject a 135 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Ultrabase C18 (Scharlau Science, Spain)

Mobile phase: MeOH:20 mM KH_2PO_4 30:70 adjusted to pH 4.0 with orthophosphoric acid

Flow rate: 1.5
Injection volume: 135
Detector: UV 246

CHROMATOGRAM

Retention time: 2.3
Limit of quantitation: 1.9 µg/mL

OTHER SUBSTANCES

Simultaneous: aspirin, caffeine, salicylic acid

KEY WORDS

tablets

REFERENCE

Gámiz-Gracia,L.; Luque de Castro,M.D. An HPLC method for the determination of vitamin B1, caffeine, acetylsalicylic acid, and the impurities of salicylic acid in a pharmaceutical preparation, *J.Liq.Chromatogr.Rel.Technol.*, **1997**, 20, 2123–2133.

SAMPLE

Matrix: formulations

Sample preparation: Pulverize tablets if necessary. Add tablets to 100 mL 5 mM pH 4.5 potassium phosphate buffer, sonicate at 75 W for 2 min, cool to room temperature, make up to 200 mL with buffer, filter (0.45 µm), inject a 10 µL aliquot of the filtrate.

HPLC VARIABLES

Column: 250 × 4.6 10 µm LiChrosorb NH2 aminopropyl

Mobile phase: MeCN:5 mM KH₂PO₄ 87:13 (Wash column with MeCN:water 10:90 at the end of the day.)

Column temperature: 25

Flow rate: 2

Injection volume: 10

Detector: UV 210

CHROMATOGRAM

Retention time: 3.5

OTHER SUBSTANCES

Simultaneous: pantothenic acid, riboflavin, niacinamide, pyridoxine

KEY WORDS

tablets

REFERENCE

Hudson,T.J.; Allen,R.J. Determination of pantothenic acid in multivitamin pharmaceutical preparations by reverse-phase high-performance liquid chromatography, *J.Pharm.Sci.*, **1984**, 73, 113–115.

SAMPLE

Matrix: formulations

Sample preparation: Tablets without iron. Grind 5 tablets to a fine powder, add 10 mL monothioglycerol and 800 mL buffer, sonicate for 30 min, add 150 mL MeOH, make up to 1 L with buffer, filter (GF/C paper), discard first few mL, remove a 10 mL aliquot, make up to 25 mL with mobile phase, inject an aliquot. Tablets with dioctyl sodium sulfosuccinate. Grind 5 tablets to a fine powder, add 10 mL 2-monothioglycerol and 1 g barium chloride, make up to 1 L with buffer, stir vigorously for 30 min, filter (GF/C paper), discard first few mL, inject an aliquot. Capsules with iron. Contents of one capsule + 5 mL 2-monothioglycerol + 2 mL glacial acetic acid + 75 mL buffer, sonicate for 5 min, make up to 100 mL with buffer, stir vigorously for 30 min, filter (GF/C paper), add 300 mg cupferron, stir for 10 min, let stand for 1 h at room temperature, filter (GF/C paper), let stand for 30 min, filter again (if necessary), discard first few mL, inject an aliquot. (Buffer was 48 mL glacial acetic acid and 10 mL triethylamine in 1

L water, adjust pH to 3.6 ± 0.05 with acetic acid or triethylamine, make up to 1.7 L with water.)

HPLC VARIABLES

Column: 100 \times 8 Radial Pak A C18 (Waters)

Mobile phase: MeOH:buffer 15:85 (Buffer was 2.20 g sodium heptanesulfonate, 100 mg EDTA, 48 mL glacial acetic acid, and 10 mL triethylamine made up to 1.7 L with water, adjust pH to 3.6 ± 0.05 with acetic acid or triethylamine.)

Flow rate: 2

Injection volume: 10

Detector: UV 280

CHROMATOGRAM

Retention time: 9

OTHER SUBSTANCES

Simultaneous: niacinamide, riboflavin, pyridoxine, ascorbic acid (UV 254)

KEY WORDS

multi-vitamin; protect from light; tablets; capsules

REFERENCE

Lam,F.-L.; Holcomb,I.J.; Fusari,S.A. Liquid chromatographic assay of ascorbic acid, niacinamide, pyridoxine, thiamine, and riboflavin in multivitamin-mineral preparations, *J.Assoc.Off.Anal.Chem.*, **1984**, 67, 1007–1011.

SAMPLE

Matrix: formulations

Sample preparation: Dilute injections with water, inject a 50 μ L aliquot. Dissolve tablets or capsule contents in water (warm if necessary), filter (0.5 μ m PTFE), inject a 50 μ L aliquot of the filtrate. (Dissolve tablets or other formulations containing proteinaceous material in water at 60°, add 5% trichloroacetic acid (to pH 4.4), filter, inject a 50 μ L aliquot.)

HPLC VARIABLES

Guard column: pellicular Corasil

Column: 10 μ m μ Bondapak C18

Mobile phase: Gradient. A was prepared by dissolving 1 g sodium dioctylsulfosuccinate in 170 mL MeOH, add 10 mL concentrated formic acid, make up to 800 mL with water, adjust pH to 2.5 with 1 M KOH, make up to 1 L. B was prepared by dissolving 1 g sodium dioctylsulfosuccinate in 450 mL MeOH, add 10 mL concentrated formic acid, make up to 800 mL with water, adjust pH to 4.6, make up to 1 L. A:B 100:0 for 19 min then 0:100 (step gradient) or A: B from 100:0 to 0:100 over 25 min (concave curve 9), maintain at 0:100 for 3 min, return to initial conditions over 2 min.

Flow rate: 1.5

Injection volume: 50

Detector: UV 254

CHROMATOGRAM

Retention time: 31 (step gradient), 24 (curve gradient)

OTHER SUBSTANCES

Simultaneous: folic acid (UV 280), niacin, niacinamide, pyridoxamine (UV 280), riboflavin, pyridoxine (UV 280), ascorbic acid (UV 280)

KEY WORDS

injections; capsules; tablets

REFERENCE

Woollard,D.C. New ion-pair reagent for the high-performance liquid chromatographic separation of B-group vitamins in pharmaceuticals, *J.Chromatogr.*, **1984**, 301, 470–476.

SAMPLE**Matrix:** formulations**Sample preparation:** Weigh out 500 mg ground tablets, extract with water, make up to 50 or 100 mL with water, filter, inject an aliquot.

HPLC VARIABLES**Column:** 250 × 4.6 Nucleosil 10 C18**Mobile phase:** MeOH:1% acetic acid 25:75**Flow rate:** 1**Injection volume:** 20**Detector:** UV 270

CHROMATOGRAM**Retention time:** 2

OTHER SUBSTANCES**Simultaneous:** menadione hydrogen sulfite, niacinamide, pyridoxine, riboflavin, ascorbic acid

KEY WORDS

tablets; multi-vitamin

REFERENCE

Sadlej-Sosnowska, N.; Blitek, D.; Wilczynska-Wojtulewicz, I. Determination of menadione sodium hydrogen sulphite and nicotinamide in multivitamin formulations by high-performance liquid chromatography, *J. Chromatogr.*, **1986**, *357*, 227–232.

SAMPLE**Matrix:** formulations

HPLC VARIABLES**Column:** 100 × 4.3 µm Hypersil BDS-C18**Mobile phase:** Gradient. MeCN:water adjusted to pH 2.1 from 0.3:99.7 to 25:75 over 11 min**Flow rate:** 0.5**Detector:** UV 220

CHROMATOGRAM**Retention time:** 2

OTHER SUBSTANCES**Simultaneous:** biotin, caffeine, citric acid, folic acid, niacinamide, niacin, pantothenic acid, pyridoxine, riboflavin, saccharin, vitamin B12, ascorbic acid

KEY WORDS

tablets

REFERENCE

Hewlett Packard Leaflet 12-5091-7351 EUS, **1993**.

SAMPLE**Matrix:** formulations**Sample preparation:** Dilute liquid multivitamin formulations, filter (0.45 µm), inject an aliquot.

HPLC VARIABLES**Column:** 250 × 4.5 µm Lichrosorb RP-8**Mobile phase:** Gradient. A was 10 mM KH₂PO₄ containing 5 mM sodium hexanesulfonate adjusted to pH 2.8 with phosphoric acid. B was MeOH. A:B from 90:10 to 71.8:28.2 over 4 min, maintain at 71.8:28.2 for 1.5 min, to 50:50 over 6.5 min, maintain at 50:50 for 5 min, return to initial conditions over 5 min**Flow rate:** 1

Injection volume: 5

Detector: UV 272

CHROMATOGRAM

Retention time: 9.90

Internal standard: theobromine (8)

Limit of detection: 0.430 ng

OTHER SUBSTANCES

Simultaneous: folic acid, niacin, niacinamide, riboflavin, pyridoxine (UV 290)

KEY WORDS

liquid multivitamins; degas solutions with helium; protect from light

REFERENCE

Blanco,D.; Sánchez,L.A.; Gutiérrez,M.D. Determination of water soluble vitamins by liquid chromatography with ordinary and narrow-bore columns, *J.Liq.Chromatogr.*, **1994**, 17, 1525–1539.

SAMPLE

Matrix: rice

Sample preparation: Heat ground rice with at least 10 volumes of 100 mM HCl at 95–100° for 30 min with frequent mixing, cool, dilute to a thiamine concentration of 200 ng/mL with 100 mM HCl, adjust the pH of a 65 mL aliquot to 4.0–4.5 with about 5 mL 2 M sodium acetate, add 5 mL 10% takadiastase in water, heat at 45–50° for 3 h, adjust pH to 3.5, make up to 100 mL with water, filter (paper) (AOAC Official Methods of Analysis, 1990, 1049), centrifuge at 3500 rpm for 20 min, inject a 25 µL aliquot of the supernatant.

HPLC VARIABLES

Column: 150 × 4 Nucleosil 5 C18

Mobile phase: 10 mM NaH₂PO₄ containing 150 mM sodium perchlorate, adjusted to pH 2.2 with perchloric acid

Column temperature: 55

Flow rate: 0.6

Injection volume: 25

Detector: F ex 375 em 435 following post-column reaction. The column effluent mixed with 0.1% potassium hexacyanoferrate(III) in 12% NaOH pumped at 0.6 mL/min and this mixture flowed through a 30 cm × 0.8 mm ID stainless steel coil at 55° to the detector.

CHROMATOGRAM

Retention time: 8.5

KEY WORDS

post-column reaction

REFERENCE

Ohta,H.; Baba,T.; Suzuki,Y.; Okada,E. High-performance liquid chromatographic analysis of thiamine in rice flour with fluorimetric post-column derivatization, *J.Chromatogr.*, **1984**, 284, 281–284.

SAMPLE

Matrix: rice

Sample preparation: Grind rice to pass 30-mesh screen. Stir 3 g ground rice and 50 mL MeOH: 100 mM HCl 40:60 with a glass rod until homogeneous, reflux for 30 min, vortex for 1 min, sonicate for 20 min, centrifuge at 3000 g for 20 min, filter (0.45 µm) the supernatant, inject a 5 µL aliquot of the filtrate.

HPLC VARIABLES

Column: 250 × 4.6 Zorbax TMS

Mobile phase: 10 mM NaH₂PO₄ containing 500 mM sodium perchlorate adjusted to pH 2.5 with 3 M perchloric acid

Column temperature: 55

Flow rate: 0.4

Injection volume: 5

Detector: F ex 375 em 435 following post-column reaction. The column effluent mixed with reagent pumped at 0.4 mL/min, the mixture flowed through a 300×0.8 stainless steel mixing coil to the detector. (Reagent was 0.1% potassium hexacyanoferrate(III) in 15% NaOH.)

CHROMATOGRAM

Retention time: 8

Limit of quantitation: 1.5 µg/g

KEY WORDS

post-column reaction

REFERENCE

Ohta,H.; Maeda,M.; Nogata,Y.; Yoza,K.-I.; Takeda,Y.; Osajima,Y. A simple determination of thiamine in rice (*Oryza sativa* L.) by high-performance liquid chromatography with post-column derivatization, *J.Liq.Chromatogr.*, **1993**, 16, 2617-2629.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250×4.6 5 µm Accubond Amino (J & W)

Mobile phase: MeCN:buffer 10:90 (Buffer was 20 mM phosphoric acid adjusted to pH 3.0 with 20 mM NaOH.)

Flow rate: 1

Detector: UV 254

CHROMATOGRAM

Retention time: 1.2

OTHER SUBSTANCES

Simultaneous: p-aminobenzoic acid, niacinamide, pyridoxal, pyridoxamine, riboflavin, pyridoxine, vitamin B12

REFERENCE

J & W Catalog, 1992-3, p. 277.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 33×4.6 3 µm Supelcosil LC-8-DB

Mobile phase: MeOH:buffer 15:85 (Buffer was 4.3 mM sodium hexanesulfonate containing 0.1% triethylamine, adjusted to pH 2.8 with phosphoric acid.)

Column temperature: 35

Flow rate: 1

Detector: UV 200

CHROMATOGRAM

Retention time: 2.2

OTHER SUBSTANCES

Simultaneous: niacin, pantothenic acid, pyridoxine, riboflavin, niacinamide, ascorbic acid

REFERENCE

Rainin Catalog, C1-94, **1994**, p. 780.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 100 × 4.6 Spheri-5 RP-8

Mobile phase: Gradient. A was 100 mM pH 4.7 acetate buffer. B was MeCN:100 mM pH 4.7 acetate buffer 25:75.

Column temperature: 26

Flow rate: 4

Detector: UV 254

CHROMATOGRAM

Retention time: 2.2

OTHER SUBSTANCES

Simultaneous: niacin, pyridoxine, riboflavin, niacinamide, ascorbic acid

REFERENCE

Rainin Catalog, C1-94, 1994, p. 7.21.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitrityline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbitol, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenpropofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocodone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-stilbene, imipramine, indomethacin, isocarboxystyl, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methylpyrrol, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nyldrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, predni-

solone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopalamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, **1994**, *18*, 233–242.

SAMPLE

Matrix: tissue

Sample preparation: Homogenize (motor-driven glass homogenizer) 100 mg tissue with 300 μ L 5% trichloroacetic acid at 4°, centrifuge at 4° at 5000 g for 1 h. Wash the supernatant with 3 volumes of water-saturated diethyl ether for 1 h. Remove an 80 μ L aliquot of the aqueous phase and add it to 50 μ L reagent, mix for 5–10 s, inject an aliquot 1 min after the addition of the reagent. (Prepare reagent by mixing 50 μ L 10 mg/mL potassium ferricyanide with 2.5 mL 15% NaOH, store in the dark, discard after 1 day.)

HPLC VARIABLES

Guard column: 4.2 \times 3.2 30–44 μ m 201RP (Vydac)

Column: 150 or 250 \times 4.6 5 μ m Ultrasphere ODS

Mobile phase: Gradient. MeOH:25 mM pH 8.4 phosphate buffer 10:90 for 1 min, to 100:0 over 3 min, maintain at 100:0 for 2 min, return to initial conditions over 2 min.

Flow rate: 1

Injection volume: 20

Detector: F ex 390 em 475

CHROMATOGRAM

Retention time: 4.2 (thiamine triphosphate), 4.5 (thiamine pyrophosphate), 5 (thiamine monophosphate), 6 (thiamine)

Limit of detection: 0.05 pmole

KEY WORDS

derivatization; rat; nerve; heart

REFERENCE

Bontemps, J.; Philippe, P.; Bettendorff, L.; Lombet, J.; Dandrifosse, G.; Schoffeniels, E.; Crommen, J. Determination of thiamine and thiamine phosphates in excitable tissues as thiochrome derivatives by reversed-phase high-performance liquid chromatography on octadecyl silica, *J. Chromatogr.*, **1984**, *307*, 283–294.

SAMPLE

Matrix: tissue

Sample preparation: Condition a 500 mg Nucleosil C18 SPE cartridge with 2 mL MeOH, 2 mL MeOH containing 5 mM sodium heptanesulfonate, and two 2 mL portions of water. Suspend 5 g homogenized tissue with 35 mL 10 mM HCl, autoclave at 121° for 30 min, add 2 mL 25 mg/mL taka-diestase (Fluka) in 2.5 M sodium acetate, add 2 mL 10 (muscle) or 20 (liver) mg/mL clara-diastase (Fluka) in water, add 2 mL 50 mg/mL papain (Merck) in water, adjust pH to 4.5, heat at 37° for 16–18 h, filter (paper), adjust pH to 6.5, filter again, make up to 50 mL with water, add 4 mL to the SPE cartridge, wash with 2 mL MeOH:water 20:80 containing 5 mM sodium heptanesulfonate, elute with 2 mL MeOH:water 50:50 containing 5 mM sodium heptanesulfonate, inject a 50 μ L aliquot.

HPLC VARIABLES

Guard column: 20 \times 4.6 10 μ m Nucleosil C18

Column: 150 \times 4.6 3 μ m Nucleosil C18

Mobile phase: MeCN:10 mM pH 3.0 KH_2PO_4 16:84 (muscle) or 15:85 (liver) containing 5 mM sodium heptane sulfonate (Wash with MeCN:water 20:80 at the end of the day, store column in MeCN.)

Column temperature: 45

Injection volume: 50

Detector: UV 254

CHROMATOGRAM

Retention time: 6.5 (muscle), 8 (liver)

Limit of detection: 500 ng/g

OTHER SUBSTANCES

Extracted: riboflavin

KEY WORDS

pig; muscle; liver; protect from light; SPE

REFERENCE

Barna,I.; Dworschák,E. Determination of thiamine (vitamin B1) and riboflavin (vitamin B2) in meat and liver by high-performance liquid chromatography, *J.Chromatogr.A*, **1994**, 668, 359–363.

SAMPLE

Matrix: yeast

Sample preparation: Condition CM-cellulose (44-105 μm CM-Cellulofine C-200, Biochemical Industry, Tokyo) by washing with 200 mM HCl, 200 mM NaOH, water, and 200 mM HCl, rinse with water until the rinses are neutral. Add a 2.5 mL aliquot to a 170×10 column, wash with 300 mM phosphoric acid, wash with water until the eluate is neutral. Heat 1 g dried yeast, 1 mL 10% HCl, and 80 mL water with frequent shaking at 80-85° for 30 min, cool, make up to 100 mL with water, centrifuge for 10 min. Remove a 4 mL aliquot of the supernatant and add it to 5 mL 200 mM pH 4.5 acetic acid/sodium acetate buffer and 1 mL 30 mg/mL Taka-diestase supernatant (Sankyo) in 5 mM HCl, heat at 45-50° for 3 h, add 2 mL of this mixture to the column at 0.5 mL/min, wash with two 10 mL portions of water at 1 mL/min, elute with two 2.5 mL portions of 300 mM phosphoric acid at 0.5 mL/min. Add 1 mL 1 $\mu\text{g/mL}$ IS and 10 mg sodium 1-octanesulfonate to the eluate, mix, inject a 200 μL aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 Capcell-pak C18 (Shiseido)

Mobile phase: MeCN:buffer 20:80 (Buffer was 20 mM pH 3.5 phosphate buffer containing 0.2% sodium 1-octanesulfonate.)

Column temperature: 40

Flow rate: 1

Injection volume: 200

Detector: UV 254

CHROMATOGRAM

Retention time: 7.5

Internal standard: phenacetin (10)

Limit of quantitation: 250 ng/mL

KEY WORDS

SPE

REFERENCE

Yamanaka,K.; Matsuoka,M.; Banno,K. Determination of thiamine in dried yeast by high-performance liquid chromatography using a clean-up column of CM-cellulose, *J.Chromatogr.A*, **1996**, 726, 237–240.

Thiamphenicol


Molecular formula: $C_{12}H_{15}Cl_2NO_3S$

Molecular weight: 356.23

CAS Registry No.: 15318-45-3

Merck Index: 9436

Lednicer No.: 2 45



SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 200.5

CHROMATOGRAM

Retention time: 7

KEY WORDS

whole blood

REFERENCE

Gaillard, Y.; Pépin, G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J. Chromatogr. A*, **1997**, 763, 149-163.

Thiamylal


Molecular formula: $C_{12}H_{18}N_2O_2S$

Molecular weight: 254.35

CAS Registry No.: 77-27-0, 337-47-3 (Na salt)

Merck Index: 9437

Lednicer No.: 1 274



SAMPLE

Matrix: bile, blood, gastric contents, tissue, urine

Sample preparation: Homogenize tissue with four volumes 100 mM pH 5 phthalate buffer. 2 mL Blood, bile, urine, stomach contents, or tissue homogenate + 2 mL 100 mM pH 5 phthalate buffer + 100 μ L 1 mg/mL phenolphthalein in MeOH, vortex for 10 s, add to a Clin Elut SPE

cartridge (Analytichem), let stand for 10 min, elute with three 5 mL portions of dichloromethane. Evaporate the eluate to dryness under a stream of air at 60°, reconstitute with 200 μ L MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m ODS (Altex)

Mobile phase: MeOH:water 5:1

Flow rate: 2

Injection volume: 20

Detector: UV 290

CHROMATOGRAM

Internal standard: phenolphthalein

KEY WORDS

SPE; liver; kidney

REFERENCE

Costantino,A.G.; Caplan,Y.H.; Levine,B.S.; Dixon,A.M.; Smialek,J.E. Thiamylal: review of the literature and report of a suicide, *J.Forensic Sci.*, **1990**, 35, 89–96.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Plasma + 100 μ L MeCN, vortex for 10 s, let stand for 10 min, vortex for 10 s, centrifuge at 12000 g for 2 min, inject a 20 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 150 \times 6 5 μ m Shim-pack CLC-ODS (Shimadzu)

Mobile phase: MeCN:water 55:45

Flow rate: 1.2

Injection volume: 20

Detector: UV 288

CHROMATOGRAM

Retention time: 4.9

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Extracted: thiopental

Noninterfering: acetaminophen, allopurinol, amikacin, amobarbital, amphotericin B, ampicillin, aspirin, barbital, caffeine, carbenicillin, chloramphenicol, chlorpromazine, cimetidine, cisplatin, cyclophosphamide, cyclosporin A, cytarabine, dactinomycin, doxorubicin, droperidol, ethosuximide, 5-fluorocytosine, 5-fluorouracil, furosemide, gentamicin, hexobarbital, ketamine, ketoconazole, 6-mercaptopurine, metharbital, methotrexate, miconazole, mizoribine, pentobarbital, phenobarbital, procainamide, secobarbital, tegafur, vancomycin

KEY WORDS

pharmacokinetics; plasma

REFERENCE

Hosotsubo,H.; Takeda,K.; Hosotsubo,K.; Yoshiya,I. Measurement of thiamylal in human plasma using reversed-phase high-performance liquid chromatography, *J.Chromatogr.*, **1989**, 487, 204–209.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Serum or plasma + 100 μ L MeOH + 1 mL 70 mM pH 6.4 phosphate buffer + 5 mL n-pentane, shake vigorously for 10 min, centrifuge at 1500 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 200 μ L mobile phase, vortex for 1 min, inject a 5 μ L aliquot.

HPLC VARIABLES**Column:** Nucleosil 5 C8**Mobile phase:** MeCN:buffer 35:65 (Buffer was 0.2 mM phosphoric acid containing 0.175 mM KH_2PO_4)**Flow rate:** 1.2**Injection volume:** 5**Detector:** UV 290

CHROMATOGRAM**Retention time:** 12**Internal standard:** thiamylal

OTHER SUBSTANCES**Extracted:** thiopental

KEY WORDSserum; plasma; cow; human; comparison with capillary electrophoresis; thiamylal is IS

REFERENCEMeier,P.; Thormann,W. Determination of thiopental in human serum and plasma by high-performance capillary electrophoresis-micellar electrokinetic chromatography, *J.Chromatogr.*, **1991**, 559, 505–513.

SAMPLE**Matrix:** blood**Sample preparation:** Mix serum with an equal volume of 1 M pH 5.0 phosphate buffer, add IS, add 3 mL ethyl acetate, rotate for 10 min, centrifuge at 3500 rpm for 5 min. Remove 2 mL of the organic layer and evaporate it to dryness, reconstitute the residue in 200 μL mobile phase, inject a 10-20 μL aliquot.

HPLC VARIABLES**Column:** 150 \times 6 Shim-pack CLC-ODS (Shimadzu)**Mobile phase:** MeOH:10 mM pH 5.0 sodium phosphate buffer 60:40**Flow rate:** 1**Injection volume:** 10-20**Detector:** UV

CHROMATOGRAM**Internal standard:** 5-(p-methylphenyl)-5-phenylhydantoin**Limit of quantitation:** 1 $\mu\text{g/mL}$

KEY WORDSserum; rat; pharmacokinetics

REFERENCENakashima,E.; Matsushita,R.; Ohshima,T.; Tsuji,A.; Ichimura,F. Quantitative relationship between structure and peritoneal membrane transport based on physiological pharmacokinetic concepts for acidic drugs, *Drug Metab.Dispos.*, **1995**, 23, 1220–1224.

SAMPLE**Matrix:** blood**Sample preparation:** Condition a 1 mL 50 mg Bond Elut C18 SPE cartridge with two 1 mL portions of MeOH and two 1 mL portions of 10 mM KH_2PO_4 . 100 μL Serum + 50 μL 25 $\mu\text{g/mL}$ n-propyl p-hydroxybenzoate in water + 500 μL 10 mM KH_2PO_4 , add to the SPE cartridge, wash with two 1 mL portions of water, elute with 300 μL MeOH. Evaporate the eluate, reconstitute in 100 μL EtOH:water 50:50, inject a 50 μL aliquot.

HPLC VARIABLES**Column:** 75 \times 4.6 TSK-gel ODS 80TM**Mobile phase:** EtOH:10 mM KH_2PO_4 15:85 containing 17 mM β -cyclodextrin (Suspend 17 mmoles β -cyclodextrin in 150 mL EtOH, make up to 1 L with 10 mM KH_2PO_4 .)**Column temperature:** 25

Flow rate: 0.9
Injection volume: 50
Detector: UV 288

CHROMATOGRAM

Retention time: 37 (S(-)), 39 (R(+))
Internal standard: n-propyl p-hydroxybenzoate (25)
Limit of detection: 50 ng/mL

OTHER SUBSTANCES

Noninterfering: muscle relaxants, benzodiazepines

KEY WORDS

serum; chiral; SPE

REFERENCE

Sueyasu,M.; Ikeda,T.; Otsubo,K.; Taniyama,T.; Aoyama,T.; Oishi,R. Enantioselective determination of thiamylal in human serum by high-performance liquid chromatography, *J.Chromatogr.B*, **1995**, 665, 133–137.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 125 × 4.6 5 µm C-18 (Perkin-Elmer)
Mobile phase: MeOH:THF:buffer 50:7:43 (Buffer was prepared from equal volumes of 22 g/L KH₂PO₄ and 18 g/L Na₂HPO₄, pH 6.6 ± 0.2.)
Flow rate: 2
Injection volume: 6
Detector: UV 240

CHROMATOGRAM

Retention time: 5.2

OTHER SUBSTANCES

Simultaneous: acetaminophen, acetazolamide, amobarbital, aspirin, barbital, butabarbital, ce-fazolin, chlorothiazide, chlorzoxazone, dicloxacillin, ethosuximide, furosemide, hydrochlorothia-zide, ibuprofen, oxacillin, pentobarbital, phenobarbital, phenylbutazone, phenytoin, salicylic acid, secobarbital, sulfamethoxazole, theophylline, thiopental, ascorbic acid
Noninterfering: ampicillin, penicillin G, valproic acid

REFERENCE

Kelner,M.; Bailey,D.N. Reversed-phase liquid-chromatographic simultaneous analysis for thiopental and pen-tobarbital in serum, *Clin.Chem.*, **1983**, 29, 1097–1100.

SAMPLE

Matrix: solutions
Sample preparation: Prepare a solution in MeCN:water 25:75, inject a 100 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 10 µm C18 (Alltech)
Mobile phase: MeCN:4 mM potassium phosphate buffer 50:50, pH 4.0
Flow rate: 1.2
Injection volume: 100
Detector: UV 290

CHROMATOGRAM

Retention time: 6.9

OTHER SUBSTANCES

Simultaneous: amobarbital, heptabarbital, hexobarbital, methohexital, pentobarbital, pheno-barbital, secobarbital, thiopental

Flow rate: 0.9
Injection volume: 50
Detector: UV 288

CHROMATOGRAM

Retention time: 37 (S(-)), 39 (R(+))
Internal standard: n-propyl p-hydroxybenzoate (25)
Limit of detection: 50 ng/mL

OTHER SUBSTANCES

Noninterfering: muscle relaxants, benzodiazepines

KEY WORDS

serum; chiral; SPE

REFERENCE

Sueyasu,M.; Ikeda,T.; Otsubo,K.; Taniyama,T.; Aoyama,T.; Oishi,R. Enantioselective determination of thiamylal in human serum by high-performance liquid chromatography, *J.Chromatogr.B*, **1995**, 665, 133–137.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 125 × 4.6 5 µm C-18 (Perkin-Elmer)
Mobile phase: MeOH:THF:buffer 50:7:43 (Buffer was prepared from equal volumes of 22 g/L KH₂PO₄ and 18 g/L Na₂HPO₄, pH 6.6 ± 0.2.)
Flow rate: 2
Injection volume: 6
Detector: UV 240

CHROMATOGRAM

Retention time: 5.2

OTHER SUBSTANCES

Simultaneous: acetaminophen, acetazolamide, amobarbital, aspirin, barbital, butabarbital, ce-fazolin, chlorothiazide, chlorzoxazone, dicloxacillin, ethosuximide, furosemide, hydrochlorothia-zide, ibuprofen, oxacillin, pentobarbital, phenobarbital, phenylbutazone, phenytoin, salicylic acid, secobarbital, sulfamethoxazole, theophylline, thiopental, ascorbic acid
Noninterfering: ampicillin, penicillin G, valproic acid

REFERENCE

Kelner,M.; Bailey,D.N. Reversed-phase liquid-chromatographic simultaneous analysis for thiopental and pen-tobarbital in serum, *Clin.Chem.*, **1983**, 29, 1097–1100.

SAMPLE

Matrix: solutions
Sample preparation: Prepare a solution in MeCN:water 25:75, inject a 100 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 10 µm C18 (Alltech)
Mobile phase: MeCN:4 mM potassium phosphate buffer 50:50, pH 4.0
Flow rate: 1.2
Injection volume: 100
Detector: UV 290

CHROMATOGRAM

Retention time: 6.9

OTHER SUBSTANCES

Simultaneous: amobarbital, heptabarbital, hexobarbital, methohexital, pentobarbital, pheno-barbital, secobarbital, thiopental

mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methyprylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripelennamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, **1994**, *18*, 233–242.


Thiethylperazine

Molecular formula: $C_{22}H_{29}N_3S_2$

Molecular weight: 399.62

CAS Registry No.: 1420-55-9, 52239-63-1 (maleate)

Merck Index: 9449



SAMPLE

Matrix: blood

Sample preparation: 3 mL Serum + 3 mL diluted Titrisol (pH 10 borate buffer, Merck) + 4 mL heptane:isoamyl alcohol 97:3, shake thoroughly for 15 s, centrifuge at 2500 g for 5 min. Remove the organic phase and add it to 1.5 mL 50 mM sulfuric acid containing 0.1% $Na_2S_2O_5$, mix for 15 s, centrifuge at 2500 g for 10 min. Remove the aqueous phase. Repeat the extraction and back extraction. Combine the aqueous phases and add them to 1.5 mL 2 M pH 9.1 glycine buffer, add 200 μ L n-hexane:isoamyl alcohol 97:3, vortex for 25 s. Remove the organic phase and evaporate it to dryness in a desiccator, reconstitute in mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 300 \times 3.9 10 μ m μ Bondapak C18

Mobile phase: MeCN:water:acetic acid 65:35:3 containing 10 mM dodecyl hydrogen sulfate

Flow rate: 2

Detector: UV 257

CHROMATOGRAM

Retention time: 13

Internal standard: thiethylperazine

OTHER SUBSTANCES

Extracted: perphenazine

Noninterfering: alimemazine, biperidine, carbamazepine, chlorpromazine, clomipramine, diazepam, dihydroergotamine, disulfiram, dixyrazine, haloperidol, levomepromazine, nitrazepam, orphenadrine, promethazine, propiomazine, thioridazine, trimipramine, vitamins

KEY WORDS

serum; thiethylperazine is IS

REFERENCE

Larsson,M.; Forsman,A. A high-performance liquid chromatographic method for the assay of perphenazine and its dealkylated metabolite in serum after therapeutic doses, *Ther.Drug Monit.*, **1983**, 5, 225–228.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a solution in mobile phase, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 100 × 3.5 µm Lichrosorb SI60

Mobile phase: MeCN:MeOH:ammonium hydroxide 250:55:13

Flow rate: 1.2

Injection volume: 20

Detector: UV 240

CHROMATOGRAM

Retention time: 3.4

OTHER SUBSTANCES

Simultaneous: N-acetylprocainamide, amitriptyline, chlorimipramine, chlorpromazine, codeine, desipramine, dimethacrine, diphenhydramine, disopyramide, doxepin, hydroquinidine, maprotiline, melitracene, mesoridazine, nortriptyline, opipramol, perazine, perphenazine, procainamide, prochlorperazine, promazine, prothipendyl, protriptyline, thioperazine, thioridazine, trifluoperazine

Noninterfering: acenocoumaron, acetaminophen, acetophenetidine, aspirin, benzodiazepines, bibenzepin, butriptyline, caffeine, chlorprothixene, clopenthixol, clothiapine, dixyrazine, droperidol, fluphenazine, haloperidol, hydroxyzine, isoniazid, methotrimeprazine, metopimazine, moperone, noxiptyline, orphenadrine, pericyazine, phenprocoumon, pipothiazine, promethazine, salicylic acid, theophylline, thiopropazate, trimeprazine, trimipramine

Interfering: butaperazine, imipramine, pipamperone, quinidine, thiothixene

REFERENCE

Edelbroek,P.M.; de Haas,E.J.M.; de Wolff,F.A. Liquid-chromatographic determination of amitriptyline and its metabolites in serum, with adsorption onto glass minimized, *Clin.Chem.*, **1982**, 28, 2143–2148.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 125 × 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 4.3

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cy-

clizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipranone, diprenorphine, dipyridamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, flupromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, proprietyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocainide, tolpropamine, tolycaine, tranlylcypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, L.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, *323*, 191-225.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.5 µm LiChrospher 100 RP-8

Mobile phase: MeCN:0.025% phosphoric acid:buffer 60:25:15 (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 9.45

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordiazepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyridamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenopropfen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxy-

chloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazin-
dol, mefenamic acid, meperidine, mephenytoin, mepivacaine, mesoridazine, metaproterenol,
metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine,
methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, met-
ronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline,
naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymet-
azoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, phen-
iramine, phenobarbital, phenol, phenolphthalein, phenolamine, phenylbutazone, phenyltolox-
amine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid,
procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone,
propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quin-
ine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, seco-
barbital, sertraline, sotalol, spironolactone, sulfapyrazone, sulindac, temazepam, terbutaline,
terfenadine, tetracaine, theophylline, thiopental, thioridazine, thiothixene, timolol, tocamide,
tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, tri-
meprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine,
zopiclone

KEY WORDS

details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, 692, 103–119.


Thimerosal

Molecular formula: C₆H₆HgNaO₂S

Molecular weight: 404.82

CAS Registry No.: 54-64-8

Merck Index: 9451



SAMPLE

Matrix: formulations

Sample preparation: Add IS to vaccine at a concentration of 40 $\mu\text{g/mL}$, centrifuge at 3400 g for 15 min, inject a 25 μL aliquot of the supernatant. (IS stock solution was prepared in mobile phase:water 1:4.)

HPLC VARIABLES

Guard column: 5 × 45 μm Hypersil C18

Column: 210 × 4.6 5 μm Hypersil C18

Mobile phase: MeOH:water:orthophosphoric acid 35:35:0.9, pH 2.5

Flow rate: 0.6

Injection volume: 25

Detector: UV 222

CHROMATOGRAM

Retention time: 10.1

Internal standard: salicylic acid (6.5)

OTHER SUBSTANCES

Simultaneous: degradation products

KEY WORDS

vaccine

REFERENCE

Tleugabulova,D.; Gonzalez Perez,I. Reversed-phase high-performance liquid chromatographic study of thimerosal stability in Cuban recombinant hepatitis B vaccine, *J.Chromatogr.A*, **1996**, 729, 219–227.

SAMPLE

Matrix: formulations

HPLC VARIABLES

Column: 300 × 3.9 10 µm LiChrosorb Si-60

Mobile phase: MeOH:water 60:40 containing 4 mM disodium citrate and 4 mM tetrabutylammonium bromide, pH 6.0

Flow rate: 1

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 6

OTHER SUBSTANCES

Simultaneous: benzalkonium chloride, domiphen bromide, xylometazoline

Interfering: inorganic salts

KEY WORDS

nasal drops

REFERENCE

Lingeman,H.; van Munster,H.A.; Beynen,J.H.; Underberg,W.J.; Hulshoff,A. High-performance liquid chromatographic analysis of basic compounds on non-modified silica gel and aluminium oxide with aqueous solvent mixtures, *J.Chromatogr.*, **1986**, 352, 261–274.

SAMPLE

Matrix: formulations

Sample preparation: 9.5 mL Contact lens solution + 0.5 mL 3 mg/mL methylparaben, inject a 20 µL aliquot.

HPLC VARIABLES

Guard column: 7 µm Nucleosil C18 pre-column

Column: 7 µm Nucleosil C18

Mobile phase: MeOH:100 mM KH₂PO₄ adjusted to pH 3.5 with phosphoric acid 60:40

Flow rate: 1

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 11.4

Internal standard: methyl paraben (7.8)

Limit of detection: 100 ng

OTHER SUBSTANCES

Simultaneous: chlorhexidine gluconate, thiosalicylic acid

KEY WORDS

stability-indicating; contact lens solutions

REFERENCE

Hu,O.Y.-P.; Wang,S.-Y.; Fang,Y.-J.; Chen,Y.-H.; King,M.-L. Simultaneous determination of thimerosal and chlorhexidine in solutions for soft contact lenses and its applications in stability studies, *J.Chromatogr.*, **1990**, 523, 321–326.

SAMPLE**Matrix:** formulations**Sample preparation:** Inject directly or extract as follows. Condition a Sep-Pak C18 SPE cartridge with MeOH and water. 10 mL Ophthalmic solution + 100 μ L concentrated phosphoric acid, add to the SPE cartridge, wash with 5 mL water, wash with 2 mL MeOH:water 20:80, elute with 2 mL MeOH, dilute the eluate to 10 mL with MeOH, inject an aliquot.

HPLC VARIABLES**Column:** 150 \times 4 5 μ m Spherisorb C18**Mobile phase:** MeOH:water 50:50 containing 2 mM tetraethylammonium perchlorate, adjusted to pH 4.8 with perchloric acid**Flow rate:** 1**Injection volume:** 20**Detector:** E, Metrohm Model 461, Metrohm Model 656 flow cell, carbon paste electrode, 0.9 V

CHROMATOGRAM**Retention time:** 5**Limit of detection:** 90 ng/mL

OTHER SUBSTANCES**Simultaneous:** degradation products, 2,2'-dithiodibenzoic acid, thiosalicylic acid

KEY WORDSophthalmic solutions; SPE

REFERENCE

del Pilar da Silva,M.; Procopio,J.R.; Hernández,L. Evaluation of the capability of different chromatographic systems for the monitoring of thimerosal and its degradation products by high-performance liquid chromatography with amperometric detection, *J.Chromatogr.A*, **1993**, 653, 267–273.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 10 μ m Spherisorb 10 ODS**Mobile phase:** MeOH:water:phosphoric acid 60:50:1**Flow rate:** 2.6**Injection volume:** 25**Detector:** UV 222

CHROMATOGRAM**Retention time:** 3**Limit of detection:** 100 ng/mL

OTHER SUBSTANCES**Simultaneous:** degradation products, 2,2'-dithiosalicylic acid, thiosalicylic acid

REFERENCE

Reader,M.J.; Lines,C.B. Decomposition of thimerosal in aqueous solution and its determination by high-performance liquid chromatography, *J.Pharm.Sci.*, **1983**, 72, 1406–1409.

SAMPLE**Matrix:** solutions**Sample preparation:** Inject an aliquot of the ophthalmic solution directly.

HPLC VARIABLES**Column:** two 250 \times 4 10 μ m Spherisorb ODS 10 columns in series**Mobile phase:** MeOH:water:phosphoric acid 60:50:1**Flow rate:** 4**Injection volume:** 50

Detector: UV 222

KEY WORDS

ophthalmic solutions

REFERENCE

Reader,M.J. Influence of isotonic agents on the stability of thimerosal in ophthalmic formulations, *J.Pharm.Sci.*, 1984, 73, 840–841.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 210 × 4.6 5 µm Spherisorb RP-18

Mobile phase: MeOH:water:phosphoric acid 65:35:0.9

Flow rate: 0.6

Injection volume: 20

Detector: UV 222

CHROMATOGRAM

Retention time: 7.6

Limit of quantitation: 5 ppm

OTHER SUBSTANCES

Simultaneous: degradation products, dithiosalicylic acid, thiosalicylic acid

REFERENCE

Caraballo,I.; Rabasco,A.M.; Fernández-Arévalo,M. Study of thimerosal degradation mechanism, *Int.J.Pharm.*, 1993, 89, 213–221.

SAMPLE

Matrix: urine

Sample preparation: Add 0.05% thymol to urine, filter, dilute 1:10 with water, inject an aliquot.

HPLC VARIABLES

Guard column: Guard-Pak C18

Column: 300 × 3.9 µBondapak C18

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 33

OTHER SUBSTANCES

Simultaneous: creatinine, thymol, uric acid

REFERENCE

Chen,Y.; Pietrzyk,R.A.; Whitson,P.A. Quantification of urinary uric acid in the presence of thymol and thimerosal by high-performance liquid chromatography, *J.Chromatogr.A*, 1997, 763, 187–192.


Thiobutabarbital

Molecular formula: C₁₀H₁₆N₂O₂S, C₁₀H₁₅N₂NaO₂S

Molecular weight: 228.32


CAS Registry No.: 2095-57-0, 947-08-0 (sodium salt)

Merck Index: 9457



SAMPLE**Matrix:** solutions**HPLC VARIABLES****Guard column:** 4 × 4 5 µm LiChroCART LiChrospher 60 RP Select B**Column:** 125 × 4 5 µm LiChroCART LiChrospher 60 RP Select B**Mobile phase:** MeCN:buffer 50:50 (Buffer was 25 mM pH 3.0 triethylammonium phosphate containing 2% MeCN.)**Flow rate:** 2**Injection volume:** 50**Detector:** UV 283**CHROMATOGRAM****Retention time:** 2.70**OTHER SUBSTANCES****Simultaneous:** thiopental**REFERENCE**Hannak,D.; Scharbert,F.; Kattermann,R. Stepwise binary gradient high-performance liquid chromatographic system for routine drug monitoring, *J.Chromatogr.A*, **1996**, 728, 307–310.

Thiocctic acid

**Molecular formula:** C₈H₁₄O₂S₂**Molecular weight:** 206.33**CAS Registry No.:** 62-46-4, 1200-22-2 (d-form), 1077-28-7 (dl-form), 1077-27-6 (l-form), 2319-84-8 (sodium salt)**Merck Index:** 9462**SAMPLE****Matrix:** solutions**Sample preparation:** Mix 200 µL of a solution in MeOH:water 10:90 with 150 µL 20 mM tetraethylammonium bromide in 100 mM pH 7.0 phosphate buffer and 100 µL 4.2 mg/mL 2-bromoacetyl-6-methoxynaphthalene in acetone, sonicate for 3 min, let stand for 22 min, add 150 µL 4.5 µg/mL IS in MeCN, sonicate at room temperature for 1 min, inject a 50 µL aliquot. (Prepare 2-bromoacetyl-6-methoxynaphthalene by stirring equimolar amounts of 2-acetyl-6-methoxynaphthalene (Janssen Chimica, Belgium) and phenyltrimethylammonium tribromide in THF at room temperature for 3 h (Phosphorus and Sulfur 1985, 25, 357), purify by column chromatography on silica gel with chloroform:petroleum ether 50:50 (mp 109–112°) (Chromatographia 1992, 33, 13).)**HPLC VARIABLES****Column:** 250 × 4.5 Hypersil 5 ODS**Mobile phase:** MeCN:MeOH:water 37.4:30.6:32**Column temperature:** 35**Flow rate:** 1.1**Injection volume:** 50**Detector:** F ex 300 em 460**CHROMATOGRAM****Retention time:** 16.5**Internal standard:** n-hexanoic acid 6-methoxynaphthacyl ester (?) (Dissolve 2 mmole n-hexanoic acid and 1 mmole 2-bromoacetyl-6-methoxynaphthalene in 10 mL anhydrous MeCN, add 500 µL triethylamine, heat to 60° for 30 min, cool, add 30 mL water, extract three times with 10 mL portions of diethyl ether. Combine the organic layers and wash them with 5% sodium bicarbonate solution, wash three times with 10 mL portions of water, dry over anhydrous sodium sulfate, evaporate to dryness under reduced pressure, recrystallize from MeOH/water

to give 6-methoxynaphthacyl ester of n-hexanoic acid (mp 79-80°) (J. Pharm. Biomed. Anal. 1993, 11, 761) (19)


KEY WORDS

derivatization

REFERENCE

Gatti,R.; Bousquet,E.; Bonazzi,D.; Cavrini,V. Determination of carboxylic acid salts in pharmaceuticals by high-performance liquid chromatography after pre-column fluorogenic labelling, *Biomed.Chromatogr.*, **1996**, 10, 19-24.

Thioguanine



Molecular formula: C₅H₅N₅S

Molecular weight: 167.19

CAS Registry No.: 154-42-7, 5580-03-0 (hemihydrate)

Merck Index: 9473

Lednicer No.: 2 464

SAMPLE

Matrix: blood

Sample preparation: Plasma. 1 mL Plasma + 1 mL cold 2 M perchloric acid, mix, centrifuge at 4° at 48000 g for 20 min. Remove a 1 mL aliquot of the supernatant and adjust to pH 10-12 with 150 µL 4 M KOH, let stand at 4° for 2 days. Adjust the pH of the supernatant to 2-3 with 150 µL 1 M HCl, centrifuge at 700 g for 10 min, inject an aliquot of the supernatant.

HPLC VARIABLES

Column: 300 × 3.9 10 µm µBondapak C18

Mobile phase: MeOH:buffer 10:90 (Buffer was 10 mM acetic acid adjusted to pH 3.5 with 5 M NaOH. Flush column daily with MeOH:water 50:50.)

Flow rate: 2

Injection volume: 100

Detector: UV 340

CHROMATOGRAM

Retention time: 4

Limit of detection: 800 nM

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Andrews,P.A.; Egorin,M.J.; May,M.E.; Bachur,N.R. Reversed-phase high-performance liquid chromatography analysis of 6-thioguanine applicable to pharmacologic studies in humans, *J.Chromatogr.*, **1982**, 227, 83-91.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 100 µL 400 mM NaOH + 1 mL 0.3% phenylmercuric acetate in ethyl acetate + 3 mL diethyl ether, shake on a tumble mixer for 10 min, centrifuge for 5 min. Remove the organic layer and add it to 500 µL 100 mM HCl, whirlmix for 2 min, centrifuge for 5 min, discard the organic layer, evaporate traces of organic solvent under a stream of nitrogen at room temperature for 15 min, add 10 µL 3 mg/mL dithioerythritol in water, inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.6 LiChrosorb 10 RP-18

Mobile phase: Isopropanol:water 3:97 containing 13.80 g/L NaH₂PO₄·H₂O, 200 µL/L 85% phosphoric acid, 60 mg/L dithioerythritol, and 500 mg/L sodium octanesulfonate, pH 3.6-3.7

Flow rate: 1.5

Detector: F ex 295 em 380 following post-column reaction. The column effluent mixed with 8 mM potassium chromate in 500 mM HCl pumped at 0.16 mL/min and with air flowing at 0.32 mL/min and the mixture flowed through a single mixing coil. The effluent from this coil mixed with 1.6% sodium metabisulfite pumped at 0.16 mL/min and this mixture flowed through a single mixing coil. The effluent from this coil mixed with 4 M ammonium hydroxide pumped at 0.23 mL/min and this mixture flowed through a double mixing coil to a debubbler. The liquid effluent from the debubbler flowed to the detector.

CHROMATOGRAM

Retention time: 8

Internal standard: 6-thioguanine

OTHER SUBSTANCES

Extracted: 6-mercaptopurine

KEY WORDS

post-column reaction; plasma; 6-thioguanine is IS

REFERENCE

Jonkers,R.E.; Oosterhuis,B.; ten Berge,R.J.M.; van Boxtel,C.J. Analysis of 6-mercaptopurine in human plasma with a high-performance liquid chromatographic method including post-column derivatization and fluorimetric detection, *J.Chromatogr.*, **1982**, 233, 249-255.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 80 µL water + 10 µL 1 M dithiothreitol, vortex for 10 s, add 2 mL MeCN, vortex for 30 s, centrifuge at 2000 g for 5 min. Remove the supernatant and add it to 2 mL dichloromethane, shake on a reciprocating shaker for 5 min, centrifuge at 2000 g for 5 min. Remove 750 µL from the top aqueous layer and evaporate it to dryness under a stream of nitrogen at 37°, reconstitute the residue in 150 µL distilled water, vortex for 1 min, inject a 15 µL aliquot.

HPLC VARIABLES

Guard column: 70 × 2.2 30-38 µm Co:Pell ODS

Column: 250 × 4.6 5 µm Ultrasphere ODS

Mobile phase: MeCN:acetic acid:water 3.5:0.2:96.3

Flow rate: 1.4

Injection volume: 15

Detector: UV 322

CHROMATOGRAM

Retention time: 6.5

Internal standard: 6-thioguanine

OTHER SUBSTANCES

Extracted: 6-mercaptopurine

Noninterfering: caffeine, cytarabine, 5-fluorouracil, prednisone, theophylline, vinblastine, vincristine

KEY WORDS

plasma; protect from light; monkey; human; 6-thioguanine is IS

REFERENCE

Narang,P.K.; Yeager,R.L.; Chatterji,D.C. Quantitation of 6-mercaptopurine in biologic fluids using high-performance liquid chromatography: a selective and novel procedure, *J.Chromatogr.*, **1982**, 230, 373-380.

SAMPLE**Matrix:** blood**Sample preparation:** Dilute 10 mL blood with 15 mL phosphate buffered saline containing 3 mM ethyleneglycoltetraacetic acid. Layer 8 mL diluted blood on 2 mL Ficoll-Hypaque (Sigma), centrifuge at 9000 g for 10 min, collect lymphocyte band, wash with PBS. 200 μ L Lymphocytes + 50 μ L 1 M sulfuric acid, heat at 100° for 45 min, cool, centrifuge. Remove 200 μ L and add to 55 μ L 1 M pH 10.1 sodium bicarbonate, add 1.5 mL ethyl acetate:dichloromethane containing 100 ng/mL sulfamethoxazole, vortex, centrifuge. Remove 175 μ L of the upper aqueous layer and add it to 25 μ L 1 M pH 10.1 sodium bicarbonate, add 25 μ L 0.5% potassium permanganate (freshly prepared), let stand for 5 min, add 5 μ L 15% hydrogen peroxide, vortex, inject a 10 μ L aliquot.

HPLC VARIABLES**Column:** 200 \times 4.6 5 μ m Hypersil ODS**Mobile phase:** MeCN:buffer 11:89 (Buffer was 10 mM pH 7 sodium phosphate containing 0.06% tetrabutylammonium chloride.)**Column temperature:** 45**Flow rate:** 1.5**Injection volume:** 10**Detector:** F ex 330 em 412 (370 nm cut-off filter)

CHROMATOGRAM**Retention time:** 2.30**Limit of detection:** 1 ng

KEY WORDS

whole blood; lymphocytes

REFERENCEErdmann,G.R.; Steury,J.C.; Carleton,B.C.; Stafford,R.J.; Bostrom,B.C.; Canafax,D.M. Reversed-phase high-performance liquid chromatographic approach to determine total lymphocyte concentrations of 6-thioguanine, methylmercaptapurine and methylthioguanine in humans, *J.Chromatogr.*, **1991**, 571, 149-156.

SAMPLE**Matrix:** blood**Sample preparation:** 250 μ L Plasma + 25 μ L 1 M aluminum perchlorate in water, let stand at room temperature for 15 min, chill in ice water for 15 min, centrifuge at 15600 g for 15 min, discard the supernatant. Suspend the precipitate in 500 μ L 50 mM aluminum perchlorate in water by stirring to break up the precipitate, vortex for 20 s, centrifuge at 15600 g for 15 min, discard the supernatant. Suspend the precipitate in 150 μ L 400 mM perchloric acid, add 5 μ L freshly prepared 200 mM aqueous sodium hydrosulfite, mix, let stand at room temperature for 30 min, chill in ice water, centrifuge at 15600 g for 15 min, inject a 50 μ L aliquot of the supernatant.

HPLC VARIABLES**Guard column:** 45 \times 2 40 μ m ODS**Column:** 150 \times 4.3 5 μ m Ultrasphere ODS**Mobile phase:** Water:85% phosphoric acid 99.32:0.68 containing 154.3 mg/L dithiothreitol**Flow rate:** 1**Injection volume:** 50**Detector:** UV 340

CHROMATOGRAM**Retention time:** 4.73**Internal standard:** 6-thioguanine**Limit of quantitation:** 3 ng/mL

OTHER SUBSTANCES**Extracted:** 6-mercaptopurine, 6-thiouric acid

KEY WORDS

plasma; pharmacokinetics; 6-thioguanine is IS

REFERENCE

Lin, K.T.; Varin, F.; Rivard, G.E.; Leclerc, J.M. Isolation of 6-mercaptopurine in human plasma by aluminum ion complexation for high-performance liquid chromatographic analysis, *J. Chromatogr.*, **1991**, 536, 349–355.

SAMPLE

Matrix: blood

Sample preparation: Add 1 volume ice-cold 8 M perchloric acid to 20 volumes plasma, mix, keep on ice for 10 min, centrifuge at 10 000 g for 15 min, remove the supernatant. Adjust the pH of the supernatant to 6–7 with 10 volumes ice-cold 4 M K_2KHPO_4 , keep on ice for 10 min, centrifuge at 10 000 g for 5 min, inject a 100 μ L aliquot of the supernatant.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Supelcosil LC-18-DB

Mobile phase: Gradient. A was 25 mM KH_2PO_4 . B was MeOH:50 mM KH_2PO_4 25:75. A:B from 98:2 to 96:4 over 5 min, to 85:15 over 3 min, to 80:20 over 2 min, to 40:60 over 10 min, maintain at 40:60 over 2 min, to 20:80 over 3 min, maintain at 20:80 for 20 min, return to initial conditions over 3 min, re-equilibrate for 12 min.

Flow rate: 1.25

Injection volume: 100

Detector: UV 342

CHROMATOGRAM

Retention time: 10

Limit of detection: 20–50 nM

OTHER SUBSTANCES

Extracted: metabolites, mercaptopurine (UV 320)

KEY WORDS

plasma

REFERENCE

Keuzenkamp-Jansen, S.W.; De Abreu, R.A.; Böklerink, J.P.M.; Trijbels, J.M.F. Determination of extracellular and intracellular thiopurines and methylthiopurines by high-performance liquid chromatography, *J. Chromatogr. B*, **1995**, 672, 53–61.

SAMPLE

Matrix: blood, CSF, urine

Sample preparation: Blood, CSF. Collect 2 mL blood in tubes containing heparin and 120 μ g dithiothreitol (DTT). Mix, centrifuge at 2000 g for 5 min, remove plasma. CSF. Collect 0.5 mL CSF in tubes containing 30 μ g DTT. Cool CSF and plasma samples on ice, add freshly prepared ice-cold 50% trichloroacetic acid equal to 10% of sample volume. Urine. Collect 2 mL urine in tubes containing 120 μ g DTT, filter (0.22 μ m), inject an aliquot.

HPLC VARIABLES

Column: two 250 \times 4.6 10 μ m Nucleosil 10 C18 columns in series

Mobile phase: Gradient. A was 25 mM pH 2.75 phosphoric acid. B was MeOH:water 50:50. C was 100 mM pH 6.6 KH_2PO_4 . A:B:C from 100:0:0 to 98:2:0 over 5 min, to 30:3.5:66.5 over 5 min, maintain at 30:3.5:66.5 for 10 min, re-equilibrate at initial conditions for 10 min.

Column temperature: 33

Flow rate: 1.7

Injection volume: 195, 500

Detector: UV 342

CHROMATOGRAM

Retention time: 13

Limit of detection: 25 nM

OTHER SUBSTANCES

Extracted: metabolites, mercaptopurine (UV 312), 6-mercaptopurine riboside, 6-thioguanosine

KEY WORDS

plasma; goat; pharmacokinetics

REFERENCE

van Baal, J.M.; van Leeuwen, M.B.; Schouten, T.J.; De Abreu, R.A. Sensitive high-performance liquid chromatographic determination of 6-mercaptopurine, 6-thioguanine, 6-mercaptopurine riboside and 6-thioguanosine in biological fluids, *J.Chromatogr.*, **1984**, 336, 422–428.

SAMPLE

Matrix: blood, urine

Sample preparation: Plasma. Filter (Amicon Model MM 302, type PM 10 Diaflo membrane) plasma while centrifuging at 8000 g for 3 min, inject a 20–100 μ L aliquot of the ultrafiltrate. Urine. Inject an aliquot directly.

HPLC VARIABLES

Column: 150 \times 4 Nucleosil 5C8

Mobile phase: 50 mM pH 7.0 citrate-phosphate buffer (Buffer was 9.61 g $\text{Na}_2\text{HPO}_4 \cdot \text{H}_2\text{O}$ and 1.33 g citric acid monohydrate in 1 L water, pH 7.0.)

Flow rate: 1.4

Injection volume: 20–100

Detector: UV 343

CHROMATOGRAM

Retention time: 5.7

Limit of quantitation: 200 ng/mL

OTHER SUBSTANCES

Extracted: 6-mercaptopurine

Simultaneous: allopurinol, azaguanine, guanine, oxipurinol, uric acid

Noninterfering: adenine, 2-amino-6-methylthiopurine, aspirin, benzbromarone, caffeine, diazepam, dihydralazine, dipyridamole, fluorouracil, hypoxanthine, methotrexate, procarbazine, propranolol, spironolactone, sulfamethoxazole, sulfapyrazone, theophylline, thiouric acid, thioxanthine, trimethoprim, xanthine

KEY WORDS

plasma; ultrafiltrate

REFERENCE

Breithaupt, H.; Goebel, G. Quantitative high pressure liquid chromatography of 6-thioguanine in biological fluids, *J.Chromatogr.Sci.*, **1981**, 19, 496–499.

Thiopental


Molecular formula: $\text{C}_{11}\text{H}_{18}\text{N}_2\text{O}_2\text{S}$

Molecular weight: 242.34

CAS Registry No.: 76-75-5, 71-73-8 (Na salt)

Merck Index: 9487

Lednicer No.: 1 274

**SAMPLE**

Matrix: blood

Sample preparation: Precipitate 100 μ L serum with 200 μ L 10 μ g/mL IS in MeCN, centrifuge at 12000 g for 5 min. Inject a 50 μ L aliquot of the supernatant.

HPLC VARIABLES

Guard column: 4 \times 4 5 μ m LiChroCART LiChrospher 60 RP Select B