

Tolazamide


Molecular formula: $C_{14}H_{21}N_3O_3S$

Molecular weight: 311.41

CAS Registry No.: 1156-19-0

Merck Index: 9644

Lednicer No.: 1 241



SAMPLE

Matrix: blood

Sample preparation: 500 μ L Serum + 500 μ L 5-(p-methylphenyl)-5-phenylhydantoin in chloroform + 500 μ L pH 4.5 sodium acetate buffer + 5 mL dichloromethane, shake horizontally for 5 min, centrifuge at 3000 g for 3 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at room temperature, reconstitute the residue in 80 μ L MeOH, vortex, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 10 μ m Lichrosorb C-18

Mobile phase: MeOH:pH 5.6 acetate buffer 52.3:47.7

Flow rate: 2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 5.5

Internal standard: 5-(p-methylphenyl)-5-phenylhydantoin (7.5)

Limit of quantitation: 1000 ng/mL

KEY WORDS

serum; pharmacokinetics

REFERENCE

Welling,P.G.; Patel,R.B.; Patel,U.R.; Gillespie,W.R.; Craig,W.A.; Albert,K.S. Bioavailability of tolazamide from tablets: comparison of in vitro and in vivo results, *J.Pharm.Sci.*, **1982**, *71*, 1259–1263.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Serum + 2 mL water + 200 μ L 1 (?) M HCl + 200 μ L 2.5 μ g/mL glibornuride in MeOH + 7 mL diethyl ether, mix, centrifuge at 2000 rpm for 5 min. Remove 6.5 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 500 μ L 2 mg/mL dinitrofluorobenzene in butyl acetate, heat at 120° for 1 h, cool, evaporate to dryness under a stream of nitrogen, reconstitute the residue in 150 μ L mobile phase, inject a 120 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Spherisorb ODS 2

Mobile phase: MeCN:0.4% aqueous phosphoric acid 75:25

Column temperature: 40

Flow rate: 1.2

Injection volume: 120

Detector: UV 360

CHROMATOGRAM

Retention time: 7.3

Internal standard: glibornuride (5.8)

Limit of detection: 40 ng/mL

OTHER SUBSTANCES

Extracted: chlorpropamide, glipizide, glyburide (glibenclamide), tolbutamide

KEY WORDS

serum; derivatization

REFERENCE

Starkey,B.J.; Mould,G.P.; Teale,J.D. The determination of sulphonylurea drugs by HPLC and its clinical application, *J.Liq.Chromatogr.*, **1989**, 12, 1889-1896.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma + 200 μ L 2 M HCl + 2 mL diethyl ether, vortex for 30 s, centrifuge at 1500 g for 5 min, freeze in dry ice for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 35-40°, reconstitute the residue in 100 μ L mobile phase, vortex for 30 s, inject a 25-50 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.1 10 μ m Versapak C18 (Alltech)

Mobile phase: MeCN:10 mM orthophosphoric acid 50:50

Flow rate: 1

Injection volume: 25-50

Detector: UV 230

CHROMATOGRAM

Retention time: 8.2

OTHER SUBSTANCES

Extracted: chlorpropamide, gliclazide, glyburide (glibenclamide), glipizide, tolbutamide

Simultaneous: sulfamethoxazole, N-acetylsulfamethoxazole

Noninterfering: trimethoprim

KEY WORDS

plasma

REFERENCE

Shenfield,G.M.; Boutagy,J.S.; Webb,C. A screening test for detecting sulfonylureas in plasma, *Ther.Drug Monit.*, **1990**, 12, 393-397.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitriptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepan, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, dil-

CHROMATOGRAM**Retention time:** 3.7**Internal standard:** clonidine**Limit of detection:** 50 ng/mL

OTHER SUBSTANCES**Noninterfering:** amikacin, ampicillin, dopamine, furosemide, gentamicin, heparin, morphine, pancuronium

KEY WORDSserum

REFERENCETodesco, L.M.; Thoma, J.J.; Barth, R.D.; Myers, N.J.; White, R.; Ward, R.M. Quantitative determination of tolazoline in human serum by high performance liquid chromatography, *Ther. Drug Monit.*, **1987**, 9, 78–84.

SAMPLE**Matrix:** blood, urine**Sample preparation:** 50–100 μ L Serum or urine + 30 μ L 1 (serum) or 100 (urine) μ g/mL naphazoline + 250 μ L buffer + 5 mL dichloromethane, shake mechanically for 20 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 35°, reconstitute the residue in 100 (serum) or 500 (urine) μ L 20 mM pH 3.75 KH_2PO_4 , vortex for 15 s, inject a 10 (urine) or 90 (serum) μ L aliquot. (Buffer was 100 mM potassium hydrogen carbonate and 100 mM potassium carbonate, pH 10.0.)

HPLC VARIABLES**Guard column:** 50 \times 2 Co:Pell ODS**Column:** Resolve C18 (Waters)**Mobile phase:** MeCN:buffer 40:60 (Buffer was 20 mM KH_2PO_4 adjusted to pH 3.75 with 85% phosphoric acid.)**Flow rate:** 1.2**Injection volume:** 10–90**Detector:** UV 210

CHROMATOGRAM**Retention time:** 7**Internal standard:** naphazoline (10)**Limit of detection:** 100 ng/mL

OTHER SUBSTANCES**Noninterfering:** acetaminophen, N-acetylprocainamide, carbamazepine, chloramphenicol, desipramine, digoxin, disopyramide, dopamine, ethosuximide, gentamicin, imipramine, lidocaine, methotrexate, phenobarbital, phenytoin, primidone, procainamide, propranolol, quinidine, salicylic acid, theophylline, valproic acid

KEY WORDSserum

REFERENCECwik, M.J.; Chiu, G.P.; Fischer, J.H.; Chow-Tung, E.; Currie, B.L. Quantitative determination of tolazoline in serum and urine, *J. Chromatogr.*, **1985**, 338, 123–130.

SAMPLE**Matrix:** formulations**Sample preparation:** Mix 5 mL nasal solution and 10 mL 500 μ g/mL tolazoline hydrochloride in MeOH:water 40:60, make up to 50 mL with MeOH:water 40:60, inject a 10 μ L aliquot.

HPLC VARIABLES**Column:** 150 \times 4.1 RSIL C18 (RSL, Eke, Belgium)**Mobile phase:** MeOH:water 40:60 containing 20 mM sodium 1-octanesulfonate and 10 mM N,N-dimethyloctylamine, pH adjusted to 3.0 with orthophosphoric acid

Column temperature: 25

Flow rate: 1

Injection volume: 10

Detector: UV 220

CHROMATOGRAM

Retention time: 2.5

Internal standard: tolazoline

OTHER SUBSTANCES

Simultaneous: degradation products, antazoline, coumazoline, lidocaine, naphazoline, oxy-metazoline, prednisolone, sulfadimidine, sulfanilamide, sulfathiazole, tenaphthoxaline, tetrahydrozoline, tramazoline, xylometazoline

KEY WORDS

nasal solutions; stability-indicating; tolazoline is IS

REFERENCE

De Schutter, J.A.; Van den Bossche, W.; De Moerloose, P. Stability-indicating analysis of tetryzoline hydrochloride in pharmaceutical formulations by reversed-phase ion-pair liquid chromatography, *J.Chromatogr.*, **1987**, 391, 303-308.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, am-triptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepan, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, difunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenoprofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, iminostilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isoxsuprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methyprylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic

acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phenacyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyrithyldione, quazepam, quinaldic acid, quinidine, quinine, ranitidine, rescinamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.

Tolbutamide


Molecular formula: $C_{12}H_{18}N_2O_3S$

Molecular weight: 270.35

CAS Registry No.: 64-77-7, 473-41-6 (Na salt)

Merck Index: 9646

Lednicer No.: 1 136



SAMPLE

Matrix: blood

Sample preparation: 2 mL Plasma + 60 ng chlorpropamide + 500 μ L 1 M HCl + 8 mL chloroform, shake on a reciprocal shaker, shake for 10 min in a reciprocal shaker, centrifuge at 2000 g for 15 min. Remove 7 mL of the lower organic layer and evaporate it to dryness under a stream of nitrogen at 45°, reconstitute the residue in 100 μ L 3 mg/mL dinitrofluorobenzene in n-butyl acetate, heat at 120° for 30 min, evaporate to dryness under a stream of nitrogen at 60°, dissolve the residue in 100 μ L mobile phase, inject a 30-70 μ L aliquot. (Recrystallize dinitrofluorobenzene from diethyl ether. Prepare solutions weekly, store at 4° in the dark.)

HPLC VARIABLES

Column: 125 \times 4.6 5 μ m C8 (Perkin-Elmer)

Mobile phase: MeCN:water 50:50 containing 0.15% phosphoric acid

Flow rate: 1.5

Injection volume: 30-70

Detector: UV 350

CHROMATOGRAM

Retention time: 4.5

Internal standard: chlorpropamide (6.2)

Limit of detection: 2 ng/mL

OTHER SUBSTANCES

Extracted: glyburide

Noninterfering: acetaminophen, aspirin, diazepam, chlordiazepoxide, quinidine, phenytoin, theophylline, phenobarbital

KEY WORDS

plasma; derivatization

REFERENCE

Zecca,L.; Trivulzio,S.; Pinelli,A.; Colombo,R.; Tofanetti,O. Determination of glibenclamide, chlorpropamide and tolbutamide in plasma by high-performance liquid chromatography with ultraviolet detection, *J.Chromatogr.*, **1985**, 339, 203–209.

SAMPLE

Matrix: blood

Sample preparation: 1 mL PLasma or urine + 10 μ L MeOH + 10 μ L concentrated HCl + 8 mL toluene, extract for 15 min, centrifuge at 1500 g for 10 min, repeat extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at room temperature, reconstitute the residue in 200 μ L mobile phase, inject a 5 μ L aliquot.

HPLC VARIABLES

Column: 300 \times 4.6 10 μ m C18 (Merck)

Mobile phase: MeOH:50 mM pH 5.0 phosphate buffer 50:50

Flow rate: 1

Injection volume: 5

Detector: UV 230

CHROMATOGRAM

Retention time: 6

Internal standard: tolbutamide

OTHER SUBSTANCES

Extracted: nimesulide

Simultaneous: acetaminophen, aspirin, doxepin, glibenclamide, salicylic acid, theophylline

Noninterfering: digoxin, flurazepam, tiadenol

Interfering: bezafibrate

KEY WORDS

plasma; tolbutamide is IS

REFERENCE

Castoldi,D.; Monzani,V.; Tofanetti,O. Simultaneous determination of nimesulide and hydroxynimesulide in human plasma and urine by high-performance liquid chromatography, *J.Chromatogr.*, **1988**, 425, 413–418.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Serum + 2 mL 2 M HCl + 1 mL 2.5 μ g/mL chlorpropamide in water + 5 mL diethyl ether, extract, repeat extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 45°, reconstitute the residue in 200 μ L DMF:water 10:90, inject a 50 μ L aliquot.

HPLC VARIABLES

Guard column: 10 μ m μ Bondapak C16

Column: 250 \times 2.5 5 μ m Apex ODS (Rayonic Scientifique)

Mobile phase: MeCN:buffer 22.5:77.5 containing 3 mL/L Pic A (tetrabutylammonium phosphate) (Buffer was 392 mL 67 mM KH_2PO_4 and 608 mL 67 mM Na_2HPO_4 , pH 7.0.)

Injection volume: 50

Detector: UV 254

CHROMATOGRAM

Retention time: 14.8

Internal standard: chlorpropamide (10.0)

Limit of detection: 100-200 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

rat; serum; pharmacokinetics

REFERENCE

St-Hilaire,S.; Belanger,P.M. Simultaneous determinations of tolbutamide and its hydroxy and carboxy metabolites in serum and urine: application to pharmacokinetic studies of tolbutamide in the rat, *J.Pharm.Sci.*, **1989**, 78, 863-866.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Serum + 2 mL water + 200 μ L 1 (?) M HCl + 200 μ L 2.5 μ g/mL glibornuride in MeOH + 7 mL diethyl ether, mix, centrifuge at 2000 rpm for 5 min. Remove 6.5 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 500 μ L 2 mg/mL dinitrofluorobenzene in butyl acetate, heat at 120° for 1 h, cool, evaporate to dryness under a stream of nitrogen, reconstitute the residue in 150 μ L mobile phase, inject a 120 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Spherisorb ODS 2

Mobile phase: MeCN:0.4% aqueous phosphoric acid 75:25

Column temperature: 40

Flow rate: 1.2

Injection volume: 120

Detector: UV 360

CHROMATOGRAM

Retention time: 4.9

Internal standard: glibornuride (5.8)

Limit of detection: 40 ng/mL

OTHER SUBSTANCES

Extracted: chlorpropamide, glipizide, glyburide (glibenclamide), tolazamide

KEY WORDS

serum; derivatization

REFERENCE

Starkey,B.J.; Mould,G.P.; Teale,J.D. The determination of sulphonylurea drugs by HPLC and its clinical application, *J.Liq.Chromatogr.*, **1989**, 12, 1889-1896.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Serum + 200 μ L 1 M HCl + 10 μ g glyburide + 5 mL toluene, shake gently for 15 min, centrifuge at 1500 g for 3 min. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 25 μ L 6 mg/mL dinitrofluorobenzene in n-butyl acetate (prepare fresh each week, store at 4° in the dark), heat at 120° for 30 min, evaporate to dryness, reconstitute with 50 μ L mobile phase, inject a 25-50 μ L aliquot. Alternatively, filter (Amicon YMT membrane, 30000 MW cutoff) 200 μ L 100 mM NaOH while centrifuging at 4°, rinse filter with 500 μ L water, filter 1 mL serum in the same unit while centrifuging at 4° at 2500 g for 1.5 h. Remove a 700 μ L aliquot of the ultrafiltrate, add 200 μ L 1 M HCl, add 10 μ g glyburide, add 5 mL toluene, shake gently for 15 min, centrifuge at 1500 g for 3 min. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 25 μ L 6 mg/mL dinitrofluorobenzene in n-butyl acetate (prepare fresh each week, store at 4° in the dark), heat at 120° for 30 min, evaporate to dryness, reconstitute with 50 μ L mobile phase, inject a 25-50 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 7 μ m LiChrosorb RP18

Mobile phase: MeCN:10 mM pH 3.5 phosphate buffer 80:20

Flow rate: 1

Injection volume: 25-50

Detector: UV 360

CHROMATOGRAM

Retention time: 5

Internal standard: glyburide (7)

Limit of detection: 2 ng/mL

KEY WORDS

derivatization; serum; ultrafiltrate; pharmacokinetics

REFERENCE

Arcelloni,C.; Fermo,I.; Calderara,A.; Pacchioni,M.; Pontiroli,A.E.; Paroni,R. Glibenclamide and tolbutamide in human serum: Rapid measurement of the free fraction, *J.Liq.Chromatogr.*, **1990**, *13*, 175–189.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma + 200 μ L 2 M HCl + 2 mL diethyl ether, vortex 30 s, centrifuge at 1500 g for 5 min, freeze in dry ice for 5 min. Decant ether layer and evaporate it to dryness under a stream of nitrogen at 35–40°. Reconstitute extract in 100 μ L mobile phase, vortex 30 s, inject 25–50 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.1 10 μ m Versapack C18

Mobile phase: MeCN:10 mM orthophosphoric acid 50:50

Flow rate: 1

Injection volume: 25–50

Detector: UV 230

CHROMATOGRAM

Retention time: 7.67

Limit of quantitation: 1000 ng/mL

OTHER SUBSTANCES

Simultaneous: glyburide, gliclazide, chlorpropamide, glipizide, tolazamide

Noninterfering: trimethoprim, sulfamethoxazole

KEY WORDS

plasma

REFERENCE

Shenfield,G.M.; Boutagy,J.S.; Webb,C. A screening test for detecting sulfonylureas in plasma, *Ther.Drug Monit.*, **1990**, *12*, 393–397.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma + 500 μ L water + 500 μ L 100 mM HCl + 3 mL dichloromethane, mix for 15 s, centrifuge. Remove an aliquot of the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 50 μ L MeOH, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 1000 \times 8 10 μ m radial pak C18 (Waters)

Mobile phase: MeOH:0.2% acetic acid 60:40 adjusted to pH 6.7 with 1 M NaOH (Wash with MeCN at 1 mL/min for 20 min at the end of each day.)

Flow rate: 0.8

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 11

Internal standard: tolbutamide

OTHER SUBSTANCES

Extracted: chlorpropamide

KEY WORDS

plasma; tolbutamide is IS

REFERENCE

Bakare,M.T.; Mustapha,A.; Abdu-Aguye,I. An improved high-performance liquid chromatographic determination of chlorpropamide in human plasma, *Chromatographia*, **1994**, 39, 107–109.

SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 229

CHROMATOGRAM

Retention time: 3.95

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfapyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; vinoxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; videsine; mexiletine; dipyrindamole; trazodone; pipamperone; pyrimethamine; benzepiril; vincristine; propranolol; acepromazine; dothiepin; dextromoramide; fenopropfen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil;

lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpi-pramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254-262.

SAMPLE

Matrix: blood, tissue, urine

Sample preparation: Plasma. 100 μ L Plasma + 300 μ L MeOH, mix, centrifuge at 3000 rpm for 5 min, filter (0.50 μ m) the supernatant, inject an aliquot of the filtrate. Tissue. 500 μ L Tissue homogenate + 2 mL MeOH, mix for 30 min, centrifuge at 4000 rpm for 20 min, filter (0.50 μ m) the supernatant, inject an aliquot of the filtrate. Urine. Filter (0.45 μ m) urine, inject an aliquot of the filtrate.

HPLC VARIABLES

Column: 150 \times 4.6 TSKgel ODS-80TM (Tosoh)

Mobile phase: MeOH:isopropanol:0.05% phosphoric acid 35:6:60

Column temperature: 50

Detector: UV 227

CHROMATOGRAM

Limit of quantitation: 4 μ g/mL

KEY WORDS

rat; pharmacokinetics; plasma; liver; small intestine; pharmacokinetics

REFERENCE

Yamao,T.; Nakagami,H.; Furuhashi,K.; Onodera,T.; Kurosaki,Y.; Nakayama,T.; Kimura,T. Pharmacokinetics of tolbutamide following intravenous and oral administrations in rats with obstructive jaundice, *Biol.Pharm.Bull.*, **1994**, *17*, 691-695.

SAMPLE

Matrix: blood, urine

Sample preparation: Caution! Benzene is a carcinogen! 500 μ L Plasma or urine + 1 mL 50 mM HCl + 3 mL benzene, shake gently for 15 min, centrifuge at 3250 g for 5 min. Remove organic layer and evaporate it to dryness under a stream of air. Dissolve residue in 50 μ L mobile phase, vortex, inject 20 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.5 5 μ m Spherisorb ODS C18

Mobile phase: MeCN:10 mM pH 3.5 phosphate buffer 35:65

Flow rate: 1.5

Injection volume: 20

Detector: UV 275

CHROMATOGRAM

Retention time: 8.2

Internal standard: tolbutamide

OTHER SUBSTANCES

Simultaneous: glipizide

KEY WORDS

plasma; tolbutamide is IS

REFERENCE

Emilsson,H. High-performance liquid chromatographic determination of glipizide in human plasma and urine, *J.Chromatogr.*, **1987**, *421*, 319-326.

SAMPLE

Matrix: blood, urine

Sample preparation: 500 μL Plasma or 200 μL urine + 100 μL 50 (plasma) or 250 (urine) $\mu\text{g}/\text{mL}$ chlorpropamide in water + 100 μL 1 M HCl + 4 mL diethyl ether, shake for 10 min, centrifuge at 2000 g for 2-3 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100 μL MeCN:water 50:50, inject a 5 μL aliquot.

HPLC VARIABLES

Column: 250 \times 4.7 μm BST C8 (BST, Budapest)

Mobile phase: MeCN:isopropanol:0.1% orthophosphoric acid 17:17:66

Flow rate: 1.2

Injection volume: 5

Detector: UV 235

CHROMATOGRAM

Retention time: 10.5

Internal standard: chlorpropamide (8.5)

Limit of detection: 50 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

plasma; pharmacokinetics

REFERENCE

Csillag,K.; Vereczkey,L.; Gachalyi,B. Simple high-performance liquid chromatographic method for the determination of tolbutamide and its metabolites in human plasma and urine using photodiode-array detection, *J.Chromatogr.*, **1989**, 490, 355-363.

SAMPLE

Matrix: microsomal incubation

Sample preparation: 1 mL Microsomal incubation + 100 μL 2 M HCl, cool rapidly on ice, extract with 8 mL hexane:chloroform:isobutyl alcohol 100:25:0.5, add 100 μL 10 $\mu\text{g}/\text{mL}$ tolbutamide to the aqueous layer, extract with 3 mL diethyl ether. Combine (?) the organic layers and evaporate them to dryness at 60° under a stream of nitrogen, reconstitute with 100 μL mobile phase, inject a 40 μL aliquot.

HPLC VARIABLES

Mobile phase: MeCN:5 mM pH 4.3 acetate buffer 25:75

Flow rate: 1.5

Injection volume: 40

Detector: UV 230

CHROMATOGRAM

Retention time: 15

Internal standard: chlorpropamide (8)

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

rat; liver

REFERENCE

Rieutord,A.; Stupans,I.; Shenfield,G.M.; Gross,A.S. Glucalazide hydroxylation by rat liver microsomes, *Xenobiotica*, **1995**, 25, 1345-1354.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 200 μ L Microsomal incubation + 25 μ L 3 M HCl + 5 μ L 1.78 μ M chlorpropamide + 1 mL diethyl ether, vortex for 2 min, centrifuge at 2000 rpm for 5 min. Remove the organic layer and evaporate it to dryness, reconstitute the residue in 50 μ L MeOH, inject a 4 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m C18 (Analytichem)

Mobile phase: MeOH:10 mM (NH₄)H₂PO₄ 55:45, pH 5.4

Flow rate: 0.7

Injection volume: 4

Detector: UV 240

CHROMATOGRAM

Retention time: 11

Internal standard: chlorpropamide (7)

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

rat; liver

REFERENCE

Ho,J.W.; Moody,D.E. Determination of tolbutamide hydroxylation in rat liver microsomes by high-performance liquid chromatography: effect of psychoactive drugs on in vitro activity, *Life Sci.*, **1993**, 52, 21–28.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 2 mL Microsomal incubation + 50 μ L concentrated HCl + 100 μ L 50 μ g/mL chlorpropamide in MeOH + 4 mL chloroform, rotate for 30 min, centrifuge at 2000 rpm for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 60°, reconstitute the residue in 400 μ L mobile phase, inject a 100 μ L aliquot.

HPLC VARIABLES

Column: Spherisorb ODS2 C18

Mobile phase: MeOH:0.05% orthophosphoric acid 50:50

Flow rate: 1

Injection volume: 100

Detector: UV 230

CHROMATOGRAM

Retention time: 12.3

Internal standard: chlorpropamide (9.5)

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

rat; liver

REFERENCE

Ashforth,E.I.L.; Carlile,D.J.; Chenery,R.; Houston,J.B. Prediction of *in vivo* disposition from *in vitro* systems: Clearance of phenytoin and tolbutamide using rat hepatic microsomal and hepatocyte data, *J.Pharmacol.Exp.Ther.*, **1995**, 274, 761–766.

SAMPLE

Matrix: microsomal incubations

Sample preparation: 1 mL Microsomal incubation + 200 μ L 1 M orthophosphoric acid + 200 μ L 2% chlorpropamide in water + 11 mL MTBE, extract for 40 min. Remove the organic layer and evaporate it to dryness under vacuum at 40°, reconstitute the residue in 25 μ L MeCN: water 25:75, inject a 100 μ L (sic) aliquot.

HPLC VARIABLES**Column:** Novapack C18 radial compression**Mobile phase:** MeCN:water 26:74**Injection volume:** 100**Detector:** UV 230**CHROMATOGRAM****Retention time:** 25.6**Internal standard:** chlorpropamide (16.5)**OTHER SUBSTANCES****Extracted:** metabolites**KEY WORDS**

human; liver

REFERENCE

Weaver,R.J.; Dickens,M.; Burke,M.D. Hydroxylation of the antimalarial drug 58C80 by CYP2C9 in human liver microsomes: Comparison with mephenytoin and tolbutamide hydroxylations, *Biochem.Pharmacol.*, **1995**, *49*, 997-1004.

SAMPLE**Matrix:** solutions**HPLC VARIABLES****Column:** 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30**Flow rate:** 2**Detector:** UV 210**OTHER SUBSTANCES**

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitrityline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fencamfamine, fenoprofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-stilbene, imipramine, indomethacin, isocarboxtyril, isocarboxazid, isoniazid, isoproterenol, isox-suprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, methapyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methypylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic

acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, oxymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendimetrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phentermine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, prednisolone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidone, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, **1994**, *18*, 233-242.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4.5 µm LiChrospher 100 RP-8 (B)

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 6.74 (A), 6.94 (B)

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordinazepam, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazinol, mefenamic acid, meperidine, mephénytoin, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozone, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfapyrazole, sulindac, temazepam, terbutaline,

terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocainide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yohimbine, zopiclone


KEY WORDS

details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, 692, 103–119.

Tolfenamic acid



Molecular formula: C₁₄H₁₂ClNO₂

Molecular weight: 261.71

CAS Registry No.: 13710-19-5

Merck Index: 9650

SAMPLE

Matrix: solutions

Sample preparation: Inject a 20 µL aliquot of a 100–500 µg/mL solution in mobile phase.

HPLC VARIABLES

Column: 100 × 4.6 5 µm Hypersil C8 MOS 100A coated with phosphatidylcholine (95% pure soybean lecithin, Epikuron, Lucas Meyer & Co.) (Coat column by recycling a 1 mM solution of phosphatidylcholine in MeOH:water 80:20 for 24 h.)

Mobile phase: MeCN:35 mM pH 7.4 sodium phosphate buffer 40:60

Flow rate: 0.5–2

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: k' 3.31

OTHER SUBSTANCES

Also analyzed: amoxicillin, antipyrine, carbamazepine, chlorpheniramine, chlorpromazine, clonidine, codeine, desipramine, diphenhydramine, dipyrindamole, ephedrine, flufenamic acid, haloperidol, hydroxyzine, imipramine, indomethacin, lidocaine, megestrol acetate, metoprolol, nabumetone, nadolol, phenobarbital, phenol, promazine, propranolol, pyrilamine, quinidine, ropinirole, testosterone, thioridazine, verapamil

Noninterfering: acetaminophen, aspirin, azathioprine, caffeine, carprofen, chlorambucil, cimetidine, fenoterol, flurbiprofen, ibuprofen, ketoprofen, ranitidine, salicylic acid, sulfamethoxazole, theophylline, thioguanine, tiaprofenic acid, trimethoprim, valproic acid


KEY WORDS

comparison with capillary electrophoresis

REFERENCE

Hanna, M.; de Biasi, V.; Bond, B.; Salter, C.; Hutt, A.J.; Camilleri, P. Estimation of the partitioning characteristics of drugs: A comparison of a large and diverse drug series utilizing chromatographic and electrophoretic methodology, *Anal. Chem.*, **1998**, 70, 2092–2099.

Tolmetin



Molecular formula: $C_{15}H_{15}NO_3$

Molecular weight: 257.29

CAS Registry No.: 26171-23-3, 64490-92-2 (sodium salt dihydrate), 35711-34-3 (sodium salt)

Merck Index: 9655

Lednicer No.: 2 234

SAMPLE

Matrix: blood

Sample preparation: 250 μ L Plasma + 125 μ L 1 M sulfuric acid + 5 mL 10 μ g/mL diphenylacetic acid in dichloromethane, vortex for 10 s, centrifuge at 500 g for 5 min. Remove 4 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 100 μ L MeOH, inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 Zorbax ODS

Mobile phase: Gradient. MeOH:100 mM pH 5 acetate buffer 22:78 for 5 min, to 53:47 over 2 min, maintain at 53:47 for 8 min.

Column temperature: 40

Flow rate: 1.5

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 13.11

Internal standard: diphenylacetic acid (13.93)

Limit of quantitation: 1000 ng/mL

OTHER SUBSTANCES

Extracted: metabolites

Simultaneous: acetaminophen, caffeine, carbamazepine, ethosuximide, phenobarbital, phenytoin, primidone, quinidine, salicylic acid, theophylline

Noninterfering: fenopropfen, ibuprofen, indomethacin, naproxen, sulindac, valproic acid

KEY WORDS

plasma

REFERENCE

Shimek, J.L.; Rao, N.G.S.; Wahba, Khalil, S.K. High performance liquid chromatographic analysis of tolmetin, indomethacin and sulindac in plasma, *J. Liq. Chromatogr.*, **1981**, *4*, 1987–2013.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 0.5 mL 1 M HCl + 10 mL dichloromethane, shake 10 min, centrifuge at 1000 g for 5 min. Remove the organic phase and evaporate it to dryness under a stream of nitrogen at 40°. Reconstitute with 200 μ L mobile phase, inject 10–30 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.5 Partisil 10 ODS-3

Mobile phase: MeOH:water:phosphoric acid 700:300:1

Flow rate: 2.5

Injection volume: 10–30

Detector: UV 220

CHROMATOGRAM

Retention time: 3.2

Internal standard: tolmetin

OTHER SUBSTANCES

Simultaneous: ibuprofen, tolmetin is IS

KEY WORDS

plasma

REFERENCE

Lockwood,G.F.; Wagner,J.G. High-performance liquid chromatographic determination of ibuprofen and its major metabolites in biological fluids, *J.Chromatogr.*, **1982**, 232, 335-343.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 90 μ L 1 M HCl + 20 μ L 1 mg/mL p-phenylphenol in MeCN, vortex for 10 s, add 6 mL benzene:tert-butyl alcohol 90:10 (Caution! Benzene is a carcinogen!), shake mechanically for 20 min, centrifuge at 2000 g for 15 min. Remove the organic layer and evaporate it to dryness under a stream of air at 40°, reconstitute the residue in 200 μ L mobile phase, sonicate for a few s, inject an aliquot.

HPLC VARIABLES

Guard column: 40 \times 4.6 30-40 μ m Perisorb RP-18

Column: 150 \times 4.6 5 μ m Ultrasphere ODS

Mobile phase: MeCN:50 mM KH_2PO_4 48:52, pH adjusted to 4.7 with 8.5% phosphoric acid

Flow rate: 0.5 for 5 min then increase to 0.7 over 1 min, after another 4 min increase to 1.0 over 1 min

Injection volume: 20

Detector: UV 313

CHROMATOGRAM

Retention time: 6.5

Internal standard: p-phenylphenol (11)

Limit of quantitation: 50 ng/mL

OTHER SUBSTANCES

Extracted: Med 5, Med 15, metabolites

Noninterfering: sulindac, naproxen, indomethacin, flufenamic acid, ibuprofen, aspirin

KEY WORDS

plasma; rat

REFERENCE

Mancinelli,A.; Bruno,G.; Cardace,G.; Morabito,E.; Marzo,A.; Arrigoni Martelli,E. High-performance liquid chromatographic evaluation of Med 15 and its metabolites Med 5 and tolmetin in rat plasma, *J.Chromatogr.*, **1991**, 553, 81-86.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 200 μ L 100 mM pH 4 sodium acetate, extract twice with 5 mL diethyl ether. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 50°, reconstitute the residue in 200 μ L water, inject an 80 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 3.9 4 μ m Novapak C18

Mobile phase: MeCN:1 mM pH 3 phosphoric acid 32:68

Flow rate: 1

Injection volume: 80

Detector: UV 313

CHROMATOGRAM

Retention time: 11

Internal standard: tolmetin

OTHER SUBSTANCES**Extracted:** ketorolac

KEY WORDS

plasma; pharmacokinetics; tolmetin is IS

REFERENCE

Flores-Murrieta,F.J.; Granados-Soto,V.; Castañeda-Hernández,G.; Herrera,J.E.; Hong,E. Comparative bio-availability of two oral formulations of ketorolac tromethamine: Dolac and Exodol, *Biopharm.Drug Dispos.*, **1994**, *15*, 129–136.

SAMPLE**Matrix:** blood

Sample preparation: 1 mL Plasma + 500 μ L pH 1.8 phosphate buffer, extract with 1-butanol/MTBE. Remove the organic layer and add it to 500 μ L pH 6.1 ammonium acetate buffer, mix, inject an aliquot of the aqueous layer.

HPLC VARIABLES**Column:** 150 \times 4.6 5 μ m Cosmosil C18**Mobile phase:** MeCN:250 mM pH 5.0 ammonium acetate buffer 20:80**Flow rate:** 1.8**Detector:** UV 258

CHROMATOGRAM**Internal standard:** tolmetin

OTHER SUBSTANCES**Extracted:** ketoprofen (UV 350)

KEY WORDS

plasma; tolmetin is IS

REFERENCE

Shah,A.K.; Wei,G.; Lanman,R.C.; Bhargava,V.O.; Weir,S.J. Percutaneous absorption of ketoprofen from different anatomical sites in man, *Pharm.Res.*, **1996**, *13*, 168–172.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 250 \times 4 ODS (Hitachi)**Mobile phase:** MeCN:50 mM phosphoric acid 40:60 adjusted to pH 5.5 with NaOH**Column temperature:** 55**Flow rate:** 0.6**Injection volume:** 20**Detector:** UV 322

OTHER SUBSTANCES

Also analyzed: carbamazepine, fenbufen, indomethacin, ketoprofen, α -naphthoquinone, naproxen

REFERENCE

Sugawara,M.; Takekuma,Y; Yamada,H.; Kobayashi,M.; Iseki,K.; Miyazaki,K. A general approach for the prediction of the intestinal absorption of drugs: regression analysis using the physicochemical properties and drug-membrane electrostatic interactions, *J.Pharm.Sci.*, **1998**, *87*, 960–966.

SAMPLE**Matrix:** solutions

HPLC VARIABLES

Column: 250 × 4 OmniPac PAX-500 (Dionex)

Mobile phase: Gradient. A was MeCN:10 mM sodium carbonate 18:82. B was MeCN:50 mM sodium carbonate 33:67. A:B from 100:0 to 0:100 over 10 min.

Flow rate: 1

Detector: UV 254

CHROMATOGRAM

Retention time: 11

OTHER SUBSTANCES

Simultaneous: aspirin, ibuprofen, naproxen, fenbufen, indomethacin, carprofen, diflunisal

REFERENCE

Slingsby, R.W.; Rey, M. Determination of pharmaceuticals by multi-phase chromatography: Combined reversed phase and ion exchange in one column, *J.Liq.Chromatogr.*, **1990**, *13*, 107–134.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 Zorbax RX

Mobile phase: Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

Column temperature: 30

Flow rate: 2

Detector: UV 210

OTHER SUBSTANCES

Also analyzed: acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, amitrityline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspirin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphetamine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clenbuterol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamorphine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, diltiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, doxapram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenpropofen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-stilbene, imipramine, indomethacin, isocarboxystyryl, isocarboxazid, isoniazid, isoproterenol, isox-suprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesis, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, meth-apyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyl dopa, methyl dopamine, methylphenidate, methylprednisolone, methyltestosterone, methylpyrrol, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nylidrin, oxazepam, oxycodone, ox-ymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendi-metrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phenter-mine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, predni-

solone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidone, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, salicylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sulfadiazine, sulfadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetracycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, tyramine, verapamil, vincamine, warfarin, yohimbine, zoxazolamine

REFERENCE

Hill, D.W.; Kind, A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J. Anal. Toxicol.*, **1994**, *18*, 233–242.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4.5 µm LiChrospher 100 RP-8 (B)

Mobile phase: MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

Flow rate: 0.6

Injection volume: 25

Detector: UV 229

CHROMATOGRAM

Retention time: 6.45 (A), 6.29 (B)

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordi-azepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazin-
dol, mefenamic acid, meperidine, mephenytoin, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl-
dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymet-
azoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, phen-
iramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltolox-
amine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, qui-
nine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, seco-
barbital, sertraline, sotalol, spironolactone, sulfinpyrazone, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocainide, tolbutamide, trazodone, triamterene, triazolam, trifluoperazine, trifluop-
romazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, yo-
himbine, zopiclone


KEY WORDS

details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, 692, 103–119.

Tolnaftate

**Molecular formula:** C₁₉H₁₇NOS**Molecular weight:** 307.42**CAS Registry No.:** 2398-96-1**Merck Index:** 9656**Lednicer No.:** 2 211**SAMPLE****Matrix:** blood, tissue

Sample preparation: Skin. Homogenize 100–200 mg skin with 3 mL 1.15% KCl, add 3 mL naphthyl phenylacetate in dichloromethane (?), add 3 mL dichloromethane, agitate for 10 min, centrifuge at 3000 rpm for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in mobile phase, inject a 10 µL aliquot. Whole blood. 200 µL Whole blood + 200 µL water, mix, add 3 mL naphthyl phenylacetate in dichloromethane (?), add 3 mL dichloromethane, agitate for 10 min, centrifuge at 3000 rpm for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in mobile phase, inject a 10 µL aliquot.

HPLC VARIABLES**Column:** 250 × 4.6 5 µm Nucleosil C18**Mobile phase:** MeOH:water 80:20**KEY WORDS**

whole blood; skin; mouse; pharmacokinetics

REFERENCE

Szeman, J.; Ueda, H.; Szejtli, J.; Fenyvesi, E.; Watanabe, Y.; Machida, Y.; Nagai, T. Enhanced percutaneous absorption of homogenized tolnaftate/β-cyclodextrin polymer ground mixture, *Drug Des. Deliv.*, **1987**, 1, 325–332.

SAMPLE**Matrix:** formulations

Sample preparation: Solutions. Measure out amount containing 10 mg tolnaftate, dilute to 100 mL with MeCN:water 80:20, dilute 10 fold with MeCN:water 80:20, filter (0.45 µm), discard the first 1–2 mL, inject a 20 µL aliquot of the filtrate. Aerosol liquids. Discharge aerosol into a container, warm at 25° until propellant dissipates, dilute a 1 mL aliquot to 100 mL with MeCN:water 80:20, filter (0.45 µm), discard the first 1–2 mL, inject a 20 µL aliquot of the filtrate. Powders. Weigh out amount containing 5 mg tolnaftate, add 50 mL MeCN:water 80:20, shake mechanically for 15 min, make up to 100 mL with MeCN:water 80:20, let stand for 15 min, dilute an aliquot 5 fold with MeCN:water 80:20, filter (0.45 µm), discard the first 1–2 mL, inject a 20 µL aliquot of the filtrate. Aerosol powders. Discharge aerosol into a container, heat at 25° under a current of air to remove propellant, weigh out amount containing 5 mg tolnaftate, add 50 mL MeCN:water 80:20, shake mechanically for 15 min, make up to 100 mL with MeCN:water 80:20, let stand for 15 min, dilute an aliquot 5 fold with MeCN:water 80:20, filter (0.45 µm), discard the first 1–2 mL, inject a 20 µL aliquot of the filtrate. Creams, gels. Weigh out amount containing 5 mg tolnaftate, add 50 mL MeCN:water 80:20, shake mechanically for 15 min, make up to 100 mL with MeCN:water 80:20, let stand for 15 min, dilute an aliquot 5 fold with MeCN:water 80:20, filter (0.45 µm), discard the first 1–2 mL, inject a 20 µL aliquot of the filtrate.

HPLC VARIABLES**Column:** 250 × 4.6 5 µm Ultrasphere ODS**Mobile phase:** MeCN:water 80:20 containing 500 µL/L phosphoric acid, pH 2.4**Flow rate:** 1**Injection volume:** 20**Detector:** UV 257

CHROMATOGRAM**Retention time:** 5

KEY WORDSliquids; powders; aerosols; solutions; cream; gels

REFERENCEThompson,R.D.; Carlson,M. Liquid chromatographic determination of tolinaftate in commercial products, *J.Assoc.Off.Anal.Chem.*, **1991**, 74, 603–607.

SAMPLE**Matrix:** solutions**Sample preparation:** Prepare a 100 µg/mL solution in MeOH, inject an aliquot.


HPLC VARIABLES**Column:** 125 × 4 5 µm LiChrospher 60 RP-Select-B (A) or 125 × 2.4 µm Superspher 100 RP-18 (B)**Mobile phase:** MeCN:pH 2 trifluoroacetic acid 70:30 (A) or Gradient. C was MeCN:water 90:10. D was MeOH:pH 2 trifluoroacetic acid buffer 15:85. C:D from 30:70 to 80:20 over 8 min (B)**Flow rate:** 1 (A) or 0.4 (B)**Detector:** UV 254 (A), UV 245 (B)

CHROMATOGRAM**Retention time:** 9.75 (A), 9.6 (B)

OTHER SUBSTANCES**Simultaneous:** cloxyquin, chlorphenesin, naftifine, sulbentine, degradation products

REFERENCEThoma,K.; Kübler,N.; Reimann,E. Untersuchung der Photostabilität von Antimykotika. 3. Mitteilung: Photostabilität lokal wirksamer Antimykotika [Photodegradation of antimycotic drugs. 3. Communication: Photodegradation of topical antimycotics], *Pharmazie*, **1997**, 52, 362–373.

Toloxatone

**Molecular formula:** C₁₁H₁₃NO₃**Molecular weight:** 207.23**CAS Registry No.:** 29218-27-7**Merck Index:** 9659

SAMPLE**Matrix:** blood**Sample preparation:** 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol:n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 × 3.9 4 µm NovaPack C18

Mobile phase: MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH₂PO₄ adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

Column temperature: 30

Flow rate: 0.8

Injection volume: 50

Detector: UV 238

CHROMATOGRAM

Retention time: 3.35

Limit of detection: <120 ng/mL

KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; tolaxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melfalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vandesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenoprofen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafennine; alpidem; loratadine; chlorpromazine; clomipramine; carpi-pramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

REFERENCE

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254–262.

SAMPLE

Matrix: blood, CSF

Sample preparation: 50 µL Plasma or CSF + 50 µL 2 µg/mL trazodone, vortex briefly, add 20 µL 4 M NaOH, vortex briefly, add 750 µL diethyl ether, vortex for 1 min, centrifuge at 2600 g for 1 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 50 µL MeOH, inject a 20 µL aliquot.

HPLC VARIABLES**Guard column:** 25 × 4 Hibar LiChroCART C8 (Merck)**Column:** 250 × 4 5 µm LiChrospher 100 CH-8 II C8**Mobile phase:** MeCN:10 mM pH 3.0 phosphate buffer 60:40 containing 20 mM tetramethylammonium chloride**Flow rate:** 1**Injection volume:** 20**Detector:** UV 240

CHROMATOGRAM**Retention time:** 4.9**Internal standard:** trazodone (6.0)**Limit of detection:** 70 ng/mL

KEY WORDSplasma; rabbit; pharmacokinetics

REFERENCE

Vistelle,R.; Lamiabile,D.; Zinzou,M. Simple high-performance liquid chromatographic method for the measurement of tolaxatone in rabbit cerebrospinal fluid and plasma, *J.Chromatogr.*, **1989**, 490, 387–394.

SAMPLE**Matrix:** blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES**Guard column:** 20 mm long Symmetry C18**Column:** 250 × 4.6 5 µm Symmetry C8 (Waters)**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.**Column temperature:** 30**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)**Injection volume:** 10–30**Detector:** UV 204

CHROMATOGRAM**Retention time:** 14.107

KEY WORDSwhole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.


Tolpropamine

Molecular formula: $C_{16}H_{23}N$

Molecular weight: 253.39

CAS Registry No.: 5632-44-0

Merck Index: 9662



SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 $\mu\text{g/mL}$ solution in MeOH, inject a 20 μL aliquot.

HPLC VARIABLES

Column: 125 \times 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 3.3

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzotamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazepine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipranone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, flupromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserine, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypropazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, naltrexone, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscipine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxylbenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thiopropazine, thioridazine, thiothixene, thonzylamine, timolol, tocinide, tolycaine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane,I.; McKinnon,A.; Flanagan,R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J.Chromatogr.*, **1985**, 323, 191-225.

Tolrestat


Molecular formula: C₁₆H₁₄F₃NO₃S

Molecular weight: 357.35

CAS Registry No.: 82964-04-3

Merck Index: 9663

Lednicer No.: 4 56



SAMPLE

Matrix: blood, tissue

Sample preparation: Serum. 1 mL Serum + 4 mL 1 M HCl + 10 mL isopropyl ether (Caution! isopropyl ether readily forms explosive peroxides!), agitate mechanically for 15 min, centrifuge at 1500 rpm for 10 min. Remove 8.5 mL of the upper organic phase and add it to 1 mL 100 mM pH 11.0 glycine buffer, agitate mechanically for 15 min, centrifuge at 1500 rpm for 10 min. Remove 800 µL of the aqueous phase and add it to 20 µL 2.5 M phosphoric acid, inject a 50-150 µL aliquot. (For increased sensitivity use 2 mL serum and 500 µL glycine buffer.) Tissue. Homogenize (all-glass, Kontes) 50-100 mg lens or nerve tissue with 3 mL isotonic saline. 2.5 mL Homogenate + 4 mL 100 mM phosphoric acid + 10 mL isopropyl ether (Caution! isopropyl ether readily forms explosive peroxides!), agitate mechanically for 30 min, centrifuge at 1500 rpm for 15 min. Remove 9.0 mL of the upper organic phase and add it to 1 mL 100 mM pH 11.0 glycine buffer, agitate mechanically for 15 min, centrifuge at 1500 rpm for 10 min. Remove 800 µL of the aqueous phase and add it to 20 µL 2.5 M phosphoric acid, inject a 50-150 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Hypersil ODS

Mobile phase: MeCN:100 mM pH 6.0 phosphate buffer 30:70

Column temperature: 50

Flow rate: 1.1

Injection volume: 50-150

Detector: UV 226

CHROMATOGRAM

Retention time: 6.4

Limit of detection: 25 ng/mL (2 mL serum), 200 ng/mL (1 mL serum), 50 ng/g (tissue)

OTHER SUBSTANCES

Extracted: rotamer

Noninterfering: acetaminophen, diazepam, glyburide, hydrochlorothiazide, indomethacin, niacin, phenobarbital, phenytoin, propoxyphene, salicylic acid, tolbutamide

Interfering: dicoumarol, phenylbutazone

KEY WORDS

rat; dog; human; serum; sciatic nerve; lens

REFERENCE

Hicks,D.R.; Kraml,M. Determination of tolrestat, a novel aldose reductase inhibitor, in serum and tissues, *Ther.Drug Monit.*, **1984**, 6, 328-333.

SAMPLE

Matrix: bulk

Sample preparation: Prepare a 100 µg/mL solution in MeCN:buffer 35:65, inject a 5 µL aliquot. (Buffer was 10 mM phosphoric acid adjusted to pH 7.0 with 10 mM ammonium hydroxide.)

HPLC VARIABLES

Column: 150 × 3.9 5 µm Resolve octadecylsilane (Waters)

Mobile phase: MeCN:THF:buffer:40% tetrabutylammonium hydroxide in water 20.5:18.5:61.5:0.3 (Buffer was 50 mM (NH₄)H₂PO₄ adjusted to pH 3.5 with 50 mM phosphoric acid.)

Flow rate: 1

Injection volume: 5

Detector: UV 229

CHROMATOGRAM

Retention time: 4.6

Limit of quantitation: 0.05% (of tolrestat)

OTHER SUBSTANCES

Simultaneous: impurities

KEY WORDS

rugged

REFERENCE

Beaulieu,N.; Lacroix,P.M.; Sears,R.W.; Lovering,E.G. Liquid chromatographic determination of tolrestat and related compounds in raw materials, *JAOAC Int.*, **1995**, *78*, 647–650.

SAMPLE

Matrix: solutions

Sample preparation: Inject buffer solution directly.

HPLC VARIABLES

Column: 150 × 4 5 µm Spherisorb ODS

Mobile phase: MeCN:50 mM pH 6.5 phosphate buffer 30:70 containing 0.8 g/L tetraethylammonium chloride (Flush column with MeOH at the end of each day.)

Flow rate: 1

Injection volume: 50

Detector: UV 254

CHROMATOGRAM

Retention time: 5.7 (rotamer A), 3.9 (rotamer B)

OTHER SUBSTANCES

Simultaneous: rotamers

KEY WORDS

buffer

REFERENCE

Lee,H.-K.; Querijero,G. Kinetics and mechanisms of thioamide rotational isomerism: *N*-thionaphthoyl-*N*-methyl glycine derivative, *J.Pharm.Sci.*, **1985**, *74*, 273–276.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a solution in mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 150 × 4.6 5 µm Spherisorb ODS

Mobile phase: MeCN:50 mM KH₂PO₄:1 M tetrabutylammonium hydroxide 35:65:0.5


Flow rate: 1.5

Injection volume: 50

Detector: UV 220

CHROMATOGRAM**Retention time:** 8 (rotamer A), 6 (rotamer B)**OTHER SUBSTANCES****Simultaneous:** degradation products, rotamers**REFERENCE**Lee,Y.J.; Lee,H.-K. Degradation kinetics of tolrestat, *J.Pharm.Sci.*, **1990**, 79, 628–633.

Tolterodine

Molecular formula: C₂₂H₃₁NO**Molecular weight:** 325.49**CAS Registry No.:** 124937-51-5**SAMPLE****Matrix:** blood**Sample preparation:** Mix plasma or serum with two volumes of acetone, centrifuge. Remove the supernatant, wash the pellets with two portions of acetone:20 mM pH 4.5 ammonium acetate 50:50, evaporate the combined supernatant and pellet extract to dryness, dissolve the residue in MeOH:20 mM pH 4.5 ammonium acetate 10:90. Inject a 100 µL aliquot.**HPLC VARIABLES****Guard column:** 20 × 4.5 Supelco PKB 100**Column:** 150 × 4.5 Supelco PKB 100**Mobile phase:** Gradient. A was MeOH. B was 20 mM pH 4.5 ammonium acetate. A:B from 10:90 to 20:80 in 5 min, from 20:80 to 45:55 in 30 min, from 45:55 to 100:0 in 5 min, maintain at 100:0 for 10 min**Flow rate:** 1**Injection volume:** 100**Detector:** UV 280**CHROMATOGRAM****Retention time:** 29**OTHER SUBSTANCES****Extracted:** metabolites**KEY WORDS**

dog; mouse; serum; plasma

REFERENCEAndersson,S.H.G.; Lindgren,A.; Postlind,H. Biotransformation of tolterodine, a new muscarinic receptor antagonist, in mice, rats, and dogs, *Drug Metab.Dispos.*, **1998**, 26, 528–535.**SAMPLE****Matrix:** microsomal incubations**Sample preparation:** Mix 250 µL or 1 mL microsomal incubation with an equal volume of acetone, centrifuge at 3200 rpm, evaporate the acetone in the supernatant under a stream of nitrogen. Inject a 100-200 µL aliquot of the remaining supernatant.**HPLC VARIABLES****Guard column:** 20 mm Supelco PKB 100**Column:** 150 × 4.5 Supelco PKB 100

Mobile phase: Gradient. A was MeOH. B was 20 mM pH 4.5 ammonium acetate. A:B from 10:90 to 20:80 over 5 min, to 45:55 over 30 min, to 100:0 over 5 min, maintain at 100:0 for 10 min

Flow rate: 1

Injection volume: 100-200

Detector: UV 280; Radioactivity, Beckman 171 radioisotope detector

CHROMATOGRAM

Retention time: 26

OTHER SUBSTANCES

Simultaneous: metabolites

KEY WORDS

pharmacokinetics; liver

REFERENCE

Postlind,H.; Danielson,Å.; Lindgren,A.; Andersson,H.,G. Tolterodine, a new muscarinic receptor antagonist, is metabolized by cytochromes P450 2D6 and 3A in human liver microsomes, *Drug Metab.Dispos.*, **1998**, 26, 289-293.

SAMPLE

Matrix: microsomal incubations

Sample preparation: Add an equal volume of acetone to the microsomal incubation, centrifuge at 3200 rpm, evaporate acetone under a stream of nitrogen. Inject a 100-200 µL aliquot of the remaining supernatant.

HPLC VARIABLES

Guard column: 20 mm long Supelco PKB 100

Column: 150 × 4.5 Supelco PKB 100

Mobile phase: Gradient. A. MeOH. B. 20 mM ammonium acetate. A:B from 10:90 to 20:80 over 5 min, to 45:55 over 30 min, to 100:0 over 5 min, maintain at 100:0 for 10 min.

Flow rate: 1

Injection volume: 100-200

Detector: UV 280

CHROMATOGRAM

Retention time: 26

OTHER SUBSTANCES

Simultaneous: metabolites

KEY WORDS

liver

REFERENCE

Postlind,H.; Danielson,Å.; Lindgren,A.; Andersson,H.,G. Tolterodine, a new muscarinic receptor antagonist, is metabolized by cytochromes P450 2D6 and 3A in human liver microsomes, *Drug Metab.Dispos.*, **1998**, 26, 289-293.

Tolycaine


Molecular formula: C₁₅H₂₂N₂O₃

Molecular weight: 278.35

CAS Registry No.: 3686-58-6, 7210-92-6 (HCl)

Merck Index: 9679

Lednicer No.: 1 17



SAMPLE**Matrix:** solutions**Sample preparation:** Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.**HPLC VARIABLES****Column:** 125 × 4.9 Spherisorb S5W silica**Mobile phase:** MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7**Flow rate:** 2**Injection volume:** 20**Detector:** E, LeCarbone, V25 glassy carbon electrode, + 1.2 V**CHROMATOGRAM****Retention time:** 1.4**OTHER SUBSTANCES**

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazepine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipipranone, diprenorphine, dipyrindamole, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, fluopromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserin, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypromazine, methylephedrine, methylegonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, piminodine, pimozone, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, propranolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenylamine, theophylline, thietilperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tranlycypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil, xylometazoline

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R. J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J. Chromatogr.*, **1985**, *323*, 191–225.


Topotecan

Molecular formula: $C_{23}H_{23}N_3O_5$

Molecular weight: 421.45

CAS Registry No.: 123948-87-8, 119413-54-6 (HCl)

Merck Index: 9687



SAMPLE

Matrix: formulations

Sample preparation: Dilute topotecan hydrochloride injection with 5% dextrose or 0.9% NaCl, inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 μ m Cosmosil 5C18-AR (Nacalai Tesque, Japan)

Mobile phase: Gradient. A was MeCN:water:trifluoroacetic acid 15:85:0.1. B was MeCN:water:trifluoroacetic acid 40:60:0.1. A:B from 100:0 to 0:100 in 20 min, from 0:100 to 100:0 in 1 min

Flow rate: 1

Injection volume: 40

Detector: UV 228

CHROMATOGRAM

Retention time: 10.0

KEY WORDS

injections

REFERENCE

Craig,S.B.; Bhatt,U.H.; Patel,K. Stability and compatibility of topotecan hydrochloride for injection with common infusion solutions and containers, *J.Pharm.Biomed.Anal.*, **1997**, 16, 199–205.

SAMPLE

Matrix: formulations

Sample preparation: Reconstitute 4-mg topotecan hydrochloride injection with water to a nominal concentration of 1 mg/mL, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 μ m Cosmosil 5C18-AR (Nacalai Tesque, Kyoto)

Mobile phase: Gradient. A was MeCN:water:trifluoroacetic acid 15:85:0.1. B was MeCN:water:trifluoroacetic acid 40:60:0.1. A:B 100:0 for 16 min, to 0:100 over 24 min, to 100:0 over 2 min, maintain at 100:0 for 18 min

Flow rate: 1

Injection volume: 20

Detector: UV 228

CHROMATOGRAM

Retention time: 10.7-12.0

OTHER SUBSTANCES

Simultaneous: degradation products

KEY WORDS

injections; stability-indicating

REFERENCE

Patel,K.; Craig,S.B.; McBride,M.G.; Palepu,N.R. Microbial inhibitory properties and stability of topotecan hydrochloride injection, *Am.J.Health-Syst.Pharm.*, **1998**, 55, 1584–1587.