

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

Vitamin A



Molecular formula: $C_{20}H_{30}O$

Molecular weight: 286.46

CAS Registry No.: 68-26-8

Merck Index: 10150

SAMPLE

Matrix: blood

Sample preparation: 2 mL Serum + 2 mL EtOH + 10 mL hexane, mix for 30 s, centrifuge at 3000 rpm for 5 min, store the hexane layer at 15°, repeat the extraction with 10 mL hexane. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute with 200 μ L isopropanol, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m TSKgel ODS-80Ts

Mobile phase: Gradient. EtOH:water 80:20 for 11.5 min then 87:13 (step gradient)

Column temperature: 40

Flow rate: 1

Injection volume: 20

Detector: F ex 340 em 460

CHROMATOGRAM

Retention time: 7.5

OTHER SUBSTANCES

Extracted: vitamin E (F ex 298 em 325)

KEY WORDS

serum

REFERENCE

Moriyama,H.; Yamasaki,H.; Masumoto,S.; Adachi,K.; Katsura,N.; Onimaru,T. Rapid determination of vitamins A and E in serum with surfactant as a diluent by column-switching high-performance liquid chromatography, *J.Chromatogr.A*, **1998**, 798, 125–130.

SAMPLE

Matrix: blood

Sample preparation: Dilute 100 μ L serum with 900 μ L 12.62 mg/mL pyrogallol in EtOH, filter (450 μ m cellulose disk), cool at 15° in the autosampler, inject a 300 μ L aliquot onto column A and elute to waste with mobile phase A. After 3 min backflush the contents of column A onto column B with mobile phase B, after another 1 min remove column A from the circuit. Elute column B with mobile phase B for another 7.5 min then elute with mobile phase C. Monitor the effluent from column B.

HPLC VARIABLES

Column: A 10 \times 4.6 13 μ m TSK BSA-80Ts; B 15 \times 3.2 5 μ m TSK ODS-80Ts + 150 \times 4.6 5 μ m TSKgel ODS-80Ts

Mobile phase: A 200 mM sodium dodecyl sulfate solution:EtOH 70:30 containing 200 mM ethylenediaminetetraacetic acid 4 sodium salt and 0.3% phosphoric acid; B EtOH:water 80:20; C EtOH:water 87:13

Column temperature: 40

Flow rate: A 1.5; B 1; C 1
Injection volume: 300
Detector: F ex 340 em 460

CHROMATOGRAM

Retention time: 7.5
Limit of detection: 1.67 IU/dL

OTHER SUBSTANCES

Extracted: vitamin E (F ex 298 em 325)

KEY WORDS

serum; column-switching

REFERENCE

Moriyama,H.; Yamasaki,H.; Masumoto,S.; Adachi,K.; Katsura,N.; Onimaru,T. Rapid determination of vitamins A and E in serum with surfactant as a diluent by column-switching high-performance liquid chromatography, *J.Chromatogr.A*, **1998**, 798, 125–130.

SAMPLE

Matrix: blood

Sample preparation: Mix 100 μ L serum with 200 μ L MeCN and 10 μ L 20 mM ascorbic acid, centrifuge at 16000 g for 5 min. Mix the supernatant with 200 μ L water, inject a 2 μ L aliquot. (Protect all solutions from light.)

HPLC VARIABLES

Column: 50 \times 0.18 3 μ m ODS-AQ (YMC, Wilmington, NC) (The separation capillary column was formed from fused-silica capillaries (Polymicro Technologies, Phoenix) by inserting a small 50 μ m I.D. capillary ca. 15 mm into a larger 180 μ m I.D. capillary and fixed by applying epoxy (No. 353ND, Epoxy Technology, Billerica MA). A glass filter paper frit (Whatman GF/A) was inserted into the larger capillary and forced against the smaller capillary with a stream of isopropanol. The stationary phase was suspended in 3 mL isopropanol and pumped into the larger capillary until a 50 mm bed was formed. The larger and smaller diameter capillaries extended no more than 100 and 16 mm from the frit, respectively.)

Mobile phase: MeCN:MeOH:water 65:2.5:32.5 containing 1% tetrabutylammonium perchlorate, adjusted to pH 5.0 with acetic acid and 174 mM sodium acetate

Flow rate: 0.004

Injection volume: 2

Detector: E, carbon-fiber working electrode +900 mV, Ag/AgCl reference electrode (details of preparation in paper)

CHROMATOGRAM

Retention time: 11
Limit of detection: 38 pg/mL
Limit of quantitation: 70 fmol

OTHER SUBSTANCES

Extracted: isotretinoin, retinaldehyde, tretinoin

KEY WORDS

cow; serum; capillary HPLC

REFERENCE

Hagen,J.J.; Washco,K.A.; Monnig,C.A. Determination of retinoids by reversed-phase capillary liquid chromatography with amperometric electrochemical detection, *J.Chromatogr.B*, **1996**, 677, 225–231.

SAMPLE

Matrix: blood

Sample preparation: Add 100 μ L 837 ng/mL retinyl acetate in EtOH and 400 μ L EtOH to 500 μ L plasma, vortex for 30 s, add 2 mL hexane, vortex for 30 s, centrifuge at 1500 rpm for 5 min. Remove the upper hexane layer, add 2 mL hexane to the lower layer, reextract. Evaporate

the combined hexane layers to dryness under nitrogen at 40°, reconstitute the residue in 200 μ L EtOH, inject a 50 μ L aliquot.

HPLC VARIABLES

Guard column: 20 \times 4.6 5 μ m 100 Å pore size Supelguard column (Supelco)

Column: 250 \times 4.6 5 μ m 100 Å pore size Suplex pKb-100 RP (Supelco)

Mobile phase: MeOH:MTBE:water 80:20:5

Column temperature: 0

Flow rate: 0.8

Injection volume: 50

Detector: UV 328

CHROMATOGRAM

Retention time: 6.7

Internal standard: retinyl acetate (7.8)

KEY WORDS

plasma

REFERENCE

Lane, J.R.; Webb, L.W.; Acuff, R.V. Concurrent liquid chromatographic separation and photodiode array detection of retinol, tocopherols, all-trans- α -carotene, all-trans- β -carotene and the mono-cis isomers of β -carotene in extracts of human plasma, *J.Chromatogr.A*, **1997**, 787, 111–118.

SAMPLE

Matrix: blood

Sample preparation: Centrifuge 200 μ L serum, add 200 μ L IS and 200 μ L EtOH, mix on orbital shaker for 5 min, add 200 μ L water and 500 μ L hexane, mix for 10 min, centrifuge at 2000 g for 10 min at 17°, remove 300 μ L upper organic layer. Re-extract with 300 μ L hexane, mix for 10 min, centrifuge at 4000 g for 10 min at 17°, remove 300 μ L upper organic layer. Combine the organic layers, and evaporate them to dryness under vacuum in 15 min. Reconstitute the residue with 300 μ L MeOH:EtOH:hexane 88:10:2, vortex for 10 min, inject a 40 μ L aliquot.

HPLC VARIABLES

Column: 100 \times 4.6 3 μ m Adsorbosphere HS C18 + 150 \times 4.6 3 μ m Adsorbosphere HS C18 in series

Mobile phase: Gradient. A was MeCN:MeOH 60:40 containing 0.05% acetic acid. B was MeCN:MeCN:dichloromethane 45.6:30.4:24 containing 0.04% acetic acid. A:B 100:0 for 7 min then 0:100 for 10.4 min (step gradient), re-equilibrate at initial conditions for 5.6 min.

Column temperature: 37

Flow rate: 0.9

Injection volume: 40

Detector: UV 325

CHROMATOGRAM

Retention time: 5.7

Internal standard: tocol (UV 292) (10.1), echinenone (UV 450) (12.8)

Limit of detection: 25 ng/mL

OTHER SUBSTANCES

Extracted: canthaxanthine (UV 473), α -carotene (UV 450), β -carotene (UV 450), β -cryptoxanthine (UV 450), lutein (UV 450), lycopene (UV 473), vitamin E (UV 292), zeaxanthin (UV 450), nonidentified carotenoids

KEY WORDS

serum

REFERENCE

Steghens, J.-P.; van Kappel, A.L.; Riboli, E.; Collombel, C. Simultaneous measurement of seven carotenoids, retinol and α -tocopherol in serum by high-performance liquid chromatography, *J.Chromatogr.B*, **1997**, 694, 71–81.

SAMPLE**Matrix:** blood**Sample preparation:** 500 μ L Plasma + 100 μ L 5% perchloric acid, mix rapidly, add 500 μ L ethyl acetate, mix for 60-90 s, centrifuge at 13000 g for 5 min, inject an aliquot of the supernatant.

HPLC VARIABLES**Column:** Reversed-phase C18 (Waters or Bio-Rad)**Mobile phase:** MeCN:1% ammonium acetate 75:25**Flow rate:** 2.5**Detector:** UV 340

CHROMATOGRAM**Retention time:** 5

OTHER SUBSTANCES**Extracted:** isotretinoin, tretinoin

KEY WORDS

plasma; pharmacokinetics; protect from light

REFERENCEDavis,T.P.; Peng,Y.-M.; Goodman,G.E.; Alberts,D.S. HPLC, MS, and pharmacokinetics of melphalan, bisantrene and 13-cis retinoic acid, *J.Chromatogr.Sci.*, **1982**, 20, 511-516.

SAMPLE**Matrix:** blood**Sample preparation:** 500 μ L Plasma + 50 μ L 5% perchloric acid, vortex for 30 s, add 500 μ L ethyl acetate, whirl for 1 min, centrifuge at 13000 g for 1 min, inject a 50 μ L aliquot of the organic layer.

HPLC VARIABLES**Guard column:** present but not specified**Column:** 250 \times 4.6 5 μ m Ultrasphere ODS + 300 \times 4 10 μ m μ Bondapak in series**Mobile phase:** MeCN:1% ammonium acetate 95:5**Flow rate:** 2.5**Injection volume:** 50**Detector:** UV 340, UV 365

CHROMATOGRAM**Retention time:** 5.3**Limit of quantitation:** 200 ng/mL

OTHER SUBSTANCES**Extracted:** retinoic acid

KEY WORDS

protect from light; plasma

REFERENCEPeng,Y.-M.; Xu,M.-J.; Alberts,D.S. Analysis and stability of retinol in plasma, *J.Natl.Cancer Inst.*, **1987**, 78, 95-99.

SAMPLE**Matrix:** blood**Sample preparation:** Add 1 mL 0.5 μ g/mL retinyl acetate, 1 μ g/mL retinyl palmitate, and 25 μ g/mL α -tocopheryl acetate in EtOH to 1 mL serum or plasma while continuously vortexing, add 3 mL hexane, vortex for 2 min, centrifuge at 2500 g for 2 min, remove the upper phase, add 2 mL hexane to the lower layer, repeat extraction. Combine the upper layers and evaporate them to dryness under a stream of nitrogen at 37°, reconstitute the residue in 200 μ L mobile phase, inject a 40 μ L aliquot,

HPLC VARIABLES**Guard column:** C18 (Waters)**Column:** 5 μ m Biophase ODS C18 (Bioanalytical Systems)**Mobile phase:** MeCN:chloroform:isopropanol:water 78:16:3.5:2.5**Flow rate:** 2**Injection volume:** 40**Detector:** UV 460 (UV 292 for tocopherol)

CHROMATOGRAM**Retention time:** 2.60**Internal standard:** retinyl acetate (3.07), retinyl palmitate (18.66), α -tocopheryl acetate (8.33)

OTHER SUBSTANCES**Extracted:** β -carotene, vitamin E (α -tocopherol), gamma-tocopherol, α -carotene, lycopene, cryptoxanthin

KEY WORDS

serum; plasma

REFERENCEKaplan,L.A.; Miller,J.A.; Stein,E.A.; Stampfer,M.J. Simultaneous, high-performance liquid chromatographic analysis of retinol, tocopherols, lycopene, and α - and β -carotene in serum and plasma, *Methods Enzymol.*, **1990**, 189, 155–167.

SAMPLE**Matrix:** blood**Sample preparation:** 250 μ L Serum + 25 μ L 80 μ g/mL tocol in EtOH + 250 μ L 20 μ g/mL BHT (butylated hydroxytoluene) in EtOH + 1.5 mL hexane, vortex for 1 min, remove 1 mL of upper layer, add 500 μ L hexane, vortex for 1 min, remove 300 μ L of upper layer. Combine the hexane extracts, evaporate to dryness under a stream of inert gas. Reconstitute in 250 μ L 20 μ g/mL BHT in EtOH, sonicate, centrifuge if necessary, inject a 25 μ L aliquot.

HPLC VARIABLES**Column:** 250 \times 4.6 5 μ m Vydac 201TP54 (wide pore, polymerically bonded C18)**Mobile phase:** Gradient. A was MeOH:n-butanol:water 75:10:15 containing 50 mM ammonium acetate, pH 5.5. B was MeOH:n-butanol:water 88:10:2 containing 50 mM ammonium acetate, pH 5.5. A:B 100:0 for 3 min, to 0:100 over 15 min, maintain at 0:100 for 17 min**Injection volume:** 25**Detector:** UV 325 for 7 min, UV 295 for 13 min, UV 450 for 14 min or E, glassy carbon electrode, Ag/AgCl reference electrode +1050 mV for retinol, +900 mV for tocol, +750 mV for α -tocopherol, +700 mV for β -carotene

CHROMATOGRAM**Retention time:** 5**Internal standard:** tocol (13)**Limit of detection:** 4.1 μ g/mL (E), 6 μ g/mL (UV)

OTHER SUBSTANCES**Extracted:** β -carotene, vitamin E (α -tocopherol), gamma-tocopherol, lutein, zeaxanthin, cryptoxanthin, α -carotene, 9-cis- β -carotene

KEY WORDS

serum

REFERENCEMacCrehan,W.A. Determination of retinol, α -tocopherol, and β -carotene in serum by liquid chromatography, *Methods Enzymol.*, **1990**, 189, 172–181.

SAMPLE**Matrix:** blood

Sample preparation: 200 μ L Serum + 100 μ L EtOH + 100 μ L α -tocopheryl acetate in EtOH, vortex for 5 s, add 500 μ L hexane, vortex for 2 min, centrifuge at 700 g for 5 min. Remove 250 μ L of the hexane layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 200 μ L mobile phase, mix for 2 min, inject a 50 μ L aliquot.

HPLC VARIABLES

Guard column: 30 \times 4.6 10 μ m Spheri-10 RP18

Column: 150 \times 4.6 5 μ m Ultrasphere ODS

Mobile phase: MeCN:dichloromethane:MeOH 70:20:10

Flow rate: 1.2

Injection volume: 50

Detector: UV 325 for 3.5 min, UV 291 for 4.5 min, UV 450 for 6 min

CHROMATOGRAM

Retention time: 2.31

Internal standard: α -tocopheryl acetate (6.30)

Limit of detection: 16 nM

OTHER SUBSTANCES

Extracted: beta carotene, vitamin E

KEY WORDS

protect from light; serum

REFERENCE

Arnaud,J.; Fortis,I.; Blachier,S.; Kia,D.; Favier,A. Simultaneous determination of retinol, α -tocopherol and β -carotene in serum by isocratic high-performance liquid chromatography, *J.Chromatogr.*, **1991**, 572, 103–116.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma or serum + 100 μ L 6 μ g/mL 9-methylanthracene in MeOH + 1.5 mL MeCN + 100 μ L 100 mM perchloric acid, flush headspace of vial with argon, vortex, centrifuge, inject a 50 μ L aliquot of the supernatant. Sonicate serum, solvents, and mobile phase under vacuum before use. Use low actinic glassware and yellow light.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Zorbax C18

Mobile phase: MeCN:0.5% acetic acid 85:15 containing 0.05% sodium hexanesulfonate

Flow rate: 2

Injection volume: 50

Detector: UV 365

CHROMATOGRAM

Retention time: 9

Internal standard: 9-methylanthracene (5)

Limit of detection: 12 ng/mL

OTHER SUBSTANCES

Simultaneous: 13-cis-retinoic acid, all-trans-retinoic acid, 4-oxo-13-cis-retinoic acid

KEY WORDS

plasma; serum

REFERENCE

Gadde,R.R.; Burton,F.W. Simple reversed-phase high-performance liquid chromatographic method for 13-cis-retinoic acid in serum, *J.Chromatogr.*, **1992**, 593, 41–46.

SAMPLE

Matrix: blood

Sample preparation: 2.5 mL Plasma + 2.5 mL 18 ng/mL IS1 and 10 ng/mL IS2 in EtOH, shake vigorously for 20 s, centrifuge at 1200 g for 5 min, add 5 mL diethyl ether, shake vigorously, centrifuge for 5 min, extract twice more with 5 mL ether. Combine ether layers, wash with 15 mL 5% NaCl, dry over sodium sulfate, evaporate to dryness under vacuum at 35°. Dissolve residue in 1-2 mL dichloromethane, filter (0.45 μ m). Evaporate to dryness under a stream of nitrogen, make up to 100 μ L with MeCN:MeOH:dichloromethane:hexane 45:10:22.5:22.5, inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: 30 \times 4.6 5 μ m Spheri-5-C18 (Brownlee)

Column: 250 \times 4.6 5 μ m Microsorb C18 (Rainin)

Mobile phase: Gradient. MeCN:MeOH:dichloromethane:hexane 85:10:2.5:2.5 for 10 min then to 45:10:22.5:22.5 over 30 min, re-equilibrate for 15 min

Flow rate: 0.7

Injection volume: 20

Detector: UV 325

CHROMATOGRAM

Retention time: 7

Internal standard: IS1 ethyl β -apo-8'-carotenate (18), IS2 (3R)-8'-apo- β -carotene-3,8'-diol (5)

OTHER SUBSTANCES

Extracted: carotenoids, β -carotene, vitamin E (α -tocopherol)

KEY WORDS

plasma; handle under yellow lights

REFERENCE

Khachik,F.; Beecher,G.R.; Goli,M.B.; Lusby,W.R.; Smith,J.C.,Jr. Separation and identification of carotenoids and their oxidation products in the extracts of human plasma, *Anal.Chem.*, **1992**, *64*, 2111-2122.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Serum or plasma + 200 μ L 25 μ g/mL tocopheryl acetate in EtOH, vortex, add 400 μ L butanol:ethyl acetate 50:50, mix for 1 min, add 20 mg sodium sulfate, vortex for 1 min, let stand at -20° for 20 min, centrifuge at 15000 g for 2 min, inject a 10 μ L aliquot of the upper organic layer.

HPLC VARIABLES

Guard column: 5 μ m C18

Column: 110 \times 4.7 5 μ m Partisphere 5 C18 (Whatman)

Mobile phase: MeOH:butanol:water 89.5:5:5.5

Column temperature: 45

Flow rate: 1.5

Injection volume: 10

Detector: UV 340 for 3 min, UV 290 for 1.5 min, UV 280 for 10.5 min, UV 450 for 7 min

CHROMATOGRAM

Retention time: 1.7

Internal standard: tocopheryl acetate (5.3)

Limit of detection: 20 ng/mL

OTHER SUBSTANCES

Extracted: α -carotene, β -carotene, lycopene, δ -tocopherol, gamma-tocopherol, vitamin E, xanthophyll

KEY WORDS

serum; plasma; protect from light

REFERENCE

Lee,B.L.; Chua,S.C.; Ong,H.Y.; Ong,C.N. High-performance liquid chromatographic method for routine determination of vitamins A and E and β -carotene in plasma, *J.Chromatogr.*, **1992**, *581*, 41-47.

SAMPLE

Matrix: blood

Sample preparation: Dilute 1 mL serum 0.5-5 times with saline. Add 1 mL 0.25 µg/mL IS in EtOH to 1 mL diluted serum dropwise while vortexing, add 1.5 mL n-heptane, vortex for 1 min, centrifuge at 3000 rpm (Labofuge) for 15 min. Remove 1.3 mL of the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue in 40 µL MeCN:THF 50:50, inject a 5 µL aliquot.

HPLC VARIABLES

Column: 200 × 2.1 5 µm ODS Hypersil

Mobile phase: MeCN:water:THF 81.3:5.7:13

Column temperature: 40

Flow rate: 0.4

Injection volume: 5

Detector: UV 326

CHROMATOGRAM

Retention time: 1.880

Internal standard: retinol acetate (2.216)

Limit of detection: 1 ng

OTHER SUBSTANCES

Extracted: vitamin E (α-tocopherol), probucol, gamma-tocopherol, lycopene, α-carotene, β-carotene, metabolites

KEY WORDS

serum

REFERENCE

Schäfer Elinder,L.; Walldius,G. Simultaneous measurement of serum probucol and lipid-soluble antioxidants, *J.Lipid Res.*, **1992**, 33, 131-137.

SAMPLE

Matrix: blood

Sample preparation: 20-500 µL Serum + 2 volumes EtOH + 1 mL ethyl acetate + 4-7 µL of a solution containing 16 mg/mL tocopheryl acetate, 2-3 µg/mL canthaxanthin, and 10 µg/mL retinoic acid, vortex for 30 s, centrifuge for 30 s, extract the pellet twice with 0.5-1 mL portions of ethyl acetate, extract the pellet with 0.5-1 mL hexane. Combine the supernatants, add 500 µL water, vortex, centrifuge. Remove the upper organic layer and evaporate it to dryness under a stream of argon, reconstitute the residue in 100 µL MeOH:dichloromethane 2:1, inject a 10-90 µL aliquot.

HPLC VARIABLES

Guard column: C18 (Upchurch)

Column: 300 × 3.9 5 µm Resolve C18 (Waters)

Mobile phase: MeCN:dichloromethane:MeOH:1-octanol 90:15:10:0.1

Flow rate: 1

Injection volume: 10-90

Detector: UV 325

CHROMATOGRAM

Retention time: 4

Internal standard: tocopheryl acetate, canthaxanthin, retinoic acid

OTHER SUBSTANCES

Extracted: beta carotene (UV 450), carotenoids (UV 450), vitamin E (UV 290)

KEY WORDS

protect from light; serum

REFERENCE

Barua,A.B.; Kostic,D.; Olson,J.A. New simplified procedures for the extraction and simultaneous high-performance liquid chromatographic analysis of retinol, tocopherols and carotenoids in human serum, *J.Chromatogr.*, **1993**, 617, 257-264.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Serum + 100 μ L retinyl hexanoate in MeOH, extract three times with 200 μ L portions of hexane. Combine the hexane layers and evaporate them to dryness under a stream of argon, reconstitute the residue in 50 μ L isopropanol:dichloromethane, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 5 μ m Resolve C18 (Waters)

Mobile phase: MeCN:dichloromethane:MeOH:n-butanol 90:15:10:0.1 containing 0.1% ammonium acetate

Flow rate: 1

Injection volume: 20

Detector: UV 300

CHROMATOGRAM

Internal standard: retinyl hexanoate

KEY WORDS

serum; comparison with capillary electrophoresis

REFERENCE

Ma,Y.; Wu,Z.; Furr,H.C.; Lammi-Keefe,C.; Craft,N.E. Fast minimicroassay of serum retinol (vitamin A) by capillary zone electrophoresis with laser-excited fluorescence detection, *J.Chromatogr.*, **1993**, 616, 31-37.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Serum or plasma + 500 μ L EtOH containing 4.27 μ M retinyl acetate and 0.31 μ M echinenone, rotamix for 30 s, add 2 mL n-hexane, rotamix for 30 s, centrifuge at 2000 g for 2 min, repeat extraction with 2 mL n-hexane. Combine the organic layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in 50 μ L THF, make up to 200 μ L with EtOH, inject a 50 μ L aliquot.

HPLC VARIABLES

Guard column: 50 \times 4.6 5 μ m Spherisorb ODS1

Column: 250 \times 4.6 5 μ m Spherisorb ODS1

Mobile phase: Gradient. A was MeCN:MeOH 20:80 containing 100 mM ammonium acetate. B was 100 mM ammonium acetate in water. A:B from 90:10 to 100:0 over 12 min, maintain at 100:0 for 10 min, re-equilibrate at initial conditions for 5 min

Flow rate: 2

Injection volume: 50

Detector: UV 325 for 7.5 min, UV 292 for 5.5 min, then UV 450

CHROMATOGRAM

Retention time: 4.37

Internal standard: retinyl acetate (5.96), echinenone (15.15)

Limit of detection: 0.35 μ M

OTHER SUBSTANCES

Extracted: β -carotene, cryptoxanthin, lutein, lycopene, vitamin E

KEY WORDS

plasma; protect from light; serum

REFERENCE

Zaman,Z.; Fielden,P.; Frost,P.G. Simultaneous determination of vitamins A and E and carotenoids in plasma by reversed-phase HPLC in elderly and younger subjects, *Clin.Chem.*, **1993**, 39, 2229-2234.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Plasma or serum + 200 μ L 850 ng/mL retinyl acetate in EtOH, mix for 1 min, add 1 mL 0.4 g/L BHT (2,6-di-tert-butyl-4-methylphenol) in n-hexane, shake on a mechanical shaker for 10 min, centrifuge at 2000 g for 5 min, remove 800 μ L of the supernatant, evaporate to dryness at 40° under a stream of nitrogen, reconstitute in 100 μ L MeCN:THF:MeOH 68:22:7, inject a 15 μ L aliquot.

HPLC VARIABLES

Guard column: 15 \times 3.2 7 μ m Lichrosorb RP18

Column: 250 \times 4.6 5 μ m Nucleosil 100-5 C18

Mobile phase: MeCN:THF:MeOH 68:22:7 made up to 100 with 1% ammonium acetate

Flow rate: 1.5

Injection volume: 15

Detector: UV 325 for 3 min, UV 450 for 1.9 min, UV 290 for 2.5 min, UV 470 for 4.6 min, UV 450 for 3 min, then UV 325 for rest of run

CHROMATOGRAM

Retention time: 2.5

Internal standard: retinyl acetate (2.7)

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Extracted: β -carotene, vitamin E (α -tocopherol), lutein, lycopene, α -carotene, zeaxanthin, trans β -carotene, δ -tocopherol

KEY WORDS

plasma; serum; protect from sunlight

REFERENCE

Bui,M.H. Simple determination of retinol, α -tocopherol and carotenoids (lutein, all-*trans*-lycopene, α - and β -carotenes) in human plasma by isocratic liquid chromatography, *J.Chromatogr.B*, **1994**, 654, 129–133.

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma + 1 mL IS in EtOH, vortex for 30 s, add 5 mL water, add 7.5 mL n-hexane, add 300 μ L 2 M HCl, rotate for 10 min, centrifuge at 1250 g for 8 min. Remove the organic layer and evaporate it at room temperature under a stream of nitrogen. Dissolve the residue in 100 μ L mobile phase, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Spherisorb S5W

Mobile phase: n-Hexane:isopropanol:acetic acid 200:0.7:0.135

Flow rate: 0.9

Injection volume: 50

Detector: UV 350

CHROMATOGRAM

Retention time: 25

Internal standard: Ro 15-1570 (22)

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Simultaneous: 13-*cis*-retinoic acid, trans-retinoic acid

KEY WORDS

plasma; normal phase

REFERENCE

Meyer,E.; Lambert,W.E.; De Leenheer,A.P. Simultaneous determination of endogenous retinoic acid isomers and retinol in human plasma by isocratic normal-phase HPLC with ultraviolet detection, *Clin.Chem.*, **1994**, 40, 48–51.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Serum + 200 μ L nonapreno- β -carotene and retinyl butyrate in EtOH, vortex for 10 s, add 1 mL hexane, vortex for 30 s, centrifuge at 1500 g for 5 min. Remove 900 μ L of the hexane layer and evaporate it to a waxy or glassy consistency (not dryness) under vacuum, dissolve in 100 μ L EtOH, add 100 μ L MeCN, vortex, filter (0.45 μ m), inject a 30 μ L aliquot.

HPLC VARIABLES

Column: 1540 \times 4.6 5 μ m Ultramex C18 (Phenomenex)

Mobile phase: MeCN:EtOH 50:50 containing 0.1 mL/L diethylamine

Column temperature: 29

Flow rate: 0.9

Injection volume: 30

Detector: UV 300

CHROMATOGRAM

Retention time: 2.66

Internal standard: nonapreno- β -carotene (9.5, UV 450), retinyl butyrate (3.5, UV 300)

Limit of detection: 17 nM

OTHER SUBSTANCES

Extracted: β -carotene, vitamin E (α -tocopherol), lutein, zeaxanthin, β -cryptoxanthin, lycopene, α -carotene, retinyl linoleate, retinyl oleate, retinyl palmitate, retinyl stearate

KEY WORDS

serum; use gold fluorescent lamps; hold sample at 4° before injection

REFERENCE

Sowell,A.L.; Huff,D.L.; Yeager,P.R.; Caudill,S.P.; Gunter,E.W. Retinol, α -tocopherol, lutein/zeaxanthin, β -cryptoxanthin, lycopene, α -carotene, trans- β -carotene, and four retinyl esters in serum determined simultaneously by reversed-phase HPLC with multiwavelength detection, *Clin.Chem.*, **1994**, *40*, 411–416.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Serum + 2 mL diluting agent (amber tube), mix for 10 s, add 5 mL n-hexane, vortex for 30 s, centrifuge at 1600 g at 4° for 5 min. Remove the organic phase and concentrate to < 1 mL under vacuum below 30°, evaporate the rest of the solvent under nitrogen, dissolve the residue in 25 μ L MeCN:MeOH 2:1, mix for 30 s, centrifuge at 8000 g for 1 min, inject a 20 μ L aliquot. (Diluting agent was MeCN:100 mM ammonium acetate 25:75, pH adjusted to 5.5 with acetic acid.)

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Chemcosorb 5-ODS-H

Mobile phase: MeCN:MeOH:100 mM ammonium acetate 46.7:23.3:30, pH adjusted to 7.0

Column temperature: 50

Flow rate: 1

Injection volume: 20

Detector: UV 340

CHROMATOGRAM

Retention time: 31.9

Limit of quantitation: 0.5 ng/mL

OTHER SUBSTANCES

Simultaneous: 13-cis-retinoic acid, all-trans-retinoic acid

KEY WORDS

serum

REFERENCE

Takeda,N.; Yamamoto,A. Simultaneous determination of 13-cis- and all-trans-retinoic acids and retinol in human serum by high-performance liquid chromatography, *J.Chromatogr.B*, **1994**, 657, 53–59.

SAMPLE

Matrix: blood

Sample preparation: 10 μ L Serum + 30 μ L 200-500 ng/mL retinyl acetate in isopropanol:dichloroethane 2:1 + 5 μ L glacial acetic acid, vortex for 30 s, centrifuge for 1 min, inject a 10-20 μ L aliquot.

HPLC VARIABLES

Guard column: C18 (Upchurch)

Column: 150 \times 4.6 Ultracarb 5 ODS30 (Phenomenex)

Mobile phase: MeCN:dichloromethane:MeOH 85:12:3 containing 0.1% ammonium acetate (dissolve ammonium acetate in MeOH first)

Flow rate: 1

Injection volume: 10-20

Detector: UV 335

CHROMATOGRAM

Retention time: 4.5

Internal standard: retinyl acetate (5.5)

Limit of detection: 10 ng/mL

Limit of quantitation: 50 ng/mL

OTHER SUBSTANCES

Extracted: retinoic acid

KEY WORDS

protect from light; serum

REFERENCE

Barua,A.B.; Kostic,D.; Barua,M.; Olson,J.A. Determination of retinol and retinoic acid in capillary blood by high-performance liquid chromatography, *J.Liq.Chromatogr.*, **1995**, 18, 1459–1471.

SAMPLE

Matrix: blood

Sample preparation: 1 mL Serum + 2.5 mL EtOH, mix for 5 min, add 5 mL n-hexane, mix vigorously, centrifuge at 2000 g for 5 min, repeat extraction with 3 mL n-hexane. Combine the n-hexane layers and evaporate them to dryness under a stream of nitrogen, reconstitute the residue in dichloromethane, inject an aliquot.

HPLC VARIABLES

Guard column: μ Bondapak C18 Guard-Pak + 50 mm long C18 ODS(4) (Shimadzu)

Column: 250 \times 4.6 5 μ m Vydac 201 TP 54 C18

Mobile phase: MeOH:MeCN 90:10 (Every 100 injections wash column with MeOH:MeCN:dichloromethane 8:1:1.)

Flow rate: 1

Detector: UV 324

CHROMATOGRAM

Retention time: 3.8, 3.85 (9-cis)

OTHER SUBSTANCES

Extracted: beta carotene (UV 451), vitamin E (UV 291)

KEY WORDS

serum

REFERENCE

Ben-Amotz, A. Simultaneous profiling and identification of carotenoids, retinols, and tocopherols by high performance liquid chromatography equipped with three-dimensional photodiode array detection, *J.Liq.Chromatogr.*, **1995**, *18*, 2813–2825.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Serum + 100 μ L retinyl acetate in MeOH, mix, extract three times with 200 μ L portions of hexane. Combine the extracts and evaporate them to dryness under a stream of argon, reconstitute with 25 μ L isopropanol:dichloromethane, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 3.9 5 μ m Resolve C18 (Waters)

Mobile phase: MeOH:water 95:5

Flow rate: 1

Injection volume: 20

Detector: UV 325

KEY WORDS

comparison with capillary electrophoresis; serum

REFERENCE

Shi, H.; Ma, Y.; Humphrey, J.N.; Craft, N.E. Determination of vitamin A in dried human blood spots by high-performance capillary electrophoresis with laser-excited fluorescence detection, *J.Chromatogr.B*, **1995**, *665*, 89–96.

SAMPLE

Matrix: blood

Sample preparation: 100 μ L Serum + 300 μ L MeCN, vortex for 1.5 min, centrifuge at 10500 g for 1 min, inject an aliquot of the supernatant.

HPLC VARIABLES

Guard column: 10 \times 4.5 5 μ m CLC-G-ODS C18 (Shimadzu)

Column: 250 \times 4.6 5 μ m CLC-ODS C18 (Shimadzu)

Mobile phase: MeCN:water 90:10

Column temperature: 37

Flow rate: 2

Injection volume: 100

Detector: UV 328

CHROMATOGRAM

Retention time: 5.5

Limit of detection: 10 ng/mL

KEY WORDS

protect from light; serum

REFERENCE

Siddiqui, F.Q.; Malik, F.; Fazli, F.R. Determination of serum retinol by reversed-phase, *J.Chromatogr.B*, **1995**, *666*, 342–346.

SAMPLE

Matrix: blood

Sample preparation: 200 μ L Serum + 200 μ L 650 ng/mL tocopheryl acetate in MeOH, vortex for 30 s, add 200 μ L n-hexane, shake for 15 min, centrifuge at 3000 rpm for 10 min. Remove 120 μ L of the hexane layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 20 μ L dichloromethane, add 100 μ L MeCN:MeOH 50:50, inject a 50 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4 5 μ m Nucleosil C18

Mobile phase: MeCN:MeOH:dichloromethane 50:45:5

Flow rate: 0.7

Injection volume: 50

Detector: F ex 325 em 480 for 5 min, then ex 295 em 330 for 10 min, then ex 325 em 480

CHROMATOGRAM

Retention time: 2.9

Internal standard: tocopheryl acetate (9.3)

OTHER SUBSTANCES

Extracted: vitamin E, γ -tocopherol, retinyl palmitate, β -carotene (UV 450), α -carotene (UV 450)

KEY WORDS

serum

REFERENCE

Yakushina,L.; Taranova,A. Rapid HPLC simultaneous determination of fat-soluble vitamins, including carotenoids, in human serum, *J.Pharm.Biomed.Anal.*, **1995**, 13, 715–718.

SAMPLE

Matrix: blood, tissue

Sample preparation: Plasma. 100 μ L Plasma + 100 μ L 0.9% NaCl + 200 μ L MeOH, vortex for 30 s, let stand for 10 min, add 400 μ L chloroform, vortex for 4 min. Remove the chloroform layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in MeOH, inject an aliquot. Tissue. Homogenize (liver with Mikro-dismembrator II in liquid nitrogen; other tissue with Ultra-Turrax) tissue with 3 mL 1% acetic acid containing 1 mg/mL ascorbic acid and 10 mM EDTA, add 2 mL MeOH, vortex for 30 s, let stand for 10 min, add 4 mL chloroform, vortex for 4 min. Remove the chloroform layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in MeOH, inject an aliquot.

HPLC VARIABLES

Column: 125 \times 4.6 3 μ m Hypersil ODS

Mobile phase: MeCN:dichloromethane:MeOH:water 70:10:15:5

Flow rate: 0.5 for 13 min, to 1 over 1 min, maintain at 1 for 10 min, to 1.5 over 1 min, maintain at 1.5 for 21 min, to 2 over 1 min, maintain at 2 for 10 min, return to 0.5 over 1 min, maintain at 0.5 for 2 min.

Injection volume: 50

Detector: UV 350

CHROMATOGRAM

Retention time: 6

Limit of detection: 10 ng/mL

OTHER SUBSTANCES

Extracted: beta carotene (UV 445), vitamin E (UV 292)

KEY WORDS

rat; protect from light; liver; plasma; lung

REFERENCE

Van Vliet,T.; Van Schaik,F.; Van Schoonhoven,J.; Schrijver,J. Determination of several retinoids, carotenoids and E vitamers by high-performance liquid chromatography. Application to plasma and tissues of rats fed a diet rich in either β -carotene or canthaxanthin, *J.Chromatogr.*, **1991**, 553, 179–186.

SAMPLE

Matrix: blood, urine

Sample preparation: Condition a Lichrolut RP-18 (Merck) SPE cartridge with 3 mL MeOH and 3 mL water. Mix 40 μ L plasma or 100 μ L urine with twice the volume of MeCN for 2 min, add 100 μ L water, centrifuge at 3500 rpm for 15 min, evaporate the supernatant under nitrogen at 45° to remove the organic solvents, add slowly to the SPE cartridge, discard the effluent,

elute with 3 mL MeOH. Evaporate the eluate to dryness under a stream of nitrogen at 45°. Reconstitute the residue with 500 μ L MeOH containing 3.5 μ g/mL IS. Inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 220 \times 4.6 5 μ m Spherisorb RP-18

Mobile phase: MeCN:MeOH 70:30

Flow rate: 1.5

Injection volume: 10

Detector: UV 290

CHROMATOGRAM

Retention time: 2.6

Internal standard: anthraquinone (1.89)

Limit of detection: 2 ng

OTHER SUBSTANCES

Extracted: vitamin E, vitamin E acetate

KEY WORDS

plasma; SPE

REFERENCE

Papadoyannis, I.N.; Tsioni, G.K.; Samanidou, V.F. Simultaneous determination of nine water and fat soluble vitamins after SPE separation and RP-HPLC analysis in pharmaceutical preparations and biological fluids, *J. Liq. Chromatogr. Rel. Technol.*, **1997**, 20, 3203–3231.

SAMPLE

Matrix: cheese

Sample preparation: 500 mg Cheese + 2 mL 60% KOH + 2 mL 95% EtOH + 1 mL 1% NaCl + 5 mL 6% pyrogallol in EtOH, flush tube with nitrogen, seal, heat at 70° for 30 min, cool in ice water, add 15 mL 1% NaCl, extract twice with 15 mL portions of n-hexane:ethyl acetate 90:10. Combine the organic layers and evaporate them to dryness, dissolve the residue in 2 mL mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Ultrasphere Si

Mobile phase: Gradient. A was n-hexane:isopropanol 99:1. B was n-hexane. A:B 50:50 for 7 min; to 90:10 over 4 min, maintain at 90:10 for 7 min, to 50:50 over 1 min, maintain at 50:50 for 4 min. (About every 100 injections recondition column with 50 mL dichloromethane, 50 mL isopropanol, and 50 mL dichloromethane.)

Flow rate: 1.5

Injection volume: 20

Detector: F ex 325 em 475 (vitamin A), UV 450 (β -carotene), and F ex 325 em 475 for 3.5 min, ex 280 em 475 for 10.5 min, ex 325 em 475 for 9 min (others)

CHROMATOGRAM

Retention time: 18

Limit of detection: 0.32 ng

OTHER SUBSTANCES

Extracted: β -carotene, vitamin E (α -tocopherol), β -tocopherol, gamma-tocopherol, δ -tocopherol, 13-cis-retinol

KEY WORDS

normal phase; cheese

REFERENCE

Panfili, G.; Manzi, P.; Pizzoferrato, L. High-performance liquid chromatographic method for the simultaneous determination of tocopherols, carotenes, and retinol and its geometric isomers in Italian cheeses, *Analyst*, **1994**, 119, 1161–1165.

SAMPLE

Matrix: culture media

Sample preparation: 100 μ L Culture media + 200 μ L ice-cold EtOH, mix thoroughly, let stand for 15 min, centrifuge at 12000 g for 15 min, inject an aliquot of the supernatant.

HPLC VARIABLES

Guard column: Whatman CO:PELL ODS guard column

Column: 100 \times 8 5 μ m Nova-Pak C18 (radial-packed)

Mobile phase: MeOH:100 mM pH 7.0 ammonium acetate 90:10

Flow rate: 1

Detector: UV 340

CHROMATOGRAM

Retention time: 12.50

OTHER SUBSTANCES

Extracted: isotretin, motretinid, acitretin, all-trans-retinoic acid, retinal, etretinate

REFERENCE

Kochhar,D.M.; Penner,J.D.; Minutella,L.M. Biotransformation of etretinate and developmental toxicity of etretin and other aromatic retinoids in teratogenesis bioassays, *Drug Metab.Dispos.*, **1989**, 17, 618-624.

SAMPLE

Matrix: flour, milk

Sample preparation: Mix 5 g powdered milk or flour with 20 mL 500 g/L NaOH in water. Warm at 30° for 3 min, add 100 mL EtOH and 2 mL 200 g/L hydroquinone in EtOH. Heat at 80° for 30 min. Cool, add 100 mL water and 50 mL diethyl ether. Shake and add 50 mL petroleum ether, shake. Remove the organic layer and reextract the aqueous layer as above. Combine the organic layers, filter (cellulose), wash twice with 100 mL water, evaporate to dryness under reduced pressure, dissolve in 1 mL MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Guard column: 20 \times 4.6 5 μ m Zorbax ODS

Column: 150 \times 4.6 5 μ m Zorbax ODS

Mobile phase: MeOH:water 92:8

Column temperature: 25 \pm 1

Flow rate: 1.5

Injection volume: 20

Detector: UV 330

CHROMATOGRAM

Retention time: 13.45

Limit of detection: 20 ng/g

Limit of quantitation: 120 ng/g

KEY WORDS

powdered milk

REFERENCE

Ake,M.; Mandrou,B.; Malan,A. Determination of vitamin A in milk and flour consumed by one- to four-year-old children in Côte d'Ivoire, *J.AOAC Int.*, **1998**, 81, 111-114.

SAMPLE

Matrix: food

Sample preparation: 500 mg Cheese + 2 mL 60% KOH + 2 mL 95% EtOH + 1 mL 1% NaCl + 5 mL 6% pyrogallol in EtOH, flush tube with nitrogen, seal, heat at 70, for 30 min, cool in ice water, add 15 mL 1% NaCl, extract twice with 15 mL portions of n-hexane:ethyl acetate 90:10 (Analyst 1994, 119, 1161). Combine the organic layers and evaporate them to dryness, dissolve the residue in MeOH:dichloromethane 90:10, inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.5 µm Supelco C18

Mobile phase: MeOH

Flow rate: 2

Detector: F ex 325 em 475

CHROMATOGRAM

Retention time: 2.4

Limit of detection: 400 pg

OTHER SUBSTANCES

Extracted: vitamin E (α-tocopherol, β and γ-tocopherol, δ-tocopherol) (F ex 280 em 325), sterols (UV 208), carotenes (UV 450)

KEY WORDS

cheese

REFERENCE

Manzi,P.; Panfili,G.; Pizzoferrato,L. Normal and reversed-phase HPLC for more complete evaluation of tocopherols, retinols, carotenes and sterols in dairy products, *Chromatographia*, **1996**, 43, 89–91.

SAMPLE

Matrix: food

Sample preparation: 500 mg Cheese + 2 mL 60% KOH + 2 mL 95% EtOH + 1 mL 1% NaCl + 5 mL 6% pyrogallol in EtOH, flush tube with nitrogen, seal, heat at 70, for 30 min, cool in ice water, add 15 mL 1% NaCl, extract twice with 15 mL portions of n-hexane:ethyl acetate 90:10 (Analyst 1994, 119, 1161). Combine the organic layers and evaporate them to dryness, dissolve the residue in n-hexane:2-propanol 99:1, inject an aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Ultrasphere Si (Beckman)

Mobile phase: Gradient. A was n-hexane:isopropanol 99:1. B was n-hexane. A:B 50:50 for 7 min; to 90:10 over 4 min, maintain at 90:10 for 7 min, to 50:50 over 1 min, maintain at 50:50 for 4 min. (About every 100 injections recondition column with 50 mL dichloromethane, 50 mL isopropanol, and 50 mL dichloromethane.)

Flow rate: 1.5

Detector: F ex 325 em 475

CHROMATOGRAM

Retention time: 15 (13-cis-retinol), 17.3 (all trans-retinols)

Limit of detection: 100 pg (13-cis-retinol), 300 pg (all trans-retinols)

OTHER SUBSTANCES

Extracted: vitamin E (α, β, γ, δ-tocopherol) (F ex 280 em 325), sterols (UV 208), carotenes (UV 450)

KEY WORDS

cheese; normal phase

REFERENCE

Manzi,P.; Panfili,G.; Pizzoferrato,L. Normal and reversed-phase HPLC for more complete evaluation of tocopherols, retinols, carotenes and sterols in dairy products, *Chromatographia*, **1996**, 43, 89–91.

SAMPLE

Matrix: food, milk, tissue

Sample preparation: 40 mL Fluid sample or 1-10 g solid sample homogenized with 30 mL water + 12 mL 60% KOH + 80 mL EtOH + 0.5 mL 1% BHT in EtOH + 0.5 g ascorbic acid, stir for 16 h under nitrogen, make up to 250 mL with water and EtOH so that the water:EtOH ratio is 50:50, add 20 mL to a 150 × 25 Kieselguhr SPE cartridge, let stand for 20 min, elute with 50 mL light petroleum, force out all the eluate with nitrogen, evaporate the eluate under reduced pressure, take up the residue in 2-50 mL isooctane, inject an aliquot.

HPLC VARIABLES

Guard column: 7 × 2.3 µm Spherisorb SW silica gel

Column: 100 × 2.3 µm Spherisorb SW silica gel

Mobile phase: n-Hexane:1-octanol 99.7:0.3

Flow rate: 0.4

Injection volume: 10

Detector: UV 325

CHROMATOGRAM

Retention time: 6.4 (11-cis), 7.1 (11,13-di-cis), 8.4 (13-cis), 9.3 (9,13-di-cis), 11.0 (9-cis), 11.6 (7-cis), 12.5 (all-trans)

Limit of detection: 140 ng/mL

KEY WORDS

normal phase; liver; sausage; sour cream; food; SPE

REFERENCE

Brinkmann,E.; Dehne,L.; Oei,H.B.; Tiebach,R.; Baltes,W. Separation of geometrical retinol isomers in food samples by using narrow-bore high-performance liquid chromatography, *J.Chromatogr.A*, **1995**, 693, 271–279.

SAMPLE

Matrix: formula, reference material

Sample preparation: 10 g Zero-control reference material (ZRM) powder + 50 g hot water, mix. 6.5 g Reconstituted ZRM or 3.5 g concentrated commercial formula + 10 mL boiling isopropanol, mix, add 7.5 g anhydrous magnesium sulfate to the ZRM and 4 g to the concentrated commercial formula, mix thoroughly, add 25 mL hexane:ethyl acetate 85:15, add 1 mL 360 µg/mL BHT, mix, homogenize (Polytron) for 1 min, filter through 60 mL coarse-porosity fritted glass filter using vacuum, wash with two 15 mL portions of hexane:ethyl acetate 85:15. Re-extract with 20 mL hexane:ethyl acetate 85:15 and 5 mL isopropyl alcohol, homogenize for 1 min, filter through 60 mL coarse-porosity fritted glass filter using vacuum, wash with two 15 mL portions of hexane:ethyl acetate 85:15. Mix the combined filtrate with 500 mg anhydrous magnesium sulfate, evaporate to dryness, add 15 mL hexane to the residue, filter (0.45 µm nylon) using vacuum, wash with three 7 mL portions of hexane, evaporate to 1 mL with nitrogen at 45°, dilute to 10 mL with hexane, inject a 50 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 mm 5 µm Lichrosorb Si 60

Mobile phase: Hexane:isopropanol 99.875:0.125

Flow rate: 1

Injection volume: 50

Detector: F ex 325 em 470

CHROMATOGRAM

Retention time: 3.4 (vitamin A palmitate)

Limit of quantitation: 300 ng/mL

OTHER SUBSTANCES

Also analyzed: vitamin E

KEY WORDS

normal phase; soy-based infant formula

REFERENCE

Chase,G.W.,Jr.; Long,A.R.; Eitenmiller,R.R. Liquid chromatographic method for analysis of all-rac- α -tocopheryl acetate and retinyl palmitate in soy-based infant formula using a zero-control reference material (ZRM) as a method development tool, *JAOAC Int.*, **1998**, 81, 577–581.

SAMPLE

Matrix: formulations

Sample preparation: Reconstitute 28 g milk based infant formula with ca. 145 g 78-80° water, mix thoroughly. Add 15 mL boiling isopropanol to 6.5 g reconstituted infant formula, mix thoroughly. Add 7.5 g magnesium sulfate, 30 mL hexane:ethyl acetate 85:15 and 1 mL BHT. Homogenize mixture for 1 min, filter through a coarse porosity glass filter by vacuum, wash the magnesium sulfate cake with two 15 mL portions of hexane:ethyl acetate 85:15. Repeat the extraction with 5 mL isopropanol and 15 mL hexane:ethyl acetate 85:15. Add 1 g magnesium sulfate to the combined filtrate and evaporate it to dryness under nitrogen. Dissolve the residue in 10 mL hexane, filter (0.45 µm), evaporate to a volume of less than 5 mL at 45°. Dilute to 10 mL with mobile phase. Inject a 50 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Lichrosorb Si 60

Mobile phase: Hexane:isopropanol 99.5:0.5

Flow rate: 0.5

Injection volume: 50

Detector: F ex 325 em 470

CHROMATOGRAM

Retention time: 2.9

Limit of detection: 187 ng/mL

KEY WORDS

infant formula; normal phase

REFERENCE

Chase, Jr., G.W.; Eitenmiller, R.R.; Long, A.R. Liquid chromatographic analysis of all-RAC- α -tocopheryl acetate, tocopherols, and retinyl palmitate in SRM 1846, *J.Liq.Chromatogr.Rel.Technol.*, **1997**, 20, 3317-3327.

SAMPLE

Matrix: formulations

Sample preparation: Tablets, capsules. Add 5 crushed tablets or the contents from 5 capsules to 10 mL DMSO, add 15 mL hexane, shake at 60° for 45 min, centrifuge at 3000 rpm for 10 min, remove hexane layer, add 15 mL hexane, vortex for 5 min at room temperature, remove hexane layer, repeat hexane extraction three more times, combine all hexane layers, filter, make up to 100 mL with hexane, dilute if necessary, inject a 100 µL aliquot. Syrup. 10 mL Syrup + 10 mL DMSO + 15 mL ether:hexane 10:90, vortex for 5 min, remove hexane layer, repeat extraction four more times, combine the extracts, evaporate with nitrogen until the ether is removed, filter, make up to 100 mL with hexane, dilute if necessary, inject a 100 µL aliquot.

HPLC VARIABLES

Column: 150 × 4.6 3 µm amino bonded phase (Chromatography Sciences Co.)

Mobile phase: Hexane:isopropanol 99:1

Flow rate: 1-2

Injection volume: 100

Detector: UV 265

CHROMATOGRAM

Retention time: 12 (13-cis), 13 (9-cis), 14.5 (all-trans)

OTHER SUBSTANCES

Simultaneous: cholecalciferol (vitamin D3), ergocalciferol (vitamin D2)

KEY WORDS

work under subdued light; tablets; capsules; syrup

REFERENCE

Beaulieu, N.; Curran, N.M.; Gagné, C.; Gravelle, M.; Lovering, E.G. Liquid chromatographic methods for vitamins A and D in multivitamin-mineral formulations, *J.Assoc.Off.Anal.Chem.*, **1989**, 72, 247-254.

SAMPLE

Matrix: formulations

Sample preparation: Weigh out a sample containing about 2500 UI of vitamin A, add 20 mL MeOH, protect from light, stir vigorously for 2-2.5 h, filter (0.45 μm), inject a 2 μL aliquot of the filtrate.

HPLC VARIABLES

Column: 150 \times 3.9 4 μm Nova-Pack C18

Mobile phase: MeCN:THF:water 55:37:8

Flow rate: 1.5

Injection volume: 2

Detector: UV 325

CHROMATOGRAM

Retention time: 1.08 (retinol), 1.25 (retinyl acetate)

KEY WORDS

protect from light; capsules; tablets

REFERENCE

Genestar,C.; Grases,F. Determination of vitamin A in pharmaceutical preparations by high-performance liquid chromatography with diode-array detection, *Chromatographia*, **1995**, *40*, 143-146.

SAMPLE

Matrix: formulations

Sample preparation: Grind tablets to a fine powder, weigh out 100-150 mg, mix with sea sand.

Extract with supercritical carbon dioxide (Dionex) in the dynamic mode at 250 atmospheres and 40° for 15 min (restrictor 60°, gaseous flow rate 190-220 mL/min), collect in 6 mL THF at 0°, make up to 10 mL with THF, inject an aliquot. [Alternatively, add 5 crushed tablets or the contents from 5 capsules to 10 mL DMSO, add 15 mL hexane, shake at 60° for 45 min, centrifuge at 3000 rpm for 10 min, remove hexane layer, add 15 mL hexane, vortex for 5 min at room temperature, remove hexane layer, repeat hexane extraction three more times, combine all hexane layers, dilute with THF, inject a 25 μL aliquot. *J.Assoc.Off.Anal.Chem.* 1989, *72*, 247.]

HPLC VARIABLES

Guard column: 4 \times 4 5 μm (Merck)

Column: 250 \times 4 5 μm Lichrospher CH-8

Mobile phase: MeOH:MeCN 75:25

Flow rate: 1.3

Injection volume: 50

Detector: UV 325

CHROMATOGRAM

Retention time: 9 (vitamin A palmitate)

OTHER SUBSTANCES

Simultaneous: vitamin E acetate (UV 280)

KEY WORDS

tablets; SFE

REFERENCE

Scalia,S.; Ruberto,G.; Bonina,F. Determination of vitamin A, vitamin E, and their esters in, *J.Pharm.Sci.*, **1995**, *84*, 433-436.

SAMPLE

Matrix: formulations

Sample preparation: Mix 100-150 mg of the formulation with celite and extract with supercritical carbon dioxide at 250 atmospheres at 40° at 190-220 mL/min with the restrictor at 100° (Dionex SFE-703), collect in 4 mL THF:MeOH 80:20 at 0°, make up to 5 mL, inject an aliquot.

HPLC VARIABLES**Guard column:** 10 μm Guard-Pak (Waters)**Column:** 300 \times 3.9 10 μm μ Bondapak C18**Mobile phase:** MeCN:MeOH 25:75**Flow rate:** 1.5**Detector:** UV 325

CHROMATOGRAM**Retention time:** 10 (vitamin A palmitate)

OTHER SUBSTANCES**Simultaneous:** vitamin E acetate (UV 280)

KEY WORDSSFE; cream; lotion; protect from light

REFERENCE

Scalia,S.; Renda,A.; Ruberto,G.; Bonina,F.; Menegatti,E. Assay of vitamin A palmitate and vitamin E acetate in cosmetic creams and lotions by supercritical fluid extraction and HPLC, *J.Pharm.Biomed.Anal.*, **1995**, *13*, 273–277.

SAMPLE**Matrix:** milk**Sample preparation:** Mix 50 mL milk with 30 mL ethanolic KOH (10:30). Saponify the mixture at 80° for 20 min. Extract twice with 10 mL n-hexane. Evaporate the extracts to dryness, reconstitute the residue with 1 mL mobile phase, inject a 5 μL aliquot.

HPLC VARIABLES**Column:** 250 \times 4.6 5 μm C18 (Alltech)**Mobile phase:** MeOH:EtOH 80:20 (A) or EtOH:water 95:5 (B)**Flow rate:** 1**Injection volume:** 5**Detector:** UV 250

CHROMATOGRAM**Retention time:** 4.45 (A), 4 (B)

OTHER SUBSTANCES**Extracted:** retinal, isotretinoin, tretinoin, vitamin D2, vitamin D3, vitamin E, vitamin K1, vitamin K2

REFERENCE

Gong,B.Y.; Ho,J.W. Simultaneous separation and detection of ten common fat-soluble vitamins in milk, *J.Liq.Chromatogr.Rel.Technol.*, **1997**, *20*, 2389–2397.

SAMPLE**Matrix:** milk**Sample preparation:** Add 500 mg ascorbic acid, 50 mL EtOH, and 10 mL 60% KOH to 25 g liquid or reconstituted powdered milk under a nitrogen stream, stir overnight at room temperature. Extract with three 50 mL portions of n-hexane and two 25 mL portions of n-hexane by shaking for 2 min each. Combine the n-hexane extracts, wash with 50 mL portions of water containing a few drops of phenolphthalein until the aqueous phase appears colorless, add 1 g butylated hydroxytoluene, filter through a Whatman No.1 filter containing 20 g anhydrous sodium sulfate, concentrate the filtrate under reduced pressure at 40°, reconstitute with 10 mL MeOH, filter (0.45 μm), inject an aliquot of the filtrate.

HPLC VARIABLES**Guard column:** Tracer Spherisorb ODS 2 C18**Column:** 250 \times 4.6 5 μm Tracer Spherisorb ODS 2 C18 (Teknokroma, Spain)**Mobile phase:** MeCN:MeOH:water 1:95:4**Injection volume:** 20

Detector: UV 323 for 14 min then UV 292

CHROMATOGRAM

Retention time: 5.5

Limit of detection: 10 ng/mL

Limit of quantitation: 20 ng/mL

OTHER SUBSTANCES

Extracted: vitamin E

REFERENCE

Albal-Hurtado,S.; Novella-Rodríguez,S.; Veciana-Nogués,M.; Mariné-Font,A. Determination of vitamins A and E in infant milk formulae by high-performance liquid chromatography, *J.Chromatogr.A*, **1997**, 778, 243–246.

SAMPLE

Matrix: milk

Sample preparation: 1 g Powdered milk or 25 mL liquid milk + alcoholic KOH, let stand overnight, extract with hexane. Remove the organic layer and evaporate it to dryness, reconstitute the residue in MeOH, filter, inject an aliquot. (Alcoholic KOH was prepared from 50 mL EtOH and 15 mL 60% KOH in water.)

HPLC VARIABLES

Guard column: 15 × 3.2 7 µm RP18 (Brownlee)

Column: 220 × 4.6 5 µm OD-224 RP18 (Brownlee)

Mobile phase: MeOH:water 99:1 containing 2.5 mM acetic acid/sodium acetate

Flow rate: 1.25

Injection volume: 10

Detector: E, EG & G PAR Model 400, MP 1304 glassy carbon series dual electrode, E1 (upstream) -1100 mV, E2 (downstream) +700 mV (Condition electrodes for 30 min at E1 -1200 mV and E2 +1500 mV at the start of each day.)

CHROMATOGRAM

Retention time: 4

Limit of detection: 0.06 ng

OTHER SUBSTANCES

Extracted: vitamin E

REFERENCE

Delgado-Zamarreño,M.M.; Sanchez Perez,A.; Gomez Perez,M.C.; Fernandez Moro,M.A.; Hernandez Mendez,J. Determination of vitamins A, E and K1 in milk by high-performance liquid chromatography with dual amperometric detection, *Analyst*, **1995**, 120, 2489–2492.

SAMPLE

Matrix: milk

Sample preparation: Dilute milk to 30% with water, mix with reagent and pass through a 5 m × 0.5 mm i.d. PTFE tube knotted reactor at 1.25 mL/min, mix with 2.5 M acetic acid, pass onto a Sep-Pak Plus C18 SPE cartridge for 5 min, wash SPE cartridge with MeOH:water 40:60 for 4 min, elute SPE cartridge with MeOH for 4 min, inject the last 100 µL of the eluate. (Reagent was 50 mL EtOH + 15 mL 60% aqueous NaOH + 5 mL 10% ascorbic acid.)

HPLC VARIABLES

Guard column: 15 × 3.2 7 µm RP18 (Brownlee)

Column: 220 × 4.6 5 µm OD-224 RP18 (Brownlee)

Mobile phase: MeOH:water 99:1 containing 2.5 mM acetic acid-sodium acetate

Flow rate: 1

Injection volume: 100

Detector: E, glassy carbon working electrode +1300 mV

CHROMATOGRAM

Retention time: 6

Limit of detection: 34.9 nM

OTHER SUBSTANCES

Extracted: vitamin E, vitamin D3

KEY WORDS

SPE; derivatization

REFERENCE

Delgado-Zamarreño, M.M.; Sanchez-Perez, A.; Gomez-Perez, M.C.; Hernandez-Mendez, J. Directly coupled sample treatment-high-performance liquid chromatography for on-line automatic determination of liposoluble vitamins in milk, *J.Chromatogr.A*, **1995**, 694, 399–406.

SAMPLE

Matrix: silicone oils

Sample preparation: Condition a 1 g Si Bond-Elut SPE cartridge with 5 mL n-hexane. Mix 1 g silicone oil with 2 mL dichloromethane, vortex for 2 min, centrifuge at 3000 g. Withdrawn the supernatant, repeat this procedure twice, filter (0.45 μ m), heat the filtrate at 50°, expose to a stream of helium for 30 min. Add 2.5 μ g retinol acetate, 2.5 μ g α -tocopherol acetate, and 25 μ g BHT. Add the mixture to the SPE cartridge, elute with 500 μ L MeOH, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Zorbax C8

Mobile phase: Gradient. A was MeCN:200 mM ammonium acetate 72:25. B was MeOH:water 95:5. A:B 100:0 for 10 min, to 0:100 over 1 min, maintain at 0:100 for 14 min

Flow rate: 2 for 10 min then 1.5

Injection volume: 20

Detector: UV 350

CHROMATOGRAM

Retention time: 4

Internal standard: retinol acetate (9.5)

Limit of detection: 51.4 ng/mL

Limit of quantitation: 171.2 ng/mL

OTHER SUBSTANCES

Extracted: cholesterol (UV 210), retinal (UV 350), retinoic acid (UV 350), α -tocopherol acetate, vitamin E (UV 210)

KEY WORDS

ophthalmic silicone oils; SPE

REFERENCE

Del Nozal, M.J.; Bernal, J.L.; Marinero, P. Simultaneous HPLC determination of cholesterol, α -tocopherol, retinol, retinal and retinoic acid in silicone oils used as vitreous substitutes in eye surgery, *J.Liq.Chromatogr.Rel.Technol.*, **1996**, 19, 1151–1167.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4.6 8 μ m Unisphere-PBD (polybutadiene on alumina) (Biotage, Charlottesville, VA)

Mobile phase: MeOH:water 92:8

Flow rate: 1

Detector: UV 330

CHROMATOGRAM

Retention time: 4

OTHER SUBSTANCES

Simultaneous: retinoic acid, retinol acetate, vitamin E

REFERENCE

Jedrejewski,P.T.; Taylor,L.T. Comparison of silica-, alumina-, and polymer-based stationary phases for reversed-phase liquid chromatography, *J.Chromatogr.Sci.*, **1995**, 33, 438–445.

SAMPLE

Matrix: tissue

Sample preparation: 50 mg Tissue + 50 μ L 560 μ g/mL vitamin K in EtOH, extract twice with 1 mL n-hexane using a sonicator, centrifuge at 2000 g for 5 min. Evaporate the supernatant to dryness, reconstitute it in 200 μ L chloroform:MeOH 25:75, inject an aliquot. (Caution! Chloroform is a carcinogen!)

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Nucleosil 120 C18

Mobile phase: MeOH:water 96.5:3.5

Column temperature: 40

Flow rate: 2

Detector: UV 325

CHROMATOGRAM

Retention time: ca. 2

Internal standard: vitamin K (9)

OTHER SUBSTANCES

Extracted: vitamin A_p, vitamin E (F ex 295 em 350)

KEY WORDS

rat; liver; placenta; brain

REFERENCE

Barbas,C.; Castro,M.; Bonet,B.; Viana,M.; Herrera,E. Simultaneous determination of vitamins A and E in rat tissues by high-performance liquid chromatography, *J.Chromatogr.A*, **1997**, 778, 415–420.

SAMPLE

Matrix: tissue

Sample preparation: 100-120 mg Tadpole embryos + 100 mL isopropanol, sonicate on ice, vortex for 1 min, centrifuge at 4000 g at 4° for 20 min, inject an aliquot of the supernatant.

HPLC VARIABLES

Guard column: 20 \times 4 10 μ m LiChrosorb RB 18

Column: 125 \times 4.6 3 μ m Spherisorb ODS II

Mobile phase: Gradient. MeOH:40 mM pH 7.3 ammonium acetate from 55:45 to 100:0 over 20 min

Flow rate: 1.6

Detector: UV 354

CHROMATOGRAM

Retention time: 15.2

OTHER SUBSTANCES

Extracted: isotretinoin, tretinoin, metabolites

KEY WORDS

handle under yellow light; tadpoles; embryos

REFERENCE

Creech Kraft,J.; Kimelman,D.; Juchau,M.R. *Xenopus Laevis*: A model system for the study of embryonic retinoid metabolism. I. Embryonic metabolism of 9-*cis*- and all-*trans*-retinals and retinols and their corresponding acid forms, *Drug Metab.Dispos.*, **1995**, 23, 72–82.

SAMPLE

Matrix: tissue

Sample preparation: 100-120 mg Frog embryos + 100 μ L isopropanol, sonicate on ice, vortex for 1 min, centrifuge at 4° at 4000 g for 20 min, inject an aliquot of the supernatant.

HPLC VARIABLES

Guard column: 20 \times 4 10 μ m LiChrosorb RB 18

Column: 125 \times 4.6 3 μ m Spherisorb ODS II

Mobile phase: Gradient. MeOH:40 mM pH 7.3 ammonium acetate from 55:45 to 100:0 over 18 min.

Flow rate: 1.6

Detector: UV 354

CHROMATOGRAM

Retention time: 15

OTHER SUBSTANCES

Extracted: retinoic acid

KEY WORDS

frog; embryo

REFERENCE

Creech Kraft,J.; Juchau,M.R. *Xenopus laevis*: A model system for the study of embryonic retinoid metabolism. III. Isomerization and metabolism of all-*trans*-retinoic acid and 9-*cis*-retinoic acid and their dysmorphogenic effects in embryos during neurulation, *Drug Metab.Dispos.*, **1995**, 23, 1058–1071.

SAMPLE

Matrix: tissue

Sample preparation: Homogenize (7 mL Tenbroeck grinder, Bioblock) 2-10 mg liver and 1.5 mL 1% pyrogallol in EtOH, shake at 4° for 30 min, add 100 μ L water, add 3 mL n-hexane, extract, centrifuge at 4° at 3000 g for 15 min, repeat extraction. Evaporate the hexane layers separately as rapidly as possible under a stream of nitrogen at 50°, reconstitute with 200 μ L 30 μ M retinyl acetate in MeOH, inject a 50 μ L aliquot. (For total retinol homogenize liver with 1.5 mL 10% KOH in 95% EtOH, heat at 60° for 30 min, cool in an ice bath, add 800 μ L water, extract twice with 1.9 mL n-hexane. Evaporate the hexane layers separately as rapidly as possible under a stream of nitrogen at 50°, reconstitute with 200 μ L MeOH, inject a 50 μ L aliquot.)

HPLC VARIABLES

Guard column: 15 \times 3.2 7 μ m RP-8 (Brownlee)

Column: 250 \times 4.6 5 μ m Supelcosil LC-8

Mobile phase: MeOH:water 94:6

Flow rate: 1.5

Injection volume: 50

Detector: UV 325

CHROMATOGRAM

Retention time: 3

Internal standard: retinyl acetate (3.7)

Limit of quantitation: 6 nmole/g

OTHER SUBSTANCES

Extracted: retinyl oleate, retinyl palmitate, retinyl stearate

KEY WORDS

liver; protect from light

REFERENCE

Got,L.; Gousson,T.; Delacoux,E. Simultaneous determination of retinyl esters and retinol in human livers by reversed-phase high-performance liquid chromatography, *J.Chromatogr.B*, **1995**, 668, 233–239.


Vitamin B12

Molecular formula: $C_{63}H_{88}CoN_{14}O_{14}P$

Molecular weight: 1355.38

CAS Registry No.: 68-19-9

Merck Index: 10152



SAMPLE

Matrix: blood, feces

Sample preparation: Feces. Condition a Sep-Pak C18 SPE cartridge with 2 mL MeCN and 6 mL water. 1 g Feces + 3 mL buffer, rotate at 4° for 24 h, centrifuge at 19000 g for 1 h. Remove 3 mL of the supernatant and add cadmium acetate to a concentration of 200 mM, let stand for 2 h, add four volumes of EtOH preheated to 80°, mix vigorously, heat at 80° for 20 min, cool in an ice bath, centrifuge at 2000 g for 10 min, remove the supernatant, mix the precipitate with two volumes of cold EtOH:water 80:20, centrifuge. Combine the supernatants and evaporate them to dryness under reduced pressure at 40°, reconstitute the residue in 2 mL water, add to the SPE cartridge, wash with 12 mL water, elute with 6 mL t-butanol:water 20:80, evaporate the eluate to dryness, reconstitute with 2 mL 1% acetic acid, add to the Amberlite XAD2 column, wash with 12 mL 1% acetic acid, wash with 12 mL MeOH:1% acetic acid 10:90, elute with 30 mL MeOH:1% acetic acid 50:50, lyophilize the eluate at -80°, reconstitute with 2 mL water, filter (0.2 μ m), lyophilize, reconstitute with 300 μ L 1% acetic acid, inject a 250 μ L aliquot. Plasma. Condition a Sep-Pak C18 SPE cartridge with 2 mL MeCN and 6 mL water. Add cadmium acetate to a concentration of 200 mM, let stand for 2 h, add four volumes of EtOH preheated to 80°, mix vigorously, heat at 80° for 20 min, cool in an ice bath, centrifuge at 2000 g for 10 min, remove the supernatant, mix the precipitate with two volumes of cold EtOH:water 80:20, centrifuge. Combine the supernatants and evaporate them to dryness under reduced pressure at 40°, reconstitute the residue in 2 mL water, add to the SPE cartridge, wash with 12 mL water, elute with 6 mL t-butanol:water 20:80, evaporate the eluate to dryness, lyophilize the eluate at -80°, reconstitute with 2 mL water, filter (0.2 μ m), lyophilize, reconstitute with 300 μ L 1% acetic acid, inject a 250 μ L aliquot (Meth. Enzymol. 1986, 123, 3). (Buffer was 100 mM pH 7.4 sodium phosphate buffer containing 5 U/mL aprotinin, 0.02 mM phenylmethylsulfonyl fluoride, 3 mM sodium azide, and 0.05% Triton X100. Prepare Amberlite XAD2 column as follows. Suspend 30 g resin in 50 mL acetone, filter, wash with 50 mL acetone, dry at 80°, suspend in 200 mL MeOH, allow to settle, discard the supernatant, repeat 2-4 times until supernatant is clear, suspend in 100 mL MeOH, fill a 330 \times 8 glass column to a bed height of 40 mm (Meth. Enzymol. 1986, 123, 3).)

HPLC VARIABLES

Column: 250 \times 5 μ m LiChrospher RP18 glass column

Mobile phase: Gradient. A was 85 mM phosphoric acid adjusted to pH 3.0 with triethanolamine. B was MeCN. A:B from 90:10 to 50:50 over 20 min.

Flow rate: 0.5

Injection volume: 250

Detector: UV 365

CHROMATOGRAM

Retention time: 17.5

KEY WORDS

protect from light; plasma; SPE

REFERENCE

Djalali,M.; Gueant,J.-L.; Lambert,D.; el Kholty,S.; Saunier,M.; Nicolas,J.-P. High-performance liquid chromatographic separation and dual competitive binding assay of corrinoids in biological material, *J.Chromatogr.*, **1990**, 529, 81–91.

SAMPLE

Matrix: blood, formulations, urine

Sample preparation: Tablets. Powder tablets, dissolve in water, inject a 10 μ L aliquot. Injections. Dilute with water, inject a 10 μ L aliquot. Plasma, urine. Condition a Lichrolut RP-18 (Merck) SPE cartridge with 3 mL MeOH and 3 mL water. Mix 40 μ L plasma or 100 μ L urine with twice the volume of MeCN for 2 min, add 100 μ L water, centrifuge at 3500 rpm for 15 min, evaporate the supernatant under nitrogen at 45° to remove the organic solvents, add slowly to the SPE cartridge, collect the eluate. Evaporate to dryness under a stream of nitrogen at 45°. Reconstitute the residue with 500 μ L MeOH containing 4.2 μ g/mL IS. Inject a 10 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.5 μ m Lichrosorb RP-18

Mobile phase: Gradient. A was MeOH. B was 50 mM ammonium acetate. A:B from 5:95 to 15:85 over 6 min, to 30:70 over 7 min, maintain at 30:70 over 7 min

Flow rate: 1

Injection volume: 10

Detector: UV 270

CHROMATOGRAM

Retention time: 17.12

Internal standard: xanthine (4.65)

Limit of detection: 5 ng

OTHER SUBSTANCES

Extracted: ascorbic acid, folic acid, niacin, niacinamide, riboflavin

KEY WORDS

plasma; SPE; tablets; injections

REFERENCE

Papadoyannis,I.N.; Tsioni,G.K.; Samanidou,V.F. Simultaneous determination of nine water and fat soluble vitamins after SPE separation and RP-HPLC analysis in pharmaceutical preparations and biological fluids, *J.Liq.Chromatogr.Rel.Technol.*, **1997**, 20, 3203–3231.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 μ L MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) μ L aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 \times 4.6 5 μ m Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10-30

Detector: UV 207.5

CHROMATOGRAM

Retention time: 3.777

KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

SAMPLE

Matrix: bulk, formulations, premix

Sample preparation: Weigh out amount containing 20-100 µg vitamin B12, extract with DMSO: water 50:50 (or 0.5% ammonium pyrrolidine dithiocarbamate and 2% citric acid in DMSO: water 50:50 for samples containing ascorbic acid with iron or copper) by shaking at 55° for 45 min, centrifuge, filter (0.45 µm), inject an aliquot.

HPLC VARIABLES

Column: 2540 × 4.6 µm Bondapak C18

Mobile phase: Gradient. MeOH:water from 15:85 to 50:50 over 15 min.

Detector: UV 546

CHROMATOGRAM

Retention time: 8

Limit of quantitation: 500 ng/mL

OTHER SUBSTANCES

Noninterfering: riboflavin

REFERENCE

Hudson,T.S.; Subramanian,S.; Allen,R.J. Determination of pantothenic acid, biotin, and vitamin B12 in nutritional products, *J.Assoc.Off.Anal.Chem.*, **1984**, 67, 994–998.

SAMPLE

Matrix: formula, milk

Sample preparation: Mix 8.0 g powdered infant milk with 10 mL water to it. Mix the diluted powder or 10.5 g liquid infant milk with 1 g solid trichloroacetic acid, shake thoroughly with magnetic stirring for 10 min, centrifuge at 1250 g for 10 min, add 3 mL 4% trichloroacetic acid to the solid residue, mix thoroughly for 10 min, centrifuge, discard the solid phase. Combine the two acid extracts and make up to 10 mL with 4% trichloroacetic acid, filter (0.45 µm), inject an aliquot of the filtrate.

HPLC VARIABLES

Guard column: 5 µm Tracer Spherisorb ODS 2 C18 (Teknokroma, Spain)

Column: 250 × 4.6 5 µm Tracer Spherisorb ODS 2 C18 (Teknokroma, Spain)

Mobile phase: MeOH:buffer 15:85 (Buffer was 5 mM octanesulfonic acid and 0.5% triethylamine, pH 3.6.)

Flow rate: 1

Injection volume: 20

Detector: UV 261 for 6 min, UV 287 for 2 min, UV 290 for 5 min, UV 282 for 3 min, UV 268 for 3.5 min, UV 361 for 20.5 min, UV 246 for 20 min

CHROMATOGRAM

Retention time: 22

Limit of quantitation: ≤300 ng/mL

OTHER SUBSTANCES

Extracted: thiamine, riboflavin, pyridoxine, folic acid, niacinamide, pyridoxal, pyridoxamine

REFERENCE

Albalá-Hurtado,S.; Veciana-Nogués,M.; Izquierdo-Pulido,M.; Mariné-Font,A. Determination of water-soluble vitamins in infant milk by high-performance liquid chromatography, *J.Chromatogr.A*, **1997**, 778, 247–253.

SAMPLE

Matrix: formulations

HPLC VARIABLES

Column: 150 × 6 5 µm Capcell Pak C8 (Shiseido, Japan)

Mobile phase: MeOH:50 mM KH₂PO₄ containing 5 mM tetra-n-butylammonium phosphate 15:85, adjusted to pH 2.6 with 5% orthophosphoric acid (After one week of use, wash the column with water and MeOH:water 70:30 at 1 mL/min for 30 min.)

Column temperature: 30

Flow rate: 1

Injection volume: 10-20

Detector: UV 215

CHROMATOGRAM

Retention time: 27

OTHER SUBSTANCES

Simultaneous: chlorpheniramine, dipotassium glycyrrizate, fumaric acid, m-hydroxybenzoic acid, p-hydroxybenzoic acid, maleic acid, neostigmine methylsulfate, pyridoxine, tetrahydrozoline

Noninterfering: chondroitin sulfate, lysozyme

KEY WORDS

ophthalmic solutions; ion-pair agents

REFERENCE

Yamato,S.; Nakajima,M.; Shimada,K. Simultaneous determination of chlorpheniramine and maleate by high-performance liquid chromatography using tetra-n-butylammonium phosphate as an ion-pair reagent, *J.Chromatogr.A*, **1996**, 731, 346–350.

SAMPLE

Matrix: formulations

Sample preparation: Condition a Sep-Pak C18 SPE cartridge with 5 mL MeCN and 10 mL water. Dissolve 20 g powder formulation in 60 mL water, add 10 g NaCl, heat at 50° until completely dissolved, let stand at room temperature for 30 min, make up to 100 mL with water, wash with 10 mL hexane for 3 min. Add the aqueous layer to the SPE cartridge, elute with 8 mL MeCN:water 50:50, evaporate the eluate to dryness under reduced pressure at 50°, reconstitute with 4 mL water, inject a 2 mL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Capcellpak C18 (Shiseido)

Mobile phase: MeCN:water 13:87

Flow rate: 0.6

Injection volume: 2000

Detector: UV 550

CHROMATOGRAM

Retention time: 13

Limit of quantitation: 2.2 ng/g

KEY WORDS

SPE; powders

REFERENCE

Iwase, H. Ultramicrodetermination of cyanocobalamin in elemental diet by solid-phase extraction and high-performance liquid chromatography with visible detection, *J. Chromatogr.*, **1992**, 590, 359–363.

SAMPLE

Matrix: formulations

HPLC VARIABLES

Column: 100 × 4 3 μm Hypersil BDS-C18

Mobile phase: Gradient. MeCN:water adjusted to pH 2.1 from 0.3:99.7 to 25:75 over 11 min

Flow rate: 0.5

Detector: UV 220

CHROMATOGRAM

Retention time: 9.5

OTHER SUBSTANCES

Simultaneous: biotin, caffeine, citric acid, folic acid, niacinamide, niacin, pantothenic acid, pyridoxine, riboflavin, saccharin, thiamine, ascorbic acid

KEY WORDS

tablets

REFERENCE

Hewlett Packard Leaflet 12-5091-7351 *EUS*, **1993**.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 × 4.6 5 μm Accubond Amino (J & W)

Mobile phase: MeCN:buffer 10 :90 (Buffer was 20 mM phosphoric acid adjusted to pH 3.0 with 20 mM NaOH.)

Flow rate: 1

Detector: UV 254

CHROMATOGRAM

Retention time: 1.9

OTHER SUBSTANCES

Simultaneous: p-aminobenzoic acid, niacinamide, pyridoxal, pyridoxamine, thiamine, riboflavin, pyridoxine

REFERENCE

J & W Catalog, 1992-3, p. 277.


Vitamin D2

Molecular formula: $C_{28}H_{44}O$

Molecular weight: 396.66

CAS Registry No.: 50-14-6

Merck Index: 10156



SAMPLE

Matrix: blood

Sample preparation: 3-5 mL Plasma + 3 volumes ether, shake horizontally at 120 oscillations/min for 5 min, let stand for 1-2 min, freeze in dry ice/acetone, repeat ether extraction, extract aqueous layer with 4 volumes dichloromethane:MeOH 75:25, shake for 3-5 min, add 1 mL MeOH, shake for 15 s. Remove the organic layer and wash it twice with 100 mM pH 10.5 phosphate buffer, combine the dichloromethane and ether extracts and evaporate them to dryness under a stream of nitrogen, chromatograph on a 155×6 Sephadex LH-20 column with hexane:chloroform:MeOH 90:10:10, discard first 1 mL eluate, collect next 2.5 mL eluate. Evaporate to dryness under a stream of nitrogen, reconstitute with 500 μ L hexane:chloroform 95:5 and chromatograph on a 145×6 Lipidex 5000 column (Packard) with hexane:chloroform 95:5, discard first 6 mL, collect next 4 mL. Evaporate to dryness under a stream of nitrogen, reconstitute the residue in 150 μ L mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 250×4.5 Zorbax Sil

Mobile phase: Dichloromethane:isopropanol 99.75:0.25

Flow rate: 2

Injection volume: 150

Detector: UV 254

CHROMATOGRAM

Retention time: 16

OTHER SUBSTANCES

Interfering: vitamin D3

KEY WORDS

plasma; normal phase; turkey; chicken; cow; pig; sheep; SPE

REFERENCE

Horst, R.L.; Littledike, E.T.; Riley, J.L.; Napoli, J.L. Quantitation of vitamin D and its metabolites and their plasma concentrations in five species of animals, *Anal. Biochem.*, **1981**, *116*, 189-203.

SAMPLE

Matrix: blood

Sample preparation: 2-3 mL Serum + 3.75 volumes hexane:isopropanol 1:2, shake for 30 min, let stand for 5 min, add 1.25 volume hexane, shake for 5 min, centrifuge at 600 g for 5 min, remove the upper organic layer, extract the aqueous layer twice with 1.25 volumes of hexane. Combine all the organic layers and evaporate them to dryness under a stream of nitrogen at 35°, reconstitute the residue in 50 μ L hexane:isopropanol 80:20 and inject on to a 200×1 column of Kieselgel Si-60, Size A (Merck) and elute with hexane:isopropanol 80:20 at 2 mL/min, monitor at UV 254 and collect the vitamin D fraction at 6.0-7.5 min, evaporate the eluate to dryness under a stream of nitrogen. Reconstitute with 10 μ L hexane:isopropanol 95:5 and inject on to a 50×4.6 Polygosil (Macherey-Nagel) + 250×4.6 7 μ m Nucleosil 50-7 silica column, elute with hexane:isopropanol 85:5 at 2 mL/min, monitor the effluent at UV 254 and

collect the vitamin D fraction at 4.5 min, evaporate the vitamin D fraction, reconstitute, inject an aliquot.

HPLC VARIABLES

Guard column: 30-40 μm pellicular C18 (Vydac)

Column: 250 \times 4.6 5 μm TP C18 (Vydac)

Mobile phase: MeOH

Flow rate: 1

Detector: UV 254

CHROMATOGRAM

Retention time: 5.5

OTHER SUBSTANCES

Extracted: vitamin D3

KEY WORDS

serum; normal phase; reverse phase

REFERENCE

Parviainen, M.T.; Savolainen, K.E.; Korhonen, P.H.; Alhava, E.M.; Visakorpi, J.K. An improved method for routine determination of vitamin D and its hydroxylated metabolites in serum from children and adults, *Clin. Chim. Acta*, **1981**, *114*, 233-247.

SAMPLE

Matrix: formula

Sample preparation: Condition a 2.8 mL 500 mg silica SPE cartridge with 4 mL dichloromethane:isopropanol 80:20 and 5 mL dichloromethane:isopropanol 99.8:0.2. 15 mL Formula + 4 mL 46 ng/mL vitamin D₃ in EtOH + 15 mL ethanolic KOH, shake at 60° for 30 min, cool to room temperature, add 15 mL water, add 60 mL hexane, shake vigorously for 1.5 min, let stand for 10 min, discard aqueous layer. Wash the hexane layer with 15 mL water, add 15 mL water and 1 drop phenolphthalein to the hexane layer, add 10% acetic acid dropwise with shaking until the aqueous layer is colorless. Filter the hexane layer through anhydrous sodium sulfate, wash through with a few mL hexane, evaporate to dryness under reduced pressure at 40°, reconstitute with 2 mL dichloromethane:isopropanol 99.8:0.2, add to the SPE cartridge, rinse flask with 1 mL dichloromethane:isopropanol 99.8:0.2, add rinse to cartridge, wash with 2 mL dichloromethane:isopropanol 99.8:0.2, elute with 7 mL dichloromethane:isopropanol 99.8:0.2. Evaporate the eluate to dryness under a stream of nitrogen at 40°, reconstitute the residue in 1 mL MeCN, inject an aliquot. (Prepare ethanolic KOH by dissolving 140 g KOH in 310 mL EtOH, add 50 mL water. Prepare fresh each day.)

HPLC VARIABLES

Column: 250 \times 4.6 5 μm C18 (not end-capped)

Mobile phase: Gradient. MeCN:MeOH:ethyl acetate 91:9:0 for 28 min, to 0:0:100 over 0.5 min, maintain at 0:0:100 for 2.5 min, return to initial conditions over 0.5 min, re-equilibrate for 2.5 min.

Column temperature: 27

Flow rate: 0.7 for 28 min, to 2.5 over 0.5 min, maintain at 2.5 for 4.5 min, return to initial conditions over 1 min

Injection volume: 250

Detector: UV 265

CHROMATOGRAM

Retention time: 19.5

Internal standard: vitamin D₃ (23)

KEY WORDS

protect from light and oxygen; SPE

REFERENCE

Sliva, M.G.; Sanders, J.K. Vitamin D in infant formula and enteral products by liquid chromatography: Collaborative study, *JAOAC Int.*, **1996**, *79*, 73-80.

SAMPLE**Matrix:** formulations**Sample preparation:** Powder tablets, weigh out powder equivalent to 200 IU vitamin D3, add 5 µg vitamin D2 and 10 mL EtOH:water 50:50. Extract with 15 mL hexane 3 times. Remove the organic layer and dry it under reduced pressure. Dissolve the residue in 1 mL MeOH:water 90:10 and inject a 20 µL aliquot.

HPLC VARIABLES**Column:** 150 × 4.5 µm TSK-gel ODS 80TM (TOSOH, Japan)**Mobile phase:** MeOH:water 90:10**Flow rate:** 1**Injection volume:** 20**Detector:** MS, Hitachi M-1000, APCI interface, drift voltage 20 V, focus voltage 120 V, vaporizer 399°, desolvation chamber 200°, multiplier voltage 2 kV

CHROMATOGRAM**Retention time:** 14.4**Internal standard:** vitamin D2**Limit of detection:** 400 pg

OTHER SUBSTANCES**Simultaneous:** vitamin D3

KEY WORDStablets; vitamin D2 is IS

REFERENCE

Adachi,T.; Nishio,M.; Yunoki,N.; Hayashi,H. Determination of vitamin D₃ and D₂ in multi-vitamin tablets by high-performance liquid chromatography/atmospheric pressure chemical ionization mass spectrometry, *Anal.Sci.*, **1994**, 10, 457-460.

SAMPLE**Matrix:** formulations**Sample preparation:** Tablets, capsules. Add 5 crushed tablets or the contents from 5 capsules to 10 mL DMSO, add 15 mL hexane, shake at 60° for 45 min, centrifuge at 3000 rpm for 10 min, remove hexane layer, add 15 mL hexane, vortex for 5 min at room temperature, remove hexane layer, repeat hexane extraction three more times, combine all hexane layers, filter, make up to 100 mL with hexane, dilute if necessary, inject a 100 µL aliquot. Syrup. 10 mL Syrup + 10 mL DMSO + 15 mL ether:hexane 10:90, vortex for 5 min, remove hexane layer, repeat extraction four more times, combine the extracts, evaporate with nitrogen until the ether is removed, filter, make up to 100 mL with hexane, dilute if necessary, inject a 100 µL aliquot.

HPLC VARIABLES**Column:** 150 × 4.6 3 µm amino bonded phase (Chromatography Sciences Co.)**Mobile phase:** Hexane:isopropanol 99:1**Flow rate:** 1-2**Injection volume:** 100**Detector:** UV 265

CHROMATOGRAM**Retention time:** 9

OTHER SUBSTANCES**Simultaneous:** vitamin A**Interfering:** cholecalciferol (vitamin D3)

KEY WORDS

work under subdued light; tablets; capsules; syrup

REFERENCE

Beaulieu,N.; Curran,N.M.; Gagné,C.; Gravelle,M.; Lovering,E.G. Liquid chromatographic methods for vitamins A and D in multivitamin-mineral formulations, *J.Assoc.Off.Anal.Chem.*, **1989**, 72, 247–254.

SAMPLE

Matrix: milk

Sample preparation: Mix 50 mL milk with 30 mL ethanolic KOH (10:30). Saponify the mixture at 80° for 20min. Extract the vitamins twice with 10 mL n-hexane. Evaporate to dryness, reconstitute the residue with 1 mL mobile phase, inject a 5 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5µm C18 (Alltech, CA, USA)

Mobile phase: MeOH:EtOH 80:20 (A) or EtOH:water 95:5 (B)

Flow rate: 1

Injection volume: 5

Detector: UV 250

CHROMATOGRAM

Retention time: 8 (A), 9 (B)

OTHER SUBSTANCES

Extracted: isotretinoin, retinal, tretinoin, vitamin A, vitamin D3, vitamin E, vitamin K1, vitamin K2

REFERENCE

Gong,B.Y.; Ho,J.W. Simultaneous separation and detection of ten common fat-soluble vitamins in milk, *J.Liq.Chromatogr.Rel.Technol.*, **1997**, 20, 2389–2397.

SAMPLE

Matrix: solutions

Sample preparation: Make up a solution in mobile phase, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 Zorbax PRO-10 C18

Mobile phase: MeOH:MeCN:hexane 95:3:2

Flow rate: 1

Injection volume: 20

Detector: UV 265

CHROMATOGRAM

Retention time: 8

OTHER SUBSTANCES

Simultaneous: cholecalciferol (vitamin D3), tachysterol

KEY WORDS

keep at 4° away from light; also preparative details

REFERENCE

Letter,W.S. Preparative isolation of vitamin D2 from previtamin D2 by recycle high-performance liquid chromatography, *J.Chromatogr.*, **1992**, 590, 169–173.