
REFERENCE

Pouliquen,H.; Fauconnet,V.; Morvan,M.-L.; Pinault,L. Determination of warfarin in the yolk and the white of hens' eggs by reversed-phase high-performance liquid chromatography, *J.Chromatogr.B*, **1997**, 702, 143–148.

SAMPLE

Matrix: microsomal incubations

Sample preparation: Add 250 μ L 1 M HCl to 500 μ L microsomal incubation, extract with 3 mL MTBE, evaporate the organic layer under nitrogen, reconstitute the residue in 100 μ L MeCN: water 50:50, inject an aliquot.

HPLC VARIABLES

Column: 150 \times 3.2 Novapak C18

Mobile phase: MeCN:buffer 35:65 (Buffer was 20 mM acetic acid adjusted to pH 4.8 with ammonium hydroxide.)

Flow rate: 1

Injection volume: 120

Detector: Radioactivity, Inus β -Ram using Inus Tru-Count scintillation fluid at a flow rate of 5 mL/min

CHROMATOGRAM

Retention time: 3.5

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

human; liver

REFERENCE

Obach,R.S. Nonspecific binding to microsomes: Impact on scale-up of in vitro intrinsic clearance to hepatic clearance as assessed through examination of warfarin, imipramine, and propranolol, *Drug Metab.Dispos.*, **1997**, 25, 1359–1369.

SAMPLE

Matrix: solutions

Sample preparation: Filter (0.45 μ m), dilute the filtrate with mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 5 μ m Hypersil ODS

Mobile phase: MeCN:10 mM pH 4.7 acetate buffer 50:50

Detector: UV 214

REFERENCE

Okimoto,K.; Rajewski,R.A.; Uekama,K.; Jona,J.A.; Stella,V.J. The interaction of charged and uncharged drugs with neutral (HP- β -CD) and anionically charged (SBE7- β -CD) β -cyclodextrins, *Pharm.Res.*, **1996**, 13, 256–264.

Xamoterol

Molecular formula: C₁₆H₂₅N₃O₅


Molecular weight: 339.39

CAS Registry No.: 81801-12-9, 73210-73-8 (hemifumarate)

Merck Index: 10189

SAMPLE

Matrix: blood



Sample preparation: Condition a Sep-Pak C18 SPE cartridge with MeOH and water. 1 mL Plasma + 2 mL pH 7.0 phosphate-citrate buffer, add to the SPE cartridge, wash with 10 mL water, wash with 5 mL MeOH:water 10:90, elute with 2 mL MeOH. Evaporate the eluate to dryness at 60°, reconstitute with 200 µL mobile phase, inject a 100 µL aliquot.

HPLC VARIABLES

Column: 150 × 4.1 5 µm Hypersil-ODS

Mobile phase: MeOH:THF:30 mM perchloric acid 8:0.6:91.4

Flow rate: 1.5

Injection volume: 100

Detector: E, Bioanalytical Systems Model LC-4A, TL-5A thin-layer glassy carbon electrode + 0.85 V, RE-1 Ag/AgCl reference electrode

CHROMATOGRAM

Retention time: 8

Limit of detection: 2 ng/mL

OTHER SUBSTANCES

Noninterfering: acetaminophen, acetazolamide, caffeine, chlordiazepoxide, cimetidine, codeine, diazepam, digoxin, docusate sodium, flurazepam, furosemide, heparin, hydrochlorothiazide, meclufenamate, methyl dopa, nitroglycerin, quinidine, simethicone, spironolactone, thioridazine, triamterene, warfarin

Interfering: aspirin, benzyl alcohol, hydralazine

KEY WORDS

plasma; SPE

REFERENCE

Davis, P.C. Determination of xamoterol in human plasma by high-performance liquid chromatography with electrochemical detection, *J. Chromatogr.*, **1987**, 417, 233–235.

SAMPLE

Matrix: blood, urine

Sample preparation: 1 mL Plasma or 100 µL urine + 40 ng prenatalol + 5 mL 0.1% disodium EDTA + 100 µL 1 M pH 6.5 ammonium acetate, add to the ion-exchange column, wash with two 10 mL portions of water, elute with 3 mL 1 M ammonium hydroxide. Evaporate the eluate to dryness under reduced pressure at 50°, reconstitute with 250 µL ice-cold 100 mM perchloric acid, centrifuge at 4° at 10000 g for 10 min, inject a 100 µL aliquot of the supernatant. (Prepare an ion-exchange column by packing 50-100 mesh Bio-Rex 70 Na⁺ cation-exchange resin (Bio-Rad) into a 40 × 10 column, wash with 3 M HCl, 3 M NaOH, 3 M acetic acid, 1 M pH 6.5 ammonium acetate, and water.)

HPLC VARIABLES

Column: 250 × 4.6 5 µm Spherisorb ODS

Mobile phase: MeCN:10 mM perchloric acid 15:85

Flow rate: 2

Injection volume: 100

Detector: F ex 190 em 320-400 (filter)

CHROMATOGRAM

Retention time: 6.9

Internal standard: prenatalol (5.8)

Limit of detection: 10 ng/mL (urine), 1 ng/mL (plasma)

OTHER SUBSTANCES

Simultaneous: albuterol, fenoterol

Noninterfering: atropine, clonidine, diazepam, nitrazepam, prazosin, quinidine

KEY WORDS

plasma; SPE; pharmacokinetics

REFERENCE

Oddie, C.J.; Jackman, G.P.; Bobik, A. Measurement of xamoterol in plasma and urine by high-performance liquid chromatography, *J. Chromatogr.*, **1984**, 308, 370–375.

SAMPLE

Matrix: bulk

Sample preparation: Dissolve 10 μ mole compound (as free base or hydrochloride) in 500 μ L MeCN, add 250 μ L 5% sodium carbonate (for hydrochlorides only), add 500 μ L 100 mM reagent in MeCN, vortex for 1 min, heat at 60° for 2 h, add 100 μ mole L-proline, heat at 60° for 30 min. Remove a 100 μ L aliquot and dilute it with mobile phase, neutralize with acetic acid, inject a 10 μ L aliquot. Prepare the reagent ((R,R)-N-(3,5-dinitrobenzoyl)-2-aminocyclohexylisothiocyanate) as follows. Add 0.7 mL carbon disulfide to 6 mL (1R,2R)-(-)-1,2-diaminocyclohexane, 12 mL water, and 12 mL EtOH, heat the oil bath to 80°, add 2.8 mL carbon disulfide dropwise (making sure that the product does not start to precipitate), when addition is complete reflux for 1 h, acidify with 500 μ L 5 M HCl, reflux for 12 h, cool, filter, wash the solid with a little cold EtOH to give trans-4,5-tetramethyleneimidazolidine-2-thione as a white fluffy solid (mp 148–150°) (Tetrahedron 1993, 49, 4419). Stir 7.97 g 3,5-dinitrobenzoyl chloride in 30 mL dichloroethane at 50°, add a solution of 6 g trans-4,5-tetramethyleneimidazolidine-2-thione in 120 mL dichloroethane containing a catalytic amount of 4-(dimethylamino)pyridine over 15 min, reflux for 2 h, remove the crystals of (R,R)-N-(3,5-dinitrobenzoyl)-2-aminocyclohexylisothiocyanate by filtration, evaporate the filtrate to dryness and dissolve the residue in 60 mL dichloroethane, reflux for 16 h to obtain more (R,R)-N-(3,5-dinitrobenzoyl)-2-aminocyclohexylisothiocyanate (mp >250°, $[\alpha]_{546} = -133^\circ$ (c = 1) in MeCN).

HPLC VARIABLES

Column: 125 \times 4 5 μ m Lichrospher 60 RP Select B

Mobile phase: MeCN:20 mM ammonium acetate 55:45

Flow rate: 1

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: k' 1.13, k' 2.27 (enantiomers)

OTHER SUBSTANCES

Also analyzed: acebutolol, alprenolol, atenolol, carazolol, carvedilol, formoterol, methamphetamine, metipranolol, metoprolol, nifenanol, nitrilo atenolol, oxprenolol, pindolol, propranolol

KEY WORDS

derivatization; chiral

REFERENCE

Kleidermigg, O.P.; Posch, K.; Lindner, W. Synthesis and application of a new isothiocyanate as a chiral derivatizing agent for the indirect resolution of chiral amino alcohols and amines, *J. Chromatogr. A*, **1996**, 729, 33–42.


Xanomeline

Molecular formula: C₁₄H₂₃N₃OS

Molecular weight: 281.42

CAS Registry No.: 131986-45-3, 141064-23-5 (oxalate), 152854-19-8 (tartrate)

Merck Index: 10190



SAMPLE

Matrix: microsomal incubations

Sample preparation: Condition a 500 mg C18 Bond-Elut SPE cartridge with MeOH and water. Mix 2 mL microsomal incubation with 4.5 mL ice-cold MeOH, centrifuge at 4000 g for 10 min, dilute the supernatant with 15 mL water. Add a 2 mL aliquot to the SPE cartridge, wash with

2 mL water, elute with 2 mL MeOH:25% aqueous ammonia 96:4, evaporate the eluate to dryness under reduced pressure, reconstitute with mobile phase, inject an aliquot.

HPLC VARIABLES

Column: 100 × 4.6 3 µm ChromSpher Si (Chrompack)

Mobile phase: n-Heptane:2-propanol:water:25% ammonia 90:10:0.075:0.075 (A) or n-heptane:2-propanol:MeOH:25% ammonia 50:50:10:1 (B)

Detector: UV 295; MS, VG TRIO 1000, particle beam interface at 50°, helium at 25-30 psi, ion source 200°, positive ionization mode, electron current 150 µA, electron energy 70 eV

OTHER SUBSTANCES

Extracted: metabolites

KEY WORDS

normal phase; rat; liver; SPE

REFERENCE

Andersen, J.V.; Hansen, K.T. Normal-phase liquid chromatography-particle-beam mass spectrometry in drug metabolism studies of the dopamine receptor antagonist Odapipam and the muscarine M1 receptor agonist Xanomeline, *Xenobiotica*, **1997**, 27, 901-912.

SAMPLE

Matrix: perfusate

Sample preparation: Inject an aliquot of perfusate directly onto the column.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Spherisorb Phenyl

Mobile phase: MeCN:buffer 40:60 (Buffer was 30 mM tetramethylammonium hydroxide adjusted to pH 3.0 with perchloric acid.)

Column temperature: 35

Flow rate: 1

Injection volume: 50

Detector: UV 295

CHROMATOGRAM

Limit of detection: 100 nM

KEY WORDS

liver

REFERENCE

Andersen, J.V.; Hansen, K.T. Normal-phase liquid chromatography-particle-beam mass spectrometry in drug metabolism studies of the dopamine receptor antagonist Odapipam and the muscarine M1 receptor agonist Xanomeline, *Xenobiotica*, **1997**, 27, 901-912.

Xipamide


Molecular formula: C₁₅H₁₅ClN₂O₄S

Molecular weight: 354.81

CAS Registry No.: 14293-44-8

Merck Index: 10212

Lednicer No.: 2 93



SAMPLE

Matrix: blood

Sample preparation: 1 mL Plasma + 10 mL dichloromethane:2-propanol 75:25, shake for 10 min. Centrifuge at 2000 g for 10 min at 4°. Remove the organic phase and evaporate it to

dryness under a stream of nitrogen at 50°. Reconstitute the residue in 200 µL mobile phase, mix for 10 s. Centrifuge at 6500 g for 10 min. Inject a 40 µL aliquot of the supernatant.

HPLC VARIABLES

Guard column: 4 × 4 5 µm LiChrospher 100 RP-18

Column: 250 × 4 5 µm Supelcosil LC-18 (Supelco)

Mobile phase: n-Propanol:buffer 5:95 (Buffer was 50 mM sodium dodecyl sulfate in 10 mM pH 5.8 sodium phosphate buffer.)

Flow rate: 1.3

Injection volume: 40

Detector: UV 225

CHROMATOGRAM

Retention time: 8.58

Internal standard: xipamide

OTHER SUBSTANCES

Extracted: albuterol, atenolol, chlorthalidone

KEY WORDS

plasma; xipamide is IS

REFERENCE

Giachetti,C.; Tenconi,A.; Canali,S.; Zanolò,G. Simultaneous determination of atenolol and chlorthalidone in plasma by high-performance liquid chromatography. Application to pharmacokinetic studies in man, *J.Chromatogr.B*, **1997**, 698, 187–194.

SAMPLE

Matrix: blood, urine

Sample preparation: Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200–350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

HPLC VARIABLES

Guard column: 20 mm long Symmetry C18

Column: 250 × 4.6 5 µm Symmetry C8 (Waters)

Mobile phase: Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

Column temperature: 30

Flow rate: 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

Injection volume: 10–30

Detector: UV 218.1

CHROMATOGRAM

Retention time: 18.823


KEY WORDS

whole blood

REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149–163.

Xylazine



Molecular formula: C₁₂H₁₆N₂S

Molecular weight: 220.34

CAS Registry No.: 7361-61-7, 23076-35-9 (HCl)

Merck Index: 10213

Lednicer No.: 2 307

SAMPLE

Matrix: blood

Sample preparation: 500 μ L Plasma + 10 μ L 100 μ g/mL doxapram in MeOH + 1 mL 50 mM pH 11 borax buffer, vortex for 5 s, add 10 mL chloroform, shake for 10 min, centrifuge at 11400 g for 10 min, filter (Whatman No. 1 PS phase-separating paper). Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 70°, reconstitute the residue in 100 μ L mobile phase, inject the whole sample.

HPLC VARIABLES

Column: 300 \times 4 10 μ m μ Bondapak C18

Mobile phase: MeOH:water:heptanesulfonic acid (Pic B7) 45:55:0.2 containing 2% glacial acetic acid

Flow rate: 2

Injection volume: 100

Detector: UV 225

CHROMATOGRAM

Retention time: 4

Internal standard: doxapram (5.5)

Limit of detection: 20 ng/mL

KEY WORDS

plasma; sheep

REFERENCE

Alvinerie,M.; Toutain,P.L. Determination of xylazine in plasma using high-performance liquid chromatography, *J.Chromatogr.*, **1981**, 222, 308–310.

SAMPLE

Matrix: blood

Sample preparation: Condition a Sep-Pak C18 SPE cartridge with water, MeOH, and 100 mM ammonium acetate. Add 200 μ L plasma to the SPE cartridge, wash with 100 mM ammonium acetate, elute with MeOH:100 mM ammonium acetate 3:1. Evaporate the eluate to dryness under reduced pressure, dissolve the residue in 200 μ L mobile phase, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 150 \times 4.6 Hitachi gel 3056 octadecylsilica

Mobile phase: MeOH:100 mM ammonium acetate 60:40

Flow rate: 1

Injection volume: 20

Detector: MS, Hitachi M1000, APCI, nebulizer 260°, vaporizer 399

CHROMATOGRAM

Retention time: 7.0

Limit of detection: 0.5-2.5 ng/mL

OTHER SUBSTANCES

Simultaneous: atipamezole, atropine, butorphanol, flumazenil, ketamine, medetomidine, midazolam

KEY WORDS

plasma; SPE; dog

REFERENCE

Kanazawa,H.; Nagata,Y.; Matsushima,Y.; Takai,N.; Uchiyama,H.; Nishimura,R.; Takeuchi,A. Liquid chromatography-mass spectrometry for the determination of medetomidine and other anaesthetics in plasma, *J.Chromatogr.*, **1993**, 631, 215–220.

SAMPLE**Matrix:** solutions**Sample preparation:** Prepare a 1-10 µg/mL solution in water, inject an aliquot.

HPLC VARIABLES**Column:** 250 × 4.6 5 µm Hypersil SCX/C18**Mobile phase:** MeCN:25 mM pH 3 Na₂HPO₄ 50:50**Injection volume:** 20**Detector:** UV 254

CHROMATOGRAM**Retention time:** k' 4.53

OTHER SUBSTANCES

Also analyzed: amitriptyline, barbital, benzoic acid, butabarbital, clomipramine, clonazepam, desipramine, diazepam, flurazepam, furosemide, imipramine, nitrazepam, phenobarbital, phenol, phenolphthalein, pindolol, propranolol, resorcinol, salicylic acid, secobarbital, terbutaline

KEY WORDS

effect of mobile phase pH on capacity factor is discussed

REFERENCE

Walshe,M.; Kelly,M.T.; Smyth,M.R.; Ritchie,H. Retention studies on mixed-mode columns in high-performance liquid chromatography, *J.Chromatogr.A*, **1995**, 708, 31–40.

SAMPLE**Matrix:** tissue

Sample preparation: Condition a Sep-Pak C18 SPE cartridge with 5 mL MeOH and 5 mL water. Homogenize kidney with a kitchen grinder. Weigh out a 5 g sample and add 20 mL MeCN with continuous gentle mixing, mix vigorously on a vibromixer at 1500 rpm for 30 s, sonicate for 2 min, centrifuge at 4000 g for 5 min. Mix 7.5 mL sample extract and 40 mL 10% NaCl and add to SPE cartridge, wash with 1 mL 10 mM sulfuric acid, wash with 2 mL air, elute with 2 mL acidic MeCN. Place eluate in a washed tube and evaporate to 300 µL at 70° under a stream of nitrogen, mix gently, add 1 mL n-hexane, mix on a vibromixer for 30 s, centrifuge at 2000 g, inject a 50 µL aliquot of the aqueous phase. (Acidic MeCN was 1 mL 50 mM sulfuric acid and 100 mL MeCN. The washed tube was prepared by rinsing with concentrated ammonia, water, and acetone and drying under a stream of nitrogen.)

HPLC VARIABLES**Guard column:** 10 × 2.1 37-50 µm Bondapak C18**Column:** 300 × 3.9 Bondapak C18**Mobile phase:** MeCN:water 55:45 containing 2.46 g/L anhydrous sodium acetate, pH adjusted to 6.5 with acetic acid**Flow rate:** 1.2**Injection volume:** 50**Detector:** UV 240

CHROMATOGRAM**Retention time:** 6.5**Limit of detection:** 4 ng/g

OTHER SUBSTANCES

Extracted: azaperol, carazolol, acepromazine, azaperone, haloperidol, propiomazine, chlorpromazine

KEY WORDS

SPE; pig; kidney

REFERENCE

Keukens,H.J.; Aerts,M.M.L. Determination of residues of carazolol and a number of tranquilizers in swine kidney by high-performance liquid chromatography with ultraviolet and fluorescence detection, *J.Chromatogr.*, **1989**, *464*, 149–161.

SAMPLE

Matrix: tissue

Sample preparation: Condition a Bond-Elut C18 SPE cartridge with 5 mL MeOH and 5 mL water. Cut pig kidney or liver into small pieces and homogenize. 5 g Homogenate + 10 mL MeCN, shake, vortex for 30 s, sonicate for 3 min, vortex for 30 s, sonicate for 3 min, centrifuge at 10000 g for 20 min. Add 7.5 mL supernatant + 40 mL 10% NaCl to the SPE cartridge at about 1 mL/min, do not allow cartridge to dry out, wash with 850 μ L 10 mM sulfuric acid, dry with air, elute with 3.5 mL acidic MeCN. Evaporate the eluate to dryness under a stream of nitrogen at 50°, reconstitute the residue in 300 μ L 10 mM sulfuric acid, vortex briefly, add 1 mL hexane, vortex for 30 s, centrifuge at 2000 g for 5 min, inject an aliquot of the aqueous layer. (Acidic MeCN was 1 mL 50 mM sulfuric acid in 100 mL MeCN.)

HPLC VARIABLES

Guard column: Hypersil 5 μ m SAS C1

Column: 250 mm long 5 μ m Hypersil SAS C1

Mobile phase: MeCN:water 50:50 containing 0.77 g/L ammonium acetate

Flow rate: 2

Detector: E, ESA Model 5100A Coulochem, first electrode +0.4 V, second electrode (which was monitored) +0.7 V, Model 5020 guard cell after pump but before injector at +0.75 V

CHROMATOGRAM

Retention time: 10

Limit of detection: 2 ng/g

OTHER SUBSTANCES

Extracted: azaperol, acepromazine, carazolol, azaperone, haloperidol, propiomazine, chlorpromazine

KEY WORDS

SPE; pig; kidney; liver

REFERENCE

Rose,M.D.; Shearer,G. Determination of tranquilisers and carazolol residues in animal tissue using high-performance liquid chromatography with electrochemical detection, *J.Chromatogr.*, **1992**, *624*, 471–477.

SAMPLE

Matrix: tissue

Sample preparation: Homogenize (Tissuemizer) 20 g kidney and 70 mL chloroform at medium-high speed for 1-1.5 min, rinse blade with 5 mL chloroform, centrifuge at 1800 rpm for 10 min, filter (glass-fiber) supernatant, repeat extraction, rinse filter with 30 mL chloroform. Combine filtrates, add 1 mL 1 M HCl, evaporate to dryness under reduced pressure at 55°, add 25 mL MeOH, evaporate to dryness, reconstitute with 25 mL petroleum ether, add to Celite column, rinse flask with three 10 mL portions of petroleum ether, add rinses to column, rinse flask with two 50 mL portions of petroleum ether:MeOH 98.5:1.5, add to column, discard all effluent, elute with two 50 mL portions of MeOH. Evaporate eluate to dryness under reduced pressure at 55°, reconstitute with 25 mL MeOH, evaporate to dryness, reconstitute with 4-10 mL mobile phase, filter (0.45 μ m), inject a 50 μ L aliquot. (Petroleum ether was water saturated. Prepare Celite column as follows. Slurry 250 g Celite 545 and 800 mL HCl:water 50:50, heat on a steam bath with occasional stirring for several hours, allow to settle, pour off liquid, add 800 mL HCl:water 50:50, slurry, heat on a steam bath for several hours, allow to settle, decant liquid, wash with water until pH is neutral, wash with 250 mL MeOH, wash with 250 mL n-hexane, wash with 250 mL petroleum ether, heat on a steam bath to remove residual solvent, dry at 105°. Blend 2 g acid-washed Celite and 0.5 mL water, add to a 300 \times 22 column, blend 3 g Celite and 1 mL 1 M HCl, add to the column.)

HPLC VARIABLES**Column:** 300 × 3.9 10 μm μBondapak phenyl**Mobile phase:** MeCN:water:2 M sodium acetate:1 M acetic acid 32:64:2:2**Flow rate:** 1**Injection volume:** 50**Detector:** UV 225**CHROMATOGRAM****Retention time:** 6**OTHER SUBSTANCES****Extracted:** metabolites**KEY WORDS**

cow; pig; kidney; SPE

REFERENCE

Holland,D.C.; Munns,R.K.; Roybal,J.E.; Hurlbut,J.A.; Long,A.R. Simultaneous determination of xylazine and its major metabolite, 2,6-dimethylaniline, in bovine and swine kidney by liquid chromatography, *J.AOAC Int.*, **1993**, 76, 720–724.

SAMPLE**Matrix:** urine

Sample preparation: Condition a 1 mL 100 mg Bond-Elut cyanopropyl SPE cartridge with two 1 mL portions of MeOH and 1 mL water, do not allow to dry. 1 mL Urine + 500 μL water + 100 μL 15 μg/mL diazepam in MeOH, add to the SPE cartridge, dry under vacuum for 3 min, elute with two 250 μL portions of MeCN:methanolic HCl 50:50, inject an aliquot of the eluate. (Prepare methanolic HCl by adding 3 mL concentrated HCl to 50 mL MeOH.)


HPLC VARIABLES**Column:** 250 × 4.6 5 μm Hypersil C18**Mobile phase:** MeCN:MeOH:buffer 30:2:50 (Buffer was 4 g/L tetramethylammonium hydroxide.)**Flow rate:** 2**Injection volume:** 20**Detector:** UV 225**CHROMATOGRAM****Retention time:** 6.5**Internal standard:** diazepam (9.4)**Limit of detection:** 10 ng/mL**KEY WORDS**

dog; SPE; pharmacokinetics

REFERENCE

Moore,C.M.; Oliver,J.S. Rapid extraction and determination of xylazine in greyhound urine using high-performance liquid chromatography, *J.Chromatogr.*, **1989**, 491, 519–524.

Xylometazoline

Molecular formula: C₁₆H₂₄N₂**Molecular weight:** 244.38**CAS Registry No.:** 526-36-3, 1218-35-5 (HCl)**Merck Index:** 10219**Lednicer No.:** 1 242

SAMPLE

Matrix: formulations

HPLC VARIABLES

Column: 300 × 3.9 10 µm LiChrosorb Si-60

Mobile phase: MeOH:water 60:40 containing 4 mM disodium citrate and 4 mM tetrabutylammonium bromide, pH 6.0

Flow rate: 1

Injection volume: 10

Detector: UV 254

CHROMATOGRAM

Retention time: 9

OTHER SUBSTANCES

Simultaneous: benzalkonium chloride, domiphen bromide, thimerosal

KEY WORDS

nasal drops

REFERENCE

Lingeman,H.; van Munster,H.A.; Beynen,J.H.; Underberg,W.J.; Hulshoff,A. High-performance liquid chromatographic analysis of basic compounds on non-modified silica gel and aluminium oxide with aqueous solvent mixtures, *J.Chromatogr.*, **1986**, 352, 261–274.

SAMPLE

Matrix: formulations

Sample preparation: Dilute nasal solution 10-fold with water, filter (0.45 µm), inject a 10 µL aliquot of the filtrate.

HPLC VARIABLES

Column: 125 × 4 5 µm Aluspher RP-select B (Merck)

Mobile phase: Gradient. MeCN:1 mM NaOH from 10:90 to 80:20 over 25 min.

Column temperature: 25

Flow rate: 1.2

Injection volume: 10

Detector: UV 224

CHROMATOGRAM

Retention time: 14.5

Limit of detection: 1 ng

OTHER SUBSTANCES

Simultaneous: ephedrine, naphazoline, oxymetazoline

KEY WORDS

nasal solutions

REFERENCE

De Orsi,D.; Gagliardi,L.; Cavazzutti,G.; Mediati,M.G.; Tonelli,D. Simultaneous determination of ephedrine and 2-imidazolines in pharmaceutical formulations by reversed-phase HPLC, *J.Liq.Chromatogr.*, **1995**, 18, 3233–3242.

SAMPLE

Matrix: solutions

Sample preparation: Prepare a 10 µg/mL solution in MeOH, inject a 20 µL aliquot.

HPLC VARIABLES

Column: 125 × 4.9 Spherisorb S5W silica

Mobile phase: MeOH containing 10 mM ammonium perchlorate and 1 mL/L 100 mM NaOH in MeOH, pH 6.7

Flow rate: 2

Injection volume: 20

Detector: E, LeCarbone, V25 glassy carbon electrode, + 1.2 V

CHROMATOGRAM

Retention time: 2.3

OTHER SUBSTANCES

Also analyzed: acebutolol, acepromazine, acetophenazine, N-acetylprocainamide, albuterol, alprenolol, amethocaine, amiodarone, amitriptyline, antazoline, atenolol, azacyclonal, bamethan, benactyzine, benperidol, benzethidine, benzocaine, benzocetamine, benzphetamine, benzquinamide, bromhexine, bromodiphenhydramine, bromperidol, brompheniramine, brompromazine, buclizine, bufotenine, bupivacaine, buprenorphine, butacaine, butethamate, chlorcyclizine, chlorpheniramine, chlorphenoxamine, chlorprenaline, chlorpromazine, chlorprothixene, cimetidine, cinchonidine, cinnarizine, clemastine, clomipramine, clonidine, cocaine, cyclazocine, cyclizine, cyclopentamine, cyproheptadine, deserpidine, desipramine, dextromoramide, dextropropoxyphene, dicyclomine, diethylcarbamazepine, diethylpropion, diethylthiambutene, dihydroergotamine, dimethindene, dimethothiazine, diphenhydramine, diphenoxylate, dipiphenone, diprenorphine, dipyrizamide, disopyramide, dothiepin, doxapram, doxepin, doxylamine, droperidol, ephedrine, ergocornine, ergocristine, ergocristinine, ergocryptine, ergometrine, ergosine, ergosinine, ergotamine, ethopropazine, etorphine, etoxeridine, fenethazine, fenfluramine, fenoterol, fentanyl, flavoxate, flupromazine, flupenthixol, fluphenazine, flurazepam, haloperidol, hydroxyzine, hyoscine, ibogaine, imipramine, indapamine, iprindole, isothipendyl, isoxsuprine, ketanserin, laudanosine, lidocaine, lofepramine, loxapine, maprotiline, mecamlamine, meclorphenoxate, meclozine, medazepam, mephentermine, mepivacaine, meptazinol, mepyramine, mesoridazine, metaraminol, methadone, methamphetamine, methapyrilene, methdilazene, methotrimeprazine, methoxamine, methoxyphenamine, methoxypropazine, methylephedrine, methylergonovine, methysergide, metoclopramide, metopimazine, metoprolol, mianserin, morazone, nadolol, nalorphine, naloxone, naphazoline, nicotine, nifedipine, nomifensine, nortriptyline, noscapine, orphenadrine, oxeladin, oxprenolol, oxymetazolin, papaverine, pargyline, pecazine, penbutolol, pentazocine, penthienate, pericyazine, perphenazine, phenadoxone, phenampromide, phenazocine, phenbutrazate, phendimetrazine, phenelzine, phenglutarimide, phenindamine, pheniramine, phenmetrazine, phenomorphan, phenoperidine, phenothiazine, phenoxybenzamine, phentolamine, phenylephrine, phenyltoloxamine, physostigmine, pimindoline, pimizole, pindolol, pipamazine, pipazethate, piperacetazine, piperidolate, pipradol, pirenzepine, piritramide, pizotifen, practolol, pramoxine, prazosin, prenylamine, prilocaine, primaquine, proadifen, procainamide, procaine, prochlorperazine, procyclidine, proheptazine, prolintane, promazine, promethazine, pronethalol, properidine, propiomazine, proprenolol, prothipendyl, protriptyline, proxymetacaine, pseudoephedrine, pyrimethamine, quinidine, quinine, ranitidine, rescinnamine, sotalol, tacrine, terazosin, terbutaline, terfenadine, thenyldiamine, theophylline, thiethylperazine, thiopropazate, thioproperazine, thioridazine, thiothixene, thonzylamine, timolol, tocanide, tolpropamine, tolycaine, tranlylcypromine, trazodone, trifluoperazine, trifluoperidol, trimeperidine, trimeprazine, trimethobenzamide, trimethoprim, trimipramine, tripeleminamine, triprolidine, tryptamine, verapamil

REFERENCE

Jane, I.; McKinnon, A.; Flanagan, R.J. High-performance liquid chromatographic analysis of basic drugs on silica columns using non-aqueous ionic eluents. II. Application of UV, fluorescence and electrochemical oxidation detection, *J.Chromatogr.*, **1985**, *323*, 191-225.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 300 × 3.9 10 μm LiChrosorb Si-60

Mobile phase: MeOH:water 60:40 containing 4 mM disodium citrate and 4 mM tetrabutylammonium bromide, pH 5.9

Flow rate: 1

Injection volume: 10

Detector: UV 254

CHROMATOGRAM**Retention time:** 7

OTHER SUBSTANCES**Simultaneous:** atropine, codeine, dansylamide, dansylcadaverine, doxorubicin, methylatropine, naphazoline, noscapine

REFERENCE

Lingeman,H.; van Munster,H.A.; Beynen,J.H.; Underberg,W.J.; Hulshoff,A. High-performance liquid chromatographic analysis of basic compounds on non-modified silica gel and aluminium oxide with aqueous solvent mixtures, *J.Chromatogr.*, **1986**, 352, 261–274.

SAMPLE**Matrix:** solutions**Sample preparation:** Prepare a solution in MeOH:water 40:60, inject a 10 μ L aliquot.

HPLC VARIABLES**Column:** 150 \times 4.1 RSIL C18 (RSL, Eke, Belgium)**Mobile phase:** MeOH:water 40:60 containing 20 mM sodium 1-octanesulfonate and 10 mM N,N-dimethyloctylamine, pH adjusted to 3.0 with orthophosphoric acid**Column temperature:** 25**Flow rate:** 1**Injection volume:** 10**Detector:** UV 220

CHROMATOGRAM**Retention time:** 40

OTHER SUBSTANCES**Simultaneous:** degradation products, antazoline, coumazoline, lidocaine, naphazoline, oxy-metazoline, prednisolone, sulfadimidine, sulfanilamide, sulfathiazole, tenaphtoxaline, tetrahydrozoline, tolazoline, tramazoline

REFERENCE

De Schutter,J.A.; Van den Bossche,W.; De Moerloose,P. Stability-indicating analysis of tetrazyline hydrochloride in pharmaceutical formulations by reversed-phase ion-pair liquid chromatography, *J.Chromatogr.*, **1987**, 391, 303–308.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 250 \times 4.6 5 μ m Supelcosil LC-DP (A) or 250 \times 4 5 μ m LiChrospher 100 RP-8 (B)**Mobile phase:** MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)**Flow rate:** 0.6**Injection volume:** 25**Detector:** UV 229

CHROMATOGRAM**Retention time:** 12.68 (A), 7.16 (B)

OTHER SUBSTANCES**Also analyzed:** acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordi-azepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene,

desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyridamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazinol, mefenamic acid, meperidine, mephenytoin, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfapyrazole, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocainide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, yohimbine, zopiclone


KEY WORDS

details of plasma extraction

REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, 692, 103–119.

Xylose



Molecular formula: $C_5H_{10}O_5$

Molecular weight: 150.13

CAS Registry No.: 58-86-6

Merck Index: 10220

SAMPLE

Matrix: beverages, juice, milk

Sample preparation: Orange juice. Dilute orange juice 100-fold with water, filter (Millipore HV, 0.45 μ m), dilute filtrate 10-fold, inject an aliquot. Beverages. Dilute soft drinks 1000-fold with water, inject an aliquot. Milk. Dilute 5 mL milk to 100 mL with mobile phase, filter (Millipore HV, 0.45 μ m), dilute filtrate 50-fold, inject an aliquot.

HPLC VARIABLES

Guard column: 30 \times 4.6 Cation H (Bio-Rad)

Column: 300 \times 3.8 9 μ m HPX 87-H Aminex (Bio-Rad)

Mobile phase: 10 mM Sulfuric acid

Column temperature: 50

Flow rate: 0.5

Injection volume: 40

Detector: E following post-column reaction, Hewlett-Packard 1049A programmable electrochemical detector, Metrohm detector cell, cuprous oxide working electrode +550 mV, glassy carbon auxiliary electrode, Ag/AgCl (3 M KCl) reference electrode. The column effluent mixed with 200 mM NaOH pumped at 0.4 mL/min, the mixture flowed through a 220 \times 0.8 single-bead string reactor packed with 0.6 mm glass beads to the detector. (Prepare cuprous oxide electrode as follows. Stir 300 mg conductive carbon cement (Gerhard Neubauer, Munster), 60 mg cuprous oxide (Fluka), and 300 μ L acetone until a thick paste forms as the acetone evaporates. Pack

conductive carbon cement into the base of a 3 mm diameter cavity carbon paste electrode base (Metrohm), allow to dry, polish with dry emery paper (grade 2/0, Oakey), remove surface layer with an acetone-soaked tissue, pack the paste into the cavity, allow to dry overnight, polish with dry emery paper (grade 2/0), 3 μm imperial micro finishing film sheet (3M), 0.3 μm imperial micro finishing film sheet (3M), and 0.05 μm alumina particles on a Buehler pad, sonicate for 2 min in water (Anal. Chim. Acta 1995, 300, 5.)

CHROMATOGRAM

Retention time: 11.25

Limit of detection: 0.9 μM

OTHER SUBSTANCES

Also analyzed: arabinose, cellobiose, dextrose, fructose, fucose, galactitol, galactose, galacturonic acid, lactose, lactulose, lyxose, maltose, mannitol, mannose, myo-inositol, raffinose, rhamnose, ribose, sorbose, sucrose

KEY WORDS

orange juice; soft drinks; post-column reaction; fruit

REFERENCE

Huang,X.; Pot,J.J.; Kok,W.T. Determination of sugars by liquid chromatography and amperometric detection with a cuprous oxide modified electrode, *Chromatographia*, **1995**, *40*, 684–689.

SAMPLE

Matrix: blood

Sample preparation: 100 μL Serum + 500 μL MeOH, shake for 1 min, centrifuge at 10000 rpm for 1 min, inject a 10 μL aliquot of the supernatant.

HPLC VARIABLES

Column: 300 \times 4 Aminex A-27

Mobile phase: 500 mM Boric acid adjusted to pH 8.7 with KOH

Flow rate: 2

Injection volume: 10

Detector: F ex 357 (low-pressure mercury lamp) em 436 (420 nm cutoff filter) following post-column reaction. The column effluent mixed with the reagent pumped at 0.5 mL/min and the mixture flowed through a 10 m \times 0.8 mm ID stainless steel coil at 150° to the detector. (Reagent was 20 g boric acid and 20 g ethanolamine in 1 L water.)

CHROMATOGRAM

Retention time: 20

Limit of quantitation: 5 nmoles

OTHER SUBSTANCES

Extracted: dextrose, fructose, galactose, maltose, ribose

KEY WORDS

post-column reaction; serum

REFERENCE

Kato,T.; Kinoshita,T. Fluorometric detection and determination of carbohydrates by high-performance liquid chromatography using ethanolamine, *Anal.Biochem.*, **1980**, *106*, 238–243.

SAMPLE

Matrix: bulk

Sample preparation: Evaporate hydrolysates of glycosaminoglycans to dryness, reconstitute in 500 μL 10% benzoic anhydride in pyridine containing 5% 4-dimethylaminopyridine, heat at 37° for 1.5 h, add 4.5 mL water, shake vigorously, pass through a Sep-Pak C18 SPE cartridge three times, wash with 10 mL pyridine:water 10:90, wash with 5 mL water, reverse the direction of flow and elute with 2.5 mL MeCN, evaporate the eluate to dryness, reconstitute with MeCN, centrifuge at 11000 g for 5 min, inject an aliquot.

HPLC VARIABLES**Guard column:** 30 × 4.6 RP-18**Column:** 250 × 4.6 Supelcosil LC-18**Mobile phase:** MeCN:water 75:25**Flow rate:** 1**Injection volume:** 20**Detector:** UV 230

CHROMATOGRAM**Retention time:** 9.4

OTHER SUBSTANCES**Simultaneous:** N-acetylgalactosamine, N-acetylglucosamine, 1,6-anhydroidose, dextrose, fucose, galactosamine, galactose, glucosamine, mannose, 1-methylfucose, 1-methylgalactose, 1-methylmannose, 1-methylxylose**Interfering:** 1-methylglucose

KEY WORDSderivatization; SPE

REFERENCE

Karamanos,N.K.; Hjerpe,A.; Tsegendis,T.; Engfeldt,B.; Antonopoulos,C.A. Determination of iduronic acid and glucuronic acid in glycosaminoglycans after stoichiometric reduction and depolymerization using high-performance liquid chromatography and ultraviolet detection, *Anal.Biochem.*, **1988**, 172, 410–419.

SAMPLE**Matrix:** carbohydrates**Sample preparation:** Mix 10 nmoles total monosaccharides with 200 μ L 2 M trifluoroacetic acid, flush with nitrogen for a few min, seal, heat at 100° for 6 h, evaporate to dryness under reduced pressure in a desiccator over NaOH pellets, reconstitute with water, inject an aliquot.

HPLC VARIABLES**Column:** 80 × 8 11 μ m Hitachi No. 2633 resin (quaternary ammonium)**Mobile phase:** Gradient. A was 250 mM pH 8.2 borate buffer. B was 400 mM pH 7.4 borate buffer. C was 600 mM pH 9.3 borate buffer. A:B:C from 100:0:0 to 0:100:0 over 15 min, maintain at 0:100:0 for 20 min, to 0:0:100 over 10 min, maintain at 0:0:100.**Column temperature:** 65**Flow rate:** 1**Injection volume:** 20**Detector:** F ex 331 em 383 following post-column reaction. The column effluent mixed with 10% 2-cyanoacetamide in water pumped at 0.25 mL/min and 600 mM pH 9.3 borate buffer pumped at 0.25 mL/min and the mixture flowed through a 10 m × 0.5 mm ID PTFE coil at 100 ± 0.5° to the detector.

CHROMATOGRAM**Retention time:** 58**Limit of detection:** 0.1-1 nmole

OTHER SUBSTANCES**Simultaneous:** arabinose, dextrose, fucose, galactose, lyxose, mannose, rhamnose, ribose

KEY WORDSpost-column reaction

REFERENCE

Honda,S.; Takahashi,M.; Kakehi,K.; Ganno,S. Rapid, automated analysis of monosaccharides by high-performance anion-exchange chromatography of borate complexes with fluorimetric detection using 2-cyanoacetamide, *Anal.Biochem.*, **1981**, 113, 130–138.

SAMPLE**Matrix:** glycoconjugates

Sample preparation: Mix 0.1-1.5 mg glycoconjugate with 200 μ L 2 M trifluoroacetic acid, flush with nitrogen for a few min, seal, heat at 100° for 6 h, evaporate to dryness under reduced pressure in a desiccator over NaOH pellets, reconstitute with 200 μ L water, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: Hitachi No. 2633 resin

Mobile phase: Gradient. A was 250 mM pH 8.2 borate buffer. B was 400 mM pH 7.4 borate buffer. C was 600 mM pH 9.3 borate buffer. A:B:C from 100:0:0 to 0:100:0 over 15 min, maintain at 0:100:0 for 20 min, to 0:0:100 over 11 min, maintain at 0:0:100.

Column temperature: 65 \pm 1

Flow rate: 1

Injection volume: 20

Detector: UV 276 following post-column reaction. The column effluent mixed with 1% 2-cyanoacetamide pumped at 0.5 mL/min and 600 mM pH 10.5 borate buffer pumped at 0.5 mL/min and the mixture flowed through a 10 m \times 0.5 mm ID PTFE coil at 100 \pm 0.2° and a 1 m \times 0.5 mm ID PTFE cooling coil to the detector.

CHROMATOGRAM

Retention time: 57

Limit of detection: 1 nmole

OTHER SUBSTANCES

Extracted: arabinose, dextrose, fucose, galactose, lyxose, mannose, rhamnose, ribose

KEY WORDS

post-column reaction

REFERENCE

Honda,S.; Takahashi,M.; Nishimura,Y.; Kakehi,K.; Ganno,S. Sensitive ultraviolet monitoring of aldoses in automated borate complex anion-exchange chromatography with 2-cyanoacetamide, *Anal.Biochem.*, **1981**, *118*, 162-167.

SAMPLE

Matrix: glycoproteins

Sample preparation: 200 μ g Glycoprotein + 100 μ L water + 100 μ L 4 M trifluoroacetic acid, heat at 100° for 6 h, cool to room temperature, evaporate to dryness under reduced pressure at 35°, add 40 μ L reagent, heat at 80° for 1 h, cool to room temperature, add 200 μ L water, add 200 μ L chloroform, vortex vigorously, centrifuge for 1 min, inject an aliquot of the upper aqueous layer. (Prepare the reagent by mixing 165 mg ethyl p-aminobenzoate, 35 mg sodium cyanoborohydride, 41 μ L glacial acetic acid, and 350 μ L glacial acetic acid.)

HPLC VARIABLES

Column: 150 \times 3.9 Pico.Tag (Waters)

Mobile phase: MeCN:MeOH:50 mM pH 4.5 sodium acetate 10:5:85

Column temperature: 45

Flow rate: 1.2

Detector: UV 254

CHROMATOGRAM

Retention time: 8

OTHER SUBSTANCES

Extracted: N-acetylgalactosamine, N-acetylglucosamine, 2-deoxyglucose, fucose, galactosamine, galactose, glucosamine, glucose, lactose, maltose, mannose

KEY WORDS

derivatization

REFERENCE

Kwon,H.; Kim,J. Determination of monosaccharides in glycoproteins by reverse-phase high-performance liquid chromatography, *Anal.Biochem.*, **1993**, *215*, 243-252.

SAMPLE**Matrix:** plants**Sample preparation:** Condition trimethylaminopropylsilica SAX and cyclohexylsilica SPE cartridges (Analytichem) with 4 mL MeOH and 4 mL water. Heat 1 g plant material with 10 mL EtOH:water 80:20 in a sealed tube at 100° for 15-30 min, evaporate the extract to dryness, reconstitute with water, pass through the SPE cartridges, inject a 50 µL aliquot of the eluate.

HPLC VARIABLES**Column:** 300 × 6.5 Sugar Pak-1 microparticulate gel, calcium form (Waters)**Mobile phase:** 100 µM Calcium EDTA**Column temperature:** 70**Flow rate:** 0.4**Injection volume:** 50**Detector:** F ex 360 em 470 following post-column reaction. The column effluent mixed with 30 mM benzamidine in 1 M KOH pumped at 1 mL/min and the mixture flowed through a 530 µL reaction coil (Varian PCR1) at 100° to the detector.

CHROMATOGRAM**Retention time:** 10.92**Limit of detection:** 15.8-62.5

OTHER SUBSTANCES**Extracted:** arabinose, dextrose, fructose, fucose, galactose, lactose, mannose**Interfering:** rhamnose

KEY WORDS

post-column reaction; SPE

REFERENCECoquet,A.; Veuthey,J.-L.; Haerdi,W. Selective post-column fluorogenic reaction with benzamidine for trace level detection of reducing saccharides in liquid chromatography, *J.Chromatogr.*, **1991**, 553, 255-263.

SAMPLE**Matrix:** solutions

HPLC VARIABLES**Column:** 160 × 4 8 µm DAX8 anion-exchange resin, sulfate form (Durrum Chemical Co.) (Regenerate resin outside the column by washing 10 g resin with 400 mL water, 400 mL 500 mM NaCl at 50°, water, 0.5 N sodium sulfate at 50° (until a negative chloride test is obtained), water, and EtOH:water 95:5. Slurry pack below 70° with mobile phase at 1.14 mL/min.)**Mobile phase:** EtOH:water 87.6:12.4**Column temperature:** 88**Flow rate:** 0.57**Injection volume:** 5-100**Detector:** UV 562 following post-column reaction. The column effluent mixed with the reagent pumped at 0.3 mL/min and the mixture flowed for 5 min through a coil of 0.3 mm ID PTFE at 100° to the detector. (Prepare reagent by mixing equal volumes of solution A and solution B, the mixture is stable for at least 1 month. Solution A is 1 g of copper sulfate pentahydrate and 3.7 g aspartic acid in 1 L water. Solution B is 38 g sodium carbonate decahydrate and 2 g sodium bicinchoninate (Pierce Chemical Co.) in 1 L water.)

CHROMATOGRAM**Retention time:** 83**Limit of detection:** <500 pmole

OTHER SUBSTANCES**Simultaneous:** arabinose, dextrose, digitose, fructose, fucose, galactose, 2-d-galactose, 6-d-glucose, gulose, lyxose, mannose, 3-O-methylglucose, rhamnose, 2-d-ribose, ribose, sorbose, tagatose

KEY WORDS

post-column reaction

REFERENCE

Mopper, K. Improved chromatographic separations on anion-exchange resins. I. Partition chromatography of sugars in ethanol, *Anal. Biochem.*, **1978**, *85*, 528–532.

SAMPLE

Matrix: solutions

Sample preparation: Inject a 20 μL aliquot of a solution in mobile phase.

HPLC VARIABLES

Column: 150 \times 6 μm Shodex RSPak DC-613 sulfonated polystyrene 55% cross-linked with divinylbenzene calcium form (Showa Denko)

Mobile phase: MeCN:water 80:20

Column temperature: 4

Flow rate: 0.5

Injection volume: 20

Detector: UV 280 following post-column reaction. The column effluent mixed with 500 mM pH 8.5 borate buffer pumped at 0.5 mL/min and 1% 2-cyanoacetamide in water pumped at 0.5 mL/min and the mixture flowed through a 5 m \times 0.5 mm ID PTFE coil at $100 \pm 1^\circ$ and a 1 m \times 0.5 mm PTFE cooling coil to the detector.

CHROMATOGRAM

Retention time: k' 2.68 (α -D-xylose), k' 2.20 (β -D-xylose)

OTHER SUBSTANCES

Also analyzed: allose, altrose, arabinose, dextrose, fucose, galactose, gulose, idose, lyxose, mannose, rhamnose

KEY WORDS

post-column reaction

REFERENCE

Honda, S.; Suzuki, S.; Kakehi, K. Improved analysis of aldose anomers by high-performance liquid chromatography on cation-exchange columns, *J. Chromatogr.*, **1984**, *291*, 317–325.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 250 \times 4 μm Shodex RSPak DC-613 sulfonated polystyrene 55% cross-linked with divinylbenzene (H^+) (Showa Denko)

Mobile phase: MeCN:water 90:10

Flow rate: 0.6

Injection volume: 20

Detector: E, Irika E-502, glassy carbon working electrode 0.30 V, Ag/AgCl reference electrode, following post-column reaction. The column effluent mixed with 200 mM pH 9.5 borate buffer pumped at 0.25 mL/min and 1.5% 2-cyanoacetamide in water pumped at 0.25 mL/min and the mixture flowed through a 10 m \times 0.5 mm ID PTFE coil at 100° and a 1 m \times 0.5 mm PTFE cooling coil to the detector.

CHROMATOGRAM

Retention time: 13

Limit of detection: 20 pmole

OTHER SUBSTANCES

Simultaneous: fucose, galactose, rhamnose

KEY WORDS

post-column reaction

REFERENCE

Honda,S.; Konishi,T.; Suzuki,S. Electrochemical detection of reducing carbohydrates in high-performance liquid chromatography after post-column derivatization with 2-cyanoacetamide, *J.Chromatogr.*, **1984**, 299, 245–251.

SAMPLE

Matrix: solutions

Sample preparation: Condition a Sep-Pak C18 SPE cartridge with 3 mL MeCN and 5 mL water. Mix 250 μ L of an aqueous solution with 225 μ L 1% dansylhydrazine in EtOH and 45 μ L 10%trichloroacetic acid in water, heat at 65° for 20 min, dilute with water to an organic solvent concentration of \leq 5%, add a 5 mL aliquot to the SPE cartridge, wash with 5 mL MeCN:water 5:95 at \leq 2 mL/min, elute with 6 mL MeCN:water 20:80 at \leq 2 mL/min (*J. Chromatogr.* 1983, 256, 27), lyophilize the eluate, reconstitute with MeCN:water 20:80, inject a 20 μ L aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 10 μ m 600 RPB C18 (Alltech)

Mobile phase: MeCN:water 20:80 containing 10 mM formic acid, 40 mM acetic acid, and 1 mM triethylamine. (After each run flush column with MeCN:MeOH 20:80 for 5 min.)

Flow rate: 1

Injection volume: 20

Detector: UV 254

CHROMATOGRAM

Retention time: 12

Limit of detection: 200-300 pmole

OTHER SUBSTANCES

Simultaneous: dextrose, fucose, galactose, lyxose, mannose

KEY WORDS

derivatization; SPE

REFERENCE

Eggert,F.M.; Jones,M. Measurement of neutral sugars in glycoproteins as dansyl derivatives by automated high-performance liquid chromatography, *J.Chromatogr.*, **1985**, 333, 123–131.

SAMPLE

Matrix: solutions

Sample preparation: Heat 100-200 pmole sample with 20 μ L 4 M trifluoroacetic acid and 20 μ L 4 M HCl in a tube sealed under vacuum at 100° for 6 h, add 500 pmole L-rhamnose, evaporate to dryness under reduced pressure at 50°, add 50 μ L 9.8% sodium bicarbonate solution (freshly prepared), add 2 μ L acetic anhydride, let stand at room temperature with occasional stirring for 30 min, add 200 μ L 100-200 mesh Dowex 50W-X2 (H⁺), check that pH is about 3. Add the mixture to a 100 \times 5 column and wash it with 5 bed volumes of water, evaporate to dryness under reduced pressure, add 5 μ L reagent, seal tube, heat at 100° for 13-15 min, add 2 μ L 20 mg/mL sodium cyanoborohydride in water (freshly prepared), reseal the tube, heat at 90° for 8 h, dilute with 20 μ L water, inject the whole amount on to a 600 \times 7.5 10 μ m TSK-GEL G2000PW column (Toyo Soda) and elute with 20 mM pH 7.5 ammonium acetate buffer at 0.5 mL/min, collect the sugar fraction at 40-55 min. Evaporate the eluate to dryness and reconstitute it with 250 μ L water, inject a 5 μ L aliquot. (Prepare reagent by mixing 500 mg 2-aminopyridine, 400 μ L concentrated HCl, and 11 mL water.)

HPLC VARIABLES

Column: two 250 \times 4.6 5 μ m Ultrasphere-ODS column in series

Mobile phase: MeCN:250 mM pH 4.0 sodium citrate buffer 1:99

Flow rate: 0.5

Injection volume: 5

Detector: F ex 320 em 400

CHROMATOGRAM

Retention time: 39

Internal standard: L-rhamnose
Limit of quantitation: 10 pmoles

OTHER SUBSTANCES

Simultaneous: N-acetyl-D-mannosamine, N-acetylgalactosamine, N-acetylglucosamine, 2-deoxy-D-ribose, dextrose, fucose, galactose, mannose, ribose

KEY WORDS

derivatization; SPE

REFERENCE

Takemoto, H.; Hase, S.; Ikenaka, T. Microquantitative analysis of neutral and amino sugars as fluorescent pyridylamino derivatives by high-performance liquid chromatography, *Anal. Biochem.*, **1985**, *145*, 245–250.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 80 × 3 Hitachi 2633

Mobile phase: 700 mM pH 8.5 Borate buffer containing 0.01% EDTA

Column temperature: 60

Flow rate: 0.7

Injection volume: 20

Detector: E, Irika E-502, glassy carbon working electrode 350 mV, Ag/AgCl reference electrode, following post-column reaction. The column effluent mixed with 100 mM ethylenediamine sulfate pumped at 0.25 mL/min and 700 mM pH 9.0 borate buffer pumped at 0.25 mL/min and the mixture flowed through a 30 m × 0.5 mm ID PTFE coil at 140° and a 10 m × 0.2 mm ID cooling coil to the detector.

CHROMATOGRAM

Retention time: 30

Limit of detection: 1 pmole

OTHER SUBSTANCES

Simultaneous: dextrose, galactose, mannose, rhamnose

KEY WORDS

post-column reaction

REFERENCE

Honda, S.; Enami, K.; Konishi, T.; Suzuki, S.; Kakehi, K. Use of ethylenediamine sulphate for post-column derivatization of reducing carbohydrates to electrochemically oxidizable compounds in high-performance liquid chromatography, *J. Chromatogr.*, **1986**, *361*, 321–329.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: HPIC-AS6 anion-exchange (Dionex)

Mobile phase: 150 mM NaOH

Column temperature: 36

Flow rate: 0.5

Injection volume: 20

Detector: UV 500 following post-column reaction. The column effluent mixed with the reagent pumped at 0.2 mL/min (?) and the mixture flowed through a knitted 10 m × 0.3 mm ID PTFE coil at 90° and then a short knitted PTFE coil at 22° to the detector. (Reagent was 2 mg/mL thymol in concentrated sulfuric acid, let stand for 30 min after preparation, discard after 48 h. (The reagent was displaced from a pressure vessel into the post-column reaction system by pumping n-heptane into the vessel.)

CHROMATOGRAM**Retention time:** 9.3**Limit of detection:** 100 ng

OTHER SUBSTANCES**Simultaneous:** arabinose, desoxyribose, dextrose, fructose, galactose, lactose, maltose, mannose, raffinose, ribose, saccharose**Noninterfering:** methyl arabinose, methyl glycoside, rhamnose, rutinose, trehalose

KEY WORDSpost-column reaction

REFERENCEEngelhardt, H.; Ohs, P. Trace analysis of sugars by HPLC and post-column derivatization, *Chromatographia*, **1987**, *23*, 657–662.

SAMPLE**Matrix:** solutions**Sample preparation:** Dissolve 50 mg sugars in 700 μ L pyridine, add 700 μ L 720 mM hydroxylamine hydrochloride in pyridine, heat at 60° for 10 min, add 250 μ L acetic anhydride, heat at 75° for 10 min, evaporate to dryness under reduced pressure, reconstitute with 3 mL chloroform. Wash the organic layer three times with 6 mL portions of water and dry it over anhydrous sodium sulfate, evaporate to dryness under reduced pressure, take up in chloroform, pass through silica gel using chloroform, evaporate the eluate to dryness, reconstitute, inject a 5 μ L aliquot.

HPLC VARIABLES**Column:** 250 \times 4.5 μ m μ Bondapak C18**Mobile phase:** Gradient. MeCN:water from 35:75 to 50:50 over 15 min.**Flow rate:** 1**Injection volume:** 5**Detector:** UV 207

CHROMATOGRAM**Retention time:** k' 3.1**Limit of detection:** 3 μ g

OTHER SUBSTANCES**Also analyzed:** allose, altrose, arabinose, dextrose, fructose, galactose, gulose, idose, lyxose, mannose, ribose, talose

KEY WORDSderivatization

REFERENCEVelasco, D.; Castells, J.; Lopez-Calahorra, F. High-performance liquid chromatographic separation of monosaccharides as their peracetylated ketoximes and aldonoitriles, *J. Chromatogr.*, **1990**, *519*, 228–236.

SAMPLE**Matrix:** solutions**Sample preparation:** 10 μ L EtOH containing sugars + 110 μ L EtOH:acetic acid 99.9:0.1 + 100 μ L Fmoc-hydrazine in MeCN, mix, heat at 65° for 3 h, cool to room temperature, dilute with EtOH or EtOH:acetic acid 99.9:0.1, inject an aliquot. (Prepare Fmoc-hydrazine as follows. Dissolve 100 mg 9-fluorenylmethylchloroformate in 25 mL MeCN, add this solution dropwise with stirring to 1 mL hydrazine hydrate (Caution! Hydrazine hydrate is a carcinogen!), stir for 30 min, evaporate under reduced pressure, use the crude product or recrystallize from EtOH or MeCN (mp 173–5°). Prepare a solution of Fmoc-hydrazine in MeCN so that the hydrazine:sugar ratio is 10:1.)

HPLC VARIABLES**Column:** 150 \times 4.6 5 μ m Zorbax ODS

Mobile phase: Gradient. A was MeCN:water 27:73 containing 80 mM acetic acid. B was MeCN:water 30:70 containing 80 mM acetic acid. A:B from 100:0 to 0:100 over 30 min.

Flow rate: 1

Detector: F ex 270 em 320

CHROMATOGRAM

Retention time: 15

Limit of detection: 0.1 pmole

OTHER SUBSTANCES

Simultaneous: fructose, galactose, lactose, maltose, mannose, ribose

KEY WORDS

derivatization

REFERENCE

Zhang,R.-E.; Cao,Y.-L.; Hearn,M.W. Synthesis and application of Fmoc-hydrazine for the quantitative determination of saccharides by reversed-phase high-performance liquid chromatography in the low and subpicomole range, *Anal.Biochem.*, **1991**, 195, 160–167.

SAMPLE

Matrix: solutions

Sample preparation: Add 55 μ L phenylisocyanate to a 1 mg/mL solution in DMF, heat at 55° for 95 min, cool, add 500 μ L MeOH, let stand for 5 min, make up to 6 mL with DMF, dilute an aliquot 10-fold with DMF, inject an aliquot.

HPLC VARIABLES

Column: 220 \times 4.6 5 μ m ODS 224 RP18 (Brownlee)

Mobile phase: MeCN:water 60:40

Flow rate: 2

Injection volume: 10

Detector: UV 240

CHROMATOGRAM

Retention time: 9.42

Limit of detection: 0.2 ng

OTHER SUBSTANCES

Simultaneous: allose, arabinose, deoxyglucose, deoxyribose, dextrose, fucose, galactose, lyxose, mannose, methylgalactoside, methylglucoside, methylmannoside, rhamnose, ribose

KEY WORDS

derivatization; more than one derivative was observed, retention time is for major derivative

REFERENCE

Rakotomanga,S.; Baillet,A.; Pellerin,F.; Baylocq-Ferrier,D. Liquid chromatographic analysis of monosaccharides with phenylisocyanate derivatization, *J.Pharm.Biomed.Anal.*, **1992**, 10, 587–591.

SAMPLE

Matrix: solutions

Sample preparation: Add 25 nmoles 3-O-methylglucose, evaporate the solution to dryness, add 50 μ L 300 mM sodium cyanoborohydride in 2 M pH 7.0 ammonium acetate (freshly prepared), heat at 105° for 4 h, add 100 μ L water, add 40 μ L 6 M formic acid, evaporate to dryness under reduced pressure, add 500 μ L MeOH, evaporate to dryness, repeat MeOH evaporation twice more, add 100 μ L EtOH:water:triethylamine 40:40:20, evaporate to dryness, add 100 μ L EtOH:triethylamine:water:phenylisothiocyanate 70:10:10:10, let stand at room temperature for 20 min, evaporate to dryness under reduced pressure, reconstitute with 20 μ L MeCN:water 60:40, add 180 μ L MeCN:5 mM pH 7.4 sodium phosphate buffer 5:95, filter, inject an aliquot.

HPLC VARIABLES

Column: 250 \times 4.6 5 μ m Microsorb C18

Mobile phase: Gradient. A was 50 mM pH 6.8 ammonium acetate. B was 100 mM pH 6.8 ammonium acetate in MeCN:MeOH:water 44:10:46. A:B 78:22 until the run is over, to 0:100 over 5 min, maintain at 0:100 for 6 min, return to initial conditions over 5 min

Column temperature: 30

Flow rate: 0.8

Detector: UV 254

CHROMATOGRAM

Retention time: 13.5

Internal standard: 3-O-methylglucose (20)

Limit of detection: 50 pmole

OTHER SUBSTANCES

Simultaneous: dextrose, fucose, galactose, mannose, ribose

KEY WORDS

derivatization

REFERENCE

Spiro,M.J.; Spiro,R.G. Monosaccharide determination of glycoconjugates by reverse-phase high-performance liquid chromatography of their phenylthiocarbamyl derivatives, *Anal.Biochem.*, **1992**, 204, 152–157.

SAMPLE

Matrix: solutions

HPLC VARIABLES

Column: 300 × 7.8 Aminex HPX-87P cation-exchange (Bio-Rad)

Mobile phase: Water

Column temperature: 85

Flow rate: 0.8

Injection volume: 20

Detector: UV 550 following post-column reaction. The column effluent mixed with the reagent pumped at 0.4 mL/min and the mixture flowed through a 20 m × 0.3 mm ID PTFE coil at 90° and a coil at 0° to the detector. (Reagent was prepared from 4000 ppm Purpald (4-amino-3-hydrazino-5-mercapto-1,2,4-triazole) in 2 M NaOH (A) and 40 mM hydrogen peroxide (B). A:B 70:30.)

CHROMATOGRAM

Retention time: 16.5

Limit of detection: 30 ng

OTHER SUBSTANCES

Simultaneous: arabinose, dextrose, fructose, galactose, mannose, ribose

KEY WORDS

post-column reaction

REFERENCE

Del Nozal,M.J.; Bernal,J.L.; Hernandez,V.; Toribio,L.; Mendez,R. Purpald (4-amino-3-hydrazino-5-mercapto-1,2,4-triazole) as a reagent for post-column derivatization of neutral monosaccharides in high pressure liquid chromatography, *J.Liq.Chromatogr.*, **1993**, 16, 1105–1116.

SAMPLE

Matrix: solutions

Sample preparation: Mix a 50 µL aliquot of a 500 µM saccharide solution in MeCN:water 30:70 with 50 µL 10 mM reagent in MeCN and 100 µL 0.5% trichloroacetic acid in MeCN, heat at 65° in the dark for 3 h. Remove a 50 µL aliquot of the reaction mixture, add 200 µL water, add 200 µL ethyl acetate, mix, centrifuge at 3000 rpm for 2 min, repeat the ethyl acetate wash twice more. Dry the aqueous layer under reduced pressure, reconstitute with 200 µL MeCN, inject an aliquot. (Synthesis of reagent, R-(+)-DBD-ProCZ, is as follows. Dissolve 0.5 g magnesium sulfate heptahydrate and 6 g NaOH in 60 mL water, throughout the reaction keep the

flask at about 20° with cold water cooling, add 15 mL 30% hydrogen peroxide, add 75 mL MeOH, add 12.1 g powdered benzoyl peroxide in one go, stir for 10 min, pour into 150 mL 20% sulfuric acid, extract three times with 50 mL portions of chloroform, determine peroxybenzoic acid concentration by iodometric titration (Tetrahedron 1967, 23, 3327). Slowly add 110 mL 1 M peroxybenzoic acid in chloroform to 7 g 2,6-difluoroaniline dissolved in 100 mL chloroform, stir at room temperature, when reaction is complete (iodometric titration) wash with 2% sodium thiosulfate, wash with 5% sodium carbonate, wash with water, dry over anhydrous sodium sulfate, evaporate to dryness under reduced pressure, recrystallize 2,6-difluoronitrosobenzene from EtOH (mp 108.5-109.5). Stir 8.5 g 2,6-difluoronitrosobenzene in 85 mL DMSO at room temperature and add a solution of 3.91 g sodium azide in 85 mL DMSO dropwise, let stand for about 1 h, add to a large volume of water, extract with ether, dry the extracts over anhydrous sodium sulfate, evaporate to dryness under reduced pressure and distil to give 4-fluoro-2,1,3-benzoxadiazole as a colorless oil (bp 83°/12 mm Hg) (J.Chem.Soc.(C) 1970, 1433). Add 11 mL chlorosulfonic acid dropwise to 3 g 4-fluoro-2,1,3-benzoxadiazole in 10 mL chloroform at 0-10° (use a calcium chloride drying tube), stir at room temperature for 1 h, reflux for 2 h, cool, slowly pour into ice water, remove the organic layer, extract the aqueous layer with chloroform, combine the organic layer, wash, dry over anhydrous magnesium sulfate, evaporate under reduced pressure, take up the residue in 5 mL benzene (Caution! Benzene is a carcinogen!), chromatograph on a 150 × 30 column of silica gel (100-200 mesh Kanto Chemical) with n-hexane:benzene 50:50, evaporate the appropriate fractions to give 4-(chlorosulfonyl)-7-fluoro-2,1,3-benzoxadiazole (CBD-F) as pale yellow needles (mp 64-66°) (Anal. Chem. 1984, 56, 2461). Stir 0.76 g CBD-F in 70 mL MeCN at 0-10° and add 1 g dimethylamine hydrochloride in 10 mL 100 mM pH 10 borax dropwise, adjust pH to 5 with 1 M HCl, concentrate to about 10 mL under reduced pressure, extract three times with 200 mL portions of diethyl ether, wash with water, dry over anhydrous magnesium sulfate, evaporate under reduced pressure, chromatograph on a 500 × 20 column of silica gel with chloroform, isolate the appropriate fraction and re-chromatograph on the same column with ethyl acetate:benzene 1:2 to give 4-(N,N-dimethylaminosulfonyl)-7-fluoro-2,1,3-benzoxadiazole (DBD-F) as white needles (mp 124-125°) (yield = 1% !). On a Merck no. 5714 60F₂₅₄ tlc plate eluted with chloroform DBD-F has R_f 0.32 and lies between two other reaction products (Analyst 1989, 114, 413). It is also reported that DBD-F can be purchased from Tokyo Kasei. Add 100 mg DBD-F in 10 mL MeCN to 47 mg R-(+)-proline in 20 mL 250 mM pH 11.5 sodium carbonate solution, stir at room temperature for 30 min, wash with ethyl acetate, adjust the pH of the aqueous layer to 1-2 with 2 M HCl, extract three times with 30 mL ethyl acetate. Combine the extracts and evaporate them under reduced pressure, recrystallize from benzene/ethyl acetate to give R-(+)-4-(N,N-dimethylaminosulfonyl)-7-(2-carboxypyrrolidin-1-yl)-2,1,3-benzoxadiazole (DBD-Pro) as yellow needles (mp 187-9° d) (Analyst 1989, 114, 1233). Suspend 55 mg (R)-(+)-DBD-Pro in 55 mL anhydrous diethyl ether at 0°, add 110 mg phosphorus pentachloride, stir at 5° for 1 h, filter quickly, evaporate to dryness under reduced pressure, dry under vacuum over phosphorus pentoxide for 12 h to give R-(+)-4-(N,N-dimethylaminosulfonyl)-7-(2-chloroformylpyrrolidin-1-yl)-2,1,3-benzoxadiazole (DBD-Pro-Cl) as yellow crystals (mp 116-17°) (Analyst 1993, 118, 759). Add 130 mg DBD-Pro-Cl dissolved in 25 mL anhydrous benzene dropwise to 100 mL MeOH containing 70 mg hydrazine hydrate, stir for 30 min at room temperature, evaporate under reduced pressure, recrystallize from ethyl acetate:MeOH 90:10 to give R-(+)-4-(2-carbazolylpyrrolidin-1-yl)-7-(N,N-dimethylaminosulfonyl)-2,1,3-benzoxadiazole (R-(+)-DBD-ProCZ) as orange crystals (mp 107-109°) (Anal. Proc. 1994, 31, 265.).

HPLC VARIABLES

Column: 150 × 4.6 5 µm Inertsil ODS-80A

Mobile phase: MeCN:water 15:85

Column temperature: 40

Flow rate: 1

Detector: F ex 450 em 540

CHROMATOGRAM

Retention time: 15.81

OTHER SUBSTANCES

Simultaneous: N-acetyl-D-glucosamine, arabinose, dextrose, galactose, mannose

KEY WORDS

derivatization

REFERENCE

Toyooka,T.; Kuze,A. Determination of saccharides labelled with a fluorescent reagent, DBD-ProCZ, by liquid chromatography, *Biomed.Chromatogr.*, **1997**, *11*, 132–136.

SAMPLE

Matrix: urine

Sample preparation: Mix acetone with urine so as to make a 63:47 acetone:urine mixture, centrifuge a 6 mL aliquot. Evaporate the acetone from the supernatant under a stream of helium at 35°, add 30 mg Dowex 50W-X8, add 30 mg Dowex 1-X8, agitate, centrifuge, inject a 1-10 µL aliquot of the supernatant.

HPLC VARIABLES

Column: 150 × 4.3 Hitachi 3013 N anion-exchange resin, phosphate form

Mobile phase: MeCN:water 83:17

Column temperature: 60

Flow rate: 1

Injection volume: 1-10

Detector: UV 530 following post-column reaction. The column effluent mixed with the reagent pumped at 1.5 mL/min, the mixture flowed through a 3 m × 0.5 mm i.d. coil of PTFE tubing at 85° and a 1 m × 0.5 mm i.d. coil of PTFE tubing at room temperature to the detector. (Reagent was 2 g/L blue tetrazolium in EtOH:water 50:50 containing 180 mM NaOH.)

CHROMATOGRAM

Retention time: 10

Limit of detection: 10 ng

OTHER SUBSTANCES

Extracted: arabinose, dextrose, fructose, fucose, galactose, lactose, ribose

KEY WORDS

post-column reaction

REFERENCE

D'Amboise,M.; Hanai,T.; Noël,D. Liquid-chromatographic measurement of urinary monosaccharides, *Clin. Chem.*, **1980**, *26*, 1348–1350.

SAMPLE

Matrix: urine

Sample preparation: 100 µL Urine + 10 µL 10% trichloroacetic acid in water + 50 µL 5% dansyl hydrazine in MeCN, heat at 65° for 20 min, cool in ice, add an equal volume of water, inject a 10-300 µL aliquot.

HPLC VARIABLES

Column: 250 × 4.6 5 µm Nucleosil ODS

Mobile phase: MeCN:80 mM acetic acid 21:79

Flow rate: 1

Injection volume: 10

Detector: F ex 360 em >470

CHROMATOGRAM

Retention time: 11.5

Limit of detection: 5-15 pmole

OTHER SUBSTANCES

Extracted: cellobiose, 2-deoxyglucose, 2-deoxyribose, dextrose, fucose, galactose, gentobiose, lactose, maltose, mannose, rhamnose, ribose

Interfering: arabinose, fructose

KEY WORDS

derivatization

REFERENCE

Mopper, K.; Johnson, L. Reversed-phase liquid chromatographic analysis of Dns-sugars. Optimization of derivatization and chromatographic procedures and applications to natural samples, *J. Chromatogr.*, **1983**, 256, 27–38.

SAMPLE

Matrix: urine

Sample preparation: 10 μ L Urine + 200 μ L reagent, heat at 65° for 16 h, cool to room temperature, inject a 5 μ L aliquot of the clear supernatant. (Prepare reagent by dissolving 5 mg Fmoc-hydrazine in 1 mL MeCN, add 10 μ L buffer. Buffer was 1.44 M formic acid containing 600 mM NaOH. Prepare Fmoc-hydrazine as follows. Dissolve 1 g 9-fluorenylmethyl chloroformate in 100 mL EtOH, add this solution dropwise with stirring to 10 mL hydrazine hydrate (Caution! Hydrazine hydrate is a carcinogen!), stir for 30 min, filter off the precipitate, wash it twice with 20 mL portions of ice-cold EtOH, dry at room temperature.)

HPLC VARIABLES

Guard column: 10 \times 4.6 3 μ m Spherisorb ODS II

Column: 125 \times 4.6 3 μ m Spherisorb ODS II

Mobile phase: Gradient. Isopropanol:isobutyl alcohol:water 6:6:88 for 13 min, to 80:0:20 (step gradient), maintain at 80:0:20 for 6 min, re-equilibrate at initial conditions.

Column temperature: 50

Injection volume: 5

Detector: F ex 270 em 315

CHROMATOGRAM

Retention time: 9.8

Limit of detection: 20 nM

OTHER SUBSTANCES

Extracted: lactulose, 3-O-methyl-D-glucose, rhamnose

KEY WORDS

derivatization

REFERENCE

Rooyakkers, D.R.; van Eijk, H.M.H.; Deutz, N.E.P. Simple and sensitive multi-sugar-probe gut permeability test by high-performance liquid chromatography with fluorescence labelling, *J. Chromatogr. A*, **1996**, 730, 99–105.


Yohimbine

Molecular formula: C₂₁H₂₆N₂O₃

Molecular weight: 354.45

CAS Registry No.: 146-48-5

Merck Index: 10236



SAMPLE

Matrix: blood

Sample preparation: 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 μ L mobile phase, centrifuge at 2800 g for 5 min, inject a 50 μ L aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

HPLC VARIABLES

Column: 300 \times 3.9 4 μ m NovaPack C18