

**HPLC VARIABLES**

**Column:** 250 × 4.6 5 µm Supelcosil LC-DP (A) or 250 × 4.5 µm LiChrospher 100 RP-8 (B)

**Mobile phase:** MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

**Flow rate:** 0.6

**Injection volume:** 25

**Detector:** UV 229

**CHROMATOGRAM**

**Retention time:** 8.15 (A), 4.37 (B)

**OTHER SUBSTANCES**

**Also analyzed:** acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordiazepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenoprofen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxychloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazinol, mefenamic acid, meperidine, mephénytoin, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, metronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymetazoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, pheniramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltoloxamine, phenytoin, pimozide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quinine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, secobarbital, sertraline, sotalol, spironolactone, sulfapyrazone, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocainide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, trifluoromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazoline, zopiclone

**KEY WORDS**

details of plasma extraction

**REFERENCE**

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J.Chromatogr.A*, **1995**, 692, 103-119.


# Zafirlukast

**Molecular formula:** C<sub>31</sub>H<sub>33</sub>N<sub>3</sub>O<sub>6</sub>S

**Molecular weight:** 575.69

**CAS Registry No.:** 107753-78-6

**Merck Index:** 10241



**SAMPLE****Matrix:** blood

**Sample preparation:** Condition a Bond Elut (LRBE) C18 SPE cartridge (Varian) with two 1 mL portions of MeOH and two 1 mL portions of water. Add 1 mL MeCN to 1 mL plasma, mix, centrifuge at 1000 g for 5 min, decant the supernatant, add it to 5 mL 10 mM pH 7.0 triethylammonium phosphate buffer. Add the diluted supernatant to the SPE cartridge, let the whole supernatant volume pass through the cartridge, elute with three separate aliquots of 500  $\mu$ L MeCN:THF:triethylamine 90:10:0.2, evaporate the solvent under nitrogen at 40°, reconstitute the sample in 250  $\mu$ L mobile phase, filter (0.45  $\mu$ m syringe filter), inject an aliquot.

**HPLC VARIABLES****Guard column:** 12.5  $\times$  4.0 5  $\mu$ m Zorbax CN**Column:** 150  $\times$  4.6 5  $\mu$ m Zorbax CN**Mobile phase:** THF:hexane:90% glacial acetic acid 30:70:0.1**Flow rate:** 0.9**Injection volume:** 150**Detector:** F ex 250 em 452**CHROMATOGRAM****Retention time:** 9-10

**Internal standard:** 4-5(-cyclopentyloxycarbonylamino-1-methylindol-3-ylmethyl)-3-methoxy-*a*-benzylsulfonylbenzamide (ICI 198 707) (10-11)


**Limit of quantitation:** 750 pg/mL**OTHER SUBSTANCES****Noninterfering:** acetaminophen, albuterol, aspirin, benzyl alcohol, caffeine, ibuprofen**KEY WORDS**

plasma; SPE

**REFERENCE**

Bui, K.H.; Kennedy, C.M.; Azumaya, C.T.; Birmingham, B.K. Determination of zafirlukast, a selective leukotriene antagonist, human plasma by normal-phase high-performance liquid chromatography with fluorescence detection, *J. Chromatogr. B*, **1997**, 696, 131-136.

# Zalcitabine

**Molecular formula:** C<sub>9</sub>H<sub>13</sub>N<sub>3</sub>O<sub>3</sub>**Molecular weight:** 211.22**CAS Registry No.:** 7481-89-2**Merck Index:** 10242**Lednicer No.:** 5 98**SAMPLE****Matrix:** blood

**Sample preparation:** Condition a Sep-Pak C18 SPE cartridge with 6 mL MeOH and 12 mL MeOH. Add 1.25 mL plasma to the SPE cartridge at 0.5 mL/min, wash with 1.5 mL water, elute with 2 mL MeOH at 0.5 mL/min. Evaporate the eluate to dryness under a stream of nitrogen, reconstitute the residue in 200  $\mu$ L water, filter (Ultrafree-MC 0.45  $\mu$ m) while centrifuging at 11900 g for 10 min, inject a 100  $\mu$ L aliquot of the filtrate.

**HPLC VARIABLES****Column:** C18**Mobile phase:** MeCN:water:heptafluorobutyric acid 6.5:93.4:0.1**Flow rate:** 2**Injection volume:** 100**Detector:** UV 288

---

**CHROMATOGRAM****Retention time:** 5.8**Internal standard:** 5-methyldeoxycytidine (UV 306) (3.8)**Limit of detection:** 30 nM

---

**KEY WORDS**SPE; pharmacokinetics; plasma

---

**REFERENCE**

Hawkins, M.E.; Poplack, D.G.; Pizzo, P.A.; Balis, F.M. High-performance liquid chromatographic method for the analysis of 2',3'-dideoxycytidine in human plasma, *J.Chromatogr.*, **1990**, 532, 442-444.

---

**SAMPLE****Matrix:** blood

**Sample preparation:** Condition a 3 mL Supelclean LC-18 SPE cartridge (Supelco) with 2 mL MeOH and 2 mL water. 1 mL Plasma + 5 ng IS, vortex, add to the SPE cartridge, wash with 2 mL water, elute with 2 mL MeOH:water 20:80. Evaporate the eluate to dryness, reconstitute the residue in 50  $\mu$ L MeOH:water 10:90, inject a 25  $\mu$ L aliquot.

---

**HPLC VARIABLES****Column:** 250  $\times$  4.6 5  $\mu$ m Ultrasphere ODS**Mobile phase:** MeOH:50 mM ammonium acetate 10:90**Flow rate:** 1**Injection volume:** 25

**Detector:** MS, VG Trio-3 triple quadrupole, thermospray interface with heated capillary at 280°, source 200°, repeller electrode 200 V, SIM, m/z 212 for zalcitabine, m/z 216 for IS

---

**CHROMATOGRAM****Retention time:** 8.9**Internal standard:** [ $^{15}\text{N}_3$ ,  $^2\text{H}_2$ ]zalcitabine**Limit of quantitation:** 0.25 ng/mL

---

**KEY WORDS**plasma; SPE; pharmacokinetics

---

**REFERENCE**

Jajoo, H.K.; Bennett, S.M.; Kornhauser, D.M. Thermospray liquid chromatographic-mass spectrometric analysis of anti-AIDS nucleosides: quantification of 2',3'-dideoxycytidine in plasma samples, *J.Chromatogr.*, **1992**, 577, 299-304.

---

**SAMPLE****Matrix:** blood

**Sample preparation:** 200  $\mu$ L Plasma + 400  $\mu$ L ethyl acetate:MeCN 50:50, vortex for 30 s. Remove a 200  $\mu$ L aliquot of the supernatant and evaporate it to dryness under a stream of nitrogen at 50°, reconstitute the residue in 200  $\mu$ L mobile phase, vortex for 10 s, inject a 10  $\mu$ L aliquot.

---

**HPLC VARIABLES****Guard column:** 4  $\times$  4 5  $\mu$ m Lichrospher 60 RP-select B**Column:** 125  $\times$  4 5  $\mu$ m Lichrospher 60 RP-select B**Mobile phase:** MeOH:pH 7.0 phosphate buffer 5:95**Flow rate:** 1**Injection volume:** 10**Detector:** UV 270

---

**OTHER SUBSTANCES****Also analyzed:** didanosine (UV 250)

---

**KEY WORDS**

plasma; rabbit; pharmacokinetics

---

**REFERENCE**

Mirchandani,H.L.; Chien,Y.W. Intestinal absorption of dideoxynucleosides: Characterization using a multiloop in situ technique, *J.Pharm.Sci.*, **1995**, *84*, 44–48.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Filter (Millipore Ultrafree-MC, 10000 molecular mass limit) 250  $\mu$ L serum while centrifuging at 17000 g for 1.5 h, inject a 50  $\mu$ L aliquot of the clear ultrafiltrate.

---

**HPLC VARIABLES**

**Column:** 150  $\times$  3.9 4  $\mu$ m Nova-Pak phenyl

**Mobile phase:** Isopropanol:20 mM pH 5 sodium citrate 2.5:97.5

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 250

---

**CHROMATOGRAM**

**Retention time:** 3.0

---

**OTHER SUBSTANCES**

**Extracted:** didanosine, zidovudine

---

**KEY WORDS**

serum; ultrafiltrate

---

**REFERENCE**

Rosell-Rovira,M.L.; Pou-Clavé,L.; Lopez-Galera,R.; Pascual-Mostaza,C. Determination of free serum didanosine by ultrafiltration and high-performance liquid chromatography, *J.Chromatogr.B*, **1996**, *675*, 89–92.

---

**SAMPLE**

**Matrix:** blood, cell suspensions, perfusate

**Sample preparation:** Centrifuge cellular suspensions at 17000 g for 5 min, inject a 25  $\mu$ L aliquot. Centrifuge perfusion fluid at 17000 g for 5 min, inject a 50  $\mu$ L aliquot. Dilute 1 mL plasma with 1 mL saturated ammonium sulfate, vortex for 30 s, centrifuge at 3000 g for 2 min, inject a 50  $\mu$ L aliquot of the supernatant.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4 5  $\mu$ m Phenyl Hypersil NC-04

**Mobile phase:** MeOH:1.4 g/L sodium acetate 3:97, adjusted to pH 6.55

**Flow rate:** 1

**Injection volume:** 25-50

**Detector:** UV 271

---

**CHROMATOGRAM**

**Retention time:** 10

**Limit of detection:** 50 ng/mL

---

**KEY WORDS**

plasma

---

**REFERENCE**

Frijus-Plessen,N.; Michaelis,H.C.; Foth,H.; Kahl,G.F. Determination of 3'-azido-3'-deoxythymidine, 2',3'-dideoxycytidine, 3'-fluoro-3'-deoxythymidine and 2',3'-dideoxyinosine in biological samples by high-performance liquid chromatography, *J.Chromatogr.*, **1990**, *534*, 101–107.

---

**SAMPLE**

**Matrix:** bulk

**Sample preparation:** Prepare a 100  $\mu$ g/mL aqueous solution, inject a 5  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250 × 4.6 5 µm Cyclobond I RSP (Advanced Separation Technologies)

**Mobile phase:** 2.5 mL/L Triethylamine in water adjusted to pH 6.5 with glacial acetic acid (Pass mobile phase through a silica column to saturate it with silica.)

**Flow rate:** 0.25

**Injection volume:** 5

**Detector:** UV 270

---

**CHROMATOGRAM**

**Retention time:** 29

**Limit of detection:** 0.05-0.1% (of zalcitabine)

---

**OTHER SUBSTANCES**

**Simultaneous:** impurities

---

**REFERENCE**

Scypinski,S.; Ross,A.J. Liquid chromatographic separation of zalcitabine and its stereoisomers, *J.Pharm.Biomed.Anal.*, **1994**, 12, 1271–1276.

---

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Inject a 15 µL aliquot.

---

**HPLC VARIABLES**

**Column:** 200 × 4.6 5 µm HP Hypersil ODS

**Mobile phase:** MeCN:20 mM pH 7.0 Na<sub>2</sub>HPO<sub>4</sub> 5:95

**Column temperature:** 37

**Flow rate:** 1

**Injection volume:** 15

**Detector:** UV 265

---

**CHROMATOGRAM**

**Retention time:** 3.89

---

**OTHER SUBSTANCES**

**Simultaneous:** degradation products

---

**REFERENCE**

Kim,D.-D.; Chien,Y.W. Transdermal delivery of dideoxynucleoside-type anti-HIV drugs. 1. Stability studies for hairless rat skin permeation, *J.Pharm.Sci.*, **1995**, 84, 1061–1066.

---

---

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Briefly vortex 100 µL of an aqueous solution with 300 µL MeCN:water:acetic acid 80:19:1, centrifuge at 15000 g for 5 min. Remove the supernatant and add it to 100 µL 1.25 mM phenacyl bromide in MeCN, heat at 80° for 45 min, evaporate to dryness under reduced pressure at room temperature, reconstitute with 60 µL water, vortex briefly, centrifuge at 15000 g for 5 min, inject a 20 µL aliquot of the supernatant.

---

**HPLC VARIABLES**

**Column:** 100 × 4.6 5 µm Prodigy ODS-2 (Phenomenex)

**Mobile phase:** MeCN:water:trifluoroacetic acid 16:84:0.1

**Column temperature:** 45

**Flow rate:** 2.2

**Injection volume:** 20

**Detector:** F ex 305 em 370

---

**CHROMATOGRAM**

**Retention time:** 2

---

**OTHER SUBSTANCES**

**Simultaneous:** (-)-cis-5-(4-amino-1,2-dihydro-2-oxo-1-pyrimidinyl)-1,3-oxathiolane-2-methanol (3-TC)), cytidine, cytosine monophosphate, deoxycytidine

**Noninterfering:** guanidine, hypoxanthine, thymine, uracil, xanthine

**Interfering:** cytosine

**KEY WORDS**

derivatization

**REFERENCE**

Eisenberg, E.J.; Cundy, K.C. High-performance liquid chromatographic determination of cytosine-containing compounds by precolumn fluorescence derivatization with phenacyl bromide: application to antiviral nucleosides and nucleotides, *J. Chromatogr. B*, **1996**, 679, 119–127.

**SAMPLE**

**Matrix:** solutions

**Sample preparation:** Inject a 15  $\mu$ L aliquot.

**HPLC VARIABLES**

**Column:** 200  $\times$  4.6 5  $\mu$ m HP Hypersil ODS

**Mobile phase:** MeCN:20 mM pH 7.0  $\text{Na}_2\text{HPO}_4$  5:95

**Column temperature:** 37

**Flow rate:** 1

**Injection volume:** 15

**Detector:** UV 265

**CHROMATOGRAM**

**Retention time:** 3.89

**REFERENCE**

Kim, D.-D.; Chien, Y.W. Transdermal delivery of dideoxynucleoside-type Anti-HIV drugs. 2. The effect of vehicle and enhancer on skin permeation, *J. Pharm. Sci.*, **1996**, 85, 214–219.


# Zearalenone

**Molecular formula:**  $\text{C}_{18}\text{H}_{22}\text{O}_5$

**Molecular weight:** 318.37

**CAS Registry No.:** 17924-92-4

**Merck Index:** 10246

**SAMPLE**

**Matrix:** solutions

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu$ m Spherisorb ODS-1

**Mobile phase:** MeOH:water 80:20

**Flow rate:** 1

**Detector:** F ex 285 em 440 following post-column reaction. The column effluent mixed with 250 mM aluminum chloride hexahydrate in MeOH:water 75:25 pumped at 0.5 mL/min and the mixture flowed through a 5 m  $\times$  0.3 mm ID PTFE coil at 50° to the detector. (At the end of each day flush system with MeOH:water 75:25.)

**CHROMATOGRAM**

**Retention time:** 6.5

**KEY WORDS**

post-column reaction

---

**REFERENCE**

Hetmanski, M.T.; Scudamore, K.A. Detection of zearalenone in cereal extracts using high-performance liquid chromatography with post-column derivatization, *J. Chromatogr.*, **1991**, 588, 47–52.

---

# Zeranol


**Molecular formula:** C<sub>18</sub>H<sub>26</sub>O<sub>5</sub>

**Molecular weight:** 322.40

**CAS Registry No.:** 26538-44-3, 55331-29-8

**Merck Index:** 10251

---



---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** Hydrolyse serum or plasma with  $\beta$ -glucuronidase and sulfatase, extract with diethyl ether. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at room temperature, reconstitute the residue in 40  $\mu$ L buffer, add 100  $\mu$ L 1.5 mg/mL dansyl chloride in acetone, shake vigorously for 30 s, heat at 100° for 5 min, inject a 20  $\mu$ L aliquot. (Prepare buffer by adjusting the pH of 4 g/L sodium bicarbonate in water to 10.5 with 5 M NaOH.)

---

**HPLC VARIABLES**

**Column:** 250  $\times$  2.6 PAH-10 C18 (Perkin-Elmer)

**Mobile phase:** Gradient. MeCN:water from 60:40 to 95:5 over 15 min (Perkin-Elmer curve 1), maintain at 95:5 for 10 min. (Flush column with MeCN at 0.1 mL/min overnight.)

**Flow rate:** 1

**Injection volume:** 20

**Detector:** F ex 335 em 540

---

**CHROMATOGRAM**

**Retention time:** 11.5

**Limit of detection:** 80 ng

---

**OTHER SUBSTANCES**

**Extracted:** diethylstilbestrol, estriol, estrone, hexestrol, zanone, zenone

---

**KEY WORDS**

derivatization; cow; sheep; plasma; serum; LOD is too high for practical detection of compounds in serum and plasma.

---

**REFERENCE**

Rhys Williams, A.T.; Winfield, S.A.; Belloli, R.C. Dns derivatization of anabolic agents with high-performance liquid chromatographic separation and fluorescence detection, *J. Chromatogr.*, **1982**, 240, 224–229.

---

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 10 mL dichloromethane, shake for 20 min, centrifuge at 4000 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in mobile phase, inject an aliquot. To measure conjugates add 1 mL residual plasma from above procedure to 1 mL pH 5.5 phosphate buffer, add 1 mL 1000 U/mL helicase (type H-2 from *Helix pomatia*, Sigma), heat at 37° overnight, add 10 mL dichloromethane, shake for 20 min, centrifuge at 4000 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in mobile phase, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4 5  $\mu$ m Spherisorb ODS-1 C18

**Mobile phase:** MeCN:MeOH:0.2% acetic acid 30:20:50

**Flow rate:** 1.2  
**Detector:** radioactivity

---

**CHROMATOGRAM**  
**Retention time:** 8

---

**OTHER SUBSTANCES**  
**Extracted:** taleranol, zeralanone

---

**KEY WORDS**  
pig; plasma; tritium-labeled

---

**REFERENCE**  
Bories,G.; Suarez,A.F. Profiling of free and conjugated [<sup>3</sup>H]zeranol metabolites in pig plasma, *J.Chromatogr.*, **1989**, 489, 191–197.

---

**SAMPLE**  
**Matrix:** filters  
**Sample preparation:** Sonicate filter with 2 mL MeOH for 1 h, add 1 mL water, filter (0.5 µm PTFE), inject an aliquot.

---

**HPLC VARIABLES**  
**Column:** C18 Radial Compression (Waters)  
**Mobile phase:** MeOH:water 60:40  
**Flow rate:** 2  
**Detector:** UV 236

---

**CHROMATOGRAM**  
**Retention time:** 10.0  
**Limit of detection:** 7 ng/mL  
**Limit of quantitation:** 20 ng/mL

---

**OTHER SUBSTANCES**  
**Simultaneous:** taleranol, zearalanone, zearalenol, zearalenone

---

**REFERENCE**  
Neumeister,C.E. Environmental sampling and analysis for zeranol, *Am.Ind.Hyg.Assoc.J.*, **1987**, 48, 919–921.

---

**SAMPLE**  
**Matrix:** solutions

---

**HPLC VARIABLES**  
**Column:** 150 × 4.6 5 µm Hypersil ODS  
**Mobile phase:** MeOH:water 60:40  
**Injection volume:** 250  
**Detector:** UV

---

**CHROMATOGRAM**  
**Retention time:** 4.8

---

**OTHER SUBSTANCES**  
**Simultaneous:** diethylstilbestrol, trenbolone, nandrolone, dienestrol, hexestrol, 17α-methyltestosterone, medroxyprogesterone

---

**REFERENCE**  
Jansen,E.H.J.M.; Both-Miedema,R.; van den Berg,R.H. Application of optimization procedures for the separation of anabolic compounds by high-performance liquid chromatography, *J.Chromatogr.*, **1989**, 489, 57–64.

---

**SAMPLE**  
**Matrix:** tissue



**Sample preparation:** Homogenize 2.5 g tissue with 10 mL acetone for 20 s, sonicate for 5 min, centrifuge at 3200 rpm. Decant the supernatant into a silanized tube. Add 8 mL acetone to the pellet and repeat the extraction. Combine the supernatants. Add to a 5 mL pipette tip containing 1.5 g alumina (80-200 mesh, Brockman activity 1) followed by an Econo-Column filled with 1.0 g AGMP-1 resin (Bio-Rad), allow to pass through by gravity. Wash with four 1 mL portions of acetone:water 95:5. Remove the alumina column, wash the ion-exchange column with 1 mL acetone:water 95:5, elute with four 1 mL portions of 10% acetic acid in acetone. Evaporate the combined eluates to dryness with nitrogen at 40°. Add 500 µL water to the residue, extract twice with 2 mL portions of ether. Combine the ether layers and evaporate them to dryness. Reconstitute the residue in mobile phase B. Inject a 20 µL aliquot.

---

#### HPLC VARIABLES

**Column:** 250 × 4.6 5 µm Supelco silica

**Mobile phase:** Gradient. A was hexane. B was MeOH:hexane:2-propanol 45:40:15. A:B from 100:0 to 60:40 over 15 min.

**Flow rate:** 2.0

**Injection volume:** 20

**Detector:** UV 280

---

#### CHROMATOGRAM

**Retention time:** 11.54

**Limit of detection:** 10 ng

---

#### OTHER SUBSTANCES

**Extracted:** diethylstilbestrol, estradiol

**Simultaneous:** estrone, zeranol, zeralanone, zeralenone

---

#### KEY WORDS

chicken; muscle; normal phase; SPE

---

#### REFERENCE

Medina, M.B.; Sherman, J.T. High performance liquid chromatographic separation of anabolic oestrogens and ultraviolet detection of 17β-oestradiol, zeranol, diethylstilboestrol or zearalenone in avian muscle tissue extracts, *Food Addit. Contam.*, **1986**, 3, 263-272.

---

#### SAMPLE

**Matrix:** tissue

**Sample preparation:** Homogenize (Tekmar Model 3412-A20) 80 g ground tissue and 160 mL MeOH at medium-high speed for two 1 min intervals, centrifuge at 1500 rpm for 10 min, remove the supernatant, rinse apparatus with two 10 mL portions of MeOH. Combine the organic layers and evaporate them to 25 mL under a stream of nitrogen, adjust pH to 5.0 ± 0.2 with 20% acetic acid, add 4 mL freshly prepared β-glucuronidase solution in water (1 mg/mL, Type B-1, 500 000 U/g, Sigma), add 5 mL chloroform, mix by swirling, heat at 37° for 12-16 h, cool, add 45 mL chloroform, shake for 10 s, pass organic layer through 30-35 g anhydrous sodium sulfate (pewashed with 30 mL chloroform), extract aqueous layer with 50 mL chloroform, pass organic layer through sodium sulfate, wash sodium sulfate with two 5 mL portions of chloroform. Combine the chloroform layers and add 10 mL 2 M NaOH, shake vigorously for 30 s, discard organic layer, wash aqueous layer with two 50 mL portions of chloroform. Adjust the pH of the aqueous layer to 10.6-10.8 with 19-20 mL 1 M sodium bicarbonate solution, extract three times with 25 mL portions of chloroform, pass extracts through 30-35 g anhydrous sodium sulfate (pewashed with chloroform), wash sodium sulfate with 10 mL chloroform, evaporate chloroform layer to dryness under reduced pressure at 40°, reconstitute with 2 mL MeOH, inject a 25 µL aliquot.

---

#### HPLC VARIABLES

**Column:** 150 × 3.9 5 µm Nova-Pak RP-C18

**Mobile phase:** MeOH:buffer 50:50 (Buffer was 90 mM sodium acetate containing 10 mg/L EDTA adjusted to pH 6.9 with acetic acid.)

**Flow rate:** 1

**Injection volume:** 25

**Detector:** E, Bioanalytical Systems LC4B-17D, glassy carbon electrode +0.90 V, Ag/AgCl reference electrode

---

**CHROMATOGRAM****Retention time:** 20**Limit of detection:** 100 pg

---

**OTHER SUBSTANCES****Extracted:** zearalanone,  $\beta$ -zearalenol, zearalenone

---

**KEY WORDS**

chicken; cow; muscle; liver

---

**REFERENCE**

Roybal, J.E.; Munns, R.K.; Morris, W.J.; Hurlbut, J.A.; Shimoda, W. Determination of zeranol/zearalenone and their metabolites in edible animal tissue by liquid chromatography with electrochemical detection and confirmation by gas chromatography/mass spectrometry, *J. Assoc. Off. Anal. Chem.*, **1988**, 71, 263–271.

---

**SAMPLE****Matrix:** tissue

**Sample preparation:** Dry pack  $60 \times 8$  mm glass columns with 250 mg Carbowpack B (200-400 mesh) and  $60 \times 4$  mm glass columns with 50 mg Amberlite CG-400 I (100-200 mesh). Wash Carbowpack column with 5 mL MeOH, 15 mL dichloromethane:MeOH 70:30, and MeOH:water 85:15. Wash Amberlite column with 3 mL 0.5 M NaOH, 8 mL dichloromethane:MeOH 70:30, 1 mL water, and 3 mL 1 M HCl. Repeat this cycle 4 times. Finally pass through 20 mL 50 mM NaOH then 1 mL water. Keep column in water. (Process converts Amberlite to OH form.) Homogenize 1 g of tissue in 5 mL MeOH, sonicate 5 min, centrifuge at 6000 rpm for 10 min. Add another 5 mL MeOH to pellet and repeat. Combine supernatants, make up to 6.8 mL with MeOH, add 1.2 mL water. Pass through Carbowpack column, wash column with 2 mL MeOH:water 85:15 then 2 mL MeOH, elute column with 8 mL dichloromethane:MeOH 70:30. Pass eluate onto Amberlite column, wash with 1 mL MeOH, 1 mL 1 M HCl, elute with 2 mL 30 mM HCl in MeCN:MeOH 20:80. Evaporate eluate to dryness with nitrogen at 40°, take up in 100  $\mu$ L MeCN:MeOH:THF:10 mM  $\text{KH}_2\text{PO}_4$  adjusted to pH 3.0 with phosphoric acid 22:8:13:57, inject 50  $\mu$ L aliquot

---

**HPLC VARIABLES****Guard column:**  $20 \times 4.6$  5  $\mu$ m Supelguard LC-18**Column:**  $250 \times 4.6$  5  $\mu$ m Supelco C18**Mobile phase:** MeCN:MeOH:THF:10 mM  $\text{KH}_2\text{PO}_4$  adjusted to pH 3.0 with phosphoric acid 21:7:12:60**Flow rate:** 1.5**Injection volume:** 50**Detector:** E, Coulochem 5100A, detector 1 0.05 V, detector 2 0.60 V

---

**CHROMATOGRAM****Retention time:** 12**Limit of detection:** 1 ng/g

---

**OTHER SUBSTANCES****Simultaneous:** taleranol, zearalenol, zearalenone

---

**KEY WORDS**

muscle; liver; chicken; ox

---

**REFERENCE**

Laganà, A.; Marino, A. General and selective isolation procedure for high-performance liquid chromatographic determination of anabolic steroids in tissues, *J. Chromatogr.*, **1991**, 588, 89–98.

---

**SAMPLE****Matrix:** tissue

**Sample preparation:** Homogenize (Waring blender) tissue at full speed for 2 min, lyophilize, grind. Extract with supercritical carbon dioxide at 60° at 400 atmospheres with a  $20 \text{ cm} \times 21 \mu\text{m}$  restrictor for 1 h, collect the extract in 1 mL MeOH cooled to 5°. Evaporate the MeOH to dryness under a stream of nitrogen, reconstitute the residue in 100  $\mu$ L MeCN:MeOH:20 mM

ammonium formate 15:15:70, inject an aliquot. Alternatively, vortex 5 g ground tissue with 10 mL 40 mM sodium acetate, adjust pH to 4.2-4.7 with glacial acetic acid, add 100  $\mu$ L  $\beta$ -glucuronidase (Sigma), heat at 37° for 8 h, add 20 mL MeCN, vortex for 30 s, centrifuge at 5000 rpm for 20 min. Remove a 30 mL aliquot of the supernatant and add it to 8 mL hexane and 2 mL dichloromethane, rotate for 3 min, centrifuge at 2000 rpm for 2 min. Remove a 15 mL aliquot of the middle layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 1 mL dichloromethane, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 50  $\times$  4.6 5  $\mu$ m Supelcosil

**Mobile phase:** Gradient. MeCN:MeOH:20 mM ammonium formate from 2.5:2.5:95 to 47.5:47.5:5 over 19 min.

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 245 or MS, Sciex TAGA 6000E tandem triple quadrupole, APCI

---

**CHROMATOGRAM**

**Retention time:** 11.3

**Limit of detection:** 100 ppb

---

**OTHER SUBSTANCES**

**Extracted:** dexamethasone, diethylstilbestrol, medroxyprogesterone, melengestrol acetate, trenbolone, triamcinolone acetonide

---

**KEY WORDS**

cow; muscle; liver; SFE

---

**REFERENCE**

Huopalahti,R.P.; Henion,J.D. Application of supercritical fluid extraction and high performance liquid chromatography/mass spectrometry for the determination of some anabolic agents directly from bovine tissue samples, *J.Liq.Chromatogr.Rel.Technol.*, **1996**, 19, 69-87.

---

**SAMPLE**

**Matrix:** urine

**Sample preparation:** 10 mL Urine + glucuronidase/sulfatase (*Helix pomatia*), incubate at 37° for 1 h, extract twice with 5 mL diethyl ether, add 225  $\mu$ L water and evaporate ether under nitrogen, add 400  $\mu$ L MeOH, inject a 250  $\mu$ L aliquot of this mixture.

---

**HPLC VARIABLES**

**Guard column:** 75  $\times$  2.1 Corasil C18

**Column:** 150  $\times$  4.6 5  $\mu$ m Hypersil ODS

**Mobile phase:** MeOH:water 60:40

**Flow rate:** 2

**Injection volume:** 250

**Detector:** UV 240

---

**CHROMATOGRAM**

**Retention time:** 6.5

**Limit of detection:** about 6 ng/mL

---

**OTHER SUBSTANCES**

**Simultaneous:** 17 $\alpha$ -methyltestosterone, 17 $\beta$ -trenbolone, trans-diethylstilbestrol, medroxyprogesterone, nandrolone

---

**KEY WORDS**

cow

---

**REFERENCE**

Jansen,E.H.; Both-Miedema,R.; van Blitterswijk,H.; Stephany,R.W. Separation and purification of several anabolics present in bovine urine by isocratic high-performance liquid chromatography, *J.Chromatogr.*, **1984**, 299, 450-455.

# Zidovudine


**Molecular formula:** C<sub>10</sub>H<sub>13</sub>N<sub>5</sub>O<sub>4</sub>

**Molecular weight:** 267.24

**CAS Registry No.:** 30516-87-1

**Merck Index:** 10252

**Lednicer No.:** 4 118



## SAMPLE

**Matrix:** blood

**Sample preparation:** Condition a 3 mL 200 mg Bond Elut C18 SPE cartridge with 3 mL MeOH and 2 mL buffer, do not allow to go dry. 1 mL Plasma + 25  $\mu$ L 40  $\mu$ g/mL 7-ethyltheophylline in water, vortex, add to the SPE cartridge at  $\leq 0.5$  mL/min, wash with 1 mL buffer, air dry for 3 min, elute with 2 mL MeOH. Evaporate the eluate to dryness under a stream of nitrogen at 50°, reconstitute the residue in 100  $\mu$ L MeCN:5% acetic acid 10:90, vortex vigorously, centrifuge at 10000 g for 5 min, inject a 20  $\mu$ L aliquot of the supernatant. (Prepare buffer by diluting 1.44 mL concentrated phosphoric acid to 1 L with water and adjusting the pH to 6.55 with concentrated ammonia.)

## HPLC VARIABLES

**Column:** 75  $\times$  4.6 3  $\mu$ m Ultrasphere ODS

**Mobile phase:** Gradient. A was 1.44 mL concentrated phosphoric acid and 4 mL n-octylamine in 1 L water, pH adjusted to 6.55 with concentrated ammonia. B was MeCN. A:B from 95:5 to 70:30 over 7 min, to 20:80 over 1.5 min, return to initial conditions over 1 min, re-equilibrate for 2.5 min.

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 266

## CHROMATOGRAM

**Retention time:** 3.6

**Internal standard:** 7-ethyltheophylline (4.2)

**Limit of detection:** 7 ng/mL

**Limit of quantitation:** 22 ng/mL

## OTHER SUBSTANCES

**Extracted:** metabolites, caffeine

## KEY WORDS

plasma; SPE

## REFERENCE

Nadal,T; Ortuño,J.; Pascual,J.A. Rapid and sensitive determination of zidovudine and zidovudine glucuronide in human plasma by ion-pair high-performance liquid chromatography, *J.Chromatogr.A*, **1996**, 721, 127-137.

## SAMPLE

**Matrix:** blood

**Sample preparation:** Filter (Millipore Ultrafree-MC, 10000 molecular mass limit) 250  $\mu$ L serum while centrifuging at 17000 g for 1.5 h, inject a 50  $\mu$ L aliquot of the clear ultrafiltrate.

## HPLC VARIABLES

**Column:** 150  $\times$  3.9 4  $\mu$ m Nova-Pak phenyl

**Mobile phase:** Isopropanol:20 mM pH 5 sodium citrate 2.5:97.5

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 250

---

**CHROMATOGRAM****Retention time:** 22.4

---

**OTHER SUBSTANCES****Extracted:** didanosine, zalcitabine

---

**KEY WORDS**serum; ultrafiltrate

---

**REFERENCE**

Rosell-Rovira, M.L.; Pou-Clavé, L.; Lopez-Galera, R.; Pascual-Mostaza, C. Determination of free serum didanosine by ultrafiltration and high-performance liquid chromatography, *J.Chromatogr.B*, **1996**, 675, 89–92.

---

**SAMPLE****Matrix:** intestinal mucosal homogenate**Sample preparation:** Homogenate mixture + 100  $\mu$ L 250 mM NaCN, mix, centrifuge at 4° at 34000 g for 10 min, filter (0.45  $\mu$ m) the supernatant, inject an aliquot of the filtrate.

---

**HPLC VARIABLES****Column:** 150  $\times$  3.9 Nova-Pak C18**Mobile phase:** MeOH:100 mM potassium phosphate 25:75**Flow rate:** 1**Detector:** UV 254

---

**KEY WORDS**rat

---

**REFERENCE**

Sinko, P.J.; Hu, P. Determining intestinal metabolism and permeability for several compounds in rats. Implications on regional bioavailability in humans, *Pharm.Res.*, **1996**, 13, 108–113.

---

**SAMPLE****Matrix:** solutions**Sample preparation:** Inject a 15  $\mu$ L aliquot.

---

**HPLC VARIABLES****Column:** 200  $\times$  4.6 5  $\mu$ m HP Hypersil ODS**Mobile phase:** MeCN:20 mM pH 7.0 Na<sub>2</sub>HPO<sub>4</sub> 20:80**Column temperature:** 37**Flow rate:** 1**Injection volume:** 15**Detector:** UV 265

---

**CHROMATOGRAM****Retention time:** 4.18


---

**REFERENCE**

Kim, D.-D.; Chien, Y.W. Transdermal delivery of dideoxynucleoside-type Anti-HIV drugs. 2. The effect of vehicle and enhancer on skin permeation, *J.Pharm.Sci.*, **1996**, 85, 214–219.

---

# Zileuton

**Molecular formula:** C<sub>11</sub>H<sub>12</sub>N<sub>2</sub>O<sub>2</sub>S**Molecular weight:** 236.29**CAS Registry No.:** 111406-87-2**Merck Index:** 10253

---

**SAMPLE****Matrix:** blood**Sample preparation:** Condition a 300  $\mu$ L 200 mg Analytichem Bond Elut C8 SPE cartridge with two 2 mL portions of MeOH and two 2 mL portions of water. 1 mL Plasma + 0.1-1  $\mu$ g/mL IS in MeOH, mix, add to the SPE cartridge, wash with 1 mL water, dry for 3 min. Elute with two 1 mL portions of diethyl ether, evaporate under a gentle stream of air or nitrogen at 40°. Reconstitute the residue in 200  $\mu$ L mobile phase or water, inject a 50  $\mu$ L aliquot.

---

**HPLC VARIABLES****Column:** 250  $\times$  4.6 5  $\mu$ m Supelcosil LC-18**Mobile phase:** MeCN:MeOH:THF:water 15:5:10:70 (or 10:2:14:74 for rat blood) containing 40 mM NaH<sub>2</sub>PO<sub>4</sub>, 7.5 mM H<sub>3</sub>PO<sub>4</sub>, and 5 mM acetohydroxamic acid, pH 4.0**Flow rate:** 1.0-1.5**Injection volume:** 50**Detector:** UV 260

---

**CHROMATOGRAM****Retention time:** 15-19**Internal standard:** A-66649 (Abbott Laboratories, USA) (19-25)**Limit of quantitation:** 10-20 ng/mL

---

**OTHER SUBSTANCES****Extracted:** metabolites

---

**KEY WORDS**

plasma; SPE; human; monkey; rat; pharmacokinetics

---

**REFERENCE**Granneman; G.R.; Breackman; R.A.; Erdman; K.A. Determination of a new 5-lipoxygenase inhibitor, zileuton, and its inactive *N*-dehydroxylated metabolite in plasma by high performance liquid chromatography, *Clin.Pharmacokinet.*, **1995**, 29, 1-8.

---

**SAMPLE****Matrix:** blood**Sample preparation:** Condition Bond Elut C8 cartridge with MeOH and water. Add 1 mL plasma or diluted plasma to the SPE cartridge, aspirate, wash with water, elute with ether, evaporate the eluate, reconstitute in mobile phase, inject an aliquot.

---

**HPLC VARIABLES****Column:** 100  $\times$  4 5  $\mu$ m AGP 100.4**Mobile phase:** 50 mM Sodium perchlorate containing 17 mM acetic acid, adjusted to pH 2.0**Detector:** F ex 260 em 321

---

**CHROMATOGRAM****Internal standard:** Abbott-66649 (Abbott Laboratories, USA)**Limit of quantitation:** 10 ng/mL (R(+)), 10 ng/mL (S(-))

---

**KEY WORDS**

SPE; plasma; chiral; pharmacokinetics

---

**REFERENCE**Wong; S.L.; Awni; W.M.; Cavanaugh; J.H.; El-Shourbagy; T.; Locke; Ch.S.; Dubé; L.M. The pharmacokinetics of zileuton 200 to 800 mg, its enantiomers, and its metabolites in normal healthy volunteers, *Clin.Pharmacokinet.*, **1995**, 29, 9-21.

---

**SAMPLE****Matrix:** urine**Sample preparation:** Directly inject a urine sample.

---

**HPLC VARIABLES****Column:** 250  $\times$  4.6 5  $\mu$ m IB-sil CN-bonded silica (Phenomenex)

**Mobile phase:** Isopropanol:25 mM sodium dodecyl sulfate containing 10 mM phosphate buffer 3:97, pH 3.0  
**Column temperature:** 50  
**Flow rate:** 1  
**Injection volume:** 50  
**Detector:** UV 262

---

#### CHROMATOGRAM

**Retention time:** 26.19  
**Limit of quantitation:** 0.08-0.10 ppm

---

#### OTHER SUBSTANCES

**Extracted:** metabolites

---

#### REFERENCE

Thomas,S.B.; Albazi,S.J. Simultaneous determination of the 5-lipoxygenase inhibitor "Zileuton" and its N-dehydroxylated metabolite in untreated rat urine by micellar liquid chromatography, *J.Liq.Chromatogr.Rel.Technol.*, **1996**, 19, 977-991.

---

## Zinostatin

**Molecular weight:** ca. 16000  
**CAS Registry No.:** 9014-02-2, 123760-07-6  
**Merck Index:** 10302

---

#### SAMPLE

**Matrix:** bulk  
**Sample preparation:** Inject 200  $\mu$ L of a 660  $\mu$ g/mL solution in 15 mM pH 5 sodium acetate buffer.

---

#### HPLC VARIABLES

**Column:** Mono Q HR 5/5 anion exchange (Pharmacia)  
**Mobile phase:** Gradient. A was 20 mM pH 5 ammonium acetate. B was 20 mM pH 5 ammonium acetate containing 1 M NaCl. A:B 100:0 to 3:97 over 2.5 min, to 75:25 over 11 min, to 0:100 over 2 min.  
**Flow rate:** 2  
**Injection volume:** 200  
**Detector:** UV 280

---

#### CHROMATOGRAM

**Retention time:** 6

---

#### OTHER SUBSTANCES

**Simultaneous:** degradation products

---

#### KEY WORDS

protect from daylight

---

#### REFERENCE

Denklau,D.; Kohnlein,W.; Luders,G.; Stellmach,J. Isolation and fast purification of neocarzinostatin by FPLC -ion exchange chromatography, *Z.Naturforsch.[C]*, **1983**, 38, 939-942.

---

#### SAMPLE

**Matrix:** solutions

---

#### HPLC VARIABLES

**Column:** TSK G-3000SW (Toyo Soda)

**Mobile phase:** 10 mM pH 7.9 ammonium bicarbonate containing 30 mM NaCl

**Flow rate:** 1

**Detector:** UV 254, UV 280

---

## REFERENCE

Maeda,H.; Ueda,M.; Morinaga,T.; Matsumoto,T. Conjugation of poly(styrene-co-maleic acid) derivatives to the antitumor protein neocarzinostatin: pronounced improvements in pharmacological properties, *J.Med.Chem.*, **1985**, *28*, 455–461.

---


# Zipeprol

**Molecular formula:** C<sub>23</sub>H<sub>32</sub>N<sub>2</sub>O<sub>3</sub>

**Molecular weight:** 384.52

**CAS Registry No.:** 34758-83-3, 34758-84-4 (2.HCl)

**Merck Index:** 10303



---

## SAMPLE

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

## HPLC VARIABLES

**Guard column:** 20 mm long Symmetry C18

**Column:** 250 × 4.6 5 µm Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 205.2

---

## CHROMATOGRAM

**Retention time:** 13.45

---

## KEY WORDS

whole blood

---

## REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, *763*, 149–163.




# Ziprasidone

**Molecular formula:**  $C_{21}H_{21}ClN_4OS$

**Molecular weight:** 412.94

**CAS Registry No.:** 146939-27-7, 138982-67-9 ( $H_2O.HCl$ )

**Merck Index:** 10304



## SAMPLE

**Matrix:** blood, feces, urine

**Sample preparation:** Serum. 5 mL Serum + 10 mL MeCN, mix, centrifuge, wash the pellet with 2 mL MeCN, combine the supernatants, concentrate, reconstitute the residue with 200  $\mu$ L mobile phase, inject an 80  $\mu$ L aliquot. Urine. Condition a Sep-Pak C18 cartridge. 10 mL Urine + pH 5.0 acetate buffer, add to the SPE cartridge, wash with water, elute with MeOH, evaporate to dryness. Reconstitute the residue with 200  $\mu$ L MeOH:ammonium acetate 20:80, inject an 80  $\mu$ L aliquot.

## HPLC VARIABLES

**Column:** 250  $\times$  4.6 5  $\mu$ m YMC basic

**Mobile phase:** Gradient. MeOH:20 mM pH 5.0 ammonium acetate 10:90 for 10 min, to 80:20 over 50 min, maintain at 80:20 for 7 min, return to initial conditions over 8 min, re-equilibrate for 10 min.

**Injection volume:** 80

**Detector:** UV; Radioactivity,  $\beta$ -RAM; MS, Perkin-Elmer Sciex API III, ion spray interface at 6000 V, collision gas argon, collision energy 25 eV

## CHROMATOGRAM

**Retention time:** 48

**Limit of quantitation:** 0.5-50 ng/mL (serum)

## OTHER SUBSTANCES

**Extracted:** metabolites

## KEY WORDS

SPE; serum; pharmacokinetics; radiolabeled

## REFERENCE

Prakash,C.; Kamel,A.; Gummerus,J.; Wilner,K. Metabolism and excretion of a new antipsychotic drug, ziprasidone, in humans, *Drug Metab.Dispos.*, **1997**, *25*, 863-872.

# Zolpidem


**Molecular formula:**  $C_{19}H_{21}N_3O$

**Molecular weight:** 307.40

**CAS Registry No.:** 82626-48-0, 99294-93-6 ((+)-tartrate (2:1))

**Merck Index:** 10321

**Lednicer No.:** 4 162



## SAMPLE

**Matrix:** blood

**Sample preparation:** Mix 2 mL plasma, 2 mL pH 9.2 carbonate buffer, and 4 mL hexane:chloromethane 4:3, centrifuge at 3000 rpm for 10 min. Remove the upper organic phase and evaporate it to dryness under a stream of air. Reconstitute the residue with mobile phase. Inject a 100  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 300 × 4.6 5 µm C18 Spherisorb ODS-2

**Mobile phase:** MeOH:THF:buffer 30:65:5 (Buffer was 10 mM potassium dihydrogen phosphate containing 0.1% triethylamine, adjusted to pH 2.6 with orthophosphoric acid.)

**Flow rate:** 0.8

**Injection volume:** 100

**Detector:** UV 305

---

**CHROMATOGRAM**

**Retention time:** 5.5

**Limit of detection:** 20 ng/mL

---

**OTHER SUBSTANCES**

**Extracted:** zopiclone

**Simultaneous:** acepromazine, alprazolam, amisulpiride, amitriptyline, clobazam, clotiazepam, cyamemazine, desipramine, diazepam, flunitrazepam, haloperidol, imipramine, levomepromazine, lormetazepam, midazolam, nitrazepam, nortriptyline, prazepam

**Interfering:** bromazepam, lorazepam, niaprazine

---

**KEY WORDS**

plasma

---

**REFERENCE**

Stanke,F.; Jourdil,N.; Lauby,V.; Bessard,G. Zopiclone and zolpidem quantification in human plasma by high performance liquid chromatography with photodiode-array detection, *J.Liq.Chromatogr.Rel.Technol.*, **1996**, *19*, 2623–2633.

---

---

**SAMPLE**

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

**HPLC VARIABLES**

**Guard column:** 20 mm long Symmetry C18

**Column:** 250 × 4.6 5 µm Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 208.7

---

**CHROMATOGRAM**

**Retention time:** 11.882

---

**KEY WORDS**


whole blood

---

**REFERENCE**

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, *763*, 149–163.

# Zonisamide



**Molecular formula:** C<sub>8</sub>H<sub>9</sub>N<sub>2</sub>O<sub>3</sub>S

**Molecular weight:** 212.23

**CAS Registry No.:** 68291-97-4

**Merck Index:** 10323

**Lednicer No.:** 4 130

## SAMPLE

**Matrix:** blood

**Sample preparation:** 100  $\mu$ L Plasma + 250  $\mu$ L 500 mM pH 4.5 NaH<sub>2</sub>PO<sub>4</sub> + 50  $\mu$ L 40 mg/mL IS in MeOH + 7 mL dichloroethane, vortex for 1 min, let stand for a few min. Remove the organic layer and filter (Whatman No. 1 paper) it, evaporate the filtrate to dryness under a stream of air in a warm water bath, reconstitute the residue in 50  $\mu$ L mobile phase, inject a 20  $\mu$ L aliquot.

## HPLC VARIABLES

**Guard column:** 10  $\times$  4.4 30-38  $\mu$ m CO:Pell ODS

**Column:** 250  $\times$  4.4 Hypersil 5 MOS

**Mobile phase:** MeCN:buffer 58:100 (Buffer was 50 mL 1 M NaOH and 58 mL 1 M acetic acid made up to 1 L with water.)

**Flow rate:** 1.5

**Injection volume:** 20

**Detector:** UV 280

## CHROMATOGRAM

**Retention time:** 4

**Internal standard:** 3-sulfamoylmethyl-6-fluoro-1,2-benzisoxazole (5)

**Limit of quantitation:** 900 ng/mL

## OTHER SUBSTANCES

**Simultaneous:** carbamazepine, carbamazepine 10,11-epoxide, carbamazepine 10,11-diol

**Noninterfering:** acetaminophen, acetazolamide, N-acetylprocainamide, ampicillin, caffeine, cefuroxime, cimetidine, chloramphenicol, clobazam, desmethylclobazam, desmethylsuximide, ethosuximide, mefenamic acid, metronidazole, naproxen, paraxanthine, pentobarbital, phenobarbital, phenytoin, primidone, procainamide, salicylic acid, sulfamethazine, sulfamethoxazole, sulthiame, theophylline, thiopental

## KEY WORDS

plasma

## REFERENCE

Berry, D.J. Determination of zonisamide (3-sulphamoylmethyl-1,2-benzisoxazole) in plasma at therapeutic concentrations by high-performance liquid chromatography, *J. Chromatogr.*, **1990**, 534, 173-181.

## SAMPLE

**Matrix:** blood

**Sample preparation:** 500  $\mu$ L Serum + 600  $\mu$ L allobarbitol in 75 mM pH 6.8 buffer, add 200 units  $\beta$ -glucuronidase (Type VII-A from *E. coli*), incubate at 37° for 30 min, add 1 mL of the sample to an Extrelut-1 SPE cartridge, after 10 min elute with 2.5 mL MTBE, dry the eluate under a stream of nitrogen, dissolve the residue in 50  $\mu$ L MeOH:water 1:1, inject a 10  $\mu$ L aliquot.

## HPLC VARIABLES

**Column:** 25  $\times$  4  $\mu$ m Superspher RP-18e (Merck)

**Mobile phase:** MeOH:11.2 mM  $\beta$ -cyclodextrin in 20 mM KH<sub>2</sub>PO<sub>4</sub> 5:95

**Flow rate:** 0.8

**Injection volume:** 10

**Detector:** UV 210

---

**CHROMATOGRAM****Retention time:** 7**Internal standard:** allobarbitol (16)**Limit of detection:** 2.1 ng/mL

---

**OTHER SUBSTANCES****Simultaneous:** phenobarbital, mephobarbital

---

**KEY WORDS**serum; SPE

---

**REFERENCE**

Eto,S.; Noda,H.; Noda,A. Simultaneous determination of antiepileptic drugs and their metabolites, including chiral compounds, via  $\beta$ -cyclodextrin inclusion complexes by a column-switching chromatographic technique, *J.Chromatogr.B*, **1994**, 658, 385–390.

---

**SAMPLE****Matrix:** blood

**Sample preparation:** Condition a 1 mL Bond-Elut SPE cartridge (cat. no. 607101) with two 1 mL portions of MeOH and two 1 mL portions of water. Add 20  $\mu$ L serum, 20  $\mu$ L 100  $\mu$ g/mL IS in MeOH:water 50:50, and 800  $\mu$ L water to the SPE cartridge, wash with 1 mL water, elute with 250  $\mu$ L MeOH, inject a 40  $\mu$ L aliquot of the eluate.

---

**HPLC VARIABLES****Column:** 300  $\times$  3.9 10  $\mu$ m  $\mu$ Bondapak C18**Mobile phase:** MeCN:MeOH:water 17:20:63**Flow rate:** 1.4**Injection volume:** 40**Detector:** UV 246

---

**CHROMATOGRAM****Retention time:** 3.9**Internal standard:** N,N-dimethylzonisamide (8.4)**Limit of detection:** 100 ng/mL

---

**OTHER SUBSTANCES****Simultaneous:** carbamazepine, carbamazepine epoxide, phenobarbital, phenytoin**Noninterfering:** valproic acid

---

**KEY WORDS**serum; SPE

---

**REFERENCE**

Furuno,K.; Oishi,R.; Gomita,Y.; Eto,K. Simple and sensitive assay of zonisamide in human serum by high-performance liquid chromatography using a solid-phase extraction technique, *J.Chromatogr.B*, **1994**, 656, 456–459.

---

**SAMPLE****Matrix:** blood

**Sample preparation:** 200  $\mu$ L Serum + 1 mL MeCN, vortex for 30 s, centrifuge at 13400 g for 2 min, inject a 20  $\mu$ L aliquot of the organic layer.

---

**HPLC VARIABLES****Column:** 75  $\times$  4.6 TSK gel ODS-80Tm (Tosoh)**Mobile phase:** MeCN:5 mM pH 4.7  $\text{KH}_2\text{PO}_4$  31:69**Flow rate:** 1**Injection volume:** 20**Detector:** UV 210

---

**KEY WORDS**

serum; comparison with capillary electrophoresis

---

**REFERENCE**

Makino,K.; Goto,Y.; Sueyasu,M.; Futagami,K.; Kataoka,Y.; Oishi,R. Micellar electrokinetic capillary chromatography for therapeutic drug monitoring of zonisamide, *J.Chromatogr.B*, **1997**, 695, 417–425.

---

**SAMPLE**

**Matrix:** blood, tissue

**Sample preparation:** Blood. 200  $\mu$ L Serum, plasma, or whole blood + 300  $\mu$ L water + 1 mL 100 mM pH 6.0 phosphate buffer + 50  $\mu$ L 100  $\mu$ g/mL dimethylzonisamide in chloroform:EtOH 10:1 + 6 mL chloroform:EtOH 10:1, shake for 10 min, centrifuge at 3500 rpm for 10 min. Remove 5 mL of the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 100  $\mu$ L EtOH, centrifuge at 10000 rpm for 1 min, inject a 20  $\mu$ L aliquot. Tissue. Homogenize (glass homogenizer) rat brain with 100 mM pH 6.0 phosphate buffer. 500  $\mu$ L Homogenate + 50  $\mu$ L 100  $\mu$ g/mL dimethylzonisamide in chloroform:EtOH 10:1 + 1 mL MeCN, let stand for 30 min at room temperature, centrifuge at 14000 rpm for 1 min. Remove the supernatant and evaporate it to 200  $\mu$ L, centrifuge at 14000 rpm, inject a 20  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 Develosil ODS-7 (Chemco, Osaka)

**Mobile phase:** MeCN:isopropanol:1% acetic acid 11:10:70

**Column temperature:** 35

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 285

---

**CHROMATOGRAM**

**Internal standard:** dimethylzonisamide

**Limit of detection:** 50 ng/g (brain), 100 ng/mL (blood)

---

**KEY WORDS**

rat; serum; whole blood; brain; human; plasma; pharmacokinetics

---

**REFERENCE**

Nishiguchi,K.; Ohnishi,N.; Iwakawa,S.; Yagi,J.; Nakayama,S.; Takada,S.; Nakamura,H.; Yokoyama,T.; Okumura,K. Pharmacokinetics of zonisamide; saturable distribution into human and rat erythrocytes and into rat brain, *J.Pharmacobiodyn.*, **1992**, 15, 409–415.

---

**SAMPLE**

**Matrix:** urine

**Sample preparation:** Inject a 100  $\mu$ L aliquot of urine.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  10 10  $\mu$ m Econosil C18

**Mobile phase:** Gradient. MeCN:MeOH:50 mM ammonium acetate from 1.6:2.4:96 to 28:42:30 over 30 min.

**Flow rate:** 4

**Injection volume:** 100

**Detector:** UV 280

---

**CHROMATOGRAM**

**Retention time:** 21

---

**OTHER SUBSTANCES**

**Extracted:** metabolites

---

**KEY WORDS**

rat; preparative

---

**REFERENCE**

Stiff,D.D.; Zemaitis,M.A. Metabolism of the anticonvulsant agent zonisamide in the rat, *Drug Metab.Dispos.*, **1990**, 18, 888–894.


# Zopiclone

**Molecular formula:** C<sub>17</sub>H<sub>17</sub>ClN<sub>6</sub>O<sub>3</sub>

**Molecular weight:** 388.81

**CAS Registry No.:** 43200-80-2

**Merck Index:** 10324



## SAMPLE

**Matrix:** blood

**Sample preparation:** Condition a 500 mg Bond Elut LRC C2 SPE cartridge with 2 mL MeOH, 2 mL MeCN:water:trifluoroacetic acid 80:20:0.1, and 5 mL Tris buffer. Mix 100  $\mu$ L 100 ng/mL IS with 2 mL 50 mM Tris buffer adjusted to pH 7.5 with concentrated HCl and 500  $\mu$ L plasma. Add sample to the SPE cartridge and draw through by vacuum at 2 mL/min. Wash with 20 mL Tris buffer, dry with vacuum for 1 min. Elute with 2 mL MeCN:water:trifluoroacetic acid 80:20:0.1, evaporate the eluate to dryness under a stream of nitrogen at 40°. Reconstitute the residue in 1 mL MeOH, vortex for 40 s, centrifuge at 3000 rpm for 15 min, decant the supernatant, evaporate to dryness. Reconstitute the residue in 200  $\mu$ L mobile phase, vortex for 20 s. Inject a 50  $\mu$ L aliquot.

## HPLC VARIABLES

**Guard column:** 4  $\mu$ m Nova-Pak phenyl Guard-Pak

**Column:** 300  $\times$  4.6 5  $\mu$ m C18 Spherisorb ODS-2

**Mobile phase:** MeOH:THF:buffer 30:65:5 (Buffer was 10 mM potassium dihydrogen phosphate containing 0.1% triethylamine, adjusted to pH 2.6 with orthophosphoric acid.)

**Flow rate:** 0.8

**Injection volume:** 100

**Detector:** UV 305

## CHROMATOGRAM

**Retention time:** 4.5

**Limit of detection:** 20 ng/mL

## OTHER SUBSTANCES

**Extracted:** zolpidem

**Simultaneous:** acepromazine, alprazolam, amisulpride, amitriptyline, clobazam, clotiazepam, cyamemazine, desipramine, diazepam, flunitrazepam, haloperidol, imipramine, levomepromazine, lormetazepam, midazolam, nitrazepam, nortriptyline, prazepam

**Interfering:** bromazepam, lorazepam, niaprazine

## KEY WORDS

plasma; SPE

## REFERENCE

Stanke,F.; Jourdil,N.; Lauby,V.; Bessard,G. Zopiclone and zolpidem quantification in human plasma by high performance liquid chromatography with photodiode-array detection, *J.Liq.Chromatogr.Rel.Technol.*, **1996**, *19*, 2623-2633.

## SAMPLE

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 100  $\mu$ L 2  $\mu$ g/mL IS in MeCN + 1 mL 10 mM pH 8 sodium phosphate buffer + 10 mL dichloromethane:isopropanol 90:10, shake gently for 10 min, centrifuge at 2000 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 37°, reconstitute the residue in 120  $\mu$ L mobile phase, inject a 100  $\mu$ L aliquot on to a 250  $\times$  4.6 5  $\mu$ m Nucleosil silica column and elute with MeCN:MeOH 95:5 at 1 mL min, collect the fraction containing zopiclone (about 7.5 min), evaporate it to dryness under a stream of nitrogen at 37°, reconstitute the residue in 120  $\mu$ L mobile phase, inject a 100  $\mu$ L aliquot.

## HPLC VARIABLES

**Guard column:** cellulose carbamate (Société Française Chromato Colonne)

**Column:** 250 × 4.6 cellulose carbamate (Société Francaise Chromato Colonne)

**Mobile phase:** Hexane:EtOH 40:60

**Column temperature:** 35

**Flow rate:** 1

**Injection volume:** 100

**Detector:** F ex 300 em 470

---

#### CHROMATOGRAM

**Retention time:** 8 (-), 10 (+)

**Internal standard:** 6-(7-chloro-2-quinolyl)-5-hydroxy-6,7-dihydro-5H-pyrrolo[3,4-b]pyrazin-7-one 4-methyl-1-piperazine carboxylate (5.5 (on achiral system))

---

#### KEY WORDS

plasma; normal phase; pharmacokinetics; chiral

---

#### REFERENCE

Fernandez,C.; Baune,B.; Gimenez,F.; Thuillier,A.; Farinotti,R. Determination of zopiclone enantiomers in plasma by liquid chromatography using a chiral cellulose carbamate column, *J.Chromatogr.*, **1991**, 572, 195–202.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** 1 mL Whole blood + 10 µL 8.5 µg/mL IS in MeCN, let stand for 5 min, add 300 µL saturated borate buffer (pH 9.2), add 3 mL n-butyl chloride, vortex for 2 min, centrifuge at 2140 g for 10 min, repeat the extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100 µL MeCN, inject a 1-30 µL aliquot.

---

#### HPLC VARIABLES

**Column:** 220 × 4.6 5 µm Spheri 5 C18

**Mobile phase:** MeCN:buffer 60:40 (Buffer was 1.15 g (NH<sub>4</sub>)H<sub>2</sub>PO<sub>4</sub> and 1 mL triethylamine in 1 L water, pH 6.8.)

**Flow rate:** 1

**Injection volume:** 10-30

**Detector:** UV 305

---

#### CHROMATOGRAM

**Retention time:** 4

**Internal standard:** 6-(6-chloro-2-quinolyl)-5-hydroxy-6,7-dihydro-5H-pyrrolo[3,4-b]pyrazin-7-one 4-methyl-1-piperazine carboxylate (29481 R.P.) (6)

**Limit of quantitation:** 4 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** metabolites

---

#### KEY WORDS

whole blood

---

#### REFERENCE

Boniface,P.J.; Martin,I.C.; Nolan,S.L.; Tan,S.T. Development of a method for the determination of zopiclone in whole blood, *J.Chromatogr.*, **1992**, 584, 199–206.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** 0.5-1 mL Plasma + 250 ng hydroquinidine + 1 mL pH 8 phosphate buffer + 6 mL dichloromethane, mix, centrifuge at 1000 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 50 µL mobile phase, inject a 20 µL aliquot.

---

#### HPLC VARIABLES

**Column:** 150 × 4.6 5 µm Ultrasphere ODS C18

**Mobile phase:** MeCN:MeOH:THF:buffer 15:5:2:78 (Buffer was 10 mM trimethylamine adjusted to pH 2.5 with phosphoric acid.)

**Flow rate:** 1

**Injection volume:** 20

**Detector:** UV 305

---

#### CHROMATOGRAM

**Retention time:** 4.6

**Internal standard:** hydroquinidine (3.3)

**Limit of detection:** 5 ng/mL

---

#### OTHER SUBSTANCES

**Simultaneous:** quinidine, quinine

**Noninterfering:** benzodiazepines, phenothiazines, tricyclic antidepressants, zolpidem

---

#### KEY WORDS

plasma; pharmacokinetics

---

#### REFERENCE

Royer-Morrot,M.J.; Rambourg,M.; Jacob,I.; Bauer,P.; Royer,R.J. Determination of zopiclone in plasma using column liquid chromatography with ultraviolet detection, *J.Chromatogr.*, **1992**, 581, 297–299.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** 2 mL Plasma + 2 mL buffer + 5 mL chloroform:isopropanol:n-heptane 60:14:26, shake gently on a horizontal agitator for 10 min, centrifuge at 2800 g for 10 min. Remove the organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, inject a 50 µL aliquot. (Buffer was saturated ammonium chloride, diluted 25% with water, adjusted to pH 9.5 with 25% diluted ammonia solution.)

---

#### HPLC VARIABLES

**Column:** 300 × 3.9 4 µm Nova-Pak C18

**Mobile phase:** MeOH:THF:buffer 65:5:30 (Buffer was 10 mM KH<sub>2</sub>PO<sub>4</sub> adjusted to pH 2.6 with orthophosphoric acid. At the end of the day wash column with water at 0.8 mL/min for 1 h and MeOH at 0.8 mL/min for 1 h.)

**Column temperature:** 30

**Flow rate:** 0.8

**Injection volume:** 200

**Detector:** UV 305

---

#### CHROMATOGRAM

**Retention time:** 4.05

**Limit of detection:** 24.8 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** alpidem, zolpidem, suriclone

**Simultaneous:** p-nitrophenol, ketotifen, tiaprofenic acid, vincristine, sultopride, pyrimethamine, nimodipine

---

#### KEY WORDS

plasma

---

#### REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. High-performance liquid chromatographic assay with diode-array detection for toxicological screening of zopiclone, zolpidem, suriclone and alpidem in human plasma, *J.Chromatogr.*, **1993**, 616, 95–103.

---

#### SAMPLE

**Matrix:** blood

**Sample preparation:** 1 mL Plasma + 50 µL 100 µg/mL chlordiazepoxide hydrochloride in water + 100 µL 70 mM pH 8 phosphate buffer + 5 mL MTBE:isooctane 75:25, vortex for 30 s, cen-



trifuge at 1800 g for 5 min. Remove the organic layer and evaporate it to dryness under reduced pressure, reconstitute the residue in 75  $\mu$ L hexane:EtOH 20:80, inject a 50  $\mu$ L aliquot.

---

**HPLC VARIABLES**

**Column:** 250  $\times$  4.6 5  $\mu$ m Chiralcel OD-H

**Mobile phase:** Hexane:EtOH 40:60

**Flow rate:** 0.6

**Injection volume:** 50

**Detector:** F ex 300 em 470

---

**CHROMATOGRAM**

**Retention time:** 19 (-), 28 (+)

**Internal standard:** chlordiazepoxide (9.5)

**Limit of detection:** 1 ng/mL

**Limit of quantitation:** 2.5 ng/mL

---

**KEY WORDS**

plasma; chiral; pharmacokinetics

---

**REFERENCE**

Foster, R.T.; Caillé, G.; Ngoc, A.H.; Lemko, C.H.; Kherani, R.; Pasutto, F.M. Stereospecific high-performance liquid chromatographic assay of zopiclone in human plasma, *J. Chromatogr. B*, **1994**, 658, 161–166.

---

**SAMPLE**

**Matrix:** blood

**Sample preparation:** 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100  $\mu$ L mobile phase, centrifuge at 2800 g for 5 min, inject a 50  $\mu$ L aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

---

**HPLC VARIABLES**

**Column:** 300  $\times$  3.9 4  $\mu$ m NovaPack C18

**Mobile phase:** MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic))  $\text{KH}_2\text{PO}_4$  adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

**Column temperature:** 30

**Flow rate:** 0.8

**Injection volume:** 50

**Detector:** UV 305

---

**CHROMATOGRAM**

**Retention time:** 4.05

**Limit of detection:** <120 ng/mL

---

**KEY WORDS**

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; toloxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; almino-

profen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vindesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; metapramine; chlordiazepoxide; oxprenolol; warfarin; clorazepate; flecainide; phenacyclidine; thiopental; fenfluramine; metipranolol; triprolidine; naproxen; buprenorphine; verapamil; buspirone; tianeptine; midazolam; bupivacaine; carbinoxamine; loprazolam; cetirizine; chlorpheniramine; moperone; cibenzoline; medifoxamine; astemizole; vinblastine; nicardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextrometamide; fenopropfen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimoziide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

## REFERENCE

Tracqui, A.; Kintz, P.; Mangin, P. Systematic toxicological analysis using HPLC/DAD, *J. Forensic Sci.*, **1995**, *40*, 254-262.

## SAMPLE

**Matrix:** blood, gastric contents

**Sample preparation:** 1 mL Whole blood or gastric contents + 50  $\mu$ L 400  $\mu$ g/mL IS in MeOH + 1 mL EtOH + 5 drops 1 M pH 9 potassium carbonate + 2 mL water + 8 mL n-hexane:MTBE 25:75, rotate for 15 min, centrifuge at 2500 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen, reconstitute the residue in 100  $\mu$ L mobile phase, inject a 50  $\mu$ L aliquot.

## HPLC VARIABLES

**Column:** 150  $\times$  4.6 5  $\mu$ m Hypersil BDS C18

**Mobile phase:** Gradient. A was MeCN:MeOH:1.5 M ammonium acetate:water 10:10:3:77. B was MeCN:MeOH:1.5 M ammonium acetate:water 40:40:3:17. A:B from 95:5 to 50:50 over 20 min.

**Flow rate:** 1

**Injection volume:** 50

**Detector:** UV 230

## CHROMATOGRAM

**Retention time:** 15.6

**Internal standard:** N-allylnormetazocine (Sigma) (9)

## OTHER SUBSTANCES

**Extracted:** pentazocine

## KEY WORDS

whole blood

## REFERENCE

Van Bocxlaer, J.; Meyer, E.; Clauwaert, K.; Lambert, W.; Piette, M.; De Leenheer, A. Analysis of zopiclone (Imovane) in postmortem specimens by GC-MS and HPLC with diode-array detection, *J. Anal. Toxicol.*, **1996**, *20*, 52-54.

## SAMPLE

**Matrix:** blood, tissue

**Sample preparation:** Blood. 1 mL Whole blood + 10  $\mu$ L 8.5  $\mu$ g/mL IS in MeCN + 300  $\mu$ L saturated pH 9 borate buffer + 3 mL n-butyl chloride, vortex for 2 min, centrifuge, repeat the extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 40°, reconstitute the residue in 100  $\mu$ L MeCN, inject a 10-30  $\mu$ L aliquot. Tissue. Stir

10 g homogenized liver tissue, 40 mL Tris buffer, and 10 mg Subtilisin-Carlsberg protease at room temperature overnight, filter through glass wool. 1 mL Homogenate + 10  $\mu$ L 8.5  $\mu$ g/mL IS in MeCN + 300  $\mu$ L saturated pH 9 borate buffer + 3 mL n-butyl chloride, vortex for 2 min, centrifuge, repeat the extraction. Combine the organic layers and evaporate them to dryness under a stream of nitrogen at 40°, reconstitute the residue in 500  $\mu$ L hexane, extract with two 1 mL portions of MeCN. Combine the MeCN layers and evaporate them to dryness, reconstitute the residue in 100  $\mu$ L MeCN, inject a 10-30  $\mu$ L aliquot. (Prepare Tris buffer by dissolving 121 g Tris in 1 L water, adjust pH to 7.0 with concentrated sulfuric acid.)

---

#### HPLC VARIABLES

**Column:** 220  $\times$  4.6 5  $\mu$ m Spheri 5 silica

**Mobile phase:** MeCN:buffer 60:40 (Prepare buffer by dissolving 1.15 g (NH<sub>4</sub>)H<sub>2</sub>PO<sub>4</sub> in 1 L water and adding 1 mL triethylamine, pH 6.8.)

**Flow rate:** 1

**Injection volume:** 10-30

**Detector:** UV 305

---

#### CHROMATOGRAM

**Internal standard:** 6-(6-chloro-2-quinolyl)-7-[(4-methyl-1-piperazinyl)carbonyloxy]-6,7-dihydro [5H]pyrrolo[3,4-b]pyrazine-5-one

**Limit of detection:** 50 ng/g

---

#### KEY WORDS

whole blood; liver

---

#### REFERENCE

Boniface,P.J.; Russell,S.G.G. Two cases of fatal zopiclone overdose, *J.Anal.Toxicol.*, **1996**, 20, 131-133.

---

#### SAMPLE

**Matrix:** blood, urine

**Sample preparation:** Condition a 1 mL BondElut C18 SPE cartridge with one column volume of 1 M HCl, with 2 column volumes of MeOH, and with 1 column volume of water. Hemolyze 1 mL blood with 500  $\mu$ L water, centrifuge. Add 250  $\mu$ L 250 ng/mL harmane hydrochloride in 200 mM NaH<sub>2</sub>PO<sub>4</sub>, and 1 mL serum or hemolysate or 500  $\mu$ L urine to the SPE cartridge at 1 mL/min, wash with 2 column volumes of water, wash with two 500  $\mu$ L aliquots of MeCN (drain completely after each wash), elute with 250  $\mu$ L MeOH:35% perchloric acid 99:1 (remove the final portion of eluate by centrifuging for 20 s), inject a 10 (urine) or 25 (serum)  $\mu$ L aliquot of the eluate.

---

#### HPLC VARIABLES

**Guard column:** 15  $\times$  3.2 7  $\mu$ m RP-18 (Applied Biosystems)

**Column:** 150  $\times$  4.6 5  $\mu$ m Ultrasphere ODS C18

**Mobile phase:** MeCN:0.1% tetramethylammonium perchlorate 17:83, adjusted to pH 3.8 with 10% perchloric acid

**Flow rate:** 1.8

**Injection volume:** 10-25

**Detector:** F ex 320 em 520, UV 305

---

#### CHROMATOGRAM

**Retention time:** 12.2

**Internal standard:** harmane (8.3)

**Limit of quantitation:** 2 ng/mL (serum, F), 10 ng/mL (urine, F)

---

#### OTHER SUBSTANCES

**Extracted:** metabolites

**Noninterfering:** acetaminophen, atenolol, barbiturates, benzodiazepines, chlordiazepoxide, fluoxetine, imipramine, metoprolol, nadolol, paroxetine, salicylic acid

---

#### KEY WORDS

serum; SPE

---

**REFERENCE**

Gupta,R.N. Simultaneous determination of zopiclone and its two major metabolites (N-oxide and N-desmethyl) in human biological fluids by column liquid chromatography after solid-phase extraction, *J.Liq. Chromatogr.Rel. Technol.*, **1996**, *19*, 699–709.

---

**SAMPLE**

**Matrix:** formulations

**Sample preparation:** Grind tablets to a fine powder, weigh out amount equivalent to 5 mg zopiclone, add 40 mL 100 mM HCl, sonicate for 15 min, cool to room temperature, make up to 50 mL with 100 mM HCl, filter (1.6  $\mu$ m glass fiber, Whatman GF/A), inject a 20  $\mu$ L aliquot of the filtrate.

---

**HPLC VARIABLES**

**Column:** 125  $\times$  4.5  $\mu$ m LiChrospher-60 RP Select B

**Mobile phase:** MeCN:THF:buffer 18:1:81 (Buffer was 3.4 g monosodium hexane sulfonate and 7.0 g  $\text{NaH}_2\text{PO}_4 \cdot 2\text{H}_2\text{O}$  in 1 L water, pH 4.55.)

**Column temperature:** 25

**Flow rate:** 1.5

**Injection volume:** 20

**Detector:** UV 303

---

**CHROMATOGRAM**

**Retention time:** 7.5

**Limit of detection:** 0.05% (of zopiclone)

---

**OTHER SUBSTANCES**

**Simultaneous:** degradation products, impurities

**Noninterfering:** excipients

---

**KEY WORDS**

protect from light; tablets; rugged; stability-indicating

---

**REFERENCE**

Bounine,J.P.; Tardif,B.; Beltran,P.; Mazzo,D.J. High-performance liquid chromatographic stability-indicating determination of zopiclone in tablets, *J.Chromatogr.A*, **1994**, *677*, 87–93.

---

**SAMPLE**

**Matrix:** solutions

---

**HPLC VARIABLES**

**Guard column:** 30  $\times$  3.2  $\mu$ m SI 100 ODS (not commercially available)

**Column:** 150  $\times$  3.2  $\mu$ m SI 100 ODS (not commercially available)

**Mobile phase:** MeCN:buffer 31.2:68.8 (Buffer was 6.66 g  $\text{KH}_2\text{PO}_4$  and 4.8 g 85% phosphoric acid in 1 L water, pH 2.3.)

**Flow rate:** 0.5–1

**Detector:** UV 212, 300

---

**CHROMATOGRAM**

**Retention time:** 1.4

**Internal standard:** 5-(4-methylphenyl)-5-phenylhydantoin (7.3)

---

**OTHER SUBSTANCES**

**Also analyzed:** aspirin, caffeine, carbamazepine, chlordiazepoxide, chlorprothixene, clonazepam, diazepam, doxylamine, ethosuximide, furosemide, haloperidol, hydrochlorothiazide, methocarbamol, methotrimeprazine, nicotine, oxazepam, procaine, promazine, propafenone, propranolol, salicylamide, temazepam, tetracaine, thiopental, triamterene, verapamil, zolpidem

---

**REFERENCE**

Below,E.; Burrmann,M. Application of HPLC equipment with rapid scan detection to the identification of drugs in toxicological analysis, *J.Liq. Chromatogr.*, **1994**, *17*, 4131–4144.

## SAMPLE

**Matrix:** solutions

## HPLC VARIABLES

**Column:** 250 × 4.6 5 μm Supelcosil LC-DP (A) or 250 × 4 5 μm LiChrospher 100 RP-8 (B)

**Mobile phase:** MeCN:0.025% phosphoric acid:buffer 25:10:5 (A) or 60:25:15 (B) (Buffer was 9 mL concentrated phosphoric acid and 10 mL triethylamine in 900 mL water, adjust pH to 3.4 with dilute phosphoric acid, make up to 1 L.)

**Flow rate: 0.6**

**Injection volume: 25**

**Detector:** UV 229

## CHROMATOGRAM

**Retention time:** 7.50 (A), 3.79 (B)

## OTHER SUBSTANCES

**Also analyzed:** acebutolol, acepromazine, acetaminophen, acetazolamide, acetophenazine, albuterol, alprazolam, amitriptyline, amobarbital, amoxapine, antipyrine, atenolol, atropine, azatadine, baclofen, benzocaine, bromocriptine, brompheniramine, brotizolam, bupivacaine, buspirone, butabarbital, butalbital, caffeine, carbamazepine, cetirizine, chlorcyclizine, chlordi-azepoxide, chlormezanone, chloroquine, chlorpheniramine, chlorpromazine, chlorpropamide, chlorprothixene, chlorthalidone, chlorzoxazone, cimetidine, cisapride, clomipramine, clonazepam, clonidine, clozapine, cocaine, codeine, colchicine, cyclizine, cyclobenzaprine, dantrolene, desipramine, diazepam, diclofenac, diflunisal, diltiazem, diphenhydramine, diphenidol, diphenoxylate, dipyrindamole, disopyramide, dobutamine, doxapram, doxepin, droperidol, encainide, ethidium bromide, ethopropazine, fenopropfen, fentanyl, flavoxate, fluoxetine, fluphenazine, flurazepam, flurbiprofen, fluvoxamine, furosemide, glutethimide, glyburide, guaifenesin, haloperidol, homatropine, hydralazine, hydrochlorothiazide, hydrocodone, hydromorphone, hydroxy-chloroquine, hydroxyzine, ibuprofen, imipramine, indomethacin, ketoconazole, ketoprofen, ketorolac, labetalol, levorphanol, lidocaine, loratadine, lorazepam, lovastatin, loxapine, mazin-dol, mefenamic acid, meperidine, mephénytoin, mepivacaine, mesoridazine, metaproterenol, metformin, methadone, methdilazine, methocarbamol, methotrexate, methotrimeprazine, methoxamine, methyl dopa, methylphenidate, metoclopramide, metolazone, metoprolol, met-ronidazole, midazolam, moclobemide, morphine, nadolol, nalbuphine, naloxone, naphazoline, naproxen, nifedipine, nizatidine, norepinephrine, nortriptyline, oxazepam, oxycodone, oxymet-azoline, paroxetine, pemoline, pentazocine, pentobarbital, pentoxifylline, perphenazine, phen-iramine, phenobarbital, phenol, phenolphthalein, phentolamine, phenylbutazone, phenyltolox-amine, phenytoin, pimizide, pindolol, piroxicam, pramoxine, prazepam, prazosin, probenecid, procainamide, procaine, prochlorperazine, procyclidine, promazine, promethazine, propafenone, propantheline, propiomazine, propofol, propranolol, protriptyline, quazepam, quinidine, quin-ine, racemethorphan, ranitidine, remoxipride, risperidone, salicylic acid, scopolamine, seco-barbital, sertraline, sotalol, spirinolactone, sulfapyrazone, sulindac, temazepam, terbutaline, terfenadine, tetracaine, theophylline, thiethylperazine, thiopental, thioridazine, thiothixene, timolol, tocainide, tolbutamide, tolmetin, trazodone, triamterene, triazolam, trifluoperazine, triflupromazine, trimeprazine, trimethoprim, trimipramine, verapamil, warfarin, xylometazo-line, yohimbine

## KEY WORDS

details of plasma extraction

## REFERENCE

Koves, E.M. Use of high-performance liquid chromatography-diode array detection in forensic toxicology, *J. Chromatogr. A*, **1995**, 692, 103–119.

## SAMPLE

**Matrix:** urine

**Sample preparation:** Dilute urine 5-fold with 50 mM pH 8.4 phosphate buffer. 1 mL Diluted urine + 100  $\mu$ L 20  $\mu$ g/mL suriclone in EtOH + 10 mL dichloromethane:isopropanol 95:5, shake gently for 5 min, centrifuge at 2000 g for 10 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 37°, reconstitute the residue in 120  $\mu$ L EtOH, inject a 100  $\mu$ L aliquot on to column A and elute to waste with mobile phase A, after 11.5 min elute the fraction containing zopiclone on to column B, after 1.5 min backflush the contents of column

B on to column C with mobile phase B, elute with mobile phase B, monitor the effluent from column C. (Chiral separation of metabolites can also be achieved by collecting the fractions containing the metabolites (20.5 to 25.5 min and 46.5 to 55.5 min) on column B then chromatographing them on column C.)

#### HPLC VARIABLES

**Column:** A 150 × 4.6 5 µm Nucleosil cyanopropyl; B 60 × 4 5 µm Kromasil silica (Informatiques & Technologies); C cellulose carbamate guard column + 250 × 4.6 5 µm cellulose carbamate (Société Française Chromato Colonne)

**Mobile phase:** A Hexane:EtOH:MeOH 80:5:15 containing 0.18% MeOH:diethylamine 99.9:0.1 and 0.05% water; B Hexane:EtOH:MeOH:diethylamine 55:30:15:1

**Column temperature:** 35 (column C only)

**Flow rate:** A 0.7; B 1.2

**Injection volume:** 100

**Detector:** F ex 305 em 470

#### CHROMATOGRAM

**Retention time:** 20 (-), 25 (+)

**Internal standard:** suriclone (18 (on column A only))

**Limit of quantitation:** 10 ng

#### OTHER SUBSTANCES

**Extracted:** metabolites

#### KEY WORDS

chiral; column-switching; pharmacokinetics; heart-cut

#### REFERENCE

Fernandez,C.; Gimenez,F.; Baune,B.; Maradeix,V.; Thuillier,A.; Farinotti,R. Determination of the enantiomers of zopiclone and its two chiral metabolites in urine using an automated coupled achiral-chiral chromatographic system, *J.Chromatogr.*, **1993**, 617, 271–278.


## Zorubicin

**Molecular formula:** C<sub>34</sub>H<sub>35</sub>N<sub>3</sub>O<sub>10</sub>

**Molecular weight:** 645.67

**CAS Registry No.:** 54083-22-6, 36508-71-1 (HCl)

**Merck Index:** 10326



#### SAMPLE

**Matrix:** blood

**Sample preparation:** 2 mL Whole blood or plasma + 2 mL buffer + 5 mL chloroform:isopropanol: n-heptane 60:14:26, shake gently horizontally for 10 min, centrifuge at 2800 g for 10 min. Remove the lower organic layer and evaporate it to dryness under vacuum at 45°, reconstitute the residue in 100 µL mobile phase, centrifuge at 2800 g for 5 min, inject a 50 µL aliquot of the supernatant. (Buffer was saturated ammonium chloride solution 25% diluted with water, adjusted to pH 9.5 with 25% ammonia solution.)

#### HPLC VARIABLES

**Column:** 300 × 3.9 4 µm NovaPack C18

**Mobile phase:** MeOH:THF:buffer 65:5:30 (Buffer was 0.68 g/L (10 mM (sic)) KH<sub>2</sub>PO<sub>4</sub> adjusted to pH 2.6 with concentrated orthophosphoric acid.) (At the end of each session wash the column with water for 1 h and MeOH for 1 h, re-equilibrate for 30 min.)

**Column temperature:** 30

**Flow rate:** 0.8

**Injection volume:** 50

**Detector:** UV 234

## CHROMATOGRAM

**Retention time:** 8.10

**Limit of detection:** <120 ng/mL


## KEY WORDS

whole blood; plasma; interferences may occur—compounds (all of which are extracted) elute in this order tenoxicam; iproniazid; methocarbamol; methotrexate; caffeine; nialamide; colchicine; cytarabine; benzoylecgonine; acetaminophen; diazoxide; dacarbazine; sulfinpyrazole; flumazenil; sulpride; morphine; atenolol; tolaxatone; terbutaline; albuterol; phenobarbital; ranitidine; tiapride; phenol; chlormezanone; aspirin; metformin; ritodrine; codeine; sultopride; amisulpride; naltrexone; lisinopril; benzocaine; nizatidine; nalorphine; mephenesin; naloxone; sotalol; carteolol; procainamide; carbamazepine; bromazepam; nalbuphine; nadolol; procarbazine; dihydralazine; omeprazole; strychnine; acebutolol; glutethimide; chlorpropamide; glipizide; triazolam; prazosin; flunitrazepam; clonazepam; metoclopramide; melphalan; estazolam; tolbutamide; ephedrine; clonidine; pindolol; clobazam; minoxidil; disopyramide; nitrazepam; dextromethorphan; tofisopam; zopiclone; debrisoquine; sulindac; alprazolam; cycloguanil; lorazepam; methaqualone; ketamine; piroxicam; metoprolol; nifedipine; quinine; mephentermine; prilocaine; pentazocine; oxazepam; tiaprofenic acid; quinidine; celiprolol; ajmaline; yohimbine; lidocaine; secobarbital; viloxazine; mepivacaine; meperidine; doxylamine; labetalol; temazepam; amodiaquine; benperidol; droperidol; hydroxychloroquine; zolpidem; ketoprofen; alminoprofen; cicletanine; moclobemide; chloroquine; cocaine; timolol; nomifensine; ticlopidine; acenocoumarol; vindesine; mexiletine; dipyridamole; trazodone; pipamperone; pyrimethamine; benazepril; vincristine; nocardipine; bisoprolol; diltiazem; glibornuride; reserpine; aconitine; nitrendipine; diazepam; mianserin; ramipril; haloperidol; tetracaine; alprenolol; aceprometazine; glibenclamide; chlorophenacinone; doxepin; nimodipine; diphenhydramine; cyclizine; histapyrodine; phenylbutazone; demexiptiline; clozapine; proguanil; trifluoperidol; medazepam; cyamemazine; bumadizone; suriclone; propranolol; acepromazine; dothiepin; dextromoramide; fenopropfen; dextropropoxyphene; loxapine; betaxolol; propafenone; promethazine; thioproperazine; methadone; amoxapine; quinupramine; opiipramol; cyproheptadine; brompheniramine; mefenidramine; protriptyline; flurbiprofen; tetrazepam; zorubicin; prazepam; alimemazine; loperamide; imipramine; desipramine; levomepromazine; hydroxyzine; niflumic acid; penbutolol; fluvoxamine; pimozide; daunorubicin; indomethacin; maprotiline; tropatenine; etodolac; fluoxetine; amitriptyline; nortriptyline; tiocloamarol; diclofenac; mefloquine; trimipramine; chlorambucil; lidoflazine; ibuprofen; floctafenine; alpidem; loratadine; chlorpromazine; clomipramine; carpipramine; thioridazine; fentiazac; clemastine; mefenamic acid; fluphenazine; prochlorperazine; penfluridol; bepridil; terfenadine; trifluoperazine

## REFERENCE

Tracqui,A.; Kintz,P.; Mangin,P. Systematic toxicological analysis using HPLC/DAD, *J.Forensic Sci.*, **1995**, *40*, 254–262.

# Zoxazolamine



**Molecular formula:** C<sub>7</sub>H<sub>5</sub>ClN<sub>2</sub>O

**Molecular weight:** 168.58

**CAS Registry No.:** 61-80-3

**Merck Index:** 10328

## SAMPLE

**Matrix:** solutions

## HPLC VARIABLES

**Column:** 250 × 4.6 Zorbax RX

**Mobile phase:** Gradient. A was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 1 L water. B was 10 mL concentrated orthophosphoric acid and 7 mL triethylamine in 200 mL water, make up to 1 L with MeCN. A:B from 100:0 to 0:100 over 30 min, maintain at 0:100 for 5 min.

**Column temperature:** 30

**Flow rate:** 2

**Detector:** UV 210

## OTHER SUBSTANCES

**Also analyzed:** acepromazine, acetaminophen, acetophenazine, albuterol, aminophylline, am-triptyline, amobarbital, amoxapine, amphetamine, amylocaine, antipyrine, aprobarbital, aspi-rin, atenolol, atropine, avermectin, barbital, benzocaine, benzoic acid, benzotropine, benzphet-amine, berberine, bibucaine, bromazepam, brompheniramine, buprenorphine, buspirone, butabarbital, butacaine, butethal, caffeine, carbamazepine, carbromal, chloramphenicol, chlor-diazepoxide, chloroquine, chlorothiazide, chloroxylenol, chlorphenesin, chlorpheniramine, chlorpromazine, chlorpropamide, chlortetracycline, cimetidine, cinchonidine, cinchonine, clen-butanol, clonazepam, clonixin, clorazepate, cocaine, codeine, colchicine, cortisone, coumarin, cyclazocine, cyclobenzaprine, cyclothiazide, cyheptamide, cymarin, danazol, danthron, dapsone, debrisoquine, desipramine, dexamethasone, dextromethorphan, dextropropoxyphene, diamor-phine, diazepam, diclofenac, diethylpropion, diethylstilbestrol, diflunisal, digitoxin, digoxin, dil-tiazem, diphenhydramine, diphenoxylate, diprenorphine, dipyrone, disulfiram, dopamine, dox-apram, doxepin, dronabinol, ephedrine, epinephrine, epinine, estradiol, estriol, estrone, ethacrynic acid, ethosuximide, etonitazene, etorphine, eugenol, famotidine, fenbendazole, fen-camfamine, fenopropfen, fenproporex, fentanyl, flubendazole, flufenamic acid, flunitrazepam, 5-fluorouracil, fluoxymesterone, fluphenazine, furosemide, gentisic acid, gitoxigenin, glipizide, glunixin, glutethimide, glybenclamide, guaiaicol, halazepam, haloperidol, hydrochlorothiazide, hydrocodone, hydrocortisone, hydromorphone, hydroxyquinoline, ibogaine, ibuprofen, imino-stilbene, imipramine, indomethacin, isocarbostyryl, isocarboxazid, isoniazid, isoproterenol, isox-suprine, ivermectin, ketamine, ketoprofen, kynurenic acid, levorphanol, lidocaine, lorazepam, lormetazepam, loxapine, mazindol, mebendazole, meclizine, meclofenamic acid, medazepam, mefenamic acid, megestrol, mepacrine, meperidine, mephentermine, mephenytoin, mephesin, mephobarbital, mepivacaine, mescaline, mesoridazine, methadone, methamphetamine, meth-apyrilene, methaqualone, methazolamide, methocarbamol, methoxamine, methsuximide, methyl salicylate, methyldopa, methyldopamine, methylphenidate, methylprednisolone, meth-yltestosterone, methyprylon, metoprolol, mibolerone, morphine, nadolol, nalorphine, naloxone, naltrexone, naphazoline, naproxen, nefopam, niacinamide, nicotine, niacin, nifedipine, niflumic acid, nitrazepam, norepinephrine, nortriptyline, noscapine, nyldrin, oxazepam, oxycodone, ox-ymorphone, oxyphenbutazone, oxytetracycline, papaverine, pargyline, pemoline, pentazocine, pentobarbital, persantine, phenacetin, phenazocine, phenazopyridine, phencyclidine, phendi-metrazine, phenelzine, pheniramine, phenobarbital, phenothiazine, phensuximide, phenter-mine, phenylbutazone, phenylephrine, phenylpropanolamine, piperocaine, prazepam, predni-solone, primidone, probenecid, progesterone, propiomazine, propranolol, propylparaben, pseudoephedrine, puromycin, pyrilamine, pyridylidone, quazepam, quinaldic acid, quinidine, quinine, ranitidine, recinnamine, reserpine, resorcinol, saccharin, albuterol, salicylamide, sal-icylic acid, scopolamine, scopoletin, secobarbital, strychnine, sulfacetamide, sufadiazine, sul-fadimethoxine, sulfaethidole, sulfamerazine, sulfamethazine, sulfamethoxazole, sulfanilamide, sulfapyridine, sulfasoxazole, sulindac, tamoxifen, temazepam, testosterone, tetracaine, tetra-cycline, tetramisole, thebaine, theobromine, theophylline, thiabendazole, thiamine, thiamylal, thiobarbituric acid, thioridazine, thiosalicylic acid, thiothixene, thymol, tolazamide, tolazoline, tobutamide, tolmetin, tranlycypromine, triamcinolone, tribenzylamine, trichloromethiazide, trifluoperazine, trihexyphenidyl, trimethoprim, tripeleminamine, triprolidine, tropacocaine, ty-ramine, verapamil, vincamine, warfarin

## REFERENCE

Hill,D.W.; Kind,A.J. Reversed-phase solvent gradient HPLC retention indexes of drugs, *J.Anal.Toxicol.*, **1994**, *18*, 233-242.




# Zuclopenthixol

**Molecular formula:** C<sub>22</sub>H<sub>25</sub>ClN<sub>2</sub>OS

**Molecular weight:** 400.97

**CAS Registry No.:** 53772-83-1

**Merck Index:** 2455



## SAMPLE

**Matrix:** blood

**Sample preparation:** 2 mL Serum + 25 ng IS + 300 µL EtOH + 100 µL 7 M NaOH + 8 mL hexane:isopropylamine 99.9:0.1, shake for 15 min, centrifuge at 2400 g for 5 min. Remove the organic layer and add it to 2 mL 100 mM HCl, shake for 15 min, centrifuge for 5 min. Discard the organic layer and add 200 µL 7 M NaOH and 4 mL hexane:isopropylamine 99.9:0.1 to the aqueous layer, shake for 15 min, centrifuge for 5 min. Remove the organic layer and evaporate it to dryness under a stream of air at 30°, reconstitute the residue in 1 mL hexane, evaporate, reconstitute with 100 µL hexane:isopropylamine 99.9:0.1, inject a 70 µL aliquot.

## HPLC VARIABLES

**Column:** 250 × 4.6 5 µm Spherisorb S5W silica

**Mobile phase:** n-Heptane:isopropanol:concentrated ammonia:water 85:15:0.4:0.2

**Flow rate:** 1

**Injection volume:** 70

**Detector:** UV 254

## CHROMATOGRAM

**Retention time:** 7.5

**Internal standard:** 4-[3-(2-chloro-7-trifluoromethylthioxanthen-9-yl)propyl]-1-piperazineethanol (Lu 9-215) (6.5)

**Limit of detection:** 0.5 ng/mL

## OTHER SUBSTANCES

**Extracted:** trans(E)-isomer, metabolites

**Simultaneous:** estazolam

**Noninterfering:** amitriptyline, biperiden, imipramine, nortriptyline, orphenadrine, procyclidine

**Interfering:** benzodiazepines

## KEY WORDS

serum; rat; dog; human; pharmacokinetics; normal phase

## REFERENCE

Aaes-Jorgensen, T. Specific high-performance liquid chromatographic method for estimation of *cis*(Z) and *trans*(E)-isomers of clopenthixol and a N-dealkyl metabolite, *J. Chromatogr.*, **1980**, *183*, 239–245.

## SAMPLE

**Matrix:** blood

**Sample preparation:** Condition a 3 mL Bond Elut Certify SPE cartridge with 2 mL MeOH and 2 mL 100 mM pH 6.0 phosphate buffer, do not allow to dry. 1 mL Whole blood + 6 mL 100 mM pH 6.0 phosphate buffer, vortex, sonicate, centrifuge, add the supernatant to the SPE cartridge, wash with water, wash with 1 mM pH 3.3 acetic acid, dry under suction, wash with 2 mL acetone:chloroform 50:50, elute with 3 mL freshly prepared ethyl acetate:ammonia 98:2. Evaporate the eluate under a stream of nitrogen at 40°, reconstitute in 50 µL MeOH, inject a 10 µL aliquot.

## HPLC VARIABLES

**Column:** 125 × 4 5 µm Asahipak ODP-50

**Mobile phase:** MeCN:50 mM ammonium acetate 85:15

**Flow rate:** 0.6

**Injection volume:** 10

**Detector:** MS, Finnigan MAT TSQ 700 tandem quadrupole, Finnigan MAT TSP-2 interface, collision gas argon 3.0 mTorr, collision offset -17.5 V, repeller 70 V, vaporizer 130-5°, source 200°, filament off, multiplier 1500 V, dynode power 15 kV, scantime 1.20 s, MS/MS factor 0, monitor 316–271. (The effluent from the column was mixed with 50 mM ammonium acetate pumped at 0.6 mL/min. The mixture flowed to the detector.)

---

**CHROMATOGRAM**

**Retention time:** 3.40

**Limit of detection:** 2 ng

---

**OTHER SUBSTANCES**

**Extracted:** chlorprothixene, flupenthixol, thiothixene

---

**KEY WORDS**

whole blood; SPE

---

**REFERENCE**

Verweij, A.M.A.; Hordijk, M.L.; Lipman, P.J.L. Quantitative liquid chromatography, thermospray/tandem mass spectrometric (LC/TSP/MS/MS) analysis of some tranquilizers of the thioxanthene group in whole-blood, *J.Liq.Chromatogr.*, **1994**, *17*, 4009–4110.

---

**SAMPLE**

**Matrix:** blood, gastric contents, tissue, urine, vitreous humor

**Sample preparation:** Homogenize tissue with 4 volumes water. Extract 3 mL Blood, gastric contents, urine, vitreous humor, or homogenized tissue with 1.5 mL saturated pH 9.5 ammonium chloride buffer and 5 mL chloroform:2-propanol:n-heptane 25:10:65 for 10 min. (Caution! Chloroform is a carcinogen!). Centrifuge at 3500 g for 10 min, evaporate the organic layer at 45°. Reconstitute with 30 µL MeOH. Centrifuge at 10000 g for 5 min, remove 20 µL of the supernatant, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 150 × 2.4 µm NovaPak C18

**Mobile phase:** MeOH:2mM pH 3 ammonium acetate buffer 90:10

**Flow rate:** 0.2

**Injection volume:** 2

**Detector:** MS, PE Sciex API 100 double quadrupole, nebulizing gas nitrogen at 1.6 mL/min, curtain gas nitrogen at 1.08 mL/min, sprayer electrode +4.5 kV, electron multiplier +2.4 kV, m/z 100–450

---

**CHROMATOGRAM**

**Retention time:** 14.95 ((Z)-cis), ((E)-trans)

---

**KEY WORDS**

liver; kidney; lung; brain; skeletal muscle; comparison with HPLC/UV

---

**REFERENCE**

Tracqui, A.; Kintz, P.; Cirimele, V.; Berthault, F.; Mangin, P.; Ludes, B. HPLC-DAD and HPLC-MS findings in a fatality involving (Z)-cis-clopenthixol (zuclopenthixol), *J.Anal.Toxicol.*, **1997**, *21*, 314–318.

---

**SAMPLE**

**Matrix:** blood, gastric contents, tissue, urine, vitreous humor

**Sample preparation:** Homogenize tissue with 4 volumes water. 3 mL Blood, gastric contents, urine, vitreous humor, or homogenized tissue + 4 µg IS, extract with 1.5 mL saturated pH 9.5 ammonium chloride buffer and 5 mL chloroform:2-propanol:n-heptane 25:10:65 for 10 min. (Caution! Chloroform is a carcinogen!). Centrifuge at 3500 g for 10 min, evaporate the organic layer at 45°. Reconstitute with 30 µL MeOH. Centrifuge at 10000 g for 5 min, remove 20 µL of the supernatant, inject an aliquot.

---

**HPLC VARIABLES**

**Column:** 300 × 3.9 4 µm NovaPak C18

**Mobile phase:** MeOH:THF:10 mM pH 2.6 KH<sub>2</sub>PO<sub>4</sub> buffer 65:5:30

**Flow rate:** 0.8  
**Injection volume:** 7  
**Detector:** UV 228

---

#### CHROMATOGRAM

**Retention time:** 8.98 ((Z)-cis), 10.21 ((E)-trans)  
**Internal standard:** prochlorperazine  
**Limit of detection:** 7 ng/mL

---

#### KEY WORDS

liver; kidney; lung; brain; skeletal muscle; comparison with HPLC/MS

---

#### REFERENCE

Tracqui,A.; Kintz,P.; Cirimele,V.; Berthault,F.; Mangin,P.; Ludes,B. HPLC-DAD and HPLC-MS findings in a fatality involving (Z)-cis-clopenthixol (zuclopenthixol), *J.Anal.Toxicol.*, **1997**, 21, 314–318.

---

#### SAMPLE

**Matrix:** blood, urine

**Sample preparation:** Condition a 1 mL Bond Elut CN with 2 mL MeCN and 2 mL water. Add 2 mL plasma or urine to the SPE cartridge, wash with 2 mL water, elute with MeCN:n-butylamine 90:10. Evaporate the eluate to dryness under a stream of nitrogen at 40°, reconstitute the residue in 200 µL mobile phase, inject a 100 µL aliquot.

---

#### HPLC VARIABLES

**Column:** 120 × 4.6 Spherisorb S5 CN

**Mobile phase:** MeCN:200 mM pH 6.5 potassium phosphate buffer:water 36:5:59 containing 6 mM dodecyl-N,N,N-trimethylammonium bromide

**Column temperature:** 40

**Flow rate:** 1

**Injection volume:** 100

**Detector:** F ex 260 em 435 following post-column photolysis with a low-pressure 8 W mercury UV light in a 5 m × 0.5 mm i.d.coil of PTFE tubing

---

#### CHROMATOGRAM

**Retention time:** 9

**Limit of detection:** 0.05 ng/mL

---

#### OTHER SUBSTANCES

**Extracted:** metabolites

---

#### KEY WORDS

protect from light; plasma; SPE; post-column photochemical derivatization

---

#### REFERENCE

Hansen,B.B.; Hansen,S.H. Determination of zuclopenthixol and its main N-dealkylated metabolite in biological fluids using high-performance liquid chromatography with post-column photochemical derivatization and fluorescence detection, *J.Chromatogr.B*, **1994**, 658, 319–325.

---

#### SAMPLE

**Matrix:** blood, urine

**Sample preparation:** Add 1 mL whole blood or urine to Toxi-Tube A (Toxi-Lab, Irvine CA), add 3 mL water, mix by gentle inversion for 5 min, centrifuge at 1500 g for 5 min. Remove the organic layer and evaporate it to dryness under a stream of nitrogen at 40°, reconstitute the residue with 50 µL MeCN:water 50:50, vortex for 10 s, centrifuge at 7500 g for 2 min, inject a 10 (urine) or 30 (blood) µL aliquot. (The detector wavelength shown is the wavelength of maximum absorbance. This will not necessarily be the optimal wavelength for the separation. Multiple wavelengths from 200-350 nm can be scanned using a diode-array detector. Otherwise, 220 nm may be a reasonable choice for initial work. Matrix may interfere.)

---

#### HPLC VARIABLES

**Guard column:** 20 mm long Symmetry C18

**Column:** 250 × 4.6 5 µm Symmetry C8 (Waters)

**Mobile phase:** Gradient. A was 50 mM pH 3.8 sodium phosphate buffer. B was MeCN. A:B 85:15 for 6.5 min, 65:35 for 18.5 min, 20:80 for 3 min (step gradient), re-equilibrate at initial conditions for 7 min.

**Column temperature:** 30

**Flow rate:** 1 for 6.5 min, to 1.5 over 18.5 min, maintain at 1.5 for 3 min (re-equilibrate at 1.5 mL/min)

**Injection volume:** 10-30

**Detector:** UV 206.4

---

#### CHROMATOGRAM

**Retention time:** 16.325

---

#### KEY WORDS

whole blood

---

#### REFERENCE

Gaillard,Y.; Pépin,G. Use of high-performance liquid chromatography with photodiode-array UV detection for the creation of a 600-compound library. Application to forensic toxicology, *J.Chromatogr.A*, **1997**, 763, 149-163.

---

#### SAMPLE

**Matrix:** solutions

**Sample preparation:** Inject 1 mL onto column A. Elute column A onto column B with mobile phase for 30 s then remove it from the circuit. Elute column B with mobile phase and monitor the effluent.

---

#### HPLC VARIABLES

**Column:** A 10 × 6 packed with 40 µm material from a Bond Elut cartridge (cat. no. 620303); B 100 × 4 3 µm Spherisorb ODS Suprapac

**Mobile phase:** MeCN:85% phosphoric acid:triethylamine:water 49.55:0.225:0.225:50

**Flow rate:** 0.65

**Injection volume:** 1000

**Detector:** UV 238

---

#### CHROMATOGRAM

**Retention time:** 3.30

---

#### OTHER SUBSTANCES

**Simultaneous:** alprazolam, amitriptyline, chlorpromazine, chlorprothixene, clomipramine, des-clomipramine, desmethylinipramine, diazepam, flunitrazepam, imipramine, levomepromazine, maprotiline, nortriptyline, promethazine, thioridazine, thioridazine sulfone, thioridazine sulf-oxide, trimipramine, zimeldine

**Noninterfering:** carbamazepine, clonazepam, lorazepam, nitrazepam, oxazepam, phenytoin

**Interfering:** fluphenazine, haloperidol, perphenazine, protriptyline

---

#### KEY WORDS

column-switching

---

#### REFERENCE

Svensson,C.; Nyberg,G.; Mårtensson,E. High-performance liquid chromatographic quantitation of amitriptyline and nortriptyline in dialysate from plasma or serum using on-line solid-phase extraction, *J.Chromatogr.*, **1988**, 432, 363-369.